



Stationary Phases & Columns

BIO-Stationary Phases

Selection Guide

Peptides	Polar	Mid & non-polar	Hydrophobic	Natural, Fatty Acids
<p>< 40AA MW: up to 5KDa</p> <p>pH: 1.5 to 8.0</p> <p>max. pH: 10</p>	<p>puriFlash® BIO 100 C18N</p>	<p>puriFlash® BIO 100 C18T</p> <p>puriFlash® BIO 100 C18XS</p>	<p>Screening Of puriFlash® BIO 100 (C18N /C18T)</p> <p>puriFlash® BIO 100 C18XS</p>	
<p>< 80AA MW: up to 10KDa</p> <p>pH: 1.5 to 8.0</p> <p>max. pH: 10</p>	<p>puriFlash® BIO 200 C18N</p>	<p>puriFlash® BIO 200 C18T</p> <p>puriFlash® BIO 200 C18XS</p>	<p>Screening Of puriFlash® BIO 200 (C18N /C18T)</p> <p>puriFlash® BIO 200 C18XS</p>	
<p>< 160AA MW: up to 20KDa</p> <p>pH: 1.5 to 8.0</p>	<p>puriFlash® BIO 200 C18N</p>	<p>puriFlash® BIO 200 C8N</p>	<p>puriFlash® BIO 200 C8N</p>	
<p>< 80AA MW: up to 100KDa</p> <p>pH: 1.5 to 8.0</p>				<p>puriFlash® BIO 300 C4AQ</p>
<p>In-Process QA/QC of Peptides Synthesis</p>	<p>In-Process QA/QC of Peptides Synthesis puriFlash® BIO CS 2.6C18N => puriFlash® BIO 100 2.5C18N</p>			


Notes:

Polar Peptides => HILIC mode using higher % of ACN 95 -to- 85%

Hydrophobic Peptides => it is useful to work with Water/ACN using a few % Formic Acid or 0.05% TFA ~ pH 2. In case your peptides have Lysine, Arginine etc. it is better to have an alkali environment in the solvent. You need real buffer and according to buffer solubility it is to suggest to switch to MeOH instead of ACN. Usually step-Gradients (Ramp Gradients) or Pseudo-Isocratic or very flat gradients lead to highest capacity.




Peptides




puriFlash® BIO C18-N

100Å - 320m²/g
 2.5, 3.5, 5, 10, 15 & 30µm
 C18 - octadecyl
 Mono-functional
 %C: 15.0
 End-capping: None
 pH stability: 1.5 to 8.0
 Use mode: Reverse
In-Process QA/QC of Peptides Synthesis. Analysis & Purification of polar Peptides with less than 40AA & mw. up to 5KDa under pseudo hilic mode with 85% -to- 95% ACN. Analysis & Purification of hydrophobic Peptides with less than 40AA & mw. up to 5KDa.




puriFlash® BIO C18-T

100Å - 320m²/g
 2.5, 3.5, 5, 10, 15 & 30µm
 C18 - octadecyl
 Tri-functional
 %C: 17.0
 End-capping: One-step
 pH stability: 1.5 to 8.0
 Use mode: Reverse
Analysis & Purification of mid & non-polar Peptides, hydrophobic Peptides with less than 40AA & mw. up to 5KDa.




puriFlash® BIO C18-XS

100Å - 320m²/g
 2.5, 3.5, 5, 10, 15 & 30µm
 C18 - octadecyl
 Mono-functional
 %C: 17.0
 End-capping: Multi-step
 pH stability: 1.0 to 10.0
 Use mode: Reverse
Analysis & Purification of mid & non-polar Peptides, hydrophobic Peptides with less than 40AA & mw. up to 5KDa under basic conditions up to pH: 10.0



puriFlash® BIO C8-N

200Å - 200m²/g
 2.5, 3.5, 5, 10, 15 & 30µm
 C8 - octadecyl
 Mono-functional
 %C: 7.0
 End-capping: None
 pH stability: 1.5 to 8.0
 Use mode: Reverse
Analysis & Purification of polar Peptides less than 160AA & mw. up to 20KDa under pseudo hilic mode with 85% -to- 95% ACN. Analysis & Purification of hydrophobic Peptides with less than 80AA & mw. up to 10KDa.



puriFlash® BIO C4-AQ

300Å - 100m²/g
 3.5, 5, 10, 15 & 30µm
 C4 - butyl
 Mono-functional
 %C: 3.0
 End-capping: Mixte
 pH stability: 1.5 to 8.0
 Use mode: Reverse
Analysis & Purification of natural Peptides, fatty acids with larger than 80AA & mw. up to 100KDa.