



REFERENCE STANDARDS FOR RESIDUE ANALYSIS 2017

2017 | Volume 1.1

lgcstandards.com/foodandenvironment
dr.ehrenstorfer@lgcgroup.com

LGC Quality | ISO Guide 34 | ISO 9001 |
ISO/IEC 17025 | ISO/IEC 17043



Since 1975, Dr. Ehrenstorfer has led the way in producing pesticide reference standards. Today, our portfolio has adapted and expanded to meet changing regulations and technology.

You'll find our latest products in the Dr. Ehrenstorfer catalogue, including:

- Our leading range of neat materials. Those that are manufactured under ISO Guide 34 are labelled in this catalogue with the symbol ‡
- New pesticide and pesticide metabolites
- Stable isotope labelled standards for analysis using mass spectrometry
- A wide range of veterinary and pharmaceutical residue standards, including marker metabolites
- Popular mixtures for EPA and other regulatory methods



Order at
lgcstandards.com/foodandenvironment

email
dr.ehrenstorfer@lgcgroup.com
or find contact details
for your local sales office.

Valentin Wagner, Head of RM Manufacturing

TECH TIPS

To help our customers, we've captured the most common questions asked about our products. We're also creating a series of Dr. Ehrenstorfer Tech Tips videos and PDFs to help you.

Our most frequently asked questions include:

1

How much material is in the bottle/ampoule?

Dr. Ehrenstorfer neat reference materials are supplied with a nominal weight, and the vial or bottle is typically overfilled with 10% more of the product than stated. However, this means the exact amount of material supplied cannot be provided. In order to use the material, it is usually practical to prepare a solution.

The preparation procedure described below can be used to calculate the amount of material present.

2

The container looks empty – is there anything inside?

Where small quantities of solid material are supplied, this can be dispersed over the inside surface of the container. Liquid may also coat the inner surface of the container which may not be visible. Therefore, please refer to the following procedure to extract the material from the container.

3

How can I prepare a solution from neat Dr. Ehrenstorfer product?

In order to recover all the material from the container the following procedure can be used. We recommend storing the vial / bottle in an upright position for at least 24 hours prior to handling.

- 1 Ensure the container is clean and dry. Using an analytical balance, weigh the container (including the lid for bottles) and record the weight.
- 2 Using a suitable solvent, carefully transfer the contents to a volumetric flask. Rinse the container (including the lid) at least three times and combine all rinses to ensure a complete transfer of material.
- 3 Dry the empty container completely, then using an analytical balance, weigh the dry, empty container and lid and record the weight.
- 4 Calculate the difference in weight between the first and second weighing. The difference in weight is the amount of material that has been transferred.
- 5 Make the solution up to volume in the volumetric flask. The concentration of the solution can then be determined.

Where larger quantities of material are supplied, you may wish to only transfer an aliquot of the material. In these cases, you may need to use a weighing boat and a spatula to weigh the material before transferring it into the volumetric flask.



To access more FAQs, visit
lgcstandards.com

QUALITY AND ACCREDITATION

1 Certification

Dr. Ehrenstorfer use a quality management system certified to **ISO 9001** and methodology accredited to **ISO/IEC 17025** to manufacture products. Traceability back to the SI unit is demonstrated by using accredited procedures measuring weight, volume and temperature.

2 Gravimetric certification

Dr. Ehrenstorfer solution standards are designed for calibration or recalibration of chromatographic or other analytical systems used for quantitative or qualitative purposes. The products may also be used as reference standards for interlaboratory studies and to validate analytical procedures.

4 ISO Guide 34

A range of materials are produced under Dr. Ehrenstorfer's **ISO Guide 34** scope of accreditation. These materials are identified in this catalogue by the symbol ‡ and lot numbers are prefixed with the letter "G". This range is constantly increasing, therefore if you are unable to find the product that you require, please contact your local sales office or email us at dr.ehrenstorfer@lgcgroup.com.

The main difference between products produced under **ISO Guide 34** scope of accreditation and other products, is that the homogeneity testing and the calculation of the measurement uncertainties are performed strictly according to **ISO Guide 34** and **ISO Guide 35**. These products are provided with a Certificate of Analysis, which fulfills the requirements of **ISO Guide 34** and **ISO Guide 31**.

All raw materials used to prepare the multicomponent mixtures are of the highest purity. After our production process, each material is checked using several of the following methods, if applicable: UV, IR spectroscopy, elemental analysis, chromatographic properties (GC/FID, GC/ECD, HPLC/DAD, TLC) and physical properties (phase, colour, odour, melting point) and Karl Fischer titration for detection of water content.

3 Measurement uncertainty

In accordance with the Guide to the Expression of Uncertainty in Measurement (GUM), our Certificates of Analysis list the expanded concentration uncertainty at the 95% confidence level, using a coverage factor of $k = 2$. The expanded uncertainty is based on the combination of uncertainties associated with the individual steps involved in the preparation of the product e.g. purity, material density and weighing technique.

5 Expiry dates

All materials are supplied with an assigned expiry date which is printed on the certificate of analysis and the product label. Expiry dates are determined by real time stability testing.

6 Storage and transportation

Every reference material is stored under controlled conditions as stated on the certificate. These conditions are for the long term storage of the material. Normally products are not shipped under controlled conditions as shipping times are generally **<72 hours** and therefore short term. One sample of each lot is kept aside to enable checks on the specific lot to be undertaken if required.

HOW TO READ YOUR CERTIFICATE OF ANALYSIS

Dr. Ehrenstorfer Certificates of Analysis provide a full description of the material to which they relate and summarise the analyses undertaken during the characterisation process. Current Certificates of Analysis for individual lots of products in the range are available on our website, www.lgcstandards.com on the relevant product pages. If you are unable to find the product that you require please contact one of our local sales offices.

Products produced under Dr. Ehrenstorfer's **ISO Guide 34** scope of accreditation are identified by a lot number pre-fixed by the letter "G".

An example of a Certificate of Analysis for a neat product is given below. For certificates of **gravimetric solutions**, the following information is included:

- **Gravimetric data**
Concentration of the product, purity and the weight of product.
- **Solvent information**
Identity, lot number and exact quantity of solvent used.
- **Traceability data**
Identification of materials used including lot numbers for any neat and solution products used.

Article Code
The unique identification code for the product is given in the certificate.

Formula
Molecular formula of the material.

Mol. weight
Molecular weight of the material.

CAS No.
A unique identifier assigned by the Chemical Abstracts Service.

Note 1 Storage conditions
These are the minimum storage requirements based on stability studies.

Note 2 Tolerance / Uncertainty
For products produced under ISO Guide 34 scope of accreditation the uncertainty is calculated in accordance with ISO Guide 34 and RUECHEM/ CITAC Guide - Quantifying Uncertainty in Analytical Measurement, Second Edition.

Lot number
Identification number for a specific batch of material.

Expiry date
Determined by real time stability testing.

Storage Temperature
Long time storage solutions. See Note 1.

Purity
% purity that has been determined for this particular lot of this product.

Tolerance / Expanded Uncertainty
See Note 2.

Analytical data
Details of the methodology used to determine the purity of this particular lot product. Analytical chromatograms are also supplied where appropriate on supplementary pages.

Certificate of Analysis
ISO Guide 34 Reference Material
Product Identification
Article Code: DRE-C16080000
Article Name: Phosphate buffer solution
Formula: C₁₂H₂₂O₁₁P₃
Mol. Weight: 276.06
CAS No.: 2388-05-6

Lot Number: G142157
Expiry Date: 06.06.2020
Storage Temperature: 4°C ± 4°C


Purity: 99.91% (g/g)
Expanded Uncertainty U₉₅: 0.08% (g/g)

The uncertainty of this standard is calculated in accordance with the ISO Guide 34 and EURACHEM/CITAC Guide - Quantifying Uncertainty in Analytical Measurement, Second Edition. The expanded uncertainty is U₉₅(g/g) = 0.08% (g/g), where k is the coverage factor at the 95% confidence level (k=2). Uncertainty U₉₅(g/g) is based on the combination of the uncertainties associated with each individual operation involved in the analysis of the product: u₁(g/g) = u₁(g/g) + u₂(g/g) + u₃(g/g) + u₄(g/g) + u₅(g/g) + u₆(g/g) + u₇(g/g) + u₈(g/g) + u₉(g/g) + u₁₀(g/g) + u₁₁(g/g) + u₁₂(g/g) + u₁₃(g/g) + u₁₄(g/g) + u₁₅(g/g) + u₁₆(g/g) + u₁₇(g/g) + u₁₈(g/g) + u₁₉(g/g) + u₂₀(g/g) + u₂₁(g/g) + u₂₂(g/g) + u₂₃(g/g) + u₂₄(g/g) + u₂₅(g/g) + u₂₆(g/g) + u₂₇(g/g) + u₂₈(g/g) + u₂₉(g/g) + u₃₀(g/g) + u₃₁(g/g) + u₃₂(g/g) + u₃₃(g/g) + u₃₄(g/g) + u₃₅(g/g) + u₃₆(g/g) + u₃₇(g/g) + u₃₈(g/g) + u₃₉(g/g) + u₄₀(g/g) + u₄₁(g/g) + u₄₂(g/g) + u₄₃(g/g) + u₄₄(g/g) + u₄₅(g/g) + u₄₆(g/g) + u₄₇(g/g) + u₄₈(g/g) + u₄₉(g/g) + u₅₀(g/g) + u₅₁(g/g) + u₅₂(g/g) + u₅₃(g/g) + u₅₄(g/g) + u₅₅(g/g) + u₅₆(g/g) + u₅₇(g/g) + u₅₈(g/g) + u₅₉(g/g) + u₆₀(g/g) + u₆₁(g/g) + u₆₂(g/g) + u₆₃(g/g) + u₆₄(g/g) + u₆₅(g/g) + u₆₆(g/g) + u₆₇(g/g) + u₆₈(g/g) + u₆₉(g/g) + u₇₀(g/g) + u₇₁(g/g) + u₇₂(g/g) + u₇₃(g/g) + u₇₄(g/g) + u₇₅(g/g) + u₇₆(g/g) + u₇₇(g/g) + u₇₈(g/g) + u₇₉(g/g) + u₈₀(g/g) + u₈₁(g/g) + u₈₂(g/g) + u₈₃(g/g) + u₈₄(g/g) + u₈₅(g/g) + u₈₆(g/g) + u₈₇(g/g) + u₈₈(g/g) + u₈₉(g/g) + u₉₀(g/g) + u₉₁(g/g) + u₉₂(g/g) + u₉₃(g/g) + u₉₄(g/g) + u₉₅(g/g) + u₉₆(g/g) + u₉₇(g/g) + u₉₈(g/g) + u₉₉(g/g) + u₁₀₀(g/g) + u₁₀₁(g/g) + u₁₀₂(g/g) + u₁₀₃(g/g) + u₁₀₄(g/g) + u₁₀₅(g/g) + u₁₀₆(g/g) + u₁₀₇(g/g) + u₁₀₈(g/g) + u₁₀₉(g/g) + u₁₁₀(g/g) + u₁₁₁(g/g) + u₁₁₂(g/g) + u₁₁₃(g/g) + u₁₁₄(g/g) + u₁₁₅(g/g) + u₁₁₆(g/g) + u₁₁₇(g/g) + u₁₁₈(g/g) + u₁₁₉(g/g) + u₁₂₀(g/g) + u₁₂₁(g/g) + u₁₂₂(g/g) + u₁₂₃(g/g) + u₁₂₄(g/g) + u₁₂₅(g/g) + u₁₂₆(g/g) + u₁₂₇(g/g) + u₁₂₈(g/g) + u₁₂₉(g/g) + u₁₃₀(g/g) + u₁₃₁(g/g) + u₁₃₂(g/g) + u₁₃₃(g/g) + u₁₃₄(g/g) + u₁₃₅(g/g) + u₁₃₆(g/g) + u₁₃₇(g/g) + u₁₃₈(g/g) + u₁₃₉(g/g) + u₁₄₀(g/g) + u₁₄₁(g/g) + u₁₄₂(g/g) + u₁₄₃(g/g) + u₁₄₄(g/g) + u₁₄₅(g/g) + u₁₄₆(g/g) + u₁₄₇(g/g) + u₁₄₈(g/g) + u₁₄₉(g/g) + u₁₅₀(g/g) + u₁₅₁(g/g) + u₁₅₂(g/g) + u₁₅₃(g/g) + u₁₅₄(g/g) + u₁₅₅(g/g) + u₁₅₆(g/g) + u₁₅₇(g/g) + u₁₅₈(g/g) + u₁₅₉(g/g) + u₁₆₀(g/g) + u₁₆₁(g/g) + u₁₆₂(g/g) + u₁₆₃(g/g) + u₁₆₄(g/g) + u₁₆₅(g/g) + u₁₆₆(g/g) + u₁₆₇(g/g) + u₁₆₈(g/g) + u₁₆₉(g/g) + u₁₇₀(g/g) + u₁₇₁(g/g) + u₁₇₂(g/g) + u₁₇₃(g/g) + u₁₇₄(g/g) + u₁₇₅(g/g) + u₁₇₆(g/g) + u₁₇₇(g/g) + u₁₇₈(g/g) + u₁₇₉(g/g) + u₁₈₀(g/g) + u₁₈₁(g/g) + u₁₈₂(g/g) + u₁₈₃(g/g) + u₁₈₄(g/g) + u₁₈₅(g/g) + u₁₈₆(g/g) + u₁₈₇(g/g) + u₁₈₈(g/g) + u₁₈₉(g/g) + u₁₉₀(g/g) + u₁₉₁(g/g) + u₁₉₂(g/g) + u₁₉₃(g/g) + u₁₉₄(g/g) + u₁₉₅(g/g) + u₁₉₆(g/g) + u₁₉₇(g/g) + u₁₉₈(g/g) + u₁₉₉(g/g) + u₂₀₀(g/g) + u₂₀₁(g/g) + u₂₀₂(g/g) + u₂₀₃(g/g) + u₂₀₄(g/g) + u₂₀₅(g/g) + u₂₀₆(g/g) + u₂₀₇(g/g) + u₂₀₈(g/g) + u₂₀₉(g/g) + u₂₁₀(g/g) + u₂₁₁(g/g) + u₂₁₂(g/g) + u₂₁₃(g/g) + u₂₁₄(g/g) + u₂₁₅(g/g) + u₂₁₆(g/g) + u₂₁₇(g/g) + u₂₁₈(g/g) + u₂₁₉(g/g) + u₂₂₀(g/g) + u₂₂₁(g/g) + u₂₂₂(g/g) + u₂₂₃(g/g) + u₂₂₄(g/g) + u₂₂₅(g/g) + u₂₂₆(g/g) + u₂₂₇(g/g) + u₂₂₈(g/g) + u₂₂₉(g/g) + u₂₃₀(g/g) + u₂₃₁(g/g) + u₂₃₂(g/g) + u₂₃₃(g/g) + u₂₃₄(g/g) + u₂₃₅(g/g) + u₂₃₆(g/g) + u₂₃₇(g/g) + u₂₃₈(g/g) + u₂₃₉(g/g) + u₂₄₀(g/g) + u₂₄₁(g/g) + u₂₄₂(g/g) + u₂₄₃(g/g) + u₂₄₄(g/g) + u₂₄₅(g/g) + u₂₄₆(g/g) + u₂₄₇(g/g) + u₂₄₈(g/g) + u₂₄₉(g/g) + u₂₅₀(g/g) + u₂₅₁(g/g) + u₂₅₂(g/g) + u₂₅₃(g/g) + u₂₅₄(g/g) + u₂₅₅(g/g) + u₂₅₆(g/g) + u₂₅₇(g/g) + u₂₅₈(g/g) + u₂₅₉(g/g) + u₂₆₀(g/g) + u₂₆₁(g/g) + u₂₆₂(g/g) + u₂₆₃(g/g) + u₂₆₄(g/g) + u₂₆₅(g/g) + u₂₆₆(g/g) + u₂₆₇(g/g) + u₂₆₈(g/g) + u₂₆₉(g/g) + u₂₇₀(g/g) + u₂₇₁(g/g) + u₂₇₂(g/g) + u₂₇₃(g/g) + u₂₇₄(g/g) + u₂₇₅(g/g) + u₂₇₆(g/g) + u₂₇₇(g/g) + u₂₇₈(g/g) + u₂₇₉(g/g) + u₂₈₀(g/g) + u₂₈₁(g/g) + u₂₈₂(g/g) + u₂₈₃(g/g) + u₂₈₄(g/g) + u₂₈₅(g/g) + u₂₈₆(g/g) + u₂₈₇(g/g) + u₂₈₈(g/g) + u₂₈₉(g/g) + u₂₉₀(g/g) + u₂₉₁(g/g) + u₂₉₂(g/g) + u₂₉₃(g/g) + u₂₉₄(g/g) + u₂₉₅(g/g) + u₂₉₆(g/g) + u₂₉₇(g/g) + u₂₉₈(g/g) + u₂₉₉(g/g) + u₃₀₀(g/g) + u₃₀₁(g/g) + u₃₀₂(g/g) + u₃₀₃(g/g) + u₃₀₄(g/g) + u₃₀₅(g/g) + u₃₀₆(g/g) + u₃₀₇(g/g) + u₃₀₈(g/g) + u₃₀₉(g/g) + u₃₁₀(g/g) + u₃₁₁(g/g) + u₃₁₂(g/g) + u₃₁₃(g/g) + u₃₁₄(g/g) + u₃₁₅(g/g) + u₃₁₆(g/g) + u₃₁₇(g/g) + u₃₁₈(g/g) + u₃₁₉(g/g) + u₃₂₀(g/g) + u₃₂₁(g/g) + u₃₂₂(g/g) + u₃₂₃(g/g) + u₃₂₄(g/g) + u₃₂₅(g/g) + u₃₂₆(g/g) + u₃₂₇(g/g) + u₃₂₈(g/g) + u₃₂₉(g/g) + u₃₃₀(g/g) + u₃₃₁(g/g) + u₃₃₂(g/g) + u₃₃₃(g/g) + u₃₃₄(g/g) + u₃₃₅(g/g) + u₃₃₆(g/g) + u₃₃₇(g/g) + u₃₃₈(g/g) + u₃₃₉(g/g) + u₃₄₀(g/g) + u₃₄₁(g/g) + u₃₄₂(g/g) + u₃₄₃(g/g) + u₃₄₄(g/g) + u₃₄₅(g/g) + u₃₄₆(g/g) + u₃₄₇(g/g) + u₃₄₈(g/g) + u₃₄₉(g/g) + u₃₅₀(g/g) + u₃₅₁(g/g) + u₃₅₂(g/g) + u₃₅₃(g/g) + u₃₅₄(g/g) + u₃₅₅(g/g) + u₃₅₆(g/g) + u₃₅₇(g/g) + u₃₅₈(g/g) + u₃₅₉(g/g) + u₃₆₀(g/g) + u₃₆₁(g/g) + u₃₆₂(g/g) + u₃₆₃(g/g) + u₃₆₄(g/g) + u₃₆₅(g/g) + u₃₆₆(g/g) + u₃₆₇(g/g) + u₃₆₈(g/g) + u₃₆₉(g/g) + u₃₇₀(g/g) + u₃₇₁(g/g) + u₃₇₂(g/g) + u₃₇₃(g/g) + u₃₇₄(g/g) + u₃₇₅(g/g) + u₃₇₆(g/g) + u₃₇₇(g/g) + u₃₇₈(g/g) + u₃₇₉(g/g) + u₃₈₀(g/g) + u₃₈₁(g/g) + u₃₈₂(g/g) + u₃₈₃(g/g) + u₃₈₄(g/g) + u₃₈₅(g/g) + u₃₈₆(g/g) + u₃₈₇(g/g) + u₃₈₈(g/g) + u₃₈₉(g/g) + u₃₉₀(g/g) + u₃₉₁(g/g) + u₃₉₂(g/g) + u₃₉₃(g/g) + u₃₉₄(g/g) + u₃₉₅(g/g) + u₃₉₆(g/g) + u₃₉₇(g/g) + u₃₉₈(g/g) + u₃₉₉(g/g) + u₄₀₀(g/g) + u₄₀₁(g/g) + u₄₀₂(g/g) + u₄₀₃(g/g) + u₄₀₄(g/g) + u₄₀₅(g/g) + u₄₀₆(g/g) + u₄₀₇(g/g) + u₄₀₈(g/g) + u₄₀₉(g/g) + u₄₁₀(g/g) + u₄₁₁(g/g) + u₄₁₂(g/g) + u₄₁₃(g/g) + u₄₁₄(g/g) + u₄₁₅(g/g) + u₄₁₆(g/g) + u₄₁₇(g/g) + u₄₁₈(g/g) + u₄₁₉(g/g) + u₄₂₀(g/g) + u₄₂₁(g/g) + u₄₂₂(g/g) + u₄₂₃(g/g) + u₄₂₄(g/g) + u₄₂₅(g/g) + u₄₂₆(g/g) + u₄₂₇(g/g) + u₄₂₈(g/g) + u₄₂₉(g/g) + u₄₃₀(g/g) + u₄₃₁(g/g) + u₄₃₂(g/g) + u₄₃₃(g/g) + u₄₃₄(g/g) + u₄₃₅(g/g) + u₄₃₆(g/g) + u₄₃₇(g/g) + u₄₃₈(g/g) + u₄₃₉(g/g) + u₄₄₀(g/g) + u₄₄₁(g/g) + u₄₄₂(g/g) + u₄₄₃(g/g) + u₄₄₄(g/g) + u₄₄₅(g/g) + u₄₄₆(g/g) + u₄₄₇(g/g) + u₄₄₈(g/g) + u₄₄₉(g/g) + u₄₅₀(g/g) + u₄₅₁(g/g) + u₄₅₂(g/g) + u₄₅₃(g/g) + u₄₅₄(g/g) + u₄₅₅(g/g) + u₄₅₆(g/g) + u₄₅₇(g/g) + u₄₅₈(g/g) + u₄₅₉(g/g) + u₄₆₀(g/g) + u₄₆₁(g/g) + u₄₆₂(g/g) + u₄₆₃(g/g) + u₄₆₄(g/g) + u₄₆₅(g/g) + u₄₆₆(g/g) + u₄₆₇(g/g) + u₄₆₈(g/g) + u₄₆₉(g/g) + u₄₇₀(g/g) + u₄₇₁(g/g) + u₄₇₂(g/g) + u₄₇₃(g/g) + u₄₇₄(g/g) + u₄₇₅(g/g) + u₄₇₆(g/g) + u₄₇₇(g/g) + u₄₇₈(g/g) + u₄₇₉(g/g) + u₄₈₀(g/g) + u₄₈₁(g/g) + u₄₈₂(g/g) + u₄₈₃(g/g) + u₄₈₄(g/g) + u₄₈₅(g/g) + u₄₈₆(g/g) + u₄₈₇(g/g) + u₄₈₈(g/g) + u₄₈₉(g/g) + u₄₉₀(g/g) + u₄₉₁(g/g) + u₄₉₂(g/g) + u₄₉₃(g/g) + u₄₉₄(g/g) + u₄₉₅(g/g) + u₄₉₆(g/g) + u₄₉₇(g/g) + u₄₉₈(g/g) + u₄₉₉(g/g) + u₅₀₀(g/g) + u₅₀₁(g/g) + u₅₀₂(g/g) + u₅₀₃(g/g) + u₅₀₄(g/g) + u₅₀₅(g/g) + u₅₀₆(g/g) + u₅₀₇(g/g) + u₅₀₈(g/g) + u₅₀₉(g/g) + u₅₁₀(g/g) + u₅₁₁(g/g) + u₅₁₂(g/g) + u₅₁₃(g/g) + u₅₁₄(g/g) + u₅₁₅(g/g) + u₅₁₆(g/g) + u₅₁₇(g/g) + u₅₁₈(g/g) + u₅₁₉(g/g) + u₅₂₀(g/g) + u₅₂₁(g/g) + u₅₂₂(g/g) + u₅₂₃(g/g) + u₅₂₄(g/g) + u₅₂₅(g/g) + u₅₂₆(g/g) + u₅₂₇(g/g) + u₅₂₈(g/g) + u₅₂₉(g/g) + u₅₃₀(g/g) + u₅₃₁(g/g) + u₅₃₂(g/g) + u₅₃₃(g/g) + u₅₃₄(g/g) + u₅₃₅(g/g) + u₅₃₆(g/g) + u₅₃₇(g/g) + u₅₃₈(g/g) + u₅₃₉(g/g) + u₅₄₀(g/g) + u₅₄₁(g/g) + u₅₄₂(g/g) + u₅₄₃(g/g) + u₅₄₄(g/g) + u₅₄₅(g/g) + u₅₄₆(g/g) + u₅₄₇(g/g) + u₅₄₈(g/g) + u₅₄₉(g/g) + u₅₅₀(g/g) + u₅₅₁(g/g) + u₅₅₂(g/g) + u₅₅₃(g/g) + u₅₅₄(g/g) + u₅₅₅(g/g) + u₅₅₆(g/g) + u₅₅₇(g/g) + u₅₅₈(g/g) + u₅₅₉(g/g) + u₅₆₀(g/g) + u₅₆₁(g/g) + u₅₆₂(g/g) + u₅₆₃(g/g) + u₅₆₄(g/g) + u₅₆₅(g/g) + u₅₆₆(g/g) + u₅₆₇(g/g) + u₅₆₈(g/g) + u₅₆₉(g/g) + u₅₇₀(g/g) + u₅₇₁(g/g) + u₅₇₂(g/g) + u₅₇₃(g/g) + u₅₇₄(g/g) + u₅₇₅(g/g) + u₅₇₆(g/g) + u₅₇₇(g/g) + u₅₇₈(g/g) + u₅₇₉(g/g) + u₅₈₀(g/g) + u₅₈₁(g/g) + u₅₈₂(g/g) + u₅₈₃(g/g) + u₅₈₄(g/g) + u₅₈₅(g/g) + u₅₈₆(g/g) + u₅₈₇(g/g) + u₅₈₈(g/g) + u₅₈₉(g/g) + u₅₉₀(g/g) + u₅₉₁(g/g) + u₅₉₂(g/g) + u₅₉₃(g/g) + u₅₉₄(g/g) + u₅₉₅(g/g) + u₅₉₆(g/g) + u₅₉₇(g/g) + u₅₉₈(g/g) + u₅₉₉(g/g) + u₆₀₀(g/g) + u₆₀₁(g/g) + u₆₀₂(g/g) + u₆₀₃(g/g) + u₆₀₄(g/g) + u₆₀₅(g/g) + u₆₀₆(g/g) + u₆₀₇(g/g) + u₆₀₈(g/g) + u₆₀₉(g/g) + u₆₁₀(g/g) + u₆₁₁(g/g) + u₆₁₂(g/g) + u₆₁₃(g/g) + u₆₁₄(g/g) + u₆₁₅(g/g) + u₆₁₆(g/g) + u₆₁₇(g/g) + u₆₁₈(g/g) + u₆₁₉(g/g) + u₆₂₀(g/g) + u₆₂₁(g/g) + u₆₂₂(g/g) + u₆₂₃(g/g) + u₆₂₄(g/g) + u₆₂₅(g/g) + u₆₂₆(g/g) + u₆₂₇(g/g) + u₆₂₈(g/g) + u₆₂₉(g/g) + u₆₃₀(g/g) + u₆₃₁(g/g) + u₆₃₂(g/g) + u₆₃₃(g/g) + u₆₃₄(g/g) + u₆₃₅(g/g) + u₆₃₆(g/g) + u₆₃₇(g/g) + u₆₃₈(g/g) + u₆₃₉(g/g) + u₆₄₀(g/g) + u₆₄₁(g/g) + u₆₄₂(g/g) + u₆₄₃(g/g) + u₆₄₄(g/g) + u₆₄₅(g/g) + u₆₄₆(g/g) + u₆₄₇(g/g) + u₆₄₈(g/g) + u₆₄₉(g/g) + u₆₅₀(g/g) + u₆₅₁(g/g) + u₆₅₂(g/g) + u₆₅₃(g/g) + u₆₅₄(g/g) + u₆₅₅(g/g) + u₆₅₆(g/g) + u₆₅₇(g/g) + u₆₅₈(g/g) + u₆₅₉(g/g) + u₆₆₀(g/g) + u₆₆₁(g/g) + u₆₆₂(g/g) + u₆₆₃(g/g) + u₆₆₄(g/g) + u₆₆₅(g/g) + u₆₆₆(g/g) + u₆₆₇(g/g) + u₆₆₈(g/g) + u₆₆₉(g/g) + u₆₇₀(g/g) + u₆₇₁(g/g) + u₆₇₂(g/g) + u₆₇₃(g/g) + u₆₇₄(g/g) + u₆₇₅(g/g) + u₆₇₆(g/g) + u₆₇₇(g/g) + u₆₇₈(g/g) + u₆₇₉(g/g) + u₆₈₀(g/g) + u₆₈₁(g/g) + u₆₈₂(g/g) + u₆₈₃(g/g) + u₆₈₄(g/g) + u₆₈₅(g/g) + u₆₈₆(g/g) + u₆₈₇(g/g) + u₆₈₈(g/g) + u₆₈₉(g/g) + u₆₉₀(g/g) + u₆₉₁(g/g) + u₆₉₂(g/g) + u₆₉₃(g/g) + u₆₉₄(g/g) + u₆₉₅(g/g) + u₆₉₆(g/g) + u₆₉₇(g/g) + u₆₉₈(g/g) + u₆₉₉(g/g) + u₇₀₀(g/g) + u₇₀₁(g/g) + u₇₀₂(g/g) + u₇₀₃(g/g) + u₇₀₄(g/g) + u₇₀₅(g/g) + u₇₀₆(g/g) + u₇₀₇(g/g) + u₇₀₈(g/g) + u₇₀₉(g/g) + u₇₁₀(g/g) + u₇₁₁(g/g) + u₇₁₂(g/g) + u₇₁₃(g/g) + u₇₁₄(g/g) + u₇₁₅(g/g) + u₇₁₆(g/g) + u₇₁₇(g/g) + u₇₁₈(g/g) + u₇₁₉(g/g) + u₇₂₀(g/g) + u₇₂₁(g/g) + u₇₂₂(g/g) + u₇₂₃(g/g) + u₇₂₄(g/g) + u₇₂₅(g/g)

ABOUT THIS CATALOGUE

Dr. Ehrenstorfer is a leading manufacturer of pesticide standards and other organic reference materials. We typically have over **7,000 different products available**, in a variety of formats, from neat/powder materials to solutions of individual compounds and mixtures of compounds.

Our product codes



Pre-fix

Each product code has a pre-fix that indicates the certification level and how it is supplied	
DRE-C Neat product	DRE-C/G/L/X/Y/Z A Product supplied in an ampoule
DRE-E Neat product supplied without certificate	DRE-G/L/X/Y/Z S Product supplied as a set or in multiple units
DRE-G/L/X/Y/Z Solution product	

Solvent code

Each product prepared in solution ends with a two letter code to indicate the solvent used. Commonly used codes are given below:

AB	Acetone/Benzene
AC	Acetone
AD	Acetone D6
AL	Acetonitrile
AM	Acetone/Methanol 1:1
BD	Benzene/Dichloromethane 1:1
CS	Carbondisulfide
CT	Cyclohexane/Toluene 1:1
CY	Cyclohexane
DA	N,N-Dimethylacetamide
DI	Dichloromethane
DM	DMF
DO	Deuteriumoxide
DS	DMSO
EA	Ethylacetate
EL	Ethanol
HA	n-Hexane/Acetone 94:6
HC	Hexane:Carbon disulfide (80:20)
HE	n-Hexane
HP	n-Heptane
IO	Iso-Octane
IP	Isopropanol
IT	Iso-Octane/Toluene
MB	Methyl-tert-butyl ether
ME	Methanol
MW	Methanol/Water
OG	Oxygenate free gasoline
PE	n-Pentane
PY	Pyridine
SP	Benzene:MeCl2:ACN (2:2:1)
TH	Toluene/Hexane
TN	Triacetin
TO	Toluene
WA	Water



CONTENTS

Key product categories

Food analysis

Pesticides and metabolites: Fungicides	1
Pesticides and metabolites: Herbicides	33
Pesticides and metabolites: Insecticides	104
Other pesticides and metabolites	172
Pesticide mixtures	195
Food related compounds	457
Pharmaceutical and Veterinary compounds and metabolites	312
Dyes and metabolites	400
Food packaging contaminants	484

Environmental testing

Pesticides and metabolites: Fungicides	1
Pesticides and metabolites: Herbicides	33
Pesticides and metabolites: Insecticides	104
Other pesticides and metabolites	172
Pesticide mixtures	195
PCB's and related compounds	214
Polycyclic aromatic hydrocarbons (PAHs)	237
Hydrocarbons and Petrochemicals	263
Flame retardants - polybrominated biphenyls	309
Pharmaceutical and Veterinary compounds and metabolites	312
Volatile organic compounds (VOCs)	383
Dyes and metabolites	400
Phenol and aromatic compounds	420
Additional organic contaminants	560
Standards for regulatory methods	495
Stable isotope labelled compounds	523



How to order on our Webshop

An up to date listing of all products is available on the LGC website lgcstandards.com.

Once registered, you can use the website to order products, check delivery times and access a range of resources for a wide variety of sectors.









LGC Standards' technical staff are available to advise on the use and suitability of a particular product. Customers requiring assistance with the use or application of a particular reference substance should contact their local office.

For our latest product offering, check out our continually updated [online catalogue](#).










PESTICIDES AND METABOLITES: FUNGICIDES












Pesticide and metabolites: Fungicides

Product code	Description			
Acibenzolar Acid (1,2,3-Benzothiadiazole-7-carboxylic acid)				
CAS 35272-27-6 DRE-C10027900	MW 180.1839 Acibenzolar acid	C ₇ H ₄ N ₂ O ₂ S	50mg	
Acibenzolar-S-methyl				
CAS 135158-54-2 DRE-C10028000 DRE-L10028000CY	MW 210.276 Acibenzolar-S-methyl Acibenzolar-S-methyl 10 µg/mL in Cyclohexane	C ₈ H ₆ N ₂ OS ₂	100mg 10ml	
Aldimorph				
CAS 91315-15-0 DRE-C10083000	MW n/a Aldimorph		20mg	No Structure
Ametoctradin				
CAS 865318-97-4 DRE-C10148900	MW 275.3925 Ametoctradin(‡)	C ₁₅ H ₂₅ N ₅	100mg	
2-Aminobenzimidazole				
CAS 934-32-7 DRE-C10170000	MW 133.1506 2-Aminobenzimidazole	C ₇ H ₇ N ₃	250mg	
Amisulbrom				
CAS 348635-87-0 DRE-C10229000	MW 466.3058 Amisulbrom(‡)	C ₁₃ H ₁₃ BrFN ₅ O ₄ S ₂	25mg	
Amobam				
CAS 3566-10-7 DRE-C10242000	MW 246.4408 Amobam	C ₄ H ₆ N ₂ S ₄ ·2H ₄ N	50mg	
Anilazine				
CAS 101-05-3 DRE-C10260000	MW 275.5218 Anilazine	C ₉ H ₅ Cl ₃ N ₄	250mg	
Azaconazole				
CAS 60207-31-0 DRE-C10339000 DRE-L10339000AL DRE-L10339000CY DRE-XA09010150ME	MW 300.1406 Azaconazole Azaconazole 10 µg/mL in Acetonitrile Azaconazole 10 µg/mL in Cyclohexane Azaconazole 100 µg/mL in Methanol(‡)	C ₁₂ H ₁₁ Cl ₂ N ₃ O ₂	100mg 10ml 10ml 1ml	









Pesticide and metabolites: Fungicides

Product code	Description			
Azoxystrobin				
CAS 131860-33-8 DRE-C10413000 DRE-L10413000CY	MW 403.3875 Azoxystrobin(±) Azoxystrobin 10 µg/mL in Cyclohexane	C ₂₂ H ₁₇ N ₃ O ₅	100mg 10ml	
Azoxystrobin (free acid)				
CAS 1185255-09-7 DRE-LA10413200AL	MW 389.3609 Azoxystrobin (free acid) 10 µg/mL in Acetonitrile (*)	C ₂₁ H ₁₅ N ₃ O ₅	1ml	
Benalaxyl				
CAS 71626-11-4 DRE-C10440000 DRE-L10440000CY DRE-XA10440000CY	MW 325.4015 Benalaxyl(±) Benalaxyl 10 µg/mL in Cyclohexane Benalaxyl 100 µg/mL in Cyclohexane	C ₂₀ H ₂₃ NO ₃	250mg 10ml 1ml	
Benalaxyl-M				
CAS 98243-83-5 DRE-C10442000	MW 325.4015 Benalaxyl-M	C ₂₀ H ₂₃ NO ₃	10mg	
Benodanil				
CAS 15310-01-7 DRE-C10480000	MW 323.1291 Benodanil	C ₁₃ H ₁₀ INO	250mg	
Benomyl				
CAS 17804-35-2 DRE-C10490000	MW 290.3177 Benomyl	C ₁₄ H ₁₈ N ₄ O ₃	250mg	
Benthiavalicarb-isopropyl				
CAS 177406-68-7 DRE-C10516000 DRE-XA10516000AL	MW 381.4649 Benthiavalicarb-isopropyl(±) Benthiavalicarb-isopropyl 100 µg/mL in Acetonitrile	C ₁₈ H ₂₄ FN ₃ O ₃ S	10mg 1ml	
Benzalkonium Chloride				
CAS 8001-54-5 DRE-C10532200	MW 227.7735 Benzalkonium chloride	C ₁₁ H ₁₈ N(C ₂ H ₄) _n .Cl	100mg	
Bethoxazin				
CAS 163269-30-5 DRE-C10576000	MW 251.3247 Bethoxazin	C ₁₁ H ₉ NO ₂ S ₂	25mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Binapacryl				
CAS 485-31-4 DRE-C10590000 DRE-L10590000CY	MW 322.3132 Binapacryl Binapacryl 10 µg/mL in Cyclohexane	C ₁₅ H ₁₈ N ₂ O ₆	250mg 10ml	
Biphenyl				
CAS 92-52-4 DRE-C10630000	MW 154.2078 Biphenyl(‡)	C ₁₂ H ₁₀	250mg	
Biphenyl D10				
CAS 1486-01-7 DRE-C10630010	MW 164.2694 Biphenyl D10	C ₁₂ ² H ₁₀	100mg	
Bitertanol				
CAS 55179-31-2 DRE-C10660000 DRE-L10660000AL	MW 337.4155 Bitertanol Bitertanol 10 µg/mL in Acetonitrile	C ₂₀ H ₂₃ N ₃ O ₂	250mg 10ml	
Bithionol				
CAS 97-18-7 DRE-C10660500	MW 356.0518 Bithionol	C ₁₂ H ₆ Cl ₄ O ₂ S	250mg	
Bixafen				
CAS 581809-46-3 DRE-C10661480	MW 414.2086 Bixafen	C ₁₈ H ₁₂ Cl ₂ F ₃ N ₃ O	100mg	
Bixafen-desmethyl				
CAS 1655498-06-8 DRE-C10661486	MW 400.182 Bixafen-desmethyl	C ₁₇ H ₁₀ Cl ₂ F ₃ N ₃ O	10mg	
Boscalid (2-Chloro-N-(4'-chlorobiphenyl-2-yl)nicotinamide)				
CAS 188425-85-6 DRE-C10663000 DRE-L10663000AL DRE-L10663000CY	MW 343.2067 Boscalid(‡) Boscalid 10 µg/mL in Acetonitrile Boscalid 10 µg/mL in Cyclohexane	C ₁₈ H ₁₂ Cl ₂ N ₂ O	100mg 10ml 10ml	
Bromuconazole				
CAS 116255-48-2 DRE-C10802200 DRE-L10802200AL DRE-L10802200CY	MW 377.0639 Bromuconazole Bromuconazole 10 µg/mL in Acetonitrile Bromuconazole 10 µg/mL in Cyclohexane	C ₁₃ H ₁₂ BrCl ₂ N ₃ O	100mg 10ml 10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Bupirimate				
CAS 41483-43-6	MW 316.4197	C ₁₃ H ₂₄ N ₄ O ₃ S		
DRE-C10850000	Bupirimate		250mg	
DRE-L10850000CY	Bupirimate 10 µg/mL in Cyclohexane		10ml	
DRE-XA10850000AL	Bupirimate 100 µg/mL in Acetonitrile		1ml	
Captafol				
CAS 2425-06-1	MW 349.061	C ₁₀ H ₈ Cl ₄ NO ₂ S		
DRE-C10950000	Captafol		250mg	
DRE-L10950000IO	Captafol 10 µg/mL in Isooctane		10ml	
DRE-XA10950000CY	Captafol 100 µg/mL in Cyclohexane		1ml	
Captan				
CAS 133-06-2	MW 300.5893	C ₉ H ₈ Cl ₃ NO ₂ S		
DRE-C10960000	Captan(‡)		250mg	
DRE-L10960000CY	Captan 10 µg/mL in Cyclohexane		10ml	
DRE-XA10960000AL	Captan 100 µg/mL in Acetonitrile		1ml	
DRE-XA10960000CY	Captan 100 µg/mL in Cyclohexane		1ml	
Captan-4,4,5,6,7,7-D6				
CAS 1330190-00-5	MW 306.6263	C ₉ ² H ₆ H ₂ Cl ₃ NO ₂ S		
DRE-XA10960100AC	Captan D6 100 µg/mL in Acetone (*)		1.1ml	
Carbendazim				
CAS 10605-21-7	MW 191.1867	C ₈ H ₈ N ₂ O ₂		
DRE-C10990000	Carbendazim(‡)		250mg	
DRE-L10990000AL	Carbendazim 10 µg/mL in Acetonitrile (*)		10ml	
DRE-L10990000ME	Carbendazim 10 µg/mL in Methanol (*)		10ml	
DRE-GS09010053AL	Carbendazim 100 µg/mL in Acetonitrile w/ 0.1% Formic Acid(‡)		5x1ml	
DRE-GH09010098AL	Carbendazim 100 µg/mL in Acetonitrile(‡)		10x1ml	
Carbendazim D3 (methyl D3)				
CAS 1255507-88-0	MW 194.2051	C ₈ ² H ₅ H ₃ N ₂ O ₂		
DRE-C10990100	Carbendazim D3		10mg	
DRE-LA10990100AL	Carbendazim D3 (methyl D3) 10µg/mL in Acetonitrile (*)		1ml	
Carbendazim D4 (ring D4)				
CAS 291765-95-2	MW 195.2113	C ₈ ² H ₄ H ₅ N ₂ O ₂		
DRE-C10990200	Carbendazim D4 (ring)		10mg	
DRE-XA10990200AC	Carbendazim D4 (ring) 100 µg/mL in Acetone (*)		1.1ml	
Carboxin				
CAS 5234-68-4	MW 235.3021	C ₁₂ H ₁₃ NO ₂ S		
DRE-C11040000	Carboxin		250mg	
DRE-L11040000CY	Carboxin 10 µg/mL in Cyclohexane		10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Carpropamid				
CAS 104030-54-8 DRE-C11045900	MW 334.6685 Carpropamid(±)	C ₁₅ H ₁₈ Cl ₃ NO	100mg	
Chloraniformethan				
CAS 20856-57-9 DRE-L11127000AC	MW 336.4297 Chloraniformethan 10 µg/mL in Acetone	C ₈ H ₇ Cl ₃ N ₂ O	10ml	
4-Chloranil				
CAS 118-75-2 DRE-C11130000	MW 245.875 4-Chloranil	C ₆ Cl ₄ O ₂	250mg	
Chloroneb				
CAS 2675-77-6 DRE-C11450000 DRE-L11450000IO	MW 207.0539 Chloroneb Chloroneb 10 µg/mL in Isooctane	C ₈ H ₈ Cl ₂ O ₂	100mg 10ml	
Chlorothalonil				
CAS 1897-45-6 DRE-C11510000 DRE-L11510000AL DRE-L11510000CY DRE-XA11510000CY DRE-LA09010256CY	MW 265.911 Chlorothalonil(±) Chlorothalonil 10 µg/mL in Acetonitrile Chlorothalonil 10 µg/mL in Cyclohexane Chlorothalonil 100 µg/mL in Cyclohexane Chlorothalonil 10 µg/mL in Cyclohexane(±)	C ₆ Cl ₄ N ₂	250mg 10ml 10ml 1ml 10ml	
Chlorothalonil-4-hydroxy (4-Hydroxychlorothalonil)				
CAS 28343-61-5 DRE-LA11510400AL	MW 247.4653 Chlorothalonil-4-hydroxy 10 µg/mL in Acetonitrile (*)	C ₆ HCl ₃ N ₂ O	1ml	
Chlozolate				
CAS 84332-86-5 DRE-C11665000 DRE-L11665000CY DRE-XA11665000CY	MW 332.1361 Chlozolate(±) Chlozolate 10 µg/mL in Cyclohexane Chlozolate 100 µg/mL in Cyclohexane	C ₁₃ H ₁₁ Cl ₂ NO ₅	10mg 10ml 1ml	
Copper(II) Hydroxide				
CAS 20427-59-2 DRE-C11698000	MW 97.5607 Copper hydroxide	CuH ₂ O ₂	250mg	
Copper Oxychloride				
CAS 1332-40-7 DRE-C11700000	MW 213.567 Copper Oxychloride	Cl·2Cu·3HO	250mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Cyazofamid				
CAS 120116-88-3	MW 324.7859	C ₁₃ H ₁₃ ClN ₄ O ₂ S		
DRE-CA11816000	Cyazofamid(±)		100mg	
DRE-L11816000CY	Cyazofamid 10 µg/mL in Cyclohexane		10ml	
Cyflufenamid				
CAS 180409-60-3	MW 412.3532	C ₂₀ H ₁₇ F ₅ N ₂ O ₂		
DRE-C11843000	Cyflufenamid(±)		50mg	
DRE-L11843000CY	Cyflufenamid 10 µg/mL in Cyclohexane		10ml	
Cymoxanil				
CAS 57966-95-7	MW 198.1793	C ₇ H ₁₀ N ₄ O ₃		
DRE-C11880000	Cymoxanil(±)		100mg	
DRE-XA11880000MB	Cymoxanil 100 µg/mL in Methyl-tert. butyl ether		1ml	
Cyproconazole				
CAS 94361-06-5	MW 291.7759	C ₁₅ H ₁₆ ClN ₃ O		
DRE-C11908000	Cyproconazole(±)		100mg	
DRE-L11908000AL	Cyproconazole 10 µg/mL in Acetonitrile		10ml	
DRE-XA11908000CY	Cyproconazole 100 µg/mL in Cyclohexane		1ml	
Cyprodinil				
CAS 121552-61-2	MW 225.289	C ₁₄ H ₁₅ N ₃		
DRE-C11909000	Cyprodinil(±)		100mg	
DRE-L11909000AL	Cyprodinil 10 µg/mL in Acetonitrile		10ml	
DRE-L11909000IO	Cyprodinil 10 µg/mL in Isooctane		10ml	
DRE-XA11909000ME	Cyprodinil 100 µg/mL in Methanol		1ml	
Cyprofuram				
CAS 69581-33-5	MW 279.7189	C ₁₄ H ₁₄ ClNO ₃		
DRE-C11910000	Cyprofuram		250mg	
Dichlofluanid				
CAS 1085-98-9	MW 333.2302	C ₉ H ₁₁ Cl ₂ FN ₂ O ₂ S ₂		
DRE-C12300000	Dichlofluanid		250mg	
DRE-L12300000IO	Dichlofluanid 10 µg/mL in Isooctane		10ml	
DRE-XA12300000IO	Dichlofluanid 100 µg/mL in Isooctane		1ml	
Dichlone				
CAS 117-80-6	MW 227.0436	C ₁₀ H ₄ Cl ₂ O ₂		
DRE-C12310000	Dichlone		250mg	
Dichlorophen				
CAS 97-23-4	MW 269.1233	C ₁₃ H ₁₀ Cl ₂ O ₂		
DRE-C12440000	Dichlorophen		250mg	

Pesticide and metabolites: Fungicides

Product code	Description			
Diclobutrazol				
CAS 75736-33-3	MW 328.2369	C ₁₅ H ₁₉ Cl ₂ N ₃ O		
DRE-C12535000	Diclobutrazol		50mg	
DRE-L12535000AC	Diclobutrazol 10 µg/mL in Acetone		10ml	
DRE-XA09010180ME	Diclobutrazol 100 µg/mL in Methanol(‡)		1ml	
Diclomezine				
CAS 62865-36-5	MW 255.1	C ₁₁ H ₈ Cl ₂ N ₂ O		
DRE-C12545000	Diclomezine		100mg	
DRE-XA12545000AC	Diclomezine 100 µg/mL in Acetone		1ml	
Dicloran (2,6-Dichloro-4-nitroaniline)				
CAS 99-30-9	MW 207.0142	C ₆ H ₄ Cl ₂ N ₂ O ₂		
DRE-C12560000	Dicloran		250mg	
DRE-L12560000CY	Dicloran 10 µg/mL in Cyclohexane		10ml	
Diethofencarb				
CAS 87130-20-9	MW 267.3208	C ₁₄ H ₂₁ NO ₄		
DRE-C12603500	Diethofencarb(‡)		100mg	
DRE-L12603500AL	Diethofencarb 10 µg/mL in Acetonitrile		10ml	
DRE-L12603500CY	Diethofencarb 10 µg/mL in Cyclohexane		10ml	
Difenoconazole				
CAS 119446-68-3	MW 406.2626	C ₁₉ H ₁₇ Cl ₂ N ₃ O ₃		
DRE-C12609000	Difenoconazole(‡)		250mg	
DRE-L12609000AL	Difenoconazole 10 µg/mL in Acetonitrile		10ml	
DRE-L12609000CY	Difenoconazole 10 µg/mL in Cyclohexane		10ml	
DRE-XA12609000CY	Difenoconazole 100 µg/mL in Cyclohexane		1ml	
Difenoconazole D6 (1,1,2,3,3,3-propyl-D6)				
CAS n/a	MW 412.2996	C ₁₉ ² H ₆ H ₁₁ Cl ₂ N ₃ O ₃		
DRE-XA12609010AL	Difenoconazole D6 (1,1,2,3,3,3-propyl-D6) 100 µg/mL in Acetonitrile		1.1ml	
Difenoconazole-alcohol				
CAS 117018-19-6	MW 350.1993	C ₁₆ H ₁₃ Cl ₂ N ₃ O ₂		
DRE-C12609050	Difenoconazole-alcohol		10mg	
Diflumetorim				
CAS 130339-07-0	MW 327.7568	C ₁₅ H ₁₆ ClF ₂ N ₃ O		
DRE-C12631500	Diflumetorim		25mg	
Dimethachlon				
CAS 24096-53-5	MW 244.0741	C ₁₀ H ₇ Cl ₂ NO ₂		
DRE-C12669950	Dimethachlon		100mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Dimethirimol				
CAS 5221-53-4 DRE-C12690000	MW 209.2881 Dimethirimol(‡)	C ₁₁ H ₁₉ N ₃ O	100mg	
Dimethomorph				
CAS 110488-70-5 DRE-C12710000 DRE-L12710000AL DRE-L12710000CY	MW 387.8567 Dimethomorph(‡) Dimethomorph 10 µg/mL in Acetonitrile Dimethomorph 10 µg/mL in Cyclohexane	C ₂₁ H ₂₂ ClNO ₄	100mg 10ml 10ml	
Z-Dimethomorph				
CAS 113210-98-3 DRE-XA12710020AL	MW 387.8567 Z-Dimethomorph 100 µg/mL in Acetonitrile	C ₂₁ H ₂₂ ClNO ₄	1ml	
Dimoxystrobin				
CAS 149961-52-4 DRE-C12775000 DRE-L12775000AL	MW 326.3896 Dimoxystrobin Dimoxystrobin 10 µg/mL in Acetonitrile	C ₁₉ H ₂₂ N ₂ O ₃	100mg 10ml	
Diniconazole				
CAS 83657-24-3 DRE-C12777000 DRE-L12777000IO	MW 326.221 Diniconazole(‡) Diniconazole 10 µg/mL in Isooctane	C ₁₅ H ₁₇ Cl ₂ N ₃ O	100mg 10ml	
Dinocap				
CAS 39300-45-3 DRE-C12800000	MW 350.3664 Dinocap	C ₉ H ₅ NO ₄ ·C ₈ H ₁₇ ·NO ₂	250mg	
Diphenylamine				
CAS 122-39-4 DRE-C12890000 DRE-L12890000CY DRE-XA12890000AL DRE-GS09010101DI	MW 169.2224 Diphenylamine Diphenylamine 10 µg/mL in Cyclohexane Diphenylamine 100 µg/mL in Acetonitrile Diphenylamine 1000 µg/mL in Dichloromethane(‡)	C ₁₂ H ₁₁ N	250mg 10ml 1ml 5x1ml	
Ditalimfos				
CAS 5131-24-8 DRE-C13000000 DRE-L13000000CY	MW 299.2826 Ditalimfos Ditalimfos 10 µg/mL in Cyclohexane	C ₁₂ H ₁₄ NO ₄ PS	50mg 10ml	
Dithianon				
CAS 3347-22-6 DRE-C13010000 DRE-L13010000CY	MW 296.3238 Dithianon Dithianon 10 µg/mL in Cyclohexane	C ₁₄ H ₄ N ₂ O ₂ S ₂	250mg 10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Dodemorph				
CAS 1593-77-7 DRE-C13070000	MW 281.4766 Dodemorph(±)	C ₁₈ H ₃₅ NO	250mg	
trans-Dodemorph				
CAS 91269-48-6 DRE-C13070030	MW 281.4766 trans-Dodemorph	C ₁₈ H ₃₅ NO	10mg	
Dodine (Dodecylguanidinium Acetate)				
CAS 2439-10-3 DRE-C13080000	MW 287.4414 Dodine	C ₁₃ H ₂₉ N ₃ ·C ₂ H ₄ O ₂	250mg	
Drazoxolon				
CAS 5707-69-7 DRE-C13090000	MW 237.6424 Drazoxolon	C ₁₀ H ₈ ClN ₃ O ₂	100mg	
Edifenphos				
CAS 17109-49-8 DRE-CA13110000 DRE-L13110000CY DRE-XA09010189ME	MW 310.3715 Edifenphos(±) Edifenphos 10 µg/mL in Cyclohexane Edifenphos 100 µg/mL in Methanol(±)	C ₁₄ H ₁₅ O ₂ PS ₂	250mg 10ml 1ml	
Enilconazole (Imazalil)				
CAS 35554-44-0 DRE-C14280000 DRE-L14280000AL DRE-L14280000CY DRE-XA14280000AL	MW 297.1798 Imazalil(±) Imazalil 10 µg/mL in Acetonitrile Imazalil 10 µg/mL in Cyclohexane Imazalil 100 µg/mL in Acetonitrile	C ₁₄ H ₁₄ Cl ₂ N ₂ O	100mg 10ml 10ml 1ml	
Enilconazole-D5 (2-propenyl-D5) (Imazalil-D5)				
CAS 1398065-91-2 DRE-XA14280100AC	MW 302.2106 Imazalil D5 100 µg/mL in in Acetone(±)	C ₁₄ ² H ₅ H ₉ Cl ₂ N ₂ O	1.1ml	
Enoxastrobin				
CAS 238410-11-2 DRE-C13167500	MW 399.8674 Enoxastrobin	C ₂₂ H ₂₂ ClNO ₄	50mg	
Epichlorohydrin				
CAS 106-89-8 DRE-CA13175000 DRE-XA13175000CY	MW 92.5242 Epichlorohydrin Epichlorohydrin 100 µg/mL in Cyclohexane (*)	C ₃ H ₅ ClO	1ml 1ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Epoxiconazole				
CAS 133855-98-8 DRE-C13185000 DRE-XA13185000AL	MW 329.756 Epoxiconazole(±) Epoxiconazole 100 µg/mL in Acetonitrile	C ₁₇ H ₁₃ ClFN ₃ O	100mg 1ml	
Etaconazole				
CAS 60207-93-4 DRE-C13215000 DRE-L13215000CY	MW 328.1938 Etaconazole Etaconazole 10 µg/mL in Cyclohexane	C ₁₄ H ₁₅ Cl ₂ N ₃ O ₂	100mg 10ml	
Ethaboxam				
CAS 162650-77-3 DRE-C13217000 DRE-L13217000ME	MW 320.433 Ethaboxam(±) Ethaboxam 10 µg/mL in Methanol (*)	C ₁₄ H ₁₆ N ₄ OS ₂	50mg 10ml	
Ethirimol				
CAS 23947-60-6 DRE-C13280000 DRE-L13280000CY	MW 209.2881 Ethirimol Ethirimol 10 µg/mL in Cyclohexane	C ₁₁ H ₁₆ N ₃ O	100mg 10ml	
Ethoxyquin				
CAS 91-53-2 DRE-CA13310000	MW 217.3068 Ethoxyquin (*)	C ₁₄ H ₁₉ NO	250mg	
Ethoxyquin hydrochloride				
CAS 3659-01-6 DRE-C13310006	MW 253.7677 Ethoxyquin hydrochloride	C ₁₄ H ₁₉ NO·ClH	100mg	
Etridiazole				
CAS 2593-15-9 DRE-C13370000 DRE-L13370000CY	MW 247.53 Etridiazole(±) Etridiazole 10 µg/mL in Cyclohexane	C ₉ H ₉ Cl ₃ N ₂ OS	100mg 10ml	
Famoxadone				
CAS 131807-57-3 DRE-C13399000 DRE-L13399000CY	MW 374.3893 Famoxadone(±) Famoxadone 10 µg/mL in Cyclohexane	C ₂₂ H ₁₈ N ₂ O ₄	100mg 10ml	
Fenamidone				
CAS 161326-34-7 DRE-C13408000 DRE-L13408000AL DRE-L13408000CY DRE-XA13408000CY	MW 311.4014 Fenamidone(±) Fenamidone 10 µg/mL in Acetonitrile Fenamidone 10 µg/mL in Cyclohexane Fenamidone 100 µg/mL in Cyclohexane	C ₁₇ H ₁₇ N ₃ OS	100mg 10ml 10ml 1ml	









Pesticide and metabolites: Fungicides

Product code	Description			
Fenaminosulf				
CAS 140-56-7 DRE-C13410000	MW 251.2381 Fenaminosulf	$C_8H_{10}Na_3O_3S \cdot Na$	100mg	
Fenarimol				
CAS 60168-88-9 DRE-C13430000 DRE-L13430000CY DRE-L13430000EA DRE-XA13430000CY	MW 331.196 Fenarimol Fenarimol 10 µg/mL in Cyclohexane Fenarimol 10 µg/mL in Ethylacetate Fenarimol 100 µg/mL in Cyclohexane	$C_{17}H_{12}Cl_2N_2O$	100mg 10ml 10ml 1ml	
Fenbuconazole				
CAS 114369-43-6 DRE-C13448500 DRE-L13448500CY	MW 336.8181 Fenbuconazole(‡) Fenbuconazole 10 µg/mL in Cyclohexane	$C_{19}H_{17}ClN_4$	100mg 10ml	
Fenbuconazole (phenyl D5)				
CAS 1398066-06-2 DRE-XA13448510AC	MW 341.8489 Fenbuconazole (phenyl D5) 100 µg/mL in Acetone	$C_{19}^2H_5H_{12}ClN_4$	1.1ml	
Fenfuram				
CAS 24691-80-3 DRE-C13470000	MW 201.2212 Fenfuram	$C_{12}H_{11}NO_2$	100mg	
Fenhexamid				
CAS 126833-17-8 DRE-C13476000 DRE-L13476000CY	MW 302.1963 Fenhexamid(‡) Fenhexamid 10 µg/mL in Cyclohexane	$C_{14}H_{17}Cl_2NO_2$	100mg 10ml	
Fenoxanil				
CAS 115852-48-7 DRE-C13498000 DRE-L13498000AL	MW 329.2216 Fenoxanil Fenoxanil 10 µg/mL in Acetonitrile	$C_{15}H_{15}Cl_2N_2O_2$	100mg 10ml	
Fenpiclonil				
CAS 74738-17-3 DRE-C13525000	MW 237.0847 Fenpiclonil	$C_{11}H_6Cl_2N_2$	100mg	
Fenpropidin				
CAS 67306-00-7 DRE-C13537000 DRE-L13537000CY DRE-XA13537000CY	MW 273.4561 Fenpropidin Fenpropidin 10 µg/mL in Cyclohexane Fenpropidin 100 µg/mL in Cyclohexane	$C_{19}H_{31}N$	100mg 10ml 1ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Fenpropidin D10 (piperidine D10)				
CAS n/a	MW 283.5178	C ₁₉ ² H ₁₀ H ₂₁ N		
DRE-XA13537100CY	Fenpropidin D10 (piperidine D10) 100 µg/mL in Cyclohexane		1.1ml	
Fenpropimorph				
CAS 67564-91-4	MW 303.4821	C ₂₀ H ₃₃ NO		
DRE-CA13540000	Fenpropimorph(±)		250mg	
DRE-L13540000AL	Fenpropimorph 10 µg/mL in Acetonitrile		10ml	
DRE-L13540000CY	Fenpropimorph 10 µg/mL in Cyclohexane		10ml	
DRE-XA13540000CY	Fenpropimorph 100 µg/mL in Cyclohexane		1ml	
Fenpyrazamine				
CAS 473798-59-3	MW 331.4325	C ₁₇ H ₂₁ N ₃ O ₂ S		
DRE-C13544000	Fenpyrazamine(±)		100mg	
Fentin Chloride (Triphenyltin chloride)				
CAS 639-58-7	MW 385.4747	C ₁₈ H ₁₅ ClSn		
DRE-C13601000	Fentin-chloride(±)		250mg	
DRE-GA09010359ME	Triphenyltin Chloride 1000 µg/mL in Methanol(±)		1ml	
Fentin Hydroxide (Triphenyltin hydroxide)				
CAS 76-87-9	MW 367.029	C ₁₈ H ₁₅ OSn		
DRE-C13602000	Fentin-hydroxide		250mg	
Fentin Oxide				
CAS 1262-21-1	MW 716.0428	C ₃₆ H ₃₀ OSn ₂		
DRE-C13603000	Fentin-oxide		250mg	
Fentin-acetate				
CAS 900-95-8	MW 409.0657	C ₂₀ H ₁₈ O ₂ Sn		
DRE-CA13600000	Fentin-acetate		250mg	
Ferbam				
CAS 14484-64-1	MW 416.4943	C ₉ H ₁₈ FeN ₂ S ₆		
DRE-C13640000	Ferbam		250mg	
Ferimzone				
CAS 89269-64-7	MW 254.3302	C ₁₅ H ₁₈ N ₄		
DRE-C13642000	Ferimzone		100mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Fluazinam				
CAS 79622-59-6	MW 465.0917	C ₁₃ H ₄ Cl ₂ F ₆ N ₄ O ₄		
DRE-C13671500	Fluazinam(±)		100mg	
DRE-L13671500CY	Fluazinam 10 µg/mL in Cyclohexane		10ml	
Fludioxonil				
CAS 131341-86-1	MW 248.185	C ₁₂ H ₆ F ₂ N ₂ O ₂		
DRE-C13705000	Fludioxonil(±)		100mg	
DRE-GA13705000AL	Fludioxonil 100 µg/mL in Acetonitrile		1ml	
DRE-GS13705000AL	Fludioxonil 100 µg/mL in Acetonitrile		5x1ml	
DRE-L13705000AC	Fludioxonil 10 µg/mL in Acetone		10ml	
DRE-L13705000AL	Fludioxonil 10 µg/mL in Acetonitrile		10ml	
DRE-XA13705000AL	Fludioxonil 100 µg/mL in Acetonitrile		1ml	
Flumorph				
CAS 211867-47-9	MW 371.4021	C ₂₁ H ₂₂ FNO ₄		
DRE-C13726000	Flumorph		25mg	
DRE-XA13726000AL	Flumorph 100 µg/mL in Acetonitrile		1ml	
Fluopicolide				
CAS 239110-15-7	MW 383.5803	C ₁₄ H ₆ Cl ₃ F ₃ N ₂ O		
DRE-C13740000	Fluopicolide(±)		100mg	
DRE-XA13740000AL	Fluopicolide 100 µg/mL in Acetonitrile		1ml	
Fluopyram				
CAS 658066-35-4	MW 396.7148	C ₁₆ H ₁₁ ClF ₆ N ₂ O		
DRE-C13743000	Fluopyram(±)		50mg	
Fluoroimide				
CAS 41205-21-4	MW 260.0487	C ₁₀ H ₄ Cl ₂ FNO ₂		
DRE-C13793000	Fluoroimide		100mg	
Fluotrimazole				
CAS 31251-03-3	MW 379.3777	C ₂₂ H ₁₆ F ₃ N ₃		
DRE-C13800000	Fluotrimazole		100mg	
Fluoxastrobin				
CAS 361377-29-9	MW 458.8269	C ₂₁ H ₁₆ ClFN ₄ O ₅		
DRE-C13801000	Fluoxastrobin(±)		100mg	
DRE-L13801000AL	Fluoxastrobin 10 µg/mL in Acetonitrile		10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Fluquinconazole				
CAS 136426-54-5	MW 376.172	C ₁₆ H ₈ Cl ₂ FN ₃ O		
DRE-C13805000	Fluquinconazole(±)		100mg	
DRE-L13805000IO	Fluquinconazole 10 µg/mL in Isooctane		10ml	
Flusilazole				
CAS 85509-19-9	MW 315.3927	C ₁₆ H ₁₅ F ₂ N ₃ Si		
DRE-C13860000	Flusilazole(±)		100mg	
DRE-L13860000AL	Flusilazole 10 µg/mL in Acetonitrile		10ml	
DRE-L13860000EA	Flusilazole 10 µg/mL in Ethylacetate		10ml	
DRE-XA13860000EA	Flusilazole 100 µg/mL in Ethylacetate		1ml	
Flusulfamide				
CAS 106917-52-6	MW 415.1719	C ₁₃ H ₇ Cl ₂ F ₃ N ₂ O ₄ S		
DRE-C13861000	Flusulfamide		100mg	
Flutolanil				
CAS 66332-96-5	MW 323.3096	C ₁₇ H ₁₆ F ₃ NO ₂		
DRE-C13863500	Flutolanil		100mg	
DRE-L13863500CY	Flutolanil 10 µg/mL in Cyclohexane		10ml	
Flutriafol				
CAS 76674-21-0	MW 301.2907	C ₁₆ H ₁₃ F ₂ N ₃ O		
DRE-C13865000	Flutriafol(±)		100mg	
DRE-L13865000AL	Flutriafol 10 µg/mL in Acetonitrile		10ml	
Flutrimazole				
CAS 119006-77-8	MW 346.3726	C ₂₂ H ₁₆ F ₂ N ₂		
DRE-C13868000	Flutrimazole		250mg	
Fluxapyroxad				
CAS 907204-31-3	MW 381.2994	C ₁₈ H ₁₂ F ₅ N ₃ O		
DRE-C13875000	Fluxapyroxad(±)		100mg	
Folpet				
CAS 133-07-3	MW 296.5576	C ₉ H ₄ Cl ₃ NO ₂ S		
DRE-C13890000	Folpet(±)		250mg	
DRE-L13890000CY	Folpet 10 µg/mL in Cyclohexane		10ml	
DRE-XA13890000CY	Folpet 100 µg/mL in Cyclohexane		1ml	
Folpet D4				
CAS 1327204-12-5	MW 300.5822	C ₉ ² H ₄ Cl ₃ NO ₂ S		
DRE-C13890100	Folpet D4(±)		10mg	
DRE-XA13890100AC	Folpet D4 100 µg/mL in Acetone(±)		1.1ml	









Pesticide and metabolites: Fungicides

Product code	Description			
Fosetyl-aluminium D15				
CAS n/a	MW 369.197	3C ₂ H ₅ HO ₃ P·Al		
DRE-CA13940010	Fosetyl-aluminium D15 (*)		10mg	
Fosetyl Aluminium Salt				
CAS 39148-24-8	MW 354.1045	3C ₂ H ₅ O ₃ P·Al		
DRE-CA13940000	Fosetyl-aluminium(‡)		250mg	
Fuberidazole				
CAS 3878-19-1	MW 184.194	C ₁₁ H ₈ N ₂ O		
DRE-C13950000	Fuberidazole(‡)		250mg	
DRE-L13950000ME	Fuberidazole 10 µg/mL in Methanol		10ml	
Furalaxyl				
CAS 57646-30-7	MW 301.3371	C ₁₇ H ₁₉ NO ₄		
DRE-C13960000	Furalaxyl		100mg	
DRE-L13960000CY	Furalaxyl 10 µg/mL in Cyclohexane		10ml	
Furalaxyl (free acid)				
CAS 118597-19-6	MW 287.3105	C ₁₆ H ₁₇ NO ₄		
DRE-C13960010	Furalaxyl (free acid)		50mg	
Furametpyr				
CAS 123572-88-3	MW 333.8126	C ₁₇ H ₂₀ ClN ₃ O ₂		
DRE-C13964000	Furametpyr		100mg	
Furfural (Furane-2-carbaldehyde, 2-Furaldehyde)				
CAS 98-01-1	MW 96.0841	C ₅ H ₄ O ₂		
DRE-C13972100	Furfural (2-Furfurol, 2-Furaldehyde)(‡)		250mg	
Furmecyclox				
CAS 60568-05-0	MW 251.3214	C ₁₄ H ₂₁ NO ₃		
DRE-C13980000	Furmecyclox		100mg	
DRE-L13980000CY	Furmecyclox 10 µg/mL in Cyclohexane		10ml	
Hexachloroacetone				
CAS 116-16-5	MW 264.7495	C ₃ Cl ₆ O		
DRE-C14155000	Hexachloroacetone		500mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Hexachlorobenzene				
CAS 118-74-1	MW 284.7822	C ₆ Cl ₆		
DRE-C14160000	Hexachlorobenzene(‡)		250mg	
DRE-L14160000AL	Hexachlorobenzene 10 µg/mL in Acetonitrile		10ml	
DRE-L14160000CY	Hexachlorobenzene 10 µg/mL in Cyclohexane		10ml	
Hexachlorobenzene 13C6				
CAS 93952-14-8	MW 290.7381	¹³ C ₆ Cl ₆		
DRE-C14160100	Hexachlorobenzene 13C6(‡)		10mg	
DRE-XA14160100AC	Hexachlorobenzene 13C6 100 µg/mL in Acetone(‡)		1ml	
Hexachlorophene				
CAS 70-30-4	MW 406.9035	C ₁₂ H ₆ Cl ₆ O ₂		
DRE-C14180000	Hexachlorophen		250mg	
DRE-GA09010397DI	Hexachlorophene Solution, 5000 µg/mL in Dichloromethane(‡)		1ml	
DRE-GS09010397DI-SS	Hexachlorophene Solution, Second Source, 5000 µg/mL in Dichloromethane (‡)		5x1ml	
Hexaconazole				
CAS 79983-71-4	MW 314.2103	C ₁₄ H ₁₇ Cl ₂ N ₃ O		
DRE-C14190000	Hexaconazole(‡)		100mg	
DRE-L14190000CY	Hexaconazole 10 µg/mL in Cyclohexane		10ml	
DRE-XA14190000AL	Hexaconazole 100 µg/mL in Acetonitrile		1ml	
Hymexazol				
CAS 10004-44-1	MW 99.088	C ₄ H ₅ NO ₂		
DRE-C14270000	Hymexazol(‡)		25mg	
Imazalil sulfate				
CAS 58594-72-2	MW 395.2582	C ₁₄ H ₁₄ Cl ₂ N ₂ O·H ₂ O ₄ S		
DRE-C14280500	Imazalil sulfate		100mg	
Imibenconazole				
CAS 86598-92-7	MW 411.7359	C ₁₇ H ₁₃ Cl ₃ N ₃ S		
DRE-C14283600	Imibenconazole(‡)		100mg	
DRE-L14283600AC	Imibenconazole 10 µg/mL in Acetone		10ml	
Imibenconazole-desbenzyl (N-(2,4-Dichlorophenyl)-1H-1,2,4-triazole-1-acetamide)				
CAS 154221-27-9	MW 271.1027	C ₁₀ H ₆ Cl ₂ N ₄ O		
DRE-C14283620	Imibenconazole-desbenzyl(‡)		10mg	
Iminoctadine triacetate				
CAS 57520-17-9	MW 535.7209	C ₁₈ H ₄₁ N ₇ ·3C ₂ H ₄ O ₂		
DRE-C14284900	Iminoctadine triacetate		10mg	
DRE-L14284900ME	Iminoctadine triacetate 10 µg/mL in Methanol (*)		10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Iminoctadine tris(albesilate)				
CAS 169202-06-6 DRE-C14284950	MW 1335.047 Iminoctadine tris(albesilate)	$C_{18}H_{41}N_7 \cdot 3C_{18}H_{30}O_3S$	10mg	
Iodocarb (IPBC)				
CAS 55406-53-6 DRE-C14335000	MW 281.0909 Iodocarb (IPBC)	$C_8H_{12}INO_2$	100mg	
Ipconazole				
CAS 125225-28-7 DRE-C14365000 DRE-XA09010204AL	MW 333.8557 Ipconazole(‡) Ipconazole 100 µg/mL in Acetonitrile(‡)	$C_{18}H_{24}ClN_3O$	100mg 1ml	
Iprobenfos				
CAS 26087-47-8 DRE-CA14368000 DRE-L14368000AL DRE-L14368000CY DRE-XA09010205AL	MW 288.3428 Iprobenfos(‡) Iprobenfos 10 µg/mL in Acetonitrile Iprobenfos 10 µg/mL in Cyclohexane Iprobenfos 100 µg/mL in Acetonitrile(‡)	$C_{13}H_{21}O_3PS$	100mg 10ml 10ml 1ml	
Iprodione				
CAS 36734-19-7 DRE-C14370000 DRE-L14370000AL DRE-L14370000CY DRE-XA14370000AL DRE-XA14370000CY DRE-LA09010266CY	MW 330.1666 Iprodione(‡) Iprodione 10 µg/mL in Acetonitrile Iprodione 10 µg/mL in Cyclohexane Iprodione 100 µg/mL in Acetonitrile Iprodione 100 µg/mL in Cyclohexane Iprodione 10 µg/mL in Cyclohexane(‡)	$C_{13}H_{13}Cl_2N_3O_3$	100mg 10ml 10ml 1ml 1ml 10ml	
Iprodione D7 (isopropyl D7)				
CAS n/a DRE-XA14370012AL	MW 337.2098 Iprodione D7 100 µg/mL in Acetonitrile	$C_{13}^2H_{17}H_6Cl_2N_3O_3$	1ml	
Iprodione des-(N-isopropylcarboxamid)				
CAS 27387-87-7 DRE-C14370040	MW 245.0621 Iprodione des-(N-isopropylcarboxamid)	$C_9H_6Cl_2N_2O_2$	10mg	
Iprodione isomer 1 (Isoiprodione)				
CAS 63637-89-8 DRE-C14370500	MW 330.1666 Iprodione isomer 1	$C_{13}H_{13}Cl_2N_3O_3$	100mg	








Pesticide and metabolites: Fungicides

Product code	Description			
Iprovalicarb				
CAS 140923-17-7	MW 320.4265	C ₁₆ H ₂₈ N ₂ O ₃		
DRE-C14371000	Iprovalicarb(‡)		100mg	
DRE-L14371000AL	Iprovalicarb 10 µg/mL in Acetonitrile		10ml	
Isoprothiolane				
CAS 50512-35-1	MW 290.3989	C ₁₂ H ₁₈ O ₄ S ₂		
DRE-C14467500	Isoprothiolane(‡)		100mg	
DRE-L14467500AL	Isoprothiolane 10 µg/mL in Acetonitrile		10ml	
DRE-L14467500CY	Isoprothiolane 10 µg/mL in Cyclohexane		10ml	
Isopyrazam				
CAS 881685-58-1	MW 359.4129	C ₂₀ H ₂₃ F ₂ N ₃ O		
DRE-C14473000	Isopyrazam		10mg	
Kasugamycin Hydrochloride				
CAS 19408-46-9	MW 415.8239	C ₁₄ H ₂₃ N ₃ O ₉ ·ClH		
DRE-C14515000	Kasugamycin hydrochloride		250mg	
Kresoxim-methyl				
CAS 143390-89-0	MW 313.3478	C ₁₈ H ₁₉ NO ₄		
DRE-C14570000	Kresoxim-methyl(‡)		100mg	
DRE-L14570000CY	Kresoxim-methyl 10 µg/mL in Cyclohexane		10ml	
Mancozeb				
CAS 8018-01-7	MW 541.0747	C ₄ H ₆ MnN ₂ S ₄ ·C ₄ H ₆ N ₂ S ₄ Zn		
DRE-C14740000	Mancozeb		250mg	
Mandipropamid				
CAS 374726-62-2	MW 411.8781	C ₂₃ H ₂₂ ClNO ₄		
DRE-C14745000	Mandipropamid(‡)		100mg	
DRE-L14745000CY	Mandipropamid 10 µg/mL in Cyclohexane		10ml	
DRE-XA14745000AL	Mandipropamid 100 µg/mL in Acetonitrile		1ml	
Maneb				
CAS 12427-38-2	MW 265.3019	C ₄ H ₆ MnN ₂ S ₄		
DRE-CA14750000	Maneb		250mg	
Mebenil				
CAS 7055-03-0	MW 211.2591	C ₁₄ H ₁₃ NO		
DRE-C14798500	Mebenil		100mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Mepanipyrim				
CAS 110235-47-7	MW 223.2731	C ₁₄ H ₁₃ N ₃		
DRE-C14867000	Mepanipyrim(±)		50mg	
DRE-L14867000CY	Mepanipyrim 10 µg/mL in Cyclohexane		10ml	
Mepanipyrim-2-hydroxypropyl				
CAS 204571-52-8	MW 243.3043	C ₁₄ H ₁₇ N ₃ O		
DRE-C14867050	Mepanipyrim-2-hydroxypropyl		10mg	
Mepronil				
CAS 55814-41-0	MW 269.3383	C ₁₇ H ₁₉ N ₂ O ₂		
DRE-C14890000	Mepronil		100mg	
DRE-L14890000CY	Mepronil 10 µg/mL in Cyclohexane		10ml	
Meptyldinocap				
CAS 131-72-6	MW 364.393	C ₁₈ H ₂₄ N ₂ O ₆		
DRE-C14895000	Meptyldinocap(±)		100mg	
Meptyldinocap-phenol				
CAS 3687-22-7	MW 296.319	C ₁₄ H ₂₀ N ₂ O ₅		
DRE-C14895050	Meptyldinocap-phenol		100mg	
Metalaxyl				
CAS 57837-19-1	MW 279.3315	C ₁₅ H ₂₁ NO ₄		
DRE-C14920000	Metalaxyl(±)		100mg	
DRE-L14920000AL	Metalaxyl 10 µg/mL in Acetonitrile		10ml	
DRE-L14920000CY	Metalaxyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA14920000TO	Metalaxyl 100 µg/mL in Toluene		1ml	
Metalaxyl-M				
CAS 70630-17-0	MW 279.3315	C ₁₅ H ₂₁ NO ₄		
DRE-C14920500	Metalaxyl-M(±)		100mg	
DRE-L14920500CY	Metalaxyl-M 10 µg/mL in Cyclohexane		10ml	
Metconazole				
CAS 125116-23-6	MW 319.8291	C ₁₇ H ₂₂ ClN ₃ O		
DRE-C14955000	Metconazole(±)		100mg	
DRE-L14955000CY	Metconazole 10 µg/mL in Cyclohexane		10ml	
DRE-XA14955000AL	Metconazole 100 µg/mL in Acetonitrile		1ml	
Methasulfocarb				
CAS 66952-49-6	MW 261.3179	C ₉ H ₁₁ NO ₄ S ₂		
DRE-C14998000	Methasulfocarb		25mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Methfuroxam				
CAS 28730-17-8 DRE-C15010000	MW 229.2744 Methfuroxam	C ₁₄ H ₁₅ NO ₂	100mg	
Metiram				
CAS 9006-42-2 DRE-C15150000	MW n/a Metiram		250mg	No Structure
(E)-Metominostrobin				
CAS 133408-50-1 DRE-C15175500 DRE-L15175500AL	MW 284.3098 (E)-Metominostrobin (E)-Metominostrobin 10 µg/mL in Acetonitrile	C ₁₆ H ₁₆ N ₂ O ₃	10mg 10ml	
(Z)-Metominostrobin				
CAS 133408-51-2 DRE-LA15175510AL	MW 284.3098 (Z)-Metominostrobin 10 µg/mL in Acetonitrile	C ₁₆ H ₁₆ N ₂ O ₃	1ml	
Metrafenone				
CAS 220899-03-6 DRE-C15190000 DRE-L15190000AL	MW 409.271 Metrafenone(‡) Metrafenone 10 µg/mL in Acetonitrile	C ₁₉ H ₂₁ BrO ₅	100mg 10ml	
Z-Mevinphos				
CAS 338-45-4 DRE-LA15222000CY	MW 224.1483 Z-Mevinphos (cis) 10 µg/mL in Cyclohexane	C ₇ H ₁₃ O ₆ P	1ml	
Myclobutanil				
CAS 88671-89-0 DRE-C15390000 DRE-L15390000CY	MW 288.7753 Myclobutanil(‡) Myclobutanil 10 µg/mL in Cyclohexane	C ₁₅ H ₁₇ ClN ₄	100mg 10ml	
Nabam				
CAS 142-59-6 DRE-C15400000	MW 256.3434 Nabam	C ₄ H ₆ N ₂ S ₄ ·2Na	100mg	
Naphthenic acids, copper salts (Copper (II) naphthenate)				
CAS 1338-02-9 DRE-C11699000	MW n/a Copper (II) naphthenate		250mg	No Structure










Pesticide and metabolites: Fungicides

Product code	Description			
Nitrothal-isopropyl				
CAS 10552-74-6 DRE-C15610000	MW 295.2879 Nitrothal-isopropyl	C ₁₄ H ₁₇ NO ₆	100mg	
Nuarimol				
CAS 63284-71-9 DRE-L15660000CY	MW 314.7414 Nuarimol 10 µg/mL in Cyclohexane	C ₁₇ H ₁₂ ClFN ₂ O	10ml	
Octhilinone				
CAS 26530-20-1 DRE-C15711500 DRE-L15711500AL	MW 213.3397 Octhilinone Octhilinone 10 µg/mL in Acetonitrile	C ₁₁ H ₁₉ NOS	100mg 10ml	
Ofurace				
CAS 58810-48-3 DRE-C15718000 DRE-L15718000CY	MW 281.7347 Ofurace(±) Ofurace 10 µg/mL in Cyclohexane	C ₁₄ H ₁₆ ClNO ₃	100mg 10ml	
Orysastrobin				
CAS 248593-16-0 DRE-C15749000 DRE-XA15749000AL	MW 391.4216 Orysastrobin Orysastrobin 100 µg/mL in Acetonitrile	C ₁₈ H ₂₅ N ₅ O ₅	100mg 1ml	
Oxadixyl				
CAS 77732-09-3 DRE-C15770000 DRE-L15770000AL DRE-L15770000CY DRE-XA15770000AC	MW 278.3037 Oxadixyl(±) Oxadixyl 10 µg/mL in Acetonitrile Oxadixyl 10 µg/mL in Cyclohexane Oxadixyl 100 µg/mL in Acetone	C ₁₄ H ₁₈ N ₂ O ₄	100mg 10ml 10ml 1ml	
Oxine Copper				
CAS 10380-28-6 DRE-C15786000	MW 351.8461 Oxine-Copper	C ₁₈ H ₁₂ CuN ₂ O ₂	250mg	
Oxolinic Acid				
CAS 14698-29-4 DRE-C15788000	MW 261.2301 Oxolinic acid(±)	C ₁₃ H ₁₁ NO ₅	250mg	
Oxpoconazole fumarate				
CAS 174212-12-5 DRE-C15789000	MW 477.9379 Oxpoconazole fumarate	C ₁₉ H ₂₄ ClN ₃ O ₂ ·C ₄ H ₄ O ₄	100mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Oxycarboxin				
CAS 5259-88-1	MW 267.3009	C ₁₂ H ₁₃ NO ₄ S		
DRE-C15790000	Oxycarboxin		100mg	
DRE-XA09010217AL	Oxycarboxin 100 µg/mL in Acetonitrile(‡)		1ml	
Oxyquinoline (8-Hydroxyquinoline)				
CAS 148-24-3	MW 145.158	C ₉ H ₇ NO		
DRE-C14249500	8-Hydroxyquinoline		250mg	
Oxyquinoline Sulfate (8-Hydroxyquinoline sulfate)				
CAS 134-31-6	MW 388.3944	2C ₉ H ₇ NO·H ₂ O ₄ S		
DRE-C14250000	8-Hydroxyquinoline sulfate		250mg	
Oxytetracycline Hydrochloride				
CAS 2058-46-0	MW 496.8949	C ₂₂ H ₂₄ N ₂ O ₉ ·ClH		
DRE-C15820000	Oxytetracycline hydrochloride(‡)		250mg	
Penconazole				
CAS 66246-88-6	MW 284.1843	C ₁₃ H ₁₅ Cl ₂ N ₃		
DRE-C15910000	Penconazole		100mg	
DRE-L15910000CY	Penconazole 10 µg/mL in Cyclohexane		10ml	
DRE-XA15910000AL	Penconazole 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010220AL	Penconazol 100 µg/mL in Acetonitrile(‡)		1ml	
Pencycuron				
CAS 66063-05-6	MW 328.8358	C ₁₉ H ₂₁ ClN ₂ O		
DRE-C15920000	Pencycuron		100mg	
Penthiopyrad				
CAS 183675-82-3	MW 359.4097	C ₁₆ H ₂₀ F ₃ N ₃ OS		
DRE-C15981760	Penthiopyrad		25mg	
Phenylmercury Chloride				
CAS 100-56-1	MW 313.1469	C ₆ H ₅ ClHg		
DRE-C16065000	Phenyl-mercury-chloride		250mg	
2-Phenylphenol				
CAS 90-43-7	MW 170.2072	C ₁₂ H ₁₀ O		
DRE-C16070000	2-Phenylphenol		250mg	
DRE-L16070000AL	2-Phenylphenol 10 µg/mL in Acetonitrile		10ml	
DRE-L16070000CY	2-Phenylphenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA16070000AC	2-Phenylphenol 100 µg/mL in Acetone		1ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Phosphonic Acid				
CAS 13598-36-2 DRE-C16144000	MW 81.9958 Phosphonic acid	H ₃ O ₃ P	250mg	
Phthalide (4,5,6,7-Tetrachlorophthalide)				
CAS 27355-22-2 DRE-C16185000 DRE-LA16185000AC	MW 271.9123 Phthalide(‡) Phthalide 10 µg/mL in Acetone (*)	C ₈ H ₂ Cl ₄ O ₂	10mg 1ml	
Picoxystrobin				
CAS 117428-22-5 DRE-C16206000	MW 367.3191 Picoxystrobin	C ₁₈ H ₁₆ F ₃ NO ₄	100mg	
Piperalin				
CAS 3478-94-2 DRE-C16220000 DRE-XA09010223AL	MW 330.2494 Piperalin Piperalin 100 µg/mL in Acetonitrile(‡)	C ₁₆ H ₂₁ Cl ₂ NO ₂	100mg 1ml	
Polycarbamate				
CAS 64440-88-6 DRE-C16282000	MW 581.6147 Polycarbamate	C ₁₀ H ₁₈ N ₄ S ₆ Zn ₂	100mg	
Polyoxin B				
CAS 19396-06-6 DRE-E16283500	MW 507.4061 Polyoxin B	C ₁₇ H ₂₅ N ₉ O ₁₃	100mg	
Probenazole				
CAS 27605-76-1 DRE-C16289000	MW 223.2484 Probenazole	C ₁₀ H ₉ NO ₃ S	100mg	
Prochloraz				
CAS 67747-09-5 DRE-C16290000 DRE-L16290000CY DRE-XA16290000AL	MW 376.6654 Prochloraz(‡) Prochloraz 10 µg/mL in Cyclohexane Prochloraz 100 µg/mL in Acetonitrile	C ₁₅ H ₁₆ Cl ₃ N ₃ O ₂	250mg 10ml 1ml	
Prochloraz D7 (propyl D7)				
CAS n/a DRE-XA16290010AL	MW 383.7086 Prochloraz D7 (propyl D7) 100 µg/mL in Acetonitrile	C ₁₅ ² H ₁₇ H ₉ Cl ₃ N ₃ O ₂	1.1ml	









Pesticide and metabolites: Fungicides

Product code	Description			
Prochloraz desimidazole-amino				
CAS 139520-94-8 DRE-XA16290100AL	MW 325.6187 Prochloraz desimidazole-amino 100 µg/mL in Acetonitrile (*)	C ₁₂ H ₁₅ Cl ₃ N ₂ O ₂	1ml	
Prochloraz desimidazole-formylamino				
CAS 139542-32-8 DRE-C16290150	MW 353.6288 Prochloraz desimidazole-formylamino	C ₁₃ H ₁₅ Cl ₃ N ₂ O ₃	10mg	
Procymidone				
CAS 32809-16-8 DRE-C16310000 DRE-L16310000CY DRE-XA16310000IO	MW 284.1379 Procymidone(‡) Procymidone 10 µg/mL in Cyclohexane Procymidone 100 µg/mL in Isooctane	C ₁₃ H ₁₁ Cl ₂ N ₂ O ₂	250mg 10ml 1ml	
Propamocarb free base				
CAS 24579-73-5 DRE-C16390000 DRE-L16390000AL	MW 188.2673 Propamocarb free base(‡) Propamocarb free base 10 µg/mL in Acetonitrile	C ₉ H ₂₀ N ₂ O ₂	100mg 10ml	
Propamocarb free base D7 (O-propyl D7)				
CAS 1398065-89-8 DRE-C16390100 DRE-XA16390100AC	MW 195.3104 Propamocarb D7 (propyl D7) Propamocarb free base D7 100 µg/mL in Acetone	C ₉ ² H ₇ ² H ₁₃ N ₂ O ₂	10mg 1.1ml	
Propamocarb hydrochloride				
CAS 25606-41-1 DRE-C16400000 DRE-L16400000ME DRE-XA16400000ME	MW 224.7282 Propamocarb hydrochloride Propamocarb hydrochloride 10 µg/mL in Methanol Propamocarb hydrochloride 100 µg/mL in Methanol	C ₉ H ₂₀ N ₂ O ₂ ·ClH	100mg 10ml 1ml	
Propiconazole				
CAS 60207-90-1 DRE-C16480000 DRE-L16480000AL DRE-L16480000CY DRE-XA16480000AL DRE-XA16480000CY	MW 342.2204 Propiconazole(‡) Propiconazole 10 µg/mL in Acetonitrile Propiconazole 10 µg/mL in Cyclohexane Propiconazole 100 µg/mL in Acetonitrile Propiconazole 100 µg/mL in Cyclohexane	C ₁₅ H ₁₇ Cl ₂ N ₃ O ₂	250mg 10ml 10ml 1ml 1ml	
Propiconazole D5 (2,2,3,3,3-propyl-D5)				
CAS 1246818-14-3 DRE-XA16480100AC	MW 347.2512 Propiconazole D5 (2,2,3,3,3-propyl-D5) 100 µg/mL in Acetone	C ₁₅ ² H ₅ ² H ₁₂ Cl ₂ N ₃ O ₂	1.1ml	
Propineb				
CAS 12071-83-9 DRE-C16490000	MW 289.7994 Propineb	C ₈ H ₈ N ₂ S ₄ Zn	250mg	









Pesticide and metabolites: Fungicides

Product code	Description			
Proquinazid				
CAS 189278-12-4 DRE-C16542000 DRE-L16542000CY	MW 372.2014 Proquinazid Proquinazid 10 µg/mL in Cyclohexane	C ₁₄ H ₁₇ IN ₂ O ₂	50mg 10ml	
Prothiocarb hydrochloride				
CAS 19622-19-6 DRE-C16550000	MW 226.7673 Prothiocarb hydrochloride	C ₈ H ₁₈ N ₂ OS·ClH	10mg	
Prothioconazole				
CAS 178928-70-6 DRE-C16555000 DRE-L16555000AL	MW 344.2594 Prothioconazole(‡) Prothioconazole 10 µg/mL in Acetonitrile (*)	C ₁₄ H ₁₅ Cl ₂ N ₃ OS	100mg 10ml	
Prothioconazole-desthio				
CAS 120983-64-4 DRE-C16555500 DRE-L16555500AL	MW 312.1944 Prothioconazole-desthio(‡) Prothioconazole-desthio 10 µg/mL in Acetonitrile	C ₁₄ H ₁₅ Cl ₂ N ₃ O	25mg 10ml	
Pyracarbolid				
CAS 24691-76-7 DRE-C16590000 DRE-XA09010230AL	MW 217.2637 Pyracarbolid Pyracarbolid 100 µg/mL in Acetonitrile(‡)	C ₁₃ H ₁₅ NO ₂	100mg 1ml	
Pyraclostrobin				
CAS 175013-18-0 DRE-C16595000 DRE-L16595000AL	MW 387.8169 Pyraclostrobin(‡) Pyraclostrobin 10 µg/mL in Acetonitrile	C ₁₉ H ₁₈ ClN ₃ O ₄	100mg 10ml	
Pyrazophos				
CAS 13457-18-6 DRE-C16610000 DRE-L16610000IO	MW 373.3645 Pyrazophos Pyrazophos 10 µg/mL in Isooctane	C ₁₄ H ₂₀ N ₃ O ₅ PS	250mg 10ml	
Pyribencarb				
CAS 799247-52-2 DRE-C16623000	MW 361.8227 Pyribencarb	C ₁₈ H ₂₀ ClN ₃ O ₃	25mg	
Pyrifenox				
CAS 88283-41-4 DRE-C16655000 DRE-L16655000IO	MW 295.1639 Pyrifenox Pyrifenox 10 µg/mL in Isooctane	C ₁₄ H ₁₂ Cl ₂ N ₂ O	100mg 10ml	










Pesticide and metabolites: Fungicides

Product code	Description			
E-Pyriphenox				
CAS 83227-22-9	MW 295.1639	C ₁₄ H ₁₂ Cl ₂ N ₂ O		
DRE-LA16655100IO	E-Pyriphenox 10 µg/mL in Isooctane		1ml	
Pyrimethanil				
CAS 53112-28-0	MW 199.2517	C ₁₂ H ₁₃ N ₃		
DRE-C16658500	Pyrimethanil(‡)		100mg	
DRE-L16658500CY	Pyrimethanil 10 µg/mL in Cyclohexane		10ml	
DRE-XA16658500ME	Pyrimethanil 100 µg/mL in Methanol		1ml	
Pyroquilon				
CAS 57369-32-1	MW 173.2111	C ₁₁ H ₁₁ NO		
DRE-C16665000	Pyroquilon		250mg	
DRE-L16665000CY	Pyroquilon 10 µg/mL in Cyclohexane		10ml	
Quinoxifen				
CAS 124495-18-7	MW 308.1345	C ₁₅ H ₈ Cl ₂ FNO		
DRE-C16715000	Quinoxifen(‡)		100mg	
DRE-L16715000AL	Quinoxifen 10 µg/mL in Acetonitrile		10ml	
DRE-L16715000CY	Quinoxifen 10 µg/mL in Cyclohexane		10ml	
Quintozene				
CAS 82-68-8	MW 295.3347	C ₆ Cl ₅ NO ₂		
DRE-C16730000	Quintozene(‡)		250mg	
DRE-GA16730000AL	Quintozene 100 µg/mL in Acetonitrile		1ml	
DRE-GS16730000AL	Quintozene 100 µg/mL in Acetonitrile		5x1ml	
DRE-L16730000CY	Quintozene 10 µg/mL in Cyclohexane		10ml	
DRE-XA16730000AL	Quintozene 100 µg/mL in Acetonitrile		1ml	
DRE-XA16730000CY	Quintozene 100 µg/mL in Cyclohexane		1ml	
Rabenzazole				
CAS 40341-04-6	MW 212.2505	C ₁₂ H ₁₂ N ₄		
DRE-C16800000	Rabenzazole		10mg	
Silthiofam				
CAS 175217-20-6	MW 267.4624	C ₁₃ H ₂₁ NOSSi		
DRE-C16947000	Silthiofam		25mg	
DRE-L16947000CY	Silthiofam 10 µg/mL in Cyclohexane		10ml	
Simeconazole				
CAS 149508-90-7	MW 293.412	C ₁₄ H ₂₀ FN ₃ OSi		
DRE-C16957000	Simeconazole(‡)		100mg	









Pesticide and metabolites: Fungicides

Product code	Description			
Spiroxamine				
CAS 118134-30-8	MW 297.476	C ₁₈ H ₃₅ NO ₂		
DRE-C16973000	Spiroxamine		100mg	
DRE-L16973000AL	Spiroxamine 10 µg/mL in Acetonitrile		10ml	
DRE-L16973000CY	Spiroxamine 10 µg/mL in Cyclohexane		10ml	
Sulfur				
CAS 7704-34-9	MW 32.065	S		
DRE-C17025000	Sulfur		250mg	
DRE-L17025000IO	Sulfur 10 µg/mL in Isooctane		10ml	
TCMTB (Busan; 2-(Thiocyanomethylthio)benzothiazole)				
CAS 21564-17-0	MW 238.3523	C ₉ H ₆ N ₂ S ₃		
DRE-C17178500	TCMTB (Busan)(‡)		10mg	
DRE-XA09010239AL	TCMTB 100 µg/mL in Acetonitrile(‡)		1ml	
Tebuconazole				
CAS 107534-96-3	MW 307.8184	C ₁₆ H ₂₂ ClN ₃ O		
DRE-C17178700	Tebuconazole(‡)		250mg	
DRE-L17178700AL	Tebuconazole 10 µg/mL in Acetonitrile		10ml	
DRE-L17178700IO	Tebuconazole 10 µg/mL in Isooctane		10ml	
DRE-XA17178700AL	Tebuconazole 100 µg/mL in Acetonitrile		1ml	
DRE-XA17178700IO	Tebuconazole 100 µg/mL in Isooctane		1ml	
DRE-GA09010371ME	Tebuconazole 100 µg/mL in Methanol(‡)		1ml	
Tebuconazole D6				
CAS n/a	MW 313.8554	C ₁₆ ² H ₆ H ₁₆ ClN ₃ O		
DRE-XA17178710AC	Tebuconazole D6 100 µg/mL in Acetone(‡)		1ml	
Tecloftalam				
CAS 76280-91-6	MW 447.9124	C ₁₄ H ₅ Cl ₆ NO ₃		
DRE-C17195000	Tecloftalam		10mg	
Tecnazene (2,3,5,6-Tetrachloronitrobenzene)				
CAS 117-18-0	MW 260.8896	C ₆ HCl ₄ NO ₂		
DRE-C17200000	Tecnazene		250mg	
DRE-L17200000CY	Tecnazene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17200000IO	Tecnazene 100 µg/mL in Isooctane		1ml	
DRE-LA09010252CY	Tecnazene 10 µg/mL in Cyclohexane(‡)		10ml	
Tetraconazole				
CAS 112281-77-3	MW 372.1456	C ₁₃ H ₁₁ Cl ₂ F ₄ N ₃ O		
DRE-C17395000	Tetraconazole(‡)		100mg	
DRE-L17395000AL	Tetraconazole 10 µg/mL in Acetonitrile		10ml	
DRE-L17395000IO	Tetraconazole 10 µg/mL in Isooctane		10ml	
DRE-XA17395000ME	Tetraconazole 100 µg/mL in Methanol		1ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Thiabendazole-5-hydroxy (5-Hydroxytiabendazole)				
CAS 948-71-0	MW 217.2471	C ₁₀ H ₇ N ₃ OS		
DRE-C17450500	Thiabendazole-5-hydroxy(‡)		10mg	
DRE-LA17450500ME	Thiabendazole-5-hydroxy 10 µg/mL in Methanol		1ml	
Thifluzamide				
CAS 130000-40-7	MW 528.0624	C ₁₃ H ₆ Br ₂ F ₆ N ₂ O ₂ S		
DRE-C17468000	Thifluzamide		100mg	
DRE-L17468000AL	Thifluzamide 10 µg/mL in Acetonitrile		10ml	
DRE-L17468000CY	Thifluzamide 10 µg/mL in Cyclohexane		10ml	
Thiophanate (Thiophanate-ethyl)				
CAS 23564-06-9	MW 370.4471	C ₁₄ H ₁₄ N ₄ O ₄ S ₂		
DRE-C17540000	Thiophanate-ethyl		100mg	
Thiophanate-methyl				
CAS 23564-05-8	MW 342.394	C ₁₂ H ₁₄ N ₄ O ₄ S ₂		
DRE-C17545000	Thiophanate-methyl(‡)		250mg	
Thioquinox				
CAS 93-75-4	MW 236.3365	C ₉ H ₄ N ₂ S ₃		
DRE-C17560000	Thioquinox		10mg	
Thiram				
CAS 137-26-8	MW 240.4329	C ₆ H ₁₂ N ₂ S ₄		
DRE-C17570000	Thiram		250mg	
DRE-L17570000CY	Thiram 10 µg/mL in Cyclohexane		10ml	
Thiram D12				
CAS 69193-86-8	MW 252.5068	C ₆ ² H ₁₂ N ₂ S ₄		
DRE-X17570100AL	Thiram D12 100 µg/mL in Acetonitrile (*)		10ml	
Thiabendazole (Thiabendazole)				
CAS 148-79-8	MW 201.2477	C ₁₀ H ₇ N ₃ S		
DRE-C17450000	Thiabendazole(‡)		250mg	
DRE-L17450000AL	Thiabendazole 10 µg/mL in Acetonitrile		10ml	
DRE-L17450000ME	Thiabendazole 10 µg/mL in Methanol		10ml	
DRE-XA17450000ME	Thiabendazole 100 µg/mL in Methanol		1ml	
Thiabendazole D6 (Thiabendazole D6)				
CAS 1262551-89-2	MW 207.2847	C ₁₀ ² H ₆ HN ₃ S		
DRE-C17450100	Thiabendazole NH D6		10mg	
DRE-XA17450100AC	Thiabendazole NH D6 100 µg/mL in Acetone		1.1ml	










Pesticide and metabolites: Fungicides

Product code	Description			
Tiadinil				
CAS 223580-51-6	MW 267.7346	C ₁₁ H ₁₀ ClN ₃ OS		
DRE-C17575400	Tiadinil		100mg	
DRE-XA17575400AL	Tiadinil 100 µg/mL in Acetonitrile		1ml	
Tolclofos-methyl				
CAS 57018-04-9	MW 301.1266	C ₉ H ₁₁ Cl ₂ O ₃ PS		
DRE-C17590000	Tolclofos-methyl(‡)		250mg	
DRE-L17590000CY	Tolclofos-methyl 10 µg/mL in Cyclohexane		10ml	
Tolnifanide				
CAS 304911-98-6	MW 354.8086	C ₁₅ H ₁₅ ClN ₂ O ₄ S		
DRE-C17591800	Tolnifanide		50mg	
Tolyfluanid				
CAS 731-27-1	MW 347.2568	C ₁₀ H ₁₃ Cl ₂ FN ₂ O ₂ S ₂		
DRE-C17600000	Tolyfluanid(‡)		250mg	
DRE-L17600000CY	Tolyfluanid 10 µg/mL in Cyclohexane		10ml	
DRE-XA17600000CY	Tolyfluanid 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010246AL	Tolyfluanid 100 µg/mL in Acetonitrile(‡)		1ml	
Triadimefon				
CAS 43121-43-3	MW 293.7487	C ₁₄ H ₁₆ ClN ₃ O ₂		
DRE-C17610000	Triadimefon(‡)		250mg	
DRE-L17610000AL	Triadimefon 10 µg/mL in Acetonitrile		10ml	
DRE-L17610000CY	Triadimefon 10 µg/mL in Cyclohexane		10ml	
DRE-LA09010264CY	Triadimefon 10 µg/mL in Cyclohexane(‡)		10ml	
Triadimenol				
CAS 55219-65-3	MW 295.7646	C ₁₄ H ₁₈ ClN ₃ O ₂		
DRE-C17620000	Triadimenol(‡)		250mg	
DRE-L17620000AL	Triadimenol 10 µg/mL in Acetonitrile		10ml	
DRE-L17620000IQ	Triadimenol 10 µg/mL in Isooctane		10ml	
DRE-XA17620000CY	Triadimenol 100 µg/mL in Cyclohexane		1ml	
erythro-Triadimenol (isomer A)				
CAS 70585-35-2	MW 295.7646	C ₁₄ H ₁₈ ClN ₃ O ₂		
DRE-C17621000	Triadimenol isomer A		100mg	
threo-Triadimenol (isomer B)				
CAS 70585-37-4	MW 295.7646	C ₁₄ H ₁₈ ClN ₃ O ₂		
DRE-XA17621100AL	Triadimenol isomer B 100 µg/mL in Acetonitrile		1ml	



Pesticide and metabolites: Fungicides

Product code	Description			
Triamphos				
CAS 1031-47-6	MW 294.2926	C ₁₂ H ₁₉ N ₆ OP		
DRE-L17640000CY	Triamphos 10 µg/mL in Cyclohexane		10ml	
DRE-XA17640000AL	Triamphos 100 µg/mL in Acetonitrile		1ml	
Triazoxide				
CAS 72459-58-6	MW 247.6405	C ₁₀ H ₆ ClN ₅ O		
DRE-C17660000	Triazoxide		100mg	
Trichlamide (N-(1-Butoxy-2,2,2-trichloroethyl)salicylamide)				
CAS 70193-21-4	MW 340.63	C ₁₃ H ₁₆ Cl ₃ NO ₃		
DRE-C17669300	Trichlamide		10mg	
DRE-L17669300CY	Trichlamide 10 µg/mL in Cyclohexane		10ml	
Triclopyricarb				
CAS 902760-40-1	MW 391.6337	C ₁₅ H ₁₃ Cl ₃ N ₂ O ₄		
DRE-C17801000	Triclopyricarb		25mg	
Tricyclazole				
CAS 41814-78-2	MW 189.237	C ₉ H ₇ N ₃ S		
DRE-C17810000	Tricyclazole(‡)		100mg	
DRE-L17810000CY	Tricyclazole 10 µg/mL in Cyclohexane		10ml	
Tridemorph				
CAS 81412-43-3	MW 143.2267	C ₇ H ₁₅ NO(CH ₂) _n		
DRE-C17820000	Tridemorph		100mg	
DRE-L17820000IO	Tridemorph 10 µg/mL in Isooctane		10ml	
Trifloxystrobin				
CAS 141517-21-7	MW 408.3711	C ₂₀ H ₁₉ F ₃ N ₂ O ₄		
DRE-C17842000	Trifloxystrobin(‡)		100mg	
DRE-L17842000CY	Trifloxystrobin 10 µg/mL in Cyclohexane		10ml	
Triflumizole				
CAS 68694-11-1	MW 345.7473	C ₁₅ H ₁₅ ClF ₃ N ₃ O		
DRE-C17844000	Triflumizole(‡)		100mg	
DRE-L17844000IO	Triflumizole 10 µg/mL in Isooctane		10ml	
Triflumizole 13C3 (O-propyl 13C3)				
CAS n/a	MW 348.7253	¹³ C ₃ C ₁₂ H ₁₅ ClF ₃ N ₃ O		
DRE-XA17844010AC	Triflumizole 13C3 100 µg/mL in Acetone		1.1ml	

Pesticide and metabolites: Fungicides

Product code	Description			
Triflumizole-amino				
CAS 131549-75-2	MW 294.7006	C ₁₂ H ₁₄ ClF ₃ N ₂ O		
DRE-C17844030	Triflumizole-amino		10mg	
DRE-L17844030AL	Triflumizole-amino 10 µg/mL in Acetonitrile		10ml	
Triforine				
CAS 26644-46-2	MW 434.9618	C ₁₀ H ₁₄ Cl ₆ N ₄ O ₂		
DRE-C17860000	Triforine		250mg	
Triforine D8 (piperazine D8)				
CAS 948595-12-8	MW 443.0111	C ₁₀ ² H ₈ H ₆ Cl ₆ N ₄ O ₂		
DRE-C17860100	Triforine D8 (piperazine D8)		10mg	
Triticonazole				
CAS 131983-72-7	MW 317.8132	C ₁₇ H ₂₀ ClN ₃ O		
DRE-C17894600	Triticonazole(±)		100mg	
DRE-L17894600IO	Triticonazole 10 µg/mL in Isooctane		10ml	
Urbacide				
CAS 2445-07-0	MW 330.389	C ₇ H ₁₅ AsN ₂ S ₄		
DRE-C17897300	Urbacid		100mg	
Validamycin				
CAS 37248-47-8	MW 497.4908	C ₂₀ H ₃₅ NO ₁₃		
DRE-C17899900	Validamycin (technical)		250mg	
Valifenalate				
CAS 283159-90-0	MW 398.8811	C ₁₉ H ₂₇ ClN ₂ O ₅		
DRE-C17899960	Valifenalate		20mg	
Vinclozolin				
CAS 50471-44-8	MW 286.1108	C ₁₂ H ₉ Cl ₂ NO ₃		
DRE-C17920000	Vinclozolin		250mg	
DRE-L17920000IO	Vinclozolin 10 µg/mL in Isooctane		10ml	
DRE-XA17920000CY	Vinclozolin 100 µg/mL in Cyclohexane		1ml	
Zineb				
CAS 12122-67-7	MW 275.7728	C ₄ H ₆ N ₂ S ₄ Zn		
DRE-C17950000	Zineb		250mg	










Pesticide and metabolites: Fungicides

Product code	Description			
Ziram				
CAS 137-30-4 DRE-C17970000	MW 305.8419 Ziram	C ₆ H ₁₂ N ₂ S ₄ Zn	250mg	
Zoxamide				
CAS 156052-68-5 DRE-C17980000	MW 336.6413 Zoxamide(‡)	C ₁₄ H ₁₆ Cl ₃ NO ₂	100mg	










PESTICIDES AND METABOLITES: HERBICIDES



Pesticide and metabolites: Herbicides

Product code	Description			
Acetamidomethylphosphonic Acid (Aminomethylphosphonic acid N-acetyl)				
CAS 57637-97-5 DRE-C10205150	MW 153.0737 Aminomethyl phosphonicacid N-acetyl(‡)	C ₃ H ₈ NO ₄ P	10mg	
Acetochlor				
CAS 34256-82-1 DRE-C10018000 DRE-L10018000IO DRE-XA10018000CY	MW 269.7671 Acetochlor(‡) Acetochlor 10 µg/mL in Isooctane Acetochlor 100 µg/mL in Cyclohexane	C ₁₄ H ₂₀ ClNO ₂	100mg 10ml 1ml	
Acetochlor D11				
CAS 1189897-44-6 DRE-XA10018100AC	MW 280.8349 Acetochlor D11 100 µg/mL in Acetone	C ₁₄ ² H ₁₁ H ₉ ClNO ₂	1.1ml	
Acetochlor Ethane Sulfonic Acid Sodium Salt				
CAS 947601-84-5 DRE-CA10018210	MW 337.3671 Acetochlor ethane sulfonic acid (ESA) sodium salt	C ₁₄ H ₂₀ NO ₅ S·Na	10mg	
Acetochlor Oxalamic Acid (OA)				
CAS 194992-44-4 DRE-CA10018400	MW 265.305 Acetochlor oxalamic acid (OA)	C ₁₄ H ₁₉ NO ₄	10mg	
Acetochlor-2-hydroxy				
CAS 60090-47-3 DRE-XA10018250AL	MW 251.3214 Acetochlor-2-hydroxy 100 µg/mL in Acetonitrile	C ₁₄ H ₂₁ NO ₃	1ml	
Acifluorfen				
CAS 50594-66-6 DRE-C10030000	MW 361.6573 Acifluorfen	C ₁₄ H ₇ ClF ₃ NO ₅	50mg	
Acifluorfen Sodium				
CAS 62476-59-9 DRE-C10030500	MW 383.6391 Acifluorfen sodium	C ₁₄ H ₆ ClF ₃ NO ₅ ·Na	250mg	
Aclonifen				
CAS 74070-46-5 DRE-C10042000 DRE-L10042000AL DRE-XA10042000AL	MW 264.6645 Aclonifen Aclonifen 10 µg/mL in Acetonitrile Aclonifen 100 µg/mL in Acetonitrile	C ₁₂ H ₉ ClN ₂ O ₃	250mg 10ml 1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Alachlor				
CAS 15972-60-8	MW 269.7671	C ₁₄ H ₂₀ ClNO ₂		
DRE-C10060000	Alachlor(±)		250mg	
DRE-L10060000AC	Alachlor 10 µg/mL in Acetone		10ml	
DRE-L10060000AL	Alachlor 10 µg/mL in Acetonitrile		10ml	
DRE-XA10060000ME	Alachlor 100 µg/mL in Methanol		1ml	
Alachlor D13 (2,6-diethylphenyl D13)				
CAS 1015856-63-9	MW 282.8472	C ₁₄ ² H ₁₃ H ₇ ClNO ₂		
DRE-C10060100	Alachlor D13 (2,6-diethylphenyl D13)		10mg	
DRE-XA10060100AC	Alachlor D13 100 µg/mL in Acetone(±)		1.1ml	
Alachlor D3 (methoxy D3)				
CAS n/a	MW 272.7856	C ₁₄ ² H ₃ H ₁₇ ClNO ₂		
DRE-XA10060001AC	Alachlor D3 (methoxy D3) 100 µg/mL in Acetone		1.1ml	
Alachlor ethane sulfonic acid (ESA) sodium salt				
CAS 140939-15-7	MW 337.3671	C ₁₄ H ₂₀ NO ₅ S·Na		
DRE-CA10060210	Alachlor ethane sulfonic acid (ESA) sodium salt		10mg	
Alachlor oxalamic acid (OA)				
CAS 171262-17-2	MW 265.305	C ₁₄ H ₁₉ NO ₄		
DRE-CA10060400	Alachlor oxalamic acid (OA)		10mg	
Alachlor-2-hydroxy				
CAS 56681-55-1	MW 285.7665	C ₁₄ H ₂₀ ClNO ₃		
DRE-XA10060250AL	Alachlor-2-hydroxy 100 µg/mL in Acetonitrile		1ml	
Allidochlor				
CAS 93-71-0	MW 173.64	C ₈ H ₁₂ ClNO		
DRE-C10110000	Allidochlor		100mg	
DRE-XA09010147ME	Allidochlor 100 µg/mL in Methanol(±)		1ml	
Alloxydim-sodium				
CAS 55635-13-7	MW 345.3659	C ₁₇ H ₂₄ NO ₅ ·Na		
DRE-C10120000	Alloxydim-sodium		250mg	
Ametryn				
CAS 834-12-8	MW 227.3298	C ₉ H ₁₇ N ₅ S		
DRE-C10150000	Ametryn		250mg	
DRE-L10150000AL	Ametryn 10 µg/mL in Acetonitrile		10ml	
DRE-XA10150000AL	Ametryn 100 µg/mL in Acetonitrile		1ml	









Pesticide and metabolites: Herbicides

Product code	Description				
Ametryn-N-nitroso					
CAS n/a	MW 256.3279	C ₉ H ₁₆ N ₆ OS			
DRE-C10150200	Ametryn-N-nitroso		50mg		
Amicarbazone					
CAS 129909-90-6	MW 241.2902	C ₁₀ H ₁₉ N ₅ O ₂			
DRE-C10155000	Amicarbazone		10mg		
Amidosulfuron					
CAS 120923-37-7	MW 369.3747	C ₉ H ₁₅ N ₅ O ₇ S ₂			
DRE-C10162000	Amidosulfuron(‡)		100mg		
Amidosulfuron D6 (dimethoxy D6)					
CAS n/a	MW 375.4117	C ₉ ² H ₆ ² H ₉ N ₅ O ₇ S ₂			
DRE-C10162100	Amidosulfuron D6		10mg		
2-Amino-4-methoxy-6-methyl-1,3,5-triazine					
CAS 1668-54-8	MW 140.1432	C ₅ H ₈ N ₄ O			
DRE-C10204000	2-Amino-4-methoxy-6-methyl-1,3,5-triazine		100mg		
Aminomethylphosphonic Acid (AMPA)					
CAS 1066-51-9	MW 111.037	CH ₆ NO ₃ P			
DRE-C10205000	Aminomethyl phosphonic acid (AMPA)(‡)		100mg		
DRE-L10205000WA	Aminomethyl phosphonic acid (AMPA) 10 µg/mL in Water		10ml		
DRE-XA10205000WA	Aminomethyl phosphonic acid (AMPA) 100 µg/mL in Water		1ml		
DRE-XA09010249WA	1-Aminomethylphosphonic Acid 100 µg/mL in Water(‡)		1ml		
(Aminomethyl) phosphonic Acid 13C 15N (AMPA)					
CAS n/a	MW 113.0231	¹³ CH ₆ ¹⁵ NO ₃ P			
DRE-XA10205100WA	Aminomethylphosphonic acid 13C 15N (AMPA) 100 µg/mL in Water(‡)		1ml		
Aminomethylphosphonic Acid-FMOC (AMPA-FMOC)					
CAS 195306-88-8	MW 333.2757	C ₁₆ H ₁₆ NO ₅ P			
DRE-C10205800	AMPA-FMOC(‡)		10mg		
Aminopyralid					
CAS 150114-71-9	MW 207.0142	C ₆ H ₄ Cl ₂ N ₂ O ₂			
DRE-C10218000	Aminopyralid		100mg		










Pesticide and metabolites: Herbicides

Product code	Description			
Amiprofos-methyl				
CAS 36001-88-4 DRE-C10228500	MW 304.3024 Amiprofos-methyl	C ₁₁ H ₁₇ N ₂ O ₄ PS	50mg	
Amitrole				
CAS 61-82-5 DRE-C10240000 DRE-L10240000ME	MW 84.08 Amitrole Amitrole 10 µg/mL in Methanol	C ₂ H ₄ N ₄	250mg 10ml	
Amitrole (1-15N)				
CAS 367498-28-0 DRE-XA10240100ME	MW 85.0734 Amitrole (1-15N) 100 µg/mL in Methanol	C ₂ H ₄ ¹⁵ NN ₃	1.1ml	
Amitrole 1-15N 5-13C				
CAS n/a DRE-XA10240110AL	MW 86.066 Amitrole 1-15N 5-13C 100µg/mL in Acetonitrile (*) (‡)	¹³ CCH ₄ ¹⁵ NN ₃	1.1ml	
Anilofos				
CAS 64249-01-0 DRE-CA10264000 DRE-L10264000IO	MW 367.8516 Anilofos Anilofos 10 µg/mL in Isooctane	C ₁₃ H ₁₉ ClNO ₃ PS ₂	100mg 10ml	
Arylpyrazol (ACE5599; 2-Chloro-N-[1-[2,6-dichloro-4-(difluoromethyl)phenyl]-4-nitro-1H-pyrazol-5-yl]propanamide)				
CAS 121750-17-2 DRE-C10300600	MW 413.5914 Arylpyrazol	C ₁₃ H ₉ Cl ₃ F ₂ N ₄ O ₃	100mg	
Asulam				
CAS 3337-71-1 DRE-C10310000 DRE-L10310000AL	MW 230.241 Asulam Asulam 10 µg/mL in Acetonitrile	C ₈ H ₁₀ N ₂ O ₄ S	250mg 10ml	
Atraton				
CAS 1610-17-9 DRE-C10320000 DRE-L10320000AL	MW 211.2642 Atraton Atraton 10 µg/mL in Acetonitrile	C ₉ H ₁₇ N ₅ O	250mg 10ml	
Atrazine				
CAS 1912-24-9 DRE-C10330000 DRE-CR10330000 DRE-L10330000AL DRE-L10330000CY DRE-XA10330000AL DRE-XA10330000TO	MW 215.6833 Atrazine(‡) Atrazine Atrazine 10 µg/mL in Acetonitrile Atrazine 10 µg/mL in Cyclohexane Atrazine 100 µg/mL in Acetonitrile Atrazine 100 µg/mL in Toluene	C ₆ H ₁₄ ClN ₅	250mg 50mg 10ml 10ml 1ml 1ml	










Pesticide and metabolites: Herbicides

Product code	Description		
DRE-GS09010036ME	Atrazine 200 µg/mL in Methanol(±)		4x1ml
DRE-YS09010007DI	Atrazine 2000 µg/mL in Dichloromethane(±)		5x1ml
DRE-YS09010011AC	Atrazine 1000 µg/mL in Acetone(±)		5x1ml
Atrazine 13C3 (ring 13C3)			
CAS 1443685-80-0	MW 218.6612	$^{13}\text{C}_3\text{H}_5\text{H}_{14}\text{ClN}_5$	
DRE-XA10330200AC	Atrazine 13C3 (ring 13C3) 100 µg/mL in Acetone		1.1ml
			
Atrazine D5 (ethyl-D5)			
CAS 163165-75-1	MW 220.7141	$\text{C}_8\text{H}_5\text{H}_9\text{ClN}_5$	
DRE-C10330100	Atrazine D5(±)		10mg
DRE-XA10330100AC	Atrazine D5 100 µg/mL in Acetone(±)		1ml
DRE-YA10330100AL	Atrazine D5 (ethylamino D5) 1000 µg/mL in Acetonitrile		1.1ml
			
Atrazine-desethyl			
CAS 6190-65-4	MW 187.6301	$\text{C}_6\text{H}_{10}\text{ClN}_5$	
DRE-C10331000	Atrazine-desethyl		250mg
DRE-L10331000AL	Atrazine-desethyl 10 µg/mL in Acetonitrile		10ml
DRE-XA10331000AL	Atrazine-desethyl 100 µg/mL in Acetonitrile		1ml
DRE-XA10331000CY	Atrazine-desethyl 100 µg/mL in Cyclohexane		1ml
			
Atrazine-desethyl D6 (dimethyl D6)			
CAS 1216649-31-8	MW 193.6671	$\text{C}_6\text{H}_6\text{H}_4\text{ClN}_5$	
DRE-XA10331100AC	Atrazine-desethyl D6 100 µg/mL in Acetone		1.1ml
			
Atrazine-desethyl-desisopropyl			
CAS 3397-62-4	MW 145.5504	$\text{C}_6\text{H}_4\text{ClN}_5$	
DRE-C10331500	Atrazine-desethyl-desisopropyl		100mg
			
Atrazine-desethyl-desisopropyl-2-hydroxy (Ammelide)			
CAS 645-92-1	MW 127.1047	$\text{C}_3\text{H}_5\text{N}_5\text{O}$	
DRE-C10331700	Atrazine-desethyl-desisopropyl-2-hydroxy (Ammelide)		250mg
			
Atrazine-desethyl-2-hydroxy			
CAS 19988-24-0	MW 169.1844	$\text{C}_6\text{H}_{11}\text{N}_5\text{O}$	
DRE-C10331300	Atrazine-desethyl-2-hydroxy		100mg
DRE-XA10331300AL	Atrazine-desethyl-2-hydroxy 100 µg/mL in Acetonitrile		1ml
			
Atrazine-desisopropyl (6-Chloro-N-ethyl-1,3,5-triazine-2,4-diamine)			
CAS 1007-28-9	MW 173.6035	$\text{C}_8\text{H}_8\text{ClN}_5$	
DRE-C10332000	Atrazine-desisopropyl(±)		250mg
DRE-L10332000AL	Atrazine-desisopropyl 10 µg/mL in Acetonitrile		10ml
DRE-XA10332000AL	Atrazine-desisopropyl 100 µg/mL in Acetonitrile		1ml
DRE-XA09010149ME	Atrazine Desisopropyl 100 µg/mL in Methanol(±)		1ml
			










Pesticide and metabolites: Herbicides

Product code	Description			
Atrazine-desisopropyl D5 (ethylamino D5)				
CAS 1189961-78-1	MW 178.6343	C ₈ ² H ₉ H ₃ ClN ₅		
DRE-C10332100	Atrazine-desisopropyl D5 (ethylamino D5)		10mg	
DRE-XA10332100AC	Atrazine-desisopropyl D5 (ethylamino D5) 100 µg/mL in Acetone		1.1ml	
Atrazine-desisopropyl-2-hydroxy				
CAS 7313-54-4	MW 155.1579	C ₈ H ₉ N ₅ O		
DRE-C10332300	Atrazine-desisopropyl-2-hydroxy		100mg	
DRE-XA10332300ME	Atrazine-desisopropyl-2-hydroxy 100 µg/mL in Methanol		1ml	
Atrazine-2-hydroxy				
CAS 2163-68-0	MW 197.2376	C ₈ H ₁₅ N ₅ O		
DRE-C10333000	Atrazine-2-hydroxy		100mg	
DRE-L10333000ME	Atrazine-2-hydroxy 10 µg/mL in Methanol		10ml	
DRE-XA10333000ME	Atrazine-2-hydroxy 100 µg/mL in Methanol		1ml	
Atrazine-2-hydroxy D5 (ethyl D5)				
CAS 1276197-25-1	MW 202.2684	C ₈ ² H ₉ H ₁₀ N ₅ O		
DRE-XA10333100ME	Atrazine-2-hydroxy D5 100 µg/mL in Methanol		1.1ml	
Atrazine-mercapturate				
CAS 138722-96-0	MW 342.4172	C ₁₃ H ₂₂ N ₅ O ₃ S		
DRE-XA10333200AL	Atrazine-mercapturate 100 µg/mL in Acetonitrile		1.1ml	
Azimsulfuron				
CAS 120162-55-2	MW 424.3951	C ₁₃ H ₁₆ N ₁₀ O ₅ S		
DRE-C10355000	Azimsulfuron(‡)		10mg	
Aziprotryne				
CAS 4658-28-0	MW 225.2741	C ₇ H ₁₁ N ₇ S		
DRE-C10380000	Aziprotryne		250mg	
Barban				
CAS 101-27-9	MW 258.1007	C ₁₁ H ₉ Cl ₂ NO ₂		
DRE-C10420000	Barban		100mg	
Beflubutamid				
CAS 113614-08-7	MW 355.3267	C ₁₈ H ₁₇ F ₄ NO ₂		
DRE-C10430000	Beflubutamid		100mg	
DRE-XA09010153AL	Beflubutamid 100 µg/mL in Acetonitrile(‡)		1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Benazolin				
CAS 3813-05-6	MW 243.6668	C ₉ H ₆ ClNO ₃ S		
DRE-C10450000	Benazolin		100mg	
DRE-XA10450000ME	Benazolin 100 µg/mL in Methanol		1ml	
Benazolin Ethyl Ester				
CAS 25059-80-7	MW 271.72	C ₁₁ H ₁₀ ClNO ₃ S		
DRE-C10450500	Benazolin-ethyl ester		100mg	
Benfluralin				
CAS 1861-40-1	MW 335.279	C ₁₃ H ₁₆ F ₃ N ₃ O ₄		
DRE-C10470000	Benfluralin(‡)		250mg	
Benfuresate				
CAS 68505-69-1	MW 256.318	C ₁₂ H ₁₆ O ₄ S		
DRE-C10476000	Benfuresate(‡)		100mg	
Bensulfuron-methyl				
CAS 83055-99-6	MW 410.4017	C ₁₆ H ₁₈ NaO ₇ S		
DRE-C10498000	Bensulfuron-methyl		100mg	
Bensulide				
CAS 741-58-2	MW 397.5134	C ₁₄ H ₂₄ NO ₄ PS ₃		
DRE-C10500000	Bensulide(‡)		250mg	
DRE-L10500000AL	Bensulide 10 µg/mL in Acetonitrile		10ml	
DRE-XA09010155ME	Bensulide 100 µg/mL in Methanol(‡)		1ml	
Bentazone				
CAS 25057-89-0	MW 240.2789	C ₁₀ H ₁₂ N ₂ O ₃ S		
DRE-C10510000	Bentazone		250mg	
DRE-L10510000AL	Bentazone 10 µg/mL in Acetonitrile		10ml	
DRE-XA10510000AL	Bentazone 100 µg/mL in Acetonitrile		1ml	
Bentazone (isopropyl-1,1,1,3,3,3) D6				
CAS n/a	MW 246.3159	C ₁₀ ² H ₆ HeN ₂ O ₃ S		
DRE-C10510100	Bentazone (isopropyl-1,1,1,3,3,3) D6		10mg	
DRE-XA10510100AC	Bentazone (isopropyl-1,1,1,3,3,3) D6 100 µg/mL in Acetone (*)		1.1ml	
Bentazone-8-hydroxy				
CAS 60374-43-8	MW 258.2942	C ₁₀ H ₁₄ N ₂ O ₄ S		
DRE-C10512000	Bentazone-8-hydroxy		5mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Bentazone-methyl (N-Methylbentazone)				
CAS 61592-45-8 DRE-C10512500 DRE-L10512500AL	MW 254.3055 Bentazone-methyl Bentazone-methyl 10 µg/mL in Acetonitrile	C ₁₁ H ₁₄ N ₂ O ₃ S	50mg 10ml	
Bentranil				
CAS 1022-46-4 DRE-C10520000	MW 223.2268 Bentranil	C ₁₄ H ₉ NO ₂	10mg	
Benzobicyclon				
CAS 156963-66-5 DRE-C10536800	MW 446.9669 Benzobicyclon	C ₂₂ H ₁₉ ClO ₄ S ₂	10mg	
Benzofenap				
CAS 82692-44-2 DRE-C10537400 DRE-L10537400CY DRE-XA10537400CY	MW 431.3118 Benzofenap Benzofenap 10 µg/mL in Cyclohexane Benzofenap 100 µg/mL in Cyclohexane	C ₂₂ H ₂₀ Cl ₂ N ₂ O ₃	10mg 10ml 1ml	
Benzoylprop Ethyl Ester (Benzoylprop-ethyl)				
CAS 22212-55-1 DRE-C10550000	MW 366.2385 Benzoylprop-ethyl	C ₁₈ H ₁₇ Cl ₂ NO ₃	100mg	
Benzthiazuron				
CAS 1929-88-0 DRE-C10560000	MW 207.2523 Benzthiazuron	C ₉ H ₉ N ₃ OS	100mg	
Bifenox				
CAS 42576-02-3 DRE-C10580000 DRE-L10580000AL DRE-L10580000CY	MW 342.131 Bifenox Bifenox 10 µg/mL in Acetonitrile Bifenox 10 µg/mL in Cyclohexane	C ₁₄ H ₉ Cl ₂ NO ₅	250mg 10ml 10ml	
Bifenox (free acid)				
CAS 53774-07-5 DRE-C10580500	MW 328.1044 Bifenox (free acid)	C ₁₃ H ₇ Cl ₂ NO ₅	25mg	
Bilanafos sodium				
CAS 71048-99-2 DRE-C10587200	MW 345.2645 Bilanafos sodium	C ₁₁ H ₂₁ N ₃ O ₆ P·Na	10mg	










Pesticide and metabolites: Herbicides

Product code	Description			
N,N'-Bis-(3,4-dichlorophenyl)urea				
CAS 4300-43-0 DRE-C10651800	MW 350.0274 N,N'-Bis-(3,4-dichlorophenyl)urea	C ₁₃ H ₈ Cl ₄ N ₂ O	100mg	
Bispyribac-sodium				
CAS 125401-92-5 DRE-CA10657000	MW 452.35 Bispyribac-sodium(±)	C ₁₉ H ₁₇ N ₄ O ₈ ·Na	100mg	
Bromacil				
CAS 314-40-9 DRE-C10670000 DRE-L10670000AL DRE-XA10670000AL	MW 261.1157 Bromacil Bromacil 10 µg/mL in Acetonitrile Bromacil 100 µg/mL in Acetonitrile	C ₉ H ₁₃ BrN ₂ O ₂	250mg 10ml 1ml	
Bromacil D3 (methyl D3)				
CAS n/a DRE-XA10670100AL	MW 264.1342 Bromacil D3 (methyl D3) 100 µg/mL in Acetonitrile	C ₉ ² H ₁₃ H ₁₀ BrN ₂ O ₂	1.1ml	
Bromobutide				
CAS 74712-19-9 DRE-C10711000 DRE-L10711000CY	MW 312.2453 Bromobutide Bromobutide 10 µg/mL in Cyclohexane	C ₁₅ H ₂₂ BrNO	100mg 10ml	
Bromofenoxim				
CAS 13181-17-4 DRE-C10730000	MW 461.0192 Bromofenoxim	C ₁₃ H ₇ Br ₂ N ₃ O ₆	250mg	
Bromoxynil				
CAS 1689-84-5 DRE-C10770000 DRE-L10770000AL DRE-XA10770000AL	MW 276.9128 Bromoxynil(±) Bromoxynil 10 µg/mL in Acetonitrile Bromoxynil 100 µg/mL in Acetonitrile	C ₇ H ₅ Br ₂ NO	250mg 10ml 1ml	
Bromoxynil Heptanoate				
CAS 56634-95-8 DRE-C10775000	MW 389.0824 Bromoxynil-heptanoate	C ₁₄ H ₁₅ Br ₂ NO ₂	50mg	
Bromoxynil Methyl Ether				
CAS 3336-39-8 DRE-C10777500	MW 290.9394 Bromoxynil-methyl ether	C ₈ H ₅ Br ₂ NO	50mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Bromoxynil Octanoate				
CAS 1689-99-2 DRE-C10780000	MW 403.109 Bromoxynil-octanoate	C ₁₅ H ₁₇ Br ₂ NO ₂	250mg	
Brompyrazon				
CAS 3042-84-0 DRE-C10800000	MW 266.094 Brompyrazon	C ₁₀ H ₈ BrN ₂ O	250mg	
Bromuron				
CAS 3408-97-7 DRE-C10803000	MW 243.1004 Bromuron	C ₉ H ₁₁ BrN ₂ O	10mg	
Butachlor				
CAS 23184-66-9 DRE-C10860000 DRE-L10860000CY	MW 311.8468 Butachlor(±) Butachlor 10 µg/mL in Cyclohexane	C ₁₇ H ₂₆ ClNO ₂	100mg 10ml	
Butachlor ethane sulfonic acid (ESA) calcium salt				
CAS n/a DRE-CA10860200	MW 752.9921 Butachlor ethane sulfonic acid (ESA) calcium salt	2C ₁₇ H ₂₆ NO ₅ S·Ca	10mg	
Butafenacil				
CAS 134605-64-4 DRE-C10860800 DRE-L10860800AL	MW 474.8149 Butafenacil Butafenacil 10 µg/mL in Acetonitrile	C ₂₀ H ₁₈ ClF ₃ N ₂ O ₆	100mg 10ml	
Butamifos				
CAS 36335-67-8 DRE-C10861000 DRE-LA10861000CY	MW 332.3556 Butamifos Butamifos 10 µg/mL in Cyclohexane	C ₁₃ H ₂₁ N ₂ O ₄ PS	100mg 1ml	
Butralin				
CAS 33629-47-9 DRE-C10910000 DRE-XA09010158AL	MW 295.3342 Butralin(±) Butralin 100 µg/mL in Acetonitrile(±)	C ₁₄ H ₂₁ N ₃ O ₄	250mg 1ml	
Butroxydim				
CAS 138164-12-2 DRE-C10910500 DRE-L10910500CY	MW 399.5231 Butroxydim Butroxydim 10 µg/mL in Cyclohexane (*)	C ₂₄ H ₃₃ NO ₄	25mg 10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Buturon				
CAS 3766-60-7	MW 236.6974	C ₁₂ H ₁₃ ClN ₂ O		
DRE-C10920000	Buturon		100mg	
DRE-XA09010159ME	Buturon 100 µg/mL in Methanol(±)		1ml	
Butylate				
CAS 2008-41-5	MW 217.3714	C ₁₁ H ₂₃ NOS		
DRE-C10930000	Butylate		250mg	
DRE-XA09010160AL	Butylate 100 µg/mL in Acetonitrile(±)		1ml	
Cacodylic Acid				
CAS 75-60-5	MW 137.9974	C ₂ H ₇ AsO ₂		
DRE-C10933000	Cacodylic acid		250mg	
Cacodylic Acid Sodium Salt Trihydrate				
CAS 6131-99-3	MW 214.0251	C ₂ H ₆ AsO ₂ ·Na·3H ₂ O		
DRE-C10933200	Cacodylic acid sodium salt trihydrate		250mg	
Cafenstrole				
CAS 125306-83-4	MW 350.4359	C ₁₆ H ₂₂ N ₄ O ₃ S		
DRE-C10934500	Cafenstrole		100mg	
DRE-XA10934500AL	Cafenstrole 100 µg/mL in Acetonitrile		1ml	
Carbetamide				
CAS 16118-49-3	MW 236.267	C ₁₂ H ₁₆ N ₂ O ₃		
DRE-C11000000	Carbetamide		250mg	
DRE-L11000000AL	Carbetamide 10 µg/mL in Acetonitrile		10ml	
DRE-XA11000000AC	Carbetamide 100 µg/mL in Acetone		1ml	
Carbetamide D5 (phenyl D5)				
CAS n/a	MW 241.2978	C ₁₂ ² H ₅ H ₁₁ N ₂ O ₃		
DRE-XA11000100AC	Carbetamide D5 (phenyl D5) 100 µg/mL in Acetone(±)		1.1ml	
Carfentrazone (free acid)				
CAS 128621-72-7	MW 384.138	C ₁₃ H ₁₀ Cl ₂ F ₃ N ₃ O ₃		
DRE-C11043100	Carfentrazone (free acid)(±)		10mg	
DRE-XA11043100AL	Carfentrazone (free acid) 100 µg/mL in Acetonitrile		1.1ml	
Carfentrazone-ethyl				
CAS 128639-02-1	MW 412.1912	C ₁₅ H ₁₄ Cl ₂ F ₃ N ₃ O ₃		
DRE-C11043000	Carfentrazone-ethyl(±)		10mg	
DRE-L11043000AL	Carfentrazone-ethyl 10 µg/mL in Acetonitrile		10ml	
DRE-L11043000CY	Carfentrazone-ethyl 10 µg/mL in Cyclohexane		10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Chlormethoxyfen				
CAS 32861-85-1 DRE-L11085000AL	MW 314.1209 Chlormethoxyfen 10 µg/mL in Acetonitrile	C ₁₃ H ₉ Cl ₂ NO ₄	10ml	
Chloramben				
CAS 133-90-4 DRE-C11110000	MW 206.0261 Chloramben	C ₇ H ₅ Cl ₂ NO ₂	100mg	
Chloramben Methyl Ester				
CAS 7286-84-2 DRE-C11110800	MW 220.0527 Chloramben methyl ester	C ₈ H ₇ Cl ₂ NO ₂	10mg	
Chloranocryl				
CAS 2164-09-2 DRE-C11140000	MW 230.0906 Chloranocryl	C ₁₀ H ₈ Cl ₂ NO	250mg	
Chlorbromuron				
CAS 13360-45-7 DRE-C11180000 DRE-XA11180000AL DRE-XA09010164AL	MW 293.5449 Chlorbromuron Chlorbromuron 100 µg/mL in Acetonitrile Chlorbromuron 100 µg/mL in Acetonitrile(‡)	C ₉ H ₁₀ BrClN ₂ O ₂	250mg 1ml 1ml	
Chlorbufam				
CAS 1967-16-4 DRE-C11190000 DRE-L11190000EA	MW 223.6556 Chlorbufam(‡) Chlorbufam 10 µg/mL in Ethylacetate	C ₁₁ H ₁₀ ClNO ₂	250mg 10ml	
Chlorfenac				
CAS 85-34-7 DRE-C11247000	MW 239.4831 Chlorfenac	C ₈ H ₅ Cl ₃ O ₂	100mg	
Chlorfenprop-methyl				
CAS 14437-17-3 DRE-C11260000	MW 233.0912 Chlorfenprop-methyl	C ₁₀ H ₁₀ Cl ₂ O ₂	100mg	
Chlorflurenol				
CAS 2464-37-1 DRE-C11300000	MW 260.6725 Chlorflurenol	C ₁₄ H ₉ ClO ₃	100mg	









Pesticide and metabolites: Herbicides

Product code	Description			
Chlorflurenol Methyl Ester				
CAS 2536-31-4 DRE-C11305000	MW 274.699 Chlorflurenol-methyl ester	C ₁₅ H ₁₁ ClO ₃	250mg	
Chloridazon				
CAS 1698-60-8 DRE-C11320000 DRE-L11320000AL DRE-XA11320000AL DRE-XA09010165ME	MW 221.643 Chloridazon Chloridazon 10 µg/mL in Acetonitrile Chloridazon 100 µg/mL in Acetonitrile Chloridazon 100 µg/mL in Methanol(‡)	C ₁₀ H ₈ ClN ₃ O	250mg 10ml 1ml 1ml	
Chloridazon D5 (phenyl D5)				
CAS 1346818-99-4 DRE-XA11320100AL	MW 226.6738 Chloridazon D5 (phenyl D5) 100 µg/mL in Acetonitrile	C ₁₀ ² H ₈ H ₃ ClN ₃ O	1.1ml	
Chloridazon-desphenyl				
CAS 6339-19-1 DRE-C11322000 DRE-L11322000AL DRE-XA11322000AL	MW 145.5471 Chloridazon-desphenyl(‡) Chloridazon-desphenyl 10 µg/mL in Acetonitrile Chloridazon-desphenyl 100 µg/mL in Acetonitrile (*)	C ₄ H ₄ ClN ₃ O	10mg 10ml 1ml	
Chloridazon-methyl-desphenyl (5-Amino-4-chloro-2-methylpyridazin-3-one)				
CAS 17254-80-7 DRE-L11322500AL DRE-XA11322500AL	MW 159.5736 Chloridazon-methyl-desphenyl 10 µg/mL in Acetonitrile Chloridazon-methyl-desphenyl 100 µg/mL in Acetonitrile	C ₅ H ₆ ClN ₃ O	10ml 1ml	
Chlorimuron-ethyl				
CAS 90982-32-4 DRE-C11325000	MW 414.8208 Chlorimuron-ethyl(‡)	C ₁₅ H ₁₅ ClN ₄ O ₆ S	100mg	
Chlornitrofen				
CAS 1836-77-7 DRE-C11342000	MW 318.5399 Chlornitrofen	C ₁₂ H ₆ Cl ₃ NO ₃	100mg	
Chlornitrofen-amino (4-(2,4,6-Trichlorophenoxy)aniline)				
CAS 26306-61-6 DRE-XA11342100CY	MW 288.557 Chlornitrofen-amino 10 µg/mL in Cyclohexane	C ₁₂ H ₆ Cl ₃ NO	1ml	
Chloroacetic Acid				
CAS 79-11-8 DRE-C11348500	MW 94.497 Chloroacetic acid(‡)	C ₂ H ₃ ClO ₂	250mg	










Pesticide and metabolites: Herbicides

Product code	Description			
1-(3-Chloro-4-methylphenyl)-3-methyl-urea				
CAS 22175-22-0 DRE-C11443000	MW 198.6494 1-(3-Chloro-4-methylphenyl)-3-methyl-urea	C ₉ H ₁₁ ClN ₂ O	100mg	
1-(3-Chloro-4-methylphenyl)urea				
CAS 13142-64-8 DRE-C11442000	MW 184.6229 1-(3-Chloro-4-methylphenyl)urea	C ₈ H ₉ ClN ₂ O	100mg	
(4-Chlorophenyl)urea				
CAS 140-38-5 DRE-C11492000	MW 170.5963 1-(4-Chlorophenyl)urea	C ₇ H ₇ ClN ₂ O	100mg	
Chlorotoluron				
CAS 15545-48-9 DRE-C11530000 DRE-L11530000AL DRE-XA11530000AL DRE-XA09010166AL	MW 212.676 Chlorotoluron(‡) Chlorotoluron 10 µg/mL in Acetonitrile Chlorotoluron 100 µg/mL in Acetonitrile Chlortoluron 100 µg/mL in Acetonitrile(‡)	C ₁₀ H ₁₃ ClN ₂ O	250mg 10ml 1ml 1ml	
Chlorotoluron D6 (N,N-dimethyl D6)				
CAS 1219803-48-1 DRE-C11530100 DRE-XA11530100AC	MW 218.713 Chlorotoluron (N,N-dimethyl) D6 Chlorotoluron D6 100 µg/mL in Acetone	C ₁₀ ² H ₆ H ₇ ClN ₂ O	5mg 1.1ml	
Chloroxuron				
CAS 1982-47-4 DRE-C11540000 DRE-L11540000AL	MW 290.7448 Chloroxuron Chloroxuron 10 µg/mL in Acetonitrile	C ₁₅ H ₁₅ ClN ₂ O ₂	250mg 10ml	
Chloroxuron-N-monodesmethyl				
CAS 28365-79-9 DRE-C11542000	MW 276.7182 Chloroxuron-N-monodesmethyl	C ₁₄ H ₁₃ ClN ₂ O ₂	10mg	
Chloroxynil				
CAS 1891-95-8 DRE-C11550000	MW 188.0108 Chloroxynil	C ₇ H ₃ Cl ₂ NO	25mg	
Chlorphonium Chloride				
CAS 115-78-6 DRE-L11565000ME	MW 397.7901 Chlorphonium-Cl 10 µg/mL in Methanol	C ₁₉ H ₃₂ Cl ₂ P·Cl	10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Chlorpropham				
CAS 101-21-3	MW 213.6608	C ₁₀ H ₁₂ ClNO ₂		
DRE-C11580000	Chlorpropham		250mg	
DRE-L11580000AL	Chlorpropham 10 µg/mL in Acetonitrile		10ml	
DRE-L11580000IO	Chlorpropham 10 µg/mL in Isooctane		10ml	
DRE-XA11580000AL	Chlorpropham 100 µg/mL in Acetonitrile		1ml	
DRE-GS09010096ME	Chlorpropham 100 µg/mL in Methanol(‡)		5x1ml	
DRE-LA09010254IO	Chlorpropham 10 µg/mL in Isooctane(‡)		10ml	
Chlorsulfuron				
CAS 64902-72-3	MW 357.7728	C ₁₂ H ₁₂ ClN ₅ O ₄ S		
DRE-C11610000	Chlorsulfuron(‡)		100mg	
Chlorthal-diacid				
CAS 2136-79-0	MW 303.9111	C ₈ H ₂ Cl ₄ O ₄		
DRE-C11619000	Chlorthal-diacid		10mg	
Chlorthal-dimethyl				
CAS 1861-32-1	MW 331.9642	C ₁₀ H ₆ Cl ₄ O ₄		
DRE-C11620000	Chlorthal-dimethyl		250mg	
DRE-GH11620000HE	Chlorthal-dimethyl 200 µg/mL in Hexane		5x1ml	
DRE-L11620000CY	Chlorthal-dimethyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA11620000CY	Chlorthal-dimethyl 100 µg/mL in Cyclohexane		1ml	
Chlorthiamid				
CAS 1918-13-4	MW 206.0923	C ₇ H ₅ Cl ₂ NS		
DRE-C11630000	Chlorthiamid		100mg	
DRE-XA09010168ME	Chlorthiamid 100 µg/mL in Methanol(‡)		1ml	
Cinidon (free acid)				
CAS 175156-71-5	MW 366.1954	C ₁₇ H ₁₃ Cl ₂ NO ₄		
DRE-XA11667410AL	Cinidon (free acid) 100 µg/mL in Acetonitrile		1ml	
Cinidon-ethyl				
CAS 142891-20-1	MW 394.2486	C ₁₉ H ₁₇ Cl ₂ NO ₄		
DRE-C11667400	Cinidon-ethyl		100mg	
Cinmethylin				
CAS 87818-31-3	MW 274.3978	C ₁₈ H ₂₆ O ₂		
DRE-C11667450	Cinmethylin		50mg	
DRE-L11667450CY	Cinmethylin 10 µg/mL in Cyclohexane		10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Cinosulfuron				
CAS 94593-91-6 DRE-C11668000	MW 413.4057 Cinosulfuron	C ₁₅ H ₁₉ N ₅ O ₇ S	100mg	
Clethodim				
CAS 99129-21-2 DRE-C11669000	MW 359.9112 Clethodim (*)	C ₁₇ H ₂₆ ClNO ₃ S	100mg	
Clodinafop free acid				
CAS 114420-56-3 DRE-C11678900 DRE-L11678900AL	MW 311.6928 Clodinafop free acid Clodinafop free acid 10 µg/mL in Acetonitrile	C ₁₄ H ₁₁ ClFNO ₄	100mg 10ml	
Clodinafop-propargyl ester				
CAS 105512-06-9 DRE-C11679000 DRE-L11679000AL DRE-XA11679000AL DRE-XA09010169AL	MW 349.7408 Clodinafop-propargyl ester Clodinafop-propargyl ester 10 µg/mL in Acetonitrile Clodinafop-propargyl ester 100 µg/mL in Acetonitrile Clodinafop-propargyl 100 µg/mL in Acetonitrile(‡)	C ₁₇ H ₁₃ ClFNO ₄	100mg 10ml 1ml 1ml	
Clomazone (Command)				
CAS 81777-89-1 DRE-C11685000 DRE-L11685000AL DRE-L11685000CY DRE-XA09010170ME	MW 239.6981 Clomazone(‡) Clomazone 10 µg/mL in Acetonitrile Clomazone 10 µg/mL in Cyclohexane Command 100 µg/mL in Methanol(‡)	C ₁₂ H ₁₄ ClNO ₂	100mg 10ml 10ml 1ml	
Clomeprop				
CAS 84496-56-0 DRE-C11686000 DRE-XA11686000AL	MW 324.2018 Clomeprop(‡) Clomeprop 100 µg/mL in Acetonitrile	C ₁₆ H ₁₅ Cl ₂ NO ₂	10mg 1ml	
Clopyralid				
CAS 1702-17-6 DRE-C11690000 DRE-L11690000AL DRE-XA11690000AL	MW 191.9995 Clopyralid Clopyralid 10 µg/mL in Acetonitrile Clopyralid 100 µg/mL in Acetonitrile	C ₈ H ₅ Cl ₂ NO ₂	250mg 10ml 1ml	
Cloransulam-methyl				
CAS 147150-35-4 DRE-C11691310	MW 429.8106 Cloransulam-methyl(‡)	C ₁₅ H ₁₃ ClFN ₅ O ₅ S	100mg	
Cumyluron				
CAS 99485-76-4 DRE-C11775000	MW 302.7986 Cumyluron(‡)	C ₁₇ H ₁₉ ClN ₂ O	50mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Cyanamide				
CAS 420-04-2	MW 42.04	CH ₂ N ₂		
DRE-CA11788000	Cyanamide		100mg	
Cyanazine				
CAS 21725-46-2	MW 240.6927	C ₈ H ₁₃ ClN ₆		
DRE-C11790000	Cyanazine		250mg	
DRE-L11790000AL	Cyanazine 10 µg/mL in Acetonitrile		10ml	
DRE-XA11790000AL	Cyanazine 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010174AL	Cyanazine 100 µg/mL in Acetonitrile(±)		1ml	
Cyanazine D5 (N-ethyl D5)				
CAS 1190003-29-2	MW 245.7235	C ₉ H ₅ H ₈ ClN ₆		
DRE-XA11790100AC	Cyanazine D5 (N-ethyl D5) 100 µg/mL in Acetone		1.1ml	
Cycloate				
CAS 1134-23-2	MW 215.3555	C ₁₁ H ₂₁ NOS		
DRE-C11820000	Cycloate		250mg	
Cyclosulfamuron				
CAS 136849-15-5	MW 421.4277	C ₁₇ H ₁₉ N ₅ O ₆ S		
DRE-C11836500	Cyclosulfamuron		100mg	
Cycloxydim				
CAS 101205-02-1	MW 325.4662	C ₁₇ H ₂₇ NO ₃ S		
DRE-C11837000	Cycloxydim		100mg	
DRE-XA09010175ME	Cycloxydim 100 µg/mL in Methanol(±)		1ml	
Cycluron				
CAS 2163-69-1	MW 198.3052	C ₁₁ H ₂₃ N ₂ O		
DRE-C11840000	Cycluron		250mg	
DRE-XA09010176AL	Cycluron 100 µg/mL in Acetonitrile(±)		1ml	
Cyhalofop				
CAS 122008-78-0	MW 301.2692	C ₁₆ H ₁₂ FNO ₄		
DRE-C11857000	Cyhalofop		25mg	
Cyhalofop-butyl				
CAS 122008-85-9	MW 357.3755	C ₂₀ H ₂₀ FNO ₄		
DRE-C11858000	Cyhalofop-butyl		50mg	
DRE-L11858000CY	Cyhalofop-butyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010177AL	Cyhalofop-butyl 100 µg/mL in Acetonitrile(±)		1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Cyprazine				
CAS 22936-86-3	MW 227.694	C ₉ H ₁₄ ClN ₅		
DRE-C11900000	Cyprazine		100mg	
DRE-XA09010178AL	Cyprazine 100 µg/mL in Acetonitrile(±)		1ml	
Cyprazine-desisopropyl				
CAS 35516-73-5	MW 185.6142	C ₈ H ₉ ClN ₅		
DRE-C11900200	Cyprazine-desisopropyl		25mg	
Cyprazine-desisopropyl-2-hydroxy				
CAS 92510-62-8	MW 167.1686	C ₈ H ₉ N ₅ O		
DRE-C11900400	Cyprazine-desisopropyl-2-hydroxy		10mg	
Cyprazine-2-hydroxy				
CAS 39095-16-4	MW 209.2483	C ₉ H ₁₅ N ₅ O		
DRE-C11900800	Cyprazine-2-hydroxy		25mg	
2,4-D (phenyl-13C6)				
CAS 150907-52-1	MW 226.9934	¹³ C ₆ C ₂ H ₆ Cl ₂ O ₃		
DRE-XA11940200AC	2,4-D 13C6 100 µg/mL in Acetone(±)		1ml	
2,4-D 2-Butoxyethyl Ester ((2,4-Dichlorophenoxy)acetic Acid 2-Butoxyethyl Ester)				
CAS 1929-73-3	MW 321.1963	C ₁₄ H ₁₈ Cl ₂ O ₄		
DRE-C11942000	2,4-D-butylglycol ester		250mg	
2,4-D Butyl Ester ((2,4-Dichlorophenoxy)acetic Acid Butyl Ester)				
CAS 94-80-4	MW 277.1438	C ₁₂ H ₁₄ Cl ₂ O ₃		
DRE-C11941000	2,4-D-1-butyl ester		250mg	
DRE-XA11941000CY	2,4-D-1-butyl ester 100 µg/mL in Cyclohexane		1ml	
2,4-D Dimethylamine Salt				
CAS 2008-39-1	MW 266.1211	C ₈ H ₆ Cl ₂ O ₃ ·C ₂ H ₇ N		
DRE-C11940500	2,4-D dimethylamine salt		250mg	
2,4-D Ethyl Ester ((2,4-Dichlorophenoxy)acetic Acid Ethyl Ester)				
CAS 533-23-3	MW 249.0906	C ₁₀ H ₁₀ Cl ₂ O ₃		
DRE-C11942500	2,4-D-ethyl ester		250mg	










Pesticide and metabolites: Herbicides

Product code	Description			
2,4-D 2-Ethylhexyl Ester				
CAS 1928-43-4 DRE-C11942700	MW 333.2501 2,4-D 2-ethylhexyl ester(‡)	C ₁₆ H ₂₂ Cl ₂ O ₃	100mg	
2,4-D Isobutyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isobutyl Ester)				
CAS 1713-15-1 DRE-C11943000 DRE-XA11943000CY	MW 277.1438 2,4-D-iso-butyl ester 2,4-D-iso-butyl ester 100 µg/mL in Cyclohexane	C ₁₂ H ₁₄ Cl ₂ O ₃	250mg 1ml	
2,4-D Isooctyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isooctyl Ester)				
CAS 25168-26-7 DRE-C11944000	MW 333.2501 2,4-D-iso-octyl ester (technical)	C ₁₆ H ₂₂ Cl ₂ O ₃	250mg	
2,4-D Isopropyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isopropyl Ester)				
CAS 94-11-1 DRE-C11944300 DRE-XA11944300CY	MW 263.1172 2,4-D-iso-propyl ester 2,4-D-iso-propyl ester 100 µg/mL in Cyclohexane	C ₁₁ H ₁₂ Cl ₂ O ₃	250mg 1ml	
2,4-D Methyl Ester ((2,4-Dichlorophenoxy)acetic Acid Methyl Ester)				
CAS 1928-38-7 DRE-C11945000 DRE-L11945000CY	MW 235.064 2,4-D-methyl ester 2,4-D-methyl ester 10 µg/mL in Cyclohexane	C ₉ H ₈ Cl ₂ O ₃	250mg 10ml	
2,4-D Sodium Salt Monohydrate				
CAS 7084-86-8 DRE-C11946000	MW 261.0346 2,4-D-sodium salt monohydrate(‡)	C ₈ H ₅ Cl ₂ O ₃ ·Na·H ₂ O	250mg	
2,3-D (2,3-Dichlorophenoxyacetic Acid)				
CAS 2976-74-1 DRE-C11939000	MW 221.0374 2,3-D	C ₈ H ₆ Cl ₂ O ₃	100mg	
2,4-D ((2,4-Dichlorophenoxy)acetic Acid)				
CAS 94-75-7 DRE-C11940000 DRE-L11940000AL DRE-XA11940000AL DRE-YS09010032AL DRE-XA09010250ME	MW 221.0374 2,4-D(‡) 2,4-D 10 µg/mL in Acetonitrile 2,4-D 100 µg/mL in Acetonitrile 2,4-Dichlorophenoxyacetic Acid 200 µg/mL in Acetonitrile(‡) 2,4-D 100 µg/mL in Methanol(‡)	C ₈ H ₆ Cl ₂ O ₃	250mg 10ml 1ml 5x1ml 1ml	
2,4-D D3 ((2,4-Dichloro-3,5,6-trideuteriophenoxy)acetic Acid)				
CAS 202480-67-9 DRE-C11940100 DRE-XA11940100AC	MW 224.0559 2,4-D D3 2,4-D D3 100 µg/mL in Acetone(‡)	C ₈ ² H ₃ H ₃ Cl ₂ O ₃	10mg 1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
2,6-D				
CAS 575-90-6 DRE-C11946600	MW 221.0374 2,6-D	C ₈ H ₆ Cl ₂ O ₃	10mg	
3,4-D ((3,4-Dichlorophenoxy)acetic Acid)				
CAS 588-22-7 DRE-C11947000	MW 221.0374 3,4-D	C ₈ H ₆ Cl ₂ O ₃	250mg	
2,4-DB Dimethylamine Salt				
CAS 2758-42-1 DRE-C11980500	MW 294.1743 2,4-DB dimethylamine salt	C ₁₀ H ₁₀ Cl ₂ O ₃ ·C ₂ H ₇ N	100mg	
2,4-DB-2-ethylhexyl ester				
CAS 1320-15-6 DRE-C11981000	MW 361.3032 2,4-DB-2-ethylhexyl ester	C ₁₈ H ₂₆ Cl ₂ O ₃	250mg	
2,4-DB Methyl Ester (4-(2,4-Dichlorophenoxy)butanoic Acid Methyl Ester)				
CAS 18625-12-2 DRE-C11982000	MW 263.1172 2,4-DB-methyl ester	C ₁₁ H ₁₂ Cl ₂ O ₃	250mg	
2,4-DB (4-(2,4-Dichlorophenoxy)butanoic Acid)				
CAS 94-82-6 DRE-C11980000 DRE-L11980000AL DRE-XA11980000AL DRE-YS09010033MB	MW 249.0906 2,4-DB 2,4-DB 10 µg/mL in Acetonitrile 2,4-DB 100 µg/mL in Acetonitrile 2,4-DB 200 µg/mL in Methyl tert-Butyl Ether(†)	C ₁₀ H ₁₀ Cl ₂ O ₃	250mg 10ml 1ml 5x1ml	
DNOC (2-Methyl-4,6-dinitrophenol)				
CAS 534-52-1 DRE-C13050000 DRE-L13050000AL DRE-XA13050000ME	MW 198.1329 DNOC DNOC 10 µg/mL in Acetonitrile DNOC 100 µg/mL in Methanol	C ₇ H ₅ N ₂ O ₅	250mg 10ml 1ml	
DNOC D5 (ring D2, methyl D3)				
CAS n/a DRE-XA13050100AC	MW 203.1637 DNOC D5 (ring D2, methyl D3) 100 µg/mL in Acetone	C ₇ H ₅ HN ₂ O ₅	1.1ml	
Daimuron (Dymron)				
CAS 42609-52-9 DRE-C11948000 DRE-LA11948000AL	MW 268.3535 Daimuron (Dymron) Daimuron (Dymron) 10 µg/mL in Acetonitrile	C ₁₇ H ₂₀ N ₂ O	100mg 1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Dicamba-desmethyl				
CAS 3401-80-7 DRE-C12260300	MW 207.0109 Dicamba-desmethyl	C ₇ H ₄ Cl ₂ O ₃	10mg	
Dicamba-5-hydroxy (5-Hydroxydicamba)				
CAS 7600-50-2 DRE-C12260500 DRE-XA12260500AL	MW 237.0368 Dicamba-5-hydroxy Dicamba-5-hydroxy 100 µg/mL in Acetonitrile	C ₈ H ₆ Cl ₂ O ₄	25mg 1ml	
Dichlobenil				
CAS 1194-65-6 DRE-C12280000 DRE-L12280000CY	MW 172.0114 Dichlobenil Dichlobenil 10 µg/mL in Cyclohexane	C ₇ H ₃ Cl ₂ N	250mg 10ml	
1-(3,4-Dichlorophenyl)-3-methoxyurea				
CAS 17356-61-5 DRE-C12471250	MW 235.0673 1-(3,4-Dichlorophenyl)-3-methoxyurea	C ₈ H ₈ Cl ₂ N ₂ O ₂	50mg	
1-(3,4-Dichlorophenyl)-3-methylurea				
CAS 3567-62-2 DRE-C12471500 DRE-L12471500AL DRE-XA12471500AL	MW 219.0679 1-(3,4-Dichlorophenyl)-3-methyl urea 1-(3,4-Dichlorophenyl)-3-methyl urea 10 µg/mL in Acetonitrile 1-(3,4-Dichlorophenyl)-3-methylurea 100 µg/mL in Acetonitrile	C ₈ H ₈ Cl ₂ N ₂ O	100mg 10ml 1ml	
1-(3,5-Dichlorophenyl)-3-methylurea				
CAS 38655-13-9 DRE-C12471510	MW 219.0679 1-(3,5-Dichlorophenyl)-3-methylurea	C ₈ H ₈ Cl ₂ N ₂ O	100mg	
1-(3,4-Dichlorophenyl)urea				
CAS 2327-02-8 DRE-C12472000 DRE-L12472000AL	MW 205.0413 1-(3,4-Dichlorophenyl)urea 1-(3,4-Dichlorophenyl)urea 10 µg/mL in Acetonitrile	C ₇ H ₆ Cl ₂ N ₂ O	100mg 10ml	
Dichlorprop (2,4-DP)				
CAS 120-36-5 DRE-C12510000 DRE-L12510000AL DRE-XA12510000AL DRE-YS09010034MB	MW 235.064 Dichlorprop Dichlorprop 10 µg/mL in Acetonitrile Dichlorprop 100 µg/mL in Acetonitrile Dichlorprop (2,4-DP) 200 µg/mL in Methyl tert-Butyl Ether(‡)	C ₈ H ₈ Cl ₂ O ₃	250mg 10ml 1ml 5x1ml	
Dichlorprop D6 (ring D3, 3,3,3-D3)				
CAS n/a DRE-C12510100 DRE-XA12510100AC	MW 241.101 Dichlorprop D6 (ring D3, 3,3,3-D3) Dichlorprop (ring-D3, 3,3,3-D3) 100 µg/mL in Acetone	C ₈ ² H ₆ H ₂ Cl ₂ O ₃	10mg 1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Dichlorprop 2-Ethylhexyl Ester				
CAS 79270-78-3 DRE-C12512000	MW 347.2767 Dichlorprop-2-ethylhexyl ester	C ₁₇ H ₂₄ Cl ₂ O ₃	250mg	
Dichlorprop-P (D-(+)-Dichlorprop)				
CAS 15165-67-0 DRE-C12510500	MW 235.064 D(+) Dichlorprop	C ₉ H ₉ Cl ₂ O ₃	250mg	
Dichlorprop-butoxyethyl ester				
CAS 53404-31-2 DRE-C12511000	MW 335.2229 Dichlorprop-butoxyethyl ester	C ₁₅ H ₂₀ Cl ₂ O ₄	50mg	
Dichlorprop-methyl Ester				
CAS 23844-57-7 DRE-C12514000 DRE-L12514000CY	MW 249.0906 Dichlorprop-methyl ester Dichlorprop-methyl ester 10 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₂ O ₃	100mg 10ml	
2,6-Dichlorprop				
CAS 25140-90-3 DRE-C12515000	MW 235.064 2,6-Dichlorprop	C ₉ H ₉ Cl ₂ O ₃	50mg	
Diclofop methyl				
CAS 51338-27-3 DRE-C12540000 DRE-L12540000CY DRE-GA09010063ME	MW 341.186 Diclofop methyl Diclofop methyl 10 µg/mL in Cyclohexane Diclofop-methyl 100 µg/mL in Methanol(‡)	C ₁₆ H ₁₄ Cl ₂ O ₄	250mg 10ml 1ml	
Diclofop-methyl-5-hydroxy (2-[4-(2,4-Dichloro-5-hydroxyphenoxy)phenoxy]propanoic Acid Methyl Ester)				
CAS 124992-48-9 DRE-C12542000	MW 357.1854 Diclofop-methyl-5-hydroxy	C ₁₆ H ₁₄ Cl ₂ O ₅	10mg	
Diclosulam				
CAS 145701-21-9 DRE-C12560300	MW 406.2196 Diclosulam	C ₁₃ H ₁₀ Cl ₂ FN ₅ O ₃ S	25mg	
Diethatyl Ethyl Ester				
CAS 38727-55-8 DRE-C12602000	MW 311.8038 Diethatyl-ethyl	C ₁₆ H ₂₂ ClNO ₃	100mg	

Pesticide and metabolites: Herbicides

Product code	Description			
Difenoxuron				
CAS 14214-32-5	MW 286.3257	C ₁₆ H ₁₈ N ₂ O ₃		
DRE-C12610000	Difenoxuron		250mg	
DRE-L12610000AL	Difenoxuron 10 µg/mL in Acetonitrile		10ml	
Difenzoquat Methyl Sulfate (1,2-Dimethyl-3,5-diphenylpyrazolium methyl sulfate)				
CAS 43222-48-6	MW 360.4274	C ₁₇ H ₁₇ N ₂ ·CH ₃ O ₄ S		
DRE-C12620000	Difenzoquat methyl sulfate		250mg	
Diflufenican				
CAS 83164-33-4	MW 394.2949	C ₁₉ H ₁₁ F ₅ N ₂ O ₂		
DRE-C12631000	Diflufenican		100mg	
DRE-L12631000AC	Diflufenican 10 µg/mL in Acetone		10ml	
DRE-L12631000AL	Diflufenican 10 µg/mL in Acetonitrile		10ml	
DRE-XA12631000AL	Diflufenican 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010181ME	Diflufenican 100 µg/mL in Methanol(‡)		1ml	
Diflufenican D3 (3-trifluoromethylphenoxy-2,4,6 D3)				
CAS 1185009-29-3	MW 397.3133	C ₁₉ ² H ₃ H ₈ F ₉ N ₂ O ₂		
DRE-XA12631001AL	Diflufenican D3 (3-trifluoromethylphenoxy-2,4,6 D3) 100 µg/mL in Acetonitrile		1.1ml	
Diflufenzopyr				
CAS 109293-97-2	MW 334.2776	C ₁₅ H ₁₂ F ₂ N ₄ O ₃		
DRE-C12631030	Diflufenzopyr		100mg	
Diflufenzopyr sodium salt				
CAS 109293-98-3	MW 356.2594	C ₁₅ H ₁₁ F ₂ N ₄ O ₃ ·Na		
DRE-C12631032	Diflufenzopyr sodium salt		100mg	
Dimefuron				
CAS 34205-21-5	MW 338.7894	C ₁₅ H ₁₉ ClN ₄ O ₃		
DRE-C12660000	Dimefuron(‡)		100mg	
DRE-L12660000AL	Dimefuron 10 µg/mL in Acetonitrile		10ml	
DRE-XA09010182ME	Dimefuron 100 µg/mL in Methanol(‡)		1ml	
Dimepiperate				
CAS 61432-55-1	MW 263.3983	C ₁₅ H ₂₁ NOS		
DRE-C12665000	Dimepiperate(‡)		100mg	
DRE-XA12665000EA	Dimepiperate 100 µg/mL in Ethylacetate		1ml	
Dimethachlor				
CAS 50563-36-5	MW 255.7405	C ₁₃ H ₁₈ ClNO ₂		
DRE-C12670000	Dimethachlor(‡)		250mg	
DRE-L12670000CY	Dimethachlor 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010183ME	Dimethachlor 100 µg/mL in Methanol(‡)		1ml	










Pesticide and metabolites: Herbicides

Product code	Description				
Dimethachlor-deschloro					
CAS 1231710-71-6 DRE-C12670100	MW 221.2955 Dimethachlor-deschloro	C ₁₃ H ₁₉ NO ₂	25mg		
Dimethachlor ethane sulfonic acid (ESA) sodium salt					
CAS 1231710-75-0 DRE-CA12670200	MW 323.3405 Dimethachlor ethane sulfonic acid (ESA) sodium salt	C ₁₃ H ₁₈ NO ₅ S·Na	10mg		
Dimethachlor oxalamic acid (OA)					
CAS 1086384-49-7 DRE-CA12670400	MW 251.2784 Dimethachlor oxalamic acid (OA)	C ₁₃ H ₁₇ NO ₄	10mg		
Dimethametryn					
CAS 22936-75-0 DRE-C12675000	MW 255.3829 Dimethametryn	C ₁₁ H ₂₁ N ₅ S	250mg		
Dimethenamid					
CAS 87674-68-8 DRE-C12677500 DRE-L12677500ME	MW 275.7948 Dimethenamid(‡) Dimethenamid 10 µg/mL in Methanol	C ₁₂ H ₁₈ ClNO ₂ S	100mg 10ml		
Dimethenamid-P					
CAS 163515-14-8 DRE-C12678000 DRE-L12678000ME	MW 275.7948 Dimethenamid-P Dimethenamid-P 10 µg/mL in Methanol	C ₁₂ H ₁₈ ClNO ₂ S	100mg 10ml		
Dimethipin					
CAS 55290-64-7 DRE-C12680000 DRE-GA09010367ME	MW 210.2712 Dimethipin Dimethipin 100 µg/mL in Methanol(‡)	C ₈ H ₁₀ O ₄ S ₂	100mg 1ml		
Dinitramine					
CAS 29091-05-2 DRE-C12780000 DRE-XA09010185ME	MW 322.2405 Dinitramine(‡) Dinitramine 100 µg/mL in Methanol(‡)	C ₁₁ H ₁₃ F ₃ N ₄ O ₄	250mg 1ml		
Dinoseb					
CAS 88-85-7 DRE-C12810000 DRE-L12810000AL DRE-XA12810000AL DRE-GS09010080MB	MW 240.2127 Dinoseb Dinoseb 10 µg/mL in Acetonitrile Dinoseb 100 µg/mL in Acetonitrile Dinoseb 200 µg/mL in Methyl tert-Butyl Ether(‡)	C ₁₀ H ₁₂ N ₂ O ₅	100mg 10ml 1ml 5x1ml		

Pesticide and metabolites: Herbicides

Product code	Description			
Dinoseb Methyl Ether				
CAS 6099-79-2 DRE-C12811500	MW 254.2393 Dinoseb-methyl ether	C ₁₁ H ₁₄ N ₂ O ₅	100mg	
Dinoseb acetate				
CAS 2813-95-8 DRE-C12811000 DRE-XA12811000CY	MW 282.2494 Dinoseb acetate Dinoseb acetate 100 µg/mL in Cyclohexane	C ₁₂ H ₁₄ N ₂ O ₆	250mg 1ml	
Dinoterb				
CAS 1420-07-1 DRE-C12830000 DRE-L12830000AL DRE-XA12830000AL	MW 240.2127 Dinoterb Dinoterb 10 µg/mL in Acetonitrile Dinoterb 100 µg/mL in Acetonitrile	C ₁₀ H ₁₂ N ₂ O ₅	100mg 10ml 1ml	
Dinoterb Acetate				
CAS 3204-27-1 DRE-C12831000	MW 282.2494 Dinoterb acetate	C ₁₂ H ₁₄ N ₂ O ₆	100mg	
Diphenamid				
CAS 957-51-7 DRE-C12880000	MW 239.3123 Diphenamid	C ₁₆ H ₁₇ NO	250mg	
Dipropetryn				
CAS 4147-51-7 DRE-C12930000 DRE-L12930000CY DRE-XA09010186ME	MW 255.3829 Dipropetryn Dipropetryn 10 µg/mL in Cyclohexane Dipropetryn 100 µg/mL in Methanol(‡)	C ₁₁ H ₂₁ N ₅ S	100mg 10ml 1ml	
Diquat Dibromide				
CAS 85-00-7 DRE-CA12960000 DRE-XA12960000WA	MW 344.0451 Diquat dibromide(‡) Diquat dibromide 100 µg/mL in Water	C ₁₂ H ₁₂ N ₂ ·2Br	250mg 1ml	
Diquat dibromide D4				
CAS n/a DRE-CA12960010	MW 348.0697 Diquat dibromide D4(‡)	C ₁₂ ² H ₄ H ₈ N ₂ ·2Br	50mg	
Disul Sodium				
CAS 136-78-7 DRE-C12974000	MW 309.099 Disul sodium	C ₈ H ₇ Cl ₂ O ₅ S·Na	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Dithiopyr				
CAS 97886-45-8 DRE-C13013700	MW 401.4151 Dithiopyr(±)	C ₁₅ H ₁₆ F ₃ NO ₂ S ₂	100mg	
Diuron				
CAS 330-54-1 DRE-C13020000 DRE-L13020000AL DRE-XA13020000AL	Diuron Diuron 10 µg/mL in Acetonitrile Diuron 100 µg/mL in Acetonitrile	C ₉ H ₆ Cl ₂ N ₂ O	250mg 10ml 1ml	
Diuron D6				
CAS 1007536-67-5 DRE-C13020100 DRE-XA13020100AC	Diuron D6 Diuron D6 100 µg/mL in Acetone	C ₉ ² H ₆ H ₄ Cl ₂ N ₂ O	10mg 1.1ml	
EPTC				
CAS 759-94-4 DRE-C13190000 DRE-L13190000CY	EPTC(±) EPTC 10 µg/mL in Cyclohexane	C ₉ H ₁₉ NOS	250mg 10ml	
Endothal				
CAS 145-73-3 DRE-C13150000	Endothal	C ₈ H ₁₀ O ₅	100mg	
Erbon				
CAS 136-25-4 DRE-L13200000CY	Erbon 10 µg/mL in Cyclohexane	C ₁₁ H ₉ Cl ₅ O ₃	10ml	
Esprocarb				
CAS 85785-20-2 DRE-C13212000 DRE-LA13212000AC DRE-XA13212000AC	Esprocarb(±) Esprocarb 10 µg/mL in Acetone Esprocarb 100 µg/mL in Acetone	C ₁₅ H ₂₃ NOS	100mg 1ml 1ml	
Ethalfuralin				
CAS 55283-68-6 DRE-C13220000 DRE-L13220000CY	Ethalfuralin(±) Ethalfuralin 10 µg/mL in Cyclohexane	C ₁₃ H ₁₄ F ₃ N ₃ O ₄	250mg 10ml	
Ethametsulfuron-methyl				
CAS 97780-06-8 DRE-C13222000	Ethametsulfuron-methyl(±)	C ₁₅ H ₁₈ N ₆ O ₆ S	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Ethidimuron				
CAS 30043-49-3 DRE-C13240000	MW 264.3252 Ethidimuron	C ₇ H ₁₂ N ₄ O ₃ S ₂	250mg	
Ethiolate				
CAS 2941-55-1 DRE-C13260000	MW 161.2651 Ethiolate	C ₇ H ₁₅ NOS	100mg	
Ethiozin				
CAS 64529-56-2 DRE-C13272000 DRE-XA09010247ME	MW 228.3145 Ethiozin Ethiozin 100 µg/mL in Methanol(‡)	C ₉ H ₁₆ N ₄ OS	100mg 1ml	
Ethofumesate				
CAS 26225-79-6 DRE-C13290000 DRE-L13290000AL DRE-XA13290000AL	MW 286.344 Ethofumesate(‡) Ethofumesate 10 µg/mL in Acetonitrile Ethofumesate 100 µg/mL in Acetonitrile	C ₁₃ H ₁₈ O ₅ S	250mg 10ml 1ml	
Ethofumesate-2-hydroxy				
CAS 26322-82-7 DRE-C13290800	MW 258.2909 Ethofumesate-2-hydroxy	C ₁₁ H ₁₄ O ₅ S	10mg	
Ethofumesate-2-keto				
CAS 26244-33-7 DRE-C13291000 DRE-LA13291000CY	MW 256.275 Ethofumesate-2-keto Ethofumesate-2-keto 10 µg/mL in Cyclohexane	C ₁₁ H ₁₂ O ₅ S	10mg 1ml	
Ethoxyfen-ethyl				
CAS 131086-42-5 DRE-L13308500CY	MW 416.7755 Ethoxyfen-ethyl 10 µg/mL in Cyclohexane	C ₁₉ H ₁₆ ClF ₃ O ₅	10ml	
Ethoxysulfuron				
CAS 126801-58-9 DRE-C13311000	MW 398.391 Ethoxysulfuron	C ₁₅ H ₁₈ N ₄ O ₇ S	100mg	
Ethylene Glycol Bis(trichloroacetate)				
CAS 2514-53-6 DRE-C13328500	MW 352.8116 Ethylene glycol bis(trichloroacetate)	C ₆ H ₄ Cl ₆ O ₄	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Etobenzanid				
CAS 79540-50-4 DRE-C13360000	MW 340.2012 Etobenzanid	C ₁₆ H ₁₅ Cl ₂ NO ₃	50mg	
Fenchlorazol-ethyl				
CAS 103112-35-2 DRE-C13457000	MW 403.4758 Fenchlorazol-ethyl	C ₁₂ H ₆ Cl ₅ N ₃ O ₂	100mg	
Fenoprop				
CAS 93-72-1 DRE-C13490000 DRE-L13490000AL DRE-XA13490000AL	MW 269.5091 Fenoprop Fenoprop 10 µg/mL in Acetonitrile Fenoprop 100 µg/mL in Acetonitrile	C ₉ H ₇ Cl ₃ O ₃	100mg 10ml 1ml	
Fenoprop-methyl ester				
CAS 4841-20-7 DRE-C13495000 DRE-L13495000EA	MW 283.5357 Fenoprop-methyl ester Fenoprop-methyl ester 10 µg/mL in Ethylacetate	C ₁₀ H ₇ Cl ₃ O ₃	100mg 10ml	
Fenoxaprop				
CAS 95617-09-7 DRE-C13499500	MW 333.7232 Fenoxaprop	C ₁₆ H ₁₂ ClNO ₅	100mg	
Fenoxaprop-P (R-enantiomer)				
CAS 113158-40-0 DRE-C13499600 DRE-L13499600AL	MW 333.7232 Fenoxaprop-P (R-enantiomer) Fenoxaprop-P (R-enantiomer) 10 µg/mL in Acetonitrile	C ₁₆ H ₁₂ ClNO ₅	100mg 10ml	
Fenoxaprop-P-ethyl (R-enantiomer)				
CAS 71283-80-2 DRE-C13510500	MW 361.7763 Fenoxaprop-P-ethyl (R-enantiomer)(±)	C ₁₈ H ₁₆ ClNO ₅	250mg	
Fenoxaprop-ethyl				
CAS 66441-23-4 DRE-C13510000 DRE-XA09010196ME	MW 361.7763 Fenoxaprop-ethyl Fenoxaprop-ethyl 100 µg/mL in Methanol(±)	C ₁₈ H ₁₆ ClNO ₅	100mg 1ml	
Fentrazamide				
CAS 158237-07-1 DRE-C13610000	MW 349.8153 Fentrazamide	C ₁₆ H ₂₀ ClN ₅ O ₂	100mg	










Pesticide and metabolites: Herbicides

Product code	Description		
Fentrazamide Metabolite 1			
CAS 98377-35-6	MW 196.5938	C ₇ H ₅ ClN ₄ O	
DRE-L13610200AL	Fentrazamide Metabolite 1 10 µg/mL in Acetonitrile		10ml
DRE-XA13610200AL	Fentrazamide Metabolite 1 100 µg/mL in Acetonitrile		1ml
			
Fenuron			
CAS 101-42-8	MW 164.2044	C ₉ H ₁₂ N ₂ O	
DRE-C13620000	Fenuron(‡)		250mg
DRE-L13620000AL	Fenuron 10 µg/mL in Acetonitrile		10ml
			
Fenuron D5 (phenyl D5)			
CAS 1219802-06-8	MW 169.2352	C ₉ ² H ₅ H ₇ N ₂ O	
DRE-XA13620010AL	Fenuron D5 (phenyl D5) 100 µg/mL in Acetonitrile		1.1ml
			
Flamprop (free acid)			
CAS 58667-63-3	MW 321.7307	C ₁₆ H ₁₃ ClFNO ₃	
DRE-C13649000	Flamprop (free acid)		100mg
			
Flamprop-M-isopropyl			
CAS 63782-90-1	MW 363.8105	C ₁₉ H ₁₉ ClFNO ₃	
DRE-C13650500	Flamprop-M-isopropyl		250mg
			
Flamprop-M-methyl			
CAS 63729-98-6	MW 335.7573	C ₁₇ H ₁₅ ClFNO ₃	
DRE-C13655200	Flamprop-M-methyl		10mg
			
Flamprop-isopropyl			
CAS 52756-22-6	MW 363.8105	C ₁₉ H ₁₉ ClFNO ₃	
DRE-C13650000	Flamprop-isopropyl		250mg
DRE-L13650000CY	Flamprop-isopropyl 10 µg/mL in Cyclohexane		10ml
			
Flamprop-methyl			
CAS 52756-25-9	MW 335.7573	C ₁₇ H ₁₅ ClFNO ₃	
DRE-C13655000	Flamprop-methyl		100mg
DRE-L13655000CY	Flamprop-methyl 10 µg/mL in Cyclohexane		10ml
			
Flazasulfuron			
CAS 104040-78-0	MW 407.3251	C ₁₃ H ₁₂ F ₃ N ₅ O ₅ S	
DRE-C13655600	Flazasulfuron(‡)		50mg
			










Pesticide and metabolites: Herbicides

Product code	Description			
Florasulam				
CAS 145701-23-1	MW 359.2838	C ₁₂ H ₈ F ₃ N ₅ O ₃ S		
DRE-C13662200	Florasulam(±)		50mg	
DRE-L13662200AL	Florasulam 10 µg/mL in Acetonitrile		10ml	
Fluazifop-P-butyl				
CAS 79241-46-6	MW 383.3616	C ₁₉ H ₂₀ F ₃ NO ₄		
DRE-C13670200	Fluazifop-P-butyl		100mg	
DRE-L13670200AC	Fluazifop-P-butyl 10 µg/mL in Acetone		10ml	
Fluazifop-butyl				
CAS 69806-50-4	MW 383.3616	C ₁₉ H ₂₀ F ₃ NO ₄		
DRE-C13670000	Fluazifop-butyl		250mg	
DRE-L13670000AL	Fluazifop-butyl 10 µg/mL in Acetonitrile		10ml	
DRE-XA13670000AL	Fluazifop-butyl 100 µg/mL in Acetonitrile		1ml	
Fluazifop-butyl D9 (n-butyl D9)				
CAS n/a	MW 392.4171	C ₁₉ ² H ₁₉ F ₃ NO ₄		
DRE-XA13670100AC	Fluazifop-butyl D9 100 µg/mL in Acetone		1.1ml	
Fluazifop-methyl				
CAS 69335-90-6	MW 341.2819	C ₁₆ H ₁₄ F ₃ NO ₄		
DRE-C13670400	Fluazifop-methyl		50mg	
Flucarbazone-sodium				
CAS 181274-17-9	MW 418.281	C ₁₂ H ₁₀ F ₃ N ₄ O ₆ S·Na		
DRE-C13685000	Flucarbazone-sodium		100mg	
Flucetosulfuron				
CAS 412928-75-7	MW 487.4594	C ₁₈ H ₂₂ FN ₅ O ₈ S		
DRE-C13687000	Flucetosulfuron		25mg	
Fluchloralin				
CAS 33245-39-5	MW 355.6975	C ₁₂ H ₁₃ ClF ₃ N ₃ O ₄		
DRE-C13690000	Fluchloralin(±)		100mg	
DRE-L13690000CY	Fluchloralin 10 µg/mL in Cyclohexane		10ml	
Flufenacet				
CAS 142459-58-3	MW 363.3305	C ₁₄ H ₁₃ F ₄ N ₃ O ₂ S		
DRE-C13711000	Flufenacet(±)		100mg	
DRE-L13711000AL	Flufenacet 10 µg/mL in Acetonitrile		10ml	
DRE-L13711000CY	Flufenacet 10 µg/mL in Cyclohexane		10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Flufenacet-alcohol				
CAS 54041-17-7 DRE-C13711010	MW 211.2328 Flufenacet-alcohol	C ₁₁ H ₁₄ FNO ₂	10mg	
Flufenacet ethane sulfonic acid ESA sodium salt				
CAS 947601-87-8 DRE-CA13711017	MW 297.2784 Flufenacet ethane sulfonic acid ESA sodium salt	C ₁₁ H ₁₃ FNO ₄ S·Na	10mg	
Flufenacet-methylsulfide				
CAS n/a DRE-C13711018	MW 241.3249 Flufenacet-methylsulfide	C ₁₂ H ₁₆ FNOS	50mg	
Flufenacet-methylsulfone				
CAS 861391-80-2 DRE-C13711019	MW 273.3237 Flufenacet-methylsulfone	C ₁₂ H ₁₆ FNO ₃ S	50mg	
Flufenacet-oxalat (OA)				
CAS 201668-31-7 DRE-CA13711020	MW 225.2163 Flufenacet-oxalat (OA)	C ₁₁ H ₁₂ FNO ₃	10mg	
Flufenacet-thioglycolate sulfoxide				
CAS 201668-33-9 DRE-C13711050	MW 301.3338 Flufenacet-thioglycolate sulfoxide	C ₁₃ H ₁₆ FNO ₄ S	10mg	
Flufenpyr-ethyl				
CAS 188489-07-8 DRE-C13713000	MW 408.732 Flufenpyr-ethyl	C ₁₆ H ₁₃ ClF ₄ N ₂ O ₄	10mg	
Flumetsulam				
CAS 98967-40-9 DRE-C13721000	MW 325.294 Flumetsulam	C ₁₂ H ₉ F ₂ N ₅ O ₂ S	100mg	
Flumiclorac-pentyl				
CAS 87546-18-7 DRE-C13724000	MW 423.8624 Flumiclorac-pentyl	C ₂₁ H ₂₃ ClFNO ₅	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Flumioxazin				
CAS 103361-09-7	MW 354.3318	C ₁₉ H ₁₅ FN ₂ O ₄		
DRE-C13725000	Flumioxazin(±)		100mg	
DRE-L13725000AL	Flumioxazin 10 µg/mL in Acetonitrile		10ml	
DRE-L13725000CY	Flumioxazin 10 µg/mL in Cyclohexane		10ml	
Fluometuron				
CAS 2164-17-2	MW 232.2023	C ₁₀ H ₁₁ F ₃ N ₂ O		
DRE-C13730000	Fluometuron		250mg	
DRE-L13730000AL	Fluometuron 10 µg/mL in Acetonitrile		10ml	
DRE-XA09010197AL	Fluometuron 100 µg/mL in Acetonitrile(±)		1ml	
Fluometuron-desmethyl (1-Methyl-3-[3-(trifluoromethyl)phenyl]urea)				
CAS 3032-40-4	MW 218.1758	C ₉ H ₉ F ₃ N ₂ O		
DRE-L13730500ME	Fluometuron-desmethyl 10 µg/mL in Methanol		10ml	
Fluorodifen				
CAS 15457-05-3	MW 328.2003	C ₁₃ H ₇ F ₃ N ₂ O ₅		
DRE-C13790000	Fluorodifen(±)		100mg	
Fluoroglycofen-ethyl				
CAS 77501-90-7	MW 447.7465	C ₁₈ H ₁₃ ClF ₃ NO ₇		
DRE-C13792000	Fluoroglycofen-ethyl		250mg	
DRE-XA09010198AL	Fluoglycofene-ethyl 100 µg/mL in Acetonitrile(±)		1ml	
Flupropanate				
CAS 756-09-2	MW 146.0404	C ₃ H ₂ F ₄ O ₂		
DRE-C13802000	Flupropanate		10mg	
Fluprimsulfuron-methyl sodium				
CAS 144740-54-5	MW 473.3164	C ₁₄ H ₁₁ F ₃ N ₅ O ₇ S·Na		
DRE-C13802500	Fluprimsulfuron-methyl sodium(±)		10mg	
Flurenol				
CAS 467-69-6	MW 226.2274	C ₁₄ H ₁₀ O ₃		
DRE-C13810000	Flurenol		250mg	
Flurenol Butyl Ester				
CAS 2314-09-2	MW 282.3337	C ₁₈ H ₁₈ O ₃		
DRE-C13812000	Flurenol-butyl ester		100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Flurenol Methyl Ester				
CAS 1216-44-0 DRE-C13815000	MW 240.254 Flurenol-methyl ester	C ₁₅ H ₁₂ O ₃	100mg	
Fluridone				
CAS 59756-60-4 DRE-C13840000 DRE-XA09010199ME	MW 329.3158 Fluridone Fluridone 100 µg/mL in Methanol(±)	C ₁₉ H ₁₄ F ₃ NO	100mg 1ml	
Flurochloridone				
CAS 61213-25-0 DRE-C13847500 DRE-L13847500CY DRE-XA13847500AL	MW 312.1151 Flurochloridone Flurochloridone 10 µg/mL in Cyclohexane Flurochloridone 100 µg/mL in Acetonitrile	C ₁₂ H ₁₀ Cl ₂ F ₃ NO	100mg 10ml 1ml	
Fluroxypyr				
CAS 69377-81-7 DRE-C13849000 DRE-L13849000AL DRE-XA13849000AL	MW 255.0306 Fluroxypyr Fluroxypyr 10 µg/mL in Acetonitrile Fluroxypyr 100 µg/mL in Acetonitrile	C ₇ H ₅ Cl ₂ FN ₂ O ₃	100mg 10ml 1ml	
Fluroxypyr-1-methylheptylester				
CAS 81406-37-3 DRE-C13850000	MW 367.2432 Fluroxypyr-1-methylheptylester	C ₁₅ H ₂₁ Cl ₂ FN ₂ O ₃	100mg	
Flurtamone				
CAS 96525-23-4 DRE-C13855000	MW 333.3045 Flurtamone	C ₁₈ H ₁₄ F ₃ NO ₂	100mg	
Fluthiacet-methyl				
CAS 117337-19-6 DRE-C13862000	MW 403.8793 Fluthiacet-methyl	C ₁₅ H ₁₅ ClF ₃ N ₃ O ₃ S ₂	100mg	
Fomesafen				
CAS 72178-02-0 DRE-C13895000 DRE-L13895000AL	MW 438.7629 Fomesafen Fomesafen 10 µg/mL in Acetonitrile	C ₁₅ H ₁₀ ClF ₃ N ₂ O ₆ S	100mg 10ml	
Foramsulfuron				
CAS 173159-57-4 DRE-C13905000	MW 452.4417 Foramsulfuron	C ₁₇ H ₂₀ N ₆ O ₇ S	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Fosamine Ammonium Salt				
CAS 25954-13-6 DRE-C13930000	MW 170.1042 Fosamine-ammonium	C ₃ H ₇ NO ₄ P·H ₄ N	250mg	
Glufosinate-FMOC				
CAS 1822429-60-6 DRE-CA14031000 DRE-L14031000AL	MW 403.3655 Glufosinate-FMOC Glufosinate-FMOC 10 µg/mL in Acetonitrile (*)	C ₂₀ H ₂₂ NO ₅ P	10mg 10ml	
Glufosinate-N-acetyl				
CAS 73634-73-8 DRE-CA14031500 DRE-L14031500ME	MW 223.1635 Glufosinate-N-acetyl Glufosinate-N-acetyl 10 µg/mL in Methanol	C ₇ H ₁₄ NO ₅ P	10mg 10ml	
Glufosinate ammonium				
CAS 77182-82-2 DRE-CA14030000 DRE-L14030000WA DRE-XA14030000WA	MW 198.1574 Glufosinate ammonium(±) Glufosinate ammonium 10 µg/mL in Water Glufosinate ammonium 100 µg/mL in Water	C ₅ H ₁₁ NO ₄ P·H ₄ N	100mg 10ml 1ml	
Glyphosate (N-(Phosphonomethyl)glycine)				
CAS 1071-83-6 DRE-C14050000 DRE-L14050000WA DRE-XA14050000WA DRE-GS09010310WA DRE-XA09010248WA	MW 169.0731 Glyphosate(±) Glyphosate 10 µg/mL in Water Glyphosate 100 µg/mL in Water Glyphosate (n-(phosphonomethyl)glycine) 100 µg/mL in Water(±) Glyphosate (n-(phosphonomethyl)glycine) 100 µg/mL in Water(±)	C ₃ H ₈ NO ₅ P	250mg 10ml 1ml 5x1ml 1ml	
Glyphosate 1,2-13C2 15N				
CAS 1185107-63-4 DRE-XA14050100WA	MW 172.0518 Glyphosate 1,2-13C2 15N 100 µg/mL in Water(±)	¹³ C ₂ H ₈ ¹⁵ NO ₅ P	1ml	
Glyphosate-FMOC				
CAS 1373205-41-4 DRE-CA14051000	MW 391.3118 Glyphosate-FMOC	C ₁₈ H ₁₈ NO ₇ P	50mg	
Glyphosate-N-methyl				
CAS 24569-83-3 DRE-C14055300	MW 183.0997 Glyphosate-N-methyl	C ₄ H ₁₀ NO ₅ P	10mg	
Glyphosate-N-nitroso mono sodium salt				
CAS 56516-71-3 DRE-C14055400 DRE-XA14055400WA	MW 220.0531 Glyphosate-N-nitroso mono sodium salt(±) Glyphosate-N-nitroso mono sodium salt 100 µg/mL in Water (*)	C ₃ H ₈ N ₂ O ₅ P·Na	10mg 1ml	









Pesticide and metabolites: Herbicides

Product code	Description			
Halosulfuron-methyl				
CAS 100784-20-1 DRE-C14059500	MW 434.8122 Halosulfuron-methyl	C ₁₃ H ₁₅ ClN ₆ O ₇ S	50mg	
Haloxypop (free acid)				
CAS 69806-34-4 DRE-C14060000 DRE-L14060000AL	MW 361.7003 Haloxypop (free acid) Haloxypop (free acid) 10 µg/mL in Acetonitrile	C ₁₉ H ₁₁ ClF ₃ NO ₄	100mg 10ml	
Haloxypop-R-methyl				
CAS 72619-32-0 DRE-C14062500 DRE-L14062500AL	MW 375.7269 Haloxypop-R-methyl(±) Haloxypop-R-methyl 10 µg/mL in Acetonitrile	C ₁₆ H ₁₃ ClF ₃ NO ₄	100mg 10ml	
Haloxypop-ethoxyethyl				
CAS 87237-48-7 DRE-C14061000 DRE-L14061000IO	MW 433.8061 Haloxypop-ethoxyethyl Haloxypop-ethoxyethyl 10 µg/mL in Isooctane	C ₁₉ H ₁₉ ClF ₃ NO ₅	100mg 10ml	
Haloxypop-methyl				
CAS 69806-40-2 DRE-C14062000 DRE-L14062000IO	MW 375.7269 Haloxypop-methyl Haloxypop-methyl 10 µg/mL in Isooctane	C ₁₆ H ₁₃ ClF ₃ NO ₄	100mg 10ml	
Hexazinone (Velpar)				
CAS 51235-04-2 DRE-C14200000 DRE-L14200000AC DRE-L14200000AL DRE-XA14200000AL DRE-XA09010202AL	MW 252.3128 Hexazinone Hexazinone 10 µg/mL in Acetone Hexazinone 10 µg/mL in Acetonitrile Hexazinone 100 µg/mL in Acetonitrile Velpar (Hexazinone) 100 µg/mL in Acetonitrile(±)	C ₁₂ H ₂₀ N ₄ O ₂	100mg 10ml 10ml 1ml 1ml	
Hexazinone D6 (N,N-dimethyl D6)				
CAS 1219804-22-4 DRE-XA14200010AL	MW 258.3498 Hexazinone D6 (N,N-dimethyl D6) 100 µg/mL in Acetonitrile	C ₁₂ ² H ₆ ¹⁴ N ₄ O ₂	1ml	
Hydroxymethyl phosphonic acid				
CAS 2617-47-2 DRE-CA14233050	MW 112.0218 Hydroxymethyl phosphonic acid	CH ₅ O ₄ P	100mg	
Imazamethabenz (free acid)				
CAS 100728-84-5 DRE-XA14281400AL	MW 274.315 Imazamethabenz (free acid) 100 µg/mL in Acetonitrile	C ₁₅ H ₁₆ N ₂ O ₃	1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Imazamethabenz-methyl				
CAS 81405-85-8	MW 576.6832	2C ₁₆ H ₂₀ N ₂ O ₃		
DRE-C14281500	Imazamethabenz-methyl		100mg	
DRE-XA14281500AL	Imazamethabenz-methyl 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010203AL	Imazamethabenz-methyl 100 µg/mL in Acetonitrile(±)		1ml	
Imazamox				
CAS 114311-32-9	MW 305.3291	C ₁₅ H ₁₉ N ₃ O ₄		
DRE-C14282000	Imazamox		100mg	
DRE-L14282000AL	Imazamox 10 µg/mL in Acetonitrile		10ml	
Imazapic				
CAS 104098-48-8	MW 275.3031	C ₁₄ H ₁₇ N ₃ O ₃		
DRE-C14282500	Imazapic		100mg	
Imazapyr				
CAS 81334-34-1	MW 261.2765	C ₁₃ H ₁₅ N ₃ O ₃		
DRE-C14283000	Imazapyr(±)		100mg	
DRE-L14283000AL	Imazapyr 10 µg/mL in Acetonitrile		10ml	
Imazaquin				
CAS 81335-37-7	MW 311.3352	C ₁₇ H ₁₇ N ₃ O ₃		
DRE-C14283300	Imazaquin		100mg	
Imazethapyr				
CAS 81335-77-5	MW 289.3297	C ₁₅ H ₁₉ N ₃ O ₃		
DRE-C14283500	Imazethapyr		100mg	
DRE-L14283500AL	Imazethapyr 10 µg/mL in Acetonitrile		10ml	
Imazosulfuron				
CAS 122548-33-8	MW 412.8082	C ₁₄ H ₁₃ ClN ₆ O ₅ S		
DRE-C14283550	Imazosulfuron		100mg	
Inabenfide				
CAS 82211-24-3	MW 338.7876	C ₁₉ H ₁₅ ClN ₂ O ₂		
DRE-C14287000	Inabenfide(±)		50mg	
Indanofan				
CAS 133220-30-1	MW 340.8002	C ₂₀ H ₁₇ ClO ₃		
DRE-C14288000	Indanofan		100mg	
DRE-L14288000CY	Indanofan 10 µg/mL in Cyclohexane		10ml	









Pesticide and metabolites: Herbicides

Product code	Description			
Indaziflam				
CAS 950782-86-2 DRE-C14288300	MW 301.3619 Indaziflam	C ₁₆ H ₂₀ FN ₅	100mg	
Iodosulfuron-methyl sodium				
CAS 144550-36-7 DRE-C14343500	MW 529.2422 Iodosulfuron-methyl sodium	C ₁₄ H ₁₃ IN ₅ O ₆ S·Na	100mg	
Ioxynil				
CAS 1689-83-4 DRE-C14350000 DRE-L14350000AL	MW 370.9138 Ioxynil Ioxynil 10 µg/mL in Acetonitrile	C ₇ H ₅ I ₂ NO	100mg 10ml	
Ioxynil-methyl				
CAS 3336-40-1 DRE-C14355000	MW 384.9403 Ioxynil-methyl	C ₈ H ₅ I ₂ NO	100mg	
Ioxynil-octanoate				
CAS 3861-47-0 DRE-C14360000 DRE-L14360000IO	MW 497.1099 Ioxynil-octanoate Ioxynil-octanoate 10 µg/mL in Isooctane	C ₁₅ H ₁₇ I ₂ NO ₂	100mg 10ml	
Ipfencarbazone				
CAS 212201-70-2 DRE-C14366000	MW 427.2322 Ipfencarbazone	C ₁₈ H ₁₄ Cl ₂ F ₂ N ₄ O ₂	25mg	
Isocarbamid				
CAS 30979-48-7 DRE-C14400000	MW 185.2236 Isocarbamid	C ₈ H ₁₅ N ₃ O ₂	100mg	
Isochloridazon				
CAS 162354-96-3 DRE-L11321000EA	MW n/a iso-Chloridazon 10 µg/mL in Ethylacetate		10ml	No Structure
Isomethiozin				
CAS 57052-04-7 DRE-C14430000	MW 268.3784 Isomethiozin	C ₁₂ H ₂₀ N ₄ OS	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Isonoruron				
CAS 28805-78-9 DRE-C14440000	MW 444.6532 Isonoruron	$((C_{13}H_{22}N_2O)c(C_{13}H_{22}N_2O)c)_{mix}$	100mg	
Isopropalin				
CAS 33820-53-0 DRE-C14460000	MW 309.3608 Isopropalin(‡)	$C_{15}H_{23}N_3O_4$	100mg	
1-(4-Isopropylphenyl)-3-methylurea				
CAS 34123-57-4 DRE-C14464800 DRE-XA14464800AL	MW 192.2575 1-(4-Isopropylphenyl)-3-methylurea 1-(4-Isopropylphenyl)-3-methylurea 100 µg/mL in Acetonitrile	$C_{11}H_{16}N_2O$	100mg 1ml	
1-(4-Isopropylphenyl)urea				
CAS 56046-17-4 DRE-C14465000 DRE-L14465000AL	MW 178.231 1-(4-Isopropylphenyl)urea(‡) 1-(4-Isopropylphenyl)urea 10 µg/mL in Acetonitrile	$C_{10}H_{14}N_2O$	100mg 10ml	
Isoproturon				
CAS 34123-59-6 DRE-C14470000 DRE-L14470000AL DRE-L14470000EA DRE-XA14470000AL	MW 206.2841 Isoproturon Isoproturon 10 µg/mL in Acetonitrile Isoproturon 10 µg/mL in Ethylacetate Isoproturon 100 µg/mL in Acetonitrile	$C_{12}H_{18}N_2O$	100mg 10ml 10ml 1ml	
Isoproturon D6				
CAS 1007461-76-8 DRE-C14470100 DRE-XA14470100AC DRE-YA14470100AL	MW 212.3211 Isoproturon D6 Isoproturon D6 100 µg/mL in Acetone(‡) Isoproturon D6 1000 µg/mL in Acetonitrile	$C_{12}^2H_6H_{12}N_2O$	10mg 1.1ml 1.1ml	
Isouron				
CAS 55861-78-4 DRE-C14479000	MW 211.2609 Isouron	$C_{10}H_{17}N_3O_2$	100mg	
Isoxaben				
CAS 82558-50-7 DRE-C14480000 DRE-XA09010206ME	MW 332.3942 Isoxaben Isoxaben 100 µg/mL in Methanol(‡)	$C_{18}H_{24}N_2O_4$	100mg 1ml	
Isoxaflutole				
CAS 141112-29-0 DRE-C14481000 DRE-L14481000CY	MW 359.3203 Isoxaflutole Isoxaflutole 10 µg/mL in Cyclohexane	$C_{15}H_{12}F_3NO_4S$	100mg 10ml	

Pesticide and metabolites: Herbicides

Product code	Description			
Karbutilate				
CAS 4849-32-5	MW 279.3348	C ₁₄ H ₂₁ N ₃ O ₃		
DRE-C14510000	Karbutilate		250mg	
DRE-GS09010095ME	Karbutilate 100 µg/mL in Methanol(±)		5x1ml	
Lactofen				
CAS 77501-63-4	MW 461.7731	C ₁₉ H ₁₅ ClF ₃ NO ₇		
DRE-C14590000	Lactofen		100mg	
Lenacil				
CAS 2164-08-1	MW 234.2942	C ₁₃ H ₁₈ N ₂ O ₂		
DRE-C14610000	Lenacil		100mg	
DRE-L14610000AL	Lenacil 10 µg/mL in Acetonitrile		10ml	
DRE-XA14610000AL	Lenacil 100 µg/mL in Acetonitrile		1ml	
Lenacil (2,2,2',2'-cyclohexyl D4)				
CAS n/a	MW 238.3189	C ₁₃ ² H ₁₄ H ₄ N ₂ O ₂		
DRE-XA14610100AL	Lenacil (cyclohexyl D4) 100 µg/mL in Acetonitrile		1.1ml	
Linuron				
CAS 330-55-2	MW 249.0939	C ₉ H ₁₀ Cl ₂ N ₂ O ₂		
DRE-C14640000	Linuron		100mg	
DRE-L14640000AL	Linuron 10 µg/mL in Acetonitrile		10ml	
DRE-XA14640000AL	Linuron 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010263AC	Linuron 100 µg/mL in Acetone(±)		1ml	
Linuron D6				
CAS 1219804-76-8	MW 255.1309	C ₉ ² H ₆ H ₄ Cl ₂ N ₂ O ₂		
DRE-C14640100	Linuron D6 (methyl D3 methoxy D3)(±)		10mg	
DRE-XA14640100AC	Linuron D6 100 µg/mL in Acetone(±)		1.1ml	
MCPA ((4-Chloro-2-methylphenoxy)acetic Acid)				
CAS 94-74-6	MW 200.619	C ₉ H ₉ ClO ₃		
DRE-C14760000	MCPA		250mg	
DRE-L14760000AL	MCPA 10 µg/mL in Acetonitrile		10ml	
DRE-XA14760000AL	MCPA 100 µg/mL in Acetonitrile		1ml	
DRE-YA09010001MB	MCPA Standard 2000 µg/mL in Methyl tert-Butyl Ether		1ml	
DRE-YS09010001MB	MCPA Standard 2000 µg/mL in Methyl tert-Butyl Ether(±)		5x1ml	
MCPA Butoxyethyl Ester				
CAS 19480-43-4	MW 300.7778	C ₁₅ H ₂₁ ClO ₄		
DRE-C14762000	MCPA-butoxyethyl ester		100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
MCPA Butyl Ester				
CAS 1713-12-8 DRE-C14763000	MW 256.7253 MCPA-1-butyl ester	$C_{13}H_{17}ClO_3$	100mg	
MCPA D3 (phenyl D3)				
CAS 352431-14-2 DRE-C14760100 DRE-XA14760100AC	MW 203.6374 MCPA D3 MCPA D3 100 µg/mL in Acetone	$C_9^2H_6H_3ClO_3$	10mg 1.1ml	
MCPA D6 (methyl-D3,phenoxy-D3)				
CAS n/a DRE-XA14760200AC	MW 206.6559 MCPA D6 100 µg/mL in Acetone	$C_9^2H_6H_3ClO_3$	1.1ml	
MCPA Ethyl Ester				
CAS 2698-38-6 DRE-C14763800	MW 228.6721 MCPA-ethyl ester	$C_{11}H_{13}ClO_3$	100mg	
MCPA-dimethylammonium				
CAS 2039-46-5 DRE-C14763500	MW 245.7026 MCPA-dimethylammonium	$C_9H_9ClO_3 \cdot C_2H_7N$	100mg	
MCPA-2-ethylhexyl ester				
CAS 29450-45-1 DRE-C14764000	MW 312.8316 MCPA-2-ethylhexyl ester	$C_{17}H_{25}ClO_3$	100mg	
MCPA-methyl ester				
CAS 2436-73-9 DRE-C14768000 DRE-L14768000CY	MW 214.6455 MCPA-methyl ester MCPA-methyl ester 10 µg/mL in Cyclohexane	$C_{10}H_{11}ClO_3$	100mg 10ml	
MCPA-sodium salt monohydrate				
CAS 3653-48-3 DRE-C14775000	MW 240.6161 MCPA-sodium salt monohydrate	$C_9H_8ClO_3 \cdot Na \cdot H_2O$	250mg	
MCPA-thioethyl				
CAS 25319-90-8 DRE-C14776000	MW 244.7377 MCPA-thioethyl	$C_{11}H_{13}ClO_2S$	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
MCPB (4-(4-Chloro-2-methylphenoxy)butanoic Acid)				
CAS 94-81-5	MW 228.6721	C ₁₁ H ₁₃ ClO ₃		
DRE-C14790000	MCPB		100mg	
DRE-L14790000AL	MCPB 10 µg/mL in Acetonitrile		10ml	
DRE-XA14790000AL	MCPB 100 µg/mL in Acetonitrile		1ml	
MCPB D6 (ring D3, methyl D3) (4-(4-Chloro-2-trideuteriomethyl-3,5,6-trideuteriophenoxy)butanoic Acid)				
CAS n/a	MW 234.7091	C ₁₁ ² H ₆ H ₇ ClO ₃		
DRE-C14790100	MCPB D6 (ring D3, methyl D3)		10mg	
DRE-XA14790100AC	MCPB D6 (ring D3, methyl D3) 100 µg/mL in Acetone		1.1ml	
MCPB Ethyl Ester				
CAS 10443-70-6	MW 256.7253	C ₁₃ H ₁₇ ClO ₃		
DRE-C14794500	MCPB-ethyl ester		100mg	
MCPB-ethylhexylester				
CAS 94232-74-3	MW 340.8848	C ₁₉ H ₂₅ ClO ₃		
DRE-C14794700	MCPB-ethylhexylester		100mg	
MCPB-methyl ester				
CAS 57153-18-1	MW 242.6987	C ₁₂ H ₁₅ ClO ₃		
DRE-C14795000	MCPB-methyl ester		100mg	
DRE-L14795000CY	MCPB-methyl ester 10 µg/mL in Cyclohexane		10ml	
Mecoprop (MCP)				
CAS 93-65-2	MW 214.6455	C ₁₀ H ₁₁ ClO ₃		
DRE-C14820000	Mecoprop		100mg	
DRE-L14820000AL	Mecoprop 10 µg/mL in Acetonitrile		10ml	
DRE-XA14820000AL	Mecoprop 100 µg/mL in Acetonitrile		1ml	
DRE-YA09010002MB	MCP (Mecoprop) 2000 µg/mL in Methyl tert-Butyl Ether		1ml	
DRE-YS09010002MB	MCP (Mecoprop) 2000 µg/mL in Methyl tert-Butyl Ether(‡)		5x1ml	
Mecoprop D3 (phenyl D3)				
CAS 352431-15-3	MW 217.664	C ₁₀ ² H ₅ H ₆ ClO ₃		
DRE-XA14820100AC	Mecoprop D3 100 µg/mL in Acetone		1.1ml	
Mecoprop D6 (ring D3, methyl D3)				
CAS 1705649-54-2	MW 220.6825	C ₁₀ ² H ₆ H ₅ ClO ₃		
DRE-C14820110	Mecoprop D6 (ring D3, methyl D3)		10mg	
DRE-XA14820110AL	Mecoprop D6 (ring D3, methyl D3) 100 µg/mL in Acetonitrile		1.1ml	
Mecoprop-2-butoxyethyl ester				
CAS 23359-62-8	MW 314.8044	C ₁₆ H ₂₃ ClO ₄		
DRE-C14825000	Mecoprop-2-butoxyethyl ester		100mg	

Pesticide and metabolites: Herbicides

Product code	Description			
Mecoprop-n/iso-butyl ester				
CAS n/a	MW 541.5037	2C ₁₄ H ₁₉ ClO ₃		
DRE-C14833500	Mecoprop-n/iso-butyl ester		100mg	
Mecoprop methyl Ester				
CAS 23844-56-6	MW 228.6721	C ₁₁ H ₁₃ ClO ₃		
DRE-C14835000	Mecoprop-methyl ester		100mg	
DRE-L14835000IO	Mecoprop-methyl ester 10 µg/mL in Isooctane		10ml	
Mecoprop-methyl ester D3 (ring D3)				
CAS n/a	MW 231.6906	C ₁₁ ² H ₃ H ₁₀ ClO ₃		
DRE-XA14835100AC	Mecoprop-methyl ester D3 (ring) 100 µg/mL in Acetone		1.1ml	
Mecoprop 2-Octyl Ester				
CAS 28473-03-2	MW 326.8582	C ₁₈ H ₂₇ ClO ₃		
DRE-C14838200	Mecoprop-2-octyl ester		100mg	
Mecoprop-2,4,4-trimethylpentylester				
CAS 217487-13-3	MW 326.8582	C ₁₈ H ₂₇ ClO ₃		
DRE-C14840000	Mecoprop-2,4,4-trimethylpentylester		100mg	
2,6-Mecoprop				
CAS 35851-12-8	MW 214.6455	C ₁₀ H ₁₁ ClO ₃		
DRE-C14843000	2,6-Mecoprop		10mg	
D(+)-Mecoprop				
CAS 16484-77-8	MW 214.6455	C ₁₀ H ₁₁ ClO ₃		
DRE-C14820200	D(+)-Mecoprop		250mg	
Mecoprop-2-ethylhexyl ester				
CAS 71526-69-7	MW 326.8582	C ₁₈ H ₂₇ ClO ₃		
DRE-C14830000	Mecoprop-2-ethylhexyl ester		100mg	
Medinoterb Acetate				
CAS 2487-01-6	MW 296.2759	C ₁₃ H ₁₆ N ₂ O ₆		
DRE-C14851000	Medinoterb acetate		100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Mefenacet				
CAS 73250-68-7	MW 298.3596	C ₁₆ H ₁₄ N ₂ O ₂ S		
DRE-C14860000	Mefenacet		100mg	
DRE-L14860000IO	Mefenacet 10 µg/mL in Isooctane		10ml	
Mefluidide				
CAS 53780-34-0	MW 310.2927	C ₁₁ H ₁₃ F ₃ N ₂ O ₃ S		
DRE-C14861000	Mefluidide		25mg	
DRE-L14861000AL	Mefluidide 10 µg/mL in Acetonitrile		10ml	
Mesosulfuron-methyl				
CAS 208465-21-8	MW 503.5067	C ₁₇ H ₂₁ N ₅ O ₉ S ₂		
DRE-C14913500	Mesosulfuron-methyl(±)		100mg	
Mesotrione				
CAS 104206-82-8	MW 339.3205	C ₁₄ H ₁₃ NO ₇ S		
DRE-C14914000	Mesotrione(±)		100mg	
Metam-sodium dihydrate				
CAS 6734-80-1	MW 165.2102	C ₂ H ₄ NS ₂ ·Na·2H ₂ O		
DRE-C14935000	Metam-sodium dihydrate (*)		250mg	
(±)-Metamifop				
CAS 256412-89-2	MW 440.8514	C ₂₃ H ₁₈ ClFN ₂ O ₄		
DRE-C14938000	Metamifop		25mg	
Metamitron				
CAS 41394-05-2	MW 202.2126	C ₁₀ H ₁₀ N ₄ O		
DRE-C14940000	Metamitron		100mg	
DRE-L14940000AL	Metamitron 10 µg/mL in Acetonitrile		10ml	
DRE-XA14940000AL	Metamitron 100 µg/mL in Acetonitrile		1ml	
Metamitron-desamino				
CAS 36993-94-9	MW 187.198	C ₁₀ H ₉ N ₃ O		
DRE-C14940300	Metamitron-desamino		10mg	
DRE-XA14940300AL	Metamitron-desamino 100 µg/mL in Acetonitrile		1ml	
Metazachlor				
CAS 67129-08-2	MW 277.7493	C ₁₄ H ₁₆ ClN ₃ O		
DRE-C14950000	Metazachlor		100mg	
DRE-L14950000AL	Metazachlor 10 µg/mL in Acetonitrile		10ml	
DRE-XA14950000AC	Metazachlor 100 µg/mL in Acetone		1ml	
DRE-XA14950000AL	Metazachlor 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010209ME	Metazachlor 100 µg/mL in Methanol(±)		1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Metazachlor Ethane Sulfonic Acid (2-[(2,6-Dimethylphenyl)(1H-pyrazol-1-ylmethyl)amino]-2-oxoethanesulfonic acid)				
CAS 172960-62-2 DRE-CA14950020	MW 323.3675 Metazachlor ethane sulfonic acid (ESA)	C ₁₄ H ₁₇ N ₃ O ₄ S	10mg	
Metazachlor oxanilic acid (OA)				
CAS 1231244-60-2 DRE-CA14950050	MW 273.2872 Metazachlor oxanilic acid (OA)	C ₁₄ H ₁₅ N ₃ O ₃	10mg	
Methabenzthiazuron				
CAS 18691-97-9 DRE-C14960000 DRE-L14960000AL DRE-XA14960000AL	MW 221.2788 Methabenzthiazuron(‡) Methabenzthiazuron 10 µg/mL in Acetonitrile Methabenzthiazuron 100 µg/mL in Acetonitrile	C ₁₀ H ₁₁ N ₃ OS	100mg 10ml 1ml	
Methoprotryne				
CAS 841-06-5 DRE-C15050000	MW 271.3823 Methoprotryne	C ₁₁ H ₂₁ N ₅ OS	100mg	
Methoxyphenone (3,3'-Dimethyl-4-methoxybenzophenone)				
CAS 41295-28-7 DRE-C15081500	MW 240.297 Methoxyphenon	C ₁₆ H ₁₆ O ₂	100mg	
Methyl N-(3,4-Dichlorophenyl)carbamate (Sweep)				
CAS 1918-18-9 DRE-C17060000	MW 220.0527 Sweep	C ₈ H ₇ Cl ₂ NO ₂	250mg	
Methylarsonic acid mono sodium salt sesquihydrat				
CAS 2163-80-6 DRE-C15083775	MW 341.9193 Methylarsonic acid mono sodium salt sesquihydrat	2CH ₄ AsO ₃ ·2Na·H ₂ O	100mg	
Methyldymron				
CAS 42609-73-4 DRE-C15086020	MW 268.3535 Methyldymron	C ₁₇ H ₂₀ N ₂ O	10mg	
Metobromuron				
CAS 3060-89-7 DRE-C15160000 DRE-L15160000AL DRE-XA15160000AL	MW 259.0998 Metobromuron(‡) Metobromuron 10 µg/mL in Acetonitrile Metobromuron 100 µg/mL in Acetonitrile	C ₉ H ₁₁ BrN ₂ O ₂	100mg 10ml 1ml	









Pesticide and metabolites: Herbicides

Product code	Description			
Metobromuron D6 (methyl D3 methoxy D3)				
CAS n/a	MW 265.1368	C ₉ ² H ₆ H ₅ BrN ₂ O ₂		
DRE-XA15160100AC	Metobromuron D6 (methyl D3 methoxy D3) 100 µg/mL in Acetone			1.1ml
				
Metolachlor				
CAS 51218-45-2	MW 283.7937	C ₁₅ H ₂₂ ClNO ₂		
DRE-C15170000	Metolachlor(‡)			100mg
DRE-L15170000AL	Metolachlor 10 µg/mL in Acetonitrile			10ml
DRE-L15170000CY	Metolachlor 10 µg/mL in Cyclohexane			10ml
DRE-XA15170000AL	Metolachlor 100 µg/mL in Acetonitrile			1ml
DRE-XA15170000CY	Metolachlor 100 µg/mL in Cyclohexane			1ml
				
Metolachlor D6 (propyl D6)				
CAS 1219803-97-0	MW 289.8307	C ₁₅ ² H ₆ H ₁₆ ClNO ₂		
DRE-CA15170100	Metolachlor D6 (propyl D6)			10mg
DRE-XA15170100AC	Metolachlor D6 100 µg/mL in Acetone			1.1ml
				
Metolachlor Ethane Sulfonic Acid (ESA) Sodium Salt				
CAS 947601-85-6	MW 351.3937	C ₁₅ H ₂₂ NO ₅ S·Na		
DRE-CA15171100	Metolachlor ethane sulfonic acid (ESA) sodium salt			10mg
				
Metolachlor deschloro				
CAS 126605-22-9	MW 249.3486	C ₁₅ H ₂₃ NO ₂		
DRE-XA15171050AL	Metolachlor deschloro 100 µg/mL in Acetonitrile			1ml
				
Metolachlor-2-ethoxy				
CAS 68544-97-8	MW 293.4012	C ₁₇ H ₂₇ NO ₃		
DRE-XA15171150AL	Metolachlor-2-ethoxy 100 µg/mL in Acetonitrile			1ml
				
Metolachlor-2-hydroxy				
CAS 131068-72-9	MW 265.348	C ₁₅ H ₂₃ NO ₃		
DRE-XA15171170AL	Metolachlor-2-hydroxy 100 µg/mL in Acetonitrile			1ml
				
Metolachlor mercapturate				
CAS 159956-64-6	MW 410.5276	C ₂₀ H ₃₀ N ₂ O ₅ S		
DRE-C15171190	Metolachlor mercapturate			100mg
				
Metolachlor oxanilic acid (OA)				
CAS 152019-73-3	MW 279.3315	C ₁₅ H ₂₁ NO ₄		
DRE-C15171200	Metolachlor oxanilic acid (OA)(‡)			10mg
				









Pesticide and metabolites: Herbicides

Product code	Description			
S-Metolachlor				
CAS 87392-12-9 DRE-C15171000 DRE-L15171000CY	MW 283.7937 S-Metolachlor(±) S-Metolachlor 10 µg/mL in Cyclohexane	C ₁₅ H ₂₂ ClNO ₂	100mg 10ml	
S-Metolachlor CGA 357704				
CAS 1217465-10-5 DRE-C15171020	MW 279.2885 S-Metolachlor CGA 357704	C ₁₄ H ₁₇ NO ₅	10mg	
S-Metolachlor CGA 368208				
CAS 1173021-76-5 DRE-C15171022	MW 279.2879 S-Metolachlor CGA 368208	C ₁₁ H ₁₄ NO ₄ S·Na	10mg	
S-Metolachlor CGA 37735				
CAS 97055-05-5 DRE-C15171024	MW 193.2423 S-Metolachlor CGA 37735	C ₁₁ H ₁₅ NO ₂	10mg	
S-Metolachlor CGA 50267				
CAS 82508-03-0 DRE-C15171026	MW 207.2689 S-Metolachlor CGA 50267	C ₁₂ H ₁₇ NO ₂	10mg	
S-Metolachlor CGA 50720				
CAS 152019-74-4 DRE-C15171028	MW 207.2258 S-Metolachlor CGA 50720	C ₁₁ H ₁₃ NO ₃	10mg	
S-Metolachlor NOA 413173				
CAS n/a DRE-C15171030	MW 373.3324 S-Metolachlor NOA 413173	C ₁₄ H ₁₇ NO ₆ S·2Na	10mg	
Metosulam				
CAS 139528-85-1 DRE-C15177500 DRE-L15177500AL	MW 418.2551 Metosulam Metosulam 10 µg/mL in Acetonitrile	C ₁₄ H ₁₃ Cl ₂ N ₅ O ₄ S	100mg 10ml	
Metoxuron				
CAS 19937-59-8 DRE-C15180000 DRE-L15180000AL DRE-XA15180000AL	MW 228.6754 Metoxuron Metoxuron 10 µg/mL in Acetonitrile Metoxuron 100 µg/mL in Acetonitrile	C ₁₀ H ₁₃ ClN ₂ O ₂	100mg 10ml 1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Metoxuron-monomethyl				
CAS 20782-57-4	MW 214.6488	C ₉ H ₁₁ ClN ₂ O ₂		
DRE-LA15185000EA	Metoxuron-monomethyl 10 µg/mL in Ethylacetate		1ml	
Metribuzin				
CAS 21087-64-9	MW 214.288	C ₈ H ₁₄ N ₄ OS		
DRE-C15200000	Metribuzin		100mg	
DRE-L15200000AL	Metribuzin 10 µg/mL in Acetonitrile		10ml	
DRE-L15200000CY	Metribuzin 10 µg/mL in Cyclohexane		10ml	
DRE-XA15200000AL	Metribuzin 100 µg/mL in Acetonitrile		1ml	
DRE-GS09010037ME	Metribuzin 200 µg/mL in Methanol(±)		4x1ml	
DRE-GS09010299AC	Metribuzin 500 µg/mL in Acetone(±)		5x1ml	
Metribuzin-desamino				
CAS 35045-02-4	MW 199.2733	C ₈ H ₁₃ N ₃ OS		
DRE-C15200300	Metribuzin-desamino		5mg	
DRE-LA15200300AL	Metribuzin-desamino 10 µg/mL in Acetonitrile		1ml	
Metribuzin-desamino-diketo				
CAS 52236-30-3	MW 169.1811	C ₇ H ₁₁ N ₃ O ₂		
DRE-C15200500	Metribuzin-desamino-diketo		5mg	
DRE-LA15200500AL	Metribuzin-desamino-diketo 10 µg/mL in Acetonitrile		1ml	
Metribuzin-diketo				
CAS 56507-37-0	MW 184.1958	C ₇ H ₁₂ N ₄ O ₂		
DRE-C15200700	Metribuzin-diketo		5mg	
DRE-LA15200700AL	Metribuzin-diketo 10 µg/mL in Acetonitrile		1ml	
Metsulfuron-methyl				
CAS 74223-64-6	MW 381.3638	C ₁₄ H ₁₅ N ₅ O ₆ S		
DRE-C15210000	Metsulfuron-methyl		100mg	
Metsulfuron-methyl D3 (triazine methoxy D3)				
CAS n/a	MW 384.3823	C ₁₄ ² H ₃ H ₁₂ N ₅ O ₆ S		
DRE-C15210100	Metsulfuron-methyl D3 (triazine methoxy D3)		10mg	
DRE-XA15210100AC	Metsulfuron-methyl D3 (triazine methoxy D3) 100 µg/mL in Acetone		1.1ml	
Molinate				
CAS 2212-67-1	MW 187.3024	C ₉ H ₁₇ NOS		
DRE-C15280000	Molinate		100mg	
DRE-L15280000CY	Molinate 10 µg/mL in Cyclohexane		10ml	
DRE-XA15280000CY	Molinate 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010212AL	Molinate 100 µg/mL in Acetonitrile(±)		1ml	

Pesticide and metabolites: Herbicides

Product code	Description			
Molinate-sulfoxide				
CAS 52236-29-0 DRE-C15281400	MW 203.3018 Molinate-sulfoxide	C ₉ H ₁₇ NO ₂ S	10mg	
Monalide				
CAS 7287-36-7 DRE-C15290000	MW 239.7411 Monalide	C ₁₃ H ₁₈ ClNO	100mg	
Monolinuron				
CAS 1746-81-2 DRE-C15310000 DRE-L15310000AL DRE-XA15310000AL DRE-LA09010255AL DRE-XA09010265HE	MW 214.6488 Monolinuron Monolinuron 10 µg/mL in Acetonitrile Monolinuron 100 µg/mL in Acetonitrile Monolinuron 10 µg/mL in Acetonitrile(‡) Monolinuron 100 µg/mL in Hexane(‡)	C ₉ H ₁₁ ClN ₂ O ₂	250mg 10ml 1ml 10ml 1ml	
Monolinuron D6 (methyl D3 methoxy D3)				
CAS n/a DRE-C15310100 DRE-XA15310100AC	MW 220.6858 Monolinuron D6 (methyl D3 methoxy D3) Monolinuron D6 (methyl D3 methoxy D3) 100 µg/mL in Acetone	C ₉ ² H ₆ ² H ₅ ClN ₂ O ₂	10mg 1.1ml	
Monuron				
CAS 150-68-5 DRE-C15320000 DRE-L15320000AL DRE-XA15320000AL DRE-XA09010213ME	MW 198.6494 Monuron Monuron 10 µg/mL in Acetonitrile Monuron 100 µg/mL in Acetonitrile Monuron 100 µg/mL in Methanol(‡)	C ₉ H ₁₁ ClN ₂ O	100mg 10ml 1ml 1ml	
Monuron D6 (dimethyl D6)				
CAS 217488-65-8 DRE-C15320100 DRE-XA15320100AC	MW 204.6864 Monuron D6 Monuron D6 100 µg/mL in Acetone	C ₉ ² H ₆ ² H ₅ ClN ₂ O	5mg 1.1ml	
Monuron Trichloroacetate (Monuron TCA)				
CAS 140-41-0 DRE-C15320200	MW 362.0366 Monuron TCA	C ₉ H ₁₁ ClN ₂ O · C ₂ HCl ₃ O ₂	100mg	
Naproanilide				
CAS 52570-16-8 DRE-LA15479300AC	MW 291.3438 Naproanilide 10 µg/mL in Acetone	C ₁₉ H ₁₇ NO ₂	1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Napropamide				
CAS 15299-99-7	MW 271.3541	C ₁₇ H ₂₁ NO ₂		
DRE-C15480000	Napropamide		100mg	
DRE-L15480000CY	Napropamide 10 µg/mL in Cyclohexane		10ml	
DRE-XA15480000CY	Napropamide 100 µg/mL in Cyclohexane		1ml	
Naptalam				
CAS 132-66-1	MW 291.3007	C ₁₈ H ₁₃ NO ₃		
DRE-C15490000	Naptalam		100mg	
Neburon				
CAS 555-37-3	MW 275.1742	C ₁₂ H ₁₆ Cl ₂ N ₂ O		
DRE-C15500000	Neburon		100mg	
DRE-L15500000AL	Neburon 10 µg/mL in Acetonitrile		10ml	
DRE-XA15500000AL	Neburon 100 µg/mL in Acetonitrile		1ml	
Nicosulfuron				
CAS 111991-09-4	MW 410.405	C ₁₅ H ₁₈ N ₆ O ₆ S		
DRE-CA15515000	Nicosulfuron(‡)		100mg	
Nitralin				
CAS 4726-14-1	MW 345.3715	C ₁₃ H ₁₈ N ₄ O ₆ S		
DRE-C15540000	Nitralin		100mg	
DRE-XA09010215AL	Nitralin 100 µg/mL in Acetonitrile(‡)		1ml	
Nitrofen				
CAS 1836-75-5	MW 284.0949	C ₁₂ H ₇ Cl ₂ NO ₃		
DRE-C15560000	Nitrofen(‡)		100mg	
DRE-L15560000AL	Nitrofen 10 µg/mL in Acetonitrile		10ml	
DRE-L15560000CY	Nitrofen 10 µg/mL in Cyclohexane		10ml	
Nitrofluorfen				
CAS 42874-01-1	MW 317.6478	C ₁₃ H ₇ ClF ₃ NO ₃		
DRE-LA15570000CY	Nitrofluorfen 10 µg/mL in Cyclohexane		1ml	
Nonanoic Acid (Pelargonic acid)				
CAS 112-05-0	MW 158.238	C ₉ H ₁₈ O ₂		
DRE-C15623100	Nonanoic acid (Pelargonic acid)(‡)		1ml	
Norflurazon				
CAS 27314-13-2	MW 303.6676	C ₁₂ H ₉ ClF ₃ N ₃ O		
DRE-C15650000	Norflurazon		100mg	
DRE-L15650000AL	Norflurazon 10 µg/mL in Acetonitrile		10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Norflurazon-desmethyl				
CAS 23576-24-1 DRE-LA15651000EA	MW 289.641 Norflurazon-desmethyl 10 µg/mL in Ethylacetate(±)	C ₁₁ H ₇ ClF ₃ N ₃ O	1ml	
Noruron				
CAS 18530-56-8 DRE-XA09010216AL	MW 222.3266 Noruron 100 µg/mL in Acetonitrile(±)	C ₁₃ H ₂₂ N ₂ O	1ml	
Oleic Acid Methyl Ester				
CAS 112-62-9 DRE-C15727060	MW 296.4879 Oleic acid methyl ester	C ₁₉ H ₃₆ O ₂	100mg	
Omethoate				
CAS 1113-02-6 DRE-C15730000 DRE-L15730000AL DRE-L15730000EA DRE-XA15730000AC	MW 213.1918 Omethoate(±) Omethoate 10 µg/mL in Acetonitrile Omethoate 10 µg/mL in Ethylacetate Omethoate 100 µg/mL in Acetone	C ₅ H ₁₂ NO ₄ PS	100mg 10ml 10ml 1ml	
Omethoate D6 (O-dimethyl D6)				
CAS 1219804-92-8 DRE-XA15730100AC	MW 219.2288 Omethoate D6 (dimethyl D6) 100 µg/mL in Acetone	C ₅ ² H ₆ H ₆ NO ₄ PS	1.1ml	
Orbencarb				
CAS 34622-58-7 DRE-CA15740000	MW 257.7795 Orbencarb	C ₁₂ H ₁₆ ClNOS	100mg	
Orthosulfamuron				
CAS 213464-77-8 DRE-C15748500	MW 424.4316 Orthosulfamuron(±)	C ₁₆ H ₂₀ NaO ₆ S	50mg	
Oryzalin				
CAS 19044-88-3 DRE-C15750000 DRE-L15750000AL	MW 346.3595 Oryzalin Oryzalin 10 µg/mL in Acetonitrile	C ₁₂ H ₁₈ N ₄ O ₆ S	100mg 10ml	
Oxadiargyl				
CAS 39807-15-3 DRE-C15758000 DRE-L15758000AL	MW 341.1893 Oxadiargyl Oxadiargyl 10 µg/mL in Acetonitrile	C ₁₅ H ₁₄ Cl ₂ N ₂ O ₃	100mg 10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Oxadiazon				
CAS 19666-30-9	MW 345.221	C ₁₅ H ₁₈ Cl ₂ N ₂ O ₃		
DRE-C15760000	Oxadiazon(‡)		100mg	
DRE-L15760000AL	Oxadiazon 10 µg/mL in Acetonitrile		10ml	
DRE-L15760000CY	Oxadiazon 10 µg/mL in Cyclohexane		10ml	
DRE-XA15760000CY	Oxadiazon 100 µg/mL in Cyclohexane		1ml	
Oxadiazon-hydroxy (5-tert-Butyl-3-(2,4-dichloro-5-hydroxyphenyl)-1,3,4-oxadiazol-2-one)				
CAS 39807-19-7	MW 303.1413	C ₁₂ H ₁₂ Cl ₂ N ₂ O ₃		
DRE-LA15760500TO	Oxadiazon-hydroxy 10 µg/mL in Toluene		1ml	
Oxasulfuron				
CAS 144651-06-9	MW 406.413	C ₁₇ H ₁₈ N ₄ O ₆ S		
DRE-C15781500	Oxasulfuron		100mg	
Oxaziclomefone				
CAS 153197-14-9	MW 376.2764	C ₂₀ H ₁₉ Cl ₂ NO ₂		
DRE-C15782000	Oxaziclomefone(‡)		10mg	
DRE-L15782000AL	Oxaziclomefone 10 µg/mL in Acetonitrile		10ml	
Oxyfluorfen				
CAS 42874-03-3	MW 361.7003	C ₁₅ H ₁₁ ClF ₃ NO ₄		
DRE-C15800000	Oxyfluorfen(‡)		100mg	
DRE-L15800000CY	Oxyfluorfen 10 µg/mL in Cyclohexane		10ml	
DRE-GA09010369ME	Oxyfluorfen 100 µg/mL in Methanol(‡)		1ml	
Paraquat Dichloride				
CAS 1910-42-5	MW 257.159	C ₁₂ H ₁₄ N ₂ ·2Cl		
DRE-C15870000	Paraquat dichloride		100mg	
DRE-L15870000WA	Paraquat dichloride 10 µg/mL in Water		10ml	
DRE-XA15870000WA	Paraquat dichloride 100 µg/mL in Water		1ml	
Paraquat Dichloride D8				
CAS 347841-45-6	MW 265.2083	C ₁₂ ² H ₆ H ₈ N ₂ ·2Cl		
DRE-CA15870100	Paraquat dichloride D8		50mg	
Paraquat dichloride D6 (dimethyl)				
CAS n/a	MW 263.1959	C ₁₂ ² H ₆ H ₈ N ₂ ·2Cl		
DRE-C15870050	Paraquat dichloride D6		50mg	
Paraquat diiodide D6 (dimethyl)				
CAS n/a	MW 446.0989	C ₁₂ ² H ₆ H ₈ N ₂ ·2I		
DRE-C15870200	Paraquat diiodide D6		50mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Pebulate				
CAS 1114-71-2	MW 203.3448	C ₁₀ H ₂₁ NOS		
DRE-C15904000	Pebulate(±)		100mg	
DRE-L15904000AL	Pebulate 10 µg/mL in Acetonitrile		10ml	
DRE-XA09010219AL	Pebulate 100 µg/mL in Acetonitrile(±)		1ml	
Pendimethalin				
CAS 40487-42-1	MW 281.3077	C ₁₃ H ₁₉ N ₃ O ₄		
DRE-C15930000	Pendimethalin(±)		100mg	
DRE-L15930000AL	Pendimethalin 10 µg/mL in Acetonitrile		10ml	
DRE-L15930000CY	Pendimethalin 10 µg/mL in Cyclohexane		10ml	
DRE-XA15930000CY	Pendimethalin 100 µg/mL in Cyclohexane		1ml	
Pendimethalin D5 (pent-3-yl (2,2,3,4,4)-D5)				
CAS 1219803-39-0	MW 286.3385	C ₁₃ ² H ₁₆ H ₁₄ N ₃ O ₄		
DRE-C15930100	Pendimethalin D5		10mg	
DRE-XA15930100AC	Pendimethalin D5 100 µg/mL in Acetone		1.1ml	
Penoxsulam				
CAS 219714-96-2	MW 483.3699	C ₁₆ H ₁₄ F ₅ N ₅ O ₅ S		
DRE-C15937000	Penoxsulam(±)		50mg	
DRE-XA15937000AL	Penoxsulam 100 µg/mL in Acetonitrile		1ml	
Pentanochlor				
CAS 2307-68-8	MW 239.7411	C ₁₃ H ₁₈ ClNO		
DRE-C15980000	Pentanochlor		25mg	
Pentoxazone				
CAS 110956-75-7	MW 353.7726	C ₁₇ H ₁₇ ClFNO ₄		
DRE-C15981770	Pentoxazone(±)		10mg	
Perfluidone				
CAS 37924-13-3	MW 379.3746	C ₁₄ H ₁₂ F ₃ NO ₄ S ₂		
DRE-C15985000	Perfluidone		25mg	
Pethoxamid				
CAS 106700-29-2	MW 295.8044	C ₁₆ H ₂₂ ClNO ₂		
DRE-C16000500	Pethoxamid		100mg	
Phenisopham				
CAS 57375-63-0	MW 342.389	C ₁₉ H ₂₂ N ₂ O ₄		
DRE-C16005000	Phenisopham		250mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Phenmedipham				
CAS 13684-63-4 DRE-C16020000	MW 300.3092 Phenmedipham	C ₁₆ H ₁₆ N ₂ O ₄	100mg	
Phenmedipham D3 (methoxy D3)				
CAS 1773497-41-8 DRE-C16020100	MW 303.3277 Phenmedipham D3	C ₁₆ ² H ₁₃ H ₁₃ N ₂ O ₄	5mg	
Phenoxyacetic Acid				
CAS 122-59-8 DRE-C16045000	MW 152.1473 Phenoxyacetic acid	C ₈ H ₈ O ₃	250mg	
Picloram				
CAS 1918-02-1 DRE-C16200000 DRE-YS09010035MB	MW 241.4592 Picloram Picloram 200 µg/mL in Methyl tert-Butyl Ether(±)	C ₈ H ₃ Cl ₃ N ₂ O ₂	250mg 5x1ml	
Picloram-2-ethylhexyl ester				
CAS 36374-99-9 DRE-C16200350	MW 353.6719 Picloram-2-ethylhexyl ester	C ₁₄ H ₁₉ Cl ₃ N ₂ O ₂	10mg	
Picolinafen				
CAS 137641-05-5 DRE-C16205000 DRE-L16205000AL DRE-L16205000CY	MW 376.3044 Picolinafen Picolinafen 10 µg/mL in Acetonitrile Picolinafen 10 µg/mL in Cyclohexane	C ₁₉ H ₁₂ F ₄ N ₂ O ₂	100mg 10ml 10ml	
Pinoxaden				
CAS 243973-20-8 DRE-C16215000	MW 400.5112 Pinoxaden(±)	C ₂₃ H ₃₂ N ₂ O ₄	25mg	
Piperophos				
CAS 24151-93-7 DRE-C16241000 DRE-L16241000IO	MW 353.4808 Piperophos Piperophos 10 µg/mL in Isooctane	C ₁₄ H ₂₈ NO ₃ PS ₂	250mg 10ml	
Pretilachlor				
CAS 51218-49-6 DRE-C16287000 DRE-L16287000CY DRE-XA09010225AL	MW 311.8468 Pretilachlor Pretilachlor 10 µg/mL in Cyclohexane Pretilachlor 100 µg/mL in Acetonitrile(±)	C ₁₇ H ₂₆ ClNO ₂	250mg 10ml 1ml	

Pesticide and metabolites: Herbicides

Product code	Description			
Primisulfuron-methyl				
CAS 86209-51-0 DRE-C16288000	MW 468.337 Primisulfuron-methyl	C ₁₅ H ₁₂ F ₄ N ₄ O ₇ S	100mg	
Procyazine				
CAS 32889-48-8 DRE-C16300000 DRE-XA09010226AL	MW 241.7205 Procyazine Procyazine 100 µg/mL in Acetonitrile(‡)	C ₁₀ H ₁₆ ClN ₅	10mg 1ml	
Prodiamine				
CAS 29091-21-2 DRE-C16320000	MW 350.2937 Prodiamine	C ₁₃ H ₁₇ F ₃ N ₄ O ₄	100mg	
Profluralin				
CAS 26399-36-0 DRE-C16340000	MW 347.2897 Profluralin	C ₁₄ H ₁₆ F ₃ N ₃ O ₄	250mg	
Profoxydim-lithium				
CAS 281664-76-4 DRE-C16341000	MW 471.9662 Profoxydim-lithium	C ₂₄ H ₃₁ ClNO ₄ S·Li	100mg	
Prometon				
CAS 1610-18-0 DRE-C16360000 DRE-L16360000AL DRE-L16360000CY	MW 225.2908 Prometon Prometon 10 µg/mL in Acetonitrile Prometon 10 µg/mL in Cyclohexane	C ₁₀ H ₁₉ N ₅ O	100mg 10ml 10ml	
Prometon D14 (di-iso-propyl D14)				
CAS 1219803-45-8 DRE-XA16360005AL	MW 239.377 Prometon D14 (di-iso-propyl D14) 100 µg/mL in Acetonitrile	C ₁₀ ² H ₁₄ H ₅ N ₅ O	1.1ml	
Prometryn				
CAS 7287-19-6 DRE-C16370000 DRE-L16370000AL DRE-L16370000CY DRE-XA16370000AL	MW 241.3564 Prometryn Prometryn 10 µg/mL in Acetonitrile Prometryn 10 µg/mL in Cyclohexane Prometryn 100 µg/mL in Acetonitrile	C ₁₀ H ₁₉ N ₅ S	250mg 10ml 10ml 1ml	
Prometryn D6 (isopropyl D6)				
CAS 1705649-52-0 DRE-XA16370100AC	MW 247.3933 Prometryn D6 100 µg/mL in Acetone	C ₁₀ ² HeH ₁₃ N ₅ S	1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Propachlor				
CAS 1918-16-7	MW 211.688	C ₁₁ H ₁₄ ClNO		
DRE-C16380000	Propachlor		250mg	
DRE-L16380000CY	Propachlor 10 µg/mL in Cyclohexane		10ml	
DRE-XA16380000CY	Propachlor 100 µg/mL in Cyclohexane		1ml	
Propachlor D5 (phenyl D5)				
CAS n/a	MW 216.7188	C ₁₁ ² H ₈ H ₉ ClNO		
DRE-XA16380010AC	Propachlor D5 (phenyl D5) 100 µg/mL in Acetone		1.1ml	
Propachlor ethane sulfonic acid (ESA) sodium salt				
CAS 947601-88-9	MW 279.2879	C ₁₁ H ₁₄ NO ₄ S·Na		
DRE-CA16380210	Propachlor ethane sulfonic acid (ESA) sodium salt		10mg	
Propachlor-2-hydroxy				
CAS 42404-06-8	MW 193.2423	C ₁₁ H ₁₃ NO ₂		
DRE-C16380180	Propachlor-2-hydroxy		100mg	
Propachlor Oxalamic Acid				
CAS 70628-36-3	MW 207.2258	C ₁₁ H ₁₃ NO ₃		
DRE-CA16380400	Propachlor oxalamic acid (OA)		10mg	
Propanil				
CAS 709-98-8	MW 218.0799	C ₉ H ₉ Cl ₂ NO		
DRE-C16410000	Propanil		250mg	
DRE-L16410000AL	Propanil 10 µg/mL in Acetonitrile		10ml	
DRE-L16410000CY	Propanil 10 µg/mL in Cyclohexane		10ml	
Propaquizafop				
CAS 111479-05-1	MW 443.8802	C ₂₂ H ₂₂ ClN ₃ O ₅		
DRE-C16425000	Propaquizafop		100mg	
DRE-L16425000CY	Propaquizafop 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010228AL	Propaquizafop 100 µg/mL in Acetonitrile(‡)		1ml	
Propazine				
CAS 139-40-2	MW 229.7098	C ₉ H ₁₆ ClN ₅		
DRE-C16440000	Propazine		250mg	
DRE-L16440000AL	Propazine 10 µg/mL in Acetonitrile		10ml	
DRE-L16440000CY	Propazine 10 µg/mL in Cyclohexane		10ml	
DRE-XA16440000AL	Propazine 100 µg/mL in Acetonitrile		1ml	
Propazine 13C3 (ring 13C3)				
CAS 446276-68-2	MW 232.6878	¹³ C ₃ C ₆ H ₁₆ ClN ₅		
DRE-XA16440200AC	Propazine 13C3 100 µg/mL in Acetone		1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Propazine D6 (isopropyl D6)				
CAS 1655498-05-7 DRE-XA16440100AC	MW 235.7468 Propazine D6 100 µg/mL in Acetone	C ₉ H ₆ H ₁₀ ClN ₃	1.1ml	
Propazine-2-hydroxy				
CAS 7374-53-0 DRE-C16445000	MW 211.2642 Propazine-2-hydroxy	C ₉ H ₁₇ N ₃ O	100mg	
Propham				
CAS 122-42-9 DRE-C16470000 DRE-L16470000AL DRE-XA16470000AL	MW 179.2157 Propham Propham 10 µg/mL in Acetonitrile Propham 100 µg/mL in Acetonitrile	C ₁₀ H ₁₃ NO ₂	250mg 10ml 1ml	
Propisochlor				
CAS 86763-47-5 DRE-C16495000 DRE-L16495000CY	MW 283.7937 Propisochlor Propisochlor 10 µg/mL in Cyclohexane	C ₁₅ H ₂₂ ClNO ₂	50mg 10ml	
Propoxycarbazone-sodium				
CAS 181274-15-7 DRE-C16500500 DRE-L16500500AL DRE-XA16500500AL	MW 420.3728 Propoxycarbazone-sodium(±) Propoxycarbazone-sodium 10 µg/mL in Acetonitrile Propoxycarbazone sodium 100 µg/mL in Acetonitrile	C ₁₅ H ₁₇ Na ₄ O ₇ S·Na	100mg 10ml 1ml	
Propyzamide				
CAS 23950-58-5 DRE-C16540000 DRE-L16540000CY DRE-XA16540000AL	MW 256.1278 Propyzamide Propyzamide 10 µg/mL in Cyclohexane Propyzamide 100 µg/mL in Acetonitrile	C ₁₂ H ₁₁ Cl ₂ NO	250mg 10ml 1ml	
Propyzamide D3 (phenyl-2,4,6 D3)				
CAS 1219805-79-4 DRE-XA16540010AL	MW 259.1463 Propyzamide D3 (phenyl-2,4,6 D3) 100 µg/mL in Acetonitrile	C ₁₂ ² H ₃ H ₈ Cl ₂ NO	1.1ml	
Prosulfocarb				
CAS 52888-80-9 DRE-C16545000 DRE-L16545000AL	MW 251.3876 Prosulfocarb Prosulfocarb 10 µg/mL in Acetonitrile	C ₁₄ H ₂₁ NOS	250mg 10ml	
Prosulfuron				
CAS 94125-34-5 DRE-C16546000	MW 419.3788 Prosulfuron(±)	C ₁₅ H ₁₆ F ₃ N ₅ O ₄ S	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Proximpham				
CAS 2828-42-4 DRE-C16575000	MW 192.2145 Proximpham	C ₁₀ H ₁₂ N ₂ O ₂	100mg	
Pyraclonil				
CAS 158353-15-2 DRE-C16593000	MW 314.7728 Pyraclonil	C ₁₅ H ₁₃ ClN ₆	10mg	
Pyraflufen (free acid)				
CAS 129630-17-7 DRE-C16597100	MW 385.1228 Pyraflufen (free acid)	C ₁₃ H ₉ Cl ₂ F ₃ N ₂ O ₄	10mg	
Pyraflufen-ethyl				
CAS 129630-19-9 DRE-C16597000 DRE-L16597000CY	MW 413.1759 Pyraflufen-ethyl(±) Pyraflufen-ethyl 10 µg/mL in Cyclohexane	C ₁₅ H ₁₃ Cl ₂ F ₃ N ₂ O ₄	100mg 10ml	
Pyrasulfotole				
CAS 365400-11-9 DRE-C16604000	MW 362.3242 Pyrasulfotole	C ₁₄ H ₁₃ F ₃ N ₂ O ₄ S	100mg	
Pyrasulfotole-desmethyl				
CAS 936621-05-5 DRE-C16604500	MW 348.2976 Pyrasulfotole-desmethyl	C ₁₃ H ₁₁ F ₃ N ₂ O ₄ S	10mg	
Pyrazolynate				
CAS 58011-68-0 DRE-C16608500	MW 439.3123 Pyrazolynate	C ₁₉ H ₁₆ Cl ₂ N ₂ O ₄ S	25mg	
Pyrazosulfuron-ethyl				
CAS 93697-74-6 DRE-C16611500	MW 414.3937 Pyrazosulfuron-ethyl	C ₁₄ H ₁₆ NaO ₇ S	100mg	
Pyrazoxyfen				
CAS 71561-11-0 DRE-C16612500	MW 403.2586 Pyrazoxyfen	C ₂₀ H ₁₆ Cl ₂ N ₂ O ₃	100mg	










Pesticide and metabolites: Herbicides

Product code	Description			
Pyribenzoxim				
CAS 168088-61-7 DRE-C16624000	MW 609.5855 Pyribenzoxim	C ₃₂ H ₂₇ N ₅ O ₈	100mg	
Pyributicarb				
CAS 88678-67-5 DRE-C16626000 DRE-L16626000CY DRE-LA16626000CY	MW 330.4445 Pyributicarb Pyributicarb 10 µg/mL in Cyclohexane Pyributicarb 10 µg/mL in Cyclohexane	C ₁₈ H ₂₂ N ₂ O ₂ S	10mg 10ml 1ml	
Pyridate				
CAS 55512-33-9 DRE-C16640000	MW 378.9161 Pyridate(±)	C ₁₉ H ₂₃ ClN ₂ O ₂ S	250mg	
Pyriftalid				
CAS 135186-78-6 DRE-C16656000	MW 318.3477 Pyriftalid	C ₁₅ H ₁₄ N ₂ O ₄ S	100mg	
(E)-Pyriminobac-methyl				
CAS 147411-69-6 DRE-C16659510 DRE-L16659510CY	MW 361.3493 (E)-Pyriminobac-methyl Pyriminobac-methyl (E) 10 µg/mL in Cyclohexane	C ₁₇ H ₁₈ N ₄ O ₆	10mg 10ml	
Pyriminobac-methyl (Z)				
CAS 147411-70-9 DRE-C16659520 DRE-LA16659520CY	MW 361.3493 Pyriminobac-methyl (Z) Pyriminobac-methyl (Z) 10 µg/mL in Cyclohexane	C ₁₇ H ₁₈ N ₄ O ₆	10mg 1ml	
Pyrithiobac sodium				
CAS 123343-16-8 DRE-C16663000	MW 348.7373 Pyrithiobac sodium	C ₁₃ H ₁₀ ClN ₂ O ₄ S·Na	100mg	
Pyroxsulam				
CAS 422556-08-9 DRE-C16667000 DRE-XA16667000AL	MW 434.3504 Pyroxsulam Pyroxsulam 100 µg/mL in Acetonitrile	C ₁₄ H ₁₃ F ₃ N ₅ O ₅ S	100mg 1ml	
Quinchlorac				
CAS 84087-01-4 DRE-C16705000 DRE-L16705000AL	MW 242.0582 Quinchlorac Quinchlorac 10 µg/mL in Acetonitrile	C ₁₀ H ₅ Cl ₂ NO ₂	100mg 10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Quinclorac-methyl ester				
CAS 84087-33-2 DRE-C16706000	MW 256.0848 Quinclorac-methyl ester	C ₁₁ H ₇ Cl ₂ NO ₂	100mg	
Quinmerac				
CAS 90717-03-6 DRE-C16708000 DRE-L16708000AL DRE-XA16708000AL	MW 221.6397 Quinmerac Quinmerac 10 µg/mL in Acetonitrile Quinmerac 100 µg/mL in Acetonitrile	C ₁₁ H ₈ ClNO ₂	100mg 10ml 1ml	
Quinmerac D4 (methyl(D3)-quinoline-2-D)				
CAS n/a DRE-XA16708100AL	MW 225.6644 Quinmerac D4 (methyl(D3)-quinoline-2-D) 100 µg/mL in Acetonitrile	C ₁₁ ² H ₈ H ₄ ClNO ₂	1.1ml	
Quinoclamine				
CAS 2797-51-5 DRE-C16709500 DRE-XA16709500AC	MW 207.6131 Quinoclamine Quinoclamine 100 µg/mL in Acetone	C ₁₀ H ₆ ClNO ₂	50mg 1ml	
Quinonamid				
CAS 27541-88-4 DRE-C16710000	MW 318.5399 Quinonamid	C ₁₂ H ₆ Cl ₃ NO ₃	100mg	
Quizalofop-P-ethyl				
CAS 100646-51-3 DRE-C16740500	MW 372.8023 Quizalofop-P-ethyl	C ₁₉ H ₁₇ ClN ₂ O ₄	100mg	
Quizalofop-P-tefuryl				
CAS 200509-41-7 DRE-C16741000 DRE-L16741000AL	MW 428.8655 Quizalofop-P-tefuryl Quizalofop-P-tefuryl 10 µg/mL in Acetonitrile	C ₂₂ H ₂₁ ClN ₂ O ₅	25mg 10ml	
Quizalofop-ethyl				
CAS 76578-14-8 DRE-C16740000 DRE-L16740000IO DRE-XA09010232AL	MW 372.8023 Quizalofop-ethyl Quizalofop-ethyl 10 µg/mL in Isooctane Quizalofop Ethyl 100 µg/mL in Acetonitrile(±)	C ₁₉ H ₁₇ ClN ₂ O ₄	50mg 10ml 1ml	
Quizalofop-ethyl D3 (3,3,3-D3)				
CAS 1398065-84-3 DRE-XA16740100AC	MW 375.8208 Quizalofop-ethyl D3 (3,3,3-D3) 100 µg/mL in Acetone	C ₁₉ ² H ₃ H ₁₄ ClN ₂ O ₄	1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Quizalofop free acid				
CAS 76578-12-6	MW 344.7491	C ₁₇ H ₁₃ ClN ₂ O ₄		
DRE-C16739990	Quizalofop free acid		50mg	
DRE-L16739990AL	Quizalofop free acid 10 µg/mL in Acetonitrile		10ml	
DRE-GS09010057DI	Quizalofop 1000 µg/mL in Dichloromethane(‡)		4x1ml	
Quizalofop free acid D3 (methyl D3)				
CAS n/a	MW 347.7676	C ₁₇ ² H ₁₃ H ₁₀ ClN ₂ O ₄		
DRE-XA16739991AC	Quizalofop free acid D3 100 µg/mL in Acetone		1.1ml	
Quizalofop-methyl				
CAS 76578-13-7	MW 358.7757	C ₁₈ H ₁₅ ClN ₂ O ₄		
DRE-C16740800	Quizalofop-methyl		20mg	
Rimsulfuron				
CAS 122931-48-0	MW 431.4441	C ₁₄ H ₁₇ N ₅ O ₇ S ₂		
DRE-C16815000	Rimsulfuron(‡)		100mg	
Rimsulfuron-desulfon				
CAS 138724-53-5	MW 367.3803	C ₁₄ H ₁₇ N ₅ O ₅ S		
DRE-C16815060	Rimsulfuron-desulfon		10mg	
Rimsulfuron metabolite 2				
CAS 151331-80-5	MW 324.3555	C ₁₃ H ₁₆ N ₄ O ₄ S		
DRE-XA16815090AL	Rimsulfuron metabolite 2 (100 µg/mL) in Acetonitrile		1ml	
Saflufenacil				
CAS 372137-35-4	MW 500.8523	C ₁₇ H ₁₇ ClF ₄ N ₄ O ₅ S		
DRE-C16901600	Saflufenacil(‡)		100mg	
Sebuthylazine				
CAS 7286-69-3	MW 229.7098	C ₉ H ₁₆ ClN ₅		
DRE-C16920000	Sebuthylazine		250mg	
DRE-L16920000AL	Sebuthylazine 10 µg/mL in Acetonitrile		10ml	
DRE-XA16920000AL	Sebuthylazine 100 µg/mL in Acetonitrile		1ml	
Sebuthylazine D5 (N-ethyl D5)				
CAS 1219805-56-7	MW 234.7406	C ₉ ² H ₅ H ₁₁ ClN ₅		
DRE-XA16920100AC	Sebuthylazine D5 100 µg/mL in Acetone		1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Sebuthylazine-desethyl (6-Chloro-N-sec-butyl-1,3,5-triazine-2,4-diamine)				
CAS 37019-18-4	MW 201.6567	C ₇ H ₁₂ ClN ₅		
DRE-C16920500	Sebuthylazine-desethyl		250mg	
DRE-XA09010235AL	Sebuthylazin-desethyl 100 µg/mL in Acetonitrile(±)		1ml	
Sebuthylazine-2-hydroxy				
CAS 33124-61-7	MW 211.2642	C ₉ H ₁₇ N ₅ O		
DRE-L16920900AL	Sebuthylazine-2-hydroxy 10 µg/mL in Acetonitrile		10ml	
Secbumenton				
CAS 26259-45-0	MW 225.2908	C ₁₀ H ₁₉ N ₅ O		
DRE-C16930000	Secbumenton		250mg	
DRE-L16930000AL	Secbumenton 10 µg/mL in Acetonitrile		10ml	
DRE-XA09010236AL	Secbumenton 100 µg/mL in Acetonitrile(±)		1ml	
Secbumenton D5 (N-ethyl D5)				
CAS 1705649-53-1	MW 230.3216	C ₁₀ ² H ₁₉ H ₁₄ N ₅ O		
DRE-C16930100	Secbumenton D5 (ethyl D5)		10mg	
DRE-XA16930100AC	Secbumenton D5 100 µg/mL in Acetone		1.1ml	
Sethoxydim				
CAS 74051-80-2	MW 327.4821	C ₁₇ H ₂₉ NO ₃ S		
DRE-C16940000	Sethoxydim		10mg	
DRE-L16940000CY	Sethoxydim 10 µg/mL in Cyclohexane (*)		10ml	
Siduron				
CAS 1982-49-6	MW 232.3214	C ₁₄ H ₂₀ N ₂ O		
DRE-C16945000	Siduron		100mg	
trans-Siduron				
CAS 19123-21-8	MW 232.3214	C ₁₄ H ₂₀ N ₂ O		
DRE-L16945500AL	trans-Siduron 10 µg/mL in Acetonitrile		10ml	
Simazine				
CAS 122-34-9	MW 201.6567	C ₇ H ₁₂ ClN ₅		
DRE-C16950000	Simazine(±)		250mg	
DRE-L16950000AL	Simazine 10 µg/mL in Acetonitrile		10ml	
DRE-L16950000EA	Simazine 10 µg/mL in Ethylacetate		10ml	
DRE-XA16950000AL	Simazine 100 µg/mL in Acetonitrile		1ml	
Simazine D10 (diethyl D5)				
CAS 220621-39-6	MW 211.7183	C ₇ ² H ₁₀ H ₂ ClN ₅		
DRE-C16950100	Simazine D10		10mg	
DRE-XA16950100AC	Simazine D10 100 µg/mL in Acetone		1.1ml	

Pesticide and metabolites: Herbicides

Product code	Description			
Simazine D5 (ethyl D5)				
CAS 220621-41-0	MW 206.6875	C ₇ H ₅ H ₇ ClN ₅		
DRE-C16950200	Simazine D5		10mg	
DRE-XA16950200AL	Simazine D5 100 µg/mL in Acetonitrile(‡)		1ml	
Simazine-2-hydroxy				
CAS 2599-11-3	MW 183.211	C ₇ H ₁₃ N ₅ O		
DRE-C16955000	Simazine-2-hydroxy(‡)		25mg	
DRE-L16955000ME	Simazine-2-hydroxy 10 µg/mL in Methanol		10ml	
Simazine mercapturate				
CAS 205764-69-8	MW 328.3906	C ₁₂ H ₂₀ N ₆ O ₃ S		
DRE-C16955500	Simazine mercapturate		100mg	
Simeton				
CAS 673-04-1	MW 197.2376	C ₈ H ₁₅ N ₅ O		
DRE-C16960000	Simeton		50mg	
Simetryn				
CAS 1014-70-6	MW 213.3032	C ₈ H ₁₅ N ₅ S		
DRE-C16970000	Simetryn		100mg	
DRE-XA16970000CY	Simetryn 100 µg/mL in Cyclohexane		1ml	
Sulcotrione				
CAS 99105-77-8	MW 328.768	C ₁₄ H ₁₃ ClO ₅ S		
DRE-C16988000	Sulcotrione		100mg	
Sulfallate				
CAS 95-06-7	MW 223.7865	C ₈ H ₁₄ ClN ₂ S ₂		
DRE-C16995000	Sulfallate		100mg	
Sulfentrazone				
CAS 122836-35-5	MW 387.1899	C ₁₁ H ₁₀ Cl ₂ F ₂ N ₄ O ₅ S		
DRE-C17000300	Sulfentrazone(‡)		100mg	
DRE-L17000300AL	Sulfentrazone 10 µg/mL in Acetonitrile		10ml	
Sulfometuron-methyl				
CAS 74222-97-2	MW 364.3763	C ₁₅ H ₁₆ N ₄ O ₅ S		
DRE-C17009000	Sulfometuron-methyl		100mg	

Pesticide and metabolites: Herbicides

Product code	Description			
Sulfosulfuron				
CAS 141776-32-1 DRE-C17009500	MW 470.4801 Sulfosulfuron(±)	C ₁₆ H ₁₈ N ₆ O ₇ S ₂	100mg	
2,4,5-T 2-Butoxyethyl Ester				
CAS 2545-59-7 DRE-C17102000	MW 355.6414 2,4,5-T-butoxyethyl ester	C ₁₄ H ₁₇ Cl ₃ O ₄	100mg	
2,4,5-T 1-Butoxy-1-methylethyl Ester				
CAS n/a DRE-C17103000	MW 369.668 2,4,5-T-2-butoxy isopropyl ester	C ₁₅ H ₁₉ Cl ₃ O ₄	100mg	
2,4,5-T Butyl Ester				
CAS 93-79-8 DRE-C17105500	MW 311.5888 2,4,5-T-1-butyl ester	C ₁₂ H ₁₃ Cl ₃ O ₃	100mg	
2,4,5-T-ester (technical)				
CAS n/a DRE-C17105000	MW 283.5357 2,4,5-T-ester (technical)	C ₉ H ₇ Cl ₃ O ₃ (CH ₂) _n	100mg	
2,4,5-T 2-Ethylhexyl Ester				
CAS 1928-47-8 DRE-C17110000	MW 367.6951 2,4,5-T-2-ethylhexyl ester	C ₁₆ H ₂₁ Cl ₃ O ₃	250mg	
2,4,5-T Isooctyl Ester				
CAS 25168-15-4 DRE-C17120000	MW 367.6951 2,4,5-T-iso-octyl ester (technical)	C ₁₆ H ₂₁ Cl ₃ O ₃	100mg	
2,4,5-T Isopropyl Ester				
CAS 93-78-7 DRE-C17121000	MW 297.5622 2,4,5-T-iso-propyl ester	C ₁₁ H ₁₁ Cl ₃ O ₃	100mg	
2,4,5-T Methyl Ester ((2,4,5-Trichlorophenoxy)acetic Acid Methyl Ester)				
CAS 1928-37-6 DRE-C17125000 DRE-L17125000CY	MW 269.5091 2,4,5-T-methyl ester 2,4,5-T-methyl ester 10 µg/mL in Cyclohexane	C ₉ H ₇ Cl ₃ O ₃	250mg 10ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Terbacil				
CAS 5902-51-2	MW 216.6647	C ₉ H ₁₃ ClN ₂ O ₂		
DRE-C17250000	Terbacil		250mg	
DRE-L17250000IO	Terbacil 10 µg/mL in Isooctane		10ml	
Terbucarb				
CAS 1918-11-2	MW 277.4018	C ₁₇ H ₂₇ NO ₂		
DRE-C17260000	Terbucarb		100mg	
Terbumeton				
CAS 33693-04-8	MW 225.2908	C ₁₀ H ₁₉ N ₅ O		
DRE-C17290000	Terbumeton		250mg	
DRE-L17290000CY	Terbumeton 10 µg/mL in Cyclohexane		10ml	
DRE-L17290000ME	Terbumeton 10 µg/mL in Methanol		10ml	
Terbumeton-desethyl (N-tert-Butyl-6-methoxy-1,3,5-triazine-2,4-diamine)				
CAS 30125-64-5	MW 197.2376	C ₈ H ₁₅ N ₅ O		
DRE-C17290500	Terbumeton-desethyl		250mg	
DRE-L17290500ME	Terbumeton-desethyl 10 µg/mL in Methanol		10ml	
Terbuthylazine				
CAS 5915-41-3	MW 229.7098	C ₈ H ₁₆ ClN ₅		
DRE-C17300000	Terbuthylazine		250mg	
DRE-L17300000AL	Terbuthylazine 10 µg/mL in Acetonitrile		10ml	
DRE-L17300000CY	Terbuthylazine 10 µg/mL in Cyclohexane		10ml	
DRE-XA17300000AC	Terbuthylazine 100 µg/mL in Acetone		1ml	
DRE-XA17300000AL	Terbuthylazine 100 µg/mL in Acetonitrile		1ml	
Terbuthylazine D5 (ethyl D5)				
CAS 222986-60-9	MW 234.7406	C ₈ ² H ₅ H ₁₁ ClN ₅		
DRE-C17300100	Terbuthylazine D5		5mg	
DRE-XA17300100AC	Terbuthylazine D5 100 µg/mL in Acetone		1.1ml	
Terbuthylazine-desethyl (N-tert-Butyl-6-chloro-1,3,5-triazine-2,4-diamine)				
CAS 30125-63-4	MW 201.6567	C ₇ H ₁₂ ClN ₅		
DRE-C17303000	Terbuthylazine-desethyl		100mg	
DRE-L17303000AL	Terbuthylazine-desethyl 10 µg/mL in Acetonitrile		10ml	
DRE-XA17303000AL	Terbuthylazine-desethyl 100 µg/mL in Acetonitrile		1ml	
DRE-XA09010243AL	Terbuthylazine-desethyl 100 µg/mL in Acetonitrile(‡)		1ml	
Terbuthylazine-desethyl D9 (tert-butyl D9)				
CAS 1219798-52-3	MW 210.7121	C ₇ ² H ₉ H ₃ ClN ₅		
DRE-C17303100	Terbuthylazine-desethyl D9 (tert-butyl D9)		10mg	
DRE-XA17303100AC	Terbuthylazine-desethyl D9 (tert-butyl D9) 100 µg/mL in Acetone		1.1ml	










Pesticide and metabolites: Herbicides

Product code	Description			
Terbutylazine-desethyl-2-hydroxy (4-Amino-6-(tert-butylamino)-1,3,5-triazin-2-ol)				
CAS 66753-06-8	MW 183.211	C ₇ H ₁₃ N ₅ O		
DRE-C17303300	Terbutylazine-desethyl-2-hydroxy(‡)		10mg	
DRE-L17303300AL	Terbutylazine-desethyl-2-hydroxy 10 µg/mL in Acetonitrile		10ml	
Terbutylazine-2-hydroxy (4-(tert-Butylamino)-6-(ethylamino)-1,3,5-triazin-2-ol)				
CAS 66753-07-9	MW 211.2642	C ₉ H ₁₇ N ₅ O		
DRE-C17305000	Terbutylazine-2-hydroxy		100mg	
DRE-L17305000AL	Terbutylazine-2-hydroxy 10 µg/mL in Acetonitrile		10ml	
Terbutryn				
CAS 886-50-0	MW 241.3564	C ₁₀ H ₁₉ N ₅ S		
DRE-C17320000	Terbutryn		250mg	
DRE-L17320000AL	Terbutryn 10 µg/mL in Acetonitrile		10ml	
DRE-L17320000CY	Terbutryn 10 µg/mL in Cyclohexane		10ml	
DRE-XA17320000AL	Terbutryn 100 µg/mL in Acetonitrile		1ml	
Terbutryn D5 (ethyl D5)				
CAS 1219804-47-3	MW 246.3872	C ₁₀ ² H ₁₅ H ₁₄ N ₅ S		
DRE-C17320100	Terbutryn D5		10mg	
DRE-XA17320100AC	Terbutryn D5 100 µg/mL in Acetone		1.1ml	
Thenylchlor				
CAS 96491-05-3	MW 323.8376	C ₁₆ H ₁₉ ClNO ₂ S		
DRE-C17445500	Thenylchlor(‡)		10mg	
DRE-LA17445500CY	Thenylchlor 10 µg/mL in Cyclohexane		1ml	
Thiazafluron				
CAS 25366-23-8	MW 240.2062	C ₆ H ₇ F ₃ N ₄ OS		
DRE-C17460000	Thiazafluron		250mg	
Thiazopyr				
CAS 117718-60-2	MW 396.3754	C ₁₆ H ₁₇ F ₅ N ₂ O ₂ S		
DRE-C17462000	Thiazopyr		100mg	
DRE-L17462000CY	Thiazopyr 10 µg/mL in Cyclohexane		10ml	
Thien carbazone-methyl				
CAS 317815-83-1	MW 390.3922	C ₁₂ H ₁₄ N ₄ O ₇ S ₂		
DRE-C17465500	Thien carbazone-methyl		100mg	
Thifensulfuron-methyl				
CAS 79277-27-3	MW 387.3915	C ₁₂ H ₁₃ N ₅ O ₆ S ₂		
DRE-C17466000	Thifensulfuron-methyl		100mg	

Pesticide and metabolites: Herbicides

Product code	Description			
Thifensulfuron-methyl D3				
CAS n/a DRE-C17466100	MW 390.41 Thifensulfuron-methyl D3	C ₁₂ H ₁₃ H ₁₀ N ₅ O ₆ S ₂	5mg	
Thiobencarb				
CAS 28249-77-6 DRE-C17470000 DRE-L17470000CY	Thiobencarb(‡) Thiobencarb 10 µg/mL in Cyclohexane	C ₁₂ H ₁₆ ClNOS	250mg 10ml	
Tiocarbazil				
CAS 36756-79-3 DRE-C17585000	Tiocarbazil	C ₁₆ H ₂₅ NOS	100mg	
Topramezone				
CAS 210631-68-8 DRE-C17602500 DRE-L17602500AL	Topramezone Topramezone 10 µg/mL in Acetonitrile	C ₁₆ H ₁₇ N ₃ O ₅ S	50mg 10ml	
Tralkoxydim				
CAS 87820-88-0 DRE-C17605000	Tralkoxydim	C ₂₀ H ₂₇ NO ₃	250mg	
Triallate				
CAS 2303-17-5 DRE-C17630000 DRE-L17630000CY DRE-XA17630000AL	Triallate Triallate 10 µg/mL in Cyclohexane Triallate 100 µg/mL in Acetonitrile	C ₁₀ H ₁₆ Cl ₃ NOS	250mg 10ml 1ml	
Triasulfuron				
CAS 82097-50-5 DRE-C17648000	Triasulfuron	C ₁₄ H ₁₆ ClN ₅ O ₅ S	100mg	
Triaziflam				
CAS 131475-57-5 DRE-C17648900	Triaziflam	C ₁₇ H ₂₄ FN ₅ O	25mg	
Tribenuron-methyl (technical)				
CAS 101200-48-0 DRE-C17662000	Tribenuron-methyl (technical)	C ₁₅ H ₁₇ N ₅ O ₆ S	100mg	

Pesticide and metabolites: Herbicides

Product code	Description			
Trichloroacetic Acid (TCA)				
CAS 76-03-9 DRE-C17683500	MW 163.3871 Trichloroacetic acid (TCA)	$C_2HCl_3O_2$	250mg	
Trichloroacetic Acid Methyl Ester (TCA-methyl ester)				
CAS 598-99-2 DRE-C17684000 DRE-YA17684000MB	MW 177.4137 Trichloroacetic acid methyl ester (TCA-methyl ester) Trichloroacetic acid methyl ester 1000 µg/mL in Methyl-tert. butyl ether	$C_3H_3Cl_3O_2$	250mg 1ml	
Trichloroacetic Acid Sodium Salt				
CAS 650-51-1 DRE-C17684500	MW 185.369 Trichloroacetic acid sodium salt	$C_2Cl_3O_2\cdot Na$	250mg	
Triclopyr				
CAS 55335-06-3 DRE-C17800000 DRE-L17800000AC DRE-XA17800000AC DRE-XA17800000AL	MW 256.4706 Triclopyr Triclopyr 10 µg/mL in Acetone Triclopyr 100 µg/mL in Acetone Triclopyr 100 µg/mL in Acetonitrile	$C_7H_4Cl_3NO_3$	250mg 10ml 1ml 1ml	
Triclopyr-butotyl (Triclopyr-2-butoxyethyl ester)				
CAS 64700-56-7 DRE-C17800500	MW 356.6294 Triclopyr-2-butoxyethyl ester	$C_{13}H_{16}Cl_3NO_4$	100mg	
Triclopyr-methyl				
CAS 60825-26-5 DRE-C17800600	MW 270.4971 Triclopyr-methyl	$C_8H_6Cl_3NO_3$	100mg	
Tridiphane (2-(3,5-Dichlorophenyl)-2-(2,2,2-trichloroethyl)oxirane)				
CAS 58138-08-2 DRE-C17823000 DRE-L17823000IO	MW 320.427 Tridiphane Tridiphane 10 µg/mL in Isooctane	$C_{10}H_7Cl_5O$	100mg 10ml	
Trietazine				
CAS 1912-26-1 DRE-C17830000 DRE-L17830000AL	MW 229.7098 Trietazine Trietazine 10 µg/mL in Acetonitrile	$C_9H_{16}ClN_5$	100mg 10ml	
Trietazine D5 (ethyl D5)				
CAS 1397243-73-0 DRE-XA17830100AC	MW 234.7406 Trietazine D5 100 µg/mL in Acetone	$C_9^2H_{16}ClN_5$	1.1ml	

Pesticide and metabolites: Herbicides

Product code	Description			
Trietazine-desethyl				
CAS 38902-68-0 DRE-C17830500	MW 201.6567 Trietazine-desethyl	C ₇ H ₁₂ ClN ₅	100mg	
Trietazine-2-hydroxy				
CAS 13532-25-7 DRE-L17830900AL	MW 211.2642 Trietazine-2-hydroxy 10 µg/mL in Acetonitrile	C ₉ H ₁₇ N ₅ O	10ml	
Trifloxysulfuron sodium				
CAS 199119-58-9 DRE-C17843000	MW 459.3329 Trifloxysulfuron sodium(±)	C ₁₄ H ₁₃ F ₃ N ₅ O ₆ S·Na	100mg	
Trifluralin				
CAS 1582-09-8 DRE-C17850000 DRE-L17850000CY	MW 335.279 Trifluralin(±) Trifluralin 10 µg/mL in Cyclohexane (*)	C ₁₃ H ₁₆ F ₃ N ₃ O ₄	250mg 10ml	
Trifluralin D14 (di-n-propyl D14)				
CAS 347841-79-6 DRE-C17850100 DRE-XA17850100AC	MW 349.3653 Trifluralin D14 (di-n-propyl D14) Trifluralin D14 (di-n-propyl D14) 100 µg/mL in Acetone (*)	C ₁₃ ² H ₁₄ H ₂ F ₃ N ₃ O ₄	10mg 1.1ml	
Triflusulfuron-methyl				
CAS 126535-15-7 DRE-C17851500	MW 492.4296 Triflusulfuron-methyl	C ₁₇ H ₁₉ F ₃ N ₆ O ₆ S	100mg	
Tritosulfuron				
CAS 142469-14-5 DRE-C17894700	MW 445.2971 Tritosulfuron	C ₁₃ H ₉ F ₆ N ₅ O ₄ S	100mg	
Uniconazole				
CAS 83657-22-1 DRE-L17897000CY	MW 291.7759 Uniconazole 10 µg/mL in Cyclohexane	C ₁₅ H ₁₆ ClN ₃ O	10ml	
Uniconazole-P				
CAS 83657-17-4 DRE-L17897100CY	MW 291.7759 Uniconazole-P 10 µg/mL in Cyclohexane	C ₁₅ H ₁₆ ClN ₃ O	10ml	

Pesticide and metabolites: Herbicides

Product code	Description

Vernolate

CAS 1929-77-7


MW 203.3448

C₁₀H₂₁NOS

[DRE-C17910000](#)

Vernolate










250mg












PESTICIDES AND METABOLITES: INSECTICIDES












Pesticide and metabolites: Insecticides

Product code	Description			
Acephate				
CAS 30560-19-1	MW 183.1659	C ₄ H ₁₀ NO ₃ PS		
DRE-C10010000	Acephate(‡)		250mg	
DRE-L10010000AL	Acephate 10 µg/mL in Acetonitrile		10ml	
DRE-XA10010000AL	Acephate 100 µg/mL in Acetonitrile		1ml	
Acequinocyl				
CAS 57960-19-7	MW 384.5085	C ₂₄ H ₃₂ O ₄		
DRE-C10010500	Acequinocyl		50mg	
Acequinocyl-hydroxy (2-Dodecyl-3-hydroxy-1,4-naphthoquinone)				
CAS 57960-31-3	MW 342.4718	C ₂₂ H ₃₀ O ₃		
DRE-C10010520	Acequinocyl-hydroxy		10mg	
Acetamiprid				
CAS 135410-20-7	MW 222.6741	C ₁₀ H ₁₁ ClN ₄		
DRE-C10013000	Acetamiprid(‡)		100mg	
DRE-L10013000EA	Acetamiprid 10 µg/mL in Ethylacetate		10ml	
Acetamiprid D3 (N-methyl D3)				
CAS 1353869-35-8	MW 225.6926	C ₁₀ ² H ₉ H ₆ ClN ₄		
DRE-C10013100	Acetamiprid D3 (N-methyl D3)		50mg	
DRE-XA10013100AC	Acetamiprid D3 (N-methyl D3) 100 µg/mL in Acetone		1.1ml	
Acetamiprid-N-desmethyl				
CAS 190604-92-3	MW 208.6476	C ₉ H ₉ ClN ₄		
DRE-C10013200	Acetamiprid-N-desmethyl(‡)		10mg	
Acethion				
CAS 919-54-0	MW 272.3219	C ₈ H ₁₇ O ₄ PS ₂		
DRE-XA10015000CY	Acethion 100 µg/mL in Cyclohexane		1ml	
Acrinathrin				
CAS 101007-06-1	MW 541.4391	C ₂₆ H ₂₁ F ₆ NO ₅		
DRE-C10044000	Acrinathrin(‡)		100mg	
DRE-L10044000CY	Acrinathrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA10044000CY	Acrinathrin 100 µg/mL in Cyclohexane		1ml	
Acrolein (2-Propenal; 2-Propen-1-one)				
CAS 107-02-8	MW 56.0633	C ₃ H ₄ O		
DRE-XA10045000AC	Acrolein (2-Propenal) 100 µg/mL in Acetone (*)		1ml	
DRE-GA09010399ME	Acrolein Solution, 1000 µg/mL in Methanol(‡)		1ml	
DRE-GA09010384ME	Acrolein Solution, 5000 µg/mL in p/t Methanol(‡)		1ml	
DRE-GA09010384WA	Acrolein Solution, 5000 µg/mL in Water(‡)		1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Acrylonitrile				
CAS 107-13-1 DRE-C10045500	MW 53.0626 Acrylonitrile	C ₃ H ₃ N	5ml	
Akton (O-[2-Chloro-1-(2,5-dichlorophenyl)vinyl] O,O-diethyl phosphorothioate)				
CAS 1757-18-2 DRE-C10050000	MW 375.6355 Akton	C ₁₂ H ₁₄ Cl ₃ O ₃ PS	100mg	
Alanycarb				
CAS 83130-01-2 DRE-C10063000 DRE-XA10063000CY	MW 399.5281 Alanycarb Alanycarb 100 µg/mL in Cyclohexane	C ₁₇ H ₂₈ N ₃ O ₄ S ₂	100mg 1ml	
Aldicarb				
CAS 116-06-3 DRE-C10070000 DRE-XA10070000AL	MW 190.2633 Aldicarb Aldicarb 100 µg/mL in Acetonitrile (*)	C ₇ H ₁₄ N ₂ O ₂ S	100mg 1ml	
Aldicarb D3 (N-methyl D3)				
CAS 1795142-83-4 DRE-C10070100	MW 193.2817 Aldicarb D3 (methyl-D3)	C ₇ ² H ₃ H ₁₁ N ₂ O ₂ S	10mg	
Aldicarb Sulfone				
CAS 1646-88-4 DRE-C10080000	MW 222.2621 Aldicarb-sulfone(‡)	C ₇ H ₁₄ N ₂ O ₄ S	100mg	
Aldicarb Sulfoxide				
CAS 1646-87-3 DRE-C10080500	MW 206.2627 Aldicarb-sulfoxide (*)	C ₇ H ₁₄ N ₂ O ₃ S	100mg	
Aldicarb-oxime				
CAS 1646-75-9 DRE-LA10079500AL	MW 133.2119 Aldicarb-oxime 10 µg/mL in Acetonitrile	C ₈ H ₁₁ NOS	1ml	
Aldicarb-sulfone D3 (N-methyl D3)				
CAS 1795135-15-7 DRE-XA10080100AC	MW 225.2805 Aldicarb-sulfone D3 100 µg/mL in Acetone (*)	C ₇ ² H ₃ H ₁₁ N ₂ O ₄ S	1.1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Aldicarb-sulfoxide-nitrile				
CAS 14668-28-1 DRE-C10080800	MW 131.1961 Aldicarb-sulfoxide-nitrile	C ₅ H ₉ NOS	10mg	
Aldrin (HHDN)				
CAS 309-00-2 DRE-C10090000 DRE-L10090000AL DRE-L10090000CY DRE-XA10090000IO	MW 364.9099 Aldrin (HHDN)(‡) Aldrin 10 µg/mL in Acetonitrile Aldrin 10 µg/mL in Cyclohexane Aldrin 100 µg/mL in Isooctane	C ₁₂ H ₈ Cl ₆	250mg 10ml 10ml 1ml	
Allethrin				
CAS 584-79-2 DRE-C10100000 DRE-L10100000CY	MW 302.4079 Allethrin Allethrin 10 µg/mL in Cyclohexane (*)	C ₁₉ H ₂₆ O ₃	100mg 10ml	
D-Allethrin				
CAS 231937-89-6 DRE-C10100500	MW 302.4079 D-Allethrin	C ₁₉ H ₂₆ O ₃	100mg	
Allyxycarb				
CAS 6392-46-7 DRE-C10141200	MW 274.3581 Allyxycarb	C ₁₆ H ₂₂ N ₂ O ₂	100mg	
Amidithion				
CAS 919-76-6 DRE-LA10160000EA	MW 273.31 Amidithion 10 µg/mL in Ethylacetate(‡)	C ₇ H ₁₆ NO ₄ PS ₂	1ml	
Aminocarb				
CAS 2032-59-9 DRE-C10190000 DRE-XA10190000AL	MW 208.2569 Aminocarb Aminocarb 100 µg/mL in Acetonitrile	C ₁₁ H ₁₆ N ₂ O ₂	250mg 1ml	
Amitraz				
CAS 33089-61-1 DRE-C10230000 DRE-L10230000CY DRE-XA10230000AL DRE-XA09010148AL	MW 293.406 Amitraz(‡) Amitraz 10 µg/mL in Cyclohexane Amitraz 100 µg/mL in Acetonitrile Amitraz 100 µg/mL in Acetonitrile(‡)	C ₁₉ H ₂₃ N ₃	250mg 10ml 1ml 1ml	
Amitraz D12 (methylphenyl D12)				
CAS n/a DRE-XA10230100AC	MW 305.48 Amitraz D12 100 µg/mL in Acetone(‡)	C ₁₉ ² H ₁₂ H ₁₁ N ₃	1.1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Avermectin B1 (Abamectine)				
CAS 71751-41-2	MW 1732.1272	((C ₄₈ H ₇₂ O ₁₄)c(C ₄₇ H ₇₀ O ₁₄)c)mix		
DRE-CA10001000	Abamectin (*)		100mg	
DRE-GS10001000AL	Abamectine, 100 µg/mL in Acetonitrile		5x1ml	
DRE-XA10001000AL	Abamectine 100 µg/mL in Acetonitrile (*)		1ml	
Azamethiphos				
CAS 35575-96-3	MW 324.6779	C ₉ H ₁₀ ClN ₂ O ₅ PS		
DRE-C10340000	Azamethiphos		250mg	
DRE-XA09010151AC	Azamethiphos 100 µg/mL in Acetone(‡)		1ml	
DRE-XA09010151ME	Azamethiphos 100 µg/mL in Methanol(‡)		1ml	
Azidithion				
CAS 78-57-9	MW 281.2955	C ₆ H ₁₂ N ₅ O ₂ PS ₂		
DRE-C10350000	Azidithion		10mg	
Azinphos-ethyl (Guthion Ethyl)				
CAS 2642-71-9	MW 345.3775	C ₁₂ H ₁₆ N ₃ O ₃ PS ₂		
DRE-C10360000	Azinphos-ethyl		250mg	
DRE-L10360000AL	Azinphos-ethyl 10 µg/mL in Acetonitrile		10ml	
DRE-L10360000IO	Azinphos-ethyl 10 µg/mL in Isooctane		10ml	
DRE-X10360000AL	Azinphos-ethyl 100 µg/mL in Acetonitrile		10ml	
DRE-XA10360000CY	Azinphos-ethyl 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010152AL	Guthion Ethyl 100 µg/mL in Acetonitrile(‡)		1ml	
DRE-XA10360100AC	Azinphos-ethyl D10 100 µg/mL in Acetone		1.1ml	
Azinphos-ethyl D10 (ethyl D10)				
CAS n/a	MW 355.4391	C ₁₂ ² H ₁₆ H ₆ N ₃ O ₃ PS ₂		
DRE-XA10360100AC	Azinphos-ethyl D10 100 µg/mL in Acetone		1.1ml	
Azinphos-methyl				
CAS 86-50-0	MW 317.3243	C ₁₀ H ₁₂ N ₃ O ₃ PS ₂		
DRE-C10365000	Azinphos-methyl(‡)		250mg	
DRE-L10365000CY	Azinphos-methyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA10365000CY	Azinphos-methyl 100 µg/mL in Cyclohexane		1ml	
Azinphos-methyl D6 (dimethyl D6)				
CAS n/a	MW 323.3613	C ₁₀ ² H ₆ H ₆ N ₃ O ₃ PS ₂		
DRE-C10365100	Azinphos-methyl D6 (dimethyl D6)		10mg	
DRE-XA10365100AC	Azinphos-methyl D6 100 µg/mL in Acetone		1.1ml	
Azinphos-methyl-oxon				
CAS 961-22-8	MW 301.2587	C ₁₀ H ₁₂ N ₃ O ₄ PS		
DRE-LA10370000CY	Azinphos-methyl-oxon 10 µg/mL in Cyclohexane (*)		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Bifenthrin alcohol metabolite ((2-Methyl-3-biphenyl)methanol)				
CAS 76350-90-8 DRE-LA10584010CY	MW 198.2604 Bifenthrin alcohol metabolite 10 µg/mL in Cyclohexane	C ₁₄ H ₁₄ O	1ml	
Bioallethrin ((RS)-3-Allyl-2-methyl-4-oxocyclopent-2-enyl (1R)-trans-2,2-dimethyl-3-(2-methylprop-1-enyl) cyclopropanecarboxylate)				
CAS 260359-57-7 DRE-C10610000	MW 302.4079 Bioallethrin	C ₁₉ H ₂₆ O ₃	250mg	
S-Bioallethrin (Esbiol)				
CAS 28434-00-6 DRE-C10611000 DRE-XA09010234AL	MW 302.4079 S-Bioallethrin(±) S-Bioallethrin (Esbiol) 100 µg/mL in Acetonitrile(±)	C ₁₉ H ₂₆ O ₃	100mg 1ml	
Bioresmethrin				
CAS 28434-01-7 DRE-C10620000	MW 338.44 Bioresmethrin(±)	C ₂₂ H ₂₆ O ₃	250mg	
Bis(2-chloroethyl) Ether				
CAS 111-44-4 DRE-CA10651500	MW 143.0117 Bis-(2-chloroethyl)ether	C ₄ H ₈ Cl ₂ O	250mg	
Bis(4-chlorophenyl)acetic Acid (4,4'-DDA)				
CAS 83-05-6 DRE-C12020000	MW 281.134 4,4'-DDA	C ₁₄ H ₁₀ Cl ₂ O ₂	100mg	
2,2-Bis(4-chlorophenyl)ethanol (DDOH)				
CAS 2642-82-2 DRE-C12070000	MW 267.1505 DDOH	C ₁₄ H ₁₂ Cl ₂ O	100mg	
Bis(4-chlorophenyl)methane (4,4'-DDM)				
CAS 101-76-8 DRE-C12051000	MW 237.1245 4,4'-DDM	C ₁₃ H ₁₀ Cl ₂	100mg	
Bromfenvinfos (Bromfenvinphos)				
CAS 33399-00-7 DRE-C10690000	MW 404.0209 Bromfenvinphos	C ₁₂ H ₁₄ BrCl ₂ O ₄ P	100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Bromfenvinfos-methyl (2-Bromo-1-(2,4-dichlorophenyl)vinyl dimethyl phosphate)				
CAS 13104-21-7 DRE-C10690500	MW 375.9678 Bromfenvinfos-methyl	C ₁₀ H ₁₀ BrCl ₂ O ₄ P	100mg	
Bromocyclen				
CAS 1715-40-8 DRE-C10726000 DRE-L10726000IO	MW 393.7473 Bromocyclen(±) Bromocyclen 10 µg/mL in Isooctane	C ₈ H ₈ BrCl ₆	100mg 10ml	
Bromophos (Bromophos-methyl)				
CAS 2104-96-3 DRE-C10745000 DRE-L10745000AL DRE-L10745000CY DRE-XA10745000CY	MW 365.9961 Bromophos-methyl(±) Bromophos-methyl 10 µg/mL in Acetonitrile Bromophos-methyl 10 µg/mL in Cyclohexane Bromophos-methyl 100 µg/mL in Cyclohexane	C ₈ H ₈ BrCl ₂ O ₃ PS	100mg 10ml 10ml 1ml	
Bromophos-ethyl				
CAS 4824-78-6 DRE-C10740000 DRE-L10740000IO DRE-XA10740000CY	MW 394.0492 Bromophos-ethyl Bromophos-ethyl 10 µg/mL in Isooctane Bromophos-ethyl 100 µg/mL in Cyclohexane	C ₁₀ H ₁₂ BrCl ₂ O ₃ PS	100mg 10ml 1ml	
Bromophos-ethyl-oxon				
CAS 4855-63-4 DRE-C10740100	MW 377.9836 Bromophos-ethyl-oxon	C ₁₀ H ₁₂ BrCl ₂ O ₄ P	25mg	
Bromophos-methyl-oxon				
CAS 4855-62-3 DRE-C10745100	MW 349.9305 Bromophos-methyl-oxon	C ₈ H ₈ BrCl ₂ O ₄ P	25mg	
Bromopropylate				
CAS 18181-80-1 DRE-C10762000 DRE-L10762000AL DRE-L10762000CY DRE-XA10762000CY DRE-LA09010257CY	MW 428.1151 Bromopropylate Bromopropylate 10 µg/mL in Acetonitrile Bromopropylate 10 µg/mL in Cyclohexane Bromopropylate 100 µg/mL in Cyclohexane Bromopropylate 10 µg/mL in Cyclohexane(±)	C ₁₇ H ₁₆ Br ₂ O ₃	250mg 10ml 10ml 1ml 10ml	
Bromopropylate D7 (isopropyl D7)				
CAS n/a DRE-XA10762100AC	MW 435.1583 Bromopropylate D7 (isopropyl D7) 100 µg/mL in Acetone	C ₁₇ ² H ₁₇ H ₈ Br ₂ O ₃	1.1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Bufencarb				
CAS 8065-36-9 DRE-LA10830000CY	MW 442.5909 Bufencarb 10 µg/mL in Cyclohexane	2C ₁₃ H ₁₉ NO ₂	1ml	
(EZ)-Buprofezin				
CAS 69327-76-0 DRE-C10854000 DRE-L10854000CY DRE-XA10854000CY	MW 305.4383 Buprofezin(‡) Buprofezin 10 µg/mL in Cyclohexane Buprofezin 100 µg/mL in Cyclohexane	C ₁₆ H ₂₃ N ₃ OS	100mg 10ml 1ml	
Butocarboxim				
CAS 34681-10-2 DRE-C10880000 DRE-L10880000AL	MW 190.2633 Butocarboxim Butocarboxim 10 µg/mL in Acetonitrile	C ₇ H ₁₄ N ₂ O ₂ S	100mg 10ml	
Butocarboxim Sulfoxide				
CAS 34681-24-8 DRE-C10890000 DRE-XA09010157ME	MW 206.2627 Butocarboxim sulfoxide Butocarboxim Sulfoxide 100 µg/mL in Methanol(‡)	C ₇ H ₁₄ N ₂ O ₃ S	100mg 1ml	
Butoxycarboxim				
CAS 34681-23-7 DRE-C10900000	MW 222.2621 Butoxycarboxim	C ₇ H ₁₄ N ₂ O ₄ S	100mg	
Cadusafos				
CAS 95465-99-9 DRE-C10934000 DRE-L10934000AL DRE-L10934000CY DRE-XA10934000IO DRE-XA09010161AL	MW 270.3922 Cadusafos(‡) Cadusafos 10 µg/mL in Acetonitrile Cadusafos 10 µg/mL in Cyclohexane Cadusafos 100 µg/mL in Isooctane Cadusafos 100 µg/mL in Acetonitrile(‡)	C ₁₀ H ₂₃ O ₂ PS ₂	100mg 10ml 10ml 1ml 1ml	
Cadusafos D5 (ethyl-D5)				
CAS 1215380-71-4 DRE-XA10934010AL	MW 275.423 Cadusafos D5 (ethyl-D5) 100 µg/mL in Acetonitrile	C ₁₀ ² H ₅ H ₁₈ O ₂ PS ₂	1.1ml	
Carbanolate (Banol)				
CAS 671-04-5 DRE-L10970000CY	MW 213.6608 Carbanolate (Banol) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₂ ClNO ₂	10ml	
Carbaryl (1-Naphthol N-methylcarbamate)				
CAS 63-25-2 DRE-C10980000 DRE-L10980000CY DRE-XA10980000CY	MW 201.2212 Carbaryl(‡) Carbaryl 10 µg/mL in Cyclohexane Carbaryl 100 µg/mL in Cyclohexane	C ₁₂ H ₁₁ NO ₂	250mg 10ml 1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Carbaryl D3 (methyl D3)				
CAS 1433961-56-8	MW 204.2397	C ₁₂ H ₉ H ₈ NO ₂		
DRE-XA10980010AC	Carbaryl D3 (methyl D3) 100 µg/mL in Acetone		1.1ml	
Carbaryl D7 (naphthyl D7)				
CAS 362049-56-7	MW 208.2644	C ₁₂ H ₇ H ₄ NO ₂		
DRE-C10980100	Carbaryl D7 (naphthyl D7)		50mg	
DRE-XA10980100AC	Carbaryl D7 (naphthyl D7) 100 µg/mL in Acetone		1.1ml	
Carbofuran				
CAS 1563-66-2	MW 221.2524	C ₁₂ H ₁₃ NO ₃		
DRE-C11010000	Carbofuran		250mg	
DRE-L11010000CY	Carbofuran 10 µg/mL in Cyclohexane		10ml	
DRE-XA11010000ME	Carbofuran 100 µg/mL in Methanol		1ml	
DRE-GS09010311AL	Carbofuran 100 µg/mL in Acetonitrile(‡)		5x1ml	
Carbofuran D3 (N-methyl D3)				
CAS 1007459-98-4	MW 224.2709	C ₁₂ H ₃ H ₁₂ NO ₃		
DRE-C11010100	Carbofuran D3 (N-methyl D3)(‡)		10mg	
DRE-XA11010100AC	Carbofuran D3 (N-methyl D3) 100 µg/mL in Acetone		1.1ml	
Carbofuran-3-hydroxy (3-Hydroxycarbofuran)				
CAS 16655-82-6	MW 237.2518	C ₁₂ H ₁₃ NO ₄		
DRE-C11011000	Carbofuran-3-hydroxy(‡)		10mg	
DRE-L11011000EA	Carbofuran-3-hydroxy 10 µg/mL in Ethylacetate		10ml	
DRE-XA11011000EA	Carbofuran-3-hydroxy 100 µg/mL in Ethylacetate		1ml	
DRE-GH09010099AL	3-Hydroxycarbofuran 100 µg/mL in Acetonitrile(‡)		10x1ml	
Carbofuran-3-keto (3-Ketocarbofuran)				
CAS 16709-30-1	MW 235.2359	C ₁₂ H ₁₃ NO ₄		
DRE-C11012000	Carbofuran-3-keto		10mg	
DRE-XA09010146AC	3-Ketocarbofuran 100 µg/mL in Acetone(‡)		1ml	
Carbofuranphenol (2,2-Dimethyl-3H-benzofuran-7-ol)				
CAS 1563-38-8	MW 164.2011	C ₁₀ H ₁₂ O ₂		
DRE-C11012100	Carbofuranphenol		250mg	
DRE-GS09010097ME	Carbofuran Phenol 100 µg/mL in Methanol(‡)		5x1ml	
Carbon Disulfide				
CAS 75-15-0	MW 76.1407	CS ₂		
DRE-CA11016000	Carbon disulfide		250mg	
DRE-L11016000ME	Carbondisulfide 10 µg/mL in Methanol		10ml	
DRE-YA11016000ME	Carbondisulfide, 5000 µg/mL in Methanol		1ml	
Carbophenothion				
CAS 786-19-6	MW 342.8653	C ₁₁ H ₁₆ ClO ₂ PS ₃		
DRE-C11020000	Carbophenothion		250mg	
DRE-L11020000IO	Carbophenothion 10 µg/mL in Isooctane		10ml	
DRE-XA11020000CY	Carbophenothion 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Carbophenothion D10 (di(ethyl D5))				
CAS n/a	MW 352.9269	C ₁₁ H ₁₀ ClO ₂ PS ₃		
DRE-XA11020100AC	Carbophenothion D10 (di(ethyl D5)) 100 µg/mL in Acetone(±)		1.1ml	
Carbophenothion-methyl				
CAS 953-17-3	MW 314.8121	C ₉ H ₁₂ ClO ₂ PS ₃		
DRE-C11020500	Carbophenothion-methyl(±)		10mg	
Carbophenothion-methyl sulfone				
CAS 62059-34-1	MW 346.8109	C ₉ H ₁₂ ClO ₄ PS ₃		
DRE-C11020900	Carbophenothion-methyl sulfone		10mg	
DRE-LA11020900CY	Carbophenothion-methyl sulfone 10 µg/mL in Cyclohexane(±)		1ml	
Carbophenothion-methyl sulfoxide				
CAS 62059-33-0	MW 330.8115	C ₉ H ₁₂ ClO ₃ PS ₃		
DRE-LA11021100CY	Carbophenothion-methyl sulfoxide 10 µg/mL in Cyclohexane		1ml	
Carbophenothion-oxon				
CAS 7173-84-4	MW 326.7997	C ₁₁ H ₁₆ ClO ₃ PS ₂		
DRE-C11020150	Carbophenothion-oxon		25mg	
Carbophenothion-sulfone				
CAS 16662-85-4	MW 374.8641	C ₁₁ H ₁₆ ClO ₄ PS ₃		
DRE-C11020210	Carbophenothion-sulfone		25mg	
Carbophenothion-sulfoxide				
CAS 17297-40-4	MW 358.8647	C ₁₁ H ₁₆ ClO ₃ PS ₃		
DRE-C11020220	Carbophenothion-sulfoxide		25mg	
Carbosulfan				
CAS 55285-14-8	MW 380.5447	C ₂₀ H ₃₂ N ₂ O ₃ S		
DRE-C11030000	Carbosulfan(±)		250mg	
DRE-L11030000IO	Carbosulfan 10 µg/mL in Isooctane		10ml	
DRE-XA11030000IO	Carbosulfan 100 µg/mL in Isooctane		1ml	
Cartap				
CAS 15263-53-3	MW 237.3429	C ₇ H ₁₅ N ₃ O ₂ S ₂		
DRE-XA09010163ME	Cartap 100 µg/mL in Methanol(±)		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Cartap Hydrochloride				
CAS 15263-52-2 DRE-C11050000	MW 273.8038 Cartap hydrochloride	C ₇ H ₁₅ N ₃ O ₂ S ₂ ·ClH	100mg	
Chinomethionat				
CAS 2439-01-2 DRE-C11080000 DRE-L11080000CY DRE-XA11080000CY	MW 234.2974 Chinomethionat Chinomethionat 10 µg/mL in Cyclohexane Chinomethionat 100 µg/mL in Cyclohexane	C ₁₀ H ₆ N ₂ OS ₂	250mg 10ml 1ml	
Chlorantraniliprole				
CAS 500008-45-7 DRE-C11145000	MW 483.1461 Chlorantraniliprole(±)	C ₁₈ H ₁₄ BrCl ₂ N ₅ O ₂	100mg	
Chlorbenside				
CAS 103-17-3 DRE-C11150000 DRE-L11150000CY	MW 269.1895 Chlorbenside Chlorbenside 10 µg/mL in Cyclohexane	C ₁₃ H ₁₀ Cl ₂ S	100mg 10ml	
Chlorbenside Sulfone				
CAS 7082-99-7 DRE-C11155000	MW 301.1883 Chlorbenside sulfone	C ₁₃ H ₁₀ Cl ₂ O ₂ S	10mg	
Chlorbicyclen				
CAS 2550-75-6 DRE-C11170000	MW 397.7679 Chlorbicyclen	C ₉ H ₆ Cl ₈	250mg	
Chlordane				
CAS 57-74-9 DRE-C11200000 DRE-GA11200000HE DRE-GH11200000HE DRE-L11200000CY DRE-YA11200000ME	MW 409.7786 Chlordane (technical) Chlordane (technical) 1000 µg/mL in n-Hexane Chlordane (technical) 1000 µg/mL in Hexane Chlordane (technical) 10 µg/mL in Cyclohexane Chlordane (technical) 2000 µg/mL in Methanol	C ₁₀ H ₆ Cl ₈	250mg 1ml 5x1ml 10ml 1ml	
cis-Chlordane (α)				
CAS 5103-71-9 DRE-C11201000 DRE-L11201000CY DRE-XA11201000CY	MW 409.7786 cis-Chlordane (α)(±) cis-Chlordane (α) 10 µg/mL in Cyclohexane cis-Chlordane (α) 100 µg/mL in Cyclohexane	C ₁₀ H ₆ Cl ₈	10mg 10ml 1ml	
trans-Chlordane (γ)				
CAS 5103-74-2 DRE-C11202000 DRE-L11202000CY DRE-XA11202000CY	MW 409.7786 trans-Chlordane (γ)(±) trans-Chlordane (γ) 10 µg/mL in Cyclohexane(±) trans-Chlordane (γ) 100 µg/mL in Cyclohexane(±)	C ₁₀ H ₆ Cl ₈	10mg 10ml 1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
oxy-Chlordane				
CAS 27304-13-8	MW 423.7622	C ₁₀ H ₄ Cl ₈ O		AND Enantiomer
DRE-L11203000CY	oxy-Chlordane 10 µg/mL in Cyclohexane		10ml	
DRE-LA11203000AL	oxy-Chlordane 10 µg/mL in Acetonitrile		1ml	
DRE-LA11203000CY	oxy-Chlordane 10 µg/mL in Cyclohexane(‡)		1ml	
DRE-XA11203000CY	oxy-Chlordane 100 µg/mL in Cyclohexane		1ml	
Chlordecone				
CAS 143-50-0	MW 490.6364	C ₁₀ Cl ₁₀ O		AND Enantiomer
DRE-C11220000	Chlordecone		100mg	
DRE-L11220000IO	Chlordecone 10 µg/mL in Isooctane		10ml	
Chlordecone-5b-hydro				
CAS 53308-47-7	MW 456.1913	C ₁₀ HCl ₉ O		AND Enantiomer
DRE-LA11220200CY	Chlordecone-5b-hydro 10 µg/mL in Cyclohexane		1ml	
Chlordene				
CAS 3734-48-3	MW 338.8726	C ₁₀ H ₆ Cl ₆		AND Enantiomer
DRE-C11230000	Chlordene		10mg	
β-Chlordene				
CAS 56534-03-3	MW 338.8726	C ₁₀ H ₆ Cl ₆		AND Enantiomer
DRE-LA11230600CY	β-Chlordene 10 µg/mL in Cyclohexane		1ml	
cis-Chlordene				
CAS 56534-02-2	MW 338.8726	C ₁₀ H ₆ Cl ₆		AND Enantiomer
DRE-LA11230700IO	cis-Chlordene 10 µg/mL in Isooctane		1ml	
oxy-Chlordene				
CAS 6058-23-7	MW 354.872	C ₁₀ H ₆ Cl ₆ O		AND Enantiomer
DRE-LA11232000CY	oxy-Chlordene 10 µg/mL in Cyclohexane		1ml	
Chlordimeform				
CAS 6164-98-3	MW 196.6766	C ₁₀ H ₁₃ ClN ₂		
DRE-C11240000	Chlordimeform free base		100mg	
DRE-L11240000CY	Chlordimeform free base 10 µg/mL in Cyclohexane		10ml	
Chlordimeform Hydrochloride				
CAS 19750-95-9	MW 233.1376	C ₁₀ H ₁₃ ClN ₂ ·ClH		
DRE-C11240100	Chlordimeform HCl		100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Chlorfenapyr				
CAS 122453-73-0 DRE-C11247500 DRE-L11247500CY	MW 407.6128 Chlorfenapyr(±) Chlorfenapyr 10 µg/mL in Cyclohexane	C ₁₅ H ₁₁ BrClF ₃ N ₂ O	100mg 10ml	
Chlorfenethol				
CAS 80-06-8 DRE-C11250000	MW 267.1505 Chlorfenethol	C ₁₄ H ₁₂ Cl ₂ O	50mg	
Chlorfenson				
CAS 80-33-1 DRE-C11270000 DRE-L11270000IO	MW 303.1611 Chlorfenson Chlorfenson 10 µg/mL in Isooctane	C ₁₂ H ₈ Cl ₂ O ₃ S	100mg 10ml	
Chlorfenvinphos				
CAS 470-90-6 DRE-C11290000 DRE-L11290000AL DRE-L11290000CY DRE-XA11290000AL DRE-XA11290000CY	MW 359.5699 Chlorfenvinphos Chlorfenvinphos 10 µg/mL in Acetonitrile Chlorfenvinphos 10 µg/mL in Cyclohexane Chlorfenvinphos 100 µg/mL in Acetonitrile Chlorfenvinphos 100 µg/mL in Cyclohexane	C ₁₂ H ₁₄ Cl ₃ O ₄ P	250mg 10ml 10ml 1ml 1ml	
Chlorfenvinphos D10 (ethyl D10)				
CAS 1346606-54-9 DRE-XA11290100AC	MW 369.6315 Chlorfenvinphos (ethyl) D10 100 µg/mL in Acetone	C ₁₂ ² H ₁₀ H ₄ Cl ₃ O ₄ P	1.1ml	
cis-Chlorfenvinphos				
CAS 18708-87-7 DRE-LA11291000CY	MW 359.5699 cis-Chlorfenvinphos 10 µg/mL in Cyclohexane(±)	C ₁₂ H ₁₄ Cl ₃ O ₄ P	1ml	
Chlorfluazuron				
CAS 71422-67-8 DRE-C11297000 DRE-L11297000CY	MW 540.6548 Chlorfluazuron(±) Chlorfluazuron 10 µg/mL in Cyclohexane	C ₂₀ H ₉ Cl ₃ F ₅ N ₃ O ₃	100mg 10ml	
Chlormephos				
CAS 24934-91-6 DRE-C11330000	MW 234.7043 Chlormephos	C ₈ H ₁₂ ClO ₂ PS ₂	100mg	
Chlorobenzilate				
CAS 510-15-6 DRE-C11390000 DRE-L11390000IO DRE-XA11390000CY	MW 325.1866 Chlorobenzilate Chlorobenzilate 10 µg/mL in Isooctane Chlorobenzilate 100 µg/mL in Cyclohexane	C ₁₆ H ₁₄ Cl ₂ O ₃	100mg 10ml 1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Chlorobenzuron				
CAS 57160-47-1 DRE-C11392500	MW 309.1474 Chlorobenzuron	C ₁₄ H ₁₀ Cl ₂ N ₂ O ₂	100mg	
Chloropropylate				
CAS 5836-10-2 DRE-C11503400 DRE-L11503400IO	MW 339.2131 Chloropropylate Chloropropylate 10 µg/mL in Isooctane	C ₁₇ H ₁₆ Cl ₂ O ₃	100mg 10ml	
Chlorphoxim				
CAS 14816-20-7 DRE-C11570000	MW 332.7429 Chlorphoxim	C ₁₂ H ₁₄ ClN ₂ O ₃ PS	100mg	
Chlorpyrifos				
CAS 2921-88-2 DRE-C11600000 DRE-CR11600000 DRE-L11600000AL DRE-L11600000CY DRE-XA11600000AL DRE-XA11600000CY DRE-GA09010366AC	MW 350.5863 Chlorpyrifos(‡) Chlorpyrifos Chlorpyrifos 10 µg/mL in Acetonitrile Chlorpyrifos 10 µg/mL in Cyclohexane Chlorpyrifos 100 µg/mL in Acetonitrile Chlorpyrifos 100 µg/mL in Cyclohexane Chlorpyrifos 100 µg/mL in Acetone(‡)	C ₉ H ₁₁ Cl ₃ NO ₃ PS	250mg 100mg 10ml 10ml 1ml 1ml 1ml	
Chlorpyrifos (diethyl-D10)				
CAS 285138-81-0 DRE-C11600100 DRE-XA11600100AC	MW 360.6479 Chlorpyrifos D10(‡) Chlorpyrifos D10 100 µg/mL in Acetone(‡)	C ₉ ² H ₁₀ HCl ₃ NO ₃ PS	25mg 1ml	
Chlorpyrifos-methyl				
CAS 5598-13-0 DRE-C11601000 DRE-L11601000CY DRE-XA11601000CY	MW 322.5331 Chlorpyrifos-methyl(‡) Chlorpyrifos-methyl 10 µg/mL in Cyclohexane Chlorpyrifos-methyl 100 µg/mL in Cyclohexane	C ₇ H ₇ Cl ₃ NO ₃ PS	250mg 10ml 1ml	
Chlorpyrifos-methyl D6				
CAS n/a DRE-C11601100 DRE-XA11601100AC	MW 328.5701 Chlorpyrifos-methyl D6 Chlorpyrifos-methyl D6 100 µg/mL in Acetone (*)	C ₇ ² H ₆ HCl ₃ NO ₃ PS	10mg 1.1ml	
Chlorpyrifos-methyl-oxon				
CAS 5598-52-7 DRE-C11601500	MW 306.4675 Chlorpyrifos-methyl-oxon	C ₇ H ₇ Cl ₃ NO ₄ P	25mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Chlorpyrifos-oxon				
CAS 5598-15-2	MW 334.5207	C ₉ H ₁₁ Cl ₃ NO ₄ P		
DRE-C11603000	Chlorpyrifos-oxon		50mg	
DRE-XA09010167ME	Chlorpyrifos-oxon 100 µg/mL in Methanol(±)		1ml	
Chlorthion				
CAS 500-28-7	MW 297.6525	C ₈ H ₉ ClNO ₃ PS		
DRE-C11640000	Chlorthion		25mg	
DRE-L11640000IO	Chlorthion 10 µg/mL in Isooctane		10ml	
Chlorthiophos				
CAS 60238-56-4	MW 1083.7343	((C ₁₁ H ₁₅ Cl ₂ O ₃ PS ₂)c(C ₁₁ H ₁₅ Cl ₂ O ₃ PS ₂)c(C ₁₁ H ₁₅ Cl ₂ O ₃ PS ₂)c)mix		
DRE-C11650000	Chlorthiophos		100mg	
Chlorthiophos-I-D10 (diethyl D10)				
CAS n/a	MW 371.3064	C ₁₁ ² H ₁₀ H ₅ Cl ₂ O ₃ PS ₂		
DRE-C11650010	Chlorthiophos D10 (diethyl D10)		10mg	
Chlorthiophos I sulfone				
CAS 25900-20-3	MW 393.2436	C ₁₁ H ₁₅ Cl ₂ O ₃ PS ₂		
DRE-LA11651000CY	Chlorthiophos sulfone 10 µg/mL in Cyclohexane		1ml	
Chromafenozide				
CAS 143807-66-3	MW 394.5066	C ₂₄ H ₃₀ N ₂ O ₃		
DRE-C11665500	Chromafenozide(±)		100mg	
DRE-XA11665500AL	Chromafenozide 100 µg/mL in Acetonitrile		1ml	
Chrysanthemic Acid				
CAS 10453-89-1	MW 168.2328	C ₁₀ H ₁₆ O ₂		
DRE-L11666000CY	Chrysanthemic acid 10 µg/mL in Cyclohexane		10ml	
Cloethocarb				
CAS 51487-69-5	MW 259.6862	C ₁₁ H ₁₄ ClNO ₄		
DRE-C11679300	Cloethocarb		50mg	
Clofentezine				
CAS 74115-24-5	MW 303.1461	C ₁₄ H ₈ Cl ₂ N ₄		
DRE-C11680000	Clofentezine		100mg	
DRE-L11680000CY	Clofentezine 10 µg/mL in Cyclohexane		10ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Clothianidin				
CAS 210880-92-5 DRE-C11691700	MW 249.678 Clothianidin(±)	C ₈ H ₈ ClN ₃ O ₂ S	100mg	
Clothianidin D3 (N'-methyl D3)				
CAS 1262776-24-8 DRE-C11691710	MW 252.6965 Clothianidin D3 (N'-methyl D3)	C ₈ H ₇ ClN ₃ O ₂ S	50mg	
Clothianidin Urea				
CAS 634192-72-6 DRE-C11691730	MW 205.6652 Clothianidin Urea	C ₈ H ₈ ClN ₃ OS	25mg	
Coumaphos				
CAS 56-72-4 DRE-C11730000 DRE-L11730000AL DRE-XA11730000AL	MW 362.7656 Coumaphos Coumaphos 10 µg/mL in Acetonitrile Coumaphos 100 µg/mL in Acetonitrile	C ₁₄ H ₁₆ ClO ₅ PS	100mg 10ml 1ml	
Coumaphos D10 (di(ethyl-D5))				
CAS 287397-86-8 DRE-XA11730010AL	MW 372.8272 Coumaphos D10 di(ethyl-D5) 100 µg/mL in Acetonitrile	C ₁₄ ² H ₁₀ H ₆ ClO ₅ PS	1.1ml	
Coumaphos alcohol metabolite				
CAS 6174-86-3 DRE-C11730030	MW 210.6138 Coumaphos alcohol metabolite	C ₁₀ H ₇ ClO ₃	100mg	
Coumaphos-oxon				
CAS 321-54-0 DRE-C11731000	MW 346.7 Coumaphos-oxon	C ₁₄ H ₁₆ ClO ₆ P	10mg	
Crotoxyphos				
CAS 7700-17-6 DRE-LA11760000IO DRE-XA09010172ME	MW 314.2708 Crotoxyphos 10 µg/mL in Isooctane (*) Crotoxyphos 100 µg/mL in Methanol(±)	C ₁₄ H ₁₈ O ₆ P	1ml 1ml	
Cruformate				
CAS 299-86-5 DRE-C11770000 DRE-XA09010173AL	MW 291.7109 Cruformate Cruformate 100 µg/mL in Acetonitrile(±)	C ₁₂ H ₁₉ ClNO ₃ P	50mg 1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Cyanofenphos				
CAS 13067-93-1	MW 303.3159	C ₁₅ H ₁₄ NO ₂ PS		
DRE-C11800000	Cyanofenphos		100mg	
DRE-L11800000IO	Cyanofenphos 10 µg/mL in Isooctane		10ml	
Cyanophos				
CAS 2636-26-2	MW 243.2194	C ₉ H ₁₀ NO ₃ PS		
DRE-C11810000	Cyanophos		25mg	
DRE-L11810000IO	Cyanophos 10 µg/mL in Isooctane		10ml	
Cycloprothrin				
CAS 63935-38-6	MW 482.3552	C ₂₆ H ₂₁ Cl ₂ NO ₄		
DRE-C11836000	Cycloprothrin		10mg	
DRE-L11836000CY	Cycloprothrin 10 µg/mL in Cyclohexane		10ml	
Cycloxaprid				
CAS 1203791-41-6	MW 322.7469	C ₁₄ H ₁₅ ClN ₄ O ₃		
DRE-C11836700	Cycloxaprid		10mg	
Cyenopyrafen				
CAS 560121-52-0	MW 393.5218	C ₂₄ H ₃₁ N ₃ O ₂		
DRE-C11841500	Cyenopyrafen		25mg	
DRE-XA11841500AL	Cyenopyrafen 100 µg/mL in Acetonitrile(±)		1ml	
Cyflumetofen				
CAS 400882-07-7	MW 447.4469	C ₂₄ H ₂₄ F ₃ NO ₄		
DRE-CA11846000	Cyflumetofen		100mg	
Cyfluthrin				
CAS 68359-37-5	MW 434.2876	C ₂₂ H ₁₈ Cl ₂ FNO ₃		
DRE-C11850000	Cyfluthrin(±)		250mg	
DRE-L11850000AL	Cyfluthrin 10 µg/mL in Acetonitrile		10ml	
DRE-L11850000CY	Cyfluthrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA11850000CY	Cyfluthrin 100 µg/mL in Cyclohexane		1ml	
β-Cyfluthrin				
CAS 1820573-27-0	MW 434.2876	C ₂₂ H ₁₈ Cl ₂ FNO ₃		
DRE-C11850200	β-Cyfluthrin		250mg	
DRE-L11850200CY	β-Cyfluthrin 10 µg/mL in Cyclohexane		10ml	
trans-Cyfluthrin D6 (2,2-dimethyl D6)				
CAS n/a	MW 440.3246	C ₂₂ ² H ₁₈ H ₁₂ Cl ₂ FNO ₃		
DRE-XA11850010AL	trans-Cyfluthrin D6 (2,2-dimethyl D6) 100 µg/mL in Acetonitrile		1.1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
γ-Cyhalothrin				
CAS 76703-62-3	MW 449.8501	C ₂₃ H ₁₉ ClF ₃ NO ₃		
DRE-C11859500	γ-Cyhalothrin(±)		100mg	
DRE-L11859500CY	γ-Cyhalothrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA11859500CY	γ-Cyhalothrin 100 µg/mL in Cyclohexane		1ml	
λ-Cyhalothrin				
CAS 91465-08-6	MW 449.8501	C ₂₃ H ₁₉ ClF ₃ NO ₃		
DRE-C11860000	λ-Cyhalothrin(±)		100mg	
DRE-L11860000CY	λ-Cyhalothrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA11860000CY	λ-Cyhalothrin 100 µg/mL in Cyclohexane		1ml	
Cyhexatin				
CAS 13121-70-5	MW 385.172	C ₁₈ H ₃₄ OSn		
DRE-C11870000	Cyhexatin		250mg	
Cypermethrin				
CAS 52315-07-8	MW 416.2972	C ₂₂ H ₁₉ Cl ₂ NO ₃		
DRE-C11890000	Cypermethrin (technical)(±)		100mg	
DRE-CR11890000	Cypermethrin (technical)		50mg	
DRE-GA11890000AL	Cypermethrin (technical) Solution in Acetonitrile		1ml	
DRE-GS11890000AL	Cypermethrin (technical) Solution in Acetonitrile		5x1ml	
DRE-L11890000AL	Cypermethrin (technical) 10 µg/mL in Acetonitrile		10ml	
DRE-L11890000IO	Cypermethrin (technical) 10 µg/mL in Isooctane		10ml	
DRE-XA11890000CY	Cypermethrin (technical) 100 µg/mL in Cyclohexane		1ml	
α-Cypermethrin				
CAS 67375-30-8	MW 416.2972	C ₂₂ H ₁₉ Cl ₂ NO ₃		
DRE-C11890100	α-Cypermethrin(±)		100mg	
DRE-L11890100CY	α-Cypermethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA11890100CY	α-Cypermethrin 100 µg/mL in Cyclohexane		1ml	
β-Cypermethrin				
CAS 1224510-29-5	MW 416.2972	C ₂₂ H ₁₉ Cl ₂ NO ₃		
DRE-C11890200	β-Cypermethrin		100mg	
DRE-L11890200CY	β-Cypermethrin 10 µg/mL in Cyclohexane		10ml	
θ-Cypermethrin				
CAS 71697-59-1	MW 416.2972	C ₂₂ H ₁₉ Cl ₂ NO ₃		
DRE-L11890300CY	θ-Cypermethrin 10 µg/mL in Cyclohexane		10ml	
trans-Cypermethrin D6 (dimethyl D6)				
CAS 82523-65-7	MW 422.3341	C ₂₂ ² H ₆ H ₁₃ Cl ₂ NO ₃		
DRE-XA11890400AC	trans-Cypermethrin D6 (dimethyl D6) 100 µg/mL in Acetone(±)		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
ζ-Cypermethrin				
CAS 1315501-18-8	MW 416.2972	C ₂₂ H ₁₉ Cl ₂ NO ₃		
DRE-C11890500	ζ-Cypermethrin		100mg	
DRE-L11890500CY	ζ-Cypermethrin 10 µg/mL in Cyclohexane		10ml	
Cyphenothrin				
CAS 39515-40-7	MW 375.4602	C ₂₄ H ₂₅ NO ₃		
DRE-C11895000	Cyphenothrin		100mg	
DRE-L11895000CY	Cyphenothrin 10 µg/mL in Cyclohexane		10ml	
Cyromazine				
CAS 66215-27-8	MW 166.1838	C ₆ H ₁₀ N ₆		
DRE-C11920000	Cyromazine		250mg	
Cyromazine D4				
CAS 1219804-19-9	MW 170.2084	C ₆ ² H ₄ H ₆ N ₆		
DRE-C11920010	Cyromazine D4 (cyclopropyl-2,2,3,3 D4)		10mg	
DRE-XA11920010EA	Cyromazine D4 (cyclopropyl-2,2,3,3 D4) 100 µg/mL in Ethylacetate (*)		1.1ml	
Cythioate				
CAS 115-93-5	MW 297.2883	C ₈ H ₁₂ NO ₅ PS ₂		
DRE-C11930000	Cythioate		10mg	
4,4'-DDD (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethane)				
CAS 72-54-8	MW 320.0412	C ₁₄ H ₁₀ Cl ₄		
DRE-C12031000	4,4'-DDD(‡)		250mg	
DRE-L12031000CY	4,4'-DDD 10 µg/mL in Cyclohexane		10ml	
DRE-XA12031000CY	4,4'-DDD 100 µg/mL in Cyclohexane		1ml	
2,4'-DDE (1,1-Dichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethene)				
CAS 3424-82-6	MW 318.0253	C ₁₄ H ₈ Cl ₄		
DRE-C12040000	2,4'-DDE(‡)		50mg	
DRE-L12040000CY	2,4'-DDE 10 µg/mL in Cyclohexane		10ml	
DRE-XA12040000CY	2,4'-DDE 100 µg/mL in Cyclohexane		1ml	
2,4'-DDE D8 (1,1-Dichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethene D8)				
CAS 1402834-57-4	MW 326.0746	C ₁₄ ² H ₈ Cl ₄		
DRE-XA12040100AC	2,4'-DDE D8 100 µg/mL in Acetone		1.1ml	
4,4'-DDE (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethene)				
CAS 72-55-9	MW 318.0253	C ₁₄ H ₈ Cl ₄		
DRE-C12041000	4,4'-DDE(‡)		100mg	
DRE-L12041000AL	4,4'-DDE 10 µg/mL in Acetonitrile		10ml	
DRE-L12041000CY	4,4'-DDE 10 µg/mL in Cyclohexane		10ml	
DRE-XA12041000CY	4,4'-DDE 100 µg/mL in Cyclohexane(‡)		1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
4,4'-DDE D8 (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethene D8)				
CAS 93952-19-3 DRE-XA12041100AC	MW 326.0746 4,4'-DDE D8 100 µg/mL in Acetone(‡)	C ₁₄ ² H ₈ Cl ₄	1ml	
4,4'-DDMU (1,1-Bis(p-chlorophenyl)-2-chloroethene)				
CAS 1022-22-6 DRE-C12061000	MW 283.5803 4,4'-DDMU	C ₁₄ H ₉ Cl ₃	100mg	
4,4'-DDNU (1,1-Bis(p-chlorophenyl)ethene)				
CAS 2642-81-1 DRE-C12062000	MW 249.1352 4,4'-DDNU	C ₁₄ H ₁₀ Cl ₂	100mg	
4,4'-DDT (ring-13C12)				
CAS 104215-84-1 DRE-XA12082200AC	MW 366.3981 4,4'-DDT 13C12 100 µg/mL in Acetone(‡)	¹³ C ₁₂ C ₂ H ₉ Cl ₅	1.1ml	
DDT (technical)				
CAS 8017-34-3 DRE-C12080000 DRE-L12080000CY	MW 708.9725 DDT (technical) DDT (technical) 10 µg/mL in Cyclohexane	((C ₁₄ H ₉ Cl ₅)c(C ₁₄ H ₉ Cl ₅)c)mix	250mg 10ml	
2,4'-DDT 13C12 (benzen 13C12)				
CAS 1396995-26-8 DRE-XA12081200AC	MW 366.3981 2,4'-DDT 13C12 100 µg/mL in Acetone	¹³ C ₁₂ C ₂ H ₉ Cl ₅	1.1ml	
2,4'-DDT (1,1,1-Trichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethane)				
CAS 789-02-6 DRE-C12081000 DRE-L12081000CY DRE-XA12081000CY	MW 354.4863 2,4'-DDT(‡) 2,4'-DDT 10 µg/mL in Cyclohexane 2,4'-DDT 100 µg/mL in Cyclohexane	C ₁₄ H ₉ Cl ₅	50mg 10ml 1ml	
2,4'-DDT D8 (1,1,1-Trichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethane D8 (ring D8))				
CAS 221899-88-3 DRE-C12081100 DRE-XA12081100AC	MW 362.5356 2,4'-DDT D8 2,4'-DDT D8 100 µg/mL in Acetone	C ₁₄ ² H ₈ HCl ₅	5mg 1.1ml	
4,4'-DDT (1,1,1-Trichloro-2,2-bis(p-chlorophenyl)ethane)				
CAS 50-29-3 DRE-C12082000 DRE-L12082000CY DRE-XA12082000CY	MW 354.4863 4,4'-DDT(‡) 4,4'-DDT 10 µg/mL in Cyclohexane 4,4'-DDT 100 µg/mL in Cyclohexane	C ₁₄ H ₉ Cl ₅	100mg 10ml 1ml	

Pesticide and metabolites: Insecticides

Product code	Description			
4,4'-DDT D8 (1,1,1-Trichloro-2,2-bis(p-chlorophenyl)ethane D8 (ring D8))				
CAS 93952-18-2	MW 362.5356	C ₁₄ H ₈ Cl ₅		
DRE-C12082100	4,4'-DDT D8(‡)		10mg	
DRE-XA12082100AC	4,4'-DDT D8 100 µg/mL in Acetone		1.1ml	
DRE-XA12082100CY	4,4'-DDT D8 100 µg/mL in Cyclohexane		1.1ml	
Deltamethrin				
CAS 52918-63-5	MW 505.1992	C ₂₂ H ₁₉ Br ₂ NO ₃		
DRE-C12120000	Deltamethrin(‡)		250mg	
DRE-L12120000AL	Deltamethrin 10 µg/mL in Acetonitrile		10ml	
DRE-L12120000CY	Deltamethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA12120000CY	Deltamethrin 100 µg/mL in Cyclohexane		1ml	
trans-Deltamethrin D6 (dimethyl D6)				
CAS n/a	MW 511.2361	C ₂₂ H ₁₆ H ₁₃ Br ₂ NO ₃		
DRE-C12120100	trans-Deltamethrin D6 (dimethyl D6)		10mg	
DRE-XA12120100AL	trans-Deltamethrin D6 (dimethyl D6) 100 µg/mL in Acetonitrile(‡)		1.1ml	
Demeton (O+S)				
CAS 8065-48-3	MW 516.6768	2C ₈ H ₁₉ O ₃ PS ₂		
DRE-C12140000	Demeton (O+S)		100mg	
DRE-L12140000CY	Demeton (O+S) 10 µg/mL in Cyclohexane		10ml	
DRE-XA12140000AL	Demeton (O+S) 100 µg/mL in Acetonitrile		1ml	
Demeton-O				
CAS 298-03-3	MW 258.3384	C ₈ H ₁₉ O ₃ PS ₂		
DRE-L12141000CY	Demeton-O 10 µg/mL in Cyclohexane (*)		10ml	
Demeton-S				
CAS 126-75-0	MW 258.3384	C ₈ H ₁₉ O ₃ PS ₂		
DRE-C12142000	Demeton-S(‡)		100mg	
DRE-L12142000IO	Demeton-S 10 µg/mL in Isooctane		10ml	
Demeton-S-methyl				
CAS 919-86-8	MW 230.2853	C ₈ H ₁₅ O ₃ PS ₂		
DRE-CA12143000	Demeton-S-methyl(‡) (*)		50mg	
DRE-XA12143000AL	Demeton-S-methyl 100 µg/mL in Acetonitrile		1ml	
Demeton-S-methyl D6 (dimethyl D6)				
CAS n/a	MW 236.3222	C ₈ ² H ₆ H ₉ O ₃ PS ₂		
DRE-XA12143100CY	Demeton-S-methyl D6 (dimethyl D6) 100 µg/mL in Cyclohexane (*)		1.1ml	
Demeton-S-methyl sulfone				
CAS 17040-19-6	MW 262.2841	C ₈ H ₁₅ O ₅ PS ₂		
DRE-C12144000	Demeton-S-methyl sulfone		100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Demeton-S-methyl sulfoxide (Oxydemeton-methyl)				
CAS 301-12-2	MW 246.2847	C ₆ H ₁₅ O ₄ PS ₂		
DRE-CA12145000	Demeton-S-methyl sulfoxide(±)		25mg	
DRE-L12145000AL	Demeton-S-methyl sulfoxide 10 µg/mL in Acetonitrile (*)		10ml	
DRE-XA12145000CT	Demeton-S-methyl sulfoxide 100 µg/mL in Cyclohexane/Toluene		1ml	
DRE-XA09010218AL	Oxydemeton-methyl 100 µg/mL in Acetonitrile(±)		1ml	
Demeton-S-methyl-sulfoxide D6 (di-O-methyl D6)				
CAS n/a	MW 252.3216	C ₆ ² H ₆ H ₉ O ₄ PS ₂		
DRE-XA12145010AC	Demeton-S-methyl-sulfoxide D6 (di-O-methyl D6) 100 µg/mL in Acetonitrile		1.1ml	
Diafenthiuron				
CAS 80060-09-9	MW 384.578	C ₂₃ H ₃₂ N ₂ OS		
DRE-C12177000	Diafenthiuron		100mg	
Dialifos				
CAS 10311-84-9	MW 393.8458	C ₁₄ H ₁₇ ClNO ₄ PS ₂		
DRE-C12180000	Dialifos(±)		100mg	
DRE-L12180000CY	Dialifos 10 µg/mL in Cyclohexane		10ml	
Diazinon				
CAS 333-41-5	MW 304.3455	C ₁₂ H ₂₁ N ₂ O ₃ PS		
DRE-C12210000	Diazinon(±)		250mg	
DRE-L12210000AL	Diazinon 10 µg/mL in Acetonitrile		10ml	
DRE-L12210000CY	Diazinon 10 µg/mL in Cyclohexane		10ml	
DRE-XA12210000AL	Diazinon 100 µg/mL in Acetonitrile		1ml	
DRE-XA12210000CY	Diazinon 100 µg/mL in Cyclohexane		1ml	
Diazinon D10 (diethyl D10)				
CAS 100155-47-3	MW 314.4071	C ₁₂ ² H ₁₀ H ₁₁ N ₂ O ₃ PS		
DRE-C12210100	Diazinon D10 (diethyl D10)(±)		10mg	
DRE-XA12210100AC	Diazinon D10 (diethyl D10) 100 µg/mL in Acetone		1.1ml	
3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropane)carboxylic				
CAS 63597-73-9	MW 297.9718	C ₈ H ₁₀ Br ₂ O ₂		
DRE-LA12244000ME	3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropane)carboxylic 10 µg/mL in Methanol		1ml	 AND Enantiomer
Dicamba-methyl ester				
CAS 6597-78-0	MW 235.064	C ₉ H ₈ Cl ₂ O ₃		
DRE-C12261000	Dicamba-methyl ester		100mg	









Pesticide and metabolites: Insecticides

Product code	Description			
Dicapthon				
CAS 2463-84-5 DRE-C12270000	MW 297.6525 Dicapthon	C ₈ H ₉ ClNO ₅ PS	100mg	
Dichlobenil				
CAS 1194-65-6 DRE-XA12280000AL DRE-XA12280000CY	MW 172.0114 Dichlobenil 100 µg/mL in Acetonitrile Dichlobenil 100 µg/mL in Cyclohexane	C ₇ H ₃ Cl ₂ N	1ml 1ml	
Dichlofenthion				
CAS 97-17-6 DRE-C12290000 DRE-L12290000CY	MW 315.1532 Dichlofenthion Dichlofenthion 10 µg/mL in Cyclohexane	C ₁₀ H ₁₃ Cl ₂ O ₃ PS	100mg 10ml	
(E)-3-(2,2-Dichlorovinyl)-2,2-di(methyl D3)-(1-cyclopropane)carboxylic acid D6				
CAS n/a DRE-XA12507510AC	MW 215.1068 (E)-3-(2,2-Dichlorovinyl)-2,2-di(methyl D3)-(1-cyclopropane)carboxylic acid D6 100 µg/mL in Acetone	C ₈ ² H ₆ H ₄ Cl ₂ O ₂	1ml	
Dichlorvos				
CAS 62-73-7 DRE-C12530000 DRE-L12530000CY DRE-XA12530000AC	MW 220.9757 Dichlorvos (*) Dichlorvos 10 µg/mL in Cyclohexane (*) Dichlorvos 100 µg/mL in Acetone (*)	C ₄ H ₇ Cl ₂ O ₄ P	250mg 10ml 1ml	
Dichlorvos D6 (dimethyl D6)				
CAS 203645-53-8 DRE-C12530100 DRE-XA12530100CY	MW 227.0127 Dichlorvos D6 (dimethyl D6) Dichlorvos D6 100 µg/mL in Cyclohexane (*)	C ₄ ² H ₆ HCl ₂ O ₄ P	10mg 1.1ml	
Diclocymet				
CAS 139920-32-4 DRE-C12536000	MW 313.2222 Diclocymet(‡)	C ₁₅ H ₁₀ Cl ₂ N ₂ O	50mg	
Diclofop (free acid)				
CAS 40843-25-2 DRE-C12539000	MW 327.1594 Diclofop (free acid)	C ₁₅ H ₁₂ Cl ₂ O ₄	100mg	
Dicofol				
CAS 115-32-2 DRE-C12570000 DRE-L12570000CY DRE-L12570000IQ DRE-XA12570000CY	MW 370.4857 Dicofol(‡) Dicofol 10 µg/mL in Cyclohexane Dicofol 10 µg/mL in Isooctane Dicofol 100 µg/mL in Cyclohexane	C ₁₄ H ₉ Cl ₅ O	100mg 10ml 10ml 1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Dicofol D8 (ring D8)				
CAS n/a	MW 378.535	C ₁₄ ² H ₈ HCl ₅ O		
DRE-XA12570100CY	Dicofol D8 (ring D8) 100 µg/mL in Cyclohexane(‡)			1ml
				
2,4'-Dicofol				
CAS 10606-46-9	MW 370.4857	C ₁₄ H ₈ Cl ₅ O		
DRE-LA12571000CY	2,4'-Dicofol 10 µg/mL in Cyclohexane			1ml
DRE-LA09010253CY	2,4'-Dicofol 10 µg/mL in Cyclohexane(‡)			10ml
				
Dicrotophos				
CAS 141-66-2	MW 237.1901	C ₈ H ₁₆ NO ₅ P		
DRE-C12580000	Dicrotophos			100mg
DRE-L12580000CY	Dicrotophos 10 µg/mL in Cyclohexane			10ml
				
Dicyclanil				
CAS 112636-83-6	MW 190.2052	C ₈ H ₁₀ N ₆		
DRE-C12583000	Dicyclanil			100mg
DRE-L12583000AL	Dicyclanil 10 µg/mL in Acetonitrile			10ml
				
Dieldrin				
CAS 60-57-1	MW 380.9093	C ₁₂ H ₈ Cl ₆ O		
DRE-C12590000	Dieldrin			50mg
DRE-L12590000CY	Dieldrin 10 µg/mL in Cyclohexane			10ml
DRE-XA12590000CY	Dieldrin 100 µg/mL in Cyclohexane			1ml
				
Dienochlor				
CAS 2227-17-0	MW 474.637	C ₁₀ Cl ₁₀		
DRE-C12600000	Dienochlor			250mg
				
Diflubenzuron				
CAS 35367-38-5	MW 310.6833	C ₁₄ H ₉ ClF ₂ N ₂ O ₂		
DRE-C12630000	Diflubenzuron(‡)			250mg
				
Dimefluthrin				
CAS 271241-14-6	MW 374.3698	C ₁₉ H ₂₂ F ₄ O ₃		
DRE-C12655000	Dimefluthrin			100mg
				
Dimefox				
CAS 115-26-4	MW 154.123	C ₄ H ₁₂ FN ₂ OP		
DRE-C12650000	Dimefox			100mg
DRE-L12650000IO	Dimefox 10 µg/mL in Isooctane			10ml
				










Pesticide and metabolites: Insecticides

Product code	Description			
Dimethoate				
CAS 60-51-5	MW 229.2574	C ₅ H ₁₂ NO ₃ PS ₂		
DRE-C12700000	Dimethoate(±)		100mg	
DRE-L12700000AL	Dimethoate 10 µg/mL in Acetonitrile (*)		10ml	
DRE-L12700000CY	Dimethoate 10 µg/mL in Cyclohexane		10ml	
DRE-XA12700000AC	Dimethoate 100 µg/mL in Acetone		1ml	
DRE-XA12700000AL	Dimethoate 100 µg/mL in Acetonitrile		1ml	
DRE-LA09010261CY	Dimethoate 10 µg/mL in Cyclohexane(±)		10ml	
DRE-XA09010262AC	Dimethoate 100 µg/mL in Acetone(±)		1ml	
Dimethoate D6 (O,O dimethyl D6)				
CAS 1219794-81-6	MW 235.2944	C ₅ ² H ₆ H ₆ NO ₃ PS ₂		
DRE-C12700100	Dimethoate D6 (O,O-dimethyl-D6)		10mg	
DRE-XA12700100AC	Dimethoate D6 (O,O dimethyl D6) 100 µg/mL in Acetone		1.1ml	
N-(2,4-Dimethylphenyl)formamide (Form-2',4'-xylidide)				
CAS 60397-77-5	MW 149.1897	C ₉ H ₁₁ NO		
DRE-C12737000	N-(2,4-Dimethylphenyl)formamide		100mg	
(3,5-Dimethylphenyl) N-Methylcarbamate (XMC)				
CAS 2655-14-3	MW 179.2157	C ₁₀ H ₁₃ NO ₂		
DRE-C17943000	XMC (3,5-xylol methylcarbamate)(±)		100mg	
DRE-L17943000CY	XMC (3,5-xylol methylcarbamate) 10 µg/mL in Cyclohexane		10ml	
N-2,4-Dimethylphenyl-N'-methylformamidine (N-Methyl-N'-(2,4-xylol)formamidine)				
CAS 33089-74-6	MW 162.2316	C ₁₀ H ₁₄ N ₂		
DRE-C12738000	N-2,4-Dimethylphenyl-N'-methylformamidine(±)		10mg	
Dimethylvinphos (Z type)				
CAS 67628-93-7	MW 331.5168	C ₁₀ H ₁₀ Cl ₃ O ₄ P		
DRE-C12765060	Dimethylvinphos (Z type)		50mg	
DRE-L12765060CY	Dimethylvinphos (Z type)10 µg/mL in Cyclohexane		10ml	
Dimetilan				
CAS 644-64-4	MW 240.259	C ₁₀ H ₁₆ N ₂ O ₃		
DRE-C12770000	Dimetilan		100mg	
DRE-L12770000CY	Dimetilan 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010184AL	Dimetilan 100 µg/mL in Acetonitrile(±)		1ml	
Dinobuton				
CAS 973-21-7	MW 326.3019	C ₁₄ H ₁₈ N ₂ O ₇		
DRE-C12790000	Dinobuton		100mg	









Pesticide and metabolites: Insecticides

Product code	Description			
Dinotefuran				
CAS 165252-70-0	MW 202.2111	C ₇ H ₁₄ N ₄ O ₃		
DRE-C12820000	Dinotefuran(±)		50mg	
DRE-L12820000AL	Dinotefuran 10 µg/mL in Acetonitrile		10ml	
Diofenolan				
CAS 63837-33-2	MW 300.349	C ₁₈ H ₂₀ O ₄		
DRE-C12840000	Diofenolan		100mg	
Dioxabenzofos (Salithion)				
CAS 3811-49-2	MW 216.194	C ₈ H ₉ O ₃ PS		
DRE-C12850000	Dioxabenzofos (Salithion)		10mg	
DRE-L12850000CY	Dioxabenzofos (Salithion) 10 µg/mL in Cyclohexane		10ml	
Dioxacarb				
CAS 6988-21-2	MW 223.2252	C ₁₁ H ₁₃ NO ₄		
DRE-C12860000	Dioxacarb(±)		100mg	
Dioxathion				
CAS 78-34-2	MW 456.5388	C ₁₂ H ₂₆ O ₆ P ₂ S ₄		
DRE-C12870000	Dioxathion(±)		100mg	
Disulfoton				
CAS 298-04-4	MW 274.404	C ₈ H ₁₉ O ₂ PS ₃		
DRE-CA12980000	Disulfoton		250mg	
DRE-L12980000CY	Disulfoton 10 µg/mL in Cyclohexane		10ml	
DRE-XA12980000AL	Disulfoton 100 µg/mL in Acetonitrile		1ml	
DRE-XA12980000CY	Disulfoton 100 µg/mL in Cyclohexane		1ml	
DRE-GA09010396DI	Disulfoton Solution, 1000 µg/mL in Dichloromethane(±)		1ml	
DRE-GA09010396DI-SS	Disulfoton Solution, Second Source, 2000 µg/mL in Dichloromethane(±)		1ml	
Disulfoton D10 (Di-ethyl D10)				
CAS n/a	MW 284.4656	C ₈ ² H ₁₀ H ₉ O ₂ PS ₃		
DRE-XA12980100AC	Disulfoton D10 (Di-ethyl D10) 100 µg/mL in Acetone		1.1ml	
DRE-XA12980100CY	Disulfoton D10 (Di-ethyl D10) 100 µg/mL in Cyclohexane		1.1ml	
Disulfoton-oxon-sulfon				
CAS 2496-91-5	MW 290.3372	C ₈ H ₁₉ O ₅ PS ₂		
DRE-CA12982200	Disulfoton-oxon-sulfon		10mg	
DRE-XA12982200AL	Disulfoton-oxon-sulfon 100 µg/mL in Acetonitrile		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Disulfoton-sulfone				
CAS 2497-06-5	MW 306.4028	C ₈ H ₁₉ O ₄ PS ₃		
DRE-CA12985000	Disulfoton-sulfone		50mg	
DRE-L12985000CY	Disulfoton-sulfone 10 µg/mL in Cyclohexane (*)		10ml	
DRE-XA09010187ME	Disulfoton Sulfone 100 µg/mL in Methanol(‡)		1ml	
Disulfoton-sulfoxide				
CAS 2497-07-6	MW 290.4034	C ₈ H ₁₉ O ₃ PS ₃		
DRE-CA12985500	Disulfoton-sulfoxide		50mg	
DRE-L12985500CY	Disulfoton-sulfoxide 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010188ME	Disulfoton sulfoxide 100 µg/mL in Methanol(‡)		1ml	
EPN				
CAS 2104-64-5	MW 323.304	C ₁₄ H ₁₄ NO ₄ PS		
DRE-C13180000	EPN		100mg	
DRE-L13180000CY	EPN 10 µg/mL in Cyclohexane		10ml	
DRE-XA13180000AL	EPN 100 µg/mL in Acetonitrile		1ml	
DRE-XA13180000CY	EPN 100 µg/mL in Cyclohexane		1ml	
Emamectin Benzoate				
CAS 155569-91-8	MW 2002.4535	C ₄₉ H ₇₅ NO ₁₃ ·C ₁₄ H ₇ NO ₁₃ ·2C ₇ H ₆ O ₂		
DRE-CA13117000	Emamectin benzoate(‡)		100mg	
Empenthrin				
CAS 54406-48-3	MW 274.3978	C ₁₈ H ₂₆ O ₂		
DRE-C13118000	Empenthrin		100mg	
Endosulfan (α + β)				
CAS 115-29-7	MW 406.9251	C ₉ H ₆ Cl ₆ O ₃ S		
DRE-C13120000	Endosulfan (α + β)(‡)		250mg	
DRE-L13120000AL	Endosulfan (α + β) 10 µg/mL in Acetonitrile		10ml	
DRE-L13120000IO	Endosulfan (α + β) 10 µg/mL in Isooctane		10ml	
DRE-XA13120000IO	Endosulfan (α + β) 100 µg/mL in Isooctane		1ml	
Endosulfan-alcohol				
CAS 2157-19-9	MW 360.8766	C ₉ H ₆ Cl ₆ O ₂		
DRE-C13130000	Endosulfan-alcohol		100mg	
Endosulfan-ether				
CAS 3369-52-6	MW 342.8613	C ₉ H ₆ Cl ₆ O		
DRE-C13131000	Endosulfan-ether		100mg	
Endosulfan-lacton				
CAS 3868-61-9	MW 356.8449	C ₉ H ₄ Cl ₆ O ₂		
DRE-C13132000	Endosulfan-lacton		100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Endosulfan-sulfate				
CAS 1031-07-8	MW 422.9245	C ₉ H ₆ Cl ₆ O ₄ S		AND Enantiomer 
DRE-C13133000	Endosulfan-sulfate(±)		100mg	
DRE-L13133000CY	Endosulfan-sulfate 10 µg/mL in Cyclohexane		10ml	
DRE-XA13133000CY	Endosulfan-sulfate 100 µg/mL in Cyclohexane		1ml	
α-Endosulfan				
CAS 959-98-8	MW 406.9251	C ₉ H ₆ Cl ₆ O ₃ S		AND Enantiomer 
DRE-C13121000	α-Endosulfan(±)		100mg	
DRE-L13121000AL	α-Endosulfan 10 µg/mL in Acetonitrile		10ml	
DRE-L13121000CY	α-Endosulfan 10 µg/mL in Cyclohexane		10ml	
DRE-XA13121000IO	α-Endosulfan 100 µg/mL in Isooctane		1ml	
α-Endosulfan D4				
CAS 203645-57-2	MW 410.9498	C ₉ ² H ₄ H ₂ Cl ₆ O ₃ S		AND Enantiomer 
DRE-C13121100	α-Endosulfan D4		10mg	
DRE-XA13121100AC	α-Endosulfan D4 100 µg/mL in Acetone		1.1ml	
β-Endosulfan				
CAS 33213-65-9	MW 406.9251	C ₉ H ₆ Cl ₆ O ₃ S		AND Enantiomer 
DRE-C13122000	β-Endosulfan(±)		100mg	
DRE-L13122000AL	β-Endosulfan 10 µg/mL in Acetonitrile		10ml	
DRE-L13122000CY	β-Endosulfan 10 µg/mL in Cyclohexane		10ml	
DRE-XA13122000IO	β-Endosulfan 100 µg/mL in Isooctane		1ml	
β-Endosulfan D4				
CAS 203716-99-8	MW 410.9498	C ₉ ² H ₄ H ₂ Cl ₆ O ₃ S		AND Enantiomer 
DRE-C13122100	β-Endosulfan D4		10mg	
DRE-XA13122100AC	β-Endosulfan D4 100 µg/mL in Acetone		1.1ml	
Endrin				
CAS 72-20-8	MW 380.9093	C ₁₂ H ₈ Cl ₆ O		AND Enantiomer 
DRE-C13160000	Endrin(±)		100mg	
DRE-L13160000AL	Endrin 10 µg/mL in Acetonitrile		10ml	
DRE-L13160000CY	Endrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA13160000IO	Endrin 100 µg/mL in Isooctane		1ml	
Endrin-aldehyde				
CAS 7421-93-4	MW 380.9093	C ₁₂ H ₈ Cl ₆ O		AND Enantiomer 
DRE-C13161000	Endrin-aldehyde		10mg	
DRE-L13161000CY	Endrin-aldehyde 10 µg/mL in Cyclohexane		10ml	
DRE-XA13161000CY	Endrin-aldehyde 100 µg/mL in Cyclohexane		1ml	
Endrin-ketone				
CAS 53494-70-5	MW 380.9093	C ₁₂ H ₈ Cl ₆ O		AND Enantiomer 
DRE-C13165000	Endrin-ketone		10mg	
DRE-L13165000CY	Endrin-ketone 10 µg/mL in Cyclohexane		10ml	
DRE-XA13165000CY	Endrin-ketone 100 µg/mL in Cyclohexane		1ml	

Pesticide and metabolites: Insecticides

Product code	Description			
Esfenvalerate				
CAS 66230-04-4	MW 419.9001	C ₂₅ H ₂₂ ClNO ₃		
DRE-C13211000	Esfenvalerate		100mg	
DRE-L13211000CY	Esfenvalerate 10 µg/mL in Cyclohexane		10ml	
DRE-XA13211000AL	Esfenvalerate 100 µg/mL in Acetonitrile		1ml	
Esfenvalerate free acid metabolite				
CAS 55332-38-2	MW 212.6727	C ₁₁ H ₁₃ ClO ₂		
DRE-LA13211020CY	Esfenvalerate free acid metabolite 10 µg/mL in Cyclohexane		1ml	
Ethiofencarb				
CAS 29973-13-5	MW 225.3073	C ₁₁ H ₁₅ NO ₂ S		
DRE-C13250000	Ethiofencarb		100mg	
DRE-L13250000CY	Ethiofencarb 10 µg/mL in Cyclohexane		10ml	
DRE-XA13250000CY	Ethiofencarb 100 µg/mL in Cyclohexane		1ml	
Ethiofencarb Sulfone				
CAS 53380-23-7	MW 257.3061	C ₁₁ H ₁₅ NO ₄ S		
DRE-C13250500	Ethiofencarb-sulfone		10mg	
DRE-XA09010190AL	Ethiofencarb Sulfone 100 µg/mL in Acetonitrile(‡)		1ml	
Ethiofencarb-phenol-sulfone				
CAS 53380-27-1	MW 200.2548	C ₉ H ₁₂ O ₃ S		
DRE-C13250200	Ethiofencarb-phenol-sulfone		50mg	
Ethiofencarb-phenol-sulfoxide				
CAS 53380-26-0	MW 184.2554	C ₉ H ₁₂ O ₂ S		
DRE-C13250300	Ethiofencarb-phenol-sulfoxide		50mg	
Ethiofencarb-sulfoxide				
CAS 53380-22-6	MW 241.3067	C ₁₁ H ₁₅ NO ₃ S		
DRE-CA13251000	Ethiofencarb-sulfoxide(‡)		10mg	
Ethion				
CAS 563-12-2	MW 384.4761	C ₉ H ₂₂ O ₄ P ₂ S ₄		
DRE-CA13270000	Ethion(‡)		250mg	
DRE-L13270000AL	Ethion 10 µg/mL in Acetonitrile		10ml	
DRE-L13270000CY	Ethion 10 µg/mL in Cyclohexane		10ml	
DRE-XA13270000CY	Ethion 100 µg/mL in Cyclohexane		1ml	
Ethion D20 (tetraethyl D20)				
CAS n/a	MW 404.5993	C ₉ ² H ₂₀ H ₂ O ₄ P ₂ S ₄		
DRE-XA13270100AC	Ethion D20 100 µg/mL in Acetone		1.1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Ethiprole				
CAS 181587-01-9 DRE-C13275000	MW 397.203 Ethiprole(‡)	C ₁₃ H ₉ Cl ₂ F ₃ N ₄ OS	100mg	
Ethoprophos (Ethoprop)				
CAS 13194-48-4 DRE-C13300000 DRE-L13300000AL DRE-XA13300000CY DRE-GA09010368AC	MW 242.339 Ethoprophos(‡) Ethoprophos 10 µg/mL in Acetonitrile Ethoprophos 100 µg/mL in Cyclohexane Ethoprop 100 µg/mL in Acetone(‡)	C ₈ H ₁₉ O ₂ PS ₂	100mg 10ml 1ml 1ml	
Ethoprophos D5 (ethyl D5)				
CAS n/a DRE-XA13300010CY	MW 247.3698 Ethoprophos D5 (ethyl D5) 100 µg/mL in Cyclohexane	C ₈ ² H ₅ H ₁₄ O ₂ PS ₂	1.1ml	
Etofenprox				
CAS 80844-07-1 DRE-C13363000 DRE-L13363000AL DRE-L13363000CY DRE-XA13363000AL	MW 376.488 Etofenprox(‡) Etofenprox 10 µg/mL in Acetonitrile Etofenprox 10 µg/mL in Cyclohexane Etofenprox 100 µg/mL in Acetonitrile	C ₂₈ H ₂₈ O ₃	100mg 10ml 10ml 1ml	
Etofenprox D5 (ethyl D5)				
CAS 1705649-55-3 DRE-XA13363010AC	MW 381.5188 Etofenprox D5 (ethyl D5) 100 µg/mL in Acetone(‡)	C ₂₈ ² H ₅ H ₂₃ O ₃	1.1ml	
Etoxazole				
CAS 153233-91-1 DRE-C13368000 DRE-L13368000CY DRE-XA13368000AL	MW 359.4096 Etoxazole(‡) Etoxazole 10 µg/mL in Cyclohexane Etoxazole 100 µg/mL in Acetonitrile	C ₂₁ H ₂₃ F ₂ NO ₂	50mg 10ml 1ml	
Etrinfos				
CAS 38260-54-7 DRE-CA13380000 DRE-L13380000CY DRE-XA09010191ME	MW 292.2917 Etrinfos (technical mixture) (*) Etrinfos 10 µg/mL in Cyclohexane Etrinfos 100 µg/mL in Methanol(‡)	C ₁₀ H ₁₇ N ₂ O ₄ PS	100mg 10ml 1ml	
Famphur				
CAS 52-85-7 DRE-C13400000 DRE-L13400000CY DRE-XA09010192ME	MW 325.3415 Famphur Famphur 10 µg/mL in Cyclohexane Famphur 100 µg/mL in Methanol(‡)	C ₁₀ H ₁₆ NO ₅ PS ₂	100mg 10ml 1ml	
Fenamiphos				
CAS 22224-92-6 DRE-C13420000 DRE-L13420000CY DRE-XA13420000AL	MW 303.3574 Fenamiphos Fenamiphos 10 µg/mL in Cyclohexane Fenamiphos 100 µg/mL in Acetonitrile(‡)	C ₁₃ H ₂₂ NO ₃ PS	250mg 10ml 1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Fenamiphos-sulfone				
CAS 31972-44-8	MW 335.3562	C ₁₃ H ₂₂ NO ₅ PS		
DRE-C13421000	Fenamiphos-sulfone		100mg	
DRE-L13421000EA	Fenamiphos sulfone 10 µg/mL in Ethylacetate		10ml	
DRE-XA09010193ME	Fenamiphos Sulfone 100 µg/mL in Methanol(±)		1ml	
Fenamiphos-sulfoxide				
CAS 31972-43-7	MW 319.3568	C ₁₃ H ₂₂ NO ₄ PS		
DRE-C13422000	Fenamiphos-sulfoxide		100mg	
DRE-L13422000MB	Fenamiphos sulfoxide 10 µg/mL in Methyl-tert. butyl ether		10ml	
DRE-XA09010194ME	Fenamiphos Sulfoxide 100 µg/mL in Methanol(±)		1ml	
Fenazaflor				
CAS 14255-88-0	MW 375.1295	C ₁₅ H ₇ Cl ₂ F ₃ N ₂ O ₂		
DRE-C13440000	Fenazaflor		100mg	
Fenazaquin				
CAS 120928-09-8	MW 306.4015	C ₂₀ H ₂₂ N ₂ O		
DRE-C13441000	Fenazaquin		100mg	
DRE-L13441000AL	Fenazaquin 10 µg/mL in Acetonitrile		10ml	
DRE-L13441000CY	Fenazaquin 10 µg/mL in Cyclohexane		10ml	
DRE-XA13441000CY	Fenazaquin 100 µg/mL in Cyclohexane		1ml	
Fenazox				
CAS 495-48-7	MW 198.2206	C ₁₂ H ₁₀ N ₂ O		
DRE-C13442000	Fenazox		250mg	
Fenbutatin oxide				
CAS 13356-08-6	MW 1052.6807	C ₆₀ H ₇₈ OSn ₂		
DRE-C13450000	Fenbutatin oxide		250mg	
DRE-L13450000CY	Fenbutatin oxide 10 µg/mL in Cyclohexane		10ml	
Fenchlorphos				
CAS 299-84-3	MW 321.5451	C ₈ H ₈ Cl ₃ O ₃ PS		
DRE-C13460000	Fenchlorphos		100mg	
DRE-L13460000CY	Fenchlorphos 10 µg/mL in Cyclohexane		10ml	
DRE-XA13460000CY	Fenchlorphos 100 µg/mL in Cyclohexane		1ml	
Fenchlorphos-oxon				
CAS 3983-45-7	MW 305.4795	C ₈ H ₈ Cl ₃ O ₄ P		
DRE-C13460500	Fenchlorphos-oxon		100mg	
DRE-XA13460500AL	Fenchlorphos-oxon 100 µg/mL in Acetonitrile		1ml	
Fenfluthrin				
CAS 75867-00-4	MW 389.1447	C ₁₅ H ₁₁ Cl ₂ F ₅ O ₂		
DRE-XA13469000AL	Fenfluthrin 100 µg/mL in Acetonitrile(±)		1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
cis-Fenfluthrin				
CAS 74337-23-8 DRE-XA13469020AL	MW 389.1447 cis-Fenfluthrin 100 µg/mL in Acetonitrile	C ₁₅ H ₁₁ Cl ₂ F ₅ O ₂	1ml	
Fenitrothion				
CAS 122-14-5 DRE-C13480000 DRE-L13480000AL DRE-L13480000CY DRE-XA13480000AL DRE-XA13480000CY	MW 277.234 Fenitrothion(±) Fenitrothion 10 µg/mL in Acetonitrile Fenitrothion 10 µg/mL in Cyclohexane Fenitrothion 100 µg/mL in Acetonitrile Fenitrothion 100 µg/mL in Cyclohexane	C ₉ H ₁₂ NO ₅ PS	250mg 10ml 10ml 1ml 1ml	
Fenitrothion D6 (O,O-dimethyl D6)				
CAS 203645-59-4 DRE-C13480100 DRE-XA13480100CY	MW 283.271 Fenitrothion D6 (O,O-dimethyl D6) Fenitrothion D6 (O,O-dimethyl D6) 100 µg/mL in Cyclohexane	C ₉ ² H ₆ H ₆ NO ₅ PS	10mg 1.1ml	
Fenitrothion-oxon				
CAS 2255-17-6 DRE-LA13481000EA	MW 261.1684 Fenitrothion-oxon 10 µg/mL in Ethylacetate	C ₉ H ₁₂ NO ₆ P	1ml	
Fenitrothion-oxon D6 (O,O-dimethyl D6)				
CAS 1185155-54-7 DRE-XA13481100AC	MW 267.2054 Fenitrothion-oxon D6 (O,O-dimethyl D6) 100 µg/mL in Acetone	C ₉ ² H ₆ H ₆ NO ₆ P	1.1ml	
Fenobucarb				
CAS 3766-81-2 DRE-C13485000 DRE-L13485000CY	MW 207.2689 Fenobucarb Fenobucarb 10 µg/mL in Cyclohexane	C ₁₂ H ₁₇ NO ₂	250mg 10ml	
Fenothiocarb				
CAS 62850-32-2 DRE-C13497500 DRE-L13497500AC DRE-L13497500TO DRE-XA13497500TO DRE-XA09010195AL	MW 253.3605 Fenothiocarb(±) Fenothiocarb 10 µg/mL in Acetone Fenothiocarb 10 µg/mL in Toluene Fenothiocarb 100 µg/mL in Toluene Fenothiocarb 100 µg/mL in Acetonitrile(±)	C ₁₃ H ₁₉ NO ₂ S	100mg 10ml 10ml 1ml 1ml	
Fenoxycarb				
CAS 72490-01-8 DRE-C13520000 DRE-L13520000AL	MW 301.3371 Fenoxycarb Fenoxycarb 10 µg/mL in Acetonitrile	C ₁₇ H ₁₉ NO ₄	250mg 10ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Fenpropathrin				
CAS 39515-41-8	MW 349.4229	C ₂₂ H ₂₃ NO ₃		
DRE-C13530000	Fenpropathrin(‡)		250mg	
DRE-L13530000AL	Fenpropathrin 10 µg/mL in Acetonitrile		10ml	
DRE-L13530000IO	Fenpropathrin 10 µg/mL in Isooctane		10ml	
DRE-XA13530000IO	Fenpropathrin 100 µg/mL in Isooctane		1ml	
(E)-Fenpyroximate				
CAS 134098-61-6	MW 421.4889	C ₂₄ H ₂₇ N ₃ O ₄		
DRE-C13545000	(E)-Fenpyroximate(‡)		100mg	
DRE-L13545000AL	(E)-Fenpyroximate 10 µg/mL in Acetonitrile		10ml	
DRE-L13545000CY	(E)-Fenpyroximate 10 µg/mL in Cyclohexane		10ml	
DRE-XA13545000CY	(E)-Fenpyroximate 100 µg/mL in Cyclohexane		1ml	
(Z)-Fenpyroximate				
CAS 149054-53-5	MW 421.4889	C ₂₄ H ₂₇ N ₃ O ₄		
DRE-XA13545200AL	(Z)-Fenpyroximate 100 µg/mL in Aceonitrile		1ml	
Fenson				
CAS 80-38-6	MW 268.7161	C ₁₂ H ₉ ClO ₃ S		
DRE-C13560000	Fenson		250mg	
DRE-L13560000CY	Fenson 10 µg/mL in Cyclohexane		10ml	
Fensulfothion				
CAS 115-90-2	MW 308.354	C ₁₁ H ₁₇ O ₄ PS ₂		
DRE-C13570000	Fensulfothion		100mg	
DRE-XA13570000AL	Fensulfothion 100 µg/mL in Acetonitrile		1ml	
DRE-LA09010260CY	Fensulfothion 10 µg/mL in Cyclohexane(‡)		10ml	
Fensulfothion-oxon				
CAS 6552-21-2	MW 292.2884	C ₁₁ H ₁₇ O ₅ PS		
DRE-C13570013	Fensulfothion-oxon		50mg	
DRE-XA13570013CY	Fensulfothion-oxon 100 µg/mL in Cyclohexane(‡)		1ml	
Fensulfothion-oxon-sulfone				
CAS 6132-17-8	MW 308.2878	C ₁₁ H ₁₇ O ₆ PS		
DRE-C13570016	Fensulfothion-oxon-sulfone		50mg	
DRE-XA13570016CY	Fensulfothion-oxon-sulfone 100 µg/mL in Cyclohexane		1ml	
Fensulfothion-sulfone				
CAS 14255-72-2	MW 324.3534	C ₁₁ H ₁₇ O ₅ PS ₂		
DRE-C13570020	Fensulfothion-sulfone(‡)		50mg	
DRE-XA13570020CY	Fensulfothion-sulfone 100 µg/mL in Cyclohexane		1ml	
Fenthion				
CAS 55-38-9	MW 278.3281	C ₁₀ H ₁₅ O ₃ PS ₂		
DRE-C13580000	Fenthion(‡)		250mg	
DRE-L13580000AL	Fenthion 10 µg/mL in Acetonitrile		10ml	
DRE-L13580000CY	Fenthion 10 µg/mL in Cyclohexane		10ml	

Pesticide and metabolites: Insecticides

Product code	Description		
Fenthion D6 (dimethoxy D6)			
CAS 1189662-83-6	MW 284.365	C ₁₀ H ₁₆ H ₉ O ₃ PS ₂	
DRE-C13580100	Fenthion D6 (O,O-dimethyl D6)		
DRE-XA13580100AC	Fenthion D6 (dimethoxy D6) 100 µg/mL in Acetone		
			10mg 1.1ml
Fenthion-oxon			
CAS 6552-12-1	MW 262.2625	C ₁₀ H ₁₅ O ₄ PS	
DRE-C13585000	Fenthion-oxon		
DRE-L13585000AL	Fenthion-oxon 10 µg/mL in Acetonitrile(‡)		
DRE-LA13585000AL	Fenthion-oxon 10 µg/mL in Acetonitrile		
DRE-XA13585000CY	Fenthion-oxon 100 µg/mL in Cyclohexane		
			10mg 10ml 1ml 1ml
Fenthion-oxon-sulfone			
CAS 14086-35-2	MW 294.2613	C ₁₀ H ₁₅ O ₆ PS	
DRE-C13585200	Fenthion-oxon-sulfone		
DRE-L13585200AL	Fenthion-oxon-sulfone 10 µg/mL in Acetonitrile		
DRE-LA13585200AL	Fenthion-oxon-sulfone 10 µg/mL in Acetonitrile(‡)		
DRE-XA13585200AL	Fenthion-oxon-sulfone 100 µg/mL in Acetonitrile		
			50mg 10ml 1ml 1ml
Fenthion-oxon-sulfoxide			
CAS 6552-13-2	MW 278.2619	C ₁₀ H ₁₅ O ₅ PS	
DRE-C13585400	Fenthion-oxon-sulfoxide		
DRE-L13585400AL	Fenthion-oxon-sulfoxide 10 µg/mL in Acetonitrile		
DRE-LA13585400AL	Fenthion-oxon-sulfoxide 10 µg/mL in Acetonitrile(‡)		
DRE-LA13585400CY	Fenthion-oxon-sulfoxide 10 µg/mL in Cyclohexane		
DRE-XA13585400AL	Fenthion-oxon-sulfoxide 100 µg/mL in Acetonitrile		50mg 10ml 1ml 1ml
Fenthion-sulfone			
CAS 3761-42-0	MW 310.3269	C ₁₀ H ₁₅ O ₅ PS ₂	
DRE-C13586000	Fenthion-sulfone(‡)		
DRE-LA13586000EA	Fenthion-sulfone 10 µg/mL in Ethylacetate		
DRE-XA13586000EA	Fenthion-sulfone 100 µg/mL in Ethylacetate(‡)		
			10mg 1ml 1ml
Fenthion-sulfoxide			
CAS 3761-41-9	MW 294.3275	C ₁₀ H ₁₅ O ₄ PS ₂	
DRE-C13586500	Fenthion-sulfoxide(‡)		
			50mg
Fenvalerate			
CAS 51630-58-1	MW 419.9001	C ₂₅ H ₂₂ ClNO ₃	
DRE-C13630000	Fenvalerate(‡)		
DRE-L13630000AL	Fenvalerate 10 µg/mL in Acetonitrile		
DRE-L13630000IO	Fenvalerate 10 µg/mL in Isooctane		
DRE-XA13630000IO	Fenvalerate 100 µg/mL in Isooctane		
			250mg 10ml 10ml 1ml
Fenvalerate D7 (isopropyl D7)			
CAS n/a	MW 426.9432	C ₂₅ H ₂₇ H ₁₅ ClNO ₃	
DRE-XA13630010IO	Fenvalerate D7(isopropyl D7) 100 µg/mL in Isooctane		
			1.1ml










Pesticide and metabolites: Insecticides

Product code	Description			
Fenvalerate free acid metabolite				
CAS 2012-74-0 DRE-C13630020	MW 212.6727 Fenvalerate free acid metabolite	C ₁₁ H ₁₃ ClO ₂	100mg	
Fipronil				
CAS 120068-37-3 DRE-C13645000 DRE-L13645000AL DRE-L13645000CY DRE-XA13645000AL	Fipronil(‡) Fipronil 10 µg/mL in Acetonitrile Fipronil 10 µg/mL in Cyclohexane Fipronil 100 µg/mL in Acetonitrile	C ₁₂ H ₄ Cl ₂ F ₆ N ₄ OS	100mg 10ml 10ml 1ml	
Fipronil Amide				
CAS 205650-69-7 DRE-XA13645200AL	Fipronil-carboxamide 100 µg/mL in Acetonitrile	C ₁₂ H ₄ Cl ₂ F ₆ N ₄ O ₂ S	1ml	
Fipronil Sulfide				
CAS 120067-83-6 DRE-C13645400 DRE-L13645400AL DRE-XA13645400AL	Fipronil-sulfide(‡) Fipronil-sulfide 10 µg/mL in Acetonitrile Fipronil-sulfide 100 µg/mL in Acetonitrile	C ₁₂ H ₄ Cl ₂ F ₆ N ₄ S	25mg 10ml 1ml	
Fipronil Sulfone				
CAS 120068-36-2 DRE-C13645500 DRE-L13645500AL DRE-XA13645500AL	Fipronil-sulfone(‡) Fipronil-sulfone 10 µg/mL in Acetonitrile Fipronil-sulfone 100 µg/mL in Acetonitrile	C ₁₂ H ₄ Cl ₂ F ₆ N ₄ O ₂ S	50mg 10ml 1ml	
Fipronil des F3				
CAS 154807-27-9 DRE-XA13645100AL	Fipronil des F3 100 µg/mL in Acetonitrile	C ₁₂ H ₇ Cl ₂ F ₃ N ₄ OS	1.1ml	
Fipronil-desulfinyl				
CAS 205650-65-3 DRE-L13645300AL DRE-XA13645300AL	Fipronil-desulfinyl 10 µg/mL in Acetonitrile Fipronil-desulfinyl 100 µg/mL in Acetonitrile	C ₁₂ H ₄ Cl ₂ F ₆ N ₄	10ml 1ml	
Flonicamid				
CAS 158062-67-0 DRE-C13662100 DRE-XA13662100AL	Flonicamid(‡) Flonicamid 100 µg/mL in Acetonitrile(‡)	C ₉ H ₆ F ₃ N ₃ O	25mg 1ml	
Flonicamid free acid				
CAS 158063-66-2 DRE-C13662120	Flonicamid free acid	C ₇ H ₄ F ₃ NO ₂	100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Flonicamid-carboxylic acid (4-(Trifluoromethyl)nicotinoyl glycine)				
CAS 207502-65-6 DRE-C13662110	MW 248.1587 Flonicamid-carboxylic acid	C ₉ H ₇ F ₃ N ₂ O ₃	25mg	
Fluacrypyrim				
CAS 229977-93-9 DRE-C13667000	MW 426.3863 Fluacrypyrim	C ₂₀ H ₁₇ F ₃ N ₂ O ₅	50mg	
Fluazifop (free acid)				
CAS 69335-91-7 DRE-C13669000 DRE-XA13669000AL	MW 327.2553 Fluazifop (free acid)(‡) Fluazifop (free acid) 100 µg/mL in Acetonitrile	C ₁₅ H ₁₂ F ₃ NO ₄	10mg 1ml	
Fluazifop-P (free acid)				
CAS 83066-88-0 DRE-C13669500	MW 327.2553 Fluazifop-P (free acid)	C ₁₅ H ₁₂ F ₃ NO ₄	10mg	
Flubendiamide				
CAS 272451-65-7 DRE-C13679000	MW 682.3901 Flubendiamide(‡)	C ₂₃ H ₂₂ F ₇ IN ₂ O ₄ S	100mg	
Flubenzimine				
CAS 37893-02-0 DRE-C13680000	MW 416.3435 Flubenzimine	C ₁₇ H ₁₀ F ₆ N ₄ S	250mg	
Flucifuron				
CAS 370-50-3 DRE-C13697000	MW 417.1332 Flucifuron	C ₁₅ H ₈ Cl ₂ F ₆ N ₂ O	10mg	
Flucycloxuron				
CAS 113036-88-7 DRE-C13698500	MW 483.8944 Flucycloxuron	C ₂₅ H ₂₀ ClF ₂ N ₃ O ₃	10mg	
Flucythrinate				
CAS 70124-77-5 DRE-C13700000 DRE-L13700000CY DRE-XA13700000IO	MW 451.4619 Flucythrinate(‡) Flucythrinate 10 µg/mL in Cyclohexane Flucythrinate 100 µg/mL in Isooctane	C ₂₆ H ₂₃ F ₂ NO ₄	100mg 10ml 1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Fluenetil (Lambrol)				
CAS 4301-50-2	MW 258.2875	C ₁₆ H ₁₅ FO ₂		
DRE-LA13710000CY	Fluenetil (Lambrol) 10 µg/mL in Cyclohexane		1ml	
Flufenoxuron				
CAS 101463-69-8	MW 488.7671	C ₂₁ H ₁₁ ClF ₆ N ₂ O ₃		
DRE-C13712000	Flufenoxuron(±)		100mg	
DRE-L13712000CY	Flufenoxuron 10 µg/mL in Cyclohexane		10ml	
Flufenzine				
CAS 162320-67-4	MW 304.682	C ₁₄ H ₇ ClF ₂ N ₄		
DRE-C13715000	Flufenzine		25mg	
DRE-XA13715000AL	Flufenzine 100 µg/mL in Acetonitrile		1ml	
Flufiprole				
CAS 704886-18-0	MW 491.2382	C ₁₆ H ₁₀ Cl ₂ F ₆ N ₄ OS		
DRE-C13717000	Flufiprole		100mg	
Flumethrin				
CAS 69770-45-2	MW 510.3836	C ₂₈ H ₂₂ Cl ₂ FNO ₃		
DRE-C13719000	Flumethrin		100mg	
DRE-L13719000CY	Flumethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA13719000AL	Flumethrin 100 µg/mL in Acetonitrile		1ml	
tau-Fluvalinate				
CAS 102851-06-9	MW 502.9127	C ₂₆ H ₂₂ ClF ₃ N ₂ O ₃		
DRE-C13870000	tau-Fluvalinate		100mg	
DRE-L13870000CY	tau-Fluvalinate 10 µg/mL in Cyclohexane		10ml	
DRE-XA13870000CY	tau-Fluvalinate 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010200AL	Tau-fluvaliante 100 µg/mL in Acetonitrile(±)		1ml	
Fonofos				
CAS 944-22-9	MW 246.3293	C ₁₀ H ₁₅ OPS ₂		
DRE-C13900000	Fonofos(±)		100mg	
DRE-L13900000CY	Fonofos 10 µg/mL in Cyclohexane		10ml	
DRE-XA13900000AL	Fonofos 100 µg/mL in Acetonitrile		1ml	
Formetanate hydrochloride				
CAS 23422-53-9	MW 257.7166	C ₁₁ H ₁₅ N ₃ O ₂ ·ClH		
DRE-C13910000	Formetanate hydrochloride		100mg	
DRE-L13910000AL	Formetanate hydrochloride 10 µg/mL in Acetonitrile		10ml	
Formothion				
CAS 2540-82-1	MW 257.2675	C ₆ H ₁₂ NO ₄ PS ₂		
DRE-CA13920000	Formothion		100mg	
DRE-L13920000EA	Formothion 10 µg/mL in Ethylacetate		10ml	
DRE-XA13920000AL	Formothion 100 µg/mL in Acetonitrile		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Fosthiazate				
CAS 98886-44-3	MW 283.3479	C ₉ H ₁₈ NO ₃ PS ₂		
DRE-CA13944500	Fosthiazate(±)		25mg	
DRE-L13944500AL	Fosthiazate 10 µg/mL in Acetonitrile		10ml	
Furathiocarb				
CAS 65907-30-4	MW 382.4744	C ₁₈ H ₂₆ N ₂ O ₅ S		
DRE-XA13970000AL	Furathiocarb 100 µg/mL in Acetonitrile		1ml	
α-HCH 13C6				
CAS 222966-66-7	MW 296.7858	¹³ C ₆ H ₆ Cl ₆		
DRE-XA14071300CY	α-HCH 13C6 100 µg/mL in Cyclohexane(±)		1.1ml	
α-HCH D6				
CAS 86194-41-4	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-C14071400	α-HCH D6(±)		10mg	
DRE-XA14071400CY	α-HCH D6 100 µg/mL in Cyclohexane(±)		1ml	
β-HCH				
CAS 319-85-7	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-C14072000	β-HCH(±)		100mg	
DRE-L14072000CY	β-HCH 10 µg/mL in Cyclohexane		10ml	
DRE-XA14072000TO	β-HCH 100 µg/mL in Toluene		1ml	
β-HCH D6				
CAS 86194-42-5	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-L14072100CY	β-HCH D6 10 µg/mL in Cyclohexane		10ml	
δ-HCH				
CAS 319-86-8	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-C14074000	δ-HCH(±)		50mg	
DRE-L14074000CY	δ-HCH 10 µg/mL in Cyclohexane(±)		10ml	
DRE-XA14074000CY	δ-HCH 100 µg/mL in Cyclohexane		1ml	
δ-HCH D6				
CAS n/a	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-XA14074100CY	δ-HCH D6 100 µg/mL in cylohexane		1.1ml	
ε-HCH				
CAS 6108-10-7	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-LA14075000CY	ε-HCH 10 µg/mL in Cyclohexane(±)		1ml	
DRE-XA14075000CY	ε-HCH 100 µg/mL in Cyclohexane(±)		1ml	
DRE-XA14075000ME	ε-HCH 100 µg/mL in Methanol(±)		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
γ-HCH D6				
CAS 60556-82-3	MW 296.8668	C ₆ H ₆ Cl ₆		
DRE-C14073100	γ-HCH D6		10mg	
DRE-XA14073100CY	γ-HCH D6 100 µg/mL in Cyclohexane		1ml	
α-HCH (alpha-HCH)				
CAS 319-84-6	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-C14071000	α-HCH(±)		100mg	
DRE-L14071000CY	α-HCH 10 µg/mL in Cyclohexane(±)		10ml	
DRE-XA14071000CY	α-HCH 100 µg/mL in Cyclohexane		1ml	
DRE-GA09010365AC	alpha-HCH 100 µg/mL in Acetone(±)		1ml	
Halfenprox				
CAS 111872-58-3	MW 477.3384	C ₂₄ H ₂₃ BrF ₂ O ₃		
DRE-C14059000	Halfenprox(±)		10mg	
DRE-L14059000CY	Halfenprox 10 µg/mL in Cyclohexane		10ml	
Halofenozide				
CAS 112226-61-6	MW 330.8087	C ₁₈ H ₁₉ ClN ₂ O ₂		
DRE-C14059200	Halofenozide		100mg	
Hepachlor-β-dihydro				
CAS 14168-01-5	MW 375.3336	C ₁₀ H ₇ Cl ₇		
DRE-C14090300	Hepachlor-β-dihydro		25mg	
Heptachlor				
CAS 76-44-8	MW 373.3177	C ₁₀ H ₅ Cl ₇		
DRE-C14090000	Heptachlor(±)		100mg	
DRE-L14090000AL	Heptachlor 10 µg/mL in Acetonitrile		10ml	
DRE-L14090000CY	Heptachlor 10 µg/mL in Cyclohexane		10ml	
DRE-XA14090000CY	Heptachlor 100 µg/mL in Cyclohexane		1ml	
(±)-cis-Heptachlorepoxyde (Heptachlor-exo-epoxide, isomer B)				
CAS 1024-57-3	MW 389.3171	C ₁₀ H ₅ Cl ₇ O		
DRE-C14101000	Heptachlor-exo-epoxide (cis-, isomer B)(±)		10mg	
DRE-L14101000CY	Heptachlor-exo-epoxide (Isomer B) 10 µg/mL in Cyclohexane		10ml	
DRE-XA14101000CY	Heptachlor-exo-epoxide (isomer B) 100 µg/mL in Cyclohexane		1ml	
DRE-XA14101000ME	Heptachlor-exo-epoxide (isomer B) 100 µg/mL in Methanol		1ml	
(±)-trans-Heptachlorepoxyde (Heptachlor-endo-epoxide, isomer A)				
CAS 28044-83-9	MW 389.3171	C ₁₀ H ₅ Cl ₇ O		
DRE-C14102000	Heptachlor-endo-epoxide (trans-, isomer A)(±)		10mg	
DRE-L14102000CY	Heptachlor-endo-epoxide (Isomer A) 10 µg/mL in Cyclohexane		10ml	
DRE-XA14102000CY	Heptachlor-endo-epoxide (Isomer A) 100 µg/mL in Cyclohexane		1ml	
DRE-XA14102000ME	Heptachlor-endo-epoxide (Isomer A) 100 µg/mL in Methanol		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Heptenophos				
CAS 23560-59-0	MW 250.6159	C ₉ H ₁₂ ClO ₄ P		
DRE-CA14130000	Heptenophos		100mg	
DRE-XA09010201AL	Heptanophos 100 µg/mL in Acetonitrile(‡)		1ml	
1,2,3,4,5,6-Hexachlorocyclohexane (HCH; BHC)				
CAS 608-73-1	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-C14070000	HCH (BHC) (technical)		100mg	
DRE-L14070000IO	HCH (BHC) (technical) 10 µg/mL in Isooctane		10ml	
Hexaflumuron				
CAS 86479-06-3	MW 461.1427	C ₁₆ H ₈ Cl ₂ F ₆ N ₂ O ₃		
DRE-C14194000	Hexaflumuron		100mg	
DRE-L14194000EA	Hexaflumuron 10 µg/mL in Ethylacetate		10ml	
Hexythiazox				
CAS 78587-05-0	MW 352.8788	C ₁₇ H ₂₁ ClN ₂ O ₂ S		
DRE-C14210000	Hexythiazox(‡)		100mg	
DRE-L14210000TO	Hexythiazox 10 µg/mL in Toluene		10ml	
DRE-XA14210000AL	Hexythiazox 100 µg/mL in Acetonitrile		1ml	
Hydramethylnon				
CAS 67485-29-4	MW 494.4753	C ₂₅ H ₂₄ F ₆ N ₄		
DRE-C14220000	Hydramethylnon		100mg	
DRE-L14220000TO	Hydramethylnon 10 µg/mL in Toluene		10ml	
Imidacloprid				
CAS 138261-41-3	MW 255.661	C ₉ H ₁₀ ClN ₅ O ₂		
DRE-C14283700	Imidacloprid(‡)		100mg	
DRE-GS14283700AL	Imidacloprid 100 µg/mL in Acetonitrile		5x1ml	
DRE-L14283700AL	Imidacloprid 10 µg/mL in Acetonitrile		10ml	
DRE-XA14283700AL	Imidacloprid 100 µg/mL in Acetonitrile(‡)		1ml	
Imidacloprid D4 (imidazolidin-4,4,5,5 D4)				
CAS 1015855-75-0	MW 259.6856	C ₉ ² H ₄ H ₆ ClN ₅ O ₂		
DRE-C14283710	Imidacloprid D4 (imidazolidin-4,4,5,5 D4)		10mg	
DRE-XA14283710AC	Imidacloprid D4 (imidazolidin-4,4,5,5 D4) 100 µg/mL in Acetone(‡)		1.1ml	
Imidacloprid Impurity 1				
CAS 105828-41-9	MW 381.2167	C ₁₅ H ₁₄ Cl ₂ N ₆ O ₂		
DRE-C14283750	Imidacloprid Impurity 1		50mg	
Imidacloprid-guanidine Hydrochloride				
CAS 127202-53-3	MW 247.1244	C ₉ H ₁₁ ClN ₄ ·ClH		
DRE-C14283715	Imidacloprid-guanidine hydrochloride		20mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Imidacloprid-olefin				
CAS 115086-54-9 DRE-C14283760	MW 254.6531 Imidacloprid-olefin	C ₉ H ₁₀ ClN ₅ O ₂	10mg	
Imidacloprid-urea				
CAS 120868-66-8 DRE-C14283780	MW 211.6482 Imidacloprid-urea	C ₉ H ₁₀ ClN ₅ O	50mg	
Imidaclothiz				
CAS 105843-36-5 DRE-C14283850	MW 261.6887 Imidaclothiz	C ₇ H ₈ ClN ₅ O ₂ S	25mg	
Imiprothrin				
CAS 72963-72-5 DRE-C14286000	MW 318.3676 Imiprothrin(‡)	C ₁₇ H ₂₂ N ₂ O ₄	50mg	
Indoxacarb				
CAS 144171-61-9 DRE-C14325500 DRE-L14325500CY	MW 527.8345 Indoxacarb(‡) Indoxacarb 10 µg/mL in Cyclohexane	C ₂₂ H ₁₇ ClF ₃ N ₃ O ₇	100mg 10ml	
Iodofenphos				
CAS 18181-70-9 DRE-C14340000 DRE-XA14340000CY	MW 412.9966 Iodofenphos(‡) Iodofenphos 100 µg/mL in Cyclohexane	C ₈ H ₈ Cl ₂ I O ₃ PS	100mg 1ml	
Isazofos				
CAS 42509-80-8 DRE-C14377000 DRE-L14377000IO	MW 313.7413 Isazofos Isazofos 10 µg/mL in Isooctane	C ₉ H ₁₇ ClN ₃ O ₃ PS	100mg 10ml	
Isobenzan (Telodrin)				
CAS 297-78-9 DRE-C14380000 DRE-L14380000CY DRE-XA14380000ME	MW 411.7515 Isobenzan (Telodrin) Isobenzan 10 µg/mL in Cyclohexane Isobenzan 100 µg/mL in Methanol	C ₉ H ₄ Cl ₆ O	10mg 10ml 1ml	
Isocarbofos				
CAS 24353-61-5 DRE-C14402000 DRE-L14402000CY DRE-XA14402000CY	MW 289.2878 Isocarbofos(‡) Isocarbofos 10 µg/mL in Cyclohexane Isocarbofos 100 µg/mL in Cyclohexane(‡)	C ₁₁ H ₁₆ NO ₄ PS	100mg 10ml 1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Isodrin				
CAS 465-73-6	MW 364.9099	C ₁₂ H ₈ Cl ₆		AND Enantiomer 
DRE-C14410000	Isodrin		100mg	
DRE-L14410000CY	Isodrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA14410000IQ	Isodrin 100 µg/mL in Isooctane		1ml	
Isofenphos				
CAS 25311-71-1	MW 345.3941	C ₁₅ H ₂₄ NO ₄ PS		
DRE-C14420000	Isofenphos(‡)		100mg	
DRE-L14420000CY	Isofenphos 10 µg/mL in Cyclohexane		10ml	
DRE-XA14420000CY	Isofenphos 100 µg/mL in Cyclohexane		1ml	
Isofenphos-des-N-isopropyl				
CAS 25205-08-7	MW 303.3144	C ₁₂ H ₁₈ NO ₄ PS		
DRE-C14422000	Isofenphos-des-N-isopropyl		100mg	
Isofenphos-methyl				
CAS 99675-03-3	MW 331.3675	C ₁₄ H ₂₂ NO ₄ PS		
DRE-C14421000	Isofenphos-methyl(‡)		50mg	
DRE-XA14421000AL	Isofenphos-methyl 100 µg/mL in Acetonitrile		1ml	
Isofenphos-oxon				
CAS 31120-85-1	MW 329.3285	C ₁₅ H ₂₄ NO ₅ P		
DRE-C14423000	Isofenphos-oxon		10mg	
Isofenphos-des-N-isopropyl-oxon (Isopropyl 2-[Amino(ethoxy)phosphoryl]oxybenzoate)				
CAS 31120-83-9	MW 287.2488	C ₁₂ H ₁₈ NO ₅ P		
DRE-C14422200	Isofenphos-des-N-isopropyl-oxon (*)		10mg	
Isomalathion				
CAS 3344-12-5	MW 330.358	C ₁₀ H ₁₈ O ₆ PS ₂		
DRE-C14710100	iso-Malathion		50mg	
Isoprobenphos (O,O-Diethyl S-benzyl Phosphorothioate)				
CAS 13286-32-3	MW 260.2896	C ₁₁ H ₁₇ O ₃ PS		
DRE-C14449500	Isoprobenphos		50mg	
Isoprocarb				
CAS 2631-40-5	MW 193.2423	C ₁₁ H ₁₅ NO ₂		
DRE-C14450000	Isoprocarb		100mg	










Pesticide and metabolites: Insecticides

Product code	Description			
Isoxathion				
CAS 18854-01-8	MW 313.3092	C ₁₃ H ₁₆ NO ₄ PS		
DRE-C14483000	Isoxathion		100mg	
DRE-L14483000CY	Isoxathion 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010207AL	Isoxathion 100 µg/mL in Acetonitrile(‡)		1ml	
Kadethrin				
CAS 58769-20-3	MW 396.4993	C ₂₃ H ₂₄ O ₄ S		
DRE-C14500000	Kadethrin		100mg	
Karanjin				
CAS 521-88-0	MW 292.2855	C ₁₈ H ₁₂ O ₄		
DRE-XA14509000AL	Karanjin 100 µg/mL in Acetonitrile		1ml	
Kelevan				
CAS 4234-79-1	MW 634.8048	C ₁₇ H ₁₂ Cl ₁₀ O ₄		
DRE-C14530000	Kelevan		100mg	
Kinoprene				
CAS 42588-37-4	MW 276.4137	C ₁₈ H ₂₈ O ₂		
DRE-C14538000	Kinoprene		100mg	
Lead Arsenate				
CAS 7784-40-9	MW 347.1271	Pb·AsHO ₄		
DRE-C14605000	Lead arsenate		100mg	
Leptophos				
CAS 21609-90-5	MW 412.0661	C ₁₃ H ₁₀ BrCl ₂ O ₂ PS		
DRE-C14620000	Leptophos		100mg	
Leptophos-oxon				
CAS 25006-32-0	MW 396.0005	C ₁₃ H ₁₀ BrCl ₂ O ₃ P		
DRE-C14625000	Leptophos-oxon		10mg	
Lindane (γ-HCH)				
CAS 58-89-9	MW 290.8298	C ₆ H ₆ Cl ₆		
DRE-C14073000	γ-HCH(‡)		250mg	
DRE-L14073000CY	γ-HCH 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-XA14073000CY	γ-HCH 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Lufenuron Anhydrous				
CAS 103055-07-8 DRE-C14650000	MW 511.1502 Lufenuron(±)	C ₁₇ H ₆ Cl ₂ F ₈ N ₂ O ₃	100mg	
Malaoxon				
CAS 1634-78-2 DRE-C14700000	MW 314.2924 Malaoxon	C ₁₀ H ₁₈ O ₇ PS	100mg	
Malathion				
CAS 121-75-5 DRE-C14710000 DRE-CR14710000 DRE-L14710000CY DRE-XA14710000CY	MW 330.358 Malathion(±) Malathion Malathion 10 µg/mL in Cyclohexane Malathion 100 µg/mL in Cyclohexane	C ₁₀ H ₁₈ O ₆ PS ₂	100mg 50mg 10ml 1ml	
Malathion D10 (diethyl D10)				
CAS 347841-48-9 DRE-C14710010	MW 340.4196 Malathion D10 (diethyl D10)	C ₁₀ ² H ₁₀ H ₈ O ₆ PS ₂	20mg	
Malathion D6 (dimethyl D6)				
CAS 1189877-72-2 DRE-XA14710020CY DRE-GH09010100AL	MW 336.395 Malathion D6 (dimethyl D6) 100 µg/mL in Cyclohexane (*) Malathion-d6 100 µg/mL in Acetonitrile(±)	C ₁₀ ² H ₆ H ₁₃ O ₆ PS ₂	1.1ml 10x1ml	
Malathion D7 (dimethyl D6,3-D1)				
CAS 352438-94-9 DRE-C14710030	MW 337.4012 Malathion D7 (dimethyl D6,3-D1)	C ₁₀ ² H ₇ H ₁₂ O ₆ PS ₂	20mg	
Malathion dicarboxylic acid				
CAS 1190-28-9 DRE-C14713000	MW 274.2517 Malathion dicarboxylic acid	C ₈ H ₁₁ O ₆ PS ₂	10mg	
Mecarbam				
CAS 2595-54-2 DRE-C14800000 DRE-L14800000AL DRE-L14800000CY	MW 329.3733 Mecarbam Mecarbam 10 µg/mL in Acetonitrile Mecarbam 10 µg/mL in Cyclohexane	C ₁₀ H ₂₀ NO ₅ PS ₂	100mg 10ml 10ml	
Mecarbam D10 (diethoxy D10)				
CAS n/a DRE-XA14800100AC	MW 339.4349 Mecarbam D10 (diethoxy D10) 100 µg/mL in Acetone	C ₁₀ ² H ₁₀ H ₁₀ NO ₅ PS ₂	1.1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Mephosfolan				
CAS 950-10-7 DRE-C14870000	MW 269.3213 Mephosfolan	$C_8H_{16}NO_3PS_2$	100mg	
Metaflumizone				
CAS 139968-49-3 DRE-C14918500 DRE-L14918500CY	MW 506.3999 Metaflumizone Metaflumizone 10 µg/mL in Cyclohexane	$C_{24}H_{16}F_6N_4O_2$	100mg 10ml	
Methacrifos				
CAS 62610-77-9 DRE-C14970000 DRE-L14970000CY DRE-XA14970000CY DRE-XA09010210AL	MW 240.2139 Methacrifos Methacrifos 10 µg/mL in Cyclohexane Methacrifos 100 µg/mL in Cyclohexane Methacrifos 100 µg/mL in Acetonitrile(‡)	$C_7H_{13}O_3PS$	100mg 10ml 1ml 1ml	
Methamidophos				
CAS 10265-92-6 DRE-C14980000 DRE-L14980000AL DRE-XA14980000EA	MW 141.1292 Methamidophos(‡) Methamidophos 10 µg/mL in Acetonitrile Methamidophos 100 µg/mL in Ethylacetate	$C_2H_5NO_2PS$	100mg 10ml 1ml	
Methamidophos (dimethyl D6)				
CAS 1219799-41-3 DRE-C14980100	MW 147.1662 Methamidophos D6 (dimethyl D6)	$C_2^2H_6H_2NO_2PS$	10mg	
Methamidophos-N-methyl (O,S-Dimethylmethylphosphoramidothioate)				
CAS 28167-49-9 DRE-C14981000	MW 155.1558 Methamidophos-N-methyl	$C_3H_{10}NO_2PS$	10mg	
Methidathion				
CAS 950-37-8 DRE-C15020000 DRE-L15020000CY DRE-XA15020000CY DRE-LA09010258CY	MW 302.3313 Methidathion Methidathion 10 µg/mL in Cyclohexane Methidathion 100 µg/mL in Cyclohexane Methidathion 10 µg/mL in Cyclohexane(‡)	$C_8H_{11}N_2O_4PS_3$	100mg 10ml 1ml 10ml	
Methiocarb (Mercaptodimethur)				
CAS 2032-65-7 DRE-C15020500 DRE-L15020500CY DRE-XA15020500CY	MW 225.3073 Methiocarb (Mercaptodimethur)(‡) Methiocarb 10 µg/mL in Cyclohexane Methiocarb 100 µg/mL in Cyclohexane	$C_{11}H_{15}NO_2S$	100mg 10ml 1ml	
Methiocarb D3 (N-methyl D3)				
CAS 1581694-94-1 DRE-XA15020501CY	MW 228.3258 Methiocarb D3 (N-methyl D3) 100 µg/mL in Cyclohexane	$C_{11}^2H_5H_{12}NO_2S$	1.1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Methiocarb-phenol-sulfone				
CAS 14763-62-3 DRE-C15020503	MW 200.2548 Methiocarb-phenol-sulfone	C ₉ H ₁₂ O ₃ S	25mg	
Methiocarb-phenol-sulfoxide				
CAS 22454-92-8 DRE-C15020506	MW 184.2554 Methiocarb-phenol-sulfoxide	C ₉ H ₁₂ O ₂ S	25mg	
Methiocarb-sulfone				
CAS 2179-25-1 DRE-C15020510 DRE-XA15020510MB	MW 257.3061 Methiocarb sulfone(‡) Methiocarb-sulfone 100 µg/mL in Methyl-tert. butyl ether	C ₁₁ H ₁₈ NO ₄ S	50mg 1ml	
Methiocarb sulfoxide				
CAS 2635-10-1 DRE-C15020520	MW 241.3067 Methiocarb sulfoxide	C ₁₁ H ₁₈ NO ₃ S	100mg	
Methomyl				
CAS 16752-77-5 DRE-C15030000 DRE-L15030000CY DRE-XA15030000CY	MW 162.2101 Methomyl(‡) Methomyl 10 µg/mL in Cyclohexane Methomyl 100 µg/mL in Cyclohexane	C ₈ H ₁₀ N ₂ O ₂ S	100mg 10ml 1ml	
Methomyl D3				
CAS 1398109-07-3 DRE-XA15030100AC	MW 165.2286 Methomyl D3 100 µg/mL in Acetone	C ₈ ² H ₉ N ₂ O ₂ S	1.1ml	
Methomyl-oxime				
CAS 13749-94-5 DRE-C15035000	MW 105.1588 Methomyl-oxime	C ₃ H ₇ NOS	100mg	
Methomyl-sulfone				
CAS 55620-24-1 DRE-LA15035015AL	MW 194.2089 Methomyl-sulfone 10 µg/mL in Acetonitrile	C ₈ H ₁₀ N ₂ O ₄ S	1ml	
Methoprene				
CAS 40596-69-8 DRE-C15045000 DRE-L15045000CY	MW 310.4715 Methoprene(‡) Methoprene 10 µg/mL in Cyclohexane	C ₁₉ H ₃₄ O ₃	100mg 10ml	










Pesticide and metabolites: Insecticides

Product code	Description			
trans-S-Methoprene-epoxide				
CAS 65277-96-5 DRE-LA15045300CY	MW 300.4567 trans-S-Methoprene-epoxide 10 µg/mL in Cyclohexane	C ₁₆ H ₂₆ O ₃ S	1ml	
S-Methoprene				
CAS 65733-16-6 DRE-C15045200	MW 310.4715 (S)-Methoprene	C ₁₉ H ₃₄ O ₃	10mg	
Methothrin				
CAS 34388-29-9 DRE-C15055000	MW 302.4079 Methothrin	C ₁₉ H ₂₆ O ₃	100mg	
Methoxychlor				
CAS 72-43-5 DRE-C15060000 DRE-L15060000AL DRE-L15060000IO DRE-XA15060000CY	MW 345.6481 Methoxychlor Methoxychlor 10 µg/mL in Acetonitrile Methoxychlor 10 µg/mL in Isooctane Methoxychlor 100 µg/mL in Cyclohexane	C ₁₆ H ₁₅ Cl ₃ O ₂	100mg 10ml 10ml 1ml	
Methoxychlor D14 (bis(4-methoxyphenyl-D7))				
CAS 1644449-82-0 DRE-C15060100 DRE-XA15060100AC	MW 359.7344 Methoxychlor D14 (bis(4-methoxyphenyl-D7)) Methoxychlor D14 100 µg/mL in Acetone	C ₁₆ ² H ₁₄ HCl ₃ O ₂	10mg 1.1ml	
2,4'-Methoxychlor				
CAS 30667-99-3 DRE-C15061000	MW 345.6481 2,4'-Methoxychlor	C ₁₆ H ₁₅ Cl ₃ O ₂	10mg	
4,4'-Methoxychlor olefin (1,1-Dichloro-2,2-bis(p-methoxyphenyl)ethene)				
CAS 2132-70-9 DRE-C15062000 DRE-L15062000CY	MW 309.1872 4,4'-Methoxychlor olefin 4,4'-Methoxychlor olefin 10 µg/mL in Cyclohexane	C ₁₆ H ₁₄ Cl ₂ O ₂	10mg 10ml	
Methoxyfenozide				
CAS 161050-58-4 DRE-C15080500 DRE-L15080500AL	MW 368.4693 Methoxyfenozide(‡) Methoxyfenozide 10 µg/mL in Acetonitrile	C ₂₂ H ₂₈ N ₂ O ₃	100mg 10ml	
3-Methylphosphinicopropionic Acid				
CAS 15090-23-0 DRE-XA15141200AL	MW 152.0856 3-Methylphosphinicopropionic acid 100 µg/mL in Acetonitrile	C ₄ H ₉ O ₄ P	1ml	










Pesticide and metabolites: Insecticides

Product code	Description				
Metofluthrin					
CAS 240494-70-6 DRE-C15167000	MW 360.3432 Metofluthrin	C ₁₈ H ₂₀ F ₄ O ₃	10mg		
Metolcarb					
CAS 1129-41-5 DRE-C15175000 DRE-L15175000CY DRE-XA09010211AL	MW 165.1891 Metolcarb Metolcarb 10 µg/mL in Cyclohexane Metolcarb 100 µg/mL in Acetonitrile(‡)	C ₉ H ₁₁ NO ₂	100mg 10ml 1ml		
Metrifonate (Trichlorphon)					
CAS 52-68-6 DRE-C17680000	MW 257.4367 Trichlorphon(‡)	C ₄ H ₈ Cl ₃ O ₄ P	250mg		
Mevinphos					
CAS 7786-34-7 DRE-C15220000 DRE-L15220000AL DRE-XA15220000AL	MW 224.1483 Mevinphos(‡) Mevinphos 10 µg/mL in Acetonitrile Mevinphos 100 µg/mL in Acetonitrile	C ₇ H ₁₃ O ₆ P	100mg 10ml 1ml		
E-Mevinphos (cis-butenic acid)					
CAS 26718-65-0 DRE-C15221000	MW 224.1483 E-Mevinphos (cis-butenic acid)	C ₇ H ₁₃ O ₆ P	50mg		
Z-Mevinphos					
CAS 338-45-4 DRE-C15222000	MW 224.1483 Z-Mevinphos (trans-butenic acid)	C ₇ H ₁₃ O ₆ P	50mg		
Mexacarbate					
CAS 315-18-4 DRE-C15230000	MW 222.2835 Mexacarbate	C ₁₂ H ₁₈ N ₂ O ₂	25mg		
Milbemectin A3					
CAS 51596-10-2 DRE-L15265020AL	MW 528.6769 Milbemectin A3 10 µg/mL in Acetonitrile (*)	C ₃₁ H ₄₄ O ₇	10ml		
Milbemectin A4					
CAS 51596-11-3 DRE-L15265040AL	MW 542.7034 Milbemectin A4 10 µg/mL in Acetonitrile (*)	C ₃₂ H ₄₆ O ₇	10ml		









Pesticide and metabolites: Insecticides

Product code	Description			
Mirex				
CAS 2385-85-5	MW 545.543	C ₁₀ Cl ₁₂		
DRE-C15270000	Mirex		100mg	
DRE-L15270000IQ	Mirex 10 µg/mL in Isooctane		10ml	
DRE-XA15270000IQ	Mirex 100 µg/mL in Isooctane		1ml	
Mitotane (2,4'-DDD)				
CAS 53-19-0	MW 320.0412	C ₁₄ H ₁₀ Cl ₄		
DRE-C12030000	2,4'-DDD		100mg	
DRE-L12030000CY	2,4'-DDD 10 µg/mL in Cyclohexane		10ml	
DRE-XA12030000CY	2,4'-DDD 100 µg/mL in Cyclohexane		1ml	
Momfluorothrin				
CAS 609346-29-4	MW 385.3527	C ₁₉ H ₁₉ F ₄ NO ₃		
DRE-C15285000	Momfluorothrin		10mg	
Monocrotophos				
CAS 6923-22-4	MW 223.1635	C ₇ H ₁₄ NO ₅ P		
DRE-C15300000	Monocrotophos		100mg	
DRE-L15300000AL	Monocrotophos 10 µg/mL in Acetonitrile		10ml	
DRE-L15300000IQ	Monocrotophos 10 µg/mL in Isooctane		10ml	
DRE-XA15300000AL	Monocrotophos 100 µg/mL in Acetonitrile		1ml	
Monosultap				
CAS 29547-00-0	MW 333.4016	C ₅ H ₁₂ NO ₆ S ₄ ·Na		
DRE-C15313000	Monosultap		100mg	
Naftalofos				
CAS 1491-41-4	MW 349.2751	C ₁₆ H ₁₆ NO ₆ P		
DRE-C15405000	Naftalofos		50mg	
Naled (Dibrom)				
CAS 300-76-5	MW 380.7837	C ₄ H ₇ Br ₂ Cl ₂ O ₄ P		
DRE-C15410000	Naled (*)		250mg	
DRE-GA09010403AL	Naled (Dibrom) Solution, 100 µg/mL in Acetonitrile(‡)		1ml	
Nereistoxin Oxalate D6 (dimethyl D6)				
CAS n/a	MW 245.3494	C ₅ ² H ₆ H ₅ NS ₂ ·C ₂ H ₂ O ₄		
DRE-C15502010	Nereistoxin oxalate D6 (dimethyl D6)		10mg	
Nereistoxin oxalate				
CAS 1631-52-3	MW 239.3124	C ₅ H ₁₁ NS ₂ ·C ₂ H ₂ O ₄		
DRE-C15502000	Nereistoxin oxalate		100mg	









Pesticide and metabolites: Insecticides

Product code	Description			
(E)-Nitenpyram				
CAS 150824-47-8	MW 270.7154	C ₁₁ H ₁₅ ClN ₄ O ₂		
DRE-C15535000	Nitenpyram		100mg	
DRE-L15535000AL	Nitenpyram 10 µg/mL in Acetonitrile		10ml	
cis-Nonachlor				
CAS 5103-73-1	MW 444.2237	C ₁₀ H ₅ Cl ₉		
DRE-C15620100	cis-Nonachlor		10mg	
DRE-L15620100CY	cis-Nonachlor 10 µg/mL in Cyclohexane		10ml	
DRE-XA15620100IO	cis-Nonachlor 100 µg/mL in Isooctane		1ml	
trans-Nonachlor				
CAS 39765-80-5	MW 444.2237	C ₁₀ H ₅ Cl ₉		
DRE-C15620200	trans-Nonachlor(‡)		10mg	
DRE-L15620200CY	trans-Nonachlor 10 µg/mL in Cyclohexane		10ml	
DRE-XA15620200CY	trans-Nonachlor 100 µg/mL in Cyclohexane		1ml	
Novaluron				
CAS 116714-46-6	MW 492.7046	C ₁₇ H ₉ ClF ₈ N ₂ O ₄		
DRE-C15653000	Novaluron(‡)		50mg	
Noviflumuron				
CAS 121451-02-3	MW 529.1407	C ₁₇ H ₇ Cl ₂ F ₉ N ₂ O ₃		
DRE-C15654000	Noviflumuron		10mg	
Oxamyl				
CAS 23135-22-0	MW 219.2614	C ₇ H ₁₃ N ₃ O ₃ S		
DRE-C15780000	Oxamyl		100mg	
DRE-L15780000MB	Oxamyl 10 µg/mL in Methyl-tert. butyl ether (*)		10ml	
Oxamyl D3 (2-methyl-D3)				
CAS n/a	MW 222.2799	C ₇ ² H ₉ H ₁₀ N ₃ O ₃ S		
DRE-XA15780100MB	Oxamyl D3 (2-methyl-D3) 100 µg/mL in Methyl-tert.butyl ether (*)		1.1ml	
Oxamyl-oxime				
CAS 30558-43-1	MW 162.2101	C ₈ H ₁₀ N ₂ O ₃ S		
DRE-C15780500	Oxamyl-oxime(‡)		10mg	
DRE-LA15780500AL	Oxamyl-oxime 10 µg/mL in Acetonitrile		1ml	
Paraoxon (Paraoxon-ethyl)				
CAS 311-45-5	MW 275.195	C ₁₀ H ₁₄ NO ₆ P		
DRE-C15850000	Paraoxon-ethyl		100mg	
DRE-L15850000CY	Paraoxon-ethyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA15850000CY	Paraoxon-ethyl 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Paraoxon-methyl				
CAS 950-35-6	MW 247.1419	C ₈ H ₁₀ NO ₅ P		
DRE-CA15860000	Paraoxon-methyl		100mg	
DRE-L15860000IQ	Paraoxon-methyl 10 µg/mL in Isooctane		10ml	
Parathion-ethyl				
CAS 56-38-2	MW 291.2606	C ₁₀ H ₁₄ NO ₅ PS		
DRE-C15880000	Parathion-ethyl(‡)		100mg	
DRE-L15880000AL	Parathion-ethyl 10 µg/mL in Acetonitrile		10ml	
DRE-L15880000CY	Parathion-ethyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA15880000AL	Parathion-ethyl 100 µg/mL in Acetonitrile		1ml	
DRE-XA15880000CY	Parathion-ethyl 100 µg/mL in Cyclohexane		1ml	
DRE-YS09010031HE	Parathion 1000 µg/mL in Hexane(‡)		5x1ml	
DRE-GA09010344ME	Parathion 1000 µg/mL in Methanol(‡)		1ml	
Parathion-ethyl D10 (diethyl D10)				
CAS 350820-04-1	MW 301.3222	C ₁₀ ² H ₁₀ H ₄ NO ₅ PS		
DRE-C15880100	Parathion-ethyl D10 (diethyl D10)		10mg	
DRE-XA15880100AC	Parathion-ethyl D10 (diethyl D10) 100 µg/mL in Acetone		1.1ml	
Parathion-methyl				
CAS 298-00-0	MW 263.2075	C ₈ H ₁₀ NO ₅ PS		
DRE-C15890000	Parathion-methyl(‡)		100mg	
DRE-L15890000CY	Parathion-methyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA15890000AL	Parathion-methyl 100 µg/mL in Acetonitrile (*)		1ml	
DRE-XA15890000TO	Parathion-methyl 100 µg/mL in Toluene		1ml	
Parathion-methyl D6 (dimethyl D6)				
CAS 96740-32-8	MW 269.2444	C ₈ ² H ₆ H ₄ NO ₅ PS		
DRE-C15890100	Parathion-methyl D6 (dimethyl D6)		25mg	
DRE-XA15890100AC	Parathion-methyl D6 (dimethyl D6) 100 µg/mL in Acetone (*)		1.1ml	
Penfluron				
CAS 35367-31-8	MW 344.2362	C ₁₅ H ₉ F ₅ N ₂ O ₂		
DRE-C15934000	Penfluron		100mg	
Pentachloroaniline				
CAS 527-20-8	MW 265.3518	C ₆ H ₂ Cl ₅ N		
DRE-L15940000CY	Pentachloroaniline 10 µg/mL in Cyclohexane		10ml	
Pentachlorobenzene				
CAS 608-93-5	MW 250.3371	C ₆ HCl ₅		
DRE-L15960000AL	Pentachlorobenzene 10 µg/mL in Acetonitrile		10ml	
DRE-L15960000CY	Pentachlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA15960000CY	Pentachlorobenzene 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Pentachlorophenol				
CAS 87-86-5	MW 266.3365	C ₆ HCl ₅ O		
DRE-C15970000	Pentachlorophenol		100mg	
DRE-L15970000AL	Pentachlorophenol 10 µg/mL in Acetonitrile		10ml	
DRE-L15970000CY	Pentachlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-L15970000ME	Pentachlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA15970000ME	Pentachlorophenol 100 µg/mL in Methanol		1ml	
DRE-YS09010008ME	Pentachlorophenol 1000 µg/mL in Methanol(‡)		5x1ml	
Pentachlorophenol 13C6				
CAS 85380-74-1	MW 272.2925	¹³ C ₆ HCl ₅ O		
DRE-C15970100	Pentachlorophenol 13C6		10mg	
DRE-XA15970100CY	Pentachlorophenol 13C6 100 µg/mL in Cyclohexane		1.1ml	
DRE-GS09010309ME	Pentachlorophenol-13C6 1000 µg/mL in Methanol(‡)		5x1ml	
Pentachlorophenol Acetate				
CAS 1441-02-7	MW 308.3732	C ₈ H ₃ Cl ₅ O ₂		
DRE-C15971000	Pentachlorophenol acetate		100mg	
DRE-L15971000CY	Pentachlorophenol acetate 10 µg/mL in Cyclohexane		10ml	
DRE-XA15971000CY	Pentachlorophenol acetate 100 µg/mL in Cyclohexane		1ml	
Pentachlorophenol Sodium Salt Hydrate				
CAS 123333-54-0	MW 306.3337	C ₆ Cl ₅ O·Na·H ₂ O		
DRE-C15973000	Pentachlorophenol sodium salt hydrate		250mg	
Permethrin				
CAS 52645-53-1	MW 391.2877	C ₂₁ H ₂₀ Cl ₂ O ₃		
DRE-C15990000	Permethrin(‡)		250mg	
DRE-L15990000CY	Permethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA15990000CY	Permethrin 100 µg/mL in Cyclohexane		1ml	
DRE-X15990000CY	Permethrin 100 µg/mL in Cyclohexane		10ml	
DRE-GA09010373AC	Total Permethrin 100 µg/mL in Acetone(‡)		1ml	
cis-Permethrin				
CAS 61949-76-6	MW 391.2877	C ₂₁ H ₂₀ Cl ₂ O ₃		
DRE-C15990100	cis-Permethrin(‡)		10mg	
DRE-L15990100AL	cis-Permethrin 10 µg/mL in Acetonitrile		10ml	
DRE-L15990100CY	cis-Permethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA15990100AL	cis-Permethrin 100 µg/mL in Acetonitrile		1ml	
DRE-XA15990100CY	cis-Permethrin 100 µg/mL in Cyclohexane		1ml	
trans-Permethrin				
CAS 61949-77-7	MW 391.2877	C ₂₁ H ₂₀ Cl ₂ O ₃		
DRE-C15990200	trans-Permethrin(‡)		10mg	
DRE-L15990200AL	trans-Permethrin 10 µg/mL in Acetonitrile		10ml	
DRE-L15990200CY	trans-Permethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA15990200CY	trans-Permethrin 100 µg/mL in Cyclohexane		1ml	
trans-Permethrin D6 (dimethyl D6)				
CAS 82523-59-9	MW 397.3247	C ₂₁ ² H ₆ H ₁₄ Cl ₂ O ₃		
DRE-XA15990201AC	trans-Permethrin D6 (dimethyl D6) 100 µg/mL in Acetone(‡)		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Permethrinic Acid (3-(2,2-Dichlorovinyl)-2,2-dimethylcyclopropanecarboxylic Acid)				
CAS 55701-05-8	MW 209.0698	C ₈ H ₁₀ Cl ₂ O ₂		
DRE-LA12507500ME	3-(2,2-Dichlorovinyl)-2,2-dimethyl(1-cyclopropane)carboxylic 10 µg/mL in Methanol		1ml	
Permethrinic Acid Methyl Ester (Methyl 3-(2,2-Dichlorovinyl)-2,2-dimethylcyclopropanecarboxylate)				
CAS 61898-95-1	MW 223.0964	C ₉ H ₁₂ Cl ₂ O ₂		
DRE-C15086000	Methyl-3-(2,2-dichlorovinyl)-2,2-dimethyl-(1-cyclopropane)carboxylate		100mg	
Perthane				
CAS 72-56-0	MW 307.2574	C ₁₈ H ₂₀ Cl ₂		
DRE-C16000000	Perthane		100mg	
DRE-L16000000CY	Perthane 10 µg/mL in Cyclohexane		10ml	
Phenkapton				
CAS 2275-14-1	MW 377.3104	C ₁₁ H ₁₅ Cl ₂ O ₂ PS ₃		
DRE-L16010000CY	Phenkapton 10 µg/mL in Cyclohexane(±)		10ml	
DRE-LA16010000CY	Phenkapton 10 µg/mL in Cyclohexane		1ml	
Phenothrin				
CAS 26002-80-2	MW 350.4507	C ₂₃ H ₂₆ O ₃		
DRE-C16040000	Phenothrin(±)		100mg	
DRE-L16040000IO	Phenothrin 10 µg/mL in Isooctane		10ml	
D-trans-Phenothrin				
CAS 26046-85-5	MW 350.4507	C ₂₃ H ₂₆ O ₃		
DRE-C16041000	D-trans-Phenothrin		100mg	
2-Phenoxybenzoic Acid				
CAS 2243-42-7	MW 214.2167	C ₁₃ H ₁₀ O ₃		
DRE-C16045100	2-Phenoxybenzoic acid		100mg	
3-Phenoxybenzoic Acid				
CAS 3739-38-6	MW 214.2167	C ₁₃ H ₁₀ O ₃		
DRE-C16045200	3-Phenoxybenzoic acid		100mg	
3-Phenoxybenzyl Alcohol				
CAS 13826-35-2	MW 200.2332	C ₁₃ H ₁₂ O ₂		
DRE-C16045300	3-Phenoxybenzyl alcohol		100mg	









Pesticide and metabolites: Insecticides

Product code	Description			
Phenthoate				
CAS 2597-03-7	MW 320.3647	C ₁₂ H ₁₇ O ₄ PS ₂		
DRE-C16050000	Phenthoate(±)		100mg	
DRE-L16050000IQ	Phenthoate 10 µg/mL in Isooctane		10ml	
DRE-XA16050000IQ	Phenthoate 100 µg/mL in Isooctane		1ml	
Phorate				
CAS 298-02-2	MW 260.3774	C ₇ H ₁₇ O ₂ PS ₃		
DRE-C16080000	Phorate		100mg	
DRE-L16080000CY	Phorate 10 µg/mL in Cyclohexane		10ml	
DRE-XA16080000CY	Phorate 100 µg/mL in Cyclohexane		1ml	
Phorate (diethyl-D10)				
CAS 1219805-45-4	MW 270.4391	C ₇ ² H ₁₀ H ₇ O ₂ PS ₃		
DRE-XA16080100AC	Phorate D10 100 µg/mL in Acetone		1.1ml	
Phorate-oxon				
CAS 2600-69-3	MW 244.3118	C ₇ H ₁₇ O ₃ PS ₂		
DRE-C16085000	Phorate-oxon		25mg	
DRE-LA16085000CY	Phorate-oxon 10 µg/mL in Cyclohexane		1ml	
DRE-XA16085000CY	Phorate-oxon 100 µg/mL in Cyclohexane		1ml	
Phorate-oxon-sulfone				
CAS 2588-06-9	MW 276.3106	C ₇ H ₁₇ O ₅ PS ₂		
DRE-C16085500	Phorate-oxon-sulfone		10mg	
DRE-XA16085500CY	Phorate-oxon-sulfone 100 µg/mL in Cyclohexane(±)		1ml	
Phorate-oxon-sulfoxide				
CAS 2588-05-8	MW 260.3112	C ₇ H ₁₇ O ₄ PS ₂		
DRE-C16086000	Phorate-oxon-sulfoxide		10mg	
DRE-XA16086000CY	Phorate-oxon-sulfoxide 100 µg/mL in Cyclohexane		1ml	
Phorate-sulfone				
CAS 2588-04-7	MW 292.3762	C ₇ H ₁₇ O ₄ PS ₃		
DRE-C16088000	Phorate-sulfone(±)		100mg	
DRE-L16088000CY	Phorate sulfone 10 µg/mL in Cyclohexane		10ml	
Phorate-sulfoxide				
CAS 2588-03-6	MW 276.3768	C ₇ H ₁₇ O ₃ PS ₃		
DRE-CA16089000	Phorate-sulfoxide(±)		100mg	
DRE-L16089000CY	Phorate-sulfoxide 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010221ME	Phorate Sulfoxide 100 µg/mL in Methanol(±)		1ml	
Phosalone				
CAS 2310-17-0	MW 367.8086	C ₁₂ H ₁₅ ClNO ₄ PS ₂		
DRE-C16100000	Phosalone(±)		100mg	
DRE-L16100000AL	Phosalone 10 µg/mL in Acetonitrile		10ml	
DRE-L16100000CY	Phosalone 10 µg/mL in Cyclohexane		10ml	
DRE-XA16100000CY	Phosalone 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Phosalone D10 (di-ethyl D5)				
CAS n/a	MW 377.8702	C ₁₂ ² H ₁₀ H ₅ ClNO ₄ PS ₂		
DRE-XA16100100AC	Phosalone D10 (di-ethyl D5) 100 µg/mL in Acetone		1.1ml	
Phosfolan				
CAS 947-02-4	MW 255.2947	C ₇ H ₁₄ NO ₃ PS ₂		
DRE-C16110000	Phosfolan(±)		100mg	
DRE-XA16110000CY	Phosfolan 100 µg/mL in Cyclohexane		1ml	
Phosfolan-methyl				
CAS 5120-23-0	MW 227.2416	C ₈ H ₁₀ NO ₃ PS ₂		
DRE-C16115000	Phosfolan-methyl		25mg	
Phosmet				
CAS 732-11-6	MW 317.321	C ₁₁ H ₁₂ NO ₄ PS ₂		
DRE-C16120000	Phosmet(±)		100mg	
DRE-L16120000CY	Phosmet 10 µg/mL in Cyclohexane		10ml	
DRE-XA16120000CY	Phosmet 100 µg/mL in Cyclohexane		1ml	
Phosmet D6 (dimethoxy D3)				
CAS 2083623-41-8	MW 323.358	C ₁₁ ² H ₆ H ₆ NO ₄ PS ₂		
DRE-XA16120100AC	Phosmet D6 100 µg/mL in Acetone		1.1ml	
Phosmet-oxon				
CAS 3735-33-9	MW 301.2554	C ₁₁ H ₁₂ NO ₅ PS		
DRE-C16125000	Phosmet-oxon		10mg	
DRE-L16125000IO	Phosmet-oxon 10 µg/mL in Isooctane (*)		10ml	
DRE-XA16125000IO	Phosmet oxon 100 µg/mL in Isooctane(±) (*)		1ml	
Phosphamidon				
CAS 13171-21-6	MW 299.6883	C ₁₀ H ₁₉ ClNO ₅ P		
DRE-C16140000	Phosphamidon		100mg	
DRE-L16140000CY	Phosphamidon 10 µg/mL in Cyclohexane		10ml	
DRE-XA16140000CY	Phosphamidon 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010222AL	Phosphamidon 100 µg/mL in Acetonitrile(±)		1ml	
Photoheptachlor				
CAS 33442-83-0	MW 373.3177	C ₁₀ H ₅ Cl ₇		
DRE-ZA16147600CY	Photoheptachlor, 1 µg/mL in Cyclohexane		1ml	
Phoxim				
CAS 14816-18-3	MW 298.2979	C ₁₂ H ₁₅ N ₂ O ₃ PS		
DRE-C16150000	Phoxim(±)		100mg	
DRE-L16150000CY	Phoxim 10 µg/mL in Cyclohexane		10ml	
DRE-XA16150000CY	Phoxim 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Piperonyl Butoxide				
CAS 51-03-6	MW 338.4385	C ₁₉ H ₃₀ O ₅		
DRE-C16240000	Piperonyl butoxide(±)		100mg	
DRE-GA16240000AL	Piperonyl butoxide 100 µg/mL in Acetonitrile		1ml	
DRE-GS16240000AL	Piperonyl butoxide 100 µg/mL in Acetonitrile		5x1ml	
DRE-L16240000AL	Piperonyl butoxide 10 µg/mL in Acetonitrile		10ml	
DRE-L16240000CY	Piperonyl butoxide 10 µg/mL in Cyclohexane		10ml	
Pirimetaphos				
CAS 31377-69-2	MW 288.2832	C ₁₁ H ₂₁ N ₄ O ₃ P		
DRE-C16249000	Pirimethaphos		50mg	
Pirimicarb				
CAS 23103-98-2	MW 238.2862	C ₁₁ H ₁₈ N ₄ O ₂		
DRE-C16250000	Pirimicarb(±)		250mg	
DRE-L16250000AL	Pirimicarb 10 µg/mL in Acetonitrile		10ml	
DRE-L16250000CY	Pirimicarb 10 µg/mL in Cyclohexane		10ml	
DRE-XA16250000AL	Pirimicarb 100 µg/mL in Acetonitrile		1ml	
DRE-XA16250000CY	Pirimicarb 100 µg/mL in Cyclohexane		1ml	
Pirimicarb D6 (dimethylcarbamate D6)				
CAS 1015854-66-6	MW 244.3232	C ₁₁ ² H ₁₆ H ₁₂ N ₄ O ₂		
DRE-C16250100	Pirimicarb D6 (dimethylcarbamate D6)		10mg	
DRE-XA16250100AC	Pirimicarb D6 100 µg/mL in Acetone		1.1ml	
Pirimicarb-desamido				
CAS 40778-16-3	MW 167.2083	C ₈ H ₁₃ N ₃ O		
DRE-LA16250800AL	Pirimicarb-desamido 10 µg/mL in Acetonitrile		1ml	
Pirimicarb-desamido-desmethyl				
CAS 78195-30-9	MW 153.1817	C ₇ H ₁₁ N ₃ O		
DRE-LA16251100AL	Pirimicarb-desamido-desmethyl 10 µg/mL in Acetonitrile		1ml	
Pirimicarb-desmethyl				
CAS 30614-22-3	MW 224.2596	C ₁₀ H ₁₆ N ₄ O ₂		
DRE-CA16251000	Pirimicarb-desmethyl(±)		10mg	
DRE-LA16251000AL	Pirimicarb-desmethyl 10 µg/mL in Acetonitrile		1ml	
Pirimicarb-desmethyl-formamido				
CAS 27218-04-8	MW 252.2697	C ₁₁ H ₁₆ N ₄ O ₃		
DRE-CA16251300	Pirimicarb-desmethyl-formamido		10mg	
DRE-LA16251300AL	Pirimicarb-desmethyl-formamido 10 µg/mL in Acetonitrile		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Pirimiphos-ethyl				
CAS 23505-41-1	MW 333.3867	C ₁₃ H ₂₄ N ₃ O ₃ PS		
DRE-C16260000	Pirimiphos-ethyl(‡)		100mg	
DRE-L16260000CY	Pirimiphos-ethyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA16260000CY	Pirimiphos-ethyl 100 µg/mL in Cyclohexane		1ml	
Pirimiphos-ethyl D10 (diethoxy D5)				
CAS n/a	MW 343.4483	C ₁₃ ² H ₁₀ H ₁₄ N ₃ O ₃ PS		
DRE-XA16260100AC	Pirimiphos-ethyl D10 100 µg/mL in Acetone		1.1ml	
Pirimiphos-ethyl-oxon				
CAS 36378-61-7	MW 317.3211	C ₁₃ H ₂₄ N ₃ O ₄ P		
DRE-XA16261000CY	Pirimiphos-ethyl-oxon 100 µg/mL in Cyclohexane		1ml	
Pirimiphos-methyl				
CAS 29232-93-7	MW 305.3336	C ₁₁ H ₂₀ N ₃ O ₃ PS		
DRE-C16270000	Pirimiphos-methyl(‡)		250mg	
DRE-L16270000CY	Pirimiphos-methyl 10 µg/mL in Cyclohexane		10ml	
DRE-XA16270000CY	Pirimiphos-methyl 100 µg/mL in Cyclohexane		1ml	
Pirimiphos-methyl D6 (dimethoxy D3)				
CAS 1793055-06-7	MW 311.3705	C ₁₁ ² H ₆ H ₁₄ N ₃ O ₃ PS		
DRE-XA16270100AC	Pirimiphos-methyl D6 100 µg/mL in Acetone		1.1ml	
Pirimiphos-methyl-N-desethyl				
CAS 67018-59-1	MW 277.2804	C ₉ H ₁₆ N ₃ O ₃ PS		
DRE-CA16270300	Pirimiphos-methyl-N-desethyl (*)		100mg	
Plifenate				
CAS 21757-82-4	MW 336.4264	C ₁₀ H ₇ Cl ₅ O ₂		
DRE-C16280000	Plifenate(‡)		25mg	
Prallethrin				
CAS 23031-36-9	MW 300.3921	C ₁₉ H ₂₄ O ₃		
DRE-C16286200	Prallethrin		100mg	
DRE-L16286200CY	Prallethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010224AL	Prallethrin 100 µg/mL in Acetonitrile(‡)		1ml	
Profenofos				
CAS 41198-08-7	MW 373.6308	C ₁₁ H ₁₅ BrClO ₃ PS		
DRE-C16330000	Profenofos(‡)		250mg	
DRE-L16330000CY	Profenofos 10 µg/mL in Cyclohexane		10ml	
DRE-XA16330000CY	Profenofos 100 µg/mL in Cyclohexane		1ml	
DRE-GA09010370ME	Profenofos 100 µg/mL in Methanol(‡)		1ml	








Pesticide and metabolites: Insecticides

Product code	Description			
Promecarb				
CAS 2631-37-0	MW 207.2689	C ₁₂ H ₁₇ NO ₂		
DRE-C16350000	Promecarb		100mg	
DRE-L16350000CY	Promecarb 10 µg/mL in Cyclohexane		10ml	
DRE-XA16350000CY	Promecarb 100 µg/mL in Cyclohexane		1ml	
DRE-XA09010227AL	Promecarb 100 µg/mL in Acetonitrile(‡)		1ml	
Propaphos				
CAS 7292-16-2	MW 304.3422	C ₁₃ H ₂₁ O ₄ PS		
DRE-C16420000	Propaphos		100mg	
DRE-L16420000IO	Propaphos 10 µg/mL in Isooctane		10ml	
Propargite				
CAS 2312-35-8	MW 350.4723	C ₁₉ H ₂₆ O ₄ S		
DRE-C16430000	Propargite(‡)		100mg	
DRE-L16430000AL	Propargite 10 µg/mL in Acetonitrile		10ml	
DRE-XA16430000IO	Propargite 100 µg/mL in Isooctane		1ml	
Propetamphos				
CAS 31218-83-4	MW 281.3089	C ₁₀ H ₂₀ NO ₄ PS		
DRE-C16460000	Propetamphos(‡)		100mg	
DRE-L16460000CY	Propetamphos 10 µg/mL in Cyclohexane		10ml	
DRE-XA16460000CY	Propetamphos 100 µg/mL in Cyclohexane		1ml	
Propoxur				
CAS 114-26-1	MW 209.2417	C ₁₁ H ₁₃ NO ₃		
DRE-C16500000	Propoxur(‡)		250mg	
DRE-L16500000CY	Propoxur 10 µg/mL in Cyclohexane		10ml	
Propoxur D3 (N-methyl D3)				
CAS 1219798-56-7	MW 212.2602	C ₁₁ ² H ₃ H ₁₂ NO ₃		
DRE-XA16500100AC	Propoxur D3 (methyl D3) 100 µg/mL in Acetone		1.1ml	
Prothiophos (Tokuthion®)				
CAS 34643-46-4	MW 345.2454	C ₁₁ H ₁₅ Cl ₂ O ₂ PS ₂		
DRE-C16560000	Prothiophos(‡)		50mg	
DRE-L16560000AL	Prothiophos 10 µg/mL in Acetonitrile		10ml	
DRE-L16560000CY	Prothiophos 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010229ME	Tokuthion® 100 µg/mL in Methanol(‡)		1ml	
Prothoate				
CAS 2275-18-5	MW 285.3638	C ₉ H ₂₀ NO ₃ PS ₂		
DRE-L16570000IO	Prothoate 10 µg/mL in Isooctane		10ml	
DRE-L16570000ME	Prothoate 10 µg/mL in Methanol		10ml	
DRE-XA16570000IO	Prothoate 100 µg/mL in Isooctane		1ml	
Pymetrozine				
CAS 123312-89-0	MW 217.2272	C ₁₀ H ₁₁ N ₅ O		
DRE-C16587000	Pymetrozine		100mg	
DRE-L16587000ME	Pymetrozine 10 µg/mL in Methanol		10ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Pyraclofos				
CAS 89784-60-1	MW 360.7961	C ₁₄ H ₁₈ ClN ₂ O ₃ PS		
DRE-C16592000	Pyraclofos(±)		50mg	
DRE-L16592000CY	Pyraclofos 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010231ME	Pyraclofos 100 µg/mL in Methanol(±)		1ml	
Pyrethrin 1				
CAS 121-21-1	MW 328.4452	C ₂₁ H ₂₈ O ₃		
DRE-XA16621000CY	Pyrethrin 1 100 µg/mL in Cyclohexane(±)		1ml	
Pyrethrins				
CAS 8003-34-7	MW 700.8999	C ₂₂ H ₂₈ O ₅ ·C ₂₁ H ₂₈ O ₃		
DRE-C16620000	Pyrethrins (technical)		100mg	
DRE-L16620000AL	Pyrethrins (technical) 10 µg/mL in Acetonitrile		10ml	
Pyridaben				
CAS 96489-71-3	MW 364.9326	C ₁₉ H ₂₅ ClN ₂ OS		
DRE-C16628000	Pyridaben(±)		25mg	
DRE-L16628000AL	Pyridaben 10 µg/mL in Acetonitrile		10ml	
DRE-L16628000CY	Pyridaben 10 µg/mL in Cyclohexane		10ml	
DRE-XA16628000AL	Pyridaben 100 µg/mL in Acetonitrile		1ml	
Pyridalyl				
CAS 179101-81-6	MW 491.1159	C ₁₈ H ₁₄ Cl ₄ F ₃ NO ₃		
DRE-C16629000	Pyridalyl(±)		100mg	
DRE-XA16629000CY	Pyridalyl 100 µg/mL in Cyclohexane		1.1ml	
Pyridaphenthion				
CAS 119-12-0	MW 340.3345	C ₁₄ H ₁₇ N ₂ O ₄ PS		
DRE-C16630000	Pyridaphenthion(±)		100mg	
DRE-L16630000CY	Pyridaphenthion 10 µg/mL in Cyclohexane		10ml	
DRE-XA16630000CY	Pyridaphenthion 100 µg/mL in Cyclohexane		1ml	
Pyrifluquinazon				
CAS 337458-27-2	MW 464.3368	C ₁₉ H ₁₅ F ₇ N ₄ O ₂		
DRE-XA16655700AL	Pyrifluquinazon 100 µg/mL in Acetonitrile (*) (±)		1ml	
Pyrimidifen				
CAS 105779-78-0	MW 377.9082	C ₂₀ H ₂₆ ClN ₃ O ₂		
DRE-C16659300	Pyrimidifen		10mg	
DRE-L16659300AL	Pyrimidifen 10 µg/mL in Acetonitrile		10ml	
Pyrimitate				
CAS 5221-49-8	MW 305.3336	C ₁₁ H ₂₀ N ₃ O ₃ PS		
DRE-LA16660000IO	Pyrimitate 10 µg/mL in Isooctane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Pyriproxyfen				
CAS 95737-68-1	MW 321.3698	C ₂₀ H ₁₉ NO ₃		
DRE-C16662500	Pyriproxyfen(±)		100mg	
DRE-L16662500AL	Pyriproxyfen 10 µg/mL in Acetonitrile		10ml	
DRE-L16662500CY	Pyriproxyfen 10 µg/mL in Cyclohexane		10ml	
DRE-XA16662500AL	Pyriproxyfen 100 µg/mL in Acetonitrile		1ml	
Quinalphos				
CAS 13593-03-8	MW 298.2979	C ₁₂ H ₁₅ N ₂ O ₃ PS		
DRE-C16700000	Quinalphos(±)		250mg	
DRE-L16700000CY	Quinalphos 10 µg/mL in Cyclohexane		10ml	
Quinthiophos				
CAS 1776-83-6	MW 329.3532	C ₁₇ H ₁₆ NO ₂ PS		
DRE-LA16720000EA	Quinthiophos 10 µg/mL in Ethylacetate(±)		1ml	
Resmethrin				
CAS 10453-86-8	MW 338.44	C ₂₂ H ₂₆ O ₃		
DRE-C16810000	Resmethrin		100mg	
DRE-L16810000CY	Resmethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010233AL	Resmethrin 100 µg/mL in Acetonitrile(±)		1ml	
Rotenone				
CAS 83-79-4	MW 394.4172	C ₂₃ H ₂₂ O ₆		
DRE-C16820000	Rotenone(±)		250mg	
Schradan				
CAS 152-16-9	MW 286.2487	C ₈ H ₂₄ N ₄ O ₃ P ₂		
DRE-C16910000	Schradan		10mg	
Silafluofen				
CAS 105024-66-6	MW 408.5805	C ₂₅ H ₂₉ FO ₂ Si		
DRE-C16946000	Silafluofen		100mg	
Spinetoram				
CAS 935545-74-7	MW n/a			
DRE-C16972770	Spinetoram(±)		100mg	No Structure
DRE-L16972770AL	Spinetoram 10 µg/mL in Acetonitrile (*)		10ml	
Spinosad				
CAS 168316-95-8	MW n/a			
DRE-C16972830	Spinosad(±)		100mg	No Structure
DRE-L16972830AL	Spinosad 10 µg/mL in Acetonitrile		10ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Spirodiclofen				
CAS 148477-71-8	MW 411.3189	C ₂₁ H ₂₄ Cl ₂ O ₄		
DRE-C16972950	Spirodiclofen		100mg	
DRE-L16972950CY	Spirodiclofen 10 µg/mL in Cyclohexane		10ml	
Spiromesifen				
CAS 283594-90-1	MW 370.4819	C ₂₃ H ₃₀ O ₄		
DRE-C16972970	Spiromesifen(±)		100mg	
DRE-L16972970CY	Spiromesifen 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010237AL	Spiromesifen 100 µg/mL in Acetonitrile(±)		1ml	
Spiromesifen-alcohol				
CAS 148476-30-6	MW 272.3389	C ₁₇ H ₂₀ O ₃		
DRE-C16972975	Spiromesifen-alcohol		10mg	
Spiromesifen-alcohol-4-hydroxymethyl				
CAS 873423-07-5	MW 288.3383	C ₁₇ H ₂₀ O ₄		
DRE-C16972978	Spiromesifen-alcohol-4-hydroxymethyl		10mg	
Spirotetramat				
CAS 203313-25-1	MW 373.4428	C ₂₁ H ₂₇ NO ₅		
DRE-C16972985	Spirotetramat(±)		100mg	
Spirotetramat-enol				
CAS 203312-38-3	MW 301.3801	C ₁₈ H ₂₃ NO ₃		
DRE-C16972990	Spirotetramat-enol(±)		10mg	
Spirotetramat-enol-glucoside				
CAS 1172614-86-6	MW 463.5207	C ₂₄ H ₃₃ NO ₈		
DRE-C16972993	Spirotetramat-enol-glucoside(±)		10mg	
DRE-LA16972993AL	Spirotetramat-enol-glucoside 10 µg/mL in Acetonitrile (*)		1ml	
Spirotetramat-keto-hydroxy				
CAS 1172134-11-0	MW 317.3795	C ₁₈ H ₂₃ NO ₄		
DRE-C16972996	Spirotetramat-keto-hydroxy(±)		10mg	
Spirotetramat-mono-hydroxy				
CAS 1172134-12-1	MW 303.396	C ₁₈ H ₂₅ NO ₃		
DRE-C16972998	Spirotetramat-mono-hydroxy(±)		10mg	
DRE-LA16972998AL	Spirotetramat-mono-hydroxy 10 µg/mL in Acetonitrile		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Strobane				
CAS 8001-50-1	MW 377.3495	C ₁₀ H ₉ Cl ₇		AND Enantiomer 
DRE-C16975000	Strobane		10mg	
DRE-L16975000CY	Strobane 10 µg/mL in Cyclohexane		10ml	
Sulcofuron-sodium Monohydrate				
CAS 1049723-02-5	MW 562.1831	C ₁₉ H ₁₁ Cl ₄ N ₂ O ₅ S·Na·H ₂ O		
DRE-C16987505	Sulcofuron-sodium monohydrate		100mg	
Sulfluramid				
CAS 4151-50-2	MW 527.198	C ₁₀ H ₆ F ₁₇ NO ₂ S		
DRE-C17004000	Sulfluramid(±)		100mg	
Sulfotep				
CAS 3689-24-5	MW 322.3189	C ₈ H ₂₀ O ₅ P ₂ S ₂		
DRE-C17010000	Sulfotep		100mg	
DRE-L17010000IO	Sulfotep 10 µg/mL in Isooctane		10ml	
DRE-XA09010238AL	Sulfotep 100 µg/mL in Acetonitrile(±)		1ml	
Sulfoxaflor				
CAS 946578-00-3	MW 277.2661	C ₁₀ H ₁₀ F ₃ N ₃ OS		
DRE-C17015000	Sulfoxaflor		10mg	
Sulprofos				
CAS 35400-43-2	MW 322.4468	C ₁₂ H ₁₉ O ₂ PS ₃		
DRE-C17030000	Sulprofos		100mg	
DRE-XA17030000CY	Sulprofos 100 µg/mL in Cyclohexane		1ml	
O,O-TEPP				
CAS 107-49-3	MW 290.1877	C ₈ H ₂₀ O ₇ P ₂		
DRE-C17240000	O,O-TEPP		50mg	
DRE-XA09010241AL	TEPP 100 µg/mL in Acetonitrile(±)		1ml	
Tebufenozide				
CAS 112410-23-8	MW 352.4699	C ₂₂ H ₂₈ N ₂ O ₂		
DRE-C17178800	Tebufenozide(±)		100mg	
DRE-L17178800AL	Tebufenozide 10 µg/mL in Acetonitrile		10ml	
DRE-XA17178800AC	Tebufenozide 100 µg/mL in Acetone		1ml	
Tebufenpyrad				
CAS 119168-77-3	MW 333.8557	C ₁₈ H ₂₄ ClN ₃ O		
DRE-C17179300	Tebufenpyrad(±)		100mg	
DRE-L17179300AL	Tebufenpyrad 10 µg/mL in Acetonitrile		10ml	
DRE-L17179300CY	Tebufenpyrad 10 µg/mL in Cyclohexane		10ml	
DRE-XA17179300CY	Tebufenpyrad 100 µg/mL in Cyclohexane		1ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Tebupirimfos				
CAS 96182-53-5	MW 318.3721	C ₁₃ H ₂₃ N ₂ O ₃ PS		
DRE-C17179700	Tebupirimfos		100mg	
DRE-L17179700CY	Tebupirimfos 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010240AL	Tebupirimfos 100 µg/mL in Acetonitrile(‡)		1ml	
Teflubenzuron				
CAS 83121-18-0	MW 381.1093	C ₁₄ H ₆ Cl ₂ F ₄ N ₂ O ₂		
DRE-C17210000	Teflubenzuron		250mg	
Tefluthrin				
CAS 79538-32-2	MW 418.7337	C ₁₇ H ₁₄ ClF ₇ O ₂		
DRE-C17213000	Tefluthrin(‡)		100mg	
DRE-L17213000AL	Tefluthrin 10 µg/mL in Acetonitrile		10ml	
DRE-L17213000CY	Tefluthrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA17213000CY	Tefluthrin 100 µg/mL in Cyclohexane		1ml	
Temephos				
CAS 3383-96-8	MW 466.4689	C ₁₆ H ₂₀ O ₆ P ₂ S ₃		
DRE-C17220000	Temephos		100mg	
DRE-L17220000CY	Temephos 10 µg/mL in Cyclohexane		10ml	
Temephos D12 (tetramethyl D12)				
CAS 1219795-39-7	MW 478.5429	C ₁₆ ² H ₁₂ H ₈ O ₆ P ₂ S ₃		
DRE-XA17220100CY	Temephos D12 (O,O,O',O'-tetramethyl D12) 100 µg/mL in Cylohexane		1.1ml	
Terbufos				
CAS 13071-79-9	MW 288.4306	C ₉ H ₂₁ O ₂ PS ₃		
DRE-C17270000	Terbufos		100mg	
DRE-L17270000AL	Terbufos 10 µg/mL in Acetonitrile		10ml	
DRE-L17270000CY	Terbufos 10 µg/mL in Cyclohexane		10ml	
DRE-LA09010259CY	Terbufos 10 µg/mL in Cyclohexane(‡)		10ml	
Terbufos-oxon				
CAS 56070-14-5	MW 272.365	C ₉ H ₂₁ O ₃ PS ₂		
DRE-LA17270100AL	Terbufos-oxon 10 µg/mL in Acetonitrile		1ml	
Terbufos-oxon-sulfone				
CAS 56070-15-6	MW 304.3638	C ₉ H ₂₁ O ₅ PS ₂		
DRE-LA17270200AL	Terbufos-oxon-sulfone 10 µg/mL in Acetonitrile		1ml	
DRE-LA17270200CY	Terbufos-oxon-sulfone 10 µg/mL in Cyclohexane (*)		1ml	
Terbufos-oxon-sulfoxide				
CAS 56165-57-2	MW 288.3644	C ₉ H ₂₁ O ₄ PS ₂		
DRE-LA17270300AL	Terbufos-oxon-sulfoxide 10 µg/mL in Acetonitrile		1ml	









Pesticide and metabolites: Insecticides

Product code	Description			
Terbufos-sulfone				
CAS 56070-16-7	MW 320.4294	C ₉ H ₂₁ O ₄ PS ₃		
DRE-C17270400	Terbufos-sulfone		50mg	
DRE-L17270400CY	Terbufos-sulfone 10 µg/mL in Cyclohexane		10ml	
Terbufos-sulfoxide				
CAS 10548-10-4	MW 304.43	C ₉ H ₂₁ O ₃ PS ₃		
DRE-C17270500	Terbufos-sulfoxide		50mg	
DRE-XA09010242AL	Terbufos-Sulfoxide 100 µg/mL in Acetonitrile(‡)		1ml	
Tetrachlorvinphos				
CAS 22248-79-9	MW 365.9618	C ₁₀ H ₆ Cl ₄ O ₄ P		
DRE-C17390000	Tetrachlorvinphos		250mg	
DRE-L17390000IO	Tetrachlorvinphos 10 µg/mL in Isooctane		10ml	
Tetradifon				
CAS 116-29-0	MW 356.0518	C ₁₂ H ₆ Cl ₄ O ₂ S		
DRE-C17400000	Tetradifon		250mg	
DRE-L17400000AL	Tetradifon 10 µg/mL in Acetonitrile		10ml	
DRE-L17400000CY	Tetradifon 10 µg/mL in Cyclohexane		10ml	
DRE-XA17400000CY	Tetradifon 100 µg/mL in Cyclohexane		1ml	
Tetramethrin				
CAS 7696-12-0	MW 331.4061	C ₁₉ H ₂₅ NO ₄		
DRE-C17410000	Tetramethrin		100mg	
DRE-L17410000CY	Tetramethrin 10 µg/mL in Cyclohexane		10ml	
DRE-XA09010244AL	Tetramethrin 100 µg/mL in Acetonitrile(‡)		1ml	
Tetrasul				
CAS 2227-13-6	MW 324.053	C ₁₂ H ₆ Cl ₄ S		
DRE-C17420000	Tetrasul		100mg	
Thiacloprid				
CAS 111988-49-9	MW 252.7233	C ₁₀ H ₉ ClN ₄ S		
DRE-C17451000	Thiacloprid(‡)		100mg	
DRE-L17451000AL	Thiacloprid 10 µg/mL in Acetonitrile		10ml	
Thiacloprid-amide				
CAS 676228-91-4	MW 270.7385	C ₁₀ H ₁₁ ClN ₅ O ₃ S		
DRE-C17451050	Thiacloprid-amide		10mg	
Thiamethoxam				
CAS 153719-23-4	MW 291.7147	C ₈ H ₁₀ ClN ₅ O ₃ S		
DRE-C17453000	Thiamethoxam(‡)		100mg	
DRE-GA17453000AL	Thiamethoxam 100 µg/mL in Acetonitrile		1ml	
DRE-GS17453000AL	Thiamethoxam 100 µg/mL in Acetonitrile		5x1ml	
DRE-L17453000AL	Thiamethoxam 10 µg/mL in Acetonitrile		10ml	










Pesticide and metabolites: Insecticides

Product code	Description			
Thiamethoxam D4 (oxadiazine D4)				
CAS 1331642-98-8 DRE-XA17453010AC	MW 295.7393 Thiamethoxam D4 (oxadiazine D4) 100 µg/mL in Acetone(±)	C ₈ H ₄ H ₆ ClN ₅ O ₃ S	1.1ml	
Thiamethoxam urea				
CAS 902493-06-5 DRE-C17453030	MW 247.7019 Thiamethoxam urea	C ₈ H ₁₀ ClN ₃ O ₂ S	10mg	
Thiocyclam hydrogenoxalate				
CAS 31895-22-4 DRE-C17480000	MW 271.3774 Thiocyclam hydrogenoxalate	C ₈ H ₁₁ NS ₃ ·C ₂ H ₂ O ₄	100mg	
Thiodicarb				
CAS 59669-26-0 DRE-C17490000 DRE-XA09010245AL	MW 354.4693 Thiodicarb Thiodicarb 100 µg/mL in Acetonitrile(±)	C ₁₀ H ₁₈ N ₄ O ₄ S ₃	100mg 1ml	
Thiodiphosphoric Acid Tetrapropyl Ester				
CAS 3244-90-4 DRE-CA10305000	MW 378.4252 Aspon	C ₁₂ H ₂₈ O ₅ P ₂ S ₂	10mg	
Thiofanox				
CAS 39196-18-4 DRE-C17500000 DRE-L17500000CY	MW 218.3164 Thiofanox Thiofanox 10 µg/mL in Cyclohexane	C ₉ H ₁₈ N ₂ O ₂ S	100mg 10ml	
Thiofanox D3 (N-methyl D3)				
CAS 74829-46-2 DRE-XA17500100AC	MW 221.3349 Thiofanox D3 (N-methyl D3) 100 µg/mL in Acetone	C ₉ ² H ₁₈ H ₁₉ N ₂ O ₂ S	1ml	
Thiofanox-sulfone				
CAS 39184-59-3 DRE-C17505000	MW 250.3152 Thiofanox-sulfone	C ₉ H ₁₈ N ₂ O ₄ S	100mg	
Thiofanox-sulfoxide				
CAS 39184-27-5 DRE-C17508000	MW 234.3158 Thiofanox-sulfoxide	C ₉ H ₁₈ N ₂ O ₃ S	100mg	







Pesticide and metabolites: Insecticides

Product code	Description			
Thiometon				
CAS 640-15-3 DRE-CA17520000	MW 246.3509 Thiometon (technical mixture)	C ₆ H ₁₅ O ₂ PS ₃	100mg	
Thiometon-sulfone				
CAS 20301-63-7 DRE-XA17521000CY	MW 278.3497 Thiometon-sulfone 100 µg/mL in Cyclohexane	C ₆ H ₁₅ O ₄ PS ₃	1ml	
Thiometon-sulfoxide				
CAS 2703-37-9 DRE-XA17522000CY	MW 262.3503 Thiometon-sulfoxide 100 µg/mL in Cyclohexane (*)	C ₆ H ₁₅ O ₃ PS ₃	1ml	
Thionazin				
CAS 297-97-2 DRE-C17530000 DRE-L17530000CY	MW 248.2392 Thionazin Thionazin 10 µg/mL in Cyclohexane	C ₈ H ₁₃ N ₂ O ₃ PS	50mg 10ml	
Thionazin D10 (diethyl D10)				
CAS 1173019-39-0 DRE-C17530100	MW 258.3008 Thionazin D10 (diethyl-d10)	C ₈ ² H ₁₀ H ₃ N ₂ O ₃ PS	10mg	
Thiosultap-sodium				
CAS 52207-48-4 DRE-C17560500	MW 355.3835 Thiosultap-sodium	C ₅ H ₁₁ NO ₆ S ₄ ·2Na	100mg	
Tolfenpyrad				
CAS 129558-76-5 DRE-C17591500 DRE-L17591500AL	MW 383.8713 Tolfenpyrad(±) Tolfenpyrad 10 µg/mL in Acetonitrile	C ₂₁ H ₂₂ ClN ₃ O ₂	100mg 10ml	
Toxaphene (Camphechlor)				
CAS 8001-35-2 DRE-C10940000 DRE-L10940000IO DRE-YA10940000HE DRE-GA09010300ME	MW n/a Camphechlor (Toxaphene) Camphechlor (Toxaphene) 10 µg/mL in Isooctane Camphechlor (Toxaphene) 4000 µg/mL in n-Hexane Toxaphene 500 µg/mL in Methanol(±)		250mg 10ml 1ml 1ml	No Structure
Tralomethrin				
CAS 66841-25-6 DRE-C17605500	MW 665.0072 Tralomethrin(±)	C ₂₂ H ₁₉ Br ₄ NO ₃	100mg	

Pesticide and metabolites: Insecticides

Product code	Description			
Transfluthrin				
CAS 118712-89-3	MW 371.1542	C ₁₅ H ₁₂ Cl ₂ F ₄ O ₂		
DRE-C17606000	Transfluthrin		250mg	
DRE-L17606000CY	Transfluthrin 10 µg/mL in Cyclohexane		10ml	
Triazamate				
CAS 112143-82-5	MW 314.4038	C ₁₃ H ₂₂ N ₄ O ₃ S		
DRE-L17648800CY	Triazamate 10 µg/mL in Cyclohexane		10ml	
DRE-XA17648800CY	Triazamate 100 µg/mL in Cyclohexane		1ml	
Triazophos				
CAS 24017-47-8	MW 313.3125	C ₁₂ H ₁₆ N ₃ O ₃ PS		
DRE-CA17650000	Triazophos(±)		100mg	
DRE-XA17650000AC	Triazophos 100 µg/mL in Acetone		1ml	
DRE-XA17650000CY	Triazophos 100 µg/mL in Cyclohexane (*)		1ml	
Tribufos				
CAS 78-48-8	MW 314.5109	C ₁₂ H ₂₇ OPS ₃		
DRE-CA17667000	Tribufos(±)		100mg	
DRE-XA17667000CY	Tribufos 100 µg/mL in Cyclohexane		1ml	
DRE-GA09010372ME	Tribufos 100 µg/mL in Methanol(±)		1ml	
DRE-XA09010179ME	Tribufos 100 µg/mL in Methanol(±)		1ml	
Trichloronate				
CAS 327-98-0	MW 333.5988	C ₁₀ H ₁₂ Cl ₃ O ₂ PS		
DRE-C17750000	Trichloronate		10mg	
DRE-L17750000CY	Trichloronate 10 µg/mL in Cyclohexane		10ml	
DRE-XA17750000CY	Trichloronate 100 µg/mL in Cyclohexane		1ml	
Trichloronitromethane				
CAS 76-06-2	MW 164.3752	CCl ₃ NO ₂		
DRE-CA11500000	Chloropicrin		100mg	
3,5,6-Trichloro-2-pyridinol				
CAS 6515-38-4	MW 198.4345	C ₅ H ₂ Cl ₃ NO		
DRE-LA17785000AL	3,5,6-Trichloro-2-pyridinol 10 µg/mL in Acetonitrile		1ml	
Trifenmorph				
CAS 1420-06-0	MW 329.4348	C ₂₃ H ₂₃ NO		
DRE-C17840000	Trifenmorph		10mg	
Triflumuron				
CAS 64628-44-0	MW 358.6997	C ₁₅ H ₁₀ ClF ₃ N ₂ O ₃		
DRE-C17844300	Triflumuron		100mg	

Pesticide and metabolites: Insecticides










Product code	Description			
2,3,5-Trimethacarb				
CAS 2655-15-4	MW 193.2423	C ₁₁ H ₁₅ NO ₂		
DRE-C17873500	2,3,5-Trimethacarb		10mg	
DRE-XA09010145ME	2,3,5-Trimethacarb 100 µg/mL in Methanol(‡)		1ml	
3,4,5-Trimethacarb				
CAS 2686-99-9	MW 193.2423	C ₁₁ H ₁₅ NO ₂		
DRE-C17874500	3,4,5-Trimethacarb		10mg	
DRE-L17874500CY	3,4,5-Trimethacarb 10 µg/mL in Cyclohexane		10ml	
Vamidothion				
CAS 2275-23-2	MW 287.3366	C ₈ H ₁₈ NO ₄ PS ₂		
DRE-C17900000	Vamidothion(‡)		100mg	
DRE-XA17900000AL	Vamidothion 100 µg/mL in Acetonitrile		1ml	
Vamidothion Sulfone				
CAS 70898-34-9	MW 319.3354	C ₈ H ₁₈ NO ₆ PS ₂		
DRE-XA17900007AL	Vamidothion-sulfon 100 µg/mL in Acetonitrile (*)		1ml	
Vamidothion Sulfoxide				
CAS 20300-00-9	MW 303.336	C ₈ H ₁₈ NO ₅ PS ₂		
DRE-C17900010	Vamidothion-sulfoxide		25mg	
Xylylcarb				
CAS 2425-10-7	MW 179.2157	C ₁₀ H ₁₃ NO ₂		
DRE-XA17946500AL	Xylylcarb 100 µg/mL in Acetonitrile		1ml	












OTHER PESTICIDES AND METABOLITES












Other pesticides and metabolites

Product code	Description			
ANTU				
CAS 86-88-4 DRE-C10290000	MW 202.2755 Antu	C ₁₁ H ₁₀ N ₂ S	250mg	
N-Acetylthiazolidine-4-carboxylic Acid				
CAS 5025-82-1 DRE-C10024100	MW 175.2056 N-Acetylthiazolidine-4-carboxylic acid	C ₆ H ₉ NO ₃ S	100mg	
Amidoflumet				
CAS 84466-05-7 DRE-C10161500	MW 317.6694 Amidoflumet	C ₉ H ₇ ClF ₃ NO ₄ S	25mg	
2-Amino-4,6-dimethoxypyrimidine (4,6-Dimethoxypyrimidin-2-ylamine)				
CAS 36315-01-2 DRE-C10201000	MW 155.1546 2-Amino-4,6-dimethoxypyrimidine	C ₆ H ₉ N ₃ O ₂	100mg	
3-Aminophenol				
CAS 591-27-5 DRE-C10211000	MW 109.1259 3-Aminophenol(‡)	C ₆ H ₇ NO	500mg	
4-Aminopyridine (Dalfampridine)				
CAS 504-24-5 DRE-C10222000	MW 94.1145 4-Aminopyridine	C ₅ H ₆ N ₂	500mg	
Ancymidol				
CAS 12771-68-5 DRE-C10250000	MW 256.2997 Ancymidol	C ₁₅ H ₁₆ N ₂ O ₂	100mg	
Anthranilic Acid Isopropylamide				
CAS 30391-89-0 DRE-C10280000	MW 178.231 Anthranilic acid isopropylamide	C ₁₀ H ₁₄ N ₂ O	100mg	
Anthraquinone				
CAS 84-65-1 DRE-C10281000 DRE-XA10281000AL	MW 208.2121 Anthraquinone(‡) Anthraquinone 100 µg/mL in Acetonitrile	C ₁₄ H ₈ O ₂	250mg 1ml	










Other pesticides and metabolites

Product code	Description			
Aramite				
CAS 140-57-8	MW 334.8587	C ₁₅ H ₂₃ ClO ₄ S		
DRE-C10300000	Aramite		10mg	
DRE-L10300000IO	Aramite 10 µg/mL in Isooctane		10ml	
DRE-GA09010391HE	Aramite Solution, 2000 µg/mL in Hexane(‡)		1ml	
Aviglycine Hydrochloride				
CAS 55720-26-8	MW 196.632	C ₆ H ₁₂ N ₂ O ₃ ·ClH		
DRE-C10334000	Aviglycine hydrochloride		5mg	
Azobenzene (1,2-Diphenyldiazene)				
CAS 103-33-3	MW 182.2212	C ₁₂ H ₁₀ N ₂		
DRE-C10390000	Azobenzene		500mg	
DRE-L10390000IO	Azobenzene 10 µg/mL in Isooctane		10ml	
Benoxacor				
CAS 98730-04-2	MW 260.1165	C ₁₁ H ₁₁ Cl ₂ NO ₂		
DRE-C10491000	Benoxacor(‡)		250mg	
DRE-L10491000AL	Benoxacor 10 µg/mL in Acetonitrile		10ml	
1,2-Benzisothiazol-3(2H)-one				
CAS 2634-33-5	MW 151.1857	C ₇ H ₅ NOS		
DRE-C10536600	1,2-Benzisothiazol-3(2H)-one(‡)		100mg	
6-Benzylaminopurine				
CAS 1214-39-7	MW 225.2492	C ₁₂ H ₁₁ N ₅		
DRE-C10569600	6-Benzylaminopurine(‡)		100mg	
Bifenazate-diazene				
CAS 149878-40-0	MW 298.3364	C ₁₇ H ₁₈ N ₂ O ₃		
DRE-C10579510	Bifenazate-diazene		25mg	
Bis-(2-chloro-1-methylethyl)ether				
CAS 108-60-1	MW 171.0649	C ₆ H ₁₂ Cl ₂ O		
DRE-C10651700	Bis-(2-chloro-1-methylethyl)ether		100mg	
DRE-YA10651700ME	Bis-(2-chloro-1-methylethyl)ether (DCIP) 2000 µg/mL in Methanol		1ml	
Brodifacoum				
CAS 56073-10-0	MW 523.4165	C ₃₁ H ₂₃ BrO ₃		
DRE-C10667500	Brodifacoum(‡)		100mg	
DRE-L10667500MB	Brodifacoum 10 µg/mL in Methyl-tert. butyl ether		10ml	










Other pesticides and metabolites

Product code	Description				
Bromadiolone					
CAS 28772-56-7 DRE-C10680000	MW 527.4052 Bromadiolone	$C_{30}H_{23}BrO_4$	100mg		
Bromethalin					
CAS 63333-35-7 DRE-XA10685000CY	MW 577.9303 Bromethalin 100 µg/mL in Cyclohexane	$C_{14}H_7Br_3F_3N_3O_4$	1ml		
4-Bromoaniline					
CAS 106-40-1 DRE-C10700000	MW 172.0225 4-Bromoaniline	C_6H_6BrN	500mg		
4-Bromo-2-chlorophenol					
CAS 3964-56-5 DRE-C10721500	MW 207.4524 4-Bromo-2-chlorophenol	C_6H_4BrClO	500mg		
Bronopol					
CAS 52-51-7 DRE-C10810000	MW 199.988 Bronopol	$C_3H_6BrNO_4$	250mg		
Buminafos					
CAS 51249-05-9 DRE-C10843000	MW 347.473 Buminafos	$C_{18}H_{38}NO_3P$	100mg		
Carvone					
CAS 99-49-0 DRE-C11052000	MW 150.2176 Carvone	$C_{10}H_{14}O$	250mg		
Cekafix					
CAS 121227-99-4 DRE-XA11066000CY	MW 407.2166 Cekafix 100 µg/mL in Cyclohexane	$C_{14}H_{16}BrO_5PS$	1ml		
Chloralose					
CAS 15879-93-3 DRE-C11100000	MW 309.5283 α-Chloralose	$C_8H_{11}Cl_3O_6$	250mg		










Other pesticides and metabolites

Product code	Description			
β-Chloralose				
CAS 16376-36-6 DRE-C11100200	MW 309.5283 β-Chloralose	C ₈ H ₁₁ Cl ₃ O ₆	50mg	
Chlormequat-chloride				
CAS 999-81-5 DRE-CA11340000 DRE-XA11340000AL	Chlormequat chloride Chlormequat-chloride 100 µg/mL in Acetonitrile	C ₅ H ₁₃ ClN·Cl	250mg 1ml	
Chlormequat-chloride 1,1,2,2-D4				
CAS n/a DRE-CA11340100 DRE-X11340100DO DRE-XA11340100DO	Chlormequat chloride 1,1,2,2-D4(‡) Chlormequat-chloride D4 100 µg/mL in Deuterium oxide Chlormequat-chloride D4 (CCC) 100 µg/mL in Deuterium oxide	C ₅ ² H ₄ H ₉ ClN·Cl	5mg 10ml 1.1ml	
N-Chloroacetyl-2,6-diethylaniline				
CAS 6967-29-9 DRE-C11349700	N-Chloroacetyl-2,6-diethylaniline	C ₁₂ H ₁₆ ClNO	100mg	
3-Chloroaniline				
CAS 108-42-9 DRE-C11351000 DRE-XA11351000ME	3-Chloroaniline 3-Chloroaniline 100 µg/mL in Methanol	C ₆ H ₆ ClN	500mg 1ml	
4-Chloroaniline				
CAS 106-47-8 DRE-C11352000 DRE-XA11352000AL	4-Chloroaniline(‡) 4-Chloroaniline 100 µg/mL in Acetonitrile	C ₆ H ₆ ClN	500mg 1ml	
6-Chloro-2-benzoxazolinone				
CAS 19932-84-4 DRE-C11392300	6-Chloro-2-benzoxazolinone	C ₇ H ₄ ClNO ₂	100mg	
2-Chloro-5-(chloromethyl)pyridine				
CAS 70258-18-3 DRE-C11397000	2-Chloro-5-(chloromethyl)pyridine	C ₆ H ₅ Cl ₂ N	100mg	
3-Chloro-5-hydroxy-6-phenylpyridazine (6-Chloro-4-hydroxy-3-phenyl-pyridazine)				
CAS 40020-01-7 DRE-C11417500	6-Chloro-4-hydroxy-3-phenyl-pyridazin(‡)	C ₁₀ H ₇ ClN ₂ O	10mg	










Other pesticides and metabolites

Product code	Description			
3-Chloro-4-methoxyaniline				
CAS 5345-54-0 DRE-C11420000	MW 157.5975 3-Chloro-4-methoxyaniline	C ₇ H ₈ ClNO	250mg	
3-Chloro-4-methylaniline				
CAS 95-74-9 DRE-C11429800	MW 141.5981 3-Chloro-4-methylaniline	C ₇ H ₈ ClN	500mg	
4-Chloro-2-methylaniline (4-Chloro-o-toluidine)				
CAS 95-69-2 DRE-C11430000 DRE-L11430000CY DRE-XA11430000ME	MW 141.5981 4-Chloro-2-methylaniline(‡) 4-Chloro-2-methylaniline 10 µg/mL in Cyclohexane 4-Chloro-2-methylaniline 100 µg/mL in Methanol	C ₇ H ₈ ClN	250mg 10ml 1ml	
5-Chloro-2-methyl-3(2H)-isothiazolone				
CAS 26172-55-4 DRE-CA11433000	MW 149.5987 5-Chloro-2-methyl-4-isothiazolin-3-one(‡)	C ₄ H ₄ ClNOS	25mg	
4-Chloro-2-methylphenol				
CAS 1570-64-5 DRE-C11440000 DRE-XA11440000AL	MW 142.5829 4-Chloro-2-methylphenol 4-Chloro-2-methylphenol 100 µg/mL in Acetonitrile	C ₇ H ₇ ClO	250mg 1ml	
6-Chloronicotinic Acid				
CAS 5326-23-8 DRE-C11452000	MW 157.5545 6-Chloronicotinic acid	C ₆ H ₄ ClNO ₂	100mg	
2-Chloro-4-nitrophenol				
CAS 619-08-9 DRE-C11454000	MW 173.5539 2-Chloro-4-nitrophenol	C ₆ H ₄ ClNO ₃	250mg	
Chlorophacinone				
CAS 3691-35-8 DRE-C11460000	MW 374.8164 Chlorophacinone	C ₂₃ H ₁₅ ClO ₃	100mg	
4-Chlorophenol				
CAS 106-48-9 DRE-C11472000 DRE-L11472000ME	MW 128.5563 4-Chlorophenol 4-Chlorophenol 10 µg/mL in Methanol	C ₆ H ₅ ClO	500mg 10ml	

Other pesticides and metabolites

Product code	Description			
(4-Chlorophenoxy)acetic Acid				
CAS 122-88-3 DRE-C11480000	MW 186.5924 4-Chlorophenoxyacetic acid	C ₈ H ₇ ClO ₃	250mg	
6-Chloro-2-picolinic Acid				
CAS 4684-94-0 DRE-C11495000	MW 157.5545 6-Chloro-2-picolinic acid	C ₆ H ₄ ClNO ₂	10mg	
6-Chloropyrimidine-2,4-diamine (2,6-Diamino-4-chloropyrimidine)				
CAS 156-83-2 DRE-C12193500	MW 144.5623 2,6-Diamino-4-chloropyrimidine	C ₄ H ₅ ClN ₄	100mg	
Cloprop				
CAS 101-10-0 DRE-C11688000	MW 200.619 Cloprop	C ₉ H ₉ ClO ₃	100mg	
Cloquintocet free acid ((5-Chloro-8-quinolinyloxy)acetic acid)				
CAS 88349-88-6 DRE-C11690500	MW 237.6391 Cloquintocet (free acid)	C ₁₁ H ₈ ClNO ₃	25mg	
Cloquintocet-1-methylhexyl ester (Cloquintocet-mexyl)				
CAS 99607-70-2 DRE-C11691000 DRE-L11691000AL DRE-XA09010171AC DRE-XA09010171ME	MW 335.8252 Cloquintocet-1-methylhexyl ester Cloquintocet-1-methylhexyl ester 10 µg/mL in Acetonitrile Cloquintocet-mexyl 100 µg/mL in Acetone(‡) Cloquintocet-mexyl 100 µg/mL in Methanol(‡)	C ₁₈ H ₂₂ ClNO ₃	100mg 10ml 1ml 1ml	
Cloxyfonac				
CAS 6386-63-6 DRE-C11692150	MW 216.6184 Cloxyfonac	C ₉ H ₉ ClO ₄	25mg	
Codlemone				
CAS 33956-49-9 DRE-C11692500	MW 182.3025 Codlemone	C ₁₂ H ₂₂ O	100mg	
Coumachlor				
CAS 81-82-3 DRE-C11710000	MW 342.773 Coumachlor	C ₁₉ H ₁₅ ClO ₄	100mg	










Other pesticides and metabolites

Product code	Description			
Coumafuryl				
CAS 117-52-2	MW 298.2901	C ₁₇ H ₁₄ O ₅		
DRE-C11720000	Coumafuryl		10mg	
DRE-XA11720000CY	Coumafuryl 100 µg/mL in Cyclohexane		1ml	
Coumatetralyl				
CAS 5836-29-3	MW 292.3285	C ₁₉ H ₁₆ O ₃		
DRE-C11740000	Coumatetralyl		250mg	
Crimidine				
CAS 535-89-7	MW 171.6274	C ₇ H ₁₀ ClN ₃		
DRE-C11750000	Crimidine		100mg	
DRE-L11750000AL	Crimidine 10 µg/mL in Acetonitrile		10ml	
Cuelure				
CAS 3572-06-3	MW 206.2378	C ₁₂ H ₁₄ O ₃		
DRE-C11773000	Cuelure		50mg	
Cyclanilide				
CAS 113136-77-9	MW 274.1001	C ₁₁ H ₉ Cl ₂ NO ₃		
DRE-C11817500	Cyclanilide		100mg	
Cyclohexanol-2-(4-tert-butyl-phenoxy)				
CAS 1942-71-8	MW 248.3606	C ₁₆ H ₂₄ O ₂		
DRE-C11825000	Cyclohexanol-2-(4-tert-butyl-phenoxy)		100mg	
Cycloheximide				
CAS 66-81-9	MW 281.3474	C ₁₅ H ₂₃ NO ₄		
DRE-C11830000	Cycloheximide		100mg	
Cyprosulfamide				
CAS 221667-31-8	MW 374.4109	C ₁₈ H ₁₈ N ₂ O ₅ S		
DRE-C11915000	Cyprosulfamide		100mg	
DEET D7 (methyl D3 phenyl D4)				
CAS 1219799-37-7	MW 198.3126	C ₁₂ ² H ₇ ² H ₁₀ NO		
DRE-XA12100010ME	DEET D7 (methyl D3 benzeneamide D4) 100 µg/mL in Methanol		1.1ml	










Other pesticides and metabolites

Product code	Description			
DMST				
CAS 66840-71-9	MW 214.2847	C ₉ H ₁₄ N ₂ O ₂ S		
DRE-C13040000	DMST(±)		50mg	
DRE-XA13040000AL	DMST 100 µg/mL in Acetonitrile		1ml	
Daminozide				
CAS 1596-84-5	MW 160.1711	C ₆ H ₁₂ N ₂ O ₃		
DRE-C11960000	Daminozide		250mg	
Daminozide D6 (dimethyl D6)				
CAS n/a	MW 166.2081	C ₆ ² H ₆ ² N ₂ O ₃		
DRE-C11960100	Daminozide D6 (dimethyl D6)		10mg	
DRE-XA11960100AL	Daminozide D6 100 µg/mL in Acetonitrile		1.1ml	
1-Decanol				
CAS 112-30-1	MW 158.2811	C ₁₀ H ₂₂ O		
DRE-C12095200	1-Decanol(±)		1ml	
4,4'-Dibromobenzophenone				
CAS 3988-03-2	MW 340.01	C ₁₃ H ₈ Br ₂ O		
DRE-C12220000	4,4'-Dibromobenzophenone		250mg	
Dichlorflurenol-methyl				
CAS 21634-96-8	MW 309.1441	C ₁₅ H ₁₀ Cl ₂ O ₃		
DRE-C12314020	Dichlorflurenol-methyl		10mg	
Dichlormid				
CAS 37764-25-3	MW 208.085	C ₈ H ₁₁ Cl ₂ NO		
DRE-C12315000	Dichlormid		100mg	
2,5-Dichloroaniline				
CAS 95-82-9	MW 162.0166	C ₆ H ₅ Cl ₂ N		
DRE-C12322500	2,5-Dichloroaniline		500mg	
DRE-XA12322500ME	2,5-Dichloroaniline 100 µg/mL in Methanol		1ml	
3,4-Dichloroaniline				
CAS 95-76-1	MW 162.0166	C ₆ H ₅ Cl ₂ N		
DRE-C12323400	3,4-Dichloroaniline		500mg	
DRE-XA12323400AL	3,4-Dichloroaniline 100 µg/mL in Acetonitrile		1ml	










Other pesticides and metabolites

Product code	Description			
3,5-Dichloroaniline				
CAS 626-43-7	MW 162.0166	C ₆ H ₅ Cl ₂ N		
DRE-C12323500	3,5-Dichloroaniline		500mg	
DRE-XA12323500ME	3,5-Dichloroaniline 100 µg/mL in Methanol		1ml	
2,6-Dichlorobenzamide				
CAS 2008-58-4	MW 190.0267	C ₇ H ₅ Cl ₂ NO		
DRE-C12355000	2,6-Dichlorobenzamide		500mg	
DRE-L12355000AL	2,6-Dichlorobenzamide 10 µg/mL in Acetonitrile		10ml	
DRE-XA12355000AL	2,6-Dichlorobenzamide 100 µg/mL in Acetonitrile		1ml	
2,6-Dichlorobenzoic Acid				
CAS 50-30-6	MW 191.0115	C ₇ H ₄ Cl ₂ O ₂		
DRE-C12402000	2,6-Dichlorobenzoic acid		250mg	
3,5-Dichlorobenzoic Acid				
CAS 51-36-5	MW 191.0115	C ₇ H ₄ Cl ₂ O ₂		
DRE-C12403100	3,5-Dichlorobenzoic acid		250mg	
2,4'-Dichlorobenzophenone				
CAS 85-29-0	MW 251.108	C ₁₃ H ₈ Cl ₂ O		
DRE-C12409000	2,4'-Dichlorobenzophenone		250mg	
4,4'-Dichlorobenzophenone				
CAS 90-98-2	MW 251.108	C ₁₃ H ₈ Cl ₂ O		
DRE-C12410000	4,4'-Dichlorobenzophenone(‡)		250mg	
DRE-L12410000CY	4,4'-Dichlorobenzophenone 10 µg/mL in Cyclohexane		10ml	
4,4'-Dichlorobenzophenone D8				
CAS 1219806-01-5	MW 259.1573	C ₁₃ ² H ₆ Cl ₂ O		
DRE-XA12410100AC	4,4'-Dichlorobenzophenone D8 100 µg/mL in Acetone		1.1ml	
2,4-Dichlorophenol 13C6				
CAS 1202864-83-2	MW 168.9573	¹³ C ₆ H ₄ Cl ₂ O		
DRE-XA12451200AC	2,4-Dichlorophenol 13C6 100 µg/mL in Acetone		1.1ml	
2,4-Dichlorophenol				
CAS 120-83-2	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12451000	2,4-Dichlorophenol(‡)		250mg	
DRE-L12451000ME	2,4-Dichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA12451000ME	2,4-Dichlorophenol 100 µg/mL in Methanol		1ml	
DRE-GA09010065ME	2,4-Dichlorophenol 100 µg/mL in Methanol(‡)		1ml	
DRE-GS09010065ME	2,4-Dichlorophenol 100 µg/mL in Methanol(‡)		5x1ml	










Other pesticides and metabolites

Product code	Description			
2,4-Dichlorophenol D3				
CAS 93951-74-7	MW 166.0198	C ₆ H ₃ HCl ₂ O		
DRE-C12451100	2,4-Dichlorophenol-3,5,6 D3		50mg	
DRE-XA12451100MB	2,4-Dichlorophenol-3,5,6 D3 100 µg/mL in Methyl-tert. butyl ether		1ml	
Dicyclopentadiene				
CAS 77-73-6	MW 132.2023	C ₁₀ H ₁₂		
DRE-C12587000	Dicyclopentadiene		250mg	
Diethyl Phosphorochloridothionate				
CAS 2524-04-1	MW 188.6128	C ₄ H ₁₀ ClO ₂ PS		
DRE-C12606700	Diethyl phosphorochloridothionate		100mg	
2-Diethylaminoethyl Hexanoate				
CAS 10369-83-2	MW 215.3324	C ₁₂ H ₂₅ NO ₂		
DRE-C14169100	2-Diethylaminoethyl Hexanoate		100mg	
2-Diethylamino-6-methyl-4-pyrimidinol				
CAS 42487-72-9	MW 181.2349	C ₈ H ₁₅ N ₃ O		
DRE-C12604750	2-Diethylamino-6-methyl-4-pyrimidinol		10mg	
2,6-Diethylaniline				
CAS 579-66-8	MW 149.2328	C ₁₀ H ₁₅ N		
DRE-C12605000	2,6-Diethylaniline		1g	
Diethyltoluamide (DEET)				
CAS 134-62-3	MW 191.2695	C ₁₂ H ₁₇ NO		
DRE-C12100000	DEET		250mg	
DRE-L12100000ME	DEET 10 µg/mL in Methanol		10ml	
DRE-XA12100000ME	DEET 100 µg/mL in Methanol		1ml	
Difenacoum				
CAS 56073-07-5	MW 444.5205	C ₃₁ H ₂₄ O ₃		
DRE-C12608000	Difenacoum(±)		10mg	
DRE-L12608000ME	Difenacoum 10 µg/mL in Methanol		10ml	
DRE-XA12608000ME	Difenacoum 100 µg/mL in Methanol		1ml	
Difethialone				
CAS 104653-34-1	MW 539.4821	C ₃₁ H ₂₃ BrO ₂ S		
DRE-L12625000AL	Difethialone 10 µg/mL in Acetonitrile (*)		10ml	










Other pesticides and metabolites

Product code	Description			
2,6-Difluorobenzoic Acid				
CAS 385-00-2 DRE-C12632010	MW 158.1023 2,6-Difluorobenzoic acid	C ₇ H ₄ F ₂ O ₂	250mg	
Diisopropylnaphthalene (mixture of isomers)				
CAS 38640-62-9 DRE-C12637690	MW 212.33 Diisopropylnaphthalene (mixture of isomers)	C ₁₆ H ₂₀	100mg	
2,6-Diisopropylnaphthalene				
CAS 24157-81-1 DRE-C12637700	MW 212.33 2,6-Diisopropylnaphthalene	C ₁₆ H ₂₀	100mg	
Dikegulac Monohydrate				
CAS 68539-16-2 DRE-C12639700	MW 292.2824 Dikegulac monohydrate	C ₁₂ H ₁₆ O ₇ ·H ₂ O	100mg	
Dikegulac Na				
CAS 52508-35-7 DRE-C12640000 DRE-L12640000AL	MW 296.249 Dikegulac Na Dikegulac Na 10 µg/mL in Acetonitrile	C ₁₂ H ₁₇ O ₇ ·Na	250mg 10ml	
2,4-Dimethylaniline				
CAS 95-68-1 DRE-C12724500 DRE-L12724500CY DRE-XA12724500AL	MW 121.1796 2,4-Dimethylaniline(‡) 2,4-Dimethylaniline 10 µg/mL in Cyclohexane 2,4-Dimethylaniline 100 µg/mL in Acetonitrile	C ₈ H ₁₁ N	500mg 10ml 1ml	
Dimethyldioctylammonium Chloride				
CAS 5538-94-3 DRE-C12726470	MW 305.9699 Dimethyldioctylammonium chloride	C ₁₈ H ₄₀ N·Cl	100mg	
Diphenadione (Diphacinone)				
CAS 82-66-6 DRE-C12878000	MW 340.3713 Diphacinone	C ₂₃ H ₁₆ O ₃	100mg	
Diphenyl Sulfone				
CAS 127-63-9 DRE-C12915000	MW 218.2716 Diphenyl sulfone	C ₁₂ H ₁₀ O ₂ S	250mg	










Other pesticides and metabolites

Product code	Description			
Dipropyl Pyridine-2,5-dicarboxylate (MGK 326)				
CAS 136-45-8 DRE-C15260000	MW 251.2784 MGK 326	C ₁₃ H ₁₇ NO ₄	100mg	
Ethephon				
CAS 16672-87-0 DRE-C13230000 DRE-L13230000AL	MW 144.494 Ethephon Ethephon 10 µg/mL in Acetonitrile	C ₂ H ₆ ClO ₃ P	250mg 10ml	
Ethephon D4 (2-Chloroethyl-1,1,2,2-D4)				
CAS 1020719-29-2 DRE-C13230100 DRE-XA13230100AC	MW 148.5186 Ethephon D4 (2-Chloroethyl-1,1,2,2-D4) Ethephon D4 (2-Chloroethyl-1,1,2,2-D4) 100 µg/mL in Acetone (*)	C ₂ ² H ₄ H ₂ ClO ₃ P	10mg 1ml	
Ethephon-hydroxy (2-Hydroxyethanephosphonic Acid)				
CAS 22987-21-9 DRE-CA13230200	MW 126.0483 Ethephon-hydroxy	C ₂ H ₇ O ₄ P	10mg	
Ethychlozate				
CAS 27512-72-7 DRE-C13315000	MW 238.6702 Ethychlozate(‡)	C ₁₁ H ₁₁ ClN ₂ O ₂	100mg	
Ethyl (1-Naphthyl)acetate (1-Naphthyl acetic acid ethyl ester)				
CAS 2122-70-5 DRE-C15465000	MW 214.2598 1-Naphthyl acetic acid-ethyl ester	C ₁₄ H ₁₄ O ₂	250mg	
Ethylene thiourea D4				
CAS 352431-28-8 DRE-C13330100 DRE-XA13330100AC	MW 106.1828 Ethylene thiourea D4 Ethylene thiourea D4 100 µg/mL in Acetone	C ₂ ² H ₄ H ₂ N ₂ S	50mg 1.1ml	
N-(2-Ethylhexyl)-5-norbornen-2,3-dicarboximide (MGK 264)				
CAS 113-48-4 DRE-C15250000	MW 275.3859 MGK 264	C ₁₇ H ₂₅ NO ₂	100mg	
2-Ethyl-6-methylaniline				
CAS 24549-06-2 DRE-C13348000	MW 135.2062 2-Ethyl-6-methylaniline	C ₉ H ₁₃ N	250mg	










Other pesticides and metabolites

Product code	Description			
2-Ethyl-6-methyl-2-chloroacetanilide				
CAS 32428-71-0 DRE-C13348010	MW 211.688 2-Ethyl-6-methyl-2-chloroacetanilide	C ₁₁ H ₁₄ ClNO	10mg	
Farnesol, mixture of isomers				
CAS 4602-84-0 DRE-C13405000	MW 222.3663 Farnesol, mixture of isomers	C ₁₅ H ₂₆ O	250mg	
Fenclorim				
CAS 3740-92-9 DRE-C13463000	MW 225.074 Fenclorim	C ₁₀ H ₆ Cl ₂ N ₂	100mg	
Fenridazon				
CAS 68254-10-4 DRE-C13552000	MW 264.6645 Fenridazon	C ₁₂ H ₉ ClN ₂ O ₃	100mg	
Fenridazon-propyl				
CAS 78778-15-1 DRE-C13553000	MW 306.7442 Fenridazon-propyl	C ₁₅ H ₁₅ ClN ₂ O ₃	25mg	
Flocoumafen				
CAS 90035-08-8 DRE-C13662000 DRE-LA13662000AL	MW 542.5444 Flocoumafen Flocoumafen 10 µg/mL in Acetonitrile	C ₃₃ H ₂₅ F ₃ O ₄	50mg 1ml	
Fluensulfone				
CAS 318290-98-1 DRE-C13710900	MW 291.6983 Fluensulfone	C ₇ H ₉ ClF ₃ NO ₂ S ₂	25mg	
Flumetralin				
CAS 62924-70-3 DRE-C13720000 DRE-L13720000CY	MW 421.7308 Flumetralin Flumetralin 10 µg/mL in Cyclohexane	C ₁₆ H ₁₂ ClF ₄ N ₃ O ₄	250mg 10ml	
Fluoroacetamide				
CAS 640-19-7 DRE-C13760000	MW 77.0577 Fluoroacetamide	C ₂ H ₄ FNO	100mg	










Other pesticides and metabolites

Product code	Description			
4-Fluoro-3-phenoxy benzoic acid				
CAS 77279-89-1 DRE-XA13797500AL	MW 232.2072 4-Fluoro-3-phenoxy benzoic acid 100 µg/mL in Acetonitrile	C ₁₃ H ₉ FO ₃	1ml	
Flurprimidol				
CAS 56425-91-3 DRE-C13851000 DRE-L13851000AL	MW 312.287 Flurprimidol Flurprimidol 10 µg/mL in Acetonitrile	C ₁₅ H ₁₅ F ₃ N ₂ O ₂	100mg 10ml	
Fluxofenim				
CAS 88485-37-4 DRE-C13880000	MW 309.6688 Fluxofenim	C ₁₂ H ₁₁ ClF ₃ NO ₃	100mg	
Forchlorfenuron				
CAS 68157-60-8 DRE-C13907000 DRE-L13907000ME	MW 247.6803 Forchlorfenuron(±) Forchlorfenuron 10 µg/mL in Methanol	C ₁₂ H ₁₀ ClN ₃ O	100mg 10ml	
Forchlorfenuron D5 (phenyl D5)				
CAS 1398065-87-6 DRE-XA13907010ME	MW 252.7111 Forchlorfenuron D5 (phenyl D5) 100 µg/mL in Methanol	C ₁₂ ² H ₁₀ H ₅ ClN ₃ O	1.1ml	
Furilazole				
CAS 121776-33-8 DRE-C13977000	MW 278.1318 Furilazole	C ₁₁ H ₁₃ Cl ₂ NO ₃	10mg	
Gibberellic Acid				
CAS 77-06-5 DRE-C14020000	MW 346.3744 Gibberellic acid(±)	C ₁₉ H ₂₂ O ₆	100mg	
Glyphosine				
CAS 2439-99-8 DRE-C14056000	MW 263.0796 Glyphosine	C ₄ H ₁₁ NO ₆ P ₂	100mg	
Icaridin (Picaridin)				
CAS 119515-38-7 DRE-C16195000	MW 229.3159 Picaridin	C ₁₂ H ₂₃ NO ₃	10mg	










Other pesticides and metabolites

Product code	Description			
Imicyafos				
CAS 140163-89-9 DRE-C14283650	MW 304.3488 Imicyafos	C ₁₁ H ₂₁ N ₄ O ₂ PS	25mg	
Imidazolidine-2-thione (Ethylene thiourea)				
CAS 96-45-7 DRE-C13330000 DRE-L13330000AL	MW 102.1581 Ethylene thiourea Ethylene thiourea 10 µg/mL in Acetonitrile	C ₃ H ₆ N ₂ S	250mg 10ml	
3-Indolylacetic Acid Ethyl Ester				
CAS 778-82-5 DRE-C14291000	MW 203.2371 3-Indolyl-acetic acid-ethyl ester	C ₁₂ H ₁₃ NO ₂	250mg	
3-Indolylacetic Acid				
CAS 87-51-4 DRE-C14290000	MW 175.184 3-Indolyl-acetic acid	C ₁₀ H ₉ NO ₂	250mg	
4-(3-Indolyl)butyric Acid				
CAS 133-32-4 DRE-C14310000	MW 203.2371 4-(3-Indolyl)-butyric acid	C ₁₂ H ₁₃ NO ₂	250mg	
3-(3-Indolyl)propionic Acid				
CAS 830-96-6 DRE-C14320000	MW 189.2105 3-(3-Indolyl)-propionic acid	C ₁₁ H ₁₁ NO ₂	250mg	
Irgarol 1051				
CAS 28159-98-0 DRE-C14374000 DRE-L14374000CY	MW 253.3671 Irgarol 1051 Irgarol 1051 10 µg/mL in Cyclohexane	C ₁₁ H ₁₉ N ₅ S	100mg 10ml	
2-Isopropoxyphenol				
CAS 4812-20-8 DRE-C14461000	MW 152.1904 2-Isopropoxyphenol	C ₉ H ₁₂ O ₂	250mg	
4-Isopropylaniline				
CAS 99-88-7 DRE-C14463000	MW 135.2062 4-Isopropylaniline	C ₉ H ₁₃ N	250mg	










Other pesticides and metabolites

Product code	Description			
2-Isopropyl-6-methyl-4-pyrimidinol				
CAS 2814-20-2 DRE-C14463800	MW 152.1937 2-Isopropyl-6-methyl-4-pyrimidinol	C ₈ H ₁₂ N ₂ O	250mg	
2-Isopropyl-naphthalene				
CAS 2027-17-0 DRE-GA09010402HD	MW 170.2503 2-Isopropyl-naphthalene Solution, 10000 µg/mL in Hexadecane(‡)	C ₁₃ H ₁₄	250ml	
Isoxadifen-ethyl				
CAS 163520-33-0 DRE-C14480500 DRE-L14480500AL DRE-L14480500CY	MW 295.3325 Isoxadifen-ethyl Isoxadifen-ethyl 10 µg/mL in Acetonitrile Isoxadifen-ethyl 10 µg/mL in Cyclohexane	C ₁₈ H ₁₇ NO ₃	100mg 10ml 10ml	
Lamprecide				
CAS 654-66-0 DRE-L14592000AL	MW 229.0886 Lamprecide 10 µg/mL in Acetonitrile	C ₇ H ₃ F ₃ NO ₃ ·Na	10ml	
Maleic Hydrazide				
CAS 123-33-1 DRE-C14730000	MW 112.0868 Maleic hydrazide	C ₄ H ₄ N ₂ O ₂	250mg	
Maleic Hydrazide D2				
CAS n/a DRE-C14730100	MW 114.0991 Maleic hydrazide D2	C ₄ ² H ₂ H ₂ N ₂ O ₂	10mg	
Mefenpyr-diethyl				
CAS 135590-91-9 DRE-C14860400	MW 373.2311 Mefenpyr-diethyl	C ₁₆ H ₁₈ Cl ₂ N ₂ O ₄	100mg	
Mepiquat Chloride (1,1-Dimethylpiperidinium chloride)				
CAS 24307-26-4 DRE-CA14880000 DRE-L14880000WA	MW 149.6616 Mepiquat chloride Mepiquat chloride 10 µg/mL in Water	C ₇ H ₁₆ N·Cl	100mg 10ml	
Mepiquat iodide D3 (methyl-d3)				
CAS 32317-85-4 DRE-X14880100DO DRE-XA14880100DO	MW 244.1316 Mepiquat Iodide D3 100 µg/ml in Deuterium oxide(‡) Mepiquat iodide D3 100 µg/mL in Deuterium oxide	C ₇ ² H ₉ H ₃ N·I	10ml 1.1ml	










Other pesticides and metabolites

Product code	Description			
Merphos (Tributylphosphoro-trithioite)				
CAS 150-50-5	MW 298.5115	C ₁₂ H ₂₇ PS ₃		
DRE-CA14910000	Merphos (*)		10mg	
DRE-XA09010208HE	Tributylphosphoro-trithioite (Merphos) 100 µg/mL in Hexane(‡)		1ml	
DRE-XA09010208ME	Tributylphosphoro-trithioite (Merphos) 100 µg/mL in Methanol(‡)		1ml	
Metaldehyde				
CAS 108-62-3	MW 176.2102	C ₈ H ₁₆ O ₄		
DRE-C14930000	Metaldehyde		250mg	
Methyl N-(3-Hydroxyphenyl)carbamate				
CAS 13683-89-1	MW 167.162	C ₈ H ₉ NO ₃		
DRE-L15088280CY	Methyl N-(3-hydroxyphenyl)-carbamate 10 µg/mL in Cyclohexane		10ml	
Methylene Dithiocyanate				
CAS 6317-18-6	MW 130.1914	C ₂ H ₂ N ₂ S ₂		
DRE-C15086030	Methylene dithiocyanate		100mg	
4-Methylimidazolidine-2-thione (Propylene thiourea)				
CAS 2122-19-2	MW 116.1847	C ₄ H ₈ N ₂ S		
DRE-C16530000	Propylene thiourea		250mg	
2-Methyl-4-isothiazolin-3-one Hydrochloride				
CAS 26172-54-3	MW 151.6145	C ₄ H ₅ NOS·ClH		
DRE-C15089050	2-Methyl-4-isothiazolin-3-on hydrochloride		100mg	
2-Methyl-4-isothiazolin-3-one				
CAS 2682-20-4	MW 115.1536	C ₄ H ₅ NOS		
DRE-C15089000	2-Methyl-4-isothiazolin-3-one(‡)		100mg	
3-Methyl-4-nitroimino-tetrahydro-1,3,5-oxadiazine				
CAS 153719-38-1	MW 160.1313	C ₆ H ₈ N ₄ O ₃		
DRE-XA15107500AC	3-Methyl-4-nitroimino-tetrahydro-1,3,5-oxadiazine 100 µg/mL in Acetone		1.1ml	
3-Methyl-4-nitrophenol				
CAS 2581-34-2	MW 153.1354	C ₇ H ₇ NO ₃		
DRE-C15110000	3-Methyl-4-nitrophenol		100mg	










Other pesticides and metabolites

Product code	Description			
3-Methylphenol (m-Cresol; Metacresol)				
CAS 108-39-4	MW 108.1378	C ₇ H ₈ O		
DRE-C15140300	3-Methylphenol (m-Cresol)(‡)		500mg	
DRE-XA15140300ME	3-Methylphenol 100 µg/mL in Methanol		1ml	
3-Methyl-quinoxaline-2-carboxylic acid				
CAS 74003-63-7	MW 188.1827	C ₁₀ H ₆ N ₂ O ₂		
DRE-C15143300	3-Methyl-quinoxaline-2-carboxylic acid(‡)		10mg	
2-(Methylthio)benzothiazole				
CAS 615-22-5	MW 181.2779	C ₈ H ₇ NS ₂		
DRE-C15144200	2-(Methylthio)benzothiazole		100mg	
Muscalure				
CAS 27519-02-4	MW 322.6113	C ₂₃ H ₄₆		
DRE-C15350000	Muscalure		100mg	
2-Naphthoxyacetic Acid				
CAS 120-23-0	MW 202.206	C ₁₂ H ₁₀ O ₃		
DRE-C15439000	2-Naphthoxy acetic acid		250mg	
1-Naphthylacetamide (1-Naphthaleneacetamide)				
CAS 86-86-2	MW 185.2218	C ₁₂ H ₁₁ NO		
DRE-C15450000	1-Naphthyl acetamide		500mg	
DRE-XA09010214ME	1-Naphthaleneacetamide 100 µg/mL in Methanol(‡)		1ml	
1-Naphthylacetic Acid Methyl Ester				
CAS 2876-78-0	MW 200.2332	C ₁₃ H ₁₂ O ₂		
DRE-C15470000	1-Naphthyl acetic acid methyl ester		500mg	
1-Naphthylacetic Acid				
CAS 86-87-3	MW 186.2066	C ₁₂ H ₁₀ O ₂		
DRE-C15460000	1-Naphthyl acetic acid(‡)		500mg	
Niclosamide				
CAS 50-65-7	MW 327.1196	C ₁₃ H ₈ Cl ₂ N ₂ O ₄		
DRE-C15510000	Niclosamid(‡)		100mg	










Other pesticides and metabolites

Product code	Description			
Nitrapyrin				
CAS 1929-82-4 DRE-C15550000	MW 230.9067 Nitrapyrin	C ₈ H ₅ Cl ₄ N	100mg	
4-Nitrophenol Sodium Salt				
CAS 824-78-2 DRE-C15594400	MW 161.0906 4-Nitrophenol sodium salt	C ₆ H ₄ NO ₃ ·Na	250mg	
4-Nitrophenol				
CAS 100-02-7 DRE-C15590400 DRE-L15590400AL	MW 139.1088 4-Nitrophenol 4-Nitrophenol 10 µg/mL in Acetonitrile	C ₆ H ₅ NO ₃	500mg 10ml	
4-Nitrophenol-2,3,5,6-D4				
CAS 93951-79-2 DRE-C15590404	MW 143.1334 4-Nitrophenol-2,3,5,6-D4	C ₆ ² H ₄ HNO ₃	100mg	
Oxabetrinil				
CAS 74782-23-3 DRE-C15755000	MW 232.2353 Oxabetrinil	C ₁₂ H ₁₂ N ₂ O ₃	100mg	
Paclobutrazol				
CAS 76738-62-0 DRE-C15840000 DRE-L15840000IO DRE-XA15840000AL	MW 293.7918 Paclobutrazol(‡) Paclobutrazol 10 µg/mL in Isooctane Paclobutrazol 100 µg/mL in Acetonitrile	C ₁₅ H ₂₀ ClN ₃ O	100mg 10ml 1ml	
Pentachloroaniline				
CAS 527-20-8 DRE-C15940000	MW 265.3518 Pentachloroaniline(‡)	C ₆ H ₂ Cl ₅ N	100mg	
Pentachlorobenzene				
CAS 608-93-5 DRE-C15960000	MW 250.3371 Pentachlorobenzene(‡)	C ₆ HCl ₅	500mg	
Pindone				
CAS 83-26-1 DRE-C16210000	MW 230.2592 Pindone	C ₁₄ H ₁₄ O ₃	250mg	










Other pesticides and metabolites

Product code	Description			
Prohexadione-calcium				
CAS 127277-53-6 DRE-C16345000	MW 462.4607 Prohexadione-calcium	$2C_{10}H_{11}O_5 \cdot Ca$	100mg	
Prohydrojasmon				
CAS 158474-72-7 DRE-C16345400 DRE-XA16345400CY	MW 254.3651 Prohydrojasmon Prohydrojasmon 100 µg/mL in Cyclohexane	$C_{15}H_{26}O_3$	100mg 1.1ml	
RH 5849				
CAS 112225-87-3 DRE-C16813000	MW 296.3636 RH 5849	$C_{18}H_{20}N_2O_2$	100mg	
S 421				
CAS 127-90-2 DRE-C16900000 DRE-L16900000CY	MW 377.7352 S 421(±) S 421 10 µg/mL in Cyclohexane	$C_8H_6Cl_6O$	250mg 10ml	
Strychnine				
CAS 57-24-9 DRE-C16980000 DRE-L16980000AL	MW 334.4116 Strychnine Strychnine 10 µg/mL in Acetonitrile	$C_{21}H_{22}N_2O_2$	250mg 10ml	
Sulfaquinoxaline				
CAS 59-40-5 DRE-C16990000	MW 300.3357 Sulfachinoxalin(±)	$C_{14}H_{12}N_4O_2S$	250mg	
2,3,5,6-Tetrachloroaniline				
CAS 3481-20-7 DRE-C17330600 DRE-L17330600CY	MW 230.9067 2,3,5,6-Tetrachloroaniline 2,3,5,6-Tetrachloroaniline 10 µg/mL in Cyclohexane	$C_6H_3Cl_4N$	100mg 10ml	
1,2,3,4-Tetrachlorobenzene				
CAS 634-66-2 DRE-C17353400 DRE-L17353400CY DRE-XA17353400ME	MW 215.8921 1,2,3,4-Tetrachlorobenzene 1,2,3,4-Tetrachlorobenzene 10 µg/mL in Cyclohexane 1,2,3,4-Tetrachlorobenzene 100 µg/mL in Methanol	$C_6H_2Cl_4$	100mg 10ml 1ml	
2,3,4,6-Tetrachlorophenol				
CAS 58-90-2 DRE-C17374600 DRE-L17374600CY	MW 231.8915 2,3,4,6-Tetrachlorophenol 2,3,4,6-Tetrachlorophenol 10 µg/mL in Cyclohexane	$C_6H_2Cl_4O$	50mg 10ml	





Other pesticides and metabolites

Product code	Description			
2,3,5,6-Tetrachlorophenol				
CAS 935-95-5	MW 231.8915	C ₆ H ₂ Cl ₄ O		
DRE-C17375600	2,3,5,6-Tetrachlorophenol(‡)		10mg	
DRE-L17375600CY	2,3,5,6-Tetrachlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA17375600CY	2,3,5,6-Tetrachlorophenol 100 µg/mL in Cyclohexane		1ml	
cis-1,2,3,6-Tetrahydrophthalimide				
CAS 1469-48-3	MW 151.1626	C ₈ H ₇ NO ₂		
DRE-C17406500	cis-1,2,3,6-Tetrahydrophthalimide		250mg	
Thidiazuron				
CAS 51707-55-2	MW 220.251	C ₉ H ₈ N ₄ OS		
DRE-C17465000	Thidiazuron		250mg	
Tralopyril				
CAS 122454-29-9	MW 349.5337	C ₁₂ H ₅ BrClF ₃ N ₂		
DRE-C17605700	Tralopyril		10mg	
1-Triacontanol (Melissyl alcohol)				
CAS 593-50-0	MW 438.8127	C ₃₀ H ₆₂ O		
DRE-XA17609200MB	1-Triacontanol (Melissyl alcohol) 100 µg/mL in Methyl-tert. butyl ether		1ml	
Triapenthenol				
CAS 76608-88-3	MW 263.3785	C ₁₅ H ₂₅ N ₃ O		
DRE-C17644000	Triapenthenol		100mg	
1H-1,2,4-Triazol				
CAS 288-88-0	MW 69.0653	C ₂ H ₃ N ₃		
DRE-C17649500	1 H-1,2,4-Triazole		250mg	
2,4,5-Trichlorophenol				
CAS 95-95-4	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17774500	2,4,5-Trichlorophenol(‡)		100mg	
DRE-L17774500AL	2,4,5-Trichlorophenol 10 µg/mL in Acetonitrile		10ml	
DRE-L17774500CY	2,4,5-Trichlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA17774500ME	2,4,5-Trichlorophenol 100 µg/mL in Methanol		1ml	
DRE-GA09010071ME	2,4,5-Trichlorophenol 100 µg/mL in Methanol(‡)		1ml	
3,5,6-Trichloro-2-pyridinol				
CAS 6515-38-4	MW 198.4345	C ₅ H ₂ Cl ₃ NO		
DRE-C17785000	3,5,6-Trichloro-2-pyridinol		10mg	

Other pesticides and metabolites

Product code	Description			
O,O,O-Triethylphosphorothioate				
CAS 126-68-1 DRE-C17837000	MW 198.2203 O,O,O-Triethylphosphorothioate	C ₆ H ₁₅ O ₃ PS	10mg	
4-(Trifluoromethoxy)aniline				
CAS 461-82-5 DRE-C17844810	MW 177.1238 4-(Trifluoromethoxy)aniline	C ₇ H ₆ F ₃ NO	100mg	
3-(Trifluoromethyl)aniline				
CAS 98-16-8 DRE-C17845000	MW 161.1244 3-Trifluoromethylaniline	C ₇ H ₆ F ₃ N	500mg	
O,O,S-Trimethyldithiophosphate				
CAS 2953-29-9 DRE-C17881500	MW 172.2061 O,O,S-Trimethyldithiophosphate(‡)	C ₃ H ₉ O ₂ PS ₂	100mg	
2,3,5-Trimethylphenol				
CAS 697-82-5 DRE-C17883500 DRE-XA17883500ME	MW 136.191 2,3,5-Trimethylphenol 2,3,5-Trimethylphenol 100 µg/mL in Methanol	C ₉ H ₁₂ O	250mg 1ml	
3,4,5-Trimethylphenol				
CAS 527-54-8 DRE-C17883800	MW 136.191 3,4,5-Trimethylphenol	C ₉ H ₁₂ O	100mg	
Trinexapac				
CAS 143294-89-7 DRE-L17888490AL	MW 224.21 Trinexapac (free acid) 10 µg/mL in Acetonitrile	C ₁₁ H ₁₂ O ₅	10ml	
Trinexapac-ethyl				
CAS 95266-40-3 DRE-C17888500	MW 252.2631 Trinexapac-ethyl	C ₁₃ H ₁₆ O ₅	100mg	
Trinexapac-methyl				
CAS 104273-71-4 DRE-C17888510	MW 238.2366 Trinexapac-methyl	C ₁₂ H ₁₄ O ₅	10mg	

Other pesticides and metabolites

Product code	Description			
Uniconazole				
CAS 83657-22-1	MW 291.7759	C ₁₅ H ₁₈ ClN ₃ O		
DRE-C17897000	Uniconazole(±)		10mg	
DRE-GA09010084AL	Uniconazole 1000 µg/mL in Acetonitrile(±)		1ml	
Verbenone				
CAS 80-57-9	MW 150.2176	C ₁₀ H ₁₄ O		
DRE-C17908000	Verbenone		100mg	
Warfarin				
CAS 81-81-2	MW 308.3279	C ₁₉ H ₁₆ O ₄		
DRE-C17940000	Warfarin(±)		250mg	
DRE-XA17940000AL	Warfarin 100 µg/mL in Acetonitrile		1ml	
(±)-Warfarin D5 (phenyl-D5)				
CAS 75472-93-4	MW 313.3587	C ₁₉ ² H ₅ H ₁₁ O ₄		
DRE-C17940100	(±)-Warfarin D5 (phenyl-D5)		10mg	
DRE-XA17940100AL	(±)-Warfarin D5 (phenyl-D5) 100 µg/mL in Acetonitrile		1.1ml	

PESTICIDE MIXTURES



Pesticide mixtures

Product code	Description	
Chlorinated Pesticide Mix 1		
DRE-YA05080100MB	Chlorinated Pesticide Mix 1 1000 µg/mL in Methyl-tert. butyl	1ml
4,4'-DDD alpha-Endosulfan delta-HCH Endrin-aldehyde Methoxychlor	4,4'-DDE alpha-HCH Dieldrin gamma-HCH	4,4'-DDT beta-Endosulfan Endosulfan-sulfate Heptachlor Aldrin beta-HCH Endrin Heptachlor-exo-epoxide (cis-isomer B)
Diquat and Paraquat Mix 1		
DRE-YA05490100WA	Diquat and Paraquat Mix 1 1000 µg/mL in Water	1ml
	Diquat dibromide hydrate Paraquat dichloride hydrate	
Organochlorine Pesticide Mix 1		
DRE-XA06080100TO	Organochlorine Pesticide Mix 1 in Toluene	1ml
	4,4'-DDD [600 µg/mL] 4,4'-DDT [600 µg/mL] alpha-Endosulfan [200 µg/mL] beta-HCH [100 µg/mL] Dieldrin [200 µg/mL] Endrin [200 µg/mL] gamma-HCH [100 µg/mL] Heptachlor-endo-epoxide (trans-isom. A) [100 µg/mL]	4,4'-DDE [200 µg/mL] Aldrin [100 µg/mL] alpha-HCH [100 µg/mL] delta-HCH [100 µg/mL] Endosulfan-sulfate [600 µg/mL] Endrin-aldehyde [600 µg/mL] Heptachlor [100 µg/mL]
Organochlorine Pesticide Mix 2		
DRE-YA08080100TH	Organochlorine Pesticide Mix 2 2000 µg/mL in Toluene/Hexane	1ml
4,4'-DDD alpha-Endosulfan delta-HCH Endrin-aldehyde Heptachlor-endo-epoxide (trans-isom. A)	4,4'-DDE alpha-HCH Dieldrin Endrin-ketone Methoxychlor	4,4'-DDT beta-Endosulfan Endosulfan-sulfate gamma-HCH Aldrin beta-HCH Endrin Heptachlor
Organochlorine Pesticide Mix 3		
DRE-LA06170100TH	Organochlorine Pesticide Mix 3 20 µg/mL in Toluene/Hexane	1ml
DRE-YA06170100TH	Organochlorine Pesticide Mix 3 1000 µg/mL in Toluene/Hexane	1ml
4,4'-DDD alpha-Endosulfan delta-HCH Endrin-aldehyde Methoxychlor	4,4'-DDE alpha-HCH Dieldrin gamma-HCH	4,4'-DDT beta-Endosulfan Endosulfan-sulfate Heptachlor Aldrin beta-HCH Endrin Heptachlor-endo-epoxide (trans-isom. A)
Organohalide Pesticide Mix 1		
DRE-XA05050100AC	Organohalide Pesticide Mix 1 100 µg/mL in Acetone	1ml
Alachlor Dieldrin Heptachlor-endo-epoxide (trans-isom. A) Methoxychlor	Aldrin Endrin Heptachlor-exo-epoxide (cis-isomer B) Simazine	Atrazine gamma-HCH Hexachlorobenzene trans-Chlordane (gamma) cis-Nonachlor Heptachlor Hexachlorocyclopentadiene trans-Nonachlor
Organohalide Pesticide Mix 2		
DRE-XA05050200ME	Organohalide Pesticide Mix 2 100 µg/mL in Methanol	1ml
	Alachlor Dieldrin gamma-HCH Heptachlor-endo-epoxide (trans-isom. A) Hexachlorobenzene	Aldrin Endrin Heptachlor Heptachlor-exo-epoxide (cis-isomer B) Methoxychlor

Pesticide mixtures

Product code	Description	
Organophosphorous Pesticide Mix 1		
DRE-XA06140100AC	Organophosphorous Pesticide Mix 1 200 µg/mL in Acetone	1ml
	Azinphos-methyl Diazinon Ethion Parathion-ethyl Demeton (O+S) Disulfoton Malathion Parathion-methyl	
Organophosphorous Pesticides Mix 4		
DRE-XA06220100CY	Organophosphorous Pesticides Mix 4 100 µg/mL in Cyclohexane	1ml
	Azinphos-methyl Demeton (O+S) Fensulfothion Phorate Sulprofos Coumaphos Disulfoton Fenthion Prothiophos Trichloronate	
Pesticide-Mix 2		
DRE-L18000002CY	Pesticide-Mix 2 10 µg/mL in Cyclohexane	10ml
	Diazinon Ethion Malathion Parathion-ethyl Parathion-methyl	
Pesticide-Mix 5		
DRE-L18000005CY	Pesticide Mix 5 10 µg/mL in Cyclohexane(±)	10ml
	2,4'-DDD 2,4'-DDT 4,4'-DDE Dieldrin 2,4'-DDE 4,4'-DDD 4,4'-DDT Endrin	
Pesticide-Mix 7		
DRE-L18000007CY	Pesticide Mix 7 10 µg/mL in Cyclohexane	10ml
	alpha-HCH beta-HCH delta-HCH gamma-HCH	
Pesticide-Mix 8		
DRE-LA18000008CY	Pesticide-Mix 8 10 µg/mL in Cyclohexane	1ml
	alpha-HCH beta-HCH delta-HCH epsilon-HCH gamma-HCH	
Pesticide-Mix 10		
DRE-Z18000010IO	Pesticide-Mix 10, 1 µg/mL in Isooctane	10ml
	2,4'-DDE 2,4'-DDE Dieldrin alpha-HCH gamma-HCH Heptachlor-endo-epoxide (trans-isom. A) Methoxychlor 4,4'-DDD 4,4'-DDE Endrin beta-HCH delta-HCH Hexachlorobenzene Mirex	
Pesticide-Mix 11		
DRE-L18000011CY	Pesticide Mix 11 10 µg/mL in Cyclohexane	10ml
	alpha-HCH beta-HCH delta-HCH gamma-HCH Hexachlorobenzene	

Pesticide mixtures

Product code	Description	
Pesticide-Mix 13		
DRE-L18000013CY	Pesticide Mix 13 10 µg/mL in Cyclohexane(‡)	10ml
DRE-L18000013TO	Pesticide Mix 13 10 µg/mL in Toluene	10ml
DRE-LA18000013CY	Pesticide-Mix 13 10 µg/mL in Cyclohexane(‡)	1ml
DRE-LA18000013TO	Pesticide-Mix 13 10 µg/mL in Toluene	1ml
2,4'-DDD	2,4'-DDE	2,4'-DDT
4,4'-DDE	4,4'-DDT	Aldrin
alpha-HCH	beta-Endosulfan	beta-HCH
delta-HCH	Dieldrin	Endrin
gamma-HCH	Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)
Hexachlorobenzene	Isodrin	Methoxychlor
oxy-Chlordane	PCB 101	PCB 138
PCB 180	PCB 28	PCB 52
		4,4'-DDD
		alpha-Endosulfan
		cis-Chlordane (alpha)
		epsilon-HCH
		Heptachlor-exo-epoxide (cis-isomer B)
		Mirex
		PCB 153
		trans-Chlordane (gamma)
Pesticide-Mix 14		
DRE-L18000014CY	Pesticide Mix 14 10 µg/mL in Cyclohexane(‡)	10ml
DRE-LA18000014CY	Pesticide-Mix 14 10 µg/mL in Cyclohexane	1ml
Aldrin	4,4'-DDD	4,4'-DDE
4,4'-DDT	Dieldrin	alpha-Endosulfan
Endrin	alpha-HCH	beta-HCH
Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)	Hexachlorobenzene
		2,4'-DDT
		beta-Endosulfan
		gamma-HCH
		Methoxychlor
Pesticide-Mix 15		
DRE-LA18000015AC	Pesticide-Mix 15 10 µg/mL in Acetone	1ml
DRE-X18000015AC	Pesticide-Mix 15 100 µg/mL in Acetone	10ml
	Atrazine	
	Atrazine-desethyl	
	Atrazine-desisopropyl	
	Propazine	
	Simazine	
	Terbutylazine	
Pesticide-Mix 17		
DRE-L18000017CY	Pesticide-Mix 17 10 µg/mL in Cyclohexane	10ml
1,2,4,5-Tetrachlorobenzene	1,2,4-Trichlorobenzene	2,4'-DDT
4,4'-DDE	4,4'-DDT	Aldrin
alpha-HCH	beta-Endosulfan	beta-HCH
Endrin	gamma-HCH	Heptachlor
Hexachlorobenzene	Methoxychlor	Pentachlorobenzene
		4,4'-DDD
		alpha-Endosulfan
		Dieldrin
		Heptachlor-endo-epoxide (trans-isom. A)
		Quintozene
Pesticide-Mix 18		
DRE-XA18000018AL	Pesticide-Mix 18 100 µg/mL in Acetonitrile	1ml
	Atrazine	
	Atrazine-desethyl	
	Atrazine-desisopropyl	
	Propazine	
	Sebutylazine	
	Simazine	
	Terbutylazine	
Pesticide-Mix 20		
DRE-L18000020CY	Pesticide-Mix 20 10 µg/mL in Cyclohexane	10ml
DRE-LA18000020CY	Pesticide-Mix 20 10 µg/mL in Cyclohexane	1ml
2,4'-DDD	2,4'-DDE	2,4'-DDT
4,4'-DDE	4,4'-DDT	Aldrin
alpha-HCH	beta-Endosulfan	beta-HCH
Endrin	gamma-HCH	Heptachlor
Hexachlorobenzene	Isobenzan	Isodrin
		4,4'-DDD
		alpha-Endosulfan
		Dieldrin
		Heptachlor-exo-epoxide (cis-isomer B)

Pesticide mixtures

Product code	Description	
Pesticide-Mix 31		
DRE-L18000031AC	Pesticide-Mix 31 10 µg/mL in Acetone	10ml
	Atrazine Atrazine-desethyl Metazachlor Metolachlor Simazine Terbuthylazine	
Pesticide-Mix 32		
DRE-L18000032CY	Pesticide-Mix 32 10 µg/mL in Cyclohexane	10ml
	2,4,5-T-methyl ester 2,4-D-methyl ester Dichlorprop-methyl ester MCPA-methyl ester MCPB-methyl ester Mecoprop-methyl ester	
Pesticide-Mix 33		
DRE-LA18000033IO	Pesticide-Mix 33 10 µg/mL in Isooctane	1ml
	4,4'-DDD alpha-Endosulfan Dieldrin Heptachlor-endo-epoxide (trans-isom. A) Isobenzan PCB 138 PCB 52	4,4'-DDE alpha-HCH Endrin Heptachlor-exo-epoxide (cis-isomer B) Isodrin PCB 153 Quintozone
	4,4'-DDT beta-HCH gamma-HCH Hexachloro-1,3-butadiene PCB 101 PCB 180	Aldrin Dichlobenil Heptachlor Hexachlorobenzene PCB 118 PCB 28
Pesticide-Mix 34		
DRE-L18000034AL	Pesticide-Mix 34 10 µg/mL in Acetonitrile	10ml
DRE-XA18000034AL	Pesticide-Mix 34 100 µg/mL in Acetonitrile	1ml
	Atrazine Chloroxuron Diuron Metamitron Metolachlor Prometryn Simazine	Atrazine-desethyl Chlorpropham Fenuron Metazachlor Metoxuron Propazine Terbuthylazine
	Atrazine-desethyl-desisopropyl Crimidine Isoproturon Methabenzthiazuron Metribuzin Propham Terbutryn	Chlorotoluron Cyanazine Linuron Metobromuron Monolinuron Sebuthylazine
Pesticide-Mix 40		
DRE-L18000040AL	Pesticide-Mix 40 10 µg/mL in Acetonitrile	10ml
DRE-LA18000040AL	Pesticide-Mix 40 10 µg/mL in Acetonitrile	1ml
	2,4-D 2,4-DB Dichlorprop MCPA MCPB Mecoprop	
Pesticide-Mix 41		
DRE-L18000041AL	Pesticide-Mix 41 10 µg/mL in Acetonitrile	10ml
	Atrazine Atrazine-desisopropyl Sebuthylazine Terbuthylazine	Atrazine-desethyl Propazine Simazine Terbuthylazine-desethyl
Pesticide-Mix 44		
DRE-L18000044AL	Pesticide-Mix 44 10 µg/mL in Acetonitrile	10ml
DRE-LA18000044AL	Pesticide-Mix 44 10 µg/mL in Acetonitrile	1ml
	Atrazine Diuron Metazachlor Metoxuron Terbuthylazine	Atrazine-desethyl Hexazinone Methabenzthiazuron Monolinuron
	Chlorotoluron Isoproturon Metobromuron Sebuthylazine	Cyanazine Linuron Metolachlor Simazine

Pesticide mixtures

Product code	Description	
Pesticide-Mix 51		
DRE-L18000051AL	Pesticide Mix 51 10 µg/mL in Acetonitrile	10ml
	Atrazine Atrazine-desisopropyl Prometryn Sebuthylazine Terbutylazine	Atrazine-desethyl Cyanazine Propazine Simazine Terbutryn
Pesticide-Mix 56		
DRE-L18000056AL	Pesticide-Mix 56 10 µg/mL in Acetonitrile	10ml
	Atrazine Atrazine-desisopropyl Diuron Propazine	Atrazine-desethyl Chlorotoluron Isoproturon Simazine
Pesticide-Mix 60		
DRE-LA18000060AC	Pesticide-Mix 60 10 µg/mL in Acetone	1ml
	Atrazine Atrazine-desisopropyl Sebuthylazine Simazine Terbutylazine-desethyl	Atrazine-desethyl Propazine Sebuthylazine-desethyl Terbutylazine
Pesticide-Mix 62		
DRE-L18000062CY	Pesticide-Mix 62 10 µg/mL in Cyclohexane	10ml
	Cypermethrin (technical) Deltamethrin Fenpropathrin Fenvalerate Permethrin Resmethrin	
Pesticide-Mix 64		
DRE-L18000064CY	Pesticide-Mix 64 10 µg/mL in Cyclohexane (*)	10ml
	Chlorpyrifos Diazinon Fenitrothion Methacrifos Pirimiphos-methyl	Chlorpyrifos-methyl Dichlorvos Malathion Phosphamidon
Pesticide-Mix 67		
DRE-L18000067EA	Pesticide-Mix 67 10 µg/mL in Ethylacetate (*)	10ml
	2,4,5-T 2,4-D Dichlorprop MCPA MCPB Mecoprop	
Pesticide-Mix 68		
DRE-XA18000068AC	Pesticide-Mix 68 100 µg/mL in Acetone	1ml
	Ametryn Prometryn Simazine Terbutryn	Atrazine Propazine Terbutylazine Trietazine
Pesticide-Mix 71		
DRE-L18000071CY	Pesticide-Mix 71 10 µg/mL in Cyclohexane(‡)	10ml
2,4'-DDD 4,4'-DDE alpha-HCH Dieldrin Heptachlor-endo-epoxide (trans-isom. A)	2,4'-DDE 4,4'-DDT beta-Endosulfan Endrin Heptachlor-exo-epoxide (cis-isomer B)	2,4'-DDT Aldrin beta-HCH gamma-HCH Hexachlorobenzene
		4,4'-DDD alpha-Endosulfan cis-Chlordane (alpha) Heptachlor Methoxychlor

Pesticide mixtures

Product code	Description			
Quintozene	Tecnazene	trans-Chlordane (gamma)		
Pesticide-Mix 95				
DRE-YA18000095EA	Pesticide-Mix 95 1000 µg/mL in Ethylacetate			1ml
	2,3,6-Trichlorobenzoic acid	2,4,5-T		
	2,4-D	2,4-DB		
	Bentazone	Bromoxynil		
	Dicamba	Dichlorprop		
	Ioxynil	MCPA		
	MCPB	Mecoprop		
	Triclopyr			
Pesticide-Mix 101				
DRE-LA18000101AL	Pesticide-Mix 101, 50 µg/mL in Acetonitrile (*)			1ml
Alachlor	Ametryn	Atrazine	Azinphos-ethyl	
Azinphos-methyl	Chlorpyrifos	Diazinon	Malathion	
Metolachlor	Molinate	Parathion-ethyl	Parathion-methyl	
Pendimethalin	Pirimicarb	Prometryn	Propazine	
Simazine	Terbuthylazine	Terbutryn	Trifluralin	
Pesticide-Mix 102				
DRE-LA18000102AL	Pesticide-Mix 102, 50 µg/mL in Acetonitrile			1ml
2,4'-DDT	4,4'-DDD	4,4'-DDT	Alachlor	
Aldrin	alpha-Endosulfan	alpha-HCH	beta-Endosulfan	
beta-HCH	Chlorpyrifos	Dieldrin	Endrin	
gamma-HCH	Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)	Hexachlorobenzene	
Metolachlor				
Pesticide-Mix 113				
DRE-X18000113AC	Pesticide-Mix 113 100 µg/mL in Acetone			10ml
	Ametryn			
	Atrazine			
	Prometryn			
	Propazine			
	Simazine			
	Terbuthylazine			
	Terbutryn			
Pesticide-Mix 114				
DRE-XA18000114IO	Pesticide-Mix 114 in Isooctane			1ml
	Bifenthrin [10 µg/mL]			
	Cyfluthrin [40 µg/mL]			
	Cypermethrin (technical) [40 µg/mL]			
	Deltamethrin [40 µg/mL]			
	Fenvalerate [30 µg/mL]			
	Permethrin [50 µg/mL]			
	Phenothrin [200 µg/mL]			
Pesticide-Mix 118				
DRE-L18000118CY	Pesticide-Mix 118 10 µg/mL in Cyclohexane(‡)			10ml
	Cyfluthrin			
	Cypermethrin (technical)			
	Deltamethrin			
	Fenvalerate			
	Permethrin			
	Phenothrin			
	Tetramethrin			
Pesticide-Mix 123				
DRE-LA18000123AL	Pesticide-Mix 123 in Acetonitrile			1ml
	Alachlor [20 µg/mL]	Ametryn [10 µg/mL]		
	Atrazine [10 µg/mL]	Atrazine-desethyl [10 µg/mL]		
	Cyanazine [10 µg/mL]	Pendimethalin [20 µg/mL]		
	Prometryn [10 µg/mL]	Propazine [10 µg/mL]		
	Simazine [10 µg/mL]	Terbuthylazine [10 µg/mL]		

(‡) ISO Guide 34

(*) Shorter expiry due to chemical nature of component(s)

Order online at www.lgcstandards.com

Pesticide mixtures

Product code	Description			
	Terbutryn [10 µg/mL]		Trifluralin [10 µg/mL]	
Pesticide-Mix 128				
DRE-L18000128AL	Pesticide-Mix 128 20 µg/mL in Acetonitrile			10ml
Clofibric acid	2,4,5-T	2,4-D	2,4-DB	
Bentazone	Bromacil	Clopyralid	Dicamba	
Dichlorprop	Fenoprop	Fluazifop (free acid)	Haloxypop (free acid)	
MCPA	MCPB	Mecoprop	Picloram	
Triclopyr				
Pesticide-Mix 129				
DRE-L18000129AL	Pesticide-Mix 129 in Acetonitrile			10ml
	Chlorotoluron [20 µg/mL]	Diuron [20 µg/mL]		
	Fenuron [20 µg/mL]	Fluometuron [20 µg/mL]		
	Isoproturon [20 µg/mL]	Linuron [20 µg/mL]		
	Metobromuron [20 µg/mL]	Metoxuron [40 µg/mL]		
	Monolinuron [20 µg/mL]	Monuron [20 µg/mL]		
Pesticide-Mix 134				
DRE-L18000134AL	Pesticide-Mix 134 10 µg/mL in Acetonitrile			10ml
	Chlorotoluron			
	Diuron			
	Isoproturon			
	Linuron			
Pesticide-Mix 157				
DRE-L18000157CY	Pesticide-Mix 157 10 µg/mL in Cyclohexane			10ml
Fenpropathrin	Fenvalerate	Folpet	HCH (BHC) (technical)	
Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)	Hexachlorobenzene	Iprodione	
Methoxychlor	Metolachlor	Nuarimol	Penconazole	
Pentachlorobenzene	Permethrin	Propachlor	Propiconazole	
tau-Fluvalinate	Tetradifon	Tetramethrin	Trifluralin	
Vinclozolin				
Pesticide-Mix 163				
DRE-YA18000163TH	Pesticide-Mix 163 2000 µg/mL in Toluene/Hexane			1ml
4,4'-DDD	4,4'-DDE	4,4'-DDT	Aldrin	
alpha-Endosulfan	alpha-HCH	beta-Endosulfan	beta-HCH	
delta-HCH	Dieldrin	Endosulfan-sulfate	Endrin	
gamma-HCH	Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)	Methoxychlor	
Pesticide-Mix 164				
DRE-L18000164CY	Pesticide-Mix 164 10 µg/mL in Cyclohexane			10ml
	2,4'-DDD			
	4,4'-DDD			
	2,4'-DDE			
	4,4'-DDE			
	2,4'-DDT			
	4,4'-DDT			
Pesticide-Mix 167				
DRE-L18000167TO	Pesticide-Mix 167 10 µg/mL in Toluene (*)			10ml
Azinphos-ethyl	Azinphos-methyl	Bromophos-ethyl	Bromophos-methyl	
Carbophenothion	Chlorfenvinphos	Chlorpyrifos	Chlorpyrifos-methyl	
Diazinon	Dichlofenthion	Dichlorvos	Ethion	
Etrinfos	Fenchlorphos	Fonofos	Malathion	
Methacrifos	Methidathion	Pirimiphos-ethyl	Pirimiphos-methyl	
Propetamphos	Sulfotep	Tetrachlorvinphos		

Pesticide mixtures

Product code	Description	
Pesticide-Mix 168		
DRE-L18000168TO	Pesticide-Mix 168 10 µg/mL in Toluene (*)	10ml
Dicrotophos Fenitrothion Methamidophos Parathion-ethyl Triazophos	Dimefox Formothion Mevinphos Parathion-methyl	Dimethoate Iodofenphos Omethoate Phosalone Disulfoton Malaoxon Paraoxon-ethyl Pyrazophos
Pesticide-Mix 172		
DRE-X18000172AL	Pesticide-Mix 172 100 µg/mL in Acetonitrile	10ml
	Ametryn Atrazine-desethyl Cyanazine Hexazinone Prometryn Sebuthylazine Terbutylazine	Atrazine Atrazine-desisopropyl Desmetryn Metribuzin Propazine Simazine Terbutryn
Pesticide-Mix 173		
DRE-X18000173EA	Pesticide-Mix 173 200 µg/mL in Ethylacetate (*)	10ml
	Acephate Chlorpyrifos Dichlorvos Fenchlorphos Formothion Paraoxon-ethyl Parathion-methyl Pirimiphos-methyl	Chlorfenvinphos Demeton-S-methyl-sulfone Dimethoate Fonofos Mevinphos Parathion-ethyl Phosalone
Pesticide-Mix 174		
DRE-X18000174EA	Pesticide-Mix 174 200 µg/mL in Ethylacetate(‡)	10ml
	Azinphos-methyl Chlorpyrifos-methyl Diazinon Fenitrothion Malathion Methidathion Phosphamidon	Bromophos-ethyl Demeton-S-methyl Ethion Malaoxon Methamidophos Paraoxon-methyl
Pesticide-Mix 192		
DRE-L18000192IO	Pesticide-Mix 192 Isooctane	10ml
	Bifenthrin [5 µg/mL] Cyfluthrin [10 µg/mL] Flucythrinate [10 µg/mL] lambda-Cyhalothrin [5 µg/mL] tau-Fluvalinate [10 µg/mL]	
Pesticide-Mix 195		
DRE-LA18000195IO	Pesticide-Mix 195 in Isooctane	1ml
	Cypermethrin (technical) [10 µg/mL] Deltamethrin [10 µg/mL] Fenpropathrin [5 µg/mL] Fenvalerate [10 µg/mL] Permethrin [10 µg/mL]	
Pesticide-Mix 235		
DRE-XA18000235AC	Pesticide-Mix 235 200 µg/mL in Acetone (*)	1ml
	Bromophos-ethyl Chlorfenvinphos Diazinon Dimethoate Ethion Malathion Paraoxon-ethyl Parathion-methyl	Bromophos-methyl Chlorpyrifos Dichlorvos Disulfoton Fenthion Mevinphos Parathion-ethyl

Pesticide mixtures

Product code	Description	
Pesticide-Mix 323		
DRE-XA18000323IO	Pesticide-Mix 323 in Isooctane(‡)	1ml
2,4'-DDD [20 µg/mL] 4,4'-DDD [20 µg/mL] alpha-Endosulfan [20 µg/mL] Chlorothalonil [40 µg/mL] Dieldrin [20 µg/mL] Heptachlor [20 µg/mL] trans-Chlordane (gamma) [20 µg/mL]	2,4'-DDE [20 µg/mL] 4,4'-DDE [20 µg/mL] alpha-HCH [20 µg/mL] cis-Chlordane (alpha) [20 µg/mL] Endosulfan-sulfate [20 µg/mL] Heptachlor-exo-epoxide [20 µg/mL]	2,4'-DDT [20 µg/mL] 4,4'-DDT [20 µg/mL] beta-Endosulfan [20 µg/mL] delta-HCH [20 µg/mL] Endrin [20 µg/mL] Hexachlorobenzene [20 µg/mL]
		2,4'-Dicofol [40 µg/mL] Aldrin [20 µg/mL] beta-HCH [20 µg/mL] Dicofol [20 µg/mL] gamma-HCH [20 µg/mL] oxy-Chlordane [20 µg/mL]
Pesticide-Mix 1037		
DRE-LA18001037CY	Pesticide-Mix 1037 10 µg/mL in Cyclohexane	1ml
2,4'-DDD 4,4'-DDE beta-HCH Dieldrin Heptachlor-endo-epoxide (trans-isom. A) trans-Chlordane (gamma)	2,4'-DDE 4,4'-DDT cis-Chlordane (alpha) Endrin Heptachlor-exo-epoxide (cis-isomer B) trans-Nonachlor	2,4'-DDT Aldrin cis-Nonachlor gamma-HCH Hexachlorobenzene
		4,4'-DDD alpha-HCH delta-HCH Heptachlor oxy-Chlordane
Pesticide-Mix 1111		
DRE-LA18001111CY	Pesticide-Mix 1111 10 µg/mL in Cyclohexane	1ml
	2,4'-DDT Aldrin Dicofol Endrin Heptachlor-endo-epoxide (trans-isom. A) Hexachlorobenzene trans-Chlordane (gamma)	4,4'-DDT cis-Chlordane (alpha) Dieldrin Heptachlor Heptachlor-exo-epoxide (cis-isomer B) oxy-Chlordane
Pesticide-Mix 1471		
DRE-LA18001471CY	Pesticide-Mix 1471 10 µg/mL in Cyclohexane	1ml
2,4'-DDT beta-HCH Endrin Heptachlor-exo-epoxide (cis-isomer B)	4,4'-DDT cis-Chlordane (alpha) gamma-HCH Hexachlorobenzene	Aldrin Dicofol Heptachlor oxy-Chlordane
		alpha-HCH Dieldrin Heptachlor-endo-epoxide (trans-isom. A) trans-Chlordane (gamma)
Pesticide-Mix 1584		
DRE-LA18001584AL	Pesticide-Mix 1584 10 µg/mL in Acetonitrile	1ml
Atrazine Cyanazine Isoproturon Metribuzin Propyzamide Triadimenol	Atrazine-desethyl Dimethoate Linuron Pirimicarb Simazine	Atrazine-desisopropyl Diuron Metamitron Prochloraz Terbutylazine
		Chloridazon Hexazinone Methabenzthiazuron Propiconazole Terbutylazine-desethyl
Pesticide-Mix 1598		
DRE-LA18001598CY	Pesticide-Mix 1598, 10 µg/mL in Cyclohexane	1ml
2,4'-DDT alpha-HCH Dicofol Heptachlor oxy-Chlordane	4,4'-DDT beta-Endosulfan Dieldrin Heptachlor-endo-epoxide (trans-isom. A) trans-Chlordane (gamma)	Aldrin beta-HCH Endrin Heptachlor-exo-epoxide (cis-isomer B)
		alpha-Endosulfan cis-Chlordane (alpha) gamma-HCH Hexachlorobenzene
Pesticide-Mix 1604		
DRE-LA18001604TO	Pesticide-Mix 1604, 10 µg/mL in Toluene	1ml
2,4'-DDD 4,4'-DDE beta-HCH gamma-HCH Hexachlorobenzene	2,4'-DDE 4,4'-DDT Chloridazon Heptachlor Tecnazene	2,4'-DDT Aldrin delta-HCH Heptachlor-endo-epoxide (trans-isom. A) Tetrachlorvinphos
		4,4'-DDD alpha-HCH Endrin Heptachlor-exo-epoxide (cis-isomer B)

Pesticide mixtures

Product code	Description	
Single Column Analytes Mix 4		
DRE-XA08080400TH	Single Column Analytes Mix 4 100 µg/mL in Toluene/Hexane	1ml
DRE-YA08080400TH	Single Column Analytes Mix 4 2000 µg/mL in Toluene/Hexane	1ml
4,4'-DDD alpha-Endosulfan cis-Chlordane (alpha) Endrin Heptachlor	4,4'-DDE alpha-HCH delta-HCH Endrin-aldehyde Heptachlor-endo-epoxide (trans-isom. A)	4,4'-DDT beta-Endosulfan Dieldrin Endrin-ketone trans-Chlordane (gamma)
		Aldrin beta-HCH Endosulfan-sulfate gamma-HCH Methoxychlor
Triazine Pesticide Mix 1		
DRE-XA06190100AC	Triazine Pesticide Mix 1 100 µg/mL in Acetone	1ml
	Ametryn Atrazine Prometryn Secbumeton Simetryn Terbutryn	Atraton Prometon Propazine Simazine Terbutylazine
Acrolein/Acrylonitrile Solution		
DRE-YA09000016WA	Acrolein/Acrylonitrile Mix 10000 µg/mL in Water(±)	1ml
	acrylonitrile acrolein	
Florida Pesticide Mix #2 25 Comp.		
DRE-GS09000249AL	Florida Pesticide Mix #2 25 Comp. 200 µg/mL in Acetonitrile(±)	1ea
acetamiprid asulam bromacil clothianidin diuron imidacloprid 3-hydroxycarbofuran	aldicarb azoxystrobin carbaryl bladex fenpyroximate (Mixture of Raceimers) Imidacloprid-olefin	aldicarb sulfone bifenazate carbofuran Diffubenzuron flupyradifurone iodosulfuron-methyl
		aldicarb sulfoxide boscalid chlorantraniliprole dinotefuran fluridone metalaxyl
Nevada Pesticide Mix		
DRE-GA09000062AL	Nevada Pesticide Mix, 22-1983 100 µg/mL in Acetonitrile(±)	1ml
DRE-GS09000062AL	Nevada Pesticide Mix, 22-1983 100 µg/mL in Acetonitrile(±)	5x1ml
acequinocyl baythroid (mixture of isomers) pentachloronitrobenzene cypermethrin (mix of isomers) fenhexamid fludioxonil	bifenazate dimethomorph spinosad (Spinosyn A & D) piperonyl butoxide flonicamid pyrethrin (mixture of isomers)	bifenthrin etoxazole thiamethoxam imidacloprid spinetoram (mixture of isomers)
		captan Systhane TM trifloxystrobin abamectin spirotetramat
Oregon Cannabis Pesticide Mix, 11-3485		
DRE-GA09000255AL	Oregon Cannabis Pesticide Mix, 11-3485 in Acetonitrile(±)	1ml
DRE-GS09000255AL	Oregon Cannabis Pesticide Mix, 11-3485 in Acetonitrile(±)	5x1ml
	acequinocyl [200 mg/L] baythroid (mixture of isomers) [100 mg/L] daminozide [100 mg/L] flonicamid [100 mg/L] oxamyl [200 mg/L] pyrethrin (mixture of isomers) [100 mg/L]	chlorfenapyr [100 mg/L] cypermethrin (mix of isomers) [100 mg/L] dichlorvos [100 mg/L] hexythiazox [100 mg/L] piperonyl butoxide [100 mg/L]
Oregon Cannabis Pesticide Mix, 16-3486		
DRE-GA09000256AL	Oregon Cannabis Pesticide Mix, 16-3486 in Acetonitrile(±)	1ml
DRE-GS09000256AL	Oregon Cannabis Pesticide Mix, 16-3486 in Acetonitrile(±)	5x1ml
acephate [40 mg/L] fenpyroximate [40 mg/L] dibrom [50 mg/L] paclobutrazol (mixt. isomers) [40 mg/L]	aldicarb [40 mg/L] kresoxim methyl [40 mg/L] propiconazol (mixt. isomers) [40 mg/L] fipronil [40 mg/L]	boscalid [40 mg/L] imidacloprid [40 mg/L] spiroxamine (mixture isomers) [40 mg/L] abamectin [50 mg/L]
		ethofenprox [40 mg/L] methomyl [40 mg/L] tebuconazol (Folicur) [40 mg/L] fludioxonil [40 mg/L]

Pesticide mixtures

Product code	Description	
Oregon Cannabis Pesticide Mix, 32-3487		
DRE-GA09000257AL	Oregon Cannabis Pesticide Mix, 32-3487 20 µg/mL in Acetonitrile(‡)	1ml
DRE-GS09000257AL	Oregon Cannabis Pesticide Mix, 32-3487 20 µg/mL in Acetonitrile(‡)	5x1ml
acetamiprid	azoxystrobin	bifenthrin
chlorpyrifos	diazinon	dimethoate
etoxazole	fenoxycarb	imazalil
metalaxyl	methiocarb	methyl parathion
Systhane TM	permethrin (mixture of isomers)	phosmet
pyridaben	spinosad (Spinosyn A & D)	spiromesifen
thiamethoxam	trifloxystrobin	spirotetramat
carbaryl	chlorantraniliprole	clofentezine
carbofuran		
ethoprophos (prophos)		
malathion		
MGK-264 - isomer b		
propoxur		
thiacloprid		
bifenazate		
prallethrin		
Oregon Cannabis Pesticide Mix, 59-2885		
DRE-GA09000244AL	Oregon Cannabis Pesticide Mix, 59-2885 in Acetonitrile(‡)	1ml
DRE-GS09000244AL	Oregon Cannabis Pesticide Mix, 59-2885 in Acetonitrile(‡)	5x1ml
abamectin [5 mg/L]	acequinocyl [20 mg/L]	aldicarb [4 mg/L]
daminozide [10 mg/L]	dichlorvos [10 mg/L]	ethofenprox [4 mg/L]
flonicamid [10 mg/L]	fludioxonil [4 mg/L]	hexythiazox [10 mg/L]
kresoxim methyl [4 mg/L]	methomyl [4 mg/L]	MGK-264 - isomer b [2 mg/L]
paclobutrazol (mixt. isomers) [4 mg/L]	piperonyl butoxide [20 mg/L]	pyrethrin (mixture of isomers) [10 mg/L]
tebuconazol (Folicur) [4 mg/L]	azoxystrobin [2 mg/L]	bifenthrin [2 mg/L]
permethrin (mixture of isomers) [2 mg/L]	phosmet [2 mg/L]	prallethrin [2 mg/L]
pyridaben [2 mg/L]	trifloxystrobin [2 mg/L]	acephate [4 mg/L]
diazinon [2 mg/L]	baythroid (mixture of isomers) [10 mg/L]	cypermethrin (mix of isomers) [10 mg/L]
malathion [2 mg/L]	methyl parathion [2 mg/L]	Systhane TM [2 mg/L]
spiromesifen [2 mg/L]	spirotetramat [2 mg/L]	thiacloprid [2 mg/L]
acetamiprid [2 mg/L]	bifenazate [2 mg/L]	boscalid [4 mg/L]
carbofuran [2 mg/L]	chlorantraniliprole [2 mg/L]	clofentezine [2 mg/L]
metalaxyl [2 mg/L]	methiocarb [2 mg/L]	dibrom [5 mg/L]
fenoxycarb [2 mg/L]	fenpyroximate [4 mg/L]	propoxur [2 mg/L]
chlorfenapyr [10 mg/L]		
fipronil [4 mg/L]		
imidacloprid [4 mg/L]		
oxamyl [10 mg/L]		
spiroxamine (mixture isomers) [4 mg/L]		
ethoprophos (prophos) [2 mg/L]		
propiconazol (mixt. isomers) [4 mg/L]		
chlorpyrifos [2 mg/L]		
dimethoate [2 mg/L]		
spinosad (Spinosyn A&D) [2 mg/L]		
thiamethoxam [2 mg/L]		
carbaryl [2 mg/L]		
imazalil [2 mg/L]		
etoxazole [2 mg/L]		
Oregon Cannabis Pesticide Mix, 59-2885, 10 µg/mL		
DRE-GS09000261AL	Oregon Cannabis Pesticide Mix, 59-2885 10 µg/mL in Acetonitrile(‡)	5x1ml
abamectin	acequinocyl	acetamiprid
bifenazate	boscalid	carbaryl
chlorantraniliprole	chlorfenapyr	clofentezine
cypermethrin (mix of isomers)	daminozide	dichlorvos
ethofenprox	etoxazole	fenoxycarb
fipronil	flonicamid	fludioxonil
imazalil	imidacloprid	kresoxim methyl
metalaxyl	methiocarb	methomyl
MGK-264 - isomer b	Systhane TM	dibrom
paclobutrazol (mixture of stereo isomers)	piperonyl butoxide	propoxur
spinosad (Spinosyn A & D)	spiromesifen	spirotetramat
tebuconazol (Folicur)	thiacloprid	thiamethoxam
bifenthrin	ethoprophos (prophos)	permethrin (mixture of isomers)
prallethrin	Propiconazol (mixture of isomers)	pyridaben
acephate	chlorpyrifos	diazinon
aldicarb		
carbofuran		
baythroid (mixture of isomers)		
dimethoate		
fenpyroximate		
hexythiazox		
malathion		
methyl parathion		
oxamyl		
pyrethrin (mixture of isomers)		
spiroxamine (mixture isomers)		
azoxystrobin		
phosmet		
trifloxystrobin		
Oregon Pest Mix 1		
DRE-GA09000231AL	Oregon Pesticide Mix #1 600 µg/mL in Acetonitrile(‡)	1ml
abamectin		
spinosad (Spinosyn A & D)		
Oregon Pesticide Mix #2		
DRE-GA09000232AL	Oregon Pesticide Mix #2 600 µg/mL in Acetonitrile(‡)	1ml
fenpyroximate	acequinocyl	acetamiprid
bifenazate	boscalid	chlorfenapyr
fludioxonil	imidacloprid	kresoxim methyl
MGK-264 - isomer a	piperonyl butoxide	spiromesifen
spiroxamine	trifloxystrobin	
azoxystrobin		
etoxazole		
metalaxyl		
spirotetramat		

Pesticide mixtures

Product code	Description	
Oregon Pesticide Mix #3, 10-3075		
DRE-GA09000233AL	Oregon Pesticide Mix #3, 10-3075 600 µg/mL in Acetonitrile(‡)	1ml
	aldicarb flonicamid methiocarb oxamyl thiacloprid	fipronil hexythiazox methomyl pyridaben thiamethoxam
Oregon Pesticide Mix #4, 12-3078		
DRE-GA09000234AL	Oregon Pesticide Mix #4, 12-3078 600 µg/mL in Acetonitrile(‡)	1ml
	carbaryl chlorantraniliprole daminozide Imazalil paclobutrazol (mixture of stereo isomers) propoxur	carbofuran clofentazine fenoxycarb (same as compd# 1871) Systhane TM Propiconazol (mixture of isomers) tebuconazol (Folicur)
Oregon Pesticide Mix #5, 7-3076		
DRE-GA09000235AL	Oregon Pesticide Mix #5, 7-3076 in Acetonitrile(‡)	1ml
	bifenthrin [600 mg/L] cypermethrin (mix of isomers) [600 mg/L] permethrin (mixture of isomers) [600 mg/L] pyrethrin I [350 mg/L]	baythroid (mixture of isomers) [600 mg/L] ethofenprox [600 mg/L] prallethrin [600 mg/L] pyrethrin II [250 mg/L]
Oregon Pesticide Mix 1		
DRE-GA09000058AL	Oregon Pesticide Mix 1, 20-1978 100 µg/mL in Acetonitrile(‡)	1ml
	abamectin aldicarb boscalid chlorfenapyr cypermethrin (mix of isomers)	acephate azoxystrobin carbaryl chlorpyrifos daminozide acequinocyl bifenazate carbofuran clofentazine dichlorvos acetamiprid bifenthrin chlorantraniliprole baythroid (mixture of isomers) diazinon
Oregon Pesticide Mix 2		
DRE-GA09000059AL	Oregon Pesticide Mix 2, 20-1979 100 µg/mL in Acetonitrile(‡)	1ml
	dimethoate fenpyroximate hexythiazox metalaxyl MGK-264 - isomer b	ethoprophos (prophos) fipronil imazalil methiocarb Systhane TM ethofenprox flonicamid imidacloprid methomyl malathion fenoxycarb fludioxonil kresoxim methyl methyl parathion etoxazole
Oregon Pesticide Mix 3		
DRE-GA09000060AL	Oregon Pesticide Mix 3, 19-1980 100 µg/mL in Acetonitrile(‡)	1ml
	phosmet spiromesifen thiamethoxam permethrin (mixture of isomers) Propiconazol (mixture of isomers)	piperonyl butoxide spirotetramat propoxur prallethrin trifloxystrobin pyridaben spiroxamine oxamyl pyrethrin (mixture of isomers) spinosad (Spinosyn A & D) tebuconazol (Folicur) paclobutrazol (mixture of stereo isomers) thiacloprid
Organochlorine Pesticide Mix		
DRE-GS09000057DI	Organochlorine Pesticide Mix 200 µg/mL in Dichloromethane(‡)	5x1ml
	a-BHC aldrin trans-chlordane dieldrin endosulfan I	b-BHC heptachlor p,p'-DDD endrin endosulfan II d-BHC heptachlor epoxide isomer B p,p'-DDE endrin aldehyde endosulfan sulfate g-BHC cis-chlordane p,p'-DDT endrin ketone methoxychlor
Pesticide Mix #6, 9-3079		
DRE-GA09000236AL	Pesticide Mix #6, 9-3079 600 µg/mL in Acetonitrile(‡)	1ml
	acephate diazinon ethoprophos (prophos) methyl parathion phosmet	chlorpyrifos dimethoate malathion dibrom dichlorvos

Pesticide mixtures

Product code	Description			
Pesticide Mix 12-3385				
DRE-GA09000300AL	Pesticide Mix 12-3385 100 µg/mL in Acetonitrile(‡)	10ml		
	alpha-cypermethrin bifenthrin danitol esfenvalerate fenvalerate (mixture of diastereoisomers) permethrin (mixture of isomers)	allethrin baythroid (mixture of isomers) deltamethrin ethofenprox lambda cyhalothrin phenothrin		
Pesticide Mix 15-3390				
DRE-GA09000301AL	Pesticide Mix 15-3390 100 µg/mL in Acetonitrile(‡)	10ml		
	thiophanate methyl pyraclostrobin acetamiprid azoxystrobin Imazalil methoxyfenozide piperonyl butoxide thiabendazole	3-hydroxycarbofuran abamectin amitraz chlorantraniliprole imidacloprid Systhane TM fenamidone		
Pesticide Mix 20-1124				
DRE-GS09000252AC	Pesticide Mix 20-1124 100 µg/mL in Acetone(‡)	5x1ml		
	dichlorvos prometon diazinon acetochlor prometryn	butylate propazine pyrimethanil chloropyrifos-methyl pirimiphos methyl	ethoprophos (prophos) propetamphos disulfoton metalaxyl metolachlor	trifluralin terbufos dimethenamide-P ametryne chlorpyrifos
Pesticide Mix A				
DRE-GS09000051DI	Pesticide Mix A in Dichloromethane(‡)	4x1ml		
	atrazine desisopropyl [1000 mg/L] methoxychlor [1000 mg/L] oxycarboxin [1000 mg/L] g-BHC [200 mg/L] metalaxyl [200 mg/L]	chlorothalonil [200 mg/L] metolachlor [200 mg/L] dieldrin [200 mg/L] diazinon [200 mg/L] parathion [200 mg/L]	bladex [1000 mg/L] simazine [200 mg/L] aldrin [200 mg/L] endosulfan I [200 mg/L] triallate [200 mg/L]	desethyl atrazine [1000 mg/L] aldicarb [4000 mg/L] a-BHC [200 mg/L] ethofumesate [200 mg/L]
Pesticide Mix B				
DRE-GS09000052DI	Pesticide Mix B in Dichloromethane(‡)	4x1ml		
	methomyl [1000 mg/L] chlorpyrifos [50 mg/L] Ethalfuralin [50 mg/L] fenoxaprop-P-ethyl [250 mg/L] terbufos [250 mg/L]	Guthion® [1000 mg/L] dimethoate [50 mg/L] ethion [1000 mg/L] linuron [250 mg/L] napropamide [250 mg/L]	bromacil [250 mg/L] disulfoton [1000 mg/L] phorate [50 mg/L] malathion [250 mg/L] vinclozolin [50 mg/L]	carboxine [1000 mg/L] diuron [1000 mg/L] thiamethoxam [250 mg/L] pyridaben [250 mg/L]
Pesticide Mix C				
DRE-GS09000053DI	Pesticide Mix C in MeCl2:Acetone:MeOH:MtBE 5:2:2:1(‡)	4x1ml		
	brominal [200 mg/L] 2,4-dichlorophenol [200 mg/L] imazethapyr [1000 mg/L] picloram [200 mg/L] imazamox [200 mg/L] 2,4-DB [200 mg/L]	4-chloro-2-methylphenol [200 mg/L] dichlorprop (2,4-DP) [200 mg/L] MCPA acid [200 mg/L] Triclopyr [200 mg/L] quinclorac [200 mg/L]	Lontrel (clopyralid) [1000 mg/L] diclofop-methyl [1000 mg/L] MCPB [1000 mg/L] Hexaconazole [200 mg/L] trifluralin [200 mg/L]	Dicamba [200 mg/L] imazamethabenz-methyl [1000 mg/L] MCPP acid [200 mg/L] fluroxypyr [200 mg/L] 2,4-D [200 mg/L]
Pesticide Mix D				
DRE-GS09000054DI	Pesticide Mix D in Dichloromethane(‡)	4x1ml		
	iprodione [1000 mg/L] clodinafop (acid free) [1000 mg/L] aminopyralid [1000 mg/L] quizalofop [1000 mg/L]	clodinafop-propargyl [1000 mg/L] fluzifop [200 mg/L] Propiconazol (mixture of isomers) [1000 mg/L] bentazon [200 mg/L]		

Pesticide mixtures

Product code	Description	
Pesticide Mix E		
DRE-GS09000139TO	Pesticide Mix E 100 µg/mL in Toluene(±)	5x1ml
acetochlor anilazine b-BHC bifenthrin bromopropylate	alachlor atrazine d-BHC Bioresmethrin butachlor	aldrin Benfluralin (Benefin) g-BHC bitertanole (mixture of isomers) butralin
		ametryne a-BHC bifenox boscalid butylate
Pesticide Mix F		
DRE-GS09000140TO	Pesticide Mix F 100 µg/mL in Toluene(±)	5x1ml
captan cis-chlordane chlorobenzilate clodinafop-propargyl gamma-cyhalothrin	carbophenothion oxychlordane chlorothalonil Command lambda cyhalothrin	carboxine trans-chlordane chlorpyrifos cloquintocet-mexyl cypermethrin (mix of isomers)
		carfentrazone-ethyl chlorfenvinphos (E/Z-mixture) chlorpyrifos-methyl baythroid (mixture of isomers) cyprodinil
Pesticide Mix G		
DRE-GS09000141TO	Pesticide Mix G 100 µg/mL in Toluene(±)	5x1ml
chlorthal-dimethyl (dacthal) o,p'-DDT diclofop-methyl dimethenamid diphenamid	o,p'-DDD p,p'-DDT Dichloran dimethenamide-P disulfoton	p,p'-DDD deltamethrin dieldrin dimethomorph endosulfan I
		p,p'-DDE diazinon Difenoconazole (mixture of isomers) dinitramine endosulfan II
Pesticide Mix H		
DRE-GS09000142TO	Pesticide Mix H 100 µg/mL in Toluene(±)	5x1ml
endosulfan sulfate esfenvalerate famoxadone fenitrothion fenvalerate (mixture of diastereoisomers)	endrin Ethfluralin fenamidone fenoxaprop-ethyl fipronil	EPN ethoprophos (prophos) fenamiphos fensulfothion fluzifop-butyl
		EPTC (s-ethyl dipropylthiocarbamate) ethofenprox fenbuconazole fenthion flucythrinate
Pesticide Mix I		
DRE-GS09000143AL	Pesticide Mix I 100 µg/mL in Acetonitrile(±)	5x1ml
flufenacet Folpet Hexaconazole isoxadifen-ethyl metalaxyl-m	flusilazole Flurilazole imazamethabenz-methyl kresoxim methyl mefenpyr-diethyl	flutolanil heptachlor indoxacarb lactofen metalaxyl
		flutriafol heptachlor epoxide isomer B ipconazole malathion metconazole
Pesticide Mix J		
DRE-GS09000144TO	Pesticide Mix J 100 µg/mL in Toluene(±)	5x1ml
methidathion metribuzin oxadiazon methyl parathion phenthoate	methoprene (mixture of isomers) Systhane TM Oxyfluorfen pendimethalin phorate	methoxychlor nitrapyrin paclobutrazol (mixture of stereo isomers) permethrin phosmet
		metolachlor norflurazon parathion phenothrin phosphamidon
Pesticide Mix K		
DRE-GS09000145TO	Pesticide Mix K 100 µg/mL in Toluene(±)	5x1ml
Pirimicarb propargite pyrimethanil simazine terbufos	procymidone Propiconazol (mixture of isomers) quinoxifen spiromesifen terbutryne	profenofos pyraflufen-ethyl pentachloronitrobenzene tebuconazol thiometon
		propachlor pyrazophos quizalofop-ethyl tefluthrin tolclofos-methyl
Pesticide Mix L		
DRE-GS09000146TO	Pesticide Mix L 100 µg/mL in Toluene(±)	5x1ml
	triadimefon triadimenol (baitan) triallate triazophos trifluralin	

Pesticide mixtures

Product code	Description			
	vinclozolin zoxamide			
Pesticide Mix M				
DRE-GS09000147AL	Pesticide Mix M 100 µg/mL in Acetonitrile(‡)			5x1ml
3-hydroxycarbofuran alanycarb aminocarb benfuracarb buprofezin	acephate aldicarb Guthion® bensulide butocarboxim	acetamiprid aldicarb sulfone azoxystrobin bentazon butocarboxim sulfoxide	acifluorfen aldicarb sulfoxide bendiocarb brominal butoxycarboxim	
Pesticide Mix N				
DRE-GS09000148AL	Pesticide Mix N 100 µg/mL in Acetonitrile(‡)			5x1ml
carbaryl chlorimuron ethyl clothianidin cyromazine dimepiperate	carbofuran chlorpropham crotoxyphos dichlofluanid dimethoate	carbosulfan chlorosulfuron cycloxdim dichlorvos dinotefuran	chlorfluazuron Cloransulam-methyl cyproconazole (diastereomers) Diflubenzuron dioxacarb	
Pesticide Mix O				
DRE-GS09000149AL	Pesticide Mix O 100 µg/mL in Acetonitrile(‡)			5x1ml
diuron ethion fenpropimorph 2 flumestulam fomesafen	esprocarb fenhexamid florasulam flumioxazin furathiocarb	ethiofencarb fenobucarb flucarbazon-sodium fluometuron halosulfuron-methyl	ethiofencarb-sulfoxide fenoxycarb fludioxonil fluthiacet-methyl haloxyfop	
Pesticide Mix P				
DRE-GS09000150AL	Pesticide Mix P 100 µg/mL in Acetonitrile(‡)			5x1ml
haloxyfop-2-ethoxyethyl flumiclorac-pentyl iprovalicarb methamidophos methoxyfenozide	haloxyfop-methyl imidacloprid isoprocarb methiocarb metobromuron	hexaflumuron iodosulfuron-methyl-sodium isoxaflutole methiocarb sulfone metsulfuron-methyl	imazalil iprodione linuron methomyl monocrotophos	
Pesticide Mix Q				
DRE-GS09000151AL	Pesticide Mix Q 100 µg/mL in Acetonitrile(‡)			5x1ml
nicosulfuron oxycarboxin pirimiphos methyl propamocarb HCl propoxycarbazon sodium	omethoate oxydemeton-methyl primisulfuron-methyl Propanil propyzamide (pronamide)	oxadixyl phenmedipham promecarb propaquizafop prothioconazole	oxamyl phoxim propamocarb propoxur pyraclostrobin	
Pesticide Mix R				
DRE-GS09000152AL	Pesticide Mix R 100 µg/mL in Acetonitrile(‡)			5x1ml
pyridate spinosad (Spinosyn A & D) terbutol thifensulfuron-methyl thiofanox-sulfone	pyroxsulam sulfentrazone thiabendazole thiobencarb thiophanate methyl	rimsulfuron tebufenozide thiacloprid thiodicarb tolylfluanid	sethoxydim tembotrione thiamethoxam thiofanox triasulfuron	
Pesticide Mix S				
DRE-GS09000153AL	Pesticide Mix S 100 µg/mL in Acetonitrile(‡)			5x1ml
	tribenuron-methyl trichlorfon tridemorph (mix of isomers) trifloxystrobin triflumizole XMC (3,5-xylyl methylcarbamate)			

Pesticide mixtures

Product code	Description	
Pesticide Mix, 11-3578		
DRE-GS09000356AL	Pesticide Mix, 11-3578 100 µg/mL in Acetonitrile(‡)	5x1.5ml
	p,p'-DDT bifenthrin endosulfan I endosulfan sulfate g-BHC o,p'-DDT	acetochlor p,p'-DDE endosulfan II a-BHC lambda cyhalothrin
Pesticide Mix, 12-3543		
DRE-GS09000375HE	Pesticide Mix, 12-3543 200 µg/mL in Hexane(‡)	5x1ml
	diazinon tetrachlorvinphos (Rabon) parathion carbophenothion hexachlorobenzene decachlorobiphenyl (BZ# 209)	malathion ethion methyl parathion chlorpyrifos mirex 2,4,5,6-tetrachloro-m-xylene
Pesticide Mix, 13-1126		
DRE-GS09000250AC	Pesticide Mix, 13-1126 100 µg/mL in Acetone(‡)	5x1ml
	resmethrin fenamiphos sulfone tetramethrin etoxazole triticonazole fenbuconazole tau-fluvalinate	fenamiphos sulfoxide bifenthrin phosmet fenamidone lambda cyhalothrin fenthion
Pesticide Mix, 13-3470		
DRE-GS09000254AC	Pesticide Mix, 13-3470 200 µg/mL in Acetonitrile(‡)	5x1ml
	Lontrel (clopyralid) triclopyr bentazon Dicamba MCPA acid imazamox imazethapyr	chloramben picloram 2,4-dichlorophenoxyacetic acid MCPP (Mecoprop) imazapyr imazaquin
Pesticide Mix, 16-3583		
DRE-GS09000357AL	Pesticide Mix, 16-3583 100 µg/mL in Acetonitrile:Water 1:1(‡)	5x1.5ml
	pyriproxyfen spinetoram (mixture of isomers) spirodiclofen streptomycin sulfate as streptomycin tebufenozide thiabendazole thiophanate methyl trifloxystrobin	sethoxydim spinosad (Spinosyn A & D) spirotetramat tebuconazol (Folicur) Abate® thiamethoxam Tolfenpyrade
Pesticide Mix, 19-3911		
DRE-GS09000394AL	Pesticide Mix, 19-3911 in Acetonitrile(‡)	5x1ml
	abamectin [5 mg/L] fenoxycarb [2 mg/L] paclobutrazol [4 mg/L] trifloxystrobin [2 mg/L] Cinerin I [3,3 mg/L]	bifenazate [2 mg/L] imazalil [2 mg/L] spinosad (Spinosyn A&D) [2 mg/L] permethrin (mixt. of isomers) [1 mg/L] Jasmolin I [3,3 mg/L]
	bifenthrin [2 mg/L] imidacloprid [4 mg/L] spiromesifen [2 mg/L] daminozide [10 mg/L]	etoxazole [2 mg/L] Systhane TM [2 mg/L] spirotetramat [2 mg/L] pyrethrin I [3,3 mg/L]
Pesticide Mix, 20-1125		
DRE-GS09000251AC	Pesticide Mix, 20-1125 100 µg/mL in Acetone(‡)	5x1ml
	flumioxazin cyprodinil fenamiphos tribufos carfentrazone-ethyl	triadimefon fipronil flutolanil Systhane (TM) hexazinon
	fosthiazate allethrin fludioxonil buprofezin piperonyl butoxide	pendimethalin methidathion oxadiazon ethion Propiconazol (mixture of isomers)

Pesticide mixtures

Product code	Description			
Pesticide Mix, 20-3435				
DRE-GA09000203EA	Pesticide Mix, 20-3435 100 µg/mL in Ethyl Acetate(‡)			1ml
chlorfenapyr	chlorthal-dimethyl (dacthal)	fenitrothion	Folpet	
Nitrofen	2-phenylphenol	methyl parathion	Perthane (TM)	
procymidone	fenchlorphos	tecnazene	iprodione	
chlorfenson	chlorpropham	dichlofluanid	fenson	
propham	tetradifon	tetrasul	cis-1,2,3,6-tetrahydrophthalimide	
Pesticide Mix, 20-3912				
DRE-GS09000395AL	Pesticide Mix, 20-3912 in Acetonitrile(‡)			5x1ml
acephate [4 mg/L]	aldicarb [4 mg/L]	azoxystrobin [2 mg/L]	carbaryl [2 mg/L]	
chlorpyrifos [2 mg/L]	baythroid (four isomers) [10 mg/L]	dichlorvos [1 mg/L]	diazinon [2 mg/L]	
dimethoate [2 mg/L]	ethoprophos (prophos) [2 mg/L]	kresoxim methyl [4 mg/L]	malathion [2 mg/L]	
metalaxyl [2 mg/L]	methiocarb [2 mg/L]	methomyl [4 mg/L]	methyl parathion [2 mg/L]	
oxamyl [10 mg/L]	prallethrin [2 mg/L]	pyridaben [2 mg/L]	cypermethrin (mix of isomers) [10 mg/L]	
Pesticide Mix, 22-3430				
DRE-GA09000202EA	Pesticide Mix, 22-3430 100 µg/mL in Ethyl Acetate(‡)			1ml
bromopropylate	chlorobenzilate	chlorothalonil	baythroid (mixture of isomers)	
2,6-Dichlorobenzonitrile (Dichlobenil)	Ethalfuralin	danitol	fenvalerate (mixture of diastereoisomers)	
lambda cyhalothrin	trifluralin	diphenylamine	kelthane ® (dicofol)	
captan	Dichloran	chlozolinate	profluralin	
tefluthrin	vinclozolin	esfenvalerate	Benfluralin (Benefin)	
Captafol	pentachloroaniline			
Pesticide Mix, 23-3579				
DRE-GS09000355AL	Pesticide Mix, 23-3579 100 µg/mL in Acetonitrile(‡)			5x1ml
2,4-D ethylhexyl ester	2,4-D-isopropyl	Guthion®	boscalid	
bromopropylate	buprofezin	cadusafos	chlorfenapyr	
chlorobenzilate	chlorpyrifos	chloropyrifos-methyl	alpha-cypermethrin	
cypermethrin (mix of isomers)	cyprodinil	diazinon	dichlorvos	
kelthane ® (dicofol)	disulfoton	ethion	fenitrothion	
danitol	fenthion	fenthion-sulfone		
Pesticide Mix, 23-3581				
DRE-GS09000359AL	Pesticide Mix, 23-3581 100 µg/mL in Acetonitrile(‡)			5x1.5ml
acephate	acetamiprid	aldicarb	abamectin	
azoxystrobin	bifenazate	carbaryl	clothianidin	
cyantraniliprole	daminozide	Difenoconazole (mixture of isomers)	Diflubenzuron	
dimethoate	diuron	endosulfan sulfate	etoxazole	
N,N'-Ethylenethiourea	famoxadone	fenazaquin	fenbuconazole	
fenbutatin oxide	fenoxycarb	fenpyroximate (Mixture of Raceimers)		
Pesticide Mix, 24-3469				
DRE-GS09000247AL	Pesticide Mix, 24-3469 200 µg/mL in Acetonitrile(‡)			5x1ml
methiocarb	methomyl	methoxyfenozide	metsulfuron-methyl	
nitenpyram	norflurazon	desmethyl norflurazon	oryzalin	
oxamyl	picoxystrobin	propoxur	pyraclostrobin	
simazine	spirodiclofen	spirotetramat	sulfometuron methyl	
sulfosulfuron	sulfoxaflor	thiacloprid	thiamethoxam	
thiazopyr	thifensulfuron-methyl	trifloxystrobin	trifloxysulfuron-sodium	
Pesticide Mix, 24-3582				
DRE-GS09000358AL	Pesticide Mix, 24-3582 100 µg/mL in Acetonitrile(‡)			5x1.5ml
fludioxonil	flufenoxuron	flusilazole	hexythiazox	
imidacloprid	isoprocarb	isoproturon	metalaxyl	
methamidophos	methiocarb	methomyl	methoxyfenozide	
monocrotophos	1-naphthol	napropamide	oxamyl	
oxytetracycline	pencycuron	pendimethalin	phenmedipham	
pymetrozine	pyraclostrobin	pyridaben	pyrimethanil	

Pesticide mixtures

Product code	Description		
Pesticide Mix, 25			
DRE-GH09000205AL	Pesticide Mix 100 µg/mL in Acetonitrile(‡)		10x1ml
lufenuron	metalaxyl	methidathion	methomyl
methoxyfenozide	novaluron	omethoate	oxamyl
phenthoate	prochloraz	profenofos	propargite
pyraclostrobin	pyrimethanil	spinetoram (mixture of isomers)	spinosad (Spinosyn A & D)
tebuconazol (Folicur)	tebufenozide	thiabendazole	thiacloprid
Tolfenpyrade	triazophos	trifloxystrobin	imidacloprid
malathion			
Pesticide Mix, 25-3580			
DRE-GS09000354AL	Pesticide Mix, 25-3580 100 µg/mL in Acetonitrile(‡)		5x1.5ml
Imazalil	isocarbofos	Isofenphos	malathion
methidathion	Systhane TM	2-phenylphenol	oxyfluorfen
parathion	methyl parathion	phenthoate	phorate
phosalone	phosmet	Pirimicarb	prochloraz
propargite	Propiconazol (mixture of isomers)	Tokuthion®	Pyridaphenthion
quinalphos	tetradifon	profenofos	triazophos
pirimiphos methyl			
Pesticide Mix, 25-3808			
DRE-GA09000385MF	Pesticide Mix, 25-3808 200 µg/mL in Methanol(‡)		5x1ml
methiocarb	desmethyl norflurazon	methomyl	oryzalin
oxamyl	propoxur	bladex	thiacloprid
thiamethoxam	thiazopir	aldicarb sulfone	chlorantraniliprole
clothianidin	Imidacloprid-olefin	aldicarb sulfoxide	bromacil
carbaryl	carbofuran	diuron	imidacloprid
metalaxyl	aldicarb	norflurazon	simazine
3-hydroxycarbofuran			
Pesticide Mix, 26			
DRE-GH09000204AL	Pesticide Mix 100 µg/mL in Acetonitrile(‡)		10x1ml
acetamiprid	azoxystrobin	bifenazate	bifenthrin
boscalid	buprofezin	carbaryl	carbofuran
chlorantraniliprole	chlorfluazuron	clothianidin	cyprodinil
dicrotophos	Difenoconazole	Diffubenzuron	dimethoate
dinotefuran	fenbuconazole	fenbutatin oxide	fipronil
flonicamid	flubendiamide	fludioxonil	flufenoxuron
hexaflumuron	imazalil		
Pesticide Mix, 26-3428			
DRE-GA09000200EA	Pesticide Mix, 26-3428 100 µg/mL in Ethyl Acetate(‡)		1ml
b-BHC	d-BHC	endosulfan I	heptachlor epoxide isomer A
o,p'-DDE	pentachloronitrobenzene	mirex	oxychlordane
a-BHC	hexachlorobenzene	g-BHC	heptachlor
aldrin	heptachlor epoxide isomer B	trans-chlordane	cis-chlordane
dieldrin	endrin	endosulfan II	endosulfan sulfate
methoxychlor	p,p'-DDT	p,p'-DDD	p,p'-DDE
o,p'-DDD	o,p'-DDT		
Pesticide Mix, 31-1123			
DRE-GS09000253AC	Pesticide Mix, 31-1123 100 µg/mL in Acetone(‡)		5x1ml
atrazine desisopropyl	desethyl atrazine	Benfluralin (Benefin)	parathion
dimethoate	simazine	atrazine	g-BHC
propyzamide (pronamide)	chlorothalonil	metribuzin	methyl parathion
alachlor	malathion	phorate	MGK-264 Mix of Isomers
endosulfan I	chlorfenapyr	fluazifop-P-butyl	endosulfan II
iprodione	danitol	Guthion®	coumaphos
permethrin (mixture of isomers)	baythroid (mixture of isomers)	cypermethrin (mix of isomers)	fenvalerate (mixture of diastereoisomers)
deltamethrin	phenothrin	prodiamine	

Pesticide mixtures

Product code	Description	
Rodenticides Mix, 6-3464		
DRE-GS09000248AL	Rodenticides Mix, 6-3464 100 µg/mL in Acetonitrile(‡) bromadiolol diphacinone chlorophacinone Brodifacoum difenacoum warfarin	5x1ml

PCBS AND RELATED COMPOUNDS



PCB's and related compounds

Product code	Description	
Aroclor		
DRE-L20101600CY	Aroclor 1016, CAS 12674-11-2, 10 µg/mL in Cyclohexane	10ml
DRE-YA20101600CY	Aroclor 1016, CAS 12674-11-2, 1000 µg/mL in Cyclohexane	1ml
DRE-X20101600IO	Aroclor 1016, CAS 12674-11-2, 100 µg/mL in Isooctane	10ml
DRE-GA20101600HE	Aroclor 1016, CAS 12674-11-2, 1000 µg/mL in n-Hexane	1ml
DRE-GH20101600HE	Aroclor 1016, CAS 12674-11-2, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20122100CY	Aroclor 1221, CAS 11104-28-2, 10 µg/mL in Cyclohexane	10ml
DRE-YA20122100CY	Aroclor 1221, CAS 11104-28-2, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20122100HE	Aroclor 1221, CAS 11104-28-2, 1000 µg/mL in n-Hexane	1ml
DRE-GH20122100HE	Aroclor 1221, CAS 11104-28-2, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20123200CY	Aroclor 1232, CAS 11141-16-5, 10 µg/mL in Cyclohexane	10ml
DRE-X20123200CY	Aroclor 1232, CAS 11141-16-5, 100 µg/mL in Cyclohexane	10ml
DRE-YA20123200CY	Aroclor 1232, CAS 11141-16-5, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20123200HE	Aroclor 1232, CAS 11141-16-5, 1000 µg/mL in n-Hexane	1ml
DRE-GH20123200HE	Aroclor 1232, CAS 11141-16-5, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20124200CY	Aroclor 1242, CAS 53469-21-9, 10 µg/mL in Cyclohexane	10ml
DRE-X20124200CY	Aroclor 1242, CAS 53469-21-9, 100 µg/mL in Cyclohexane	10ml
DRE-YA20124200CY	Aroclor 1242, CAS 53469-21-9, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20124200HE	Aroclor 1242, CAS 53469-21-9, 1000 µg/mL in n-Hexane	1ml
DRE-GH20124200HE	Aroclor 1242, CAS 53469-21-9, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20124800CY	Aroclor 1248, CAS 12672-29-6, 10 µg/mL in Cyclohexane	10ml
DRE-X20124800CY	Aroclor 1248, CAS 12672-29-6, 100 µg/mL in Cyclohexane	10ml
DRE-YA20124800CY	Aroclor 1248, CAS 12672-29-6, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20124800HE	Aroclor 1248, CAS 12672-29-6, 1000 µg/mL in n-Hexane	1ml
DRE-GH20124800HE	Aroclor 1248, CAS 12672-29-6, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20125400CY	Aroclor 1254, CAS 11097-69-1, 10 µg/mL in Cyclohexane	10ml
DRE-X20125400CY	Aroclor 1254, CAS 11097-69-1, 100 µg/mL in Cyclohexane	10ml
DRE-YA20125400CY	Aroclor 1254, CAS 11097-69-1, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20125400HE	Aroclor 1254, CAS 11097-69-1, 1000 µg/mL in n-Hexane	1ml
DRE-GH20125400HE	Aroclor 1254, CAS 11097-69-1, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20126000CY	Aroclor 1260, CAS 11096-82-5, 10 µg/mL in Cyclohexane	10ml
DRE-X20126000CY	Aroclor 1260, CAS 11096-82-5, 100 µg/mL in Cyclohexane	10ml
DRE-YA20126000CY	Aroclor 1260, CAS 11096-82-5, 1000 µg/mL in Cyclohexane	1ml
DRE-GA20126000HE	Aroclor 1260, CAS 11096-82-5, 1000 µg/mL in n-Hexane	1ml
DRE-GH20126000HE	Aroclor 1260, CAS 11096-82-5, 1000 µg/mL in n-Hexane	5x1ml
DRE-GA20126200HE	Aroclor 1262, CAS 37324-23-5, 1000 µg/mL in n-Hexane	1ml
DRE-GH20126200HE	Aroclor 1262, CAS 37324-23-5, 1000 µg/mL in n-Hexane	5x1ml
DRE-GA20126800HE	Aroclor 1268, CAS 11100-14-4, 1000 µg/mL in n-Hexane	1ml
DRE-GH20126800HE	Aroclor 1268, CAS 11100-14-4, 1000 µg/mL in n-Hexane	5x1ml
DRE-L20206000CY	Aroclor 5060, CAS 73928-90-2, 10 µg/mL in Cyclohexane	10ml
DRE-L20243200CY	Aroclor 5432, CAS 63496-31-1, 10 µg/mL in Cyclohexane	10ml
DRE-LA20244200CY	Aroclor 5442, CAS 12642-23-8, 10 µg/mL in Cyclohexane	1ml
DRE-L20246000CY	Aroclor 5460, CAS 11126-42-4, 10 µg/mL in Cyclohexane	10ml


Aroclor 1016

CAS 12674-11-2	MW n/a		
DRE-GA09010301ME	Aroclor 1016 1000 µg/mL in Methanol(±)	1ml	No Structure








Aroclor 1221

CAS 11104-28-2	MW n/a		
DRE-GA09010302ME	Aroclor 1221 1000 µg/mL in Methanol(±)	1ml	No Structure










PCB's and related compounds

Product code	Description		
Aroclor 1232			
CAS 11141-16-5	MW n/a		
DRE-GA09010303ME	Aroclor 1232 1000 µg/mL in Methanol(±)	1ml	No Structure
Aroclor 1242			
CAS 53469-21-9	MW n/a		
DRE-GA09010304ME	Aroclor 1242 1000 µg/mL in Methanol(±)	1ml	No Structure
Aroclor 1248			
CAS 12672-29-6	MW n/a		
DRE-GA09010305ME	Aroclor 1248 1000 µg/mL in Methanol(±)	1ml	No Structure
Aroclor 1260			
CAS 11096-82-5	MW n/a		
DRE-GA09010307ME	Aroclor 1260 1000 µg/mL in Methanol(±)	1ml	No Structure
Aroclor-Mix 1242,1254,1260 1:1:1			
DRE-L20258000CY	Aroclor-Mix 1242,1254,1260 1:1:1 10 µg/mL in Cyclohexane		10ml
	Aroclor 1242		
	Aroclor 1254		
	Aroclor 1260		
Chlorinated Terphenyl Mix 1			
DRE-LA20399995HE	Chlorinated Terphenyl Mix 1 10 µg/mL in n-Hexane		1ml
	2,2",4,4",5,5"-Hexachloro-p-terphenyl	3,3',3",4,4"-Pentachloro-m-terphenyl	
	3,3",4,4",5,5"-Hexachloro-p-terphenyl	3,3",4,4"-Tetrachloro-o-terphenyl	
	3,3",4,4"-Tetrachloro-p-terphenyl	3,3",5,5"-Tetrachloro-p-terphenyl	
	3,3"-Dichloro-o-terphenyl	3,3"-Dichloro-p-terphenyl	
	3',4,4"-Trichloro-m-terphenyl		
2'-Chloro-2,3,4,5,6-pentadeuterio-1,1'-biphenyl			
CAS 51624-35-2	MW 193.6837	C ₁₂ ² H ₅ H ₄ Cl	
DRE-XA20000101IO	2-Chlorobiphenyl D5 100 µg/mL in Isooctane	1.1ml	
Clophen A 30			
CAS 55600-34-5	MW n/a		
DRE-X20303000CY	Clophen A 30 100 µg/mL in Cyclohexane	10ml	No Structure
Clophen A 40			
CAS 52306-32-8	MW n/a		
DRE-X20304000CY	Clophen A 40 100 µg/mL in Cyclohexane	10ml	No Structure










PCB's and related compounds

Product code	Description			
Clophen A 50				
CAS 8068-44-8	MW n/a			No Structure
DRE-X20305000CY	Clophen A 50 100 µg/mL in Cyclohexane		10ml	
Clophen A 60				
CAS 11096-99-4	MW n/a			No Structure
DRE-X20306000CY	Clophen A 60 100 µg/mL in Cyclohexane		10ml	
3,3''-Dichloro-o-terphenyl				
CAS 2007947-35-3	MW 299.1939	C ₁₈ H ₁₂ Cl ₂		
DRE-LA20310332HE	3,3''-Dichloro-o-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3''-Dichloro-p-terphenyl				
CAS 159392-36-6	MW 299.1939	C ₁₈ H ₁₂ Cl ₂		
DRE-LA20310334HE	3,3''-Dichloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	
4,4''-Dichloro-o-terphenyl				
CAS 21711-56-8	MW 299.1939	C ₁₈ H ₁₂ Cl ₂		
DRE-LA20310442HE	4,4''-Dichloro-o-terphenyl 10 µg/mL in n-Hexane		1ml	
4,4''-Dichloro-p-terphenyl				
CAS 21711-52-4	MW 299.1939	C ₁₈ H ₁₂ Cl ₂		
DRE-LA20310444HE	4,4''-Dichloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	
2,2'',3',4,4'',5,5''-Heptachloro-m-terphenyl				
CAS n/a	MW 471.4192	C ₁₈ H ₇ Cl ₇		
DRE-LA20388553HE	2,2'',3',4,4'',5,5''-Heptachloro-m-terphenyl 10 µg/mL in n-Hexane		1ml	
2,2'',4,4'',5,5''-Hexachloro-p-terphenyl				
CAS n/a	MW 436.9741	C ₁₈ H ₆ Cl ₆		
DRE-LA20387554HE	2,2'',4,4'',5,5''-Hexachloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3'',4,4'',5,5''-Hexachloro-p-terphenyl				
CAS 1312994-27-6	MW 436.9741	C ₁₈ H ₆ Cl ₆		
DRE-LA20387654HE	3,3'',4,4'',5,5''-Hexachloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	










PCB's and related compounds

Product code	Description			
PCB 1 (2-Chlorobiphenyl)				
CAS 2051-60-7	MW 188.6529	C ₁₂ H ₉ Cl		
DRE-C20000100	PCB No. 1		50mg	
DRE-L20000100IO	PCB No. 1 10 µg/mL in Isooctane		10ml	
PCB 2 (3-Chlorobiphenyl)				
CAS 2051-61-8	MW 188.6529	C ₁₂ H ₉ Cl		
DRE-C20000200	PCB No. 2		50mg	
DRE-L20000200IO	PCB No. 2 10 µg/mL in Isooctane		10ml	
PCB 3 (4-Chlorobiphenyl)				
CAS 2051-62-9	MW 188.6529	C ₁₂ H ₉ Cl		
DRE-C20000300	PCB No. 3		50mg	
DRE-L20000300IO	PCB No. 3 10 µg/mL in Isooctane		10ml	
PCB 4 (2,2'-Dichlorobiphenyl)				
CAS 13029-08-8	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20000400	PCB No. 4		25mg	
DRE-L20000400IO	PCB No. 4 10 µg/mL in Isooctane		10ml	
PCB 5 (2,3-Dichlorobiphenyl)				
CAS 16605-91-7	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20000500	PCB No. 5		50mg	
DRE-L20000500IO	PCB No. 5 10 µg/mL in Isooctane		10ml	
PCB 7 (2,4-Dichlorobiphenyl)				
CAS 33284-50-3	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20000700	PCB No. 7		25mg	
PCB 8 (2,4'-Dichlorobiphenyl)				
CAS 34883-43-7	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20000800	PCB No. 8		25mg	
DRE-L20000800IO	PCB No. 8 10 µg/mL in Isooctane		10ml	
PCB 9 (2,5-Dichlorobiphenyl)				
CAS 34883-39-1	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20000900	PCB No. 9		50mg	
PCB 10 (2,6-Dichlorobiphenyl)				
CAS 33146-45-1	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20001000	PCB No. 10		25mg	
DRE-L20001000IO	PCB No. 10 10 µg/mL in Isooctane		10ml	










PCB's and related compounds

Product code	Description			
PCB 11 (3,3'-Dichlorobiphenyl)				
CAS 2050-67-1	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20001100	PCB No. 11		25mg	
DRE-L20001100IO	PCB No. 11 10 µg/mL in Isooctane		10ml	
PCB 13 (3,4'-Dichlorobiphenyl)				
CAS 2974-90-5	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20001300	PCB No. 13		5mg	
PCB 14 (3,5-Dichlorobiphenyl)				
CAS 34883-41-5	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20001400	PCB No. 14		50mg	
DRE-L20001400IO	PCB No. 14 10 µg/mL in Isooctane		10ml	
PCB 15 (4,4'-Dichlorobiphenyl)				
CAS 2050-68-2	MW 223.0979	C ₁₂ H ₈ Cl ₂		
DRE-C20001500	PCB No. 15		10mg	
PCB 16 (2,2',3-Trichlorobiphenyl)				
CAS 38444-78-9	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-E20001600	PCB No. 16		5mg	
PCB 18 (2,2',5-Trichlorobiphenyl)				
CAS 37680-65-2	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20001800	PCB No. 18(‡)		25mg	
DRE-L20001800IO	PCB No. 18 10 µg/mL in Isooctane		10ml	
PCB 19 (2,2',6-Trichlorobiphenyl)				
CAS 38444-73-4	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-E20001900	PCB No. 19		5mg	
DRE-L20001900IO	PCB No. 19 10 µg/mL in Isooctane		10ml	
PCB 20 (2,3,3'-Trichlorobiphenyl)				
CAS 38444-84-7	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20002000	PCB No. 20		10mg	
DRE-L20002000IO	PCB No. 20 10 µg/mL in Isooctane		10ml	
PCB 21 (2,3,4-Trichlorobiphenyl)				
CAS 55702-46-0	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20002100	PCB No. 21		10mg	
DRE-L20002100IO	PCB No. 21 10 µg/mL in Isooctane		10ml	










PCB's and related compounds

Product code	Description			
PCB 22 (2,3,4'-Trichlorobiphenyl)				
CAS 38444-85-8	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-L20002200IO	PCB No. 22 10 µg/mL in Isooctane		10ml	
PCB 23 (2,3,5-Trichlorobiphenyl)				
CAS 55720-44-0	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-L20002300IO	PCB No. 23 10 µg/mL in Isooctane		10ml	
PCB 26 (2,3',5-Trichlorobiphenyl)				
CAS 38444-81-4	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-L20002600IO	PCB No. 26 10 µg/mL in Isooctane		10ml	
PCB 28 (2,4,4'-Trichlorobiphenyl)				
CAS 7012-37-5	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20002800	PCB No. 28(±)		10mg	
DRE-L20002800IO	PCB No. 28 10 µg/mL in Isooctane		10ml	
PCB 29 (2,4,5-Trichlorobiphenyl)				
CAS 15862-07-4	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20002900	PCB No. 29		10mg	
DRE-L20002900IO	PCB No. 29 10 µg/mL in Isooctane		10ml	
PCB 30 (2,4,6-Trichlorobiphenyl)				
CAS 35693-92-6	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20003000	PCB No. 30		25mg	
DRE-L20003000CY	PCB No. 30 10 µg/mL in Cyclohexane(±)		10ml	
DRE-L20003000IO	PCB No. 30 10 µg/mL in Isooctane(±)		10ml	
DRE-X20003000IO	PCB No. 30 100 µg/mL in Isooctane(±)		10ml	
DRE-XA20003000IO	PCB No. 30 100 µg/mL in Isooctane		1ml	
PCB 31 (2,4',5-Trichlorobiphenyl)				
CAS 16606-02-3	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20003100	PCB No. 31(±)		25mg	
DRE-L20003100CY	PCB No. 31 10 µg/mL in Cyclohexane		10ml	
DRE-L20003100IO	PCB No. 31 10 µg/mL in Isooctane		10ml	
PCB 33 (2,3',4'-Trichlorobiphenyl)				
CAS 38444-86-9	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20003300	PCB No. 33		10mg	
DRE-L20003300IO	PCB No. 33 10 µg/mL in Isooctane		10ml	
PCB 34 (2,3',5'-Trichlorobiphenyl)				
CAS 37680-68-5	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20003400	PCB No. 34		5mg	
DRE-L20003400IO	PCB No. 34 10 µg/mL in Isooctane(±)		10ml	










PCB's and related compounds

Product code	Description			
PCB 35 (3,3',4-Trichlorobiphenyl)				
CAS 37680-69-6	MW 257.543	C ₁₂ H ₇ Cl ₃		
DRE-C20003500	PCB No. 35		5mg	
DRE-L20003500IO	PCB No. 35 10 µg/mL in Isooctane		10ml	
PCB 37F (3'-Fluoro-3,4,4'-trichlorobiphenyl)				
CAS 1191034-39-5	MW 275.5334	C ₁₂ H ₆ Cl ₃ F		
DRE-XA15901037IO	PCB 37F (3'-Fluoro-3,4,4'-trichlorobiphenyl) 100 µg/mL in Isooctane		1ml	
PCB 40 (2,2',3,3'-Tetrachlorobiphenyl)				
CAS 38444-93-8	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20004000	PCB No. 40		25mg	
DRE-L20004000IO	PCB No. 40 10 µg/mL in Isooctane		10ml	
PCB 43 (2,2',3,5-Tetrachlorobiphenyl)				
CAS 70362-46-8	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20004300IO	PCB No. 43 10 µg/mL in Isooctane		10ml	
PCB 44 (2,2',3,5'-Tetrachlorobiphenyl)				
CAS 41464-39-5	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20004400	PCB No. 44		25mg	
DRE-L20004400IO	PCB No. 44 10 µg/mL in Isooctane		10ml	
PCB 46 (2,2',3,6'-Tetrachlorobiphenyl)				
CAS 41464-47-5	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20004600IO	PCB No. 46 10 µg/mL in Isooctane		10ml	
PCB 47 (2,2',4,4'-Tetrachlorobiphenyl)				
CAS 2437-79-8	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20004700	PCB No. 47		25mg	
DRE-L20004700IO	PCB No. 47 10 µg/mL in Isooctane		10ml	
PCB 48 (2,2',4,5-Tetrachlorobiphenyl)				
CAS 70362-47-9	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20004800	PCB No. 48		5mg	
DRE-L20004800IO	PCB No. 48 10 µg/mL in Isooctane		10ml	
PCB 49 (2,2',4,5'-Tetrachlorobiphenyl)				
CAS 41464-40-8	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20004900	PCB No. 49		20mg	
DRE-L20004900IO	PCB No. 49 10 µg/mL in Isooctane		10ml	










PCB's and related compounds

Product code	Description			
PCB 50 (2,2',4,6-Tetrachlorobiphenyl)				
CAS 62796-65-0	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20005000	PCB No. 50		5mg	
DRE-L20005000IO	PCB No. 50 10 µg/mL in Isooctane		10ml	
PCB 51 (2,2',4,6'-Tetrachlorobiphenyl)				
CAS 68194-04-7	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20005100IO	PCB No. 51 10 µg/mL in Isooctane		10ml	
PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)				
CAS 35693-99-3	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20005200	PCB No. 52(‡)		10mg	
DRE-L20005200IO	PCB No. 52 10 µg/mL in Isooctane(‡)		10ml	
PCB 53 (2,2',5,6'-Tetrachlorobiphenyl)				
CAS 41464-41-9	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20005300	PCB No. 53		25mg	
DRE-L20005300IO	PCB No. 53 10 µg/mL in Isooctane		10ml	
PCB 54 (2,2',6,6'-Tetrachlorobiphenyl)				
CAS 15968-05-5	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20005400	PCB No. 54		25mg	
DRE-L20005400IO	PCB No. 54 10 µg/mL in Isooctane		10ml	
PCB 55 (2,3,3',4-Tetrachlorobiphenyl)				
CAS 74338-24-2	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20005500	PCB No. 55		5mg	
DRE-L20005500IO	PCB No. 55 10 µg/mL in Isooctane		10ml	
PCB 56 (2,3,3',4'-Tetrachlorobiphenyl)				
CAS 41464-43-1	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20005600	PCB No. 56		5mg	
DRE-L20005600IO	PCB No. 56 10 µg/mL in Isooctane		10ml	
PCB 61 (2,3,4,5-Tetrachlorobiphenyl)				
CAS 33284-53-6	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20006100	PCB No. 61		10mg	
DRE-L20006100IO	PCB No. 61 10 µg/mL in Isooctane		10ml	
PCB 62 (2,3,4,6-Tetrachlorobiphenyl)				
CAS 54230-22-7	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20006200	PCB No. 62		5mg	










PCB's and related compounds

Product code	Description			
PCB 65 (2,3,5,6-Tetrachlorobiphenyl)				
CAS 33284-54-7	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20006500	PCB No. 65		25mg	
DRE-L20006500IO	PCB No. 65 10 µg/mL in Isooctane		10ml	
PCB 66 (2,3',4,4'-Tetrachlorobiphenyl)				
CAS 32598-10-0	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20006600	PCB No. 66		20mg	
DRE-L20006600IO	PCB No. 66 10 µg/mL in Isooctane		10ml	
PCB 68 (2,3',4,5'-Tetrachlorobiphenyl)				
CAS 73575-52-7	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20006800IO	PCB No. 68 10 µg/mL in Isooctane		10ml	
PCB 69 (2,3',4,6-Tetrachlorobiphenyl)				
CAS 60233-24-1	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20006900	PCB No. 69		5mg	
DRE-L20006900IO	PCB No. 69 10 µg/mL in Isooctane		10ml	
PCB 70 (2,3',4',5-Tetrachlorobiphenyl)				
CAS 32598-11-1	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20007000	PCB No. 70		10mg	
DRE-L20007000IO	PCB No. 70 10 µg/mL in Isooctane		10ml	
PCB 73 (2,3',5',6-Tetrachlorobiphenyl)				
CAS 74338-23-1	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20007300	PCB No. 73		5mg	
DRE-L20007300IO	PCB No. 73 10 µg/mL in Isooctane		10ml	
PCB 74 (2,4,4',5-Tetrachlorobiphenyl)				
CAS 32690-93-0	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20007400IO	PCB No. 74 10 µg/mL in Isooctane		10ml	
PCB 75 (2,4,4',6-Tetrachlorobiphenyl)				
CAS 32598-12-2	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20007500	PCB No. 75		5mg	
PCB 77 (3,3',4,4'-Tetrachlorobiphenyl)				
CAS 32598-13-3	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20007700	PCB No. 77		20mg	
DRE-L20007700IO	PCB No. 77 10 µg/mL in Isooctane		10ml	

PCB's and related compounds

Product code	Description			
PCB 78 (3,3',4,5-Tetrachlorobiphenyl)				
CAS 70362-49-1	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-E20007800	PCB No. 78		5mg	
DRE-L20007800IO	PCB No. 78 10 µg/mL in Isooctane		10ml	
PCB 80 (3,3',5,5'-Tetrachlorobiphenyl)				
CAS 33284-52-5	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-L20008000IO	PCB No. 80 10 µg/mL in Isooctane		10ml	
PCB 81 (3,4,4',5-Tetrachlorobiphenyl)				
CAS 70362-50-4	MW 291.988	C ₁₂ H ₆ Cl ₄		
DRE-C20008100	PCB No. 81		10mg	
DRE-L20008100IO	PCB No. 81 10 µg/mL in Isooctane		10ml	
PCB 82 (2,2',3,3',4-Pentachlorobiphenyl)				
CAS 52663-62-4	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20008200	PCB No. 82		5mg	
PCB 85 (2,2',3,4,4'-Pentachlorobiphenyl)				
CAS 65510-45-4	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20008500	PCB No. 85		10mg	
PCB 87 (2,2',3,4,5'-Pentachlorobiphenyl)				
CAS 38380-02-8	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20008700	PCB No. 87		10mg	
DRE-L20008700IO	PCB No. 87 10 µg/mL in Isooctane		10ml	
PCB 89 (2,2',3,4,6'-Pentachlorobiphenyl)				
CAS 73575-57-2	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20008900	PCB No. 89		5mg	
DRE-L20008900IO	PCB No. 89 10 µg/mL in Isooctane		10ml	
PCB 90 (2,2',3,4',5-Pentachlorobiphenyl)				
CAS 68194-07-0	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20009000	PCB No. 90		5mg	
DRE-L20009000IO	PCB No. 90 10 µg/mL in Isooctane		10ml	
PCB 91 (2,2',3,4',6-Pentachlorobiphenyl)				
CAS 68194-05-8	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20009100	PCB No. 91		5mg	










PCB's and related compounds

Product code	Description			
PCB 95 (2,2',3,5',6-Pentachlorobiphenyl)				
CAS 38379-99-6	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20009500	PCB No. 95		5mg	
DRE-L20009500IO	PCB No. 95 10 µg/mL in Isooctane		10ml	
PCB 97 (2,2',3,4',5'-Pentachlorobiphenyl)				
CAS 41464-51-1	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20009700	PCB No. 97		10mg	
PCB 99 (2,2',4,4',5-Pentachlorobiphenyl)				
CAS 38380-01-7	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20009900	PCB No. 99		5mg	
DRE-L20009900IO	PCB No. 99 10 µg/mL in Isooctane		10ml	
PCB 100 (2,2',4,4',6-Pentachlorobiphenyl)				
CAS 39485-83-1	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-L20010000IO	PCB No. 100 10 µg/mL in Isooctane		10ml	
PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)				
CAS 37680-73-2	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20010100	PCB No. 101(‡)		10mg	
DRE-L20010100IO	PCB No. 101 10 µg/mL in Isooctane		10ml	
PCB 103 (2,2',4,5',6-Pentachlorobiphenyl)				
CAS 60145-21-3	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20010300	PCB No. 103		10mg	
DRE-L20010300IO	PCB No. 103 10 µg/mL in Isooctane		10ml	
PCB 104 (2,2',4,6,6'-Pentachlorobiphenyl)				
CAS 56558-16-8	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20010400	PCB No. 104		5mg	
PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl)				
CAS 32598-14-4	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20010500	PCB No. 105		5mg	
DRE-L20010500IO	PCB No. 105 10 µg/mL in Isooctane		10ml	
PCB 106 (2,3,3',4,5-Pentachlorobiphenyl)				
CAS 70424-69-0	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20010600	PCB No. 106		5mg	
DRE-L20010600IO	PCB No. 106 10 µg/mL in Isooctane		10ml	










PCB's and related compounds

Product code	Description			
PCB 107 (2,3,3',4',5-Pentachlorobiphenyl)				
CAS 70424-68-9	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-L20010700IO	PCB No. 107 10 µg/mL in Isooctane		10ml	
PCB 108 (2,3,3',4,5'-Pentachlorobiphenyl)				
CAS 70362-41-3	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20010800	PCB No. 108		5mg	
PCB 110 (2,3,3',4',6-Pentachlorobiphenyl)				
CAS 38380-03-9	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20011000	PCB No. 110		5mg	
DRE-L20011000IO	PCB No. 110 10 µg/mL in Isooctane		10ml	
PCB 112 (2,3,3',5,6-Pentachlorobiphenyl)				
CAS 74472-36-9	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20011200	PCB No. 112		5mg	
DRE-L20011200IO	PCB No. 112 10 µg/mL in Isooctane		10ml	
PCB 113 (2,3,3',5',6-Pentachlorobiphenyl)				
CAS 68194-10-5	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-L20011300IO	PCB No. 113 10 µg/mL in Isooctane		10ml	
PCB 114 (2,3,4,4',5-Pentachlorobiphenyl)				
CAS 74472-37-0	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20011400	PCB No. 114		5mg	
DRE-L20011400IO	PCB No. 114 10 µg/mL in Isooctane		10ml	
PCB 116 (2,3,4,5,6-Pentachlorobiphenyl)				
CAS 18259-05-7	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20011600	PCB No. 116		10mg	
PCB 118 (2,3',4,4',5-Pentachlorobiphenyl)				
CAS 31508-00-6	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20011800	PCB No. 118(‡)		10mg	
DRE-L20011800IO	PCB No. 118 10 µg/mL in Isooctane(‡)		10ml	
PCB 119 (2,3',4,4',6-Pentachlorobiphenyl)				
CAS 56558-17-9	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20011900	PCB No. 119		5mg	
DRE-L20011900IO	PCB No. 119 10 µg/mL in Isooctane		10ml	

PCB's and related compounds

Product code	Description			
PCB 121 (2,3',4,5',6-Pentachlorobiphenyl)				
CAS 56558-18-0	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-E20012100	PCB No. 121		5mg	
DRE-L20012100IO	PCB No. 121 10 µg/mL in Isooctane		10ml	
PCB 123 (2,3',4,4',5'-Pentachlorobiphenyl)				
CAS 65510-44-3	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20012300	PCB No. 123		5mg	
DRE-L20012300IO	PCB No. 123 10 µg/mL in Isooctane		10ml	
PCB 125 (2,3',4',5',6-Pentachlorobiphenyl)				
CAS 74472-39-2	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-L20012500IO	PCB No. 125 10 µg/mL in Isooctane		10ml	
PCB 126 (3,3',4,4',5-Pentachlorobiphenyl)				
CAS 57465-28-8	MW 326.4331	C ₁₂ H ₅ Cl ₅		
DRE-C20012600	PCB No. 126		10mg	
DRE-L20012600IO	PCB No. 126 10 µg/mL in Isooctane		10ml	
PCB 128 (2,2',3,3',4,4'-Hexachlorobiphenyl)				
CAS 38380-07-3	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20012800	PCB No. 128		25mg	
DRE-L20012800IO	PCB No. 128 10 µg/mL in Isooctane		10ml	
PCB 129 (2,2',3,3',4,5-Hexachlorobiphenyl)				
CAS 55215-18-4	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20012900	PCB No. 129		5mg	
PCB 131 (2,2',3,3',4,6-Hexachlorobiphenyl)				
CAS 61798-70-7	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20013100	PCB No. 131		5mg	
PCB 132 (2,2',3,3',4,6'-Hexachlorobiphenyl)				
CAS 38380-05-1	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20013200	PCB No. 132		5mg	
DRE-L20013200IO	PCB No. 132 10 µg/mL in Isooctane		10ml	
PCB 136 (2,2',3,3',6,6'-Hexachlorobiphenyl)				
CAS 38411-22-2	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20013600	PCB No. 136		20mg	
DRE-L20013600IO	PCB No. 136 10 µg/mL in Isooctane(‡)		10ml	










PCB's and related compounds

Product code	Description			
PCB 137 (2,2',3,4,4',5-Hexachlorobiphenyl)				
CAS 35694-06-5	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-L20013700IO	PCB No. 137 10 µg/mL in Isooctane		10ml	
PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)				
CAS 35065-28-2	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20013800	PCB No. 138(‡)		10mg	
DRE-L20013800IO	PCB No. 138 10 µg/mL in Isooctane		10ml	
PCB 141 (2,2',3,4,5,5'-Hexachlorobiphenyl)				
CAS 52712-04-6	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20014100	PCB No. 141		5mg	
DRE-L20014100IO	PCB No. 141 10 µg/mL in Isooctane		10ml	
PCB 142 (2,2',3,4,5,6-Hexachlorobiphenyl)				
CAS 41411-61-4	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-L20014200IO	PCB No. 142 10 µg/mL in Isooctane		10ml	
PCB 143 (2,2',3,4,5,6'-Hexachlorobiphenyl)				
CAS 68194-15-0	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20014300	PCB No. 143		5mg	
DRE-E20014300	PCB No. 143		5mg	
DRE-L20014300IO	PCB No. 143 10 µg/mL in Isooctane		10ml	
PCB 145 (2,2',3,4,6,6'-Hexachlorobiphenyl)				
CAS 74472-40-5	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20014500	PCB No. 145		5mg	
DRE-L20014500IO	PCB No. 145 10 µg/mL in Isooctane		10ml	
PCB 146 (2,2',3,4',5,5'-Hexachlorobiphenyl)				
CAS 51908-16-8	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20014600	PCB No. 146		5mg	
DRE-L20014600IO	PCB No. 146 10 µg/mL in Isooctane		10ml	
PCB 147 (2,2',3,4',5,6-Hexachlorobiphenyl)				
CAS 68194-13-8	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20014700	PCB No. 147		5mg	
PCB 149 (2,2',3,4',5',6-Hexachlorobiphenyl)				
CAS 38380-04-0	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20014900	PCB No. 149		5mg	
DRE-L20014900IO	PCB No. 149 10 µg/mL in Isooctane		10ml	

PCB's and related compounds

Product code	Description			
PCB 150 (2,2',3,4',6,6'-Hexachlorobiphenyl)				
CAS 68194-08-1	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20015000	PCB No. 150		5mg	
PCB 151 (2,2',3,5,5',6-Hexachlorobiphenyl)				
CAS 52663-63-5	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20015100	PCB No. 151		5mg	
DRE-L20015100IO	PCB No. 151 10 µg/mL in Isooctane		10ml	
PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)				
CAS 35065-27-1	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20015300	PCB No. 153(‡)		10mg	
DRE-L20015300IO	PCB No. 153 10 µg/mL in Isooctane		10ml	
PCB 155 (2,2',4,4',6,6'-Hexachlorobiphenyl)				
CAS 33979-03-2	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20015500	PCB No. 155		10mg	
DRE-L20015500IO	PCB No. 155 10 µg/mL in Isooctane		10ml	
PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl)				
CAS 38380-08-4	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20015600	PCB No. 156		10mg	
DRE-L20015600IO	PCB No. 156 10 µg/mL in Isooctane		10ml	
PCB 157 (2,3,3',4,4',5'-Hexachlorobiphenyl)				
CAS 69782-90-7	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20015700	PCB No. 157		10mg	
DRE-L20015700IO	PCB No. 157 10 µg/mL in Isooctane		10ml	
PCB 159 (2,3,3',4,5,5'-Hexachlorobiphenyl)				
CAS 39635-35-3	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20015900	PCB No. 159		5mg	
PCB 160 (2,3,3',4,5,6-Hexachlorobiphenyl)				
CAS 41411-62-5	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20016000	PCB No. 160		5mg	
DRE-L20016000IO	PCB No. 160 10 µg/mL in Isooctane		10ml	
PCB 161 (2,3,3',4,5',6-Hexachlorobiphenyl)				
CAS 74472-43-8	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-L20016100IO	PCB No. 161 10 µg/mL in Isooctane		10ml	

PCB's and related compounds

Product code	Description			
PCB 163 (2,3,3',4',5,6-Hexachlorobiphenyl)				
CAS 74472-44-9	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20016300	PCB No. 163		5mg	
DRE-L20016300IO	PCB No. 163 10 µg/mL in Isooctane		10ml	
PCB 164 (2,3,3',4',5',6-Hexachlorobiphenyl)				
CAS 74472-45-0	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20016400	PCB No. 164		5mg	
DRE-L20016400IO	PCB No. 164 10 µg/mL in Isooctane		10ml	
PCB 166 (2,3,4,4',5,6-Hexachlorobiphenyl)				
CAS 41411-63-6	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-E20016600	PCB No. 166		5mg	
DRE-L20016600IO	PCB No. 166 10 µg/mL in Isooctane		10ml	
PCB 167 (2,3',4,4',5,5'-Hexachlorobiphenyl)				
CAS 52663-72-6	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20016700	PCB No. 167		10mg	
DRE-L20016700IO	PCB No. 167 10 µg/mL in Isooctane		10ml	
PCB 168 (2,3',4,4',5',6-Hexachlorobiphenyl)				
CAS 59291-65-5	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-L20016800IO	PCB No. 168 10 µg/mL in Isooctane		10ml	
PCB 169 (3,3',4,4',5,5'-Hexachlorobiphenyl)				
CAS 32774-16-6	MW 360.8782	C ₁₂ H ₄ Cl ₆		
DRE-C20016900	PCB No. 169		5mg	
DRE-L20016900IO	PCB No. 169 10 µg/mL in Isooctane		10ml	
PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl)				
CAS 35065-30-6	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-C20017000	PCB No. 170		5mg	
DRE-L20017000IO	PCB No. 170 10 µg/mL in Isooctane		10ml	
PCB 171 (2,2',3,3',4,4',6-Heptachlorobiphenyl)				
CAS 52663-71-5	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20017100IO	PCB No. 171 10 µg/mL in Isooctane		10ml	
PCB 172 (2,2',3,3',4,5,5'-Heptachlorobiphenyl)				
CAS 52663-74-8	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20017200IO	PCB No. 172 10 µg/mL in Isooctane		10ml	










PCB's and related compounds

Product code	Description			
PCB 173 (2,2',3,3',4,5,6-Heptachlorobiphenyl)				
CAS 68194-16-1	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20017300	PCB No. 173		5mg	
DRE-L20017300IO	PCB No. 173 10 µg/mL in Isooctane		10ml	
PCB 174 (2,2',3,3',4,5,6'-Heptachlorobiphenyl)				
CAS 38411-25-5	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20017400IO	PCB No. 174 10 µg/mL in Isooctane		10ml	
PCB 175 (2,2',3,3',4,5',6-Heptachlorobiphenyl)				
CAS 40186-70-7	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20017500	PCB No. 175		5mg	
PCB 177 (2,2',3,3',4,5',6'-Heptachlorobiphenyl)				
CAS 52663-70-4	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20017700	PCB No. 177		5mg	
DRE-L20017700IO	PCB No. 177 10 µg/mL in Isooctane		10ml	
PCB 178 (2,2',3,3',5,5',6-Heptachlorobiphenyl)				
CAS 52663-67-9	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20017800IO	PCB No. 178 10 µg/mL in Isooctane		10ml	
PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)				
CAS 35065-29-3	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-C20018000	PCB No. 180		5mg	
DRE-L20018000IO	PCB No. 180 10 µg/mL in Isooctane(‡)		10ml	
PCB 181 (2,2',3,4,4',5,6-Heptachlorobiphenyl)				
CAS 74472-47-2	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20018100IO	PCB No. 181 10 µg/mL in Isooctane		10ml	
PCB 183 (2,2',3,4,4',5',6-Heptachlorobiphenyl)				
CAS 52663-69-1	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20018300	PCB No. 183		5mg	
DRE-L20018300IO	PCB No. 183 10 µg/mL in Isooctane		10ml	
PCB 185 (2,2',3,4,5,5',6-Heptachlorobiphenyl)				
CAS 52712-05-7	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20018500	PCB No. 185		5mg	
DRE-L20018500IO	PCB No. 185 10 µg/mL in Isooctane		10ml	









PCB's and related compounds

Product code	Description			
PCB 187 (2,2'',3,4'',5,5'',6-Heptachlorobiphenyl)				
CAS 52663-68-0	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-C20018700	PCB No. 187		10mg	
DRE-L20018700IO	PCB No. 187 10 µg/mL in Isooctane		10ml	
PCB 189 (2,3,3',4,4',5,5'-Heptachlorobiphenyl)				
CAS 39635-31-9	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-C20018900	PCB No. 189		5mg	
DRE-L20018900IO	PCB No. 189 10 µg/mL in Isooctane(‡)		10ml	
PCB 190 (2,3,3',4,4',5,6-Heptachlorobiphenyl)				
CAS 41411-64-7	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20019000IO	PCB No. 190 10 µg/mL in Isooctane		10ml	
PCB 192 (2,3,3',4,5,5',6-Heptachlorobiphenyl)				
CAS 74472-51-8	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-L20019200IO	PCB No. 192 10 µg/mL in Isooctane		10ml	
PCB 193 (2,3,3',4',5,5',6-Heptachlorobiphenyl)				
CAS 69782-91-8	MW 395.3232	C ₁₂ H ₃ Cl ₇		
DRE-E20019300	PCB No. 193		5mg	
DRE-L20019300IO	PCB No. 193 10 µg/mL in Isooctane		10ml	
PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)				
CAS 35694-08-7	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-C20019400	PCB No. 194(‡)		5mg	
DRE-L20019400IO	PCB No. 194 10 µg/mL in Isooctane		10ml	
PCB 195 (2,2',3,3',4,4',5,6-Octachlorobiphenyl)				
CAS 52663-78-2	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-C20019500	PCB No. 195		5mg	
DRE-L20019500IO	PCB No. 195 10 µg/mL in Isooctane		10ml	
PCB 196 (2,2',3,3',4,4',5,6'-Octachlorobiphenyl)				
CAS 42740-50-1	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-L20019600IO	PCB No. 196 10 µg/mL in Isooctane		10ml	
PCB 198 (2,2',3,3',4,5,5',6-Octachlorobiphenyl)				
CAS 68194-17-2	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-C20019800	PCB No. 198		5mg	
DRE-L20019800IO	PCB No. 198 10 µg/mL in Isooctane(‡)		10ml	

PCB's and related compounds

Product code	Description			
PCB 199 (2,2',3,3',4,5,5',6'-Octachlorobiphenyl)				
CAS 52663-75-9	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-L20019900IO	PCB No. 199 10 µg/mL in Isooctane		10ml	
PCB 200 (2,2',3,3',4,5,6,6'-Octachlorobiphenyl)				
CAS 52663-73-7	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-E20020000	PCB No. 200		5mg	
DRE-L20020000IO	PCB No. 200 10 µg/mL in Isooctane		10ml	
PCB 202 (2,2',3,3',5,5',6,6'-Octachlorobiphenyl)				
CAS 2136-99-4	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-L20020200IO	PCB No. 202 10 µg/mL in Isooctane		10ml	
PCB 203 (2,2',3,4,4',5,5',6-Octachlorobiphenyl)				
CAS 52663-76-0	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-L20020300IO	PCB No. 203 10 µg/mL in Isooctane		10ml	
PCB 204 (2,2',3,4,4',5,6,6'-Octachlorobiphenyl)				
CAS 74472-52-9	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-E20020400	PCB No. 204		5mg	
DRE-L20020400IO	PCB No. 204 10 µg/mL in Isooctane		10ml	
PCB 205 (2,3,3',4,4',5,5',6-Octachlorobiphenyl)				
CAS 74472-53-0	MW 429.7683	C ₁₂ H ₂ Cl ₈		
DRE-L20020500IO	PCB No. 205 10 µg/mL in Isooctane		10ml	
PCB 206 (2,2',3,3',4,4',5,5',6-Nonachlorobiphenyl)				
CAS 40186-72-9	MW 464.2133	C ₁₂ HCl ₉		
DRE-E20020600	PCB No. 206		5mg	
DRE-L20020600IO	PCB No. 206 10 µg/mL in Isooctane		10ml	
PCB 207 (2,2',3,3',4,4',5,6,6'-Nonachlorobiphenyl)				
CAS 52663-79-3	MW 464.2133	C ₁₂ HCl ₉		
DRE-E20020700	PCB No. 207		5mg	
DRE-L20020700IO	PCB No. 207 10 µg/mL in Isooctane		10ml	
PCB 208 (2,2',3,3',4,5,5',6,6'-Nonachlorobiphenyl)				
CAS 52663-77-1	MW 464.2133	C ₁₂ HCl ₉		
DRE-L20020800IO	PCB No. 208 10 µg/mL in Isooctane		10ml	

PCB's and related compounds

Product code	Description			
PCB 209 (Decachlorobiphenyl)				
CAS 2051-24-3	MW 498.6584	C ₁₂ Cl ₁₀		
DRE-C20020900	PCB No. 209(±)		10mg	
DRE-L20020900AL	PCB No. 209 10 µg/mL in Acetonitrile		10ml	
DRE-L20020900CY	PCB No. 209 10 µg/mL in Cyclohexane(±)		10ml	
DRE-L20020900IO	PCB No. 209 10 µg/mL in Isooctane		10ml	
DRE-X20020900CY	PCB No. 209 100 µg/mL in Cyclohexane(±)		10ml	
DRE-X20020900IO	PCB No. 209 100 µg/mL in Isooctane(±)		10ml	
DRE-XA20020900IO	PCB No. 209 100 µg/mL in Isooctane		1ml	
3,3',3'',4,4''-Pentachloro-m-terphenyl				
CAS 1064187-31-0	MW 402.5291	C ₁₈ H ₉ Cl ₅		
DRE-LA20386443HE	3,3',3'',4,4''-Pentachloro-m-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3',3'',5,5''-Pentachloro-m-terphenyl				
CAS n/a	MW 402.5291	C ₁₈ H ₉ Cl ₅		
DRE-LA20386553HE	3,3',3'',5,5''-Pentachloro-m-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3'',4,4''-Tetrachloro-o-terphenyl				
CAS n/a	MW 368.084	C ₁₈ H ₁₀ Cl ₄		
DRE-LA20383442HE	3,3'',4,4''-Tetrachloro-o-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3'',4,4''-Tetrachloro-p-terphenyl				
CAS 1957180-49-2	MW 368.084	C ₁₈ H ₁₀ Cl ₄		
DRE-LA20383444HE	3,3'',4,4''-Tetrachloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3'',5,5''-Tetrachloro-o-terphenyl				
CAS 566197-75-9	MW 368.084	C ₁₈ H ₁₀ Cl ₄		
DRE-LA20383552HE	3,3'',5,5''-Tetrachloro-o-terphenyl 10 µg/mL in n-Hexane		1ml	
3,3'',5,5''-Tetrachloro-p-terphenyl				
CAS 500729-81-7	MW 368.084	C ₁₈ H ₁₀ Cl ₄		
DRE-LA20383554HE	3,3'',5,5''-Tetrachloro-p-terphenyl 10 µg/mL in n-Hexane		1ml	
3',4,4''-Trichloro-m-terphenyl				
CAS n/a	MW 333.6389	C ₁₈ H ₁₁ Cl ₃		
DRE-LA20313443HE	3',4,4''-Trichloro-m-terphenyl 10 µg/mL in n-Hexane		1ml	

PCB's and related compounds

Product code	Description		
Ugilec 141			
CAS 111483-93-3	MW n/a		
DRE-L20434100TO	Ugilec 141 10 µg/mL in Toluene	10ml	No Structure
DRE-X20434100TO	Ugilec 141 100 µg/mL in Toluene	10ml	
Aroclor 1016 + 1260 Mix			
DRE-YS09000049HE	Aroclor 1016 + 1260 Mix 1000 µg/mL in Hexane(‡)		5x1ml
	Aroclor 1016		
	Aroclor 1260		
PCB-Mix 1			
DRE-L20030100AL	PCB-Mix 1 10 µg/mL in Acetonitrile		10ml
DRE-L20030100CY	PCB-Mix 1 10 µg/mL in Cyclohexane(‡)		10ml
DRE-L20030100IO	PCB-Mix 1 10 µg/mL in Isooctane(‡)		10ml
DRE-LA20030100CY	PCB-Mix 1 10 µg/mL in Cyclohexane		1ml
	PCB 28 (2,4,4'-Trichlorobiphenyl)		
	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)		
	PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)		
	PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)		
	PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)		
	PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)		
PCB-Mix 2			
DRE-L20030200CY	PCB-Mix 2 10 µg/mL in Cyclohexane		10ml
	PCB 18 (2,2',5-Trichlorobiphenyl)		
	PCB 28 (2,4,4'-Trichlorobiphenyl)		
	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)		
	PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)		
	PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)		
	PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)		
	PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)		
PCB-Mix 3			
DRE-L20030300AL	PCB-Mix 3 10 µg/mL in Acetonitrile(‡)		10ml
DRE-LA20030300CY	PCB-Mix 3 10 µg/mL in Cyclohexane(‡)		1ml
DRE-L20030300CY	PCB-Mix 3 10 µg/mL in Cyclohexane(‡)		10ml
DRE-LA20030300IO	PCB-Mix 3 10 µg/mL in Isooctane		1ml
DRE-L20030300IO	PCB-Mix 3 10 µg/mL in Isooctane(‡)		10ml
DRE-X20030300IO	PCB-Mix 3 100 µg/mL in Isooctane(‡)		10ml
DRE-X20030300ME	PCB-Mix 3 100 µg/mL in Methanol		10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl)		
	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)		
	PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)		
	PCB 118 (2,3',4,4',5-Pentachlorobiphenyl)		
	PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)		
	PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)		
	PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)		
PCB-Mix 4			
DRE-L20030400AL	PCB-Mix 4 10 µg/mL in Acetonitrile		10ml
DRE-L20030400IO	PCB-Mix 4 10 µg/mL in Isooctane		10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl)		
	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)		
	PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)		
	PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)		
	PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)		
	PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)		
	PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)		

PCB's and related compounds

Product code	Description	
PCB-Mix 7		
DRE-LA20030700IO	PCB-Mix 7 10 µg/mL in Isooctane	1ml
	PCB 8 (2,4'-Dichlorobiphenyl) PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 44 (2,2',3,5'-Tetrachlorobiphenyl) PCB 70 (2,3',4',5'-Tetrachlorobiphenyl) PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl) PCB 151 (2,2',3,5,5',6'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl) PCB 195 (2,2',3,3',4,4',5,6-Octachlorobiphenyl)	PCB 18 (2,2',5-Trichlorobiphenyl) PCB 31 (2,4',5-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)
PCB-Mix 8		
DRE-L20030800IO	PCB-Mix 8 10 µg/mL in Isooctane	10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 31 (2,4',5-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	
PCB-Mix 12		
DRE-L20031200IO	PCB-Mix 12 10 µg/mL in Isooctane	10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	PCB 31 (2,4',5-Trichlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 209 (2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl)
PCB-Mix 19		
DRE-L20031900IO	PCB-Mix 19 10 µg/mL in Isooctane(‡)	10ml
DRE-LA20031900IO	PCB-Mix 19 10 µg/mL in Isooctane(‡)	1ml
	PCB 18 (2,2',5-Trichlorobiphenyl) PCB 31 (2,4',5-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 149 (2,2',3,4',5',6-Hexachlorobiphenyl) PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl) PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 44 (2,2',3,5'-Tetrachlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl) PCB 209 (2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl)
PCB-Mix 20		
DRE-L20032000IO	PCB-Mix 20 10 µg/mL in Isooctane(‡)	10ml
DRE-LA20032000IO	PCB-Mix 20 10 µg/mL in Isooctane	1ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 128 (2,2',3,3',4,4'-Hexachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 169 (3,3',4,4',5,5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	PCB 31 (2,4',5-Trichlorobiphenyl) PCB 77 (3,3',4,4'-Tetrachlorobiphenyl) PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl) PCB 126 (3,3',4,4',5-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl) PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl)
PCB-Mix 21		
DRE-L20032100CY	PCB-Mix 21 10 µg/mL in Cyclohexane(‡)	10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)	PCB 31 (2,4',5-Trichlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)
PCB-Mix 24		
DRE-L20032400IO	PCB-Mix 24 10 µg/mL in Isooctane	10ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl)










PCB's and related compounds

Product code	Description	
	PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	
PCB-Mix 26		
DRE-L20032600CY	PCB-Mix 26 in Cyclohexane	10ml
	PCB 30 (2,4,6-Trichlorobiphenyl) [300 µg/mL] PCB 209 (2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl) [100 µg/mL]	
PCB-Mix 32		
DRE-LA20033200AC	PCB-Mix 32 10 µg/mL in Acetone(‡)	1.1ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)
PCB-Mix 37		
DRE-LA20033700IO	PCB-Mix 37 10 µg/mL in Isooctane	1ml
	PCB 28 (2,4,4'-Trichlorobiphenyl) PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl) PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl) PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl) PCB 183 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	PCB 52 (2,2',5,5'-Tetrachlorobiphenyl) PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl) PCB 146 (2,2',3,4',5,5'-Hexachlorobiphenyl) PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl) PCB 187 (2,2',3,4',5,5',6-Heptachlorobiphenyl)
	PCB 95 (2,2',3,5',6-Pentachlorobiphenyl) PCB 110 (2,3,3',4',6-Pentachlorobiphenyl) PCB 149 (2,2',3,4',5,6-Hexachlorobiphenyl) PCB 177 (2,2',3,3',4',5,6-Heptachlorobiphenyl)	PCB 99 (2,2',4,4',5-Pentachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 151 (2,2',3,5,5',6-Hexachlorobiphenyl) PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)
PCB-Mix 41		
DRE-LA20034100IO	PCB-Mix 41 10 µg/mL in Isooctane(‡)	1ml
	PCB 77 (3,3',4,4'-Tetrachlorobiphenyl) PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl) PCB 118 (2,3',4,4',5-Pentachlorobiphenyl) PCB 126 (3,3',4,4',5-Pentachlorobiphenyl) PCB 157 (2,3,3',4,4',5'-Hexachlorobiphenyl) PCB 169 (3,3',4,4',5,5'-Hexachlorobiphenyl)	PCB 81 (3,4,4',5-Tetrachlorobiphenyl) PCB 114 (2,3,4,4',5-Pentachlorobiphenyl) PCB 123 (2',3,4,4',5-Pentachlorobiphenyl) PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl) PCB 167 (2,3',4,4',5,5'-Hexachlorobiphenyl) PCB 189 (2,3,3',4,4',5,5'-Heptachlorobiphenyl)
PCB Mix, 7-3429		
DRE-GA09000201EA	PCB Mix, 7-3429 100 µg/mL in Ethyl Acetate(‡)	1ml
DRE-GX09000201AC	EPA 1664 LCS Mix in PFA Tubes 2000 µg/mL in Acetone(‡)	20x10ml
DRE-GX09000201AC-S	EPA 1664 LCS Mix in PFA Tubes 2000 µg/mL in Acetone(‡)	4x10ml
	n-hexadecane (C16) stearic acid	
Pesticide/PCB Surrogate Mix		
DRE-GS09000055AC	Pesticide/PCB Surrogate Mix 200 µg/mL in Acetone(‡)	10x1ml
	decachlorobiphenyl (BZ# 209) 2,4,5,6-tetrachloro-m-xylene	










POLYCYCLIC AROMATIC HYDROCARBONS (PAHS)











Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Acenaphthene				
CAS 83-32-9	MW 154.2078	C ₁₂ H ₁₀		
DRE-C20505000	Acenaphthene(‡)		100mg	
DRE-L20505000AL	Acenaphthene 10 µg/mL in Acetonitrile		10ml	
DRE-XA20505000AL	Acenaphthene 100 µg/mL in Acetonitrile		1ml	
Acenaphthene D10				
CAS 15067-26-2	MW 164.2694	C ₁₂ ² H ₁₀		
DRE-C20505100	Acenaphthene D10		100mg	
DRE-L20505100CY	Acenaphthene D10 10 µg/mL in Cyclohexane		10ml	
DRE-YA20505100TO	Acenaphthene D10 2000 µg/mL in Toluene		1ml	
Acenaphthylene				
CAS 208-96-8	MW 152.1919	C ₁₂ H ₈		
DRE-C20510000	Acenaphthylene(‡)		100mg	
DRE-L20510000AL	Acenaphthylene 10 µg/mL in Acetonitrile		10ml	
DRE-L20510000CY	Acenaphthylene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20510000AL	Acenaphthylene 100 µg/mL in Acetonitrile		1ml	
Acenaphthylene D8				
CAS 93951-97-4	MW 160.2412	C ₁₂ ² H ₈		
DRE-C20510100	Acenaphthylene D8		100mg	
DRE-L20510100CY	Acenaphthylene D8 10 µg/mL in Cyclohexane		10ml	
2-Acetamidofluorene				
CAS 53-96-3	MW 223.2698	C ₁₅ H ₁₃ NO		
DRE-C10012000	2-Acetamidofluorene		100mg	
Acridine				
CAS 260-94-6	MW 179.2173	C ₁₃ H ₉ N		
DRE-C20511000	Acridine		10mg	
2-Aminoanthracene				
CAS 613-13-8	MW 193.2438	C ₁₄ H ₁₁ N		
DRE-L20982200CY	2-Aminoanthracene 10 µg/mL in Cyclohexane		10ml	
1-Aminoanthraquinone				
CAS 82-45-1	MW 223.2268	C ₁₄ H ₉ NO ₂		
DRE-L20982800CY	1-Aminoanthraquinone 10 µg/mL in Cyclohexane		10ml	
6-Aminochrysene				
CAS 2642-98-0	MW 243.3025	C ₁₈ H ₁₃ N		
DRE-L20984600AL	6-Aminochrysene 10 µg/mL in Acetonitrile		10ml	









Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
2-Amino-6-ethoxynaphthalene				
CAS 293733-21-8 DRE-C10202350	MW 187.2377 2-Amino-6-ethoxynaphthalene	C ₁₂ H ₁₃ NO	10mg	
2-Aminofluorene				
CAS 153-78-6 DRE-L20984900AL	MW 181.2331 2-Aminofluorene 10 µg/mL in Acetonitrile	C ₁₃ H ₁₁ N	10ml	
2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline				
CAS 76180-96-6 DRE-L10204950AL	MW 198.2239 2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline 10 µg/mL in Acetonitrile	C ₁₁ H ₁₀ N ₄	10ml	
1-Aminonaphthalene				
CAS 134-32-7 DRE-C10206350	MW 143.1852 1-Aminonaphthalene	C ₁₀ H ₉ N	50mg	
1-Aminonaphthalene D7				
CAS 78832-53-8 DRE-XA10206351ME	MW 150.2283 1-Aminonaphthalene D7 100 µg/mL in Methanol	C ₁₀ ² H ₇ H ₂ N	1ml	
2-Aminonaphthalene				
CAS 91-59-8 DRE-C10206355 DRE-L10206355AL DRE-L10206355ME DRE-YS09010005DI	MW 143.1852 2-Aminonaphthalene(‡) 2-Aminonaphthalene 10 µg/mL in Acetonitrile 2-Aminonaphthalene 10 µg/mL in Methanol 2-Aminonaphthalene 2000 µg/mL in Dichloromethane(‡)	C ₁₀ H ₉ N	10mg 10ml 10ml 5x1ml	
2-Aminonaphthalene D7				
CAS 93951-94-1 DRE-XA10206356ME	MW 150.2283 2-Aminonaphthalene D7 100 µg/mL in Methanol	C ₁₀ ² H ₇ H ₂ N	1ml	
4-Aminoquinaldine				
CAS 6628-04-2 DRE-C10225000	MW 158.1998 4-Aminoquinaldine	C ₁₀ H ₁₀ N ₂	100mg	
Anthanthrene				
CAS 191-26-4 DRE-E20515000 DRE-L20515000AL DRE-L20515000CY DRE-XA20515000AL	MW 276.3307 Anthanthrene Anthanthrene 10 µg/mL in Acetonitrile Anthanthrene 10 µg/mL in Cyclohexane Anthanthrene 100 µg/mL in Acetonitrile	C ₂₂ H ₁₂	10mg 10ml 10ml 1ml	









Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Anthracene				
CAS 120-12-7	MW 178.2292	C ₁₄ H ₁₀		
DRE-C20520000	Anthracene(±)		50mg	
DRE-L20520000AL	Anthracene 10 µg/mL in Acetonitrile		10ml	
DRE-L20520000CY	Anthracene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20520000AL	Anthracene 100 µg/mL in Acetonitrile		1ml	
Anthracene D10				
CAS 1719-06-8	MW 188.2908	C ₁₄ ² H ₁₀		
DRE-C20520100	Anthracene D10		100mg	
DRE-L20520100CY	Anthracene D10 10 µg/mL in Cyclohexane		10ml	
DRE-XA20520100CY	Anthracene D10 100 µg/mL in Cyclohexane		1ml	
DRE-YA20520100MB	Anthracene D10 2000 µg/mL in Methyl-tert. butyl ether		1ml	
Benz(c)acridine				
CAS 225-51-4	MW 229.2759	C ₁₇ H ₁₁ N		
DRE-E20538200	Benz(c)acridine		10mg	
Benz[a]anthracene				
CAS 56-55-3	MW 228.2879	C ₁₈ H ₁₂		
DRE-C20545000	Benzo(a)anthracene		20mg	
DRE-L20545000AL	Benzo(a)anthracene 10 µg/mL in Acetonitrile		10ml	
DRE-L20545000CY	Benzo(a)anthracene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20545000AL	Benzo(a)anthracene 100 µg/mL in Acetonitrile		1ml	
Benzimidazole				
CAS 51-17-2	MW 118.1359	C ₇ H ₆ N ₂		
DRE-C10536500	Benzimidazole		250mg	
Benzo[a]anthracene D12				
CAS 1718-53-2	MW 240.3618	C ₁₈ ² H ₁₂		
DRE-C20545100	Benzo(a)anthracene D12		100mg	
DRE-L20545100AC	Benzo(a)anthracene D12 10 µg/mL in Acetone		10ml	
DRE-L20545100CY	Benzo(a)anthracene D12 10 µg/mL in Cyclohexane(±)		10ml	
Benzo[a]fluoranthene				
CAS 203-33-8	MW 252.3093	C ₂₀ H ₁₂		
DRE-L20560000AL	Benzo(a)fluoranthene 10 µg/mL in Acetonitrile		10ml	
DRE-L20560000CY	Benzo(a)fluoranthene 10 µg/mL in Cyclohexane		10ml	
11H-Benzo[a]fluorene (1,2-Benzofluorene)				
CAS 238-84-6	MW 216.2772	C ₁₇ H ₁₂		
DRE-L20585000AL	Benzo(a)fluorene (1,2-Benzofluorene) 10 µg/mL in Acetonitrile		10ml	









Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Benzo[a]pyrene				
CAS 50-32-8	MW 252.3093	C ₂₀ H ₁₂		
DRE-C20635000	Benzo(a)pyrene(‡)		10mg	
DRE-L20635000AL	Benzo(a)pyrene 10 µg/mL in Acetonitrile(‡)		10ml	
DRE-XA20635000AL	Benzo(a)pyrene 100 µg/mL in Acetonitrile(‡)		1ml	
DRE-XA20635000CY	Benzo(a)pyrene 100 µg/mL in Cyclohexane		1ml	
Benzo[a]pyrene D12				
CAS 63466-71-7	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20635100	Benzo(a)pyrene D12(‡)		10mg	
DRE-L20635100CY	Benzo(a)pyrene D12 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-LA20635100AL	Benzo(a)pyrene D12 10 µg/mL in Acetonitrile		1ml	
DRE-LA20635100CY	Benzo(a)pyrene D12 10 µg/mL in Cyclohexane		1ml	
DRE-XA20635100CY	Benzo(a)pyrene D12 100 µg/mL in Cyclohexane		1ml	
Benzo[b]chrysene				
CAS 214-17-5	MW 278.3466	C ₂₂ H ₁₄		
DRE-C20550000	Benzo(b)chrysene(‡)		10mg	
DRE-L20550000AL	Benzo(b)chrysene 10 µg/mL in Acetonitrile(‡)		10ml	
DRE-L20550000CY	Benzo(b)chrysene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20550000TO	Benzo(b)chrysene 100 µg/mL in Toluene		1ml	
Benzo[b]fluoranthene				
CAS 205-99-2	MW 252.3093	C ₂₀ H ₁₂		
DRE-C20565000	Benzo(b)fluoranthene(‡)		10mg	
DRE-L20565000AL	Benzo(b)fluoranthene 10 µg/mL in Acetonitrile		10ml	
DRE-L20565000CY	Benzo(b)fluoranthene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20565000AL	Benzo(b)fluoranthene 100 µg/mL in Acetonitrile		1ml	
Benzo[b]fluoranthene D12				
CAS 93951-98-5	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20565100	Benzo(b)fluoranthene D12		10mg	
DRE-LA20565100CY	Benzo(b)fluoranthene D12 10 µg/mL in Cyclohexane		1ml	
11H-Benzo[b]fluorene				
CAS 243-17-4	MW 216.2772	C ₁₇ H ₁₂		
DRE-C20590000	Benzo(b)fluorene		10mg	
DRE-L20590000AL	Benzo(b)fluorene (2,3-Benzfluorene) 10 µg/mL in Acetonitrile		10ml	
Benzo(b)furan				
CAS 271-89-6	MW 118.1326	C ₈ H ₆ O		
DRE-C20591500	Benzo(b)furan		25mg	
Benzo(b)naphtho(1,2-d)thiophene				
CAS 205-43-6	MW 234.3156	C ₁₆ H ₁₀ S		
DRE-L20595000CY	Benzo(b)naphtho(1,2-d)thiophene 10 µg/mL in Cyclohexane		10ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Benzo[b]naphtho[2,1-d]thiophene				
CAS 239-35-0	MW 234.3156	C ₁₆ H ₁₀ S		
DRE-C20600000	Benzo(b)naphtho(2,1-d)thiophene		10mg	
DRE-L20600000CY	Benzo(b)naphtho(2,1-d)thiophene 10 µg/mL in Cyclohexane		10ml	
7H-Benzo[c]fluorene				
CAS 205-12-9	MW 216.2772	C ₁₇ H ₁₂		
DRE-C20590400	7H-Benzo(c)fluorene		10mg	
DRE-L20590400CY	Benzo(c)fluorene 10 µg/mL in Cyclohexane		10ml	
Benzo[c]phenanthrene				
CAS 195-19-7	MW 228.2879	C ₁₈ H ₁₂		
DRE-C20631500	Benzo(c)phenanthrene		10mg	
DRE-E20631500	Benzo(c)phenanthrene		10mg	
DRE-L20631500CY	Benzo(c)phenanthrene 10 µg/mL in Cyclohexane		10ml	
Benzo[e]pyrene				
CAS 192-97-2	MW 252.3093	C ₂₀ H ₁₂		
DRE-C20645000	Benzo(e)pyrene(†)		10mg	
DRE-L20645000AL	Benzo(e)pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20645000CY	Benzo(e)pyrene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20645000AL	Benzo(e)pyrene 100 µg/mL in Acetonitrile		1ml	
DRE-XA20645000CY	Benzo(e)pyrene 100 µg/mL in Cyclohexane		1ml	
Benzo(e)pyrene D12				
CAS 205440-82-0	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-XA20645010CY	Benzo(e)pyrene D12 100 µg/mL in Cyclohexane		1ml	
Benzo[g,h,i]fluoranthene				
CAS 203-12-3	MW 226.272	C ₁₈ H ₁₀		
DRE-L20570000CY	Benzo(g,h,i)fluoranthene 10 µg/mL in Cyclohexane		10ml	
Benzo[g,h,i]perylene				
CAS 191-24-2	MW 276.3307	C ₂₂ H ₁₂		
DRE-C20630000	Benzo(g,h,i)perylene(†)		10mg	
DRE-L20630000AL	Benzo(g,h,i)perylene 10 µg/mL in Acetonitrile		10ml	
DRE-L20630000CY	Benzo(g,h,i)perylene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20630000AL	Benzo(g,h,i)perylene 100 µg/mL in Acetonitrile		1ml	
Benzo[g,h,i]perylene D12				
CAS 93951-66-7	MW 288.4046	C ₂₂ ² H ₁₂		
DRE-C20630200	Benzo(g,h,i)perylene D12		10mg	
DRE-LA20630200CY	Benzo(g,h,i)perylene D12 10 µg/mL in Cyclohexane		1ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description		
Benzo[j]fluoranthene			
CAS 205-82-3	MW 252.3093	C ₂₀ H ₁₂	
DRE-C20575000	Benzo(j)fluoranthene(‡)		
DRE-L20575000AL	Benzo(j)fluoranthene 10 µg/mL in Acetonitrile		
DRE-L20575000CY	Benzo(j)fluoranthene 10 µg/mL in Cyclohexane		
DRE-GA09010072DI	Benzo(j)fluoranthene 2000 µg/mL in Dichloromethane(‡)		
DRE-GS09010072DI	Benzo(j)fluoranthene 2000 µg/mL in Dichloromethane(‡)		
			10mg 10ml 10ml 1ml 5x1ml
Benzo[k]fluoranthene			
CAS 207-08-9	MW 252.3093	C ₂₀ H ₁₂	
DRE-C20580000	Benzo(k)fluoranthene(‡)		
DRE-L20580000AL	Benzo(k)fluoranthene 10 µg/mL in Acetonitrile		
DRE-L20580000CY	Benzo(k)fluoranthene 10 µg/mL in Cyclohexane		
DRE-XA20580000AL	Benzo(k)fluoranthene 100 µg/mL in Acetonitrile		
DRE-XA20580000CY	Benzo(k)fluoranthene 100 µg/mL in Cyclohexane		
			10mg 10ml 10ml 1ml 1ml
Benzo[k]fluoranthene D12			
CAS 93952-01-3	MW 264.3832	C ₂₀ ² H ₁₂	
DRE-C20580200	Benzo(k)fluoranthene D12		
DRE-LA20580200CY	Benzo(k)fluoranthene D12 10 µg/mL in Cyclohexane		
			10mg 1ml
1H-Benzotriazole			
CAS 95-14-7	MW 119.124	C ₆ H ₅ N ₃	
DRE-C10539500	1H-Benzotriazole		
			100mg
1,1'-Binaphthyl			
CAS 604-53-5	MW 254.3252	C ₂₀ H ₁₄	
DRE-C20655000	1,1'-Binaphthyl		
			10mg
2,2'-Binaphthyl			
CAS 612-78-2	MW 254.3252	C ₂₀ H ₁₄	
DRE-C20660000	2,2'-Binaphthyl		
DRE-L20660000AL	2,2'-Binaphthyl 10 µg/mL in Acetonitrile		
DRE-L20660000CY	2,2'-Binaphthyl 10 µg/mL in Cyclohexane		
			10mg 10ml 10ml
Carbazole			
CAS 86-74-8	MW 167.2066	C ₁₂ H ₉ N	
DRE-C10985000	Carbazole		
			100mg
Carbazole D8			
CAS 38537-24-5	MW 175.2559	C ₁₂ ² H ₈ HN	
DRE-XA10985100AC	Carbazole D8 100 µg/mL in Acetone		
			1.1ml










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
1-Chloronaphthalene				
CAS 90-13-1	MW 162.6156	C ₁₀ H ₇ Cl		
DRE-L20425100AL	1-Chloronaphthalene 10 µg/mL in Acetonitrile		10ml	
DRE-L20425100IO	1-Chloronaphthalene 10 µg/mL in Isooctane		10ml	
2-Chloronaphthalene				
CAS 91-58-7	MW 162.6156	C ₁₀ H ₇ Cl		
DRE-L20425200AL	2-Chloronaphthalene 10 µg/mL in Acetonitrile		10ml	
DRE-L20425200IO	2-Chloronaphthalene 10 µg/mL in Isooctane		10ml	
Chrysene				
CAS 218-01-9	MW 228.2879	C ₁₈ H ₁₂		
DRE-C20670000	Chrysene(‡)		25mg	
DRE-L20670000AL	Chrysene 10 µg/mL in Acetonitrile		10ml	
DRE-L20670000CY	Chrysene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20670000AL	Chrysene 100 µg/mL in Acetonitrile		1ml	
Chrysene D12				
CAS 1719-03-5	MW 240.3618	C ₁₈ ² H ₁₂		
DRE-C20670100	Chrysene D12		100mg	
DRE-L20670100AL	Chrysene D12 10 µg/mL in Acetonitrile		10ml	
DRE-L20670100CY	Chrysene D12 10 µg/mL in Cyclohexane		10ml	
Coronene				
CAS 191-07-1	MW 300.3521	C ₂₄ H ₁₂		
DRE-C20675000	Coronene		5mg	
DRE-L20675000AL	Coronene 10 µg/mL in Acetonitrile		10ml	
DRE-L20675000CY	Coronene 10 µg/mL in Cyclohexane		10ml	
Cyclopenta(c,d)pyrene				
CAS 27208-37-3	MW 226.272	C ₁₈ H ₁₀		
DRE-LA20680000AL	Cyclopenta(c,d)pyrene 10 µg/mL in Acetonitrile(‡)		1ml	
DRE-LA20680000CY	Cyclopenta(c,d)pyrene 10 µg/mL in Cyclohexane		1ml	
Decacyclene				
CAS 191-48-0	MW 450.5281	C ₃₆ H ₁₈		
DRE-C20685000	Decacyclene		25mg	
Dibenzo[a,c]anthracene				
CAS 215-58-7	MW 278.3466	C ₂₂ H ₁₄		
DRE-C20695000	Dibenzo(a,c)anthracene		10mg	
DRE-L20695000CY	Dibenzo(a,c)anthracene 10 µg/mL in Cyclohexane		10ml	
Dibenzo(a,e)fluoranthene				
CAS 5385-75-1	MW 302.368	C ₂₄ H ₁₄		
DRE-E20707000	Dibenzo(a,e)fluoranthene		10mg	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Dibenzo(a,e)pyrene				
CAS 192-65-4	MW 302.368	C ₂₄ H ₁₄		
DRE-C20715000	Dibenzo(a,e)pyrene		10mg	
DRE-L20715000AL	Dibenzo(a,e)pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20715000CY	Dibenzo(a,e)pyrene 10 µg/mL in Cyclohexane		10ml	
Dibenzo(a,h)acridine				
CAS 226-36-8	MW 279.3346	C ₂₁ H ₁₃ N		
DRE-L20694200CY	Dibenzo(a,h)acridine 10 µg/mL in Cyclohexane		10ml	
Dibenzo(a,h)anthracene				
CAS 53-70-3	MW 278.3466	C ₂₂ H ₁₄		
DRE-C20700000	Dibenzo(a,h)anthracene		10mg	
DRE-L20700000AL	Dibenzo(a,h)anthracene 10 µg/mL in Acetonitrile(‡)		10ml	
DRE-L20700000CY	Dibenzo(a,h)anthracene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20700000AL	Dibenzo(a,h)anthracene 100 µg/mL in Acetonitrile		1ml	
Dibenzo[a,h]anthracene D14				
CAS 13250-98-1	MW 292.4328	C ₂₂ ² H ₁₄		
DRE-C20700200	Dibenzo(a,h)anthracene D14		10mg	
DRE-L20700200CY	Dibenzo(a,h)anthracene D14 10 µg/mL in Cyclohexane(‡)		10ml	
Dibenzo[a,h]pyrene				
CAS 189-64-0	MW 302.368	C ₂₄ H ₁₄		
DRE-C20717000	Dibenzo[a,h]pyrene(‡)		10mg	
DRE-L20717000AL	Dibenzo(a,h)pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20717000CY	Dibenzo(a,h)pyrene 10 µg/mL in Cyclohexane		10ml	
Dibenzo(a,i)pyrene				
CAS 189-55-9	MW 302.368	C ₂₄ H ₁₄		
DRE-L20720000AL	Dibenzo(a,i)pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20720000IO	Dibenzo(a,i)pyrene 10 µg/mL in Isooctane		10ml	
Dibenzo(a,j)acridine				
CAS 224-42-0	MW 279.3346	C ₂₁ H ₁₃ N		
DRE-E20694600	Dibenzo(a,j)acridine		10mg	
Dibenzo(a,j)anthracene				
CAS 224-41-9	MW 278.3466	C ₂₂ H ₁₄		
DRE-L20705000AL	Dibenzo(a,j)anthracene 10 µg/mL in Acetonitrile		10ml	
Dibenzo(a,l)pyrene				
CAS 191-30-0	MW 302.368	C ₂₄ H ₁₄		
DRE-C20725000	Dibenzo(a,l)pyrene(‡)		10mg	
DRE-L20725000AL	Dibenzo(a,l)pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20725000CY	Dibenzo(a,l)pyrene 10 µg/mL in Cyclohexane		10ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Dibenzofuran				
CAS 132-64-9 DRE-C20710000	MW 168.1913 Dibenzofuran	C ₁₂ H ₈ O	10mg	
Dibenzofuran D8				
CAS 93952-04-6 DRE-L20710100CY	MW 176.2406 Dibenzofuran D8 10 µg/mL in Cyclohexane	C ₁₂ ² H ₆ O	10ml	
Dibenzothiophene				
CAS 132-65-0 DRE-C20727000 DRE-L20727000AL	MW 184.2569 Dibenzothiophene Dibenzothiophene 10 µg/mL in Acetonitrile	C ₁₂ H ₆ S	250mg 10ml	
Dibenzothiophene D8				
CAS 33262-29-2 DRE-L20727100CY	MW 192.3062 Dibenzothiophene D8 10 µg/mL in Cyclohexane	C ₁₂ ² H ₆ S	10ml	
1,4-Dibromonaphthalene				
CAS 83-53-4 DRE-C20431400	MW 285.9626 1,4-Dibromonaphthalene	C ₁₀ H ₆ Br ₂	100mg	
1,4-Dichloronaphthalene				
CAS 1825-31-6 DRE-L20421400AL DRE-L20421400IO	MW 197.0606 1,4-Dichloronaphthalene 10 µg/mL in Acetonitrile 1,4-Dichloronaphthalene 10 µg/mL in Isooctane	C ₁₀ H ₆ Cl ₂	10ml 10ml	
1,5-Dichloronaphthalene				
CAS 1825-30-5 DRE-LA20421500IO	MW 197.0606 1,5-Dichloronaphthalene 10 µg/mL in Isooctane	C ₁₀ H ₆ Cl ₂	1ml	
1,8-Dichloronaphthalene				
CAS 2050-74-0 DRE-LA20421800AL DRE-LA20421800IO	MW 197.0606 1,8-Dichloronaphthalene 10 µg/mL in Acetonitrile 1,8-Dichloronaphthalene 10 µg/mL in Isooctane	C ₁₀ H ₆ Cl ₂	1ml 1ml	
9,10-Dihydroanthracene				
CAS 613-31-0 DRE-C20730000	MW 180.2451 9,10-Dihydroanthracene	C ₁₄ H ₁₂	100mg	









Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description		
1,4-Dihydroxynaphthalene			
CAS 571-60-8 DRE-C12634840	MW 160.1693 1,4-Dihydroxynaphthalene	C ₁₀ H ₈ O ₂	100mg
			
7,12-Dimethylbenzo[a]anthracene			
CAS 57-97-6 DRE-C20745000	MW 256.341 7,12-Dimethylbenzo(a)anthracene	C ₂₀ H ₁₆	10mg
			
Dimethylnaphthalene (technical mixture)			
CAS 28804-88-8 DRE-L20780000CY	MW 156.2237 Dimethylnaphthalene (tech. Mixture of Isomers) 10 µg/mL in Cyclohexane	C ₁₀ H ₆ ·2CH ₃	10ml
			
1,2-Dimethylnaphthalene			
CAS 573-98-8 DRE-C20750000	MW 156.2237 1,2-Dimethylnaphthalene	C ₁₂ H ₁₂	50mg
			
1,3-Dimethylnaphthalene			
CAS 575-41-7 DRE-C20755000	MW 156.2237 1,3-Dimethylnaphthalene	C ₁₂ H ₁₂	50mg
			
1,4-Dimethylnaphthalene			
CAS 571-58-4 DRE-C20760000	MW 156.2237 1,4-Dimethylnaphthalene	C ₁₂ H ₁₂	50mg
			
1,5-Dimethylnaphthalene			
CAS 571-61-9 DRE-C20762000	MW 156.2237 1,5-Dimethylnaphthalene	C ₁₂ H ₁₂	50mg
			
1,6-Dimethylnaphthalene			
CAS 575-43-9 DRE-C20765000	MW 156.2237 1,6-Dimethylnaphthalene	C ₁₂ H ₁₂	10mg
			
1,7-Dimethylnaphthalene			
CAS 575-37-1 DRE-C20767000	MW 156.2237 1,7-Dimethylnaphthalene	C ₁₂ H ₁₂	10mg
			




Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
1,8-Dimethylnaphthalene				
CAS 569-41-5 DRE-LA20770000CY	MW 156.2237 1,8-Dimethylnaphthalene 10 µg/mL in Cyclohexane	C ₁₂ H ₁₂	1ml	
2,3-Dimethylnaphthalene				
CAS 581-40-8 DRE-C20772000	MW 156.2237 2,3-Dimethylnaphthalene	C ₁₂ H ₁₂	10mg	
2,6-Dimethylnaphthalene				
CAS 581-42-0 DRE-C20775000	MW 156.2237 2,6-Dimethylnaphthalene	C ₁₂ H ₁₂	10mg	
2,7-Dimethylnaphthalene				
CAS 582-16-1 DRE-C20775500	MW 156.2237 2,7-Dimethylnaphthalene	C ₁₂ H ₁₂	10mg	
3,6-Dimethylphenanthrene				
CAS 1576-67-6 DRE-C20785000 DRE-L20785000CY	MW 206.2824 3,6-Dimethylphenanthrene 3,6-Dimethylphenanthrene 10 µg/mL in Cyclohexane	C ₁₆ H ₁₄	10mg 10ml	
1,3-Dinitronaphthalene				
CAS 606-37-1 DRE-L20974300CY	MW 218.1656 1,3-Dinitronaphthalene 10 µg/mL in Cyclohexane	C ₁₀ H ₆ N ₂ O ₄	10ml	
1,5-Dinitronaphthalene				
CAS 605-71-0 DRE-L20974500CY	MW 218.1656 1,5-Dinitronaphthalene 10 µg/mL in Cyclohexane	C ₁₀ H ₆ N ₂ O ₄	10ml	
1,8-Dinitronaphthalene				
CAS 602-38-0 DRE-C20974800	MW 218.1656 1,8-Dinitronaphthalene	C ₁₀ H ₆ N ₂ O ₄	10mg	
1,3-Dinitropyrene				
CAS 75321-20-9 DRE-XA20975300TO	MW 292.2457 1,3-Dinitropyrene 100 µg/mL in Toluene	C ₁₆ H ₈ N ₂ O ₄	1ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
1,6-Dinitropyrene				
CAS 42397-64-8 DRE-L20975600TO	MW 292.2457 1,6-Dinitropyrene 10 µg/mL in Toluene	C ₁₆ H ₈ N ₂ O ₄	10ml	
1,8-Dinitropyrene				
CAS 42397-65-9 DRE-LA20975800AL	MW 292.2457 1,8-Dinitropyrene 10 µg/mL in Acetonitrile	C ₁₆ H ₈ N ₂ O ₄	1ml	
1-Ethynaphthalene				
CAS 1127-76-0 DRE-C20793100	MW 156.2237 1-Ethynaphthalene	C ₁₂ H ₁₂	10mg	
2-Ethynaphthalene				
CAS 939-27-5 DRE-C20793200	MW 156.2237 2-Ethynaphthalene	C ₁₂ H ₁₂	50mg	
Fluoranthene				
CAS 206-44-0 DRE-C20795000 DRE-L20795000AL DRE-L20795000CY DRE-XA20795000AL	MW 202.2506 Fluoranthene Fluoranthene 10 µg/mL in Acetonitrile Fluoranthene 10 µg/mL in Cyclohexane Fluoranthene 100 µg/mL in Acetonitrile	C ₁₆ H ₁₀	25mg 10ml 10ml 1ml	
Fluoranthene-D10				
CAS 93951-69-0 DRE-C20795100 DRE-L20795100AC DRE-L20795100ME DRE-XA20795100AL	MW 212.3122 Fluoranthene D10 Fluoranthene D10 10 µg/mL in Acetone Fluoranthene D10 10 µg/mL in Methanol Fluoranthene D10 100 µg/mL in Acetonitrile	C ₁₆ ² H ₁₀	50mg 10ml 10ml 1ml	
Fluorene				
CAS 86-73-7 DRE-C20800000 DRE-L20800000AL DRE-L20800000CY DRE-XA20800000AL	MW 166.2185 Fluorene(†) Fluorene 10 µg/mL in Acetonitrile Fluorene 10 µg/mL in Cyclohexane Fluorene 100 µg/mL in Acetonitrile	C ₁₃ H ₁₀	25mg 10ml 10ml 1ml	
Fluorene-D10				
CAS 81103-79-9 DRE-C20800200 DRE-L20800200CY	MW 176.2801 Fluorene D10 Fluorene D10 10 µg/mL in Cyclohexane	C ₁₃ ² H ₁₀	100mg 10ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
9-Fluorenone				
CAS 486-25-9 DRE-C20805000	MW 180.202 9-Fluorenone	C ₁₃ H ₈ O	250mg	
1-Fluoronaphthalene				
CAS 321-38-0 DRE-C13794000 DRE-YA13794000MB	MW 146.161 1-Fluoronaphthalene 1-Fluoronaphthalene 2000 µg/mL in Methyl-tert. butyl ether	C ₁₀ H ₇ F	500mg 1ml	
Halowax 1000				
CAS 58718-66-4 DRE-L20410000CY	MW n/a Halowax 1000 10 µg/mL in Cyclohexane		10ml	No Structure
Halowax 1001				
CAS 58718-67-5 DRE-L20410100CY	MW n/a Halowax 1001 10 µg/mL in Cyclohexane		10ml	No Structure
Halowax 1013				
CAS 12616-35-2 DRE-L20411300CY	MW n/a Halowax 1013 10 µg/mL in Cyclohexane		10ml	No Structure
Halowax 1014				
CAS 12616-36-3 DRE-L20411400CY	MW n/a Halowax 1014 10 µg/mL in Cyclohexane		10ml	No Structure
Halowax 1051				
CAS 70776-03-3 DRE-L20415100CY	MW n/a Halowax 1051 10 µg/mL in Cyclohexane		10ml	No Structure
Halowax 1099				
CAS 39450-05-0 DRE-L20419900CY	MW n/a Halowax 1099 10 µg/mL in Cyclohexane		10ml	No Structure
6-Hydroxychrysene				
CAS 37515-51-8 DRE-E20990600	MW 244.2873 6-Hydroxy-chrysene	C ₁₈ H ₁₂ O	10mg	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
1-Hydroxyphenanthrene				
CAS 2433-56-9	MW 194.2286	C ₁₄ H ₁₀ O		
DRE-L20992100AL	1-Hydroxy-phenanthrene 10 µg/mL in Acetonitrile (*)		10ml	
2-Hydroxyphenanthrene				
CAS 605-55-0	MW 194.2286	C ₁₄ H ₁₀ O		
DRE-E20992200	2-Hydroxy-phenanthrene		10mg	
DRE-L20992200AL	2-Hydroxy-phenanthrene 10 µg/mL in Acetonitrile		10ml	
3-Hydroxyphenanthrene				
CAS 605-87-8	MW 194.2286	C ₁₄ H ₁₀ O		
DRE-E20992300	3-Hydroxy-phenanthrene		10mg	
DRE-L20992300AL	3-Hydroxy-phenanthrene 10 µg/mL in Acetonitrile		10ml	
4-Hydroxy-phenanthrene				
CAS 7651-86-7	MW 194.2286	C ₁₄ H ₁₀ O		
DRE-L20992400AL	4-Hydroxy-phenanthrene 10 µg/mL in Acetonitrile		10ml	
9-Hydroxyphenanthrene				
CAS 484-17-3	MW 194.2286	C ₁₄ H ₁₀ O		
DRE-C20992900	9-Hydroxy-phenanthrene		10mg	
1-Hydroxypyrene				
CAS 5315-79-7	MW 218.25	C ₁₆ H ₁₀ O		
DRE-C20994100	1-Hydroxy-pyrene		10mg	
Indeno(1,2,3-c,d)fluoranthene				
CAS 193-43-1	MW 276.3307	C ₂₂ H ₁₂		
DRE-E20825000	Indeno(1,2,3-c,d)fluoranthene		10mg	
DRE-L20825000AL	Indeno(1,2,3-c,d)fluoranthene 10 µg/mL in Acetonitrile		10ml	
DRE-L20825000CY	Indeno(1,2,3-c,d)fluoranthene 10 µg/mL in Cyclohexane		10ml	
Indeno(1,2,3-c,d)pyrene				
CAS 193-39-5	MW 276.3307	C ₂₂ H ₁₂		
DRE-C20830000	Indeno(1,2,3-c,d)pyrene		10mg	
DRE-L20830000AL	Indeno(1,2,3-c,d)pyrene 10 µg/mL in Acetonitrile(‡)		10ml	
DRE-L20830000CY	Indeno(1,2,3-c,d)pyrene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20830000AL	Indeno(1,2,3-c,d)pyrene 100 µg/mL in Acetonitrile		1ml	
DRE-XA20830000CY	Indeno(1,2,3-c,d)pyrene 100 µg/mL in Cyclohexane		1ml	
Indeno(1,2,3-c,d)pyrene D12				
CAS 203578-33-0	MW 288.4046	C ₂₂ ² H ₁₂		
DRE-LA20830200CY	Indeno(1,2,3-c,d)pyrene D12 10 µg/mL in Cyclohexane		1ml	









Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Indole				
CAS 120-72-9 DRE-C20831000	MW 117.1479 Indole	C ₈ H ₇ N	10mg	
Isoquinoline				
CAS 119-65-3 DRE-C20833000	MW 129.1586 Isochinolin	C ₉ H ₇ N	25mg	
1-Methylantracene				
CAS 610-48-0 DRE-C20834900 DRE-L20834900CY	MW 192.2558 1-Methylantracene 1-Methylantracene 10 µg/mL in Cyclohexane	C ₁₅ H ₁₂	10mg 10ml	
2-Methylantracene				
CAS 613-12-7 DRE-C20835000	MW 192.2558 2-Methylantracene	C ₁₅ H ₁₂	10mg	
9-Methylantracene				
CAS 779-02-2 DRE-C20840000	MW 192.2558 9-Methylantracene	C ₁₅ H ₁₂	10mg	
2-Methylbenzofuran				
CAS 4265-25-2 DRE-C15083785	MW 132.1592 2-Methylbenzofuran	C ₉ H ₈ O	100mg	
5-Methylbenzothiophene				
CAS 14315-14-1 DRE-C20846000	MW 148.2248 5-Methylbenzo(b)thiophene	C ₉ H ₆ S	100mg	
3-Methylcholanthrene				
CAS 56-49-5 DRE-C20850000 DRE-L20850000CY	MW 268.3517 3-Methylcholanthrene 3-Methylcholanthrene 10 µg/mL in Cyclohexane	C ₂₁ H ₁₆	10mg 10ml	
1-Methylchrysene				
CAS 3351-28-8 DRE-E20865000	MW 242.3145 1-Methylchrysene	C ₁₉ H ₁₄	10mg	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
2-Methylchrysene				
CAS 3351-32-4	MW 242.3145	C ₁₉ H ₁₄		
DRE-E20870000	2-Methylchrysene		10mg	
DRE-L20870000AL	2-Methylchrysene 10 µg/mL in Acetonitrile		10ml	
DRE-L20870000CY	2-Methylchrysene 10 µg/mL in Cyclohexane		10ml	
4-Methylchrysene				
CAS 3351-30-2	MW 242.3145	C ₁₉ H ₁₄		
DRE-E20880000	4-Methylchrysene		10mg	
5-Methylchrysene				
CAS 3697-24-3	MW 242.3145	C ₁₉ H ₁₄		
DRE-L20885000AL	5-Methylchrysene 10 µg/mL in Acetonitrile		10ml	
DRE-L20885000CY	5-Methylchrysene 10 µg/mL in Cyclohexane		10ml	
6-Methylchrysene				
CAS 1705-85-7	MW 242.3145	C ₁₉ H ₁₄		
DRE-C20890000	6-Methylchrysene		10mg	
DRE-L20890000AL	6-Methylchrysene 10 µg/mL in Acetonitrile		10ml	
DRE-YS09010018DI	6-Methylchrysene 1000 µg/mL in Dichloromethane(‡)		5x1ml	
1-Methylfluoranthene				
CAS 25889-60-5	MW 216.2772	C ₁₇ H ₁₂		
DRE-L20892500AL	1-Methylfluoranthene 10 µg/mL in Acetonitrile		10ml	
DRE-L20892500CY	1-Methylfluoranthene 10 µg/mL in Cyclohexane		10ml	
2-Methylfluoranthene				
CAS 33543-31-6	MW 216.2772	C ₁₇ H ₁₂		
DRE-L20892600AL	2-Methylfluoranthene 10 µg/mL in Acetonitrile(‡)		10ml	
DRE-L20892600CY	2-Methylfluoranthene 10 µg/mL in Cyclohexane		10ml	
3-Methylindole				
CAS 83-34-1	MW 131.1745	C ₉ H ₉ N		
DRE-C20893000	3-Methylindole		10mg	
1-Methylnaphthalene				
CAS 90-12-0	MW 142.1971	C ₁₁ H ₁₀		
DRE-C20895000	1-Methylnaphthalene(‡)		50mg	
DRE-GA20895000ME	1-Methylnaphthalene 1000 µg/mL in Methanol		1ml	
DRE-GH20895000ME	1-Methylnaphthalene 1000 µg/mL in Methanol		5x1ml	
DRE-L20895000AL	1-Methylnaphthalene 10 µg/mL in Acetonitrile		10ml	
DRE-L20895000CY	1-Methylnaphthalene 10 µg/mL in Cyclohexane		10ml	
1-Methylnaphthalene D10				
CAS 38072-94-5	MW 152.2587	C ₁₁ ² H ₁₀		
DRE-C20895100	1-Methylnaphthalene D10		10mg	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
2-Methylnaphthalene				
CAS 91-57-6	MW 142.1971	C ₁₁ H ₁₀		
DRE-C20895200	2-Methylnaphthalene		50mg	
DRE-GA20895200DI	2-Methylnaphthalene 1000 µg/mL in Dichloromethane		1ml	
DRE-GH20895200DI	2-Methylnaphthalene 1000 µg/mL in Dichloromethane		5x1ml	
DRE-L20895200AL	2-Methylnaphthalene 10 µg/mL in Acetonitrile		10ml	
DRE-L20895200CY	2-Methylnaphthalene 10 µg/mL in Cyclohexane		10ml	
1-Methylphenanthrene				
CAS 832-69-9	MW 192.2558	C ₁₅ H ₁₂		
DRE-L20900000AL	1-Methylphenanthrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20900000CY	1-Methylphenanthrene 10 µg/mL in Cyclohexane		10ml	
2-Methylphenanthrene				
CAS 2531-84-2	MW 192.2558	C ₁₅ H ₁₂		
DRE-L20900100CY	2-Methylphenanthrene 10 µg/mL in Cyclohexane		10ml	
3-Methylphenanthrene				
CAS 832-71-3	MW 192.2558	C ₁₅ H ₁₂		
DRE-C20900200	3-Methylphenanthrene		10mg	
9-Methylphenanthrene				
CAS 883-20-5	MW 192.2558	C ₁₅ H ₁₂		
DRE-E20900400	9-Methylphenanthrene		10mg	
1-Methylpyrene				
CAS 2381-21-7	MW 216.2772	C ₁₇ H ₁₂		
DRE-C20901000	1-Methylpyrene		10mg	
DRE-L20901000CY	1-Methylpyrene 10 µg/mL in Cyclohexane		10ml	
2-Methylquinoline				
CAS 91-63-4	MW 143.1852	C ₁₀ H ₉ N		
DRE-C20848500	2-Methylquinoline		250mg	
Naphthalene				
CAS 91-20-3	MW 128.1705	C ₁₀ H ₈		
DRE-C20905000	Naphthalene		100mg	
DRE-L20905000AL	Naphthalene 10 µg/mL in Acetonitrile		10ml	
DRE-L20905000CY	Naphthalene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20905000AL	Naphthalene 100 µg/mL in Acetonitrile		1ml	










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Naphthalene D8				
CAS 1146-65-2	MW 136.2198	C ₁₀ H ₈		
DRE-C20905100	Naphthalene D8(‡)		100mg	
DRE-L20905100CY	Naphthalene D8 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-YA20905100MB	Naphthalene D8 2000 µg/mL in Methyl-tert. butyl ether		1.1ml	
Naphthalene-1-aldehyde				
CAS 66-77-3	MW 156.1806	C ₁₁ H ₈ O		
DRE-C15419800	Naphthalene-1-aldehyde		100mg	
5-Nitroacenaphthene				
CAS 602-87-9	MW 199.2054	C ₁₂ H ₉ NO ₂		
DRE-C20961800	5-Nitroacenaphthene		10mg	
7-Nitrobenzo[a]anthracene				
CAS 20268-51-3	MW 273.2854	C ₁₈ H ₁₁ NO ₂		
DRE-E20962600	7-Nitrobenz(a)anthracene		10mg	
6-Nitrobenzo[a]pyrene				
CAS 63041-90-7	MW 297.3068	C ₂₀ H ₁₁ NO ₂		
DRE-E20962800	6-Nitrobenz(a)pyrene		10mg	
3-Nitrofluoranthene				
CAS 892-21-7	MW 247.2482	C ₁₆ H ₉ NO ₂		
DRE-E20964700	3-Nitrofluoranthene		10mg	
2-Nitrofluorene				
CAS 607-57-8	MW 211.2161	C ₁₃ H ₉ NO ₂		
DRE-L20965000CY	2-Nitrofluorene 10 µg/mL in Cyclohexane		10ml	
1-Nitronaphthalene				
CAS 86-57-7	MW 173.1681	C ₁₀ H ₇ NO ₂		
DRE-L20965100CY	1-Nitronaphthalene 10 µg/mL in Cyclohexane		10ml	
2-Nitronaphthalene				
CAS 581-89-5	MW 173.1681	C ₁₀ H ₇ NO ₂		
DRE-E20965200	2-Nitronaphthalene		10mg	
DRE-L20965200CY	2-Nitronaphthalene 10 µg/mL in Cyclohexane		10ml	








Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
3-Nitroperylene				
CAS 20589-63-3 DRE-E20965700	MW 297.3068 3-Nitroperylene	C ₂₀ H ₁₁ NO ₂	10mg	
3-Nitrophenanthrene				
CAS 17024-19-0 DRE-L20966300CY	MW 223.2268 3-Nitrophenanthrene 10 µg/mL in Cyclohexane	C ₁₄ H ₉ NO ₂	10ml	
9-Nitrophenanthrene				
CAS 954-46-1 DRE-L20966600CY	MW 223.2268 9-Nitrophenanthrene 10 µg/mL in Cyclohexane	C ₁₄ H ₉ NO ₂	10ml	
1-Nitropyrene				
CAS 5522-43-0 DRE-C20967100	MW 247.2482 1-Nitropyrene	C ₁₆ H ₉ NO ₂	10mg	
4-Nitropyrene				
CAS 57835-92-4 DRE-E20967400	MW 247.2482 4-Nitropyrene	C ₁₆ H ₉ NO ₂	10mg	
4-Nitroquinoline N-Oxide				
CAS 56-57-5 DRE-C15558000	MW 190.1555 4-Nitroquinoline N-Oxide	C ₉ H ₆ N ₂ O ₃	100mg	
Octachloronaphthalene				
CAS 2234-13-1 DRE-C20425800 DRE-L20425800AL DRE-L20425800CY DRE-L20425800IO	MW 403.731 Octachloronaphthalene Octachloronaphthalene 10 µg/mL in Acetonitrile Octachloronaphthalene 10 µg/mL in Cyclohexane Octachloronaphthalene 10 µg/mL in Isooctane	C ₁₀ Cl ₈	5mg 10ml 10ml 10ml	
Octafluoronaphthalene				
CAS 313-72-4 DRE-C15710700 DRE-XA15710700AL	MW 272.0942 Octafluoronaphthalene Octafluoronaphthalene 100 µg/mL in Acetonitrile	C ₁₀ F ₈	100mg 1ml	
Perylene				
CAS 198-55-0 DRE-C20915000 DRE-L20915000AL DRE-L20915000CY	MW 252.3093 Perylene Perylene 10 µg/mL in Acetonitrile Perylene 10 µg/mL in Cyclohexane	C ₂₀ H ₁₂	10mg 10ml 10ml	

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
Perylene-d12				
CAS 1520-96-3	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20915100	Perylene D12(‡)		100mg	
DRE-L20915100CY	Perylene D12 10 µg/mL in Cyclohexane		10ml	
DRE-YA20915100TO	Perylene D12 2000 µg/mL in Toluene		1ml	
Phenanthrene				
CAS 85-01-8	MW 178.2292	C ₁₄ H ₁₀		
DRE-C20920000	Phenanthrene		50mg	
DRE-L20920000AL	Phenanthrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20920000CY	Phenanthrene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20920000AL	Phenanthrene 100 µg/mL in Acetonitrile		1ml	
Phenanthrene D10				
CAS 1517-22-2	MW 188.2908	C ₁₄ ² H ₁₀		
DRE-C20920100	Phenanthrene D10(‡)		100mg	
DRE-L20920100AC	Phenanthrene D10 10 µg/mL in Acetone		10ml	
DRE-L20920100CY	Phenanthrene D10 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-YA20920100MB	Phenanthrene D10 2000 µg/mL in Methyl-tert. butyl ether		1ml	
Phenazine				
CAS 92-82-0	MW 180.2053	C ₁₂ H ₈ N ₂		
DRE-C20921500	Phenazine		25mg	
9-Phenylanthracene				
CAS 602-55-1	MW 254.3252	C ₂₀ H ₁₄		
DRE-C20922500	9-Phenylanthracene		25mg	
1-Phenylnaphthalene				
CAS 605-02-7	MW 204.2665	C ₁₆ H ₁₂		
DRE-L20923000AL	1-Phenylnaphthalene 10 µg/mL in Acetonitrile		10ml	
2-Phenylnaphthalene				
CAS 612-94-2	MW 204.2665	C ₁₆ H ₁₂		
DRE-L20923100CY	2-Phenylnaphthalene 10 µg/mL in Cyclohexane		10ml	
Picene				
CAS 213-46-7	MW 278.3466	C ₂₂ H ₁₄		
DRE-L20925000CY	Picene 10 µg/mL in Cyclohexane		10ml	
Pyrene				
CAS 129-00-0	MW 202.2506	C ₁₆ H ₁₀		
DRE-C20930000	Pyrene(‡)		50mg	
DRE-L20930000AL	Pyrene 10 µg/mL in Acetonitrile		10ml	
DRE-L20930000CY	Pyrene 10 µg/mL in Cyclohexane		10ml	
DRE-XA20930000AL	Pyrene 100 µg/mL in Acetonitrile		1ml	

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description				
Pyrene-d10					
CAS 1718-52-1	MW 212.3122	C ₁₆ ² H ₁₀			
DRE-C20930100	Pyrene D10(±)				100mg
DRE-L20930100CY	Pyrene D10 10 µg/mL in Cyclohexane				10ml
DRE-XA20930100AC	Pyrene D10 100 µg/mL in Acetone				1ml
DRE-XA20930100AL	Pyrene D10 100 µg/mL in Acetonitrile				1ml
Retene					
CAS 483-65-8	MW 234.3355	C ₁₈ H ₁₈			
DRE-L16812000CY	Retene 10 µg/mL in Cyclohexane				10ml
1,2,3,4-Tetrachloronaphthalene					
CAS 20020-02-4	MW 265.9508	C ₁₀ H ₄ Cl ₄			
DRE-L17360000IO	1,2,3,4-Tetrachloronaphthalene 10 µg/mL in Isooctane				10ml
DRE-XA17360000CY	1,2,3,4-Tetrachloronaphthalene 100 µg/mL in Cyclohexane				1ml
1,2,3,4-Tetrahydronaphthalene					
CAS 119-64-2	MW 132.2023	C ₁₀ H ₁₂			
DRE-C20940000	1,2,3,4-Tetrahydronaphthalene				50mg
2,3,5-Trimethylnaphthalene					
CAS 2245-38-7	MW 170.2503	C ₁₃ H ₁₄			
DRE-L20943000CY	2,3,5-Trimethylnaphthalene 10 µg/mL in Cyclohexane				10ml
2,4,7-Trinitro-9-fluorenone					
CAS 129-79-3	MW 315.1947	C ₁₃ H ₅ N ₃ O ₇			
DRE-L20977000CY	2,4,7-Trinitro-9-fluorenone 10 µg/mL in Cyclohexane				10ml
Triphenylene					
CAS 217-59-4	MW 228.2879	C ₁₈ H ₁₂			
DRE-C20945000	Triphenylene				25mg
DRE-L20945000AL	Triphenylene 10 µg/mL in Acetonitrile				10ml
DRE-L20945000CY	Triphenylene 10 µg/mL in Cyclohexane				10ml
Deuterated Mix, 3-2135					
DRE-GS09000271TO	Deuterated Mix, 3-2135 in Toluene(±)			5x1ml	
	1-aminonaphthalene-d7 [50 mg/L] 2-aminonaphthalene-d7 [50 mg/L] 4-aminobiphenyl-d9 [25 mg/L]				
Deuterated PAH Mix, 4-1998					
DRE-GS09000189TO	Deuterated PAH Mix, 4-1998 10 µg/mL in Toluene(±)			5x1ml	
	benzo[a]anthracene-d12 benzo(a)pyrene-d12 chrysene-d12 benzo(b)fluoranthene-d12				

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description	
EPH MA Aromatics Mix		
DRE-YS09000044DI	EPH MA Aromatics Mix 1000 µg/mL in Dichloromethane(‡)	5x1ml
acenaphthene benzo[a]pyrene chrysene indeno[1,2,3-cd]pyrene pyrene	acenaphthylene benzo[b]fluoranthene dibenz[a,h]anthracene 2-methylnaphthalene	anthracene benzo[ghi]perylene fluoranthene naphthalene
		benzo[a]anthracene benzo[k]fluoranthene fluorene phenanthrene
Internal Standards Mix 26		
DRE-YA08272600BD	Internal Standards Mix 26, 4000 µg/mL in Benzene/Dichloromethane 1:1	1ml
	1,4-Dichlorobenzene D4 Acenaphthene D10 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10	
Internal Standards Mix 33		
DRE-Y08273300TO	Internal Standards Mix 33, 2000 µg/mL in Toluene(‡)	10ml
DRE-YA08273300TO	Internal Standards Mix 33 2000 µg/mL in Toluene(‡)	1ml
	1,4-Dichlorobenzene D4 Acenaphthene D10 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10	
Internal Standards Mix 37		
DRE-LA08273700IO	Internal Standards Mix 37 15 µg/mL in Isooctane	1ml
	Acenaphthene D10 Benzo[g,h,i]perylene D12 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10 Pyrene D10	
PAH-Mix 1		
DRE-L20950001AL	PAH-Mix 1, 2-10 µg/mL in Acetonitrile	10ml
DRE-L20950001CY	PAH-Mix 1, 2-10 µg/mL in Cyclohexane	10ml
	Benzo(a)pyrene [2 µg/mL] Benzo(ghi)perylene [2 µg/mL] Benzo(k)fluoranthene [2 µg/mL] Benzo[b]fluoranthene [2 µg/mL] Fluoranthene [10 µg/mL] Indeno(1,2,3-c,d)pyrene [2 µg/mL]	
PAH-Mix 3		
DRE-L20950003AL	PAH-Mix 3 in Acetonitrile	10ml
DRE-L20950003CY	PAH-Mix 3 in Cyclohexane	10ml
	Benzo(a)pyrene [20 µg/mL] Benzo(g,h,i)perylene [20 µg/mL] Benzo(k)fluoranthene [20 µg/mL] Benzo[b]fluoranthene [20 µg/mL] Fluoranthene [50 µg/mL] Indeno(1,2,3-c,d)pyrene [40 µg/mL]	
PAH-Mix 4		
DRE-GS09000177TO	PAH Mix 4-142 10 µg/mL in Toluene(‡)	5x1ml
	benzo[a]anthracene chrysene benzo[b]fluoranthene benzo[a]pyrene	

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description	
PAH-Mix 8		
DRE-L20950008AL	PAH-Mix 8 in Acetonitrile	10ml
	7,12-Dimethylbenzo(a)anthracene [10 µg/mL] Benzo(a)anthracene [10 µg/mL] Benzo(ghi)perylene [10 µg/mL] Benzo(b)fluoranthene [10 µg/mL] Dibenz[a,h]anthracene [10 µg/mL] Fluorene [10 µg/mL] Perylene [10 µg/mL]	Anthracene [50 µg/mL] Benzo(a)pyrene [10 µg/mL] Benzo(k)fluoranthene [10 µg/mL] Chrysene [50 µg/mL] Fluoranthene [50 µg/mL] Indeno(1,2,3-c,d)pyrene [10 µg/mL] Pyrene [50 µg/mL]
PAH-Mix 9		
DRE-L20950009AL	PAH-Mix 9 10 µg/mL in Acetonitrile(‡)	10ml
DRE-LS20950009AL	PAH-Mix 9 10 µg/mL in Acetonitrile	5x1ml
DRE-L20950009CY	PAH-Mix 9 10 µg/mL in Cyclohexane(‡)	10ml
DRE-LS20950009CY	PAH-Mix 9 10 µg/mL in Cyclohexane	5x1ml
DRE-XA20950009AL	PAH-Mix 9 100 µg/mL in Acetonitrile(‡)	1ml
DRE-X20950009AL	PAH-Mix 9 100 µg/mL in Acetonitrile(‡)	10ml
DRE-XA20950009CY	PAH-Mix 9 100 µg/mL in Cyclohexane(‡)	1ml
DRE-X20950009CY	PAH-Mix 9 100 µg/mL in Cyclohexane(‡)	10ml
Acenaphthene Benzo(a)pyrene Chrysene Indeno(1,2,3-c,d)pyrene	Acenaphthylene Benzo(g,h,i)perylene Dibenz[a,h]anthracene Naphthalene	Anthracene Benzo(k)fluoranthene Fluoranthene Phenanthrene
		Benz[a]anthracene Benzo(b)fluoranthene Fluorene Pyrene
PAH-Mix 9 deuterated		
DRE-L20950902CY	PAH-Mix 9 deuterated 10 µg/mL in Cyclohexane(‡)	10ml
DRE-XA20950902CY	PAH-Mix 9 deuterated 100 µg/mL in Cyclohexane(‡)	1ml
Acenaphthene D10 Benzo(a)pyrene D12 Chrysene D12 Indeno(1,2,3-c,d)pyrene D12	Acenaphthylene D8 Benzo(k)fluoranthene D12 Dibenz[a,h]anthracene D14 Naphthalene D8	Anthracene D10 Benzo(b)fluoranthene D12 Fluoranthene D10 Phenanthrene D10
		Benz[a]anthracene D12 Benzo(g,h,i)perylene D12 Fluorene D10 Pyrene D10
PAH-Mix 13		
DRE-L20950013AL	PAH-Mix 13 in Acetonitrile	10ml
Acenaphthene [100 µg/mL] Benzo(a)pyrene [10 µg/mL] Chrysene [10 µg/mL] Indeno(1,2,3-c,d)pyrene [10 µg/mL]	Acenaphthylene [100 µg/mL] Benzo(g,h,i)perylene [10 µg/mL] Dibenz[a,h]anthracene [10 µg/mL] Naphthalene [100 µg/mL]	Anthracene [10 µg/mL] Benzo(k)fluoranthene [10 µg/mL] Fluoranthene [10 µg/mL] Phenanthrene [10 µg/mL]
		Benz[a]anthracene [10 µg/mL] Benzo(b)fluoranthene [10 µg/mL] Fluorene [10 µg/mL] Pyrene [10 µg/mL]
PAH-Mix 14		
DRE-L20950014AL	PAH-Mix 14 10 µg/mL in Acetonitrile(‡)	10ml
DRE-L20950014CY	PAH-Mix 14 10 µg/mL in Cyclohexane	10ml
DRE-YA20950014AB	PAH-Mix 14 2000 µg/mL in Acetone/Benzene	1ml
1-Methylnaphthalene Anthracene Benzo(k)fluoranthene Fluoranthene Phenanthrene	2-Methylnaphthalene Benz[a]anthracene Benzo(b)fluoranthene Fluorene Pyrene	Acenaphthene Benzo(a)pyrene Chrysene Indeno(1,2,3-c,d)pyrene
		Acenaphthylene Benzo(g,h,i)perylene Dibenz[a,h]anthracene Naphthalene
PAH-Mix 18		
DRE-L20950018AL	PAH-Mix 18 10 µg/mL in Acetonitrile	10ml
Acenaphthene Benzo(a)pyrene Chrysene Indeno(1,2,3-c,d)pyrene Pyrene	Acenaphthylene Benzo(g,h,i)perylene Dibenz[a,h]anthracene Naphthalene	Anthracene Benzo(k)fluoranthene Fluoranthene Perylene
		Benz[a]anthracene Benzo(b)fluoranthene Fluorene Phenanthrene

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description	
PAH-Mix 20		
DRE-L20950020AL	PAH-Mix 20 10 µg/mL in Acetonitrile	10ml
	Benzo(a)pyrene Benzo(g,h,i)perylene Benzo(k)fluoranthene Benzo(b)fluoranthene Fluoranthene Indeno(1,2,3-c,d)pyrene	
PAH-Mix 24 deuterated		
DRE-LA20950024HE	PAH-Mix 24 deuterated 10 µg/mL in n-Hexane	1ml
	Acenaphthene D10 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10	
PAH-Mix 25		
DRE-YA20950025AB	PAH-Mix 25 2000 µg/mL in Acetone/Benzene(‡)	1ml
Acenaphthene Benzo(a)pyrene Chrysene Indeno(1,2,3-c,d)pyrene	Acenaphthylene Benzo(g,h,i)perylene Dibenz[a,h]anthracene Naphthalene	Anthracene Benzo(k)fluoranthene Fluoranthene Phenanthrene
		Benz[a]anthracene Benzo[b]fluoranthene Fluorene Pyrene
PAH-Mix 27		
DRE-LA20950027AL	PAH-Mix 27 in Acetonitrile	1ml
	Benzo(a)pyrene [25 µg/mL] Benzo(g,h,i)perylene [125 µg/mL] Benzo(k)fluoranthene [25 µg/mL] Benzo(b)fluoranthene [125 µg/mL] Fluoranthene [200 µg/mL] Indeno(1,2,3-c,d)pyrene [125 µg/mL]	
PAH-Mix 31 deuterated		
DRE-YA20950031TO	PAH-Mix 31 deuterated 1000 µg/mL in Toluene(‡)	1ml
	Acenaphthene D10 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10	
PAH-Mix 39		
DRE-X20950039AL	PAH-Mix 39 in Acetonitrile(‡)	10ml
Acenaphthene [50 µg/mL] Benzo(a)pyrene [25 µg/mL] Chrysene [25 µg/mL] Indeno(1,2,3-c,d)pyrene [100 µg/mL]	Acenaphthylene [100 µg/mL] Benzo(ghi)perylene [50 µg/mL] Dibenz[a,h]anthracene [50 µg/mL] Naphthalene [100 µg/mL]	Anthracene [10 µg/mL] Benzo(k)fluoranthene [10 µg/mL] Fluoranthene [50 µg/mL] Phenanthrene [50 µg/mL]
		Benz[a]anthracene [25 µg/mL] Benzo[b]fluoranthene [25 µg/mL] Fluorene [25 µg/mL] Pyrene [50 µg/mL]
PAH-Mix 45		
DRE-L20950045AL	PAH-Mix 45 10 µg/mL in Acetonitrile	10ml
DRE-L20950045CY	PAH-Mix 45 10 µg/mL in Cyclohexane	10ml
Acenaphthene Benzo(a)pyrene Benzo(b)fluoranthene Fluorene Phenanthrene	Acenaphthylene Benzo(e)pyrene Chrysene Indeno(1,2,3-c,d)pyrene Pyrene	Anthracene Benzo(g,h,i)perylene Dibenz[a,h]anthracene Naphthalene
		Benz[a]anthracene Benzo(k)fluoranthene Fluoranthene Perylene

Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description	
PAH-Mix 50		
DRE-XA05500100AL	PAH-Mix 50 in Acetonitrile	1ml
Acenaphthene [1000 µg/mL] Benzo(a)pyrene [5 µg/mL] Chrysene [50 µg/mL] Indeno(1,2,3-c,d)pyrene [10 µg/mL]	Acenaphthylene [1000 µg/mL] Benzo(ghi)perylene [5 µg/mL] Dibenz[a,h]anthracene [10 µg/mL] Naphthalene [1000 µg/mL]	Anthracene [50 µg/mL] Benzo(k)fluoranthene [1 µg/mL] Fluoranthene [50 µg/mL] Phenanthrene [50 µg/mL]
		Benz[a]anthracene [1 µg/mL] Benzo[b]fluoranthene [1 µg/mL] Fluorene [100 µg/mL] Pyrene [50 µg/mL]
PAH-Mix 61		
DRE-XA06100100AM	PAH Mix 61 in Methanol/Acetone 1:1	1ml
Acenaphthene [1000 µg/mL] Benzo(a)pyrene [100 µg/mL] Chrysene [100 µg/mL] Indeno(1,2,3-c,d)pyrene [100 µg/mL]	Acenaphthylene [2000 µg/mL] Benzo(g,h,i)perylene [200 µg/mL] Dibenz[a,h]anthracene [200 µg/mL] Naphthalene [1000 µg/mL]	Anthracene [100 µg/mL] Benzo(k)fluoranthene [100 µg/mL] Fluoranthene [200 µg/mL] Phenanthrene [100 µg/mL]
		Benz[a]anthracene [100 µg/mL] Benzo[b]fluoranthene [200 µg/mL] Fluorene [200 µg/mL] Pyrene [100 µg/mL]
PAH-Mix 63		
DRE-YA06100300TO	PAH Mix 63, 1000 µg/mL in Toluene(†)	1ml
	Acenaphthene Anthracene Benzo(a)pyrene Benzo(k)fluoranthene Chrysene Indeno(1,2,3-c,d)pyrene Phenanthrene	Acenaphthylene Benz[a]anthracene Benzo(g,h,i)perylene Benzo[b]fluoranthene Dibenz[a,h]anthracene Naphthalene Pyrene
PAH-Mix 64		
DRE-YA06100400BD	PAH Mix 64, 2000 µg/mL in Benzene/Dichloromethane 1:1	1ml
Acenaphthene Benzo(a)pyrene Carbazole Fluorene Pyrene	Acenaphthylene Benzo(g,h,i)perylene Chrysene Indeno(1,2,3-c,d)pyrene	Anthracene Benzo(k)fluoranthene Dibenz[a,h]anthracene Naphthalene
		Benz[a]anthracene Benzo[b]fluoranthene Fluoranthene Phenanthrene
PAH-Mix 77		
DRE-LA20950077TO	PAH-Mix 77 10 µg/mL in Toluene	1ml
	Acenaphthylene D8 Benzo(a)pyrene D12 Pyrene D10	
PAH-Mix 158		
DRE-LA20950158AL	PAH-Mix 158 10 µg/mL in Acetonitrile	1ml
Acenaphthene Benzo[b]fluoranthene Chrysene Indeno[1,2,3-c,d]pyrene Phenanthrene	Acenaphthylene Benzo[k]fluoranthene Dibenzo(a,h)anthracene 2-Methylfluoranthene Pyrene	Anthracene Benzo[g,h,i]perylene Fluoranthene 2-Methylnaphthalene
		Benz[a]anthracene Benzo[a]pyrene Fluorene Naphthalene
PAH-Mix 170		
DRE-LA20950170AL	PAH-Mix 170 10 µg/mL in Acetonitrile(‡)	1ml
	5-Methylchrysene Benzo(a)pyrene Benzo(j)fluoranthene Benzo[b]fluoranthene Cyclopenta(c,d)pyrene Dibenz[a,i]pyrene Dibenzo[a,e]pyrene Indeno(1,2,3-c,d)pyrene	Benz[a]anthracene Benzo(g,h,i)perylene Benzo(k)fluoranthene Chrysene Dibenz[a,h]anthracene Dibenz[a,i]pyrene Dibenzo[a,h]pyrene










Polycyclic aromatic hydrocarbons (PAHs)

Product code	Description			
PAH-Mix 183				
DRE-LA20950183CY	PAH-Mix 183 10 µg/mL in Cyclohexane(‡)			1ml
5-Methylchrysene	7H-Benzo(c)fluorene	Benz[a]anthracene	Benzo(a)pyrene	
Benzo(g,h,i)perylene	Benzo(j)fluoranthene	Benzo(k)fluoranthene	Benzo[b]fluoranthene	
Chrysene	Cyclopenta(c,d)pyrene	Dibenz[a,h]anthracene	Dibenz[a,i]pyrene	
Dibenz[a,l]pyrene	Dibenzo[a,e]pyrene	Dibenzo[a,h]pyrene	Indeno(1,2,3-c,d)pyrene	
PAH-Mix 197				
DRE-LS20950197CY	PAH-Mix 197, 10 µg/mL in Cyclohexane			3x10ml
Acenaphthene	Acenaphthylene	Anthracene	Benz[a]anthracene	
Benzo[a]pyrene	Benzo[b]fluoranthene	Benzo(g,h,i)perylene	Benzo[j]fluoranthene	
Benzo[k]fluoranthene	Chrysene	Dibenz[a,h]anthracene	Fluoranthene	
Fluorene	Indeno[1,2,3-c,d]pyrene	Naphthalene	Phenanthrene	
Pyrene				
PAH-Mix 525				
DRE-XA05250100AC	PAH Mix 525, 100 µg/mL in Acetone			1ml
	Acenaphthylene	Anthracene		
	Benz[a]anthracene	Benzo(a)pyrene		
	Benzo(g,h,i)perylene	Benzo(k)fluoranthene		
	Benzo[b]fluoranthene	Chrysene		
	Dibenz[a,h]anthracene	Fluorene		
	Indeno(1,2,3-c,d)pyrene	Phenanthrene		
	Pyrene			
PAH Mix A				
DRE-GA09000163DI	PAH Mix A 2000 µg/mL in Dichloromethane(‡)			1.5ml
DRE-GS09000163DI	PAH Mix A 2000 µg/mL in Dichloromethane(‡)			5x1.5ml
perylene	quinoline	acridine	benzo[k]fluoranthene	
1-methylnaphthalene	2-methylnaphthalene	acenaphthene	acenaphthylene	
anthracene	fluorene	naphthalene	phenanthrene	
benzo[a]anthracene	benzo[a]pyrene	chrysene	fluoranthene	
indeno[1,2,3-cd]pyrene	pyrene	benzo[b]fluoranthene	benzo[ghi]perylene	
dibenz[a,h]anthracene	benzo[e]pyrene			
PAH Mix, 25-6591				
DRE-GS09000390DI	PAH Mix, 25-6591 1000 µg/mL in Dichloromethane(‡)			5x1ml
benzo[k]fluoranthene	1-methylnaphthalene	2-methylnaphthalene	acenaphthene	
acenaphthylene	anthracene	fluorene	naphthalene	
phenanthrene	benzo[a]anthracene	benzo[a]pyrene	chrysene	
fluoranthene	indeno[1,2,3-cd]pyrene	pyrene	benzo[b]fluoranthene	
benzo[ghi]perylene	dibenz[a,h]anthracene	acridine	benzo(j)fluoranthene	
benzo[e]pyrene	perylene	quinoline	2-chloronaphthalene	
1-chloronaphthalene				
PAH Mix, 3-2011				
DRE-GA09000188TO	PAH Mix, 3-2011 10 µg/mL in Toluene(‡)			1ml
	triphenylene			
	benzo(j)fluoranthene			
	benzo[k]fluoranthene			







HYDROCARBONS AND PETROCHEMICALS









Hydrocarbons and Petrochemicals

Product code	Description			
Acetaldehyde-2,4-dinitrophenylhydrazone (DNPH)				
CAS 1019-57-4	MW 224.1735	C ₉ H ₈ N ₄ O ₄		
DRE-C10011300	Acetaldehyde-2,4-dinitrophenylhydrazone (DNPH)		100mg	
DRE-YA10011300AL	Acetaldehyde-2,4-dinitrophenylhydrazone (DNPH) 1000 µg/mL in Acetonitrile		1ml	
Acetamide				
CAS 60-35-5	MW 59.0672	C ₂ H ₅ NO		
DRE-C10011800	Acetamide		250mg	
Acetoacetic Acid Methyl Ester				
CAS 105-45-3	MW 116.1152	C ₅ H ₈ O ₃		
DRE-C10017400	Acetoacetic acid methyl ester		5ml	
Acetone				
CAS 67-64-1	MW 58.0791	C ₃ H ₆ O		
DRE-C10019000	Acetone		5ml	
Acetone-2,4-dinitrophenylhydrazone				
CAS 1567-89-1	MW 238.2001	C ₉ H ₁₀ N ₄ O ₄		
DRE-C10019500	Acetone-2,4-dinitrophenylhydrazone(‡)		25mg	
DRE-YA10019500AL	Acetone-2,4-dinitrophenylhydrazone (DNPH) 1000 µg/mL in Acetonitrile		1ml	
Acetophenone				
CAS 98-86-2	MW 120.1485	C ₈ H ₈ O		
DRE-CA10022000	Acetophenone		1ml	
Acetylacetamide (Acetoacetamide)				
CAS 5977-14-0	MW 101.1039	C ₄ H ₇ NO ₂		
DRE-C10017000	Acetoacetamide		250mg	
Acrolein-2,4-dinitrophenylhydrazone (DNPH)				
CAS 888-54-0	MW 236.1842	C ₉ H ₈ N ₄ O ₄		
DRE-CA10045200	Acrolein-2,4-dinitrophenylhydrazone (DNPH)		25mg	
DRE-YA10045200AL	Acrolein-2,4-dinitrophenylhydrazone (DNPH) 1000 µg/mL in Acetonitrile		1ml	
Acrylamide				
CAS 79-06-1	MW 71.0779	C ₃ H ₅ NO		
DRE-C10045300	Acrylamide(‡)		250mg	








Hydrocarbons and Petrochemicals

Product code	Description			
Allyl alcohol				
CAS 107-18-6	MW 58.0791	C ₃ H ₆ O		
DRE-C10130000	Allyl alcohol		0.5g	
DRE-GA09010377ME	Allyl Alcohol (SS) 500 µg/mL in p/t Methanol(‡)		1ml	
3-Amino-1-propanol				
CAS 156-87-6	MW 75.1097	C ₃ H ₉ NO		
DRE-C10217000	3-Amino-1-propanol		250mg	
tert-Amyl Methyl Ether				
CAS 994-05-8	MW 102.1748	C ₆ H ₁₄ O		
DRE-C10246440	tert-Amyl-methyl ether		1ml	
DRE-L10246440ME	Methyl-tert-amylether 10 µg/mL in Methanol		10ml	
tert-Amyl-ethyl Ether				
CAS 919-94-8	MW 116.2013	C ₇ H ₁₆ O		
DRE-CA10246400	tert-Amyl-ethyl ether		1ml	
Benzene				
CAS 71-43-2	MW 78.1118	C ₆ H ₆		
DRE-C10535000	Benzene(‡)		1ml	
DRE-L10535000ME	Benzene 10 µg/mL in Methanol		10ml	
DRE-XA10535000ME	Benzene 100 µg/mL in Methanol		1ml	
Benzene D6				
CAS 1076-43-3	MW 84.1488	C ₆ ² H ₆		
DRE-C10535200	Benzene D6		2ml	
DRE-YA10535200ME	Benzene D6 2000 µg/mL in Methanol		1ml	
Benzene 40-60/80°				
CAS 64742-49-0	MW n/a			No Structure
DRE-YA03011000ME	Benzene 40-60° 2500 µg/mL in Methanol		1ml	
DRE-YA03012000ME	Benzene 40-80° 2500 µg/mL in Methanol		1ml	
Benzylamine				
CAS 100-46-9	MW 107.1531	C ₇ H ₉ N		
DRE-C10569200	Benzylamine		250mg	
1,3-Butadiene				
CAS 106-99-0	MW 54.0904	C ₄ H ₆		
DRE-YA10860500ME	1,3-Butadiene 200 µg/mL in Methanol (*)		1ml	
DRE-YS09010022ME	1,3-Butadiene 1000 µg/mL in p/t Methanol(‡)		5x1ml	










Hydrocarbons and Petrochemicals

Product code	Description			
1,4-Butanediamine (1,4-Diaminobutane)				
CAS 110-60-1 DRE-C12193000	MW 88.1515 1,4-Diaminobutane	C ₄ H ₁₂ N ₂	250mg	
Butane-1,3-diol (Butylene Glycol)				
CAS 107-88-0 DRE-C10861300	MW 90.121 1,3-Butanediol	C ₄ H ₁₀ O ₂	250mg	
1,2-Butanediol				
CAS 584-03-2 DRE-C10861290	MW 90.121 1,2-Butanediol	C ₄ H ₁₀ O ₂	250mg	
1,4-Butanediol				
CAS 110-63-4 DRE-C10861310	MW 90.121 1,4-Butanediol	C ₄ H ₁₀ O ₂	1ml	
2,3-Butanediol				
CAS 513-85-9 DRE-C10861320	MW 90.121 2,3-Butanediol	C ₄ H ₁₀ O ₂	250mg	
1-Butanethiol				
CAS 109-79-5 DRE-CA10861410	MW 90.1872 1-Butanethiol	C ₄ H ₁₀ S	1ml	
1-Butanol				
CAS 71-36-3 DRE-C10861500 DRE-GA09010375ME	MW 74.1216 1-Butanol(‡) 1-Butanol (SS) 500 µg/mL in Methanol(‡)	C ₄ H ₁₀ O	1ml 1ml	
2-Butanol				
CAS 78-92-2 DRE-C10861600	MW 74.1216 2-Butanol	C ₄ H ₁₀ O	5ml	
tert.-Butanol				
CAS 75-65-0 DRE-C10861700	MW 74.1216 tert.-Butanol(‡)	C ₄ H ₁₀ O	1ml	










Hydrocarbons and Petrochemicals

Product code	Description			
2-Butanone				
CAS 78-93-3	MW 72.1057	C ₄ H ₈ O		
DRE-C10862000	2-Butanone(‡)		5ml	
DRE-YA10862000ME	2-Butanone 2000 µg/mL in Methanol		1ml	
2-Butanone-2,4-dinitrophenylhydrazone				
CAS 958-60-1	MW 252.2267	C ₁₀ H ₁₂ N ₄ O ₄		
DRE-C10862050	2-Butanone-2,4-dinitrophenylhydrazone		100mg	
2-(2-Butoxyethoxy)ethyl Acetate (Diethylene glycol monobutylether acetate)				
CAS 124-17-4	MW 204.2634	C ₁₀ H ₂₀ O ₄		
DRE-C12606000	Diethylene glycol-monobutylether acetate		1g	
1-Butoxy-2-propanol				
CAS 5131-66-8	MW 132.2007	C ₇ H ₁₆ O ₂		
DRE-C10900100	1-Butoxy-2-propanol		100mg	
Butyl Acetate				
CAS 123-86-4	MW 116.1583	C ₈ H ₁₆ O ₂		
DRE-C10929000	n-Butylacetate		1ml	
tert-Butyl Ethyl Ether				
CAS 637-92-3	MW 102.1748	C ₆ H ₁₄ O		
DRE-CA13321300	Ethyl-tert-butylether		100mg	
Butyl Hexanoate (n-Caproic acid n-butyl ester)				
CAS 626-82-4	MW 172.2646	C ₁₀ H ₂₀ O ₂		
DRE-C10948031	n-Caproic acid-n-butyl ester		100mg	
Butyl Methacrylate (Methacrylic acid butyl ester)				
CAS 97-88-1	MW 142.1956	C ₈ H ₁₄ O ₂		
DRE-CA14971720	Methacrylic acid butyl ester		1ml	
1-Butylamine				
CAS 109-73-9	MW 73.1368	C ₄ H ₁₁ N		
DRE-C10929100	1-Butylamine		250mg	

Hydrocarbons and Petrochemicals

Product code	Description			
2-Butylamine				
CAS 13952-84-6 DRE-CA10929200	MW 73.1368 2-Butylamine	C ₄ H ₁₁ N	250mg	
n-Butylbenzene				
CAS 104-51-8 DRE-C10930900	MW 134.2182 n-Butylbenzene	C ₁₀ H ₁₄	1ml	
sec-Butylbenzene				
CAS 135-98-8 DRE-CA10931000	MW 134.2182 sec-Butylbenzene	C ₁₀ H ₁₄	1ml	
tert-Butylbenzene				
CAS 98-06-6 DRE-C10931100	MW 134.2182 tert.-Butylbenzene	C ₁₀ H ₁₄	1ml	
2-sec-Butyl-1-(decyloxy)-4-tritylbenzene				
CAS 1404190-37-9 DRE-C10931165	MW 532.7978 2-sec-Butyl-1-(decyloxy)-4-tritylbenzene	C ₃₉ H ₄₆ O	25mg	
2-Butyl-2-ethyl-1,3-propanediol				
CAS 115-84-4 DRE-C10931196	MW 160.2539 2-Butyl-2-ethyl-1,3-propanediol	C ₉ H ₂₀ O ₂	1g	
3-(4-tert-Butylphenyl)isobutyraldehyde				
CAS 80-54-6 DRE-CA10931620	MW 204.308 3-(4-tert-Butylphenyl)iso-butyraldehyde	C ₁₄ H ₂₀ O	250mg	
4-tert-Butyltoluene				
CAS 98-51-1 DRE-C10931730	MW 148.2447 4-tert.-Butyltoluene	C ₁₁ H ₁₆	100mg	
Butyraldehyde (Butanal)				
CAS 123-72-8 DRE-C10931740	MW 72.1057 Butyraldehyde (Butanal)	C ₄ H ₈ O	250mg	







Hydrocarbons and Petrochemicals

Product code	Description			
Butyraldehyde-2,4-dinitrophenylhydrazone				
CAS 1527-98-6 DRE-C10931745	MW 252.2267 Butyraldehyde-2,4-dinitrophenylhydrazone	C ₁₀ H ₁₂ N ₄ O ₄	100mg	
n-Butyrophenone				
CAS 495-40-9 DRE-C10931900	MW 148.2017 n-Butyrophenone	C ₁₀ H ₁₂ O	250mg	
1-Chlorodecane				
CAS 1002-69-3 DRE-C11398000	MW 176.7267 1-Chlorodecane	C ₁₀ H ₂₁ Cl	250mg	
1-Chlorooctadecane				
CAS 3386-33-2 DRE-C11457400 DRE-YS09010015DI	MW 288.9394 1-Chlorooctadecane 1-Chlorooctadecane 10000 µg/mL in Dichloromethane(±)	C ₁₈ H ₃₇ Cl	250mg 5x1ml	
1-Chlorooctane				
CAS 111-85-3 DRE-C11457450	MW 148.6736 1-Chlorooctane	C ₈ H ₁₇ Cl	250mg	
3-Chloropropane-1,2-diol				
CAS 96-24-2 DRE-C11502630	MW 110.5395 3-Chloro-1,2-propanediol(±)	C ₃ H ₇ ClO ₂	250mg	
2-Chloropropionic Acid Ethyl Ester				
CAS 535-13-7 DRE-C11502800	MW 136.5768 2-Chloropropionic acid-ethyl ester	C ₅ H ₉ ClO ₂	100mg	
1-Chlorotetradecane				
CAS 2425-54-9 DRE-C11509500	MW 232.8331 1-Chlorotetradecane	C ₁₄ H ₂₉ Cl	250mg	
Crotonaldehyde (Butenal)				
CAS 123-73-9 DRE-C11755000	MW 70.0898 Crotonaldehyde (Butenal)	C ₄ H ₆ O	5ml	









Hydrocarbons and Petrochemicals

Product code	Description			
Crotonaldehyde-2,4-dinitrophenylhydrazone (DNPH)				
CAS 1527-96-4 DRE-CA11755200	MW 250.2108 Crotonaldehyde-2,4-dinitrophenylhydrazone (DNPH)	C ₁₀ H ₁₀ N ₄ O ₄	100mg	
Cyclohexane				
CAS 110-82-7 DRE-C11824500 DRE-GA11824500ME DRE-GH11824500ME	MW 84.1595 Cyclohexane Cyclohexane 2000 µg/mL in Methanol Cyclohexane 2000 µg/mL in Methanol	C ₆ H ₁₂	5ml 1ml 5x1ml	
Cyclohexanone				
CAS 108-94-1 DRE-C11825500 DRE-GA09010400ME	MW 98.143 Cyclohexanone Cyclohexanone Solution, 5000 µg/mL in Methanol(‡)	C ₆ H ₁₀ O	1ml 1ml	
Cyclohexanone-2,4-dinitrophenylhydrazone				
CAS 1589-62-4 DRE-C11825510	MW 278.264 Cyclohexanone-2,4-dinitrophenylhydrazone	C ₁₂ H ₁₄ N ₄ O ₄	100mg	
Cyclohexylamine (Cyclohexanamine)				
CAS 108-91-8 DRE-C11830500 DRE-GA09010092ME	MW 99.1741 Cyclohexylamine(‡) Cyclohexylamine 1000 µg/mL in Methanol(‡)	C ₆ H ₁₃ N	250mg 1ml	
N-Cyclohexyldimethylamine				
CAS 98-94-2 DRE-C11830510	MW 127.2273 N-Cyclohexyldimethylamine	C ₈ H ₁₇ N	250mg	
Cyclopentadecanone				
CAS 502-72-7 DRE-C11832000	MW 224.3822 Cyclopentadecanone	C ₁₅ H ₂₈ O	100mg	
Cyclopentane				
CAS 287-92-3 DRE-C11833000	MW 70.1329 Cyclopentane	C ₅ H ₁₀	5ml	
Cyclopentanol				
CAS 96-41-3 DRE-C11833520	MW 86.1323 Cyclopentanol	C ₅ H ₁₀ O	1g	









Hydrocarbons and Petrochemicals

Product code	Description			
Cyclopentanone				
CAS 120-92-3 DRE-C11833500	MW 84.1164 Cyclopentanone	C ₅ H ₈ O	1ml	
Decahydronaphthalene				
CAS 91-17-8 DRE-C12093900	MW 138.2499 Decahydronaphthalene	C ₁₀ H ₁₈	100mg	
cis-Decahydronaphthalene				
CAS 493-01-6 DRE-C12094000	MW 138.2499 cis-Decahydronaphthalene	C ₁₀ H ₁₈	250mg	AND Enantiomer 
trans-Decahydronaphthalene				
CAS 493-02-7 DRE-C12094100	MW 138.2499 trans-Decahydronaphthalene	C ₁₀ H ₁₈	250mg	AND Enantiomer 
Decanal-2,4-dinitrophenylhydrazone				
CAS 1527-95-3 DRE-XA12094810AL	MW 336.3862 Decanal-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitrile	C ₁₆ H ₂₄ N ₄ O ₄	1ml	
n-Decane				
CAS 124-18-5 DRE-C12095000 DRE-GA09010087ME	MW 142.2817 n-Decane n-Decane (C10) 1000 µg/mL in p/t Methanol(†)	C ₁₀ H ₂₂	1ml 1ml	
1,10-Decanediol				
CAS 112-47-0 DRE-C12095010	MW 174.2805 1,10-Decanediol	C ₁₀ H ₂₂ O ₂	250mg	
1-Decene				
CAS 872-05-9 DRE-C12096000	MW 140.2658 1-Decene	C ₁₀ H ₂₀	1ml	
1-Decylamine				
CAS 2016-57-1 DRE-C12098500	MW 157.2963 1-Decylamine	C ₁₀ H ₂₃ N	250mg	









Hydrocarbons and Petrochemicals

Product code	Description			
Decylbenzene				
CAS 104-72-3 DRE-C12098800	MW 218.3776 Decylbenzene	C ₁₆ H ₂₆	100mg	
1,5-Diaminopentane				
CAS 462-94-2 DRE-CA12196100	MW 102.1781 1,5-Diaminopentane	C ₅ H ₁₄ N ₂	250mg	
Dibenzylamine				
CAS 103-49-1 DRE-C12214500	MW 197.2756 Dibenzylamine	C ₁₄ H ₁₉ N	250mg	
Dibutyl Sulfide				
CAS 544-40-1 DRE-CA12256500	MW 146.2935 Dibutylsulfide	C ₈ H ₁₈ S	250mg	
Dibutylamine				
CAS 111-92-2 DRE-C12250000	MW 129.2432 Dibutylamine	C ₈ H ₁₉ N	1g	
Di-n-butylether				
CAS 142-96-1 DRE-C12252000	MW 130.2279 Di-n-butylether	C ₈ H ₁₈ O	5ml	
Dicyclohexylamine (N-Cyclohexylcyclohexanamine)				
CAS 101-83-7 DRE-C12585000	MW 181.3177 Dicyclohexylamine	C ₁₂ H ₂₃ N	250mg	
Diesel Fuel				
CAS 68334-30-5 DRE-C03009000	MW n/a Diesel Oil without additives		1g	
DRE-XA03009000MB	Diesel Oil (without additives, DIN H53) (EN9377-2) 100 µg/mL in Methyl-tert. butyl ether		1ml	
DRE-YA03004000ME	Regular Diesel Fuel 2500 µg/mL in Methanol		1ml	
DRE-CA03009020	Set Diesel Oil and Mineral Oil (DIN H53) (EN9377-2)			
DRE-CA03009000	Diesel Oil (without additives, DIN H53) (EN9377-2)		3x1ml	
DRE-CA03009011	Mineral Oil Heavy (without additives, DIN H53) (EN9377-2)		3x1ml	







Hydrocarbons and Petrochemicals

Product code	Description			
Bio Diesel Fuel				
CAS 870530-78-2 DRE-CA03005600	MW n/a Bio Diesel Fuel		1ml	No Structure
1,1-Diethoxyethane (Acetaldehyde diethylacetal)				
CAS 105-57-7 DRE-C10011200	MW 118.1742 Acetaldehyde diethylacetal(†)	C ₆ H ₁₄ O ₂	1ml	
Diethylamine				
CAS 109-89-7 DRE-CA12604500	MW 73.1368 Diethylamine	C ₄ H ₁₁ N	1ml	
N,N-Diethylaniline				
CAS 91-66-7 DRE-C12604800	MW 149.2328 N,N-Diethylaniline	C ₁₀ H ₁₅ N	250mg	
1,2-Diethylbenzene				
CAS 135-01-3 DRE-C12605400 DRE-XA12605400ME	MW 134.2182 1,2-Diethylbenzene 1,2-Diethylbenzene 100 µg/mL in Methanol	C ₁₀ H ₁₄	100mg 1ml	
1,3-Diethylbenzene				
CAS 141-93-5 DRE-C12605500 DRE-XA12605500ME	MW 134.2182 1,3-Diethylbenzene 1,3-Diethylbenzene 100 µg/mL in Methanol	C ₁₀ H ₁₄	100mg 1ml	
1,4-Diethylbenzene				
CAS 105-05-5 DRE-CA12605600 DRE-XA12605600ME	MW 134.2182 1,4-Diethylbenzene 1,4-Diethylbenzene 100 µg/mL in Methanol	C ₁₀ H ₁₄	.5ml 1ml	
Diethylene Glycol Monobutyl Ether				
CAS 112-34-5 DRE-C12605900	MW 162.2267 Diethylene glycol-monobutyl ether	C ₈ H ₁₈ O ₃	250mg	
Diethylene Glycol Monohexyl Ether				
CAS 112-59-4 DRE-C12606300	MW 190.2799 Diethylene glycol-monohexylether	C ₁₀ H ₂₂ O ₃	250mg	








Hydrocarbons and Petrochemicals

Product code	Description			
Diethylene Glycol Monomethyl Ether				
CAS 111-77-3 DRE-C12606400	MW 120.147 Diethylene glycol-monomethylether	C ₅ H ₁₂ O ₃	250mg	
Diethylene glycol monoethyl ether				
CAS 111-90-0 DRE-C12606200	MW 134.1736 Diethylene glycol-monoethylether	C ₆ H ₁₄ O ₃	250mg	
Diethyleneglycol (2,2'-Oxydiethanol)				
CAS 111-46-6 DRE-C12605780	MW 106.1204 Diethyleneglycol(†)	C ₄ H ₁₀ O ₃	1ml	
Diethylether (Ether)				
CAS 60-29-7 DRE-C12606500	MW 74.1216 Diethylether DRUG PRECURSOR	C ₄ H ₁₀ O	5ml	
Diglyme (Bis(2-methoxyethyl)ether; Diethylene glycol dimethyl ether)				
CAS 111-96-6 DRE-CA10653800	MW 134.1736 Bis(2-methoxyethyl)ether (Diethylene glycol dimethyl ether)	C ₆ H ₁₄ O ₃	1ml	
1,3-Diisopropylbenzene				
CAS 99-62-7 DRE-C12637300	MW 162.2713 1,3-Diisopropylbenzene	C ₁₂ H ₁₈	250mg	
1,4-Diisopropylbenzene				
CAS 100-18-5 DRE-C12637400	MW 162.2713 1,4-Diisopropylbenzene	C ₁₂ H ₁₈	1ml	
Diisopropylether				
CAS 108-20-3 DRE-C12637420	MW 102.1748 Diisopropylether	C ₆ H ₁₄ O	1ml	
Dimethyl sulfoxide				
CAS 67-68-5 DRE-C12745000	MW 78.1334 Dimethylsulfoxide	C ₂ H ₆ OS	5ml	

Hydrocarbons and Petrochemicals

Product code	Description				
1,2-Dimethylbenzene (o-Xylene)					
CAS 95-47-6	MW 106.165	C ₈ H ₁₀			
DRE-C17945000	o-Xylene (1,2-Dimethylbenzene)(‡)		1ml		
DRE-L17945000ME	o-Xylene 10 µg/mL in Methanol		10ml		
DRE-XA17945000ME	o-Xylene 100 µg/mL in Methanol		1ml		
1,3-Dimethylbenzene (m-Xylene)					
CAS 108-38-3	MW 106.165	C ₈ H ₁₀			
DRE-C17945100	m-Xylene (1,3-Dimethylbenzene)(‡)		1ml		
DRE-L17945100ME	m-Xylene 10 µg/mL in Methanol		10ml		
DRE-XA17945100ME	m-Xylene 100 µg/mL in Methanol		1ml		
1,4-Dimethylbenzene (p-Xylene)					
CAS 106-42-3	MW 106.165	C ₈ H ₁₀			
DRE-C17945200	p-Xylene (1,4-Dimethylbenzene)(‡)		1ml		
DRE-L17945200ME	p-Xylene 10 µg/mL in Methanol		10ml		
DRE-YA17945200ME	p-Xylene, 5000 µg/mL in Methanol		1ml		
2,2-Dimethylbutane					
CAS 75-83-2	MW 86.1754	C ₆ H ₁₄			
DRE-CA12726100	2,2-Dimethylbutane		100mg		
DRE-GS09010150TN	2,2-Dimethylbutane 1000 µg/mL in Triacetin(‡)		5x1ml		
2,3-Dimethylbutane					
CAS 79-29-8	MW 86.1754	C ₆ H ₁₄			
DRE-C12726110	2,3-Dimethylbutane(‡)		100mg		
1,3-Dimethylcyclohexane					
CAS 591-21-9	MW 112.2126	C ₈ H ₁₆			
DRE-CA12726430	1,3-Dimethylcyclohexane (cis/trans mixture)		100mg		
1,4-Dimethylcyclohexane					
CAS 589-90-2	MW 112.2126	C ₈ H ₁₆			
DRE-C12726460	1,4-Dimethylcyclohexane		100mg		
N,N-Dimethylformamide					
CAS 68-12-2	MW 73.0938	C ₃ H ₇ NO			
DRE-C12727000	N,N-Dimethylformamide		5ml		
2,6-Dimethyl-4-heptanone					
CAS 108-83-8	MW 142.2386	C ₉ H ₁₈ O			
DRE-C12727300	2,6-Dimethyl-4-heptanone		250mg		










Hydrocarbons and Petrochemicals

Product code	Description			
2,2-Dimethylpentane				
CAS 590-35-2 DRE-C12728200	MW 100.2019 2,2-Dimethylpentane	C ₇ H ₁₆	100mg	
2,3-Dimethylpentane				
CAS 565-59-3 DRE-C12728300	MW 100.2019 2,3-Dimethylpentane	C ₇ H ₁₆	1ml	
2,4-Dimethylpentane				
CAS 108-08-7 DRE-C12728400 DRE-CA12728400	MW 100.2019 2,4-Dimethylpentane 2,4-Dimethylpentane	C ₇ H ₁₆	1g .5ml	
2,2-Dimethylpropane (Neopentane)				
CAS 463-82-1 DRE-GS09010149TN	MW 72.1488 2,2-Dimethylpropane 1000 µg/mL in Triacetin(‡)	C ₅ H ₁₂	5x1ml	
1,4-Dioxane-d8 (Octadeuterodioxane)				
CAS 17647-74-4 DRE-GA09010386ME	MW 96.1544 1,4-Dioxane-d8 Solution, 10000 µg/mL in p/t Methanol(‡)	C ₆ H ₈ O ₂	1ml	
1,4-Dioxane				
CAS 123-91-1 DRE-C12865000 DRE-GA09010345ME	MW 88.1051 1,4-Dioxane 1,4-Dioxane 1000 µg/mL in p/t Methanol(‡)	C ₄ H ₈ O ₂	5ml 1ml	
Diphenyl Sulfide				
CAS 139-66-2 DRE-CA12910000	MW 186.2728 Diphenylsulfide	C ₁₂ H ₁₀ S	500mg	
Dipropylamine				
CAS 142-84-7 DRE-C12938000	MW 101.19 Dipropylamine	C ₆ H ₁₅ N	250mg	
Dipropylene Glycol Monomethyl Ether				
CAS 34590-94-8 DRE-C12938500	MW 148.2001 Dipropyleneglycol-monomethylether	C ₇ H ₁₆ O ₃	250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
Dipropyleneglycol				
CAS 25265-71-8 DRE-C12938480	MW 402.5207 Dipropyleneglycol	$((C_6H_{14}O_3)_c(C_6H_{14}O_3)_c(C_6H_{14}O_3)_c)_{mix}$	100mg	
n-Docosane				
CAS 629-97-0 DRE-C13057700	MW 310.6006 n-Docosane	$C_{22}H_{46}$	0.5g	
1-Docosanol (Behenyl alcohol)				
CAS 661-19-8 DRE-C13058000	MW 326.6 1-Docosanol (Behenyl alcohol)	$C_{22}H_{46}O$	250mg	
1-Docosene				
CAS 1599-67-3 DRE-C13058050	MW 308.5848 1-Docosene	$C_{22}H_{44}$	100mg	
Dodecanal				
CAS 112-54-9 DRE-C13059800	MW 184.3184 Dodecanal	$C_{12}H_{24}O$	250mg	
n-Dodecane				
CAS 112-40-3 DRE-C13060000	MW 170.3348 n-Dodecane(†)	$C_{12}H_{26}$	5ml	
1-Dodecanol				
CAS 112-53-8 DRE-C13061000	MW 186.3342 1-Dodecanol	$C_{12}H_{26}O$	250mg	
2-Dodecanol				
CAS 10203-28-8 DRE-C13061010	MW 186.3342 2-Dodecanol	$C_{12}H_{26}O$	250mg	
1-Dodecene				
CAS 112-41-4 DRE-C13061500	MW 168.319 1-Dodecene	$C_{12}H_{24}$	250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
n-Dotriacontane				
CAS 544-85-4 DRE-C13084000	MW 450.8664 n-Dotriacontane	C ₃₂ H ₆₆	250mg	
n-Eicosane				
CAS 112-95-8 DRE-C13112700	MW 282.5475 n-Eicosane	C ₂₀ H ₄₂	250mg	
1-Eicosanol (Arachidyl alcohol)				
CAS 629-96-9 DRE-C13113000	MW 298.5469 1-Eicosanol (Arachidyl alcohol)	C ₂₀ H ₄₂ O	100mg	
1-Eicosene				
CAS 3452-07-1 DRE-C13113300	MW 280.5316 1-Eicosene	C ₂₀ H ₄₀	100mg	
Ethanol				
CAS 64-17-5 DRE-C13223000 DRE-YA13223000ME	MW 46.0684 Ethanol(±) Ethanol 2000 µg/mL in Methanol	C ₂ H ₆ O	1ml 1ml	
2-Ethoxyethanol (Ethyleneglycol monoethyl ether)				
CAS 110-80-5 DRE-C13328100	MW 90.121 Ethyleneglycol-monoethyl ether	C ₄ H ₁₀ O ₂	250mg	
1-Ethoxy-2-propanol				
CAS 1569-02-4 DRE-C13309000	MW 104.1476 1-Ethoxy-2-propanol	C ₅ H ₁₂ O ₂	250mg	
Ethyl acetate				
CAS 141-78-6 DRE-C13319000	MW 88.1051 Ethylacetate	C ₄ H ₈ O ₂	5ml	
N-Ethylaniline				
CAS 103-69-5 DRE-C13319500	MW 121.1796 N-Ethylaniline	C ₈ H ₁₁ N	100mg	







Hydrocarbons and Petrochemicals

Product code	Description			
Ethylbenzene				
CAS 100-41-4	MW 106.165	C ₈ H ₁₀		
DRE-C13320000	Ethylbenzene(‡)		1ml	
DRE-L13320000CY	Ethylbenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA13320000CY	Ethylbenzene 100 µg/mL in Cyclohexane		1ml	
Ethylbenzene-D10				
CAS 25837-05-2	MW 116.2266	C ₈ ² H ₁₀		
DRE-YA13320100ME	Ethylbenzene D10 2000 µg/mL in Methanol		1ml	
2-Ethyl-1-butanol				
CAS 97-95-0	MW 102.1748	C ₆ H ₁₄ O		
DRE-C13321000	2-Ethyl-1-butanol		250mg	
2-Ethyl-1-butene				
CAS 760-21-4	MW 84.1595	C ₆ H ₁₂		
DRE-CA13321200	2-Ethyl-1-butene		100mg	
Ethylcyclohexane				
CAS 1678-91-7	MW 112.2126	C ₈ H ₁₆		
DRE-C13324000	Ethylcyclohexane		250mg	
Ethylene Glycol				
CAS 107-21-1	MW 62.0678	C ₂ H ₆ O ₂		
DRE-C13327000	Ethyleneglycol(‡)		1ml	
DRE-GA13327000ME	Ethyleneglycol 1000 µg/mL in Methanol		1ml	
DRE-GH13327000ME	Ethyleneglycol 1000 µg/mL in Methanol		5x1ml	
Ethylenediamine				
CAS 107-15-3	MW 60.0983	C ₂ H ₈ N ₂		
DRE-C13326500	Ethylenediamine		1ml	
Ethyleneglycol Monobutyl Ether				
CAS 111-76-2	MW 118.1742	C ₆ H ₁₄ O ₂		
DRE-C13328000	Ethyleneglycol-monobutyl ether(‡)		1g	
2-Ethyl-1,3-hexandiol				
CAS 94-96-2	MW 146.2273	C ₈ H ₁₈ O ₂		
DRE-C13340000	2-Ethyl-1,3-hexandiol		250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
2-Ethylhexanoic Acid				
CAS 149-57-5 DRE-C13340100	MW 144.2114 2-Ethylhexanoic acid	C ₈ H ₁₆ O ₂	1ml	
2-Ethyl-1-hexanol				
CAS 104-76-7 DRE-C13340200	MW 130.2279 2-Ethyl-1-hexanol	C ₈ H ₁₈ O	250mg	
2-Ethyltoluene				
CAS 611-14-3 DRE-CA13356200 DRE-XA13356200ME	MW 120.1916 2-Ethyltoluene 2-Ethyltoluene 100 µg/mL in Methanol	C ₉ H ₁₂	250mg 1ml	
3-Ethyltoluene				
CAS 620-14-4 DRE-CA13356300 DRE-XA13356300ME	MW 120.1916 3-Ethyltoluene 3-Ethyltoluene 100 µg/mL in Methanol	C ₉ H ₁₂	100mg 1ml	
4-Ethyltoluene				
CAS 622-96-8 DRE-C13356400 DRE-XA13356400ME	MW 120.1916 4-Ethyltoluene 4-Ethyltoluene 100 µg/mL in Methanol	C ₉ H ₁₂	100mg 1ml	
Formaldehyde-2,4-dinitrophenylhydrazone (DNPH)				
CAS 1081-15-8 DRE-XA13909200AL DRE-YA13909200AL	MW 210.1469 Formaldehyde-2,4-dinitrophenylhydrazone (DNPH) 100 µg/mL in Acetonitrile Formaldehyde-2,4-dinitrophenylhydrazone (DNPH) 1000 µg/mL in Acetonitrile	C ₇ H ₆ N ₄ O ₄	1ml 1ml	
Furan				
CAS 110-00-9 DRE-CA13965000 DRE-XA13965000AL	MW 68.074 Furan Furan 100 µg/mL in Acetonitrile	C ₄ H ₄ O	1g 1ml	
Furan-d4				
CAS 6142-90-1 DRE-XA13965010AL	MW 72.0986 Furan D4 100µg/mL in Acetonitrile	C ₄ ² H ₄ O	1ml	
Furfural-2,4-dinitrophenylhydrazone				
CAS 2074-02-4 DRE-C13972120	MW 276.205 Furfural-2,4-dinitrophenylhydrazone	C ₁₁ H ₈ N ₄ O ₅	100mg	










Hydrocarbons and Petrochemicals

Product code	Description		
Gasoline			
CAS 8032-32-4	MW n/a		No Structure
DRE-YA03001100ME	Regular Unleaded Gasoline, 2500 µg/mL in Methanol	1ml	
DRE-YA03001200ME	Regular Unleaded Gasoline unweathered 5000 µg/mL in Methanol	1ml	
DRE-YA03003100ME	Unleaded Super Gasoline Plus (unweathered) 5000 µg/mL in Methanol	1ml	
Natural Gasoline			
CAS 8006-61-9	MW n/a		No Structure
DRE-YS09010020ME	Gasoline (Unleaded) 20000 µg/mL in p/t Methanol(‡)	5x1ml	
Gasoline with Certified Content of BTEX			
DRE-YA03001700ME	Regular Unleaded Gasoline 10000 µg/mL in Methanol, with Certified Content of BTEX		1ml
	Benzene Ethylbenzene o-Xylene p-Xylene Toluene		
Glutaraldehyd-bis(2,4-dinitrophenylhydrazone)			
CAS 5085-07-4	MW 460.3577	C ₁₇ H ₁₆ N ₈ O ₈	1ml
DRE-XA14034480AL	Glutaraldehyd-bis(2,4-dinitrophenylhydrazone) 100 µg/mL in Acetonitrile		
			
Glycerol			
CAS 56-81-5	MW 92.0938	C ₃ H ₈ O ₃	1ml
DRE-C14036500	Glycerol		
			
n-Heneicosane			
CAS 629-94-7	MW 296.5741	C ₂₁ H ₄₄	20mg
DRE-C14085000	n-Heneicosane		
			
n-Hentriacontane			
CAS 630-04-6	MW 436.8399	C ₃₁ H ₆₄	20mg
DRE-C14085500	n-Hentriacontane		
			
n-Heptacosane			
CAS 593-49-7	MW 380.7335	C ₂₇ H ₅₆	20mg
DRE-C14122300	n-Heptacosane		
			
n-Heptadecane			
CAS 629-78-7	MW 240.4677	C ₁₇ H ₃₆	1ml
DRE-C14122500	n-Heptadecane		
			










Hydrocarbons and Petrochemicals

Product code	Description			
1-Heptadecanol				
CAS 1454-85-9 DRE-C14122700	MW 256.4671 1-Heptadecanol	C ₁₇ H ₃₆ O	100mg	
2,2,4,4,6,8,8-Heptamethylnonane				
CAS 4390-04-9 DRE-C14123500	MW 226.4412 2,2,4,4,6,8,8-Heptamethylnonane	C ₁₆ H ₃₄	250mg	
Heptanal				
CAS 111-71-7 DRE-CA14125000	MW 114.1855 1-Heptanal	C ₇ H ₁₄ O	250mg	
Heptanal-2,4-dinitrophenylhydrazone				
CAS 2074-05-7 DRE-XA14125010AL	MW 294.3064 Heptanal-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitrile	C ₁₃ H ₁₈ N ₄ O ₄	1ml	
n-Heptane				
CAS 142-82-5 DRE-C14126000	MW 100.2019 n-Heptane	C ₇ H ₁₆	5ml	
1-Heptanethiol				
CAS 1639-09-4 DRE-C14126600	MW 132.2669 1-Heptanethiol	C ₇ H ₁₆ S	250mg	
2-Heptanol				
CAS 543-49-7 DRE-C14127200	MW 116.2013 2-Heptanol	C ₇ H ₁₆ O	1ml	
1-Heptanol				
CAS 111-70-6 DRE-C14127100	MW 116.2013 1-Heptanol	C ₇ H ₁₆ O	1ml	
2-Heptanone				
CAS 110-43-0 DRE-C14127600	MW 114.1855 2-Heptanone	C ₇ H ₁₄ O	1ml	










Hydrocarbons and Petrochemicals

Product code	Description			
3-Heptanone				
CAS 106-35-4 DRE-C14127700	MW 114.1855 3-Heptanone	C ₇ H ₁₄ O	1ml	
4-Heptanone				
CAS 123-19-3 DRE-C14127800	MW 114.1855 4-Heptanone	C ₇ H ₁₄ O	1ml	
n-Heptatriacontane				
CAS 7194-84-5 DRE-C14128000	MW 520.9993 n-Heptatriacontane	C ₃₇ H ₇₆	100mg	
trans-2-Heptene				
CAS 14686-13-6 DRE-C14128300	MW 98.1861 trans-2-Heptene	C ₇ H ₁₄	100mg	
1-Heptylamine				
CAS 111-68-2 DRE-C14135000	MW 115.2166 Heptylamine	C ₇ H ₁₇ N	250mg	
n-Hexacontane				
CAS 7667-80-3 DRE-C14191000	MW 843.6107 n-Hexacontane	C ₆₀ H ₁₂₂	10mg	
n-Hexacosane				
CAS 630-01-3 DRE-C14191200	MW 366.707 n-Hexacosane	C ₂₆ H ₅₄	20mg	
1-Hexacosanol (Ceretyl alcohol)				
CAS 506-52-5 DRE-L14191400MB	MW 382.7064 1-Hexacosanol (Ceretyl alcohol) 10 µg/mL in Methyl-tert. butyl ether	C ₂₆ H ₅₄ O	10ml	
n-Hexadecane				
CAS 544-76-3 DRE-C14191500 DRE-L14191500AC	MW 226.4412 n-Hexadecane n-Hexadecane 10 µg/mL in Acetone	C ₁₆ H ₃₄	100mg 10ml	








Hydrocarbons and Petrochemicals

Product code	Description			
1-Hexadecanol (Cetyl alcohol)				
CAS 36653-82-4 DRE-C14192500	MW 242.4406 1-Hexadecanol (Cetyl alcohol)	C ₁₆ H ₃₄ O	250mg	
1-Hexadecene				
CAS 629-73-2 DRE-L14192800CY	MW 224.4253 1-Hexadecene 10 µg/mL in Cyclohexane	C ₁₆ H ₃₂	10ml	
Hexamethylbenzene				
CAS 87-85-4 DRE-C14194400	MW 162.2713 Hexamethylbenzene	C ₁₂ H ₁₈	250mg	
Hexanal (Capronaldehyde)				
CAS 66-25-1 DRE-CA14195000	MW 100.1589 Hexanal (Capronaldehyde)	C ₆ H ₁₂ O	1ml	
Hexanal-2,4-dinitrophenylhydrazone				
CAS 1527-97-5 DRE-XA14195010AL	MW 280.2798 Hexanal-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitril	C ₁₂ H ₁₆ N ₄ O ₄	1ml	
n-Hexane				
CAS 110-54-3 DRE-C14195500	MW 86.1754 n-Hexane	C ₆ H ₁₄	5ml	
2,5-Hexanediol				
CAS 2935-44-6 DRE-C14195700	MW 118.1742 2,5-Hexanediol	C ₆ H ₁₄ O ₂	250mg	
1-Hexanol				
CAS 111-27-3 DRE-C14196700	MW 102.1748 1-Hexanol	C ₆ H ₁₄ O	1ml	
2-Hexanol				
CAS 626-93-7 DRE-C14196800	MW 102.1748 2-Hexanol	C ₆ H ₁₄ O	1ml	







Hydrocarbons and Petrochemicals

Product code	Description			
3-Hexanol				
CAS 623-37-0 DRE-C14196900	MW 102.1748 3-Hexanol	C ₆ H ₁₄ O	1ml	
2-Hexanone				
CAS 591-78-6 DRE-C14197500	MW 100.1589 2-Hexanone	C ₆ H ₁₂ O	1ml	
n-Hexatriacontane				
CAS 630-06-8 DRE-C14199000	MW 506.9728 n-Hexatriacontane	C ₃₆ H ₇₄	100mg	
1-Hexene				
CAS 592-41-6 DRE-CA14202100	MW 84.1595 1-Hexene	C ₆ H ₁₂	250mg	
1-Hexylamine				
CAS 111-26-2 DRE-C14206000	MW 101.19 Hexylamine	C ₆ H ₁₅ N	250mg	
n-Hexylbenzene				
CAS 1077-16-3 DRE-C14207000	MW 162.2713 n-Hexylbenzene	C ₁₂ H ₁₈	100mg	
4-Hydroxy-4-methyl-2-pentanone				
CAS 123-42-2 DRE-C14233000	MW 116.1583 4-Hydroxy-4-methyl-2-pentanone	C ₆ H ₁₂ O ₂	5ml	
Indane				
CAS 496-11-7 DRE-C14287500 DRE-L14287500ME	MW 118.1757 Indan Indan 10 µg/mL in Methanol	C ₉ H ₁₀	100mg 10ml	
Indene				
CAS 95-13-6 DRE-C14288500 DRE-L14288500ME	MW 116.1598 Indene Indene 10 µg/mL in Methanol	C ₉ H ₈	100mg 10ml	









Hydrocarbons and Petrochemicals

Product code	Description			
Isobutylbenzene				
CAS 538-93-2 DRE-C14394600	MW 134.2182 Isobutylbenzene	C ₁₀ H ₁₄	1ml	
Isobutyraldehyde				
CAS 78-84-2 DRE-CA14395000	MW 72.1057 Isobutyraldehyde	C ₄ H ₈ O	1ml	
Isodecanol-ethoxylate (technical)				
CAS 61827-42-7 DRE-C14405500	MW 202.3336 Isodecanol-ethoxylate (technical)	C ₁₀ H ₂₂ O(C ₂ H ₄ O) _n	100mg	
Isododecane				
CAS 31807-55-3 DRE-CA14408000	MW 170.3348 Isododecane	C ₁₂ H ₂₆	100mg	
Isophorone				
CAS 78-59-1 DRE-C14446000	MW 138.2069 Isophorone(‡)	C ₉ H ₁₄ O	1ml	
Isoprene				
CAS 78-79-5 DRE-C14449000	MW 68.117 Isoprene	C ₅ H ₈	1ml	
Isopropylbenzene (Cumene)				
CAS 98-82-8 DRE-C14463500 DRE-L14463500ME DRE-XA14463500ME	MW 120.1916 Isopropylbenzene (Cumene) Isopropylbenzene 10 µg/mL in Methanol Isopropylbenzene 100 µg/mL in Methanol	C ₉ H ₁₂	500mg 10ml 1ml	
Isopropylcyclohexane				
CAS 696-29-7 DRE-C14463600	MW 126.2392 Isopropylcyclohexane	C ₉ H ₁₈	250mg	
2-Isopropyltoluene (o-Cymol)				
CAS 527-84-4 DRE-CA14465480	MW 134.2182 2-Isopropyltoluene (o-Cymol)	C ₁₀ H ₁₄	10mg	






Hydrocarbons and Petrochemicals

Product code	Description			
3-Isopropyltoluene (m-Cymol)				
CAS 535-77-3 DRE-CA14465490	MW 134.2182 3-Isopropyltoluene (m-Cymol)	C ₁₀ H ₁₄	10mg	
4-Isopropyltoluene (p-Cymol)				
CAS 99-87-6 DRE-CA14465500	MW 134.2182 4-Isopropyltoluene	C ₁₀ H ₁₄	1ml	
Isovaleraldehyd-2,4-dinitrophenylhydrazone				
CAS 2256-01-1 DRE-C14479450	MW 266.2533 Isovaleraldehyd-2,4-dinitrophenylhydrazone	C ₁₁ H ₁₄ N ₄ O ₄	100mg	
Malonitrile				
CAS 109-77-3 DRE-C14733500	MW 66.0614 Malonitrile	C ₃ H ₂ N ₂	250mg	
Methacrylaldehyd-2,4-dinitrophenylhydrazone				
CAS 5077-73-6 DRE-XA14971400AL	MW 250.2108 Methacrylaldehyd-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitril	C ₁₀ H ₁₀ N ₄ O ₄	1ml	
Methacrylic Acid Ethyl Ester				
CAS 97-63-2 DRE-CA14971740	MW 114.1424 Methacrylic acid-ethyl ester	C ₆ H ₁₀ O ₂	1ml	
Methacrylonitrile				
CAS 126-98-7 DRE-CA14972000	MW 67.0892 Methacrylonitrile	C ₄ H ₅ N	1ml	
Methane				
CAS 74-82-8 DRE-GA09010081DS	MW 16.0425 Methane 500 µg/mL in Triacetin(‡)	CH ₄	1ml	
Methanol				
CAS 67-56-1 DRE-C14995000	MW 32.0419 Methanol(‡)	CH ₄ O	1ml	









Hydrocarbons and Petrochemicals

Product code	Description			
2-Methoxyethanol (Ethyleneglycol monomethyl ether)				
CAS 109-86-4 DRE-C13328200 DRE-GA09010378ME	MW 76.0944 Ethyleneglycol-monomethyl ether 2-Methoxyethanol (SS) 500 µg/mL in p/t Methanol(±)	C ₃ H ₈ O ₂	250mg 1ml	
2-Methoxyethyl Acetate				
CAS 110-49-6 DRE-CA15077000	MW 118.1311 2-Methoxyethyl acetate	C ₅ H ₁₀ O ₃	1ml	
1-Methoxy-2-propanol				
CAS 107-98-2 DRE-C15083000	MW 90.121 1-Methoxy-2-propanol	C ₄ H ₁₀ O ₂	500mg	
Methyl Methacrylate (Methacrylic acid methyl ester)				
CAS 80-62-6 DRE-CA14971770 DRE-XA14971770CY	MW 100.1158 Methacrylic acid-methyl ester Methacrylic acid-methyl ester 100 µg/mL in Cyclohexane	C ₅ H ₈ O ₂	250mg 1ml	
2-Methyl-1-butanol				
CAS 137-32-6 DRE-C15084400	MW 88.1482 2-Methyl-1-butanol	C ₅ H ₁₂ O	1ml	
2-Methyl-2-butanol				
CAS 75-85-4 DRE-C15084410	MW 88.1482 2-Methyl-2-butanol	C ₅ H ₁₂ O	1ml	
3-Methyl-1-butanol				
CAS 123-51-3 DRE-C15084440	MW 88.1482 3-Methyl-1-butanol	C ₅ H ₁₂ O	1ml	
3-Methyl-2-butanone				
CAS 563-80-4 DRE-C15084470	MW 86.1323 3-Methyl-2-butanone	C ₅ H ₁₀ O	250mg	
3-Methylbutylamine				
CAS 107-85-7 DRE-C15084500	MW 87.1634 3-Methylbutylamine	C ₅ H ₁₃ N	250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
Methyl-tert-butylether				
CAS 1634-04-4	MW 88.1482	C ₅ H ₁₂ O		
DRE-C15084600	Methyl-tert-butylether		5ml	
DRE-L15084600ME	Methyl-tert-butylether 10 µg/mL in Methanol		10ml	
DRE-XA15084600ME	Methyl-tert-butylether 100 µg/mL in Methanol		1ml	
3-Methylbutyraldehyde				
CAS 590-86-3	MW 86.1323	C ₅ H ₁₀ O		
DRE-C15084700	3-Methylbutyraldehyde		250mg	
2-Methylbutyric Acid Methyl Ester				
CAS 868-57-5	MW 116.1583	C ₆ H ₁₂ O ₂		
DRE-CA15084850	2-Methylbutyric acid-methyl ester		1ml	
Methylcyclohexane				
CAS 108-87-2	MW 98.1861	C ₇ H ₁₄		
DRE-C15085000	Methylcyclohexane		1ml	
DRE-YS09010024ME	Methylcyclohexane 1000 µg/mL in p/t Methanol(‡)		5x1ml	
2-Methylcyclohexanone				
CAS 583-60-8	MW 112.1696	C ₇ H ₁₂ O		
DRE-C15085010	2-Methylcyclohexanone		250mg	
3-Methylcyclohexanone				
CAS 591-24-2	MW 112.1696	C ₇ H ₁₂ O		
DRE-C15085020	3-Methylcyclohexanone		250mg	
4-Methylcyclohexanone				
CAS 589-92-4	MW 112.1696	C ₇ H ₁₂ O		
DRE-C15085030	4-Methylcyclohexanone		250mg	
3-Methylcyclopentadecanone				
CAS 541-91-3	MW 238.4088	C ₁₆ H ₃₀ O		
DRE-XA15085033CY	3-Methylcyclopentadecanone 100 µg/mL in Cyclohexane		1ml	
Methylcyclopentane				
CAS 96-37-7	MW 84.1595	C ₆ H ₁₂		
DRE-CA15085040	Methylcyclopentane		1ml	








Hydrocarbons and Petrochemicals

Product code	Description			
3-Methylfuran				
CAS 930-27-8 DRE-XA15086070ME	MW 82.1005 3-Methylfuran 100 µg/mL in Methanol	C ₅ H ₆ O	1ml	
3-Methylheptane				
CAS 589-81-1 DRE-C15087800	MW 114.2285 3-Methylheptane	C ₈ H ₁₈	100mg	
2-Methyl-3-heptanone				
CAS 13019-20-0 DRE-C15087900	MW 128.212 2-Methyl-3-heptanone	C ₈ H ₁₆ O	250mg	
5-Methyl-2-hexanone				
CAS 110-12-3 DRE-C15088200	MW 114.1855 5-Methyl-2-hexanone	C ₇ H ₁₄ O	250mg	
2-Methylpentane				
CAS 107-83-5 DRE-CA15121200	MW 86.1754 2-Methylpentane	C ₆ H ₁₄	1ml	
3-Methylpentane				
CAS 96-14-0 DRE-CA15121300	MW 86.1754 3-Methylpentane	C ₆ H ₁₄	1ml	
2-Methyl-1-pentanol				
CAS 105-30-6 DRE-C15121800	MW 102.1748 2-Methyl-1-pentanol	C ₆ H ₁₄ O	250mg	
3-Methyl-3-pentanol				
CAS 77-74-7 DRE-C15122400	MW 102.1748 3-Methyl-3-pentanol	C ₆ H ₁₄ O	250mg	
4-Methyl-pentan-2-ol (2-Methyl-4-pentanol)				
CAS 108-11-2 DRE-C15122100	MW 102.1748 2-Methyl-4-pentanol(‡)	C ₆ H ₁₄ O	1ml	










Hydrocarbons and Petrochemicals

Product code	Description			
4-Methyl-2-pentanone-2,4-dinitrophenylhydrazone				
CAS 1655-42-1 DRE-YA15123010AL	MW 280.2798 4-Methyl-2-pentanone-2,4-dinitrophenylhydrazone 1000 µg/mL in Acetonitrile	C ₁₂ H ₁₆ N ₄ O ₄	1ml	
4-Methyl-1-pentene				
CAS 691-37-2 DRE-CA15124400	MW 84.1595 4-Methyl-1-pentene	C ₆ H ₁₂	1ml	
cis-4-Methyl-2-pentene				
CAS 691-38-3 DRE-CA15124700	MW 84.1595 cis-4-Methyl-2-pentene	C ₆ H ₁₂	100mg	
2-(4'-Methylphenyl)propionic Acid 3,3,3-D3				
CAS n/a DRE-C15140801	MW 167.2196 2-(4'-Methylphenyl)propionic acid D3	C ₁₀ ² H ₁₃ O ₂	10mg	
2-(4'-Methylphenyl)propionic Acid ((2RS)-2-(4-Methylphenyl)propanoic acid)				
CAS 938-94-3 DRE-C15140800	MW 164.2011 2-(4'-Methylphenyl)propionic acid	C ₁₀ H ₁₂ O ₂	100mg	
2-Methylpropanol (2-Methyl-1-propanol)				
CAS 78-83-1 DRE-C15142000 DRE-L15142000ME	MW 74.1216 2-Methyl-1-propanol(‡) 2-Methyl-1-propanol 10 µg/mL in Methanol	C ₄ H ₁₀ O	1g 10ml	
α-Methylstyrene				
CAS 98-83-9 DRE-CA15143500	MW 118.1757 α-Methylstyrene	C ₉ H ₁₀	1ml	
trans-β-Methylstyrene				
CAS 873-66-5 DRE-CA15143600	MW 118.1757 trans-β-Methylstyrene	C ₉ H ₁₀	1ml	
4-Methylvaleric Acid				
CAS 646-07-1 DRE-CA15147400	MW 116.1583 4-Methylvaleric acid	C ₈ H ₁₂ O ₂	250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
Mineral Oil				
CAS 8042-47-5	MW n/a			No Structure
DRE-CA03009010	Mineral Oil (without additives, DIN H53) (EN9377-2)		1ml	
DRE-CA03009011	Mineral Oil Heavy (without additives, DIN H53) (EN9377-2)		1ml	
Naphthenic acids				
CAS 1338-24-5	MW n/a			No Structure
DRE-C15423000	Naphthenic acid (technical mixture)		250mg	
n-Nonacosane				
CAS 630-03-5	MW 408.7867	C ₂₉ H ₆₀		
DRE-C15621000	n-Nonacosane(‡)		50mg	
n-Nonadecane				
CAS 629-92-5	MW 268.5209	C ₁₉ H ₄₀		
DRE-C15622000	n-Nonadecane		250mg	
Nonanal-2,4-dinitrophenylhydrazone				
CAS 2348-19-8	MW 322.3596	C ₁₅ H ₂₂ N ₄ O ₄		
DRE-XA15622910AL	Nonanal-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitrile		1ml	
n-Nonane				
CAS 111-84-2	MW 128.2551	C ₉ H ₂₀		
DRE-C15623000	n-Nonane		5ml	
DRE-L15623000CY	n-Nonane 10 µg/mL in Cyclohexane		10ml	
DRE-GA09010086ME	n-Nonane (C9) 1000 µg/mL in p/t Methanol(‡)		1ml	
1-Nonanol				
CAS 143-08-8	MW 144.2545	C ₉ H ₂₀ O		
DRE-C15623200	1-Nonanol		1ml	
2-Nonanone				
CAS 821-55-6	MW 142.2386	C ₉ H ₁₈ O		
DRE-C15623220	2-Nonanone		1ml	
n-Nonatriacontane				
CAS 7194-86-7	MW 549.0525	C ₃₉ H ₈₀		
DRE-C15624500	n-Nonatriacontane		20mg	










Hydrocarbons and Petrochemicals

Product code	Description			
n-Nonylbenzene (Phenylnonane)				
CAS 1081-77-2 DRE-C15628000	MW 204.3511 n-Nonylbenzene (Phenylnonane)	C ₁₉ H ₂₄	100mg	
n-Octacosane				
CAS 630-02-4 DRE-C15710050 DRE-GA09010088ME DRE-GA09010088ME-1	MW 394.7601 n-Octacosane n-Octacosane (C28) 1000 µg/mL in p/t Methanol(±) Octacosane 1000 µg/mL in Dichloromethane(±)	C ₂₈ H ₅₈	250mg 1ml 1ml	
1-Octacosanol (Montanyl alcohol)				
CAS 557-61-9 DRE-L15710090MB	MW 410.7595 1-Octacosanol (Montanyl alcohol) 10 µg/mL in Methyl-tert. butyl ether	C ₂₈ H ₅₈ O	10ml	
n-Octadecane				
CAS 593-45-3 DRE-C15710100	MW 254.4943 n-Octadecane	C ₁₈ H ₃₈	0.5g	
1-Octadecanol (Stearyl alcohol)				
CAS 112-92-5 DRE-C15710300	MW 270.4937 1-Octadecanol (Stearyl alcohol)	C ₁₈ H ₃₈ O	250mg	
1-Octadecene				
CAS 112-88-9 DRE-C15710400	MW 252.4784 1-Octadecene	C ₁₈ H ₃₆	250mg	
Octanal				
CAS 124-13-0 DRE-C15710900	MW 128.212 Octanal	C ₈ H ₁₆ O	250mg	
Octanal-2,4-dinitrophenylhydrazone				
CAS 1726-77-8 DRE-XA15710910AL	MW 308.333 Octanal-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitrile	C ₁₄ H ₂₀ N ₄ O ₄	1ml	
n-Octane				
CAS 111-65-9 DRE-C15711000	MW 114.2285 n-Octane	C ₈ H ₁₈	5ml	










Hydrocarbons and Petrochemicals

Product code	Description			
1-Octanol				
CAS 111-87-5 DRE-C15711100	MW 130.2279 1-Octanol	C ₈ H ₁₈ O	1ml	
2-Octanol				
CAS 123-96-6 DRE-C15711200	MW 130.2279 2-Octanol	C ₈ H ₁₈ O	1ml	
2-Octanone				
CAS 111-13-7 DRE-C15711250	MW 128.212 2-Octanone	C ₈ H ₁₆ O	1ml	
3-Octanone				
CAS 106-68-3 DRE-C15711260	MW 128.212 3-Octanone	C ₈ H ₁₆ O	250mg	
n-Octatriacontane				
CAS 7194-85-6 DRE-C15711280	MW 535.0259 n-Octatriacontane	C ₃₈ H ₇₈	20mg	
1-Octene				
CAS 111-66-0 DRE-CA15711300	MW 112.2126 1-Octene	C ₈ H ₁₆	1ml	
trans-2-Octene				
CAS 13389-42-9 DRE-CA15711400	MW 112.2126 trans-2-Octene	C ₈ H ₁₆	100mg	
1-Octylamine				
CAS 111-86-4 DRE-C15711700	MW 129.2432 Octylamine	C ₈ H ₁₉ N	1ml	
Octylbenzene				
CAS 2189-60-8 DRE-C15711800	MW 190.3245 Octylbenzene	C ₁₄ H ₂₂	100mg	










Hydrocarbons and Petrochemicals

Product code	Description			
Paromomycin Sulfate				
CAS 1263-89-4 DRE-C15893500	MW 713.707 Paromomycin sulfate	C ₂₃ H ₄₅ N ₅ O ₁₄ ·H ₂ O ₄ S	100mg	
n-Pentacosane				
CAS 629-99-2 DRE-C15973500 DRE-GA09010093DI	MW 352.6804 n-Pentacosane Pentacosane 10000 µg/mL in Dichloromethane(±)	C ₂₅ H ₅₂	20mg 1ml	
n-Pentadecane				
CAS 629-62-9 DRE-C15973700	MW 212.4146 n-Pentadecane	C ₁₅ H ₃₂	1ml	
1-Pentadecanol				
CAS 629-76-5 DRE-C15973800	MW 228.414 1-Pentadecanol	C ₁₅ H ₃₂ O	250mg	
2-Pentadecanone				
CAS 2345-28-0 DRE-C15973802	MW 226.3981 2-Pentadecanone	C ₁₅ H ₃₀ O	100mg	
Pentamethylbenzene				
CAS 700-12-9 DRE-C15975000	MW 148.2447 Pentamethylbenzene	C ₁₁ H ₁₆	500mg	
2,2,4,6,6-Pentamethylheptane				
CAS 13475-82-6 DRE-CA15975500	MW 170.3348 2,2,4,6,6-Pentamethylheptane	C ₁₂ H ₂₆	100mg	
Pentanal-2,4-dinitrophenylhydrazone				
CAS 2057-84-3 DRE-C15977400	MW 266.2533 Pentanal-2,4-dinitrophenylhydrazone	C ₁₁ H ₁₄ N ₄ O ₄	100mg	
n-Pentane				
CAS 109-66-0 DRE-C15977500 DRE-YA09010026ME DRE-YS09010026ME	MW 72.1488 n-Pentane n-Pentane 1000 µg/mL in p/t Methanol(±) Pentane 1000 µg/mL in p/t Methanol(±)	C ₅ H ₁₂	5ml 1ml 5x1ml	





Hydrocarbons and Petrochemicals

Product code	Description			
1,5-Pentanediol				
CAS 111-29-5 DRE-C15977700	MW 104.1476 1,5-Pentanediol	C ₅ H ₁₂ O ₂	250mg	
1-Pentanol				
CAS 71-41-0 DRE-C15981100	MW 88.1482 1-Pentanol	C ₅ H ₁₂ O	5ml	
2-Pentanol				
CAS 6032-29-7 DRE-C15981200	MW 88.1482 2-Pentanol	C ₅ H ₁₂ O	5ml	
3-Pentanol				
CAS 584-02-1 DRE-C15981300	MW 88.1482 3-Pentanol(†)	C ₅ H ₁₂ O	1ml	
2-Pentanone				
CAS 107-87-9 DRE-C15981500	MW 86.1323 2-Pentanone	C ₅ H ₁₀ O	1ml	
3-Pentanone				
CAS 96-22-0 DRE-C15981600	MW 86.1323 3-Pentanone	C ₅ H ₁₀ O	1ml	
n-Pentatriacontane				
CAS 630-07-9 DRE-C15981700	MW 492.9462 n-Pentatriacontane	C ₃₅ H ₇₂	20mg	
n-Pentylbenzene (n-Amylbenzene)				
CAS 538-68-1 DRE-C15983100	MW 148.2447 n-Pentylbenzene (n-Amylbenzene)	C ₁₁ H ₁₆	250mg	
sec-Pentylbenzene (sec-Amylbenzene)				
CAS 2719-52-0 DRE-C15983200	MW 148.2447 sec-Pentylbenzene (sec-Amylbenzene)	C ₁₁ H ₁₆	250mg	







Hydrocarbons and Petrochemicals

Product code	Description			
Phenoxyethanol (Ethyleneglycol monophenyl ether)				
CAS 122-99-6 DRE-C13328300	MW 138.1638 Ethyleneglycol-monophenyl ether	C ₈ H ₁₀ O ₂	250mg	
1-Phenoxy-2-propanol				
CAS 770-35-4 DRE-C16045600	MW 152.1904 1-Phenoxy-2-propanol	C ₉ H ₁₂ O ₂	250mg	
4-Phenyl-1-cyclohexene				
CAS 4994-16-5 DRE-XA16057000CY	MW 158.2396 4-Phenyl-1-cyclohexene 100 µg/mL in Cyclohexane(†)	C ₁₂ H ₁₄	1ml	
2-Phenylethylamine				
CAS 64-04-0 DRE-CA16058500	MW 121.1796 2-Phenylethylamine	C ₈ H ₁₁ N	1ml	
2-Phenyl-2-propanol				
CAS 617-94-7 DRE-C16073000	MW 136.191 2-Phenyl-2-propanol	C ₉ H ₁₂ O	25mg	
Phytane				
CAS 638-36-8 DRE-L16193200IO	MW 282.5475 Phytane 10 µg/mL in Isooctane	C ₂₀ H ₄₂	10ml	
Pristane				
CAS 1921-70-6 DRE-C16288100	MW 268.5209 Pristane	C ₁₉ H ₄₀	100mg	
Propanal (Propionaldehyde)				
CAS 123-38-6 DRE-C16492500	MW 58.0791 Propionaldehyde (Propanal)	C ₃ H ₆ O	250mg	
1,3-Propanediol				
CAS 504-63-2 DRE-C16405300	MW 76.0944 1,3-Propanediol	C ₃ H ₈ O ₂	250mg	









Hydrocarbons and Petrochemicals

Product code	Description			
1-Propanethiol				
CAS 107-03-9 DRE-CA16406100	MW 76.1606 1-Propanethiol	C ₃ H ₈ S	250mg	
2-Propanethiol				
CAS 75-33-2 DRE-CA16406200	MW 76.1606 2-Propanethiol	C ₃ H ₈ S	250mg	
1-Propanol				
CAS 71-23-8 DRE-C16415100	MW 60.095 1-Propanol(‡)	C ₃ H ₈ O	1ml	
2-Propanol (Isopropyl alcohol)				
CAS 67-63-0 DRE-C16415200	MW 60.095 2-Propanol(‡)	C ₃ H ₈ O	1ml	
Propargyl Alcohol				
CAS 107-19-7 DRE-C16433000	MW 56.0633 Propargyl alcohol	C ₃ H ₄ O	250mg	
Propionaldehyde-2,4-dinitrophenylhydrazone				
CAS 725-00-8 DRE-C16492505	MW 238.2001 Propionaldehyde-2,4-dinitrophenylhydrazone	C ₉ H ₁₀ N ₄ O ₄	100mg	
Propionitrile				
CAS 107-12-0 DRE-C16494000	MW 55.0785 Propionitrile	C ₃ H ₅ N	1ml	
1-Propylamine (n-Propylamine)				
CAS 107-10-8 DRE-C16508000	MW 59.1103 n-Propylamine	C ₃ H ₉ N	250mg	
n-Propylbenzene				
CAS 103-65-1 DRE-C16519500 DRE-XA16519500ME	MW 120.1916 n-Propylbenzene n-Propylbenzene 100 µg/mL in Methanol	C ₉ H ₁₂	1ml 1ml	










Hydrocarbons and Petrochemicals

Product code	Description			
Propylcyclohexane				
CAS 1678-92-8 DRE-C16525000	MW 126.2392 Propylcyclohexane	C ₉ H ₁₈	250mg	
Propylene Glycol (Propane-1,2-diol)				
CAS 57-55-6 DRE-C16405200	MW 76.0944 1,2-Propanediol(‡)	C ₃ H ₈ O ₂	1ml	
Propylene Glycol 1-Methyl Ether 2-Acetate				
CAS 108-65-6 DRE-C16527500	MW 132.1577 Propylene glycol 1-methyl ether 2-acetate	C ₆ H ₁₂ O ₃	250mg	
Propylene Oxide				
CAS 75-56-9 DRE-C16528000	MW 58.0791 Propylene oxide	C ₃ H ₆ O	1ml	
1,2-Propylene carbonate				
CAS 108-32-7 DRE-C16527000	MW 102.0886 1,2-Propylene carbonate	C ₄ H ₆ O ₃	5ml	
Pyridine				
CAS 110-86-1 DRE-C16646000	MW 79.0999 Pyridine(‡)	C ₅ H ₅ N	5ml	
Pyridine-d5				
CAS 7291-22-7 DRE-C16646100	MW 84.1307 Pyridine D5	C ₅ ² H ₅ N	1ml	
Pyrrolidine				
CAS 123-75-1 DRE-C16675000	MW 71.121 Pyrrolidine	C ₄ H ₉ N	1ml	
Squalane				
CAS 111-01-3 DRE-C16973500	MW 422.8133 Squalane	C ₃₀ H ₆₂	100mg	










Hydrocarbons and Petrochemicals

Product code	Description			
Styrene				
CAS 100-42-5	MW 104.1491	C ₈ H ₈		
DRE-C16982000	Styrene		1ml	
DRE-XA16982000ME	Styrene 100 µg/mL in Methanol		1ml	
DRE-YA16982000ME	Styrene 5000 µg/mL in Methanol		1ml	
Styrene D8				
CAS 19361-62-7	MW 112.1984	C ₈ ² H ₈		
DRE-C16982010	Styrene D8		100mg	
Sulfolan				
CAS 126-33-0	MW 120.1701	C ₄ H ₈ O ₂ S		
DRE-C17008500	Sulfolan		5ml	
m-Terphenyl				
CAS 92-06-8	MW 230.3038	C ₁₈ H ₁₄		
DRE-C20934900	m-Terphenyl		100mg	
o-Terphenyl				
CAS 84-15-1	MW 230.3038	C ₁₈ H ₁₄		
DRE-C20934800	o-Terphenyl		100mg	
DRE-YS09010012DI	o-Terphenyl 10000 µg/mL in Dichloromethane(‡)		5x1ml	
p-Terphenyl				
CAS 92-94-4	MW 230.3038	C ₁₈ H ₁₄		
DRE-C20935000	p-Terphenyl		100mg	
p-Terphenyl D14				
CAS 1718-51-0	MW 244.39	C ₁₈ ² H ₁₄		
DRE-C20935300	p-Terphenyl D14		10mg	
DRE-GA09010298DI	p-Terphenyl-d14 500 µg/mL in Dichloromethane(‡)		1ml	
n-Tetracontane				
CAS 4181-95-7	MW 563.0791	C ₄₀ H ₈₂		
DRE-C17395500	n-Tetracontane(‡)		100mg	
n-Tetracosane				
CAS 646-31-1	MW 338.6538	C ₂₄ H ₅₀		
DRE-C17395800	n-Tetracosane		250mg	










Hydrocarbons and Petrochemicals

Product code	Description			
n-Tetradecane				
CAS 629-59-4	MW 198.388	C ₁₄ H ₃₀		
DRE-C17396500	n-Tetradecane		5ml	
DRE-L17396500HE	n-Tetradecane 10 µg/mL in n-Hexane		10ml	
1-Tetradecanol (Myristyl alcohol)				
CAS 112-72-1	MW 214.3874	C ₁₄ H ₃₀ O		
DRE-C17396800	1-Tetradecanol (Myristyl alcohol)		250mg	
1-Tetradecene				
CAS 1120-36-1	MW 196.3721	C ₁₄ H ₂₈		
DRE-C17397100	1-Tetradecene		1ml	
Tetradecylamine				
CAS 2016-42-4	MW 213.4026	C ₁₄ H ₃₁ N		
DRE-C17397500	Tetradecylamine		250mg	
2-Tetradecylcyclobutanone				
CAS 35493-47-1	MW 266.462	C ₁₈ H ₃₄ O		
DRE-C17397700	2-Tetradecylcyclobutanone		10mg	
Tetrahydrofuran (THF)				
CAS 109-99-9	MW 72.1057	C ₄ H ₈ O		
DRE-C17405700	Tetrahydrofuran (THF)		5ml	
1,2,3,4-Tetramethylbenzene				
CAS 488-23-3	MW 134.2182	C ₁₀ H ₁₄		
DRE-CA17412300	1,2,3,4-Tetramethylbenzene		100mg	
DRE-XA17412300ME	1,2,3,4-Tetramethylbenzene 100 µg/mL in Methanol		1ml	
1,2,3,5-Tetramethylbenzene (Isodurene)				
CAS 527-53-7	MW 134.2182	C ₁₀ H ₁₄		
DRE-C17412400	1,2,3,5-Tetramethylbenzene (Isodurene)		250mg	
1,2,4,5-Tetramethylbenzene (Durene)				
CAS 95-93-2	MW 134.2182	C ₁₀ H ₁₄		
DRE-C17412500	1,2,4,5-Tetramethylbenzene (Durene)		250mg	
DRE-L17412500CY	1,2,4,5-Tetramethylbenzene (Durene) 10 µg/mL in Cyclohexane		10ml	







Hydrocarbons and Petrochemicals

Product code	Description			
n-Tetratetracontane				
CAS 7098-22-8	MW 619.1854	C ₄₄ H ₉₀		
DRE-C17430000	n-Tetratetracontane		20mg	
DRE-YS09010017CS	Tetratetracontane (C44) 1000 µg/mL in Carbon Disulfide(±)		5x1ml	
n-Tetratriacontane				
CAS 14167-59-0	MW 478.9196	C ₃₄ H ₇₀		
DRE-C17430500	n-Tetratriacontane		20mg	
Thiophene				
CAS 110-02-1	MW 84.1396	C ₄ H ₄ S		
DRE-CA17547000	Thiophene		1ml	
Toluene (Methylbenzene)				
CAS 108-88-3	MW 92.1384	C ₇ H ₈		
DRE-C17594000	Toluene (Methylbenzene) DRUG PRECURSOR		5ml	
DRE-L17594000ME	Toluene 10 µg/mL in Methanol		10ml	
DRE-XA17594000ME	Toluene 100 µg/mL in Methanol		1ml	
Toluene D8				
CAS 2037-26-5	MW 100.1877	C ₇ ² H ₈		
DRE-C17594100	Toluene D8		.5ml	
n-Triacontane				
CAS 638-68-6	MW 422.8133	C ₃₀ H ₆₂		
DRE-C17609000	n-Triacontane		100mg	
Tributylamine				
CAS 102-82-9	MW 185.3495	C ₁₂ H ₂₇ N		
DRE-C17667500	Tributylamine		250mg	
n-Tricosane				
CAS 638-67-5	MW 324.6272	C ₂₃ H ₄₈		
DRE-C17805000	n-Tricosane		100mg	
1-Tricosanol				
CAS 3133-01-5	MW 340.6266	C ₂₃ H ₄₈ O		
DRE-L17806000MB	1-Tricosanol 10 µg/mL in Methyl-tert. butyl ether		10ml	

Hydrocarbons and Petrochemicals

Product code	Description			
n-Tridecane				
CAS 629-50-5 DRE-C17818000	MW 184.3614 n-Tridecane	C ₁₃ H ₂₈	1ml	
1-Tridecanol				
CAS 112-70-9 DRE-C17818300	MW 200.3608 1-Tridecanol	C ₁₃ H ₂₈ O	250mg	
2-Tridecanone				
CAS 593-08-8 DRE-C17818500	MW 198.3449 2-Tridecanone	C ₁₃ H ₂₆ O	250mg	
Triethylamine				
CAS 121-44-8 DRE-CA17832000	MW 101.19 Triethylamine	C ₆ H ₁₅ N	250mg	
1,3,5-Triethylbenzene				
CAS 102-25-0 DRE-C17833500	MW 162.2713 1,3,5-Triethylbenzene	C ₁₂ H ₁₈	100mg	
1,2,3-Trimethylbenzene				
CAS 526-73-8 DRE-C17880600 DRE-L17880600ME DRE-XA17880600ME	MW 120.1916 1,2,3-Trimethylbenzene 1,2,3-Trimethylbenzene 10 µg/mL in Methanol 1,2,3-Trimethylbenzene 100 µg/mL in Methanol	C ₉ H ₁₂	250mg 10ml 1ml	
1,2,4-Trimethylbenzene				
CAS 95-63-6 DRE-CA17880800 DRE-XA17880800ME	MW 120.1916 1,2,4-Trimethylbenzene 1,2,4-Trimethylbenzene 100 µg/mL in Methanol	C ₉ H ₁₂	1ml 1ml	
1,3,5-Trimethylbenzene (Mesitylene)				
CAS 108-67-8 DRE-C17881000 DRE-XA17881000ME	MW 120.1916 1,3,5-Trimethylbenzene (Mesitylene)(‡) 1,3,5-Trimethylbenzene 100 µg/mL in Methanol	C ₉ H ₁₂	250mg 1ml	
2,2,4-Trimethylpentane (Isooctane)				
CAS 540-84-1 DRE-C17883000 DRE-YS09010023ME	MW 114.2285 2,2,4-Trimethylpentane Isooctane 1000 µg/mL in p/t Methanol(‡)	C ₈ H ₁₈	5ml 5x1ml	

Hydrocarbons and Petrochemicals

Product code	Description			
2,4,4-Trimethyl-1-pentene				
CAS 107-39-1	MW 112.2126	C ₈ H ₁₆		
DRE-C17883200	2,4,4-Trimethyl-1-pentene		1ml	
1,3,5-Triphenylbenzene				
CAS 612-71-5	MW 306.3997	C ₂₄ H ₁₈		
DRE-C20500000	1,3,5-Triphenylbenzene		100mg	
DRE-L20500000AL	1,3,5-Triphenylbenzene 10 µg/mL in Acetonitrile		10ml	
n-Tritriacontane				
CAS 630-05-7	MW 464.893	C ₃₃ H ₆₈		
DRE-C17894800	n-Tritriacontane		20mg	
Undecanal				
CAS 112-44-7	MW 170.2918	C ₁₁ H ₂₂ O		
DRE-C17896200	Undecanal		250mg	
n-Undecane				
CAS 1120-21-4	MW 156.3083	C ₁₁ H ₂₄		
DRE-C17896300	n-Undecane		1ml	
2-Undecanone				
CAS 112-12-9	MW 170.2918	C ₁₁ H ₂₂ O		
DRE-C17896600	2-Undecanone		1ml	
1-Undecene				
CAS 821-95-4	MW 154.2924	C ₁₁ H ₂₂		
DRE-C17896700	1-Undecene		100mg	
Vinyl Acetate				
CAS 108-05-4	MW 86.0892	C ₄ H ₆ O ₂		
DRE-C17922500	Vinylacetate		1ml	
DRE-GA09010383ME	Vinyl Acetate Solution, 2000 µg/mL in p/t Methanol(‡)		1ml	
4-Vinyl-1-cyclohexene				
CAS 100-40-3	MW 108.1809	C ₈ H ₁₂		
DRE-CA17923100	4-Vinyl-1-cyclohexene		1ml	
DRE-YA17923100ME	4-Vinyl-1-cyclohexene 2000 µg/mL in Methanol		1ml	

Hydrocarbons and Petrochemicals

Product code	Description	
AEOS Residual Solvents Mix #1		
DRE-GS09000383DA	AEOS Residual Solvents Mix #1 4000 µg/mL in N,N-Dimethylacetamide(‡)	5x1ml
	isobutane butane (C4)	
AEOS residual solvents mix #2		
DRE-GS09000382DA	AEOS residual solvents Mix #2 4000 µg/mL in N,N-Dimethylacetamide(‡)	5x1ml
	benzene toluene o-xylene heptane (C7) 2,3-dimethylpentane 3,3-dimethylpentane 2-methylhexane 2,2,3-trimethylbutane	n-hexane (C6) m-xylene p-xylene 2,2-dimethylpentane 2,4-dimethylpentane 3-ethylpentane 3-methylhexane
Alkanes-Mix 10		
DRE-YA03010010TO	Alkanes-Mix 10 500 µg/mL in Toluene(‡)	1ml
	n-Decane n-Eicosane n-Heptadecane n-Nonadecane n-Pentadecane n-Tetatriacontane n-Tritriacontane	n-Docosane n-Heneicosane n-Hexacosane n-Octacosane n-Pentatriacontane n-Triacontane n-Undecane
	n-Dodecane n-Hentriacontane n-Hexadecane n-Octadecane n-Tetracosane n-Tricosane	n-Dotriacontane n-Heptacosane n-Nonacosane n-Pentacosane n-Tetradecane n-Tridecane
Alkanes-Mix 12		
DRE-XA03010012TO	Alkanes-Mix 12 100 µg/mL in Toluene	1ml
	n-Decane n-Eicosane n-Octacosane n-Tetracontane n-Triacontane	n-Docosane n-Hexacosane n-Octadecane n-Tetracosane
	n-Dodecane n-Hexadecane n-Octane n-Tetradecane	n-Dotriacontane n-Hexatriacontane n-Octatriacontane n-Tetatriacontane
Cannabis Residual Solvents Mix		
DRE-GA09000138TN	Cannabis Residual Solvents Mix 1000 µg/mL in Triacetin(‡)	1ml
DRE-GS09000138TN	Cannabis Residual Solvents Mix 1000 µg/mL in Triacetin(‡)	5x1ml
	butane (C4) 2-methylbutane 1-pentanol isopropyl alcohol 1,2-dimethoxyethane acetone acetonitrile N,N-dimethylacetamide 2-methylpentane heptane (C7) o-xylene	isobutane 2,2-dimethylbutane 1-propanol ethanol 1,4-dioxane 2-butanone (MEK) isopropylbenzene N,N-dimethylformamide 3-methylpentane benzene m-xylene
	n-propane 2,3-dimethylbutane 2-butanol ethylene glycol ethyl ether ethyl acetate methylene chloride pyridine n-hexane (C6) toluene p-xylene	n-pentane (C5) 1-butanol 2-ethoxyethanol methanol tetrahydrofuran (THF) isopropyl acetate dimethyl sulfoxide (DMSO) tetramethylene sulfone cyclohexane ethylbenzene
Cannabis Solvent Mix, 27-1625		
DRE-GS09000367TN	Cannabis Solvent Mix, 27-1625 1000 µg/mL in Triacetin(‡)	5x1ml
	methanol 2,2-dimethylbutane methylene chloride ethyl acetate isopropyl acetate toluene p-xylene	2-methylbutane acetone 2-methylpentane tetrahydrofuran (THF) heptane (C7) ethylene glycol m-xylene
	n-pentane (C5) acetonitrile 3-methylpentane cyclohexane 1,4-dioxane ethylbenzene cumene	ethyl ether 2,3-dimethylbutane n-hexane (C6) benzene 2-ethoxyethanol o-xylene

Hydrocarbons and Petrochemicals

Product code	Description	
Cannabis Solvent Mix, 4-1624		
DRE-GS09000366TN	Cannabis Solvent Mix, 4-1624 1000 µg/mL in Triacetin(‡)	5x1ml
	n-propane isobutane butane (C4) 2,2-dimethylpropane	
Cannabis Spiking Mix, 3-1437		
DRE-GA09000227DS	Cannabis Spiking Mix, 3-1437 100 µg/mL in Dimethyl Sulfoxide(‡)	1ml
	butane (C4) isobutane n-propane	
DRO Mixture 1		
DRE-YA03010001HE	DRO Mixture 1 1000 µg/mL in n-Hexane	1ml
	n-Decane n-Eicosane n-Heptadecane n-Octadecane n-Tetradecane n-Undecane	n-Docosane n-Heneicosane n-Nonadecane n-Pentacosane n-Tridecane
EPH CT Aliphatics Mix		
DRE-YS09000040DI	EPH CT Aliphatics Mix, 15-2 1000 µg/mL in Dichloromethane(‡)	5x1ml
	nonane dodecane hexadecane eicosane tetracosane octacosane dotriacontane hexatriacontane	decane tetradecane octadecane docosane hexacosane triacontane tetatriacontane
EPH Extraction Surrogate Mix		
DRE-YS09000041AC	EPH Extraction Surrogate Mix 2000 µg/mL in Acetone(‡)	5x1ml
	1-chlorooctadecane o-terphenyl	
EPH MA Aliphatics Mix		
DRE-YS09000043HE	EPH MA Aliphatics Mix 1000 µg/mL in Hexane(‡)	5x1ml
	decane (C10) dodecane (C12) hexacosane (C26) hexatriacontane (C36) octacosane (C28) tetradecane (C14) tetracosane (C24)	docosane (C22) eicosane (C20) n-hexadecane (C16) nonane (C9) triacontane (C30) octadecane (C18) nonadecane (C19)
EPH NJ Aliphatics Mix		
DRE-YS09000046CY	EPH NJ Aliphatics Mix 1000 µg/mL in Cyclohexane(‡)	5x1ml
	decane (C10) n-eicosane (C20) hexatriacontane (C36) octatriacontane (C38) tetatriacontane (C34)	n-docosane (C22) n-heneicosane (C21) nonane (C9) tetracontane (C40) triacontane (C30)
		n-dodecane (C12) hexacosane (C26) octacosane (C28) n-tetracosane (C24)
		dotriacontane (C32) n-hexadecane (C16) n-octadecane (C18) n-tetradecane (C14)
EPH NJ Rev. 2 Aliphatics Mix		
DRE-YS09000045HC	EPH NJ Rev. 2 Aliphatics Mix 2000 µg/mL in Hexane:Carbon Disulfide 80:20(‡)	5x1ml
	decane (C10) n-eicosane (C20) hexatriacontane (C36) octacosane (C28) n-tetracosane (C24)	n-docosane (C22) n-heneicosane (C21) 2-methylnaphthalene n-octadecane (C18) n-tetradecane (C14)
		n-dodecane (C12) hexacosane (C26) naphthalene octatriacontane (C38) tetatriacontane (C34)
		dotriacontane (C32) n-hexadecane (C16) nonane (C9) tetracontane (C40) triacontane (C30)

(‡) ISO Guide 34

(*) Shorter expiry due to chemical nature of component(s)

Order online at www.lgcstandards.com

Hydrocarbons and Petrochemicals

Product code	Description	
Florida TRPH Mix		
DRE-GA09000228CD	Florida TRPH Mix 2000 µg/mL in Carbon Disulfide(‡)	10ml
octane (C8) n-hexadecane (C16) tetracosane (C24) dotriacontane (C32) tetracontane (C40)	decane (C10) n-octadecane (C18) hexacosane (C26) tetratriacontane (C34)	dodecane (C12) eicosane (C20) octacosane (C28) hexatriacontane (C36)
		n-tetradecane (C14) docosane (C22) triacontane (C30) octatriacontane (C38)
Florida TRPH Standard		
DRE-YA09000020CD	Florida TRPH Mix 2000 µg/mL in Carbon Disulfide(‡)	1ml
DRE-YS09000020CD	Florida TRPH Mix 2000 µg/mL in Carbon Disulfide(‡)	5x1ml
octane (C8) n-hexadecane (C16) tetracosane (C24) dotriacontane (C32) tetracontane (C40)	decane (C10) octadecane (C18) hexacosane (C26) tetratriacontane (C34)	dodecane (C12) eicosane (C20) octacosane (C28) hexatriacontane (C36)
		tetradecane (C14) docosane (C22) triacontane (C30) octatriacontane (C38)
Hawaii Solvent Mix #1		
DRE-GS09000245AL	Hawaii Solvent Mix #1 10000 µg/mL in Acetonitrile(‡)	5x1ml
	n-hexane (C6) toluene p-xylene isobutane	benzene m-xylene o-xylene butane (C4)
Hawaii Solvent Mix #2		
DRE-GS09000246AL	Hawaii Solvent Mix #2 10000 µg/mL in Acetonitrile(‡)	5x1ml
DRE-GS09000246DS	Hawaii Solvent Mix #2 10000 µg/mL in Dimethyl Sulfoxide(‡)	5x1ml
	heptane (C7) 3-methylhexane 2,3-dimethylpentane 3,3-dimethylpentane 2,2,3-trimethylbutane	2-methylhexane 2,2-dimethylpentane 2,4-dimethylpentane 3-ethylpentane
Ketones Mix		
DRE-GS09000064DM	Ketones Mix 10000 µg/mL in Dimethyl Formamide(‡)	5x1ml
	2-butanone (MEK) acetone 2-hexanone	
Low Level Oxygenates Mix		
DRE-GH09000168HE	Low-level Oxygenate Mix 10 µg/mL in Hexane(‡)	10x2ml
propyl ether Acetone 2-butanone (MEK) Ethanol isopropyl ether 1-propanol	Isobutyraldehyde allyl alcohol tert-butyl ethyl ether (ETBE) ethyl ether Methanol propionaldehyde	Valeraldehyde tert-amyl methyl ether (TAME) Butyraldehyde isobutyl alcohol methyl t-butyl ether
		Acetaldehyde 1-butanol dimethyl ether isopropyl alcohol 2-methyl-2-propanol
Oregon Cannabis Solvent Mix, 35-2987		
DRE-GA09000238TN	Oregon Cannabis Solvent Mix, 35-2987 1000 µg/mL in Triacetin(‡)	1ml
DRE-GS09000238TN	Oregon Cannabis Solvent Mix, 35-2987 1000 µg/mL in Triacetin(‡)	5x1ml
DRE-GS09000238TN-C	Oregon Cannabis Solvent Mix, 35-2987 1000 µg/mL in Triacetin(‡)	5x1ml
DRE-GS09000238TN-SS	Oregon Cannabis Solvent Mix, SECOND SOURCE 35-2987 1000 µg/mL in Triacetin(‡)	5x1ml
butane (C4) 2-methylbutane 3-methylpentane benzene m-xylene isopropylbenzene tetrahydrofuran (THF) isopropyl alcohol	isobutane 2,3-dimethylbutane n-hexane (C6) toluene p-xylene methylene chloride ethyl ether acetone	ethylene oxide 2,2-dimethylbutane cyclohexane ethylbenzene 1,4-dioxane ethanol 2-butanol methanol
		n-propane 2-methylpentane heptane (C7) o-xylene acetonitrile ethyl acetate 2-ethoxyethanol isopropyl acetate

Hydrocarbons and Petrochemicals

Product code	Description	
n-pentane (C5)	ethylene glycol	2,2-dimethylpropane

Residual Solvent Gases Spiking Mix

DRE-GH09000206DS	Residual Solvent Gases Spiking Mix 100 µg/mL in Dimethyl Sulfoxide(‡)	10x1.5ml
	acetylene butane (C4) 2-Methylpropene n-pentane (C5)	

Residual Solvents - FET Mix

DRE-GA09000241DS	Residual Solvents - FET Mix 100 µg/mL in Dimethyl Sulfoxide(‡)	1ml
	acetone ethanol methanol butane (C4) n-hexane (C6)	acetonitrile isopropyl alcohol n-propane isobutane heptane (C7)

Residual Solvents Gases Spiking Mix

DRE-GS09000187DS	Residual Solvents Gases Spiking Mix 100 µg/mL in Dimethyl Sulfoxide(‡)	5x1ml
	butane (C4) isobutane n-propane	

Residual Solvents Mix

DRE-GS09000177DS	Residual Solvents Mix 50 µg/mL in Dimethyl Sulfoxide(‡)	5x1ml
DRE-GS09000178DS	Residual Solvents Mix 500 µg/mL in Dimethyl Sulfoxide(‡)	5x1ml
DRE-GS09000179DS	Residual Solvents Mix 2500 µg/mL in Dimethyl Sulfoxide(‡)	5x1ml
2-methylbutane cyclohexane ethylbenzene isopropyl alcohol 3-methylpentane toluene chloroform ethylene glycol	acetone ethanol heptane (C7) methanol n-pentane (C5) o-xylene 2,2-dimethylbutane	benzene ethyl ether n-hexane (C6) methylene chloride 1-pentanol m-xylene 2,3-dimethylbutane
		2-butanone (MEK) ethyl acetate isooctane 2-methylpentane 1-propanol p-xylene 1,1,1,2-Tetrafluoroethane

Texas TPH Mix

DRE-GA09000169PE	Texas TPH Mix 20000 µg/mL in Pentane(‡)	1ml
	gasoline, mixed grades composite diesel fuel #2	

USP Class 3 Mix B

DRE-GA09000170DM	USP Class 3 Mix B 1000 µg/mL in Dimethyl Formamide(‡)	1ml
DRE-GS09000170DM	USP Class 3 Mix B 1000 µg/mL in Dimethyl Formamide(‡)	5x1ml
ethanol acetone dimethyl sulfoxide (DMSO) isobutyl alcohol propyl acetate isobutyl acetate formic acid	n-pentane (C5) ethyl formate 2-butanol isopropyl acetate 3-methyl-1-butanol butyl acetate 2-pentanol	ethyl ether methyl acetate 2-butanone (MEK) 1-butanol 4-methyl-2-pentanone (MIBK) anisole 1-propanol
		isopropyl alcohol methyl t-butyl ether ethyl acetate heptane (C7) 1-pentanol acetic acid

VPH Calibration Mix

DRE-YS09000048ME	VPH Calibration Mix 50 µg/mL in p/t Methanol(‡)	5x1ml
benzene ethylbenzene naphthalene o-xylene	butylcyclohexane isooctane toluene p-xylene	decane (C10) methyl t-butyl ether 1,2,4-trimethylbenzene nonane
		2,5-dibromotoluene 2-methylpentane m-xylene n-pentane










Hydrocarbons and Petrochemicals

Product code	Description	
Volatile BTEX Mix		
DRE-GA09000162ME	Volatile BTEX Mix 2000 µg/mL in p/t Methanol(‡) benzene toluene ethylbenzene o-xylene m-xylene p-xylene	1ml










FLAME RETARDANTS - POLYBROMINATED BIPHENYLS





Flame retardants - polybrominated biphenyls

Product code	Description			
PBB 1 (2-Bromobiphenyl)				
CAS 2052-07-5 DRE-C21000100	MW 233.1039 PBB-No. 1	C ₁₂ H ₉ Br	50mg	
PBB 2 (3-Bromobiphenyl)				
CAS 2113-57-7 DRE-C21000200	MW 233.1039 PBB-No. 2	C ₁₂ H ₉ Br	50mg	
PBB 3 (4-Bromobiphenyl)				
CAS 92-66-0 DRE-C21000300 DRE-L21000300CY	MW 233.1039 PBB-No. 3 PBB-No. 3 10 µg/mL in Cyclohexane	C ₁₂ H ₉ Br	50mg 10ml	
PBB 4 (2,2'-Dibromobiphenyl)				
CAS 13029-09-9 DRE-C21000400	MW 311.9999 PBB No. 4	C ₁₂ H ₈ Br ₂	10mg	
PBB 9 (2,5-Dibromobiphenyl)				
CAS 57422-77-2 DRE-E21000900	MW 311.9999 PBB-No. 9	C ₁₂ H ₈ Br ₂	15mg	
PBB 15 (4,4'-Dibromobiphenyl)				
CAS 92-86-4 DRE-C21001500 DRE-L21001500IO	MW 311.9999 PBB-No. 15 PBB-No. 15 10 µg/mL in Isooctane	C ₁₂ H ₈ Br ₂	10mg 10ml	
PBB 18 (2,2',5-Tribromobiphenyl)				
CAS 59080-34-1 DRE-E21001800 DRE-L21001800IO	MW 390.896 PBB-No. 18 PBB-No. 18 10 µg/mL in Isooctane	C ₁₂ H ₇ Br ₃	10mg 10ml	
PBB 26 (2,3',5-Tribromobiphenyl)				
CAS 59080-35-2 DRE-E21002600	MW 390.896 PBB-No. 26	C ₁₂ H ₇ Br ₃	10mg	
PBB 29 (2,4,5-Tribromobiphenyl)				
CAS 115245-07-3 DRE-C21002900 DRE-L21002900AL	MW 390.896 PBB-No. 29 PBB-No. 29 10 µg/mL in Acetonitrile	C ₁₂ H ₇ Br ₃	10mg 10ml	

Flame retardants - polybrominated biphenyls

Product code	Description			
PBB 30 (2,4,6-Tribromobiphenyl)				
CAS 59080-33-0	MW 390.896	C ₁₂ H ₇ Br ₃		
DRE-E21003000	PBB-No. 30		15mg	
DRE-L21003000CY	PBB-No. 30 10 µg/mL in Cyclohexane		10ml	
PBB 31 (2,4',5-Tribromobiphenyl)				
CAS 59080-36-3	MW 390.896	C ₁₂ H ₇ Br ₃		
DRE-E21003100	PBB-No. 31		10mg	
PBB 49 (2,2',4,5'-Tetrabromobiphenyl)				
CAS 60044-24-8	MW 469.792	C ₁₂ H ₆ Br ₄		
DRE-C21004900	PBB-No. 49		5mg	
DRE-L21004900AL	PBB-No. 49 10 µg/mL in Acetonitrile		10ml	
PBB 52 (2,2',5,5'-Tetrabromobiphenyl)				
CAS 59080-37-4	MW 469.792	C ₁₂ H ₆ Br ₄		
DRE-C21005200	PBB-No. 52		10mg	
DRE-L21005200AL	PBB-No. 52 10 µg/mL in Acetonitrile		10ml	
PBB 53 (2,2',5,6'-Tetrabromobiphenyl)				
CAS 60044-25-9	MW 469.792	C ₁₂ H ₆ Br ₄		
DRE-E21005300	PBB-No. 53		5mg	
PBB 80 (3,3',5,5'-Tetrabromobiphenyl)				
CAS 16400-50-3	MW 469.792	C ₁₂ H ₆ Br ₄		
DRE-L21008000CY	PBB-No. 80 10 µg/mL in Cyclohexane		10ml	
PBB 101 (2,2',4,5,5'-Pentabromobiphenyl)				
CAS 67888-96-4	MW 548.6881	C ₁₂ H ₅ Br ₅		
DRE-LA21010100CY	PBB-No. 101 10 µg/mL in Cyclohexane		1ml	
PBB 103 (2,2',4,5',6-Pentabromobiphenyl)				
CAS 59080-39-6	MW 548.6881	C ₁₂ H ₅ Br ₅		
DRE-E21010300	PBB-No. 103		5mg	
PBB 153 (2,2',4,4',5,5'-Hexabromobiphenyl)				
CAS 59080-40-9	MW 627.5842	C ₁₂ H ₄ Br ₆		
DRE-LA21015300CY	PBB-No. 153 10 µg/mL in Cyclohexane		1ml	










Flame retardants - polybrominated biphenyls

Product code	Description			
PBB 155 (2,2',4,4',6,6'-Hexabromobiphenyl)				
CAS 59261-08-4	MW 627.5842	C ₁₂ H ₄ Br ₆		5mg 
DRE-E21015500	PBB-No. 155			
PBB 209 (Decabromobiphenyl)				
CAS 13654-09-6	MW 943.1684	C ₁₂ Br ₁₀		25mg 10ml 
DRE-E21020900	PBB-No. 209			
DRE-L21020900CY	PBB-No. 209 10 µg/mL in Cyclohexane			
PBB-Mix 5				
DRE-LA21030500CY	PBB-Mix 5 10 µg/mL in Cyclohexane(‡)			1ml
PBB 1 (2-Bromobiphenyl)	PBB 2 (3-Bromobiphenyl)	PBB 3 (4-Bromobiphenyl)	PBB 4 (2,2'-Dibromobiphenyl)	
PBB 7 (2,4-Dibromobiphenyl)	PBB 9 (2,5-Dibromobiphenyl)	PBB 10 (2,6-Dibromobiphenyl)	PBB 15 (4,4'-Dibromobiphenyl)	
PBB 18 (2,2',5-Tribromobiphenyl)	PBB 26 (2,3',5-Tribromobiphenyl)	PBB 29 (2,4,5-Tribromobiphenyl)	PBB 30 (2,4,6-Tribromobiphenyl)	
PBB 31 (2,4',5-Tribromobiphenyl)	PBB 38 (3,4,5-Tribromobiphenyl)	PBB 49 (2,2',4,5'-Tetrabromobiphenyl)	PBB 52 (2,2',5,5'-Tetrabromobiphenyl)	
PBB 53 (2,2',5,6'-Tetrabromobiphenyl)	PBB 80 (3,3',5,5'-Tetrabromobiphenyl)	PBB 101 (2,2',4,5,5'-Pentabromobiphenyl)	PBB 103 (2,2',4,5',6-Pentabromobiphenyl)	
PBB 153 (2,2',4,4',5,5'-Hexabromobiphenyl)	PBB 155 (2,2',4,4',6,6'-Hexabromobiphenyl)	PBB 209 (Decabromobiphenyl)		










PHARMACEUTICAL AND VETERINARY COMPOUNDS AND METABOLITES












Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Acepromazine Hydrochloride				
CAS 973-12-6 DRE-C10010290	MW 362.9167 Acepromazine hydrochloride	C ₁₉ H ₂₂ N ₂ OS·ClH	10mg	
Acepromazine Maleate				
CAS 3598-37-6 DRE-C10010300	MW 442.5279 Acepromazine maleate salt(‡)	C ₁₉ H ₂₂ N ₂ OS·C ₄ H ₄ O ₄	100mg	
Acepromazine-d6 Maleate				
CAS 1331655-50-5 DRE-C10010320	MW 448.5649 Acepromazine D6 maleate	C ₁₉ ² H ₆ H ₁₆ N ₂ OS·C ₄ H ₄ O ₄	10mg	
4-Acetylaminophenazone				
CAS 83-15-8 DRE-C10011900	MW 245.2771 4-Acetylaminophenazone	C ₁₃ H ₁₅ N ₃ O ₂	100mg	
Acetylsalicylic Acid				
CAS 50-78-2 DRE-C10024000	MW 180.1574 Acetylsalicylic acid	C ₉ H ₈ O ₄	250mg	
Acetylsulfamethoxazole (N-[4-[(5-Methylisoxazol-3-yl)sulphamoyl]phenyl]acetamide)				
CAS 21312-10-7 DRE-C10024050	MW 295.3143 Acetylsulfamethoxazole(‡)	C ₁₂ H ₁₃ N ₃ O ₄ S	100mg	
Acetylsulfamethoxazole D4				
CAS 1215530-54-3 DRE-C10024051	MW 299.339 Acetylsulfamethoxazole D4	C ₁₂ ² H ₄ H ₉ N ₃ O ₄ S	10mg	
Aciclovir (Acyclovir)				
CAS 59277-89-3 DRE-C10045600	MW 225.2046 Acyclovir	C ₈ H ₁₁ N ₅ O ₃	100mg	
Acriflavine				
CAS 8048-52-0 DRE-C10043000	MW 468.9806 Acriflavine	C ₁₄ H ₁₄ N ₃ ·C ₁₃ H ₁₁ N ₃ ·Cl	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Adonit				
CAS 488-81-3 DRE-C10046500	MW 152.1458 Adonit	C ₅ H ₁₂ O ₅	100mg	
Adrenaline (L-Adrenaline; (-)-Epinephrine)				
CAS 51-43-4 DRE-C13177900	MW 183.2044 (-)-Epinephrine	C ₉ H ₁₃ NO ₃	0.1g	
DL-Adrenaline Hydrochloride ((±)-Epinephrine hydrochloride)				
CAS 329-63-5 DRE-C13178000	MW 219.6654 (±)-Epinephrine hydrochloride	C ₉ H ₁₃ NO ₃ ·ClH	100mg	
Agaric Acid				
CAS 666-99-9 DRE-C10047700	MW 416.5488 Agaric acid	C ₂₂ H ₄₀ O ₇	250mg	
Aklomide				
CAS 3011-89-0 DRE-C10049000	MW 200.5792 Aklomide	C ₇ H ₅ ClN ₂ O ₃	100mg	
Albendazole				
CAS 54965-21-8 DRE-C10065000	MW 265.3314 Albendazole	C ₁₂ H ₁₅ N ₃ O ₂ S	100mg	
Albendazole Sulfone (Methyl [5-Propylsulphonyl]-1H-benzimidazol-2-yl]carbamate)				
CAS 75184-71-3 DRE-C10065300	MW 297.3302 Albendazole sulfone(‡)	C ₁₂ H ₁₅ N ₃ O ₄ S	10mg	
Albendazole Sulfoxide (Methyl [5-Propylsulphiny]-1H-benzimidazol-2-yl]carbamate)				
CAS 54029-12-8 DRE-C10065400	MW 281.3308 Albendazole sulfoxide(‡)	C ₁₂ H ₁₅ N ₃ O ₃ S	10mg	
Albendazole-2-amino (5-(Propylsulfany)-1H-benzimidazol-2-amine)				
CAS 80983-36-4 DRE-C10065020	MW 207.2953 Albendazole-2-amino	C ₁₀ H ₁₃ N ₃ S	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Albendazole-2-aminosulfone (5-(Propylsulfonyl) -1H-benzimidazol-2-amine)				
CAS 80983-34-2 DRE-C10065200	MW 239.2941 Albendazole-2-aminosulfone	C ₁₀ H ₁₃ N ₃ O ₂ S	50mg	
Albendazole-2-aminosulfoxide				
CAS 80983-35-3 DRE-C10065250	MW 223.2947 Albendazole-2-aminosulfoxide	C ₁₀ H ₁₃ N ₃ OS	50mg	
Allopurinol				
CAS 315-30-0 DRE-C10118000	MW 136.1115 Allopurinol	C ₅ H ₄ N ₄ O	100mg	
Altrenogest				
CAS 850-52-2 DRE-C10144000	MW 310.4299 Altrenogest	C ₂₁ H ₂₆ O ₂	100mg	
Amantadine Hydrochloride				
CAS 665-66-7 DRE-C10146000	MW 187.7096 Amantadine Hydrochloride	C ₁₀ H ₁₇ N·ClH	100mg	
Ambroxol Hydrochloride				
CAS 23828-92-4 DRE-C10148800	MW 414.5638 Ambroxol hydrochloride	C ₁₃ H ₁₆ Br ₂ N ₂ O·ClH	100mg	
Amidate (O,O-Dimethylphosphoramidothioate)				
CAS 17321-47-0 DRE-CA10158000	MW 141.1292 Amidate	C ₂ H ₆ NO ₂ PS	100mg	
Amikacin hydrate				
CAS 1257517-67-1 DRE-C10164000	MW 603.6178 Amikacin hydrate	C ₂₂ H ₄₃ N ₅ O ₁₃ ·H ₂ O	100mg	
Aminotriazole (2-Acetamido-5-nitrothiazole)				
CAS 140-40-9 DRE-C10012500	MW 187.1765 2-Acetamido-5-nitrothiazole	C ₅ H ₅ N ₃ O ₃ S	100mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
2-Aminoflubendazole ((2-Amino-1H-benzimidazol-5-yl)(4-fluorophenyl)methanone)				
CAS 82050-13-3	MW 255.2471	C ₁₄ H ₁₀ FN ₃ O		
DRE-C10202370	2-Aminoflubendazole		10mg	
DRE-L10202370ME	2-Aminoflubendazole 10 µg/mL in Methanol		10ml	
1-Aminohydantoin (2,4,5-13C3)				
CAS 957509-31-8	MW 118.0687	¹³ C ₃ H ₅ N ₃ O ₂		
DRE-XA10203190AL	1-Aminohydantoin (2,4,5-13C3) 100 µg/mL in Acetonitrile		1.1ml	
1-Aminohydantoin Hydrochloride				
CAS 2827-56-7	MW 151.5516	C ₃ H ₅ N ₃ O ₂ ·ClH		
DRE-C10203200	1-Aminohydantoin hydrochloride(±)		100mg	
3-Amino-5-morpholinomethyl-1,3-oxazolidin-2-one (AMOZ)				
CAS 43056-63-9	MW 201.223	C ₈ H ₁₅ N ₃ O ₃		
DRE-C10206300	3-Amino-5-morpholinomethyl-2-oxazolidinone (AMOZ)(±)		50mg	
DRE-L10206300AL	3-Amino-5-morpholinomethyl-2-oxazolidinone 10 µg/mL in Acetonitrile		10ml	
3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5)				
CAS 1017793-94-0	MW 206.2538	C ₈ ² H ₁₅ H ₁₀ N ₃ O ₃		
DRE-C10206310	3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5)(±)		10mg	
DRE-XA10206310AL	3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5) 100 µg/mL in Acetonitrile		1ml	
3-Amino-2-oxazolidinone (AOZ)				
CAS 80-65-9	MW 102.0919	C ₃ H ₆ N ₂ O ₂		
DRE-C10209000	3-Amino-2-oxazolidinone AOZ(±)		50mg	
DRE-L10209000AL	3-Amino-2-oxazolidinone AOZ 10 µg/mL in Acetonitrile (*)		10ml	
3-Amino-2-oxazolidinone Sulfate				
CAS 32957-26-9	MW 200.1704	C ₃ H ₆ N ₂ O ₂ ·H ₂ O ₄ S		
DRE-C10209200	3-Amino-2-oxazolidinone sulfate		10mg	
3-Amino-2-oxazolidinone D4 (AOZ D4)				
CAS 1188331-23-8	MW 106.1166	C ₃ ² H ₄ H ₂ N ₂ O ₂		
DRE-C10209010	3-Amino-2-oxazolidinone (AOZ) D4(±)		10mg	
Aminophenazone (4,4-Dimethylaminophenazone)				
CAS 58-15-1	MW 231.2936	C ₁₃ H ₁₇ N ₃ O		
DRE-C10209600	Aminophenazone		250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Androstanolone (5α-Androstan-17β-ol-3-one)				
CAS 521-18-6 DRE-C10255010	MW 290.4403 5α-Androstan-17β-ol-3-one	C ₁₉ H ₃₀ O ₂	100mg	
Androstenedione (4-Androstene-3,17-dione)				
CAS 63-05-8 DRE-C10255030	MW 286.4085 4-Androstene-3,17-dione	C ₁₉ H ₂₈ O ₂	100mg	
Androsterone				
CAS 53-41-8 DRE-C10255040	MW 290.4403 Androsterone	C ₁₉ H ₃₀ O ₂	100mg	
8'-Apoaldehyde (Apocarotenal)				
CAS 1107-26-2 DRE-CA10290900	MW 416.638 8'-Apoaldehyde (Apocarotenal)	C ₃₀ H ₄₆ O	50mg	
Apramycin Sulfate				
CAS 65710-07-8 DRE-C10293000	MW 637.6556 Apramycin sulfate	C ₂₁ H ₄₁ N ₅ O ₁₁ ·H ₂ O ₄ S	250mg	
4-Arsanilic Acid				
CAS 98-50-0 DRE-C10300250	MW 217.0542 4-Arsanilic acid	C ₆ H ₆ AsNO ₃	100mg	
β-Asarone				
CAS 5273-86-9 DRE-C10301000	MW 208.2536 β-Asarone	C ₁₂ H ₁₆ O ₃	25mg	
Atenolol				
CAS 29122-68-7 DRE-C10313000	MW 266.3361 Atenolol(±)	C ₁₄ H ₂₂ N ₂ O ₃	100mg	
Atorvastatin Calcium				
CAS 134523-03-8 DRE-C10318000	MW 1155.3417 Atorvastatin calcium	2C ₃₃ H ₃₄ FN ₂ O ₅ ·Ca	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Atropine				
CAS 51-55-8 DRE-C10333500	MW 289.3694 Atropine	C ₁₇ H ₂₃ NO ₃	250mg	
Azaperol				
CAS 2804-05-9 DRE-C10340500 DRE-XA10340500ME	MW 329.4118 Azaperol Azaperol 100 µg/mL in Methanol	C ₁₉ H ₂₄ FN ₃ O	10mg 1ml	
Azaperone				
CAS 1649-18-9 DRE-C10340510	MW 327.3959 Azaperone(±)	C ₁₉ H ₂₂ FN ₃ O	50mg	
Azaperone D4				
CAS 1173021-72-1 DRE-C10340512	MW 331.4205 Azaperone D4	C ₁₉ ² H ₄ H ₁₈ FN ₃ O	10mg	
Azithromycin dihydrate				
CAS 117772-70-0 DRE-C10386000	MW 785.015 Azithromycin dihydrate	C ₃₈ H ₇₂ N ₂ O ₁₂ ·2H ₂ O	50mg	
Bacitracin				
CAS 1405-87-4 DRE-C10418000	MW n/a Bacitracin		250mg	No Structure
Bambuterol hydrochloride				
CAS 81732-46-9 DRE-C10418500	MW 403.9009 Bambuterol hydrochloride	C ₁₈ H ₂₉ N ₃ O ₅ ·ClH	10mg	
Beclometasone (9-Chloro-11β,17,21-trihydroxy-16β-methylpregna-1,4-diene-3,20-dione)				
CAS 4419-39-0 DRE-C10429000	MW 408.9157 Beclometasone (*)	C ₂₂ H ₂₉ ClO ₅	100mg	
Benzylpenicillin Potassium (Penicilline G potassium salt)				
CAS 113-98-4 DRE-C15935000	MW 372.4805 Penicilline G potassium salt(±)	C ₁₆ H ₁₇ N ₂ O ₄ ·S·K	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Berberine Hemisulfate				
CAS 316-41-6 DRE-C10574000	MW 768.785 Berberine sulfate	$2C_{20}H_{18}NO_4 \cdot O_4S$	100mg	
Betamethasone				
CAS 378-44-9 DRE-C10575000	MW 392.4611 Betamethasone(‡)	$C_{22}H_{28}FO_5$	100mg	
Betamethasone Acetate				
CAS 987-24-6 DRE-C10575005	MW 434.4977 Betamethasone Acetate	$C_{24}H_{31}FO_6$	50mg	
Betamethasone 17,21-Dipropionate				
CAS 5593-20-4 DRE-C10575010	MW 504.5876 Betamethasone-17,21-dipropionate	$C_{28}H_{37}FO_7$	100mg	
Betamethasone 17-Valerate				
CAS 2152-44-5 DRE-C10575020	MW 476.5775 Betamethasone 17a-valerate	$C_{27}H_{37}FO_6$	100mg	
Bezafibrate				
CAS 41859-67-0 DRE-C10578000	MW 361.8194 Bezafibrate	$C_{19}H_{20}ClNO_4$	100mg	
Bezafibrate D6 (dimethyl D6)				
CAS 1219802-74-0 DRE-C10578010	MW 367.8564 Bezafibrate D6 (dimethyl D6)	$C_{19}^2H_6H_{14}ClNO_4$	10mg	
Bisacodyl				
CAS 603-50-9 DRE-C10645000	MW 361.3906 Bisacodyl	$C_{22}H_{19}NO_4$	100mg	
N,N'-Bis-(4-nitrophenyl)urea				
CAS 587-90-6 DRE-C15598590	MW 302.2423 N,N'-Bis-(4-nitrophenyl)urea	$C_{13}H_{10}N_4O_5$	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Bisotrizole				
CAS 103597-45-1 DRE-C10653900	MW 658.8747 Bisotrizole	C ₄₁ H ₅₀ N ₈ O ₂	100mg	
Bisoprolol Fumarate				
CAS 104344-23-2 DRE-C10654000	MW 766.9582 Bisoprolol fumarate	2C ₁₈ H ₃₁ NO ₄ ·C ₄ H ₄ O ₄	100mg	
Boldenone				
CAS 846-48-0 DRE-C10662000	MW 286.4085 Boldenone	C ₁₉ H ₂₆ O ₂	10mg	
Boldenone 17-Acetate				
CAS 2363-59-9 DRE-XA10662010AL	MW 328.4452 Boldenon 17-acetate 100 µg/mL in Acetonitrile	C ₂₁ H ₂₈ O ₃	1ml	
Brombuterol Hydrochloride				
CAS 21912-49-2 DRE-C10683000	MW 402.5531 Brombuterol hydrochloride	C ₁₂ H ₁₈ BrN ₂ O·ClH	10mg	
Bromchlorbuterol Hydrochloride (1-(4-Amino-3-bromo-5-chlorophenyl)-2-[(1,1-dimethylethyl)amino]-ethanol hydrochloride)				
CAS 78982-84-0 DRE-C10683500	MW 358.1021 Bromchlorbuterol hydrochloride	C ₁₂ H ₁₈ BrClN ₂ O·ClH	10mg	
Brompheniramine Maleate				
CAS 980-71-2 DRE-C10781600	MW 435.3116 Brompheniramine maleate	C ₁₆ H ₁₉ BrN ₂ ·C ₄ H ₄ O ₄	100mg	
Budesonide				
CAS 51333-22-3 DRE-C10825000	MW 430.5339 Budesonide	C ₂₅ H ₃₄ O ₆	100mg	
Bufexamac				
CAS 2438-72-4 DRE-C10830500	MW 223.2683 Bufexamac	C ₁₂ H ₁₇ NO ₃	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Cadaverine dihydrochloride				
CAS 1476-39-7 DRE-C10933500	MW 175.0999 Cadaverine dihydrochloride	$C_5H_{14}N_2 \cdot 2ClH$	100mg	
Caffeic Acid (3,4-Dihydroxycinnamic Acid)				
CAS 331-39-5 DRE-C10934700	MW 180.1574 Caffeic acid	$C_9H_8O_4$	100mg	
Caffeine				
CAS 58-08-2 DRE-C11693000 DRE-CR11693000 DRE-L11693000ME	MW 194.1906 Caffeine(†) Caffeine Caffeine 10 µg/mL in Methanol	$C_8H_{10}N_4O_2$	250mg 250mg 10ml	
Cambendazol				
CAS 26097-80-3 DRE-C10937000	MW 302.3516 Cambendazol	$C_{14}H_{14}N_4O_2S$	10mg	
Carazolol				
CAS 57775-29-8 DRE-C10968000	MW 298.3795 Carazolol	$C_{18}H_{22}N_2O_2$	10mg	
Carazolol D7				
CAS 1173021-02-7 DRE-C10968010	MW 305.4226 Carazolol D7	$C_{18}^2H_{17}H_{15}N_2O_2$	10mg	
Carbadox				
CAS 6804-07-5 DRE-C10968300	MW 262.2215 Carbadox	$C_{11}H_{10}N_4O_4$	100mg	
Carbamazepine				
CAS 298-46-4 DRE-C10968500	MW 236.2686 Carbamazepin(†)	$C_{15}H_{12}N_2O$	100mg	
Carbazochrome (Adenochrome semicarbazone)				
CAS 69-81-8 DRE-C10045850	MW 236.2273 Adenochrome semicarbazone	$C_{10}H_{12}N_4O_3$	100mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Carbenicillin Sodium				
CAS 4800-94-6 DRE-C10993000	MW 422.3633 Carbenicillin disodium	C ₁₇ H ₁₆ N ₂ O ₆ S·2Na	100mg	
Carprofen				
CAS 53716-49-7 DRE-C11045850	MW 273.7143 Carprofen	C ₁₅ H ₁₂ ClNO ₂	100mg	
Cefadroxil				
CAS 50370-12-2 DRE-C11063000	MW 363.3883 Cefadroxil	C ₁₆ H ₁₇ N ₃ O ₅ S	100mg	
Cefalexin Monohydrate				
CAS 23325-78-2 DRE-C11064000	MW 365.4042 Cefalexin monohydrate(±)	C ₁₆ H ₁₇ N ₃ O ₄ S·H ₂ O	250mg	
Cefazolin Sodium				
CAS 27164-46-1 DRE-C11064100	MW 476.489 Cefazolin sodium salt(±)	C ₁₄ H ₁₃ N ₆ O ₄ S ₃ ·Na	100mg	
Cefoperazone Sodium				
CAS 62893-20-3 DRE-C11064300	MW 667.6492 Cefoperazone sodium salt	C ₂₅ H ₂₆ N ₉ O ₈ S ₂ ·Na	100mg	
Cefotaxime Sodium				
CAS 64485-93-4 DRE-C11064400	MW 477.4473 Cefotaxime sodium salt	C ₁₆ H ₁₆ N ₅ O ₇ S ₂ ·Na	100mg	
Cefpirome Sulfate				
CAS 98753-19-6 DRE-C11064500	MW 612.6558 Cefpirome sulfate	C ₂₂ H ₂₂ N ₆ O ₅ S ₂ ·H ₂ O ₄ S	100mg	
Cefquinome Sulfate				
CAS 118443-89-3 DRE-C11064700	MW 626.6823 Cefquinome sulfate(±)	C ₂₃ H ₂₄ N ₆ O ₅ S ₂ ·H ₂ O ₄ S	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Cefradine				
CAS 38821-53-3 DRE-C11064800	MW 349.4048 Cefradine	C ₁₆ H ₁₉ N ₃ O ₄ S	100mg	
Ceftiofur				
CAS 80370-57-6 DRE-C11065000	MW 523.5626 Ceftiofur(‡)	C ₁₉ H ₁₇ N ₅ O ₇ S ₃	100mg	
Ceftiofur Hydrochloride				
CAS 103980-44-5 DRE-C11065020	MW 560.0235 Ceftiofur hydrochloride	C ₁₉ H ₁₇ N ₅ O ₇ S ₃ ·ClH	100mg	
Cefuroxime Sodium				
CAS 56238-63-2 DRE-C11065300	MW 446.3671 Cefuroxime sodium salt	C ₁₆ H ₁₅ N ₄ O ₈ S·Na	100mg	
Cetrimide				
CAS 1119-97-7 DRE-C11070000	MW 336.3943 Cetrimide	C ₁₇ H ₃₈ N·Br	250mg	
Cetylpyridinium Chloride Monohydrate				
CAS 6004-24-6 DRE-C11075000	MW 358.0014 Cetylpyridinium chloride monohydrate	C ₂₁ H ₃₈ N·Cl·H ₂ O	250mg	
Chloral Hydrate				
CAS 302-17-0 DRE-C11098000 DRE-YA11098000AL	MW 165.403 Chloralhydrate Chloralhydrate 1000 µg/mL in Acetonitrile	C ₂ H ₃ Cl ₃ O ₂	250mg 1ml	
Chloramphenicol				
CAS 56-75-7 DRE-C11120000 DRE-L11120000AL DRE-XA11120000EA	MW 323.1294 Chloramphenicol(‡) Chloramphenicol 10 µg/mL in Acetonitrile (*) Chloramphenicol 100 µg/mL in Ethylacetate	C ₁₁ H ₁₂ Cl ₂ N ₂ O ₅	250mg 10ml 1ml	
Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2R)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)				
CAS 202480-68-0 DRE-C11120100 DRE-XA11120100AL	MW 328.1602 Chloramphenicol D5 (ring D4, benzyl D) Chloramphenicol D5 100 µg/mL in Acetonitrile(‡)	C ₁₁ ² H ₉ ² H ₇ Cl ₂ N ₂ O ₅	10mg 1ml	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description				
erythro-Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2S)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)					
CAS 1426174-26-6 DRE-XA11120110AL	MW 328.1602 (Erythro)-Chloramphenicol D5 (ring D4, benzyl D) 100 µg/mL in Acetonitrile	C ₁₁ ² H ₉ H7Cl ₂ N ₂ O ₅	1.1ml		
Chloramphenicol Palmitate					
CAS 530-43-8 DRE-C11121000	MW 561.5382 Chloramphenicol palmitat	C ₂₇ H ₄₂ Cl ₂ N ₂ O ₆	250mg		
Chlorhexidine					
CAS 55-56-1 DRE-C11310000	MW 505.4466 Chlorhexidine (free base)	C ₂₂ H ₃₀ Cl ₂ N ₁₀	250mg		
Chlorhexidine Diacetate					
CAS 56-95-1 DRE-C11310100	MW 625.5505 Chlorhexidine diacetate salt	C ₂₂ H ₃₀ Cl ₂ N ₁₀ ·2C ₂ H ₄ O ₂	250mg		
Chlorhexidine Dihydrochloride					
CAS 3697-42-5 DRE-C11310200	MW 578.3685 Chlorhexidine dihydrochloride	C ₂₂ H ₃₀ Cl ₂ N ₁₀ ·2ClH	250mg		
Chlormadinone 17-Acetate					
CAS 302-22-7 DRE-C11418500	MW 404.927 Chloromadinon 17-acetate	C ₂₃ H ₂₉ ClO ₄	100mg		
2-Chloronicotinic Acid					
CAS 2942-59-8 DRE-C11451900	MW 157.5545 2-Chloronicotinic acid	C ₆ H ₄ ClNO ₂	250mg		
Chloroquine Phosphate					
CAS 50-63-5 DRE-CA11506000	MW 515.8625 Chloroquine diphosphate	C ₁₈ H ₂₆ ClN ₃ ·2H ₃ O ₄ P	100mg		
Chlorothiazide					
CAS 58-94-6 DRE-C11510700	MW 295.7232 Chlorothiazide	C ₇ H ₆ ClN ₃ O ₄ S ₂	250mg		

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Chlorphenamine Maleate (Chlorpheniramine Maleate)				
CAS 113-92-8 DRE-C11555000	MW 390.8606 Chlorpheniramine maleate(±)	$C_{16}H_{19}ClN_2 \cdot C_4H_4O_4$	100mg	
Chlorphenesin				
CAS 104-29-0 DRE-C11553000	MW 202.6348 Chlorphenesin	$C_9H_{11}ClO_3$	250mg	
Chlorpromazine Hydrochloride				
CAS 69-09-0 DRE-C11575000	MW 355.3251 Chlorpromazine hydrochloride(±)	$C_{17}H_{19}ClN_2S \cdot ClH$	250mg	
Chlorpropamide				
CAS 94-20-2 DRE-C11577000	MW 276.7398 Chlorpropamide	$C_{10}H_{13}ClN_2O_3S$	250mg	
Chlorquinaldol				
CAS 72-80-0 DRE-C11505800	MW 228.0747 Chloroquinaldol	$C_{10}H_7Cl_2NO$	100mg	
Chlortetracycline Hydrochloride				
CAS 64-72-2 DRE-C11509100	MW 515.3406 Chlorotetracycline hydrochloride	$C_{22}H_{23}ClN_2O_8 \cdot ClH$	250mg	
Ciclesonide				
CAS 126544-47-6 DRE-C11666300	MW 540.6876 Ciclesonide	$C_{32}H_{44}O_7$	10mg	
Ciclosporin (Ciclosporin A; Cyclosporin A)				
CAS 59865-13-3 DRE-C11836300	MW 1202.6112 Cyclosporin A	$C_{62}H_{111}N_9O_{12}$	100mg	
Cimaterol				
CAS 54239-37-1 DRE-C11666350	MW 219.2829 Cimaterol	$C_{12}H_{17}N_3O$	10mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Cimaterol D7 (isopropyl D7)				
CAS 1228182-44-2 DRE-C11666352	MW 226.326 Cimaterol D7	C ₁₂ H ₁₇ H ₁₀ N ₃ O	10mg	
Cimbuterol				
CAS 54239-39-3 DRE-C11666400	MW 233.3095 Cimbuterol	C ₁₃ H ₁₉ N ₃ O	10mg	
Cinoxacin				
CAS 28657-80-9 DRE-C11668100	MW 262.2182 Cinoxacin(±)	C ₁₂ H ₁₀ N ₂ O ₅	100mg	
Ciprofloxacin D8 Hydrochloride				
CAS 1216659-54-9 DRE-C11668501 DRE-XA11668501WA	MW 375.8518 Ciprofloxacin D8 hydrochloride Ciprofloxacin D8 hydrochloride 100 µg/mL in Water	C ₁₇ ² H ₁₈ FN ₃ O ₃ ·ClH	10mg 1.1ml	
Ciprofloxacin Hydrochloride				
CAS 93107-08-5 DRE-C11668500	MW 367.8025 Ciprofloxacin hydrochloride(±)	C ₁₇ H ₁₈ FN ₃ O ₃ ·ClH	100mg	
Ciprofloxacin Lactate				
CAS 97867-33-9 DRE-C11668505	MW 421.4195 Ciprofloxacin lactate	C ₁₇ H ₁₈ FN ₃ O ₃ ·C ₃ H ₆ O ₃	100mg	
Citicoline Sodium				
CAS 33818-15-4 DRE-C11668507	MW 510.3058 Citicoline sodium	C ₁₄ H ₂₅ NaO ₁₁ P ₂ ·Na	100mg	
Clarithromycin				
CAS 81103-11-9 DRE-C11668540	MW 747.9534 Clarithromycin	C ₃₈ H ₆₉ NO ₁₃	100mg	
Clenbuterol Hydrochloride				
CAS 21898-19-1 DRE-C11668550	MW 313.6511 Clenbuterol hydrochloride(±)	C ₁₂ H ₁₈ Cl ₂ N ₂ O·ClH	100mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Clenbuterol-hydroxymethyl				
CAS 38339-18-3 DRE-C11668565	MW 293.1895 Clenbuterol-hydroxymethyl	C ₁₂ H ₁₈ Cl ₂ N ₂ O ₂	10mg	
(±)-Clenbuterol				
CAS 37148-27-9 DRE-C11668560	MW 277.1901 (±)-Clenbuterol	C ₁₂ H ₁₈ Cl ₂ N ₂ O	25mg	
(±)-Clenbuterol D9 (trimethyl D9)				
CAS 129138-58-5 DRE-C11668561 DRE-XA11668561AL	MW 286.2456 (±)-Clenbuterol D9 (trimethyl D9)(‡) (±)-Clenbuterol D9 (trimethyl D9) 100 µg/mL in Acetonitrile (*)	C ₁₂ ² H ₁₈ H ₉ Cl ₂ N ₂ O	25mg 1ml	
Clenpenterol Hydrochloride				
CAS 37158-47-7 DRE-C11668700	MW 327.6776 Clenpenterol hydrochloride	C ₁₃ H ₂₀ Cl ₂ N ₂ O·ClH	10mg	
Clenproperol				
CAS 38339-11-6 DRE-C11668740	MW 263.1635 Clenproperol(‡)	C ₁₁ H ₁₆ Cl ₂ N ₂ O	10mg	
Clenproperole D7				
CAS 1173021-09-4 DRE-C11668742	MW 270.2067 Clenproperole D7	C ₁₁ ² H ₁₇ H ₉ Cl ₂ N ₂ O	10mg	
Climbazole				
CAS 38083-17-9 DRE-C11670000 DRE-L11670000CY	MW 292.7607 Climbazole(‡) Climbazol 10 µg/mL in Cyclohexane	C ₁₅ H ₁₇ ClN ₂ O ₂	100mg 10ml	
Clindamycin Hydrochloride				
CAS 21462-39-5 DRE-C11670100	MW 461.444 Clindamycin hydrochloride	C ₁₈ H ₃₂ ClN ₂ O ₅ S·ClH	50mg	
Clobetasol Propionate				
CAS 25122-46-7 DRE-C11678000	MW 466.97 Clobetasol propionate	C ₂₅ H ₃₂ ClFO ₅	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Clofibrate (Clofibric acid ethyl ester)				
CAS 637-07-0	MW 242.6987	C ₁₂ H ₁₅ ClO ₃		
DRE-C11682000	Clofibrate (Clofibric acid ethyl ester)		100mg	
Clofibric Acid (2-(4-Chlorophenoxy)-2-methylpropionic acid)				
CAS 882-09-7	MW 214.6455	C ₁₀ H ₁₁ ClO ₃		
DRE-C11484000	2-(4-Chlorophenoxy)-2-methylpropionic acid (Clofibric acid)		250mg	
DRE-L11484000AL	2-(4-Chlorophenoxy)-2-methylpropionic acid (Clofibric acid) 10 µg/mL in Acetonitrile		10ml	
Clomifene Citrate (Clomiphene Citrate)				
CAS 50-41-9	MW 598.0831	C ₂₆ H ₂₈ ClNO·C ₆ H ₈ O ₇		
DRE-C11686500	Clomiphene Citrate salt		250mg	
Clopidol				
CAS 2971-90-6	MW 192.0426	C ₇ H ₇ Cl ₂ NO		
DRE-C11687000	Clopidol(±)		100mg	
Clorprenaline Hydrochloride				
CAS 6933-90-0	MW 250.1648	C ₁₁ H ₁₆ ClNO·ClH		
DRE-C11687510	Clorprenaline Hydrochloride		100mg	
Clorsulon				
CAS 60200-06-8	MW 380.6558	C ₈ H ₉ Cl ₃ N ₃ O ₄ S ₂		
DRE-C11691400	Clorsulon		100mg	
Closantel				
CAS 57808-65-8	MW 663.0737	C ₂₂ H ₁₄ Cl ₂ I ₂ N ₂ O ₂		
DRE-C11691500	Closantel		100mg	
Clostebol				
CAS 1093-58-9	MW 322.8695	C ₁₉ H ₂₇ ClO ₂		
DRE-C11691600	Clostebol		100mg	
Clostebol Acetate				
CAS 855-19-6	MW 364.9062	C ₂₁ H ₂₉ ClO ₃		
DRE-C11691620	Clostebol acetate		100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Clotrimazole				
CAS 23593-75-1 DRE-C11692000	MW 344.8368 Clotrimazole(‡)	C ₂₂ H ₁₇ ClN ₂	100mg	
Cloxacillin Sodium Monohydrate				
CAS 7081-44-9 DRE-C11692100	MW 475.8784 Cloxacillin sodium salt monohydrate(‡)	C ₁₉ H ₁₇ ClN ₃ O ₅ S·Na·H ₂ O	250mg	
Colistin sulfate				
CAS 1264-72-8 DRE-C11693500	MW n/a Colistin sulfate		100mg	No Structure
Corticosterone				
CAS 50-22-6 DRE-C11705000	MW 346.4605 Corticosterone	C ₂₁ H ₃₀ O ₄	100mg	
Cortisone				
CAS 53-06-5 DRE-C11705400	MW 360.444 Cortisone	C ₂₁ H ₂₈ O ₅	500mg	
Cortisone Acetate				
CAS 50-04-4 DRE-C11705500	MW 402.4807 Cortison acetate	C ₂₃ H ₃₀ O ₆	100mg	
Coumarin				
CAS 91-64-5 DRE-C11735000	MW 146.1427 Coumarin	C ₉ H ₆ O ₂	250mg	
Coumarin 5,6,7,8-D4				
CAS 185056-83-1 DRE-XA11735010AC	MW 150.1674 Coumarin 5,6,7,8-D4 100 µg/mL in Acetone	C ₉ ² H ₄ H ₂ O ₂	1.1ml	
Curcumin				
CAS 458-37-7 DRE-C11780000	MW 368.3799 Curcumin	C ₂₁ H ₂₀ O ₆	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Cyclophosphamide Monohydrate				
CAS 6055-19-2 DRE-C11833600	MW 279.1012 Cyclophosphamide monohydrate	$C_7H_{15}Cl_2N_2O_2P \cdot H_2O$	100mg	
Cymiazole				
CAS 61676-87-7 DRE-C11875000 DRE-L11875000AL	MW 218.318 Cymiazole Cymiazole 10 µg/mL in Acetonitrile	$C_{12}H_{14}N_2S$	100mg 10ml	
Danofloxacin Mesilate				
CAS 119478-55-6 DRE-C11960500	MW 453.4845 Danofloxacin mesylate(‡)	$C_{19}H_{20}FN_3O_3 \cdot CH_4O_3S$	100mg	
Dantron (Danthron)				
CAS 117-10-2 DRE-C11961000	MW 240.2109 Danthron	$C_{14}H_8O_4$	100mg	
Dapsone				
CAS 80-08-0 DRE-C11963000	MW 248.3009 Dapson	$C_{12}H_{12}N_2O_2S$	250mg	
Decoquinat				
CAS 18507-89-6 DRE-C12097000	MW 417.5384 Decoquinat(‡)	$C_{24}H_{35}NO_5$	100mg	
Decoquinat D5 (7-ethoxy D5)				
CAS 1228182-55-5 DRE-C12097010	MW 422.5692 Decoquinat D5	$C_{24}^2H_{35}NO_5$	10mg	
Demeclocycline Hydrochloride				
CAS 64-73-3 DRE-C12128000	MW 501.314 Demeclocycline hydrochloride	$C_{21}H_{27}ClN_2O_8 \cdot ClH$	100mg	
Denatonium Benzoate (Bitrex (TM))				
CAS 3734-33-6 DRE-C10661000	MW 446.5812 Bitrex (TM)	$C_{21}H_{29}N_2O \cdot C_7H_5O_2$	100mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Desoxycortone (21-Hydroxyprogesterone)				
CAS 64-85-7 DRE-C14241050	MW 330.4611 21-Hydroxyprogesterone	C ₂₁ H ₃₀ O ₃	100mg	
Dexamethasone				
CAS 50-02-2 DRE-C12170400	MW 392.4611 Dexamethasone(‡)	C ₂₂ H ₂₈ FO ₅	100mg	
Dexamethasone 21-Acetate				
CAS 1177-87-3 DRE-C12170410	MW 434.4977 Dexamethasone 21-acetate	C ₂₄ H ₃₁ FO ₆	100mg	
Diatrizoate Sodium Dihydrate				
CAS 1185748-84-8 DRE-C12206000	MW 671.926 Diatrizoate sodium dihydrate	C ₁₁ H ₉ I ₃ N ₂ O ₄ ·Na·2H ₂ O	100mg	
Diatrizoic Acid				
CAS 117-96-4 DRE-C12207000	MW 613.9136 Diatrizoic acid	C ₁₁ H ₉ I ₃ N ₂ O ₄	100mg	
Diatrizoic Acid Meglumine Salt				
CAS 131-49-7 DRE-C12207010	MW 809.1272 Diatrizoic acid meglumine salt	C ₁₁ H ₉ I ₃ N ₂ O ₄ ·C ₇ H ₁₇ NO ₅	100mg	
Diaveridine				
CAS 5355-16-8 DRE-C12208000 DRE-L12208000AL	MW 260.2917 Diaveridine Diaveridine 10 µg/mL in Acetonitrile	C ₁₃ H ₁₆ N ₄ O ₂	100mg 10ml	
2-[(2,6-Dichlorophenyl)amino]benzaldehyde				
CAS 22121-58-0 DRE-C12537100	MW 266.1227 Diclofenac Impurity B	C ₁₃ H ₉ Cl ₂ NO	50mg	
(1RS)-1-(2,4-Dichlorophenyl)-2-(1H-imidazol-1-yl)ethanol				
CAS 24155-42-8 DRE-C12471000	MW 257.1159 1-(2,4-Dichlorophenyl)-2-imidazol-1-yl ethanol	C ₁₁ H ₁₀ Cl ₂ N ₂ O	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Diclazuril				
CAS 101831-37-2 DRE-C12533000	MW 407.638 Diclazuril(‡)	C ₁₇ H ₉ Cl ₃ N ₄ O ₂	100mg	
Diclazuril-methyl				
CAS 103337-71-9 DRE-C12533200	MW 421.6645 Diclazuril-methyl	C ₁₈ H ₁₁ Cl ₃ N ₄ O ₂	10mg	
Diclofenac				
CAS 15307-86-5 DRE-C12537000	MW 296.1486 Diclofenac acid(‡)	C ₁₄ H ₁₁ Cl ₂ NO ₂	100mg	
Diclofenac D4 (phenyl D4)				
CAS 153466-65-0 DRE-XA12537010AC	MW 300.1733 Diclofenac D4 acid (phenyl D4) 100 µg/mL in Acetonitrile (*)	C ₁₄ ² H ₄ H ₇ Cl ₂ NO ₂	1.1ml	
Diclofenac Methyl Ester				
CAS 15307-78-5 DRE-C12537400	MW 310.1752 Diclofenac methyl ester	C ₁₅ H ₁₃ Cl ₂ NO ₂	10mg	
Diclofenac Sodium				
CAS 15307-79-6 DRE-C12537500	MW 318.1305 Diclofenac sodium	C ₁₄ H ₁₀ Cl ₂ NO ₂ ·Na	100mg	
Dicloxacinil Sodium Monohydrate				
CAS 13412-64-1 DRE-C12560500	MW 510.3235 Dicloxacinil sodium hydrate	C ₁₉ H ₁₆ Cl ₂ N ₃ O ₅ S·Na·H ₂ O	100mg	
Dienestrol				
CAS 84-17-3 DRE-C12598000	MW 266.3343 Dienestrol	C ₁₈ H ₁₈ O ₂	25mg	
Dienestrol Diacetate				
CAS 84-19-5 DRE-C12598020	MW 350.4077 Dienestrol diacetate	C ₂₂ H ₂₂ O ₄	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Diethyl Fumarate				
CAS 623-91-6 DRE-C12606600	MW 172.1785 Diethyl fumarate	$C_8H_{12}O_4$	250mg	
Diethylcarbamazine Citrate				
CAS 1642-54-2 DRE-C12605700	MW 391.4168 Diethylcarbamazine citrate	$C_{10}H_{21}N_3O \cdot C_6H_8O_7$	100mg	
Diethylstilbestrol				
CAS 56-53-1 DRE-C12607000	MW 268.3502 Diethylstilbestrol(‡)	$C_{18}H_{20}O_2$	250mg	
Difloxacin Hydrochloride				
CAS 91296-86-5 DRE-C12627000	MW 435.8516 Difloxacin hydrochloride(‡)	$C_{21}H_{19}F_2N_3O_3 \cdot ClH$	100mg	
Digitoxin				
CAS 71-63-6 DRE-C12633800	MW 764.9391 Digitoxin	$C_{41}H_{64}O_{13}$	100mg	
Digoxin				
CAS 20830-75-5 DRE-C12633850	MW 780.9385 Digoxin	$C_{41}H_{64}O_{14}$	100mg	
Dihydrostreptomycin Sulfate				
CAS 5490-27-7 DRE-C12635300	MW 1461.4153 Dihydrostreptomycin sesquisulfate	$2C_{21}H_{41}N_7O_{12} \cdot 3H_2O_4S$	100mg	
Diltiazem Hydrochloride				
CAS 33286-22-5 DRE-C12645000	MW 450.9788 Diltiazem hydrochloride	$C_{22}H_{26}N_2O_4S \cdot ClH$	100mg	
Dimethicone 1000				
DRE-C12679100	Dimethicone 1000		250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Dimeticone (Dimethicone)				
CAS 9006-65-9 DRE-C12679000	MW 236.5315 Dimethicone 350	$C_6H_{18}OSi_2(C_2H_6OSi)_n$	250mg	
Dimetridazole				
CAS 551-92-8 DRE-C12772000	MW 141.128 Dimetridazole(‡)	$C_5H_7N_3O_2$	250mg	
Dimetridazole D3 (N-methyl D3)				
CAS 64678-69-9 DRE-C12772010 DRE-XA12772010AC	MW 144.1465 Dimetridazole D3 Dimetridazole D3 100 µg/mL in Acetone	$C_5^2H_3H_4N_3O_2$	10mg 1.1ml	
Dimetridazole-2-hydroxy D3				
CAS 1015855-78-3 DRE-C12772051	MW 160.1459 Dimetridazole-2-hydroxy D3	$C_5^2H_3H_4N_3O_3$	10mg	
Dimetridazole-2-hydroxy ((1-Methyl-5-nitroimidazol-2-yl)methanol)				
CAS 936-05-0 DRE-C12772050 DRE-XA12772050AC	MW 157.1274 Dimetridazole-2-hydroxy Dimetridazole-2-hydroxy 100 µg/mL in Acetone	$C_5H_7N_3O_3$	10mg 1.1ml	
Diminazene Aceturate				
CAS 908-54-3 DRE-C12773000	MW 515.5224 Diminazene aceturate(‡)	$C_{14}H_{15}N_7 \cdot 2C_4H_7NO_3$	100mg	
Diphenhydramine Hydrochloride				
CAS 147-24-0 DRE-C12881000	MW 291.8157 Diphenhydramine hydrochloride	$C_{17}H_{21}NO \cdot ClH$	100mg	
Disulfiram				
CAS 97-77-8 DRE-C12975000	MW 296.5392 Disulfiram	$C_{10}H_{20}N_2S_4$	250mg	
Dithranol				
CAS 1143-38-0 DRE-C13013900	MW 226.2274 Dithranol	$C_{14}H_{10}O_3$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Dopamine Hydrochloride				
CAS 62-31-7 DRE-C13082000	MW 189.6394 Dopamine hydrochloride	$C_8H_{11}NO_2 \cdot ClH$	100mg	
Doramectin				
CAS 117704-25-3 DRE-CA13083000	MW 899.1142 Doramectin (*)	$C_{50}H_{74}O_{14}$	100mg	
Doxepin Hydrochloride				
CAS 1229-29-4 DRE-C13084200	MW 315.8371 Doxepin hydrochlorid	$C_{19}H_{21}NO \cdot ClH$	100mg	
Doxycycline Hyclate				
CAS 24390-14-5 DRE-C13084280	MW 1025.8747 Doxycycline hyclate	$2C_{22}H_{24}N_2O_8 \cdot C_2H_6O \cdot 2ClH \cdot H_2O$	100mg	
Econazole Nitrate				
CAS 24169-02-6 DRE-C13101100	MW 444.6963 Econazol nitrate	$C_{18}H_{15}Cl_3N_2O \cdot HNO_3$	100mg	
Embutramide				
CAS 15687-14-6 DRE-C13117200	MW 293.4012 Embutramide	$C_{17}H_{27}NO_3$	25mg	
Enalapril Maleate				
CAS 76095-16-4 DRE-C13119000	MW 492.5189 Enalapril maleate	$C_{20}H_{28}N_2O_5 \cdot C_4H_4O_4$	100mg	
Enoxacin				
CAS 74011-58-8 DRE-C13166950	MW 320.3189 Enoxacin(†)	$C_{15}H_{17}FN_4O_3$	100mg	
Enoxacin Sesquihydrate				
CAS 84294-96-2 DRE-C13167000	MW 694.6836 Enoxacin sesquihydrate	$2C_{15}H_{17}FN_4O_3 \cdot 3H_2O$	100mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Enrofloxacin				
CAS 93106-60-6 DRE-C13170000	MW 359.3947 Enrofloxacin(±)	C ₁₉ H ₂₂ FN ₃ O ₃	100mg	
Enrofloxacin D5 Hydrochloride (ethyl d5)				
CAS n/a DRE-C13170100	MW 400.8864 Enrofloxacin D5 hydrochloride(±)	C ₁₉ ² H ₅ H ₁₇ FN ₃ O ₃ ·ClH	10mg	
4-Epianhydrotetracycline hydrochloride				
CAS 4465-65-0 DRE-C13174500	MW 462.8802 4-Epianhydrotetracycline hydrochloride	C ₂₂ H ₂₂ N ₂ O ₇ ·ClH	10mg	
Epigallocatechin-3-gallate				
CAS 989-51-5 DRE-C13176500	MW 458.3717 Epigallocatechin-3-gallate	C ₂₂ H ₁₈ O ₁₁	100mg	
Eprinomectin				
CAS 123997-26-2 DRE-CA13187000 DRE-XA13187000AL	MW 1814.231 Eprinomectin(±) Eprinomectin 100 µg/mL in Acetonitrile	C ₅₀ H ₇₅ NO ₁₄ ·C ₄₉ H ₇₃ NO ₁₄	100mg 1ml	
Equilenin				
CAS 517-09-9 DRE-XA13193000AL	MW 266.3343 Equilenin 100 µg/mL in Acetonitrile	C ₁₈ H ₁₈ O ₂	1ml	
Equilin				
CAS 474-86-2 DRE-C13193050	MW 268.3502 Equilin	C ₁₈ H ₂₀ O ₂	100mg	
Erythromycin				
CAS 114-07-8 DRE-C13203490	MW 733.9268 Erythromycin (mixture A,B,C)	C ₃₇ H ₆₇ NO ₁₃	100mg	
Erythromycin Ethylsuccinate				
CAS 1264-62-6 DRE-C13203520	MW 862.0527 Erythromycin ethyl succinate	C ₄₃ H ₇₅ NO ₁₆	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Erythromycin Stearate				
CAS 643-22-1 DRE-C13203550	MW 1018.404 Erythromycin stearate	$C_{37}H_{67}NO_{13} \cdot C_{18}H_{36}O_2$	100mg	
Esomeprazole Magnesium Trihydrate				
CAS 217087-09-7 DRE-C13211700	MW 767.1671 Esomeprazole magnesium trihydrate	$2C_{17}H_{18}N_3O_3S \cdot Mg \cdot 3H_2O$	50mg	
Estradiol (17β-Estradiol)				
CAS 50-28-2 DRE-C13213100 DRE-XA13213100AL	MW 272.382 17-β-Estradiol 17-β-Estradiol hemihydrate 100 µg/mL in Acetonitrile	$C_{18}H_{24}O_2$	250mg 1ml	
Estradiol Benzoate (17β-Estradiol 3-benzoate)				
CAS 50-50-0 DRE-C13213120	MW 376.488 17-β-Estradiol 3-benzoate	$C_{26}H_{26}O_3$	100mg	
Estradiol Diacetate (17β-Estradiol 3,17-diacetate)				
CAS 3434-88-6 DRE-C13213130	MW 356.4553 17-β-Estradiol 3,17-diacetate	$C_{22}H_{26}O_4$	100mg	
Estradiol Valerate (17β-Estradiol 17-valerate)				
CAS 979-32-8 DRE-C13213180	MW 356.4984 17-β-Estradiol 17-valerate	$C_{23}H_{32}O_3$	100mg	
17-epi-Estradiol (17α-Estradiol)				
CAS 57-91-0 DRE-C13213000 DRE-XA13213000AL	MW 272.382 17-α-Estradiol 17-α-Estradiol 100 µg/mL in Acetonitrile	$C_{18}H_{24}O_2$	100mg 1ml	
Estriol				
CAS 50-27-1 DRE-C13213200 DRE-XA13213200AL	MW 288.3814 Estriol(±) Estriol 100 µg/mL in Acetonitrile	$C_{18}H_{24}O_3$	100mg 1ml	
Estrone				
CAS 53-16-7 DRE-C13213230 DRE-XA13213230AL	MW 270.3661 Estrone(±) Estrone 100 µg/mL in Acetonitrile	$C_{18}H_{22}O_2$	100mg 1ml	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Ethinylestradiol (17α-Ethinylestradiol)				
CAS 57-63-6	MW 296.4034	C ₂₀ H ₂₄ O ₂		
DRE-C13245100	17α-Ethinylestradiol(‡)		250mg	
DRE-XA13245100AL	17α-Ethinylestradiol 100 µg/mL in Acetonitrile		1ml	
Ethisterone				
CAS 434-03-7	MW 312.4458	C ₂₁ H ₂₈ O ₂		
DRE-C13283000	Ethisterone		100mg	
Ethopabate				
CAS 59-06-3	MW 237.2518	C ₁₂ H ₁₅ NO ₄		
DRE-C13295000	Ethopabate		100mg	
Etiocholan-3α-ol-17-one				
CAS 53-42-9	MW 290.4403	C ₁₉ H ₃₀ O ₂		
DRE-C13356500	Etiocholan-3α-ol-17-one		10mg	
Etoricoxib				
CAS 202409-33-4	MW 358.8419	C ₁₈ H ₁₅ ClN ₂ O ₂ S		
DRE-C13365000	Etoricoxib		10mg	
Febantel				
CAS 58306-30-2	MW 446.4769	C ₂₀ H ₂₂ N ₄ O ₆ S		
DRE-C13407000	Febantel(‡)		100mg	
Fenbendazole				
CAS 43210-67-9	MW 299.3476	C ₁₅ H ₁₃ N ₃ O ₂ S		
DRE-C13446000	Fenbendazole(‡)		100mg	
Fenbendazole Sulphone (Oxfendazole-sulfone)				
CAS 54029-20-8	MW 331.3464	C ₁₅ H ₁₃ N ₃ O ₄ S		
DRE-C15783010	Oxfendazole-sulfone		25mg	
Fenfluramine Hydrochloride				
CAS 404-82-0	MW 267.7183	C ₁₂ H ₁₆ F ₃ N·ClH		
DRE-XA13468000AL	Fenfluramine HCl 100 µg/mL in Acetonitrile(‡)		1ml	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Fenofibrate				
CAS 49562-28-9 DRE-C13486000	MW 360.8313 Fenofibrate	$C_{20}H_{21}ClO_4$	100mg	
Fenoprofen Calcium Dihydrate				
CAS 71720-56-4 DRE-C13488000	MW 558.6324 Fenoprofen calcium dihydrate	$2C_{15}H_{13}O_3 \cdot Ca \cdot 2H_2O$	100mg	
Fenoterol Hydrobromide				
CAS 1944-12-3 DRE-XA13497000AL	MW 384.2649 Fenoterol hydrobromide 100 µg/mL in Acetonitrile(±)	$C_{17}H_{21}NO_4 \cdot BrH$	1ml	
Finasteride				
CAS 98319-26-7 DRE-C13644200	MW 372.5441 Finasteride(±)	$C_{23}H_{36}N_2O_2$	100mg	
Flavomycin				
CAS 11015-37-5 DRE-E13655300	MW 1583.5674 Flavomycin	$C_{69}H_{107}N_4O_{35}P$	250mg	
Fleroxacin				
CAS 79660-72-3 DRE-C13658000	MW 369.3383 Fleroxacin(±)	$C_{17}H_{18}F_3N_3O_3$	100mg	
Florfenicol				
CAS 73231-34-2 DRE-C13665000	MW 358.2133 Florfenicol(±)	$C_{12}H_{14}Cl_2FNO_4S$	250mg	
Fluazuron				
CAS 86811-58-7 DRE-C13672000	MW 506.2097 Fluazuron(±)	$C_{20}H_{10}Cl_2F_5N_3O_3$	100mg	
Flubendazole				
CAS 31430-15-6 DRE-C13678000	MW 313.2832 Flubendazole(±)	$C_{16}H_{12}FN_3O_3$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Flubendazole D3 (methyl D3)				
CAS 1173021-08-3 DRE-C13678010	MW 316.3017 Flubendazole D3	C ₁₆ ² H ₁₃ H ₉ FN ₃ O ₃	10mg	
Flucloxacillin Sodium				
CAS 1847-24-1 DRE-C13696000	MW 475.8536 Flucloxacillin sodium	C ₁₉ H ₁₆ ClFN ₃ O ₅ S·Na	250mg	
Fluconazole				
CAS 86386-73-4 DRE-C13697500	MW 306.2708 Fluconazole(‡)	C ₁₃ H ₁₂ F ₂ N ₆ O	100mg	
Fluconazole D4 (bismethylene D4)				
CAS 1124197-58-5 DRE-XA13697510AC	MW 310.2954 Fluconazole D4 (bismethylene D4) 100 µg/mL in Acetone	C ₁₃ ² H ₄ H ₈ F ₂ N ₆ O	1.1ml	
Flumequine				
CAS 42835-25-6 DRE-C13718000	MW 261.2484 Flumequin(‡)	C ₁₄ H ₁₂ FNO ₃	250mg	
Flumetasone				
CAS 2135-17-3 DRE-C13718500	MW 410.4515 Flumethasone	C ₂₂ H ₂₆ F ₂ O ₅	100mg	
Flunixin				
CAS 38677-85-9 DRE-C13726900	MW 296.2445 Flunixin	C ₁₄ H ₁₁ F ₃ N ₂ O ₂	10mg	
Flunixin D3 (methyl D3)				
CAS 1015856-60-6 DRE-C13727010	MW 299.263 Flunixin D3 (methyl D3)	C ₁₄ ² H ₃ H ₈ F ₃ N ₂ O ₂	10mg	
Flunixin Meglumine				
CAS 42461-84-7 DRE-C13727000	MW 491.4581 Flunixin meglumine	C ₁₄ H ₁₁ F ₃ N ₂ O ₂ ·C ₇ H ₁₇ NO ₅	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Fluoxetine Hydrochloride				
CAS 56296-78-7 DRE-C13801500	MW 345.7871 Fluoxetine hydrochloride(±)	C ₁₇ H ₁₈ F ₃ NO·ClH	10mg	
Flurandrenolide				
CAS 1524-88-5 DRE-C13807000	MW 436.5136 Flurandrenolide (*)	C ₂₄ H ₃₃ FO ₆	250mg	
Fluticasone Propionate				
CAS 80474-14-2 DRE-C13863000	MW 500.5708 Fluticasone propionate	C ₂₈ H ₃₁ F ₃ O ₅ S	10mg	
Folic Acid				
CAS 59-30-3 DRE-C13888000	MW 441.3975 Folic acid(±)	C ₁₉ H ₁₉ N ₇ O ₆	250mg	
4-Formylaminophenazone				
CAS 1672-58-8 DRE-C13924000	MW 231.2505 4-Formylaminoantipyrine	C ₁₂ H ₁₃ N ₃ O ₂	10mg	
Fosfomycin Calcium				
CAS 26016-98-8 DRE-C13941000	MW 176.1212 Fosfomycin calcium	C ₃ H ₅ O ₄ P·Ca	100mg	
Fosinopril				
CAS 98048-97-6 DRE-C13942000	MW 563.6625 Fosinopril	C ₃₀ H ₄₆ NO ₇ P	100mg	
Furaltadone hydrochloride				
CAS 3759-92-0 DRE-C13963000	MW 360.7503 Furaltadone hydrochloride	C ₁₃ H ₁₆ N ₄ O ₆ ·ClH	100mg	
Furazolidone				
CAS 67-45-8 DRE-C13970200	MW 225.1583 Furazolidone(±)	C ₆ H ₇ N ₃ O ₅	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Furosemide				
CAS 54-31-9 DRE-C13985000	MW 330.7441 Furosemid	C ₁₂ H ₁₁ ClN ₂ O ₅ S	250mg	
Gemfibrozil				
CAS 25812-30-0 DRE-C13999000	MW 250.3334 Gemfibrozil	C ₁₅ H ₂₂ O ₃	250mg	
Genistein				
CAS 446-72-0 DRE-C13999800	MW 270.2369 Genistein	C ₁₅ H ₁₀ O ₅	250mg	
Genite				
CAS 97-16-5 DRE-C14000000	MW 303.1611 Genite	C ₁₂ H ₈ Cl ₂ O ₃ S	250mg	
Gentamicin Sulfate (Gentamycin sulfate)				
CAS 1405-41-0 DRE-C14000200	MW 1488.785 Gentamycin-2,5-sulfate hydrate	C ₂₁ H ₄₃ N ₅ O ₇ ·C ₂₀ H ₄₁ N ₅ O ₇ ·C ₁₉ H ₃₉ N ₅ O ₇ ·H ₂ O ₄ S	250mg	
Geraniol				
CAS 106-24-1 DRE-C14010000 DRE-GA09010076IP	MW 154.2493 Geraniol Geraniol 1000 µg/mL in Isopropanol(‡)	C ₁₀ H ₁₈ O	250mg 1ml	
Glibenclamide				
CAS 10238-21-8 DRE-C14025000	MW 494.0035 Glibenclamide	C ₂₃ H ₂₈ ClN ₃ O ₅ S	250mg	
Glycidol				
CAS 556-52-5 DRE-CA14036950	MW 74.0785 Glycidol	C ₃ H ₆ O ₂	1ml	
Glycolic Acid				
CAS 79-14-1 DRE-C14037500	MW 76.0514 Glycolic acid	C ₂ H ₄ O ₃	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Glycopyrronium Bromide				
CAS 51186-83-5 DRE-C14038050	MW 398.3345 Glycopyrronium bromide	C ₁₉ H ₂₈ NO ₃ ·Br	100mg	AND Enantiomer 
Gossypol				
CAS 303-45-7 DRE-C14056200	MW 518.5544 Gossypol(±)	C ₃₀ H ₃₀ O ₈	100mg	
Griseofulvin				
CAS 126-07-8 DRE-C14056500	MW 352.7663 Griseofulvin(±)	C ₁₇ H ₁₇ ClO ₆	250mg	
Guaifenesin				
CAS 93-14-1 DRE-C14056900	MW 198.2158 Guaifenesin	C ₁₀ H ₁₄ O ₄	250mg	
Haloxon (O,O-Bis(2-chloroethyl) O-(3-chloro-4-methyl-7-coumarinyl) phosphate)				
CAS 321-55-1 DRE-C14059800	MW 415.5901 Haloxon	C ₁₄ H ₁₄ Cl ₃ O ₆ P	10mg	
Hexestrol				
CAS 84-16-2 DRE-C14202800	MW 270.3661 Hexestrol(±)	C ₁₈ H ₂₂ O ₂	100mg	
Histamine Dihydrochloride				
CAS 56-92-8 DRE-C14213050	MW 184.0669 Histamine dihydrochloride(±)	C ₅ H ₉ N ₃ ·2ClH	250mg	
Hydrochlorothiazide				
CAS 58-93-5 DRE-C14223500	MW 297.7391 Hydrochlorothiazide(±)	C ₇ H ₈ ClN ₃ O ₄ S ₂	100mg	
Hydrocortisone (Cortisol)				
CAS 50-23-7 DRE-C14224000	MW 362.4599 Hydrocortisone(±)	C ₂₁ H ₃₀ O ₅	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Hydrocortisone Acetate				
CAS 50-03-3 DRE-C14224020	MW 404.4966 Hydrocortisone acetate	C ₂₃ H ₃₂ O ₆	250mg	
Hydrocortisone 17-Butyrate				
CAS 13609-67-1 DRE-C14224030	MW 432.5497 Hydrocortisone 17-Butyrate	C ₂₅ H ₃₆ O ₆	50mg	
Hydroquinone Monobenzylether (4-(Benzyloxy)phenol)				
CAS 103-16-2 DRE-C10572800	MW 200.2332 4-(Benzyloxy)phenol	C ₁₃ H ₁₂ O ₂	100mg	
1-Hydroxyibuprofen (Ibuprofen-1-hydroxy)				
CAS 53949-53-4 DRE-C14278150	MW 222.2802 Ibuprofen-1-hydroxy	C ₁₃ H ₁₈ O ₃	10mg	
2-Hydroxyibuprofen (Ibuprofen-2-hydroxy)				
CAS 51146-55-5 DRE-C14278160	MW 222.2802 Ibuprofen-2-hydroxy	C ₁₃ H ₁₈ O ₃	10mg	
Hydroxyprogesterone Acetate (17α-Hydroxyprogesterone 17-acetate)				
CAS 302-23-8 DRE-C14241010	MW 372.4978 17-α-Hydroxyprogesterone 17-acetate	C ₂₃ H ₃₂ O ₄	100mg	
17α-Hydroxyprogesterone				
CAS 68-96-2 DRE-C14241000	MW 330.4611 17-α-Hydroxyprogesterone	C ₂₁ H ₃₀ O ₃	100mg	
Hygromycin B				
CAS 31282-04-9 DRE-C14260000	MW 527.5201 Hygromycin B	C ₂₀ H ₃₇ N ₃ O ₁₃	100mg	
Hyoscine Hydrobromide Trihydrate (Scopolamine hydrobromide trihydrate)				
CAS 6533-68-2 DRE-C16915000	MW 438.3107 Scopolamine hydrobromide trihydrate	C ₁₇ H ₂₁ NO ₄ ·BrH·3H ₂ O	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Ibuprofen				
CAS 15687-27-1 DRE-C14278000	MW 206.2808 Ibuprofen(±)	C ₁₃ H ₁₈ O ₂	250mg	
Ifosfamide				
CAS 3778-73-2 DRE-C14279000	MW 261.086 Ifosfamide	C ₇ H ₁₅ Cl ₂ N ₂ O ₂ P	100mg	
Imidazole				
CAS 288-32-4 DRE-C14283950	MW 68.0773 Imidazole	C ₃ H ₄ N ₂	250mg	
Imidocarb				
CAS 27885-92-3 DRE-C14284500	MW 348.4017 Imidocarb(±)	C ₁₉ H ₂₀ N ₆ O	50mg	
Imidocarb Dipropionate				
CAS 55750-06-6 DRE-C14284520	MW 496.5588 Imidocarb dipropionate(±)	C ₁₉ H ₂₀ N ₆ O · 2C ₃ H ₆ O ₂	50mg	
Iminostilbene (5H-Dibenzo[b,f]azepine)				
CAS 256-96-2 DRE-C14285500	MW 193.2438 Iminostilbene	C ₁₄ H ₁₁ N	100mg	
Indometacin (Indomethacine)				
CAS 53-86-1 DRE-C14325000	MW 357.7876 Indomethacine	C ₁₉ H ₁₆ ClNO ₄	250mg	
Inosine				
CAS 58-63-9 DRE-C14328000	MW 268.2261 Inosine	C ₁₀ H ₁₂ N ₄ O ₅	100mg	
5'-Inosinic Acid Disodium Salt Hydrate				
CAS 352195-40-5 DRE-C14328020	MW 410.1849 5'-Inosinic acid disodium salt hydrate	C ₁₀ H ₁₁ N ₄ O ₈ P · 2Na · H ₂ O	100mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Iodipamide				
CAS 606-17-7 DRE-C14329000	MW 1139.7618 Iodipamide	C ₂₀ H ₁₄ I ₆ N ₂ O ₆	100mg	
Iodixanol				
CAS 92339-11-2 DRE-C14330000	MW 1550.1819 Iodixanol	C ₃₅ H ₄₄ I ₆ Na ₆ O ₁₅	100mg	
Iohexol				
CAS 66108-95-0 DRE-C14347000	MW 821.1379 Iohexol	C ₁₉ H ₂₆ I ₃ N ₃ O ₉	100mg	
Iomeprol				
CAS 78649-41-9 DRE-C14348000	MW 777.0853 Iomeprol	C ₁₇ H ₂₂ I ₃ N ₃ O ₈	100mg	
Iopamidol				
CAS 60166-93-0 DRE-C14348400	MW 777.0853 Iopamidol(±)	C ₁₇ H ₂₂ I ₃ N ₃ O ₈	100mg	
Iopanoic Acid				
CAS 96-83-3 DRE-C14348500	MW 570.9319 Iopanoic acid	C ₁₁ H ₁₂ I ₃ NO ₂	100mg	
Iopromide				
CAS 73334-07-3 DRE-C14348600	MW 791.1119 Iopromid	C ₁₈ H ₂₄ I ₃ N ₃ O ₈	100mg	
Iothalamate Meglumine (Iotalamic acid meglumine salt)				
CAS 13087-53-1 DRE-C14349000	MW 809.1272 Iotalamic acid meglumine salt	C ₁₁ H ₉ I ₃ N ₂ O ₄ ·C ₇ H ₁₇ NO ₅	100mg	
Ipronidazole				
CAS 14885-29-1 DRE-C14370700	MW 169.1811 Ipronidazole(±)	C ₇ H ₁₁ N ₃ O ₂	50mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Iprnidazole D3 (N-methyl D3)				
CAS 1015855-83-0 DRE-C14370701	MW 172.1996 Iprnidazole D3	C ₇ ² H ₉ H ₈ N ₃ O ₂	10mg	
Iprnidazole-hydroxy				
CAS 35175-14-5 DRE-C14370720	MW 185.1805 Iprnidazole-hydroxy	C ₇ H ₁₁ N ₃ O ₃	10mg	
Iprnidazole-hydroxy D3				
CAS 1156508-86-9 DRE-C14370721	MW 188.199 Iprnidazole-hydroxy D3	C ₇ ² H ₉ H ₈ N ₃ O ₃	10mg	
Itraconazole				
CAS 84625-61-6 DRE-C14485000	MW 705.6334 Itraconazole	C ₃₅ H ₃₈ Cl ₂ N ₈ O ₄	100mg	
Ivermectin				
CAS 70288-86-7 DRE-CA14488000 DRE-XA14488000AL	MW 1736.1589 Ivermectine Ivermectine 100 µg/mL in Acetonitrile	C ₄₈ H ₇₄ O ₁₄ ·C ₄₇ H ₇₂ O ₁₄	100mg 1ml	
Josamycin				
CAS 16846-24-5 DRE-LA14495000AL	MW 827.995 Josamycin 10 µg/mL in Acetonitrile (*)	C ₄₂ H ₆₉ NO ₁₅	1ml	
Kanamycin Sulfate				
CAS 25389-94-0 DRE-C14505520	MW 582.5771 Kanamycin sulfate	C ₁₈ H ₃₈ N ₄ O ₁₁ ·H ₂ O ₄ S	100mg	
Ketamine Hydrochloride				
CAS 1867-66-9 DRE-C14531000	MW 274.1862 Ketamine Hydrochloride	C ₁₃ H ₁₆ ClNO·ClH	100mg	
Ketoconazole				
CAS 65277-42-1 DRE-C14532000	MW 531.4309 Ketoconazole	C ₂₆ H ₂₈ Cl ₂ N ₄ O ₄	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Ketoprofen				
CAS 22071-15-4 DRE-C14532100	MW 254.2806 Ketoprofen(±)	C ₁₆ H ₁₄ O ₃	100mg	
(±)-Ketoprofen D3 (methyl D3)				
CAS 159490-55-8 DRE-XA14532110AL	MW 257.299 (±)-Ketoprofen D3 (propionic D3 acid) 100 µg/mL in Acetonitrile	C ₁₆ ² H ₁₃ H ₁₁ O ₃	1ml	
Lactic acid				
CAS 50-21-5 DRE-C14582000	MW 90.0779 Lactic acid	C ₃ H ₆ O ₃	100mg	
Lansoprazole				
CAS 103577-45-3 DRE-C14592500	MW 369.3615 Lansoprazole	C ₁₆ H ₁₄ F ₃ N ₃ O ₂ S	250mg	
Lasalocid A sodium salt				
CAS 25999-20-6 DRE-XA14593000AL	MW 612.7696 Lasalocid A sodium salt 100 µg/mL in Acetonitrile(±)	C ₃₄ H ₅₃ O ₈ ·Na	1ml	
Leucomycin				
CAS 1392-21-8 DRE-C14629600 DRE-XA14629600AL	MW n/a Leucomycin Leucomycin 100 µg/mL in Acetonitrile		250mg 1ml	No Structure
Levamisole				
CAS 14769-73-4 DRE-C14629690	MW 204.2914 Levamisole	C ₁₁ H ₁₂ N ₂ S	100mg	
Levamisole Hydrochloride				
CAS 16595-80-5 DRE-C14629700	MW 240.7523 Levamisole hydrochloride	C ₁₁ H ₁₂ N ₂ S·ClH	250mg	
Levocarnitine (L-Carnitin)				
CAS 541-15-1 DRE-C11045500	MW 161.1989 L-Carnitin	C ₇ H ₁₅ NO ₃	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Levosulpiride ((S)-(-)-Sulpiride)				
CAS 23672-07-3 DRE-C17027000	MW 341.4258 (S)-(-)-Sulpiride	C ₁₅ H ₂₃ N ₃ O ₄ S	100mg	
Levothyroxine (L-Thyroxine)				
CAS 51-48-9 DRE-C17575300	MW 776.87 L-Thyroxine	C ₁₅ H ₁₁ I ₄ NO ₄	100mg	
Lidocaine (Lignocaine)				
CAS 137-58-6 DRE-C14629790	MW 234.3373 Lignocaine base	C ₁₄ H ₂₂ N ₂ O	100mg	
Lidocaine Hydrochloride Monohydrate				
CAS 6108-05-0 DRE-C14629800	MW 288.8135 Lidocaine hydrochloride monohydrate	C ₁₄ H ₂₂ N ₂ O·ClH·H ₂ O	250mg	
Lincomycin Hydrochloride Monohydrate				
CAS 7179-49-9 DRE-C14635000	MW 461.0136 Lincomycin hydrochloride monohydrate	C ₁₈ H ₃₄ N ₂ O ₆ S·ClH·H ₂ O	250mg	
Lomefloxacin Hydrochloride				
CAS 98079-52-8 DRE-C14646000	MW 387.8088 Lomefloxacin hydrochloride(±)	C ₁₇ H ₁₉ F ₂ N ₃ O ₃ ·ClH	100mg	
Loperamide Hydrochloride				
CAS 34552-83-5 DRE-C14647000	MW 513.4985 Loperamide hydrochloride	C ₂₉ H ₃₃ ClN ₂ O ₂ ·ClH	100mg	
Mabuterol Hydrochloride				
CAS 54240-36-7 DRE-C14660000	MW 347.204 Mabuterol hydrochloride(±)	C ₁₃ H ₁₈ ClF ₃ N ₂ O·ClH	50mg	
Maduramicin ammonium salt				
CAS 84878-61-5 DRE-CA14670000	MW 934.1584 Maduramicin ammonium salt(±)	C ₄₇ H ₈₀ O ₁₇ ·H ₃ N	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Mapenterol hydrochloride				
CAS 54238-51-6 DRE-C14754000	MW 361.2305 Mapenterol hydrochloride	$C_{14}H_{20}ClF_3N_2O \cdot ClH$	5mg	
Marbofloxacin				
CAS 115550-35-1 DRE-C14755000	MW 362.3556 Marbofloxacin	$C_{17}H_{18}FN_4O_4$	100mg	
Mebendazole				
CAS 31431-39-7 DRE-C14798000	MW 295.2927 Mebendazole(±)	$C_{16}H_{13}N_3O_3$	100mg	
Mebendazole D3 (methyl D3)				
CAS 1173021-87-8 DRE-C14798010	MW 298.3112 Mebendazole D3	$C_{16}^2H_{13}H_10N_3O_3$	10mg	
Mebendazole-amine ((2-Amino-1H-benzimidazol-5-yl)phenylmethanone)				
CAS 52329-60-9 DRE-C14798015	MW 237.2566 Mebendazole-amine	$C_{14}H_{11}N_3O$	10mg	
Mebendazole-5-hydroxy				
CAS 60254-95-7 DRE-C14798020	MW 297.3086 Mebendazole-5-hydroxy	$C_{16}H_{13}N_3O_3$	50mg	
Mecillinam				
CAS 32887-01-7 DRE-C14803000	MW 325.4264 Mecillinam	$C_{15}H_{23}N_3O_3S$	100mg	
Meclofenamic Acid				
CAS 644-62-2 DRE-C14804500	MW 296.1486 Meclofenamic acid	$C_{14}H_{11}Cl_2NO_2$	10mg	
Meclofenoxate Hydrochloride				
CAS 3685-84-5 DRE-C14805000	MW 294.1743 Meclofenoxat hydrochloride	$C_{12}H_{16}ClNO_3 \cdot ClH$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Medroxyprogesterone				
CAS 520-85-4 DRE-C14852000	MW 344.4877 Medroxyprogesterone(‡)	C ₂₂ H ₃₂ O ₃	100mg	
Medroxyprogesterone 17-Acetate				
CAS 71-58-9 DRE-C14852015	MW 386.5244 Medroxyprogesterone-17-acetate	C ₂₄ H ₃₄ O ₄	100mg	
Mefenamic Acid				
CAS 61-68-7 DRE-C14860200	MW 241.2851 Mefenamic acid	C ₁₅ H ₁₃ NO ₂	250mg	
Megestrol 17-Acetate				
CAS 595-33-5 DRE-C14861300	MW 384.5085 Megestrol-17-acetate	C ₂₄ H ₃₂ O ₄	100mg	
Melatonin				
CAS 73-31-4 DRE-C14861500	MW 232.2783 Melatonin	C ₁₃ H ₁₆ N ₂ O ₂	250mg	
Melengestrol Acetate				
CAS 2919-66-6 DRE-C14861700	MW 396.5192 Melengestrol acetate	C ₂₅ H ₃₂ O ₄	100mg	
Meloxicam				
CAS 71125-38-7 DRE-C14862500	MW 351.4007 Meloxicam	C ₁₄ H ₁₃ N ₃ O ₄ S ₂	100mg	
Meloxicam D3 (2-methyl D3)				
CAS 942047-63-4 DRE-C14862510	MW 354.4192 Meloxicam D3 (2-methyl D3)	C ₁₄ ² H ₁₃ H ₁₀ N ₃ O ₄ S ₂	50mg	
Menbutone				
CAS 3562-99-0 DRE-C14864000	MW 258.2693 Menbutone	C ₁₅ H ₁₄ O ₄	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Mepyramine Maleate				
CAS 59-33-6 DRE-C14896000	MW 401.4562 Mepyramine maleate	$C_{17}H_{23}N_3O \cdot C_4H_4O_4$	250mg	
Mestranol				
CAS 72-33-3 DRE-C14915000	MW 310.4299 Mestranol	$C_{21}H_{26}O_2$	100mg	
Metacycline Hydrochloride				
CAS 3963-95-9 DRE-C14917000	MW 478.8796 Metacycline hydrochloride	$C_{22}H_{22}N_2O_8 \cdot ClH$	100mg	
Metamizole Sodium				
CAS 68-89-3 DRE-C14942000	MW 333.3386 Metamizol sodium	$C_{13}H_{16}N_2O_4S \cdot Na$	250mg	
D,L-Metanephrine Hydrochloride				
CAS 881-95-8 DRE-CA14946000	MW 233.6919 D,L-Metanephrine hydrochloride	$C_{10}H_{15}NO_3 \cdot ClH$	50mg	
Methapyrilene (2-[[2-(Dimethylamino)ethyl]-2-thenylamino]pyridine)				
CAS 91-80-5 DRE-GA09010392D DRE-GS09010392D	MW 261.3858 Methapyrilene Solution, 2000 µg/mL in Dichloromethane(‡) Methapyrilene Solution, 2000 µg/mL in Dichloromethane(‡)	$C_{14}H_{19}N_3S$	1ml 5x1ml	
Methapyrilene Hydrochloride				
CAS 135-23-9 DRE-C14996000	MW 297.8467 Methapyrilene hydrochloride	$C_{14}H_{19}N_3S \cdot ClH$	100mg	
Methazole				
CAS 20354-26-1 DRE-C15000000	MW 261.0615 Methazole	$C_9H_6Cl_2N_2O_3$	10mg	
D-Methorphan Hydrobromide				
CAS 125-69-9 DRE-C15053000	MW 352.3091 D-Methorphan hydrobromide	$C_{18}H_{25}NO \cdot BrH$	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
4-Methylimidazol				
CAS 822-36-6 DRE-C15088300	MW 82.1038 4-Methylimidazol	C ₄ H ₆ N ₂	100mg	
Methylprednisolone Acetate				
CAS 53-36-1 DRE-C15141720	MW 416.5073 Methylprednisolone Acetate	C ₂₄ H ₃₂ O ₆	25mg	
6α-Methylprednisolone				
CAS 83-43-2 DRE-C15141700 DRE-XA15141700AL	MW 374.4706 6-α-Methylprednisolone(±) 6-α-Methylprednisolone 100 µg/mL in Acetonitrile	C ₂₂ H ₃₀ O ₅	50mg 1ml	
N-Methylpyrrolidone (1-Methyl-2-pyrrolidone)				
CAS 872-50-4 DRE-C15143000	MW 99.1311 1-Methyl-2-pyrrolidon(±)	C ₅ H ₉ NO	5ml	
17α-Methyltestosterone				
CAS 58-18-4 DRE-C15144000 DRE-XA15144000AL	MW 302.451 17-α-Methyltestosterone(±) 17-α-Methyltestosterone 100 µg/mL in Acetonitrile	C ₂₀ H ₃₀ O ₂	100mg 1ml	
5-Methyl-1,3,4-thiadiazol-2-thiol				
CAS 29490-19-5 DRE-C15144100	MW 132.2073 5-Methyl-1,3,4-thiadiazol-2-thiol	C ₃ H ₄ N ₂ S ₂	100mg	
Metoclopramide Hydrochloride				
CAS 7232-21-5 DRE-C15165000	MW 336.2573 Metoclopramide hydrochloride	C ₁₄ H ₂₂ ClN ₃ O ₂ ·ClH	100mg	
Metoprolol Tartrate				
CAS 56392-17-7 DRE-C15176020	MW 684.8146 Metoprolol tartrat	2C ₁₅ H ₂₅ NO ₃ ·C ₄ H ₆ O ₆	250mg	
Metronidazole				
CAS 443-48-1 DRE-C15201000	MW 171.154 Metronidazole(±)	C ₆ H ₉ N ₃ O ₃	250mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Metronidazole D4 (ethylene D4)				
CAS 1261392-47-5 DRE-C15201001	MW 175.1786 Metronidazole D4 (ethylene D4)(‡)	C ₆ H ₄ H ₅ N ₃ O ₃	10mg	
Metronidazole-hydroxy D2				
CAS n/a DRE-C15201301	MW 189.1657 Metronidazole-hydroxy D2	C ₆ H ₂ H ₇ N ₃ O ₄	10mg	
Metronidazole-hydroxy (1-(2-Hydroxyethyl)-2-hydroxymethyl-5-nitroimidazole)				
CAS 4812-40-2 DRE-C15201300	MW 187.1534 Metronidazole-hydroxy(‡)	C ₆ H ₉ N ₃ O ₄	10mg	
Minocycline Hydrochloride Dihydrate				
CAS 128420-71-3 DRE-C15269000	MW 529.9679 Minocycline Hydrochloride Dihydrate (*)	C ₂₃ H ₂₇ N ₃ O ₇ ·ClH·2H ₂ O	250mg	
Mometasone Furoate				
CAS 83919-23-7 DRE-C15284000	MW 521.4295 Mometasone Furoate	C ₂₇ H ₃₀ Cl ₂ O ₆	250mg	
Monensin Sodium				
CAS 22373-78-0 DRE-C15291000	MW 692.8527 Monensin sodium salt	C ₃₆ H ₆₁ O ₁₁ ·Na	100mg	
Morantel Tartrate Monohydrate				
CAS 26155-31-7 DRE-C15327000	MW 370.4207 Morantel tartrate monohydrat	C ₁₂ H ₁₆ N ₂ S·C ₄ H ₆ O ₆	100mg	
Moxidectin				
CAS 113507-06-5 DRE-CA15335000	MW 639.8186 Moxidectin	C ₃₇ H ₅₃ NO ₈	25mg	
Mupirocin				
CAS 12650-69-0 DRE-C15345000	MW 500.6222 Mupirocin	C ₂₆ H ₄₄ O ₉	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
2-NP-AHD				
CAS 623145-57-3 DRE-C15654450	MW 248.1949 2-NP-AHD	C ₁₀ H ₈ N ₄ O ₄	10mg	
2-NP-AMAZ				
CAS 183193-59-1 DRE-C15654480	MW 334.3272 2-NP-AMAZ	C ₁₅ H ₁₈ N ₄ O ₅	10mg	
2-NP-AMAZ D5				
CAS 1173097-59-0 DRE-C15654481	MW 339.358 2-NP-AMAZ D5	C ₁₅ ² H ₈ H ₁₃ N ₄ O ₅	10mg	
2-NP-AOZ				
CAS 19687-73-1 DRE-C15654500	MW 235.1962 2-NP-AOZ(‡)	C ₁₀ H ₉ N ₃ O ₄	10mg	
2-NP-AOZ D4				
CAS 1007478-57-0 DRE-C15654501	MW 239.2208 2-NP-AOZ D4(‡)	C ₁₀ ² H ₄ H ₅ N ₃ O ₄	10mg	
2-NP-SCA				
CAS 16004-43-6 DRE-C15654520	MW 208.1741 2-NP-SCA	C ₈ H ₈ N ₄ O ₃	10mg	
Nadifloxacin				
CAS 124858-35-1 DRE-C15400500	MW 360.3794 Nadifloxacin	C ₁₉ H ₂₁ FN ₂ O ₄	100mg	
Nadolol				
CAS 42200-33-9 DRE-C15401000	MW 309.4006 Nadolol	C ₁₇ H ₂₇ NO ₄	100mg	
Nafcillin Sodium Monohydrate				
CAS 7177-50-6 DRE-C15402500	MW 454.4719 Nafcillin sodium salt monohydrate(‡)	C ₂₁ H ₂₁ N ₂ O ₅ S·Na·H ₂ O	100mg	









Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Nalidixic Acid				
CAS 389-08-2 DRE-C15412000	MW 232.2353 Nalidixic acid	C ₁₂ H ₁₂ N ₂ O ₃	100mg	
Nandrolone (19-Nortestosterone)				
CAS 434-22-0 DRE-C15652000	MW 274.3978 19-Nortestosterone	C ₁₈ H ₂₆ O ₂	100mg	
Nandrolone Decanoate (19-Nortestosterone 17-decanoate)				
CAS 360-70-3 DRE-C15652005	MW 428.6472 19-Nortestosterone 17-decanoate	C ₂₈ H ₄₄ O ₃	100mg	
Naproxen				
CAS 22204-53-1 DRE-C15483500	MW 230.2592 Naproxen(±)	C ₁₄ H ₁₄ O ₃	250mg	
Natamycin (Pimaricin)				
CAS 7681-93-8 DRE-C16208400	MW 665.7252 Pimaricin(±)	C ₃₃ H ₄₇ NO ₁₃	50mg	
Neomycin Sulfate				
CAS 1405-10-3 DRE-C15500900	MW n/a Neomycin sulfate		100mg	No Structure
Neostigmine Metilsulfate				
CAS 51-60-5 DRE-C15500960	MW 334.3886 Neostigmine metilsulfate	C ₁₂ H ₁₉ N ₂ O ₂ ·CH ₃ O ₄ S	250mg	
Nequinat				
CAS 13997-19-8 DRE-C15501300	MW 365.4223 Nequinat	C ₂₂ H ₂₃ NO ₄	10mg	
Nicarbazin				
CAS 330-95-0 DRE-C15508000	MW 426.3828 Nicarbazine(±)	C ₁₃ H ₁₀ N ₄ O ₅ ·C ₆ H ₈ N ₂ O	100mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Nicotine				
CAS 54-11-5	MW 162.2316	C ₁₀ H ₁₄ N ₂		
DRE-CA15520000	Nicotine(±)		500mg	
DRE-CR15520000	Nicotine		0.5g	
DRE-L15520000ME	Nicotine 10 µg/mL in Methanol		10ml	
DRE-XA15520000ME	Nicotine 100 µg/mL in Methanol		1ml	
Niflumic Acid				
CAS 4394-00-7	MW 282.218	C ₁₃ H ₉ F ₃ N ₂ O ₂		
DRE-C15522500	Niflumic acid		100mg	
Nifuroxazide				
CAS 965-52-6	MW 275.217	C ₁₂ H ₉ N ₃ O ₅		
DRE-C15523000	Nifuroxazide		100mg	
Nifuroxime (5-Nitrofurfural oxime)				
CAS 555-15-7	MW 156.0963	C ₅ H ₄ N ₂ O ₄		
DRE-C15523030	Nifuroxime		100mg	
Nifurpirinol				
CAS 13411-16-0	MW 246.2188	C ₁₂ H ₁₀ N ₂ O ₄		
DRE-C15523050	Nifurpirinol		50mg	
Nifursol				
CAS 16915-70-1	MW 365.2121	C ₁₂ H ₇ N ₅ O ₉		
DRE-C15523100	Nifursol		100mg	
Nimorazole				
CAS 6506-37-2	MW 226.2325	C ₉ H ₁₄ N ₄ O ₃		
DRE-C15527000	Nimorazole		10mg	
Nitarsons				
CAS 98-72-6	MW 247.0371	C ₆ H ₆ AsNO ₅		
DRE-C15533000	Nitarsons		100mg	
Nitrofural (Nitrofurazone)				
CAS 59-87-0	MW 198.1362	C ₆ H ₆ N ₄ O ₄		
DRE-C15571000	Nitrofurazone(±)		250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Nitrofurantoin				
CAS 67-20-9 DRE-C15570900	MW 238.157 Nitrofurantoin(±)	C ₈ H ₆ N ₄ O ₅	250mg	
Nitrovin Hydrochloride				
CAS 2315-20-0 DRE-C15616000	MW 396.7426 Nitrovin hydrochloride	C ₁₄ H ₁₂ N ₆ O ₆ ·ClH	100mg	
Nitroxinil (Nitroxynil)				
CAS 1689-89-0 DRE-C15617000	MW 290.0148 Nitroxynil(±)	C ₇ H ₃ IN ₂ O ₃	100mg	
Norethisterone Acetate				
CAS 51-98-9 DRE-C15645000	MW 340.4559 Norethisterone acetate	C ₂₂ H ₂₈ O ₃	100mg	
Norfloxacin				
CAS 70458-96-7 DRE-C15648000	MW 319.3308 Norfloxacin(±)	C ₁₆ H ₁₈ FN ₃ O ₃	100mg	
Norfloxacin D5 (ethyl D5)				
CAS 1015856-57-1 DRE-C15648010	MW 324.3616 Norfloxacin D5(±)	C ₁₆ ² H ₁₈ H ₁₃ FN ₃ O ₃	10mg	
19-Nortestosterone 17-Propionate				
CAS 7207-92-3 DRE-C15652010	MW 330.4611 19-Nortestosterone 17-propionate	C ₂₁ H ₃₀ O ₃	100mg	
Novobiocin Sodium				
CAS 1476-53-5 DRE-C15654200	MW 634.6062 Novobiocin sodium	C ₃₁ H ₃₅ N ₂ O ₁₁ ·Na	250mg	
Nystatin				
CAS 1400-61-9 DRE-C15661000	MW n/a Nystatin		250mg	No Structure










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Octisalate (Salicylic acid 2-ethyl-1-hexyl ester)				
CAS 118-60-5 DRE-C16903520	MW 250.3334 Salicylic acid-2-ethyl-1-hexyl ester(±)	C ₁₅ H ₂₂ O ₃	250mg	
Octopamine Hydrochloride				
CAS 770-05-8 DRE-CA15711650	MW 189.6394 (±)-Octopamine hydrochloride(±)	C ₈ H ₁₁ NO ₂ ·ClH	100mg	
Ofloxacin				
CAS 82419-36-1 DRE-C15717000	MW 361.3675 Ofloxacin(±)	C ₁₈ H ₂₀ FN ₃ O ₄	100mg	
Ofloxacin D3 (N-methyl D3)				
CAS 1173147-91-5 DRE-XA15717005AL	MW 364.386 Ofloxacin D3 100 µg/mL in Acetonitrile	C ₁₈ ² H ₃ H ₁₇ FN ₃ O ₄	1.1ml	
Ofloxacin D3 Hydrochloride (N-methyl D3)				
CAS 1173021-78-7 DRE-C15717010	MW 400.8469 Ofloxacin D3 hydrochloride	C ₁₈ ² H ₃ H ₁₇ FN ₃ O ₄ ·ClH	10mg	
Olaquinox				
CAS 23696-28-8 DRE-C15724000	MW 263.2493 Olaquinox(±)	C ₁₂ H ₁₃ N ₃ O ₄	100mg	
Oleandomycin phosphate, dihydrate				
CAS 7060-74-4 DRE-XA15726000AL	MW 785.8535 Oleandomycin phosphate 100 µg/mL in Acetonitrile	C ₃₅ H ₆₁ NO ₁₂ ·H ₃ O ₄ P	1ml	
Omeprazole				
CAS 73590-58-6 DRE-C15729000	MW 345.4161 Omeprazole	C ₁₇ H ₁₉ N ₃ O ₃ S	100mg	
Orbifloxacin				
CAS 113617-63-3 DRE-C15742000	MW 395.3756 Orbifloxacin(±)	C ₁₉ H ₂₀ F ₃ N ₃ O ₃	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Orciprenaline Sulfate (Metaproterenol hemisulfate salt)				
CAS 5874-97-5 DRE-C14947500	MW 520.5936 Metaproterenol hemisulfate salt	$2C_{11}H_{17}NO_3 \cdot H_2O_4S$	100mg	
Ormetoprim				
CAS 6981-18-6 DRE-C15745000	MW 274.3183 Ormetoprim	$C_{14}H_{18}N_4O_2$	250mg	
Ornidazole				
CAS 16773-42-5 DRE-C15746000	MW 219.6256 Ornidazole(‡)	$C_7H_{10}ClN_3O_3$	100mg	
Oseltamivir				
CAS 196618-13-0 DRE-C15751000	MW 312.4045 Oseltamivir (free base)	$C_{16}H_{28}N_2O_4$	10mg	
Oxacillin Sodium Monohydrate				
CAS 7240-38-2 DRE-C15755100	MW 441.4334 Oxacillin sodium salt hydrate	$C_{19}H_{18}N_3O_5 \cdot Na \cdot H_2O$	100mg	
Oxandrolone				
CAS 53-39-4 DRE-C15780800	MW 306.4397 Oxandrolone	$C_{19}H_{30}O_3$	10mg	
Oxcarbazepine				
CAS 28721-07-5 DRE-C15782500	MW 252.268 Oxcarbazepine	$C_{15}H_{12}N_2O_2$	50mg	
Oxfendazole				
CAS 53716-50-0 DRE-C15783000	MW 315.347 Oxfendazole(‡)	$C_{15}H_{13}N_3O_3S$	100mg	
Oxibendazole				
CAS 20559-55-1 DRE-C15784000	MW 249.2658 Oxibendazole	$C_{12}H_{15}N_3O_3$	100mg	








Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Oxyclozanide				
CAS 2277-92-1 DRE-C15793000	MW 401.4566 Oxyclozanide	C ₁₃ H ₆ Cl ₃ NO ₃	100mg	
Oxytetracycline Dihydrate				
CAS 6153-64-6 DRE-C15819990	MW 496.4645 Oxytetracycline dihydrate	C ₂₂ H ₂₄ N ₂ O ₉ ·2H ₂ O	250mg	
Padimate O				
CAS 21245-02-3 DRE-C15841000	MW 277.4018 Padimate O	C ₁₇ H ₂₇ NO ₂	100mg	
Paracetamol				
CAS 103-90-2 DRE-C15846000	MW 151.1626 Paracetamol(†)	C ₈ H ₉ NO ₂	250mg	
Paraldehyde				
CAS 123-63-7 DRE-C15848200	MW 132.1577 Paraldehyde	C ₆ H ₁₂ O ₃	5ml	
Parbendazole				
CAS 14255-87-9 DRE-C15892000	MW 247.293 Parbendazole	C ₁₃ H ₁₇ N ₃ O ₂	25mg	
Pefloxacin Mesilate Dihydrate (Pefloxacin methanesulfonate dihydrate)				
CAS 149676-40-4 DRE-C15905000	MW 465.4936 Pefloxacin methane sulfonate dihydrate	C ₁₇ H ₂₀ FN ₃ O ₃ ·CH ₄ O ₃ S·2H ₂ O	100mg	
(S)-Penbutolol ((S)-Penbutolol)				
CAS 38363-40-5 DRE-C15908910	MW 291.4284 (S)-Penbutolol	C ₁₈ H ₂₉ NO ₂	10mg	
Penicillin G Benzathine Tetrahydrate				
CAS 41372-02-5 DRE-C10532500	MW 981.1848 Benzathine penicilline G tetrahydrate	C ₁₆ H ₂₀ N ₂ ·2C ₁₆ H ₁₈ N ₂ O ₄ S·4H ₂ O	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Pentamidine Diisetonate				
CAS 140-64-7 DRE-C15976500	MW 592.6827 Pentamidine isethionate salt	$C_{19}H_{24}N_4O_2 \cdot 2C_2H_6O_4S$	100mg	
Pentoxifylline				
CAS 6493-05-6 DRE-C15981800	MW 278.307 Pentoxifylline	$C_{13}H_{18}N_4O_3$	100mg	
Pergolide Mesilate				
CAS 66104-23-2 DRE-C15989500	MW 410.5938 Pergolide mesilate	$C_{19}H_{26}N_2S \cdot CH_4O_3S$	100mg	
Phenacetin				
CAS 62-44-2 DRE-C16003000	MW 179.2157 Phenacetin	$C_{10}H_{13}NO_2$	250mg	
Phenazone				
CAS 60-80-0 DRE-C16003100	MW 188.2258 Phenazone	$C_{11}H_{12}N_2O$	250mg	
Phenothiazine				
CAS 92-84-2 DRE-C16030000	MW 199.2716 Phenothiazine	$C_{12}H_9NS$	250mg	
Phenoxymethylpenicillin Potassium (Penicillin V potassium salt)				
CAS 132-98-9 DRE-C15935010	MW 388.4799 Penicilline V potassium salt	$C_{16}H_{17}N_2O_5S \cdot K$	100mg	
Phenylacetoneitrile (Benzylcyanide)				
CAS 140-29-4 DRE-E10572200	MW 117.1479 Benzylcyanide	C_8H_7N	250mg	
Phenylbutazone				
CAS 50-33-9 DRE-C16056500	MW 308.3743 Phenylbutazone	$C_{19}H_{20}N_2O_2$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description		
Phenytoin			
DRE-C12895000	5,5-Diphenylhydanthoin		250mg
DRE-C16077000	Phenytoin		250mg
N4-Phthalylsulfathiazole			
CAS 85-73-4	MW 403.4322	C ₁₇ H ₁₃ N ₃ O ₅ S ₂	
DRE-C16190550	N4-Phthalylsulfathiazole		250mg
			
Pindolol			
CAS 13523-86-9	MW 248.3208	C ₁₄ H ₂₀ N ₂ O ₂	
DRE-C16209000	Pindolol		100mg
			
Piperazine			
CAS 110-85-0	MW 86.1356	C ₄ H ₁₀ N ₂	
DRE-C16220400	Piperazine		250mg
			
Piperazine Citrate Hydrate			
CAS 41372-10-5	MW 660.6691	2C ₆ H ₁₀ O ₇ ·3C ₄ H ₁₀ N ₂ ·H ₂ O	
DRE-C16220500	Piperazine citrate hydrate		250mg
			
Piracetam			
CAS 7491-74-9	MW 142.1558	C ₆ H ₁₀ N ₂ O ₂	
DRE-C16245000	Piracetam		100mg
			
Piroxicam			
CAS 36322-90-4	MW 331.3464	C ₁₅ H ₁₃ N ₃ O ₄ S	
DRE-C16278000	Piroxicam		250mg
			
Polymyxin B sulfate			
CAS 1405-20-5	MW n/a		
DRE-C16283000	Polymyxin-B-sulfate		100mg
No Structure			
Ponazuril			
CAS 69004-04-2	MW 457.3805	C ₁₈ H ₁₄ F ₃ N ₃ O ₆ S	
DRE-C16283700	Ponazuril		10mg
			










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Potassium Clavulanate (Clavulanic acid potassium salt)				
CAS 61177-45-5 DRE-C11668545	MW 237.2511 Clavulanic acid potassium salt	$C_8H_8NO_5 \cdot K$	100mg	
Praziquantel				
CAS 55268-74-1 DRE-C16286300	MW 312.4061 Praziquantel	$C_{19}H_{24}N_2O_2$	250mg	
Prednisolone				
CAS 50-24-8 DRE-C16286500	MW 360.444 Prednisolone	$C_{21}H_{28}O_5$	250mg	
Prednisolone Acetate				
CAS 52-21-1 DRE-C16286510	MW 402.4807 Prednisolone acetate	$C_{23}H_{30}O_6$	250mg	
Prednisone				
CAS 53-03-2 DRE-C16286550	MW 358.4281 Prednison(±)	$C_{21}H_{26}O_5$	100mg	
Prednisone Acetate				
CAS 125-10-0 DRE-C16286560	MW 400.4648 Prednisone Acetate	$C_{23}H_{28}O_6$	50mg	
Pregnenolone				
CAS 145-13-1 DRE-C16286700	MW 316.4776 Pregnenolone	$C_{21}H_{32}O_2$	100mg	
Primidone				
CAS 125-33-7 DRE-C16287500	MW 218.2518 Primidone	$C_{12}H_{14}N_2O_2$	100mg	
Procaine Hydrochloride				
CAS 51-05-8 DRE-C16289500	MW 272.771 Procain hydrochloride	$C_{13}H_{20}N_2O_2 \cdot ClH$	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Prochlorperazine maleate				
CAS 84-02-6 DRE-C16291000	MW 606.087 Prochlorperazine maleate (*)	$C_{20}H_{24}ClN_3S \cdot 2C_4H_4O_4$	100mg	
Progesterone				
CAS 57-83-0 DRE-C16342000	MW 314.4617 Progesterone(‡)	$C_{21}H_{30}O_2$	250mg	
Promazine Hydrochloride				
CAS 53-60-1 DRE-C16345700	MW 320.88 Promazine hydrochloride	$C_{17}H_{20}N_2S \cdot ClH$	100mg	
Promethazine Hydrochloride				
CAS 58-33-3 DRE-C16355000	MW 320.88 Promethazine hydrochloride	$C_{17}H_{20}N_2S \cdot ClH$	250mg	
Propionylpromazine hydrochloride				
CAS 7681-67-6 DRE-C16494500	MW 376.9433 Propionylpromazine hydrochloride	$C_{20}H_{24}N_2OS \cdot ClH$	10mg	
Propiophenone				
CAS 93-55-0 DRE-C16487000	MW 134.1751 Propiophenone	$C_9H_{10}O$	250mg	
Propranolol Hydrochloride				
CAS 318-98-9 DRE-C16501000	MW 295.8044 Propranolol hydrochloride	$C_{16}H_{21}NO_2 \cdot ClH$	100mg	
Propyphenazone				
CAS 479-92-5 DRE-C16535000	MW 230.3055 Propyphenazone	$C_{14}H_{18}N_2O$	250mg	
Putrescine Dihydrochloride				
CAS 333-93-7 DRE-C16584000	MW 161.0734 Putrescine dihydrochloride	$C_4H_{12}N_2 \cdot 2ClH$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Pyranocumarin				
CAS 518-20-7 DRE-C16600000	MW 322.3545 Pyranocumarin	C ₂₀ H ₁₈ O ₄	100mg	
Pyrantel Embonate (Pyrantel Pamoate)				
CAS 22204-24-6 DRE-C16600100	MW 594.6768 Pyrantel pamoate	C ₂₃ H ₁₈ O ₆ ·C ₁₁ H ₁₄ N ₂ S	250mg	
Pyrazinamide				
CAS 98-96-4 DRE-C16608000	MW 123.1127 Pyrazinamid	C ₅ H ₅ N ₃ O	100mg	
Pyrimethamine				
CAS 58-14-0 DRE-C16658000	MW 248.7114 Pyrimethamine	C ₁₂ H ₁₃ ClN ₄	250mg	
Pyrrolidone (Pyrrolidin-2-one)				
CAS 616-45-5 DRE-C16676000	MW 85.1045 2-Pyrrolidinone	C ₄ H ₇ NO	100mg	
Quinine				
CAS 130-95-0 DRE-C16706900	MW 324.4168 Quinine	C ₂₀ H ₂₄ N ₂ O ₂	100mg	
Quinine Hydrochloride Dihydrate				
CAS 6119-47-7 DRE-C16707000	MW 396.9083 Quinine monohydrochloride dihydrate	C ₂₀ H ₂₄ N ₂ O ₂ ·ClH·2H ₂ O	100mg	
Quinocetone (technical product)				
CAS 81810-66-4 DRE-C16709000	MW 306.3154 Quinocetone (technical product)	C ₁₈ H ₁₄ N ₂ O ₃	100mg	
2-Quinoxalinecarboxylic Acid				
CAS 879-65-2 DRE-C16713000	MW 174.1561 2-Quinoxalinecarboxylic acid(‡)	C ₉ H ₆ N ₂ O ₂	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Rabeprazole Sodium				
CAS 117976-90-6 DRE-C16804100	MW 381.4245 Rabeprazole sodium	$C_{18}H_{20}N_3O_3S \cdot Na$	100mg	
Ractopamine Hydrochloride				
CAS 90274-24-1 DRE-C16805000	MW 337.8411 Ractopamine hydrochloride(±)	$C_{18}H_{23}NO_3 \cdot ClH$	100mg	
Rafoxanide				
CAS 22662-39-1 DRE-C16805200	MW 626.0105 Rafoxanide	$C_{19}H_{11}Cl_2I_2NO_3$	100mg	
Resorcinol (1,3-Dihydroxybenzene)				
CAS 108-46-3 DRE-C16811200 DRE-L16811200ME DRE-YS09010006MB	MW 110.1106 Resorcinol (Resorcin)(±) Resorcinol (Resorcin) 10 µg/mL in Methanol Resorcinol 100000 µg/mL in Methyl tert-Butyl Ether(±)	$C_6H_6O_2$	250mg 10ml 5x1ml	
Ribavirin				
CAS 36791-04-5 DRE-C16813570	MW 244.2047 Ribavirin(±)	$C_8H_{12}N_4O_5$	50mg	
Ribostamycin Sulfate Salt				
CAS 53797-35-6 DRE-C16813800	MW 552.5511 Ribostamycin sulfate salt	$C_{17}H_{34}N_4O_{10} \cdot H_2O_4S$	100mg	
Rifampicin				
CAS 13292-46-1 DRE-C16814500	MW 822.9402 Rifampicin	$C_{43}H_{58}N_4O_{12}$	100mg	
Rifamycin SV Sodium Salt Hexahydrate				
CAS n/a DRE-C16814600	MW 827.8421 Rifamycin SV sodium salt hexahydrate	$C_{37}H_{46}NO_{12} \cdot Na \cdot 6H_2O$	100mg	
Rifaximin				
CAS 80621-81-4 DRE-C16814700	MW 785.8785 Rifaximin	$C_{43}H_{51}N_5O_{11}$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Robenidine Hydrochloride				
CAS 25875-50-7 DRE-C16815400	MW 370.6642 Robenidine hydrochloride	$C_{15}H_{13}Cl_2N_5 \cdot ClH$	100mg	
Ronidazole				
CAS 7681-76-7 DRE-C16815500	MW 200.1521 Ronidazol(±)	$C_6H_8N_4O_4$	100mg	
Ronidazole D3				
CAS 1015855-87-4 DRE-C16815501	MW 203.1706 Ronidazole D3	$C_6^2H_8H_8N_4O_4$	10mg	
Roxarsone				
CAS 121-19-7 DRE-C16820200	MW 263.0365 Roxarsone	$C_6H_6AsNO_6$	250mg	
Roxithromycin				
CAS 80214-83-1 DRE-C16860000	MW 837.0465 Roxithromycin	$C_{41}H_{76}N_2O_{15}$	100mg	
Rutoside Trihydrate (Rutin Hydrate)				
CAS 250249-75-3 DRE-C16880000	MW 664.5633 Rutin hydrate	$C_{27}H_{30}O_{16} \cdot 3H_2O$	500mg	
Salbutamol				
CAS 18559-94-9 DRE-C16903000 DRE-XA16903000AL	MW 239.3107 Salbutamol free base(±) Salbutamol free base 100 µg/mL in Acetonitrile(±)	$C_{13}H_{21}NO_3$	100mg 1ml	
Salbutamol-D3				
CAS 1219798-60-3 DRE-XA16903001AL	MW 242.3292 Salbutamol D3 (3-hydroxymethyl-D2;a D1) 100 µg/mL in Acetonitrile(±)	$C_{13}^2H_3H_{18}NO_3$	1ml	
Salbutamol Sulfate				
CAS 51022-70-9 DRE-C16903010	MW 576.7 Salbutamol sulfate(±)	$2C_{13}H_{21}NO_3 \cdot H_2O_4S$	100mg	








Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Salicylic Acid				
CAS 69-72-7 DRE-C16903500	MW 138.1207 Salicylic acid(±)	C ₇ H ₆ O ₃	250mg	
Salicylic Acid Benzyl Ester				
CAS 118-58-1 DRE-C16903510	MW 228.2433 Salicylic acid benzyl ester	C ₁₄ H ₁₂ O ₃	250mg	
Salinomycin sodium salt				
CAS 55721-31-8 DRE-C16904500	MW 772.9804 Salinomycin sodium salt	C ₄₂ H ₆₉ O ₁₁ ·Na	100mg	
Salmeterol				
CAS 89365-50-4 DRE-C16904900	MW 415.5656 Salmeterol	C ₂₅ H ₃₇ NO ₄	10mg	
Salmeterol Xinafoate				
CAS 94749-08-3 DRE-C16904920	MW 603.745 Salmeterol xinafoate (*)	C ₂₅ H ₃₇ NO ₄ ·C ₁₁ H ₈ O ₃	250mg	
Sarafloxacin Hydrochloride				
CAS 91296-87-6 DRE-C16908000	MW 421.825 Sarafloxacin hydrochloride(±)	C ₂₀ H ₁₇ F ₂ N ₃ O ₃ ·ClH	100mg	
Sarafloxacin D8 hydrochloride trihydrate				
CAS n/a DRE-C16908001	MW 483.9202 Sarafloxacin D8 hydrochloride trihydrate	C ₂₀ ² H ₈ H ₉ F ₂ N ₃ O ₃ ·ClH·3H ₂ O	10mg	
Secnidazole				
CAS 3366-95-8 DRE-C16931100	MW 185.1805 Secnidazole(±)	C ₇ H ₁₁ N ₃ O ₃	10mg	
Selamectin				
CAS 220119-17-5 DRE-CA16931500	MW 769.9604 Selamectin	C ₄₃ H ₆₃ NO ₁₁	10mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Semicarbazide Hydrochloride				
CAS 563-41-7 DRE-C16933500	MW 111.5308 Semicarbazide hydrochloride(±)	CH ₅ N ₃ O·ClH	100mg	
Sirolimus				
CAS 53123-88-9 DRE-CA16970700	MW 914.1719 Sirolimus	C ₅₁ H ₇₉ NO ₁₃	10mg	
Sotalol Hydrochloride				
CAS 959-24-0 DRE-C16972630	MW 308.8247 Sotalol hydrochloride	C ₁₂ H ₂₀ N ₂ O ₃ S·ClH	100mg	
Sotalol hydrochloride D6 (isopropyl-1,1,1,3,3,3-D6)				
CAS 1246820-85-8 DRE-XA16972631AL	MW 314.8617 Sotalol hydrochloride D6 100 µg/mL in Acetonitrile	C ₁₂ ² H ₆ H ₁₄ N ₂ O ₃ S·ClH	1.1ml	
Sparfloxacin				
CAS 110871-86-8 DRE-C16972650	MW 392.3998 Sparfloxacin	C ₁₉ H ₂₂ F ₂ N ₄ O ₃	100mg	
Spectinomycin Dihydrochloride Pentahydrate				
CAS 22189-32-8 DRE-C16972700	MW 495.3478 Spectinomycin dihydrochloride pentahydrate	C ₁₄ H ₂₄ N ₂ O ₇ ·2ClH·5H ₂ O	250mg	
Spectinomycin Sulfate Tetrahydrate				
CAS 64058-48-6 DRE-C16972720	MW 502.4892 Spectinomycin sulfate hydrate	C ₁₄ H ₂₄ N ₂ O ₇ ·4H ₂ O·H ₂ O ₄ S	100mg	
Spermidine trihydrochloride				
CAS 334-50-9 DRE-C16972738	MW 254.6287 Spermidine trihydrochloride	C ₇ H ₁₉ N ₃ ·3ClH	100mg	
Spermine tetrahydrochloride				
CAS 306-67-2 DRE-C16972742	MW 348.184 Spermine tetrahydrochloride	C ₁₀ H ₂₆ N ₄ ·4ClH	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Spiramycin				
CAS 8025-81-8	MW n/a			
DRE-C16972900	Spiramycin		100mg	No Structure
Spiramycin adipate base				
CAS 68880-55-7	MW n/a			
DRE-C16972905	Spiramycin adipate base		100mg	No Structure
Squalene				
CAS 111-02-4	MW 410.718	C ₃₀ H ₅₀		
DRE-C16973700	Squalene		250mg	
Stanozolol				
CAS 10418-03-8	MW 328.4916	C ₂₁ H ₃₂ N ₂ O		
DRE-C16974000	Stanozolol		100mg	
Streptomycin Sulfate				
CAS 3810-74-0	MW 1457.3836	2C ₂₁ H ₃₉ N ₇ O ₁₂ ·3H ₂ O ₄ S		
DRE-CA16974900	Streptomycin sulfate		250mg	
Succinylsulfathiazole				
CAS 116-43-8	MW 355.3894	C ₁₃ H ₁₃ N ₃ O ₅ S ₂		
DRE-C16985600	Succinylsulfathiazole		250mg	
Sulbactam				
CAS 68373-14-8	MW 233.2416	C ₈ H ₁₁ NO ₅ S		
DRE-C16986600	Sulbactam(‡)		100mg	
Sulfabenzamide				
CAS 127-71-9	MW 276.311	C ₁₃ H ₁₂ N ₂ O ₃ S		
DRE-C16988800	Sulfabenzamide(‡)		250mg	
Sulfacetamide				
CAS 144-80-9	MW 214.2416	C ₈ H ₁₀ N ₂ O ₃ S		
DRE-C16988850	Sulfacetamide(‡)		250mg	
DRE-XA16988850AL	Sulfacetamide 100 µg/mL in Acetonitrile		1ml	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Sulfacetamide Sodium Monohydrate				
CAS 6209-17-2 DRE-C16988853	MW 254.2387 Sulfacetamide sodium salt monohydrate	$C_8H_9N_2O_3S \cdot Na \cdot H_2O$	250mg	
Sulfachloropyrazine Sodium				
CAS 23307-72-4 DRE-C16990045	MW 306.7039 Sulfachloropyrazine sodium	$C_{10}H_6ClN_4O_2S \cdot Na$	100mg	
Sulfachloropyrazine sodium monohydrate				
CAS 73398-14-8 DRE-C16990050	MW 325.7271 Sulfachloropyrazine sodium monohydrate	$C_{10}H_6ClN_4O_2S \cdot HNa \cdot H_2O$	25mg	
Sulfachloropyridazine 13C6 (benzene 13C6)				
CAS n/a DRE-XA16990102AL	MW 290.678 Sulfachloropyridazine 13C6 100 µg/mL in Acetonitrile	$^{13}C_6C_8H_9ClN_4O_2S$	1ml	
Sulfachlorpyridazine				
CAS 80-32-0 DRE-C16990100 DRE-XA16990100AL	MW 284.7221 Sulfachlorpyridazine(‡) Sulfachlorpyridazine 100 µg/mL in Acetonitrile	$C_{10}H_6ClN_4O_2S$	250mg 1ml	
Sulfaclozine Sodium Monohydrate				
CAS 1392129-96-2 DRE-C16990300	MW 324.7192 Sulfaclozine sodium	$C_{10}H_6ClN_4O_2S \cdot Na \cdot H_2O$	100mg	
Sulfadiazine				
CAS 68-35-9 DRE-C16990500	MW 250.277 Sulfadiazin(‡)	$C_{10}H_{10}N_4O_2S$	100mg	
Sulfadimethoxine				
CAS 122-11-2 DRE-C16990550	MW 310.329 Sulfadimethoxine(‡)	$C_{12}H_{14}N_4O_4S$	250mg	
Sulfadimethoxine D6 (2,6-dimethoxy D6)				
CAS 73068-02-7 DRE-C16990551	MW 316.3659 Sulfadimethoxine D6 (2,6-dimethoxy D6)	$C_{12}^2H_6^2H_8N_4O_4S$	10mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Sulfadimidine (Sulfamethazine)				
CAS 57-68-1	MW 278.3302	C ₁₂ H ₁₄ N ₄ O ₂ S		
DRE-C16996500	Sulfamethazine(‡)		250mg	
DRE-XA16996500AL	Sulfamethazine 100 µg/mL in Acetonitrile		1ml	
Sulfadoxine				
CAS 2447-57-6	MW 310.329	C ₁₂ H ₁₄ N ₄ O ₄ S		
DRE-C16990600	Sulfadoxin(‡)		100mg	
Sulfafurazole (Sulfisoxazole)				
CAS 127-69-5	MW 267.3042	C ₁₁ H ₁₃ N ₃ O ₃ S		
DRE-C17000450	Sulfisoxazole(‡)		250mg	
Sulfaguanidine Monohydrate				
CAS 6190-55-2	MW 232.2602	C ₇ H ₁₀ N ₄ O ₂ S·H ₂ O		
DRE-C16990680	Sulfaguanidine monohydrate		100mg	
Sulfamerazine				
CAS 127-79-7	MW 264.3036	C ₁₁ H ₁₂ N ₄ O ₂ S		
DRE-C16995100	Sulfamerazine(‡)		250mg	
Sulfamerazine D4				
CAS 1020719-84-9	MW 268.3282	C ₁₁ ² H ₁₀ H ₈ N ₄ O ₂ S		
DRE-C16995110	Sulfamerazine D4		10mg	
Sulfameter sodium salt				
CAS 18179-67-4	MW 302.2848	C ₁₁ H ₁₁ N ₄ O ₃ S·Na		
DRE-C16995155	Sulfameter sodium salt		250mg	
Sulfamethazine-N4-acetyl				
CAS 100-90-3	MW 320.3668	C ₁₄ H ₁₆ N ₄ O ₃ S		
DRE-XA16996510AL	Sulfamethazine-N4-acetyl 100 µg/mL in Acetonitrile		1ml	
Sulfamethizole				
CAS 144-82-1	MW 270.3313	C ₉ H ₁₀ N ₄ O ₂ S ₂		
DRE-C16998000	Sulfamethizol(‡)		250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Sulfamethoxazole				
CAS 723-46-6 DRE-C16998100	MW 253.2776 Sulfamethoxazole(‡)	C ₁₀ H ₁₁ N ₃ O ₃ S	250mg	
Sulfamethoxazole D4 (benzene D4)				
CAS 1020719-86-1 DRE-C16998110 DRE-XA16998110AL	MW 257.3023 Sulfamethoxazole D4 (benzene D4) Sulfamethoxazole D4 (benzene D4) 100 µg/mL in Acetonitrile	C ₁₀ ² H ₄ H ₇ N ₃ O ₃ S	10mg 1ml	
Sulfamethoxypyridazine				
CAS 80-35-3 DRE-C16998150 DRE-XA16998150AL	MW 280.303 Sulfamethoxypyridazine(‡) Sulfamethoxypyridazine 100 µg/mL in Acetonitrile	C ₁₁ H ₁₂ N ₄ O ₃ S	250mg 1ml	
Sulfametoxydiazine (Sulfameter)				
CAS 651-06-9 DRE-C16995150	MW 280.303 Sulfameter	C ₁₁ H ₁₂ N ₄ O ₃ S	250mg	
Sulfamide				
CAS 7803-58-9 DRE-C16998170	MW 96.109 Sulfamide	H ₄ N ₂ O ₂ S	100mg	
Sulfamonomethoxine				
CAS 1220-83-3 DRE-C16998175	MW 280.303 Sulfamonomethoxine	C ₁₁ H ₁₂ N ₄ O ₃ S	100mg	
Sulfamonomethoxine sodium				
CAS 38006-08-5 DRE-C16998180	MW 302.2848 Sulfamonomethoxine sodium(‡)	C ₁₁ H ₁₁ N ₄ O ₃ S·Na	100mg	
Sulfamoxole				
CAS 729-99-7 DRE-C16998200	MW 267.3042 Sulfamoxol(‡)	C ₁₁ H ₁₃ N ₃ O ₃ S	100mg	
Sulfanilamide				
CAS 63-74-1 DRE-C17000000	MW 172.2049 Sulfanilamide(‡)	C ₆ H ₈ N ₂ O ₂ S	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Sulfanilamide-N-acetyl				
CAS 121-61-9 DRE-C17000010	MW 214.2416 Sulfanilamide-N-acetyl	C ₈ H ₁₀ N ₂ O ₃ S	100mg	
Sulfanitran				
CAS 122-16-7 DRE-C17000050	MW 335.3351 Sulfanitran	C ₁₄ H ₁₃ N ₃ O ₅ S	250mg	
Sulfaphenazole				
CAS 526-08-9 DRE-C17000080	MW 314.3623 Sulfaphenazole(‡)	C ₁₅ H ₁₄ N ₄ O ₂ S	100mg	
Sulfapyridine				
CAS 144-83-2 DRE-C17000100	MW 249.2889 Sulfapyridine(‡)	C ₁₁ H ₁₁ N ₃ O ₂ S	250mg	
Sulfathiazole				
CAS 72-14-0 DRE-C17000200 DRE-XA17000200AL	MW 255.3167 Sulfathiazole(‡) Sulfathiazole 100 µg/mL in Acetonitrile	C ₉ H ₉ N ₃ O ₂ S ₂	250mg 1ml	
Sulfatroxazole				
CAS 23256-23-7 DRE-C17000250	MW 267.3042 Sulfatroxazole	C ₁₁ H ₁₃ N ₃ O ₃ S	50mg	
Sulfisomidine Sodium				
CAS 2462-17-1 DRE-C17000405	MW 300.312 Sulfisomidin sodium	C ₁₂ H ₁₃ N ₄ O ₂ S·Na	250mg	
Sulfisozole sodium				
CAS 37514-39-9 DRE-C17002000	MW 261.2329 Sulfisozole sodium	C ₉ H ₈ N ₃ O ₃ S·Na	100mg	
Tadalafil				
CAS 171596-29-5 DRE-C17133000	MW 389.404 Tadalafil(‡)	C ₂₂ H ₁₉ N ₅ O ₄	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Tamoxifen Citrate				
CAS 54965-24-1 DRE-C17137000	MW 563.6381 Tamoxifen citrate	$C_{26}H_{29}NO \cdot C_6H_8O_7$	250mg	
Tenoxicam				
CAS 59804-37-4 DRE-C17235000	MW 337.3741 Tenoxicam	$C_{13}H_{11}N_3O_4S_2$	10mg	
Terbinafine Hydrochloride				
CAS 78628-80-5 DRE-C17255000	MW 327.8908 Terbinafine hydrochloride	$C_{21}H_{25}N \cdot ClH$	100mg	
Terbutaline Sulfate				
CAS 23031-32-5 DRE-C17295000	MW 548.6468 Terbutalin sulfate(±)	$2C_{12}H_{19}NO_3 \cdot H_2O_4S$	100mg	
Testosterone				
CAS 58-22-0 DRE-C17322500 DRE-XA17322500AL	MW 288.4244 Testosterone Testosterone 100 µg/mL in Acetonitrile	$C_{19}H_{28}O_2$	250mg 1ml	
Testosterone 17-Acetate				
CAS 1045-69-8 DRE-C17322510	MW 330.4611 Testosterone 17-acetate	$C_{21}H_{30}O_3$	100mg	
Testosterone 17-Propionate				
CAS 57-85-2 DRE-C17322540	MW 344.4877 Testosterone 17-propionate	$C_{22}H_{32}O_3$	100mg	
Tetracaine				
CAS 94-24-6 DRE-C17329000	MW 264.3633 Tetracaine	$C_{15}H_{24}N_2O_2$	250mg	
Tetracaine Hydrochloride				
CAS 136-47-0 DRE-C17329007	MW 300.8242 Tetracaine hydrochlorid	$C_{15}H_{24}N_2O_2 \cdot ClH$	100mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Tetracycline Hydrochloride				
CAS 64-75-5 DRE-C17396150	MW 480.8955 Tetracycline hydrochloride(±)	C ₂₂ H ₂₄ N ₂ O ₈ ·ClH	250mg	
Tetramisole Hydrochloride				
CAS 5086-74-8 DRE-C17414500	MW 240.7523 Tetramisole hydrochloride	C ₁₁ H ₁₂ N ₂ S·ClH	100mg	
Tetramisole D5 hydrochloride				
CAS 1173021-85-6 DRE-C17414501	MW 245.7831 Tetramisole D5 hydrochloride	C ₁₁ ² H ₈ H ₇ N ₂ S·ClH	10mg	
Thiamazole (Methimazole)				
CAS 60-56-0 DRE-C15020200	MW 114.1688 Methimazol	C ₄ H ₆ N ₂ S	250mg	
Thiamphenicol				
CAS 15318-45-3 DRE-C17457000	MW 356.2222 Thiamphenicol(±)	C ₁₂ H ₁₅ Cl ₂ NO ₅ S	100mg	
Thiomersal				
CAS 54-64-8 DRE-C17515000	MW 404.8113 Thimerosal	C ₉ H ₉ HgO ₂ S·Na	100mg	
2-Thiouracil				
CAS 141-90-2 DRE-C17561500	MW 128.1524 2-Thiouracil	C ₄ H ₄ N ₂ OS	250mg	
Tiamulin				
CAS 55297-95-5 DRE-C17575790	MW 493.7421 Tiamulin	C ₂₈ H ₄₇ NO ₄ S	100mg	
Tiamulin Hydrogen Fumarate				
CAS 55297-96-6 DRE-C17575800	MW 609.8142 Tiamulin fumarate	C ₂₈ H ₄₇ NO ₄ S·C ₄ H ₄ O ₄	100mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Tiletamine Hydrochloride				
CAS 14176-50-2 DRE-C17581000	MW 259.7954 Tiletamine Hydrochloride	C ₁₂ H ₁₇ NOS·ClH	100mg	
Tilmicosin				
CAS 108050-54-0 DRE-C17582000	MW 869.133 Tilmicosin(±)	C ₄₆ H ₈₀ N ₂ O ₁₃	100mg	
Timolol Maleate				
CAS 26921-17-5 DRE-C17583000	MW 432.4918 Timolol hydrogenmaleate	C ₁₃ H ₂₄ N ₄ O ₃ S·C ₄ H ₄ O ₄	100mg	
Tinidazole				
CAS 19387-91-8 DRE-C17584000	MW 247.2715 Tinidazole	C ₈ H ₁₃ N ₃ O ₄ S	10mg	
Tiopronin				
CAS 1953-02-2 DRE-C17587000	MW 163.1949 Tiopronin	C ₅ H ₉ NO ₃ S	100mg	
Tobramycin				
CAS 32986-56-4 DRE-CA17588000	MW 467.5145 Tobramycin	C ₁₈ H ₃₇ N ₉ O ₉	100mg	
Tolbutamide				
CAS 64-77-7 DRE-C17589900	MW 270.3479 Tolbutamide	C ₁₂ H ₁₈ N ₂ O ₃ S	250mg	
Tolfenamic Acid				
CAS 13710-19-5 DRE-C17591000	MW 261.7036 Tolfenamic acid	C ₁₄ H ₁₂ ClNO ₂	100mg	
Tolnaftate				
CAS 2398-96-1 DRE-C17591700	MW 307.4094 Tolnaftate	C ₁₉ H ₁₇ NOS	250mg	










Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Toltrazuril				
CAS 69004-03-1 DRE-C17592000	MW 425.3817 Toltrazuril	C ₁₈ H ₁₄ F ₃ N ₃ O ₄ S	100mg	
Toltrazuril-sulfoxide				
CAS 69004-15-5 DRE-C17592040	MW 441.3811 Toltrazuril-sulfoxide	C ₁₈ H ₁₄ F ₃ N ₃ O ₅ S	10mg	
Trenbolone				
CAS 10161-33-8 DRE-CA17607700	MW 270.3661 Trenbolone	C ₁₈ H ₂₂ O ₂	100mg	
Trenbolone Acetate				
CAS 10161-34-9 DRE-CA17607710	MW 312.4028 Trenbolone-acetate (*)	C ₂₀ H ₂₄ O ₃	100mg	
Tretinoin				
CAS 302-79-4 DRE-C17608000	MW 300.4351 Tretinoin	C ₂₀ H ₂₈ O ₂	250mg	
Triamcinolone Acetonide				
CAS 76-25-5 DRE-C17635000	MW 434.4977 Triamcinolone acetonide	C ₂₄ H ₃₁ FO ₆	250mg	
Triamcinolone Acetonide 21-Acetate				
CAS 3870-07-3 DRE-C17635100	MW 476.5344 Triamcinolone Acetonide 21-Acetate	C ₂₆ H ₃₃ FO ₇	50mg	
Triamcinolone 16,21-Diacetate				
CAS 67-78-7 DRE-C17636000	MW 478.5072 Triamcinolone Diacetate	C ₂₈ H ₃₁ FO ₈	50mg	
Tricaine Methanesulfonate				
CAS 886-86-2 DRE-C17669050	MW 261.2948 Tricaine Methanesulfonate	C ₉ H ₁₁ NO ₂ ·CH ₄ O ₃ S	100mg	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Triclabendazole				
CAS 68786-66-3 DRE-C17795000	MW 359.6581 Triclabendazole	C ₁₄ H ₉ Cl ₃ N ₂ OS	100mg	
Triclabendazole Sulfone				
CAS 106791-37-1 DRE-C17795010	MW 391.6569 Triclabendazole-sulfone	C ₁₄ H ₉ Cl ₃ N ₂ O ₃ S	50mg	
Triclabendazole Sulfoxide				
CAS 100648-13-3 DRE-C17795020	MW 375.6575 Triclabendazole-sulfoxide	C ₁₄ H ₉ Cl ₃ N ₂ O ₂ S	50mg	
Triflupromazine Hydrochloride				
CAS 1098-60-8 DRE-C17848000	MW 388.878 Triflupromazine hydrochloride	C ₁₈ H ₁₉ F ₃ N ₂ S·ClH	100mg	
Trihexyphenidyl Hydrochloride (Benzhexol Hydrochloride)				
CAS 52-49-3 DRE-C10535800	MW 337.9272 Benzhexol hydrochloride	C ₂₀ H ₃₁ NO·ClH	100mg	
Trimethoprim				
CAS 738-70-5 DRE-C17875000 DRE-XA17875000AL	MW 290.3177 Trimethoprim(‡) Trimethoprim 100 µg/mL in Acetonitrile	C ₁₄ H ₁₈ N ₄ O ₃	250mg 1ml	
Trimethoprim D3 (4-methoxy D3)				
CAS 1189923-38-3 DRE-C17875010 DRE-XA17875010AL	MW 293.3362 Trimethoprim D3 (4-methoxy D3) Trimethoprim D3 (4-methoxy D3) 100 µg/mL in Acetonitrile	C ₁₄ ² H ₁₈ H ₁₉ N ₄ O ₃	10mg 1ml	
Tripelennamine Hydrochloride				
CAS 154-69-8 DRE-C17892500	MW 291.819 Tripelennamine hydrochloride	C ₁₆ H ₂₁ N ₃ ·ClH	100mg	
Troleandomycin (Oleandomycin triacetate)				
CAS 2751-09-9 DRE-C15726200 DRE-XA15726200AL	MW 813.9684 Oleandomycin triacetate Oleandomycin triacetate 100 µg/mL in Acetonitrile	C ₄₁ H ₆₇ NO ₁₅	25mg 1ml	

Pharmaceutical and Veterinary compounds and metabolites

Product code	Description			
Trometamol				
CAS 77-86-1 DRE-C17894900	MW 121.135 Trometamol	$C_4H_{11}NO_3$	100mg	
Tulobuterol Hydrochloride				
CAS 56776-01-3 DRE-C17895400	MW 264.1914 Tulobuterol hydrochloride	$C_{12}H_{18}ClNO \cdot ClH$	50mg	
Tylosin Tartrate				
CAS 74610-55-2 DRE-C17895600	MW 1066.1869 Tylosin tartrate	$C_{46}H_{77}NO_{17} \cdot C_4H_6O_6$	250mg	
Valnemulin hydrochloride				
CAS 133868-46-9 DRE-C17899970	MW 601.2809 Valnemulin hydrochloride	$C_{31}H_{52}N_2O_5S \cdot ClH$	25mg	
Vancomycin Hydrochloride				
CAS 1404-93-9 DRE-C17900500	MW 1485.7145 Vancomycin hydrochloride	$C_{66}H_{72}Cl_2N_9O_{24} \cdot ClH$	100mg	
Vardenafil Hydrochloride				
CAS 224789-15-5 DRE-C17900700	MW 561.5249 Vardenafil hydrochloride(±)	$C_{23}H_{32}N_6O_4S \cdot 2ClH$	100mg	
Xylazine Hydrochloride				
CAS 23076-35-9 DRE-C17943500	MW 256.7948 Xylazine hydrochloride(±)	$C_{12}H_{16}N_2S \cdot ClH$	100mg	
Yohimbine Hydrochloride				
CAS 65-19-0 DRE-C17947100	MW 390.9037 Yohimbine Hydrochloride	$C_{21}H_{26}N_2O_3 \cdot ClH$	250mg	
Zeranol (racemic)				
CAS 55331-29-8 DRE-C17948000	MW 322.396 Zeranol (mixture of isomers)(±)	$C_{18}H_{26}O_5$	10mg	






Pharmaceutical and Veterinary compounds and metabolites

Product code	Description		
DNPH-Mix 1			
DRE-XA18001607AL	DNPH-Mix 1, 100 µg/mL in Acetonitrile (*)	1ml	
Acetone-DNPH	m-Tolualdehyde-DNPH	2-Butanone-DNPH	Crotonaldehyde-DNPH
Cyclohexanone-DNPH	2,5-Dimethylbenzaldehyde-DNPH	Glutaraldehyde-bis(DNPH)	Hexanal-DNPH
Methacrylaldehyde-DNPH	Nonanal-DNPH	Octanal-DNPH	Propionaldehyde-DNPH
o-Tolualdehyde-DNPH	p-Tolualdehyde-DNPH	Acetaldehyde-DNPH	Benzaldehyde-DNPH
Butyraldehyde-DNPH	Decanal-DNPH	Formaldehyde-DNPH	Heptanal-DNPH
Isovaleraldehyde-DNPH	4-Methyl-2-pentanone-DNPH	Pentanal-DNPH	Acrolein-DNPH






VOLATILE ORGANIC COMPOUNDS (VOCs)







Volatile organic compounds (VOCs)

Product code	Description			
Allylchloride (3-Chloro-1-propene)				
CAS 107-05-1	MW 76.5248	C ₃ H ₅ Cl		
DRE-CA10135000	Allylchloride (3-Chloro-1-propene)		250mg	
Bis(chloromethyl) Ether				
CAS 542-88-1	MW 114.9586	C ₂ H ₄ Cl ₂ O		
DRE-GA09010398MB	Bis(chloromethyl)ether Solution, 1000 µg/mL in MtBE(‡)		1.5ml	
DRE-GA09010398MB-SS	Bis(chloromethyl)ether Solution, Second Source, 2000 µg/mL in MtBE(‡)		1ml	
Bromochloroacetonitrile				
CAS 83463-62-1	MW 154.393	C ₂ HBrClN		
DRE-CA10715000	Bromochloroacetonitrile		100mg	
Bromochlorodifluoromethane				
CAS 353-59-3	MW 165.3645	CBrClF ₂		
DRE-XA10720500ME	Bromochlorodifluoromethane 100 µg/mL in Methanol		1ml	
1-Bromo-2-chloroethane				
CAS 107-04-0	MW 143.4102	C ₂ H ₄ BrCl		
DRE-CA10720700	1-Bromo-2-chloroethane		250mg	
DRE-GA09010038ME	1-Bromo-2-Chloroethane 1000 µg/mL in p/t Methanol(‡)		1ml	
DRE-GS09010038ME	1-Bromo-2-Chloroethane 1000 µg/mL in p/t Methanol(‡)		5x1ml	
Bromochloromethane				
CAS 74-97-5	MW 129.3836	CH ₂ BrCl		
DRE-CA10720800	Bromochloromethane		1ml	
DRE-GA10720800ME	Bromochloromethane		1ml	
DRE-GH10720800ME	Bromochloromethane 1000 mµg/mL in Methanol		5x1ml	
DRE-XA10720800ME	Bromochloromethane 100 µg/mL in Methanol		1ml	
1-Bromo-3-chloropropane				
CAS 109-70-6	MW 157.4367	C ₃ H ₆ BrCl		
DRE-CA10722000	1-Bromo-3-chloropropane		500mg	
2-Bromo-1-chloropropane				
CAS 3017-95-6	MW 157.4367	C ₃ H ₆ BrCl		
DRE-YA10722200ME	2-Bromo-1-chloropropane 2000 µg/mL in Methanol		1ml	
Bromodichloromethane				
CAS 75-27-4	MW 163.8286	CHBrCl ₂		
DRE-C10726700	Bromodichloromethane		1ml	
DRE-XA10726700ME	Bromodichloromethane 100 µg/mL in Methanol (*)		1ml	






Volatile organic compounds (VOCs)

Product code	Description			
Bromoethane				
CAS 74-96-4 DRE-CA10728000	MW 108.9651 Bromoethane	C ₂ H ₅ Br	500mg	
2-Bromoethanol				
CAS 540-51-2 DRE-C10728500	MW 124.9645 2-Bromoethanol	C ₂ H ₅ BrO	500mg	
1-Bromo-2-fluorobenzene				
CAS 1072-85-1 DRE-CA10730500	MW 174.9984 1-Bromo-2-fluorobenzene	C ₆ H ₄ BrF	.5ml	
Bromoform (Tribromomethane)				
CAS 75-25-2 DRE-C17665500 DRE-XA17665500ME	MW 252.7306 Tribromomethane Tribromomethane (Bromoform) 100 µg/mL in Methanol (*)	CHBr ₃	1g 1ml	
1-Bromopropane				
CAS 106-94-5 DRE-C10759500	MW 122.9917 1-Bromopropane(±)	C ₃ H ₇ Br	250mg	
2-Bromopropane				
CAS 75-26-3 DRE-CA10759600	MW 122.9917 2-Bromopropane	C ₃ H ₇ Br	250mg	
1-Bromopropan-2-ol				
CAS 19686-73-8 DRE-C10760000	MW 138.9911 1-Bromopropanol-2	C ₃ H ₇ BrO	250mg	
Bromotrichloromethane				
CAS 75-62-7 DRE-CA10765000	MW 198.2737 Bromotrichloromethane	CBrl ₃	500mg	
Carbontetrachloride (Tetrachloromethane)				
CAS 56-23-5 DRE-CA17359500 DRE-XA17359500ME	MW 153.8227 Tetrachloromethane Tetrachloromethane 100 µg/mL in Methanol	CCl ₄	1ml 1ml	






Volatile organic compounds (VOCs)

Product code	Description			
Chloroacetaldehyde				
CAS 107-20-0 DRE-C11347000	MW 78.4976 Chloroacetaldehyde	C ₂ H ₃ ClO	250mg	
1-Chlorobutane				
CAS 109-69-3 DRE-CA11395000	MW 92.5673 1-Chlorobutane	C ₄ H ₉ Cl	1ml	
Chlorobutanol Hemihydrate (Acetone chloroform)				
CAS 6001-64-5 DRE-C10020000	MW 372.9288 Acetonchloroform	2C ₄ H ₇ Cl ₃ O·H ₂ O	1g	
1-Chloro-1,1-difluoroethane				
CAS 75-68-3 DRE-XA11404000ME DRE-YS09010030ME	MW 100.495 1-Chloro-1,1-difluoroethane 100 µg/mL in Methanol 1-Chloro-1,1-difluoroethane 1000 µg/mL in p/t Methanol(‡)	C ₂ H ₃ ClF ₂	1ml 5x1ml	
Chlorodifluoromethane (Freon 22)				
CAS 75-45-6 DRE-XA11404400ME DRE-YS09010028ME	MW 86.4684 Chlorodifluoromethane 100 µg/mL in Methanol Chlorodifluoromethane (Freon 22) 1000 µg/mL in p/t Methanol(‡)	CHClF ₂	1ml 5x1ml	
Chloroethane				
CAS 75-00-3 DRE-XA11409000ME	MW 64.5141 Chloroethane 100 µg/mL in Methanol	C ₂ H ₅ Cl	1ml	
2-Chloroethyl-vinylether				
CAS 110-75-8 DRE-C11410400 DRE-GA09010380ME	MW 106.5508 2-Chloroethyl-vinylether 2-Chloroethylvinyl ether, 2000 µg/mL in p/t Methanol(‡)	C ₄ H ₇ ClO	100mg 1ml	
Chloroform				
CAS 67-66-3 DRE-C17739500 DRE-L17739500ME DRE-XA17739500ME DRE-YA17739500ME	MW 119.3776 Trichloromethane (Chloroform)(‡) Trichloromethane (Chloroform) 10 µg/mL in Methanol Trichloromethane (Chloroform) 100 µg/mL in Methanol Trichloromethane (Chloroform) 5000 µg/mL in Methanol	CHCl ₃	5ml 10ml 1ml 1ml	
1-Chlorohexane (Hexylchloride)				
CAS 544-10-5 DRE-CA11416000	MW 120.6204 1-Chlorohexane (Hexylchloride)	C ₆ H ₁₃ Cl	1ml	







Volatile organic compounds (VOCs)

Product code	Description			
Chloromethane (Methylchloride)				
CAS 74-87-3	MW 50.4875	CH ₃ Cl		
DRE-XA11419000ME	Chloromethane (Methylchloride) 100 µg/mL in Methanol (*)		1ml	
Chloromethyl Methyl Ether				
CAS 107-30-2	MW 80.5135	C ₂ H ₅ ClO		
DRE-CA11435000	Chloromethyl-methylether		100mg	
Chloropentafluoroethane				
CAS 76-15-3	MW 154.4664	C ₂ ClF ₅		
DRE-XA11458000ME	Chloropentafluoroethane 100 µg/mL in Methanol		1ml	
1-Chloropropane (Propylchloride)				
CAS 540-54-5	MW 78.5407	C ₃ H ₇ Cl		
DRE-CA11502500	1-Chloropropane (Propylchloride)		1ml	
2-Chloro-1,1,1-trifluoroethane				
CAS 75-88-7	MW 118.4855	C ₂ H ₂ ClF ₃		
DRE-XA11534000ME	2-Chloro-1,1,1-trifluoroethane 100 µg/mL in Methanol		1ml	
Dibromoacetonitrile				
CAS 3252-43-5	MW 198.844	C ₂ HBr ₂ N		
DRE-C12216500	Dibromoacetonitrile		250mg	
Dibromochloromethane				
CAS 124-48-1	MW 208.2796	CHBr ₂ Cl		
DRE-CA12234700	Dibromochloromethane		.5ml	
DRE-XA12234700ME	Dibromochloromethane 100 µg/mL in Methanol (*)		1ml	
1,2-Dibromo-3-chloropropane				
CAS 96-12-8	MW 236.3328	C ₃ H ₅ Br ₂ Cl		
DRE-CA12235000	1,2-Dibromo-3-chloropropane(‡)		250mg	
DRE-L12235000ME	1,2-Dibromo-3-chloropropane 10 µg/mL in Methanol		10ml	
DRE-YA12235000ME	1,2-Dibromo-3-chloropropane 2000 µg/mL in Methanol		1ml	
1,2-Dibromo-1,1-dichloroethane				
CAS 75-81-0	MW 256.7513	C ₂ H ₂ Br ₂ Cl ₂		
DRE-XA12236500ME	1,2-Dibromo-1,1-dichloroethane 100 µg/mL in Methanol		1ml	







Volatile organic compounds (VOCs)

Product code	Description			
Dibromodichloromethane				
CAS 594-18-3 DRE-C12237000	MW 242.7247 Dibromodichloromethane	CBr ₂ Cl ₂	100mg	
1,2-Dibromoethane				
CAS 106-93-4 DRE-XA12240000ME DRE-YA12240000ME	MW 187.8612 1,2-Dibromoethane 100 µg/mL in Methanol 1,2-Dibromoethane 1000 µg/mL in Methanol	C ₂ H ₄ Br ₂	1ml 1ml	
1,2-Dibromoethene				
CAS 540-49-8 DRE-CA12240200	MW 185.8453 1,2-Dibromoethene	C ₂ H ₂ Br ₂	250mg	
Dibromomethane				
CAS 74-95-3 DRE-C12240500	MW 173.8346 Dibromomethane	CH ₂ Br ₂	250mg	
Dichloroacetonitrile				
CAS 3018-12-0 DRE-C12321000 DRE-XA12321000CY	MW 109.942 Dichloroacetonitrile Dichloroacetonitrile 100 µg/mL in Cyclohexane	C ₂ HCl ₂ N	1ml 1ml	
1,2-Dichlorobenzene				
CAS 95-50-1 DRE-C12370000	MW 147.002 1,2-Dichlorobenzene(‡)	C ₆ H ₄ Cl ₂	1g	
1,3-Dichlorobenzene				
CAS 541-73-1 DRE-C12371000	MW 147.002 1,3-Dichlorobenzene	C ₆ H ₄ Cl ₂	1g	
1,4-Dichlorobenzene				
CAS 106-46-7 DRE-C12372000	MW 147.002 1,4-Dichlorobenzene	C ₆ H ₄ Cl ₂	1g	
1,4-Dichlorobenzene D4				
CAS 3855-82-1 DRE-C12372100	MW 151.0266 1,4-Dichlorobenzene D4	C ₆ ² H ₄ Cl ₂	100mg	









Volatile organic compounds (VOCs)

Product code	Description			
1,4-Dichlorobutane				
CAS 110-56-5 DRE-C12420500	MW 127.0123 1,4-Dichlorobutane	C ₄ H ₈ Cl ₂	250mg	
cis-1,4-Dichloro-2-butene				
CAS 1476-11-5 DRE-C12420600	MW 124.9964 cis-1,4-Dichloro-2-butene	C ₄ H ₆ Cl ₂	250mg	
trans-1,4-Dichloro-2-butene				
CAS 110-57-6 DRE-C12420700	MW 124.9964 trans-1,4-Dichloro-2-butene	C ₄ H ₆ Cl ₂	250mg	
1,1-Dichloroethane				
CAS 75-34-3 DRE-C12422000 DRE-XA12422000ME	MW 98.9592 1,1-Dichloroethane 1,1-Dichloroethane 100 µg/mL in Methanol	C ₂ H ₄ Cl ₂	.5ml 1ml	
1,2-Dichloroethane				
CAS 107-06-2 DRE-C12422200 DRE-L12422200ME DRE-XA12422200ME DRE-YA12422200ME	MW 98.9592 1,2-Dichloroethane 1,2-Dichloroethane 10 µg/mL in Methanol 1,2-Dichloroethane 100 µg/mL in Methanol 1,2-Dichloroethane 1000 µg/mL in Methanol	C ₂ H ₄ Cl ₂	1ml 10ml 1ml 1ml	
1,2-Dichloroethane D4				
CAS 17060-07-0 DRE-C12422300 DRE-YA12422300ME	MW 102.9838 1,2-Dichloroethane D4 1,2-Dichloroethane D4 2000 µg/mL in Methanol	C ₂ H ₄ Cl ₂	100mg 1ml	
2,2-Dichloroethanol				
CAS 598-38-9 DRE-C12422350	MW 114.9586 2,2-Dichloroethanol	C ₂ H ₄ Cl ₂ O	100mg	
1,1-Dichloroethene				
CAS 75-35-4 DRE-C12422400 DRE-L12422400ME DRE-XA12422400ME	MW 96.9433 1,1-Dichloroethene 1,1-Dichloroethene 10 µg/mL in Methanol 1,1-Dichloroethene 100 µg/mL in Methanol	C ₂ H ₂ Cl ₂	1ml 10ml 1ml	
cis-1,2-Dichloroethene				
CAS 156-59-2 DRE-CA12422500 DRE-L12422500ME DRE-XA12422500ME	MW 96.9433 cis-1,2-Dichloroethene cis-1,2-Dichloroethene 10 µg/mL in Methanol cis-1,2-Dichloroethene 100 µg/mL in Methanol	C ₂ H ₂ Cl ₂	100mg 10ml 1ml	






Volatile organic compounds (VOCs)

Product code	Description			
trans-1,2-Dichloroethene				
CAS 156-60-5	MW 96.9433	C ₂ H ₂ Cl ₂		
DRE-C12422600	trans-1,2-Dichloroethene		100mg	
DRE-L12422600ME	trans-1,2-Dichloroethene 10 µg/mL in Methanol		10ml	
DRE-XA12422600ME	trans-1,2-Dichloroethene 100 µg/mL in Methanol		1ml	
1,1-Dichloro-1-fluoroethane				
CAS 1717-00-6	MW 116.9496	C ₂ H ₃ Cl ₂ F		
DRE-XA12422800ME	1,1-Dichloro-1-fluoroethane 100 µg/mL in Methanol		1ml	
Dichlorofluoromethane				
CAS 75-43-4	MW 102.923	CHCl ₂ F		
DRE-XA12423100ME	Dichlorofluoromethane 100 µg/mL in Methanol		1ml	
1,1-Dichloropropane				
CAS 78-99-9	MW 112.9857	C ₃ H ₆ Cl ₂		
DRE-XA12479900ME	1,1-Dichloropropane 100 µg/mL in Methanol		1ml	
1,2-Dichloropropane				
CAS 78-87-5	MW 112.9857	C ₃ H ₆ Cl ₂		
DRE-CA12480000	1,2-Dichloropropane		1ml	
DRE-GA12480000ME	1,2-Dichloropropane 10000 µg/mL in Methanol		1ml	
DRE-GH12480000ME	1,2-Dichloropropane 10000 µg/mL in Methanol		5x1ml	
DRE-XA12480000CY	1,2-Dichloropropane 100 µg/mL in Cyclohexane		1ml	
1,3-Dichloropropane				
CAS 142-28-9	MW 112.9857	C ₃ H ₆ Cl ₂		
DRE-CA12481000	1,3-Dichloropropane		250mg	
2,2-Dichloropropane				
CAS 594-20-7	MW 112.9857	C ₃ H ₆ Cl ₂		
DRE-CA12481200	2,2-Dichloropropane		250mg	
DRE-XA12481200ME	2,2-Dichloropropane 100 µg/mL in Methanol		1ml	
1,3-Dichloropropan-2-ol				
CAS 96-23-1	MW 128.9851	C ₃ H ₆ Cl ₂ O		
DRE-C12481600	1,3-Dichloropropan-2-ol		250mg	
2,3-Dichloro-1-propanol				
CAS 616-23-9	MW 128.9851	C ₃ H ₆ Cl ₂ O		
DRE-C12482000	2,3-Dichloro-1-propanol		.5ml	





Volatile organic compounds (VOCs)

Product code	Description			
1,1-Dichloropropan-2-one				
CAS 513-88-2 DRE-C12482500	MW 126.9693 1,1-Dichloropropanone-2	C ₃ H ₄ Cl ₂ O	1ml	
1,1-Dichloro-1-propene				
CAS 563-58-6 DRE-XA12489500ME	MW 110.9699 1,1-Dichloro-1-propene 100 µg/mL in Methanol	C ₃ H ₄ Cl ₂	1ml	
1,3-Dichloropropene				
CAS 542-75-6 DRE-CA12490000	MW 110.9699 1,3-Dichloropropene (cis + trans)	C ₃ H ₄ Cl ₂	250mg	
2,3-Dichloro-1-propene				
CAS 78-88-6 DRE-CA12490500	MW 110.9699 2,3-Dichloro-1-propene	C ₃ H ₄ Cl ₂	250mg	
cis-1,3-Dichloropropene				
CAS 10061-01-5 DRE-CA12489800 DRE-L12489800ME DRE-XA12489800ME	MW 110.9699 cis-1,3-Dichloropropene cis-1,3-Dichloropropene (technical) 10 µg/mL in Methanol cis-1,3-Dichloropropene (technical) 100 µg/mL in Methanol	C ₃ H ₄ Cl ₂	100mg 10ml 1ml	
trans-1,3-Dichloropropene				
CAS 10061-02-6 DRE-C12489900	MW 110.9699 trans-1,3-Dichloropropene(‡)	C ₃ H ₄ Cl ₂	250mg	
1,2-Dichlorotetrafluoroethane (CFC-114)				
CAS 76-14-2 DRE-XA12504000ME DRE-GA09010382ME DRE-GS09010382ME	MW 170.921 1,2-Dichlorotetrafluoroethane 100 µg/mL in Methanol 1,2-Dichlorotetrafluoroethane (CFC-114) Solution, 2000 µg/mL in p/t Methanol (‡) 1,2-Dichlorotetrafluoroethane (CFC-114) Solution, 2000 µg/mL in p/t Methanol (‡)	C ₂ Cl ₂ F ₄	1ml 1ml 5x1ml	
1,2-Dichlorotrifluoroethane				
CAS 354-23-4 DRE-GS09010147ME	MW 152.9305 1,2-Dichlorotrifluoroethane 100 µg/mL in p/t Methanol(‡)	C ₂ HCl ₂ F ₃	5x1ml	








Volatile organic compounds (VOCs)

Product code	Description			
Fluorotrichloromethane (Trichlorofluoromethane)				
CAS 75-69-4	MW 137.3681	CCl ₃ F		
DRE-C13798500	Fluorotrichloromethane (Trichlorofluoromethane)		250mg	
DRE-XA13798500ME	Fluorotrichloromethane (Trichlorofluoromethane) 100 µg/mL in Methanol		1ml	
DRE-YA13798500ME	Fluorotrichloromethane (Trichlorofluoromethane), 5000 µg/mL in Methanol		1ml	
Hexachlorobutadiene				
CAS 87-68-3	MW 260.7608	C ₄ Cl ₆		
DRE-C14170000	Hexachlorobutadiene		250mg	
DRE-L14170000CY	Hexachlorobutadiene 10 µg/mL in Cyclohexane		10ml	
DRE-L14170000ME	Hexachlorobutadiene 10 µg/mL in Methanol		10ml	
DRE-XA14170000CY	Hexachlorobutadiene 100 µg/mL in Cyclohexane		1ml	
DRE-XA14170000ME	Hexachlorobutadiene 100 µg/mL in Methanol		1ml	
Hexachloro-1,3-butadiene 13C4				
CAS 93951-70-3	MW 264.7314	¹³ C ₄ Cl ₆		
DRE-XA14170100AC	Hexachloro-1,3-butadiene 13C4 100 µg/mL in Acetone(‡)		1.1ml	
Hexachloroethane				
CAS 67-72-1	MW 236.7394	C ₂ Cl ₆		
DRE-CA14172000	Hexachloroethane		250mg	
DRE-L14172000IO	Hexachloroethane 10 µg/mL in Isooctane		10ml	
DRE-XA14172000ME	Hexachloroethane 100 µg/mL in Methanol		1ml	
Hexachloropropene (Perchloropropene)				
CAS 1888-71-7	MW 248.7501	C ₃ Cl ₆		
DRE-C14183000	Hexachloropropene (Perchloropropene)		250mg	
Iodomethane				
CAS 74-88-4	MW 141.939	CH ₃ I		
DRE-C14343000	Iodomethane		250mg	
Methylene chloride (Dichloromethane)				
CAS 75-09-2	MW 84.9326	CH ₂ Cl ₂		
DRE-C12424500	Dichloromethane (Methylenechloride)		1ml	
DRE-XA12424500ME	Dichloromethane (Methylenechloride) 100 µg/mL in Methanol		1ml	
DRE-YA12424500ME	Dichloromethane (Methylenechloride) 1000 µg/mL in Methanol		1ml	
Norflurane (1,1,1,2-Tetrafluoroethane)				
CAS 811-97-2	MW 102.0309	C ₂ H ₂ F ₄		
DRE-XA17404000ME	1,1,1,2-Tetrafluoroethane 100 µg/mL in Methanol		1ml	
Pentachloroethane				
CAS 76-01-7	MW 202.2943	C ₂ HCl ₅		
DRE-CA15965000	Pentachloroethane		250mg	
DRE-GA09010385ME	Pentachloroethane Solution, 2000 µg/mL in p/t Methanol(‡)		1ml	

Volatile organic compounds (VOCs)

Product code	Description			
1,1,2,2-Tetrabromoethane				
CAS 79-27-6	MW 345.6533	C ₂ H ₂ Br ₄		
DRE-CA17325000	1,1,2,2-Tetrabromoethane		250mg	
Tetrabromomethane				
CAS 558-13-4	MW 331.6267	CBr ₄		
DRE-C17326000	Tetrabromomethane		250mg	
1,1,1,2-Tetrachloroethane				
CAS 630-20-6	MW 167.8493	C ₂ H ₂ Cl ₄		
DRE-CA17358000	1,1,1,2-Tetrachloroethane		1ml	
DRE-XA17358000ME	1,1,1,2-Tetrachloroethane 100 µg/mL in Methanol		1ml	
1,1,2,2-Tetrachloroethane				
CAS 79-34-5	MW 167.8493	C ₂ H ₂ Cl ₄		
DRE-CA17358100	1,1,2,2-Tetrachloroethane		1ml	
DRE-L17358100ME	1,1,2,2-Tetrachloroethane 10 µg/mL in Methanol		10ml	
DRE-XA17358100ME	1,1,2,2-Tetrachloroethane 100 µg/mL in Methanol		1ml	
Tetrachloroethene				
CAS 127-18-4	MW 165.8334	C ₂ Cl ₄		
DRE-C17358300	Tetrachloroethene		5ml	
DRE-XA17358300ME	Tetrachloroethene 100 µg/mL in Methanol (*)		1ml	
DRE-Y17358300ME	Tetrachloroethene 1000 µg/mL in Methanol		10ml	
Trichloroacetonitrile				
CAS 545-06-2	MW 144.3871	C ₂ Cl ₃ N		
DRE-CA17688000	Trichloroacetonitrile		250mg	
1,1,1-Trichloroethane				
CAS 71-55-6	MW 133.4042	C ₂ H ₃ Cl ₃		
DRE-CA17738300	1,1,1-Trichloroethane(‡)		.5ml	
DRE-L17738300ME	1,1,1-Trichloroethane 10 µg/mL in Methanol		10ml	
DRE-XA17738300ME	1,1,1-Trichloroethane 100 µg/mL in Methanol		1ml	
1,1,2-Trichloroethane				
CAS 79-00-5	MW 133.4042	C ₂ H ₃ Cl ₃		
DRE-C17738500	1,1,2-Trichloroethane(‡)		1ml	
DRE-L17738500ME	1,1,2-Trichloroethane 10 µg/mL in Methanol		10ml	
DRE-XA17738500ME	1,1,2-Trichloroethane 100 µg/mL in Methanol		1ml	
2,2,2-Trichloroethanol				
CAS 115-20-8	MW 149.4036	C ₂ H ₃ Cl ₃ O		
DRE-CA17739000	2,2,2-Trichloroethanol		250mg	

Volatile organic compounds (VOCs)

Product code	Description			
Trichloroethene				
CAS 79-01-6	MW 131.3883	C ₂ HCl ₃		
DRE-CA17739300	Trichloroethene		1ml	
DRE-L17739300ME	Trichloroethene 10 µg/mL in Methanol		10ml	
DRE-XA17739300ME	Trichloroethene 100 µg/mL in Methanol		1ml	
DRE-Y17739300ME	Trichloroethene 1000 µg/mL in Methanol		10ml	
DRE-YA17739300ME	Trichloroethene 1000 µg/mL in Methanol		1ml	
1,2,3-Trichloropropane				
CAS 96-18-4	MW 147.4308	C ₃ H ₅ Cl ₃		
DRE-CA17780000	1,2,3-Trichloropropane		1ml	
DRE-L17780000ME	1,2,3-Trichloropropane 10 µg/mL in Methanol		10ml	
DRE-XA17780000ME	1,2,3-Trichloropropane 100 µg/mL in Methanol		1ml	
DRE-GA09010312ME	1,2,3-Trichloropropane 1000 µg/mL in p/t Methanol(±)		1ml	
1,1,1-Trichloro-2,2,2-trifluoroethane				
CAS 354-58-5	MW 187.3756	C ₂ Cl ₃ F ₃		
DRE-L17788200ME	1,1,1-Trichlorotrifluoroethane 10 µg/mL in Methanol		10ml	
1,1,2-Trichloro-1,2,2-trifluoroethane				
CAS 76-13-1	MW 187.3756	C ₂ Cl ₃ F ₃		
DRE-CA17788300	1,1,2-Trichlorotrifluoroethane		250mg	
DRE-L17788300ME	1,1,2-Trichlorotrifluoroethane 10 µg/mL in Methanol		10ml	
DRE-XA17788300ME	1,1,2-Trichlorotrifluoroethane 100 µg/mL in Methanol		1ml	
DRE-YS09010021ME	1,1,2-trichloro-1,2,2-trifluoroethane (Freon 113) 2000 µg/mL in p/t Methanol(±)		5x1ml	
1,1,1-Trifluoroethane				
CAS 420-46-2	MW 84.0404	C ₂ H ₃ F ₃		
DRE-GS09010082ME	1,1,1-Trifluoroethane 100 µg/mL in p/t Methanol(±)		5x1ml	
2,2,2-Trifluoroethanol				
CAS 75-89-8	MW 100.0398	C ₂ H ₃ F ₃ O		
DRE-C17844600	2,2,2-Trifluoroethanol		250mg	
Vinyl Chloride				
CAS 75-01-4	MW 62.4982	C ₂ H ₃ Cl		
DRE-Y17923000ME	Vinyl chloride 1000 µg/mL in Methanol (*)		10ml	
Alcohols Mix, 2-1284				
DRE-YS09000033ME	VOC Alcohol Mix 40000 µg/mL in p/t Methanol(±)			5x1ml
	ethanol			
	isopropyl alcohol			

Volatile organic compounds (VOCs)

Product code	Description	
Aromatic VOC Mix 1		
DRE-YA08020100ME	Aromatic VOC Mix 1 2000 µg/mL in Methanol	1ml
	1,2-Dichlorobenzene 1,4-Dichlorobenzene Chlorobenzene m-Xylene p-Xylene	1,3-Dichlorobenzene Benzene Ethylbenzene o-Xylene Toluene
Aromatic VOC Mix 3		
DRE-YA08020300ME	Aromatic VOC Mix 3 2000 µg/mL in Methanol (*)	1ml
	1,2-Dichlorobenzene 1,4-Dichlorobenzene Chlorobenzene m-Xylene p-Xylene Toluene	1,3-Dichlorobenzene Benzene Ethylbenzene o-Xylene Styrene
Chlorinated VOC Mix		
DRE-GS09000175HP	Chlorinated VOC Mix 5 µg/mL in Heptane(‡)	5x1ml
DRE-GS09000176HP	Chlorinated VOC Mix 1000 µg/mL in Heptane(‡)	5x1ml
	carbon tetrachloride tetrachloroethylene 1,1,1-trichloroethane trichloroethylene trichlorofluoromethane 1,1,2-trichloro-1,2,2-trifluoroethane (Freon 113)	
Chlorinated VOC Mix, 6-3621		
DRE-YS09000034HP	Chlorinated VOC Mix 5 µg/mL in Heptane(‡)	5x1ml
	carbon tetrachloride tetrachloroethylene 1,1,1-trichloroethane 1,1,2-trichloroethane trichlorofluoromethane 1,1,2-trichloro-1,2,2-trifluoroethane (Freon 113)	
PVOC Mixture 3 (Wisconsin)		
DRE-YA03032300ME	PVOC Mixture 3 (Wisconsin) 1000 µg/mL in Methanol	1ml
	1,2,4-Trimethylbenzene Benzene Methyl-tert-butylether Naphthalene p-Xylene	1,3,5-Trimethylbenzene Ethylbenzene m-Xylene o-Xylene Toluene
Purgeable Halocarbons Mix 2		
DRE-XA06010200ME	Purgeable Halocarbons Mix 2 200 µg/mL in Methanol (*)	1ml
	Bromodichloromethane Dibromochloromethane Tribromomethane Trichloromethane	
Purgeable Halocarbons Mix 5		
DRE-YA06010500ME	Purgeable Halocarbons Mix 5 2000 µg/mL in Methanol (*)	1ml
1,1,1-Trichloroethane 1,1-Dichloroethene 1,3-Dichlorobenzene cis-1,3-Dichloropropene Tetrachloromethane Trichloroethene	1,1,2,2-Tetrachloroethane 1,2-Dichlorobenzene 1,4-Dichlorobenzene Dibromochloromethane trans-1,2-Dichloroethene Trichloromethane	1,1,2-Trichloroethane 1,2-Dichloroethane Bromodichloromethane Dichloromethane trans-1,3-Dichloropropene
		1,1-Dichloroethane 1,2-Dichloropropane Chlorobenzene Tetrachloroethene Tribromomethane

Volatile organic compounds (VOCs)

Product code	Description			
Semi Volatile Mix				
DRE-XA09000164DI	Semi Volatile Mix 1000 µg/mL in Dichloromethane			1ml
	n-nitrosodiethylamine N-nitrosopyrrolidine 1,2,4,5-tetrachlorobenzene parathion benzaldehyde	n-nitrosodi-n-butylamine pentachlorobenzene caprolactam benzoic acid		
Semi-Volatile Mix 2				
DRE-YS09000037DI	Semi-Volatile Mix 2 in Dichloromethane(‡)			5x1.5ml
	benzoic acid [2000 mg/L] 2,6-dichlorophenol [2000 mg/L] (L)-a-terpineol [500 mg/L] 2,3,4,6-Tetrachlorophenol [2000 mg/L]			
Semi Volatile Mix B				
DRE-GS09000166DI	Semi-Volatile Mix B 1000 µg/mL in Dichloromethane(‡)			5x1ml
	2,6-dichlorophenol benzoic acid 3-methylcholanthrene 1,4-dioxane			
Semi Volatile Mix C				
DRE-GS09000197AC	Semi Volatile Mix C 100 µg/mL in Acetone(‡)			5x1ml
	benzoic acid famphur n-octadecane (C18) thionazine (zinophos)	hexachlorocyclopentadiene kepone tetraethyl dithiopyrophosphate phorate	benzaldehyde methyl parathion 1,4-dioxane disulfoton	dimethoate decane (C10) O,O,O-triethylphosphorothioate parathion
Semi Volatile Mix D				
DRE-GH09000198DI	Semi Volatile Mix D in Dichloromethane(‡)			10x1ml
	1,2,3-trimethylbenzene [200 mg/L] benzo[e]pyrene [200 mg/L] dibenzothiophene [200 mg/L] perylene [200 mg/L] acenaphthylene [200 mg/L] benzo[k]fluoranthene [200 mg/L] dibenz[a,h]anthracene [200 mg/L] naphthalene [200 mg/L] carbazole [200 mg/L] ethylbenzene [200 mg/L] p-xylene [100 mg/L] isopropylbenzene [200 mg/L] methyl t-butyl ether [200 mg/L] heptadecane (C17) [200 mg/L] retene [200 mg/L]	1,2,4-trimethylbenzene [200 mg/L] biphenyl [200 mg/L] indene [200 mg/L] phenol [200 mg/L] anthracene [200 mg/L] benzo[ghi]perylene [200 mg/L] fluoranthene [200 mg/L] phenanthrene [200 mg/L] heptane (C7) [200 mg/L] toluene [200 mg/L] n-butylbenzene [200 mg/L] 4-isopropyltoluene [200 mg/L] isooctane [200 mg/L] pristane [200 mg/L] 1,2-benzodiphenylene sulfide [200 mg/L]	1,3,5-trimethylbenzene [200 mg/L] cis-decalin [200 mg/L] 1-benzothiophene [200 mg/L] trans-decalin [200 mg/L] benzo[a]anthracene [200 mg/L] benzo[a]pyrene [200 mg/L] fluorene [200 mg/L] pyrene [200 mg/L] octane (C8) [200 mg/L] o-xylene [200 mg/L] sec-butylbenzene [200 mg/L] n-propylbenzene [200 mg/L] methylcyclohexane [200 mg/L] (methyl-CP)Mn(I) tricarbonyl [200 mg/L] coronene [200 mg/L]	1-methylnaphthalene [200 mg/L] dibenzofuran [200 mg/L] n-octadecane (C18) [200 mg/L] acenaphthene [200 mg/L] benzo[b]fluoranthene [200 mg/L] chrysene [200 mg/L] indeno[1,2,3-cd]pyrene [200 mg/L] 2-methylnaphthalene [200 mg/L] benzene [200 mg/L] m-xylene [100 mg/L] tert-butylbenzene [200 mg/L] styrene [200 mg/L] phytane [200 mg/L] indane [200 mg/L]
Semivolatile Mix, 10-2206				
DRE-GA09000378ME	Semivolatile Mix, 10-2206 1000 µg/mL in Methanol(‡)			1ml
	dimethyltin dichloride di-n-propyltin dichloride di-n-butyltin dichloride tetra-n-butyltin triphenyltin chloride	methyltin trichloride n-butyltin trichloride tri-n-butyltin chloride di-n-octyltin dichloride d-(-)-quinic acid		
Semivolatile Mix, 9				
DRE-GA09000229DI	Semivolatile Mix 1000 µg/mL in Dichloromethane(‡)			1ml
	n-nitrosodiethylamine N-nitrosopyrrolidine 1,2,4,5-tetrachlorobenzene parathion benzaldehyde	n-nitrosodi-n-butylamine pentachlorobenzene caprolactam benzoic acid		

Volatile organic compounds (VOCs)

Product code	Description	
Trihalomethane Solution		
DRE-GS09000167ME	Trihalomethane Mix 200 µg/mL in p/t Methanol(±)	5x1ml
	bromodichloromethane bromoform chloroform dibromochloromethane	
VOC Mix A		
DRE-GS09000159ME	VOC Mix A in p/t Methanol(±)	5x1ml
	a-methylstyrene [2000 mg/L] trans-crotonaldehyde [10000 mg/L] epichlorohydrin [10000 mg/L] isooctane [2000 mg/L] 1,4-diethylbenzene [2000 mg/L]	benzyl chloride [2000 mg/L] dicyclopentadiene [2000 mg/L] ethyl acetate [5000 mg/L] butyl acetate [5000 mg/L]
VOC Mix B		
DRE-GS09000171DS	VOC Mix B 1000 µg/mL in Dimethyl Sulfoxide(±)	5x1ml
	1-heptanol butyraldehyde Methyl Heptanoate octane (C8) propionaldehyde	1-octen-3-ol crotonaldehyde Methyl Nonanoate Nonylaldehyde valeraldehyde
		trans,trans-2,4-decadienal decylaldehyde C8:0 methyl octanate n-pentane (C5)
		1-pentanol hexanal heptane (C7) 1-octanal
VOC Mix, 4 components		
DRE-GS09000172	VOC Mix, 4 components(±)	5x1ml
	isopropylbenzene [78000 mg/L] 1,2,3-trimethylbenzene [250000 mg/L] 1,2,4-trimethylbenzene [250000 mg/L] 1,3,5-trimethylbenzene [250000 mg/L]	
VOC Mix, 10 components		
DRE-GA09000393	VOC Mix, 10 components 0.01 Wt %	500ml
DRE-GS09000331SP	VOC Mix, 10 components 0.01 Wt %(±)	5x2ml
	a-methylstyrene [0,01 wt%] ethylbenzene [0,01 wt%] toluene [0,01 wt%] tert-butylbenzene [0,01 wt%] 1,3-diisopropylbenzene [0,01 wt%]	benzene [0,01 wt%] n-propylbenzene [0,01 wt%] sec-butylbenzene [0,01 wt%] 4-isopropyltoluene [0,01 wt%] isopropylbenzene [99,91 wt%]
VOC Mix, 29 components		
DRE-GS09000063DM	VOC Mix, 29 components 10000 µg/mL in Dimethyl Formamide(±)	5x1ml
	1,1,1-trichloroethane 1,2-dichloroethane benzene chloroform ethyl acetate methanol toluene p-xylene	1,1,2-trichloroethane 1,2-dimethoxyethane butyl acetate cyclohexane heptane (C7) methyl t-butyl ether trichloroethylene
		1,1-dichloroethane 1-butanol carbon tetrachloride methylene chloride n-hexane (C6) 1-methyl-2-pyrrolidinone o-xylene
		1,1-dichloroethylene 2-methoxyethanol chlorobenzene ethanol isopropyl alcohol 1,2,3,4-tetrahydronaphthalene m-xylene
VOC Mix, 31		
DRE-GS09000363DM	VOC Mix 10000 µg/mL in Dimethyl Formamide(±)	5x1ml
	1,1,1-trichloroethane 1,2-dichloroethane benzene chloroform ethyl acetate methanol toluene p-xylene	1,1,2-trichloroethane 1,2-dimethoxyethane butyl acetate cyclohexane heptane (C7) methyl t-butyl ether trichloroethylene 2-methylpentane
		1,1-dichloroethane 1-butanol carbon tetrachloride methylene chloride n-hexane (C6) 1-methyl-2-pyrrolidinone o-xylene 3-methylpentane
		1,1-dichloroethylene 2-methoxyethanol chlorobenzene ethanol isopropyl alcohol 1,2,3,4-tetrahydronaphthalene m-xylene

Volatile organic compounds (VOCs)

Product code	Description	
VOC Mix, 62 components		
DRE-GA09000154ME	VOC Mix, 62 components 2000 µg/mL in p/t Methanol(‡)	1.3ml
benzene	ethylbenzene	m-xylene
isopropylbenzene	n-propylbenzene	o-xylene
sec-butylbenzene	tert-butylbenzene	1,2,4-trimethylbenzene
n-butylbenzene	naphthalene	4-isopropyltoluene
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene
4-chlorotoluene	chlorobenzene	1,2,3-trichlorobenzene
bromobenzene	bromochloromethane	carbon tetrachloride
methylene chloride	bromodichloromethane	bromoform
dibromochloromethane	cis-1,2-dichloroethylene	trans-1,2-dichloroethylene
1,1-dichloroethane	1,1,1-trichloroethane	2,2-dichloropropane
1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane	1,1,2-trichloroethane
trichloroethylene	1,2-dibromo-3-chloropropane	1,2-dibromoethane
1,1-dichloropropylene	1,2,3-trichloropropane	1,2-dichloropropane
cis-1,3-dichloropropylene	1,3-dichloropropane	trichlorofluoromethane
vinyl chloride	chloromethane	chloroethane
carbon disulfide	methyl t-butyl ether	
	toluene	
	p-xylene	
	1,3,5-trimethylbenzene	
	styrene	
	2-chlorotoluene	
	1,2,4-trichlorobenzene	
	dibromomethane	
	chloroform	
	1,1-dichloroethylene	
	tetrachloroethylene	
	1,2-dichloroethane	
	hexachlorobutadiene	
	trans-1,3-dichloropropylene	
	bromomethane	
	dichlorodifluoromethane	
VOC-Mix 1		
DRE-XA05000001ME	VOC-Mix 1 100 µg/mL in Methanol (*)	1ml
	Bromodichloromethane	
	Dibromochloromethane	
	Tribromomethane	
	Trichloromethane	
VOC-Mix 2		
DRE-YA05000002ME	VOC-Mix 2 2000 µg/mL in Methanol (*)	1ml
	Bromomethane	
	Chloroethane	
	Chloromethane	
	Dichlorodifluoromethane	
	Fluorotrichloromethane	
	Vinyl Chloride	
VOC-Mix 7		
DRE-YA05000007ME	VOC-Mix 7 2000 µg/mL in Methanol (*)	1ml
	1,1,1-Trichloroethane	1,1-Dichloroethene
	1,2-Dichloroethane	1,4-Dichlorobenzene
	Benzene	Bromodichloromethane
	Dibromochloromethane	Tetrachloromethane
	Tribromomethane	Trichloroethene
	Trichloromethane	
VOC-Mix 8		
DRE-YA05000008ME	VOC-Mix 8 2000 µg/mL in Methanol (*)	1ml
	1,2-Dichlorobenzene	1,2-Dichloropropane
	Chlorobenzene	cis-1,2-Dichloroethene
	Ethylbenzene	m-Xylene
	o-Xylene	p-Xylene
	Styrene	Tetrachloroethene
	Toluene	trans-1,2-Dichloroethene
VOC-Mix 9		
DRE-YA05000009AC	VOC-Mix 9 1000 µg/mL in Acetone	1ml
	1,1,1-Trichloroethane	Benzene
	cis-1,2-Dichloroethene	Dichloromethane
	m-Xylene	o-Xylene
	p-Xylene	Tetrachloroethene
	Tetrachloromethane	Toluene
	Trichloroethene	Trichloromethane

Volatile organic compounds (VOCs)

Product code	Description		
VOC-Mix 15			
DRE-XA05000015ME	VOC-Mix 15 200 µg/mL in Methanol (*)		1ml
1,1,1,2-Tetrachloroethane	1,1,1-Trichloroethane	1,1,2,2-Tetrachloroethane	1,1,2-Trichloroethane
1,1-Dichloro-1-propene	1,1-Dichloroethane	1,1-Dichloroethene	1,2,3-Trichlorobenzene
1,2,3-Trichloropropane	1,2,4-Trichlorobenzene	1,2,4-Trimethylbenzene	1,2-Dibromo-3-chloropropane
1,2-Dibromoethane	1,2-Dichlorobenzene	1,2-Dichloroethane	1,2-Dichloropropane
1,3,5-Trimethylbenzene	1,3-Dichlorobenzene	1,3-Dichloropropane	1,3-Dichloropropene (cis + trans)
1,4-Dichlorobenzene	2,2-Dichloropropane	2-Chlorotoluene	4-Chlorotoluene
4-Isopropyltoluene	Benzene	Bromobenzene	Bromochloromethane
Bromodichloromethane	Chlorobenzene	cis-1,2-Dichloroethene	Dibromochloromethane
Dibromomethane	Dichloromethane	Ethylbenzene	Hexachloro-1,3-butadiene
Isopropylbenzene	m-Xylene	Naphthalene	n-Butylbenzene
n-Propylbenzene	o-Xylene	p-Xylene	sec-Butylbenzene
Styrene	tert-Butylbenzene	Tetrachloroethene	Tetrachloromethane
Toluene	trans-1,2-Dichloroethene	Tribromomethane	Trichloroethene
Trichloromethane			
VOC-Mix 20			
DRE-XA05000020ME	VOC-Mix 20 200 µg/mL in Methanol (*)		1ml
1,1,1,2-Tetrachloroethane	1,1,1-Trichloroethane	1,1,2,2-Tetrachloroethane	1,1,2-Trichloroethane
1,1-Dichloro-1-propene	1,1-Dichloroethane	1,1-Dichloroethene	1,2,3-Trichlorobenzene
1,2,3-Trichloropropane	1,2,4-Trichlorobenzene	1,2,4-Trimethylbenzene	1,2-Dibromo-3-chloropropane
1,2-Dibromoethane	1,2-Dichlorobenzene	1,2-Dichloroethane	1,2-Dichloropropane
1,3,5-Trimethyl benzene	1,3-Dichlorobenzene	1,3-Dichloropropane	1,3-Dichloropropene (cis + trans)
1,4-Dichlorobenzene	2,2-Dichloropropane	2-Chlorotoluene	4-Chlorotoluene
4-Isopropyltoluene	Benzene	Bromobenzene	Bromochloromethane
Bromodichloromethane	Bromomethane	Chlorobenzene	Chloroethane
Chloromethane	cis-1,2-Dichloroethene	Dibromochloromethane	Dibromomethane
Dichlorodifluoromethane	Dichloromethane	Ethylbenzene	Fluorotrichloromethane
Hexachloro-1,3-butadiene	Isopropylbenzene	m-Xylene	Naphthalene
n-Butylbenzene	n-Propylbenzene	o-Xylene	p-Xylene
sec-Butylbenzene	Styrene	tert-Butylbenzene	Tetrachloroethene
Tetrachloromethane	Toluene	trans-1,2-Dichloroethene	Tribromomethane
Trichloroethene	Trichloromethane		
VOC-Mix 21			
DRE-XA05000021ME	VOC-Mix 21 200 µg/mL in Methanol (*)		1ml
	1,1,1-Trichloroethane	1,2-Dichloroethane	
	Dibromochloromethane	Tetrachloroethene	
	Tetrachloromethane	Tribromomethane	
	Trichloroethene	Trichloromethane	
VOC-Mix 23			
DRE-XA05000023ME	VOC-Mix 23 in Methanol (*)		1ml
	Bromodichloromethane [50 µg/mL]		
	Dibromochloromethane [50 µg/mL]		
	Tetrachloroethene [20 µg/mL]		
	Tetrachloromethane [6 µg/mL]		
	Tribromomethane [50 µg/mL]		
	Trichloroethene [60 µg/mL]		
	Trichloromethane [50 µg/mL]		
VOC-Mix 61			
DRE-YA05000061ME	VOC-Mix 61 in Methanol (*)		1ml
	1,1,1-Trichloroethane [1000 µg/mL]	1,2-Dichloroethane [10000 µg/mL]	
	Bromodichloromethane [1000 µg/mL]	cis-1,2-Dichloroethene [10000 µg/mL]	
	Dibromochloromethane [1000 µg/mL]	Dichloromethane [5000 µg/mL]	
	Tetrachloroethene [1000 µg/mL]	Tetrachloromethane [1000 µg/mL]	
	Tribromomethane [1000 µg/mL]	Trichloroethene [1000 µg/mL]	
	Trichloromethane [1000 µg/mL]		










Volatile organic compounds (VOCs)

Product code	Description	
Volatile Mix, 10-1246		
DRE-GA09000352IP	Volatile Mix, 10-1246 1000 µg/mL in Isopropanol(‡)	1ml
	α-humulene mycrene b-pinene (-)-caryophyllene oxide (-)-isocaryophyllene linalool a-pinene terpinolene α-bisabolol limonene	
Volatile Mix, 15-977		
DRE-YA09000035DS	VOC Mix in Dimethyl Sulfoxide(‡)	1ml
DRE-YA09000035DS-E	VOC Mix in Dimethyl Sulfoxide(‡)	10x1ml
	n-hexane (C6) [1000 mg/L] heptane (C7) [1000 mg/L] ethanol [1000 mg/L] acetonitrile [1000 mg/L] toluene [1000 mg/L] carbon tetrachloride [1000 mg/L] o-xylene [1000 mg/L] p-xylene [500 mg/L] n-pentane (C5) [1000 mg/L] isopropyl alcohol [1000 mg/L] acetone [1000 mg/L] tetrahydrofuran [1000 mg/L] chloroform [1000 mg/L] benzene [1000 mg/L] m-xylene [500 mg/L]	
Volatile Mix, 2-1620		
DRE-GA09000365ME	Volatile Mix, 2-1620 1000 µg/mL in p/t Methanol(‡)	5ml
DRE-GS09000380ME	Volatile Mix, 2-1620 10000 µg/mL in p/t Methanol(‡)	5x5ml
	1,1,1,2-Tetrafluoroethane 1,1-difluoroethane	










DYES AND METABOLITES











Dyes and metabolites

Product code	Description			
Acid Blue 1				
CAS 129-17-9 DRE-C10028100	MW 566.6646 Acid Blue 1	$C_{27}H_{31}N_2O_6S_2 \cdot Na$	100mg	
Acid Blue 3 calcium salt				
CAS 3536-49-0 DRE-C10028130	MW 1159.4265 Acid Blue 3 calcium salt	$2C_{27}H_{31}N_2O_7S_2 \cdot Ca$	100mg	
Acid Blue 62				
CAS 4368-56-3 DRE-C10028155	MW 422.43 Acid Blue 62	$C_{20}H_{19}N_2O_5S \cdot Na$	100mg	
Acid Green 1 (technical product)				
CAS 19381-50-1 DRE-C10028200	MW 878.4605 Acid Green 1 (technical product)	$C_{30}H_{15}FeN_3O_{15}S_3 \cdot 3Na$	100mg	
Acid Green 16				
CAS 12768-78-4 DRE-C10028400	MW 616.7233 Acid Green 16	$C_{31}H_{33}N_2O_6S_2 \cdot Na$	100mg	
Acid Green 50				
CAS 3087-16-9 DRE-C10028450	MW 576.6164 Acid Green 50	$C_{27}H_{25}N_2O_7S_2 \cdot Na$	100mg	
Acid Orange 6 (4-[2-(2,4-Dihydroxyphenyl)diazenyl]benzenesulfonic Acid Sodium Salt)				
CAS 547-57-9 DRE-C10028600	MW 316.265 Acid Orange 6	$C_{12}H_9N_2O_5S \cdot Na$	100mg	
Acid Red 26				
CAS 3761-53-3 DRE-C10028800	MW 480.4225 Acid Red 26	$C_{18}H_{14}N_2O_7S_2 \cdot 2Na$	100mg	
Acid Red 2G (E 128)				
CAS 3734-67-6 DRE-C10028700	MW 509.4207 Acid Red 2G (E 128)	$C_{18}H_{13}N_3O_8S_2 \cdot 2Na$	100mg	










Dyes and metabolites

Product code	Description			
Acid Red 44				
CAS 2766-77-0 DRE-C10028812	MW 502.428 Acid Red 44	C ₂₀ H ₁₂ N ₂ O ₇ S ₂ ·2Na	25mg	
Acid Red 52				
CAS 3520-42-1 DRE-C10028820	MW 580.6481 Acid Red 52	C ₂₇ H ₂₈ N ₂ O ₇ S ₂ ·Na	100mg	
Acid Red 73				
CAS 5413-75-2 DRE-C10028830	MW 556.4787 Acid Red 73	C ₂₂ H ₁₄ N ₄ O ₇ S ₂ ·2Na	100mg	
Acid Red 87 (Eosin Y)				
CAS 17372-87-1 DRE-C10028840	MW 691.8542 Acid Red 87	C ₂₀ H ₆ Br ₄ O ₅ ·2Na	100mg	
Acid Red 95				
CAS 33239-19-9 DRE-C10028843	MW 628.063 Acid Red 95	C ₂₀ H ₆ I ₂ O ₅ ·2Na	100mg	
Acid Violet 17				
CAS 4129-84-4 DRE-C10028880	MW 761.9243 Acid Violet 17	C ₄₁ H ₄₄ N ₃ O ₆ S ₂ ·Na	10mg	
Acid Violet 49				
CAS 1694-09-3 DRE-C10028900	MW 733.8712 Acid Violet 49	C ₃₉ H ₄₀ N ₃ O ₆ S ₂ ·Na	100mg	
Acid Violet 9				
CAS 6252-76-2 DRE-C10028870	MW 612.6269 Acid Violet 9	C ₃₄ H ₂₅ N ₂ O ₆ S·Na	100mg	
Acid Yellow 11				
CAS 6359-82-6 DRE-C10029020	MW 380.3536 Acid Yellow 11	C ₁₆ H ₁₃ N ₄ O ₄ S·Na	100mg	










Dyes and metabolites

Product code	Description			
Acid Yellow 23 Aluminium lake				
CAS 12225-21-7	MW n/a			No Structure
DRE-E10029060	Acid Yellow 23 Aluminium lake		100mg	
Acid Yellow 36				
CAS 587-98-4	MW 375.3768	C ₁₈ H ₁₄ N ₃ O ₃ S·Na		
DRE-C10029500	Acid Yellow 36		100mg	
Acid mordant Black 1				
CAS 3618-58-4	MW 461.3799	C ₂₀ H ₁₂ N ₃ O ₇ S·Na		
DRE-C10028500	Acid mordant Black 1		100mg	
Allura Red AC (E 129)				
CAS 25956-17-6	MW 496.4219	C ₁₈ H ₁₄ N ₂ O ₈ S ₂ ·2Na		
DRE-C10125000	Allura Red AC (E 129)		100mg	
Amaranth				
CAS 915-67-3	MW 604.473	C ₂₀ H ₁₁ N ₂ O ₁₀ S ₃ ·3Na		
DRE-C10148500	Amaranth(†)		250mg	
4-Aminoazobenzene				
CAS 60-09-3	MW 197.2358	C ₁₂ H ₁₁ N ₃		
DRE-C10167500	4-Aminoazobenzene		100mg	
DRE-XA10167500AL	4-Aminoazobenzene 100 µg/mL in Acetonitrile		1ml	
4-Amino-2,3-dimethylazobenzene				
CAS 97-56-3	MW 225.289	C ₁₄ H ₁₅ N ₃		
DRE-C10202000	4-Amino-2,3-dimethylazobenzene		250mg	
DRE-XA10202000CY	4-Amino-2,3-dimethylazobenzene 100 µg/mL in Cyclohexane		1ml	
2-Amino-3-nitrotoluene (2-Methyl-6-nitroaniline)				
CAS 570-24-1	MW 152.1506	C ₇ H ₈ N ₂ O ₂		
DRE-C10207400	2-Amino-3-nitrotoluene (2-Methyl-6-nitroaniline)		250mg	
2-Amino-4-nitrotoluene (2-Methyl-5-nitroaniline)				
CAS 99-55-8	MW 152.1506	C ₇ H ₈ N ₂ O ₂		
DRE-C10207500	2-Amino-4-nitrotoluene (2-Methyl-5-nitroaniline)		250mg	









Dyes and metabolites

Product code	Description			
2-Amino-5-nitrotoluene (2-Methyl-4-nitroaniline)				
CAS 99-52-5 DRE-C10207600	MW 152.1506 2-Amino-5-nitrotoluene (2-Methyl-4-nitroaniline)	$C_7H_8N_2O_2$	250mg	
2-Amino-6-nitrotoluene (2-Methyl-3-nitroaniline)				
CAS 603-83-8 DRE-C10207700	MW 152.1506 2-Amino-6-nitrotoluene (2-Methyl-3-nitroaniline)	$C_7H_8N_2O_2$	250mg	
3-Amino-4-nitrotoluene (5-Methyl-2-nitroaniline)				
CAS 578-46-1 DRE-C10207900	MW 152.1506 3-Amino-4-nitrotoluene (5-Methyl-2-nitroaniline)	$C_7H_8N_2O_2$	100mg	
4-Amino-2-nitrotoluene (4-Methyl-3-nitroaniline)				
CAS 119-32-4 DRE-C10208200	MW 152.1506 4-Amino-2-nitrotoluene (4-Methyl-3-nitroaniline)	$C_7H_8N_2O_2$	250mg	
4-Amino-3-nitrotoluene (4-Methyl-2-nitroaniline)				
CAS 89-62-3 DRE-C10208300	MW 152.1506 4-Amino-3-nitrotoluene (4-Methyl-2-nitroaniline)	$C_7H_8N_2O_2$	250mg	
4-Aminophenylether (4,4'-Oxydianiline)				
CAS 101-80-4 DRE-C10215000	MW 200.2365 4-Aminophenylether (4,4'-Oxydianiline)	$C_{12}H_{12}N_2O$	250mg	
Aniline				
CAS 62-53-3 DRE-CA10262500 DRE-L10262500CY DRE-XA10262500CY	MW 93.1265 Aniline Aniline 10 µg/mL in Cyclohexane Aniline 100 µg/mL in Cyclohexane	C_6H_7N	1ml 10ml 1ml	
Aniline D5				
CAS 4165-61-1 DRE-C10262600 DRE-YA10262600MB	MW 98.1573 Aniline D5 Aniline D5 2000 µg/mL in Methyl-tert. butyl ether	$C_6^2H_5H_2N$	100mg 1.1ml	
Azorubine (E122)				
CAS 3567-69-9 DRE-C10411500	MW 502.428 Azorubin (E122)	$C_{20}H_{12}N_2O_7S_2 \cdot 2Na$	100mg	










Dyes and metabolites

Product code	Description			
Basic Blue 7				
CAS 2390-60-5 DRE-C10424000	MW 514.1438 Basic Blue 7	C ₃₃ H ₄₀ N ₃ ·Cl	100mg	
Basic Green 1				
CAS 633-03-4 DRE-C10424500	MW 482.6349 Basic Green 1	C ₂₇ H ₃₃ N ₂ ·HO ₃ S	100mg	
Basic Red 1				
CAS 989-38-8 DRE-C10424950	MW 479.0103 Basic Red 1	C ₂₈ H ₃₀ N ₂ O ₃ ·ClH	100mg	
Basic Red 2				
CAS 477-73-6 DRE-C10424960	MW 350.8447 Basic Red 2	C ₂₀ H ₁₉ N ₄ ·Cl	100mg	
Basic Red 9				
CAS 569-61-9 DRE-C10425000	MW 323.8193 Basic Red 9	C ₁₉ H ₁₇ N ₃ ·ClH	100mg	
Basic Violet 1				
CAS 603-47-4 DRE-C10427100	MW 393.9522 Basic Violet 1	C ₂₄ H ₂₇ N ₃ ·ClH	100mg	
Basic Violet 3 (Methylrosanilinium Chloride)				
CAS 548-62-9 DRE-C10427500	MW 407.9788 Basic Violet 3(±)	C ₂₅ H ₃₀ N ₃ ·Cl	100mg	
Basic Violet 3 D6				
CAS 1266676-01-0 DRE-C10427505	MW 414.0158 Basic Violet 3 D6	C ₂₅ ² H ₆ H ₂₄ N ₃ ·Cl	10mg	
4,4'-Benzidine				
CAS 92-87-5 DRE-C10536000 DRE-L10536000AL DRE-XA10536000AL	MW 184.2371 4,4'-Benzidine 4,4'-Benzidine 10 µg/mL in Acetonitrile 4,4'-Benzidine 100 µg/mL in Acetonitrile	C ₁₂ H ₁₂ N ₂	100mg 10ml 1ml	










Dyes and metabolites

Product code	Description			
Brilliant Black BN				
CAS 2519-30-4 DRE-C10665000	MW 867.6788 Brilliant Black BN	C ₂₈ H ₁₇ N ₅ O ₁₄ S ₄ ·4Na	250mg	
Brilliant Blue FCF (E133)				
CAS 3844-45-9 DRE-C10665200	MW 792.8484 Brilliant Blue FCF (E133)	C ₃₇ H ₃₄ N ₂ O ₉ S ₃ ·2Na	250mg	
Carminic acid (E120)				
CAS 1260-17-9 DRE-C11044000	MW 492.3864 Carminic acid (E120)(‡)	C ₂₂ H ₂₀ O ₁₃	100mg	
Chocolate Brown				
CAS n/a DRE-C11665150	MW n/a Chocolate Brown		100mg	No Structure
4,4'-Diaminobiphenylmethane (Bis-(4-aminophenyl)methane)				
CAS 101-77-9 DRE-C10648000	MW 198.2637 Bis-(4-aminophenyl)methane	C ₁₃ H ₁₄ N ₂	250mg	
4,4'-Diamino-3,3'-dimethyldiphenylmethane				
CAS 838-88-0 DRE-C12194800	MW 226.3168 4,4'-Diamino-3,3'-dimethyldiphenyl methane(‡)	C ₁₅ H ₁₆ N ₂	100mg	
4,4'-Diaminodiphenyl Sulfide (4-Aminophenylthioether)				
CAS 139-65-1 DRE-C10215500 DRE-L10215500AL	MW 216.3021 4-Aminophenylthioether 4-Aminophenylthioether(4,4'-Thiodianiline) 10 µg/mL in Acetonitrile	C ₁₂ H ₁₂ N ₂ S	250mg 10ml	
4',5'-Dibromofluorescein				
CAS 596-03-2 DRE-C12240230	MW 490.0984 4',5'-Dibromofluorescein (technical)	C ₂₀ H ₁₀ Br ₂ O ₅	100mg	
3,3'-Dichlorobenzidine Dihydrochloride				
CAS 612-83-9 DRE-C12378000	MW 326.0491 3,3'-Dichlorobenzidine dihydrochloride	C ₁₂ H ₁₀ Cl ₂ N ₂ ·2ClH	100mg	










Dyes and metabolites

Product code	Description			
3,3'-Dichlorobenzidine				
CAS 91-94-1	MW 253.1272	C ₁₂ H ₁₀ Cl ₂ N ₂		
DRE-C12377900	3,3'-Dichlorobenzidine(‡)		100mg	
DRE-L12377900AL	3,3'-Dichlorobenzidine 10 µg/mL in Acetonitrile		10ml	
3,3'-Dichlorobenzidine D6 (ring D6)				
CAS 93951-91-8	MW 259.1642	C ₁₂ ² H ₆ H ₄ Cl ₂ N ₂		
DRE-XA12377910AL	3,3'-Dichlorobenzidine D6 (ring D6) 100 µg/mL in Acetonitrile		1.1ml	
4-(Diethylamino)azobenzene				
CAS 2481-94-9	MW 253.3422	C ₁₆ H ₁₈ N ₂		
DRE-C12604600	4-(Diethylamino)azobenzene		100mg	
3,3'-Dimethoxybenzidine				
CAS 119-90-4	MW 244.289	C ₁₄ H ₁₆ N ₂ O ₂		
DRE-C12721000	3,3'-Dimethoxybenzidine(‡)		100mg	
3,3'-Dimethylbenzidine (o-Tolidine)				
CAS 119-93-7	MW 212.2902	C ₁₄ H ₁₆ N ₂		
DRE-C12726000	3,3'-Dimethylbenzidine (o-Tolidine)		100mg	
Direct Black 38				
CAS 1937-37-7	MW 781.7279	C ₃₄ H ₂₅ N ₉ O ₇ S ₂ ·2Na		
DRE-C12965000	Direct Black 38		100mg	
Direct Blue 86 (technical product)				
CAS 1330-38-7	MW 780.1591	C ₃₂ H ₁₄ CuN ₈ ·2Na·2O ₃ S		
DRE-C12965300	Direct Blue 86 (technical product)		100mg	
Direct Brown 95				
CAS 16071-86-6	MW 760.1	C ₃₁ H ₁₈ CuN ₆ O ₉ S·2Na		
DRE-C12965350	Direct Brown 95 (technical mixture)		100mg	
Direct Red 28				
CAS 573-58-0	MW 696.6632	C ₃₂ H ₂₂ N ₆ O ₆ S ₂ ·2Na		
DRE-C12965400	Direct Red 28		100mg	










Dyes and metabolites

Product code	Description			
Disperse Black 2				
CAS 6232-57-1 DRE-E12972005	MW 256.303 Disperse Black 2	$C_{14}H_{16}N_4O$	100mg	
Disperse Blue 1				
CAS 2475-45-8 DRE-C12972010	MW 268.2707 Disperse Blue 1	$C_{14}H_{12}N_4O_2$	100mg	
Disperse Blue 102				
CAS 69766-79-6 DRE-C12972027	MW 365.4075 Disperse Blue 102	$C_{15}H_{19}N_5O_4S$	100mg	
Disperse Blue 106				
CAS 68516-81-4 DRE-C12972030	MW 335.3815 Disperse Blue 106	$C_{14}H_{17}N_5O_3S$	100mg	
Disperse Blue 124				
CAS 15141-18-1 DRE-C12972040	MW 377.4182 Disperse Blue 124	$C_{16}H_{19}N_5O_4S$	25mg	
Disperse Blue 26				
CAS 3860-63-7 DRE-C12972019	MW 298.2934 Disperse Blue 26 (Technical product, mixture of congeners)	$C_{16}H_{14}N_2O_4$	100mg	
Disperse Blue 291				
CAS 56548-64-2 DRE-C12972047	MW 509.3106 Disperse Blue 291	$C_{19}H_{21}BrN_6O_6$	100mg	
Disperse Blue 3				
CAS 2475-46-9 DRE-C12972013	MW 296.3205 Disperse Blue 3	$C_{17}H_{16}N_2O_3$	100mg	
Disperse Blue 7				
CAS 3179-90-6 DRE-C12972017	MW 358.3453 Disperse Blue 7	$C_{18}H_{18}N_2O_6$	100mg	










Dyes and metabolites

Product code	Description			
Disperse Brown 1				
CAS 23355-64-8 DRE-C12972070	MW 433.6737 Disperse Brown 1	$C_{16}H_{15}Cl_3N_4O_4$	100mg	
Disperse Orange 1				
CAS 2581-69-3 DRE-C12972101	MW 318.3294 Disperse Orange 1(‡)	$C_{18}H_{14}N_4O_2$	25mg	
Disperse Orange 11				
CAS 82-28-0 DRE-C12972111	MW 237.2533 Disperse Orange 11(‡)	$C_{15}H_{11}NO_2$	100mg	
Disperse Orange 13				
CAS 6253-10-7 DRE-C12972113	MW 352.3886 Disperse Orange 13	$C_{22}H_{16}N_4O$	100mg	
Disperse Orange 149				
CAS 85136-74-9 DRE-C12972150	MW 458.5123 Disperse Orange 149(‡)	$C_{25}H_{26}NaO_3$	100mg	
Disperse Orange 3				
CAS 730-40-5 DRE-C12972110	MW 242.2334 Disperse Orange 3(‡)	$C_{12}H_{10}N_4O_2$	100mg	
Disperse Orange 37				
CAS 13301-61-6 DRE-C12972120	MW 392.2393 Disperse Orange 37(‡)	$C_{17}H_{15}Cl_2N_5O_2$	25mg	
Disperse Orange 61				
CAS 55281-26-0 DRE-C12972135	MW 481.1413 Disperse Orange 61	$C_{17}H_{15}Br_2N_5O_2$	100mg	
Disperse Red 1				
CAS 2872-52-8 DRE-C12972210	MW 314.3391 Disperse Red 1	$C_{16}H_{18}N_4O_3$	100mg	

Dyes and metabolites

Product code	Description			
Disperse Red 11				
CAS 2872-48-2 DRE-C12972221	MW 268.2674 Disperse Red 11(±)	C ₁₅ H ₁₂ N ₂ O ₃	100mg	
Disperse Red 13				
CAS 3180-81-2 DRE-C12972223	MW 348.7842 Disperse Red 13	C ₁₆ H ₁₇ ClN ₄ O ₃	100mg	
Disperse Red 15				
CAS 116-85-8 DRE-C12972225	MW 239.2262 Disperse Red 15	C ₁₄ H ₉ NO ₃	100mg	
Disperse Red 151				
CAS 70210-08-1 DRE-C12972264	MW 531.5829 Disperse Red 151	C ₂₇ H ₂₅ N ₅ O ₅ S	10mg	
Disperse Red 151 Metabolite 1				
CAS 58104-55-5 DRE-C12972265	MW 489.5462 Disperse Red 151 Metabolite 1	C ₂₅ H ₂₃ N ₅ O ₄ S	10mg	
Disperse Red 17				
CAS 3179-89-3 DRE-C12972227	MW 344.3651 Disperse Red 17(±)	C ₁₇ H ₂₀ N ₄ O ₄	100mg	
Disperse Red 7				
CAS 4540-00-5 DRE-C12972217	MW 364.7836 Disperse Red 7	C ₁₆ H ₁₇ ClN ₄ O ₄	10mg	
Disperse Violet 1				
CAS 128-95-0 DRE-C12972287	MW 238.2414 Disperse Violet 1	C ₁₄ H ₁₀ N ₂ O ₂	100mg	
Disperse Violet 93				
CAS 52697-38-8 DRE-C12972290	MW 479.2847 Disperse Violet 93	C ₁₈ H ₁₉ BrN ₆ O ₅	100mg	










Dyes and metabolites

Product code	Description			
Disperse Violet 93 (Cl derivative)				
CAS 66557-45-7 DRE-C12972295	MW 434.8337 Disperse Violet 93 (Cl derivative)	$C_{18}H_{19}ClN_6O_5$	100mg	
Disperse Yellow 1				
CAS 119-15-3 DRE-C12972308	MW 275.217 Disperse Yellow 1(‡)	$C_{12}H_9N_3O_5$	100mg	
Disperse Yellow 23				
CAS 6250-23-3 DRE-C12972323	MW 302.33 Disperse Yellow 23(‡)	$C_{18}H_{14}N_4O$	100mg	
Disperse Yellow 3				
CAS 2832-40-8 DRE-C12972310	MW 269.2985 Disperse Yellow 3	$C_{15}H_{15}N_3O_2$	100mg	
Disperse Yellow 3 D3 (acetyl D3)				
CAS 947601-96-9 DRE-XA12972311AL	MW 272.317 Disperse Yellow 3 D3 100 µg/mL in Acetonitrile	$C_{15}^2H_{15}H_{12}N_3O_2$	1.1ml	
Disperse Yellow 39				
CAS 12236-29-2 DRE-C12972339	MW 264.3217 Disperse Yellow 39(‡)	$C_{17}H_{16}N_2O$	10mg	
Disperse Yellow 39 surrogate				
CAS 56208-37-8 DRE-C12972338	MW 290.316 Disperse Yellow 39 surrogate	$C_{18}H_{14}N_2O_2$	100mg	
Disperse Yellow 49				
CAS 6858-49-7 DRE-C12972349	MW 374.4357 Disperse Yellow 49(‡)	$C_{22}H_{22}N_4O_2$	100mg	
Disperse Yellow 54				
CAS 7576-65-0 DRE-C12972364	MW 289.2848 Disperse Yellow 54	$C_{18}H_{11}NO_3$	25mg	










Dyes and metabolites

Product code	Description			
Disperse Yellow 56-methyl				
CAS 73287-67-9 DRE-C12972367	MW 383.4027 Disperse Yellow 56-methyl (technical)	C ₂₂ H ₁₇ N ₅ O ₂	100mg	
Disperse Yellow 7				
CAS 6300-37-4 DRE-C12972317	MW 316.3565 Disperse Yellow 7	C ₁₉ H ₁₆ N ₄ O	100mg	
Disperse Yellow 9				
CAS 6373-73-5 DRE-C12972319 DRE-L12972319AL	MW 274.2322 Disperse Yellow 9 Disperse Yellow 9, 10 µg/mL in Acetonitrile (*)	C ₁₂ H ₁₀ N ₄ O ₄	25mg 10ml	
Erythrosine B Disodium Salt (E127)				
CAS 16423-68-0 DRE-C13205000	MW 879.8561 Erythrosine B (E127)	C ₂₀ H ₁₄ O ₅ ·2Na	250mg	
Eulan WA new				
CAS 55069-01-7 DRE-E13397500	MW 974.4102 Eulan WA New (technical mixture)	((C ₁₃ H ₆ Cl ₇ NO ₃ S)(C ₁₃ H ₇ Cl ₆ NO ₃ S)c)mix	10mg	
Fast Green FCF				
CAS 2353-45-9 DRE-C13406000	MW 808.8478 Fast Green FCF	C ₃₇ H ₃₄ N ₂ O ₁₀ S ₂ ·2Na	100mg	
Fluorescein Sodium				
CAS 518-47-8 DRE-C13750000	MW 376.2699 Fluorescein disodium salt	C ₂₀ H ₁₀ O ₅ ·2Na	100mg	
Fuchsin (Basic Violet 14 HCl)				
CAS 632-99-5 DRE-C10428000	MW 337.8459 Basic Violet 14 HCl	C ₂₀ H ₁₉ N ₃ ·ClH	100mg	
Indigotine (Indigotindisulfonate Sodium; E132)				
CAS 860-22-0 DRE-C14289000	MW 466.3529 Indigotine (E132)	C ₁₆ H ₈ N ₂ O ₈ S ₂ ·2Na	250mg	








Dyes and metabolites

Product code	Description			
Iron Oxide Pigment Black				
CAS 1317-61-9 DRE-C14375000	MW 231.5326 Iron Oxide Pigment Black	FeO·Fe ₂ O ₃	100mg	
Iron Pigment Red (E172)				
CAS 1309-37-1 DRE-C14375200	MW 159.6882 Iron Pigment Red (E172)	Fe ₂ O ₃	100mg	
Iron Pigment Yellow (E172)				
CAS 51274-00-1 DRE-C14375400	MW 88.8517 Iron Pigment Yellow	FeHO ₂	100mg	
Leucocrystal Violet				
CAS 603-48-5 DRE-C14629400	MW 373.5337 Leucocrystal Violet(‡)	C ₂₅ H ₃₁ N ₃	100mg	
Leucocrystal Violet D6				
CAS 1173023-92-1 DRE-C14629401	MW 379.5707 Leucocrystal violet D6	C ₂₅ ² H ₆ H ₂₅ N ₃	10mg	
Leucomalachite Green				
CAS 129-73-7 DRE-C14629500 DRE-L14629500CY	MW 330.4659 Leucomalachite green(‡) Leucomalachite green 10 µg/mL in Cyclohexane	C ₂₃ H ₂₆ N ₂	100mg 10ml	
Leucomalachite Green D6 (phenylmethin D6)				
CAS 1173021-13-0 DRE-C14629510	MW 336.5029 Leucomalachite green D6(‡)	C ₂₃ ² H ₆ H ₂₀ N ₂	10mg	
Malachite Green D5 Oxalate Salt				
CAS n/a DRE-C14680005	MW 937.0664 Malachite green D5 oxalate salt	2C ₂₃ ² H ₅ H ₂₀ N ₂ ·C ₂ H ₂ O ₄ ·2C ₂ HO ₄	10mg	
Malachite Green D5 Picrate				
CAS 1258668-21-1 DRE-C14680010 DRE-XA14680010AC	MW 562.5848 Malachite green D5 picrate(‡) Malachite green D5 picrate 100 µg/mL in Acetone	C ₂₃ ² H ₅ H ₂₀ N ₂ ·C ₆ H ₂ N ₃ O ₇	10mg 1.1ml	









Dyes and metabolites

Product code	Description			
Malachite green oxalate salt				
CAS 2437-29-8 DRE-C14680000	MW 927.0048 Malachite green oxalate salt	$2C_{23}H_{25}N_2 \cdot C_2H_2O_4 \cdot 2C_2HO_4$	250mg	
2-Methoxy-5-methylaniline				
CAS 120-71-8 DRE-C15081000	MW 137.179 2-Methoxy-5-methylaniline	$C_8H_{11}NO$	250mg	
4-Methoxy-1,3-phenylenediamine (2,4-Diaminoanisole)				
CAS 615-05-4 DRE-CA15081900	MW 138.1671 4-Methoxy-1,3-phenylenediamine (2,4-Diaminoanisole) (*)	$C_7H_{10}N_2O$	100mg	
4,4'-Methylene-bis(2-chloroaniline)				
CAS 101-14-4 DRE-C15087500 DRE-L15087500AL	MW 267.1538 4,4'-Methylene-bis(2-chloroaniline) 4,4'-Methylene-bis(2-chloroaniline) 10 µg/mL in Acetonitrile	$C_{13}H_{12}Cl_2N_2$	100mg 10ml	
Navy Blue 018112				
CAS 118685-33-9 DRE-C15492000	MW 947.1333 Navy Blue 018112	$C_{39}H_{23}ClCrN_7O_{12}S \cdot 2Na$	100mg	
New Red				
CAS 220658-76-4 DRE-C15506000	MW 611.4657 New Red(±)	$C_{18}H_{12}N_3O_{11}S_3 \cdot 3Na$	100mg	
Orange 2 Sodium Salt				
CAS 633-96-5 DRE-C15735000	MW 350.3243 Orange 2 sodium salt(±)	$C_{16}H_{11}N_2O_4S \cdot Na$	250mg	
Para Red				
CAS 6410-10-2 DRE-C15875000	MW 293.2768 Para Red	$C_{16}H_{11}N_3O_3$	100mg	
Para Red D6				
CAS 1014689-16-7 DRE-C15875100 DRE-XA15875100AC	MW 299.3138 Para Red D6 Pararot D6 100 µg/mL in Acetone	$C_{16}^2H_6H_5N_3O_3$	10mg 1.1ml	









Dyes and metabolites

Product code	Description			
Patent Blue 5 (E131)				
CAS 20262-76-4 DRE-C15895900	MW 582.664 Patent Blue 5 (E131)	C ₂₇ H ₃₁ N ₂ O ₇ S ₂ ·Na	250mg	
Phloxine B				
CAS 18472-87-2 DRE-C16078000	MW 829.6344 Phloxine B	C ₂₀ H ₂ Br ₄ Cl ₄ O ₅ ·2Na	100mg	
Phthalimide				
CAS 85-41-6 DRE-C16190000	MW 147.1308 Phthalimide	C ₈ H ₅ NO ₂	500mg	
Pigment Blue 15				
CAS 147-14-8 DRE-C16206500	MW 576.069 Pigment Blue 15	C ₃₂ H ₁₆ CuN ₈	100mg	
Pigment Green 7				
CAS 1328-53-6 DRE-C16206700	MW n/a Pigment Green 7		100mg	No Structure
Pigment Orange 5				
CAS 3468-63-1 DRE-C16207000	MW 338.2744 Pigment Orange 5	C ₁₆ H ₁₀ N ₄ O ₅	100mg	
Pigment Red 4				
CAS 2814-77-9 DRE-C16207470	MW 327.7219 Pigment Red 4	C ₁₆ H ₁₀ ClN ₃ O ₃	100mg	
Pigment Red 64				
CAS 6371-76-2 DRE-C16207500	MW 622.6397 Pigment Red 64	2C ₁₇ H ₁₁ N ₂ O ₃ ·Ca	100mg	
Pigment Violet 3				
CAS 1325-82-2 DRE-C16208000	MW n/a Pigment Violet 3		100mg	No Structure










Dyes and metabolites

Product code	Description			
Ponceau 4RC (E124)				
CAS 2611-82-7 DRE-C16284000	MW 604.473 Ponceau 4RC (E124)	C ₂₀ H ₁₁ N ₂ O ₁₀ S ₃ ·3Na	100mg	
Ponceau SX				
CAS 4548-53-2 DRE-C16284500	MW 480.4225 Ponceau SX	C ₁₈ H ₁₄ N ₂ O ₇ S ₂ ·2Na	100mg	
Quinoline				
CAS 91-22-5 DRE-C16709600 DRE-GA09010376ME	MW 129.1586 Quinoline Quinoline (SS) 500 µg/mL in Methanol(‡)	C ₉ H ₇ N	100mg 1ml	
Quinoline Yellow (E104)				
CAS 8004-92-0 DRE-C16709700	MW n/a Quinoline Yellow (E104)		250mg	No Structure
Quinoline-d7 (Quinoline-2,3,4,5,6,7,8-D7)				
CAS 34071-94-8 DRE-XA16709601AC	MW 136.2017 Quinoline-2,3,4,5,6,7,8-D7 100 µg/mL in Acetone	C ₉ ² H ₇ N	1.1ml	
Reactive Blue 19				
CAS 2580-78-1 DRE-C16809010	MW 626.5438 Reactive Blue 19	C ₂₂ H ₁₆ N ₂ O ₁₁ S ₃ ·2Na	100mg	
Rhodamine B Chloride				
CAS 81-88-9 DRE-C16813550	MW 479.0103 Rhodamine B chloride(‡)	C ₂₈ H ₃₁ N ₂ O ₃ ·Cl	250mg	
Solvent Green 3				
CAS 128-80-3 DRE-C16971100	MW 418.4865 Solvent Green 3	C ₂₈ H ₂₂ N ₂ O ₂	100mg	
Solvent Green 7				
CAS 6358-69-6 DRE-C16971110	MW 524.3851 Solvent Green 7	C ₁₆ H ₇ O ₁₀ S ₃ ·3Na	100mg	








Dyes and metabolites

Product code	Description			
Solvent Red 49				
CAS 509-34-2 DRE-C16971300	MW 442.5494 Solvent Red 49	$C_{28}H_{30}N_2O_3$	100mg	
Solvent Violet 13				
CAS 81-48-1 DRE-C16971330	MW 329.3487 Solvent Violet 13	$C_{21}H_{15}NO_3$	100mg	
Solvent Yellow 33				
CAS 8003-22-3 DRE-C16971350	MW n/a Solvent Yellow 33		100mg	No Structure
Solvent Yellow 34 Hydrochloride				
CAS 2465-27-2 DRE-C16971353	MW 303.8297 Solvent Yellow 34 Hydrochloride(‡)	$C_{17}H_{21}N_3 \cdot ClH$	100mg	
Sudan 1				
CAS 842-07-9 DRE-C16986101	MW 248.2793 Sudan 1(‡)	$C_{16}H_{12}N_2O$	100mg	
Sudan 1 D5 (phenyl D5)				
CAS 752211-63-5 DRE-C16986105 DRE-XA16986105AC	MW 253.3101 Sudan 1 D5 Sudan 1 D5 100 µg/mL in Acetone	$C_{16}^2H_{10}H_7N_2O$	10mg 1.1ml	
Sudan 2				
CAS 3118-97-6 DRE-C16986102	MW 276.3324 Sudan 2(‡)	$C_{18}H_{16}N_2O$	100mg	
Sudan 3				
CAS 85-86-9 DRE-C16986103	MW 352.3886 Sudan 3(‡)	$C_{22}H_{16}N_4O$	100mg	
Sudan 4				
CAS 85-83-6 DRE-C16986104	MW 380.4418 Sudan 4(‡)	$C_{24}H_{20}N_4O$	100mg	

Dyes and metabolites

Product code	Description			
Sudan 4 D6 (naphthyl D6)				
CAS 1014689-18-9	MW 386.4788	C ₂₄ H ₁₆ H ₁₄ N ₄ O		
DRE-C16986108	Sudan 4 D6(±)		10mg	
DRE-XA16986108AC	Sudan 4 D6 100 µg/mL in Acetone		1.1ml	
Sudan Black B				
CAS 4197-25-5	MW 456.5411	C ₂₉ H ₂₄ N ₆		
DRE-C16986110	Sudan Black B		100mg	
Sudan Blue 2				
CAS 17354-14-2	MW 350.454	C ₂₂ H ₂₆ N ₂ O ₂		
DRE-C16986113	Sudan Blue 2		100mg	
Sudan Orange G				
CAS 2051-85-6	MW 214.22	C ₁₂ H ₁₀ N ₂ O ₂		
DRE-C16986115	Sudan Orange G		100mg	
Sudan Red 7B				
CAS 6368-72-5	MW 379.457	C ₂₄ H ₂₁ N ₅		
DRE-C16986120	Sudan Red 7B		100mg	
Sudan Red B				
CAS 3176-79-2	MW 380.4418	C ₂₄ H ₂₀ N ₄ O		
DRE-C16986122	Sudan Red B		100mg	
Sudan Red G				
CAS 1229-55-6	MW 278.3053	C ₁₇ H ₁₄ N ₂ O ₂		
DRE-C16986127	Sudan Red G		100mg	
Sudan Yellow				
CAS 60-11-7	MW 225.289	C ₁₄ H ₁₅ N ₃		
DRE-C16986150	Sudan Yellow		250mg	
Sunset Yellow (E110)				
CAS 2783-94-0	MW 452.3693	C ₁₆ H ₁₀ N ₂ O ₇ S ₂ ·2Na		
DRE-C17048000	Sunset Yellow (E110)		50mg	

Dyes and metabolites

Product code	Description			
Tartrazine (E102)				
CAS 1934-21-0 DRE-C17138000	MW 534.3634 Tartrazine (E102)(‡)	C ₁₆ H ₉ NaO ₉ S ₂ ·3Na	250mg	
3,3',4,4'-Tetrachloroazobenzene				
CAS 14047-09-7 DRE-C17340000	MW 320.0014 3,3',4,4'-Tetrachloroazobenzene	C ₁₂ H ₆ Cl ₄ N ₂	10mg	
3,3',4,4'-Tetrachloroazoxybenzene				
CAS 21232-47-3 DRE-C17341000	MW 336.0008 3,3',4,4'-Tetrachloroazoxybenzene	C ₁₂ H ₆ Cl ₄ N ₂ O	100mg	
3,3',5,5'-Tetramethylbenzidine				
CAS 54827-17-7 DRE-C17413000 DRE-L17413000AL	MW 240.3434 3,3',5,5'-Tetramethylbenzidine 3,3',5,5'-Tetramethylbenzidine 10 µg/mL in Acetonitrile	C ₁₆ H ₂₀ N ₂	100mg 10ml	
Toluidine Red				
CAS 2425-85-6 DRE-C17597000	MW 307.3034 Toluidine Red	C ₁₇ H ₁₃ NaO ₃	25mg	
VAT Red 1				
CAS 2379-74-0 DRE-C17904000	MW 393.3068 VAT Red 1	C ₁₈ H ₁₀ Cl ₂ O ₂ S ₂	50mg	
Victoria Blue B				
CAS 2580-56-5 DRE-C17915000 DRE-E17915000	MW 506.0803 Victoria Blue B Victoria Blue B	C ₃₃ H ₃₂ N ₃ ·Cl	100mg 100mg	
Azodyes-Mix 1				
DRE-L18000079AL	Azodyes-Mix 1 10 µg/mL in Acetonitrile (*)			10ml
DRE-LA18000079AL	Azodyes-Mix 1 10 µg/mL in Acetonitrile (*)			1ml
2,4,5-Trimethylaniline	2,4-Diaminotoluene	2-Amino-4-nitrotoluene	2-Aminonaphthalene	
2-Methoxy-5-methylaniline	3,3'-Dichlorobenzidine	3,3'-Dimethoxybenzidine	3,3'-Dimethylbenzidine (o-Tolidine)	
4,4'-Benzidine	4,4'-Methylene-bis(o-toluidine)	4,4'-Methylene-bis(2-chloroaniline)	4-Amino-2',3-dimethylazobenzene	
4-Aminobiphenyl	4-Aminophenylether	4-Aminophenylthioether	4-Chloro-2-methylaniline	
4-Chloroaniline	4-Methoxy-1,3-phenylenediamine	Bis-(4-aminophenyl)methane	o-Toluidine	
Azodyes-Mix 6				
DRE-L18000376AL	Azodyes-Mix 6 10 µg/mL in Acetonitrile (*)			10ml
DRE-LA18000376AL	Azodyes-Mix 6 10 µg/mL in Acetonitrile (*)			1ml
2,4,5-Trimethylaniline	2,4-Diaminotoluene	2,4-Dimethylaniline	2,6-Dimethylaniline	
2-Amino-4-nitrotoluene	2-Aminonaphthalene	2-Anisidine	2-Methoxy-5-methylaniline	
3,3'-Dichlorobenzidine	3,3'-Dimethoxybenzidine	3,3'-Dimethylbenzidine (o-Tolidine)	4,4'-Benzidine	

Dyes and metabolites

Product code	Description			
4,4'-Methylene-bis(o-toluidine)	4,4'-Methylene-bis(2-chloroaniline)	4-Amino-2',3-dimethylazobenzene	4-Aminoazobenzene	
4-Aminobiphenyl	4-Aminophenylether	4-Aminophenylthioether	4-Chloro-2-methylaniline	
4-Chloroaniline	4-Methoxy-1,3-phenyldiamine	Bis-(4-aminophenyl)methane	o-Toluidine	

Azodyes-Mix 9

DRE-XA18000466AL	Azodyes-Mix 9 100 µg/mL in Acetonitrile (*)			1ml
2,4,5-Trimethylaniline	2,4-Diaminotoluene	2,4-Dimethylaniline	2,6-Dimethylaniline	
2-Amino-4-nitrotoluene	2-Aminonaphthalene	2-Anisidine	2-Methoxy-5-methylaniline	
3,3'-Dichlorobenzidine	3,3'-Dimethoxybenzidine	3,3'-Dimethylbenzidine (o-Tolidine)	4,4'-Benzidine	
4,4'-Methylene-bis(o-toluidine)	4,4'-Methylene-bis(2-chloroaniline)	4-Amino-2',3-dimethylazobenzene	4-Aminobiphenyl	
4-Aminophenylether	4-Aminophenylthioether	4-Chloro-2-methylaniline	4-Chloroaniline	
4-Methoxy-1,3-phenyldiamine	Bis-(4-aminophenyl)methane	o-Toluidine		

Azodyes-Mix 16

DRE-XA18000856AL	Azodyes-Mix 16 50 µg/mL in Acetonitrile (*)			1ml
2,4,5-Trimethylaniline	2,4-Diaminotoluene	2-Amino-4-nitrotoluene	2-Aminonaphthalene	
2-Anisidine	2-Methoxy-5-methylaniline	3,3'-Dichlorobenzidine	3,3'-Dimethoxybenzidine	
3,3'-Dimethylbenzidine (o-Tolidine)	4,4'-Benzidine	4,4'-Methylene-bis(o-toluidine)	4,4'-Methylene-bis(2-chloroaniline)	
4-Amino-2',3-dimethylazobenzene	4-Aminoazobenzene	4-Aminobiphenyl	4-Aminophenylether	
4-Aminophenylthioether	4-Chloro-2-methylaniline	4-Chloroaniline	4-Methoxy-1,3-phenyldiamine	
Bis-(4-aminophenyl)methane	o-Toluidine			

Benzidines Mix 2

DRE-YS09000047DI	Benzidines Mix 2 2000 µg/mL in Dichloromethane(‡)			5x1ml
	3,3'-dichlorobenzidine			
	o-tolidine			
	benzidine			










Benzidines Standard Solution

DRE-AY09000026ME	Benzidines Mix 1 2000 µg/mL in Methanol(‡)	5ml
DRE-SY09000026ME	Benzidines Mix 1 2000 µg/mL in Methanol(‡)	5x5ml
DRE-YA09000026ME	Benzidines Mix 1 2000 µg/mL in Methanol(‡)	1ml
DRE-YS09000026ME	Benzidines Mix 1 2000 µg/mL in Methanol	5x1ml
	3,3'-dichlorobenzidine	
	benzidine	










PHENOL AND AROMATIC COMPOUNDS











Phenol and aromatic compounds

Product code	Description			
2-Amino-4-chlorophenol				
CAS 95-85-2 DRE-C10199500	MW 143.5709 2-Amino-4-chlorophenol	C ₆ H ₆ ClNO	250mg	
2-Amino-4,6-dinitrotoluene (2-Methyl-3,5-dinitroaniline)				
CAS 35572-78-2 DRE-C10202200	MW 197.1482 2-Amino-4,6-dinitrotoluene	C ₇ H ₇ N ₃ O ₄	10mg	
4-Amino-2,6-dinitrotoluene (4-Methyl-3,5-dinitroaniline)				
CAS 19406-51-0 DRE-C10202300 DRE-LA10202300ME	MW 197.1482 4-Amino-2,6-dinitrotoluene 4-Amino-2,6-dinitrotoluene 10 µg/mL in Methanol	C ₇ H ₇ N ₃ O ₄	10mg 1ml	
4-Amino-3-fluorophenol				
CAS 399-95-1 DRE-C10202400	MW 127.1163 4-Amino-3-fluorophenol	C ₆ H ₆ FNO	100mg	
4-Amino-3-methylphenol				
CAS 2835-99-6 DRE-C10204955	MW 123.1525 4-Amino-3-methylphenol	C ₇ H ₉ NO	100mg	
2-Amino-4-nitrotoluene (2-Methyl-5-nitroaniline)				
CAS 99-55-8 DRE-L10207500AL	MW 152.1506 2-Amino-4-nitrotoluene (2-Methyl-5-nitroaniline) 10 µg/mL in Acetonitrile	C ₇ H ₈ N ₂ O ₂	10ml	
3-Amino-2-nitrotoluene (3-Methyl-2-nitroaniline)				
CAS 601-87-6 DRE-L10207800AL	MW 152.1506 3-Amino-2-nitrotoluene (3-Methyl-2-nitroaniline) 10 µg/mL in Acetonitrile	C ₇ H ₈ N ₂ O ₂	10ml	
2-Aminophenol				
CAS 95-55-6 DRE-C10210000	MW 109.1259 2-Aminophenol	C ₆ H ₇ NO	500mg	
4-Aminophenol				
CAS 123-30-8 DRE-C10212000	MW 109.1259 4-Aminophenol(‡)	C ₆ H ₇ NO	500mg	










Phenol and aromatic compounds

Product code	Description			
4-tert-Amylphenol				
CAS 80-46-6 DRE-C10247000	MW 164.2441 4-tert-Amylphenol	C ₁₁ H ₁₆ O	250mg	
Benzaldehyd-2,4-dinitrophenylhydrazone				
CAS 1157-84-2 DRE-C10532010	MW 286.2429 Benzaldehyd-2,4-dinitrophenylhydrazone	C ₁₃ H ₁₀ N ₄ O ₄	100mg	
Benzaldehyde				
CAS 100-52-7 DRE-CA10532000 DRE-YS09010014DI	MW 106.1219 Benzaldehyde Benzaldehyde 2000 µg/mL in Dichloromethane(±)	C ₇ H ₆ O	1ml 5x1ml	
Benzene-1,2-diamine (1,2-Phenylenediamine)				
CAS 95-54-5 DRE-CA16057800	MW 108.1411 1,2-Phenylenediamine(±)	C ₆ H ₈ N ₂	100mg	
Benzene-1,2-diol (Catechol; 1,2-Dihydroxybenzene)				
CAS 120-80-9 DRE-C11060000 DRE-YS09010019DI	MW 110.1106 Catechol 1,2-Dihydroxybenzene (Catechol) 1000 µg/mL in Dichloromethane(±)	C ₆ H ₆ O ₂	500mg 5x1ml	
Benzylchloride (α-Chlorotoluene)				
CAS 100-44-7 DRE-C11519000	MW 126.5835 α-Chlorotoluene (Benzylchloride)	C ₇ H ₇ Cl	1g	
2-Benzyl-4-chlorophenol				
CAS 120-32-1 DRE-C10572000	MW 218.6788 2-Benzyl-4-chlorophenol	C ₁₃ H ₁₁ ClO	100mg	
Bromobenzene				
CAS 108-86-1 DRE-CA10710000 DRE-XA10710000ME	MW 157.0079 Bromobenzene Bromobenzene 100 µg/mL in Methanol	C ₆ H ₅ Br	1ml 1ml	
4-Bromo-3-chloroaniline				
CAS 21402-26-6 DRE-C10716000	MW 206.4676 4-Bromo-3-chloroaniline	C ₆ H ₅ BrClN	100mg	










Phenol and aromatic compounds

Product code	Description			
4-Bromofluorobenzene (BFB)				
CAS 460-00-4	MW 174.9984	C ₆ H ₄ BrF		
DRE-C10731000	4-Bromofluorobenzene		250mg	
DRE-L10731000ME	4-Bromofluorobenzene 10 µg/mL in Methanol		10ml	
DRE-XA10731000ME	4-Bromofluorobenzene 100 µg/mL in Methanol		1ml	
DRE-YA10731000ME	4-Bromofluorobenzene 2000 µg/mL in Methanol		1ml	
DRE-GA09010387ME	4-Bromofluorobenzene (BFB) Solution, 2000 µg/mL in p/t Methanol(‡)		1ml	
DRE-GA09010388ME	4-Bromofluorobenzene (BFB) Solution, 2500 µg/mL in p/t Methanol(‡)		1ml	
DRE-GA09010389ME	4-Bromofluorobenzene (BFB) Solution, 10000 µg/mL in p/t Methanol(‡)		1ml	
1-Bromo-2-nitrobenzene				
CAS 577-19-5	MW 202.0055	C ₆ H ₄ BrNO ₂		
DRE-C10735500	1-Bromo-2-nitrobenzene		100mg	
DRE-GA10735500ME	1-Bromo-2-nitrobenzene 1000 µg/mL in Methanol		1ml	
DRE-GH10735500ME	1-Bromo-2-nitrobenzene 1000 µg/mL in Methanol		5x1ml	
2-Bromophenol				
CAS 95-56-7	MW 173.0073	C ₆ H ₅ BrO		
DRE-C10736200	2-Bromophenol		250mg	
DRE-XA10736200ME	2-Bromophenol 100 µg/mL in Methanol		1ml	
3-Bromophenol				
CAS 591-20-8	MW 173.0073	C ₆ H ₅ BrO		
DRE-C10736300	3-Bromophenol		250mg	
4-tert-Butyl-2-chlorophenol				
CAS 98-28-2	MW 184.6626	C ₁₀ H ₁₃ ClO		
DRE-C10931125	4-tert-Butyl-2-chlorophenol		100mg	
6-tert-Butyl-2,4-dimethylphenol (2-tert-Butyl-4,6-dimethylphenol)				
CAS 1879-09-0	MW 178.2707	C ₁₂ H ₁₈ O		
DRE-C10931188	2-tert-Butyl-4,6-dimethylphenol		100mg	
3-tert-Butyl-4-hydroxyanisole (2-tert-Butyl-4-methoxyphenol; BHA)				
CAS 121-00-6	MW 180.2435	C ₁₁ H ₁₆ O ₂		
DRE-C10931270	2-tert-Butyl-4-methoxyphenol (BHA)		250mg	
2-tert-Butylphenol				
CAS 88-18-6	MW 150.2176	C ₁₀ H ₁₄ O		
DRE-C10931400	2-tert-Butylphenol		100mg	










Phenol and aromatic compounds

Product code	Description			
4-Butylphenol				
CAS 1638-22-8 DRE-C10931350	MW 150.2176 4-n-Butylphenol	C ₁₀ H ₁₄ O	250mg	
4-tert-Butylphenol				
CAS 98-54-4 DRE-C10931600	MW 150.2176 4-tert.-Butylphenol	C ₁₀ H ₁₄ O	100mg	
2-Chloroaniline				
CAS 95-51-2 DRE-CA11350000 DRE-GA11350000DI DRE-GH11350000DI DRE-XA11350000AL	MW 127.5715 2-Chloroaniline 2-Chloroaniline 1000 µg/mL in Dichloromethane 2-Chloroaniline 1000 µg/mL in Dichloromethane 2-Chloroaniline 100 µg/mL in Acetonitrile	C ₆ H ₆ ClN	500mg 1ml 5x1ml 1ml	
3-Chloroaniline				
CAS 108-42-9 DRE-GA11351000TO DRE-GH11351000TO	MW 127.5715 3-Chloroaniline 1000 µg/mL in Toluene 3-Chloroaniline 1000 µg/mL in Toluene	C ₆ H ₆ ClN	1ml 5x1ml	
2-Chloroanisole				
CAS 766-51-8 DRE-C11360000	MW 142.5829 2-Chloroanisole	C ₇ H ₇ ClO	500mg	
3-Chloroanisole				
CAS 2845-89-8 DRE-C11361000	MW 142.5829 3-Chloroanisole	C ₇ H ₇ ClO	500mg	
4-Chloroanisole				
CAS 623-12-1 DRE-C11362000	MW 142.5829 4-Chloroanisole	C ₇ H ₇ ClO	500mg	
Chlorobenzene				
CAS 108-90-7 DRE-C11380000 DRE-L11380000IO DRE-XA11380000ME	MW 112.5569 Chlorobenzene(‡) Chlorobenzene 10 µg/mL in Isooctane Chlorobenzene 100 µg/mL in Methanol	C ₆ H ₅ Cl	1ml 10ml 1ml	
Chlorobenzene D5				
CAS 3114-55-4 DRE-C11380100	MW 117.5877 Chlorobenzene D5	C ₆ ² H ₅ Cl	100mg	










Phenol and aromatic compounds

Product code	Description			
4-Chlorobenzotrifluoride (4-Chloro- α,α,α -trifluorotoluene)				
CAS 98-56-6 DRE-C11536400	MW 180.5549 4-Chloro- α,α,α -trifluorotoluene	$C_7H_4ClF_3$	100mg	
4-Chlorocatechol				
CAS 2138-22-9 DRE-C11395300	MW 144.5557 4-Chlorocatechol	$C_6H_5ClO_2$	100mg	
6-Chloro-2,4-dinitroaniline				
CAS 3531-19-9 DRE-C11405200	MW 217.5667 6-Chloro-2,4-dinitroaniline	$C_6H_4ClN_2O_4$	100mg	
1-Chloro-2,4-dinitrobenzene				
CAS 97-00-7 DRE-C11405500 DRE-XA11405500ME	MW 202.552 1-Chloro-2,4-dinitrobenzene 1-Chloro-2,4-dinitrobenzene 100 μ g/mL in Methanol	$C_6H_3ClN_2O_4$	500mg 1ml	
3-Chloro-4-fluoroaniline				
CAS 367-21-5 DRE-C11415200	MW 145.562 3-Chloro-4-fluoroaniline	C_6H_6ClFN	250mg	
4-Chloroguaiacol				
CAS 16766-30-6 DRE-C11415800	MW 158.5823 4-Chloroguaiacol	$C_7H_7ClO_2$	100mg	
4-Chloro-2-isopropyl-5-methylphenol				
CAS 89-68-9 DRE-C11418000	MW 184.6626 4-Chloro-2-isopropyl-5-methylphenol	$C_{10}H_{13}ClO$	250mg	
2-Chloro-4-methylaniline				
CAS 615-65-6 DRE-C11429400	MW 141.5981 2-Chloro-4-methylaniline	C_7H_8ClN	250mg	
2-Chloro-5-methylaniline				
CAS 95-81-8 DRE-C11429500	MW 141.5981 2-Chloro-5-methylaniline	C_7H_8ClN	250mg	










Phenol and aromatic compounds

Product code	Description			
2-Chloro-6-methylaniline				
CAS 87-63-8 DRE-C11429600	MW 141.5981 2-Chloro-6-methylaniline	C ₇ H ₈ ClN	250mg	
3-Chloro-2-methylaniline				
CAS 87-60-5 DRE-C11429700	MW 141.5981 3-Chloro-2-methylaniline	C ₇ H ₈ ClN	250mg	
5-Chloro-2-methylaniline				
CAS 95-79-4 DRE-C11430500	MW 141.5981 5-Chloro-2-methylaniline	C ₇ H ₈ ClN	250mg	
2-Chloro-4-methylphenol				
CAS 6640-27-3 DRE-C11439200	MW 142.5829 2-Chloro-4-methylphenol	C ₇ H ₇ ClO	100mg	
2-Chloro-5-methylphenol				
CAS 615-74-7 DRE-C11439300	MW 142.5829 2-Chloro-5-methylphenol	C ₇ H ₇ ClO	250mg	
2-Chloro-6-methylphenol (6-Chloro-o-cresol)				
CAS 87-64-9 DRE-C11439400	MW 142.5829 2-Chloro-6-methylphenol	C ₇ H ₇ ClO	100mg	
4-Chloro-3-methylphenol (Chlorocresol)				
CAS 59-50-7 DRE-C11440300 DRE-XA11440300ME	MW 142.5829 4-Chloro-3-methylphenol 4-Chloro-3-methylphenol 100 µg/mL in Methanol	C ₇ H ₇ ClO	250mg 1ml	
2-Chloro-4-nitroaniline				
CAS 121-87-9 DRE-C11452800	MW 172.5691 2-Chloro-4-nitroaniline	C ₆ H ₅ ClN ₂ O ₂	500mg	
4-Chloro-2-nitroaniline				
CAS 89-63-4 DRE-C11453000 DRE-L11453000CY	MW 172.5691 4-Chloro-2-nitroaniline 4-Chloro-2-nitroaniline 10 µg/mL in Cyclohexane	C ₆ H ₅ ClN ₂ O ₂	500mg 10ml	










Phenol and aromatic compounds

Product code	Description				
1-Chloro-2-nitrobenzene					
CAS 88-73-3	MW 157.5545	C ₆ H ₄ ClNO ₂			
DRE-C11453500	1-Chloro-2-nitrobenzene		250mg		
DRE-L11453500ME	1-Chloro-2-nitrobenzene 10 µg/mL in Methanol		10ml		
DRE-XA11453500ME	1-Chloro-2-nitrobenzene 100 µg/mL in Methanol		1ml		
1-Chloro-3-nitrobenzene					
CAS 121-73-3	MW 157.5545	C ₆ H ₄ ClNO ₂			
DRE-C11453600	1-Chloro-3-nitrobenzene		250mg		
DRE-XA11453600ME	1-Chloro-3-nitrobenzene 100 µg/mL in Methanol		1ml		
1-Chloro-4-nitrobenzene					
CAS 100-00-5	MW 157.5545	C ₆ H ₄ ClNO ₂			
DRE-C11453700	1-Chloro-4-nitrobenzene(†)		250mg		
DRE-XA11453700ME	1-Chloro-4-nitrobenzene 100 µg/mL in Methanol		1ml		
2-Chloro-4-nitrotoluene					
CAS 121-86-8	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456100	2-Chloro-4-nitrotoluene		250mg		
2-Chloro-5-nitrotoluene					
CAS 13290-74-9	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456200	2-Chloro-5-nitrotoluene		100mg		
2-Chloro-6-nitrotoluene					
CAS 83-42-1	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456300	2-Chloro-6-nitrotoluene		250mg		
3-Chloro-4-nitrotoluene					
CAS 38939-88-7	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456500	3-Chloro-4-nitrotoluene		100mg		
4-Chloro-2-nitrotoluene					
CAS 89-59-8	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456800	4-Chloro-2-nitrotoluene		250mg		
4-Chloro-3-nitrotoluene					
CAS 89-60-1	MW 171.581	C ₇ H ₆ ClNO ₂			
DRE-C11456900	4-Chloro-3-nitrotoluene		250mg		










Phenol and aromatic compounds

Product code	Description			
5-Chloro-2-nitrotoluene				
CAS 5367-28-2 DRE-C11456950	MW 171.581 5-Chloro-2-nitrotoluene	C ₇ H ₅ ClNO ₂	100mg	
2-Chlorophenol-3,4,5,6-D4				
CAS 93951-73-6 DRE-C11470100 DRE-XA11470100AC	MW 132.5809 2-Chlorophenol-3,4,5,6-D4 2-Chlorophenol D4 100 µg/mL in Acetone	C ₆ H ₄ HClO	10mg 1.1ml	
2-Chlorophenol				
CAS 95-57-8 DRE-C11470000 DRE-L11470000ME DRE-XA11470000ME	MW 128.5563 2-Chlorophenol 2-Chlorophenol 10 µg/mL in Methanol 2-Chlorophenol 100 µg/mL in Methanol	C ₆ H ₅ ClO	500mg 10ml 1ml	
3-Chlorophenol				
CAS 108-43-0 DRE-C11471000 DRE-XA11471000ME	MW 128.5563 3-Chlorophenol(‡) 3-Chlorophenol 100 µg/mL in Methanol	C ₆ H ₅ ClO	500mg 1ml	
4-Chlorophenol				
CAS 106-48-9 DRE-XA11472000ME	MW 128.5563 4-Chlorophenol 100 µg/mL in Methanol	C ₆ H ₅ ClO	1ml	
(2-Chlorophenyl)trichloromethane (α,α,α,2-Tetrachlorotoluene)				
CAS 2136-89-2 DRE-C17379800	MW 229.9187 α,α,α,2-Tetrachlorotoluene	C ₇ H ₄ Cl ₄	250mg	
2-Chlorotoluene				
CAS 95-49-8 DRE-C11520000	MW 126.5835 2-Chlorotoluene(‡)	C ₇ H ₇ Cl	500mg	
3-Chlorotoluene				
CAS 108-41-8 DRE-C11521000 DRE-XA11521000ME	MW 126.5835 3-Chlorotoluene 3-Chlorotoluene 100 µg/mL in Methanol	C ₇ H ₇ Cl	500mg 1ml	
4-Chlorotoluene				
CAS 106-43-4 DRE-C11522000 DRE-XA11522000ME	MW 126.5835 4-Chlorotoluene 4-Chlorotoluene 100 µg/mL in Methanol	C ₇ H ₇ Cl	500mg 1ml	










Phenol and aromatic compounds

Product code	Description			
Chloroxylenol (4-Chloro-3,5-dimethylphenol)				
CAS 88-04-0 DRE-C11405100	MW 156.6095 4-Chloro-3,5-dimethylphenol	C ₈ H ₉ ClO	250mg	
2-Cyclohexyl-4,6-dinitrophenol				
CAS 131-89-5 DRE-C11830600	MW 266.25 2-Cyclohexyl-4,6-dinitrophenol	C ₁₂ H ₁₄ N ₂ O ₅	20mg	
2,4-Diaminoanisole Dihydrochloride				
CAS 614-94-8 DRE-C12192000	MW 211.089 2,4-Diaminoanisole dihydrochloride	C ₇ H ₁₀ N ₂ O · 2ClH	100mg	
1,4-Diamino-2-nitrobenzene				
CAS 5307-14-2 DRE-C12195500	MW 153.1387 1,4-Diamino-2-nitrobenzene(†)	C ₆ H ₇ N ₃ O ₂	100mg	
2,6-Diamino-4-nitrotoluene				
CAS 59229-75-3 DRE-LA12195800AL	MW 167.1653 2,6-Diamino-4-nitrotoluene 10 µg/mL in Acetonitrile	C ₇ H ₉ N ₃ O ₂	1ml	
2,6-Diaminopyridine				
CAS 141-86-6 DRE-CA12197300	MW 109.1292 2,6-Diaminopyridine	C ₅ H ₇ N ₃	100mg	
2,5-Diaminotoluene sulfate (2-Methyl-p-phenylenediamine sulfate)				
CAS 615-50-9 DRE-C12198500	MW 220.2462 2,5-Diaminotoluene sulfate (2-Methyl-p-phenylenediamine sulfate)	C ₇ H ₁₀ N ₂ · H ₂ O ₄ S	250mg	
2,3-Diaminotoluene				
CAS 2687-25-4 DRE-C12197500	MW 122.1677 2,3-Diaminotoluene	C ₇ H ₁₀ N ₂	100mg	
2,4-Diaminotoluene				
CAS 95-80-7 DRE-C12197600	MW 122.1677 2,4-Diaminotoluene	C ₇ H ₁₀ N ₂	250mg	










Phenol and aromatic compounds

Product code	Description			
2,6-Diaminotoluene				
CAS 823-40-5 DRE-C12197800	MW 122.1677 2,6-Diaminotoluene	C ₇ H ₁₀ N ₂	250mg	
3,4-Diaminotoluene				
CAS 496-72-0 DRE-C12197900	MW 122.1677 3,4-Diaminotoluene	C ₇ H ₁₀ N ₂	250mg	
1,2-Dibromobenzene				
CAS 583-53-9 DRE-C12216710	MW 235.904 1,2-Dibromobenzene	C ₆ H ₄ Br ₂	100mg	
1,3-Dibromobenzene				
CAS 108-36-1 DRE-C12216711	MW 235.904 1,3-Dibromobenzene	C ₆ H ₄ Br ₂	100mg	
1,4-Dibromobenzene				
CAS 106-37-6 DRE-C12216712	MW 235.904 1,4-Dibromobenzene	C ₆ H ₄ Br ₂	100mg	
2,4-Dibromophenol				
CAS 615-58-7 DRE-C12241000 DRE-L12241000ME DRE-XA12241000ME	MW 251.9034 2,4-Dibromophenol 2,4-Dibromophenol 10 µg/mL in Methanol 2,4-Dibromophenol 100 µg/mL in Methanol	C ₆ H ₄ Br ₂ O	250mg 10ml 1ml	
2,6-Dibromophenol				
CAS 608-33-3 DRE-C12241200	MW 251.9034 2,6-Dibromophenol	C ₆ H ₄ Br ₂ O	100mg	
1,2-Dibromopropane				
CAS 78-75-1 DRE-GA12241900HE DRE-GH12241900HE	MW 201.8877 1,2-Dibromopropane 10000 µg/mL in n-Hexane 1,2-Dibromopropane 10000 µg/mL in Hexane(‡)	C ₃ H ₆ Br ₂	1ml 5x1ml	
2,5-Di-tert-butylhydroquinone				
CAS 88-58-4 DRE-C12252500	MW 222.3233 2,5-Di-tert-butylhydroquinone	C ₁₄ H ₂₂ O ₂	250mg	










Phenol and aromatic compounds

Product code	Description			
2,6-Di-tert-butyl-4-methoxyphenol				
CAS 489-01-0 DRE-C12253400	MW 236.3499 2,6-Di-tert-butyl-4-methoxyphenol	C ₁₉ H ₂₄ O ₂	100mg	
4,6-Di-tert-butyl-2-methylphenol				
CAS 616-55-7 DRE-C12253525	MW 220.3505 4,6-Di-tert-butyl-2-methylphenol	C ₁₅ H ₂₄ O	100mg	
2,4-Di-tert-butylphenol				
CAS 96-76-4 DRE-C12254600	MW 206.3239 2,4-Di-tert-butylphenol	C ₁₄ H ₂₂ O	100mg	
2,6-Di-tert-butylphenol				
CAS 128-39-2 DRE-C12254700	MW 206.3239 2,6-Di-tert-butylphenol	C ₁₄ H ₂₂ O	100mg	
3,5-Di-tert-butylphenol				
CAS 1138-52-9 DRE-C12254800	MW 206.3239 3,5-Di-tert-butylphenol	C ₁₄ H ₂₂ O	10mg	
2,3-Dichloroaniline				
CAS 608-27-5 DRE-C12322300 DRE-XA12322300ME	MW 162.0166 2,3-Dichloroaniline 2,3-Dichloroaniline 100 µg/mL in Methanol	C ₆ H ₅ Cl ₂ N	500mg 1ml	
2,4-Dichloroaniline				
CAS 554-00-7 DRE-C12322400 DRE-XA12322400ME	MW 162.0166 2,4-Dichloroaniline 2,4-Dichloroaniline 100 µg/mL in Methanol	C ₆ H ₅ Cl ₂ N	500mg 1ml	
2,6-Dichloroaniline				
CAS 608-31-1 DRE-C12322600 DRE-XA12322600ME	MW 162.0166 2,6-Dichloroaniline 2,6-Dichloroaniline 100 µg/mL in Methanol	C ₆ H ₅ Cl ₂ N	500mg 1ml	
2,3-Dichloroanisole				
CAS 1984-59-4 DRE-C12332300	MW 177.0279 2,3-Dichloroanisole	C ₇ H ₆ Cl ₂ O	500mg	










Phenol and aromatic compounds

Product code	Description			
2,4-Dichloroanisole				
CAS 553-82-2 DRE-C12332400	MW 177.0279 2,4-Dichloroanisole	C ₇ H ₆ Cl ₂ O	100mg	
2,6-Dichloroanisole				
CAS 1984-65-2 DRE-C12332600	MW 177.0279 2,6-Dichloroanisole	C ₇ H ₆ Cl ₂ O	500mg	
3,4-Dichloroanisole				
CAS 36404-30-5 DRE-C12333400	MW 177.0279 3,4-Dichloroanisole	C ₇ H ₆ Cl ₂ O	100mg	
3,5-Dichloroanisole				
CAS 33719-74-3 DRE-C12333500	MW 177.0279 3,5-Dichloroanisole	C ₇ H ₆ Cl ₂ O	500mg	
1,2-Dichlorobenzene				
CAS 95-50-1 DRE-L12370000CY DRE-XA12370000ME	MW 147.002 1,2-Dichlorobenzene 10 µg/mL in Cyclohexane 1,2-Dichlorobenzene 100 µg/mL in Methanol	C ₆ H ₄ Cl ₂	10ml 1ml	
1,2-Dichlorobenzene D4				
CAS 2199-69-1 DRE-CA12370100 DRE-YA12370100ME	MW 151.0266 1,2-Dichlorobenzene D4 1,2-Dichlorobenzene D4 2000 µg/mL in Methanol	C ₆ ² H ₄ Cl ₂	100mg 1ml	
1,3-Dichlorobenzene				
CAS 541-73-1 DRE-L12371000CY DRE-XA12371000ME	MW 147.002 1,3-Dichlorobenzene 10 µg/mL in Cyclohexane 1,3-Dichlorobenzene 100 µg/mL in Methanol	C ₆ H ₄ Cl ₂	10ml 1ml	
1,4-Dichlorobenzene				
CAS 106-46-7 DRE-L12372000IO DRE-XA12372000ME DRE-YA12372000ME	MW 147.002 1,4-Dichlorobenzene 10 µg/mL in Isooctane 1,4-Dichlorobenzene 100 µg/mL in Methanol 1,4-Dichlorobenzene 1000 µg/mL in Methanol	C ₆ H ₄ Cl ₂	10ml 1ml 1ml	
1,4-Dichlorobenzene D4				
CAS 3855-82-1 DRE-YA12372100ME	MW 151.0266 1,4-Dichlorobenzene D4 2000 µg/mL in Methanol	C ₆ ² H ₄ Cl ₂	1ml	










Phenol and aromatic compounds

Product code	Description			
3,5-Dichlorocatechol				
CAS 13673-92-2 DRE-C12420805	MW 179.0008 3,5-Dichlorocatechol	$C_6H_4Cl_2O_2$	100mg	
4,5-Dichlorocatechol				
CAS 3428-24-8 DRE-C12420820	MW 179.0008 4,5-Dichlorocatechol	$C_6H_4Cl_2O_2$	100mg	
2,4-Dichloro-6-methylphenol				
CAS 1570-65-6 DRE-L12427000ME	MW 177.0279 2,4-Dichloro-6-methylphenol 10 µg/mL in Methanol	$C_7H_6Cl_2O$	10ml	
2,6-Dichloro-4-methylphenol				
CAS 2432-12-4 DRE-C12427100	MW 177.0279 2,6-Dichloro-4-methylphenol	$C_7H_6Cl_2O$	100mg	
2,3-Dichloronitrobenzene				
CAS 3209-22-1 DRE-C12432300	MW 191.9995 2,3-Dichloronitrobenzene	$C_6H_3Cl_2NO_2$	250mg	
2,4-Dichloronitrobenzene				
CAS 611-06-3 DRE-C12432400	MW 191.9995 2,4-Dichloronitrobenzene	$C_6H_3Cl_2NO_2$	250mg	
2,5-Dichloronitrobenzene				
CAS 89-61-2 DRE-C12432500	MW 191.9995 2,5-Dichloronitrobenzene	$C_6H_3Cl_2NO_2$	250mg	
3,4-Dichloronitrobenzene				
CAS 99-54-7 DRE-C12433400	MW 191.9995 3,4-Dichloronitrobenzene	$C_6H_3Cl_2NO_2$	250mg	
3,5-Dichloronitrobenzene				
CAS 618-62-2 DRE-C12433500	MW 191.9995 3,5-Dichloronitrobenzene	$C_6H_3Cl_2NO_2$	250mg	










Phenol and aromatic compounds

Product code	Description			
2,3-Dichlorophenol				
CAS 576-24-9	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12450000	2,3-Dichlorophenol		250mg	
DRE-L12450000ME	2,3-Dichlorophenol 10 µg/mL in Methanol		10ml	
2,5-Dichlorophenol				
CAS 583-78-8	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12452000	2,5-Dichlorophenol(‡)		250mg	
DRE-L12452000ME	2,5-Dichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA12452000ME	2,5-Dichlorophenol 100 µg/mL in Methanol		1ml	
DRE-GA09010066ME	2,5-Dichlorophenol 100 µg/mL in Methanol(‡)		1ml	
2,6-Dichlorophenol				
CAS 87-65-0	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12453000	2,6-Dichlorophenol		250mg	
DRE-L12453000ME	2,6-Dichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA12453000ME	2,6-Dichlorophenol 100 µg/mL in Methanol		1ml	
DRE-YS09010009ME	2,6-Dichlorophenol 1000 µg/mL in Methanol(‡)		5x1ml	
3,4-Dichlorophenol				
CAS 95-77-2	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12453400	3,4-Dichlorophenol		250mg	
DRE-L12453400ME	3,4-Dichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA12453400ME	3,4-Dichlorophenol 100 µg/mL in Methanol		1ml	
3,5-Dichlorophenol				
CAS 591-35-5	MW 163.0014	C ₆ H ₄ Cl ₂ O		
DRE-C12453500	3,5-Dichlorophenol		250mg	
DRE-L12453500ME	3,5-Dichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA12453500ME	3,5-Dichlorophenol 100 µg/mL in Methanol		1ml	
2,3-Dichlorotoluene				
CAS 32768-54-0	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-C12505600	2,3-Dichlorotoluene		250mg	
DRE-XA12505600ME	2,3-Dichlorotoluene 100 µg/mL in Methanol		1ml	
2,4-Dichlorotoluene				
CAS 95-73-8	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-C12505700	2,4-Dichlorotoluene		250mg	
2,5-Dichlorotoluene				
CAS 19398-61-9	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-C12505800	2,5-Dichlorotoluene(‡)		250mg	
DRE-XA12505800ME	2,5-Dichlorotoluene 100 µg/mL in Methanol		1ml	
2,6-Dichlorotoluene				
CAS 118-69-4	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-C12505900	2,6-Dichlorotoluene		250mg	
DRE-XA12505900ME	2,6-Dichlorotoluene 100 µg/mL in Methanol		1ml	










Phenol and aromatic compounds

Product code	Description			
3,4-Dichlorotoluene				
CAS 95-75-0	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-C12506000	3,4-Dichlorotoluene		250mg	
DRE-XA12506000ME	3,4-Dichlorotoluene 100 µg/mL in Methanol		1ml	
α,α-Dichlorotoluene				
CAS 98-87-3	MW 161.0285	C ₇ H ₆ Cl ₂		
DRE-CA12505000	α,α-Dichlorotoluene		250mg	
Dichloroxylene (2,4-Dichloro-3,5-dimethylphenol)				
CAS 133-53-9	MW 191.0545	C ₈ H ₆ Cl ₂ O		
DRE-C12421500	2,4-Dichloro-3,5-dimethylphenol		100mg	
4-N,N-Diethyl-2-methyl-p-phenyldiamine monohydrochloride				
CAS 2051-79-8	MW 214.735	C ₁₁ H ₁₆ N ₂ ·ClH		
DRE-CA12606650	4-N,N-Diethyl-2-methyl-p-phenyldiamine monohydrochloride		250mg	
1,4-Difluorobenzene				
CAS 540-36-3	MW 114.0928	C ₆ H ₄ F ₂		
DRE-CA12632000	1,4-Difluorobenzene		100mg	
DRE-YA12632000ME	1,4-Difluorobenzene 2000 µg/mL in Methanol		1ml	
2,6-Dimethylaniline				
CAS 87-62-7	MW 121.1796	C ₈ H ₁₁ N		
DRE-C12725000	2,6-Dimethylaniline(‡)		1g	
DRE-XA12725000IO	2,6-Dimethylaniline 100 µg/mL in Isooctane		1ml	
N,N-Dimethylaniline				
CAS 121-69-7	MW 121.1796	C ₈ H ₁₁ N		
DRE-C12724300	N,N-Dimethylaniline		250mg	
2,5-Dimethylbenzaldehyde-2,4-dinitrophenylhydrazone				
CAS 152477-96-8	MW 314.2961	C ₁₅ H ₁₄ N ₄ O ₄		
DRE-XA12725510AL	2,5-Dimethylbenzaldehyd-2,4-dinitrophenylhydrazone 100 µg/mL in Acetonitrile		1ml	
1,2-Dimethyl-3-nitrobenzene				
CAS 83-41-0	MW 151.1626	C ₈ H ₉ NO ₂		
DRE-C12727900	1,2-Dimethyl-3-nitrobenzene		250mg	










Phenol and aromatic compounds

Product code	Description			
1,2-Dimethyl-4-nitrobenzene				
CAS 99-51-4	MW 151.1626	C ₈ H ₉ NO ₂		
DRE-C12727920	1,2-Dimethyl-4-nitrobenzene		250mg	
1,3-Dimethyl-2-nitrobenzene (2-Nitro-m-xylene)				
CAS 81-20-9	MW 151.1626	C ₈ H ₉ NO ₂		
DRE-C12728000	1,3-Dimethyl-2-nitrobenzene (2-Nitro-m-xylene)		250mg	
DRE-XA12728000ME	1,3-Dimethyl-2-nitrobenzene (2-Nitro-m-xylene) 100 µg/mL in Methanol		1ml	
DRE-GH09010380EA	1,3-Dimethyl-2-nitrobenzene 250 µg/mL in Ethyl Acetate(‡)		10x1ml	
2,3-Dimethylphenol				
CAS 526-75-0	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12730000	2,3-Dimethylphenol(‡)		500mg	
DRE-XA12730000ME	2,3-Dimethylphenol 100 µg/mL in Methanol		1ml	
2,4-Dimethylphenol				
CAS 105-67-9	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12731000	2,4-Dimethylphenol		500mg	
DRE-XA12731000ME	2,4-Dimethylphenol 100 µg/mL in Methanol		1ml	
2,5-Dimethylphenol				
CAS 95-87-4	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12732000	2,5-Dimethylphenol		500mg	
DRE-XA12732000ME	2,5-Dimethylphenol 100 µg/mL in Methanol		1ml	
2,6-Dimethylphenol				
CAS 576-26-1	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12733000	2,6-Dimethylphenol		500mg	
DRE-XA12733000ME	2,6-Dimethylphenol 100 µg/mL in Methanol		1ml	
3,4-Dimethylphenol				
CAS 95-65-8	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12734000	3,4-Dimethylphenol		500mg	
DRE-XA12734000ME	3,4-Dimethylphenol 100 µg/mL in Methanol		1ml	
3,5-Dimethylphenol				
CAS 108-68-9	MW 122.1644	C ₈ H ₁₀ O		
DRE-C12735000	3,5-Dimethylphenol		500mg	
DRE-XA12735000ME	3,5-Dimethylphenol 100 µg/mL in Methanol		1ml	
N,N-Dimethyl-p-phenyldiamine				
CAS 99-98-9	MW 136.1943	C ₈ H ₁₂ N ₂		
DRE-CA12736000	N,N-Dimethyl-p-phenyldiamine		100mg	









Phenol and aromatic compounds

Product code	Description			
2,4-Dimethyl-3,5,6-trideuteriophenol				
CAS 93951-75-8	MW 125.1829	C ₆ ² H ₃ H ₇ O		
DRE-C12731100	2,4-Dimethylphenol-3,5,6-D3		100mg	
DRE-XA12731100AC	2,4-Dimethylphenol-3,5,6 D3 100 µg/mL in Acetone		1.1ml	
Dinitolmide (3,5-Dinitro-o-toluamide)				
CAS 148-01-6	MW 225.1583	C ₈ H ₇ N ₃ O ₅		
DRE-C12786050	3,5-Dinitro-o-toluamide		250mg	
2,4-Dinitroaniline				
CAS 97-02-9	MW 183.1216	C ₆ H ₅ N ₃ O ₄		
DRE-C12782000	2,4-Dinitroaniline		100mg	
2,6-Dinitroaniline				
CAS 606-22-4	MW 183.1216	C ₆ H ₅ N ₃ O ₄		
DRE-C12782100	2,6-Dinitroaniline		100mg	
3,5-Dinitroaniline				
CAS 618-87-1	MW 183.1216	C ₆ H ₅ N ₃ O ₄		
DRE-C12782200	3,5-Dinitroaniline		100mg	
2,4-Dinitroanisole				
CAS 119-27-7	MW 198.1329	C ₇ H ₆ N ₂ O ₅		
DRE-C12782400	2,4-Dinitroanisole		100mg	
1,2-Dinitrobenzene				
CAS 528-29-0	MW 168.107	C ₆ H ₄ N ₂ O ₄		
DRE-C12783000	1,2-Dinitrobenzene(±)		250mg	
DRE-GA12783000ME	1,2-Dinitrobenzene 1000 µg/mL in Methanol		1ml	
DRE-GH12783000ME	1,2-Dinitrobenzene 1000 µg/mL in Methanol		5x1ml	
DRE-L12783000ME	1,2-Dinitrobenzene 10 µg/mL in Methanol		10ml	
1,3-Dinitrobenzene				
CAS 99-65-0	MW 168.107	C ₆ H ₄ N ₂ O ₄		
DRE-C12783100	1,3-Dinitrobenzene		250mg	
DRE-L12783100ME	1,3-Dinitrobenzene 10 µg/mL in Methanol		10ml	
1,4-Dinitrobenzene				
CAS 100-25-4	MW 168.107	C ₆ H ₄ N ₂ O ₄		
DRE-C12783200	1,4-Dinitrobenzene		250mg	
DRE-L12783200ME	1,4-Dinitrobenzene 10 µg/mL in Methanol		10ml	










Phenol and aromatic compounds

Product code	Description			
2,4-Dinitrobenzoic Acid				
CAS 610-30-0 DRE-C12783400	MW 212.1165 2,4-Dinitrobenzoic acid	C ₇ H ₄ N ₂ O ₆	250mg	
3,5-Dinitrobenzoic Acid				
CAS 99-34-3 DRE-C12783600	MW 212.1165 3,5-Dinitrobenzoic acid	C ₇ H ₄ N ₂ O ₆	250mg	
2,4-Dinitrodiphenylamine				
CAS 961-68-2 DRE-C12783702	MW 259.2176 2,4-Dinitrodiphenylamine	C ₁₂ H ₉ N ₃ O ₄	100mg	
2,6-Dinitro-4-methylphenol				
CAS 609-93-8 DRE-C12784000	MW 198.1329 2,6-Dinitro-4-methylphenol (moistened)	C ₇ H ₆ N ₂ O ₅	250mg	
2,6-Dinitrophenol				
CAS 573-56-8 DRE-L12785400AL	MW 184.1064 2,6-Dinitrophenol 10 µg/mL in Acetonitrile	C ₆ H ₄ N ₂ O ₅	10ml	
2,3-Dinitrotoluene				
CAS 602-01-7 DRE-C12786100	MW 182.1335 2,3-Dinitrotoluene	C ₇ H ₆ N ₂ O ₄	250mg	
2,4-Dinitrotoluene				
CAS 121-14-2 DRE-C12786200 DRE-XA12786200AL	MW 182.1335 2,4-Dinitrotoluene 2,4-Dinitrotoluene 100 µg/mL in Acetonitrile	C ₇ H ₆ N ₂ O ₄	250mg 1ml	
2,6-Dinitrotoluene				
CAS 606-20-2 DRE-C12786400 DRE-L12786400AL	MW 182.1335 2,6-Dinitrotoluene 2,6-Dinitrotoluene 10 µg/mL in Acetonitrile	C ₇ H ₆ N ₂ O ₄	250mg 10ml	
3,4-Dinitrotoluene				
CAS 610-39-9 DRE-C12786500 DRE-XA12786500AL	MW 182.1335 3,4-Dinitrotoluene 3,4-Dinitrotoluene 100 µg/mL in Acetonitrile	C ₇ H ₆ N ₂ O ₄	100mg 1ml	










Phenol and aromatic compounds

Product code	Description			
3,5-Dinitrotoluene				
CAS 618-85-9 DRE-XA12786800AL	MW 182.1335 3,5-Dinitrotoluene 100 µg/mL in Acetonitrile	C ₇ H ₅ N ₂ O ₄	1ml	
4-Dodecylphenol, mixture of isomers				
CAS 27193-86-8 DRE-C13066000	MW n/a 4-Dodecylphenol, mixture of isomers		250mg	No Structure
2-Ethylphenol				
CAS 90-00-6 DRE-C13350000 DRE-XA13350000ME	MW 122.1644 2-Ethylphenol(‡) 2-Ethylphenol 100 µg/mL in Methanol	C ₈ H ₁₀ O	500mg 1ml	
3-Ethylphenol				
CAS 620-17-7 DRE-C13351000 DRE-XA13351000ME	MW 122.1644 3-Ethylphenol 3-Ethylphenol 100 µg/mL in Methanol	C ₈ H ₁₀ O	500mg 1ml	
4-Ethylphenol				
CAS 123-07-9 DRE-C13352000 DRE-L13352000AL	MW 122.1644 4-Ethylphenol 4-Ethylphenol 10 µg/mL in Acetonitrile	C ₈ H ₁₀ O	500mg 10ml	
Fluorobenzene				
CAS 462-06-6 DRE-C13781000 DRE-L13781000ME DRE-YA13781000ME	MW 96.1023 Fluorobenzene(‡) Fluorobenzene 10 µg/mL in Methanol Fluorobenzene 2000 µg/mL in Methanol(‡)	C ₆ H ₅ F	1ml 10ml 1ml	
2-Fluorophenol				
CAS 367-12-4 DRE-CA13797000 DRE-YA13797000MB	MW 112.1017 2-Fluorophenol 2-Fluorophenol 2000 µg/mL in Methyl-tert. butyl ether	C ₆ H ₅ FO	100mg 1ml	
Glyceryl Trinitrate (Nitroglycerin)				
CAS 55-63-0 DRE-XA15586000AL	MW 227.0865 Nitroglycerin 100 µg/mL in Acetonitrile (*)	C ₃ H ₅ N ₃ O ₉	1ml	
Hexabromobenzene				
CAS 87-82-1 DRE-C14150000 DRE-L14150000CY DRE-XA14150000CY	MW 551.4882 Hexabromobenzene Hexabromobenzene 10 µg/mL in Cyclohexane Hexabromobenzene 100 µg/mL in Cyclohexane	C ₆ Br ₆	250mg 10ml 1ml	










Phenol and aromatic compounds

Product code	Description			
Hexachlorobenzene				
CAS 118-74-1	MW 284.7822	C ₆ Cl ₆		
DRE-XA14160000IO	Hexachlorobenzene 100 µg/mL in Isooctane		1ml	
DRE-XA14160000ME	Hexachlorobenzene 100 µg/mL in Methanol		1ml	
Hexogen (Hexahydro-1,3,5-trinitro-1,3,5-triazine)				
CAS 121-82-4	MW 222.1163	C ₃ H ₆ N ₆ O ₆		
DRE-LA14204000AL	Hexogen 10 µg/mL in Acetonitrile		1ml	
4-Hexylphenol				
CAS 2446-69-7	MW 178.2707	C ₁₂ H ₁₈ O		
DRE-C14209000	4-n-Hexylphenol		100mg	
Hydroquinone (Benzene-1,4-diol)				
CAS 123-31-9	MW 110.1106	C ₆ H ₆ O ₂		
DRE-C14223000	Hydrochinone(±)		250mg	
4-Hydroxybenzonitrile				
CAS 767-00-0	MW 119.1207	C ₇ H ₅ NO		
DRE-C14230000	4-Hydroxybenzonitrile		250mg	
2-Hydroxydiphenylmethane				
CAS 28994-41-4	MW 184.2338	C ₁₃ H ₁₂ O		
DRE-C14231500	2-Hydroxydiphenylmethane		500mg	
2-Hydroxy-4-nitroaniline				
CAS 121-88-0	MW 154.1234	C ₆ H ₆ N ₂ O ₃		
DRE-C14234100	2-Hydroxy-4-nitroaniline		250mg	
2-Hydroxy-5-nitroaniline				
CAS 99-57-0	MW 154.1234	C ₆ H ₆ N ₂ O ₃		
DRE-C14234200	2-Hydroxy-5-nitroaniline		250mg	
N-Isopropylaniline				
CAS 768-52-5	MW 135.2062	C ₉ H ₁₃ N		
DRE-C14463100	N-Isopropylaniline		250mg	










Phenol and aromatic compounds

Product code	Description			
4-Isopropyl-3-methylphenol				
CAS 3228-02-2 DRE-C14463650	MW 150.2176 4-Isopropyl-3-methylphenol	C ₁₀ H ₁₄ O	100mg	
2-Isopropylphenol				
CAS 88-69-7 DRE-C14463900	MW 136.191 2-Isopropylphenol	C ₉ H ₁₂ O	250mg	
3-Isopropylphenol				
CAS 618-45-1 DRE-C14464000	MW 136.191 3-Isopropylphenol	C ₉ H ₁₂ O	250mg	
4-Isopropylphenol				
CAS 99-89-8 DRE-C14464100	MW 136.191 4-Isopropylphenol	C ₉ H ₁₂ O	250mg	
Malonoben (Tyrphostin A9)				
CAS 10537-47-0 DRE-C17896100	MW 282.3801 Tyrphostin A9	C ₁₈ H ₂₂ N ₂ O	25mg	
3-Methylaniline (m-Toluidine)				
CAS 108-44-1 DRE-C17594900 DRE-XA17594900CY	MW 107.1531 m-Toluidine m-Toluidine 100 µg/mL in Cyclohexane	C ₇ H ₉ N	250mg 1ml	
4-Methylaniline (p-Toluidine)				
CAS 106-49-0 DRE-C17595000 DRE-XA17595000ME	MW 107.1531 p-Toluidine p-Toluidine 100 µg/mL in Methanol	C ₇ H ₉ N	250mg 1ml	
N-Methylaniline				
CAS 100-61-8 DRE-C15083770	MW 107.1531 N-Methylaniline	C ₇ H ₉ N	250mg	
3-Methylbenzaldehyde (m-Tolualdehyde)				
CAS 620-23-5 DRE-CA17593000	MW 120.1485 m-Tolualdehyde	C ₈ H ₈ O	250mg	










Phenol and aromatic compounds

Product code	Description			
2-Methyl-3-nitrophenol				
CAS 5460-31-1 DRE-C15108500	MW 153.1354 2-Methyl-3-nitrophenol	C ₇ H ₇ NO ₃	100mg	
3-Methyl-2-nitrophenol				
CAS 4920-77-8 DRE-C15109000	MW 153.1354 3-Methyl-2-nitrophenol	C ₇ H ₇ NO ₃	100mg	
4-Methyl-2-nitrophenol				
CAS 119-33-5 DRE-C15110800	MW 153.1354 4-Methyl-2-nitrophenol	C ₇ H ₇ NO ₃	100mg	
5-Methyl-2-nitrophenol				
CAS 700-38-9 DRE-C15112000	MW 153.1354 5-Methyl-2-nitrophenol	C ₇ H ₇ NO ₃	100mg	
2-Methylphenol (o-Cresol)				
CAS 95-48-7 DRE-C15140200 DRE-XA15140200ME	MW 108.1378 2-Methylphenol (o-Cresol)(‡) 2-Methylphenol 100 µg/mL in Methanol	C ₇ H ₈ O	500mg 1ml	
4-Methylphenol (p-Cresol)				
CAS 106-44-5 DRE-C15140400 DRE-L15140400ME DRE-XA15140400CY	MW 108.1378 4-Methylphenol (p-Cresol)(‡) 4-Methylphenol 10 µg/mL in Methanol 4-Methylphenol 100 µg/mL in Cyclohexane	C ₇ H ₈ O	500mg 10ml 1ml	
2-Nitroaniline				
CAS 88-74-4 DRE-C15554200	MW 138.124 2-Nitroaniline	C ₆ H ₆ N ₂ O ₂	250mg	
3-Nitroaniline				
CAS 99-09-2 DRE-C15554300	MW 138.124 3-Nitroaniline	C ₆ H ₆ N ₂ O ₂	250mg	
4-Nitroaniline				
CAS 100-01-6 DRE-C15554400	MW 138.124 4-Nitroaniline	C ₆ H ₆ N ₂ O ₂	250mg	









Phenol and aromatic compounds

Product code	Description			
Nitrobenzene				
CAS 98-95-3	MW 123.1094	C ₆ H ₅ NO ₂		
DRE-C15557000	Nitrobenzene		250mg	
DRE-L15557000ME	Nitrobenzene 10 µg/mL in Methanol		10ml	
DRE-XA15557000ME	Nitrobenzene 100 µg/mL in Methanol		1ml	
Nitrobenzene D5				
CAS 4165-60-0	MW 128.1402	C ₆ ² H ₅ NO ₂		
DRE-C15557100	Nitrobenzene D5		1g	
DRE-XA15557100AC	Nitrobenzene D5 100 µg/mL in Acetone		1.1ml	
DRE-YA15557100MB	Nitrobenzene D5 2000 µg/mL in Methyl-tert. butyl ether		1ml	
2-Nitrodiphenylamine				
CAS 119-75-5	MW 214.22	C ₁₂ H ₁₀ N ₂ O ₂		
DRE-C15559200	2-Nitrodiphenylamine		100mg	
4-Nitrodiphenylamine				
CAS 836-30-6	MW 214.22	C ₁₂ H ₁₀ N ₂ O ₂		
DRE-C15559400	4-Nitrodiphenylamine		100mg	
Nitroguanidine				
CAS 556-88-7	MW 104.0681	CH ₄ N ₄ O ₂		
DRE-XA15588000AL	Nitroguanidine 100 µg/mL in Acetonitrile		1ml	
4-Nitroimidazole				
CAS 3034-38-6	MW 113.0748	C ₃ H ₃ N ₃ O ₂		
DRE-C15588300	4-Nitroimidazole		100mg	
Nitromide				
CAS 121-81-3	MW 211.1317	C ₇ H ₅ N ₃ O ₅		
DRE-C15588500	Nitromid		250mg	
Nitropenta (Pentaerythritol tetranitrate)				
CAS 78-11-5	MW 316.1366	C ₅ H ₈ N ₄ O ₁₂		
DRE-LA15589000AL	Nitropenta (Pentaerythritol tetranitrate) 10 µg/mL in Acetonitrile		1ml	
2-Nitrophenol				
CAS 88-75-5	MW 139.1088	C ₆ H ₅ NO ₃		
DRE-C15590200	2-Nitrophenol		500mg	
DRE-XA15590200ME	2-Nitrophenol 100 µg/mL in Methanol		1ml	










Phenol and aromatic compounds

Product code	Description			
3-Nitrophenol				
CAS 554-84-7 DRE-C15590300	MW 139.1088 3-Nitrophenol	C ₆ H ₅ NO ₃	250mg	
4-Nitrophenol				
CAS 100-02-7 DRE-XA15590400ME	MW 139.1088 4-Nitrophenol 100 µg/mL in Methanol	C ₆ H ₅ NO ₃	1ml	
4-Nitrophenol-2,3,5,6-D4				
CAS 93951-79-2 DRE-XA15590404AC	MW 143.1334 4-Nitrophenol-2,3,5,6-D4 100 µg/mL in Acetone	C ₆ ² H ₄ HNO ₃	1.1ml	
2-Nitrotoluene				
CAS 88-72-2 DRE-C15615200 DRE-L15615200AL DRE-XA15615200AL	MW 137.136 2-Nitrotoluene 2-Nitrotoluene 10 µg/mL in Acetonitrile 2-Nitrotoluene 100 µg/mL in Acetonitrile	C ₇ H ₇ NO ₂	250mg 10ml 1ml	
2-Nitrotoluene D7				
CAS 84344-04-7 DRE-C15615205	MW 144.1791 2-Nitrotoluene D7	C ₇ ² H ₇ NO ₂	50mg	
3-Nitrotoluene				
CAS 99-08-1 DRE-C15615300	MW 137.136 3-Nitrotoluene	C ₇ H ₇ NO ₂	250mg	
4-Nitrotoluene				
CAS 99-99-0 DRE-C15615400 DRE-L15615400AL	MW 137.136 4-Nitrotoluene 4-Nitrotoluene 10 µg/mL in Acetonitrile	C ₇ H ₇ NO ₂	250mg 10ml	
4-Nitrotoluene-2-sulfonic Acid				
CAS 121-03-9 DRE-LA15615410AL	MW 217.1992 4-Nitrotoluene-2-sulfonic acid 10 µg/mL in Acetonitrile	C ₇ H ₇ NO ₅ S	1ml	
Octogen				
CAS 2691-41-0 DRE-LA15711600AL	MW 296.1551 Octogen 10 µg/mL in Acetonitrile	C ₈ H ₈ N ₈ O ₈	1ml	










Phenol and aromatic compounds

Product code	Description			
2,3,4,5,6-Pentabromoethylbenzene				
CAS 85-22-3 DRE-C15938800	MW 500.6453 2,3,4,5,6-Pentabromoethylbenzene	C ₈ H ₅ Br ₅	100mg	
Pentabromophenol				
CAS 608-71-9 DRE-C15939300	MW 488.5915 Pentabromophenol	C ₆ HBr ₅ O	100mg	
Pentachloroaniline				
CAS 527-20-8 DRE-XA15940000CY	MW 265.3518 Pentachloroaniline 100 µg/mL in Cyclohexane	C ₆ H ₂ Cl ₅ N	1ml	
Pentachloroanisole				
CAS 1825-21-4 DRE-C15950000 DRE-L15950000CY DRE-GA09010330HE	MW 280.3631 Pentachloroanisole Pentachloroanisole 10 µg/mL in Cyclohexane Pentachloroanisole 200 µg/mL in Hexane(±)	C ₇ H ₃ Cl ₅ O	100mg 10ml 1ml	
Pentachloroanisole 13C6 (ring 13C)				
CAS n/a DRE-XA15950010AC	MW 286.319 Pentachloroanisole 13C6 100 µg/mL in Acetone	¹³ C ₆ H ₃ Cl ₅ O	1.1ml	
Pentachlorothiophenol				
CAS 133-49-3 DRE-C15973180	MW 282.4021 Pentachlorothiophenol	C ₆ HCl ₅ S	100mg	
2,3,4,5,6-Pentachlorotoluene				
CAS 877-11-2 DRE-L15973200CY	MW 264.3637 2,3,4,5,6-Pentachlorotoluene 10 µg/mL in Cyclohexane	C ₇ H ₃ Cl ₅	10ml	
Pentafluorophenol				
CAS 771-61-9 DRE-C15974300	MW 184.0636 Pentafluorophenol	C ₆ HF ₅ O	500mg	
4-Pentylphenol (4-n-Amylphenol)				
CAS 14938-35-3 DRE-C10246800	MW 164.2441 4-n-Amylphenol	C ₁₁ H ₁₆ O	100mg	










Phenol and aromatic compounds

Product code	Description			
Phenol				
CAS 108-95-2	MW 94.1112	C ₆ H ₆ O		
DRE-C16025000	Phenol(‡)		1g	
DRE-L16025000ME	Phenol 10 µg/mL in Methanol		10ml	
DRE-XA16025000ME	Phenol 100 µg/mL in Methanol		1ml	
Phenol D5 (2,3,4,5,6-Pentadeuteriophenol)				
CAS 4165-62-2	MW 99.142	C ₆ ² H ₅ HO		
DRE-C16025100	Phenol D5		1g	
DRE-XA16025100AC	Phenol-D5 100 µg/mL in Acetone		1.1ml	
DRE-XA16025100ME	Phenol-D5 100 µg/mL in Methanol		1.1ml	
Phenol D6				
CAS 13127-88-3	MW 100.1482	C ₆ ² H ₆ O		
DRE-C16025200	Phenol D6		1g	
1,3-Phenylenediamine				
CAS 108-45-2	MW 108.1411	C ₆ H ₈ N ₂		
DRE-CA16057900	1,3-Phenylenediamine		100mg	
1,4-Phenylenediamine (1,4-Diaminobenzene)				
CAS 106-50-3	MW 108.1411	C ₆ H ₈ N ₂		
DRE-C16058000	1,4-Phenylenediamine (1,4-Diaminobenzene)(‡)		100mg	
DRE-YA16058000AL	1,4-Phenylenediamine (1,4-Diaminobenzene) 2000 µg/mL in Acetonitrile (*)		1ml	
3-Phenylphenol				
CAS 580-51-8	MW 170.2072	C ₁₂ H ₁₀ O		
DRE-C16070100	3-Phenylphenol		100mg	
4-Phenylphenol				
CAS 92-69-3	MW 170.2072	C ₁₂ H ₁₀ O		
DRE-C16070200	4-Phenylphenol		250mg	
4-Propylphenol				
CAS 645-56-7	MW 136.191	C ₉ H ₁₂ O		
DRE-C16530240	4-n-Propylphenol		100mg	
1,2,4,5-Tetrabromobenzene				
CAS 636-28-2	MW 393.6961	C ₆ H ₂ Br ₄		
DRE-C17324500	1,2,4,5-Tetrabromobenzene		100mg	










Phenol and aromatic compounds

Product code	Description			
2,3,4,5-Tetrachloroaniline				
CAS 634-83-3	MW 230.9067	C ₆ H ₃ Cl ₄ N		
DRE-C17330400	2,3,4,5-Tetrachloroaniline		10mg	
DRE-L17330400CY	2,3,4,5-Tetrachloroaniline 10 µg/mL in Cyclohexane		10ml	
2,3,4,5-Tetrachloroanisole				
CAS 938-86-3	MW 245.9181	C ₇ H ₄ Cl ₄ O		
DRE-L17333100IO	2,3,4,5-Tetrachloroanisole 10000 µg/mL in iso-Octane		10ml	
2,3,5,6-Tetrachloroanisole				
CAS 6936-40-9	MW 245.9181	C ₇ H ₄ Cl ₄ O		
DRE-L17333300IO	2,3,5,6-Tetrachloroanisole 10 µg/mL in Isooctane		10ml	
1,2,3,5-Tetrachlorobenzene				
CAS 634-90-2	MW 215.8921	C ₆ H ₂ Cl ₄		
DRE-C17353500	1,2,3,5-Tetrachlorobenzene		100mg	
DRE-L17353500CY	1,2,3,5-Tetrachlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17353500ME	1,2,3,5-Tetrachlorobenzene 100 µg/mL in Methanol		1ml	
1,2,4,5-Tetrachlorobenzene				
CAS 95-94-3	MW 215.8921	C ₆ H ₂ Cl ₄		
DRE-C17354500	1,2,4,5-Tetrachlorobenzene		100mg	
DRE-L17354500CY	1,2,4,5-Tetrachlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17354500ME	1,2,4,5-Tetrachlorobenzene 100 µg/mL in Methanol		1ml	
1,2,4,5-Tetrachlorobenzene-13C6				
CAS 85380-73-0	MW 221.848	¹³ C ₆ H ₂ Cl ₄		
DRE-XA17354501AL	1,2,4,5-Tetrachlorobenzene 13C6 100 µg/mL in Acetonitrile		1ml	
2,3,4,5-Tetrachloronitrobenzene				
CAS 879-39-0	MW 260.8896	C ₆ HCl ₄ NO ₂		
DRE-C17364500	2,3,4,5-Tetrachloronitrobenzene		100mg	
2,3,4,6-Tetrachlorophenol Acetate				
CAS 5435-60-9	MW 273.9282	C ₈ H ₄ Cl ₄ O ₂		
DRE-L17376100IO	2,3,4,6-Tetrachlorophenol acetate 10 µg/mL in Isooctane		10ml	
2,3,5,6-Tetrachlorophenol Acetate				
CAS 61925-90-4	MW 273.9282	C ₈ H ₄ Cl ₄ O ₂		
DRE-L17376200IO	2,3,5,6-Tetrachlorophenol acetate 10 µg/mL in Isooctane (*)		10ml	










Phenol and aromatic compounds

Product code	Description			
2,3,4,5-Tetrachlorophenol				
CAS 4901-51-3	MW 231.8915	C ₆ H ₂ Cl ₄ O		
DRE-C17374500	2,3,4,5-Tetrachlorophenol		10mg	
DRE-L17374500CY	2,3,4,5-Tetrachlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA17374500ME	2,3,4,5-Tetrachlorophenol 100 µg/mL in Methanol		1ml	
2,3,4,6-Tetrachlorophenol				
CAS 58-90-2	MW 231.8915	C ₆ H ₂ Cl ₄ O		
DRE-XA17374600ME	2,3,4,6-Tetrachlorophenol 100 µg/mL in Methanol		1ml	
2,3,5,6-Tetrachloropyridine				
CAS 2402-79-1	MW 216.8801	C ₅ HCl ₄ N		
DRE-C17376300	2,3,5,6-Tetrachloropyridine		100mg	
α,α,2,6-Tetrachlorotoluene				
CAS 81-19-6	MW 229.9187	C ₈ H ₄ Cl ₄		
DRE-C17381000	α,α,2,6-Tetrachlorotoluene		1g	
DRE-L17381000CY	α,α,2,6-Tetrachlorotoluene 10 µg/mL in Cyclohexane		10ml	
α,α,α,4-Tetrachlorotoluene				
CAS 5216-25-1	MW 229.9187	C ₈ H ₄ Cl ₄		
DRE-C17380000	α,α,α,4-Tetrachlorotoluene		1g	
2,4,5,6-Tetrachloro-m-xylene				
CAS 877-09-8	MW 243.9452	C ₈ H ₆ Cl ₄		
DRE-C17382500	2,4,5,6-Tetrachloro-m-xylene		100mg	
DRE-L17382500CY	2,4,5,6-Tetrachloro-m-xylene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17382500ME	2,4,5,6-Tetrachloro-m-xylene 100 µg/mL in Methanol		1ml	
Thiophenol				
CAS 108-98-5	MW 110.1768	C ₆ H ₆ S		
DRE-CA17548000	Thiophenol		1ml	
Thymol				
CAS 89-83-8	MW 150.2176	C ₁₀ H ₁₄ O		
DRE-C17575200	Thymol		250mg	
m-Tolualdehyd-2,4-dinitrophenylhydrazone				
CAS 2880-05-9	MW 300.2695	C ₁₄ H ₁₂ N ₄ O ₄		
DRE-C17593005	m-Tolualdehyd-2,4-dinitrophenylhydrazone		100mg	










Phenol and aromatic compounds

Product code	Description			
o-Tolualdehyd-2,4-dinitrophenylhydrazone				
CAS 1773-44-0 DRE-C17593015	MW 300.2695 o-Tolualdehyd-2,4-dinitrophenylhydrazone	C ₁₄ H ₁₂ N ₄ O ₄	100mg	
p-Tolualdehyd-2,4-dinitrophenylhydrazone				
CAS 2571-00-8 DRE-C17593025	MW 300.2695 p-Tolualdehyd-2,4-dinitrophenylhydrazone	C ₁₄ H ₁₂ N ₄ O ₄	100mg	
o-Toluidine				
CAS 95-53-4 DRE-C17594800 DRE-L17594800CY DRE-XA17594800CY DRE-GA09010374ME	MW 107.1531 o-Toluidine o-Toluidine 10 µg/mL in Cyclohexane o-Toluidine 100 µg/mL in Cyclohexane o-Toluidine (SS) 500 µg/mL in Methanol(‡)	C ₇ H ₉ N	1ml 10ml 1ml 1ml	
2,4,6-Tribromoanisole				
CAS 607-99-8 DRE-C17664000 DRE-L17664000IO	MW 344.826 2,4,6-Tribromoanisole 2,4,6-Tribromoanisole 10 µg/mL in Isooctane	C ₇ H ₅ Br ₃ O	100mg 10ml	
1,2,4-Tribromobenzene				
CAS 615-54-3 DRE-C17664900	MW 314.8 1,2,4-Tribromobenzene	C ₆ H ₃ Br ₃	100mg	
1,3,5-Tribromobenzene				
CAS 626-39-1 DRE-C17665000	MW 314.8 1,3,5-Tribromobenzene	C ₆ H ₃ Br ₃	250mg	
2,4,6-Tribromophenol				
CAS 118-79-6 DRE-C17666000 DRE-L17666000CY DRE-L17666000ME DRE-XA17666000ME	MW 330.7994 2,4,6-Tribromophenol(‡) 2,4,6-Tribromophenol 10 µg/mL in Cyclohexane 2,4,6-Tribromophenol 10 µg/mL in Methanol 2,4,6-Tribromophenol 100 µg/mL in Methanol	C ₆ H ₃ Br ₃ O	250mg 10ml 10ml 1ml	
2,4,6-Tri-tert-butylphenol				
CAS 732-26-3 DRE-C17667700	MW 262.4302 2,4,6-Tri-tert-butylphenol	C ₁₈ H ₃₀ O	250mg	
2,3,4-Trichloroaniline				
CAS 634-67-3 DRE-C17700000	MW 196.4617 2,3,4-Trichloroaniline	C ₆ H ₄ Cl ₃ N	100mg	









Phenol and aromatic compounds

Product code	Description			
2,4,5-Trichloroaniline				
CAS 636-30-6	MW 196.4617	C ₆ H ₃ Cl ₃ N		
DRE-C17700200	2,4,5-Trichloroaniline		250mg	
DRE-L17700200CY	2,4,5-Trichloroaniline 10 µg/mL in Cyclohexane		10ml	
2,4,6-Trichloroaniline				
CAS 634-93-5	MW 196.4617	C ₆ H ₃ Cl ₃ N		
DRE-C17700600	2,4,6-Trichloroaniline		250mg	
3,4,5-Trichloroaniline				
CAS 634-91-3	MW 196.4617	C ₆ H ₃ Cl ₃ N		
DRE-C17700300	3,4,5-Trichloroaniline		100mg	
2,4,6-Trichloroanisole D5				
CAS 352439-08-8	MW 216.5038	C ₇ H ₅ Cl ₃ O		
DRE-C17714601	2,4,6-Trichloroanisole D5		25mg	
DRE-XA17714601AL	2,4,6-Trichloroanisole D5 100 µg/mL in Acetonitrile		1ml	
2,3,4-Trichloroanisole				
CAS 54135-80-7	MW 211.473	C ₇ H ₅ Cl ₃ O		
DRE-C17713400	2,3,4-Trichloroanisole		100mg	
DRE-L17713400IO	2,3,4-Trichloroanisole 10 µg/mL in Isooctane		10ml	
2,3,6-Trichloroanisole				
CAS 50375-10-5	MW 211.473	C ₇ H ₅ Cl ₃ O		
DRE-C17713600	2,3,6-Trichloroanisole		100mg	
DRE-L17713600CY	2,3,6-Trichloroanisole 10 µg/mL in Cyclohexane		10ml	
2,4,5-Trichloroanisole				
CAS 6130-75-2	MW 211.473	C ₇ H ₅ Cl ₃ O		
DRE-L17714500IO	2,4,5-Trichloroanisole 10 µg/mL in Isooctane		10ml	
2,4,6-Trichloroanisole				
CAS 87-40-1	MW 211.473	C ₇ H ₅ Cl ₃ O		
DRE-C17714600	2,4,6-Trichloroanisole		100mg	
DRE-L17714600IO	2,4,6-Trichloroanisole 10 µg/mL in Isooctane		10ml	
DRE-XA17714600ME	2,4,6-Trichloroanisole 100 µg/mL in Methanol		1ml	
1,2,3-Trichlorobenzene				
CAS 87-61-6	MW 181.447	C ₆ H ₃ Cl ₃		
DRE-C17722300	1,2,3-Trichlorobenzene		1g	
DRE-GA17722300ME	1,2,3-Trichlorobenzene 1000 µg/mL in Methanol		1ml	
DRE-GH17722300ME	1,2,3-Trichlorobenzene 1000 µg/mL in Methanol		5x1ml	
DRE-L17722300AL	1,2,3-Trichlorobenzene 10 µg/mL in Acetonitrile		10ml	
DRE-L17722300CY	1,2,3-Trichlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17722300IO	1,2,3-Trichlorobenzene 100 µg/mL in Isooctane		1ml	










Phenol and aromatic compounds

Product code	Description			
1,2,4-Trichlorobenzene				
CAS 120-82-1	MW 181.447	C ₆ H ₃ Cl ₃		
DRE-C17722400	1,2,4-Trichlorobenzene		1g	
DRE-L17722400CY	1,2,4-Trichlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17722400HE	1,2,4-Trichlorobenzene 100 µg/mL in n-Hexane		1ml	
1,3,5-Trichlorobenzene				
CAS 108-70-3	MW 181.447	C ₆ H ₃ Cl ₃		
DRE-C17723500	1,3,5-Trichlorobenzene		250mg	
DRE-L17723500AL	1,3,5-Trichlorobenzene 10 µg/mL in Acetonitrile		10ml	
DRE-L17723500CY	1,3,5-Trichlorobenzene 10 µg/mL in Cyclohexane		10ml	
DRE-XA17723500ME	1,3,5-Trichlorobenzene 100 µg/mL in Methanol		1ml	
DRE-YS09010027ME	1,3,5-Trichlorobenzene 1000 µg/mL in p/t Methanol(‡)		5x1ml	
1,3,5-Trichlorobenzene D3				
CAS 1198-60-3	MW 184.4655	C ₆ ² H ₃ Cl ₃		
DRE-XA17723600AC	1,3,5-Trichlorobenzene D3 100 µg/mL in Acetone		1.1ml	
1,3,5-Trichloro-2-nitrobenzene				
CAS 18708-70-8	MW 226.4446	C ₆ H ₂ Cl ₃ NO ₂		
DRE-C17761300	1,3,5-Trichloro-2-nitrobenzene		250mg	
2,3,4-Trichloronitrobenzene				
CAS 17700-09-3	MW 226.4446	C ₆ H ₂ Cl ₃ NO ₂		
DRE-C17763400	2,3,4-Trichloronitrobenzene		250mg	
2,4,5-Trichloronitrobenzene				
CAS 89-69-0	MW 226.4446	C ₆ H ₂ Cl ₃ NO ₂		
DRE-C17762400	2,4,5-Trichloronitrobenzene		250mg	
2,3,4-Trichlorophenol				
CAS 15950-66-0	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17773400	2,3,4-Trichlorophenol(‡)		25mg	
DRE-L17773400CY	2,3,4-Trichlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA17773400ME	2,3,4-Trichlorophenol 100 µg/mL in Methanol		1ml	
2,3,5-Trichlorophenol				
CAS 933-78-8	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17773500	2,3,5-Trichlorophenol(‡)		100mg	
DRE-L17773500IO	2,3,5-Trichlorophenol 10 µg/mL in Isooctane		10ml	
DRE-GA09010070ME	2,3,5-Trichlorophenol 100 µg/mL in Methanol(‡)		1ml	
2,3,6-Trichlorophenol				
CAS 933-75-5	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17773600	2,3,6-Trichlorophenol(‡)		100mg	
DRE-L17773600ME	2,3,6-Trichlorophenol 10 µg/mL in Methanol		10ml	




Phenol and aromatic compounds

Product code	Description			
2,4,6-Trichlorophenol				
CAS 88-06-2	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17774600	2,4,6-Trichlorophenol(±)		100mg	
DRE-L17774600AL	2,4,6-Trichlorophenol 10 µg/mL in Acetonitrile		10ml	
DRE-L17774600CY	2,4,6-Trichlorophenol 10 µg/mL in Cyclohexane		10ml	
DRE-L17774600ME	2,4,6-Trichlorophenol 10 µg/mL in Methanol		10ml	
DRE-XA17774600ME	2,4,6-Trichlorophenol 100 µg/mL in Methanol		1ml	
DRE-YA17774600ME	2,4,6-Trichlorophenol 5000 µg/mL in Methanol		1ml	
3,4,5-Trichlorophenol				
CAS 609-19-8	MW 197.4464	C ₆ H ₃ Cl ₃ O		
DRE-C17774800	3,4,5-Trichlorophenol(±)		25mg	
DRE-L17774800AL	3,4,5-Trichlorophenol 10 µg/mL in Acetonitrile		10ml	
DRE-L17774800IO	3,4,5-Trichlorophenol 10 µg/mL in Isooctane		10ml	
DRE-XA17774800AL	3,4,5-Trichlorophenol 100 µg/mL in Acetonitrile		1ml	
2,3,6-Trichlorotoluene				
CAS 2077-46-5	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-LA17787900ME	2,3,6-Trichlorotoluene 10 µg/mL in Methanol		1ml	
2,4,5-Trichlorotoluene				
CAS 6639-30-1	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-LA17788000HE	2,4,5-Trichlorotoluene 10 µg/mL in n-Hexane		1ml	
α,2,4-Trichlorotoluene				
CAS 94-99-5	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-C17787400	α,2,4-Trichlorotoluene		250mg	
α,2,6-Trichlorotoluene				
CAS 2014-83-7	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-C17787600	α,2,6-Trichlorotoluene		250mg	
α,3,4-Trichlorotoluene				
CAS 102-47-6	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-C17787800	α,3,4-Trichlorotoluene		1ml	
α,α,α-Trichlorotoluene				
CAS 98-07-7	MW 195.4736	C ₇ H ₅ Cl ₃		
DRE-C17787000	α,α,α-Trichlorotoluene (Benzotrichloride)		250mg	

Phenol and aromatic compounds

Product code	Description			
2-Trifluoromethylaniline				
CAS 88-17-5 DRE-C17844900	MW 161.1244 2-Trifluoromethylaniline	C ₇ H ₆ F ₃ N	250mg	
α,α,α-Trifluorotoluene (Trifluoromethylbenzene)				
CAS 98-08-8 DRE-C17846000 DRE-YA17846000ME	α,α,α-Trifluorotoluene α,α,α-Trifluorotoluene 2000 µg/mL in Methanol	C ₇ H ₅ F ₃	250mg 1ml	
1,2,4-Trimethoxybenzene				
CAS 135-77-3 DRE-C17876420	1,2,4-Trimethoxybenzene	C ₉ H ₁₂ O ₃	100mg	
1,3,5-Trimethoxybenzene				
CAS 621-23-8 DRE-C17876430	1,3,5-Trimethoxybenzene	C ₉ H ₁₂ O ₃	100mg	
2,4,5-Trimethylaniline				
CAS 137-17-7 DRE-C17878000 DRE-LA17878000AL	2,4,5-Trimethylaniline(±) 2,4,5-Trimethylaniline 10 µg/mL in Acetonitrile	C ₉ H ₁₁ N	10mg 1ml	
2,3,6-Trimethylphenol				
CAS 2416-94-6 DRE-C17883600	2,3,6-Trimethylphenol	C ₉ H ₁₂ O	250mg	
2,4,6-Trimethylphenol				
CAS 527-60-6 DRE-C17883700	2,4,6-Trimethylphenol	C ₉ H ₁₂ O	250mg	
1,3,5-Trinitrobenzene				
CAS 99-35-4 DRE-L17889500ME	1,3,5-Trinitrobenzene 10 µg/mL in Methanol	C ₆ H ₃ N ₃ O ₆	10ml	
2,4,6-Trinitrophenol (Picric Acid)				
CAS 88-89-1 DRE-L17890500AL DRE-XA17890500AL	2,4,6-Trinitrophenol (Picric acid) 10 µg/mL in Acetonitrile 2,4,6-Trinitrophenol (Picric acid) 100 µg/mL in Acetonitrile	C ₆ H ₃ N ₃ O ₇	10ml 1ml	

Phenol and aromatic compounds

Product code	Description			
2,4,5-Trinitrotoluene				
CAS 610-25-3	MW 227.1311	C ₇ H ₅ N ₃ O ₆		
DRE-L17891100CY	2,4,5-Trinitrotoluene 10 µg/mL in Cyclohexane		10ml	
2,4,6-Trinitrotoluene (TNT)				
CAS 118-96-7	MW 227.1311	C ₇ H ₅ N ₃ O ₆		
DRE-LA17891200CY	2,4,6-Trinitrotoluene (TNT) 10 µg/mL in Cyclohexane		1ml	
DRE-XA17891200CY	2,4,6-Trinitrotoluene (TNT) 100 µg/mL in Cyclohexane		1ml	
Tyramine Hydrochloride				
CAS 60-19-5	MW 173.64	C ₈ H ₁₁ NO·ClH		
DRE-C17895800	Tyramine hydrochloride		100mg	
Acid Surrogate Solution				
DRE-YS09000039ME	Acid Surrogate Mix 10000 µg/mL in Methanol(‡)			5x5ml
	phenol-d5 2,4,6-tribromophenol 2-fluorophenol			
Aromatic Comp.-Mix 16				
DRE-LA19070016ME	Aromatic Comp.-Mix 16 in Methanol (*)			1ml
	2,3,4-Trichloroaniline [10 µg/mL] 2,6-Dimethylaniline [10 µg/mL] 3,5-Dichloroaniline [10 µg/mL] 4-Bromoaniline [10 µg/mL] o-Toluidine [10 µg/mL]	2,4,5-Trichloroaniline [10 µg/mL] 2-Chloroaniline [10 µg/mL] 3-Chloro-4-methylaniline [10 µg/mL] 4-Chloro-2-methylaniline [10 µg/mL]	2,4-Dichloroaniline [10 µg/mL] 2-Trifluoromethylaniline [20 µg/mL] 3-Chloroaniline [10 µg/mL] 4-Chloroaniline [10 µg/mL]	2,6-Dichloroaniline [10 µg/mL] 3,4-Dichloroaniline [10 µg/mL] 3-Trifluoromethylaniline [10 µg/mL] N,N-Dimethylaniline [10 µg/mL]
Mix of Aromatic Hydrocarbons 2				
DRE-L19070002ME	Mix of Aromatic Hydrocarbons 2 10 µg/mL in Methanol			10ml
DRE-X19070002ME	Mix of Aromatic Hydrocarbons 2 100 µg/mL in Methanol			10ml
	Benzene Ethylbenzene o-Xylene p-Xylene Toluene			
Aromatic Hydrocarbons Mix 11				
DRE-YA04000100ME	Aromatic Hydrocarbons Mix 11 2000 µg/mL in Methanol(‡)			1ml
	Benzene Ethylbenzene m-Xylene o-Xylene p-Xylene Toluene			
Mix of Aromatic Hydrocarbons 13 (DIN 38407-F19)				
DRE-YA19070013HE	Mix of Aromatic Hydrocarbons 13 (DIN 38407-F19) 10000 µg/mL			1ml
	1,2,3-Trimethylbenzene 1,3,5-Trimethylbenzene Ethylbenzene o-Xylene Toluene	1,2,4-Trimethylbenzene Benzene m-Xylene p-Xylene		

Phenol and aromatic compounds

Product code	Description		
Chlorinated Aromatics-Mix 1			
DRE-LA19060001HE	Chlorinated Aromatics-Mix 1 10 µg/mL in n-Hexane	1ml	
	1,2,3,4-Tetrachlorobenzene	1,2,3,5-Tetrachlorobenzene	
	1,2,4,5-Tetrachlorobenzene	1,2,4-Trichlorobenzene	
	1,2-Dichlorobenzene	1,3,5-Trichlorobenzene	
	1,3-Dichlorobenzene	Hexachlorobenzene	
	Pentachlorobenzene		
Chlorinated Aromatics-Mix 2			
DRE-LA19060002CY	Chlorinated Aromatics-Mix 2 10 µg/mL in Cyclohexane	1ml	
	1,2,3,4-Tetrachlorobenzene		
	1,2,3,5-Tetrachlorobenzene		
	1,2,3-Trichlorobenzene		
	1,2,4,5-Tetrachlorobenzene		
	1,2,4-Trichlorobenzene		
	1,3,5-Trichlorobenzene		
Nitroaromate-Explosive Mix 3			
DRE-LA08330300AL	Nitroaromate-Explosive Mix 3 10 µg/mL in Acetonitrile	1ml	
1,2-Dinitrobenzene	1,3,5-Trinitrobenzene	1,3-Dinitrobenzene	2,4,6-Trinitrophenol (Picric acid)
2,4,6-Trinitrotoluene (TNT)	2,4-Dinitrotoluene	2,6-Dinitrotoluene	2-Amino-4,6-dinitrotoluene
2-Amino-4-nitrotoluene	2-Amino-6-nitrotoluene	2-Nitrotoluene	3-Nitrotoluene
4-Amino-2,6-dinitrotoluene	4-Nitrotoluene	Hexogen	Nitrobenzene
Nitroguanidine	Nitropenta (Pentaerythritol tetranitrate)	Octogen	Tetryl
Nitroaromate-Nitroamine-Mix 2			
DRE-XA08330200AL	Nitroaromate-Nitroamine-Mix 2 100 µg/mL in Acetonitrile	1ml	
	2,6-Dinitrotoluene		
	2-Amino-4,6-dinitrotoluene		
	2-Nitrotoluene		
	3-Nitrotoluene		
	4-Amino-2,6-dinitrotoluene		
	4-Nitrotoluene		
	Tetryl		
Nitroaromate-Nitroamine-Mix 4 (DIN38407-F17)			
DRE-LA08330400AL	Nitroaromate-Nitroamine-Mix 4 (DIN38407-F17) 10 µg/mL in Acetonitrile	1ml	
1,2-Dinitrobenzene	1,3-Dinitrobenzene	1,4-Dinitrobenzene	2,3-Dinitrotoluene
2,4,6-Trinitrotoluene (TNT)	2,4-Dinitrotoluene	2,6-Dinitrotoluene	2-Amino-4,6-dinitrotoluene
2-Amino-4-nitrotoluene	2-Amino-6-nitrotoluene	2-Nitrotoluene	3,4-Dinitrotoluene
3-Nitrotoluene	4-Amino-2,6-dinitrotoluene	4-Amino-2-nitrotoluene	4-Nitrotoluene
Nitrobenzene			
Phenol-Mix 1			
DRE-L19000001ME	Phenol-Mix 1 50 µg/mL in Methanol(‡)	10ml	
	2,3,5-Trimethylphenol	2,3,6-Trimethylphenol	
	2,3-Dimethylphenol	2,4,6-Trimethylphenol	
	2,4-Dimethylphenol	2,5-Dimethylphenol	
	2,6-Dimethylphenol	2-Methylphenol	
	3,4,5-Trimethylphenol	3,4-Dimethylphenol	
	3,5-Dimethylphenol	3-Methylphenol	
	4-Methylphenol	Phenol	
Phenol-Mix 2			
DRE-L19000002AL	Phenol-Mix 2 10 µg/mL in Acetonitrile	10ml	
	2,3,4,6-Tetrachlorophenol		
	2,4,6-Trichlorophenol		
	2,4-Dichlorophenol		
	2-Chlorophenol		
	4-Chloro-2-methylphenol		
	Pentachlorophenol		

Phenol and aromatic compounds

Product code	Description	
Phenol-Mix 3		
DRE-X19000003AL	Phenol-Mix 3 100 µg/mL in Acetonitrile	10ml
	2,4,6-Trichlorophenol 2,4-Dimethylphenol 2-Chlorophenol 4-Chloro-3-methylphenol DNOC Phenol	2,4-Dichlorophenol 2,4-Dinitrophenol 2-Nitrophenol 4-Nitrophenol Pentachlorophenol
Phenol-Mix 5		
DRE-YA19000005AL	Phenol-Mix 5 2500 µg/mL in Acetonitrile	1ml
	2,3,4,6-Tetrachlorophenol 2,4-Dichlorophenol 2,4-Dinitrophenol 2-Chlorophenol 2-Nitrophenol 4-Chloro-3-methylphenol DNOC Phenol	2,4,6-Trichlorophenol 2,4-Dimethylphenol 2-Bromophenol 2-Methylphenol 3-Methylphenol 4-Nitrophenol Pentachlorophenol
Phenol-Mix 10		
DRE-L19000010AL	Phenol-Mix 10 50 µg/mL in Acetonitrile	10ml
DRE-LA19000010AL	Phenol-Mix 10 50 µg/mL in Acetonitrile	1ml
	2,3,4,5-Tetrachlorophenol 2,3,5-Trichlorophenol 2,4,6-Trichlorophenol 2-Chlorophenol 3-Chlorophenol	2,3,4,6-Tetrachlorophenol 2,3,6-Trichlorophenol 2,4-Dichlorophenol 3,4,5-Trichlorophenol 4-Chlorophenol
	2,3,4-Trichlorophenol 2,3-Dichlorophenol 2,5-Dichlorophenol 3,4-Dichlorophenol Pentachlorophenol	2,3,5,6-Tetrachlorophenol 2,4,5-Trichlorophenol 2,6-Dichlorophenol 3,5-Dichlorophenol
Phenol-Mix 15		
DRE-YA04001500ME	Phenols Mix 15 2000 µg/mL in Methanol	1ml
	2,3,4,6-Tetrachlorophenol 2,4-Dichlorophenol 2,4-Dinitrophenol 2-Chlorophenol 4-Chloro-3-methylphenol DNOC Phenol	2,4,6-Trichlorophenol 2,4-Dimethylphenol 2,6-Dichlorophenol 2-Nitrophenol 4-Nitrophenol Pentachlorophenol
Phenol-Mix 16		
DRE-YA19000016AC	Phenol-Mix 16 2000 µg/mL in Acetone	1ml
	2,4,5-Trichlorophenol 2,4-Dichlorophenol 2,5-Dimethylphenol 2-Methylphenol 3-Methylphenol 4-Chlorophenol Pentachlorophenol	2,4,6-Trichlorophenol 2,4-Dimethylphenol 2-Chlorophenol 3,5-Dimethylphenol 4-Chloro-3-methylphenol 4-Methylphenol Phenol
Phenol-Mix 17		
DRE-YA19000017ME	Phenol-Mix 17 1000 µg/mL in Methanol	1ml
	2,4,6-Trichlorophenol 2-Methylphenol 3-Methylphenol Pentachlorophenol	2-Chlorophenol 3-Chlorophenol 4-Methylphenol Phenol
Phenol-Mix 18		
DRE-LA19000018AL	Phenol-Mix 18 10 µg/mL in Acetonitrile	1ml
	2,3,4-Trichlorophenol 2,4,5-Trichlorophenol 2,6-Dichlorophenol 3,5-Dichlorophenol	2,3,5-Trichlorophenol 2,4,6-Trichlorophenol 2-Chlorophenol 3-Chlorophenol
	2,3,6-Trichlorophenol 2,4-Dichlorophenol 3,4,5-Trichlorophenol 4-Chlorophenol	2,3-Dichlorophenol 2,5-Dichlorophenol 3,4-Dichlorophenol Pentachlorophenol

Phenol and aromatic compounds

Product code	Description	
Phenol-Mix 19		
DRE-YA08271900IP	Phenol-Mix 19 2000 µg/mL in Isopropanol	1ml
2,3,4,6-Tetrachlorophenol	2,4,5-Trichlorophenol	2,4,6-Trichlorophenol
2,4-Dimethylphenol	2,4-Dinitrophenol	2,6-Dichlorophenol
2-Methylphenol	2-Nitrophenol	3-Methylphenol
4-Methylphenol	4-Nitrophenol	Dinoseb
Pentachlorophenol	Phenol	2,4-Dichlorophenol
		2-Chlorophenol
		4-Chloro-3-methylphenol
		DNOC
Phenol-Mix 27		
DRE-L19000027ME	Phenol-Mix 27 50 µg/mL in Methanol	10ml
	2,3-Dimethylphenol	2,5-Dimethylphenol
	2,6-Dimethylphenol	2-Methylphenol
	3,4-Dimethylphenol	3,5-Dimethylphenol
	3-Methylphenol	4-Methylphenol
	Phenol	
Surrogate Standard Acid Mix 18		
DRE-YA08271800ME	Surrogate Standard Acid Mix 18 2000 µg/mL in Methanol	1ml
	2,4,6-Tribromophenol	
	2-Fluorophenol	
	Phenol-2,3,4,5,6 D5	









FOOD RELATED COMPOUNDS












Food related compounds

Product code	Description			
Acesulfame Potassium				
CAS 55589-62-3 DRE-C10010800	MW 201.2422 Acesulfame potassium(‡)	C ₄ H ₄ NO ₄ S·K	250mg	
Acetic Acid 2-Ethoxyethyl Ester				
CAS 111-15-9 DRE-C10016000	MW 132.1577 Acetic acid-2-ethoxyethyl ester	C ₆ H ₁₂ O ₃	1ml	
Acetic Acid sec-Butyl Ester				
CAS 105-46-4 DRE-C10015710	MW 116.1583 Acetic acid-sec.-butyl ester	C ₆ H ₁₂ O ₂	250mg	
Acetic Acid tert-Butyl Ester				
CAS 540-88-5 DRE-CA10015720	MW 116.1583 Acetic acid-tert.-butyl ester	C ₆ H ₁₂ O ₂	1ml	
Acetic Acid Isopropyl Ester				
CAS 108-21-4 DRE-C10016200	MW 102.1317 Acetic acid-isopropyl ester	C ₅ H ₁₀ O ₂	1ml	
Acetic Acid Linalyl Ester				
CAS 115-95-7 DRE-C10016230	MW 196.286 Acetic acid-linalyl ester	C ₁₂ H ₂₀ O ₂	250mg	
Acetic Acid 3-Methylbutyl Ester				
CAS 123-92-2 DRE-C10016270	MW 130.1849 Acetic acid-3-methyl butyl ester	C ₇ H ₁₄ O ₂	250mg	
Acetic Acid 2-Pentyl Ester				
CAS 626-38-0 DRE-C10016400	MW 130.1849 Acetic acid-2-pentyl ester	C ₇ H ₁₄ O ₂	250mg	
Acetic Acid Phenyl Ester				
CAS 122-79-2 DRE-C10016450	MW 136.1479 Acetic acid-phenyl ester	C ₈ H ₈ O ₂	250mg	

Food related compounds

Product code	Description			
Acetic acid				
CAS 64-19-7 DRE-C10015500	MW 60.052 Acetic acid	C ₂ H ₄ O ₂	5ml	
L-Arabinitol (L-(-)-Arabit)				
CAS 7643-75-6 DRE-C10298000	MW 152.1458 L(-)-Arabit	C ₅ H ₁₂ O ₅	100mg	
L(+)-Arabinose				
CAS 5328-37-0 DRE-C10297500	MW 150.1299 L(+)-Arabinose	C ₅ H ₁₀ O ₅	250mg	
D(+)-Arabitol				
CAS 488-82-4 DRE-C10297900	MW 152.1458 D(+)-Arabit(±)	C ₅ H ₁₂ O ₅	100mg	
L-Ascorbic acid 2-monophosphate tri-cyclohexylammonium salt				
CAS 65718-26-5 DRE-C10303700	MW 519.5653 L-Ascorbic acid 2-monophosphate tri-cyclohexylammonium salt(±)	C ₁₈ H ₃₃ N ₃ C ₆ H ₉ O ₉ P	100mg	
L-(+)-Ascorbic Acid				
CAS 50-81-7 DRE-C10303000	MW 176.1241 L-(+)-Ascorbic acid(±)	C ₆ H ₈ O ₆	250mg	
Ascorbic acid calcium salt				
CAS 5743-27-1 DRE-C10303100	MW 390.3104 Ascorbic acid calcium salt	2C ₆ H ₇ O ₆ ·Ca	250mg	
Ascorbyl Palmitate				
CAS 137-66-6 DRE-C10303930	MW 414.5329 Ascorbyl palmitate	C ₂₂ H ₃₈ O ₇	250mg	
Aspartame				
CAS 22839-47-0 DRE-C10304940	MW 294.3031 Aspartame(±)	C ₁₄ H ₁₈ N ₂ O ₅	250mg	










Food related compounds

Product code	Description			
Astaxanthin				
CAS 472-61-7 DRE-CA10307000	MW 596.8385 Astaxanthin (*)	C ₄₀ H ₅₂ O ₄	100mg	
Azelaic Acid (1,7-Heptanedicarboxylic acid)				
CAS 123-99-9 DRE-C14126200	MW 188.2209 1,7-Heptanedicarboxylic acid	C ₉ H ₁₆ O ₄	250mg	
Benzoic Acid				
CAS 65-85-0 DRE-C10537500 DRE-XA10537500AC DRE-YS09010010DI	MW 122.1213 Benzoic acid(‡) Benzoic acid 100 µg/mL in Acetone Benzoic Acid 2000 µg/mL in Dichloromethane(‡)	C ₇ H ₆ O ₂	250mg 1ml 5x1ml	
Benzyl Benzoate				
CAS 120-51-4 DRE-C10537700 DRE-YA10537700HE	MW 212.2439 Benzoic acid-benzyl ester (Benzyl benzoate)(‡) Benzoic acid-benzyl ester 5000 µg/mL in n-Hexane	C ₁₄ H ₁₂ O ₂	500mg 1ml	
Betacarotene				
CAS 7235-40-7 DRE-CA11045800	MW 536.8726 β-Carotene	C ₄₀ H ₅₆	250mg	
D-(+)-Biotin				
CAS 58-85-5 DRE-C10625000	MW 244.3106 D(+)-Biotin	C ₁₀ H ₁₆ N ₂ O ₃ S	250mg	
Butyl Parahydroxybenzoate (4-Hydroxybenzoic acid n-butyl ester)				
CAS 94-26-8 DRE-C14228780	MW 194.2271 4-Hydroxybenzoic acid-n-butyl ester(‡)	C ₁₁ H ₁₄ O ₃	100mg	
tert-Butylhydroquinone				
CAS 1948-33-0 DRE-C10931200	MW 166.217 tert-Butylhydroquinone	C ₁₀ H ₁₄ O ₂	.25g	
Butylhydroxyanisole				
CAS 25013-16-5 DRE-C10931220 DRE-GA09010379ME	MW 360.4871 tert-Butyl-4-hydroxyanisole (mixture of 2- and 3-isomer)(‡) Butylated Hydroxyanisole (SS) 500 µg/mL in Methanol(‡)	2C ₁₁ H ₁₆ O ₂	100mg 1ml	










Food related compounds

Product code	Description			
Butylhydroxytoluene				
CAS 128-37-0 DRE-C12253500	MW 220.3505 2,6-Di-tert.-butyl-4-methylphenol(‡)	C ₁₅ H ₂₄ O	250mg	
Butylhydroxytoluene-d21 (BHT-d21; 2,6-Bis[1,1-di(methyl-d3)ethyl-2,2,2-d3]-4-methylphen-3,5-d2-ol-d)				
CAS 64502-99-4 DRE-GS09010395ME	d21-2,6-di(tert-butyl)-4-methylphenol (BHT-d21), 1000 µg/mL in Methanol(‡)	C ₁₅ ² H ₂₁ H ₃ O	5x1ml	
Butyric Acid				
CAS 107-92-6 DRE-C10931750	MW 88.1051 Butyric acid (Butanoic acid)	C ₄ H ₈ O ₂	1ml	
Butyric Acid Butyl Ester				
CAS 109-21-7 DRE-C10931760	MW 144.2114 Butyric acid-butyl ester	C ₈ H ₁₆ O ₂	250mg	
Butyric Acid Ethyl Ester				
CAS 105-54-4 DRE-CA10931770	MW 116.1583 Butyric acid-ethyl ester	C ₆ H ₁₂ O ₂	250mg	
Butyric Acid Methyl Ester				
CAS 623-42-7 DRE-CA10931780	MW 102.1317 Butyric acid-methyl ester	C ₅ H ₁₀ O ₂	1ml	
Butyric Acid Propyl Ester				
CAS 105-66-8 DRE-C10931790	MW 130.1849 Butyric acid-propyl ester	C ₇ H ₁₄ O ₂	250mg	
Calcifediol Monohydrate (Vitamin D3 25-hydroxy monohydrate)				
CAS 63283-36-3 DRE-CA17924100	MW 418.6523 Vitamin D3 25-hydroxy monohydrate	C ₂₇ H ₄₄ O ₂ ·H ₂ O	50mg	
Canthaxanthine				
CAS 514-78-3 DRE-CA10947000	MW 564.8397 Canthaxanthine (*)	C ₄₀ H ₅₂ O ₂	150mg	










Food related compounds

Product code	Description			
Caprylic Acid (Octanoic acid)				
CAS 124-07-2 DRE-CA15711050	MW 144.2114 Octanoic acid	$C_8H_{16}O_2$	250mg	
Capsaicin				
CAS 404-86-4 DRE-C10949000	MW 305.4119 Capsaicin	$C_{18}H_{27}NO_3$	100mg	
D(+)-Cellobiose				
CAS 528-50-7 DRE-C11067000	MW 342.2965 D(+)-Cellobiose	$C_{12}H_{22}O_{11}$	250mg	
Chlorogenic Acid				
CAS 327-97-9 DRE-C11415750	MW 354.3087 Chlorogenic acid	$C_{16}H_{18}O_9$	100mg	
Cholecalciferol (Vitamin D3)				
CAS 67-97-0 DRE-CA17924000	MW 384.6377 Vitamin D3(‡)	$C_{27}H_{44}O$	100mg	
Cholesterol				
CAS 57-88-5 DRE-C11665400	MW 386.6535 Cholesterol	$C_{27}H_{46}O$	250mg	
trans-Cinnamyl Alcohol (3-Phenyl-2-propen-1-ol)				
CAS 4407-36-7 DRE-C11667500	MW 134.1751 trans-Cinnamyl alcohol (3-Phenyl-2-propen-1-ol)	$C_9H_{10}O$	250mg	
Cyanocobalamin (Vitamin B12)				
CAS 68-19-9 DRE-C11798500	MW 1355.3652 Cyanocobalamin (Vitamin B12)	$C_{63}H_{88}CoN_{14}O_{14}P$	50mg	
1,10-Decanedicarboxylic Acid				
CAS 693-23-2 DRE-C12095020	MW 230.3007 1,10-Decanedicarboxylic acid	$C_{12}H_{22}O_4$	250mg	










Food related compounds

Product code	Description			
Decanoic Acid (Capric acid)				
CAS 334-48-5 DRE-C12095050	MW 172.2646 Decanoic acid (Capric acid)	$C_{10}H_{20}O_2$	250mg	
Decanoic Acid Ethyl Ester				
CAS 110-38-3 DRE-C12095070	MW 200.3178 Decanoic acid ethyl ester	$C_{12}H_{24}O_2$	250mg	
Decanoic Acid Methyl Ester				
CAS 110-42-9 DRE-C12095090	MW 186.2912 Decanoic acid methyl ester	$C_{11}H_{22}O_2$	1ml	
Dexpanthenol (D-Panthenol)				
CAS 81-13-0 DRE-C15844500	MW 205.2515 D-Panthenol	$C_9H_{19}NO_4$	500mg	
7,8-Didehydrocholesterol (Provitamin D3)				
CAS 434-16-2 DRE-C12115000	MW 384.6377 7-Dehydrocholesterol (Provitamin D3)	$C_{27}H_{44}O$	100mg	
Didodecyl 3,3'-thiodipropionate (3,3'-Thiodipropionic acid didodecyl ester)				
CAS 123-28-4 DRE-C17492300	MW 514.8441 3,3'-Thiodipropionic acid didodecyl ester	$C_{30}H_{58}O_4S$	250mg	
Dimethyl Fumarate (Fumaric acid dimethyl ester)				
CAS 624-49-7 DRE-C13955800	MW 144.1253 Fumaric acid dimethyl ester(†)	$C_8H_8O_4$	250mg	
Dimethylacetamide (Acetic acid dimethylamide)				
CAS 127-19-5 DRE-C10015900	MW 87.1204 Acetic acid-dimethylamide	C_4H_9NO	1ml	
Docosanoic Acid (Behenic acid)				
CAS 112-85-6 DRE-C13057800	MW 340.5836 Docosanoic acid (Behenic acid)	$C_{22}H_{44}O_2$	100mg	










Food related compounds

Product code	Description			
Docosanoic Acid Ethyl Ester				
CAS 5908-87-2 DRE-C13057820	MW 368.6367 Docosanoic acid ethyl ester	$C_{24}H_{48}O_2$	100mg	
Docosanoic Acid Methyl Ester				
CAS 929-77-1 DRE-C13057840	MW 354.6101 Docosanoic acid-methyl ester	$C_{23}H_{46}O_2$	100mg	
Dodecanoic Acid (Lauric acid)				
CAS 143-07-7 DRE-C13060400	MW 200.3178 Dodecanoic acid (Lauric acid)	$C_{12}H_{24}O_2$	250mg	
Dodecanoic Acid Ethyl Ester				
CAS 106-33-2 DRE-C13060500	MW 228.3709 Dodecanoic acid ethyl ester	$C_{14}H_{28}O_2$	250mg	
Dodecanoic Acid Methyl Ester				
CAS 111-82-0 DRE-C13060600	MW 214.3443 Dodecanoic acid-methyl ester	$C_{13}H_{26}O_2$	1ml	
11,14-Eicosadienoic Acid Methyl Ester				
CAS 2463-02-7 DRE-L13112550ME	MW 322.5252 11,14-Eicosadienoic acid methyl ester 10 µg/mL in Methanol	$C_{21}H_{38}O_2$	10ml	
Eicosanoic Acid (Arachidic acid)				
CAS 506-30-9 DRE-C13112800	MW 312.5304 Eicosanoic acid (Arachidic acid)	$C_{20}H_{40}O_2$	100mg	
Eicosanoic Acid Methyl Ester				
CAS 1120-28-1 DRE-C13112840	MW 326.557 Eicosanoic acid-methyl ester	$C_{21}H_{42}O_2$	100mg	
cis-8-Eicosenoic Acid				
CAS 76261-96-6 DRE-L13113500ME	MW 310.5145 cis-8-Eicosenoic acid 10 µg/mL in Methanol	$C_{20}H_{38}O_2$	10ml	










Food related compounds

Product code	Description			
Elaidic Acid				
CAS 112-79-8 DRE-C13115000	MW 282.4614 Elaidic acid	$C_{18}H_{34}O_2$	100mg	
Elaidic Acid Methyl Ester				
CAS 1937-62-8 DRE-C13115060	MW 296.4879 Elaidic acid methyl ester	$C_{19}H_{36}O_2$	100mg	
Ergocalciferol (Vitamin D2)				
CAS 50-14-6 DRE-C17923900	MW 396.6484 Vitamin D2 (*)	$C_{28}H_{44}O$	250mg	
Ergosterol (Provitamin D2)				
CAS 57-87-4 DRE-C13201500	MW 396.6484 Ergosterol (Provitamin D2)	$C_{28}H_{44}O$	100mg	
Erucic acid				
CAS 112-86-7 DRE-C13203000	MW 338.5677 Erucic acid	$C_{22}H_{42}O_2$	100mg	
Ethyl Parahydroxybenzoate (4-Hydroxybenzoic acid ethyl ester)				
CAS 120-47-8 DRE-C14228800	MW 166.1739 4-Hydroxybenzoic acid ethyl ester(‡)	$C_9H_{10}O_3$	250mg	
Ethyl Pentanoate (Valeric acid ethyl ester)				
CAS 539-82-2 DRE-CA17899700	MW 130.1849 Valeric acid ethyl ester (Pentanoic acid ethyl ester)	$C_7H_{14}O_2$	250mg	
Formic Acid (Methanoic acid)				
CAS 64-18-6 DRE-C13913000	MW 46.0254 Formic acid (Methanoic acid)	CH_2O_2	1ml	
Formic Acid Butyl Ester				
CAS 592-84-7 DRE-C13913300	MW 102.1317 Formic acid butyl ester	$C_5H_{10}O_2$	250mg	










Food related compounds

Product code	Description			
Formic acid-tert-butyl Ester				
CAS 762-75-4 DRE-C13913305	MW 102.1317 Formic acid-tert-butyl ester	C ₅ H ₁₀ O ₂	100mg	
Formic Acid Ethyl Ester				
CAS 109-94-4 DRE-C13913400	MW 74.0785 Formic acid ethyl ester	C ₃ H ₆ O ₂	250mg	
Formic Acid Propyl Ester				
CAS 110-74-7 DRE-CA13913800	MW 88.1051 Formic acid-propyl ester	C ₄ H ₈ O ₂	1ml	
Fructose (D-(-)-Fructose)				
CAS 57-48-7 DRE-C13947500	MW 180.1559 D(-)-Fructose(‡)	C ₆ H ₁₂ O ₆	250mg	
Fumaric Acid				
CAS 110-17-8 DRE-C13955500	MW 116.0722 Fumaric acid	C ₄ H ₄ O ₄	250mg	
Galactose (D-(+)-Galactopyranose)				
CAS 59-23-4 DRE-C13996500	MW 180.1559 D(+)-Galactose	C ₆ H ₁₂ O ₆	250mg	
Gallic Acid Propyl Ester (Propyl Gallate)				
CAS 121-79-9 DRE-C13998300	MW 212.1993 Gallic acid-propyl ester	C ₁₀ H ₁₂ O ₅	250mg	
Gallic Acid				
CAS 149-91-7 DRE-C13998280	MW 170.1195 Gallic acid(‡)	C ₇ H ₆ O ₅	250mg	
Glucose Anhydrous (D-(+)-Glucopyranose)				
CAS 50-99-7 DRE-C14027000	MW 180.1559 D(+)-Glucose(‡)	C ₆ H ₁₂ O ₆	250mg	










Food related compounds

Product code	Description			
Glutamine (L-Glutamine)				
CAS 56-85-9 DRE-C14034450	MW 146.1445 L-Glutamine	C ₅ H ₁₀ N ₂ O ₃	100mg	
Glutaric Acid (Pentanedioic acid)				
CAS 110-94-1 DRE-C14034500	MW 132.1146 Glutaric acid (Pentanedioic acid)	C ₅ H ₈ O ₄	250mg	
Glutaric Acid Dimethyl Ester				
CAS 1119-40-0 DRE-C14034900	MW 160.1678 Glutaric acid dimethyl ester	C ₇ H ₁₂ O ₄	250mg	
Glycine				
CAS 56-40-6 DRE-C14037000	MW 75.0666 Glycine	C ₂ H ₅ NO ₂	250mg	
Guanosine-5'-monophosphate				
CAS 85-32-5 DRE-C14057150	MW 363.2206 Guanosine-5'-monophosphate	C ₁₀ H ₁₄ N ₅ O ₈ P	100mg	
Heneicosanoic Acid				
CAS 2363-71-5 DRE-C14085100	MW 326.557 Heneicosanoic acid	C ₂₁ H ₄₂ O ₂	100mg	
Heneicosanoic Acid Methyl Ester				
CAS 6064-90-0 DRE-C14085150	MW 340.5836 Heneicosanoic acid methyl ester	C ₂₂ H ₄₄ O ₂	100mg	
Heneicosanoic acid-ethyl ester				
CAS 28898-67-1 DRE-C14085130	MW 354.6101 Heneicosanoic acid-ethyl ester	C ₂₃ H ₄₆ O ₂	100mg	
1-Heneicosanol				
CAS 15594-90-8 DRE-C14085200	MW 312.5735 1-Heneicosanol	C ₂₁ H ₄₄ O	100mg	










Food related compounds

Product code	Description			
Heptacosanoic Acid				
CAS 7138-40-1 DRE-L14122400CY	MW 410.7165 Heptacosanoic acid 10 µg/mL in Cyclohexane	C ₂₇ H ₅₄ O ₂	10ml	
Heptadecanoic Acid				
CAS 506-12-7 DRE-C14122600	MW 270.4507 Heptadecanoic acid	C ₁₇ H ₃₄ O ₂	100mg	
Heptadecanoic Acid Methyl Ester				
CAS 1731-92-6 DRE-C14122650 DRE-L14122650CY	MW 284.4772 Heptadecanoic acid-methyl ester Heptadecanoic acid methyl ester 10 µg/mL in Cyclohexane	C ₁₈ H ₃₆ O ₂	100mg 10ml	
cis-10-Heptadecenoic Acid				
CAS 29743-97-3 DRE-L14122800ME	MW 268.4348 cis-10-Heptadecenoic acid 10 µg/mL in Methanol	C ₁₇ H ₃₂ O ₂	10ml	
Heptanoic Acid (Heptylic acid)				
CAS 111-14-8 DRE-C14126700	MW 130.1849 Heptanoic acid (Heptylic acid)	C ₇ H ₁₄ O ₂	1ml	
Heptanoic Acid Ethyl Ester				
CAS 106-30-9 DRE-C14126800	MW 158.238 Heptanoic acid ethyl ester	C ₉ H ₁₈ O ₂	250mg	
Heptanoic Acid Methyl Ester				
CAS 106-73-0 DRE-C14126900	MW 144.2114 Heptanoic acid methyl ester	C ₈ H ₁₆ O ₂	1ml	
Hesperidin				
CAS 520-26-3 DRE-C14140000	MW 610.5606 Hesperidin	C ₂₈ H ₃₄ O ₁₅	250mg	
Hexacosanoic Acid (Cerotic acid)				
CAS 506-46-7 DRE-C14191300	MW 396.6899 Hexacosanoic acid (Cerotic acid)	C ₂₆ H ₅₂ O ₂	20mg	










Food related compounds

Product code	Description			
Hexacosanoic Acid Methyl Ester				
CAS 5802-82-4	MW 410.7165	C ₂₇ H ₅₄ O ₂		
DRE-L14191330MB	Hexacosanoic acid (Cerotic acid) methyl ester 10 µg/mL in Methyl-tert. butyl ether		10ml	
Hexadecanoic Acid Ethyl Ester				
CAS 628-97-7	MW 284.4772	C ₁₈ H ₃₆ O ₂		
DRE-C14192000	Hexadecanoic acid ethyl ester		100mg	
Hexadecanoic Acid Methyl Ester				
CAS 112-39-0	MW 270.4507	C ₁₇ H ₃₄ O ₂		
DRE-C14192100	Hexadecanoic acid methyl ester		250mg	
1,6-Hexanedicarboxylic Acid (Octanedioic acid, Suberic acid)				
CAS 505-48-6	MW 174.1944	C ₈ H ₁₄ O ₄		
DRE-C14195530	1,6-Hexanedicarboxylic acid (Octanedioic acid,Suberic acid)		250mg	
Hexanoic Acid (Caproic acid)				
CAS 142-62-1	MW 116.1583	C ₆ H ₁₂ O ₂		
DRE-C14196000	Hexanoic acid (Caproic acid)(±)		250mg	
Hexanoic Acid Ethyl Ester				
CAS 123-66-0	MW 144.2114	C ₈ H ₁₆ O ₂		
DRE-C14196200	Hexanoic acid ethyl ester		250mg	
Hexanoic Acid Methyl Ester				
CAS 106-70-7	MW 130.1849	C ₇ H ₁₄ O ₂		
DRE-C14196400	Hexanoic acid methyl ester		1ml	
4-Hydroxybenzoic Acid n-Heptyl Ester				
CAS 1085-12-7	MW 236.3068	C ₁₄ H ₂₀ O ₃		
DRE-C14228820	4-Hydroxybenzoic acid-n-heptyl ester		100mg	
6-Hydroxynicotinic Acid (Nicotinic acid-6-hydroxy)				
CAS 5006-66-6	MW 139.1088	C ₆ H ₅ NO ₃		
DRE-C15521080	Nicotinic acid-6-hydroxy		100mg	










Food related compounds

Product code	Description			
D(-)-Isoascorbic Acid				
CAS 89-65-6 DRE-C14379000	MW 176.1241 D(-)-Isoascorbic acid	C ₆ H ₈ O ₆	250mg	
Isobutyl Parahydroxybenzoate (4-Hydroxybenzoic acid isobutyl ester)				
CAS 4247-02-3 DRE-C14228900	MW 194.2271 4-Hydroxybenzoic acid-iso-butyl ester(±)	C ₁₁ H ₁₄ O ₃	100mg	
Isobutyl acetate (Acetic acid isobutyl ester)				
CAS 110-19-0 DRE-C10016100	MW 116.1583 Acetic acid isobutyl ester	C ₆ H ₁₂ O ₂	1ml	
Isoeugenol				
CAS 97-54-1 DRE-C14415000	MW 164.2011 Isoeugenol	C ₁₀ H ₁₂ O ₂	250mg	
D-Lactitol Monohydrate				
CAS 81025-04-9 DRE-C14588000	MW 362.3276 D-Lactitol monohydrate	C ₁₂ H ₂₄ O ₁₁ ·H ₂ O	100mg	
α-D-(+)-Lactose				
CAS 63-42-3 DRE-C14590800	MW 342.2965 α-D-(+)-Lactose(±)	C ₁₂ H ₂₂ O ₁₁	250mg	
Limonene				
CAS 138-86-3 DRE-GA09010047IP	MW 136.234 Limonene 1000 µg/mL in Isopropanol(±)	C ₁₀ H ₁₆	1ml	
R-(+)-Limonene				
CAS 5989-27-5 DRE-C14634100	MW 136.234 (R)-(+)-Limonene	C ₁₀ H ₁₆	250mg	
S-(-)-Limonene				
CAS 5989-54-8 DRE-C14634200	MW 136.234 (S)-(-)-Limonene	C ₁₀ H ₁₆	250mg	









Food related compounds

Product code	Description			
Linalol (Linalool)				
CAS 78-70-6	MW 154.2493	C ₁₀ H ₁₈ O		
DRE-C14634500	Linalool		250mg	
DRE-GA09010048IP	Linalool 1000 µg/mL in Isopropanol(±)		1ml	
Linoelaidic Methyl Ester				
CAS 2566-97-4	MW 294.4721	C ₁₉ H ₃₄ O ₂		
DRE-L14635550ME	Linoelaidic methyl ester 10 µg/mL in Methanol		10ml	
Linoleic Acid				
CAS 60-33-3	MW 280.4455	C ₁₈ H ₃₂ O ₂		
DRE-C14635400	Linoleic acid		100mg	
Linoleic Acid Methyl Ester				
CAS 112-63-0	MW 294.4721	C ₁₉ H ₃₄ O ₂		
DRE-C14635450	Linoleic acid methyl ester(±)		100mg	
Linolelaidic Acid				
CAS 506-21-8	MW 280.4455	C ₁₈ H ₃₂ O ₂		
DRE-L14635500ME	Linolelaidic 10 µg/mL in Methanol		10ml	
Linolenic Acid Methyl Ester				
CAS 301-00-8	MW 292.4562	C ₁₉ H ₃₂ O ₂		
DRE-C14635900	Linolenic acid methyl ester(±)		100mg	
DRE-L14635900ME	Linolenic acid methyl ester 10 µg/mL in Methanol		10ml	
Linolenic acid				
CAS 463-40-1	MW 278.4296	C ₁₈ H ₃₀ O ₂		
DRE-C14635600	Linolenic acid(±)		100mg	
DRE-L14635600ME	Linolenic acid 10 µg/mL in Methanol		10ml	
Malonic Acid (Propanedioic acid)				
CAS 141-82-2	MW 104.0615	C ₃ H ₄ O ₄		
DRE-C14731000	Malonic acid (Propanedioic acid)		250mg	
Malonic Acid Dimethyl Ester				
CAS 108-59-8	MW 132.1146	C ₅ H ₈ O ₄		
DRE-C14731500	Malonic acid dimethyl ester		250mg	










Food related compounds

Product code	Description			
D-(+)-Maltose Monohydrate				
CAS 6363-53-7 DRE-C14734700	MW 360.3118 D-(+)-Maltose hydrate	C ₁₂ H ₂₂ O ₁₁ ·H ₂ O	250mg	
Maltotriose Hydrate				
CAS 207511-08-8 DRE-C14735000	MW 522.4524 Maltotriose hydrate	C ₁₈ H ₃₂ O ₁₆ ·H ₂ O	100mg	
Mannitol (D-Mannit)				
CAS 69-65-8 DRE-C14752000	MW 182.1718 D-Mannit	C ₆ H ₁₄ O ₆	250mg	
D(+)-Mannose				
CAS 3458-28-4 DRE-C14752300	MW 180.1559 D(+)-Mannose	C ₆ H ₁₂ O ₆	250mg	
Melezitose Hydrate				
CAS 10030-67-8 DRE-C14861800	MW 522.4524 Melezitose hydrate	C ₁₈ H ₃₂ O ₁₆ ·H ₂ O	250mg	
α-D(+)-Melibiose Hydrate				
CAS 585-99-9 DRE-C14862000	MW 360.3118 α-D(+)-Melibiose Hydrate	C ₁₂ H ₂₂ O ₁₁ ·H ₂ O	250mg	
Menadione (Vitamin K3)				
CAS 58-27-5 DRE-C14863000	MW 172.18 Menadione (Vitamin K3)	C ₁₁ H ₈ O ₂	250mg	
Menadione Sodium Bisulfite				
CAS 130-37-0 DRE-C14863200	MW 276.2409 Menadione sodium bisulfite	C ₁₁ H ₉ O ₅ S·Na	250mg	
Methyl Acetate (Acetic acid methyl ester)				
CAS 79-20-9 DRE-C10016300	MW 74.0785 Acetic acid-methyl ester(‡)	C ₃ H ₆ O ₂	1ml	










Food related compounds

Product code	Description			
Methyl Erucate (Erucic acid methyl ester)				
CAS 1120-34-9	MW 352.5943	C ₂₃ H ₄₄ O ₂		
DRE-C13203060	Erucic acid methyl ester		100mg	
DRE-L13203060ME	Erucic acid methyl ester 10 µg/mL in Methanol		10ml	
Methyl Formate (Formic acid methyl ester)				
CAS 107-31-3	MW 60.052	C ₂ H ₄ O ₂		
DRE-CA13913600	Formic acid-methyl ester		1ml	
Methyl Parahydroxybenzoate (4-Hydroxybenzoic acid methyl ester)				
CAS 99-76-3	MW 152.1473	C ₉ H ₈ O ₃		
DRE-C14229000	4-Hydroxybenzoic acid methyl ester(±)		250mg	
Methyl Stearate (Octadecanoic acid methyl ester)				
CAS 112-61-8	MW 298.5038	C ₁₉ H ₃₈ O ₂		
DRE-C15710180	Octadecanoic acid methyl ester		1g	
Methyl Tricosanoate (Tricosanoic acid methyl ester)				
CAS 2433-97-8	MW 368.6367	C ₂₄ H ₄₈ O ₂		
DRE-C17805550	Tricosanoic acid methyl ester		100mg	
DRE-L17805550CY	Tricosanoic acid methyl ester 10 µg/mL in Cyclohexane		10ml	
Monascus Red				
CAS 874807-57-5	MW n/a			
DRE-E15290500	Monascus Red		100mg	No Structure
Myristoleic Acid				
CAS 544-64-9	MW 226.355	C ₁₄ H ₂₆ O ₂		
DRE-L15392000CY	Myristoleic acid 10 µg/mL in Cyclohexane		10ml	
Myristoleic acid methyl ester				
CAS 56219-06-8	MW 240.3816	C ₁₅ H ₂₈ O ₂		
DRE-L15392050CY	Myristoleic acid methyl ester 10 µg/mL in Cyclohexane		10ml	
Naringin				
CAS 10236-47-2	MW 580.5346	C ₂₇ H ₃₂ O ₁₄		
DRE-C15495000	Naringin		250mg	










Food related compounds

Product code	Description			
Nervonic Acid Methyl Ester				
CAS 2733-88-2 DRE-L15504040CY	MW 380.6474 Nervonic acid methyl ester 10 µg/mL in Cyclohexane	C ₂₅ H ₄₈ O ₂	10ml	
Nicotinamide				
CAS 98-92-0 DRE-C15519500	MW 122.1246 Nicotinamide(±)	C ₆ H ₆ N ₂ O	250mg	
Nicotinic Acid				
CAS 59-67-6 DRE-C15521000	MW 123.1094 Nicotinic acid(±)	C ₆ H ₅ NO ₂	250mg	
Nonadecanoic Acid (Nonadecyclic acid)				
CAS 646-30-0 DRE-C15622300	MW 298.5038 Nonadecanoic acid (Nonadecyclic acid)	C ₁₉ H ₃₈ O ₂	250mg	
Nonadecanoic Acid Ethyl Ester				
CAS 18281-04-4 DRE-L15622330MB	MW 326.557 Nonadecanoic acid ethyl ester 10 µg/mL in Methyl-tert. butyl ether	C ₂₁ H ₄₂ O ₂	10ml	
Nonadecanoic Acid Methyl Ester				
CAS 1731-94-8 DRE-C15622360 DRE-L15622360ME	MW 312.5304 Nonadecanoic acid methyl ester Nonadecanoic acid methyl ester 10 µg/mL in Methanol	C ₂₀ H ₄₀ O ₂	250mg 10ml	
Nonanoic Acid Ethyl Ester				
CAS 123-29-5 DRE-C15623130	MW 186.2912 Nonanoic acid ethyl ester	C ₁₁ H ₂₂ O ₂	250mg	
Nonanoic Acid Methyl Ester				
CAS 1731-84-6 DRE-C15623150	MW 172.2646 Nonanoic acid methyl ester	C ₁₀ H ₂₀ O ₂	1ml	
Nonivamide (N-Vanillylnonanamide)				
CAS 2444-46-4 DRE-C17900600	MW 293.4012 N-Vanillylnonanamide	C ₁₇ H ₂₇ NO ₃	100mg	










Food related compounds

Product code	Description			
Octacosanoic Acid (Montanic acid)				
CAS 506-48-9 DRE-C15710060	MW 424.743 Octacosanoic acid (Montanic acid)	$C_{28}H_{56}O_2$	50mg	
Octadecanedioic Acid				
CAS 871-70-5 DRE-C15710110	MW 314.4602 Octadecanedioic acid	$C_{18}H_{34}O_4$	100mg	
Octadecanoic Acid Ethyl Ester				
CAS 111-61-5 DRE-C15710150	MW 312.5304 Octadecanoic acid ethyl ester	$C_{20}H_{40}O_2$	100mg	
trans-11-Octadecanoic Acid Methyl Ester				
CAS 6198-58-9 DRE-L15710470CY	MW 296.4879 trans-11-Octadecanoic acid methyl ester 10 µg/mL in Cyclohexane	$C_{19}H_{36}O_2$	10ml	
trans-11-Octadecanoic Acid (Vaccenic acid)				
CAS 693-72-1 DRE-L15710430CY	MW 282.4614 trans-11-Octadecanoic acid (Vaccenic acid) 10 µg/mL in Cyclohexane	$C_{18}H_{34}O_2$	10ml	
cis-11-Octadecenoic acid methyl ester				
CAS 1937-63-9 DRE-C15710460	MW 296.4879 cis-11-Octadecenoic acid methyl ester	$C_{19}H_{36}O_2$	10mg	
1,8-Octanedicarboxylic Acid Dimethyl Ester				
CAS 106-79-6 DRE-C15711030	MW 230.3007 1,8-Octanedicarboxylic acid-bis-methyl ester	$C_{12}H_{22}O_4$	250mg	
Octanoic Acid Ethyl Ester				
CAS 106-32-1 DRE-C15711070	MW 172.2646 Octanoic acid ethyl ester	$C_{10}H_{20}O_2$	250mg	
Octanoic Acid Methyl Ester				
CAS 111-11-5 DRE-C15711090	MW 158.238 Octanoic acid methyl ester	$C_9H_{18}O_2$	1ml	










Food related compounds

Product code	Description			
Orotic Acid				
CAS 65-86-1 DRE-C15747000	MW 156.0963 Orotic acid	$C_5H_4N_2O_4$	100mg	
Orotic acid monohydrate				
CAS 50887-69-9 DRE-C15747008	MW 174.1115 Orotic acid monohydrate	$C_5H_4N_2O_4 \cdot H_2O$	100mg	
Oxalic Acid (Ethanedioic acid)				
CAS 144-62-7 DRE-C15775000	MW 90.0349 Oxalic acid (Ethanedioic acid)	$C_2H_2O_4$	250mg	
Palmitelaidic acid methyl ester				
CAS 10030-74-7 DRE-C15843000	MW 268.4348 Palmitelaidic acid methyl ester	$C_{17}H_{32}O_2$	10mg	
Palmitic Acid (Hexadecanoic acid)				
CAS 57-10-3 DRE-C14191800	MW 256.4241 Hexadecanoic acid (Palmitic acid)(±)	$C_{16}H_{32}O_2$	100mg	
Palmitic Acid n-Butyl Ester				
CAS 111-06-8 DRE-C15843100	MW 312.5304 Palmitic acid-n-butyl ester	$C_{20}H_{40}O_2$	100mg	
Palmitoleic Acid				
CAS 373-49-9 DRE-CA15843200	MW 254.4082 Palmitoleic acid	$C_{16}H_{30}O_2$	100mg	
DL-Pantothenic Acid Calcium Salt				
CAS 6381-63-1 DRE-CA15845000	MW 476.5321 Pantothenic acid calcium salt	$2C_9H_{16}NO_5 \cdot Ca$	250mg	
D-Pantothenic Acid Sodium Salt				
CAS 867-81-2 DRE-C15845100	MW 241.2168 D-Pantothenic acid sodium salt	$C_9H_{16}NO_5 \cdot Na$	250mg	










Food related compounds

Product code	Description			
Pentacosanoic Acid				
CAS 506-38-7 DRE-C15973630	MW 382.6633 Pentacosanoic acid	$C_{25}H_{50}O_2$	100mg	
1-Pentacosanol				
CAS 26040-98-2 DRE-L15973680MB	MW 368.6798 1-Pentacosanol 10 µg/mL in Methyl-tert. butyl ether	$C_{25}H_{52}O$	10ml	
Pentadecanoic Acid (Pentadecylic acid)				
CAS 1002-84-2 DRE-C15973740	MW 242.3975 Pentadecanoic acid (Pentadecylic acid)	$C_{15}H_{30}O_2$	100mg	
Pentadecanoic Acid Ethyl Ester				
CAS 41114-00-5 DRE-C15973760	MW 270.4507 Pentadecanoic acid ethyl ester	$C_{17}H_{34}O_2$	100mg	
Pentadecanoic Acid Methyl Ester				
CAS 7132-64-1 DRE-C15973770	MW 256.4241 Pentadecanoic acid methyl ester	$C_{16}H_{32}O_2$	100mg	
1,5-Pentanedicarboxylic Acid				
CAS 111-16-0 DRE-C15977600	MW 160.1678 1,5-Pentanedicarboxylic acid	$C_7H_{12}O_4$	250mg	
Phytomenadione (Vitamin K1)				
CAS n/a DRE-C17924400	MW 450.6957 Vitamin K1(‡)	$C_{31}H_{46}O_2$	250mg	
Potassium Sorbate				
CAS 24634-61-5 DRE-C16972000	MW 150.2169 trans-trans-Sorbic acid potassium salt(‡)	$C_6H_7O_2 \cdot K$	250mg	
Propanoic Acid Ethyl Ester				
CAS 105-37-3 DRE-CA16493500	MW 102.1317 Propionic acid-ethyl ester	$C_5H_{10}O_2$	250mg	










Food related compounds

Product code	Description			
Propionic Acid (Propanoic acid)				
CAS 79-09-4 DRE-C16493000	MW 74.0785 Propionic acid (Propanoic acid)	C ₃ H ₆ O ₂	1ml	
Propionic Acid Butyl Ester				
CAS 590-01-2 DRE-C16493400	MW 130.1849 Propionic acid butyl ester	C ₇ H ₁₄ O ₂	250mg	
Propionic Acid Calcium Salt				
CAS 4075-81-4 DRE-C16493100	MW 186.2192 Propionic acid calcium salt	2C ₃ H ₅ O ₂ ·Ca	250mg	
Propionic Acid Methyl Ester				
CAS 554-12-1 DRE-C16493600	MW 88.1051 Propionic acid methyl ester	C ₄ H ₈ O ₂	1ml	
Propyl Acetate (Acetic acid n-propyl ester)				
CAS 109-60-4 DRE-C10016500	MW 102.1317 Acetic acid n-propyl ester(‡)	C ₅ H ₁₀ O ₂	1g	
Propyl Parahydroxybenzoate (4-Hydroxybenzoic acid propyl ester)				
CAS 94-13-3 DRE-C14229200	MW 180.2005 4-Hydroxybenzoic acid propyl ester(‡)	C ₁₀ H ₁₂ O ₃	250mg	
Pyridoxine Hydrochloride				
CAS 58-56-0 DRE-C16652000	MW 205.6388 Pyridoxin hydrochloride (Vitamin B6-HCl)(‡)	C ₈ H ₁₁ NO ₃ ·ClH	250mg	
D-(+)-Raffinose Pentahydrate				
CAS 17629-30-0 DRE-C16806000	MW 594.5135 D-(+)-Raffinose pentahydrate	C ₁₈ H ₃₂ O ₁₆ ·5H ₂ O	250mg	
Retinol Palmitate (Vitamin A palmitate)				
CAS 79-81-2 DRE-CA17923840	MW 524.8604 Vitamin A palmitate	C ₃₆ H ₆₀ O ₂	100mg	

Food related compounds

Product code	Description			
Retinyl Acetate (Vitamin A acetate)				
CAS 127-47-9 DRE-CA17923820	MW 328.4883 Vitamin A acetate(‡)	C ₂₂ H ₃₂ O ₂	100mg	
L(+)-Rhamnose Hydrate				
CAS 10030-85-0 DRE-C16813500	MW 182.1718 L(+)-Rhamnose Hydrate	C ₆ H ₁₂ O ₅ ·H ₂ O	250mg	
(-)-Riboflavin (Vitamin B2)				
CAS 83-88-5 DRE-C16813600	MW 376.3639 (-)-Riboflavin (Vitamin B2)(‡)	C ₁₇ H ₂₀ N ₄ O ₆	250mg	
Riboflavin Sodium Phosphate				
CAS 130-40-5 DRE-C16813610	MW 478.3256 Riboflavine-5'-phosphate sodium	C ₁₇ H ₂₀ N ₄ O ₉ P·Na	250mg	
D-(-)-Ribose				
CAS 50-69-1 DRE-C16813700	MW 150.1299 D(-)-Ribose	C ₅ H ₁₀ O ₅	250mg	
Ricinoleic acid methyl ester				
CAS 141-24-2 DRE-C16814070	MW 312.4873 Ricinoleic acid methyl ester	C ₁₉ H ₃₆ O ₃	50mg	
Saccharin				
CAS 81-07-2 DRE-C16901000	MW 183.1845 Saccharin(‡)	C ₇ H ₅ NO ₃ S	250mg	
Saccharin Sodium Salt Hydrate				
CAS 128-44-9 DRE-C16901010	MW 205.1663 Saccharin sodium salt hydrate(‡)	C ₇ H ₄ NO ₃ S·Na	250mg	
Sodium Ascorbate (Ascorbic acid sodium salt; Vitamin C sodium salt)				
CAS 134-03-2 DRE-C10303900	MW 198.1059 Ascorbic acid sodium salt (Vitamin C sodium salt)	C ₆ H ₇ O ₆ ·Na	250mg	










Food related compounds

Product code	Description			
Sodium Benzoate (Benzoic acid sodium salt)				
CAS 532-32-1 DRE-C10538000	MW 144.1032 Benzoic acid sodium salt	$C_7H_5O_2 \cdot Na$	250mg	
Sodium Cyclamate (N-Cyclohexylsulfamic acid sodium salt)				
CAS 139-05-9 DRE-C11830800	MW 201.2192 N-Cyclohexylsulfamic acid sodium salt(±)	$C_6H_{12}NO_3S \cdot Na$	250mg	
Sodium Propionate (Propionic acid sodium salt)				
CAS 137-40-6 DRE-C16493300	MW 96.0604 Propionic acid sodium salt	$C_3H_5O_2 \cdot Na$	250mg	
Sorbic Acid				
CAS 110-44-1 DRE-C16971500	MW 112.1265 Sorbic acid(±)	$C_6H_8O_2$	250mg	
D-Sorbitol				
CAS 50-70-4 DRE-C16972500	MW 182.1718 D-Sorbit(±)	$C_6H_{14}O_6$	250mg	
L-(-)-Sorbitose				
CAS 87-79-6 DRE-C16972600	MW 180.1559 L-(-)-Sorbitose	$C_6H_{12}O_6$	250mg	
Stearic Acid (Octadecanoic acid)				
CAS 57-11-4 DRE-C15710120	MW 284.4772 Octadecanoic acid (Stearic acid)	$C_{18}H_{36}O_2$	250mg	
Succinic Acid (Butanedioic acid)				
CAS 110-15-6 DRE-C16985000	MW 118.088 Succinic acid (Butanedioic acid)	$C_4H_6O_4$	250mg	
Succinic Acid Dibutyl Ester				
CAS 141-03-7 DRE-C16985200	MW 230.3007 Succinic acid di-n-butyl ester	$C_{12}H_{22}O_4$	100mg	










Food related compounds

Product code	Description			
Succinic Acid Diethyl Ester				
CAS 123-25-1 DRE-C16985300	MW 174.1944 Succinic acid diethyl ester	C ₈ H ₁₄ O ₄	250mg	
Succinic Acid Dimethyl Ester				
CAS 106-65-0 DRE-C16985500	MW 146.1412 Succinic acid dimethyl ester	C ₆ H ₁₀ O ₄	250mg	
Sucralose				
CAS 56038-13-2 DRE-C16985800	MW 397.6335 Sucralose(±)	C ₁₂ H ₁₉ Cl ₃ O ₈	100mg	
Sucrose (D-(+)-Saccharose)				
CAS 57-50-1 DRE-C16901100	MW 342.2965 D-(+)-Saccharose(±)	C ₁₂ H ₂₂ O ₁₁	250mg	
Tetracosanoic Acid (Lignoceric acid)				
CAS 557-59-5 DRE-C17396000 DRE-L17396000CY	MW 368.6367 Tetracosanoic acid (Lignoceric acid) Tetracosanoic acid (Lignoceric acid) 10 µg/mL in Cyclohexane	C ₂₄ H ₄₈ O ₂	100mg 10ml	
Tetracosanoic Acid Methyl Ester				
CAS 2442-49-1 DRE-C17396040	MW 382.6633 Tetracosanoic acid methyl ester	C ₂₅ H ₅₀ O ₂	100mg	
1-Tetracosanol (Lignoceryl alcohol)				
CAS 506-51-4 DRE-C17396100 DRE-L17396100MB	MW 354.6532 1-Tetracosanol (Lignoceryl alcohol) 1-Tetracosanol (Lignoceryl alcohol) 10 µg/mL in Methyl-tert. butyl ether	C ₂₄ H ₅₀ O	100mg 10ml	
Tetradecanoic Acid (Myristic acid)				
CAS 544-63-8 DRE-C17396700	MW 228.3709 Tetradecanoic acid (Myristic acid)	C ₁₄ H ₂₈ O ₂	250mg	
Tetradecanoic Acid Ethyl Ester				
CAS 124-06-1 DRE-C17396730	MW 256.4241 Tetradecanoic acid ethyl ester	C ₁₆ H ₃₂ O ₂	250mg	







Food related compounds

Product code	Description			
Tetradecanoic Acid Methyl Ester				
CAS 124-10-7 DRE-C17396770	MW 242.3975 Tetradecanoic acid methyl ester	C ₁₅ H ₃₀ O ₂	1ml	
Thiamine Hydrochloride				
CAS 67-03-8 DRE-C17455000	MW 337.2685 Thiamine hydrochloride(‡)	C ₁₂ H ₁₇ N ₄ OS·Cl·ClH	250mg	
Thiamine Nitrate				
CAS 532-43-4 DRE-C17455500	MW 327.3595 Thiamine mononitrate	C ₁₂ H ₁₇ N ₄ OS·NO ₃	250mg	
Thioxothiamine (Thiothiamine)				
CAS 299-35-4 DRE-C17561000	MW 296.4116 Thiothiamine	C ₁₂ H ₁₆ N ₄ OS ₂	100mg	
DL-α-Tocopherol Acetate (Vitamin E acetate)				
CAS 52225-20-4 DRE-CA17924320	MW 472.7428 DL-α-Tocopherylacetate (Vitamin E acetate)(‡)	C ₃₁ H ₅₂ O ₃	500mg	
all-rac-α-Tocopherol (Vitamin E)				
CAS 10191-41-0 DRE-CA17924300	MW 430.7061 DL-α-Tocopherol (Vitamin E) (*)	C ₂₉ H ₅₀ O ₂	500mg	
D-(+)-Trehalose Dihydrate				
CAS 6138-23-4 DRE-C17607500	MW 378.327 D(+)-Trehalose dihydrate	C ₁₂ H ₂₂ O ₁₁ ·2H ₂ O	250mg	
Triacetin				
CAS 102-76-1 DRE-C17608500	MW 218.2039 Triacetin	C ₉ H ₁₄ O ₆	250mg	
Triacontanoic Acid Methyl Ester				
CAS 629-83-4 DRE-L17609150CY	MW 466.8228 Triacontanoic acid methyl ester 10 µg/mL in Cyclohexane	C ₃₁ H ₆₂ O ₂	10ml	

Food related compounds

Product code	Description			
Tricosanoic Acid				
CAS 2433-96-7 DRE-C17805500	MW 354.6101 Tricosanoic acid	C ₂₃ H ₄₆ O ₂	100mg	
Tricosanoic acid ethyl ester				
CAS 18281-07-7 DRE-C17805530	MW 382.6633 Tricosanoic acid ethyl ester	C ₂₅ H ₅₀ O ₂	100mg	
Tridecanal				
CAS 10486-19-8 DRE-C17817900	MW 198.3449 Tridecanal (technical)	C ₁₃ H ₂₆ O	100mg	
Tridecanoic Acid (Tridecylic acid)				
CAS 638-53-9 DRE-C17818100	MW 214.3443 Tridecanoic acid (Tridecylic acid)	C ₁₃ H ₂₆ O ₂	250mg	
Tridecanoic Acid Methyl Ester				
CAS 1731-88-0 DRE-C17818150	MW 228.3709 Tridecanoic acid methyl ester	C ₁₄ H ₂₈ O ₂	100mg	
D-(+)-Turanose				
CAS 547-25-1 DRE-C17895500	MW 342.2965 D-(+)-Turanose	C ₁₂ H ₂₂ O ₁₁	250mg	
Undecanoic Acid (Undecylic acid)				
CAS 112-37-8 DRE-C17896400	MW 186.2912 Undecanoic acid (Undecylic acid)	C ₁₁ H ₂₂ O ₂	250mg	
Undecanoic Acid Ethyl Ester				
CAS 627-90-7 DRE-C17896420	MW 214.3443 Undecanoic acid ethyl ester	C ₁₃ H ₂₆ O ₂	250mg	
Undecanoic Acid Methyl Ester				
CAS 1731-86-8 DRE-C17896450	MW 200.3178 Undecanoic acid methyl ester	C ₁₂ H ₂₄ O ₂	1ml	

Food related compounds










Product code	Description			
Valeric Acid Methyl Ester				
CAS 624-24-8 DRE-C17899750	MW 116.1583 Valeric acid methyl ester	C ₆ H ₁₂ O ₂	1ml	
n-Valeric Acid (n-Pentanoic Acid)				
CAS 109-52-4 DRE-C17899500	MW 102.1317 n-Valeric acid (n-Pentanoic acid)	C ₅ H ₁₀ O ₂	1ml	
Xylitol				
CAS 87-99-0 DRE-C17945500	MW 152.1458 Xylit(†)	C ₅ H ₁₂ O ₅	250mg	
D-(+)-Xylose				
CAS 58-86-6 DRE-C17946000	MW 150.1299 D(+)-Xylose	C ₅ H ₁₀ O ₅	250mg	
Zeaxanthin				
CAS 144-68-3 DRE-CA17947500	MW 568.8714 Zeaxanthin (*)	C ₄₀ H ₅₆ O ₂	50mg	
Zinc Picolinate				
CAS 17949-65-4 DRE-C17949520	MW 309.6119 Zinc picolinate	2C ₆ H ₄ NO ₂ ·Zn	100mg	











FOOD PACKAGING CONTAMINANTS









Food packaging contaminants

Product code	Description			
Adipic Acid Bis(2-ethylhexyl) Ester				
CAS 103-23-1	MW 370.5665	C ₂₂ H ₄₂ O ₄		
DRE-C10046000	Adipic acid-bis-2-ethylhexylester			250mg
DRE-XA10046000AC	Adipic acid-bis-2-ethylhexylester 100 µg/mL in Acetone			1ml
				
Adipic Acid Dibutyl Ester				
CAS 105-99-7	MW 258.3538	C ₁₄ H ₂₆ O ₄		
DRE-C10046100	Adipic acid-di-n-butyl ester			250mg
				
Adipic Acid Dimethyl Ester				
CAS 627-93-0	MW 174.1944	C ₈ H ₁₄ O ₄		
DRE-C10046200	Adipic acid-dimethyl ester			250mg
				
Adipic acid				
CAS 124-04-9	MW 146.1412	C ₆ H ₁₀ O ₄		
DRE-C10045900	Adipic acid (Hexandioic acid)			250mg
				
Cresyl Diphenyl Phosphate				
CAS 26444-49-5	MW 340.3096	C ₁₉ H ₁₇ O ₄ P		
DRE-C11749000	Cresyl Diphenyl Phosphate			250mg
				
Dibutyl Phthalate				
CAS 84-74-2	MW 278.3435	C ₁₆ H ₂₂ O ₄		
DRE-C16171000	Phthalic acid, bis-butyl ester(‡)			500mg
				
Di-n-butylphosphate				
CAS 107-66-4	MW 210.2078	C ₈ H ₁₉ O ₄ P		
DRE-C12256000	Di-n-butylphosphate			100mg
				
Diethyl Phthalate				
CAS 84-66-2	MW 222.2372	C ₁₂ H ₁₄ O ₄		
DRE-C16172000	Phthalic acid, bis-ethyl ester(‡)			500mg
				
Diethylhexylphthalate (Bis-(2-ethylhexyl) Phthalate)				
CAS 117-81-7	MW 390.5561	C ₂₄ H ₃₈ O ₄		
DRE-C16173000	Phthalic acid, bis-2-ethylhexyl ester(‡)			500mg
DRE-CR16173000	Phthalic acid, bis-2-ethylhexyl ester			100mg
DRE-L16173000CY	Phthalic acid, bis-2-ethylhexyl ester 10 µg/mL in Cyclohexane			10ml
				










Food packaging contaminants

Product code	Description			
2-Ethylhexyl Diphenyl Phosphate				
CAS 1241-94-7 DRE-C13342300	MW 362.3997 2-Ethylhexyl diphenyl phosphate (technical)	C ₂₀ H ₂₇ O ₄ P	250mg	
Isophthalic Acid				
CAS 121-91-5 DRE-C14447000	MW 166.1308 Isophthalic acid	C ₈ H ₆ O ₄	1g	
Isophthalic acid, bis-n-octyl ester				
CAS 4654-18-6 DRE-C14447400	MW 390.5561 Isophthalic acid, bis-n-octyl ester	C ₂₄ H ₃₈ O ₄	250mg	
Isophthalic Acid Dimethyl Ester				
CAS 1459-93-4 DRE-C14447300	MW 194.184 Isophthalic acid, bis-methyl ester	C ₁₀ H ₁₀ O ₄	1g	
Isophthalic Acid Diphenyl Ester				
CAS 744-45-6 DRE-C14447500	MW 318.3228 Isophthalic acid, bis-phenyl ester	C ₂₀ H ₁₄ O ₄	250mg	
Phthalic Acid Benzyl Butyl Ester				
CAS 85-68-7 DRE-C16168000 DRE-YA16168000ME	MW 312.3597 Phthalic acid, benzylbutyl ester(±) Phthalic acid, benzylbutyl ester 1000 µg/mL in Methanol	C ₁₉ H ₂₀ O ₄	250mg 1ml	
Phthalic Acid Benzyl Butyl Ester (3,4,5,6)-D4				
CAS 93951-88-3 DRE-C16168010 DRE-XA16168010CY	MW 316.3843 Phthalic acid, benzylbutyl ester (3,4,5,6) D4 Phthalic acid, benzylbutyl ester (3,4,5,6) D4 100 µg/mL in Cyclohexane	C ₁₉ ² H ₁₄ H ₁₆ O ₄	10mg 1ml	
Phthalic Acid Bis-dodecyl Ester				
CAS 2432-90-8 DRE-C16171850	MW 502.7688 Phthalic acid, bis-dodecyl ester	C ₃₂ H ₅₄ O ₄	250mg	
Phthalic Acid Bis-hexyl Ester (Mixture of Isomers)				
CAS 68515-50-4 DRE-C16173580	MW n/a Phthalic acid, bis-iso-hexyl ester (Mixture of Isomers)		250mg	No Structure










Food packaging contaminants

Product code	Description			
Phthalic Acid Bis-methylglycol Ester				
CAS 117-82-8 DRE-C16174400	MW 282.2891 Phthalic acid, bis-methylglycol ester	C ₁₄ H ₁₈ O ₆	250mg	
Phthalic acid, bis-C6-C10-alkyl ester				
CAS 68515-51-5 DRE-C16171100	MW n/a Phthalic acid, bis-C6-C10-alkyl ester		250mg	No Structure
Phthalic acid, bis-C6-C8-branched alkyl esters C7-rich				
CAS 71888-89-6 DRE-C16171830	MW 222.2372 Phthalic acid, bis-C6-C8-branched alkyl esters C7-rich	C ₁₀ H ₁₀ O ₄ (CH ₂) _n (CH ₂) _n	100mg	
Phthalic Acid Bis(2-butoxyethyl) Ester				
CAS 117-83-9 DRE-C16170500	MW 366.4486 Phthalic acid, bis-2-n-butoxyethyl ester	C ₂₀ H ₃₀ O ₆	100mg	
Phthalic Acid Bis(2-ethoxyethyl) Ester				
CAS 605-54-9 DRE-C16171900	MW 310.3423 Phthalic acid, bis-2-ethoxyethyl ester	C ₁₆ H ₂₂ O ₆	100mg	
Phthalic Acid Bis(2-ethylhexyl) Ester (3,4,5,6)-D4				
CAS 93951-87-2 DRE-C16173010	MW 394.5808 Phthalic acid, bis-2-ethylhexyl ester D4	C ₂₄ ² H ₃₄ H ₃₄ O ₄	25mg	
DRE-CR16173010	Phthalic acid, bis-2-ethylhexyl ester D4		25mg	
DRE-XA16173010CY	Phthalic acid, bis-2-ethylhexyl ester D4 100 µg/mL in Cyclohexane		1ml	
Phthalic acid, bis-iso-decyl ester				
CAS 26761-40-0 DRE-C16173550 DRE-L16173550CY	MW n/a Phthalic acid, bis-isodecyl ester(‡) Phthalic acid, bis-isodecyl ester 10 µg/mL in Cyclohexane		100mg 10ml	No Structure
Phthalic acid, bis-iso-heptyl ester				
CAS 41451-28-9 DRE-C16173570	MW n/a Phthalic acid, bis-iso-heptyl ester		100mg	No Structure
Phthalic Acid, Bis-isononyl Ester				
CAS 28553-12-0 DRE-C16173600 DRE-L16173600CY	MW 222.2372 Phthalic acid, bis-isononyl ester Phthalic acid, bis-isononyl ester 10 µg/mL in Cyclohexane	C ₁₀ H ₁₀ O ₄ (CH ₂) _n (CH ₂) _n	250mg 10ml	










Food packaging contaminants

Product code	Description			
Phthalic Acid Bis(isopentyl) Ester				
CAS 605-50-5 DRE-C16173680	MW 306.3966 Phthalic acid, bis-iso-pentyl ester(‡)	C ₁₈ H ₂₆ O ₄	10mg	
Phthalic Acid Bis(isopentyl) Ester (3,4,5,6)-D4				
CAS 1346597-80-5 DRE-C16173685	MW 310.4213 Phthalic acid, bis-iso-pentyl ester D4	C ₁₈ ² H ₁₄ H ₂₂ O ₄	10mg	
Phthalic acid, bis-2-methyloctyl ester				
CAS 70857-56-6 DRE-C16174600	MW 418.6093 Phthalic acid, bis-2-methyloctyl ester	C ₂₆ H ₄₂ O ₄	10mg	
Phthalic Acid Bis(4-methyl-2-pentyl) Ester				
CAS 84-63-9 DRE-C16174700	MW 334.4498 Phthalic acid, bis-4-methyl-2-pentyl ester(‡)	C ₂₀ H ₃₀ O ₄	100mg	
Phthalic Acid Bis(2-propylheptyl) Ester				
CAS 53306-54-0 DRE-C16177040	MW 446.6624 Phthalic acid, bis-2-propylheptyl ester	C ₂₈ H ₄₆ O ₄	100mg	
Phthalic Acid Bis(undecyl) Ester				
CAS 3648-20-2 DRE-C16177150	MW 474.7156 Phthalic acid, bis-n-undecyl ester(‡)	C ₃₀ H ₅₀ O ₄	100mg	
Phthalic Acid Butyl 2-Ethylhexyl Ester				
CAS 85-69-8 DRE-C16177250	MW 334.4498 Phthalic acid, butyl(2-ethylhexyl) ester	C ₂₀ H ₃₀ O ₄	25mg	
Phthalic Acid Butyl Octyl Ester				
CAS 84-78-6 DRE-C16177300	MW 334.4498 Phthalic acid, butyloctyl ester	C ₂₀ H ₃₀ O ₄	250mg	
Phthalic Acid Diallyl Ester				
CAS 131-17-9 DRE-C16169000	MW 246.2586 Phthalic acid, bis-allyl ester(‡)	C ₁₄ H ₁₄ O ₄	250mg	










Food packaging contaminants

Product code	Description			
Phthalic Acid Dibenzyl Ester				
CAS 523-31-9 DRE-C16170000	MW 346.3759 Phthalic acid, bis-benzylester	$C_{22}H_{18}O_4$	100mg	
Phthalic Acid Dibutyl Ester (3,4,5,6)-D4				
CAS 93952-11-5 DRE-C16171010	MW 282.3681 Phthalic acid, bis-butyl ester D4(‡)	$C_{16}^2H_4H_{18}O_4$	10mg	
Phthalic Acid Dicyclohexyl Ester				
CAS 84-61-7 DRE-C16171500	MW 330.418 Phthalic acid, bis-cyclohexyl ester(‡)	$C_{20}H_{26}O_4$	250mg	
Phthalic Acid Dicyclohexyl Ester D4				
CAS 358731-25-6 DRE-C16171510	MW 334.4427 Phthalic acid, bis-cyclohexyl ester D4	$C_{20}^2H_4H_{22}O_4$	10mg	
Phthalic Acid Didecyl Ester				
CAS 84-77-5 DRE-C16171810	MW 446.6624 Phthalic acid, bis-n-decyl ester	$C_{28}H_{46}O_4$	250mg	
Phthalic Acid Diethyl Ester (3,4,5,6)-D4				
CAS 93952-12-6 DRE-C16172010	MW 226.2618 Phthalic acid, bis-ethyl ester D4(‡)	$C_{12}^2H_4H_{10}O_4$	10mg	
Phthalic Acid Diheptyl Ester				
CAS 3648-21-3 DRE-C16173100	MW 362.503 Phthalic acid, bis-n-heptyl ester	$C_{22}H_{34}O_4$	100mg	
Phthalic Acid Dihexyl Ester				
CAS 84-75-3 DRE-C16173200 DRE-L16173200CY	MW 334.4498 Phthalic acid, bis-hexyl ester(‡) Phthalic acid, bis-hexyl ester 10 µg/mL in Cyclohexane	$C_{20}H_{30}O_4$	100mg 10ml	
Phthalic Acid Di-n-hexyl Ester (3,4,5,6)-D4				
CAS 1015854-55-3 DRE-C16173210	MW 338.4744 Phthalic acid, bis-hexyl ester D4	$C_{20}^2H_4H_{26}O_4$	10mg	







Food packaging contaminants

Product code	Description			
Phthalic Acid Diisobutyl Ester				
CAS 84-69-5	MW 278.3435	C ₁₆ H ₂₂ O ₄		
DRE-C16173500	Phthalic acid, bis-iso-butyl ester(±)		250mg	
DRE-L16173500CY	Phthalic acid, bis-iso-butyl ester 10 µg/mL in Cyclohexane		10ml	
Phthalic Acid Diisobutyl Ester (3,4,5,6)-D4				
CAS 358730-88-8	MW 282.3681	C ₁₆ ² H ₄ H ₁₈ O ₄		
DRE-C16173510	Phthalic acid, bis-iso-butyl ester D4		10mg	
Phthalic Acid Di-iso-Octyl Ester				
CAS 27554-26-3	MW 390.5561	C ₂₄ H ₃₈ O ₄		
DRE-C16173650	Phthalic acid, bis-iso-octyl ester (technical)		250mg	
Phthalic Acid Diisopropyl Ester				
CAS 605-45-8	MW 250.2903	C ₁₄ H ₁₈ O ₄		
DRE-C16173700	Phthalic acid, bis-isopropylester		100mg	
Phthalic Acid Dimethyl Ester				
CAS 131-11-3	MW 194.184	C ₁₀ H ₁₀ O ₄		
DRE-C16174000	Phthalic acid, bis-methyl ester(±)		500mg	
Phthalic Acid Dimethyl Ester D4				
CAS 93951-89-4	MW 198.2086	C ₁₀ ² H ₄ H ₆ O ₄		
DRE-C16174010	Phthalic acid, bis-methyl ester D4(±)		10mg	
Phthalic Acid Dinonyl Ester				
CAS 84-76-4	MW 418.6093	C ₂₆ H ₄₂ O ₄		
DRE-C16174800	Phthalic acid, bis-nonyl ester(±)		250mg	
Phthalic Acid Di-n-nonyl Ester (3,4,5,6)-D4				
CAS 1202865-43-7	MW 422.6339	C ₂₆ ² H ₄ H ₃₈ O ₄		
DRE-XA16174810CY	Phthalic acid, bis-n-nonyl ester D4 100 µg/mL in Cyclohexane(±)		1ml	
Phthalic Acid Di-n-octyl Ester				
CAS 117-84-0	MW 390.5561	C ₂₄ H ₃₈ O ₄		
DRE-C16175000	Phthalic acid, bis-n-octyl ester(±)		250mg	










Food packaging contaminants

Product code	Description			
Phthalic Acid Di-n-octyl Ester (3,4,5,6)-D4				
CAS 93952-13-7 DRE-C16175010	MW 394.5808 Phthalic acid, bis-n-octyl ester-3,4,5,6-D4	$C_{24}^2H_{34}O_4$	100mg	
Phthalic Acid Dipentyl Ester				
CAS 131-18-0 DRE-C16175500	MW 306.3966 Phthalic acid, bis-n-pentyl ester(‡)	$C_{18}H_{26}O_4$	250mg	
Phthalic Acid Di-n-pentyl Ester (3,4,5,6)-D4				
CAS 358730-89-9 DRE-C16175510	MW 310.4213 Phthalic acid, bis-n-pentyl ester D4	$C_{18}^2H_{22}O_4$	10mg	
Phthalic Acid Diphenyl Ester				
CAS 84-62-8 DRE-C16176000	MW 318.3228 Phthalic acid, bis-phenyl ester	$C_{20}H_{14}O_4$	500mg	
Phthalic Acid Dipropyl Ester				
CAS 131-16-8 DRE-C16177000	MW 250.2903 Phthalic acid, bis-propyl ester(‡)	$C_{14}H_{18}O_4$	250mg	
Phthalic Acid Dipropyl Ester (3,4,5,6-D4)				
CAS 358731-29-0 DRE-C16177010	MW 254.315 Phthalic acid, bis-propyl ester D4	$C_{14}^2H_{14}O_4$	10mg	
Phthalic Acid 2-Ethylhexyl Hexyl Ester				
CAS 75673-16-4 DRE-C16178500	MW 362.503 Phthalic acid, hexyl-2-ethylhexyl ester	$C_{22}H_{34}O_4$	100mg	
Phthalic Acid Monobutyl Ester				
CAS 131-70-4 DRE-C16178700	MW 222.2372 Phthalic acid, mono-n-butyl ester	$C_{12}H_{14}O_4$	100mg	
Phthalic Acid Mono-2-ethylhexyl Ester				
CAS 4376-20-9 DRE-C16178900	MW 278.3435 Phthalic acid, mono-2-ethylhexyl ester(‡)	$C_{16}H_{22}O_4$	100mg	








Food packaging contaminants

Product code	Description			
Phthalic acid, mono-2-ethyl-5-oxohexyl ester				
CAS 40321-98-0 DRE-C16178960	MW 292.327 Phthalic acid, mono-2-ethyl-5-oxohexyl ester	C ₁₆ H ₂₀ O ₅	10mg	
Phthalic Acid Monomethyl Ester				
CAS 4376-18-5 DRE-C16179000	MW 180.1574 Phthalic acid, mono-methyl ester	C ₉ H ₈ O ₄	250mg	
Phthalic Acid Octyl Decyl Ester				
CAS 119-07-3 DRE-C16179070	MW 418.6093 Phthalic acid, octyldecyl ester	C ₂₆ H ₄₂ O ₄	100mg	
Phthalic Acid n-Pentyl Isopentyl Ester (3,4,5,6)-D4 (mixture of isomers)				
CAS DRE-C16179105	MW n/a Phthalic acid, n-pentyl-isopentyl ester D4 (mixture of isomers)		10mg	No Structure
Phthalic Acid n-Pentyl Isopentyl Ester (mixture of isomers)				
CAS 84777-06-0 DRE-C16179100	MW n/a Phthalic acid, n-pentyl-isopentyl ester (mixture of isomers)(‡)		100mg	No Structure
Phthalic Acid n-Pentyl Isopentyl Ester				
CAS 776297-69-9 DRE-C16179120	MW 306.3966 Phthalic acid, n-pentyl-isopentyl ester	C ₁₈ H ₂₆ O ₄	100mg	
Phthalic acid, mixed decyl-hexyl-octyl diester				
CAS 68648-93-1 DRE-C16178600	MW n/a Phthalic acid, mixed decyl-hexyl-octyl diester		250mg	No Structure
Phthalic Anhydride				
CAS 85-44-9 DRE-C16183000	MW 148.1156 Phthalic anhydride	C ₈ H ₄ O ₃	250mg	
3,3,4,4-Tetradeuterioadipic Acid Bis(2-ethylhexyl) Ester				
CAS n/a DRE-C10046010 DRE-XA10046010AC	MW 374.5911 Adipic acid-bis-2-ethylhexyl ester D4 Adipic acid-bis-2-ethylhexyl ester D4 100 µg/mL in Acetone	C ₂₂ ² H ₄₄ H ₃₈ O ₄	10mg 1.1ml	

Food packaging contaminants

Product code	Description			
Tri-n-butyl phosphate				
CAS 126-73-8	MW 266.3141	C ₁₂ H ₂₇ O ₄ P		
DRE-C17668000	Tributylphosphate		250mg	
DRE-L17668000CY	Tributylphosphate 10 µg/mL in Cyclohexane		10ml	
Tri-m-cresyl Phosphate				
CAS 563-04-2	MW 368.3628	C ₂₁ H ₂₁ O ₄ P		
DRE-C17800820	Tri-m-cresyl phosphate		100mg	
Tri-o-cresyl Phosphate				
CAS 78-30-8	MW 368.3628	C ₂₁ H ₂₁ O ₄ P		
DRE-C17800810	Tri-o-cresyl phosphate(±)		250mg	
Tri-p-cresyl Phosphate				
CAS 78-32-0	MW 368.3628	C ₂₁ H ₂₁ O ₄ P		
DRE-C17800830	Tri-p-cresyl phosphate(±)		100mg	
Tricresylphosphate				
CAS 1330-78-5	MW 1105.0884	3C ₂₁ H ₂₁ O ₄ P		
DRE-C17800800	Tricresylphosphate		250mg	
Triethylphosphate				
CAS 78-40-0	MW 182.1547	C ₆ H ₁₅ O ₄ P		
DRE-C17835500	Triethylphosphate		250mg	
2,2,4-Trimethyl-1,3-pentanediol Diisobutyrate (TXIB)				
CAS 6846-50-0	MW 286.407	C ₁₆ H ₃₀ O ₄		
DRE-C17883150	2,2,4-Trimethyl-1,3-pentanediol diisobutyrate (TXIB)		1ml	
Triphenylphosphate				
CAS 115-86-6	MW 326.2831	C ₁₈ H ₁₅ O ₄ P		
DRE-C17893000	Triphenylphosphate		250mg	
DRE-L17893000EA	Triphenylphosphate 10 µg/mL in Ethylacetate		10ml	
DRE-XA17893000MB	Triphenylphosphate, 500 µg/ml in Methyl-tert. butyl ether		1ml	
Tripropyl Phosphate				
CAS 513-08-6	MW 224.2344	C ₉ H ₂₁ O ₄ P		
DRE-C17893830	Tri-n-propylphosphate		100mg	


Food packaging contaminants

Product code	Description			
Tris-(2-butoxyethyl)phosphate				
CAS 78-51-3 DRE-C17894100	MW 398.4718 Tris-(2-butoxyethyl)phosphate	C ₁₈ H ₃₉ O ₇ P	250mg	
Tris(2-chloro-1-(chloromethyl)ethyl) Phosphate				
CAS 13674-87-8 DRE-C17894320	Tris(2-chloro-1-(chloromethyl)ethyl) phosphate(±)	C ₉ H ₁₅ Cl ₅ O ₄ P	250mg	
Tris-(2-chloroethyl)phosphate				
CAS 115-96-8 DRE-C17894300 DRE-L17894300CY	Tris-(2-chloroethyl)phosphate Tris-(2-chloroethyl)phosphate 10 µg/mL in Cyclohexane	C ₆ H ₁₂ Cl ₃ O ₄ P	250mg 10ml	
Tris-(2-chloropropyl)phosphate				
CAS 13674-84-5 DRE-C17894330	Tris-(2-chloropropyl)phosphate(±)	C ₉ H ₁₈ Cl ₃ O ₄ P	250mg	
Tris(2,3-dibromopropyl)phosphate				
CAS 126-72-7 DRE-C17894340	Tris(2,3-dibromopropyl)phosphate	C ₉ H ₁₅ Br ₃ O ₄ P	100mg	
Tris(2,4-di-tert-butylphenyl)phosphite				
CAS 31570-04-4 DRE-C17894350	Tris(2,4-di-tert-butylphenyl)phosphite	C ₄₂ H ₆₃ O ₃ P	250mg	
Tris-(2-ethylhexyl)phosphate				
CAS 78-42-2 DRE-C17894400	Tris-(2-ethylhexyl)phosphate	C ₂₄ H ₅₁ O ₄ P	250mg	
Diisononyl Phthalate and Diisodecyl Phthalate Mix				
DRE-GS09000391IO	Diisononyl Phthalate and Diisodecyl Phthalate Mix 1000 µg/mL in Isooctane(±)			5x1ml
	diisononyl phthalate (DINP : mix of isomers) diisodecyl phthalate (mix of isomers)			
Phthalate Esters Mix 1				
DRE-XA08060100IO	Phthalate Esters Mix 1 200 µg/mL in Isooctane			1ml
DRE-YA08060100IO	Phthalate Esters Mix 1 2000 µg/mL in Isooctane			1ml
	Phthalic acid, benzylbutyl ester Phthalic acid, bis-2-ethylhexyl ester Phthalic acid, bis-butyl ester Phthalic acid, bis-ethyl ester Phthalic acid, bis-methyl ester Phthalic acid, bis-n-octyl ester			






Food packaging contaminants

Product code	Description	
Phthalate Esters - Analytes Mix 3		
DRE-YA08060300HE	Phthalate Esters - Analytes Mix 3 1000 µg/mL in n-Hexane(‡)	1ml
Benzoic acid-benzyl ester	Phthalic acid, benzylbutyl ester	Phthalic acid, bis-2-ethoxyethyl ester
Phthalic acid, bis-2-n-butoxyethyl ester	Phthalic acid, bis-4-methyl-2-pentyl ester	Phthalic acid, bis-butyl ester
Phthalic acid, bis-ethyl ester	Phthalic acid, bis-hexyl ester	Phthalic acid, bis-iso-butyl ester
Phthalic acid, bis-methylglycol ester	Phthalic acid, bis-n-octyl ester	Phthalic acid, bis-nonyl ester
Phthalic acid, hexyl-2-ethylhexyl ester		Phthalic acid, bis-n-pentyl ester
Phthalate Esters Mix 10		
DRE-YA08061000HE	Phthalate Esters Mix 10 each 1000 µg/mL in Hexane	1ml
Phthalic acid, benzylbutyl ester	Phthalic acid, bis-2-ethoxyethyl ester	Phthalic acid, bis-2-ethylhexyl ester
Phthalic acid, bis-4-methyl-2-pentyl ester	Phthalic acid, bis-butyl ester	Phthalic acid, bis-cyclohexyl ester
Phthalic acid, bis-hexyl ester	Phthalic acid, bis-iso-butyl ester	Phthalic acid, bis-methyl ester
Phthalic acid, bis-n-octyl ester	Phthalic acid, bis-nonyl ester	Phthalic acid, bis-n-pentyl ester
		Phthalic acid, bis-phenyl ester
Phthalate Mix, 10-3856		
DRE-GA09000392ME	Phthalate Mix, 10-3856 in Methanol(‡)	5ml
	diisononyl phthalate (DINP : NOT mix of isomers) [10 mg/L]	diisodecyl phthalate (mix of isomers) [10 mg/L]
	di-n-butyl phthalate [5 mg/L]	butyl benzyl phthalate [5 mg/L]
	bis(2-ethylhexyl)phthalate [5 mg/L]	di-n-octyl phthalate [5 mg/L]
	di-n-hexyl phthalate [5 mg/L]	diisobutylphthalate [5 mg/L]
	diamyl phthalate [5 mg/L]	dicyclohexyl phthalate [5 mg/L]
Phthalate and Adipate Esters Mix 1		
DRE-XA05060100ME	Phthalate and Adipate Esters Mix 1 100 µg/mL in Methanol	1ml
	Adipic acid-bis-2-ethylhexyl ester	
	Phthalic acid, benzylbutyl ester	
	Phthalic acid, bis-2-ethylhexyl ester	
	Phthalic acid, bis-butyl ester	
	Phthalic acid, bis-ethyl ester	
	Phthalic acid, bis-methyl ester	
	Phthalic acid, bis-n-octyl ester	

STANDARDS FOR REGULATORY METHODS



Standards for regulatory methods

Product code	Description			
ASTM D3606 2-Butanol				
CAS 78-92-2	MW 74.1216	C ₄ H ₁₀ O		
DRE-GA09010059	ASTM D3606 2-Butanol IS(±)		2ml	
DRE-GS09010059	ASTM D3606 2-Butanol IS(±)		5x2ml	
2-Bromobutanoic Acid (2-Bromobutyric Acid)				
CAS 80-58-0	MW 167.0012	C ₄ H ₇ BrO ₂		
DRE-GA09010313MB	EPA Method 552.3 2-Bromobutanoic Acid 10000 µg/mL in Methyl tert-Butyl Ether(±)		1ml	
Chloroprene (2-Chloro-1,3-butadiene; β-Chloroprene)				
CAS 126-99-8	MW 88.5355	C ₄ H ₅ Cl		
DRE-GA09010381ME	Chloroprene Solution, 5000 µg/mL in p/t Methanol(±)		1ml	
Decafluorotriphenylphosphine (DFTPP)				
CAS 5074-71-5	MW 442.1901	C ₁₈ H ₅ F ₁₀ P		
DRE-GA09010393AC	Decafluorotriphenylphosphine (DFTPP) SV Tuning Solution, 1000 µg/mL in Acetone(±)		1ml	
DRE-GS09010393AC	Decafluorotriphenylphosphine (DFTPP) SV Tuning Solution, 1000 µg/mL in Acetone(±)		5x1ml	
DRE-GA09010394AC	Decafluorotriphenylphosphine (DFTPP) SV Tuning Solution, 2500 µg/mL in Acetone(±)		1ml	
DRE-GS09010394AC	Decafluorotriphenylphosphine (DFTPP) SV Tuning Solution, 2500 µg/mL in Acetone(±)		5x1ml	
1,2-Dimethoxyethane				
CAS 110-71-4	MW 90.121	C ₄ H ₁₀ O ₂		
DRE-GA09010060	ASTM D4815/D5599 1,2-Dimethoxyethane IS(±)		2ml	
DRE-GS09010060	ASTM D4815/D5599 1,2-Dimethoxyethane IS(±)		5x2ml	
525.2 GC/MS Performance Check Solution				
DRE-GS09000334AC	525.2 GC/MS Performance Check Mix 1000 µg/mL in Acetone(±)			5x1ml
	decafluorotriphenylphosphine (DFTPP) p,p'-DDT endrin			
525.3 Internal Standard Solution				
DRE-GS09000347AC	525.3 Internal Standard Mix 500 µg/mL in Acetone(±)			5x1ml
	acenaphthene-d10 chrysene-d12 phenanthrene-d10			
8270 BN Solution 14-3				
DRE-GS09000207DI	8270 BN Mix 14-3 2000 µg/mL in Dichloromethane(±)			5x1ml
	2-chloronaphthalene	1,2-dichlorobenzene		
	1,3-dichlorobenzene	1,4-dichlorobenzene		
	hexachlorobenzene	hexachlorobutadiene		
	hexachlorocyclopentadiene	hexachloroethane		
	1,2,4-trichlorobenzene	2,4-dinitrotoluene		
	2,6-dinitrotoluene	isophorone		
	nitrobenzene	azobenzene		

Standards for regulatory methods

Product code	Description		
8270 Spiking Mix, 129			
DRE-GS09000239AC	8270 Spiking Mix in Acetone(‡)		
			5x25ml
benzo[k]fluoranthene [50 mg/L]	1-methylnaphthalene [50 mg/L]	2-methylnaphthalene [50 mg/L]	acenaphthene [50 mg/L]
acenaphthylene [50 mg/L]	anthracene [50 mg/L]	fluorene [50 mg/L]	naphthalene [50 mg/L]
phenanthrene [50 mg/L]	benzo[a]anthracene [50 mg/L]	benzo[a]pyrene [50 mg/L]	chrysene [50 mg/L]
fluoranthene [50 mg/L]	indeno[1,2,3-cd]pyrene [50 mg/L]	pyrene [50 mg/L]	benzo[b]fluoranthene [50 mg/L]
benzo[ghi]perylene [50 mg/L]	dibenz[a,h]anthracene [50 mg/L]	phenacetin [50 mg/L]	7,12-dimethylbenz[a]anthrac. [50 mg/L]
3-methylcholanthrene [50 mg/L]	1-naphthylamine [50 mg/L]	5-nitro-o-toluidine [50 mg/L]	4-aminobiphenyl [50 mg/L]
2-acetylaminofluorene [50 mg/L]	p-(dimethylamino)azobenzene [50 mg/L]	bis(2-chloroethoxy)methane [50 mg/L]	bis(2-chloroethyl)ether [50 mg/L]
bis(1-chloroprop-2-yl) ether [50 mg/L]	4-bromophenyl phenyl ether [50 mg/L]	4-chlorodiphenyl ether [50 mg/L]	bis(2-ethylhexyl)phthalate [50 mg/L]
butyl benzyl phthalate [50 mg/L]	diethyl phthalate [50 mg/L]	dimethyl phthalate [50 mg/L]	di-n-butyl phthalate [50 mg/L]
di-n-octyl phthalate [50 mg/L]	2-chloronaphthalene [50 mg/L]	1,2-dichlorobenzene [50 mg/L]	1,3-dichlorobenzene [50 mg/L]
1,4-dichlorobenzene [50 mg/L]	hexachlorobenzene [50 mg/L]	hexachlorobutadiene [50 mg/L]	hexachlorocyclopentadiene [50 mg/L]
hexachloroethane [50 mg/L]	1,2,4-trichlorobenzene [50 mg/L]	2,6-dichlorophenol [50 mg/L]	hexachloropropene [50 mg/L]
1,2,4,5-tetrachlorobenzene [50 mg/L]	pentachlorobenzene [50 mg/L]	pentachloronitrobenzene [50 mg/L]	aniline [50 mg/L]
4-chloroaniline [50 mg/L]	2-nitroaniline [50 mg/L]	3-nitroaniline [50 mg/L]	4-nitroaniline [50 mg/L]
pyridine [50 mg/L]	2-chlorophenol [50 mg/L]	2,4-dimethylphenol [50 mg/L]	pentachlorophenol [50 mg/L]
4-nitrophenol [50 mg/L]	2,4-dichlorophenol [50 mg/L]	4-chloro-3-methylphenol [50 mg/L]	2-methyl-4,6-dinitrophenol [50 mg/L]
2-nitrophenol [50 mg/L]	2,4-dinitrophenol [50 mg/L]	2,4,6-trichlorophenol [50 mg/L]	phenol [50 mg/L]
2-methylphenol [50 mg/L]	4-methylphenol [50 mg/L]	2,4,5-trichlorophenol [50 mg/L]	benzyl alcohol [50 mg/L]
dibenzofuran [50 mg/L]	acetophenone [50 mg/L]	2,4-dinitrotoluene [50 mg/L]	2,6-dinitrotoluene [50 mg/L]
isophorone [50 mg/L]	nitrobenzene [50 mg/L]	diphenylamine [50 mg/L]	N-nitrosodi-n-propylamine [50 mg/L]
N-nitrosomethylethylamine [50 mg/L]	N-nitrosopyrrolidine [50 mg/L]	n-nitrosodiethylamine [50 mg/L]	N-nitrosomorpholine [50 mg/L]
n-nitrosodi-n-butylamine [50 mg/L]	pentachloroethane [50 mg/L]	dinoseb [50 mg/L]	isodrin [50 mg/L]
caprolactam [50 mg/L]	kepone [50 mg/L]	diphenyl ether [50 mg/L]	1,4-dinitrobenzene [50 mg/L]
2,6-dimethylphenol [50 mg/L]	6-methylchrysene [50 mg/L]	methapyrilene hydrochloride [50 mg/L]	dibenz(a,h)acridine [50 mg/L]
aramite [50 mg/L]	atrazine [50 mg/L]	benzidine [100 mg/L]	biphenyl [50 mg/L]
chlorobenzilate [50 mg/L]	di-allate (mixture of isomers) [100 mg/L]	dibenz(a,j)acridine [50 mg/L]	n-nitrosodiphenylamine [50 mg/L]
propyzamide (pronamide) [50 mg/L]	1-chloronaphthalene [50 mg/L]	1,3-dinitrobenzene [50 mg/L]	2-naphthylamine [50 mg/L]
2-picoline [50 mg/L]	2,3,4,6-Tetrachlorophenol [50 mg/L]	carbazole [50 mg/L]	benzoic acid [50 mg/L]
N-nitrosodimethylamine [50 mg/L]	azobenzene [50 mg/L]	N-nitrosopiperidine [50 mg/L]	quinoline [50 mg/L]
octachlorostyrene [50 mg/L]	p-phenylenediamine [500 mg/L]	1,4-naphthoquinone [100 mg/L]	4,4'-methylenedianiline [50 mg/L]
methyl methanesulfonate [50 mg/L]	isosafole [100 mg/L]	4-nitroquinoline-1-oxide [50 mg/L]	3,3'-dimethylbenzidine [50 mg/L]
1,3,5-trinitrobenzene [50 mg/L]	ethyl methanesulfonate [50 mg/L]	safrole [50 mg/L]	o-toluidine [50 mg/L]
3,3'-dichlorobenzidine [50 mg/L]			

8270 Spiking Mix, 130			
DRE-XS09000028AC	8270 Spiking Mix in Acetone		
			10x25ml
a,a-dimethylphenethylamine [50 mg/L]	acenaphthene [50 mg/L]	acenaphthylene [50 mg/L]	acetophenone [50 mg/L]
aniline [50 mg/L]	anthracene [50 mg/L]	aramite [50 mg/L]	atrazine [50 mg/L]
azobenzene [50 mg/L]	benzidine [100 mg/L]	benzo[a]anthracene [50 mg/L]	benzo[a]pyrene [50 mg/L]
benzo[b]fluoranthene [50 mg/L]	benzo[ghi]perylene [50 mg/L]	benzo[k]fluoranthene [50 mg/L]	benzoic acid [50 mg/L]
benzyl alcohol [50 mg/L]	butyl benzyl phthalate [50 mg/L]	biphenyl [50 mg/L]	bis(2-chloroethoxy)methane [50 mg/L]
bis(2-chloroethyl)ether [50 mg/L]	bis(1-chloroprop-2-yl) ether [50 mg/L]	diocetyl phthalate [50 mg/L]	caprolactam [50 mg/L]
carbazole [50 mg/L]	chlorobenzilate [50 mg/L]	chrysene [50 mg/L]	di-n-butyl phthalate [50 mg/L]
di-n-octyl phthalate [50 mg/L]	di-allate (mixture of isomers) [100 mg/L]	dibenz(a,h)acridine [50 mg/L]	dibenz[a,h]anthracene [50 mg/L]
dibenz(a,j)acridine [50 mg/L]	dibenzofuran [50 mg/L]	diethyl phthalate [50 mg/L]	dimethyl phthalate [50 mg/L]
dinoseb [50 mg/L]	diphenylamine [50 mg/L]	ethyl methanesulfonate [50 mg/L]	fluoranthene [50 mg/L]
fluorene [50 mg/L]	hexachlorobenzene [50 mg/L]	hexachlorobutadiene [50 mg/L]	hexachlorocyclopentadiene [50 mg/L]
hexachloroethane [50 mg/L]	hexachloropropene [50 mg/L]	indeno[1,2,3-cd]pyrene [50 mg/L]	isodrin [50 mg/L]
isophorone [50 mg/L]	isosafole [100 mg/L]	kepone [50 mg/L]	methapyrilene hydrochloride [50 mg/L]
methyl methanesulfonate [50 mg/L]	4-nitroquinoline-1-oxide [50 mg/L]	N-nitrosomethylethylamine [50 mg/L]	n-nitrosodi-n-butylamine [50 mg/L]
N-nitrosodi-n-propylamine [50 mg/L]	n-nitrosodiethylamine [50 mg/L]	N-nitrosodimethylamine [50 mg/L]	n-nitrosodiphenylamine [50 mg/L]
N-nitrosomorpholine [50 mg/L]	N-nitrosopiperidine [50 mg/L]	N-nitrosopyrrolidine [50 mg/L]	naphthalene [50 mg/L]
nitrobenzene [50 mg/L]	o-toluidine [50 mg/L]	octachlorostyrene [50 mg/L]	p-(dimethylamino)azobenzene [50 mg/L]
pentachlorobenzene [50 mg/L]	pentachloroethane [50 mg/L]	pentachloronitrobenzene [50 mg/L]	pentachlorophenol [50 mg/L]
phenacetin [50 mg/L]	phenanthrene [50 mg/L]	phenol [50 mg/L]	diphenyl ether [50 mg/L]
propyzamide (pronamide) [50 mg/L]	pyrene [50 mg/L]	pyridine [50 mg/L]	quinoline [50 mg/L]
safrole [50 mg/L]	1-chloronaphthalene [50 mg/L]	1-methylnaphthalene [50 mg/L]	1-naphthylamine [50 mg/L]
1,2-dichlorobenzene [50 mg/L]	1,2,4-trichlorobenzene [50 mg/L]	1,2,4,5-tetrachlorobenzene [50 mg/L]	1,3-dichlorobenzene [50 mg/L]
1,3-dinitrobenzene [50 mg/L]	1,3,5-trinitrobenzene [50 mg/L]	1,4-dichlorobenzene [50 mg/L]	1,4-dinitrobenzene [50 mg/L]
1,4-naphthoquinone [100 mg/L]	p-phenylenediamine [500 mg/L]	2-acetylaminofluorene [50 mg/L]	2-chloronaphthalene [50 mg/L]
2-chlorophenol [50 mg/L]	2-methyl-4,6-dinitrophenol [50 mg/L]	2-methylnaphthalene [50 mg/L]	2-methylphenol [50 mg/L]
2-naphthylamine [50 mg/L]	2-nitroaniline [50 mg/L]	2-nitrophenol [50 mg/L]	2-picoline [50 mg/L]
2,3,4,6-tetrachlorophenol [50 mg/L]	2,4-dichlorophenol [50 mg/L]	2,4-dimethylphenol [50 mg/L]	2,4-dinitrophenol [50 mg/L]
2,4-dinitrotoluene [50 mg/L]	2,4,5-trichlorophenol [50 mg/L]	2,4,6-trichlorophenol [50 mg/L]	2,6-dichlorophenol [50 mg/L]
2,6-dimethylphenol [50 mg/L]	2,6-dinitrotoluene [50 mg/L]	3-methylcholanthrene [50 mg/L]	3-nitroaniline [50 mg/L]
3,3'-dichlorobenzidine [50 mg/L]	o-toluidine [50 mg/L]	4-aminobiphenyl [50 mg/L]	4-bromophenyl phenyl ether [50 mg/L]
4-chloro-3-methylphenol [50 mg/L]	4-chloroaniline [50 mg/L]	4-chlorodiphenyl ether [50 mg/L]	4-methylphenol [50 mg/L]
4-nitroaniline [50 mg/L]	4-nitrophenol [50 mg/L]	4,4'-methylenedianiline [50 mg/L]	5-nitro-o-toluidine [50 mg/L]
6-methylchrysene [50 mg/L]	7,12-dimethylbenz[a]anthrac. [50 mg/L]		

Standards for regulatory methods

Product code	Description	
ASTM D1319 Aromatics and Olefins by FIA Mixture		
DRE-GA09000379IO	ASTM D1319 Aromatics and Olefins by FIA Mixture, 20% LV Toluene and 5% LV 1-Hexene in Isooctane(‡) toluene [20 wt%] 1-hexene [5 wt%] isooctane [75 wt%]	1ea
ASTM D3606 Benzene in Gasoline Std. #1 w/ 10% EtOH & IS		
DRE-GA09000101IO	ASTM D3606 Benzene in Gasoline Std #1 w/ 10% EtOH & IS in Isooctane(‡) benzene [50 mL/L] toluene [200 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #1 w/ IS		
DRE-GA09000093IO	ASTM D3606 Benzene in Gasoline Std. #1 w/ IS in Isooctane(‡) benzene [50 mL/L] toluene [200 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #2 w/ 10% EtOH & IS		
DRE-GA09000102IO	ASTM D3606 Benzene in Gasoline Std #2 w/ 10% EtOH & IS in Isooctane(‡) benzene [25 mL/L] toluene [150 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #2 w/ IS		
DRE-GA09000094IO	ASTM D3606 Benzene in Gasoline Std. #2 w/ IS in Isooctane(‡) benzene [25 mL/L] toluene [150 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #3 w/ 10% EtOH & IS		
DRE-GA09000103IO	ASTM D3606 Benzene in Gasoline Std #3 w/ 10% EtOH & IS in Isooctane(‡) benzene [12,5 mL/L] toluene [100 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #3 w/ IS		
DRE-GA09000095IO	ASTM D3606 Benzene in Gasoline Std. #3 w/ IS in Isooctane(‡) benzene [10 mL/L] toluene [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #4 w/ 10% EtOH & IS		
DRE-GA09000104IO	ASTM D3606 Benzene in Gasoline Std #4 w/ 10% EtOH & IS in Isooctane(‡) benzene [6,5 mL/L] toluene [50 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #4 w/ IS		
DRE-GA09000096IO	ASTM D3606 Benzene in Gasoline Std. #4 w/ IS in Isooctane(‡) benzene [6 mL/L] toluene [50 mL/L] 2-butanol (IS) [40 mL/L]	2ml

Standards for regulatory methods

Product code	Description	
ASTM D3606 Benzene in Gasoline Std. #5 w/ 10% EtOH & IS		
DRE-GA09000105IO	ASTM D3606 Benzene in Gasoline Std #5 w/ 10% EtOH & IS in Isooctane(‡) benzene [3,3 mL/L] toluene [25 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #5 w/ IS		
DRE-GA09000097IO	ASTM D3606 Benzene in Gasoline Std. #5 w/ IS in Isooctane(‡) benzene [3 mL/L] toluene [25 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #6 w/ 10% EtOH & IS		
DRE-GA09000106IO	ASTM D3606 Benzene in Gasoline Std #6 w/ 10% EtOH & IS in Isooctane(‡) benzene [1,2 mL/L] toluene [10 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #6 w/ IS		
DRE-GA09000098IO	ASTM D3606 Benzene in Gasoline Std. #6 w/ IS in Isooctane(‡) benzene [1 mL/L] toluene [10 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #7 w/ 10% EtOH & IS		
DRE-GA09000107IO	ASTM D3606 Benzene in Gasoline Std #7 w/ 10% EtOH & IS in Isooctane(‡) benzene [0,6 mL/L] toluene [5 mL/L] ethanol [100 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. #7 w/ IS		
DRE-GA09000099IO	ASTM D3606 Benzene in Gasoline Std. #7 w/ IS in Isooctane(‡) benzene [0,6 mL/L] toluene [5 mL/L] 2-butanol (IS) [40 mL/L]	2ml
ASTM D3606 Benzene in Gasoline Std. Kit w/ 10% EtOH & IS		
DRE-GK09000108IO	ASTM D3606 Benzene in Gasoline Kit w/ 10% EtOH & IS in Isooctane(‡)	1ea
DRE-GA09000101IO	ASTM D3606 Benzene in Gasoline Std. #1 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000102IO	ASTM D3606 Benzene in Gasoline Std. #2 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000103IO	ASTM D3606 Benzene in Gasoline Std. #3 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000104IO	ASTM D3606 Benzene in Gasoline Std. #4 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000105IO	ASTM D3606 Benzene in Gasoline Std. #5 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000106IO	ASTM D3606 Benzene in Gasoline Std. #6 w/ 10% EtOH & IS (‡)	1x2ml
DRE-GA09000107IO	ASTM D3606 Benzene in Gasoline Std. #7 w/ 10% EtOH & IS (‡)	1x2ml
ASTM D3606 Benzene in Gasoline Std. Kit w/ IS		
DRE-GK09000100IO	ASTM D3606 Benzene in Gasoline Kit with IS in Isooctane(‡)	1ea
DRE-GA09000093IO	ASTM D3606 Benzene in Gasoline Std. #1 w/ IS (‡)	1x2ml
DRE-GA09000094IO	ASTM D3606 Benzene in Gasoline Std. #2 w/ IS (‡)	1x2ml
DRE-GA09000095IO	ASTM D3606 Benzene in Gasoline Std. #3 w/ IS (‡)	1x2ml
DRE-GA09000096IO	ASTM D3606 Benzene in Gasoline Std. #4 w/ IS (‡)	1x2ml
DRE-GA09000097IO	ASTM D3606 Benzene in Gasoline Std. #5 w/ IS (‡)	1x2ml
DRE-GA09000098IO	ASTM D3606 Benzene in Gasoline Std. #6 w/ IS (‡)	1x2ml
DRE-GA09000099IO	ASTM D3606 Benzene in Gasoline Std. #7 w/ IS (‡)	1x2ml

Standards for regulatory methods

Product code	Description	
ASTM D3606 Check Standard A w/ 4% 2-Butanol		
DRE-GA09000109IO	ASTM Method D3606 Check Standard A w/ 4% 2-Butanol in Isooctane(‡)	10x2ml
	benzene [1 vol%] toluene [5 vol%] ethanol [10 vol%] 2-butanol (IS) [4 vol%]	
ASTM D3606 Check Standard B w/ 4% 2-Butanol		
DRE-GA09000110IO	ASTM Method D3606 Check Standard with 4% 2-Butanol in Isooctane(‡)	10x2ml
	benzene [1 vol%] toluene [5 vol%] ethanol [0 vol%] 2-butanol (IS) [4 vol%]	
ASTM D4815 Oxygenates in Gasoline Cal. Std. #1 w/ IS		
DRE-GA09000111OG	ASTM 4815 Oxygenates in Gasoline Calibration Std. #1 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [1,25 wt%] methyl t-butyl ether [16 wt%] 1-butanol [0,25 wt%] 1-propanol [0,25 wt%] tert-amyl alcohol [0,25 wt%] oxygenate free gasoline, unleaded, premium [70 wt%]	ethanol [1,25 wt%] tert-butyl ethyl ether (ETBE) [5 wt%] isopropyl alcohol [0,5 wt%] isobutyl alcohol [0,25 wt%] 1,2-dimethoxyethane (IS) [5 wt%]
ASTM D4815 Oxygenates in Gasoline Cal. Std. #2 w/ IS		
DRE-GA09000112OG	ASTM 4815 Oxygenates in Gasoline Calibration Std. #2 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [5 wt%] 2-methyl-2-propanol [9 wt%] 1-butanol [0,5 wt%] isopropyl alcohol [0,25 wt%] isobutyl alcohol [0,5 wt%] tert-amyl alcohol [1,5 wt%] oxygenate free gasoline, unleaded, premium [69,75 wt%]	ethanol [0,25 wt%] methyl t-butyl ether [5 wt%] isopropyl ether [0,25 wt%] 1-propanol [0,5 wt%] 2-butanol [2,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%]
ASTM D4815 Oxygenates in Gasoline Cal. Std. #3 w/ IS		
DRE-GA09000113OG	ASTM 4815 Oxygenates in Gasoline Calibration Std. #3 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [9 wt%] 2-methyl-2-propanol [5 wt%] tert-butyl ethyl ether (ETBE) [9 wt%] tert-amyl methyl ether (TAME) [1,25 wt%] tert-amyl alcohol [0,75 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [70 wt%]	
ASTM D4815 Oxygenates in Gasoline Cal. Std. #4 w/ IS		
DRE-GA09000114OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #4 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [0,5 wt%] ethanol [5 wt%] methyl t-butyl ether [1 wt%] tert-butyl ethyl ether (ETBE) [15 wt%] isopropyl ether [3,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [70 wt%]	
ASTM D4815 Oxygenates in Gasoline Cal. Std. #5 w/ IS		
DRE-GA09000115OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #5 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	ethanol [10 wt%] methyl t-butyl ether [10 wt%] tert-amyl methyl ether (TAME) [0,5 wt%] 2-butanol [0,25 wt%] oxygenate free gasoline, unleaded, premium [71 wt%]	2-methyl-2-propanol [0,5 wt%] tert-butyl ethyl ether (ETBE) [1,25 wt%] isopropyl ether [1,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%]

Standards for regulatory methods

Product code	Description	
ASTM D4815 Oxygenates in Gasoline Cal. Std. #6 w/ IS		
DRE-GA09000116OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #6 w/ IS in Oxygenate Free Gasoline(‡) methanol [15 wt%] methyl t-butyl ether [7 wt%] isopropyl alcohol [2,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [70,5 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. #7 w/ IS		
DRE-GA09000117OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #7 w/ IS in Oxygenate Free Gasoline(‡) ethanol [15 wt%] tert-amyl methyl ether (TAME) [10 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [70 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. #8 w/ IS		
DRE-GA09000118OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #8 w/ IS in Oxygenate Free Gasoline(‡) 2-methyl-2-propanol [15 wt%] 1-butanol [1,25 wt%] tert-amyl methyl ether (TAME) [5 wt%] isopropyl alcohol [1,5 wt%] 1-propanol [2,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [69,75 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. #9 w/ IS		
DRE-GA09000119OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #9 w/ IS in Oxygenate Free Gasoline(‡) tert-amyl methyl ether (TAME) [15 wt%] 1-butanol [5 wt%] isobutyl alcohol [2,5 wt%] tert-amyl alcohol [0,5 wt%] oxygenate free gasoline, unleaded, premium [69,5 wt%] isopropyl ether [0,5 wt%] 1-propanol [1,5 wt%] 2-butanol [0,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. #10 w/ IS		
DRE-GA09000120OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #10 w/ IS in Oxygenate Free Gasoline (‡) tert-butyl ethyl ether (ETBE) [12 wt%] isopropyl alcohol [5 wt%] 1-propanol [5 wt%] isobutyl alcohol [1,5 wt%] 2-butanol [1,5 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate free gasoline, unleaded, premium [70 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. #11 w/ IS		
DRE-GA09000121OG	ASTM Method D4815 Oxygenates in Gasoline Calibration Mix #11 w/ IS in Oxygenate Free Gasoline (‡) 2-methyl-2-propanol [1.25 wt%] 1-butanol [2.5 wt%] 2-butanol [5 wt%] 1,2-dimethoxyethane (IS) [5 wt%] isopropyl ether [5 wt%] isobutyl alcohol [5 wt%] tert-amyl alcohol [2.5 wt%] oxygenate free gasoline, unleaded, premium [73.75 wt%]	2ml
ASTM D4815 Oxygenates in Gasoline Cal. Std. Kit w/ IS		
DRE-GK09000122OG	ASTM D4815 Oxygenates in Gasoline Calibration Kit w/ IS in Oxygenate Free Gasoline(‡)	1ea
DRE-GA09000111OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #1 w/ IS (‡)	1x2ml
DRE-GA09000112OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #2 w/ IS (‡)	1x2ml
DRE-GA09000113OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #3 w/ IS (‡)	1x2ml
DRE-GA09000114OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #4 w/ IS (‡)	1x2ml
DRE-GA09000115OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #5 w/ IS (‡)	1x2ml
DRE-GA09000116OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #6 w/ IS (‡)	1x2ml
DRE-GA09000117OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #7 w/ IS (‡)	1x2ml

Standards for regulatory methods

Product code	Description	
DRE-GA09000118OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #8 w/ IS (‡)	1x2ml
DRE-GA09000119OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #9 w/ IS (‡)	1x2ml
DRE-GA09000120OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #10 w/ IS (‡)	1x2ml
DRE-GA09000121OG	ASTM D4815 Oxygenates in Gasoline Cal. Std. #11 w/ IS (‡)	1x2ml
ASTM D4815 Quantitative Peak Mix		
DRE-GS09000186	ASTM D4815 Quantitative Peak Mix(‡)	5x1ml
tert-amyl methyl ether (TAME) [7,3 wt%] 2-methyl-2-propanol [7,3 wt%] 1,2-dimethoxyethane [6 wt%] methylcyclopentane [4 wt%]	benzene [5 wt%] tert-butyl ethyl ether (ETBE) [4 wt%] ethanol [7,3 wt%] isobutyl alcohol [7,3 wt%]	1-butanol [7,3 wt%] methyl t-butyl ether [4 wt%] methanol [7,3 wt%] 1-propanol [7,3 wt%]
		2-butanol [7,3 wt%] isopropyl ether [4 wt%] tert-amyl alcohol [7,3 wt%] isopropyl alcohol [7,3 wt%]
ASTM D4815 Valve Timing Mix		
DRE-GA09000135	ASTM D4815 Valve Timing Mix(‡)	1ml
DRE-GS09000135	ASTM D4815 Valve Timing Mix(‡)	5x1ml
	tert-butyl ethyl ether (ETBE) [10 wt%] isopropyl ether [10 wt%] n-hexane (C6) [60 wt%] methylcyclopentane [10 wt%] methyl t-butyl ether [10 wt%]	
ASTM D4815/D5599 10% Ethanol with 5% 1,2-Dimethoxyethane Internal Standard (Added Above)		
DRE-SP09000237OG	ASTM D4815/D5599 10% Ethanol with 5% 1,2-Dimethoxyethane Internal Standard Mix(‡)	60x2ml
	ethanol [10 wt%] 1,2-dimethoxyethane (IS) [5 wt%] oxygenate-free RFA gasoline [90 wt%]	
ASTM D5501 96% Ethanol QC Check		
DRE-GA09000173EL	ASTM D5501 96% Ethanol QC Check Mix in Ethanol(‡)	2ml
DRE-GH09000173EL	ASTM D5501 96% Ethanol QC Check Mix in Ethanol(‡)	10x2ml
	ethanol [960000 mg/Kg] methanol [1000 mg/Kg] heptane (C7) [39000 mg/Kg]	
ASTM D5501 Ethanol Standard 1		
DRE-GA09000180	ASTM D5501 Ethanol Standard Mix 1(‡)	2ml
	isooctane [694000 mg/Kg] heptane (C7) [100000 mg/Kg] ethanol [200000 mg/Kg] methanol [6000 mg/Kg]	
ASTM D5501 Ethanol Standard 2		
DRE-GA09000181	ASTM D5501 Ethanol Standard Mix 2(‡)	2ml
	isooctane [395000 mg/Kg] heptane (C7) [100000 mg/Kg] ethanol [500000 mg/Kg] methanol [5000 mg/Kg]	
ASTM D5501 Ethanol Standard 3		
DRE-GA09000182	ASTM D5501 Ethanol Standard Mix 3(‡)	2ml
	isooctane [148000 mg/Kg] heptane (C7) [100000 mg/Kg] ethanol [750000 mg/Kg] methanol [3000 mg/Kg]	
ASTM D5501 Ethanol Standard 4		
DRE-GA09000183	ASTM D5501 Ethanol Standard Mix 4(‡)	2ml
	isooctane [58000 mg/Kg] heptane (C7) [40000 mg/Kg] ethanol [900000 mg/Kg] methanol [2000 mg/Kg]	

Standards for regulatory methods

Product code	Description	
ASTM D5501 Ethanol Standard 5		
DRE-GA09000184	ASTM D5501 Ethanol Standard Mix 5(±) ethanol [994000 mg/Kg] methanol [1000 mg/Kg] heptane (C7) [5000 mg/Kg]	2ml
ASTM D5501 Ethanol Standard Kit		
DRE-GK09000174	ASTM D5501 Ethanol Standard Kit(±)	1ea
DRE-GA09000180	ASTM D5501 Ethanol Standard 1(±)	1x2ml
DRE-GA09000181	ASTM D5501 Ethanol Standard 2(±)	1x2ml
DRE-GA09000182	ASTM D5501 Ethanol Standard 3(±)	1x2ml
DRE-GA09000183	ASTM D5501 Ethanol Standard 4(±)	1x2ml
DRE-GA09000184	ASTM D5501 Ethanol Standard 5(±)	1x2ml
ASTM D5501 Ethanol in Fuel Std. #1		
DRE-GA09000086HP	ASTM D5501 Ethanol in Fuel Std #1 in Heptane(±) methanol [3000 mg/Kg] ethanol [100000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. #2		
DRE-GA09000087HP	ASTM D5501 Ethanol in Fuel Std #2 in Heptane(±) methanol [3000 mg/Kg] ethanol [150000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. #3		
DRE-GA09000088HP	ASTM D5501 Ethanol in Fuel Std #3 in Heptane(±) methanol [3000 mg/Kg] ethanol [200000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. #4		
DRE-GA09000089HP	ASTM D5501 Ethanol in Fuel Std #4 in Heptane(±) methanol [3000 mg/Kg] ethanol [500000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. #5		
DRE-GA09000090HP	ASTM D5501 Ethanol in Fuel Std. #5 in Heptane(±) methanol [3000 mg/Kg] ethanol [750000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. #6		
DRE-GA09000091HP	ASTM D5501 Ethanol in Fuel Std. #6 in Heptane(±) methanol [3000 mg/Kg] ethanol [850000 mg/Kg]	2ml
ASTM D5501 Ethanol in Fuel Std. Kit		
DRE-GK09000092HP	ASTM D5501 Ethanol in Fuel Calibration Kit in Heptane(±)	1ea
DRE-GA09000086HP	ASTM D5501 Ethanol in Fuel Std. #1 (±)	1x2ml
DRE-GA09000087HP	ASTM D5501 Ethanol in Fuel Std. #2 (±)	1x2ml
DRE-GA09000088HP	ASTM D5501 Ethanol in Fuel Std. #3 (±)	1x2ml
DRE-GA09000089HP	ASTM D5501 Ethanol in Fuel Std. #4 (±)	1x2ml
DRE-GA09000090HP	ASTM D5501 Ethanol in Fuel Std. #5 (±)	1x2ml
DRE-GA09000091HP	ASTM D5501 Ethanol in Fuel Std. #6 (±)	1x2ml

Standards for regulatory methods

Product code	Description	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #1 w/ IS		
DRE-GA09000123OG	ASTM Methods D4815 & D5599 OXYCAL in Gasoline Calibration Mix #1 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [13500 mg/Kg] methyl t-butyl ether [155000 mg/Kg] tert-amyl methyl ether (TAME) [20000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] ethanol [13500 mg/Kg] tert-butyl ethyl ether (ETBE) [10000 mg/Kg] isopropyl ether [20000 mg/Kg] oxygenate free gasoline, unleaded, premium [768000 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #2 w/ IS		
DRE-GA09000124OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #2 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [50000 mg/Kg] 2-methyl-2-propanol [125000 mg/Kg] tert-amyl methyl ether (TAME) [10000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] ethanol [2500 mg/Kg] methyl t-butyl ether [50000 mg/Kg] isopropyl ether [2500 mg/Kg] oxygenate free gasoline, unleaded, premium [760000 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #3 w/ IS		
DRE-GA09000125OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #3 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [90000 mg/Kg] 2-methyl-2-propanol [50000 mg/Kg] methyl t-butyl ether [110000 mg/Kg] tert-amyl methyl ether (TAME) [5000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [745000 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #4 w/ IS		
DRE-GA09000126OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #4 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [2500 mg/Kg] ethanol [50000 mg/Kg] 2-methyl-2-propanol [90000 mg/Kg] tert-butyl ethyl ether (ETBE) [90000 mg/Kg] isopropyl ether [10000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [757500 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #5 w/ IS		
DRE-GA09000127OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #5 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	ethanol [90000 mg/Kg] 2-methyl-2-propanol [2500 mg/Kg] tert-butyl ethyl ether (ETBE) [125000 mg/Kg] tert-amyl methyl ether (TAME) [15000 mg/Kg] isopropyl ether [5000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [762500 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #6 w/ IS		
DRE-GA09000128OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #6 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	methanol [125000 mg/Kg] 2-methyl-2-propanol [12500 mg/Kg] methyl t-butyl ether [2500 mg/Kg] tert-amyl methyl ether (TAME) [1000 mg/Kg] isopropyl ether [50000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [809000 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #7 w/ IS		
DRE-GA09000129OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #7 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	ethanol [125000 mg/Kg] methyl t-butyl ether [12500 mg/Kg] tert-butyl ethyl ether (ETBE) [50000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [812500 mg/Kg]	

Standards for regulatory methods

Product code	Description	
ASTM D5599 Oxygenates in Gasoline Cal. Std. #8 w/ IS		
DRE-GA09000130OG	ASTM Methods D4815 & D5599 OXYCAL in Gasoline Calibration Mix #8 w/ IS in Oxygenate Free Gasoline(‡)	2ml
	tert-butyl ethyl ether (ETBE) [180000 mg/Kg] 1,2-dimethoxyethane (IS) [50000 mg/Kg] oxygenate free gasoline, unleaded, premium [820000 mg/Kg]	
ASTM D5599 Oxygenates in Gasoline Cal. Std. Kit w/ IS		
DRE-GK09000131OG	ASTM D5599 Oxygenates in Gasoline Calibration Kit w/ IS in Oxygenate Free Gasoline(‡)	1ea
DRE-GA09000123OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #1 w/ IS (‡)	1x2ml
DRE-GA09000124OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #2 w/ IS (‡)	1x2ml
DRE-GA09000125OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #3 w/ IS (‡)	1x2ml
DRE-GA09000126OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #4 w/ IS (‡)	1x2ml
DRE-GA09000127OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #5 w/ IS (‡)	1x2ml
DRE-GA09000128OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #6 w/ IS (‡)	1x2ml
DRE-GA09000129OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #7 w/ IS (‡)	1x2ml
DRE-GA09000130OG	ASTM D5599 Oxygenates in Gasoline Cal. Std. #8 w/ IS (‡)	1x2ml
ASTM D5769 3 Internal Standards Mix		
DRE-GA09000136	ASTM Method D5769, 3 Internal Standard Mix(‡)	5ml
DRE-GS09000136	ASTM D5769 3 Internal Standard Mix(‡)	5x5ml
	benzene-d6 [40 wt%] ethylbenzene-d10 [40 wt%] naphthalene-d8 [20 wt%]	
ASTM D5769 4 Internal Standards Mix		
DRE-GA09000137	ASTM D5769 4 Internal Standard Mix(‡)	10ml
DRE-GS09000137	ASTM D5769 4 Internal Standard Mix(‡)	5x10ml
	benzene-d6 [16,66 wt%] ethylbenzene-d10 [16,66 wt%] naphthalene-d8 [8,772 wt%] toluene-d8 [57,895 wt%]	
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #1 w/ IS		
DRE-GA09000065IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #1 with IS(‡)	1ml
	benzene [5,2 wt%] o-xylene [6,25 wt%] 3-ethyltoluene [3,13 wt%] 1,2,4-trimethylbenzene [5,2 wt%] n-butylbenzene [3,13 wt%] naphthalene [2,07 wt%] benzene-d6 (IS) [1,9 wt%] toluene [19,8 wt%] p-xylene [6,25 wt%] 4-ethyltoluene [3,13 wt%] 1,2,3-trimethylbenzene [3,13 wt%] 1,2-diethylbenzene [3,13 wt%] 1-methylnaphthalene [2,08 wt%] ethylbenzene-d10 (IS) [1,9 wt%] ethylbenzene [5,2 wt%] isopropylbenzene [3,13 wt%] 1,3,5-trimethylbenzene [3,13 wt%] indane [3,13 wt%] 1,2,4,5-tetramethylbenzene [3,13 wt%] 2-methylnaphthalene [2,07 wt%] naphthalene-d8 (IS) [1 wt%] m-xylene [6,25 wt%] n-propylbenzene [3,13 wt%] 2-ethyltoluene [3,13 wt%] 1,4-diethylbenzene [3,13 wt%] 1,2,3,5-Tetramethylbenzene [2,07 wt%] isooctane [0 wt%] toluene-d8 (IS) [6,6 wt%]	
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #1 w/o IS		
DRE-GA09000072IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #1 w/o IS in Isooctane(‡)	5ml
DRE-GA09000079IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #1 w/o IS in Isooctane(‡)	10ml
	benzene [5,2 wt%] o-xylene [6,25 wt%] 3-ethyltoluene [3,13 wt%] 1,2,4-trimethylbenzene [5,2 wt%] n-butylbenzene [3,13 wt%] naphthalene [2,07 wt%] toluene [19,8 wt%] p-xylene [6,25 wt%] 4-ethyltoluene [3,13 wt%] 1,2,3-trimethylbenzene [3,13 wt%] 1,2-diethylbenzene [3,13 wt%] 1-methylnaphthalene [2,08 wt%] ethylbenzene [5,2 wt%] isopropylbenzene [3,13 wt%] 1,3,5-trimethylbenzene [3,13 wt%] indane [3,13 wt%] 1,2,4,5-tetramethylbenzene [3,13 wt%] 2-methylnaphthalene [2,07 wt%] m-xylene [6,25 wt%] n-propylbenzene [3,13 wt%] 2-ethyltoluene [3,13 wt%] 1,4-diethylbenzene [3,13 wt%] 1,2,3,5-Tetramethylbenzene [2,07 wt%] isooctane [0 wt%]	
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #2 w/ IS		
DRE-GA09000066IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #2 with IS(‡)	1ml
	benzene [4,16 wt%] o-xylene [5 wt%] 3-ethyltoluene [2,504 wt%] 1,2,4-trimethylbenzene [4,16 wt%] n-butylbenzene [2,504 wt%] toluene [15,84 wt%] p-xylene [5 wt%] 4-ethyltoluene [2,504 wt%] 1,2,3-trimethylbenzene [2,504 wt%] 1,2-diethylbenzene [2,504 wt%] ethylbenzene [4,16 wt%] isopropylbenzene [2,504 wt%] 1,3,5-trimethylbenzene [2,504 wt%] indane [2,504 wt%] 1,2,4,5-tetramethylbenzene [2,504 wt%] m-xylene [5 wt%] n-propylbenzene [2,504 wt%] 2-ethyltoluene [2,504 wt%] 1,4-diethylbenzene [2,504 wt%] 1,2,3,5-Tetramethylbenzene [1,656 wt%]	

Standards for regulatory methods

Product code	Description		
naphthalene [1,656 wt%] benzene-d6 (IS) [1,9 wt%]	1-methylnaphthalene [1,664 wt%] ethylbenzene-d10 (IS) [1,9 wt%]	2-methylnaphthalene [1,656 wt%] naphthalene-d8 (IS) [1 wt%]	isooctane [20 wt%] toluene-d8 (IS) [6,6 wt%]

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #2 w/o IS

DRE-GA09000073IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #2 w/o IS in Isooctane(‡)			5ml
DRE-GA09000080IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #2 w/o IS in Isooctane(‡)			10ml
benzene [4,16 wt%]	toluene [15,84 wt%]	ethylbenzene [4,16 wt%]	m-xylene [5 wt%]	
o-xylene [5 wt%]	p-xylene [5 wt%]	isopropylbenzene [2,5 wt%]	n-propylbenzene [2,5 wt%]	
3-ethyltoluene [2,5 wt%]	4-ethyltoluene [2,5 wt%]	1,3,5-trimethylbenzene [2,5 wt%]	2-ethyltoluene [2,5 wt%]	
1,2,4-trimethylbenzene [4,16 wt%]	1,2,3-trimethylbenzene [2,5 wt%]	indane [2,5 wt%]	1,4-diethylbenzene [2,5 wt%]	
n-butylbenzene [2,5 wt%]	1,2-diethylbenzene [2,5 wt%]	1,2,4,5-tetramethylbenzene [1,656 wt%]	1,2,3,5-Tetramethylbenzene [1,656 wt%]	
naphthalene [1,664 wt%]	1-methylnaphthalene [1,656 wt%]	2-methylnaphthalene [1,656 wt%]	isooctane [20 wt%]	

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #3 w/ IS

DRE-GA09000067IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #3 with IS(‡)			1ml
benzene [2,6 wt%]	toluene [9,9 wt%]	ethylbenzene [2,6 wt%]	m-xylene [3,125 wt%]	
o-xylene [3,125 wt%]	p-xylene [3,125 wt%]	isopropylbenzene [1,565 wt%]	n-propylbenzene [1,565 wt%]	
3-ethyltoluene [1,565 wt%]	4-ethyltoluene [1,565 wt%]	1,3,5-trimethylbenzene [1,565 wt%]	2-ethyltoluene [1,565 wt%]	
1,2,4-trimethylbenzene [2,6 wt%]	1,2,3-trimethylbenzene [1,565 wt%]	indane [1,565 wt%]	1,4-diethylbenzene [1,565 wt%]	
n-butylbenzene [1,565 wt%]	1,2-diethylbenzene [1,565 wt%]	1,2,4,5-tetramethylbenzene [1,565 wt%]	1,2,3,5-Tetramethylbenzene [1,035 wt%]	
naphthalene [1,035 wt%]	1-methylnaphthalene [1,04 wt%]	2-methylnaphthalene [1,035 wt%]	isooctane [50 wt%]	
benzene-d6 (IS) [1,9 wt%]	ethylbenzene-d10 (IS) [1,9 wt%]	naphthalene-d8 (IS) [1 wt%]	toluene-d8 (IS) [6,6 wt%]	

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #3 w/o IS

DRE-GA09000074IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #3 w/o IS in Isooctane(‡)			5ml
DRE-GA09000081IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #3 w/o IS in Isooctane(‡)			10ml
benzene [2,6 wt%]	toluene [9,9 wt%]	ethylbenzene [2,6 wt%]	m-xylene [3,125 wt%]	
o-xylene [3,125 wt%]	p-xylene [3,125 wt%]	isopropylbenzene [1,565 wt%]	n-propylbenzene [1,565 wt%]	
3-ethyltoluene [1,565 wt%]	4-ethyltoluene [1,565 wt%]	1,3,5-trimethylbenzene [1,565 wt%]	2-ethyltoluene [1,565 wt%]	
1,2,4-trimethylbenzene [2,6 wt%]	1,2,3-trimethylbenzene [1,565 wt%]	indane [1,565 wt%]	1,4-diethylbenzene [1,565 wt%]	
n-butylbenzene [1,565 wt%]	1,2-diethylbenzene [1,565 wt%]	1,2,4,5-tetramethylbenzene [1,565 wt%]	1,2,3,5-Tetramethylbenzene [1,035 wt%]	
naphthalene [1,035 wt%]	1-methylnaphthalene [1,04 wt%]	2-methylnaphthalene [1,035 wt%]	isooctane [50 wt%]	

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #4 w/ IS

DRE-GA09000068IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #4 with IS(±)			1ml
benzene [1,56 wt%]	toluene [5,94 wt%]	ethylbenzene [1,56 wt%]	m-xylene [1,875 wt%]	
o-xylene [1,875 wt%]	p-xylene [1,875 wt%]	isopropylbenzene [0,939 wt%]	n-propylbenzene [0,939 wt%]	
3-ethyltoluene [0,939 wt%]	4-ethyltoluene [0,939 wt%]	1,3,5-trimethylbenzene [0,939 wt%]	2-ethyltoluene [0,939 wt%]	
1,2,4-trimethylbenzene [1,56 wt%]	1,2,3-trimethylbenzene [0,939 wt%]	indane [0,939 wt%]	1,4-diethylbenzene [0,939 wt%]	
n-butylbenzene [0,939 wt%]	1,2-diethylbenzene [0,939 wt%]	1,2,4,5-tetramethylbenzene [0,939 wt%]	1,2,3,5-Tetramethylbenzene [0,621 wt%]	
naphthalene [0,621 wt%]	1-methylnaphthalene [0,624 wt%]	2-methylnaphthalene [0,621 wt%]	isooctane [70 wt%]	
benzene-d6 (IS) [1,9 wt%]	ethylbenzene-d10 (IS) [1,9 wt%]	naphthalene-d8 (IS) [1 wt%]	toluene-d8 (IS) [6,6 wt%]	

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #4 w/o IS

DRE-GA09000075IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #4 w/o IS in Isooctane(‡)	5ml
DRE-GA09000082IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #4 w/o IS in Isooctane(‡)	10ml
benzene [1,56 wt%] o-xylene [1,875 wt%] 3-ethyltoluene [0,939 wt%] 1,2,4-trimethylbenzene [1,56 wt%] n-butylbenzene [0,939 wt%] naphthalene [0,621 wt%]	toluene [5,94 wt%] p-xylene [1,875 wt%] 4-ethyltoluene [0,939 wt%] 1,2,3-trimethylbenzene [0,939 wt%] 1,2-diethylbenzene [0,939 wt%] 1-methylnaphthalene [0,624 wt%]	ethylbenzene [1,56 wt%] isopropylbenzene [0,939 wt%] 1,3,5-trimethylbenzene [0,939 wt%] indane [0,939 wt%] 1,2,4,5-tetramethylbenzene [0,939 wt%] 2-methylnaphthalene [0,621 wt%] m-xylene [1,875 wt%] n-propylbenzene [0,939 wt%] 2-ethyltoluene [0,939 wt%] 1,4-diethylbenzene [0,939 wt%] 1,2,3,5-Tetramethylbenzene [0,621 wt%] isooctane [70 wt%]

ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #5 w/ IS

DRE-GA09000069IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #5 with IS(‡)			1ml
benzene [0,78 wt%]	toluene [2,97 wt%]	ethylbenzene [0,78 wt%]	m-xylene [0,94 wt%]	
o-xylene [0,94 wt%]	p-xylene [0,94 wt%]	isopropylbenzene [0,47 wt%]	n-propylbenzene [0,47 wt%]	
3-ethyltoluene [0,47 wt%]	4-ethyltoluene [0,47 wt%]	1,3,5-trimethylbenzene [0,47 wt%]	2-ethyltoluene [0,47 wt%]	
1,2,4-trimethylbenzene [0,78 wt%]	1,2,3-trimethylbenzene [0,47 wt%]	indane [0,47 wt%]	1,4-diethylbenzene [0,47 wt%]	
n-butylbenzene [0,47 wt%]	1,2-diethylbenzene [0,47 wt%]	1,2,4,5-tetramethylbenzene [0,47 wt%]	1,2,3,5-tetramethylbenzene [0,31 wt%]	
naphthalene [0,31 wt%]	1-methylnaphthalene [0,31 wt%]	2-methylnaphthalene [0,31 wt%]	isooctane [85 wt%]	
benzene-d6 (IS) [1,9 wt%]	ethylbenzene-d10 (IS) [1,9 wt%]	naphthalene-d8 (IS) [1 wt%]	toluene-d8 (IS) [6,6 wt%]	

Standards for regulatory methods

Product code	Description	
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #5 w/o IS		
DRE-GA09000076IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #5 w/o IS in Isooctane(‡)	5ml
DRE-GA09000083IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #5 w/o IS in Isooctane(‡)	10ml
benzene [0,78 wt%] o-xylene [0,94 wt%] 3-ethyltoluene [0,47 wt%] 1,2,4-trimethylbenzene [0,78 wt%] n-butylbenzene [0,47 wt%] naphthalene [0,31 wt%]	toluene [2,97 wt%] p-xylene [0,94 wt%] 4-ethyltoluene [0,47 wt%] 1,2,3-trimethylbenzene [0,47 wt%] 1,2-diethylbenzene [0,47 wt%] 1-methylnaphthalene [0,31 wt%]	ethylbenzene [0,78 wt%] isopropylbenzene [0,47 wt%] 1,3,5-trimethylbenzene [0,47 wt%] indane [0,47 wt%] 1,2,4,5-tetramethylbenzene [0,47 wt%] 2-methylnaphthalene [0,31 wt%]
		m-xylene [0,94 wt%] n-propylbenzene [0,47 wt%] 2-ethyltoluene [0,47 wt%] 1,4-diethylbenzene [0,47 wt%] 1,2,3,5-tetramethylbenzene [0,31 wt%] isooctane [85 wt%]
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #6 w/ IS		
DRE-GA09000070IO	ASTM Method D5769 Aromatics in Finished Gasoline Calibration Std #6 with IS(‡)	1ml
benzene [0,39 wt%] o-xylene [0,47 wt%] 3-ethyltoluene [0,23 wt%] 1,2,4-trimethylbenzene [0,39 wt%] n-butylbenzene [0,23 wt%] naphthalene [0,16 wt%] benzene-d6 (IS) [1,9 wt%]	toluene [1,49 wt%] p-xylene [0,47 wt%] 4-ethyltoluene [0,23 wt%] 1,2,3-trimethylbenzene [0,23 wt%] 1,2-diethylbenzene [0,23 wt%] 1-methylnaphthalene [0,16 wt%] ethylbenzene-d10 (IS) [1,9 wt%]	ethylbenzene [0,39 wt%] isopropylbenzene [0,23 wt%] 1,3,5-trimethylbenzene [0,23 wt%] indane [0,23 wt%] 1,2,4,5-tetramethylbenzene [0,23 wt%] 2-methylnaphthalene [0,16 wt%] naphthalene-d8 (IS) [1 wt%]
		m-xylene [0,47 wt%] n-propylbenzene [0,23 wt%] 2-ethyltoluene [0,23 wt%] 1,4-diethylbenzene [0,23 wt%] 1,2,3,5-tetramethylbenzene [0,16 wt%] isooctane [92,5 wt%] toluene-d8 (IS) [6,6 wt%]
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #6 w/o IS		
DRE-GA09000077IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #6 w/o IS in Isooctane(‡)	5ml
DRE-GA09000084IO	ASTM D5769 Aromatics in Finished Gasoline Cal Std #6 w/o IS in Isooctane(‡)	10ml
benzene [0,39 wt%] o-xylene [0,47 wt%] 3-ethyltoluene [0,23 wt%] 1,2,4-trimethylbenzene [0,39 wt%] n-butylbenzene [0,23 wt%] naphthalene [0,16 wt%]	toluene [1,49 wt%] p-xylene [0,47 wt%] 4-ethyltoluene [0,23 wt%] 1,2,3-trimethylbenzene [0,23 wt%] 1,2-diethylbenzene [0,23 wt%] 1-methylnaphthalene [0,16 wt%]	ethylbenzene [0,39 wt%] isopropylbenzene [0,23 wt%] 1,3,5-trimethylbenzene [0,23 wt%] indane [0,23 wt%] 1,2,4,5-tetramethylbenzene [0,23 wt%] 2-methylnaphthalene [0,16 wt%]
		m-xylene [0,47 wt%] n-propylbenzene [0,23 wt%] 2-ethyltoluene [0,23 wt%] 1,4-diethylbenzene [0,23 wt%] 1,2,3,5-tetramethylbenzene [0,16 wt%] isooctane [92,5 wt%]
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. Kit w/ IS		
DRE-GK09000071IO	ASTM D5769 Aromatics in Finished Gasoline Calibration Kit with IS in Isooctane(‡)	1ea
DRE-GA09000065IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #1 w/ IS (‡)	1x1ml
DRE-GA09000066IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #2 w/ IS (‡)	1x1ml
DRE-GA09000067IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #3 w/ IS (‡)	1x1ml
DRE-GA09000068IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #4 w/ IS (‡)	1x1ml
DRE-GA09000069IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #5 w/ IS (‡)	1x1ml
DRE-GA09000070IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #6 w/ IS (‡)	1x1ml
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. Kit w/o IS, 10mL		
DRE-GK09000085IO	ASTM D5769 Aromatics in Finished Gasoline Calibration Kit w/o IS in Isooctane(‡)	1ea
DRE-GA09000079IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #1 w/o IS (‡)	1x10ml
DRE-GA09000080IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #2 w/o IS (‡)	1x10ml
DRE-GA09000081IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #3 w/o IS (‡)	1x10ml
DRE-GA09000082IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #4 w/o IS (‡)	1x10ml
DRE-GA09000083IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #5 w/o IS (‡)	1x10ml
DRE-GA09000084IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #6 w/o IS (‡)	1x10ml
ASTM D5769 Aromatics in Finished Gasoline Cal. Std. Kit w/o IS, 5mL		
DRE-GK09000078IO	ASTM D5769 Aromatics in Finished Gasoline Calibration Kit w/o IS in Isooctane(‡)	1ea
DRE-GA09000072IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #1 w/o IS (‡)	1x5ml
DRE-GA09000073IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #2 w/o IS (‡)	1x5ml
DRE-GA09000074IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #3 w/o IS (‡)	1x5ml
DRE-GA09000075IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #4 w/o IS (‡)	1x5ml
DRE-GA09000076IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #5 w/o IS (‡)	1x5ml
DRE-GA09000077IO	ASTM D5769 Aromatics in Finished Gasoline Cal. Std. #6 w/o IS (‡)	1x5ml

Standards for regulatory methods

Product code	Description	
ASTM D5769 Quality Control Reference Material w/ 3 IS		
DRE-GA09000132IO	ASTM D5769 Quality Control Reference Material w/ 3 IS in Isooctane(‡)	10x2ml
1,2,4,5-tetramethylbenzene [20 wt%] n-dodecane (C12) [5 wt%] n-hexane (C6) [12 wt%] octane (C8) [17 wt%] o-xylene [3 wt%]	benzene [1 wt%] ethylbenzene [3 wt%] naphthalene [1 wt%] toluene [9 wt%]	benzene-d6 (IS) [2 wt%] ethylbenzene-d10 (IS) [2 wt%] naphthalene-d8 (IS) [1 wt%] 1,2,4-trimethylbenzene [3 wt%]
		n-decane (C10) [12 wt%] heptane (C7) [17 wt%] isooctane [12 wt%] m-xylene [3 wt%]
ASTM D5769 Quality Control Reference Material w/ 4 IS		
DRE-GA09000134IO	ASTM D5769 Quality Control Reference Material w/ 4 IS in Isooctane(‡)	10x2ml
1,2,4,5-tetramethylbenzene [2 wt%] n-dodecane (C12) [5 wt%] n-hexane (C6) [12 wt%] octane (C8) [17 wt%] m-xylene [3 wt%]	benzene [1 wt%] ethylbenzene [3 wt%] naphthalene [1 wt%] toluene [9 wt%] o-xylene [3 wt%]	benzene-d6 (IS) [2 wt%] ethylbenzene-d10 (IS) [2 wt%] naphthalene-d8 (IS) [1 wt%] toluene d8 (IS) [7 wt%]
		n-decane (C10) [12 wt%] heptane (C7) [17 wt%] isooctane [12 wt%] 1,2,4-trimethylbenzene [3 wt%]
ASTM D5769 Quality Control Reference Material w/o IS		
DRE-GA09000133IO	ASTM Method D5769 Quality Control Reference Material w/o IS in Isooctane(‡)	10x2ml
1,2,4,5-tetramethylbenzene [2 wt%] n-decane (C10) [12 wt%] ethylbenzene [3 wt%] n-hexane (C6) [12 wt%] isooctane [12 wt%] Toluene [9 wt%] m-xylene [3 wt%]		benzene [1 wt%] n-dodecane (C12) [5 wt%] heptane (C7) [17 wt%] naphthalene [1 wt%] Octane (C8) [17 wt%] 1,2,4-trimethylbenzene [3 wt%] o-xylene [3 wt%]
Appendix IX Volatile Mix, 11-1283		
DRE-YA09000032ME	EPA App. IX VOC Mix in p/t Methanol(‡)	1ml
DRE-YS09000032ME	EPA App. IX VOC Mix in p/t Methanol(‡)	5x1ml
	acetone[nitrile [10000 mg/L] 1-butanol [20000 mg/L] ethyl methacrylate [2000 mg/L] isobutyl alcohol [20000 mg/L] methyl methacrylate [2000 mg/L] propionitrile [10000 mg/L]	allyl chloride [2000 mg/L] chloroprene [2000 mg/L] hexachloroethane [2000 mg/L] methyl acrylonitrile [10000 mg/L] pentachloroethane [2000 mg/L]
Base/Neutral Surrogate Solution		
DRE-AY09000008DI	Base/Neutral Surrogate Solution 5000 µg/mL in Dichloromethane(‡)	5ml
DRE-YA09000008DI	EPA 8270 B/N Surrogate Mix 5000 µg/mL in Dichloromethane(‡)	1ml
DRE-SY09000008DI	EPA 8270 B/N Surrogate Mix 5000 µg/mL in Dichloromethane(‡)	5x1ml
	nitrobenzene-d5 2-fluorobiphenyl p-terphenyl-d14	
CLP Acids Mix, 17-52		
DRE-AY09000023DI	CLP Acids Mix 1000 µg/mL in Dichloromethane(‡)	5ml
DRE-SY09000023DI	CLP Acids Mix 1000 µg/mL in Dichloromethane(‡)	5x5ml
DRE-YA09000023DI	CLP Acids Mix 1000 µg/mL in Dichloromethane(‡)	1ml
DRE-YS09000023DI	CLP Acids Mix 1000 µg/mL in Dichloromethane(‡)	5x1ml
benzoic acid 2,4-dimethylphenol 4-methylphenol 2,4,5-trichlorophenol phenol	4-chloro-3-methylphenol 2-methyl-4,6-dinitrophenol 2-nitrophenol 2,4,6-trichlorophenol	2-chlorophenol 2,4-dinitrophenol 4-nitrophenol 2,6-dichlorophenol
		2,4-dichlorophenol 2-methylphenol pentachlorophenol 2,3,4,6-Tetrachlorophenol
EPA 1625 Hydrocarbons Mix		
DRE-GA09000160DI	EPA 1625 Hydrocarbons Mix 4000 µg/mL in Dichloromethane(‡)	1ml
n-decane (C10) n-dodecane (C12) hexacosane (C26) octacosane (C28) n-tetracosane (C24) triacontane (C30)		n-docosane (C22) n-eicosane (C20) n-hexadecane (C16) n-octadecane (C18) n-tetradecane (C14)

Standards for regulatory methods

Product code	Description	
EPA 502/524.2 VOC Solution		
DRE-XA09000001ME	EPA VOC Mix 1 200 µg/mL in p/t Methanol(‡)	1ml
DRE-XS09000001ME	EPA VOC Mix 1 200 µg/mL in p/t Methanol(‡)	5x1ml
benzene	ethylbenzene	m-xylene
isopropylbenzene	n-propylbenzene	o-xylene
sec-butylbenzene	tert-butylbenzene	1,2,4-trimethylbenzene
n-butylbenzene	naphthalene	4-isopropyltoluene
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene
4-chlorotoluene	chlorobenzene	1,2,3-trichlorobenzene
bromobenzene	bromochloromethane	carbon tetrachloride
methylene chloride	bromodichloromethane	bromoform
dibromochloromethane	cis-1,2-dichloroethylene	trans-1,2-dichloroethylene
1,1-dichloroethane	1,1,1-trichloroethane	2,2-dichloropropane
1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane	1,1,2-trichloroethane
trichloroethylene	1,2-dibromo-3-chloropropane	1,2-dibromoethane
1,1-dichloropropylene	1,2,3-trichloropropane	1,2-dichloropropane
cis-1,3-dichloropropylene	1,3-dichloropropane	bromomethane
chloroethane	dichlorodifluoromethane	vinyl chloride
		toluene
		p-xylene
		1,3,5-trimethylbenzene
		styrene
		2-chlorotoluene
		1,2,4-trichlorobenzene
		dibromomethane
		chloroform
		1,1-dichloroethylene
		tetrachloroethylene
		1,2-dichloroethane
		hexachlorobutadiene
		trans-1,3-dichloropropylene
		chloromethane
		trichlorofluoromethane
EPA 515 Herbicide Mix		
DRE-YS09000050AC	EPA 515 Herbicide Mix in Acetone(‡)	5x1ml
MCPP acid [10000 mg/L]	MCPA acid [10000 mg/L]	3,5-dichlorobenzoic acid [100 mg/L]
pentachlorophenol [100 mg/L]	2,4,5-TP (Silvex) [100 mg/L]	2,4,5-T [100 mg/L]
tetrachloroterephthalic acid [100 mg/L]	4-nitrophenol [100 mg/L]	dichloroprop (2,4-DP) [100 mg/L]
2,4-D [100 mg/L]	acifluorfen [100 mg/L]	dalapon [100 mg/L]
bentazon [100 mg/L]	chloramben [100 mg/L]	Dicamba [100 mg/L]
		picloram [100 mg/L]
		dinoseb [100 mg/L]
		2,4-DB [100 mg/L]
EPA 524 Internal/Surrogate Standard Solution		
DRE-YA09000002ME	EPA 524 IS/SS Mix 2000 µg/mL in p/t Methanol(‡)	1ml
DRE-YS09000002ME	EPA 524 IS/SS Mix 2000 µg/mL in p/t Methanol(‡)	5x1ml
	4-bromofluorobenzene (BFB)	
	1,2-dichloroethane-d4	
	fluorobenzene	
EPA 524.2 Drinking Water VOA Mix		
DRE-YA09000018ME	EPA VOC Mix 2 2000 µg/mL in p/t Methanol(‡)	1ml
DRE-YS09000018ME	EPA VOC Mix 2 2000 µg/mL in p/t Methanol(‡)	5x1ml
benzene	ethylbenzene	m-xylene
isopropylbenzene	n-propylbenzene	o-xylene
sec-butylbenzene	tert-butylbenzene	1,2,4-trimethylbenzene
n-butylbenzene	naphthalene	4-isopropyltoluene
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene
4-chlorotoluene	chlorobenzene	1,2,3-trichlorobenzene
bromobenzene	bromochloromethane	carbon tetrachloride
methylene chloride	bromodichloromethane	bromoform
dibromochloromethane	cis-1,2-dichloroethylene	trans-1,2-dichloroethylene
1,1-dichloroethane	1,1,1-trichloroethane	2,2-dichloropropane
1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane	1,1,2-trichloroethane
trichloroethylene	1,2-dibromo-3-chloropropane	1,2-dibromoethane
1,1-dichloropropylene	1,2,3-trichloropropane	1,2-dichloropropane
cis-1,3-dichloropropylene	1,3-dichloropropane	2-nitropropane
ethyl methacrylate	hexachloroethane	methyl methacrylate
acrylonitrile	iodomethane	carbon disulfide
methyl acrylonitrile	nitrobenzene	pentachloroethane
1-chlorobutane	ethyl ether	methyl t-butyl ether
methyl acrylate		toluene
		p-xylene
		1,3,5-trimethylbenzene
		styrene
		2-chlorotoluene
		1,2,4-trichlorobenzene
		dibromomethane
		chloroform
		1,1-dichloroethylene
		tetrachloroethylene
		1,2-dichloroethane
		hexachlorobutadiene
		trans-1,3-dichloropropylene
		allyl chloride
		tetrahydrofuran
		trans-1,4-dichloro-2-butene
		chloroacetonitrile
		propionitrile
EPA 610/ 8100 PAH Mix		
DRE-GA09000161BD	EPA 610/ 8100 PAH Mix 2000 µg/mL in MeCl ₂ :Benzene 1:1(‡)	1ml
acenaphthene	acenaphthylene	anthracene
benzo[b]fluoranthene	benzo[k]fluoranthene	benzo[ghi]perylene
chrysene	fluoranthene	fluorene
naphthalene	phenanthrene	pyrene
		benzo[a]anthracene
		benzo[a]pyrene
		indeno[1,2,3-cd]pyrene
		dibenz[a,h]anthracene

Standards for regulatory methods

Product code	Description		
EPA 624 Purgeable Calibration Mix 9			
DRE-YA06240900ME	EPA 624 Purgeable Calibration Mix 9, 2000 µg/mL in Methanol (*)	1ml	
Benzene	Bromodichloromethane	Bromoform	Carbontetrachloride
Chlorobenzene	Chloroform	Dibromochloromethane	1,2-Dichlorobenzene
1,3-Dichlorobenzene	1,4-Dichlorobenzene	1,1-Dichloroethane	1,2-Dichloroethane
1,1-Dichloroethene	1,2-Dichloroethene	1,2-Dichloropropane	1,3-Dichloropropene
Ethylbenzene	Methylene chloride	1,1,2,2-Tetrachloroethane	Tetrachloroethene
Toluene	1,1,1-Trichloroethane	1,1,2-Trichloroethane	Trichloroethene
EPA 625 Base/Neutral Mix			
DRE-GH09000202DI	EPA 625 Base/Neutral Mix 200 µg/mL in Dichloromethane(‡)	5x1ml	
benzidine	benzo[a]pyrene		
benzo[ghi]perylene	benzo[k]fluoranthene		
4-chlorodiphenyl ether	di-n-octyl phthalate		
indeno[1,2,3-cd]pyrene	N-nitrosodimethylamine		
EPA 8010 Halogenated VOC Mix 2			
DRE-YA08010200ME	EPA 8010 Halogenated VOC Mix 2 2000 µg/mL in Methanol (*) (‡)	1ml	
1,1,1,2-Tetrachloroethane	1,1,1-Trichloroethane		
1,1,2,2-Tetrachloroethane	1,1,2-Trichloroethane		
1,1-Dichloroethane	1,2,3-Trichloropropane		
1,2-Dichloroethane	1,2-Dichloropropane		
Dibromomethane	Tetrachloroethene		
Tetrachloromethane	Tribromomethane		
Trichloroethene	Trichloromethane		
EPA 8140/8141 Organophosphorous Pesticide Solution A			
DRE-LA09000015HA	EPA 8140/8141 Organophosphorous Pesticide Mix 40 µg/mL in Hexane:Acetone: 95:5(‡)	1ml	
DRE-LS09000015HA	EPA 8140/8141 Organophosphorous Pesticide Mix 40 µg/mL in Hexane:Acetone: 95:5(‡)	5x1ml	
chlorpyrifos	trichloronate	sulprofos	coumaphos
Demeton O&S	fensulfothion	fenthion	fenchlorphos
Tokuthion®	tributylphosphoro-trithioite (Merphos)	tetrachlorvinphos (Rabon)	dibrom
diazinon	Guthion®	dichlorvos	disulfoton
ethoprophos (prophos)	methyl parathion	phosdrin™ (mevinphos)	phorate
EPA 8260 Additional Compounds			
DRE-YA09000012ME	EPA VOC Additional CompoundsMix 2000 µg/mL in p/t Methanol(‡)	1ml	
DRE-YS09000012ME	EPA VOC Additional CompoundsMix 2000 µg/mL in p/t Methanol(‡)	5x1ml	
acetone	2-butanone (MEK)		
4-methyl-2-pentanone (MIBK)	2-hexanone		
2-chloroethylvinyl ether	iodomethane		
carbon disulfide	vinyl acetate		
EPA 8260 Internal Standard Solution			
DRE-YA09000011ME	EPA 8260 Internal Standard Mix 2000 µg/mL in p/t Methanol(‡)	1ml	
DRE-YS09000011ME	EPA 8260 Internal Standard Mix 2000 µg/mL in p/t Methanol(‡)	5x1ml	
chlorobenzene-d5			
1,4-dichlorobenzene-d4			
fluorobenzene			
EPA 8260 Internal/Surrogate Solution			
DRE-YA09000014ME	EPA 8260 IS/SS Mix 2500 µg/mL in p/t Methanol(‡)	1ml	
DRE-YS09000014ME	EPA 8260 IS/SS Mix 2500 µg/mL in p/t Methanol(‡)	5x1ml	
1,4-dichlorobenzene-d4			
chlorobenzene-d5			
fluorobenzene			
dibromofluoromethane			
toluene-d8			
4-bromofluorobenzene (BFB)			
1,2-dichloroethane-d4			

Standards for regulatory methods

Product code	Description		
EPA 8260 VOC Liquids			
DRE-YA09000010ME	EPA VOC Mix 3 2000 µg/mL in p/t Methanol(±)	1ml	
DRE-YS09000010ME	EPA VOC Mix 3 2000 µg/mL in p/t Methanol(±)	5x1ml	
benzene	ethylbenzene	m-xylene	
isopropylbenzene	n-propylbenzene	o-xylene	
sec-butylbenzene	tert-butylbenzene	1,2,4-trimethylbenzene	
n-butylbenzene	naphthalene	4-isopropyltoluene	
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene	
4-chlorotoluene	chlorobenzene	1,2,3-trichlorobenzene	
bromobenzene	bromochloromethane	carbon tetrachloride	
methylene chloride	bromodichloromethane	bromoform	
dibromochloromethane	cis-1,2-dichloroethylene	trans-1,2-dichloroethylene	
1,1-dichloroethane	1,1,1-trichloroethane	2,2-dichloropropane	
1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane	1,1,2-trichloroethane	
trichloroethylene	1,2-dibromo-3-chloropropane	1,2-dibromoethane	
1,1-dichloropropylene	1,2,3-trichloropropane	1,2-dichloropropane	
cis-1,3-dichloropropylene	1,3-dichloropropane		
		toluene	
		p-xylene	
		1,3,5-trimethylbenzene	
		styrene	
		2-chlorotoluene	
		1,2,4-trichlorobenzene	
		dibromomethane	
		chloroform	
		1,1-dichloroethylene	
		tetrachloroethylene	
		1,2-dichloroethane	
		hexachlorobutadiene	
		trans-1,3-dichloropropylene	
EPA 8260 VOC Solution			
DRE-YA09000013ME	EPA VOC Mix 1 2000 µg/mL in p/t Methanol(±)	1ml	
DRE-YS09000013ME	EPA VOC Mix 1 2000 µg/mL in p/t Methanol(±)	5x1ml	
benzene	ethylbenzene	m-xylene	
isopropylbenzene	n-propylbenzene	o-xylene	
sec-butylbenzene	tert-butylbenzene	1,2,4-trimethylbenzene	
n-butylbenzene	naphthalene	4-isopropyltoluene	
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene	
4-chlorotoluene	chlorobenzene	1,2,3-trichlorobenzene	
bromobenzene	bromochloromethane	carbon tetrachloride	
methylene chloride	bromodichloromethane	bromoform	
dibromochloromethane	cis-1,2-dichloroethylene	trans-1,2-dichloroethylene	
1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane	1,1,2-trichloroethane	
trichloroethylene	1,2-dibromo-3-chloropropane	1,2-dibromoethane	
1,1-dichloropropylene	1,2,3-trichloropropane	1,2-dichloropropane	
cis-1,3-dichloropropylene	1,3-dichloropropane	trichlorofluoromethane	
chloromethane	chloroethane	dichlorodifluoromethane	
1,1-dichloroethane	1,1,1-trichloroethane	2,2-dichloropropane	
		toluene	
		p-xylene	
		1,3,5-trimethylbenzene	
		styrene	
		2-chlorotoluene	
		1,2,4-trichlorobenzene	
		dibromomethane	
		chloroform	
		1,1-dichloroethylene	
		1,2-dichloroethane	
		hexachlorobutadiene	
		trans-1,3-dichloropropylene	
		bromomethane	
		vinyl chloride	
		tetrachloroethylene	
EPA 8270 Acid Surrogate Mix			
DRE-YA09000006ME	EPA 8270 Acid Surrogate Mix 2000 µg/mL in Methanol(±)	1ml	
DRE-YS09000006ME	EPA 8270 Acid Surrogate Mix 2000 µg/mL in Methanol(±)	5x1ml	
DRE-YA090000025ME	EPA 8270 Acid Surrogate Mix 10000 µg/mL in Methanol(±)	1ml	
DRE-AY090000025ME	EPA 8270 Acid Surrogate Mix 10000 µg/mL in Methanol(±)	5ml	
DRE-SY090000025ME	EPA 8270 Acid Surrogate Mix 10000 µg/mL in Methanol(±)	5x5ml	
	2-fluorophenol		
	2,4,6-tribromophenol		
	phenol-d6		
EPA 8270 Base/Neutral Mix 10			
DRE-YA09000022DI	EPA 8270 B/N Mix 1 1000 µg/mL in Dichloromethane(±)	1ml	
DRE-YS09000022DI	EPA 8270 B/N Mix 1 1000 µg/mL in Dichloromethane(±)	5x1ml	
DRE-AY09000022DI	EPA 8270 B/N Mix 1 1000 µg/mL in Dichloromethane(±)	5ml	
DRE-SY09000022DI	EPA 8270 B/N Mix 1 1000 µg/mL in Dichloromethane(±)	5x5ml	
	aniline	4-chloroaniline	
	2-nitroaniline	3-nitroaniline	
	4-nitroaniline	pyridine	
	carbazole	benzyl alcohol	
	dibenzofuran	2-methylnaphthalene	
EPA 8270 Base/Neutral Solution, 44-2			
DRE-YA09000024SP	EPA 8270 B/N Mix 2 1000 µg/mL in Benz:MeCl2:ACCN 4:4:2(±)	1ml	
DRE-YS09000024SP	EPA 8270 B/N Mix 2 1000 µg/mL in Benz:MeCl2:ACCN 4:4:2(±)	5x1ml	
DRE-AY09000024SP	EPA 8270 B/N Mix 2 1000 µg/mL in Benz:MeCl2:ACCN 4:4:2(±)	5ml	
DRE-SY09000024SP	EPA 8270 B/N Mix 2 1000 µg/mL in Benz:MeCl2:ACCN 4:4:2(±)	5x5ml	
bis(2-chloroethoxy)methane	bis(2-chloroethyl)ether	bis(2-chloro-1-methylethyl) ether	4-bromophenyl phenyl ether
4-chlorodiphenyl ether	bis(2-ethylhexyl)phthalate	butyl benzyl phthalate	diethyl phthalate

Standards for regulatory methods

Product code	Description			
	dimethyl phthalate	di-n-butyl phthalate	di-n-octyl phthalate	2-chloronaphthalene
	1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene	hexachlorobenzene
	hexachlorobutadiene	hexachlorocyclopentadiene	hexachloroethane	1,2,4-trichlorobenzene
	N-nitrosodi-n-propylamine	N-nitrosodiphenylamine	N-nitrosodimethylamine	azobenzene
	nitrobenzene	isophorone	2,6-dinitrotoluene	2,4-dinitrotoluene
	acenaphthene	acenaphthylene	anthracene	benzo[a]anthracene
	benzo[b]fluoranthene	benzo[k]fluoranthene	benzo[ghi]perylene	benzo[a]pyrene
	chrysene	dibenz[a,h]anthracene	fluoranthene	fluorene
	indeno[1,2,3-cd]pyrene	naphthalene	phenanthrene	pyrene

EPA 8270 Base/Neutral Surrogate Solution

DRE-YA09000007DI	EPA 8270 B/N Surrogate Mix 1000 µg/mL in Dichloromethane(‡)	1ml
DRE-YS09000007DI	EPA 8270 B/N Surrogate Mix 1000 µg/mL in Dichloromethane(‡)	5x1ml
	nitrobenzene-d5	
	2-fluorobiphenyl	
	p-terphenyl-d14	

EPA 8270 Internal Standard Mix

DRE-YS09000038DI	EPA 8270 Internal Standard Mix 2000 µg/mL in Dichloromethane(‡)	5x1ml
	acenaphthene-d10	
	chrysene-d12	
	1,4-dichlorobenzene-d4	
	naphthalene-d8	
	perylene-d12	
	phenanthrene-d10	

EPA 8270 Internal Standard Solution

DRE-YA09000005DI	EPA 8270 Internal Standard Mix 4000 µg/mL in Dichloromethane(‡)	1ml
DRE-YS09000005DI	EPA 8270 Internal Standard Mix 4000 µg/mL in Dichloromethane(‡)	5x1ml
	acenaphthene-d10	
	chrysene-d12	
	1,4-dichlorobenzene-d4	
	naphthalene-d8	
	perylene-d12	
	phenanthrene-d10	

EPA 8270 LCS Solution, 78-1

DRE-XA09000004AD	EPA 8270 LCS Mix 100 µg/mL in Acetone:MeCl2 83.5:16.5 (‡)			25ml
	2-chlorophenol	2,4-dimethylphenol	pentachlorophenol	4-nitrophenol
	2,4-dichlorophenol	4-chloro-3-methylphenol	2-methyl-4,6-dinitrophenol	2-nitrophenol
	2,4-dinitrophenol	2,4,6-trichlorophenol	phenol	2-methylphenol
	3-methylphenol	4-methylphenol	2,4,5-trichlorophenol	2,4-dinitrotoluene
	2,6-dinitrotoluene	isophorone	nitrobenzene	aniline
	4-chloroaniline	2-nitroaniline	3-nitroaniline	4-nitroaniline
	pyridine	carbazole	bis(2-chloroethoxy)methane	bis(2-chloroethyl)ether
	bis(2-chloro-1-methylethyl) ether	4-bromophenyl phenyl ether	4-chlorodiphenyl ether	bis(2-ethylhexyl)phthalate
	butyl benzyl phthalate	diethyl phthalate	dimethyl phthalate	di-n-butyl phthalate
	di-n-octyl phthalate	2-chloronaphthalene	1,2-dichlorobenzene	1,3-dichlorobenzene
	1,4-dichlorobenzene	hexachlorobenzene	hexachlorobutadiene	hexachlorocyclopentadiene
	hexachloroethane	1,2,4-trichlorobenzene	benzo[k]fluoranthene	1-methylnaphthalene
	2-methylnaphthalene	acenaphthene	acenaphthylene	anthracene
	fluorene	naphthalene	phenanthrene	benzo[a]anthracene
	benzo[a]pyrene	chrysene	fluoranthene	indeno[1,2,3-cd]pyrene
	pyrene	benzo[b]fluoranthene	benzo[ghi]perylene	dibenz[a,h]anthracene
	benzyl alcohol	bis(2-ethylhexyl)adipate	dibenzofuran	1,2-dinitrobenzene
	1,3-dinitrobenzene	1,4-dinitrobenzene	azobenzene	2,3,4,6-Tetrachlorophenol
	2,3,5,6-tetrachlorophenol	N-nitrosodimethylamine	3,3'-dichlorobenzidine	benzoic acid
	N-nitrosodiphenylamine	N-nitrosodi-n-propylamine		

EPA Method 525.2 Organochlorine Pesticide Mix, 27-3612

DRE-GA09000336AC	EPA Method 525.2 Organochlorine Pesticide Mix, 27-3612 500 µg/mL in Acetone(‡)			1ml
	aldrin	a-BHC	b-BHC	d-BHC
	g-BHC	cis-chlordane	trans-chlordane	p,p'-DDD
	p,p'-DDE	p,p'-DDT	dieldrin	endosulfan I
	endosulfan II	endosulfan sulfate	endrin	endrin aldehyde
	endrin ketone	heptachlor	methoxychlor	chlorobenzilate
	chloroneb	chlorothalonil	chlorthal-dimethyl (dacthal)	heptachlor epoxide isomer B
	trans-nonachlor			

Standards for regulatory methods

Product code	Description		
EPA Method 525.2 Organophosphorus Pesticide Mix			
DRE-GA09000337AC	EPA Method 525.2 Organophosphorus Pesticide Mix 500 µg/mL in Acetone(‡)1ml		
	chlorpyrifos disulfoton sulfone ethoprophos (prophos) phosdrinTM (mevinphos)	dichlorvos disulfoton sulfoxide methyl paraoxon tetrachlorvinphos (Rabon)	
EPA Method 525.2 PCB Congeners Mix, 8-1046			
DRE-GA09000340AC	EPA Method 525.2 PCB Congeners Mix, 8-1046 500 µg/mL in Acetone(‡)1ml		
	2-chlorobiphenyl (BZ# 1) 2,2',4,4'-tetrachlorobiphenyl (BZ# 47) 2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171) 2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 201)	2,3-dichlorobiphenyl (BZ# 5) 2,4,5-trichlorobiphenyl (BZ# 29) 2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154) 2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	
EPA Method 525.2 Semi-Volatile Mix, 29-5			
DRE-GA09000338AC	EPA Method 525.2 Semi-Volatile Mix, 29-5 in Acetone(‡)1ml		
	acenaphthylene [1000 mg/L] benzo[a]pyrene [1000 mg/L] butyl benzyl phthalate [1000 mg/L] dibenz[a,h]anthracene [1000 mg/L] 2,4-dinitrotoluene [1000 mg/L] fluorene [1000 mg/L] isophorone [1000 mg/L] pyrene [1000 mg/L]	acetochlor [1000 mg/L] benzo[b]fluoranthene [1000 mg/L] bis(2-ethylhexyl)adipate [1000 mg/L] diethyl phthalate [1000 mg/L] 2,6-dinitrotoluene [1000 mg/L] hexachlorobenzene [1000 mg/L] naphthalene [1000 mg/L]	anthracene [1000 mg/L] benzo[ghi]perylene [1000 mg/L] bis(2-ethylhexyl)phthalate [1000 mg/L] dimethyl phthalate [1000 mg/L] di-n-octyl phthalate [1000 mg/L] hexachlorocyclopentadiene [1000 mg/L] pentachlorophenol [4000 mg/L] benzo[a]anthracene [1000 mg/L] benzo[k]fluoranthene [1000 mg/L] chrysene [1000 mg/L] di-n-butyl phthalate [1000 mg/L] fluoranthene [1000 mg/L] indeno[1,2,3-cd]pyrene [1000 mg/L] phenanthrene [1000 mg/L]
EPA Method 525.3 Organochlorine Pesticide Mix, 29-3612			
DRE-GS09000341AC	EPA Method 525.3 Organochlorine Pesticide Mix, 29-3612 500 µg/mL in Acetone(‡)5x1ml		
	acetochlor chlorobenzilate p,p'-DDD endosulfan I a-BHC heptachlor methoxychlor	aldrin chloroneb p,p'-DDE endosulfan II b-BHC heptachlor epoxide isomer B trans-nonachlor	cis-chlordane chlorothalonil p,p'-DDT endosulfan sulfate d-BHC hexachlorobenzene pentachlorophenol trans-chlordane chlorthal-dimethyl (dacthal) dieldrin endrin g-BHC hexachlorocyclopentadiene
EPA Method 525.3 Organonitrogen Pesticide Mix, 21-3615			
DRE-GS09000344AC	EPA Method 525.3 Organonitrogen Pesticide Mix, 21-3615 500 µg/mL in Acetone(‡)5x1ml		
	molinate oxyfluorfen propyzamide (pronamide) simazine terbutryne vinclozolin	napropamide pebulate propachlor tebuconazol (Folicur) triadimefon	Nitrofen prometon propazine tebuthiuron trifluralin norflurazon prometryn simetryn terbacil vernolate
EPA Method 525.3 Organonitrogen Pesticide Mix, 22-3614			
DRE-GS09000343AC	EPA Method 525.3 Organonitrogen Pesticide Mix, 22-3614 500 µg/mL in Acetone(‡)5x1ml		
	alachlor bromacil chlorpropham 2,4-dinitrotoluene etridiazole metolachlor	ametryne butachlor bladex 2,6-dinitrotoluene fenarimol MGK-264 - isomer a	atrazine 2,6-di-tert-butyl-4-methylphenol (BHT) N,N-diethyl-m-toluamide EPTC (s-ethyl dipropylthiocarbamate) hexazinon
EPA Method 525.3 Organophosphate Pesticide Mix, 16-3613			
DRE-GS09000342AC	EPA Method 525.3 Organophosphate Pesticide Mix, 16-3613 500 µg/mL in Acetone(‡)5x1ml		
	chlorfenvinphos (E/Z mixture) dimethipin parathion phosphamidon	chlorpyrifos disulfoton methyl parathion profenofos	dichlorvos ethion phosdrinTM (mevinphos) tetrachlorvinphos (ISO) diisopropyl methylphosphonate ethoprophos (prophos) phorate tribufos

Standards for regulatory methods

Product code	Description	
EPA Method 525.3 PAH Mix, 20-2158		
DRE-GS09000345AC	EPA Method 525.3 PAH Mix, 20-2158 500 µg/mL in Acetone(‡)	5x1ml
acenaphthylene benzo[b]fluoranthene chrysene bis(2-ethylhexyl)phthalate indeno[1,2,3-cd]pyrene	anthracene benzo[ghi]perylene dibenz[a,h]anthracene diethyl phthalate isophorone	benzo[a]anthracene benzo[k]fluoranthene di-n-butyl phthalate dimethyl phthalate phenanthrene
		benzo[a]pyrene butyl benzyl phthalate bis(2-ethylhexyl)adipate fluorene pyrene
EPA Method 525.3 PCB Mix, 14-3616		
DRE-GS09000346AC	EPA Method 525.3 PCB Mix, 14-3616 500 µg/mL in Acetone(‡)	5x1ml
	4-chlorobiphenyl (BZ# 3) 2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180) 2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149) 2,3,3',4',6-pentachlorobiphenyl (BZ# 110) 2,2',3,5'-tetrachlorobiphenyl (BZ# 44) 2,3',4',5-tetrachlorobiphenyl (BZ# 70) 2,4,4'-trichlorobiphenyl (BZ# 28)	2,4'-dichlorobiphenyl (BZ# 8) 2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138) 2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153) 2,3',4,4',5-pentachlorobiphenyl (BZ# 118) 2,2',5,5'-tetrachlorobiphenyl (BZ# 52) 2,2',5-trichlorobiphenyl (BZ# 18) 2-chlorobiphenyl (BZ# 1)
EPA Method 525.3 Surrogate Mix, 3-2159		
DRE-GS09000348AC	EPA Method 525.3 Surrogate Mix, 3-2159 500 µg/mL in Acetone(‡)	5x1ml
	benzo(a)pyrene-d12 2-nitro-m-xylene triphenyl phosphate	
EPA Method 525.5		
DRE-YA09000036EA	EPA 525.5 Mix 500 µg/mL in Ethyl Acetate(‡)	1ml
alachlor bis(2-ethylhexyl)adipate heptachlor methoxychlor simazine 2,4-dinitrotoluene bis(2-ethylhexyl)phthalate	aldrin butachlor heptachlor epoxide isomer B metolachlor trans-nonachlor 2,6-dinitrotoluene	atrazine dieldrin hexachlorobenzene metribuzin hexachlorocyclopentadiene cis-chlordane
		benzo[a]pyrene endrin g-BHC propachlor bladex trans-chlordane
EPA Method 530 Internal Standard Mix, 2-20		
DRE-GA09000351AC	EPA Method 530 Internal Standard Mix, 2-20 500 µg/mL in Acetone(‡)	1ml
	acenaphthene-d10 phenanthrene-d10	
EPA Method 530 Mix 1, 2-1864		
DRE-GA09000349ME	EPA Method 530 Mix 1, 2-1864 100 µg/mL in p/t Methanol(‡)	1ml
	quinoline o-toluidine	
EPA Method 530 Mix 2, 2-1865		
DRE-GA09000350ME	EPA Method 530 Mix 2, 2-1865 100 µg/mL in Methanol(‡)	1ml
	butylated hydroxyanisole dimethipin	
EPA Method 538 Mix, 11-2173		
DRE-GA09000363ME	EPA Method 538 Mix, 11-2173 in Methanol(‡)	1ml
	acephate [1000 mg/L] aldicarb sulfoxide [1000 mg/L] diisopropyl methylphosphonate [1000 mg/L] fenamiphos sulfoxide [1000 mg/L] oxydemeton-methyl [1000 mg/L] thiofanox [1000 mg/L]	aldicarb [1000 mg/L] dicrotophos [1000 mg/L] fenamiphos sulfone [1000 mg/L] methamidophos [1000 mg/L] quinoline [20000 mg/L]

Standards for regulatory methods

Product code	Description	
EPA Method 541 Mix		
DRE-GS09000360ME	EPA Method 541 Mix 2000 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000360ME-SS	EPA Method 541 Mix 2000 µg/mL in p/t Methanol (Second Source)(‡)	1ml
	allyl alcohol 1-butanol 1,4-dioxane 2-methoxyethanol	
EPA Method 8260 Acetate Mix, 7 components		
DRE-GA09000415ME	EPA Method 8260 Acetate Mix, 2000 µg/mL in p/t Methanol(‡)	1ml
	n-amyl acetate butyl acetate ethyl acetate isopropyl acetate methyl acetate propyl acetate vinyl acetate	
EPA Method 8260 Calibration Check Compounds, 6 components		
DRE-GA09000449ME	EPA Method 8260 Calibration Check Compounds, 2000 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000449ME	EPA Method 8260 Calibration Check Compounds, 2000 µg/mL in p/t Methanol(‡)	5x1ml
	chloroform 1,1-dichloroethene 1,2-dichloropropane ethylbenzene toluene vinyl chloride	
EPA Method 8260 Gases Mix, 6 components		
DRE-YA09000009ME	EPA VOC Gases Mix 2000 µg/mL in p/t Methanol(‡)	1ml
DRE-YS09000009ME	EPA VOC Gases Mix 2000 µg/mL in p/t Methanol(‡)	5x1ml
DRE-GA09000413ME	EPA Method 8260 Gases Mix, 200 µg/mL in p/t Methanol(‡)	1ml
	bromomethane chloroethane chloromethane dichlorodifluoromethane trichlorofluoromethane vinyl chloride	
EPA Method 8260 Internal Standard Mix, 4 components		
DRE-GA09000416ME	EPA Method 8260 Internal Standard Mix, 2500 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000416ME	EPA Method 8260 Internal Standard Mix, 2500 µg/mL in p/t Methanol(‡)	5x1ml
	chlorobenzene-d5 1,4-dichlorobenzene-d4 1,4-difluorobenzene pentafluorobenzene	
EPA Method 8260 Matrix Spike Mix, 5 components		
DRE-GA09000418ME	EPA Method 8260 Matrix Spike Mix, 2500 µg/mL in p/t Methanol(‡)	1ml
	benzene chlorobenzene 1,1-dichloroethene toluene trichloroethene	
EPA Method 8260 Oxygenates Mix, 6 components		
DRE-GA09000414ME	EPA Method 8260 Oxygenates Mix, 2000/10000 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000414ME	EPA Method 8260 Oxygenates Mix, 2000/10000 µg/mL in p/t Methanol(‡)	5x1ml
	4,4-dimethyl-3-oxahexane (TAEE) [2000 µg/mL] tert-amyl methyl ether (TAME) [2000 µg/mL] tert-butyl ethyl ether (ETBE) [2000 µg/mL]	

Standards for regulatory methods

Product code	Description
	isopropyl ether [2000 µg/mL]
	2-methyl-2-propanol [10000 µg/mL]
	methyl t-butyl ether [2000 µg/mL]

EPA Method 8260 Surrogates Mix, 4 components

DRE-GA09000417ME	EPA Method 8260 Surrogates Mix, 2500 µg/mL in p/t Methanol(±)	1ml
	4-bromofluorobenzene (BFB)	
	dibromofluoromethane	
	1,2-dichloroethane-d4	
	toluene-d8	

EPA Method 8260 System Performance Check Compounds, 5 components

DRE-GA09000448ME	EPA Method 8260 System Performance Check Compounds, 2000 µg/mL in p/t Methanol(±)	1ml
DRE-GS09000448ME	EPA Method 8260 System Performance Check Compounds, 2000 µg/mL in p/t Methanol(±)	5x1ml
	bromoform	
	chlorobenzene	
	chloromethane	
	1,1-dichloroethane	
	1,1,2,2-tetrachloroethane	

EPA Method 8260 UltiMix Calibration Mix, 75 components

DRE-GA09000412ME	UltiMix 8260 Calibration Mix, 2000 µg/mL in p/t Methanol(±)	1ml	
DRE-GS09000412ME	EPA Method 8260 UltiMix Calibration Mix, 2000 µg/mL in p/t Methanol(±)	5x1ml	
acetonitrile	acrylonitrile	allyl chloride	benzene
bromobenzene	bromochloromethane	bromodichloromethane	bromoform
n-butylbenzene	sec-butylbenzene	tert-butylbenzene	carbon disulfide
carbon tetrachloride	chlorobenzene	2-chloroethanol	chloroform
chloroprene	2-chlorotoluene	4-chlorotoluene	cis-1,2-dichloroethylene
dibromochloromethane	1,2-dibromo-3-chloropropane	1,2-dibromoethane	dibromomethane
1,2-dichlorobenzene	1,3-dichlorobenzene	1,4-dichlorobenzene	cis-1,4-dichloro-2-butene
trans-1,4-dichloro-2-butene	1,1-dichloroethane	1,2-dichloroethane	1,1-dichloroethylene
trans-1,2-dichloroethylene	1,2-dichloropropane	1,3-dichloropropane	2,2-dichloropropane
1,1-dichloropropylene	cis-1,3-dichloropropylene	trans-1,3-dichloropropylene	1,4-dioxane
ethyl ether	ethyl methacrylate	ethylbenzene	hexachlorobutadiene
iodomethane	isobutyl alcohol	isopropylbenzene	4-isopropyltoluene
methyl acrylate	methyl acrylonitrile	methyl methacrylate	methylene chloride
naphthalene	nitrobenzene	2-nitropropane	propionitrile
n-propylbenzene	styrene	1,1,1,2-tetrachloroethane	1,1,2,2-tetrachloroethane
tetrachloroethylene	tetrahydrofuran	toluene	1,2,3-trichlorobenzene
1,2,4-trichlorobenzene	1,1,1-trichloroethane	1,1,2-trichloroethane	trichloroethylene
1,2,3-trichloropropane	1,1,2-trichloro-1,2,2-trifluoroethane	1,2,4-trimethylbenzene	1,3,5-trimethylbenzene
m-xylene	o-xylene	p-xylene	

EPA Method 8270 Acid Calibration Check Mix, 6 components

DRE-GA09000435DI	EPA Method 8270 Acid Calibration Check Mix, 2000 µg/mL in Dichloromethane(±)	1ml
	4-chloro-3-methylphenol	
	2,4-dichlorophenol	
	2-nitrophenol	
	pentachlorophenol	
	phenol	
	2,4,6-trichlorophenol	

EPA Method 8270 Appendix IX Mix 1 Minus Phentermine, 16 components

DRE-GS09000444DI	EPA Method 8270 Appendix IX Mix 1 Minus Phentermine, 12000 µg/mL in Dichloromethane(±)			5x1ml
DRE-GA09000444DI	EPA Method 8270 Appendix IX Mix 1 Minus Phentermine, 2000 µg/mL in Dichloromethane(±)			1ml
2-acetylaminofluorene	4-aminobiphenyl	p-(dimethylamino)azobenzene	isosafole	
1-naphthylamine	2-naphthylamine	5-nitro-o-toluidine	N-nitrosodiethylamine	
N-nitrosodi-n-butylamine	N-nitrosomethylethylamine	N-nitrosomorpholine	N-nitrosopiperidine	
N-nitrosopyrrolidine	p-phenylenediamine	2-picoline	o-toluidine	

Standards for regulatory methods

Product code	Description	
EPA Method 8270 Appendix IX Mix 1 with Phentermine, 17 components		
DRE-GS09000443DI	EPA Method 8270 Appendix IX Mix 1 with Phentermine, 2000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000443DI	EPA Method 8270 Appendix IX Mix 1 with Phentermine, 2000 µg/mL in Dichloromethane(‡)	1ml
2-acetylaminofluorene	4-aminobiphenyl	p-(dimethylamino)azobenzene
isosafrrole	1-naphthylamine	2-naphthylamine
N-nitrosodiethylamine	N-nitrosodi-n-butylamine	N-nitrosomethylethylamine
N-nitrosopiperidine	N-nitrosopyrrolidine	p-phenylenediamine
o-toluidine		phentermine
		5-nitro-o-toluidine
		N-nitrosomorpholine
		2-picoline
EPA Method 8270 Appendix IX Mix 2, 31 components		
DRE-GA09000445DI	EPA Method 8270 Appendix IX Mix 2, 1000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000445DI	EPA Method 8270 Appendix IX Mix 2, 1000 µg/mL in Dichloromethane(‡)	5x1ml
acetophenone	aramite	benzaldehyde
caprolactam	chlorobenzilate	1-chloronaphthalene
dibenz(a,j)acridine	2,6-dichlorophenol	7,12-dimethylbenz[a]anthracene
diphenyl ether	ethyl methacrylate	ethyl methanesulfonate
isodrin	isosafrrole	kepone
methyl methanesulfonate	1,4-naphthoquinone	4-nitroquinoline-1-oxide
pentachloroethane	pentachloronitrobenzene	phenacetin
safrole	1,2,4,5-tetrachlorobenzene	1,3,5-trinitrobenzene
		biphenyl
		di-allate (mixture of isomers)
		1,4-dioxane
		hexachloropropene
		3-methylcholanthrene
		pentachlorobenzene
		propylamide (pronamide)
EPA Method 8270 B/N Calibration Check Mix, 7 components		
DRE-GA09000434DI	EPA Method 8270 B/N Calibration Check Mix, 2000 µg/mL in Dichloromethane(‡)	1ml
	acenaphthene	
	benzo[a]pyrene	
	1,4-dichlorobenzene	
	di-n-octyl phthalate	
	diphenylamine	
	fluoranthene	
	hexachlorobutadiene	
EPA Method 8270 BNA Surrogates, 7 components		
DRE-GA09000431DI	EPA Method 8270 BNA Surrogates Mix, 1000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000431DI	EPA Method 8270 BNA Surrogates Mix, 1000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000432DI	EPA Method 8270 BNA Surrogates Mix, 5000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000432DI	EPA Method 8270 BNA Surrogates Mix, 5000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000433DI	EPA Method 8270 BNA Surrogates Mix, 5000 µg/mL in Dichloromethane(‡)	5ml
	2-fluorobiphenyl	
	2-fluorophenol	
	nitrobenzene-d5	
	phenol-d5	
	pyrene-d10	
	p-terphenyl-d14	
	2,4,6-tribromophenol	
EPA Method 8270 Benzidines Mix, 3 components		
DRE-GA09000426ME	EPA Method 8270 Benzidines Mix, 2000 µg/mL in Methanol(‡)	1ml
DRE-GS09000426ME	EPA Method 8270 Benzidines Mix, 2000 µg/mL in Methanol(‡)	5x1ml
	benzidine	
	3,3'-dichlorobenzidine	
	3,3'-dimethylbenzidine	
EPA Method 8270 Calibration Mix, 76 components		
DRE-YA09000003DI	EPA 8270 Calibration Mix in Dichloromethane(‡)	1ml
DRE-YS09000003DI	EPA 8270 Calibration Mix in Dichloromethane(‡)	5x1ml
DRE-GA09000419BD	EPA Method 8270 UltiMix Calibration Mix, 500/1000 µg/mL in Benzene/Dichloromethane 1:1(‡)	1ml
acenaphthene [1000 µg/mL]	acenaphthylene [1000 µg/mL]	aniline [1000 µg/mL]
azobenzene [1000 µg/mL]	benzo[a]anthracene [1000 µg/mL]	benzo[b]fluoranthene [1000 µg/mL]
benzo[ghi]perylene [1000 µg/mL]	benzo[a]pyrene [1000 µg/mL]	benzyl alcohol [1000 µg/mL]
bis(2-chloroethyl)ether [1000 µg/mL]	bis(1-chloroprop-2-yl) ether [1000 µg/mL]	4-bromophenyl phenyl ether [1000 µg/mL]
carbazole [1000 µg/mL]	4-chloroaniline [1000 µg/mL]	4-chlorodiphenyl ether [1000 µg/mL]
2-chloronaphthalene [1000 µg/mL]	2-chlorophenol [1000 µg/mL]	chrysene [1000 µg/mL]
dibenzofuran [1000 µg/mL]	di-n-butyl phthalate [1000 µg/mL]	1,2-dichlorobenzene [1000 µg/mL]
1,4-dichlorobenzene [1000 µg/mL]	2,4-dichlorophenol [1000 µg/mL]	diethyl phthalate [1000 µg/mL]
		anthracene [1000 µg/mL]
		benzo[k]fluoranthene [1000 µg/mL]
		bis(2-chloroethoxy)methan [1000 µg/mL]
		butyl benzyl phthalate [1000 µg/mL]
		4-chloro-3-methylphenol [1000 µg/mL]
		dibenz[a,h]anthracene [1000 µg/mL]
		1,3-dichlorobenzene [1000 µg/mL]
		2,4-dimethylphenol [1000 µg/mL]

(‡) ISO Guide 34

(*) Shorter expiry due to chemical nature of component(s)

Order online at www.lgcstandards.com

Standards for regulatory methods

Product code	Description			
	dimethyl phthalate [1000 µg/mL]	1,2-dinitrobenzene [1000 µg/mL]	1,3-dinitrobenzene [1000 µg/mL]	1,4-dinitrobenzene [1000 µg/mL]
	2,4-dinitrophenol [1000 µg/mL]	2,4-dinitrotoluene [1000 µg/mL]	2,6-dinitrotoluene [1000 µg/mL]	di-n-octyl phthalate [1000 µg/mL]
	diphenylamine [1000 µg/mL]	2,3,5,6-tetrachlorophenol [1000 µg/mL]	fluoranthene [1000 µg/mL]	fluorene [1000 µg/mL]
	hexachlorobenzene [1000 µg/mL]	hexachlorobutadiene [1000 µg/mL]	hexachlorocyclopentadiene [1000 µg/mL]	hexachloroethane [1000 µg/mL]
	indeno[1,2,3-cd]pyrene [1000 µg/mL]	isophorone [1000 µg/mL]	2-methyl-4,6-dinitrophenol [1000 µg/mL]	1-methylnaphthalene [1000 µg/mL]
	2-methylnaphthalene [1000 µg/mL]	2-methylphenol [1000 µg/mL]	3-methylphenol [500 µg/mL]	4-methylphenol [500 µg/mL]
	naphthalene [1000 µg/mL]	2-nitroaniline [1000 µg/mL]	3-nitroaniline [1000 µg/mL]	4-nitroaniline [1000 µg/mL]
	nitrobenzene [1000 µg/mL]	2-nitrophenol [1000 µg/mL]	4-nitrophenol [1000 µg/mL]	N-nitrosodimethylamine [1000 µg/mL]
	N-nitrosodi-n-propylamine [1000 µg/mL]	pentachlorophenol [1000 µg/mL]	phenanthrene [1000 µg/mL]	phenol [1000 µg/mL]
	pyrene [1000 µg/mL]	pyridine [1000 µg/mL]	2,3,4,6-tetrachlorophenol [1000 µg/mL]	1,2,4-trichlorobenzene [1000 µg/mL]
	2,4,5-trichlorophenol [1000 µg/mL]	2,4,6-trichlorophenol [1000 µg/mL]	dioctyl phthalate [1000 µg/mL]	bis(2-ethylhexyl) adipate [1000 µg/mL]

EPA Method 8270 Calibration by Class Mix 1, 22 components

DRE-GA09000437DI	EPA Method 8270 Calibration by Class Mix 1, 2000 µg/mL in Dichloromethane(‡)			1ml
	benzoic acid	4-chloro-3-methylphenol	2-chlorophenol	2,4-dichlorophenol
	2,6-dichlorophenol	2,4-dimethylphenol	2,4-dinitrophenol	dinoseb
	2-fluorophenol	hexachlorophene	2-methyl-4,6-dinitrophenol	2-methylphenol
	3-methylphenol	4-methylphenol	2-nitrophenol	4-nitrophenol
	pentachlorophenol	phenol	2,3,4,6-tetrachlorophenol	2,4,6-tribromophenol
	2,4,5-trichlorophenol	2,4,6-trichlorophenol		

EPA Method 8270 Calibration by Class Mix 2, 11 components

DRE-GA09000438DI	EPA Method 8270 Calibration by Class Mix 2, 2000 µg/mL in Dichloromethane(‡)			1ml
	aniline	benzidine		
	4-chloroaniline	3,3'-dichlorobenzidine		
	diphenylamine	2-nitroaniline		
	3-nitroaniline	4-nitroaniline		
	N-nitrosodimethylamine	N-nitrosodi-n-propylamine		
	pyridine			

EPA Method 8270 Calibration by Class Mix 3, 23 components

DRE-GA09000439DI	EPA Method 8270 Calibration by Class Mix 3, 2000 µg/mL in Dichloromethane(‡)			1ml
	aramite	bis(2-chloroethoxy)methane	bis(2-chloroethyl)ether	Bis(2-chloro-1-methylethyl) ether
	chlorobenzilate	4-chlorodiphenyl ether	2-chloronaphthalene	1,2-dichlorobenzene
	1,3-dichlorobenzene	1,4-dichlorobenzene	1,3-dinitrobenzene	hexachlorobenzene
	hexachlorobutadiene	hexachlorocyclopentadiene	hexachloroethane	hexachloropropene
	isodrin	kepone	pentachlorobenzene	pentachloronitrobenzene
	1,2,4,5-tetrachlorobenzene	1,2,4-trichlorobenzene	4-bromophenylphenyl ether	

EPA Method 8270 Calibration by Class Mix 4, 22 components

DRE-GA09000440DI	EPA Method 8270 Calibration by Class Mix 4, 2000 µg/mL in Dichloromethane(‡)			1ml
	acetophenone	azobenzene	benzyl alcohol	butyl benzyl phthalate
	dibenzofuran	di-n-butyl phthalate	diethyl phthalate	dimethyl phthalate
	2,4-dinitrotoluene	2,6-dinitrotoluene	di-n-octyl phthalate	ethyl methanesulfonate
	isophorone	isofarole	methyl methanesulfonate	1,4-naphthoquinone
	nitrobenzene	4-nitroquinoline-1-oxide	phenacetin	safole
	1,3,5-trinitrobenzene	dioctyl phthalate		

EPA Method 8270 Calibration by Class Mix 5, 19 components

DRE-GA09000441DI	EPA Method 8270 Calibration by Class Mix 5, 2000 µg/mL in Dichloromethane(‡)			1ml
	acenaphthene	acenaphthylene	anthracene	benzo[a]anthracene
	benzo[b]fluoranthene	benzo[k]fluoranthene	benzo[ghi]perylene	benzo[a]pyrene
	chrysene	dibenz[a,h]anthracene	fluoranthene	fluorene
	indeno[1,2,3-cd]pyrene	3-methylcholanthrene	1-methylnaphthalene	2-methylnaphthalene
	naphthalene	phenanthrene	pyrene	

EPA Method 8270 Calibration by Class Mix 6, 10 components

DRE-GA09000442DI	EPA Method 8270 Calibration by Class Mix 6, 2000 µg/mL in Dichloromethane(‡)			1ml
	di-allate (mixture of isomers)	dimethoate		
	disulfoton	famphur		
	methyl parathion	parathion		
	phorate	propyzamide (pronamide)		
	thionazine (zinophos)	O,O,O-triethylphosphorothioate		

Standards for regulatory methods

Product code	Description	
EPA Method 8270 Color Changing Surrogates, 6 components		
DRE-GA09000446ME	EPA Method 8270 Color Changing Surrogates, 100/200 µg/mL in Methanol(‡)	25ml
DRE-GA09000447ME	EPA Method 8270 Color Changing Surrogates, 100/150 µg/mL in Methanol(‡)	25ml
	2-fluorobiphenyl [100 µg/mL] 2-fluorophenol [150 µg/mL] nitrobenzene-d5 [100 µg/mL] phenol-d5 [150 µg/mL] p-terphenyl-d14 [100 µg/mL] 2,4,6-tribromophenol [150 µg/mL]	
EPA Method 8270 Color Changing Surrogates, 7 components		
DRE-GA09000430ME	EPA Method 8270 Color Changing Surrogates, 100 µg/mL in Methanol(‡)	25ml
	2-fluorobiphenyl 2-fluorophenol nitrobenzene-d5 phenol-d5 pyrene-d10 p-terphenyl-d14 2,4,6-tribromophenol	
EPA Method 8270 GC/MS Tuning Mix, 4 components		
DRE-GA09000425ME	EPA Method 8270 GC/MS Tuning Mix, 500 µg/mL in Methanol(‡)	1ml
DRE-GS09000425ME	EPA Method 8270 GC/MS Tuning Mix, 500 µg/mL in Methanol(‡)	5x1ml
	benzidine p,p'-DDT decafluorotriphenylphosphine (DFTPP) pentachlorophenol	
EPA Method 8270 Internal Standard Mix, 7 components		
DRE-GA09000428DI	EPA Method 8270 Internal Standard Mix, 2000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000428DI	EPA Method 8270 Internal Standard Mix, 2000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000429DI	EPA Method 8270 Internal Standard Mix, 4000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000429DI	EPA Method 8270 Internal Standard Mix, 4000 µg/mL in Dichloromethane(‡)	5x1ml
	acenaphthene-d10 chrysene-d12 1,4-dichlorobenzene-d4 1,4-dioxane-d8 naphthalene-d8 perylene-d12 phenanthrene-d10	
EPA Method 8270 Matrix Spike Mix, 11 components		
DRE-GA09000422DI	EPA Method 8270 Matrix Spike Mix, 2000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000422DI	EPA Method 8270 Matrix Spike Mix, 2000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000423DI	EPA Method 8270 Matrix Spike Mix, 5000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000423DI	EPA Method 8270 Matrix Spike Mix, 5000 µg/mL in Dichloromethane(‡)	5x1ml
DRE-GA09000424DI	EPA Method 8270 Matrix Spike Mix, 10000 µg/mL in Dichloromethane(‡)	1ml
DRE-GS09000424DI	EPA Method 8270 Matrix Spike Mix, 10000 µg/mL in Dichloromethane(‡)	5x1ml
	acenaphthene 2-chlorophenol 2,4-dinitrotoluene N-nitrosodi-n-propylamine phenol 1,2,4-trichlorobenzene	4-chloro-3-methylphenol 1,4-dichlorobenzene 4-nitrophenol pentachlorophenol pyrene
EPA Method 8270 Matrix Spike Mix, 76 components		
DRE-GA09000420MD	EPA Method 8270 UltiMix Matrix Spike Mix, 100/200 µg/mL in Methanol/Dichloromethane 80:20(‡)	5ml
DRE-GA09000421MD	EPA Method 8270 UltiMix Matrix Spike Mix, 100/200 µg/mL in Methanol/Dichloromethane 80:20(‡)	2x5ml
acenaphthene [200 µg/mL] azobenzene [200 µg/mL] benzo[ghi]perylene [200 µg/mL] bis(2-chloroethyl)ether [200 µg/mL] carbazole [200 µg/mL]	acenaphthylene [200 µg/mL] benzo[a]anthracene [200 µg/mL] benzo[a]pyrene [200 µg/mL] bis(1-chloroprop-2-yl) ether [200 mg/L] 4-chloroaniline [200 µg/mL]	aniline [200 µg/mL] benzo[b]fluoranthene [200 µg/mL] benzyl alcohol [200 µg/mL] 4-bromophenyl phenyl ether [200 µg/mL] 4-chlorodiphenyl ether [200 µg/mL]
	anthracene [200 µg/mL] benzo[k]fluoranthene [200 µg/mL] bis(2-chloroethoxy)methane [200 µg/mL] butyl benzyl phthalate [200 µg/mL] 4-chloro-3-methylphenol [200 µg/mL]	

Standards for regulatory methods

Product code	Description			
	2-chloronaphthalene [200 µg/mL]	2-chlorophenol [200 µg/mL]	chrysene [200 µg/mL]	dibenz[a,h]anthracene [200 µg/mL]
	dibenzofuran [200 µg/mL]	di-n-butyl phthalate [200 µg/mL]	1,2-dichlorobenzene [200 µg/mL]	1,3-dichlorobenzene [200 µg/mL]
	1,4-dichlorobenzene [200 µg/mL]	2,4-dichlorophenol [200 µg/mL]	diethyl phthalate [200 µg/mL]	2,4-dimethylphenol [200 µg/mL]
	dimethyl phthalate [200 µg/mL]	1,2-dinitrobenzene [200 µg/mL]	1,3-dinitrobenzene [200 µg/mL]	1,4-dinitrobenzene [200 µg/mL]
	2,4-dinitrophenol [200 µg/mL]	2,4-dinitrotoluene [200 µg/mL]	2,6-dinitrotoluene [200 µg/mL]	di-n-octyl phthalate [200 µg/mL]
	diphenylamine [200 µg/mL]	2,3,5,6-tetrachlorophenol [200 µg/mL]	fluoranthene [200 µg/mL]	fluorene [200 µg/mL]
	hexachlorobenzene [200 µg/mL]	hexachlorobutadiene [200 µg/mL]	hexachlorocyclopentadiene [200 µg/mL]	hexachloroethane [200 µg/mL]
	indeno[1,2,3-cd]pyrene [200 µg/mL]	isophorone [200 µg/mL]	2-methyl-4,6-dinitrophenol [200 µg/mL]	1-methylnaphthalene [200 µg/mL]
	2-methylnaphthalene [200 µg/mL]	2-methylphenol [200 µg/mL]	3-methylphenol [100 µg/mL]	4-methylphenol [100 µg/mL]
	naphthalene [200 µg/mL]	2-nitroaniline [200 µg/mL]	3-nitroaniline [200 µg/mL]	4-nitroaniline [200 µg/mL]
	nitrobenzene [200 µg/mL]	2-nitrophenol [200 µg/mL]	4-nitrophenol [200 µg/mL]	N-nitrosodimethylamine [200 µg/mL]
	N-nitrosodi-n-propylamine [200 µg/mL]	pentachlorophenol [200 µg/mL]	phenanthrene [200 µg/mL]	phenol [200 µg/mL]
	pyrene [200 µg/mL]	pyridine [200 µg/mL]	2,3,4,6-tetrachlorophenol [200 µg/mL]	1,2,4-trichlorobenzene [200 µg/mL]
	2,4,5-trichlorophenol [200 µg/mL]	2,4,6-trichlorophenol [200 µg/mL]	diethyl phthalate [200 µg/mL]	bis(2-ethylhexyl) adipate [200 µg/mL]

EPA Method 8270 Organochlorine Pesticides Mix, 21 components

DRE-GA09000427MB	EPA Method 8270 Organochlorine Pesticides Mix, 1000 µg/mL in Methyl-tert-butyl ether(±)			1ml
DRE-GS09000427MB	EPA Method 8270 Organochlorine Pesticides Mix, 1000 µg/mL in Methyl-tert-butyl ether(±)			5x1ml
aldrin	α-BHC	β-BHC	δ-BHC	
γ-BHC	cis-chlordane	p,p'-DDD	p,p'-DDE	
p,p'-DDT	dieldrin	endosulfan I	endosulfan II	
endosulfan sulfate	endrin	endrin aldehyde	endrin ketone	
heptachlor	heptachlor epoxide isomer B	hydroquinone	methoxychlor	
trans-chlordane				

EPA Method 8270 SV System Performance Check Mix, 4 components

DRE-GA09000436DI	EPA Method 8270 SV System Performance Check Mix, 2000 µg/mL in Dichloromethane(‡)	1ml
	2,4-dinitrophenol	
	hexachlorocyclopentadiene	
	4-nitrophenol	
	N-nitrosodi-n-propylamine	

Ketones Solution

DRE-YA09000017MW	Ketones Mix 5000 µg/mL in p/t MeOH:Water 9:1(‡)	1ml
DRE-YS09000017MW	Ketones Mix 5000 µg/mL in p/t MeOH:Water 9:1(‡)	5x1ml
DRE-GA09000017MW	Ketones Mix, Max. Diff. from Nom.:1.5% 5000 µg/mL in p/t MeOH:Water 9:1(‡)	1.3ml
	acetone	
	2-butanone (MEK)	
	4-methyl-2-pentanone (MIBK)	
	2-hexanone	

Method 525.2 Internal Standard Solution

DRE-GA09000332AC	Method 525.2 Internal Standard Mix 1000 µg/mL in Acetone(‡)	1ml
	acenaphthene-d10	
	chrysene-d12	
	phenanthrene-d10	

Method 525.2 Nitrogen/Phosphorus Pesticide Mx, 6-13

DRE-GA09000335AC	Method 525.2 Nitrogen/Phosphorus Pesticide Mix, 6-13 1000 µg/mL in Acetone(‡)	1ml
	carboxine	
	diazinon	
	disulfoton	
	fenamiphos	
	tributylphosphoro-trithioite (Merphos)	
	terbufos	

Method 525.2 Organonitrogen Pesticide Mix, 36-1

DRE-GA09000339AC	Method 525.2 Organonitrogen Pesticide Mix, 36-1 500 µg/mL in Acetone(‡)			1ml
alachlor	ametryne	atrazine	bromacil	
chlorpropham	bladex	Cycloate	diphenamid	
EPTC (s-ethyl dipropylthiocarbamate)	etridiazole	Fenarimol	fluridone	
hexazinon	metolachlor	MGK-264 Mix of Isomers (70%a,30%b)	molinate	
napropamide	norflurazon	pebulate	prometryn	
propyzamide (pronamide)	propachlor	simazine	simetryn	
tebuthiuron	terbacil	terbutryne	triadimefon	
tricyclazole	trifluralin	vernolate	butachlor	

Standards for regulatory methods

Product code	Description			
	atraton	prometon	propazine	butylate
Method 525.2 Surrogate Solution				
DRE-GA09000333AC	Method 525.2 Surrogate Solution Mix 1000 µg/mL in Acetone(‡)			1ml
	1,3-dimethyl-2-nitrobenzene perylene-d12 pyrene-d10 triphenyl phosphate			
Method 552 Haloacetic Acids HAA9 Mix, 9-6601				
DRE-GS09000399MB	Method 552 Haloacetic Acids HAA9 Mix, 9-6601 in Methyl tert-Butyl Ether(‡)			5x1ml
	trichloroacetic acid [500 mg/L] bromochloroacetic acid [300 mg/L] dichloroacetic acid [200 mg/L] chloroacetic acid [2000 mg/L] bromodichloroacetic acid [500 mg/L]		bromoacetic acid [300 mg/L] dibromoacetic acid [300 mg/L] tribromoacetic acid [2000 mg/L] chlorodibromoacetic acid [300 mg/L]	
Method 552 Halogenated Acetic Mix				
DRE-GA09000362MB	Method 552 Halogenated Acetic Mix in Methyl tert-Butyl Ether(‡)			1ml
DRE-GA09000362ME	EPA 504.1 Mix 200 µg/mL in p/t Methanol(‡)			1ml
	1,2-dibromo-3-chloropropane 1,2-dibromoethane 1,2,3-trichloropropane			
Method 552 Halogenated Acetic Mix, 10-1				
DRE-GA09000361MB	Method 552 Halogenated Acetic Mix, 10-1 100 µg/mL in Methyl tert-Butyl Ether(‡)			1ml
	trichloroacetic acid bromochloroacetic acid dichloroacetic acid bromodichloroacetic acid tribromoacetic acid		dibromoacetic acid chloroacetic acid bromoacetic acid chlorodibromoacetic acid dalapon	
Method 8270B - Organophosphorous Pesticide Mix, 9-2				
DRE-GA09000364DI	Method 8270B - Organophosphorous Pesticide Mix, 9-2 2000 µg/mL in Dichloromethane(‡)			1ml
	dimethoate famphur methyl parathion tetraethyl dithiopyrophosphate O,O,O-triethylphosphorothioate		disulfoton parathion phorate thionazine (zinophos)	
NJDEP EPH Aromatics Calibration Standard				
DRE-YA09000021DI	EPH NJ Aromatics Mix 2000 µg/mL in Dichloromethane(‡)			1ml
DRE-YS09000021DI	EPH NJ Aromatics Mix 2000 µg/mL in Dichloromethane(‡)			5x1ml
	acenaphthene benzo[a]pyrene chrysene indeno[1,2,3-cd]pyrene pyrene	acenaphthylene benzo[b]fluoranthene dibenz[a,h]anthracene 2-methylnaphthalene 1,2,3-trimethylbenzene	anthracene benzo[ghi]perylene fluoranthene naphthalene	benzo[a]anthracene benzo[k]fluoranthene fluorene phenanthrene
SVOA Calibration Mixture, Method 8270				
DRE-GA09000397DI	SVOA Calibration Mix, Method 8270 50 µg/mL in Dichloromethane(‡)			5ml
	nitrobenzene-d5 4-chloroaniline pyridine 1,4-dichlorobenzene hexachloroethane bis(2-chloro-1-methylethyl) ether butyl benzyl phthalate di-n-octyl phthalate azobenzene chrysene benzo[b]fluoranthene 2,4-dimethylphenol 4-chloro-3-methylphenol	2-fluorobiphenyl 2-nitroaniline 2-chloronaphthalene hexachlorobenzene 1,2,4-trichlorobenzene 4-bromophenyl phenyl ether diethyl phthalate 2,4-dinitrotoluene N-nitrosodimethylamine fluoranthene benzo[ghi]perylene pentachlorophenol 2-methyl-4,6-dinitrophenol	p-terphenyl-d14 3-nitroaniline 1,2-dichlorobenzene hexachlorobutadiene bis(2-chloroethoxy)methane 4-chlorophenylphenyl ether dimethyl phthalate 2,6-dinitrotoluene benzo[a]anthracene indeno[1,2,3-cd]pyrene dibenz[a,h]anthracene 4-nitrophenol 2-nitrophenol	aniline 4-nitroaniline 1,3-dichlorobenzene hexachlorocyclopentadiene bis(2-chloroethyl)ether bis(2-ethylhexyl)phthalate di-n-butyl phthalate isophorone benzo[a]pyrene pyrene 2-chlorophenol 2,4-dichlorophenol 2,4-dinitrophenol

(‡) ISO Guide 34

(*) Shorter expiry due to chemical nature of component(s)

Order online at www.lgcstandards.com

Standards for regulatory methods

Product code	Description			
2,4,6-trichlorophenol	phenol	anthracene	phenanthrene	
fluorene	acenaphthene	acenaphthylene	naphthalene	
2-methylnaphthalene	3,3'-dichlorobenzidine	benzidine	2-methylphenol	
4-methylphenol	2,4,5-trichlorophenol	benzo[k]fluoranthene	benzoic acid	
benzyl alcohol	dibenzofuran	phenol-d5	2,4,6-tribromophenol	
2-fluorophenol	carbazole	nitrobenzene	n-nitrosodiphenylamine	
N-nitrosodi-n-propylamine				

Semi-Volatile Mix 1, 18 components

DRE-YA09000019DI	Semi-Volatile Mix 1 2000 µg/mL in Dichloromethane(‡)	1ml
DRE-YS09000019DI	Semi-Volatile Mix 1 2000 µg/mL in Dichloromethane(‡)	5x1ml

hexachlorocyclopentadiene	7,12-dimethylbenz[a]anthracene	acetophenone	2,4-dinitrophenol
isopropylbenzene	1,4-dioxane	1-methylnaphthalene	caprolactam
atrazine	biphenyl	2,6-dimethylnaphthalene	n-decane (C10)
n-octadecane (C18)	2,3-dichloroaniline	benzaldehyde	indene
quinoline	(L)-α-terpineol		

UCMR 4 Method 525.3 Internal Standard Solution

DRE-GA09000263AC	UCMR 4 Method 525.3 Internal Standard Solution 500 µg/mL in Acetone(‡)	1ml
	acenaphthene-d10	
	chrysene-d12	
	phenanthrene-d10	

UCMR 4 Method 525.3 Pesticide Mix, 9-2423

DRE-GA09000262ME	UCMR 4 Method 525.3 Pesticide Mix, 9-2423 1000 µg/mL in Methanol(‡)	1ml
	a-BHC	chlorpyrifos
	dimethipin	ethoprophos (prophos)
	Oxyfluorfen	profenofos
	tebuconazol (Folicur)	permethrin (mixture of isomers)
	tribufos	

UCMR 4 Method 525.3 Surrogate Mix, 3-1709

DRE-GA09000264AC	UCMR 4 Method 525.3 Surrogate Mix, 3-1709 500 µg/mL in Acetone(‡)	1ml
	1,3-dimethyl-2-nitrobenzene	
	triphenyl phosphate	
	benzo(a)pyrene-d12	

UCMR 4 Method 530 Internal Standard Mix

DRE-GA09000267ME	UCMR 4 Method 530 Internal Standard Mix 500 µg/mL in Methanol(‡)	1ml
	acenaphthene-d10	
	phenanthrene-d10	

UCMR 4 Method 530 Mix

DRE-GA09000265ME	UCMR 4 Method 530 Mix 1000 µg/mL in Methanol(‡)	1ml
	butylated hydroxyanisole	
	o-toluidine	
	quinoline	

UCMR 4 Method 530 Surrogate Mix

DRE-GA09000266ME	UCMR 4 Method 530 Surrogate Mix 500 µg/mL in Methanol(‡)	1ml
	o-toluidine-d9	
	quinoline-d7	


UCMR 4 Method 541 Alcohol Mix

DRE-GA09000268ME	UCMR 4 Method 541 Alcohol Mix 1000 µg/mL in p/t Methanol(‡)	1ml
	1-butanol	
	2-methoxyethanol	
	allyl alcohol	

Standards for regulatory methods

Product code	Description	
UCMR 4 Method 541 Internal Standard Mix		
DRE-GA09000269ME	UCMR 4 Method 541 Internal Standard Mix 1000 µg/mL in p/t Methanol(‡) 1,4-dioxane-d8 chlorobenzene-d5	1ml
UCMR 4 Method 541 Surrogate Mix		
DRE-GA09000270ME	UCMR 4 Method 541 Surrogate Mix 1000 µg/mL in p/t Methanol(‡) allyl alcohol-d6 1-butanol-d10	1ml










STABLE ISOTOPE LABELLED COMPOUNDS












Stable isotope labelled compounds

Product code	Description			
AHTN (Tonalide) (6-acetyl D3)				
CAS 1396967-82-0	MW 261.4169	C ₁₈ ² H ₁₃ H ₂₃ O		
DRE-XA10048600IO	AHTN (Tonalide) D3 100 µg/mL in Isooctane		1.1ml	
Acenaphthene D10				
CAS 15067-26-2	MW 164.2694	C ₁₂ ² H ₁₀		
DRE-C20505100	Acenaphthene D10		100mg	
DRE-L20505100CY	Acenaphthene D10 10 µg/mL in Cyclohexane		10ml	
DRE-YA20505100TO	Acenaphthene D10 2000 µg/mL in Toluene		1ml	
Acenaphthylene D8				
CAS 93951-97-4	MW 160.2412	C ₁₂ ² H ₈		
DRE-C20510100	Acenaphthylene D8		100mg	
DRE-L20510100CY	Acenaphthylene D8 10 µg/mL in Cyclohexane		10ml	
Acepromazine-d6 Maleate				
CAS 1331655-50-5	MW 448.5649	C ₁₉ ² H ₁₆ H ₁₆ N ₂ OS·C ₄ H ₄ O ₄		
DRE-C10010320	Acepromazine D6 maleate		10mg	
Acetamidrid D3 (N-methyl D3)				
CAS 1353869-35-8	MW 225.6926	C ₁₀ ² H ₉ H ₉ ClN ₄		
DRE-C10013100	Acetamidrid D3 (N-methyl D3)		50mg	
DRE-XA10013100AC	Acetamidrid D3 (N-methyl D3) 100 µg/mL in Acetone		1.1ml	
Acetochlor D11				
CAS 1189897-44-6	MW 280.8349	C ₁₄ ² H ₁₁ H ₉ ClNO ₂		
DRE-XA10018100AC	Acetochlor D11 100 µg/mL in Acetone		1.1ml	
Acetylsulfamethoxazole D4				
CAS 1215530-54-3	MW 299.339	C ₁₂ ² H ₄ H ₉ N ₃ O ₄ S		
DRE-C10024051	Acetylsulfamethoxazole D4		10mg	
Alachlor D13 (2,6-diethylphenyl D13)				
CAS 1015856-63-9	MW 282.8472	C ₁₄ ² H ₁₃ H ₇ ClNO ₂		
DRE-C10060100	Alachlor D13 (2,6-diethylphenyl D13)		10mg	
DRE-XA10060100AC	Alachlor D13 100 µg/mL in Acetone(±)		1.1ml	
Alachlor D3 (methoxy D3)				
CAS n/a	MW 272.7856	C ₁₄ ² H ₉ H ₁₇ ClNO ₂		
DRE-XA10060001AC	Alachlor D3 (methoxy D3) 100 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Aldicarb D3 (N-methyl D3)				
CAS 1795142-83-4 DRE-C10070100	MW 193.2817 Aldicarb D3 (methyl-D3)	C ₇ ² H ₉ H ₁₁ N ₂ O ₂ S	10mg	
Aldicarb-sulfone D3 (N-methyl D3)				
CAS 1795135-15-7 DRE-XA10080100AC	MW 225.2805 Aldicarb-sulfone D3 100 µg/mL in Acetone (*)	C ₇ ² H ₉ H ₉ N ₂ O ₄ S	1.1ml	
Amidosulfuron D6 (dimethoxy D6)				
CAS n/a DRE-C10162100	MW 375.4117 Amidosulfuron D6	C ₉ ² H ₆ H ₉ N ₅ O ₇ S ₂	10mg	
4-Aminobiphenyl D9				
CAS 344298-96-0 DRE-XA10173041AC	MW 178.2779 4-Aminobiphenyl D9 100 µg/mL in Acetone(‡)	C ₁₂ ² H ₉ H ₉ N	1ml	
1-Aminohydantoin (2,4,5-13C3)				
CAS 957509-31-8 DRE-XA10203190AL	MW 118.0687 1-Aminohydantoin (2,4,5-13C3) 100 µg/mL in Acetonitrile	¹³ C ₃ H ₅ N ₃ O ₂	1.1ml	
(Aminomethyl) phosphonic Acid 13C 15N (AMPA)				
CAS n/a DRE-XA10205100WA	MW 113.0231 Aminomethylphosphonic acid 13C 15N (AMPA) 100 µg/mL in Water(‡)	¹³ CH ₆ ¹⁵ NO ₃ P	1ml	
3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5)				
CAS 1017793-94-0 DRE-C10206310	MW 206.2538 3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5)(‡)	C ₈ ² H ₉ H ₁₀ N ₃ O ₃	10mg	
DRE-XA10206310AL	3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5) 100 µg/mL in Acetonitrile		1ml	
1-Aminonaphthalene D7				
CAS 78832-53-8 DRE-XA10206351ME	MW 150.2283 1-Aminonaphthalene D7 100 µg/mL in Methanol	C ₁₀ ² H ₇ H ₇ N	1ml	
2-Aminonaphthalene D7				
CAS 93951-94-1 DRE-XA10206356ME	MW 150.2283 2-Aminonaphthalene D7 100 µg/mL in Methanol	C ₁₀ ² H ₇ H ₇ N	1ml	










Stable isotope labelled compounds

Product code	Description			
3-Amino-2-oxazolidinone D4 (AOZ D4)				
CAS 1188331-23-8 DRE-C10209010	MW 106.1166 3-Amino-2-oxazolidinone (AOZ) D4(±)	C ₃ ² H ₄ H ₂ N ₂ O ₂	10mg	
Amitraz D12 (methylphenyl D12)				
CAS n/a DRE-XA10230100AC	MW 305.48 Amitraz D12 100 µg/mL in Acetone(±)	C ₁₉ ² H ₁₂ H ₁₁ N ₃	1.1ml	
Amitrole (1-15N)				
CAS 367498-28-0 DRE-XA10240100ME	MW 85.0734 Amitrole (1-15N) 100 µg/mL in Methanol	C ₂ H ₄ ¹⁵ NN ₃	1.1ml	
Amitrole 1-15N 5-13C				
CAS n/a DRE-XA10240110AL	MW 86.066 Amitrole 1-15N 5-13C 100µg/mL in Acetonitrile (*) (±)	¹³ CCH ₄ ¹⁵ NN ₃	1.1ml	
Aniline D5				
CAS 4165-61-1 DRE-C10262600 DRE-YA10262600MB	MW 98.1573 Aniline D5 Aniline D5 2000 µg/mL in Methyl-tert. butyl ether	C ₆ ² H ₅ H ₂ N	100mg 1.1ml	
Anthracene D10				
CAS 1719-06-8 DRE-C20520100 DRE-L20520100CY DRE-XA20520100CY DRE-YA20520100MB	MW 188.2908 Anthracene D10 Anthracene D10 10 µg/mL in Cyclohexane Anthracene D10 100 µg/mL in Cyclohexane Anthracene D10 2000 µg/mL in Methyl-tert. butyl ether	C ₁₄ ² H ₁₀	100mg 10ml 1ml 1ml	
Atrazine 13C3 (ring 13C3)				
CAS 1443685-80-0 DRE-XA10330200AC	MW 218.6612 Atrazine 13C3 (ring 13C3) 100 µg/mL in Acetone	¹³ C ₃ C ₅ H ₄ ClN ₅	1.1ml	
Atrazine D5 (ethyl-D5)				
CAS 163165-75-1 DRE-C10330100 DRE-XA10330100AC DRE-YA10330100AL	MW 220.7141 Atrazine D5(±) Atrazine D5 100 µg/mL in Acetone(±) Atrazine D5 (ethylamino D5) 1000 µg/mL in Acetonitrile	C ₆ ² H ₅ H ₄ ClN ₅	10mg 1ml 1.1ml	
Atrazine-desethyl D6 (dimethyl D6)				
CAS 1216649-31-8 DRE-XA10331100AC	MW 193.6671 Atrazine-desethyl D6 100 µg/mL in Acetone	C ₆ ² H ₆ H ₄ ClN ₅	1.1ml	










Stable isotope labelled compounds

Product code	Description			
Atrazine-desisopropyl D5 (ethylamino D5)				
CAS 1189961-78-1	MW 178.6343	C ₅ ² H ₅ H ₃ CIN ₃		
DRE-C10332100	Atrazine-desisopropyl D5 (ethylamino D5)		10mg	
DRE-XA10332100AC	Atrazine-desisopropyl D5 (ethylamino D5) 100 µg/mL in Acetone		1.1ml	
Atrazine-2-hydroxy D5 (ethyl D5)				
CAS 1276197-25-1	MW 202.2684	C ₅ ² H ₅ H ₁₀ N ₃ O		
DRE-XA10333100ME	Atrazine-2-hydroxy D5 100 µg/mL in Methanol		1.1ml	
Azaperone D4				
CAS 1173021-72-1	MW 331.4205	C ₁₉ ² H ₁₄ H ₁₈ FN ₃ O		
DRE-C10340512	Azaperone D4		10mg	
Azinphos-ethyl D10 (ethyl D10)				
CAS n/a	MW 355.4391	C ₁₂ ² H ₁₀ H ₆ N ₃ O ₃ PS ₂		
DRE-XA10360100AC	Azinphos-ethyl D10 100 µg/mL in Acetone		1.1ml	
Azinphos-methyl D6 (dimethyl D6)				
CAS n/a	MW 323.3613	C ₁₀ ² H ₆ H ₆ N ₃ O ₃ PS ₂		
DRE-C10365100	Azinphos-methyl D6 (dimethyl D6)		10mg	
DRE-XA10365100AC	Azinphos-methyl D6 100 µg/mL in Acetone		1.1ml	
Basic Violet 3 D6				
CAS 1266676-01-0	MW 414.0158	C ₂₅ ² H ₆ H ₂₄ N ₃ ⁺ Cl ⁻		
DRE-C10427505	Basic Violet 3 D6		10mg	
Bentazone (isopropyl-1,1,1,3,3,3) D6				
CAS n/a	MW 246.3159	C ₁₀ ² H ₆ H ₆ N ₂ O ₃ S		
DRE-C10510100	Bentazone (isopropyl-1,1,1,3,3,3) D6		10mg	
DRE-XA10510100AC	Bentazone (isopropyl-1,1,1,3,3,3) D6 100 µg/mL in Acetone (*)		1.1ml	
Benzene D6				
CAS 1076-43-3	MW 84.1488	C ₆ ² H ₆		
DRE-C10535200	Benzene D6		2ml	
DRE-YA10535200ME	Benzene D6 2000 µg/mL in Methanol		1ml	
Benzo[a]anthracene D12				
CAS 1718-53-2	MW 240.3618	C ₁₈ ² H ₁₂		
DRE-C20545100	Benzo(a)anthracene D12		100mg	
DRE-L20545100AC	Benzo(a)anthracene D12 10 µg/mL in Acetone		10ml	
DRE-L20545100CY	Benzo(a)anthracene D12 10 µg/mL in Cyclohexane(‡)		10ml	









Stable isotope labelled compounds

Product code	Description			
Benzo[a]pyrene D12				
CAS 63466-71-7	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20635100	Benzo(a)pyrene D12(±)		10mg	
DRE-L20635100CY	Benzo(a)pyrene D12 10 µg/mL in Cyclohexane(±)		10ml	
DRE-LA20635100AL	Benzo(a)pyrene D12 10 µg/mL in Acetonitrile		1ml	
DRE-LA20635100CY	Benzo(a)pyrene D12 10 µg/mL in Cyclohexane		1ml	
DRE-XA20635100CY	Benzo(a)pyrene D12 100 µg/mL in Cyclohexane		1ml	
Benzo[b]fluoranthene D12				
CAS 93951-98-5	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20565100	Benzo(b)fluoranthene D12		10mg	
DRE-LA20565100CY	Benzo(b)fluoranthene D12 10 µg/mL in Cyclohexane		1ml	
Benzo(e)pyrene D12				
CAS 205440-82-0	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-XA20645010CY	Benzo(e)pyrene D12 100 µg/mL in Cyclohexane		1ml	
Benzo[g,h,i]perylene D12				
CAS 93951-66-7	MW 288.4046	C ₂₂ ² H ₁₂		
DRE-C20630200	Benzo(g,h,i)perylene D12		10mg	
DRE-LA20630200CY	Benzo(g,h,i)perylene D12 10 µg/mL in Cyclohexane		1ml	
Benzo[k]fluoranthene D12				
CAS 93952-01-3	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20580200	Benzo(k)fluoranthene D12		10mg	
DRE-LA20580200CY	Benzo(k)fluoranthene D12 10 µg/mL in Cyclohexane		1ml	
Bezafibrate D6 (dimethyl D6)				
CAS 1219802-74-0	MW 367.8564	C ₁₉ ² H ₆ H ₁₄ ClNO ₄		
DRE-C10578010	Bezafibrate D6 (dimethyl D6)		10mg	
Biphenyl D10				
CAS 1486-01-7	MW 164.2694	C ₁₂ ² H ₁₀		
DRE-C10630010	Biphenyl D10		100mg	
DRE-LA10630010AC	Biphenyl D10 10 µg/mL in Acetone		1ml	
N,N'-Bis-(4-nitrophenyl)urea D8				
CAS 1156508-87-0	MW 310.2916	C ₁₃ ² H ₆ H ₂ N ₄ O ₅		
DRE-C15598600	N,N'-Bis-(4-nitrophenyl)urea D8		10mg	
Bisphenol A D14				
CAS 120155-79-5	MW 242.3726	C ₁₅ ² H ₁₄ H ₂ O ₂		
DRE-XA10655503AL	Bisphenol A D14 100 µg/mL in Acetonitrile		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Bisphenol A D16				
CAS 96210-87-6 DRE-C10655501	MW 244.3849 Bisphenol A D16	$C_{15}^2H_{16}O_2$	50mg	
Bromacil D3 (methyl D3)				
CAS n/a DRE-XA10670100AL	MW 264.1342 Bromacil D3 (methyl D3) 100 µg/mL in Acetonitrile	$C_9^2H_3H_{10}BrN_2O_2$	1.1ml	
Bromopropylate D7 (isopropyl D7)				
CAS n/a DRE-XA10762100AC	MW 435.1583 Bromopropylate D7 (isopropyl D7) 100 µg/mL in Acetone	$C_{17}^2H_7H_8Br_2O_3$	1.1ml	
Butylhydroxytoluene-d21 (BHT-d21; 2,6-Bis[1,1-di(methyl-d3)ethyl-2,2,2-d3]-4-methylphen-3,5-d2-ol-d)				
CAS 64502-99-4 DRE-GS09010395ME	MW 241.4799 d21-2,6-di(tert-butyl)-4-methylphenol (BHT-d21), 1000 µg/mL in Methanol(±)	$C_{15}^2H_{21}H_3O$	5x1ml	
Cadusafos D5 (ethyl-D5)				
CAS 1215380-71-4 DRE-XA10934010AL	MW 275.423 Cadusafos D5 (ethyl-D5) 100 µg/mL in Acetonitrile	$C_{16}^2H_5H_{18}O_2PS_2$	1.1ml	
Captan-4,4,5,6,7,7-D6				
CAS 1330190-00-5 DRE-XA10960100AC	MW 306.6263 Captan D6 100 µg/mL in Acetone (*)	$C_9^2H_6H_2Cl_3NO_2S$	1.1ml	
Carazolol D7				
CAS 1173021-02-7 DRE-C10968010	MW 305.4226 Carazolol D7	$C_{18}^2H_7H_{19}N_2O_2$	10mg	
Carbaryl D3 (methyl D3)				
CAS 1433961-56-8 DRE-XA10980010AC	MW 204.2397 Carbaryl D3 (methyl D3) 100 µg/mL in Acetone	$C_{12}^2H_3H_8NO_2$	1.1ml	
Carbaryl D7 (naphthyl D7)				
CAS 362049-56-7 DRE-C10980100 DRE-XA10980100AC	MW 208.2644 Carbaryl D7 (naphthyl D7) Carbaryl D7 (naphthyl D7) 100 µg/mL in Acetone	$C_{12}^2H_7H_4NO_2$	50mg 1.1ml	

Stable isotope labelled compounds

Product code	Description			
Carbazole D8				
CAS 38537-24-5	MW 175.2559	C ₁₂ ² H ₈ HN		
DRE-XA10985100AC	Carbazole D8 100 µg/mL in Acetone		1.1ml	
Carbendazim D3 (methyl D3)				
CAS 1255507-88-0	MW 194.2051	C ₉ ² H ₃ H ₆ N ₃ O ₂		
DRE-C10990100	Carbendazim D3		10mg	
DRE-LA10990100AL	Carbendazim D3 (methyl D3) 10µg/mL in Acetonitrile (*)		1ml	
Carbendazim D4 (ring D4)				
CAS 291765-95-2	MW 195.2113	C ₉ ² H ₄ H ₆ N ₃ O ₂		
DRE-C10990200	Carbendazim D4 (ring)		10mg	
DRE-XA10990200AC	Carbendazim D4 (ring) 100 µg/mL in Acetone (*)		1.1ml	
Carbetamide D5 (phenyl D5)				
CAS n/a	MW 241.2978	C ₁₂ ² H ₅ H ₁₁ N ₂ O ₃		
DRE-XA11000100AC	Carbetamide D5 (phenyl D5) 100 µg/mL in Acetone(‡)		1.1ml	
Carbofuran D3 (N-methyl D3)				
CAS 1007459-98-4	MW 224.2709	C ₁₂ ² H ₃ H ₇ N ₂ NO ₃		
DRE-C11010100	Carbofuran D3 (N-methyl D3)(‡)		10mg	
DRE-XA11010100AC	Carbofuran D3 (N-methyl D3) 100 µg/mL in Acetone		1.1ml	
Carbophenothion D10 (di(ethyl D5))				
CAS n/a	MW 352.9269	C ₁₁ ² H ₁₀ H ₆ ClO ₂ PS ₃		
DRE-XA11020100AC	Carbophenothion D10 (di(ethyl D5)) 100 µg/mL in Acetone(‡)		1.1ml	
Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2R)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)				
CAS 202480-68-0	MW 328.1602	C ₁₁ ² H ₅ H ₇ Cl ₂ N ₂ O ₅		
DRE-C11120100	Chloramphenicol D5 (ring D4, benzyl D)		10mg	
DRE-XA11120100AL	Chloramphenicol D5 100 µg/mL in Acetonitrile(‡)		1ml	
erythro-Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2S)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)				
CAS 1426174-26-6	MW 328.1602	C ₁₁ ² H ₅ H ₇ Cl ₂ N ₂ O ₅		
DRE-XA1120110AL	(Erythro)-Chloramphenicol D5 (ring D4, benzyl D) 100 µg/mL in Acetonitrile		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Chlorfenvinphos D10 (ethyl D10)				
CAS 1346606-54-9 DRE-XA11290100AC	MW 369.6315 Chlorfenvinphos (ethyl) D10 100 µg/mL in Acetone	$C_{12}^2H_{10}H_4Cl_3O_4P$	1.1ml	
Chloridazon D5 (phenyl D5)				
CAS 1346818-99-4 DRE-XA11320100AL	MW 226.6738 Chloridazon D5 (phenyl D5) 100 µg/mL in Acetonitrile	$C_{10}^2H_5H_3ClN_3O$	1.1ml	
Chlormequat-chloride 1,1,2,2-D4				
CAS n/a DRE-CA11340100 DRE-X11340100DO DRE-XA11340100DO	MW 162.0941 Chlormequat chloride 1,1,2,2-D4(±) Chlormequat-chloride D4 100 µg/mL in Deuterium oxide Chlormequat-chloride D4 (CCC) 100 µg/mL in Deuterium oxide	$C_5^2H_4H_3ClN.Cl$	5mg 10ml 1.1ml	
Chlorobenzene D5				
CAS 3114-55-4 DRE-C11380100	MW 117.5877 Chlorobenzene D5	$C_6^2H_5Cl$	100mg	
2'-Chloro-2,3,4,5,6-pentadeuterio-1,1'-biphenyl				
CAS 51624-35-2 DRE-XA20000101IO	MW 193.6837 2-Chlorobiphenyl D5 100 µg/mL in Isooctane	$C_{12}^2H_5H_4Cl$	1.1ml	
2-Chlorophenol-3,4,5,6-D4				
CAS 93951-73-6 DRE-C11470100 DRE-XA11470100AC	MW 132.5809 2-Chlorophenol-3,4,5,6-D4 2-Chlorophenol D4 100 µg/mL in Acetone	$C_6^2H_4HClO$	10mg 1.1ml	
Chlorotoluron D6 (N,N-dimethyl D6)				
CAS 1219803-48-1 DRE-C11530100 DRE-XA11530100AC	MW 218.713 Chlorotoluron (N,N-dimethyl) D6 Chlorotoluron D6 100 µg/mL in Acetone	$C_{10}^2H_6H_7ClN_2O$	5mg 1.1ml	
Chlorpyrifos (diethyl-D10)				
CAS 285138-81-0 DRE-C11600100 DRE-XA11600100AC	MW 360.6479 Chlorpyrifos D10(±) Chlorpyrifos D10 100 µg/mL in Acetone(±)	$C_9^2H_{10}HCl_3NO_3PS$	25mg 1ml	
Chlorpyrifos-methyl D6				
CAS n/a DRE-C11601100 DRE-XA11601100AC	MW 328.5701 Chlorpyrifos-methyl D6 Chlorpyrifos-methyl D6 100 µg/mL in Acetone (*)	$C_7^2H_6HCl_3NO_3PS$	10mg 1.1ml	










Stable isotope labelled compounds

Product code	Description			
Chlorthiophos-I-D10 (diethyl D10)				
CAS n/a DRE-C11650010	MW 371.3064 Chlorthiophos D10 (diethyl D10)	$C_{11}^2H_{10}H_5Cl_2O_3PS_2$	10mg	
Chrysene D12				
CAS 1719-03-5 DRE-C20670100 DRE-L20670100AL DRE-L20670100CY	MW 240.3618 Chrysene D12 Chrysene D12 10 µg/mL in Acetonitrile Chrysene D12 10 µg/mL in Cyclohexane	$C_{18}^2H_{12}$	100mg 10ml 10ml	
Cimaterol D7 (isopropyl D7)				
CAS 1228182-44-2 DRE-C11666352	MW 226.326 Cimaterol D7	$C_{12}^2H_7H_{10}N_3O$	10mg	
Ciprofloxacin D8 Hydrochloride				
CAS 1216659-54-9 DRE-C11668501 DRE-XA11668501WA	MW 375.8518 Ciprofloxacin D8 hydrochloride Ciprofloxacin D8 hydrochloride 100 µg/mL in Water	$C_{17}^2H_8H_{10}FN_3O_3 \cdot ClH$	10mg 1.1ml	
(±)-Clenbuterol D9 (trimethyl D9)				
CAS 129138-58-5 DRE-C11668561 DRE-XA11668561AL	MW 286.2456 (±)-Clenbuterol D9 (trimethyl D9)(‡) (±)-Clenbuterol D9 (trimethyl D9) 100 µg/mL in Acetonitrile (*)	$C_{12}^2H_8H_9Cl_2N_2O$	25mg 1ml	
Clenproperole D7				
CAS 1173021-09-4 DRE-C11668742	MW 270.2067 Clenproperole D7	$C_{11}^2H_7H_9Cl_2N_2O$	10mg	
Clothianidin D3 (N'-methyl D3)				
CAS 1262776-24-8 DRE-C11691710	MW 252.6965 Clothianidin D3 (N'-methyl D3)	$C_8^2H_3H_5ClN_5O_2S$	50mg	
Coumaphos D10 (di(ethyl-D5))				
CAS 287397-86-8 DRE-XA11730010AL	MW 372.8272 Coumaphos D10 di(ethyl-D5) 100 µg/mL in Acetonitrile	$C_{14}^2H_{10}H_6ClO_3PS$	1.1ml	
Coumarin 5,6,7,8-D4				
CAS 185056-83-1 DRE-XA11735010AC	MW 150.1674 Coumarin 5,6,7,8-D4 100 µg/mL in Acetone	$C_9^2H_4H_2O_2$	1.1ml	










Stable isotope labelled compounds

Product code	Description			
Cyanazine D5 (N-ethyl D5)				
CAS 1190003-29-2	MW 245.7235	C ₉ ² H ₅ H ₆ ClN ₆		
DRE-XA11790100AC	Cyanazine D5 (N-ethyl D5) 100 µg/mL in Acetone		1.1ml	
trans-Cyfluthrin D6 (2,2-dimethyl D6)				
CAS n/a	MW 440.3246	C ₂₂ ² H ₆ H ₁₂ Cl ₂ FNO ₃		
DRE-XA11850010AL	trans-Cyfluthrin D6 (2,2-dimethyl D6) 100 µg/mL in Acetonitrile		1.1ml	
trans-Cypermethrin D6 (dimethyl D6)				
CAS 82523-65-7	MW 422.3341	C ₂₂ ² H ₆ H ₁₂ Cl ₂ NO ₃		
DRE-XA11890400AC	trans-Cypermethrin D6 (dimethyl D6) 100 µg/mL in Acetone(‡)		1ml	
Cyromazine D4				
CAS 1219804-19-9	MW 170.2084	C ₆ ² H ₄ H ₆ N ₆		
DRE-C11920010	Cyromazine D4 (cyclopropyl-2,2,3,3 D4)		10mg	
DRE-XA11920010EA	Cyromazine D4 (cyclopropyl-2,2,3,3 D4) 100 µg/mL in Ethylacetate (*)		1.1ml	
2,4-D (phenyl-13C6)				
CAS 150907-52-1	MW 226.9934	¹³ C ₆ C ₂ H ₆ Cl ₂ O ₃		
DRE-XA11940200AC	2,4-D 13C6 100 µg/mL in Acetone(‡)		1ml	
2,4-D D3 ((2,4-Dichloro-3,5,6-trideuteriophenoxy)acetic Acid)				
CAS 202480-67-9	MW 224.0559	C ₈ ² H ₃ H ₃ Cl ₂ O ₃		
DRE-C11940100	2,4-D D3		10mg	
DRE-XA11940100AC	2,4-D D3 100 µg/mL in Acetone(‡)		1.1ml	
2,4'-DDE D8 (1,1-Dichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethene D8)				
CAS 1402834-57-4	MW 326.0746	C ₁₄ ² H ₆ Cl ₄		
DRE-XA12040100AC	2,4'-DDE D8 100 µg/mL in Acetone		1.1ml	
4,4'-DDE D8 (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethene D8)				
CAS 93952-19-3	MW 326.0746	C ₁₄ ² H ₆ Cl ₄		
DRE-XA12041100AC	4,4'-DDE D8 100 µg/mL in Acetone(‡)		1ml	
4,4'-DDT (ring-13C12)				
CAS 104215-84-1	MW 366.3981	¹³ C ₁₂ C ₂ H ₆ Cl ₅		
DRE-XA12082200AC	4,4'-DDT 13C12 100 µg/mL in Acetone(‡)		1.1ml	










Stable isotope labelled compounds

Product code	Description			
2,4'-DDT 13C12 (benzen 13C12)				
CAS 1396995-26-8 DRE-XA12081200AC	MW 366.3981 2,4'-DDT 13C12 100 µg/mL in Acetone	¹³ C ₁₂ H ₉ Cl ₅	1.1ml	
2,4'-DDT D8 (1,1,1-Trichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl)ethane D8 (ring D8))				
CAS 221899-88-3 DRE-C12081100 DRE-XA12081100AC	MW 362.5356 2,4'-DDT D8 2,4'-DDT D8 100 µg/mL in Acetone	C ₁₄ H ₆ HCl ₅	5mg 1.1ml	
4,4'-DDT D8 (1,1,1-Trichloro-2,2-bis(p-chlorophenyl)ethane D8 (ring D8))				
CAS 93952-18-2 DRE-C12082100 DRE-XA12082100AC DRE-XA12082100CY	MW 362.5356 4,4'-DDT D8(‡) 4,4'-DDT D8 100 µg/mL in Acetone 4,4'-DDT D8 100 µg/mL n Cyclohexane	C ₁₄ H ₆ HCl ₅	10mg 1.1ml 1.1ml	
DEET D7 (methyl D3 phenyl D4)				
CAS 1219799-37-7 DRE-XA12100010ME	MW 198.3126 DEET D7 (methyl D3 benzeneamide D4) 100 µg/mL in Methanol	C ₁₂ H ₇ H ₁₀ NO	1.1ml	
DNOC D5 (ring D2, methyl D3)				
CAS n/a DRE-XA13050100AC	MW 203.1637 DNOC D5 (ring D2, methyl D3) 100 µg/mL in Acetone	C ₇ H ₅ HN ₂ O ₅	1.1ml	
Daminozide D6 (dimethyl D6)				
CAS n/a DRE-C11960100 DRE-XA11960100AL	MW 166.2081 Daminozide D6 (dimethyl D6) Daminozide D6 100 µg/mL in Acetonitrile	C ₆ H ₆ H ₆ N ₂ O ₃	10mg 1.1ml	
Decoquinat D5 (7-ethoxy D5)				
CAS 1228182-55-5 DRE-C12097010	MW 422.5692 Decoquinat D5	C ₂₄ H ₂₅ H ₃₀ NO ₅	10mg	
trans-Deltamethrin D6 (dimethyl D6)				
CAS n/a DRE-C12120100 DRE-XA12120100AL	MW 511.2361 trans-Deltamethrin D6 (dimethyl D6) trans-Deltamethrin D6 (dimethyl D6) 100 µg/mL in Acetonitrile(‡)	C ₂₂ H ₁₆ H ₁₃ Br ₂ NO ₃	10mg 1.1ml	
Demeton-S-methyl D6 (dimethyl D6)				
CAS n/a DRE-XA12143100CY	MW 236.3222 Demeton-S-methyl D6 (dimethyl D6) 100 µg/mL in Cyclohexane (*)	C ₆ H ₆ H ₉ O ₃ PS ₂	1.1ml	










Stable isotope labelled compounds

Product code	Description			
Demeton-S-methyl-sulfoxide D6 (di-0-methyl D6)				
CAS n/a	MW 252.3216	C ₆ ² H ₆ H ₉ O ₄ PS ₂		
DRE-XA12145010AC	Demeton-S-methyl-sulfoxide D6 (di-0-methyl D6) 100 µg/mL in Acetonitrile		1.1ml	
Diazinon D10 (diethyl D10)				
CAS 100155-47-3	MW 314.4071	C ₁₂ ² H ₁₀ H ₁₁ N ₂ O ₃ PS		
DRE-C12210100	Diazinon D10 (diethyl D10)(‡)		10mg	
DRE-XA12210100AC	Diazinon D10 (diethyl D10) 100 µg/mL in Acetone		1.1ml	
Dibenzo[a,h]anthracene D14				
CAS 13250-98-1	MW 292.4328	C ₂₂ ² H ₁₄		
DRE-C20700200	Dibenzo(a,h)anthracene D14		10mg	
DRE-L20700200CY	Dibenzo(a,h)anthracene D14 10 µg/mL in Cyclohexane(‡)		10ml	
Dibenzofuran D8				
CAS 93952-04-6	MW 176.2406	C ₁₂ ² H ₆ O		
DRE-L20710100CY	Dibenzofuran D8 10 µg/mL in Cyclohexane		10ml	
Dibenzothiophene D8				
CAS 33262-29-2	MW 192.3062	C ₁₂ ² H ₆ S		
DRE-L20727100CY	Dibenzothiophene D8 10 µg/mL in Cyclohexane		10ml	
Dicamba 13C6 (ring 13C6)				
CAS 1173023-06-7	MW 226.9934	¹³ C ₆ C ₂ H ₆ Cl ₂ O ₃		
DRE-XA12260005AL	Dicamba 13C6 100 µg/mL in Acetonitrile		1.1ml	
Dicamba D3 (methoxy D3)				
CAS 349553-95-3	MW 224.0559	C ₈ ² H ₃ H ₃ Cl ₂ O ₃		
DRE-C12260100	Dicamba D3 (methoxy D3)		10mg	
DRE-XA12260100AC	Dicamba D3 (methoxy D3) 100 µg/mL in Acetone		1.1ml	
Dicamba D4 (phenyl D1 methoxy D3)				
CAS n/a	MW 225.0621	C ₈ ² H ₄ H ₂ Cl ₂ O ₃		
DRE-XA12260110AC	Dicamba D4 (phenyl D1 methoxy D3) 100 µg/mL in Acetone		1.1ml	
1,2-Dichlorobenzene D4				
CAS 2199-69-1	MW 151.0266	C ₆ ² H ₄ Cl ₂		
DRE-CA12370100	1,2-Dichlorobenzene D4		100mg	
DRE-YA12370100ME	1,2-Dichlorobenzene D4 2000 µg/mL in Methanol		1ml	










Stable isotope labelled compounds

Product code	Description			
1,4-Dichlorobenzene D4				
CAS 3855-82-1	MW 151.0266	C ₆ ² H ₄ Cl ₂		
DRE-YA12372100ME	1,4-Dichlorobenzene D4 2000 µg/mL in Methanol		1ml	
DRE-C12372100	1,4-Dichlorobenzene D4		100mg	
3,3'-Dichlorobenzidine D6 (ring D6)				
CAS 93951-91-8	MW 259.1642	C ₁₂ ² H ₆ H ₄ Cl ₂ N ₂		
DRE-XA12377910AL	3,3'-Dichlorobenzidine D6 (ring D6) 100 µg/mL in Acetonitrile		1.1ml	
4,4'-Dichlorobenzophenone D8				
CAS 1219806-01-5	MW 259.1573	C ₁₃ ² H ₆ Cl ₂ O		
DRE-XA12410100AC	4,4'-Dichlorobenzophenone D8 100 µg/mL in Acetone		1.1ml	
1,2-Dichloroethane D4				
CAS 17060-07-0	MW 102.9838	C ₂ ² H ₄ Cl ₂		
DRE-C12422300	1,2-Dichloroethane D4		100mg	
DRE-YA12422300ME	1,2-Dichloroethane D4 2000 µg/mL in Methanol		1ml	
2,4-Dichlorophenol 13C6				
CAS 1202864-83-2	MW 168.9573	¹³ C ₆ H ₄ Cl ₂ O		
DRE-XA12451200AC	2,4-Dichlorophenol 13C6 100 µg/mL in Acetone		1.1ml	
2,4-Dichlorophenol D3				
CAS 93951-74-7	MW 166.0198	C ₆ ² H ₃ HCl ₂ O		
DRE-C12451100	2,4-Dichlorophenol-3,5,6 D3		50mg	
DRE-XA12451100MB	2,4-Dichlorophenol-3,5,6 D3 100 µg/mL in Methyl-tert. butyl ether		1ml	
(E)-3-(2,2-Dichlorovinyl)-2,2-di(methyl D3)-(1-cyclopropane)carboxylic acid D6				
CAS n/a	MW 215.1068	C ₈ ² H ₆ H ₄ Cl ₂ O ₂		
DRE-XA12507510AC	(E)-3-(2,2-Dichlorovinyl)-2,2-di(methyl D3)-(1-cyclopropane)carboxylic acid D6 100 µg/mL in Acetone		1ml	
Dichlorprop D6 (ring D3, 3,3,3-D3)				
CAS n/a	MW 241.101	C ₉ ² H ₆ H ₂ Cl ₂ O ₃		
DRE-C12510100	Dichlorprop D6 (ring D3, 3,3,3-D3)		10mg	
DRE-XA12510100AC	Dichlorprop (ring-D3, 3,3,3-D3) 100 µg/mL in Acetone		1.1ml	
Dichlorvos D6 (dimethyl D6)				
CAS 203645-53-8	MW 227.0127	C ₈ ² H ₆ HCl ₂ O ₄ P		
DRE-C12530100	Dichlorvos D6 (dimethyl D6)		10mg	
DRE-XA12530100CY	Dichlorvos D6 100 µg/mL in Cyclohexane (*)		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Diclofenac D4 (phenyl D4)				
CAS 153466-65-0 DRE-XA12537010AC	MW 300.1733 Diclofenac D4 acid (phenyl D4) 100 µg/mL in Acetonitrile (*)	C ₁₄ ² H ₄ H ₇ Cl ₂ NO ₂	1.1ml	
Dicofol D8 (ring D8)				
CAS n/a DRE-XA12570100CY	MW 378.535 Dicofol D8 (ring D8) 100 µg/mL in Cyclohexane(±)	C ₁₄ ² H ₆ HCl ₅ O	1ml	
Difenoconazole D6 (1,1,2,3,3,3-propyl-D6)				
CAS n/a DRE-XA12609010AL	MW 412.2996 Difenoconazole D6 (1,1,2,3,3,3-propyl-D6) 100 µg/mL in Acetonitrile	C ₁₉ ² H ₆ H ₁₁ Cl ₂ N ₃ O ₃	1.1ml	
Diflufenican D3 (3-trifluoromethylphenoxy-2,4,6 D3)				
CAS 1185009-29-3 DRE-XA12631001AL	MW 397.3133 Diflufenican D3 (3-trifluoromethylphenoxy-2,4,6 D3) 100 µg/mL in Acetonitrile	C ₁₉ ² H ₅ H ₈ F ₉ N ₂ O ₂	1.1ml	
Dimethoate D6 (O,O dimethyl D6)				
CAS 1219794-81-6 DRE-C12700100 DRE-XA12700100AC	MW 235.2944 Dimethoate D6 (O,O-dimethyl-D6) Dimethoate D6 (O,O dimethyl D6) 100 µg/mL in Acetone	C ₅ ² H ₆ H ₆ NO ₃ PS ₂	10mg 1.1ml	
O,O-Dimethylphosphorothioic acid sodium salt D6				
CAS n/a DRE-XA12738610AC	MW 170.1327 O,O-Dimethylphosphorothioic acid sodium salt D6 100 µg/mL in Acetone (*)	C ₂ ² H ₆ O ₃ PS·Na	1.1ml	
2,4-Dimethyl-3,5,6-trideuteriophenol				
CAS 93951-75-8 DRE-C12731100 DRE-XA12731100AC	MW 125.1829 2,4-Dimethylphenol-3,5,6-D3 2,4-Dimethylphenol-3,5,6 D3 100 µg/mL in Acetone	C ₈ ² H ₃ H ₇ O	100mg 1.1ml	
Dimetridazole D3 (N-methyl D3)				
CAS 64678-69-9 DRE-C12772010 DRE-XA12772010AC	MW 144.1465 Dimetridazole D3 Dimetridazole D3 100 µg/mL in Acetone	C ₅ ² H ₃ H ₄ N ₃ O ₂	10mg 1.1ml	
Dimetridazole-2-hydroxy D3				
CAS 1015855-78-3 DRE-C12772051	MW 160.1459 Dimetridazole-2-hydroxy D3	C ₅ ² H ₃ H ₄ N ₃ O ₃	10mg	










Stable isotope labelled compounds

Product code	Description			
1,4-Dioxane-d8 (Octadeuterodioxane)				
CAS 17647-74-4 DRE-GA09010386ME	MW 96.1544	C ₆ ² H ₈ O ₂	1ml	
Diquat dibromide D4				
CAS n/a DRE-CA12960010	MW 348.0697	C ₁₂ ² H ₄ H ₈ N ₂ ·2Br	50mg	
Disperse Yellow 3 D3 (acetyl D3)				
CAS 947601-96-9 DRE-XA12972311AL	MW 272.317	C ₁₅ ² H ₃ H ₁₂ N ₃ O ₂	1.1ml	
Disulfoton D10 (Di-ethyl D10)				
CAS n/a DRE-XA12980100AC	MW 284.4656	C ₈ ² H ₁₀ H ₉ O ₂ PS ₃	1.1ml	
DRE-XA12980100CY	Disulfoton D10 (Di-ethyl D10) 100 µg/mL in Cyclohexane		1.1ml	
Diuron D6				
CAS 1007536-67-5 DRE-C13020100	MW 239.1315	C ₈ ² H ₆ H ₄ Cl ₂ N ₂ O	10mg	
DRE-XA13020100AC	Diuron D6 100 µg/mL in Acetone		1.1ml	
α-Endosulfan D4				
CAS 203645-57-2 DRE-C13121100	MW 410.9498	C ₉ ² H ₄ H ₂ Cl ₆ O ₃ S	10mg	
DRE-XA13121100AC	α-Endosulfan D4 100 µg/mL in Acetone		1.1ml	
β-Endosulfan D4				
CAS 203716-99-8 DRE-C13122100	MW 410.9498	C ₉ ² H ₄ H ₂ Cl ₆ O ₃ S	10mg	
DRE-XA13122100AC	β-Endosulfan D4 100 µg/mL in Acetone		1.1ml	
Enilconazole-D5 (2-propenyl-D5) (Imazalil-D5)				
CAS 1398065-91-2 DRE-XA14280100AC	MW 302.2106	C ₁₄ ² H ₅ H ₉ Cl ₂ N ₂ O	1.1ml	
Enrofloxacin D5 Hydrochloride (ethyl d5)				
CAS n/a DRE-C13170100	MW 400.8864	C ₁₉ ² H ₅ H ₁₇ FN ₃ O ₃ ·ClH	10mg	










Stable isotope labelled compounds

Product code	Description			
Ethephon D4 (2-Chloroethyl-1,1,2,2-D4)				
CAS 1020719-29-2	MW 148.5186	C ₂ ² H ₄ H ₂ ClO ₃ P		
DRE-C13230100	Ethephon D4 (2-Chloroethyl-1,1,2,2-D4)		10mg	
DRE-XA13230100AC	Ethephon D4 (2-Chloroethyl-1,1,2,2-D4) 100 µg/mL in Acetone (*)		1ml	
Ethion D20 (tetraethyl D20)				
CAS n/a	MW 404.5993	C ₈ ² H ₂₀ H ₂ O ₄ P ₂ S ₄		
DRE-XA13270100AC	Ethion D20 100 µg/mL in Acetone		1.1ml	
Ethoprophos D5 (ethyl D5)				
CAS n/a	MW 247.3698	C ₆ ² H ₅ H ₁₄ O ₂ PS ₂		
DRE-XA13300010CY	Ethoprophos D5 (ethyl D5) 100 µg/mL in Cyclohexane		1.1ml	
Ethylbenzene-D10				
CAS 25837-05-2	MW 116.2266	C ₈ ² H ₁₀		
DRE-YA13320100ME	Ethylbenzene D10 2000 µg/mL in Methanol		1ml	
Ethylene thiourea D4				
CAS 352431-28-8	MW 106.1828	C ₃ ² H ₄ H ₂ N ₂ S		
DRE-C13330100	Ethylene thiourea D4		50mg	
DRE-XA13330100AC	Ethylene thiourea D4 100 µg/mL in Acetone		1.1ml	
Etofenprox D5 (ethyl D5)				
CAS 1705649-55-3	MW 381.5188	C ₂₅ ² H ₁₅ H ₂₃ O ₃		
DRE-XA13363010AC	Etofenprox D5 (ethyl D5) 100 µg/mL in Acetone(‡)		1.1ml	
Fenbuconazole (phenyl D5)				
CAS 1398066-06-2	MW 341.8489	C ₁₉ ² H ₅ H ₁₂ ClN ₄		
DRE-XA13448510AC	Fenbuconazole (phenyl D5) 100 µg/mL in Acetone		1.1ml	
Fenitrothion D6 (O,O-dimethyl D6)				
CAS 203645-59-4	MW 283.271	C ₉ ² H ₆ H ₆ NO ₅ PS		
DRE-C13480100	Fenitrothion D6 (O,O-dimethyl D6)		10mg	
DRE-XA13480100CY	Fenitrothion D6 (O,O-dimethyl D6) 100 µg/mL in Cyclohexane		1.1ml	
Fenitrothion-oxon D6 (O,O-dimethyl D6)				
CAS 1185155-54-7	MW 267.2054	C ₉ ² H ₆ H ₆ NO ₆ P		
DRE-XA13481100AC	Fenitrothion-oxon D6 (O,O-dimethyl D6) 100 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Fenpropidin D10 (piperidine D10)				
CAS n/a	MW 283.5178	C ₁₉ ² H ₁₀ H ₂₁ N		
DRE-XA13537100CY	Fenpropidin D10 (piperidine D10) 100 µg/mL in Cyclohexane		1.1ml	
Fenthion D6 (dimethoxy D6)				
CAS 1189662-83-6	MW 284.365	C ₁₀ ² H ₆ H ₉ O ₃ PS ₂		
DRE-C13580100	Fenthion D6 (O,O-dimethyl D6)		10mg	
DRE-XA13580100AC	Fenthion D6 (dimethoxy D6) 100 µg/mL in Acetone		1.1ml	
Fenuron D5 (phenyl D5)				
CAS 1219802-06-8	MW 169.2352	C ₉ ² H ₅ H ₇ N ₂ O		
DRE-XA13620010AL	Fenuron D5 (phenyl D5) 100 µg/mL in Acetonitrile		1.1ml	
Fenvalerate D7 (isopropyl D7)				
CAS n/a	MW 426.9432	C ₂₅ ² H ₁₇ H ₁₅ ClNO ₃		
DRE-XA13630010IO	Fenvalerate D7(isopropyl D7) 100 µg/mL in Isooctane		1.1ml	
Fluazifop-butyl D9 (n-butyl D9)				
CAS n/a	MW 392.4171	C ₁₉ ² H ₉ H ₁₁ F ₃ NO ₄		
DRE-XA13670100AC	Fluazifop-butyl D9 100 µg/mL in Acetone		1.1ml	
Flubendazole D3 (methyl D3)				
CAS 1173021-08-3	MW 316.3017	C ₁₆ ² H ₅ H ₉ FN ₃ O ₃		
DRE-C13678010	Flubendazole D3		10mg	
Fluconazole D4 (bismethylene D4)				
CAS 1124197-58-5	MW 310.2954	C ₁₃ ² H ₄ H ₆ F ₂ N ₆ O		
DRE-XA13697510AC	Fluconazole D4 (bismethylene D4) 100 µg/mL in Acetone		1.1ml	
Flunixin D3 (methyl D3)				
CAS 1015856-60-6	MW 299.263	C ₁₄ ² H ₅ H ₆ F ₃ N ₂ O ₂		
DRE-C13727010	Flunixin D3 (methyl D3)		10mg	
Fluoranthene-D10				
CAS 93951-69-0	MW 212.3122	C ₁₆ ² H ₁₀		
DRE-C20795100	Fluoranthene D10		50mg	
DRE-L20795100AC	Fluoranthene D10 10 µg/mL in Acetone		10ml	
DRE-L20795100ME	Fluoranthene D10 10 µg/mL in Methanol		10ml	
DRE-XA20795100AL	Fluoranthene D10 100 µg/mL in Acetonitrile		1ml	










Stable isotope labelled compounds

Product code	Description			
Fluorene-D10				
CAS 81103-79-9	MW 176.2801	C ₁₃ ² H ₁₀		
DRE-C20800200	Fluorene D10		100mg	
DRE-L20800200CY	Fluorene D10 10 µg/mL in Cyclohexane		10ml	
Folpet D4				
CAS 1327204-12-5	MW 300.5822	C ₉ ² H ₄ Cl ₃ NO ₂ S		
DRE-C13890100	Folpet D4(‡)		10mg	
DRE-XA13890100AC	Folpet D4 100 µg/mL in Acetone(‡)		1.1ml	
Forchlorfenuron D5 (phenyl D5)				
CAS 1398065-87-6	MW 252.7111	C ₁₂ ² H ₅ H ₉ ClN ₃ O		
DRE-XA13907010ME	Forchlorfenuron D5 (phenyl D5) 100 µg/mL in Methanol		1.1ml	
Fosetyl-aluminium D15				
CAS n/a	MW 369.197	3C ₂ ² H ₅ HO ₃ P·Al		
DRE-CA13940010	Fosetyl-aluminium D15 (*)		10mg	
Furan-d4				
CAS 6142-90-1	MW 72.0986	C ₄ ² H ₄ O		
DRE-XA13965010AL	Furan D4 100µg/mL in Acetonitrile		1ml	
Glyphosate 1,2-13C2 15N				
CAS 1185107-63-4	MW 172.0518	¹³ C ₂ CH ₈ ¹⁵ NO ₅ P		
DRE-XA14050100WA	Glyphosate 1,2-13C2 15N 100 µg/mL in Water(‡)		1ml	
α-HCH 13C6				
CAS 222966-66-7	MW 296.7858	¹³ C ₆ H ₆ Cl ₆		
DRE-XA14071300CY	α-HCH 13C6 100 µg/mL in Cyclohexane(‡)		1.1ml	
α-HCH D6				
CAS 86194-41-4	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-C14071400	α-HCH D6(‡)		10mg	
DRE-XA14071400CY	α-HCH D6 100 µg/mL in Cyclohexane(‡)		1ml	
β-HCH D6				
CAS 86194-42-5	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-L14072100CY	β-HCH D6 10 µg/mL in Cyclohexane		10ml	










Stable isotope labelled compounds

Product code	Description			
δ-HCH D6				
CAS n/a	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-XA14074100CY	δ-HCH D6 100 µg/mL in cyclohexane		1.1ml	AND Enantiomer 
γ-HCH D6				
CAS 60556-82-3	MW 296.8668	C ₆ ² H ₆ Cl ₆		
DRE-C14073100	γ-HCH D6		10mg	AND Enantiomer 
DRE-XA14073100CY	γ-HCH D6 100 µg/mL in Cyclohexane		1ml	
Hexachlorobenzene 13C6				
CAS 93952-14-8	MW 290.7381	¹³ C ₆ Cl ₆		
DRE-C14160100	Hexachlorobenzene 13C6(‡)		10mg	
DRE-XA14160100AC	Hexachlorobenzene 13C6 100 µg/mL in Acetone(‡)		1ml	
Hexachloro-1,3-butadiene 13C4				
CAS 93951-70-3	MW 264.7314	¹³ C ₄ Cl ₆		
DRE-XA14170100AC	Hexachloro-1,3-butadiene 13C4 100 µg/mL in Acetone(‡)		1.1ml	
Hexazinone D6 (N,N-dimethyl D6)				
CAS 1219804-22-4	MW 258.3498	C ₁₂ ² H ₆ H ₁₄ N ₄ O ₂		
DRE-XA14200010AL	Hexazinone D6 (N,N-dimethyl D6) 100 µg/mL in Acetonitrile		1ml	
Imidacloprid D4 (imidazolidin-4,4,5,5 D4)				
CAS 1015855-75-0	MW 259.6856	C ₉ ² H ₄ H ₆ ClN ₅ O ₂		
DRE-C14283710	Imidacloprid D4 (imidazolidin-4,4,5,5 D4)		10mg	
DRE-XA14283710AC	Imidacloprid D4 (imidazolidin-4,4,5,5 D4) 100 µg/mL in Acetone(‡)		1.1ml	
Indeno(1,2,3-c,d)pyrene D12				
CAS 203578-33-0	MW 288.4046	C ₂₂ ² H ₁₂		
DRE-LA20830200CY	Indeno(1,2,3-c,d)pyrene D12 10 µg/mL in Cyclohexane		1ml	
Iprodione D7 (isopropyl D7)				
CAS n/a	MW 337.2098	C ₁₃ ² H ₇ H ₆ Cl ₂ N ₃ O ₃		
DRE-XA14370012AL	Iprodione D7 100 µg/mL in Acetonitrile		1ml	
Ipronidazole D3 (N-methyl D3)				
CAS 1015855-83-0	MW 172.1996	C ₇ ² H ₃ H ₆ N ₃ O ₂		
DRE-C14370701	Ipronidazole D3		10mg	










Stable isotope labelled compounds

Product code	Description			
Ipronidazole-hydroxy D3				
CAS 1156508-86-9 DRE-C14370721	MW 188.199 Ipronidazole-hydroxy D3	C ₇ ² H ₅ H ₈ N ₃ O ₃	10mg	
Isoproturon D6				
CAS 1007461-76-8 DRE-C14470100 DRE-XA14470100AC DRE-YA14470100AL	MW 212.3211 Isoproturon D6 Isoproturon D6 100 µg/mL in Acetone(‡) Isoproturon D6 1000 µg/mL in Acetonitrile	C ₁₂ ² H ₆ H ₁₂ N ₂ O	10mg 1.1ml 1.1ml	
(±)-Ketoprofen D3 (methyl D3)				
CAS 159490-55-8 DRE-XA14532110AL	MW 257.299 (±)-Ketoprofen D3 (propionic D3 acid) 100 µg/mL in Acetonitrile	C ₁₆ ² H ₃ H ₁₁ O ₃	1ml	
Lenacil (2,2,2',2'-cyclohexyl D4)				
CAS n/a DRE-XA14610100AL	MW 238.3189 Lenacil (cyclohexyl D4) 100 µg/mL in Acetonitrile	C ₁₃ ² H ₄ H ₁₄ N ₂ O ₂	1.1ml	
Leucocrystal Violet D6				
CAS 1173023-92-1 DRE-C14629401	MW 379.5707 Leucocrystal violet D6	C ₂₅ ² H ₆ H ₂₅ N ₃	10mg	
Leucomalachite Green D6 (phenylmethin D6)				
CAS 1173021-13-0 DRE-C14629510	MW 336.5029 Leucomalachite green D6(‡)	C ₂₃ ² H ₆ H ₂₀ N ₂	10mg	
Linuron D6				
CAS 1219804-76-8 DRE-C14640100 DRE-XA14640100AC	MW 255.1309 Linuron D6 (methyl D3 methoxy D3)(‡) Linuron D6 100 µg/mL in Acetone(‡)	C ₉ ² H ₆ H ₄ Cl ₂ N ₂ O ₂	10mg 1.1ml	
MCPA D3 (phenyl D3)				
CAS 352431-14-2 DRE-C14760100 DRE-XA14760100AC	MW 203.6374 MCPA D3 MCPA D3 100 µg/mL in Acetone	C ₉ ² H ₃ H ₆ ClO ₃	10mg 1.1ml	
MCPA D6 (methyl-D3,phenoxy-D3)				
CAS n/a DRE-XA14760200AC	MW 206.6559 MCPA D6 100 µg/mL in Acetone	C ₉ ² H ₆ H ₃ ClO ₃	1.1ml	










Stable isotope labelled compounds

Product code	Description			
MCPB D6 (ring D3, methyl D3) (4-(4-Chloro-2-trideuteriomethyl-3,5,6-trideuteriophenoxy)butanoic Acid)				
CAS n/a	MW 234.7091	C ₁₁ ² H ₆ H ₇ ClO ₃		
DRE-C14790100	MCPB D6 (ring D3, methyl D3)		10mg	
DRE-XA14790100AC	MCPB D6 (ring D3, methyl D3) 100 µg/mL in Acetone		1.1ml	
Malachite Green D5 Oxalate Salt				
CAS n/a	MW 937.0664	2C ₂₃ ² H ₅ H ₂₀ N ₂ ·C ₂ H ₂ O ₄ ·2C ₂ H ₂ O ₄		
DRE-C14680005	Malachite green D5 oxalate salt		10mg	
Malachite Green D5 Picrate				
CAS 1258668-21-1	MW 562.5848	C ₂₃ ² H ₅ H ₂₀ N ₂ ·C ₆ H ₂ N ₃ O ₇		
DRE-C14680010	Malachite green D5 picrate(±)		10mg	
DRE-XA14680010AC	Malachite green D5 picrate 100 µg/mL in Acetone		1.1ml	
Malathion D10 (diethyl D10)				
CAS 347841-48-9	MW 340.4196	C ₁₀ ² H ₁₀ H ₉ O ₆ PS ₂		
DRE-C14710010	Malathion D10 (diethyl D10)		20mg	
Malathion D6 (dimethyl D6)				
CAS 1189877-72-2	MW 336.395	C ₁₀ ² H ₆ H ₁₃ O ₆ PS ₂		
DRE-XA14710020CY	Malathion D6 (dimethyl D6) 100 µg/mL in Cyclohexane (*)		1.1ml	
DRE-GH09010100AL	Malathion-d6 100 µg/mL in Acetonitrile(±)		10x1ml	
Malathion D7 (dimethyl D6,3-D1)				
CAS 352438-94-9	MW 337.4012	C ₁₀ ² H ₇ H ₁₂ O ₆ PS ₂		
DRE-C14710030	Malathion D7 (dimethyl D6,3-D1)		20mg	
Maleic Hydrazide D2				
CAS n/a	MW 114.0991	C ₄ ² H ₂ H ₂ N ₂ O ₂		
DRE-C14730100	Maleic hydrazide D2		10mg	
Mebendazole D3 (methyl D3)				
CAS 1173021-87-8	MW 298.3112	C ₁₆ ² H ₃ H ₁₀ N ₃ O ₃		
DRE-C14798010	Mebendazole D3		10mg	
Mecarbam D10 (diethoxy D10)				
CAS n/a	MW 339.4349	C ₁₀ ² H ₁₀ H ₁₀ NO ₅ PS ₂		
DRE-XA14800100AC	Mecarbam D10 (diethoxy D10) 100 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Mecoprop D3 (phenyl D3)				
CAS 352431-15-3 DRE-XA14820100AC	MW 217.664 Mecoprop D3 100 µg/mL in Acetone	C ₁₀ ² H ₈ H ₈ ClO ₃	1.1ml	
Mecoprop D6 (ring D3, methyl D3)				
CAS 1705649-54-2 DRE-C14820110 DRE-XA14820110AL	MW 220.6825 Mecoprop D6 (ring D3, methyl D3) Mecoprop D6 (ring D3, methyl D3) 100 µg/mL in Acetonitrile	C ₁₀ ² H ₆ H ₈ ClO ₃	10mg 1.1ml	
Mecoprop-methyl ester D3 (ring D3)				
CAS n/a DRE-XA14835100AC	MW 231.6906 Mecoprop-methyl ester D3 (ring) 100 µg/mL in Acetone	C ₁₁ ² H ₈ H ₁₀ ClO ₃	1.1ml	
Melamine Triamine-15N3				
CAS 287476-11-3 DRE-L14861401AL	MW 129.1002 Melamine triamine-15N3 10 µg/mL in Acetonitrile (*)	C ₃ H ₆ ¹⁵ N ₃	10ml	
Meloxicam D3 (2-methyl D3)				
CAS 942047-63-4 DRE-C14862510	MW 354.4192 Meloxicam D3 (2-methyl D3)	C ₁₄ ² H ₈ H ₁₀ N ₃ O ₄ S ₂	50mg	
Mepiquat iodide D3 (methyl-d3)				
CAS 32317-85-4 DRE-X14880100DO DRE-XA14880100DO	MW 244.1316 Mepiquat Iodide D3 100 µg/ml in Deuterium oxide(‡) Mepiquat iodide D3 100 µg/mL in Deuterium oxide	C ₇ ² H ₈ H ₁₃ N-I	10ml 1.1ml	
Methamidophos (dimethyl D6)				
CAS 1219799-41-3 DRE-C14980100	MW 147.1662 Methamidophos D6 (dimethyl D6)	C ₂ ² H ₆ H ₂ NO ₂ PS	10mg	
Methiocarb D3 (N-methyl D3)				
CAS 1581694-94-1 DRE-XA15020501CY	MW 228.3258 Methiocarb D3 (N-methyl D3) 100 µg/mL in Cyclohexane	C ₁₁ ² H ₈ H ₁₂ NO ₂ S	1.1ml	
Methomyl D3				
CAS 1398109-07-3 DRE-XA15030100AC	MW 165.2286 Methomyl D3 100 µg/mL in Acetone	C ₈ ² H ₈ H ₇ N ₂ O ₂ S	1.1ml	









Stable isotope labelled compounds

Product code	Description			
Methoxychlor D14 (bis(4-methoxyphenyl-D7))				
CAS 1644449-82-0	MW 359.7344	C ₁₆ ² H ₁₄ Cl ₂ O ₂		
DRE-C15060100	Methoxychlor D14 (bis(4-methoxyphenyl-D7))		10mg	
DRE-XA15060100AC	Methoxychlor D14 100 µg/mL in Acetone		1.1ml	
1-Methylnaphthalene D10				
CAS 38072-94-5	MW 152.2587	C ₁₁ ² H ₁₀		
DRE-C20895100	1-Methylnaphthalene D10		10mg	
2-(4'-Methylphenyl)propionic Acid 3,3,3-D3				
CAS n/a	MW 167.2196	C ₁₀ ² H ₈ H ₃ O ₂		
DRE-C15140801	2-(4'-Methylphenyl)propionic acid D3		10mg	
Metobromuron D6 (methyl D3 methoxy D3)				
CAS n/a	MW 265.1368	C ₉ ² H ₆ H ₃ BrN ₂ O ₂		
DRE-XA15160100AC	Metobromuron D6 (methyl D3 methoxy D3) 100 µg/mL in Acetone		1.1ml	
Metolachlor D6 (propyl D6)				
CAS 1219803-97-0	MW 289.8307	C ₁₅ ² H ₆ H ₁₆ ClNO ₂		
DRE-CA15170100	Metolachlor D6 (propyl D6)		10mg	
DRE-XA15170100AC	Metolachlor D6 100 µg/mL in Acetone		1.1ml	
Metronidazole D4 (ethylene D4)				
CAS 1261392-47-5	MW 175.1786	C ₆ ² H ₄ H ₃ N ₃ O ₃		
DRE-C15201001	Metronidazole D4 (ethylene D4)(‡)		10mg	
Metronidazole-hydroxy D2				
CAS n/a	MW 189.1657	C ₆ ² H ₂ H ₇ N ₃ O ₄		
DRE-C15201301	Metronidazole-hydroxy D2		10mg	
Metsulfuron-methyl D3 (triazine methoxy D3)				
CAS n/a	MW 384.3823	C ₁₄ ² H ₃ H ₁₂ N ₅ O ₆ S		
DRE-C15210100	Metsulfuron-methyl D3 (triazine methoxy D3)		10mg	
DRE-XA15210100AC	Metsulfuron-methyl D3 (triazine methoxy D3) 100 µg/mL in Acetone		1.1ml	
Monolinuron D6 (methyl D3 methoxy D3)				
CAS n/a	MW 220.6858	C ₉ ² H ₆ H ₃ ClN ₂ O ₂		
DRE-C15310100	Monolinuron D6 (methyl D3 methoxy D3)		10mg	
DRE-XA15310100AC	Monolinuron D6 (methyl D3 methoxy D3) 100 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Monuron D6 (dimethyl D6)				
CAS 217488-65-8	MW 204.6864	C ₉ ² H ₆ H ₅ ClN ₂ O		
DRE-C15320100	Monuron D6		5mg	
DRE-XA15320100AC	Monuron D6 100 µg/mL in Acetone		1.1ml	
Musk Xylene D15				
CAS 877119-10-3	MW 312.3564	C ₁₂ ² H ₁₅ N ₃ O ₆		
DRE-XA15360100AC	Musk xylene D15 100 µg/mL in Acetone		1.1ml	
2-NP-AMZO D5				
CAS 1173097-59-0	MW 339.358	C ₁₅ ² H ₈ H ₁₃ N ₄ O ₅		
DRE-C15654481	2-NP-AMZO D5		10mg	
2-NP-AOZ D4				
CAS 1007478-57-0	MW 239.2208	C ₁₆ ² H ₄ H ₅ N ₃ O ₄		
DRE-C15654501	2-NP-AOZ D4(‡)		10mg	
Naphthalene D8				
CAS 1146-65-2	MW 136.2198	C ₁₀ ² H ₈		
DRE-C20905100	Naphthalene D8(‡)		100mg	
DRE-L20905100CY	Naphthalene D8 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-YA20905100MB	Naphthalene D8 2000 µg/mL in Methyl-tert. butyl ether		1.1ml	
Nereistoxin Oxalate D6 (dimethyl D6)				
CAS n/a	MW 245.3494	C ₅₀ ² H ₆ H ₅ NS ₂ ·C ₂ H ₂ O ₄		
DRE-C15502010	Nereistoxin oxalate D6 (dimethyl D6)		10mg	
Nitrobenzene D5				
CAS 4165-60-0	MW 128.1402	C ₆ ² H ₅ NO ₂		
DRE-C15557100	Nitrobenzene D5		1g	
DRE-XA15557100AC	Nitrobenzene D5 100 µg/mL in Acetone		1.1ml	
DRE-YA15557100MB	Nitrobenzene D5 2000 µg/mL in Methyl-tert. butyl ether		1ml	
4-Nitrophenol-2,3,5,6-D4				
CAS 93951-79-2	MW 143.1334	C ₆ ² H ₄ HNO ₃		
DRE-C15590404	4-Nitrophenol-2,3,5,6-D4		100mg	
DRE-XA15590404AC	4-Nitrophenol-2,3,5,6-D4 100 µg/mL in Acetone		1.1ml	
(S)-N-Nitrosoanabasine D4 (pyridinyl D4)				
CAS 1426174-40-4	MW 195.2544	C ₁₆ ² H ₄ H ₉ N ₃ O		
DRE-LA15600501AC	N-Nitrosoanabasine 2,4,5,6 D4 10 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
(S)-N-Nitrosoanatabine D4 (pyridinyl D4)				
CAS 1426174-82-4	MW 193.2385	C ₁₆ ² H ₄ H ₇ N ₃ O		
DRE-LA15601001AC	N-Nitrosoanatabine 2,4,5,6 D4 10 µg/mL in Acetone		1.1ml	
N-Nitroso-dimethylamine D6				
CAS 17829-05-9	MW 80.1188	C ₂ ² H ₆ N ₂ O		
DRE-C15604010	N-Nitroso-dimethylamine D6		10mg	
DRE-XA15604010AC	N-Nitroso-dimethylamine D6 100 µg/mL in Acetone		1.1ml	
N-Nitroso-di-n-propylamine D14				
CAS 93951-96-3	MW 144.2744	C ₆ ² H ₁₄ N ₂ O		
DRE-XA15605010AC	N-Nitroso-di-n-propylamine D14 100 µg/mL in Acetone		1.1ml	
N-Nitrosonornicotine-ketone D4 (pyridyl D4)				
CAS 764661-24-7	MW 211.2538	C ₁₆ ² H ₄ H ₉ N ₃ O ₂		
DRE-LA15606301AC	N-Nitrosonornicotine-ketone D4 (pyridyl D4) 10 µg/mL in Acetone		1.1ml	
2-Nitrotoluene D7				
CAS 84344-04-7	MW 144.1791	C ₇ ² H ₇ NO ₂		
DRE-C15615205	2-Nitrotoluene D7		50mg	
4-Nonylphenol-mono-ethoxylate D2				
CAS n/a	MW n/a			
DRE-LA15631017AC	4-Nonylphenol-mono-ethoxylate D2 10 µg/mL in Acetone		1.1ml	No Structure
DRE-XA15631017CY	4-Nonylphenol-mono-ethoxylate D2 100 µg/mL in Cyclohexane		1.1ml	
4-n-Nonylphenol D8 (ring D4-ethylD4)				
CAS n/a	MW 228.3998	C ₁₈ ² H ₈ H ₁₆ O		
DRE-XA15630010AC	4-n-Nonylphenol D8 100 µg/mL in Acetone(‡)		1ml	
Norfloxacin D5 (ethyl D5)				
CAS 1015856-57-1	MW 324.3616	C ₁₈ ² H ₅ H ₁₃ FN ₃ O ₃		
DRE-C15648010	Norfloxacin D5(‡)		10mg	
Ofloxacin D3 (N-methyl D3)				
CAS 1173147-91-5	MW 364.386	C ₁₈ ² H ₃ H ₁₇ FN ₃ O ₄		
DRE-XA15717005AL	Ofloxacin D3 100 µg/mL in Acetonitrile		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Ofloxacin D3 Hydrochloride (N-methyl D3)				
CAS 1173021-78-7 DRE-C15717010	MW 400.8469 Ofloxacin D3 hydrochloride	$C_{18}^2H_3H_{17}FN_3O_4 \cdot ClH$	10mg	
Omethoate D6 (O-dimethyl D6)				
CAS 1219804-92-8 DRE-XA15730100AC	MW 219.2288 Omethoate D6 (dimethyl D6) 100 µg/mL in Acetone	$C_5^2H_3H_6NO_4PS$	1.1ml	
Oxamyl D3 (2-methyl-D3)				
CAS n/a DRE-XA15780100MB	MW 222.2799 Oxamyl D3 (2-methyl-D3) 100 µg/mL in Methyl-tert.butyl ether (*)	$C_7^2H_3H_{10}N_3O_3S$	1.1ml	
Para Red D6				
CAS 1014689-16-7 DRE-C15875100 DRE-XA15875100AC	MW 299.3138 Para Red D6 Pararot D6 100 µg/mL in Acetone	$C_{16}^2H_6H_5N_3O_3$	10mg 1.1ml	
Paraquat Dichloride D8				
CAS 347841-45-6 DRE-CA15870100	MW 265.2083 Paraquat dichloride D8	$C_{12}^2H_6H_6N_2 \cdot 2Cl$	50mg	
Paraquat dichloride D6 (dimethyl)				
CAS n/a DRE-C15870050	MW 263.1959 Paraquat dichloride D6	$C_{12}^2H_6H_6N_2 \cdot 2Cl$	50mg	
Paraquat diiodide D6 (dimethyl)				
CAS n/a DRE-C15870200	MW 446.0989 Paraquat diiodide D6	$C_{12}^2H_6H_6N_2 \cdot 2I$	50mg	
Parathion-ethyl D10 (diethyl D10)				
CAS 350820-04-1 DRE-C15880100 DRE-XA15880100AC	MW 301.3222 Parathion-ethyl D10 (diethyl D10) Parathion-ethyl D10 (diethyl D10) 100 µg/mL in Acetone	$C_{10}^2H_{10}H_4NO_5PS$	10mg 1.1ml	
Parathion-methyl D6 (dimethyl D6)				
CAS 96740-32-8 DRE-C15890100 DRE-XA15890100AC	MW 269.2444 Parathion-methyl D6 (dimethyl D6) Parathion-methyl D6 (dimethyl D6) 100 µg/mL in Acetone (*)	$C_8^2H_6H_4NO_5PS$	25mg 1.1ml	

Stable isotope labelled compounds

Product code	Description			
Pendimethalin D5 (pent-3-yl (2,2,3,4,4)-D5)				
CAS 1219803-39-0	MW 286.3385	C ₁₃ ² H ₉ H ₁₄ N ₃ O ₄		
DRE-C15930100	Pendimethalin D5		10mg	
DRE-XA15930100AC	Pendimethalin D5 100 µg/mL in Acetone		1.1ml	
Pentachloroanisole 13C6 (ring 13C)				
CAS n/a	MW 286.319	¹³ C ₆ CH ₃ Cl ₅ O		
DRE-XA15950010AC	Pentachloroanisole 13C6 100 µg/mL in Acetone		1.1ml	
Pentachlorophenol 13C6				
CAS 85380-74-1	MW 272.2925	¹³ C ₆ HCl ₅ O		
DRE-C15970100	Pentachlorophenol 13C6		10mg	
DRE-XA15970100CY	Pentachlorophenol 13C6 100 µg/mL in Cyclohexane		1.1ml	
DRE-GS09010309ME	Pentachlorophenol-13C6 1000 µg/mL in Methanol(‡)		5x1ml	
trans-Permethrin D6 (dimethyl D6)				
CAS 82523-59-9	MW 397.3247	C ₂₁ ² H ₆ H ₁₄ Cl ₂ O ₃		
DRE-XA15990201AC	trans-Permethrin D6 (dimethyl D6) 100 µg/mL in Acetone(‡)		1ml	
Perylene-d12				
CAS 1520-96-3	MW 264.3832	C ₂₀ ² H ₁₂		
DRE-C20915100	Perylene D12(‡)		100mg	
DRE-L20915100CY	Perylene D12 10 µg/mL in Cyclohexane		10ml	
DRE-YA20915100TO	Perylene D12 2000 µg/mL in Toluene		1ml	
Phenanthrene D10				
CAS 1517-22-2	MW 188.2908	C ₁₄ ² H ₁₀		
DRE-C20920100	Phenanthrene D10(‡)		100mg	
DRE-L20920100AC	Phenanthrene D10 10 µg/mL in Acetone		10ml	
DRE-L20920100CY	Phenanthrene D10 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-YA20920100MB	Phenanthrene D10 2000 µg/mL in Methyl-tert. butyl ether		1ml	
Phenmedipham D3 (methoxy D3)				
CAS 1773497-41-8	MW 303.3277	C ₁₆ ² H ₃ H ₁₃ N ₂ O ₄		
DRE-C16020100	Phenmedipham D3		5mg	
Phenol D5 (2,3,4,5,6-Pentadeuteriophenol)				
CAS 4165-62-2	MW 99.142	C ₆ ² H ₅ HO		
DRE-C16025100	Phenol D5		1g	
DRE-XA16025100AC	Phenol-D5 100 µg/mL in Acetone		1.1ml	
DRE-XA16025100ME	Phenol-D5 100 µg/mL in Methanol		1.1ml	
Phenol D6				
CAS 13127-88-3	MW 100.1482	C ₆ ² H ₆ O		
DRE-C16025200	Phenol D6		1g	










Stable isotope labelled compounds

Product code	Description			
Phorate (diethyl-D10)				
CAS 1219805-45-4 DRE-XA16080100AC	MW 270.4391 Phorate D10 100 µg/mL in Acetone	C ₇ ² H ₁₀ H ₇ O ₂ PS ₃	1.1ml	
Phosalone D10 (di-ethyl D5)				
CAS n/a DRE-XA16100100AC	MW 377.8702 Phosalone D10 (di-ethyl D5) 100 µg/mL in Acetone	C ₁₂ ² H ₁₀ H ₅ ClNO ₄ PS ₂	1.1ml	
Phosmet D6 (dimethoxy D3)				
CAS 2083623-41-8 DRE-XA16120100AC	MW 323.358 Phosmet D6 100 µg/mL in Acetone	C ₁₁ ² H ₆ H ₈ NO ₄ PS ₂	1.1ml	
Phthalic Acid Benzyl Butyl Ester (3,4,5,6)-D4				
CAS 93951-88-3 DRE-C16168010 DRE-XA16168010CY	MW 316.3843 Phthalic acid, benzylbutyl ester (3,4,5,6) D4 Phthalic acid, benzylbutyl ester (3,4,5,6) D4 100 µg/mL in Cyclohexane	C ₁₉ ² H ₁₄ H ₁₆ O ₄	10mg 1ml	
Phthalic Acid Bis(2-ethylhexyl) Ester (3,4,5,6)-D4				
CAS 93951-87-2 DRE-C16173010 DRE-CR16173010 DRE-XA16173010CY	MW 394.5808 Phthalic acid, bis-2-ethylhexyl ester D4 Phthalic acid, bis-2-ethylhexyl ester D4 Phthalic acid, bis-2-ethylhexyl ester D4 100 µg/mL in Cyclohexane	C ₂₄ ² H ₁₄ H ₃₄ O ₄	25mg 25mg 1ml	
Phthalic Acid Bis(isopentyl) Ester (3,4,5,6)-D4				
CAS 1346597-80-5 DRE-C16173685	MW 310.4213 Phthalic acid, bis-iso-pentyl ester D4	C ₁₈ ² H ₁₄ H ₂₂ O ₄	10mg	
Phthalic Acid Dibutyl Ester (3,4,5,6)-D4				
CAS 93952-11-5 DRE-C16171010	MW 282.3681 Phthalic acid, bis-butyl ester D4(‡)	C ₁₆ ² H ₁₄ H ₁₈ O ₄	10mg	
Phthalic Acid Dicyclohexyl Ester D4				
CAS 358731-25-6 DRE-C16171510	MW 334.4427 Phthalic acid, bis-cyclohexyl ester D4	C ₂₀ ² H ₁₄ H ₂₂ O ₄	10mg	
Phthalic Acid Diethyl Ester (3,4,5,6)-D4				
CAS 93952-12-6 DRE-C16172010	MW 226.2618 Phthalic acid, bis-ethyl ester D4(‡)	C ₁₂ ² H ₁₄ H ₁₀ O ₄	10mg	










Stable isotope labelled compounds

Product code	Description			
Phthalic Acid Di-n-hexyl Ester (3,4,5,6)-D4				
CAS 1015854-55-3 DRE-C16173210	MW 338.4744 Phthalic acid, bis-hexyl ester D4	C ₂₆ ² H ₄₄ H ₂₀ O ₄	10mg	
Phthalic Acid Diisobutyl Ester (3,4,5,6)-D4				
CAS 358730-88-8 DRE-C16173510	MW 282.3681 Phthalic acid, bis-iso-butyl ester D4	C ₁₆ ² H ₄ H ₁₈ O ₄	10mg	
Phthalic Acid Dimethyl Ester D4				
CAS 93951-89-4 DRE-C16174010	MW 198.2086 Phthalic acid, bis-methyl ester D4(‡)	C ₁₀ ² H ₄ H ₆ O ₄	10mg	
Phthalic Acid Di-n-nonyl Ester (3,4,5,6)-D4				
CAS 1202865-43-7 DRE-XA16174810CY	MW 422.6339 Phthalic acid, bis-n-nonyl ester D4 100 µg/mL in Cyclohexane(‡)	C ₂₆ ² H ₄₄ H ₃₈ O ₄	1ml	
Phthalic Acid Di-n-octyl Ester (3,4,5,6)-D4				
CAS 93952-13-7 DRE-C16175010	MW 394.5808 Phthalic acid, bis-n-octyl ester-3,4,5,6-D4	C ₂₄ ² H ₄₄ H ₃₄ O ₄	100mg	
Phthalic Acid Di-n-pentyl Ester (3,4,5,6)-D4				
CAS 358730-89-9 DRE-C16175510	MW 310.4213 Phthalic acid, bis-n-pentyl ester D4	C ₁₈ ² H ₄ H ₂₂ O ₄	10mg	
Phthalic Acid Di-n-propyl Ester (3,4,5,6)-D4				
CAS 358731-29-0 DRE-C16177010	MW 254.315 Phthalic acid, bis-propyl ester D4	C ₁₄ ² H ₄ H ₁₄ O ₄	10mg	
Phthalic Acid n-Pentyl Isopentyl Ester (3,4,5,6)-D4 (mixture of isomers)				
CAS DRE-C16179105	MW n/a Phthalic acid, n-pentyl-isopentyl ester D4 (mixture of isomers)		10mg	No Structure
Pirimicarb D6 (dimethylcarbamate D6)				
CAS 1015854-66-6 DRE-C16250100 DRE-XA16250100AC	MW 244.3232 Pirimicarb D6 (dimethylcarbamate D6) Pirimicarb D6 100 µg/mL in Acetone	C ₁₁ ² H ₆ H ₁₂ N ₄ O ₂	10mg 1.1ml	

Stable isotope labelled compounds

Product code	Description			
Pirimiphos-ethyl D10 (diethoxy D5)				
CAS n/a	MW 343.4483	C ₁₃ ² H ₁₀ H ₁₄ N ₃ O ₃ PS		
DRE-XA16260100AC	Pirimiphos-ethyl D10 100 µg/mL in Acetone		1.1ml	
Pirimiphos-methyl D6 (dimethoxy D3)				
CAS 1793055-06-7	MW 311.3705	C ₁₁ ² H ₆ H ₁₄ N ₃ O ₃ PS		
DRE-XA16270100AC	Pirimiphos-methyl D6 100 µg/mL in Acetone		1.1ml	
Prochloraz D7 (propyl D7)				
CAS n/a	MW 383.7086	C ₁₅ ² H ₇ H ₉ Cl ₃ N ₃ O ₂		
DRE-XA16290010AL	Prochloraz D7 (propyl D7) 100 µg/mL in Acetonitrile		1.1ml	
Prometon D14 (di-iso-propyl D14)				
CAS 1219803-45-8	MW 239.377	C ₁₀ ² H ₁₄ H ₅ N ₅ O		
DRE-XA16360005AL	Prometon D14 (di-iso-propyl D14) 100 µg/mL in Acetonitrile		1.1ml	
Prometryn D6 (isopropyl D6)				
CAS 1705649-52-0	MW 247.3933	C ₁₀ ² H ₆ H ₁₃ N ₅ S		
DRE-XA16370100AC	Prometryn D6 100 µg/mL in Acetone		1.1ml	
Propachlor D5 (phenyl D5)				
CAS n/a	MW 216.7188	C ₁₁ ² H ₅ H ₉ ClNO		
DRE-XA16380010AC	Propachlor D5 (phenyl D5) 100 µg/mL in Acetone		1.1ml	
Propamocarb free base D7 (O-propyl D7)				
CAS 1398065-89-8	MW 195.3104	C ₉ ² H ₇ H ₁₃ N ₂ O ₂		
DRE-C16390100	Propamocarb D7 (propyl D7)		10mg	
DRE-XA16390100AC	Propamocarb free base D7 100 µg/mL in Acetone		1.1ml	
Propazine 13C3 (ring 13C3)				
CAS 446276-68-2	MW 232.6878	¹³ C ₃ C ₆ H ₆ ClN ₅		
DRE-XA16440200AC	Propazine 13C3 100 µg/mL in Acetone		1.1ml	
Propazine D6 (isopropyl D6)				
CAS 1655498-05-7	MW 235.7468	C ₉ ² H ₆ H ₁₀ ClN ₅		
DRE-XA16440100AC	Propazine D6 100 µg/mL in Acetone		1.1ml	










Stable isotope labelled compounds

Product code	Description			
Propiconazole D5 (2,2,3,3,3-propyl-D5)				
CAS 1246818-14-3 DRE-XA16480100AC	MW 347.2512 Propiconazole D5 (2,2,3,3,3-propyl-D5) 100 µg/mL in Acetone	C ₁₅ ² H ₅ H ₁₂ Cl ₂ N ₃ O ₂	1.1ml	
Propoxur D3 (N-methyl D3)				
CAS 1219798-56-7 DRE-XA16500100AC	MW 212.2602 Propoxur D3 (methyl D3) 100 µg/mL in Acetone	C ₁₁ ² H ₅ H ₁₂ NO ₃	1.1ml	
Propyzamide D3 (phenyl-2,4,6 D3)				
CAS 1219805-79-4 DRE-XA16540010AL	MW 259.1463 Propyzamide D3 (phenyl-2,4,6 D3) 100 µg/mL in Acetonitrile	C ₁₂ ² H ₅ H ₈ Cl ₂ NO	1.1ml	
Pyrene-d10				
CAS 1718-52-1 DRE-C20930100 DRE-L20930100CY DRE-XA20930100AC DRE-XA20930100AL	MW 212.3122 Pyrene D10(±) Pyrene D10 10 µg/mL in Cyclohexane Pyrene D10 100 µg/mL in Acetone Pyrene D10 100 µg/mL in Acetonitrile	C ₁₆ ² H ₁₀	100mg 10ml 1ml 1ml	
Pyridine-d5				
CAS 7291-22-7 DRE-C16646100	MW 84.1307 Pyridine D5	C ₅ ² H ₅ N	1ml	
Quinmerac D4 (methyl(D3)-quinoline-2-D)				
CAS n/a DRE-XA16708100AL	MW 225.6644 Quinmerac D4 (methyl(D3)-quinoline-2-D) 100 µg/mL in Acetonitrile	C ₁₁ ² H ₄ H ₄ ClNO ₂	1.1ml	
Quinoline-d7 (Quinoline-2,3,4,5,6,7,8-D7)				
CAS 34071-94-8 DRE-XA16709601AC	MW 136.2017 Quinoline-2,3,4,5,6,7,8-D7 100 µg/mL in Acetone	C ₉ ² H ₇ N	1.1ml	
Quizalofop-ethyl D3 (3,3,3-D3)				
CAS 1398065-84-3 DRE-XA16740100AC	MW 375.8208 Quizalofop-ethyl D3 (3,3,3-D3) 100 µg/mL in Acetone	C ₁₉ ² H ₅ H ₁₄ ClN ₂ O ₄	1.1ml	
Quizalofop free acid D3 (methyl D3)				
CAS n/a DRE-XA16739991AC	MW 347.7676 Quizalofop free acid D3 100 µg/mL in Acetone	C ₁₇ ² H ₅ H ₁₀ ClN ₂ O ₄	1.1ml	










Stable isotope labelled compounds

Product code	Description			
Ronidazole D3				
CAS 1015855-87-4	MW 203.1706	C ₆ ² H ₅ H ₅ N ₄ O ₄		
DRE-C16815501	Ronidazole D3		10mg	
Salbutamol-D3				
CAS 1219798-60-3	MW 242.3292	C ₁₃ ² H ₈ H ₁₈ NO ₃		
DRE-XA16903001AL	Salbutamol D3 (3-hydroxymethyl-D2;a D1) 100 µg/mL in Acetonitrile(‡)		1ml	
Sarafloxacin D8 hydrochloride trihydrate				
CAS n/a	MW 483.9202	C ₂₀ ² H ₈ H ₉ F ₂ N ₃ O ₃ ·ClH·3H ₂ O		
DRE-C16908001	Sarafloxacin D8 hydrochloride trihydrate		10mg	
Sebuthylazine D5 (N-ethyl D5)				
CAS 1219805-56-7	MW 234.7406	C ₉ ² H ₅ H ₁₁ ClN ₅		
DRE-XA16920100AC	Sebuthylazine D5 100 µg/mL in Acetone		1.1ml	
Secbumeton D5 (N-ethyl D5)				
CAS 1705649-53-1	MW 230.3216	C ₁₀ ² H ₅ H ₁₄ N ₅ O		
DRE-C16930100	Secbumeton D5 (ethyl D5)		10mg	
DRE-XA16930100AC	Secbumeton D5 100 µg/mL in Acetone		1.1ml	
Simazine D10 (diethyl D5)				
CAS 220621-39-6	MW 211.7183	C ₇ ² H ₁₀ H ₂ ClN ₅		
DRE-C16950100	Simazine D10		10mg	
DRE-XA16950100AC	Simazine D10 100 µg/mL in Acetone		1.1ml	
Simazine D5 (ethyl D5)				
CAS 220621-41-0	MW 206.6875	C ₇ ² H ₉ H ₇ ClN ₅		
DRE-C16950200	Simazine D5		10mg	
DRE-XA16950200AL	Simazine D5 100 µg/mL in Acetonitrile(‡)		1ml	
Sotalol hydrochloride D6 (isopropyl-1,1,1,3,3,3-D6)				
CAS 1246820-85-8	MW 314.8617	C ₁₂ ² H ₆ H ₁₄ N ₂ O ₃ S·ClH		
DRE-XA16972631AL	Sotalol hydrochloride D6 100 µg/mL in Acetonitrile		1.1ml	
Styrene D8				
CAS 19361-62-7	MW 112.1984	C ₈ ² H ₈		
DRE-C16982010	Styrene D8		100mg	










Stable isotope labelled compounds

Product code	Description			
Sudan 1 D5 (phenyl D5)				
CAS 752211-63-5	MW 253.3101	C ₁₆ ² H ₅ H ₇ N ₂ O		
DRE-C16986105	Sudan 1 D5		10mg	
DRE-XA16986105AC	Sudan 1 D5 100 µg/mL in Acetone		1.1ml	
Sudan 4 D6 (naphthyl D6)				
CAS 1014689-18-9	MW 386.4788	C ₂₄ ² H ₆ H ₁₄ N ₄ O		
DRE-C16986108	Sudan 4 D6(‡)		10mg	
DRE-XA16986108AC	Sudan 4 D6 100 µg/mL in Acetone		1.1ml	
Sulfachloropyridazine 13C6 (benzene 13C6)				
CAS n/a	MW 290.678	¹³ C ₆ C ₄ H ₉ ClN ₄ O ₂ S		
DRE-XA16990102AL	Sulfachloropyridazine 13C6 100 µg/mL in Acetonitrile		1ml	
Sulfadimethoxine D6 (2,6-dimethoxy D6)				
CAS 73068-02-7	MW 316.3659	C ₁₂ ² H ₆ H ₈ N ₄ O ₄ S		
DRE-C16990551	Sulfadimethoxine D6 (2,6-dimethoxy D6)		10mg	
Sulfamerazine D4				
CAS 1020719-84-9	MW 268.3282	C ₁₁ ² H ₄ H ₈ N ₄ O ₂ S		
DRE-C16995110	Sulfamerazine D4		10mg	
Sulfamethoxazole D4 (benzene D4)				
CAS 1020719-86-1	MW 257.3023	C ₁₀ ² H ₄ H ₇ N ₃ O ₃ S		
DRE-C16998110	Sulfamethoxazole D4 (benzene D4)		10mg	
DRE-XA16998110AL	Sulfamethoxazole D4 (benzene D4) 100 µg/mL in Acetonitrile		1ml	
2,4,5-T D4				
CAS 358731-37-0	MW 259.5071	C ₈ ² H ₄ HCl ₃ O ₃		
DRE-XA17100100AC	2,4,5-T D4 100 µg/mL in Acetone		1.1ml	
Tebuconazole D6				
CAS n/a	MW 313.8554	C ₁₆ ² H ₆ H ₁₆ ClN ₃ O		
DRE-XA17178710AC	Tebuconazole D6 100 µg/mL in Acetone(‡)		1ml	
Temephos D12 (tetramethyl D12)				
CAS 1219795-39-7	MW 478.5429	C ₁₆ ² H ₁₂ H ₈ O ₆ P ₂ S ₃		
DRE-XA17220100CY	Temephos D12 (O,O,O',O'-tetramethyl D12) 100 µg/mL in Cyclohexane		1.1ml	








Stable isotope labelled compounds

Product code	Description			
Terbuthylazine D5 (ethyl D5)				
CAS 222986-60-9	MW 234.7406	C ₉ ² H ₅ H ₁₁ ClN ₅		
DRE-C17300100	Terbuthylazine D5		5mg	
DRE-XA17300100AC	Terbuthylazine D5 100 µg/mL in Acetone		1.1ml	
Terbuthylazine-desethyl D9 (tert-butyl D9)				
CAS 1219798-52-3	MW 210.7121	C ₇ ² H ₅ H ₃ ClN ₅		
DRE-C17303100	Terbuthylazine-desethyl D9 (tert-butyl D9)		10mg	
DRE-XA17303100AC	Terbuthylazine-desethyl D9 (tert-butyl D9) 100 µg/mL in Acetone		1.1ml	
Terbutryn D5 (ethyl D5)				
CAS 1219804-47-3	MW 246.3872	C ₁₀ ² H ₅ H ₁₄ N ₅ S		
DRE-C17320100	Terbutryn D5		10mg	
DRE-XA17320100AC	Terbutryn D5 100 µg/mL in Acetone		1.1ml	
p-Terphenyl D14				
CAS 1718-51-0	MW 244.39	C ₁₈ ² H ₁₄		
DRE-C20935300	p-Terphenyl D14		10mg	
DRE-GA09010298DI	p-Terphenyl-d14 500 µg/mL in Dichloromethane(‡)		1ml	
Tetrabromobisphenol A D10 (dimethyl D3, bisphenol-3,5-D2)				
CAS n/a	MW 553.9322	C ₁₅ ² H ₁₀ H ₂ Br ₄ O ₂		
DRE-XA17324701AL	Tetrabromobisphenol A D10 100 µg/mL in Acetonitrile		1.1ml	
1,2,4,5-Tetrachlorobenzene-13C6				
CAS 85380-73-0	MW 221.848	¹³ C ₆ H ₂ Cl ₄		
DRE-XA17354501AL	1,2,4,5-Tetrachlorobenzene 13C6 100 µg/mL in Acetonitrile		1ml	
3,3,4,4-Tetradeuterioadipic Acid Bis(2-ethylhexyl) Ester				
CAS n/a	MW 374.5911	C ₂₂ ² H ₄ H ₃₈ O ₄		
DRE-C10046010	Adipic acid-bis-2-ethylhexyl ester D4		10mg	
DRE-XA10046010AC	Adipic acid-bis-2-ethylhexyl ester D4 100 µg/mL in Acetone		1.1ml	
Tetramisole D5 hydrochloride				
CAS 1173021-85-6	MW 245.7831	C ₁₁ ² H ₅ H ₇ N ₂ S·ClH		
DRE-C17414501	Tetramisole D5 hydrochloride		10mg	
Thiamethoxam D4 (oxadiazine D4)				
CAS 1331642-98-8	MW 295.7393	C ₈ ² H ₄ H ₆ ClN ₅ O ₃ S		
DRE-XA17453010AC	Thiamethoxam D4 (oxadiazine D4) 100 µg/mL in Acetone(‡)		1.1ml	

Stable isotope labelled compounds

Product code	Description			
Thiophenuron-methyl D3				
CAS n/a DRE-C17466100	MW 390.41 Thiophenuron-methyl D3	C ₁₂ ² H ₃ H ₁₀ N ₅ O ₆ S ₂	5mg	
Thiofanox D3 (N-methyl D3)				
CAS 74829-46-2 DRE-XA17500100AC	MW 221.3349 Thiofanox D3 (N-methyl D3) 100 µg/mL in Acetone	C ₈ ² H ₃ H ₁₅ N ₂ O ₂ S	1ml	
Thionazin D10 (diethyl D10)				
CAS 1173019-39-0 DRE-C17530100	MW 258.3008 Thionazin D10 (diethyl-d10)	C ₈ ² H ₁₀ H ₃ N ₂ O ₃ PS	10mg	
Thiram D12				
CAS 69193-86-8 DRE-X17570100AL	MW 252.5068 Thiram D12 100 µg/mL in Acetonitrile (*)	C ₆ ² H ₁₂ N ₂ S ₄	10ml	
Tiabendazole D6 (Thiabendazole D6)				
CAS 1262551-89-2 DRE-C17450100 DRE-XA17450100AC	MW 207.2847 Thiabendazole NH D6 Thiabendazole NH D6 100 µg/mL in Acetone	C ₁₀ ² H ₆ HN ₃ S	10mg 1.1ml	
Toluene D8				
CAS 2037-26-5 DRE-C17594100	MW 100.1877 Toluene D8	C ₇ ² H ₈	.5ml	
2,4,6-Trichloroanisole D5				
CAS 352439-08-8 DRE-C17714601 DRE-XA17714601AL	MW 216.5038 2,4,6-Trichloroanisole D5 2,4,6-Trichloroanisole D5 100 µg/mL in Acetonitrile	C ₇ ² H ₅ Cl ₃ O	25mg 1ml	
1,3,5-Trichlorobenzene D3				
CAS 1198-60-3 DRE-XA17723600AC	MW 184.4655 1,3,5-Trichlorobenzene D3 100 µg/mL in Acetone	C ₆ ² H ₃ Cl ₃	1.1ml	
Triclosan D3 (2,4-dichlorophenoxy D3)				
CAS 1020719-98-5 DRE-XA17803010CY	MW 292.5603 Triclosan D3 (2,4-dichlorophenoxy D3) 100 µg/mL in Cyclohexane	C ₁₂ ² H ₅ H ₄ Cl ₂ O ₂	1.1ml	

Stable isotope labelled compounds

Product code	Description			
Triclosan methyl D3 (methoxy D3)				
CAS 1020720-00-6	MW 306.5868	C ₁₃ ² H ₅ H ₆ Cl ₃ O ₂		
DRE-XA17803310AC	Triclosan methyl D3 (methoxy D3) 100 µg/mL in Acetone		1.1ml	
Trietazine D5 (ethyl D5)				
CAS 1397243-73-0	MW 234.7406	C ₉ ² H ₅ H ₁₁ ClN ₅		
DRE-XA17830100AC	Trietazine D5 100 µg/mL in Acetone		1.1ml	
Triflumizole 13C3 (O-propyl 13C3)				
CAS n/a	MW 348.7253	¹³ C ₃ C ₁₂ H ₁₅ ClF ₃ N ₃ O		
DRE-XA17844010AC	Triflumizole 13C3 100 µg/mL in Acetone		1.1ml	
Trifluralin D14 (di-n-propyl D14)				
CAS 347841-79-6	MW 349.3653	C ₁₃ ² H ₁₄ H ₂ F ₃ N ₃ O ₄		
DRE-C17850100	Trifluralin D14 (di-n-propyl D14)		10mg	
DRE-XA17850100AC	Trifluralin D14 (di-n-propyl D14) 100 µg/mL in Acetone (*)		1.1ml	
Triforine D8 (piperazine D8)				
CAS 948595-12-8	MW 443.0111	C ₁₀ ² H ₆ H ₆ Cl ₆ N ₄ O ₂		
DRE-C17860100	Triforine D8 (piperazine D8)		10mg	
Trimethoprim D3 (4-methoxy D3)				
CAS 1189923-38-3	MW 293.3362	C ₁₄ ² H ₅ H ₁₅ N ₄ O ₃		
DRE-C178750100	Trimethoprim D3 (4-methoxy D3)		10mg	
DRE-XA17875010AL	Trimethoprim D3 (4-methoxy D3) 100 µg/mL in Acetonitrile		1ml	
(±)-Warfarin D5 (phenyl-D5)				
CAS 75472-93-4	MW 313.3587	C ₁₉ ² H ₅ H ₁₁ O ₄		
DRE-C17940100	(±)-Warfarin D5 (phenyl-D5)		10mg	
DRE-XA17940100AL	(±)-Warfarin D5 (phenyl-D5) 100 µg/mL in Acetonitrile		1.1ml	
PAH-Mix 9 deuterated				
DRE-L20950902CY	PAH-Mix 9 deuterated 10 µg/mL in Cyclohexane(‡)		10ml	
DRE-XA20950902CY	PAH-Mix 9 deuterated 100 µg/mL in Cyclohexane(‡)		1ml	
Acenaphthene D10	Acenaphthylene D8	Anthracene D10	Benz[a]anthracene D12	
Benzo(a)pyrene D12	Benzo(k)fluoranthene D12	Benzo[b]fluoranthene D12	Benzo[g,h,i]perylene D12	
Chrysene D12	Dibenz[a,h]anthracene D14	Fluoranthene D10	Fluorene D10	
Indeno(1,2,3-c,d)pyrene D12	Naphthalene D8	Phenanthrene D10	Pyrene D10	
PAH-Mix 24 deuterated				
DRE-LA20950024HE	PAH-Mix 24 deuterated 10 µg/mL in n-Hexane		1ml	
	Acenaphthene D10			
	Chrysene D12			
	Naphthalene D8			
	Perylene D12			
	Phenanthrene D10			

(†) ISO Guide 34

(*) Shorter expiry due to chemical nature of component(s)

Order online at www.lgcstandards.com

Stable isotope labelled compounds










Product code	Description	
PAH-Mix 31 deuterated		
DRE-YA20950031TO	PAH-Mix 31 deuterated 1000 µg/mL in Toluene(‡) Acenaphthene D10 Chrysene D12 Naphthalene D8 Perylene D12 Phenanthrene D10	1ml

The background of the entire page is a high-speed photograph of water bubbles. The bubbles are of various sizes, from tiny specks to larger, more defined spheres, and are scattered across the white background. Some bubbles are in sharp focus, showing their internal structure and the way light reflects off their surfaces. The overall effect is a sense of movement and freshness.










ADDITIONAL ORGANIC CONTAMINANTS



Additional organic contaminants

Product code	Description			
AHMI (Phantolide)				
CAS 15323-35-0 DRE-LA10048000CY	MW 244.3719 AHMI (Phantolide) 10 µg/mL in Cyclohexane	C ₁₇ H ₂₄ O	1ml	
AHTN (Tonalide)				
CAS 1506-02-1 DRE-LA10048500CY	MW 258.3984 AHTN (Tonalide) 10 µg/mL in Cyclohexane	C ₁₈ H ₂₆ O	1ml	
AHTN (Tonalide) (6-acetyl D3)				
CAS 1396967-82-0 DRE-XA10048600IO	MW 261.4169 AHTN (Tonalide) D3 100 µg/mL in Isooctane	C ₁₈ ² H ₃ H ₂₃ O	1.1ml	
ATII (Traseolid)				
CAS 68140-48-7 DRE-LA10316000CY	MW 258.3984 ATII (Traseolid) 10 µg/mL in Cyclohexane	C ₁₈ H ₂₆ O	1ml	
Acetonitrile				
CAS 75-05-8 DRE-C10021000	MW 41.0519 Acetonitrile	C ₂ H ₃ N	5ml	
2-Acetylfuran				
CAS 1192-62-7 DRE-C10023800	MW 110.1106 2-Acetylfuran	C ₆ H ₆ O ₂	1ml	
Acrylic Acid Butyl Ester				
CAS 141-32-2 DRE-CA10045350	MW 128.169 Acrylic acid butyl ester	C ₇ H ₁₂ O ₂	1ml	
Acrylic Acid Ethyl Ester (Ethyl Acrylate)				
CAS 140-88-5 DRE-C10045380 DRE-YS09010025ME	MW 100.1158 Acrylic acid-ethyl ester Ethyl Acrylate 1000 µg/mL in p/t Methanol(±)	C ₅ H ₈ O ₂	1ml 5x1ml	
Acrylic Acid 2-Ethylhexyl Ester				
CAS 103-11-7 DRE-CA10045390	MW 184.2753 Acrylic acid-2-ethylhexyl ester	C ₁₁ H ₂₀ O ₂	1ml	










Additional organic contaminants

Product code	Description			
Acrylic Acid Methyl Ester				
CAS 96-33-3 DRE-CA10045400	MW 86.0892 Acrylic acid methyl ester	C ₄ H ₆ O ₂	1ml	
DL-Alanine				
CAS 302-72-7 DRE-C10062950	MW 89.0932 DL-Alanine	C ₃ H ₇ NO ₂	100mg	
4-Allylanisole (Estragole)				
CAS 140-67-0 DRE-C10131000	MW 148.2017 4-Allylanisole	C ₁₀ H ₁₂ O	1ml	
4-Aminobenzoic Acid				
CAS 150-13-0 DRE-C10171400	MW 137.136 4-Aminobenzoic acid	C ₇ H ₇ NO ₂	1g	
2-Aminobiphenyl				
CAS 90-41-5 DRE-C10173020	MW 169.2224 2-Aminobiphenyl	C ₁₂ H ₁₁ N	100mg	
3-Aminobiphenyl				
CAS 2243-47-2 DRE-C10173030	MW 169.2224 3-Aminobiphenyl	C ₁₂ H ₁₁ N	100mg	
4-Aminobiphenyl				
CAS 92-67-1 DRE-C10173040 DRE-L10173040CY	MW 169.2224 4-Aminobiphenyl 4-Aminobiphenyl 10 µg/mL in Cyclohexane	C ₁₂ H ₁₁ N	100mg 10ml	
4-Aminobiphenyl D9				
CAS 344298-96-0 DRE-XA10173041AC	MW 178.2779 4-Aminobiphenyl D9 100 µg/mL in Acetone(‡)	C ₁₂ ² H ₉ H ₂ N	1ml	
3-Amino-2-methyl-5-nitrobenzamide				
CAS 3572-44-9 DRE-C10204953	MW 195.1754 3-Amino-2-methyl-5-nitrobenzamide (*)	C ₈ H ₉ N ₃ O ₃	10mg	










Additional organic contaminants

Product code	Description			
5-Amino-2-methyl-3-nitrobenzamide				
CAS 90223-31-7 DRE-C10204954	MW 195.1754 5-Amino-2-methyl-3-nitrobenzamide (*)	C ₈ H ₉ N ₃ O ₃	10mg	
3-Aminopyridine				
CAS 462-08-8 DRE-C10221000	MW 94.1145 3-Aminopyridine	C ₅ H ₆ N ₂	500mg	
4-Aminosalicylic Acid				
CAS 65-49-6 DRE-C10227010	MW 153.1354 4-Aminosalicylic acid	C ₇ H ₇ NO ₃	100mg	
4-Aminotoluene-3-sulfonic Acid				
CAS 88-44-8 DRE-C10228000	MW 187.2163 4-Aminotoluene-3-sulfonic acid	C ₇ H ₉ NO ₃ S	100mg	
Ammelide				
CAS 645-93-2 DRE-C10241000	MW 128.0895 Ammelide	C ₃ H ₄ N ₄ O ₂	100mg	
α-Amylcinnamal				
CAS 122-40-7 DRE-C10246000	MW 202.2921 α-Amylcinnamal	C ₁₄ H ₁₈ O	250mg	
α-Amylcinnamyl Alcohol				
CAS 101-85-9 DRE-C10246200	MW 204.308 α-Amylcinnamyl alcohol	C ₁₄ H ₂₀ O	100mg	
rac-Anatabine				
CAS 2743-90-0 DRE-CA10249000	MW 160.2157 (R,S)-Anatabine (*)	C ₁₀ H ₁₂ N ₂	10mg	
Aniline-2-sulfonic Acid				
CAS 88-21-1 DRE-C10262900	MW 173.1897 Aniline-2-sulfonic acid	C ₆ H ₇ NO ₃ S	250mg	










Additional organic contaminants

Product code	Description			
Anise Alcohol				
CAS 105-13-5 DRE-C10265500	MW 138.1638 Anise alcohol	C ₈ H ₁₀ O ₂	250mg	
p-Anisic acid (4-Methoxybenzoic acid)				
CAS 100-09-4 DRE-C10265900	MW 152.1473 p-Anisic acid (4-Methoxybenzoic acid)	C ₈ H ₈ O ₃	250mg	
2-Anisidine (2-Methoxyaniline)				
CAS 90-04-0 DRE-C10266000 DRE-XA10266000AL	MW 123.1525 2-Anisidine (2-Methoxyaniline) 2-Anisidine (2-Methoxyaniline) 100 µg/mL in Acetonitrile	C ₇ H ₉ NO	250mg 1ml	
3-Anisidine (3-Methoxyaniline)				
CAS 536-90-3 DRE-C10266010	MW 123.1525 3-Anisidine (3-Methoxyaniline)	C ₇ H ₉ NO	250mg	
4-Anisidine (4-Methoxyaniline)				
CAS 104-94-9 DRE-C10266020	MW 123.1525 4-Anisidine (4-Methoxyaniline)	C ₇ H ₉ NO	250mg	
DL-Arginine Hydrochloride				
CAS 32042-43-6 DRE-C10300200	MW 210.6619 DL-Arginine hydrochloride	C ₆ H ₁₄ N ₄ O ₂ ·ClH	100mg	
DL-Asparagine Hydrate				
CAS 69833-18-7 DRE-C10304920	MW 150.1332 DL-Asparagine hydrate	C ₄ H ₈ N ₂ O ₃ ·H ₂ O	100mg	
DL-Aspartic Acid				
CAS 617-45-8 DRE-C10304950	MW 133.1027 DL-Aspartic acid	C ₄ H ₇ NO ₄	100mg	
Benzamide				
CAS 55-21-0 DRE-C10532250	MW 121.1366 Benzamide	C ₇ H ₇ NO	250mg	

Additional organic contaminants

Product code	Description			
Benzanilide				
CAS 93-98-1 DRE-C10532300	MW 197.2325 Benzanilide	C ₁₃ H ₁₁ NO	250mg	
Benzene-1,2-dicarboxylic Acid (Phthalic Acid)				
CAS 88-99-3 DRE-C16167500	MW 166.1308 Phthalic acid	C ₈ H ₆ O ₄	250mg	
Benzene-1,2,3-triol (Pyrogallol)				
CAS 87-66-1 DRE-C10535600	MW 126.11 Benzene-1,2,3-triol	C ₆ H ₆ O ₃	100mg	
Benzil				
CAS 134-81-6 DRE-C10536300	MW 210.228 Benzil	C ₁₄ H ₁₀ O ₂	1g	
Benzocaine (4-Aminobenzoic acid ethyl ester)				
CAS 94-09-7 DRE-C10171450	MW 165.1891 4-Aminobenzoic acid ethyl ester (Benzocaine)	C ₉ H ₁₁ NO ₂	100mg	
Benzoguanamine (6-Phenyl-1,3,5-triazine-2,4-diamine)				
CAS 91-76-9 DRE-C10537450	MW 187.2013 Benzoguanamine	C ₉ H ₉ N ₅	100mg	
Benzonitrile				
CAS 100-47-0 DRE-C10538500	MW 103.1213 Benzonitrile	C ₇ H ₅ N	250mg	
Benzophenone				
CAS 119-61-9 DRE-C10539000	MW 182.2179 Benzophenone	C ₁₃ H ₁₀ O	250mg	
1,4-Benzoquinone				
CAS 106-51-4 DRE-C10537000	MW 108.0948 1,4-Benzochinone(†)	C ₆ H ₄ O ₂	250mg	










Additional organic contaminants

Product code	Description			
5-Benzoyl-4-hydroxy-2-methoxybenzenesulfonic Acid				
CAS 4065-45-6 DRE-C10545000	MW 308.3065 5-Benzoyl-4-hydroxy-2-methoxybenzenesulfonic acid	C ₁₄ H ₁₂ O ₆ S	100mg	
Benzyl Alcohol				
CAS 100-51-6 DRE-C10569000	MW 108.1378 Benzylalcohol	C ₇ H ₈ O	1ml	
Biphenyl				
CAS 92-52-4 DRE-L10630000CY DRE-XA10630000ME	MW 154.2078 Biphenyl 10 µg/mL in Cyclohexane Biphenyl 100 µg/mL in Methanol	C ₁₂ H ₁₀	10ml 1ml	
Biphenyl D10				
CAS 1486-01-7 DRE-LA10630010AC	MW 164.2694 Biphenyl D10 10 µg/mL in Acetone	C ₁₂ ² H ₁₀	1ml	
2,2'-Bipyridyl				
CAS 366-18-7 DRE-C10640000	MW 156.1839 2,2'-Bipyridyl	C ₁₀ H ₈ N ₂	500mg	
4,4'-Bipyridyl (4,4'-Bipyridine)				
CAS 553-26-4 DRE-C10642000	MW 156.1839 4,4'-Bipyridyl	C ₁₀ H ₈ N ₂	500mg	
Alpha-Bisabolol (α-Bisabolol)				
CAS 515-69-5 DRE-GA09010039IP	MW 222.3663 Alpha-Bisabolol 1000 µg/mL in Isopropanol(±)	C ₁₅ H ₂₆ O	1ml	
Bis-(2-chloroethoxy)-methane				
CAS 111-91-1 DRE-XA10651000CY	MW 173.0377 Bis-(2-chloroethoxy)-methane 100 µg/mL in Cyclohexane	C ₆ H ₁₀ Cl ₂ O ₂	1ml	
Bis(4-dimethylaminophenyl)methane				
CAS 101-61-1 DRE-C10651900	MW 254.37 Bis-(4-dimethylaminophenyl)methane	C ₁₇ H ₂₂ N ₂	250mg	










Additional organic contaminants

Product code	Description			
Bis(2-ethylhexyl) phosphate				
CAS 298-07-7 DRE-C10652000	MW 322.4205 Bis(2-ethylhexyl) phosphate	$C_{16}H_{35}O_4P$	250mg	
Bis(4-glycidyloxyphenyl)methane				
CAS 2095-03-6 DRE-C10653400	MW 312.3597 Bis(4-glycidyloxyphenyl)methane(‡)	$C_{19}H_{20}O_4$	250mg	
2,2-Bis-(4-glycidyloxyphenyl)propane				
CAS 1675-54-3 DRE-C10653500	MW 340.4129 2,2-Bis-(4-glycidyloxyphenyl)propane	$C_{21}H_{24}O_4$	100mg	
1,3-Bis(hydroxymethyl)-2-imidazolidinone				
CAS 136-84-5 DRE-C10653530	MW 146.1445 1,3-Bis(hydroxymethyl)-2-imidazolidinone	$C_5H_{10}N_2O_3$	50mg	
N,N'-Bis-(4-nitrophenyl)urea D8				
CAS 1156508-87-0 DRE-C15598600	MW 310.2916 N,N'-Bis-(4-nitrophenyl)urea D8	$C_{13}^2H_8H_2N_4O_5$	10mg	
Bisphenol A				
CAS 80-05-7 DRE-C10655500	MW 228.2863 Bisphenol A(‡)	$C_{15}H_{16}O_2$	250mg	
Bisphenol A D14				
CAS 120155-79-5 DRE-XA10655503AL	MW 242.3726 Bisphenol A D14 100 µg/mL in Acetonitrile	$C_{15}^2H_{14}H_2O_2$	1.1ml	
Bisphenol A D16				
CAS 96210-87-6 DRE-C10655501	MW 244.3849 Bisphenol A D16	$C_{15}^2H_{16}O_2$	50mg	
Bisphenol A Bis(2,3-dihydroxypropyl) Ether				
CAS 5581-32-8 DRE-C10655620	MW 376.4434 Bisphenol A Bis(2,3-dihydroxypropyl) Ether	$C_{21}H_{28}O_6$	50mg	










Additional organic contaminants

Product code	Description			
Bisphenol B				
CAS 77-40-7 DRE-C10655670	MW 242.3129 Bisphenol B	C ₁₆ H ₁₈ O ₂	100mg	
Bisphenol F				
CAS 620-92-8 DRE-C10655800	MW 200.2332 Bisphenol F	C ₁₃ H ₁₂ O ₂	100mg	
Bisphenol S				
CAS 80-09-1 DRE-C10655940	MW 250.2704 Bisphenol S	C ₁₂ H ₁₀ O ₄ S	100mg	
Bromoacetic Acid				
CAS 79-08-3 DRE-C10697000 DRE-YA10697000MB	MW 138.948 Bromoacetic acid Bromoacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ H ₃ BrO ₂	1g 1ml	
Bromoacetic Acid Methyl Ester				
CAS 96-32-2 DRE-YA10697200MB	MW 152.9746 Bromoacetic acid methyl ester 1000 µg/mL in Methyl-tert. butyl ether	C ₃ H ₅ BrO ₂	1ml	
Bromochloroacetic Acid				
CAS 5589-96-8 DRE-YA10713000MB	MW 173.3931 Bromochloroacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ H ₂ BrClO ₂	1ml	
Bromochloroacetic Acid Methyl Ester				
CAS 20428-74-4 DRE-YA10713200MB	MW 187.4197 Bromochloroacetic acid-methyl ester 1000 µg/mL in Methyl-tert. butyl ether	C ₃ H ₄ BrClO ₂	1ml	
Bromochlorophene				
CAS 15435-29-7 DRE-C10721200	MW 426.9154 Bromochlorophene	C ₁₃ H ₈ Br ₂ Cl ₂ O ₂	50mg	
Bromodichloroacetic acid				
CAS 71133-14-7 DRE-C10726400	MW 207.8381 Bromodichloroacetic acid	C ₂ HBrCl ₂ O ₂	10mg	







Additional organic contaminants

Product code	Description			
4-Bromo-3,5-dimethylphenyl-N-methylcarbamate (BDMC)				
CAS 672-99-1 DRE-C10727000	MW 258.1118 4-Bromo-3,5-dimethylphenyl-N-methylcarbamate (BDMC)	C ₁₀ H ₁₂ BrNO ₂	100mg	
4-Bromodiphenyl Ether (PBDE 3)				
CAS 101-55-3 DRE-C15898003 DRE-XA15898003ME	MW 249.1033 PBDE 3 PBDE 3 100 µg/mL in Methanol	C ₁₂ H ₉ BrO	250mg 1ml	
2-Bromopropionic Acid				
CAS 598-72-1 DRE-C10760700 DRE-YA10760700MB	MW 152.9746 2-Bromopropionic acid 2-Bromopropionic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₃ H ₅ BrO ₂	250mg 1ml	
Bromothalonil (1,2-Dibromo-2,4-dicyanobutane)				
CAS 35691-65-7 DRE-C10764000	MW 265.9332 Bromothalonil (*)	C ₆ H ₆ Br ₂ N ₂	100mg	
Butanetriol				
CAS 42890-76-6 DRE-AY09010003PY DRE-SY09010003PY	MW 106.1204 Butanetriol 1000 µg/mL in Pyridine Butanetriol 1000 µg/mL in Pyridine	C ₄ H ₁₀ O ₃	5ml 5x5ml	
3-[N-n-Butyl-N-acetyl]aminopropionic acid ethyl ester				
CAS 52304-36-6 DRE-C10929010	MW 215.2893 3-[N-n-Butyl-N-acetyl]aminopropionic acid ethyl ester	C ₁₁ H ₂₁ NO ₃	100mg	
N-n-Butylbenzenesulfonamide				
CAS 3622-84-2 DRE-C10931120	MW 213.2966 N-n-Butylbenzenesulfonamide	C ₁₀ H ₁₃ NO ₂ S	250mg	
Butyltin Trichloride				
CAS 1118-46-3 DRE-CA10931700 DRE-GA09010356ME	MW 282.1833 n-Butyl-tin trichloride n-Butyltin Trichloride 1000 µg/mL in Methanol(‡)	C ₄ H ₉ Cl ₃ Sn	100mg 1ml	
Calflo E (Calcium Silicate)				
CAS 1344-95-2 DRE-E10935000	MW 116.1617 Calflo E	Ca·O ₃ Si	250g	

Additional organic contaminants

Product code	Description			
Camphene				
CAS 79-92-5 DRE-C10945000 DRE-GA09010073IP	MW 136.234 D,L-Camphene Camphene 1000 µg/mL in Isopropanol(±)	C ₁₀ H ₁₆	250mg 1ml	AND Enantiomer 
(1S)-(+)-3-Carene				
CAS 498-15-7 DRE-CA11042000 DRE-GA09010074IP	MW 136.234 (1S)-(+)-3-Carene (+)-3-Carene 1000 µg/mL in Isopropanol(±)	C ₁₀ H ₁₆	100mg 1ml	
3-Carene				
CAS 13466-78-9 DRE-C11041900	MW 136.234 3-Carene	C ₁₀ H ₁₆	1g	
Caryophyllene Oxide (β-Caryophyllene Epoxide)				
CAS 1139-30-6 DRE-GA09010046IP	MW 220.3505 Caryophyllene Oxide 1000 µg/mL in Isopropanol(±)	C ₁₅ H ₂₄ O	1ml	
Beta-Caryophyllene (Isocaryophyllene)				
CAS 118-65-0 DRE-GA09010043IP	MW 204.3511 Beta-Caryophyllene 1000 µg/mL in Isopropanol(±)	C ₁₅ H ₂₄	1ml	
Celestolide (ABDI)				
CAS 13171-00-1 DRE-LA10045800CY	MW 244.3719 ABDI (Celestolid) 10 µg/mL in Cyclohexane	C ₁₇ H ₂₄ O	1ml	
4-Chloroacetanilide				
CAS 539-03-7 DRE-C11348000	MW 169.6082 4-Chloroacetanilide	C ₈ H ₈ ClNO	500mg	
Chloroacetic Acid				
CAS 79-11-8 DRE-YA11348500MB	MW 94.497 Chloroacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ H ₃ ClO ₂	1ml	
Chloroacetic Acid Methyl Ester				
CAS 96-34-4 DRE-C11348700 DRE-YA11348700MB	MW 108.5236 Chloroacetic acid methyl ester Chloroacetic acid-methyl ester 1000 µg/mL in Methyl-tert. butyl ether	C ₃ H ₅ ClO ₂	250mg 1ml	

Additional organic contaminants


Product code	Description			
1-cis-3-Chloroallyl-3,5,7-triaza-1-azonia-adamantane Chloride				
CAS 51229-78-8	MW 251.1561	C ₉ H ₁₆ ClN ₄ ·Cl		
DRE-C11349900	1-cis-3-Chloroallyl-3,5,7-triaza-1-azonia-adamantane chloride	100mg		
2-Chlorobenzoic Acid				
CAS 118-91-2	MW 156.5664	C ₇ H ₅ ClO ₂		
DRE-C11390500	2-Chlorobenzoic acid	250mg		
3-Chlorobenzoic Acid				
CAS 535-80-8	MW 156.5664	C ₇ H ₅ ClO ₂		
DRE-C11390800	3-Chlorobenzoic acid	250mg		
4-Chlorobenzoic Acid				
CAS 74-11-3	MW 156.5664	C ₇ H ₅ ClO ₂		
DRE-C11391000	4-Chlorobenzoic acid	250mg		
2-Chlorobenzonitrile				
CAS 873-32-5	MW 137.5664	C ₇ H ₄ ClN		
DRE-C11392000	2-Chlorobenzonitrile	250mg		
2-Chloro-4-methylsulfonylbenzoic acid				
CAS 53250-83-2	MW 234.6568	C ₈ H ₇ ClO ₄ S		
DRE-C11447000	2-Chloro-4-methylsulfonylbenzoic acid	100mg		
Chloroparaffins				
DRE-LA11457510CY	Chloroparaffin C10 44.82% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457512CY	Chloroparaffin C10 50.18% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457514CY	Chloroparaffin C10 55.00% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457516CY	Chloroparaffin C10 60.09% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457518CY	Chloroparaffin C10 65.02% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457520CY	Chloroparaffin C11 45.50% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457522CY	Chloroparaffin C11 50.21% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457524CY	Chloroparaffin C11 55.20% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457526CY	Chloroparaffin C11 60.53% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457528CY	Chloroparaffin C11 65.25% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457530CY	Chloroparaffin C12 45.32% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457532CY	Chloroparaffin C12 50.18% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457534CY	Chloroparaffin C12 55.00% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457536CY	Chloroparaffin C12 65.08% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457538CY	Chloroparaffin C12 69.98% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457540CY	Chloroparaffin C13 44.90% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457542CY	Chloroparaffin C13 50.23% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457544CY	Chloroparaffin C13 55.03% Cl 10 µg/mL in Cyclohexane	1ml		
DRE-LA11457546CY	Chloroparaffin C13 59.98% Cl 10 µg/mL in Cyclohexane	1ml		

Additional organic contaminants

Product code	Description	
DRE-LA11457548CY	Chloroparaffin C13 65.18% Cl 10 µg/mL in Cyclohexane	1ml
DRE-X23227200CY	Chloroparaffin C22 72.1% Cl 100 µg/mL in Cyclohexane	10ml
DRE-X23105100CY	Chloroparaffin C10-C13 51.5% Cl 100 µg/mL in Cyclohexane(‡)	10ml
DRE-X23105500CY	Chloroparaffin C10-C13 55.5% Cl 100 µg/mL in Cyclohexane(‡)	10ml
DRE-Y23105500CY	Chloroparaffin C10-C13 55,5% Cl 1000 µg/mL in Cyclohexane	10ml
DRE-YS23105500CY	Chloroparaffin C10-C13 55,5% Cl 1000 µg/mL in Cyclohexane	2x10ml
DRE-X23106300CY	Chloroparaffin C10-C13 63% Cl 100 µg/mL in Cyclohexane(‡)	10ml
DRE-Y23106300CY	Chloroparaffin C10-C13 63% Cl 1000 µg/mL in Cyclohexane	10ml
DRE-YS23106300CY	Chloroparaffin C10-C13 63% Cl 1000 µg/mL in Cyclohexane	2x10ml
DRE-X23144200CY	Chloroparaffin C14-C17 42% Cl 100 µg/mL in Cyclohexane	10ml
DRE-X23145200CY	Chloroparaffin C14-C17 52% Cl 100 µg/mL in Cyclohexane(‡)	10ml
DRE-Y23145200CY	Chloroparaffin C14-C17 52% Cl 1000 µg/mL in Cyclohexane	10ml
DRE-YS23145200CY	Chloroparaffin C14-C17 52% Cl 1000 µg/mL in Cyclohexane	3x10ml
DRE-X23145700CY	Chloroparaffin C14-C17 57% Cl 100 µg/mL in Cyclohexane(‡)	10ml
DRE-Y23145700CY	Chloroparaffin C14-C17 57%Cl 1000 µg/mL in Cyclohexane	10ml
DRE-YS23145700CY	Chloroparaffin C14-C17 57%Cl 1000 µg/mL in Cyclohexane	2x10ml
DRE-X23183600CY	Chloroparaffin C18-C20 36 % Cl 100 µg/mL in Cyclohexane	10ml
DRE-X23184900CY	Chloroparaffin C18-C20 49% Cl 100 µg/mL in Cyclohexane	10ml


3-Chlorophenol Acetate

CAS 13031-39-5	MW 170.593	C ₈ H ₇ ClO ₂	
DRE-C11473100	3-Chlorophenol acetate		20mg




4-Chlorophenol Acetate

CAS 876-27-7	MW 170.593	C ₈ H ₇ ClO ₂	
DRE-C11473200	4-Chlorophenol acetate		20mg




2-Chlorophenoxyacetic Acid

CAS 614-61-9	MW 186.5924	C ₈ H ₇ ClO ₃	
DRE-C11479000	2-Chlorophenoxyacetic acid		250mg




3-Chlorophenoxyacetic Acid

CAS 588-32-9	MW 186.5924	C ₈ H ₇ ClO ₃	
DRE-C11479500	3-Chlorophenoxyacetic acid		250mg




2-(4-Chlorophenoxy)propionic Acid

CAS 3307-39-9	MW 200.619	C ₉ H ₉ ClO ₃	
DRE-C11487000	2-(4-Chlorophenoxy) propionic acid		250mg
DRE-XA11487000AL	2-(4-Chlorophenoxy) propionic acid 100 µg/mL in Acetonitrile		1ml












4-Chlorophenyl Phenyl Ether










CAS 7005-72-3	MW 204.6523	C ₁₂ H ₉ ClO	
DRE-L11489700ME	4-Chlorophenyl-phenyl ether 10 µg/mL in Methanol		10ml











Additional organic contaminants

Product code	Description			
4-Chlorophenylacetic Acid				
CAS 1878-66-6 DRE-C11489000	MW 170.593 4-Chlorophenyl acetic acid	$C_8H_7ClO_2$	500mg	
1-(2-Chlorophenyl)-2-thiourea				
CAS 5344-82-1 DRE-C11491500	MW 186.6619 1-(2-Chlorophenyl)-2-thiourea	$C_7H_7ClN_2S$	100mg	
2-Chloropyridine				
CAS 109-09-1 DRE-C11503700	MW 113.545 2-Chloropyridine	C_5H_4ClN	250mg	
4-Cholesten-3-one				
CAS 601-57-0 DRE-C11665300	MW 384.6377 4-Cholesten-3-one	$C_{27}H_{44}O$	100mg	
Cinnamaldehyde				
CAS 104-55-2 DRE-C11667480	MW 132.1592 Cinnamaldehyde	C_9H_8O	250mg	
Cinnamic Acid Benzyl Ester				
CAS 103-41-3 DRE-C11667490	MW 238.2812 Cinnamic acid-benzyl ester	$C_{16}H_{14}O_2$	250mg	
Citral				
CAS 5392-40-5 DRE-C11668510	MW 152.2334 Citral	$C_{10}H_{16}O$	250mg	
Citric Acid Monohydrate				
CAS 5949-29-1 DRE-C11668515	MW 210.1388 Citric acid monohydrate	$C_6H_8O_7 \cdot H_2O$	100mg	
(±)-β-Citronellol				
CAS 106-22-9 DRE-C11668525	MW 156.2652 (±)-β-Citronellol	$C_{10}H_{20}O$	250mg	










Additional organic contaminants

Product code	Description			
Creatine				
CAS 57-00-1 DRE-C11748500	MW 131.1332 Creatine	C ₄ H ₉ N ₃ O ₂	250mg	
Cyanuric Acid				
CAS 108-80-5 DRE-C11815000	MW 129.0742 Cyanuric acid	C ₃ H ₃ N ₃ O ₃	500mg	
Cyclohexylisocyanate				
CAS 3173-53-3 DRE-CA11831000	MW 125.1683 Cyclohexylisocyanate	C ₇ H ₁₁ NO	1g	
DL-Cysteine				
CAS 3374-22-9 DRE-C11925000	MW 121.1582 DL-Cysteine	C ₃ H ₇ NO ₂ S	100mg	
DMSA				
CAS 4710-17-2 DRE-C13030000	MW 200.2581 DMSA(‡)	C ₆ H ₁₂ N ₂ O ₂ S	100mg	
DPMI (Cashmeran)				
CAS 33704-61-9 DRE-LA13085000CY	MW 206.3239 DPMI (Cashmeran) 10 µg/mL in Cyclohexane	C ₁₄ H ₂₂ O	1ml	
Deanol (N,N-Dimethylethanolamine)				
CAS 108-01-0 DRE-CA12726800	MW 89.1362 N,N-Dimethylethanolamine	C ₄ H ₁₁ NO	1ml	
Decabromodiphenyl Ether (PBDE 209)				
CAS 1163-19-5 DRE-C15898209 DRE-XA15898209AC	MW 959.1678 PBDE 209 PBDE 209 (Decabromodiphenylether) 100 µg/mL in Acetone	C ₁₂ Br ₁₀ O	100mg 1ml	
Decafluorobiphenyl				
CAS 434-90-2 DRE-C12092500 DRE-YA12092500MB	MW 334.1124 Decafluorobiphenyl Decafluorobiphenyl 2000 µg/mL in Methyl-tert. butyl ether	C ₁₂ F ₁₀	100mg 1ml	










Additional organic contaminants

Product code	Description			
Decafluorotriphenylphosphine (DFTPP)				
CAS 5074-71-5 DRE-YA12093000AC	MW 442.1901 Decafluorotriphenylphosphine 1000 µg/mL in Acetone(±)	C ₁₈ H ₅ F ₁₀ P	1ml	
3,5-Diacetyl-1,4-dihydro-2,6-lutidine				
CAS 1079-95-4 DRE-CA12175500	MW 193.2423 3,5-Diacetyl-1,4-dihydro-2,6-lutidine	C ₁₁ H ₁₃ NO ₂	100mg	
trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic Acid Monohydrate				
CAS 125572-95-4 DRE-C12194000	MW 364.3484 trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic acid mono	C ₁₄ H ₂₂ N ₂ O ₈ ·H ₂ O	250mg	
1,2-Diaminopropane-N,N,N',N'-tetraacetic Acid				
CAS 4408-81-5 DRE-C12197000	MW 306.2692 1,2-Diaminopropane-N,N,N',N'-tetraacetic acid	C ₁₁ H ₁₈ N ₂ O ₈	100mg	
Diaryl-parl-phenylene-diamines				
CAS 68953-84-4 DRE-C12205000	MW n/a Diaryl-parl-phenylene-diamines		100mg	No Structure
Dibromoacetic Acid				
CAS 631-64-1 DRE-C12216000 DRE-YA12216000MB	MW 217.8441 Dibromoacetic acid Dibromoacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ H ₂ Br ₂ O ₂	250mg 1ml	
Dibromochloroacetic acid-methyl Ester				
CAS 20428-75-5 DRE-CA12231000	MW 266.3157 Dibromochloroacetic acid-methyl ester	C ₃ H ₃ Br ₂ ClO ₂	100mg	
4,4'-Dibromooctafluorobiphenyl				
CAS 10386-84-2 DRE-C12240600 DRE-XA12240600CY DRE-YA12240600MB	MW 455.9236 4,4'-Dibromooctafluorobiphenyl 4,4'-Dibromooctafluorobiphenyl 100 µg/mL in Cyclohexane 4,4'-Dibromooctafluorobiphenyl 2000 µg/mL in Methyl-tert. butyl ether	C ₁₂ Br ₂ F ₈	100mg 1ml 1ml	
Dibutyldiphenyltin				
CAS 6452-61-5 DRE-C12251200	MW 387.1463 Dibutyldiphenyltin	C ₂₀ H ₂₈ Sn	100mg	










Additional organic contaminants

Product code	Description			
2,6-Di-tert-butyl-4-hydroxymethylphenol				
CAS 88-26-6 DRE-C12253150	MW 236.3499 2,6-Di-tert-butyl-4-hydroxymethylphenol	C ₁₅ H ₂₄ O ₂	100mg	
2-(3,5-Di-tert-butyl-2-hydroxyphenyl)-2H-benzotriazole				
CAS 3846-71-7 DRE-C12253300	MW 323.432 2-(3,5-Di-tert-butyl-2-hydroxyphenyl)-2H-benzotriazole	C ₂₀ H ₂₅ N ₃ O	100mg	
Dibutyltin Dichloride				
CAS 683-18-1 DRE-C12258000 DRE-L12258000CY DRE-GA09010354ME	MW 303.8445 Dibutyltin-dichloride Dibutyltin-dichloride 10 µg/mL in Cyclohexane Di-n-butyltin Dichloride 1000 µg/mL in Methanol(±)	C ₈ H ₁₈ Cl ₂ Sn	250mg 10ml 1ml	
Dibutyltin Oxide				
CAS 818-08-6 DRE-C12258300	MW 248.9379 Dibutyltin-oxide	C ₈ H ₁₈ OSn	250mg	
Dichloroacetic Acid				
CAS 79-43-6 DRE-C12320000 DRE-YA12320000MB	MW 128.9421 Dichloroacetic acid Dichloroacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ H ₂ Cl ₂ O ₂	1g 1ml	
Dichloroacetic Acid Methyl Ester				
CAS 116-54-1 DRE-C12320200 DRE-YA12320200MB	MW 142.9687 Dichloroacetic acid methyl ester Dichloroacetic acid-methyl ester 1000 µg/mL in Methyl-tert. butyl ether	C ₃ H ₄ Cl ₂ O ₂	250mg 1ml	
4-(Dichloroacetyl)-1-oxa-4-azaspiro[4,5]decane				
CAS 71526-07-3 DRE-C12321500	MW 252.1376 4-(Dichloroacetyl)-1-oxa-4-azaspiro[4,5]decane	C ₁₀ H ₁₅ Cl ₂ NO ₂	100mg	
2,4-Dichlorobenzoic Acid				
CAS 50-84-0 DRE-C12400000	MW 191.0115 2,4-Dichlorobenzoic acid	C ₇ H ₄ Cl ₂ O ₂	250mg	
2-(4,6-Dichloro-2-methylphenoxy)propionic Acid				
CAS 20021-12-9 DRE-C12427500	MW 249.0906 2-(4,6-Dichloro-2-methyl-phenoxy)-propionic acid	C ₁₀ H ₉ Cl ₂ O ₃	10mg	










Additional organic contaminants

Product code	Description			
2,4-Dichlorophenyl Acetic Acid (2,4-DCAA)				
CAS 19719-28-9	MW 205.038	C ₈ H ₆ Cl ₂ O ₂		
DRE-C12468000	2,4-Dichlorophenyl acetic acid		100mg	
DRE-XA12468000MB	2,4-Dichlorophenyl acetic acid 100 µg/mL in Methyl-tert. butyl ether		1ml	
DRE-YS09010016AC	2,4-Dichlorophenylacetic Acid (2,4-DCAA) 1000 µg/mL in Acetone(‡)		5x1ml	
2,4-Dichlorophenylacetic Acid Methyl Ester (DCAA Methyl Ester)				
CAS 55954-23-9	MW 219.0646	C ₉ H ₈ Cl ₂ O ₂		
DRE-GS09010148MB	DCAA Methyl Ester 100 µg/mL in Methyl tert-Butyl Ether(‡)		5x1ml	
2,3-Dichloropropionanilide				
CAS 17839-22-4	MW 218.0799	C ₉ H ₉ Cl ₂ NO		
DRE-C12500000	2,3-Dichloropropionanilide		100mg	
Diethylaminohydroxybenzoyl hexyl benzoate				
CAS 302776-68-7	MW 397.5072	C ₂₄ H ₃₁ NO ₄		
DRE-C12604700	Diethylaminohydroxybenzoyl hexyl benzoate(‡)		100mg	
N,N-Diethyldithiocarbamate Sodium Salt Trihydrate				
CAS 20624-25-3	MW 225.3052	C ₆ H ₁₀ NS ₂ ·Na·3H ₂ O		
DRE-C12605750	N,N-Diethyldithiocarbamate sodium salt trihydrate		100mg	
Diethyleneglycoldinitrate (DEGN)				
CAS 693-21-0	MW 196.1155	C ₄ H ₈ N ₂ O ₇		
DRE-LA12605800AL	Diethyleneglycoldinitrate (DEGN) 10 µg/mL in Acetonitrile		1ml	
Diethylenetriaminepentaacetic Acid (DTPA)				
CAS 67-43-6	MW 393.3465	C ₁₄ H ₂₃ N ₃ O ₁₀		
DRE-E13095000	DTPA		250mg	
N,N-Diethyl-p-phenyldiamine sulfate salt				
CAS 6283-63-2	MW 262.3259	C ₁₀ H ₁₆ N ₂ ·H ₂ O ₄ S		
DRE-C12606680	N,N-Diethyl-p-phenyldiamine sulfate salt		100mg	
2,2'-Difluorobiphenyl (PFB 4)				
CAS 388-82-9	MW 190.1887	C ₁₂ H ₈ F ₂		
DRE-YA12632500MB	2,2'-Difluorobiphenyl (PFB 4) 2000 µg/mL in Methyl-tert. butyl ether		1ml	










Additional organic contaminants

Product code	Description			
4,4'-Difluorobiphenyl				
CAS 398-23-2	MW 190.1887	C ₁₂ H ₈ F ₂		
DRE-C12632015	4,4'-Difluorobiphenyl		100mg	
DRE-YA12632015AC	4,4'-Difluorobiphenyl (PFB 15) 2000 µg/mL in Acetone		1ml	
Di-n-heptyltin-dichloride				
CAS 74340-12-8	MW 388.004	C ₁₄ H ₃₀ Cl ₂ Sn		
DRE-C12634100	Di-n-heptyltin-dichloride		100mg	
DRE-L12634100MB	Di-n-heptyltin-dichloride 10 µg/mL in Methyl-tert. butyl ether (*)		10ml	
Dihydrosafrol				
CAS 94-58-6	MW 164.2011	C ₁₀ H ₁₂ O ₂		
DRE-L12635000ME	Dihydrosafrol 10 µg/mL in Methanol		10ml	
4,4'-Dihydroxybiphenyl				
CAS 92-88-6	MW 186.2066	C ₁₂ H ₁₀ O ₂		
DRE-C12635500	4,4'-Dihydroxybiphenyl		250mg	
Dimethyl Disulfide				
CAS 624-92-0	MW 94.199	C ₂ H ₆ S ₂		
DRE-CA12726480	Dimethyl disulfide		1ml	
Dimethyl Terephthalate (Terephthalic acid-bis-methyl ester)				
CAS 120-61-6	MW 194.184	C ₁₀ H ₁₀ O ₄		
DRE-C17322000	Terephthalic acid-bis-methyl ester		250mg	
Dimethyl trisulfide				
CAS 3658-80-8	MW 126.264	C ₂ H ₆ S ₃		
DRE-CA12755000	Dimethyl trisulfide		1ml	
1,3-Dimethyl Urea				
CAS 96-31-1	MW 88.1084	C ₃ H ₈ N ₂ O		
DRE-C12760000	1,3-Dimethyl urea		250mg	
N,N-Dimethyldithiocarbamate Nickel Salt				
CAS 15521-65-0	MW 299.1263	C ₆ H ₁₂ N ₂ NiS ₄		
DRE-C12726495	N,N-Dimethyldithiocarbamate nickel salt		100mg	










Additional organic contaminants

Product code	Description			
N,N-Dimethyldithiocarbamate Sodium Salt Hydrate				
CAS 207233-95-2 DRE-C12726500	MW 161.2215 N,N-Dimethyldithiocarbamate sodium salt hydrate	$C_3H_6NS_2 \cdot Na \cdot H_2O$	250mg	
O,O-Dimethylphosphorodithioic acid ammonium salt				
CAS 1066-97-3 DRE-C12738500	MW 175.2101 O,O-Dimethylphosphorodithioic acid ammonium salt	$C_2H_6O_2PS_2 \cdot H_4N^+$	10mg	
O,O-Dimethylphosphorothioic acid sodium salt				
CAS 23754-87-2 DRE-C12738600	MW 164.0958 O,O-Dimethylphosphorothioic acid sodium salt	$C_2H_6O_3PS \cdot Na$	10mg	
O,O-Dimethylphosphorothioic acid sodium salt D6				
CAS n/a DRE-XA12738610AC	MW 170.1327 O,O-Dimethylphosphorothioic acid sodium salt D6 100 µg/mL in Acetone (*)	$C_2^2H_6O_3PS \cdot Na$	1.1ml	
2,2-Dimethylpropionic Acid				
CAS 75-98-9 DRE-C12740000	MW 102.1317 2,2-Dimethylpropionic acid	$C_5H_{10}O_2$	250mg	
N,N-Dimethylpyridin-4-amine (4-Dimethylaminopyridine)				
CAS 1122-58-3 DRE-C12723300	MW 122.1677 4-Dimethylaminopyridine	$C_7H_{10}N_2$	100mg	
N,N-Dimethylsulfamide				
CAS 3984-14-3 DRE-C12743000	MW 124.1621 N,N-Dimethylsulfamide	$C_2H_6N_2O_2S$	100mg	
Dimethyltin Dichloride				
CAS 753-73-1 DRE-C12750000 DRE-GA09010363ME	MW 219.685 Dimethyltin-dichloride Dimethyltin Dichloride 1000 µg/mL in Methanol(‡)	$C_2H_6Cl_2Sn$	250mg 1ml	
Di-n-Octyltin Dichloride				
CAS 3542-36-7 DRE-C12837000 DRE-GA09010360ME	MW 416.0572 Di-n-Octyltin-dichloride Di-n-octyltin Dichloride 1000 µg/mL in Methanol(‡)	$C_{16}H_{34}Cl_2Sn$	100mg 1ml	










Additional organic contaminants

Product code	Description			
Diphenylmercury				
CAS 587-85-9 DRE-C12900000	MW 354.7978 Diphenylmercury	C ₁₂ H ₁₀ Hg	250mg	
Diphenylmethan-4,4'-diisocyanate				
CAS 101-68-8 DRE-CA12905000	MW 250.2521 Diphenylmethan-4,4'-diisocyanate	C ₁₅ H ₁₀ N ₂ O ₂	250mg	
Diphenylmethanol				
CAS 91-01-0 DRE-C12906000	MW 184.2338 Diphenylmethanol	C ₁₃ H ₁₂ O	250mg	
Diphenyltin Dichloride				
CAS 1135-99-5 DRE-C12921000 DRE-GA09010357ME	MW 343.8238 Diphenyl-tin-dichloride Diphenyltin Dichloride 1000 µg/mL in Methanol(±)	C ₁₂ H ₁₀ Cl ₂ Sn	100mg 1ml	
N,N'-Diphenylurea (Carbanalide)				
CAS 102-07-8 DRE-C12921500 DRE-XA09010162AL	MW 212.2472 N,N'-Diphenylurea Carbanalide 100 µg/mL in Acetonitrile(±)	C ₁₃ H ₁₂ N ₂ O	250mg 1ml	
Di-n-Propyltin Dichloride				
CAS 867-36-7 DRE-L12950000IO DRE-GA09010361ME	MW 275.7914 Di-n-propyl-tin-dichloride 10 µg/mL in Isooctane Di-n-propyltin-dichloride 1000 µg/mL in Methanol(±)	C ₆ H ₁₄ Cl ₂ Sn	10ml 1ml	
N,N'-Divinylethyleneurea				
CAS 13811-50-2 DRE-C13025000	MW 138.1671 N,N'-Divinylethyleneurea	C ₇ H ₁₀ N ₂ O	100mg	
EDTA-tetrabutyl ester				
CAS 14531-12-5 DRE-C13111000	MW 516.6679 EDTA-tetrabutyl ester	C ₂₆ H ₄₈ N ₂ O ₈	50mg	
Edetic Acid (EDTA)				
CAS 60-00-4 DRE-E13110800	MW 292.2426 EDTA	C ₁₀ H ₁₆ N ₂ O ₈	250mg	

Additional organic contaminants

Product code	Description			
Ethoxyacetic Acid Ethyl Ester				
CAS 817-95-8 DRE-C13307000	MW 132.1577 Ethoxyacetic acid-ethyl ester	$C_6H_{12}O_3$	100mg	
2-Ethylbutyric Acid				
CAS 88-09-5 DRE-C13321500	MW 116.1583 2-Ethylbutyric acid	$C_8H_{16}O_2$	250mg	
4-(1-Ethyl-1,4-dimethylpentyl)phenol				
CAS 142731-63-3 DRE-C13325500	MW 220.3505 4-(1-Ethyl-1,4-dimethylpentyl)-phenol	$C_{18}H_{24}O$	10mg	
Ethylene Oxide				
CAS 75-21-8 DRE-GA09010401TN DRE-GS09010401TN	MW 44.0526 Ethylene Oxide Solution, 1000 µg/mL in Triacetin(‡) Ethylene Oxide Solution, 1000 µg/mL in Triacetin(‡)	C_2H_4O	1ml 5x1ml	
Ethylenediamine-N,N'-diacetic Acid				
CAS 5657-17-0 DRE-C13326700	MW 176.1705 Ethylenediamine-N,N'-diacetic acid	$C_6H_{12}N_2O_4$	100mg	
2-Ethylhexylnitrate				
CAS 27247-96-7 DRE-C13342330	MW 175.2255 2-Ethylhexylnitrate	$C_8H_{17}NO_3$	250mg	
Ethylmercury Chloride				
CAS 107-27-7 DRE-C13345000	MW 265.1041 Ethylmercury-chloride	C_2H_5ClHg	100mg	
Ethylmethanesulfonate				
CAS 62-50-0 DRE-C13346500	MW 124.1588 Ethylmethanesulfonate	$C_3H_8O_3S$	100mg	
4-(2-Ethyl-6-methylphenyl)-5-methyl-3-morpholinone				
CAS 120375-14-6 DRE-C13349000	MW 233.3062 4-(2-Ethyl-6-methylphenyl)-5-methyl-3-morpholinone	$C_{14}H_{19}NO_2$	10mg	










Additional organic contaminants

Product code	Description			
N-Ethyl-N-phenylurethane				
CAS 1013-75-8 DRE-C13352500	MW 193.2423 N-Ethyl-N-phenylurethane	C ₁₁ H ₁₅ NO ₂	100mg	
Etrofol (CPMC)				
CAS 3942-54-9 DRE-C13390000	MW 185.6076 Etrofol (CPMC)	C ₈ H ₉ ClNO ₂	10mg	
Eucalyptol (Cineole)				
CAS 470-82-6 DRE-GA09010075IP	MW 154.2493 Eucalyptol 1000 µg/mL in p/t Methanol(±)	C ₁₀ H ₁₈ O	1ml	
Eugenol				
CAS 97-53-0 DRE-C13395000	MW 164.2011 Eugenol	C ₁₀ H ₁₂ O ₂	100mg	
Fluoroacetic Acid Sodium Salt				
CAS 62-74-8 DRE-L13772000WA	MW 100.0243 Fluoroacetic acid sodium salt 10 µg/mL in Water (*)	C ₂ H ₂ FO ₂ ·Na	10ml	
2-Fluorobiphenyl				
CAS 321-60-8 DRE-C13785000 DRE-YA13785000CY	MW 172.1983 2-Fluorobiphenyl 2-Fluorobiphenyl 2000 µg/mL in Cyclohexane	C ₁₂ H ₉ F	100mg 1ml	
Formaldehyde-2,4-dinitrophenylhydrazone (DNPH)				
CAS 1081-15-8 DRE-C13909200	MW 210.1469 Formaldehyde-2,4-dinitrophenylhydrazone (DNPH)	C ₇ H ₆ N ₄ O ₄	100mg	
2,5-Furandicarboxylic acid				
CAS 3238-40-2 DRE-C13966000	MW 156.093 2,5-Furandicarboxylic acid	C ₆ H ₄ O ₅	100mg	
(±)-Geosmin				
CAS 16423-19-1 DRE-C14005000 DRE-LA14005000ME DRE-XA14005000ME	MW 182.3025 (±)-Geosmin Geosmin 10 µg/mL in Methanol Geosmin 100 µg/mL in Methanol	C ₁₂ H ₂₂ O	10mg 1ml 1ml	










Additional organic contaminants

Product code	Description			
DL-Glutamic Acid				
CAS 617-65-2 DRE-C14034400	MW 147.1293 DL-Glutamic acid	C ₅ H ₉ NO ₄	100mg	
Guanidine Hydrochloride				
CAS 50-01-1 DRE-C14057000	MW 95.5314 Guanidine hydrochloride	CH ₅ N ₃ ·ClH	250mg	
HHCB (Galaxolide)				
CAS 1222-05-5 DRE-C14213000 DRE-LA14213000CY	MW 258.3984 HHCB (Galaxolide) HHCB (Galaxolide) 10 µg/mL in Cyclohexane	C ₁₈ H ₂₆ O	20mg 1ml	
1,2,4,5,6,9,10-Heptachlorodecane (CP-7)				
CAS 890302-88-2 DRE-LA14110000CY	MW 383.3971 1,2,4,5,6,9,10-Heptachlorodecane (CP-7) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₅ Cl ₇	1ml	
1,2,5,5,6,9,10-Heptachlorodecane (CP-8)				
CAS 1278446-65-3 DRE-LA14110010CY	MW 383.3971 1,2,5,5,6,9,10-Heptachlorodecane (CP-8) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₅ Cl ₇	1ml	
4-Heptylphenol				
CAS 1987-50-4 DRE-C14136500	MW 192.2973 4-Heptylphenol	C ₁₃ H ₂₀ O	100mg	
n-Heptyltin-trichloride				
CAS 59344-47-7 DRE-L14137000MB DRE-XA14137000MB	MW 324.263 n-Heptyltin-trichloride 10 µg/mL in Methyl-tert. butyl ether (*) n-Heptyltin-trichloride 100 µg/mL in Methyl-tert. butyl ether (*)	C ₇ H ₁₅ Cl ₃ Sn	10ml 1ml	
1,2,5,6,9,10-Hexabromocyclododecane				
CAS 3194-55-6 DRE-C14152000	MW 641.6953 1,2,5,6,9,10-Hexabromocyclododecane	C ₁₂ H ₁₈ Br ₆	250mg	
Hexachlorocyclopentadiene				
CAS 77-47-4 DRE-C14171000 DRE-L14171000IO DRE-XA14171000IO	MW 272.7715 Hexachlorocyclopentadiene Hexachlorocyclopentadiene 10 µg/mL in Isooctane Hexachlorocyclopentadiene 100 µg/mL in Isooctane	C ₅ Cl ₆	100mg 10ml 1ml	










Additional organic contaminants

Product code	Description			
1,2,4,5,9,10-Hexachlorodecane (CP-6)				
CAS 890302-87-1 DRE-LA14171400CY	MW 348.952 1,2,4,5,9,10-Hexachlorodecane (CP-6) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₆ Cl ₆	1ml	
1,2,5,6,9,10-Hexachlorodecane (CP-4)				
CAS 189350-94-5 DRE-LA14171500CY	MW 348.952 1,2,5,6,9,10-Hexachlorodecane (CP-4) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₆ Cl ₆	1ml	
1,4,5,6,7,7-Hexachloro-N-(ethylmercury)-5-norbornene-2.3-dicarboximide				
CAS 2597-93-5 DRE-C14175000	MW 598.4868 1,4,5,6,7,7-Hexachloro-N-(ethylmercury)-5-norbornene-2.3-dicarboximide	C ₁₁ H ₇ Cl ₆ HgNO ₂	100mg	
Hexamethylene-1,6-diisocyanate				
CAS 822-06-0 DRE-CA14194500	MW 168.1931 Hexamethylenediisocyanate	C ₈ H ₁₂ N ₂ O ₂	250mg	
Hexamethylphosphoramide (Hexametapol)				
CAS 680-31-9 DRE-C14194700	MW 179.2004 Hexamethylphosphoramide (Hexametapol)	C ₆ H ₁₈ N ₃ OP	5ml	
cis-3-Hexen-1-ol				
CAS 928-96-1 DRE-CA14202400	MW 100.1589 cis-3-Hexen-1-ol	C ₆ H ₁₂ O	100mg	
cis-3-Hexenyl acetate				
CAS 3681-71-8 DRE-CA14202500	MW 142.1956 cis-3-Hexenyl acetate	C ₈ H ₁₄ O ₂	100mg	
Hexylcinnamal				
CAS 101-86-0 DRE-C14208000	MW 216.3187 Hexylcinnamal	C ₁₅ H ₂₀ O	250mg	
DL-Histidine				
CAS 4998-57-6 DRE-C14213100	MW 155.1546 DL-Histidine	C ₆ H ₉ N ₃ O ₂	100mg	










Additional organic contaminants

Product code	Description			
Homosalate (3,3,5-Trimethylcyclohexyl salicylate)				
CAS 118-56-9 DRE-C17881300	MW 262.3441 3,3,5-Trimethylcyclohexyl salicylate	C ₁₆ H ₂₂ O ₃	250mg	
7-Hydroxycitronellal (7-Hydroxy-3,7-dimethyloctanal)				
CAS 107-75-5 DRE-CA14230300	MW 172.2646 7-Hydroxycitronellal	C ₁₀ H ₂₀ O ₂	250mg	
4-Hydroxycoumarin				
CAS 1076-38-6 DRE-C14230600	MW 162.1421 4-Hydroxycoumarin	C ₉ H ₆ O ₃	250mg	
7-Hydroxycoumarin				
CAS 93-35-6 DRE-C14230700	MW 162.1421 7-Hydroxycoumarin	C ₉ H ₆ O ₃	250mg	
N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide				
CAS 24448-09-7 DRE-C14231570	MW 557.224 N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide	C ₁₁ H ₈ F ₁₇ NO ₃ S	25mg	
Hydroxyisohexyl-3-cyclohexene carboxaldehyde				
CAS 31906-04-4 DRE-C14232000	MW 210.3126 Hydroxyisohexyl-3-cyclohexene carboxaldehyde	C ₁₃ H ₂₂ O ₂	100mg	
2-Hydroxymethylfuran (2-Furfuryl alcohol)				
CAS 98-00-0 DRE-C13972300	MW 98.0999 2-Furfuryl alcohol	C ₅ H ₆ O ₂	250mg	
2-Hydroxy-2-methylpropiophenone				
CAS 7473-98-5 DRE-C14233500	MW 164.2011 2-Hydroxy-2-methylpropiophenone	C ₁₀ H ₁₂ O ₂	250mg	
2-Hydroxy-1,4-naphthoquinone				
CAS 83-72-7 DRE-C14233800	MW 174.1528 2-Hydroxy-1,4-naphthoquinone	C ₁₀ H ₆ O ₃	100mg	










Additional organic contaminants

Product code	Description			
2-Hydroxy-5-(trifluoromethyl)pyridine				
CAS 33252-63-0 DRE-C14253000	MW 163.0973 2-Hydroxy-5-(trifluoromethyl)pyridine	C ₆ H ₄ F ₃ NO	100mg	
2-Imidazolidinone				
CAS 120-93-4 DRE-C14284000	MW 86.0925 2-Imidazolidinone	C ₃ H ₆ N ₂ O	250mg	
Imino bis(methylphosphonic Acid)				
CAS 17261-34-6 DRE-C14284800	MW 205.0435 Imino-bis(methylphosphonic acid)	C ₂ H ₉ NO ₆ P ₂	100mg	
Iminodiacetic Acid				
CAS 142-73-4 DRE-C14285000	MW 133.1027 Iminodiacetic acid	C ₄ H ₇ NO ₄	250mg	
3-Indolylacetonitrile				
CAS 771-51-7 DRE-C14300000	MW 156.1839 3-Indolyl-acetonitrile	C ₁₀ H ₈ N ₂	250mg	
Isobenzofuranone				
CAS 87-41-2 DRE-C14382000	MW 134.132 Isobenzofuranone	C ₈ H ₆ O ₂	100mg	
Isobornyl Acetate				
CAS 125-12-2 DRE-C14385000	MW 196.286 Isobornyl acetate	C ₁₂ H ₂₀ O ₂	250mg	
2-Isobutyl-3-methoxypyrazine				
CAS 24683-00-9 DRE-C14394700 DRE-XA14394700ME	MW 166.2203 2-Isobutyl-3-methoxypyrazine 2-Isobutyl-3-methoxypyrazine 100 µg/mL in Methanol	C ₉ H ₁₄ N ₂ O	100mg 1ml	
Isobutyric Acid				
CAS 79-31-2 DRE-C14395500	MW 88.1051 Isobutyric acid	C ₄ H ₈ O ₂	250mg	










Additional organic contaminants

Product code	Description			
Isobutyric Acid Methyl Ester				
CAS 547-63-7 DRE-C14396000	MW 102.1317 Isobutyric acid methyl ester	C ₅ H ₁₀ O ₂	250mg	
DL-Isoleucine				
CAS 443-79-8 DRE-C14429000	MW 131.1729 DL-Isoleucine	C ₆ H ₁₃ NO ₂	100mg	
α-Isomethylionone				
CAS 127-51-5 DRE-C14433000	MW 206.3239 α-Isomethylionone(‡)	C ₁₄ H ₂₂ O	100mg	
Isononylphenol-ethoxylate (technical mixture)				
CAS 37205-87-1 DRE-C14438000	MW 248.4036 Isononylphenol-ethoxylate (technical mixture)	C ₁₇ H ₂₆ O	100mg	
4-iso-Octylphenol-di-ethoxylate				
CAS 2315-61-9 DRE-C15712802 DRE-LA15712802AC	MW 294.429 4-iso-Octylphenol-di-ethoxylate 4-iso-Octylphenol-di-ethoxylate 10 µg/mL in Acetone	C ₁₈ H ₃₀ O ₃	10mg 1ml	
Isopropyl Myristate (Myristic acid isopropyl ester)				
CAS 110-27-0 DRE-C15391990	MW 270.4507 Myristic acid-isopropyl ester	C ₁₇ H ₃₄ O ₂	250mg	
2-Isopropyl-3-methoxypyrazine				
CAS 25773-40-4 DRE-C14463700 DRE-XA14463700ME	MW 152.1937 2-Isopropyl-3-methoxypyrazine 2-Isopropyl-3-methoxypyrazine 100 µg/mL in Methanol	C ₈ H ₁₂ N ₂ O	100mg 1ml	
Isopropylnitrate				
CAS 1712-64-7 DRE-C14463850	MW 105.0926 Isopropylnitrate	C ₃ H ₇ NO ₃	250mg	
2-Isopropylthioxantone (ITX)				
CAS 5495-84-1 DRE-C14465200 DRE-XA14465200AL	MW 254.3468 2-Isopropylthioxantone ITX 2-Isopropylthioxantone ITX 100 µg/mL in Acetonitrile	C ₁₆ H ₁₄ OS	100mg 1ml	










Additional organic contaminants

Product code	Description			
Isovaleric Acid				
CAS 503-74-2 DRE-C14479470	MW 102.1317 Isovaleric acid	C ₅ H ₁₀ O ₂	250mg	
2-Isovaleryl-1.3-indanedione				
CAS 83-28-3 DRE-C14479600	MW 230.2592 2-Isovaleryl-1.3-indanedione	C ₁₄ H ₁₄ O ₃	20mg	
Kinetin (6-Furfurylaminopurine)				
CAS 525-79-1 DRE-C13975000	MW 215.2114 6-Furfurylaminopurine	C ₁₀ H ₉ N ₅ O	100mg	
DL-Leucine				
CAS 328-39-2 DRE-C14629300	MW 131.1729 DL-Leucine	C ₆ H ₁₃ NO ₂	100mg	
Levomenthol ((-)-Menthol)				
CAS 2216-51-5 DRE-C14866020	MW 156.2652 (-)-Menthol	C ₁₀ H ₂₀ O	250mg	
DL-Lysine				
CAS 70-54-2 DRE-C14655000	MW 146.1876 DL-Lysine	C ₆ H ₁₄ N ₂ O ₂	100mg	
Maleic Acid				
CAS 110-16-7 DRE-C14726000	MW 116.0722 Maleic acid	C ₄ H ₄ O ₄	250mg	
Maleic Acid Dimethyl Ester				
CAS 624-48-6 DRE-C14727500	MW 144.1253 Maleic acid dimethyl ester	C ₆ H ₈ O ₄	250mg	
DL-Malic acid				
CAS 6915-15-7 DRE-C14730500	MW 134.0874 DL-Malic acid(‡)	C ₄ H ₆ O ₅	100mg	










Additional organic contaminants

Product code	Description			
Maltol				
CAS 118-71-8 DRE-C14734300	MW 126.11 Maltol	C ₆ H ₆ O ₃	100mg	
Melamine				
CAS 108-78-1 DRE-C14861400	MW 126.1199 Melamine(‡)	C ₃ H ₆ N ₆	250mg	
Melamine Triamine-15N3				
CAS 287476-11-3 DRE-L14861401AL	MW 129.1002 Melamine triamine-15N3 10 µg/mL in Acetonitrile (*)	C ₃ H ₆ ¹⁵ N ₃	10ml	
Menthol				
CAS 1490-04-6 DRE-C14866000	MW 156.2652 Menthol	C ₁₀ H ₂₀ O	250mg	
Mesalazine (5-Aminosalicylic acid)				
CAS 89-57-6 DRE-C10227020	MW 153.1354 5-Aminosalicylic acid	C ₇ H ₇ NO ₃	100mg	
Methacrylamide				
CAS 79-39-0 DRE-C14971500	MW 85.1045 Methacrylamide	C ₄ H ₇ NO	250mg	
dl-Methionine				
CAS 59-51-8 DRE-C15021000	MW 149.2113 DL-Methionine	C ₅ H ₁₁ NO ₂ S	100mg	
2-Methoxy-3,5-dimethylpyrazine				
CAS 92508-08-2 DRE-C15069000	MW 138.1671 2-Methoxy-3,5-dimethylpyrazine	C ₇ H ₁₀ N ₂ O	50mg	
3-Methoxy-2,5-dimethylpyrazine				
CAS 19846-22-1 DRE-C15069100	MW 138.1671 3-Methoxy-2,5-dimethylpyrazine	C ₇ H ₁₀ N ₂ O	10mg	










Additional organic contaminants

Product code	Description			
4-Methoxyphenylacetone (4'-Methoxyacetophenone)				
CAS 100-06-1 DRE-C15059000	MW 150.1745 4'-Methoxyacetophenone	C ₉ H ₁₀ O ₂	100mg	
Methyl Anthranilate (2-Aminobenzoic acid methyl ester)				
CAS 134-20-3 DRE-C10171470	MW 151.1626 2-Aminobenzoic acid methyl ester	C ₉ H ₉ NO ₂	100mg	
Methyl Salicylate				
CAS 119-36-8 DRE-C15143400	MW 152.1473 Methyl salicylate	C ₉ H ₈ O ₃	100mg	
N-Methylacetamide				
CAS 79-16-3 DRE-C15083300	MW 73.0938 N-Methylacetamide	C ₃ H ₇ NO	250mg	
N-Methylbenzenesulfonamide				
CAS 5183-78-8 DRE-C15083780	MW 171.2169 N-Methylbenzenesulfonamide	C ₇ H ₉ NO ₂ S	250mg	
5-Methyl-1H-benzotriazole				
CAS 136-85-6 DRE-C15083795	MW 133.1506 5-Methyl-1H-benzotriazole	C ₇ H ₇ N ₃	100mg	
3-(4'-Methyl)benzylidene-bornan-2-one				
CAS 36861-47-9 DRE-C15083800	MW 254.3667 3-(4'-Methyl)benzylidene-bornan-2-one	C ₁₈ H ₂₂ O	250mg	
2-Methylbutyric Acid				
CAS 116-53-0 DRE-C15084820	MW 102.1317 2-Methylbutyric acid	C ₅ H ₁₀ O ₂	250mg	
(Methylcyclopentadienyl)manganese Tricarbonyl				
CAS 12108-13-3 DRE-C15085036	MW 218.0881 (Methylcyclopentadienyl)manganese-tricarbonyl	C ₆ H ₇ ·C ₃ H ₃ MnO ₃	100mg	










Additional organic contaminants

Product code	Description			
N-Methyldiethanolamine				
CAS 105-59-9	MW 119.1622	C ₅ H ₁₃ NO ₂		
DRE-YA09010091OH	N-Methyldiethanolamine 1000 µg/mL in 15 mM Ammonium Hydroxide			1ml
				
Methyleugenol (4-Allyl-1,2-dimethoxybenzene)				
CAS 93-15-2	MW 178.2277	C ₁₁ H ₁₄ O ₂		
DRE-C13395200	Eugenol methyl ether			100mg
				
5-Methyl-2-furfural				
CAS 620-02-0	MW 110.1106	C ₆ H ₆ O ₂		
DRE-C15087700	5-Methyl-2-furfural			250mg
				
2-Methylisoborneol				
CAS 2371-42-8	MW 168.2759	C ₁₁ H ₂₀ O		
DRE-LA15088400ME	2-Methylisoborneol 10 µg/mL in Methanol			1ml
DRE-XA15088400ME	2-Methylisoborneol 100 µg/mL in Methanol			1ml
				
Methyl-mercury-chloride				
CAS 115-09-3	MW 251.0775	CH ₃ ClHg		
DRE-C15100000	Methyl-mercury-chloride			100mg
				
Methylmethanesulfonate				
CAS 66-27-3	MW 110.1322	C ₂ H ₆ O ₃ S		
DRE-C15100700	Methylmethanesulfonate			100mg
				
1-Methyl-3-nitroguanidine				
CAS 4245-76-5	MW 118.0946	C ₂ H ₅ N ₄ O ₂		
DRE-C15106500	1-Methyl-3-nitroguanidine			100mg
				
2-Methyl-4-nitroimidazole (2-Methyl-5-nitroimidazole)				
CAS 696-23-1	MW 127.1014	C ₅ H ₅ N ₃ O ₂		
DRE-C15107000	2-Methyl-5-nitroimidazole			100mg
				
Methyl-pentachlorophenyl sulfide				
CAS 1825-19-0	MW 296.4287	C ₇ H ₃ Cl ₅ S		
DRE-C15120000	Methyl-pentachlorophenyl sulfide			100mg
DRE-L15120000CY	Methyl-pentachlorophenyl sulfide 10 µg/mL in Cyclohexane			10ml
DRE-XA15120000CY	Methyl-pentachlorophenyl sulfide 100 µg/mL in Cyclohexane			1ml
				










Additional organic contaminants

Product code	Description			
(2RS)-2-Methylpentanoic Acid (2-Methylvaleric acid)				
CAS 97-61-0 DRE-C15147200	MW 116.1583 2-Methylvaleric acid	C ₆ H ₁₂ O ₂	250mg	
4-Methylpentan-2-one (MIBK)				
CAS 108-10-1 DRE-C15123000 DRE-XA15123000ME	MW 100.1589 4-Methyl-2-pentanone (MIBK) 4-Methyl-2-pentanone (MIBK) 100 µg/mL in Methanol	C ₆ H ₁₂ O	1ml 1ml	
4-Methylpent-3-en-2-one (Mesityloxide)				
CAS 141-79-7 DRE-C14913000	MW 98.143 Mesityloxide	C ₆ H ₁₀ O	250mg	
N-Methylperfluorooctanesulfonamide				
CAS 31506-32-8 DRE-C15115500	MW 513.1714 N-Methylperfluorooctanesulfonamide	C ₉ H ₄ F ₁₇ NO ₂ S	50mg	
2-Methylpyridine (2-Picolin)				
CAS 109-06-8 DRE-C16201500	MW 93.1265 2-Picolin (2-Methylpyridine)	C ₆ H ₇ N	250mg	
4-(Methylsulfonyl)-2-nitrobenzoic acid				
CAS 110964-79-9 DRE-C15143800	MW 245.2093 4-(Methylsulfonyl)-2-nitrobenzoic acid	C ₈ H ₇ NO ₆ S	100mg	
Methylthiouracil (6-Methyl-2-thiouracil)				
CAS 56-04-2 DRE-C15144300	MW 142.1789 6-Methyl-2-thiouracil	C ₅ H ₆ N ₂ OS	100mg	
Methyltin Trichloride (Trichloromethylstannane)				
CAS 993-16-8 DRE-GA09010362ME	MW 240.1035 Methyltin Trichloride 1000 µg/mL in Methanol(±)	CH ₃ Cl ₃ Sn	1ml	
Morpholine				
CAS 110-91-8 DRE-C15330000	MW 87.1204 Morpholine	C ₄ H ₉ NO	1ml	

Additional organic contaminants

Product code	Description			
Musk NN (Ethylene Brassylate)				
CAS 105-95-3	MW 270.3645	C ₁₅ H ₂₆ O ₄		
DRE-LA15359500CY	Musk NN 10 µg/mL in Cyclohexane		1ml	
Musk Xylene D15				
CAS 877119-10-3	MW 312.3564	C ₁₂ ² H ₁₅ N ₃ O ₆		
DRE-XA15360100AC	Musk xylene D15 100 µg/mL in Acetone		1.1ml	
Musk ambrette (Ambrettolide)				
CAS 83-66-9	MW 268.2658	C ₁₂ H ₁₆ N ₂ O ₅		
DRE-C15358500	Musk ambrette (technical)		10mg	
DRE-LA15358500CY	Musk ambrette (Ambrettolide) 10 µg/mL in Cyclohexane		1ml	
Musk ketone				
CAS 81-14-1	MW 294.3031	C ₁₄ H ₁₆ N ₂ O ₅		
DRE-C15359000	Musk ketone		100mg	
DRE-LA15359000CY	Musk ketone 10 µg/mL in Cyclohexane		1ml	
Musk moskene (1,1,3,3,5-Pentamethyl-4,6-dinitro-indane)				
CAS 116-66-5	MW 278.3037	C ₁₄ H ₁₆ N ₂ O ₄		
DRE-LA15359300CY	Musk moskene 10 µg/mL in Cyclohexane		1ml	
Musk tibeten (1-tert-Butyl-3,4,5-trimethyl-2,6-dinitrobenzene)				
CAS 145-39-1	MW 266.293	C ₁₃ H ₁₆ N ₂ O ₄		
DRE-L15359700CY	Musk tibeten 10 µg/mL in Cyclohexane		10ml	
Musk xylene				
CAS 81-15-2	MW 297.264	C ₁₂ H ₁₆ N ₃ O ₆		
DRE-C15360000	Musk xylene		100mg	
DRE-LA15360000CY	Musk xylene 10 µg/mL in Cyclohexane		1ml	
DRE-XA15360000CY	Musk xylene 100 µg/mL in Cyclohexane		1ml	
Mycrene (β-Myrcene)				
CAS 123-35-3	MW 136.234	C ₁₀ H ₁₆		
DRE-GA09010044IP	Mycrene 1000 µg/mL in Isopropanol(±)		1ml	
Myristoleic acid methyl ester				
CAS 56219-06-8	MW 240.3816	C ₁₅ H ₂₈ O ₂		
DRE-C15392050	Myristoleic acid methyl ester		50mg	










Additional organic contaminants

Product code	Description			
1-Naphthol				
CAS 90-15-3 DRE-C15430000	MW 144.1699 1-Naphthol(±)	C ₁₀ H ₈ O	500mg	
2-Naphthol (β-Naphthol)				
CAS 135-19-3 DRE-C15430500	MW 144.1699 2-Naphthol	C ₁₀ H ₈ O	500mg	
1,4-Naphthoquinone				
CAS 130-15-4 DRE-C15425000	MW 158.1534 1,4-Naphthoquinone	C ₁₀ H ₆ O ₂	250mg	
Nerol				
CAS 106-25-2 DRE-GA09010078IP	MW 154.2493 Nerol 1000 µg/mL in Isopropanol(±)	C ₁₀ H ₁₈ O	1ml	
Nitrilotriacetic Acid (NTA)				
CAS 139-13-9 DRE-E15655000	MW 191.1388 NTA (Nitrilotriacetic acid)	C ₆ H ₉ NO ₆	100mg	
2-Nitrobenzaldehyde				
CAS 552-89-6 DRE-C15556500	MW 151.1195 2-Nitrobenzaldehyde	C ₇ H ₅ NO ₃	100mg	
6-Nitro-1H-benzimidazole (5-Nitrobenzimidazole)				
CAS 94-52-0 DRE-C15557200	MW 163.1335 5-Nitrobenzimidazole	C ₇ H ₅ N ₃ O ₂	10mg	
2-Nitrobenzoic Acid				
CAS 552-16-9 DRE-C15557400	MW 167.1189 2-Nitrobenzoic acid	C ₇ H ₅ NO ₄	100mg	
4-Nitrobenzoic Acid				
CAS 62-23-7 DRE-C15557600	MW 167.1189 4-Nitrobenzoic acid	C ₇ H ₅ NO ₄	100mg	









Additional organic contaminants

Product code	Description			
2-Nitrobiphenyl				
CAS 86-00-0 DRE-L20963200CY	MW 199.2054 2-Nitrobiphenyl 10 µg/mL in Cyclohexane	C ₁₂ H ₉ NO ₂	10ml	
4-Nitrobiphenyl				
CAS 92-93-3 DRE-L20963400CY	MW 199.2054 4-Nitrobiphenyl 10 µg/mL in Cyclohexane	C ₁₂ H ₉ NO ₂	10ml	
5-Nitroguaiacol Sodium Salt				
CAS 67233-85-6 DRE-C15587000	MW 191.1166 5-Nitroguaiacol sodium salt	C ₇ H ₆ NO ₄ ·Na	100mg	
2-Nitropropane				
CAS 79-46-9 DRE-CA15599200	MW 89.0932 2-Nitropropane	C ₃ H ₇ NO ₂	1ml	
N-Nitroso-diphenylamine				
CAS 86-30-6 DRE-C15604500	MW 198.2206 N-Nitroso-diphenylamine	C ₁₂ H ₁₀ N ₂ O	100mg	
N-Nitroso-methyl-ethylamine				
CAS 10595-95-6 DRE-XA15605500ME	MW 88.1084 N-Nitroso-methyl-ethylamine 100 µg/mL in Methanol	C ₃ H ₈ N ₂ O	1ml	
(S)-N-Nitrosoanabasine D4 (pyridinyl D4)				
CAS 1426174-40-4 DRE-LA15600501AC	MW 195.2544 N-Nitrosoanabasine 2,4,5,6 D4 10 µg/mL in Acetone	C ₁₀ ² H ₈ H ₉ N ₃ O	1.1ml	
(RS)-N-Nitrosoanabasine				
CAS 37620-20-5 DRE-C15600500	MW 191.2297 N-Nitrosoanabasine	C ₁₀ H ₁₃ N ₃ O	10mg	
(S)-N-Nitrosoanatabine D4 (pyridinyl D4)				
CAS 1426174-82-4 DRE-LA15601001AC	MW 193.2385 N-Nitrosoanatabine 2,4,5,6 D4 10 µg/mL in Acetone	C ₁₀ ² H ₈ H ₇ N ₃ O	1.1ml	







Additional organic contaminants

Product code	Description			
(RS)-N-Nitrosoanatabine				
CAS 887407-16-1 DRE-C15601000	MW 189.2138 N-Nitrosoanatabine	C ₁₀ H ₁₁ N ₃ O	5mg	
N-Nitroso-di-n-butylamine				
CAS 924-16-3 DRE-C15602500 DRE-L15602500ME DRE-YA15602500ME	MW 158.2413 N-Nitroso-di-n-butylamine(±) N-Nitroso-di-n-butylamine 10 µg/mL in Methanol N-Nitroso-di-n-butylamine 1000 µg/mL in Methanol	C ₈ H ₁₈ N ₂ O	100mg 10ml 1ml	
N-Nitrosodiethanolamine				
CAS 1116-54-7 DRE-C15603000 DRE-L15603000AL	MW 134.1338 N-Nitroso-diethanolamine N-Nitroso-diethanolamine 10 µg/mL in Acetonitrile (*)	C ₄ H ₁₀ N ₂ O ₃	100mg 10ml	
N-Nitroso-diethylamine				
CAS 55-18-5 DRE-C15603500 DRE-L15603500ME DRE-XA15603500ME	MW 102.135 N-Nitroso-diethylamine N-Nitroso-diethylamine 10 µg/mL in Methanol N-Nitroso-diethylamine 100 µg/mL in Methanol	C ₆ H ₁₂ N ₂ O	100mg 10ml 1ml	
N-Nitrosodiisopropylamine				
CAS 601-77-4 DRE-C15604700 DRE-L15604700ME DRE-XA15604700ME	MW 130.1882 N-Nitroso-di-iso-propylamine(±) N-Nitroso-di-iso-propylamine 10 µg/mL in Methanol N-Nitroso-di-iso-propylamine 100 µg/mL in Methanol	C ₈ H ₁₄ N ₂ O	50mg 10ml 1ml	
N-Nitroso-dimethylamine				
CAS 62-75-9 DRE-C15604000 DRE-L15604000ME DRE-XA15604000ME	MW 74.0818 N-Nitroso-dimethylamine N-Nitroso-dimethylamine 10 µg/mL in Methanol N-Nitroso-dimethylamine 100 µg/mL in Methanol	C ₂ H ₆ N ₂ O	100mg 10ml 1ml	
N-Nitroso-dimethylamine D6				
CAS 17829-05-9 DRE-C15604010 DRE-XA15604010AC	MW 80.1188 N-Nitroso-dimethylamine D6 N-Nitroso-dimethylamine D6 100 µg/mL in Acetone	C ₂ ² H ₆ N ₂ O	10mg 1.1ml	
N-Nitroso-di-n-propylamine				
CAS 621-64-7 DRE-C15605000	MW 130.1882 N-Nitroso-di-n-propylamine(±)	C ₈ H ₁₄ N ₂ O	100mg	
N-Nitroso-di-n-propylamine D14				
CAS 93951-96-3 DRE-XA15605010AC	MW 144.2744 N-Nitroso-di-n-propylamine D14 100 µg/mL in Acetone	C ₆ ² H ₁₄ N ₂ O	1.1ml	









Additional organic contaminants

Product code	Description			
4-Nitrosomorpholine (N-Nitrosomorpholine)				
CAS 59-89-2	MW 116.1185	C ₄ H ₈ N ₂ O ₂		
DRE-C15606000	N-Nitrosomorpholine		100mg	
DRE-L15606000ME	N-Nitrosomorpholine 10 µg/mL in Methanol		10ml	
N-Nitrosornicotine ketone				
CAS 64091-91-4	MW 207.2291	C ₁₀ H ₁₃ N ₃ O ₂		
DRE-C15606300	N-Nitrosornicotine ketone		25mg	
N-Nitrosornicotine-ketone D4 (pyridyl D4)				
CAS 764661-24-7	MW 211.2538	C ₁₀ ² H ₄ H ₉ N ₃ O ₂		
DRE-LA15606301AC	N-Nitrosornicotine-ketone D4 (pyridyl D4) 10 µg/mL in Acetone		1.1ml	
N-Nitrosornicotine				
CAS 80508-23-2	MW 177.2031	C ₉ H ₁₁ N ₃ O		
DRE-C15606200	N-Nitrosornicotine		25mg	
N-Nitrosopiperidine				
CAS 100-75-4	MW 114.1457	C ₈ H ₁₀ N ₂ O		
DRE-C15606500	N-Nitrosopiperidine		25mg	
DRE-L15606500ME	N-Nitrosopiperidine 10 µg/mL in Methanol		10ml	
N-Nitrosopyrrolidine				
CAS 930-55-2	MW 100.1191	C ₄ H ₈ N ₂ O		
DRE-C15607000	N-Nitrosopyrrolidine		100mg	
DRE-L15607000ME	N-Nitrosopyrrolidine 10 µg/mL in Methanol		10ml	
DRE-XA15607000ME	N-Nitrosopyrrolidine 100 µg/mL in Methanol		1ml	
1,2,3,4,5,6,7,8,9-Nonachlorodecane (CP-10)				
CAS 890302-90-6	MW 452.2872	C ₁₀ H ₁₃ Cl ₉		
DRE-ZA15620500CY	1,2,3,4,5,6,7,8,9-Nonachlorodecane (CP-10) 1 µg/mL in Cyclohexane		1ml	
(2E,6Z)-Nona-2,6-dienal				
CAS 557-48-2	MW 138.2069	C ₉ H ₁₄ O		
DRE-CA15622600	trans-2, cis-6-Nonadienal		100mg	
Nonylphenol (technical)				
CAS 25154-52-3	MW n/a			No Structure
DRE-C15629000	Nonylphenol (technical)(‡)		250mg	










Additional organic contaminants

Product code	Description			
4-Nonylphenol Diethoxylate				
CAS 20427-84-3	MW 308.4556	C ₁₉ H ₃₂ O ₃		
DRE-C15631010	4-Nonylphenol-di-ethoxylate		10mg	
DRE-LA15631010AC	4-Nonylphenol-di-ethoxylate 10 µg/mL in Acetone		1ml	
4-Nonylphenol-ethoxylate (penta-) (technical) (branched)				
DRE-C15631018	4-Nonylphenol-ethoxylate (penta-) (technical) (branched)		250mg	
4-Nonylphenol-ethoxylate (technical)				
CAS 68412-54-4	MW n/a			
DRE-C15631000	4-Nonylphenol-ethoxylate (technical)		250mg	No Structure
DRE-L15631000ME	4-Nonylphenol-ethoxylate (technical) 10 µg/mL in Methanol		10ml	
4-Nonylphenol Monoethoxylate				
CAS 104-35-8	MW 264.403	C ₁₇ H ₂₈ O ₂		
DRE-C15631015	4-Nonylphenol-mono-ethoxylate(±)		25mg	
DRE-LA15631015AC	4-Nonylphenol-mono-ethoxylate 10 µg/mL in Acetone(±)		1ml	
DRE-LA15631015ME	4-Nonylphenol-mono-ethoxylate 10 µg/mL in Methanol		1ml	
4-Nonylphenol-mono-ethoxylate D2				
CAS n/a	MW n/a			
DRE-LA15631017AC	4-Nonylphenol-mono-ethoxylate D2 10 µg/mL in Acetone		1.1ml	No Structure
DRE-XA15631017CY	4-Nonylphenol-mono-ethoxylate D2 100 µg/mL in Cyclohexane		1.1ml	
4-Nonylphenol Triethoxylate				
CAS 51437-95-7	MW 352.5081	C ₂₁ H ₃₆ O ₄		
DRE-C15631019	4-Nonylphenol-tri-ethoxylate		10mg	
DRE-LA15631019AC	4-Nonylphenol-tri-ethoxylate 10 µg/mL in Acetone		1ml	
4-n-Nonylphenol				
CAS 104-40-5	MW 220.3505	C ₁₅ H ₂₄ O		
DRE-C15630000	4-n-Nonylphenol		100mg	
DRE-L15630000AL	4-n-Nonylphenol 10 µg/mL in Acetonitrile		10ml	
DRE-L15630000CY	4-n-Nonylphenol 10 µg/mL in Cyclohexane		10ml	
DRE-XA15630000CY	4-n-Nonylphenol 100 µg/mL in Cyclohexane		1ml	
4-n-Nonylphenol D8 (ring D4-ethylD4)				
CAS n/a	MW 228.3998	C ₁₅ ² H ₁₆ H ₁₆ O		
DRE-XA15630010AC	4-n-Nonylphenol D8 100 µg/mL in Acetone(±)		1ml	
4-Nonylphenoxyacetic Acid				
CAS 3115-49-9	MW 278.3865	C ₁₇ H ₂₆ O ₃		
DRE-LA15631020AC	4-Nonylphenoxy-acetic acid 10 µg/mL in Acetone		1ml	









Additional organic contaminants

Product code	Description			
Octabromodiphenylether				
CAS 32536-52-0 DRE-L15699000CY	MW 801.3757 Octabromodiphenylether 10 µg/mL in Cyclohexane	C ₁₂ H ₂ Br ₈ O	10ml	
2,3,4,5,6,7,8,9-Octachlorodecane (CP-9)				
CAS 214327-66-9 DRE-ZA15705000CY	MW 417.8422 2,3,4,5,6,7,8,9-Octachlorodecane (CP-9) 1 µg/mL in Cyclohexane	C ₁₀ H ₁₄ Cl ₈	1ml	
Octachlorostyrene				
CAS 29082-74-4 DRE-C15710000 DRE-L15710000CY DRE-XA15710000CY	MW 379.7096 Octachlorostyrene Octachlorostyrene 10 µg/mL in Cyclohexane Octachlorostyrene 100 µg/mL in Cyclohexane	C ₈ Cl ₈	20mg 10ml 1ml	
4-n-Octylphenol-di-ethoxylate				
CAS 51437-90-2 DRE-C15712803 DRE-LA15712803AC	MW 294.429 4-n-Octylphenol-di-ethoxylate 4-n-Octylphenol-di-ethoxylate 10 µg/mL in Acetone	C ₁₈ H ₃₀ O ₃	10mg 1ml	
4-Octylphenol-ethoxylate				
CAS 9036-19-5 DRE-C15712808	MW n/a 4-Octylphenol-ethoxylate (penta-) (technical)		250mg	No Structure
4-Octylphenol Polyethoxylate				
CAS 26636-32-8 DRE-C15712800 DRE-L15712800CY	MW 250.3764 4-Octylphenol-ethoxylate (technical) 4-Octylphenol-ethoxylate (technical) 10 µg/mL in Cyclohexane	C ₁₄ H ₂₂ O(C ₂ H ₄ O) _n	250mg 10ml	
4-Octylphenol				
CAS 1806-26-4 DRE-C15712000 DRE-L15712000IO	MW 206.3239 4-n-Octylphenol 4-n-Octylphenol 10 µg/mL in Isooctane	C ₁₄ H ₂₂ O	250mg 10ml	
4-tert-Octylphenol				
CAS 140-66-9 DRE-C15712100	MW 206.3239 4-tert-Octylphenol	C ₁₄ H ₂₂ O	250mg	
n-Octyltin Trichloride				
CAS 3091-25-6 DRE-C15715000 DRE-GA09010358DI	MW 338.2896 n-Octyltin-trichloride n-Octyltin-trichloride 1000 µg/mL in Dichloromethane(‡)	C ₈ H ₁₇ Cl ₃ Sn	100mg 5x1ml	








Additional organic contaminants

Product code	Description			
2-Octynoic Acid Methyl Ester				
CAS 111-12-6 DRE-C15715700	MW 154.2063 2-Octynoic acid-methyl ester	C ₉ H ₁₄ O ₂	250mg	
Oleic acid				
CAS 112-80-1 DRE-C15727000	MW 282.4614 Oleic acid	C ₁₈ H ₃₄ O ₂	0.5g	
Oxacycloheptadec-8-en-2-one				
CAS 123-69-3 DRE-C15756500 DRE-XA15756500CY	MW 252.3923 Oxacycloheptadec-8-en-2-one Oxacycloheptadec-8-en-2-one 100 µg/mL in Cyclohexane	C ₁₆ H ₂₈ O ₂	100mg 1ml	
Oxacyclohexadecan-2-one				
CAS 106-02-5 DRE-C15757000	MW 240.3816 Oxacyclohexadecan-2-one	C ₁₅ H ₂₈ O ₂	100mg	
Oxacyclohexadecen-2-one (Globalide)				
CAS 34902-57-3 DRE-C15757500	MW 238.3657 Oxacyclohexadecen-2-one (Globalide)	C ₁₅ H ₂₆ O ₂	100mg	
Oxybenzone (2-Hydroxy-4-methoxybenzophenone; Neo Heliopan BB)				
CAS 131-57-7 DRE-L14232500CY	MW 228.2433 2-Hydroxy-4-methoxybenzophenone (Neo Heliopan BB) 10 µg/mL in Cyclohexane	C ₁₄ H ₁₂ O ₃	10ml	
PBDE 15 (4,4'-Dibromodiphenyl Ether)				
CAS 2050-47-7 DRE-XA15898015ME	MW 327.9993 PBDE 15 100 µg/mL in Methanol	C ₁₂ H ₆ Br ₂ O	1ml	
PBDE 47 (2,2',4,4'-Tetrabromodiphenyl Ether)				
CAS 5436-43-1 DRE-C15898047 DRE-XA15898047ME	MW 485.7914 PBDE 47 (2,2',4,4'-Tetrabromodiphenyl ether) PBDE 47 (2,2',4,4'-Tetrabromodiphenyl ether) 100 µg/mL in Methanol	C ₁₂ H ₆ Br ₄ O	10mg 1ml	
Paraformaldehyde				
CAS 30525-89-4 DRE-C15848000	MW 30.026 Paraformaldehyde	(CH ₂ O) _n	1g	







Additional organic contaminants

Product code	Description			
Parsol 1789 (Avobenzone; 1-[4-(1,1-Dimethylethyl)phenyl]-3-(4-methoxyphenyl)-1,3-propanedione)				
CAS 70356-09-1 DRE-C15894520	MW 310.3869 Parsol 1789	C ₂₀ H ₂₂ O ₃	250mg	
Parsol 340 (Octocrylene; 2-Ethylhexyl 2-cyano-3,3-diphenylacrylate)				
CAS 6197-30-4 DRE-C15894510	MW 361.4767 Parsol 340	C ₂₄ H ₂₇ NO ₂	250mg	
Parsol HS (Ensulizole; 2-Phenyl-1H-benzimidazol-5-sulfonic acid)				
CAS 27503-81-7 DRE-C15894540	MW 274.2951 Parsol HS	C ₁₃ H ₁₀ N ₂ O ₃ S	250mg	
Parsol MCX (Octinoxate; 2-Ethylhexyl p-methoxycinnamate)				
CAS 5466-77-3 DRE-C13342000	MW 290.3972 2-Ethylhexyl p-methoxycinnamate	C ₁₈ H ₂₆ O ₃	100mg	
Parsol SLX				
CAS 207574-74-1 DRE-C15894570	MW n/a Parsol SLX		250mg	No Structure
Pentabromodiphenylether				
CAS 32534-81-9 DRE-C15938500 DRE-L15938500CY	MW 564.6875 Pentabromodiphenylether, technical Pentabromodiphenylether 10 µg/mL in Cyclohexane	C ₁₂ O ₅ Br ₅ H	10mg 10ml	
1,2,5,6,9-Pentachlorodecane (CP-3)				
CAS 175801-37-3 DRE-LA15963500CY	MW 314.507 1,2,5,6,9-Pentachlorodecane (CP-3) 10 µg/mL in Cyclohexane	C ₁₀ H ₁₇ Cl ₅	1ml	
Perfluorobutanesulfonic Acid				
CAS 375-73-5 DRE-C15986515	MW 300.0996 Perfluorobutanesulfonic acid	C ₄ HF ₉ O ₃ S	100mg	
Perfluorobutanoic Acid				
CAS 375-22-4 DRE-C15986520	MW 214.0384 Perfluorobutanoic acid	C ₄ HF ₇ O ₂	100mg	










Additional organic contaminants

Product code	Description			
Perfluorododecanoic Acid				
CAS 307-55-1 DRE-C15986620	MW 614.0984	C ₁₂ HF ₂₃ O ₂	50mg	
Perfluoroheptanesulfonic Acid				
CAS 375-92-8 DRE-C15986880	MW 450.1221	C ₇ HF ₁₉ O ₃ S	50mg	
Perfluoroheptanoic Acid				
CAS 375-85-9 DRE-C15986890	MW 364.0609	C ₇ HF ₁₃ O ₂	100mg	
Perfluorohexanesulfonic Acid				
CAS 355-46-4 DRE-C15986900	MW 400.1146	C ₆ HF ₁₃ O ₃ S	50mg	
Perfluorohexanoic Acid				
CAS 307-24-4 DRE-C15986910	MW 314.0534	C ₆ HF ₁₁ O ₂	100mg	
Perfluorononanoic Acid				
CAS 375-95-1 DRE-C15987000	MW 464.0759	C ₉ HF ₁₇ O ₂	100mg	
Perfluorooctane				
CAS 307-34-6 DRE-C15987100	MW 438.0569	C ₈ F ₁₈	100mg	
Perfluorooctane Sulfonic Acid Potassium Salt (PFOS)				
CAS 2795-39-3 DRE-C15987122	MW 538.22	C ₈ F ₁₇ O ₃ S·K	100mg	
Perfluorooctane sulfonamide PFOSA				
CAS 754-91-6 DRE-C15987110	MW 499.1448	C ₈ H ₂ F ₁₇ NO ₂ S	100mg	










Additional organic contaminants

Product code	Description			
Perfluorooctane sulfonic acid				
CAS 1763-23-1	MW 500.1296	C ₈ HF ₁₇ O ₃ S		
DRE-C15987120	Perfluorooctane sulfonic acid		10mg	
DRE-XA15987120ME	Perfluorooctane sulfonic acid 100 µg/mL in Methanol		1ml	
Perfluorooctane-1-sulfonyl Fluoride				
CAS 307-35-7	MW 502.1207	C ₈ F ₁₈ O ₂ S		
DRE-C15987130	Perfluorooctane sulfonyl fluoride		100mg	
Perfluorooctanoic Acid				
CAS 335-67-1	MW 414.0684	C ₈ HF ₁₅ O ₂		
DRE-C15987150	Perfluorooctanoic acid		100mg	
Perfluorooctanoic Acid Ammonium Salt (PFOA; POAA)				
CAS 3825-26-1	MW 431.0989	C ₈ HF ₁₅ O ₂ ·H ₃ N		
DRE-C15987152	Perfluorooctanoic acid ammonium salt PFOA POAA		100mg	
Perfluoropentanoic acid				
CAS 2706-90-3	MW 264.0459	C ₅ HF ₉ O ₂		
DRE-C15987200	Perfluoropentanoic acid		100mg	
Perfluorotributylamine (PFTBA)				
CAS 311-89-7	MW 671.092	C ₁₂ F ₂₇ N		
DRE-C15987500	Perfluorotributylamine (PFTBA)		100mg	
DRE-GA09010390ME	Perfluorotributylamine (PFTBA) MS Tuning Compound, 1000 µg/mL in p/t Methanol(‡)		1ml	
Perfluorotridecanoic Acid				
CAS 72629-94-8	MW 664.1059	C ₁₃ HF ₂₅ O ₂		
DRE-C15988000	Perfluorotridecanoic acid		100mg	
Perfluoroundecanoic Acid				
CAS 2058-94-8	MW 564.0909	C ₁₁ HF ₂₁ O ₂		
DRE-C15989000	Perfluoroundecanoic acid		100mg	
4-Phenoxybutyric Acid				
CAS 6303-58-8	MW 180.2005	C ₁₀ H ₁₂ O ₃		
DRE-C16045340	4-Phenoxybutyric acid		100mg	








Additional organic contaminants

Product code	Description			
2-Phenoxypropionic Acid				
CAS 940-31-8 DRE-C16045800	MW 166.1739 2-Phenoxypropionic acid	C ₉ H ₁₀ O ₃	500mg	
DL-Phenylalanine				
CAS 150-30-1 DRE-C16055000	MW 165.1891 DL-Phenylalanine	C ₉ H ₁₁ NO ₂	100mg	
Phenylisocyanate				
CAS 103-71-9 DRE-CA16059500	MW 119.1207 Phenylisocyanate	C ₇ H ₅ NO	1g	
N-Phenylthiourea				
CAS 103-85-5 DRE-C16075000	MW 152.2168 N-Phenylthiourea	C ₇ H ₉ N ₂ S	250mg	
N-(Phosphonomethyl)iminodiacetic acid monohydrate				
CAS 5994-61-6 DRE-C16145000	MW 245.1244 N-(Phosphonomethyl)iminodiacetic acid monohydrate	C ₅ H ₁₀ NO ₇ P·H ₂ O	250mg	
α-Pinene				
CAS 80-56-8 DRE-C16211000 DRE-GA09010041IP DRE-GS09010041IP	MW 136.234 α-Pinene α-Pinene 1000 µg/mL in Isopropanol(‡) α-Pinene 1000 µg/mL in Isopropanol(‡)	C ₁₀ H ₁₆	100mg 1ml 5x1ml	
β-Pinene				
CAS 127-91-3 DRE-C16211020 DRE-GA09010045IP	MW 136.234 β-Pinene β-Pinene 1000 µg/mL in Isopropanol(‡)	C ₁₀ H ₁₆	100mg 1ml	
Poly(hexamethylenediguanide) Hydrochloride				
CAS 32289-58-0 DRE-C16282500	MW 219.715 Polyhexamethylenediguanide hydrochloride	(C ₈ H ₁₇ N ₅) _n ·ClH	100mg	
Potassium Cyanate				
CAS 590-28-3 DRE-C16285300	MW 81.1151 Potassium cyanate	CNO·K	100mg	

Additional organic contaminants

Product code	Description			
Proline (L-(-)-Proline)				
CAS 147-85-3 DRE-C16345600	MW 115.1305 L(-)-Proline	C ₅ H ₉ NO ₂	100mg	
n-Propylcarbamate				
CAS 627-12-3 DRE-C16522500	MW 103.1198 n-Propylcarbamate (internal standard)(‡)	C ₄ H ₉ NO ₂	100mg	
Propylthiouracil (6-Propyl-2-thiouracil)				
CAS 51-52-5 DRE-C16530500	MW 170.2321 6-Propyl-2-thiouracil	C ₇ H ₁₀ N ₂ OS	250mg	
Propylthiourea				
CAS 927-67-3 DRE-C16531000	MW 118.2006 Propylthiourea	C ₄ H ₁₀ N ₂ S	100mg	
Scopoletin				
CAS 92-61-5 DRE-C16916000	MW 192.1681 Scopoletin	C ₁₀ H ₆ O ₄	50mg	
DL-Serine				
CAS 302-84-1 DRE-C16935700	MW 105.0926 DL-Serine	C ₃ H ₇ NO ₃	100mg	
Sodium Arsenate Tetrahydrate				
CAS 24992-20-9 DRE-C16970900	MW 257.9678 Sodium Arsenate Tetrahydrate	AsHO ₃	100mg	
Stearic Acid n-Butyl Ester				
CAS 123-95-5 DRE-C16974300	MW 340.5836 Stearic acid-n-butyl ester	C ₂₂ H ₄₄ O ₂	100mg	
Stearyl Stearate				
CAS 2778-96-3 DRE-C16974500	MW 536.9557 Stearyl stearate	C ₃₆ H ₇₂ O ₂	100mg	










Additional organic contaminants

Product code	Description			
Syringaldehyde				
CAS 134-96-3 DRE-C17080000	MW 182.1733 Syringaldehyde	C ₉ H ₁₀ O ₄	100mg	
O,S-TEPP				
CAS 645-78-3 DRE-C17240100	MW 306.2533 O,S-TEPP(‡)	C ₈ H ₂₀ O ₆ P ₂ S	50mg	
Tallowamine polyethoxylated				
CAS 61791-26-2 DRE-E17136000	MW n/a Tallowamine polyethoxylated		250mg	No Structure
DL-Tartaric acid				
CAS 133-37-9 DRE-C17137800	MW 150.0868 DL-Tartaric acid(‡)	C ₄ H ₆ O ₆	100mg	AND Enantiomer 
Terephthalic Acid				
CAS 100-21-0 DRE-C17321500	MW 166.1308 Terephthalic acid(‡)	C ₈ H ₆ O ₄	250mg	
Terephthalic Acid Bis(2-ethylhexyl) Ester				
CAS 6422-86-2 DRE-C17321900	MW 390.5561 Terephthalic acid, bis(2-ethylhexyl) ester	C ₂₄ H ₃₈ O ₄	250mg	
Terpene Mix 1 + 2				
DRE-KA09000238ME	Terpene Mix 1 (DRE-GA09000272ME) and 2 (DRE-GA09000273ME), 100 µg/mL in p/t Methanol(‡)			1ea
	DRE-GA09000272ME			1x1ml
	DRE-GA09000273ME			1x1ml
α-Terpinene				
CAS 99-86-5 DRE-C17322320	MW 136.234 α-Terpinene	C ₁₀ H ₁₆	100mg	
γ-Terpinene				
CAS 99-85-4 DRE-C17322325 DRE-GA09010077IP	MW 136.234 γ-Terpinene γ-Terpinene 1000 µg/mL in Isopropanol(‡)	C ₁₀ H ₁₆	100mg 1ml	










Additional organic contaminants

Product code	Description			
1-Terpinen-4-ol				
CAS 562-74-3 DRE-C17322300	MW 154.2493 1-Terpinen-4-ol	C ₁₀ H ₁₈ O	100mg	
α-Terpineol				
CAS 98-55-5 DRE-YS09010013AC	MW 154.2493 α-Terpineol 1000 µg/mL in Acetone(‡)	C ₁₀ H ₁₈ O	5x1ml	
Terpinolene (δ-Terpinene)				
CAS 586-62-9 DRE-GA09010042IP	MW 136.234 Terpinolene 1000 µg/mL in Isopropanol(‡)	C ₁₀ H ₁₆	1ml	
Tetrabromobisphenol A D10 (dimethyl D3, bisphenol-3,5-D2)				
CAS n/a DRE-XA17324701AL	MW 553.9322 Tetrabromobisphenol A D10 100 µg/mL in Acetonitrile	C ₁₅ ² H ₁₀ H ₂ Br ₄ O ₂	1.1ml	
Tetrabromobisphenol A Diallyl Ether				
CAS 25327-89-3 DRE-C17324750	MW 623.9983 Tetrabromobisphenol A diallyl ether	C ₂₁ H ₂₀ Br ₄ O ₂	100mg	
Tetrabromobisphenol A-bis(2,3-dibromopropyl Ether)				
CAS 21850-44-2 DRE-C17324800	MW 943.6143 Tetrabromobisphenol A dibromopropyl ether	C ₂₁ H ₂₀ Br ₈ O ₂	100mg	
3,3',5,5'-Tetrabromobisphenol A				
CAS 79-94-7 DRE-C17324700	MW 543.8706 Tetrabromobisphenol A	C ₁₅ H ₁₂ Br ₄ O ₂	250mg	
Tetrabutyltin				
CAS 1461-25-2 DRE-C17328000 DRE-L17328000CY DRE-GA09010353ME	MW 347.167 Tetrabutyltin Tetrabutyltin 10 µg/mL in Cyclohexane Terbutyltin 1000 µg/mL in Methanol(‡)	C ₁₆ H ₃₆ Sn	250mg 10ml 1ml	
2,3,4,6-Tetrachloroanisole				
CAS 938-22-7 DRE-X17333150HA	MW 245.9181 2,3,4,6-Tetrachloroanisole, 100 µg/mL in n-Hexane/Acetone (9:1)	C ₇ H ₄ Cl ₄ O	10ml	










Additional organic contaminants

Product code	Description			
1,2,9,10-Tetrachlorodecane (CP-2)				
CAS 205646-11-3	MW 280.0619	C ₁₀ H ₁₈ Cl ₄		
DRE-LA17356600CY	1,2,9,10-Tetrachlorodecane (CP-2) 10 µg/mL in Cyclohexane		1ml	
2,5,6,9-Tetrachlorodecane (CP-1)				
CAS n/a	MW 280.0619	C ₁₀ H ₁₈ Cl ₄		
DRE-LA17356500CY	2,5,6,9-Tetrachlorodecane (CP-1) 10 µg/mL in Cyclohexane		1ml	
2,3,7,8-Tetrachlorodibenzo-p-dioxine				
CAS 1746-01-6	MW 321.971	C ₁₂ H ₄ Cl ₄ O ₂		
DRE-ZA17356700TO	2,3,7,8-Tetrachlorodibenzo-p-dioxine 0.001 µg/mL in Toluene		1ml	
Tetrachloroguaiacol				
CAS 2539-17-5	MW 261.9175	C ₇ H ₄ Cl ₄ O ₂		
DRE-C17358700	Tetrachloroguaiacol(±)		100mg	
Tetrachlorohydroquinone				
CAS 87-87-6	MW 247.8909	C ₆ H ₂ Cl ₄ O ₂		
DRE-C17358900	Tetrachlorohydroquinone		250mg	
2,4,7,9-Tetramethyl-5-decyne-4,7-diol				
CAS 126-86-3	MW 226.355	C ₁₄ H ₂₆ O ₂		
DRE-C17413600	2,4,7,9-Tetramethyl-5-decyne-4,7-diol (mix of (±) and meso)		100mg	
Tetra-n-pentyltin				
CAS 3765-65-9	MW 403.2734	C ₂₀ H ₄₄ Sn		
DRE-C17415100	Tetra-n-pentyltin		100mg	
Tetraphenyltin				
CAS 595-90-4	MW 427.1256	C ₂₄ H ₂₀ Sn		
DRE-C17415600	Tetraphenyl tin		100mg	
Tetra-n-propyltin				
CAS 2176-98-9	MW 291.0607	C ₁₂ H ₂₈ Sn		
DRE-C17415700	Tetra-n-propyltin		100mg	










Additional organic contaminants

Product code	Description			
2,2'-Thiodiacetic Acid				
CAS 123-93-3 DRE-C17486500	MW 150.153 2,2'-Thiodiacetic acid	C ₄ H ₆ O ₄ S	250mg	
DL-Threonine				
CAS 80-68-2 DRE-C17575000	MW 119.1192 DL-Threonine	C ₄ H ₉ NO ₃	100mg	AND Enantiomer 
2,4-Toluenediisocyanate				
CAS 584-84-9 DRE-CA17594500	MW 174.1561 2,4-Toluenediisocyanate (2,4-TDI)	C ₉ H ₈ N ₂ O ₂	1ml	
2,6-Toluenediisocyanate				
CAS 91-08-7 DRE-CA17594600	MW 174.1561 2,6-Toluenediisocyanate	C ₉ H ₈ N ₂ O ₂	100mg	
p-Toluenesulfonamide				
CAS 70-55-3 DRE-C17594700	MW 171.2169 p-Toluenesulfonamide	C ₇ H ₉ NO ₂ S	100mg	
Tolyltriazole				
CAS 29385-43-1 DRE-C17600200	MW 266.3012 Tolyltriazole	2C ₇ H ₇ N ₃	100mg	
Toxaphene Parlar 11				
CAS 165820-10-0 DRE-ZA22001100CY	MW 342.9044 Toxaphene Parlar-No. 11, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₆	1ml	AND Enantiomer 
Toxaphene Parlar 26				
CAS 142534-71-2 DRE-ZA22002600CY	MW 413.8104 Toxaphene Parlar-No. 26, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₈	1ml	AND Enantiomer 
Toxaphene Parlar 32				
CAS 51775-36-1 DRE-ZA22003200CY	MW 379.3653 Toxaphene Parlar-No. 32, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₁ Cl ₇	1ml	AND Enantiomer 










Additional organic contaminants

Product code	Description			
Toxaphene Parlar 40				
CAS 166021-27-8 DRE-ZA22004000CY	MW 413.8104 Toxaphene Parlar-No. 40, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₈	1ml	AND Enantiomer 
Toxaphene Parlar 42				
CAS 187348-02-3 DRE-ZA22004200CY	MW 413.8104 Toxaphene Parlar-No. 42, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₈	1ml	Mixed 
Toxaphene Parlar 44				
CAS 165820-17-7 DRE-ZA22004400CY	MW 413.8104 Toxaphene Parlar-No. 44, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₁₀ Cl ₈	1ml	AND Enantiomer 
Toxaphene Parlar 50				
CAS 66860-80-8 DRE-ZA22005000CY	MW 448.2555 Toxaphene Parlar-No. 50, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₈ Cl ₉	1ml	AND Enantiomer 
Toxaphene Parlar 62				
CAS 154159-06-5 DRE-ZA22006200CY	MW 448.2555 Toxaphene Parlar-No. 62, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₈ Cl ₉	1ml	AND Enantiomer 
Toxaphene Parlar 69				
CAS 151183-19-6 DRE-ZA22006900CY	MW 482.7005 Toxaphene Parlar-No. 69, ca.1 µg/mL in Cyclohexane	C ₁₀ H ₆ Cl ₁₀	1ml	AND Enantiomer 
Tribromoacetic Acid				
CAS 75-96-7 DRE-C17663000 DRE-YA17663000MB	MW 296.7401 Tribromoacetic acid Tribromoacetic acid 1000 µg/mL in Methyl-tert. butyl ether	C ₂ HBr ₃ O ₂	100mg 1ml	
3,5,4'-Tribromosalicylanilide (3,5-Dibromo-N-(4-bromophenyl)-2-hydroxy-benzamide)				
CAS 87-10-5 DRE-C17666500	MW 449.9201 3,5,4'-Tribromosalicylanilide	C ₁₃ H ₈ Br ₃ NO ₂	100mg	
Tributylethyltin				
CAS 19411-60-0 DRE-C17667600	MW 319.1139 Tributylethyltin	C ₁₄ H ₃₂ Sn	50mg	










Additional organic contaminants

Product code	Description			
Tributyltin Chloride (TBTC)				
CAS 1461-22-9	MW 325.5058	C ₁₂ H ₂₇ ClSn		
DRE-C17150000	TBTC (Tributyltin-chloride)		250mg	
DRE-L17150000IO	TBTC 10 µg/mL in Isooctane		10ml	
DRE-GA09010355ME	Tri-n-butyltin Chloride 1000 µg/mL in Methanol(‡)		1ml	
Tributyltin Oxide (TBTO)				
CAS 56-35-9	MW 596.105	C ₂₄ H ₅₄ OSn ₂		
DRE-C17160000	TBTO		250mg	
Tricaprin				
CAS 621-71-6	MW 554.8418	C ₃₃ H ₆₂ O ₆		
DRE-AY09010004PY	Tricaprin 8000 µg/mL in Pyridine		5ml	
DRE-SY09010004PY	Tricaprin 8000 µg/mL in Pyridine		5x5ml	
α,2,4-Trichloroacetophenone				
CAS 4252-78-2	MW 223.4837	C ₈ H ₅ Cl ₃ O		
DRE-C17690000	α-2,4-Trichloroacetophenon		250mg	
2,3,6-Trichlorobenzoic Acid				
CAS 50-31-7	MW 225.4565	C ₇ H ₃ Cl ₃ O ₂		
DRE-C17730000	2,3,6-Trichlorobenzoic acid		100mg	
2,4,6-Trichlorobenzoic Acid				
CAS 50-43-1	MW 225.4565	C ₇ H ₃ Cl ₃ O ₂		
DRE-C17730400	2,4,6-Trichlorobenzoic acid		100mg	
2,3,4-Trichlorophenol Acetate				
CAS 61925-89-1	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-C17775100	2,3,4-Trichlorophenol acetate		25mg	
2,3,5-Trichlorophenol Acetate				
CAS 61925-88-0	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-L17775200IO	2,3,5-Trichlorophenol acetate 10 µg/mL in Isooctane		10ml	
2,3,6-Trichlorophenol Acetate				
CAS 61925-87-9	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-L17775300IO	2,3,6-Trichlorophenol acetate 10 µg/mL in Isooctane		10ml	







Additional organic contaminants

Product code	Description			
2,4,5-Trichlorophenol Acetate				
CAS 5393-75-9	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-C17775400	2,4,5-Trichlorophenol acetate		20mg	
DRE-L17775400IO	2,4,5-Trichlorophenol acetate 10 µg/mL in Isooctane		10ml	
2,4,6-Trichlorophenol Acetate				
CAS 23399-90-8	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-C17775500	2,4,6-Trichlorophenol acetate		20mg	
DRE-L17775500IO	2,4,6-Trichlorophenol acetate 10 µg/mL in Isooctane		10ml	
3,4,5-Trichlorophenol Acetate				
CAS 59190-61-3	MW 239.4831	C ₈ H ₅ Cl ₃ O ₂		
DRE-C17775600	3,4,5-Trichlorophenol acetate		25mg	
2,4,5-Trichlorophenol Ethanolamine Salt				
CAS 56961-45-6	MW 258.5295	C ₈ H ₅ Cl ₃ O·C ₂ H ₇ NO		
DRE-C17776000	2,4,5-Trichlorophenol ethanolamine salt		100mg	
2,4,6-Trichlorophenoxyacetic Acid				
CAS 575-89-3	MW 255.4825	C ₈ H ₅ Cl ₃ O ₃		
DRE-C17777000	2,4,6-Trichlorophenoxyacetic acid		100mg	
4-(2,4,5-Trichlorophenoxy)butanoic Acid (2,4,5-TB)				
CAS 93-80-1	MW 283.5357	C ₁₀ H ₅ Cl ₃ O ₃		
DRE-XA17140000AL	2,4,5-TB 100 µg/mL in Acetonitrile		1ml	
Triclocarban				
CAS 101-20-2	MW 315.5824	C ₁₃ H ₉ Cl ₃ N ₂ O		
DRE-C17797000	Triclocarban		250mg	
Triclosan				
CAS 3380-34-5	MW 289.5418	C ₁₂ H ₇ Cl ₃ O ₂		
DRE-C17803000	Triclosan		250mg	
DRE-L17803000CY	Triclosan (Irgasan) 10 µg/mL in Cyclohexane		10ml	
Triclosan D3 (2,4-dichlorophenoxy D3)				
CAS 1020719-98-5	MW 292.5603	C ₁₂ ² H ₃ H ₄ Cl ₂ O ₂		
DRE-XA17803010CY	Triclosan D3 (2,4-dichlorophenoxy D3) 100 µg/mL in Cyclohexane		1.1ml	










Additional organic contaminants

Product code	Description			
Triclosan Methyl ether				
CAS 4640-01-1 DRE-C17803300	MW 303.5684 Triclosan methyl ether	$C_{13}H_9Cl_3O_2$	50mg	
Triclosan methyl D3 (methoxy D3)				
CAS 1020720-00-6 DRE-XA17803310AC	MW 306.5868 Triclosan methyl D3 (methoxy D3) 100 µg/mL in Acetone	$C_{13}^2H_9H_6Cl_3O_2$	1.1ml	
Tricyclohexyltin Chloride				
CAS 3091-32-5 DRE-C17814000	MW 403.6176 Tricyclohexyltin chloride	$C_{18}H_{33}ClSn$	250mg	
Triethyllead Chloride				
CAS 1067-14-7 DRE-C17834500	MW 329.8363 Triethyl-lead chloride	$C_6H_{15}ClPb$	100mg	
Triethyltin Chloride				
CAS 994-31-0 DRE-L17838500CY	MW 241.3463 Triethyltin-chloride 10 µg/mL in Cyclohexane	$C_6H_{15}ClSn$	10ml	
Trifluoroacetic acid				
CAS 76-05-1 DRE-C17844400	MW 114.0233 Trifluoroacetic acid	$C_2HF_3O_2$	1ml	
2,3,5-Triiodobenzoic Acid				
CAS 88-82-4 DRE-C17870000	MW 499.8109 2,3,5-Triiodobenzoic acid	$C_7H_3I_3O_2$	250mg	
Trimethyllead Chloride				
CAS 1520-78-1 DRE-C17882500	MW 287.7566 Trimethyl lead chloride	C_3H_9ClPb	100mg	
1,1,1-Trimethylolpropane Phosphate				
CAS 1005-93-2 DRE-C17882700	MW 178.1229 1,1,1-Trimethylolpropane Phosphate	$C_6H_{11}O_4P$	25mg	

Additional organic contaminants

Product code	Description			
Trimethylphosphate				
CAS 512-56-1 DRE-C17884500	MW 140.0749 Trimethylphosphate	C ₃ H ₉ O ₄ P	250mg	
O,O,O-Trimethylphosphorthioate (O,O,O-Trimethylthiophosphate)				
CAS 152-18-1 DRE-C17886000	MW 156.1405 O,O,O-Trimethylthiophosphate	C ₃ H ₉ O ₃ PS	100mg	
O,O,S-Trimethylphosphorthioate (O,O,S-Trimethylthiophosphate)				
CAS 152-20-5 DRE-C17886050	MW 156.1405 O,O,S-Trimethylthiophosphate	C ₃ H ₉ O ₃ PS	100mg	
Trimethylsulfonium Iodide				
CAS 2181-42-2 DRE-C17885000	MW 204.073 Trimethylsulfonium iodide	C ₃ H ₉ S ⁺ I ⁻	250mg	
Trimethyltin Chloride				
CAS 1066-45-1 DRE-C17888000	MW 199.2666 Trimethyltin-chloride	C ₃ H ₉ ClSn	250mg	
Trioctyltin Chloride				
CAS 2587-76-0 DRE-C17892000	MW 493.8247 Trioctyltin-chloride	C ₂₄ H ₅₁ ClSn	25mg	
Triphenylphosphine				
CAS 603-35-0 DRE-C17893300	MW 262.2855 Triphenylphosphine	C ₁₈ H ₁₅ P	250mg	
Tripropylene Glycol Monomethyl Ether (Mixture of Isomers)				
CAS 25498-49-1 DRE-C17893800	MW n/a Tripropylene glycol-monomethylether		250mg	No Structure
Tripropyltin Chloride				
CAS 2279-76-7 DRE-C17893850	MW 283.426 Tri-n-propyltin chloride	C ₉ H ₂₁ ClSn	100mg	

Additional organic contaminants

Product code	Description			
Tryptophan				
CAS 54-12-6 DRE-C17895000	MW 204.2252 Tryptophan	C ₁₁ H ₁₂ N ₂ O ₂	500mg	
L-Tyrosine disodium salt				
CAS 69847-45-6 DRE-C17896010	MW 225.1522 L-Tyrosine disodium salt	C ₉ H ₉ NO ₃ ·2Na	250mg	
DL-Tyrosine				
CAS 556-03-6 DRE-C17896000	MW 181.1885 DL-Tyrosine	C ₉ H ₁₁ NO ₃	100mg	
Urethan (Ethyl carbamate)				
CAS 51-79-6 DRE-C13322000	MW 89.0932 Ethyl carbamate(‡)	C ₃ H ₇ NO ₂	250mg	
Urotropine (Methenamine)				
CAS 100-97-0 DRE-C17897500	MW 140.1863 Urotropine(‡)	C ₆ H ₁₂ N ₄	250mg	
(+) -Valencene (Valencene sesquiterpene)				
CAS 4630-07-3 DRE-GA09010079IP	MW 204.3511 (+) -Valencene 1000 µg/mL in Isopropanol(‡)	C ₁₅ H ₂₄	1ml	
DL-Valine				
CAS 516-06-3 DRE-C17899950	MW 117.1463 DL-Valine	C ₅ H ₁₁ NO ₂	100mg	
Vinyl Bromide				
CAS 593-60-2 DRE-YS09010029ME	MW 106.9492 Vinyl Bromide 1000 µg/mL in p/t Methanol(‡)	C ₂ H ₃ Br	5x1ml	
Zinc Phosphate				
CAS 7779-90-0 DRE-C17949500	MW 386.1697 Zinc phosphate	2O ₄ P·3Zn	100mg	

Additional organic contaminants

Product code	Description		
Zinc neodecanoate			
CAS 27253-29-8	MW n/a		No Structure
DRE-C17949400	Zinc neodecanoate	100mg	
Chlorinated Acids Mix 1			
DRE-XA05550100AL	Chlorinated Acids Mix 1, 100 µg/mL in Acetonitrile		1ml
	2,4-D	Acifluorfen	
	Bentazone	Chloramben	
	Dicamba	Dichlorprop	
	Fenoprop	Picloram	
Chlorinated Aromatics-Mix 7			
DRE-XA19060007ME	Chlorinated Aromatics-Mix 7 100 µg/mL in Methanol		1ml
	1,2,3,4-Tetrachlorobenzene	1,2,3,5-Tetrachlorobenzene	
	1,2,3-Trichlorobenzene	1,2,4,5-Tetrachlorobenzene	
	1,2,4-Trichlorobenzene	1,2-Dichlorobenzene	
	1,3,5-Trichlorobenzene	1,3-Dichlorobenzene	
	1,4-Dichlorobenzene	Chlorobenzene	
	Hexachlorobenzene	Pentachlorobenzene	
Chlorinated Organic and Desinfecant by-products Mix 1			
DRE-XA05510100AC	Chlorinated Organic and Desinfecant by-products Mix 1 100 µg/mL in Acetone (*)		1ml
	1,1,1-Trichloroethane	1,1-Dichloropropanone-2	1,2-Dibromo-3-chloropropane
	Bromochloroacetonitrile	Bromodichloromethane	Chloropicrin
	Dibromochloromethane	Dichloroacetonitrile	Tetrachloroethene
	Tribromomethane	Trichloroacetonitrile	Trichloroethene
			1,2-Dibromoethane
			Dibromoacetonitrile
			Tetrachloromethane
			Trichloromethane
Chloroparaffin C10 Mix 1			
DRE-LA22102000CY	Chloroparaffin C10 Mix 1 in Cyclohexane		1ml
	2,5,6,9-Tetrachlorodecane CP-1 [10 µg/mL]		
	1,2,9,10-Tetrachlorodecane CP-2 [10 µg/mL]		
	1,2,5,6,9-Pentachlorodecane CP-3 [10 µg/mL]		
	1,2,5,6,9,10-Hexachlorodecane CP-4 [10 µg/mL]		
	1,2,5,6,9,10-Hexachlorodecane CP-5 (Mix of 2 Diastereoisomers) [10 µg/mL]		
Chloroparaffin C10 Mix 2			
DRE-LA22102200CY	Chloroparaffin C10 Mix 2 in Cyclohexane		1ml
	1,2,5,6,9-Pentachlorodecane CP-3 [19 µg/mL]		
	1,2,4,5,9,10-Hexachlorodecane CP-6 [3.7 µg/mL]		
	1,2,4,5,6,9,10-Heptachlorodecane CP-7 [2.5 µg/mL]		
	2,3,4,5,6,7,8,9-Octachlorodecane CP-9 [0.5 µg/mL]		
	1,2,3,4,5,6,7,8,9-Nonachlorodecane CP-10 [0.9 µg/mL]		
Chloroparaffin C10-C12 Mix 5			
DRE-ZA22102105HP	Chloroparaffin C10-C12 Mix 5 in n-Heptane		1ml
DRE-ZS22102105HP	Chloroparaffin C10-C12 Mix 5 in n-Heptane		3x1ml
	Chloroparaffin C10 50.18%Cl [1 µg/mL]		
	Chloroparaffin C10 55.00%Cl [1 µg/mL]		
	Chloroparaffin C11 45.50%Cl [1 µg/mL]		
	Chloroparaffin C11 50.21%Cl [1 µg/mL]		
	Chloroparaffin C12 45.32%Cl [2 µg/mL]		
	Chloroparaffin C12 50.18%Cl [2 µg/mL]		
	Chloroparaffin C12 55.00%Cl [2 µg/mL]		

Additional organic contaminants

Product code	Description	
Chloroparaffin C10-C13 Mix 1		
DRE-LA22102101CY	Chloroparaffin C10-C13 Mix 1 in Cyclohexane	1ml
	Chloroparaffin C10 44.82%Cl [0.5 µg/mL] Chloroparaffin C10 50.18%Cl [0.5 µg/mL] Chloroparaffin C11 45.50%Cl [1.2 µg/mL] Chloroparaffin C11 50.21%Cl [2.6 µg/mL] Chloroparaffin C12 45.32%Cl [1.0 µg/mL] Chloroparaffin C12 50.18%Cl [2.4 µg/mL] Chloroparaffin C13 50.23%Cl [1.8 µg/mL]	
Chloroparaffin C10-C13 Mix 2		
DRE-LA22102102CY	Chloroparaffin C10-C13 Mix 2 in Cyclohexane	1ml
	Chloroparaffin C10 55.00%Cl [1.0 µg/mL] Chloroparaffin C11 55.20%Cl [2.5 µg/mL] Chloroparaffin C12 50.18%Cl [0.5 µg/mL] Chloroparaffin C12 65.08%Cl [0.2 µg/mL] Chloroparaffin C13 59.98%Cl [0.4 µg/mL]	Chloroparaffin C11 50.21%Cl [0.5 µg/mL] Chloroparaffin C11 60.53%Cl [1.4 µg/mL] Chloroparaffin C12 55.00%Cl [2.5 µg/mL] Chloroparaffin C13 55.03%Cl [1.0 µg/mL]
Chloroparaffin C10-C13 Mix 3		
DRE-LA22102103CY	Chloroparaffin C10-C13 Mix 3 in Cyclohexane	1ml
	Chloroparaffin C10 65.02%Cl [2.0 µg/mL] Chloroparaffin C11 65.25%Cl [3.2 µg/mL] Chloroparaffin C12 69.98%Cl [3.1 µg/mL] Chloroparaffin C13 65.18%Cl [1.7 µg/mL]	
Chloroparaffin C10-C13 Mix 4		
DRE-LA22102104CY	Chloroparaffin C10-C13 Mix 4 in Cyclohexane	1ml
	Chloroparaffin C10 50.18%Cl [0.5 µg/mL] Chloroparaffin C11 50.21%Cl [0.5 µg/mL] Chloroparaffin C11 60.53%Cl [1.9 µg/mL] Chloroparaffin C12 55.00%Cl [2.5 µg/mL] Chloroparaffin C13 55.03%Cl [1.0 µg/mL]	Chloroparaffin C10 55.00%Cl [0.5 µg/mL] Chloroparaffin C11 55.20%Cl [2.0 µg/mL] Chloroparaffin C12 50.18%Cl [0.5 µg/mL] Chloroparaffin C12 65.08%Cl [0.2 µg/mL] Chloroparaffin C13 59.98%Cl [0.4 µg/mL]
Chloroparaffin C10-C13 Mix 6 in n-Heptane		
DRE-ZA22102106HP	Chloroparaffin C10-C13 Mix 6 in n-Heptane	1ml
DRE-ZS22102106HP	Chloroparaffin C10-C13 Mix 6 in n-Heptane	3x1ml
	Chloroparaffin C10 60.09%Cl [0.5 µg/mL] Chloroparaffin C11 55.20%Cl [0.6 µg/mL] Chloroparaffin C11 65.25%Cl [3.0 µg/mL] Chloroparaffin C12 55.00%Cl [2.0 µg/mL] Chloroparaffin C13 59.98%Cl [0.1 µg/mL]	Chloroparaffin C10 65.02%Cl [1.1 µg/mL] Chloroparaffin C11 60.53%Cl [1.0 µg/mL] Chloroparaffin C12 50.18%Cl [0.8 µg/mL] Chloroparaffin C12 65.08%Cl [0.9 µg/mL]
Chloroparaffin C10-C13 Mix 7 in n-Heptane		
DRE-ZA22102107HP	Chloroparaffin C10-C13 Mix 7 in n-Heptane	1ml
DRE-ZS22102107HP	Chloroparaffin C10-C13 Mix 7 in n-Heptane	3x1ml
	Chloroparaffin C10 65.02%Cl [0.3 µg/mL] Chloroparaffin C11 60.53%Cl [0.5 µg/mL] Chloroparaffin C11 65.25%Cl [0.7 µg/mL] Chloroparaffin C12 65.08%Cl [1.0 µg/mL] Chloroparaffin C12 69.98%Cl [0.8 µg/mL] Chloroparaffin C13 59.98%Cl [0.7 µg/mL] Chloroparaffin C13 65.18%Cl [6.0 µg/mL]	
DNPH-Mix 2		
DRE-YA18001611AL	DNPH-Mix 2, 215-700 µg/mL in Acetonitrile	1ml
Pentanal-DNPH [310µg/mL] Acrolein-DNPH [420µg/mL] Crotonaldehyde-DNPH [360µg/mL] Formaldehyde-DNPH [700µg/mL] Hexanal-DNPH [280µg/mL] Nonanal-DNPH [225µg/mL] m-Tolualdehyde-DNPH [250µg/mL]	p-Tolualdehyde-DNPH [250µg/mL] Benzaldehyde-DNPH [270µg/mL] Cyclohexanone-DNPH [285µg/mL] Furfural-DNPH [290µg/mL] Isovaleraldehyde-DNPH [310µg/mL] Octanal-DNPH [240µg/mL]	Acetaldehyde-DNPH [510µg/mL] 2-Butanone-DNPH [350µg/mL] Decanal-DNPH [215µg/mL] Glutaraldehyde-bis(DNPH) [460µg/mL] Methacrylaldehyde-DNPH [360µg/mL] Propionaldehyde-DNPH [410µg/mL]
	Acetone-DNPH [410µg/mL] Butyraldehyde-DNPH [350µg/mL] 2,5-Dimethylbenzald.-DNPH [230µg/mL] Heptanal-DNPH [260µg/mL] 4-Methyl-2-pentanone-DNPH [255µg/mL] o-Tolualdehyde-DNPH [250µg/mL]	

Additional organic contaminants

Product code	Description	
Drinking Water Odor Mix 1		
DRE-XA05520700ME	Drinking Water Odor Mix 1 100 µg/mL in Methanol	1ml
	(±)-Geosmin 2-Methylisoborneol	
EPH Fractionation Surrogate Mix		
DRE-YS09000042HE	EPH Fractionation Surrogate Mix 4000 µg/mL in Tol.:Hex. 1:1(±)	5x1ml
	2-bromonaphthalene 2-fluorobiphenyl	
Fruit Juice Organic Acid Mix		
DRE-GS09000056WA	Fruit Juice Organic Acid Mix in Water(±)	5x1ml
	DL-tartaric acid [2000 mg/L] d-(-)-quinic acid [2000 mg/L] d-malic acid [2000 mg/L] citric acid [2000 mg/L] fumaric acid [25,8 mg/L]	
Furfural Mix 1		
DRE-GS09000387AC	Furfural Mix 1 5000 µg/mL in Acetone(±)	5x1ml
	furfural furfuryl alcohol phenol	
Furfural Mix 2		
DRE-GS09000388DI	Furfural Mix 2 5000 µg/mL in Dichloromethane(±)	5x1ml
	furfural furfuryl alcohol furoic acid phenol	
Internal Standard Solution Mix16		
DRE-YA05000016ME	Internal Standard Solution Mix16 2000 µg/mL in Methanol	1ml
	2-Bromo-1-chloropropane Fluorobenzene	
Internal Standards Mix 25		
DRE-XA05250600AC	Internal Standards Mix 25 500 µg/mL in Acetone	1ml
	Acenaphthene D10 Chrysene D12 Perylene D12 Phenanthrene D10	
Mix 1 of Musk a. Polycy. Musk Comp.		
DRE-LA19020100CY	Mix 1 of Musk a. Polycy. Musk Comp. 10 µg/mL in Cyclohexane	1ml
	ADBI (Celestolide) AHTN (Tonalid®) DPMI (Cashmeran®) Musk ambrette Musk xylene	AHMI (Phantolid®) ATII (Traseolide) HHCB (Galaxolid®) Musk NN Musk-ketone
PFAA Mix		
DRE-GS09000218MW	PFAA Mix 100 µg/mL in Methanol:Water 96%:4%(±)	5x1ml
	perfluoro-1-butanesulfonic acid perfluorododecanoic acid perfluorohexanoic acid perfluorooctanesulfonic acid perfluorotetradecanoic acid perfluoroundecanoic acid	perfluorodecanoic acid perfluoroheptanoic acid perfluorononanoic acid pentadecafluorooctanoic acid hydrate perfluorotridecanoic acid

Additional organic contaminants

Product code	Description	
Purgeable Aromatics for Gas.Ident.Mix 3		
DRE-XA06020300ME	Purgeable Arom. for Gas.Ident.Mix 3 200 µg/mL in Methanol	1ml
	1,2-Dichlorobenzene 1,4-Dichlorobenzene Chlorobenzene Methyl-tert-butylether o-Xylene Toluene	1,3-Dichlorobenzene Benzene Ethylbenzene m-Xylene p-Xylene
Purgeable Internal Standards Mix 4		
DRE-YA08260400ME	Purgeable Internal Standards Mix 4 2000 µg/mL in Methanol	1ml
	1,4-Dichlorobenzene-D4 1,4-Difluorobenzene Chlorobenzene D5 Pentafluorobenzene	
SV System Performance Check Mix		
DRE-GA09000384DI	SV System Performance Check Mix 1000 µg/mL in Dichloromethane(‡)	1ml
	2,4-dinitrophenol 4-nitrophenol 2-methyl-4,6-dinitrophenol 2-nitrophenol	
Semivolatile Mix, 10		
DRE-GS09000381DI	Semivolatile Mix 1000 µg/mL in Dichloromethane(‡)	5x1ml
	p-phenylenediamine 2,6-dimethylaniline a,a-dichlorotoluene 2-methoxy-5-methylaniline triethylamine	2,4-dimethylaniline acrylamide p-anisidine phthalic anhydride carbazole
Surrogate Standard Base Neutrals Mix 19		
DRE-YA04001900DI	Surrogate Standard Base Neutrals Mix 19 1000 µg/mL in Dichloromethane	1ml
	2-Fluorobiphenyl Nitrobenzene D5 p-Terphenyl D14	
Surrogate Standard Mix 13		
DRE-YA08241300ME	Surrogate Standard Mix 13 2000 µg/mL in Methanol	1ml
	1,2-Dichloroethane D4 4-Bromofluorobenzene Toluene D8	
Terpene Mix 1		
DRE-GA09000272ME	Terpene Mix 1 100 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000272ME	Terpene Mix 1 100 µg/mL in p/t Methanol(‡)	5x1ml
	(-)-α-Bisabolol (technical grade) (+)-fenchone cineole DL-menthol 3,7-dimethyl-1,3,6-octatriene cis-nerolidol	(-)-caryophyllene oxide (+)-3-carene mycrene isoborneol α-phellandrene
	(-)-Isopulegol camphene farnesene, mixture of isomers linalool (-)-trans-caryophyllene	cedrol camphor geranyl acetate nerol valencene
Terpene Mix 2		
DRE-GA09000273ME	Terpene Mix 2 100 µg/mL in p/t Methanol(‡)	1ml
DRE-GS09000273ME	Terpene Mix 2 100 µg/mL in p/t Methanol(‡)	5x1ml
	(-)-borneol (-)-camphor α-humulene g-terpinene sabinene trans-nerolidol	Borneol (+)-fenchol alpha-terpinene geraniol sabinene hydrate
	(R)-(+)-pulegone (+)-limonene α-pinene (-)-Guaial terpineol, mixed isomers	(+)-camphor α-cedrene β-pinene (-)-fenchone α-terpinolene

Additional organic contaminants

Product code	Description	
Terpene Mix A		
DRE-GA09000155IP	Terpene Mix A 1000 µg/mL in Isopropanol(‡)	1ml
	camphene (+)-3-carene cineole geraniol g-terpinene nerol valencene	
Toxaphene Mix 1		
DRE-ZA22100100CY	Toxaphene Mix 1, 1 µg/mL in Cyclohexane	1ml
	Toxaphene Parlar-No. 26 Toxaphene Parlar-No. 32 Toxaphene Parlar-No. 50 Toxaphene Parlar-No. 62 Toxaphene Parlar-No. 69	
Volatile Aromatic and Unsaturated Organic Compounds Mix 1		
DRE-XA05030100ME	Volatile Aromatic and Unsaturated Organic Compounds Mix 1 200 µg/mL in Methanol	1ml
1,2,3-Trichlorobenzene 1,3,5-Trimethylbenzene 4-Chlorotoluene Ethylbenzene n-Propylbenzene tert-Butylbenzene	1,2,4-Trichlorobenzene 1,3-Dichlorobenzene 4-Isopropyltoluene Hexachloro-1,3-butadiene o-Xylene Tetrachloroethene	1,2,4-Trimethylbenzene 1,4-Dichlorobenzene Benzene m-Xylene p-Xylene Toluene
		1,2-Dichlorobenzene 2-Chlorotoluene Chlorobenzene Naphthalene Styrene Trichloroethene

Index

Analyte name	CAS No#	Page No#
Acenaphthene	83-32-9	237
Acenaphthene D10	15067-26-2	237, 523
Acenaphthylene	208-96-8	237
Acenaphthylene D8	93951-97-4	237, 523
Acephate	30560-19-1	104
Acepromazine Hydrochloride	973-12-6	312
Acepromazine Maleate	3598-37-6	312
Acepromazine-d6 Maleate	1331655-50-5	312, 523
Acequinocyl	57960-19-7	104
Acequinocyl-hydroxy (2-Dodecyl-3-hydroxy-1,4-naphthoquinone)	57960-31-3	104
Acesulfame Potassium	55589-62-3	457
Acetaldehyde-2,4-dinitrophenylhydrazone (DNPH)	1019-57-4	263
Acetamide	60-35-5	263
2-Acetamidofluorene	53-96-3	237
Acetamidomethylphosphonic Acid (Aminomethylphosphonic acid N-acetyl)	57637-97-5	33
Acetamiprid	135410-20-7	104
Acetamiprid D3 (N-methyl D3)	1353869-35-8	104, 523
Acetamiprid-N-desmethyl	190604-92-3	104
Acethion	919-54-0	104
Acetic acid	64-19-7	458
Acetic Acid 2-Ethoxyethyl Ester	111-15-9	457
Acetic Acid sec-Butyl Ester	105-46-4	457
Acetic Acid tert-Butyl Ester	540-88-5	457
Acetic Acid Isopropyl Ester	108-21-4	457
Acetic Acid Linalyl Ester	115-95-7	457
Acetic Acid 3-Methylbutyl Ester	123-92-2	457
Acetic Acid 2-Pentyl Ester	626-38-0	457
Acetic Acid Phenyl Ester	122-79-2	457
Acetoacetic Acid Methyl Ester	105-45-3	263
Acetochlor	34256-82-1	33
Acetochlor D11	1189897-44-6	33, 523
Acetochlor Ethane Sulfonic Acid Sodium Salt	947601-84-5	33

Analyte name	CAS No#	Page No#
Acetochlor Oxalamic Acid (OA)	194992-44-4	33
Acetochlor-2-hydroxy	60090-47-3	33
Acetone	67-64-1	263
Acetone-2,4-dinitrophenylhydrazone	1567-89-1	263
Acetonitrile	75-05-8	560
Acetophenone	98-86-2	263
Acetylacetamide (Acetoacetamide)	5977-14-0	263
4-Acetylamino phenazone	83-15-8	312
2-Acetylfuran	1192-62-7	560
Acetylsalicylic Acid	50-78-2	312
Acetylsulfamethoxazole (N-[4-[(5-Methylisoxazol-3-yl)sulphamoyl]phenyl]acetamide)	21312-10-7	312
Acetylsulfamethoxazole D4	1215530-54-3	312, 523
N-Acetylthiazolidine-4-carboxylic Acid	5025-82-1	172
Acibenzolar Acid (1,2,3-Benzothiadiazole-7-carboxylic acid)	35272-27-6	1
Acibenzolar-S-methyl	135158-54-2	1
Aciclovir (Acyclovir)	59277-89-3	312
Acid Blue 1	129-17-9	400
Acid Blue 3 calcium salt	3536-49-0	400
Acid Blue 62	4368-56-3	400
Acid Green 1 (technical product)	19381-50-1	400
Acid Green 16	12768-78-4	400
Acid Green 50	3087-16-9	400
Acid mordant Black 1	3618-58-4	402
Acid Orange 6 (4-[2-(2,4-Dihydroxyphenyl) diazenyl]benzenesulfonic Acid Sodium Salt)	547-57-9	400
Acid Red 26	3761-53-3	400
Acid Red 2G (E 128)	3734-67-6	400
Acid Red 44	2766-77-0	401
Acid Red 52	3520-42-1	401
Acid Red 73	5413-75-2	401
Acid Red 87 (Eosin Y)	17372-87-1	401

Index

Analyte name	CAS No#	Page No#
Acid Red 95	33239-19-9	401
Acid Violet 17	4129-84-4	401
Acid Violet 49	1694-09-3	401
Acid Violet 9	6252-76-2	401
Acid Yellow 11	6359-82-6	401
Acid Yellow 23 Aluminium lake	12225-21-7	402
Acid Yellow 36	587-98-4	402
Acifluorfen	50594-66-6	33
Acifluorfen Sodium	62476-59-9	33
Aclonifen	74070-46-5	33
Acridine	260-94-6	237
Acriflavine	8048-52-0	312
Acrinathrin	101007-06-1	104
Acrolein (2-Propenal; 2-Propen-1-one)	107-02-8	104
Acrolein-2,4-dinitrophenylhydrazone (DNPH)	888-54-0	263
Acrylamide	79-06-1	263
Acrylic Acid Butyl Ester	141-32-2	560
Acrylic Acid Ethyl Ester (Ethyl Acrylate)	140-88-5	560
Acrylic Acid 2-Ethylhexyl Ester	103-11-7	560
Acrylic Acid Methyl Ester	96-33-3	561
Acrylonitrile	107-13-1	105
Adipic acid	124-04-9	484
Adipic Acid Bis(2-ethylhexyl) Ester	103-23-1	484
Adipic Acid Dibutyl Ester	105-99-7	484
Adipic Acid Dimethyl Ester	627-93-0	484
Adonit	488-81-3	313
Adrenaline (L-Adrenaline; (-)-Epinephrine)	51-43-4	313
DL-Adrenaline Hydrochloride ((±)-Epinephrine hydrochloride)	329-63-5	313
Agaric Acid	666-99-9	313
AHMI (Phantolide)	15323-35-0	560
AHTN (Tonalide)	1506-02-1	560
AHTN (Tonalide) (6-acetyl D3)	1396967-82-0	523, 560
Aklomide	3011-89-0	313

Analyte name	CAS No#	Page No#
Akton (O-[2-Chloro-1-(2,5-dichlorophenyl)vinyl] O,O-diethyl phosphorothioate)	1757-18-2	105
Alachlor	15972-60-8	34
Alachlor D13 (2,6-diethylphenyl D13)	1015856-63-9	34, 523
Alachlor D3 (methoxy D3)	n/a	34, 523
Alachlor ethane sulfonic acid (ESA) sodium salt	140939-15-7	34
Alachlor oxalamic acid (OA)	171262-17-2	34
Alachlor-2-hydroxy	56681-55-1	34
DL-Alanine	302-72-7	561
Alanycarb	83130-01-2	105
Albendazole	54965-21-8	313
Albendazole Sulfone (Methyl [5-Propylsulphonyl]-1H-benzimidazol-2-yl]carbamate)	75184-71-3	313
Albendazole Sulfoxide (Methyl [5-Propylsulphonyl]-1H-benzimidazol-2-yl]carbamate)	54029-12-8	313
Albendazole-2-amino (5-(Propylsulfonyl)-1H-benzimidazol-2-amine)	80983-36-4	313
Albendazole-2-aminosulfone (5-(Propylsulfonyl)-1H-benzimidazol-2-amine)	80983-34-2	314
Albendazole-2-aminosulfoxide	80983-35-3	314
Aldicarb	116-06-3	105
Aldicarb D3 (N-methyl D3)	1795142-83-4	105, 524
Aldicarb-oxime	1646-75-9	105
Aldicarb Sulfone	1646-88-4	105
Aldicarb-sulfone D3 (N-methyl D3)	1795135-15-7	105, 524
Aldicarb Sulfoxide	1646-87-3	105
Aldicarb-sulfoxide-nitrile	14668-28-1	106
Aldimorph	91315-15-0	1
Aldrin (HHDN)	309-00-2	106
Allethrin	584-79-2	106
D-Allethrin	231937-89-6	106
Allidochlor	93-71-0	34
Allopurinol	315-30-0	314
Alloxydim-sodium	55635-13-7	34

Index

Analyte name	CAS No#	Page No#
Allura Red AC (E 129)	25956-17-6	402
Allyl alcohol	107-18-6	264
4-Allylanisole (Estragole)	140-67-0	561
Allylchloride (3-Chloro-1-propene)	107-05-1	383
Allyxycarb	6392-46-7	106
Altrenogest	850-52-2	314
Amantadine Hydrochloride	665-66-7	314
Amaranth	915-67-3	402
Ambroxol Hydrochloride	23828-92-4	314
Ametoctradin	865318-97-4	1
Ametryn	834-12-8	34
Ametryn-N-nitroso	n/a	35
Amicarbazone	129909-90-6	35
Amidate (O,O-Dimethylphosphoramidothioate)	17321-47-0	314
Amidithion	919-76-6	106
Amidoflumet	84466-05-7	172
Amidosulfuron	120923-37-7	35
Amidosulfuron D6 (dimethoxy D6)	n/a	35, 524
Amikacin hydrate	1257517-67-1	314
Aminitrozole (2-Acetamido-5-nitrothiazole)	140-40-9	314
2-Aminoanthracene	613-13-8	237
1-Aminoanthraquione	82-45-1	237
4-Aminoazobenzene	60-09-3	402
2-Aminobenzimidazole	934-32-7	1
4-Aminobenzoic Acid	150-13-0	561
2-Aminobiphenyl	90-41-5	561
3-Aminobiphenyl	2243-47-2	561
4-Aminobiphenyl	92-67-1	561
4-Aminobiphenyl D9	344298-96-0	524, 561
Aminocarb	2032-59-9	106
2-Amino-4-chlorophenol	95-85-2	420
6-Aminochrysene	2642-98-0	237
2-Amino-4,6-dimethoxypyrimidine (4,6-Dimethoxypyrimidin-2-ylamine)	36315-01-2	172
4-Amino-2,3-dimethylazobenzene	97-56-3	402

Analyte name	CAS No#	Page No#
2-Amino-4,6-dinitrotoluene (2-Methyl-3,5-dinitroaniline)	35572-78-2	420
4-Amino-2,6-dinitrotoluene (4-Methyl-3,5-dinitroaniline)	19406-51-0	420
2-Amino-6-ethoxynaphthalene	293733-21-8	238
2-Aminoflubendazole ((2-Amino-1H-benzimidazol-5-yl) (4-fluorophenyl)methanone)	82050-13-3	315
2-Aminofluorene	153-78-6	238
4-Amino-3-fluorophenol	399-95-1	420
1-Aminohydantoin (2,4,5-13C3)	957509-31-8	315, 524
1-Aminohydantoin Hydrochloride	2827-56-7	315
2-Amino-4-methoxy-6-methyl-1,3,5-triazine	1668-54-8	35
2-Amino-3-methyl-3H-imidazo(4,5-f)quinoline	76180-96-6	238
3-Amino-2-methyl-5-nitrobenzamide	3572-44-9	561
5-Amino-2-methyl-3-nitrobenzamide	90223-31-7	562
4-Amino-3-methylphenol	2835-99-6	420
Aminomethylphosphonic Acid (AMPA)	1066-51-9	35
(Aminomethyl) phosphonic Acid 13C 15N (AMPA)	n/a	35, 524
Aminomethylphosphonic Acid-FMOC (AMPA-FMOC)	195306-88-8	35
3-Amino-5-morpholinomethyl-1,3-oxazolidin-2-one (AMOZ)	43056-63-9	315
3-Amino-5-morpholinomethyl-2-oxazolidinone D5 (AMOZ D5)	1017793-94-0	315, 524
1-Aminonaphthalene	134-32-7	238
1-Aminonaphthalene D7	78832-53-8	238, 524
2-Aminonaphthalene	91-59-8	238
2-Aminonaphthalene D7	93951-94-1	238, 524
2-Amino-3-nitrotoluene (2-Methyl-6-nitroaniline)	570-24-1	402
2-Amino-4-nitrotoluene (2-Methyl-5-nitroaniline)	99-55-8	402, 420
2-Amino-5-nitrotoluene (2-Methyl-4-nitroaniline)	99-52-5	403

Index

Analyte name	CAS No#	Page No#
2-Amino-6-nitrotoluene (2-Methyl-3-nitroaniline)	603-83-8	403
3-Amino-2-nitrotoluene (3-Methyl-2-nitroaniline)	601-87-6	420
3-Amino-4-nitrotoluene (5-Methyl-2-nitroaniline)	578-46-1	403
4-Amino-2-nitrotoluene (4-Methyl-3-nitroaniline)	119-32-4	403
4-Amino-3-nitrotoluene (4-Methyl-2-nitroaniline)	89-62-3	403
3-Amino-2-oxazolidinone (AOZ)	80-65-9	315
3-Amino-2-oxazolidinone Sulfate	32957-26-9	315
3-Amino-2-oxazolidinone D4 (AOZ D4)	1188331-23-8	315, 525
Aminophenazone (4,4-Dimethylaminophenazone)	58-15-1	315
2-Aminophenol	95-55-6	420
3-Aminophenol	591-27-5	172
4-Aminophenol	123-30-8	420
4-Aminophenylether (4,4'-Oxydianiline)	101-80-4	403
3-Amino-1-propanol	156-87-6	264
Aminopyralid	150114-71-9	35
2-Aminopyridine (2-Pyridylamine)	504-29-0	316
3-Aminopyridine	462-08-8	562
4-Aminopyridine (Dalfampridine)	504-24-5	172
4-Aminoquinaldine	6628-04-2	238
4-Aminosalicylic Acid	65-49-6	562
4-Aminotoluene-3-sulfonic Acid	88-44-8	562
Amiprofos-methyl	36001-88-4	36
Amisulbrom	348635-87-0	1
Amitraz	33089-61-1	106
Amitraz D12 (methylphenyl D12)	n/a	106, 525
Amitrole	61-82-5	36
Amitrole (1-15N)	367498-28-0	36, 525
Amitrole 1-15N 5-13C	n/a	36, 525
Amlodipine Besilate	111470-99-6	316
Ammelide	645-93-2	562
Amobam	3566-10-7	1
Amoxicillin Trihydrate	61336-70-7	316

Analyte name	CAS No#	Page No#
Amphotericin B	1397-89-3	316
Ampicillin Trihydrate	7177-48-2	316
Amprolium Hydrochloride	137-88-2	316
Ampyrone	83-07-8	316
D-Amygdalin	29883-15-6	316
tert-Amyl Methyl Ether	994-05-8	264
α-Amylcinnamal	122-40-7	562
α-Amylcinnamyl Alcohol	101-85-9	562
tert-Amyl-ethyl Ether	919-94-8	264
4-tert-Amylphenol	80-46-6	421
rac-Anatabine	2743-90-0	562
Ancymidol	12771-68-5	172
5α-Androstane	438-22-2	316
Androstanolone (5α-Androstan-17β-ol-3-one)	521-18-6	317
Androstenedione (4-Androstene-3,17-dione)	63-05-8	317
Androsterone	53-41-8	317
Anilazine	101-05-3	1
Aniline	62-53-3	403
Aniline D5	4165-61-1	403, 525
Aniline-2-sulfonic Acid	88-21-1	562
Anilofos	64249-01-0	36
Anise Alcohol	105-13-5	563
p-Anisic acid (4-Methoxybenzoic acid)	100-09-4	563
2-Anisidine (2-Methoxyaniline)	90-04-0	563
3-Anisidine (3-Methoxyaniline)	536-90-3	563
4-Anisidine (4-Methoxyaniline)	104-94-9	563
Anthanthrene	191-26-4	238
Anthracene	120-12-7	239
Anthracene D10	1719-06-8	239, 525
Anthranilic Acid Isopropylamide	30391-89-0	172
Anthraquinone	84-65-1	172
ANTU	86-88-4	172
8'-Apoaldehyde (Apocarotenal)	1107-26-2	317
Apramycin Sulfate	65710-07-8	317
L-Arabinitol (L-(-)-Arabit)	7643-75-6	458

Index

Analyte name	CAS No#	Page No#
L(+)-Arabinose	5328-37-0	458
D(+)-Arabitol	488-82-4	458
Aramite	140-57-8	173
DL-Arginine Hydrochloride	32042-43-6	563
Aroclor 1016	12674-11-2	214
Aroclor 1221	11104-28-2	214
Aroclor 1232	11141-16-5	215
Aroclor 1242	53469-21-9	215
Aroclor 1248	12672-29-6	215
Aroclor 1260	11096-82-5	215
4-Arsanilic Acid	98-50-0	317
Arylpyrazol (ACE5599; 2-Chloro-N-[1-[2,6-dichloro-4-(difluoromethyl)phenyl]-4-nitro-1H-pyrazol-5-yl]propanamide)	121750-17-2	36
β-Asarone	5273-86-9	317
Ascorbic acid calcium salt	5743-27-1	458
L-Ascorbic acid 2-monophosphate tri-cyclohexylammonium salt	65718-26-5	458
L-(+)-Ascorbic Acid	50-81-7	458
Ascorbyl Palmitate	137-66-6	458
DL-Asparagine Hydrate	69833-18-7	563
Aspartame	22839-47-0	458
DL-Aspartic Acid	617-45-8	563
Astaxanthin	472-61-7	459
ASTM D3606 2-Butanol	78-92-2	265, 495
Asulam	3337-71-1	36
Atenolol	29122-68-7	317
ATII (Traseolid)	68140-48-7	560
Atorvastatin Calcium	134523-03-8	317
Atraton	1610-17-9	36
Atrazine	1912-24-9	36, 37
Atrazine 13C3 (ring 13C3)	1443685-80-0	37, 525
Atrazine D5 (ethyl-D5)	163165-75-1	37, 525
Atrazine-desethyl	6190-65-4	37
Atrazine-desethyl D6 (dimethyl D6)	1216649-31-8	37, 525
Atrazine-desethyl-desisopropyl	3397-62-4	37

Analyte name	CAS No#	Page No#
Atrazine-desethyl-desisopropyl-2-hydroxy (Ammelide)	645-92-1	37
Atrazine-desethyl-2-hydroxy	19988-24-0	37
Atrazine-desisopropyl (6-Chloro-N-ethyl-1,3,5-triazine-2,4-diamine)	1007-28-9	37
Atrazine-desisopropyl D5 (ethylamino D5)	1189961-78-1	38, 526
Atrazine-desisopropyl-2-hydroxy	7313-54-4	38
Atrazine-2-hydroxy	2163-68-0	38
Atrazine-2-hydroxy D5 (ethyl D5)	1276197-25-1	38, 526
Atrazine-mercapturate	138722-96-0	38
Atropine	51-55-8	318
Avermectin B1 (Abamectine)	71751-41-2	107
Aviglycine Hydrochloride	55720-26-8	173
Azaconazole	60207-31-0	1
Azamethiphos	35575-96-3	107
Azaperol	2804-05-9	318
Azaperone	1649-18-9	318
Azaperone D4	1173021-72-1	318, 526
Azelaic Acid (1,7-Heptanedicarboxylic acid)	123-99-9	459
Azidithion	78-57-9	107
Azimsulfuron	120162-55-2	38
Azinphos-ethyl (Guthion Ethyl)	2642-71-9	107
Azinphos-ethyl D10 (ethyl D10)	n/a	107, 526
Azinphos-methyl	86-50-0	107
Azinphos-methyl D6 (dimethyl D6)	n/a	107, 526
Azinphos-methyl-oxon	961-22-8	107
Aziprotryne	4658-28-0	38
Azithromycin dihydrate	117772-70-0	318
Azobenzene (1,2-Diphenyldiazene)	103-33-3	173
Azocyclotin	41083-11-8	108
Azorubine (E122)	3567-69-9	403
Azoxystrobin	131860-33-8	2
Azoxystrobin (free acid)	1185255-09-7	2
Bacitracin	1405-87-4	318
Bambuterol hydrochloride	81732-46-9	318

Index

Analyte name	CAS No#	Page No#
Barban	101-27-9	38
Basic Blue 7	2390-60-5	404
Basic Green 1	633-03-4	404
Basic Red 1	989-38-8	404
Basic Red 2	477-73-6	404
Basic Red 9	569-61-9	404
Basic Violet 1	603-47-4	404
Basic Violet 3 (Methylrosanilinium Chloride)	548-62-9	404
Basic Violet 3 D6	1266676-01-0	404, 526
Beclometasone (9-Chloro-11β,17,21-trihydroxy-16β-methylpregna-1,4-diene-3,20-dione)	4419-39-0	318
Beflubutamid	113614-08-7	38
Benalaxyl	71626-11-4	2
Benalaxyl-M	98243-83-5	2
Benazolin	3813-05-6	39
Benazolin Ethyl Ester	25059-80-7	39
Bendiocarb	22781-23-3	108
Bendiocarb-phenol	22961-82-6	108
Benfluralin	1861-40-1	39
Benfuracarb	82560-54-1	108
Benfuresate	68505-69-1	39
Benodanil	15310-01-7	2
Benomyl	17804-35-2	2
Benoxacor	98730-04-2	173
Bensulfuron-methyl	83055-99-6	39
Bensulide	741-58-2	39
Bensultap	17606-31-4	108
Bentazone	25057-89-0	39
Bentazone (isopropyl-1,1,1,3,3,3) D6	n/a	39, 526
Bentazone-8-hydroxy	60374-43-8	39
Bentazone-methyl (N-Methylbentazone)	61592-45-8	40
Benthiavalicarb-isopropyl	177406-68-7	2
Bentranil	1022-46-4	40
Benz(c)acridine	225-51-4	239

Analyte name	CAS No#	Page No#
Benz[a]anthracene	56-55-3	239
Benzaldehyd-2,4-dinitrophenylhydrazone	1157-84-2	421
Benzaldehyde	100-52-7	421
Benzalkonium Chloride	8001-54-5	2
Benzamide	55-21-0	563
Benzanilide	93-98-1	564
Benzene	71-43-2	264
Benzene D6	1076-43-3	264, 526
Benzene-1,2-diamine (1,2-Phenylenediamine)	95-54-5	421
Benzene-1,2-dicarboxylic Acid (Phthalic Acid)	88-99-3	564
Benzene-1,2-diol (Catechol; 1,2-Dihydroxybenzene)	120-80-9	421
Benzene-1,2,3-triol (Pyrogallol)	87-66-1	564
4,4'-Benzidine	92-87-5	404
Benzil	134-81-6	564
Benzimidazole	51-17-2	239
Benzine 40-60/80°	64742-49-0	264
1,2-Benzisothiazol-3(2H)-one	2634-33-5	173
Benzo[a]anthracene D12	1718-53-2	239, 526
Benzo[a]fluoranthene	203-33-8	239
11H-Benzo[a]fluorene (1,2-Benzofluorene)	238-84-6	239
Benzo[a]pyrene	50-32-8	240
Benzo[a]pyrene D12	63466-71-7	240, 527
Benzo[b]chrysene	214-17-5	240
Benzo[b]fluoranthene	205-99-2	240
Benzo[b]fluoranthene D12	93951-98-5	240, 527
11H-Benzo[b]fluorene	243-17-4	240
Benzo(b)furan	271-89-6	240
Benzo(b)naphtho(1,2-d)thiophene	205-43-6	240
Benzo[b]naphtho[2,1-d]thiophene	239-35-0	241
7H-Benzo[c]fluorene	205-12-9	241
Benzo[c]phenanthrene	195-19-7	241
Benzo[e]pyrene	192-97-2	241
Benzo(e)pyrene D12	205440-82-0	241, 527
Benzo[g,h,i]fluoranthene	203-12-3	241

Index

Analyte name	CAS No#	Page No#
Benzo[g,h,i]perylene	191-24-2	241
Benzo[g,h,i]perylene D12	93951-66-7	241, 527
Benzo[j]fluoranthene	205-82-3	242
Benzo[k]fluoranthene	207-08-9	242
Benzo[k]fluoranthene D12	93952-01-3	242, 527
Benzobicyclon	156963-66-5	40
Benzocaine (4-Aminobenzoic acid ethyl ester)	94-09-7	564
Benzofenap	82692-44-2	40
Benzoguanamine (6-Phenyl-1,3,5-triazine-2,4-diamine)	91-76-9	564
Benzoic Acid	65-85-0	459
Benzonitrile	100-47-0	564
Benzophenone	119-61-9	564
1,4-Benzoquinone	106-51-4	564
1H-Benzotriazole	95-14-7	242
Benzoximate	29104-30-1	108
5-Benzoyl-4-hydroxy-2-methoxybenzenesulfonic Acid	4065-45-6	565
Benzoylprop Ethyl Ester (Benzoylprop-ethyl)	22212-55-1	40
Benzthiazuron	1929-88-0	40
Benzyl Alcohol	100-51-6	565
Benzyl Benzoate	120-51-4	459
Benzylamine	100-46-9	264
6-Benzylaminopurine	1214-39-7	173
Benzylchloride (α-Chlorotoluene)	100-44-7	421
2-Benzyl-4-chlorophenol	120-32-1	421
Benzylpenicillin Potassium (Penicilline G potassium salt)	113-98-4	318
Berberine Hemisulfate	316-41-6	319
Betacarotene	7235-40-7	459
Betamethasone	378-44-9	319
Betamethasone Acetate	987-24-6	319
Betamethasone 17,21-Dipropionate	5593-20-4	319
Betamethasone 17-Valerate	2152-44-5	319
Bethoxazin	163269-30-5	2
Bezafibrate	41859-67-0	319

Analyte name	CAS No#	Page No#
Bezafibrate D6 (dimethyl D6)	1219802-74-0	319, 527
Bifenazate	149877-41-8	108
Bifenazate-diazene	149878-40-0	173
Bifenox	42576-02-3	40
Bifenox (free acid)	53774-07-5	40
Bifenthrin	82657-04-3	108
Bifenthrin alcohol metabolite ((2-Methyl-3-biphenyl)methanol)	76350-90-8	109
Bilanafos sodium	71048-99-2	40
Binapacryl	485-31-4	3
1,1'-Binaphthyl	604-53-5	242
2,2'-Binaphthyl	612-78-2	242
Bioallethrin ((RS)-3-Allyl-2-methyl-4-oxocyclopent-2-enyl (1R)-trans-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropanecarboxylate)	260359-57-7	109
S-Bioallethrin (Esbiol)	28434-00-6	109
Bioresmethrin	28434-01-7	109
D-(+)-Biotin	58-85-5	459
Biphenyl	92-52-4	3, 565
Biphenyl D10	1486-01-7	3, 527, 527, 565
2,2'-Bipyridyl	366-18-7	565
4,4'-Bipyridyl (4,4'-Bipyridine)	553-26-4	565
Alpha-Bisabolol (α-Bisabolol)	515-69-5	565
Bisacodyl	603-50-9	319
Bis-(2-chloroethoxy)-methane	111-91-1	565
Bis(2-chloroethyl) Ether	111-44-4	109
Bis(chloromethyl) Ether	542-88-1	383
Bis-(2-chloro-1-methylethyl)ether	108-60-1	173
Bis(4-chlorophenyl)acetic Acid (4,4'-DDA)	83-05-6	109
2,2-Bis(4-chlorophenyl)ethanol (DDOH)	2642-82-2	109
Bis(4-chlorophenyl)methane (4,4'-DDM)	101-76-8	109
N,N'-Bis-(3,4-dichlorophenyl)urea	4300-43-0	41
Bis(4-dimethylaminophenyl) methane	101-61-1	565
Bis(2-ethylhexyl) phosphate	298-07-7	566

Index

Analyte name	CAS No#	Page No#
Bis(4-glycidyloxyphenyl)methane	2095-03-6	566
2,2-Bis-(4-glycidyloxyphenyl) propane	1675-54-3	566
1,3-Bis(hydroxymethyl)-2-imidazolidinone	136-84-5	566
N,N'-Bis-(4-nitrophenyl)urea D8	1156508-87-0	527, 566
N,N'-Bis-(4-nitrophenyl)urea	587-90-6	319
Bisotrizole	103597-45-1	320
Bisoprolol Fumarate	104344-23-2	320
Bisphenol A	80-05-7	566
Bisphenol A Bis(2,3-dihydroxypropyl) Ether	5581-32-8	566
Bisphenol A D14	120155-79-5	527, 566
Bisphenol A D16	96210-87-6	528, 566
Bisphenol B	77-40-7	567
Bisphenol F	620-92-8	567
Bisphenol S	80-09-1	567
Bispyribac-sodium	125401-92-5	41
Bitertanol	55179-31-2	3
Bithionol	97-18-7	3
Bixafen	581809-46-3	3
Bixafen-desmethyl	1655498-06-8	3
Boldenone	846-48-0	320
Boldenone 17-Acetate	2363-59-9	320
Boscalid (2-Chloro-N-(4'-chlorobiphenyl-2-yl)nicotinamide)	188425-85-6	3
Brilliant Black BN	2519-30-4	405
Brilliant Blue FCF (E133)	3844-45-9	405
Brodifacoum	56073-10-0	173
Bromacil	314-40-9	41
Bromacil D3 (methyl D3)	n/a	41, 528
Bromadiolone	28772-56-7	174
Brombuterol Hydrochloride	21912-49-2	320
Bromchlorbuterol Hydrochloride (1-(4-Amino-3-bromo-5-chlorophenyl)-2-[(1,1-dimethylethyl)amino]-ethanol hydrochloride)	78982-84-0	320
Bromethalin	63333-35-7	174

Analyte name	CAS No#	Page No#
Bromfenvinfos (Bromfenvinphos)	33399-00-7	109
Bromfenvinfos-methyl (2-Bromo-1-(2,4-dichlorophenyl)vinyl dimethyl phosphate)	13104-21-7	110
Bromoacetic Acid	79-08-3	567
Bromoacetic Acid Methyl Ester	96-32-2	567
4-Bromoaniline	106-40-1	174
Bromobenzene	108-86-1	421
2-Bromobutanoic Acid (2-Bromobutyric Acid)	80-58-0	495
Bromobutide	74712-19-9	41
Bromochloroacetic Acid	5589-96-8	567
Bromochloroacetic Acid Methyl Ester	20428-74-4	567
Bromochloroacetonitrile	83463-62-1	383
4-Bromo-3-chloroaniline	21402-26-6	421
Bromochlorodifluoromethane	353-59-3	383
1-Bromo-2-chloroethane	107-04-0	383
Bromochloromethane	74-97-5	383
Bromochlorophene	15435-29-7	567
4-Bromo-2-chlorophenol	3964-56-5	174
1-Bromo-3-chloropropane	109-70-6	383
2-Bromo-1-chloropropane	3017-95-6	383
Bromocyclen	1715-40-8	110
Bromodichloroacetic acid	71133-14-7	567
Bromodichloromethane	75-27-4	383
4-Bromo-3,5-dimethylphenyl-N-methylcarbamate (BDMC)	672-99-1	568
4-Bromodiphenyl Ether (PBDE 3)	101-55-3	568
Bromoethane	74-96-4	384
2-Bromoethanol	540-51-2	384
Bromofenoxim	13181-17-4	41
1-Bromo-2-fluorobenzene	1072-85-1	384
4-Bromofluorobenzene (BFB)	460-00-4	422
Bromoform (Tribromomethane)	75-25-2	384
1-Bromo-2-nitrobenzene	577-19-5	422
2-Bromophenol	95-56-7	422
3-Bromophenol	591-20-8	422
Bromophos (Bromophos-methyl)	2104-96-3	110

Index

Analyte name	CAS No#	Page No#
Bromophos-ethyl	4824-78-6	110
Bromophos-ethyl-oxon	4855-63-4	110
Bromophos-methyl-oxon	4855-62-3	110
1-Bromopropane	106-94-5	384
2-Bromopropane	75-26-3	384
1-Bromopropan-2-ol	19686-73-8	384
2-Bromopropionic Acid	598-72-1	568
Bromopropylate	18181-80-1	110
Bromopropylate D7 (isopropyl D7)	n/a	110, 528
Bromothalonil (1,2-Dibromo-2,4-dicyanobutane)	35691-65-7	568
Bromotrichloromethane	75-62-7	384
Bromoxynil	1689-84-5	41
Bromoxynil Heptanoate	56634-95-8	41
Bromoxynil Methyl Ether	3336-39-8	41
Bromoxynil Octanoate	1689-99-2	42
Brompheniramine Maleate	980-71-2	320
Brompyrazon	3042-84-0	42
Bromuconazole	116255-48-2	3
Bromuron	3408-97-7	42
Bronopol	52-51-7	174
Budesonide	51333-22-3	320
Bufencarb	8065-36-9	111
Bufexamac	2438-72-4	320
Buminafos	51249-05-9	174
Bupirimate	41483-43-6	4
(EZ)-Buprofezin	69327-76-0	111
Butachlor	23184-66-9	42
Butachlor ethane sulfonic acid (ESA) calcium salt	n/a	42
1,3-Butadiene	106-99-0	264
Butafenacil	134605-64-4	42
Butamifos	36335-67-8	42
1,4-Butanediamine (1,4-Diaminobutane)	110-60-1	265
1,2-Butanediol	584-03-2	265
Butane-1,3-diol (Butylene Glycol)	107-88-0	265
1,4-Butanediol	110-63-4	265

Analyte name	CAS No#	Page No#
2,3-Butanediol	513-85-9	265
1-Butanethiol	109-79-5	265
Butanetriol	42890-76-6	568
1-Butanol	71-36-3	265
tert.-Butanol	75-65-0	265
2-Butanone	78-93-3	266
2-Butanone-2,4-dinitrophenylhydrazone	958-60-1	266
Butocarboxim	34681-10-2	111
Butocarboxim Sulfoxide	34681-24-8	111
Butoxycarboxim	34681-23-7	111
2-(2-Butoxyethoxy)ethyl Acetate (Diethylene glycol monobutylether acetate)	124-17-4	266
1-Butoxy-2-propanol	5131-66-8	266
Butralin	33629-47-9	42
Butroxydim	138164-12-2	42
Buturon	3766-60-7	43
Butyl Acetate	123-86-4	266
tert-Butyl Ethyl Ether	637-92-3	266
Butyl Hexanoate (n-Caproic acid n-butyl ester)	626-82-4	266
Butyl Methacrylate (Methacrylic acid butyl ester)	97-88-1	266
Butyl Parahydroxybenzoate (4-Hydroxybenzoic acid n-butyl ester)	94-26-8	459
3-[N-n-Butyl-N-acetyl] aminopropionic acid ethyl ester	52304-36-6	568
1-Butylamine	109-73-9	266
2-Butylamine	13952-84-6	267
Butylate	2008-41-5	43
n-Butylbenzene	104-51-8	267
sec-Butylbenzene	135-98-8	267
tert-Butylbenzene	98-06-6	267
N-n-Butylbenzenesulfonamide	3622-84-2	568
4-tert-Butyl-2-chlorophenol	98-28-2	422
2-sec-Butyl-1-(decyloxy)-4-tritylbenzene	1404190-37-9	267
6-tert-Butyl-2,4-dimethylphenol (2-tert-Butyl-4,6-dimethylphenol)	1879-09-0	422

Index

Analyte name	CAS No#	Page No#
2-Butyl-2-ethyl-1,3-propanediol	115-84-4	267
tert-Butylhydroquinone	1948-33-0	459
Butylhydroxyanisole	25013-16-5	459
3-tert-Butyl-4-hydroxyanisole (2-tert-Butyl-4-methoxyphenol; BHA)	121-00-6	422
Butylhydroxytoluene	128-37-0	460
Butylhydroxytoluene-d21 (BHT-d21; 2,6-Bis[1,1-di(methyl-d3)ethyl-2,2,2-d3]-4-methylphen-3,5-d2-ol-d)	64502-99-4	460, 528
2-tert-Butylphenol	88-18-6	422
4-Butylphenol	1638-22-8	423
4-tert-Butylphenol	98-54-4	423
3-(4-tert-Butyl(phenyl) isobutylaldehyde	80-54-6	267
Butyltin Trichloride	1118-46-3	568
4-tert-Butyltoluene	98-51-1	267
Butyraldehyde (Butanal)	123-72-8	267
Butyraldehyde-2,4-dinitrophenylhydrazone	1527-98-6	268
Butyric Acid	107-92-6	460
Butyric Acid Butyl Ester	109-21-7	460
Butyric Acid Ethyl Ester	105-54-4	460
Butyric Acid Methyl Ester	623-42-7	460
Butyric Acid Propyl Ester	105-66-8	460
n-Butyrophenone	495-40-9	268
Cacodylic Acid	75-60-5	43
Cacodylic Acid Sodium Salt Trihydrate	6131-99-3	43
Cadaverine dihydrochloride	1476-39-7	321
Cadusafos	95465-99-9	111
Cadusafos D5 (ethyl-D5)	1215380-71-4	111, 528
Cafenstrole	125306-83-4	43
Caffeic Acid (3,4-Dihydroxycinnamic Acid)	331-39-5	321
Caffeine	58-08-2	321
Calcifediol Monohydrate (Vitamin D3 25-hydroxy monohydrate)	63283-36-3	460
Calflo E (Calcium Silicate)	1344-95-2	568

Analyte name	CAS No#	Page No#
Cambendazol	26097-80-3	321
Camphene	79-92-5	569
Canthaxanthine	514-78-3	460
Caprylic Acid (Octanoic acid)	124-07-2	461
Capsaicin	404-86-4	461
Captafol	2425-06-1	4
Captan	133-06-2	4
Captan-4,4,5,6,7,7-D6	1330190-00-5	4, 528
Carazolol	57775-29-8	321
Carazolol D7	1173021-02-7	321, 528
Carbadox	6804-07-5	321
Carbamazepine	298-46-4	321
Carbanolate (Banol)	671-04-5	111
Carbaryl (1-Naphthol N-methylcarbamate)	63-25-2	111
Carbaryl D3 (methyl D3)	1433961-56-8	112, 528
Carbaryl D7 (naphthyl D7)	362049-56-7	112, 528
Carbazochrome (Adenochrome semicarbazone)	69-81-8	321
Carbazole	86-74-8	242
Carbazole D8	38537-24-5	242, 529
Carbendazim	10605-21-7	4
Carbendazim D3 (methyl D3)	1255507-88-0	4, 529
Carbendazim D4 (ring D4)	291765-95-2	4, 529
Carbenicillin Sodium	4800-94-6	322
Carbetamide	16118-49-3	43
Carbetamide D5 (phenyl D5)	n/a	43, 529
Carbofuran	1563-66-2	112
Carbofuran D3 (N-methyl D3)	1007459-98-4	112, 529
Carbofuran-3-hydroxy (3-Hydroxycarbofuran)	16655-82-6	112
Carbofuran-3-keto (3-Ketocarbofuran)	16709-30-1	112
Carbofuranphenol (2,2-Dimethyl-3H-benzofuran-7-ol)	1563-38-8	112
Carbon Disulfide	75-15-0	112

Index

Analyte name	CAS No#	Page No#
Carbontetrachloride (Tetrachloromethane)	56-23-5	384
Carbophenothion	786-19-6	112
Carbophenothion D10 (di(ethyl D5))	n/a	113, 529
Carbophenothion-methyl	953-17-3	113
Carbophenothion-methyl sulfone	62059-34-1	113
Carbophenothion-methyl sulfoxide	62059-33-0	113
Carbophenothion-oxon	7173-84-4	113
Carbophenothion-sulfone	16662-85-4	113
Carbophenothion-sulfoxide	17297-40-4	113
Carbosulfan	55285-14-8	113
Carboxin	5234-68-4	4
(1S)-(+)-3-Carene	498-15-7	569
3-Carene	13466-78-9	569
Carfentrazone (free acid)	128621-72-7	43
Carfentrazone-ethyl	128639-02-1	43
Carminic acid (E120)	1260-17-9	405
Carprofen	53716-49-7	322
Carpropamid	104030-54-8	5
Cartap	15263-53-3	113
Cartap Hydrochloride	15263-52-2	114
Carvone	99-49-0	174
Caryophyllene Oxide (β-Caryophyllene Epoxide)	1139-30-6	569
Beta-Caryophyllene (Isocaryophyllene)	118-65-0	569
Cefadroxil	50370-12-2	322
Cefalexin Monohydrate	23325-78-2	322
Cefazolin Sodium	27164-46-1	322
Cefoperazone Sodium	62893-20-3	322
Cefotaxime Sodium	64485-93-4	322
Cefpirome Sulfate	98753-19-6	322
Cefquinome Sulfate	118443-89-3	322
Cefradine	38821-53-3	323
Ceftiofur	80370-57-6	323
Ceftiofur Hydrochloride	103980-44-5	323

Analyte name	CAS No#	Page No#
Cefuroxime Sodium	56238-63-2	323
Cekafix	121227-99-4	174
Celestolide (ABDI)	13171-00-1	569
D(+)-Cellobiose	528-50-7	461
Cetrimide	1119-97-7	323
Cetylpyridinium Chloride Monohydrate	6004-24-6	323
Chinomethionat	2439-01-2	114
Chlomethoxyfen	32861-85-1	44
Chloral Hydrate	302-17-0	323
Chloralose	15879-93-3	174
β-Chloralose	16376-36-6	175
Chloramben	133-90-4	44
Chloramben Methyl Ester	7286-84-2	44
Chloramphenicol	56-75-7	323
Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2R)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)	202480-68-0	323, 529
erythro-Chloramphenicol D5 (ring D4, benzyl D) (2,2-Dichloro-N-[(1R,2S)-1,3-dihydroxy-1-(4-nitrophenyl)propan-2-yl]acetamide D5)	1426174-26-6	324, 529
Chloramphenicol Palmitate	530-43-8	324
Chloraniformethan	20856-57-9	5
4-Chloranil	118-75-2	5
Chloranocryl	2164-09-2	44
Chlorantraniliprole	500008-45-7	114
Chlorbenside	103-17-3	114
Chlorbenside Sulfone	7082-99-7	114
Chlorbicyclen	2550-75-6	114
Chlorbromuron	13360-45-7	44
Chlorbufam	1967-16-4	44
Chlordane	57-74-9	114
cis-Chlordane (α)	5103-71-9	114
trans-Chlordane (γ)	5103-74-2	114
oxy-Chlordane	27304-13-8	115
Chlordecone	143-50-0	115

Index

Analyte name	CAS No#	Page No#
Chlordecone-5b-hydro	53308-47-7	115
Chlordene	3734-48-3	115
β-Chlordene	56534-03-3	115
cis-Chlordene	56534-02-2	115
oxy-Chlordene	6058-23-7	115
Chlordimeform	6164-98-3	115
Chlordimeform Hydrochloride	19750-95-9	115
Chlorfenac	85-34-7	44
Chlorfenapyr	122453-73-0	116
Chlorfenethol	80-06-8	116
Chlorfenprop-methyl	14437-17-3	44
Chlorfenson	80-33-1	116
Chlorfenvinphos	470-90-6	116
Chlorfenvinphos D10 (ethyl D10)	1346606-54-9	116, 530
cis-Chlorfenvinphos	18708-87-7	116
Chlorfluazuron	71422-67-8	116
Chlorflurenol	2464-37-1	44
Chlorflurenol Methyl Ester	2536-31-4	45
Chlorhexidine	55-56-1	324
Chlorhexidine Diacetate	56-95-1	324
Chlorhexidine Dihydrochloride	3697-42-5	324
Chloridazon	1698-60-8	45
Chloridazon D5 (phenyl D5)	1346818-99-4	45, 530
Chloridazon-desphenyl	6339-19-1	45
Chloridazon-methyl-desphenyl (5-Amino-4-chloro-2-methylpyridazin-3-one)	17254-80-7	45
Chlorimuron-ethyl	90982-32-4	45
Chlormadinone 17-Acetate	302-22-7	324
Chlormephos	24934-91-6	116
Chlormequat-chloride	999-81-5	175
Chlormequat-chloride 1,1,2,2-D4	n/a	175, 530
Chlornitrofen	1836-77-7	45
Chlornitrofen-amino (4-(2,4,6-Trichlorophenoxy)aniline)	26306-61-6	45
Chloroacetaldehyde	107-20-0	385
4-Chloroacetanilide	539-03-7	569

Analyte name	CAS No#	Page No#
Chloroacetic Acid	79-11-8	45, 569
Chloroacetic Acid Methyl Ester	96-34-4	569
N-Chloroacetyl-2,6-diethylaniline	6967-29-9	175
1-cis-3-Chloroallyl-3,5,7-triaza-1-azonia-adamantane Chloride	51229-78-8	570
2-Chloroaniline	95-51-2	423
3-Chloroaniline	108-42-9	175, 423
4-Chloroaniline	106-47-8	175
2-Chloroanisole	766-51-8	423
3-Chloroanisole	2845-89-8	423
4-Chloroanisole	623-12-1	423
Chlorobenzene	108-90-7	423
Chlorobenzene D5	3114-55-4	423, 530
Chlorobenzilate	510-15-6	116
2-Chlorobenzoic Acid	118-91-2	570
3-Chlorobenzoic Acid	535-80-8	570
4-Chlorobenzoic Acid	74-11-3	570
2-Chlorobenzonitrile	873-32-5	570
4-Chlorobenzotrifluoride (4-Chloro-α,α,α-trifluorotoluene)	98-56-6	424
6-Chloro-2-benzoxazolinone	19932-84-4	175
Chlorobenzuron	57160-47-1	117
1-Chlorobutane	109-69-3	385
Chlorobutanol Hemihydrate (Acetone chloroform)	6001-64-5	385
4-Chlorocatechol	2138-22-9	424
2-Chloro-5-(chloromethyl)pyridine	70258-18-3	175
1-Chlorodecane	1002-69-3	268
1-Chloro-1,1-difluoroethane	75-68-3	385
Chlorodifluoromethane (Freon 22)	75-45-6	385
6-Chloro-2,4-dinitroaniline	3531-19-9	424
1-Chloro-2,4-dinitrobenzene	97-00-7	424
Chloroethane	75-00-3	385
2-Chloroethyl-vinylether	110-75-8	385
3-Chloro-4-fluoroaniline	367-21-5	424
Chloroform	67-66-3	385
Chlorogenic Acid	327-97-9	461
4-Chloroguaiacol	16766-30-6	424

Index

Analyte name	CAS No#	Page No#
1-Chlorohexane (Hexylchloride)	544-10-5	385
3-Chloro-5-hydroxy-6-phenylpyridazine (6-Chloro-4-hydroxy-3-phenyl-pyridazine)	40020-01-7	175
4-Chloro-2-isopropyl-5-methylphenol	89-68-9	424
Chloromethane (Methylchloride)	74-87-3	386
3-Chloro-4-methoxyaniline	5345-54-0	176
Chloromethyl Methyl Ether	107-30-2	386
2-Chloro-4-methylaniline	615-65-6	424
2-Chloro-5-methylaniline	95-81-8	424
2-Chloro-6-methylaniline	87-63-8	425
3-Chloro-2-methylaniline	87-60-5	425
3-Chloro-4-methylaniline	95-74-9	176
4-Chloro-2-methylaniline (4-Chloro-o-toluidine)	95-69-2	176
5-Chloro-2-methylaniline	95-79-4	425
5-Chloro-2-methyl-3(2H)-isothiazolone	26172-55-4	176
2-Chloro-4-methylphenol	6640-27-3	425
2-Chloro-5-methylphenol	615-74-7	425
2-Chloro-6-methylphenol (6-Chloro-o-cresol)	87-64-9	425
4-Chloro-2-methylphenol	1570-64-5	176
4-Chloro-3-methylphenol (Chlorocresol)	59-50-7	425
1-(3-Chloro-4-methylphenyl)-3-methyl-urea	22175-22-0	46
1-(3-Chloro-4-methylphenyl)urea	13142-64-8	46
2-Chloro-4-methylsulfonylbenzoic acid	53250-83-2	570
1-Chloronaphthalene	90-13-1	243
2-Chloronaphthalene	91-58-7	243
Chloroneb	2675-77-6	5
2-Chloronicotinic Acid	2942-59-8	324
6-Chloronicotinic Acid	5326-23-8	176
2-Chloro-4-nitroaniline	121-87-9	425
4-Chloro-2-nitroaniline	89-63-4	425
1-Chloro-2-nitrobenzene	88-73-3	426
1-Chloro-3-nitrobenzene	121-73-3	426
1-Chloro-4-nitrobenzene	100-00-5	426

Analyte name	CAS No#	Page No#
2-Chloro-4-nitrophenol	619-08-9	176
2-Chloro-4-nitrotoluene	121-86-8	426
2-Chloro-5-nitrotoluene	13290-74-9	426
2-Chloro-6-nitrotoluene	83-42-1	426
3-Chloro-4-nitrotoluene	38939-88-7	426
4-Chloro-2-nitrotoluene	89-59-8	426
4-Chloro-3-nitrotoluene	89-60-1	426
5-Chloro-2-nitrotoluene	5367-28-2	427
1-Chlorooctadecane	3386-33-2	268
1-Chlorooctane	111-85-3	268
2'-Chloro-2,3,4,5,6-pentadeuterio-1,1'-biphenyl	51624-35-2	215, 530
Chloropentafluoroethane	76-15-3	386
Chlorophacinone	3691-35-8	176
3-Chlorophenol Acetate	13031-39-5	571
4-Chlorophenol Acetate	876-27-7	571
2-Chlorophenol-3,4,5,6-D4	93951-73-6	427, 530
2-Chlorophenol	95-57-8	427
3-Chlorophenol	108-43-0	427
4-Chlorophenol	106-48-9	176, 427
(4-Chlorophenoxy)acetic Acid	122-88-3	177
2-Chlorophenoxyacetic Acid	614-61-9	571
3-Chlorophenoxyacetic Acid	588-32-9	571
2-(4-Chlorophenoxy)propionic Acid	3307-39-9	571
4-Chlorophenyl Phenyl Ether	7005-72-3	571
4-Chlorophenylacetic Acid	1878-66-6	572
1-(2-Chlorophenyl)-2-thiourea	5344-82-1	572
(2-Chlorophenyl)trichlormethane (α,α,α,2-Tetrachlorotoluene)	2136-89-2	427
(4-Chlorophenyl)urea	140-38-5	46
6-Chloro-2-picolinic Acid	4684-94-0	177
Chloroprene (2-Chloro-1,3-butadiene; β-Chloroprene)	126-99-8	495
1-Chloropropane (Propylchloride)	540-54-5	386
3-Chloropropane-1,2-diol	96-24-2	268
2-Chloropropionic Acid Ethyl Ester	535-13-7	268
Chloropropylate	5836-10-2	117
2-Chloropyridine	109-09-1	572

Index

Analyte name	CAS No#	Page No#
6-Chloropyrimidine-2,4-diamine (2,6-Diamino-4-chloropyrimidine)	156-83-2	177
Chloroquine Phosphate	50-63-5	324
1-Chlorotetradecane	2425-54-9	268
Chlorothalonil	1897-45-6	5
Chlorothalonil-4-hydroxy (4-Hydroxychlorothalonil)	28343-61-5	5
Chlorothiazide	58-94-6	324
2-Chlorotoluene	95-49-8	427
3-Chlorotoluene	108-41-8	427
4-Chlorotoluene	106-43-4	427
Chlorotoluron	15545-48-9	46
Chlorotoluron D6 (N,N-dimethyl D6)	1219803-48-1	46, 530
2-Chloro-1,1,1-trifluoroethane	75-88-7	386
Chloroxuron	1982-47-4	46
Chloroxuron-N-monodesmethyl	28365-79-9	46
Chloroxylenol (4-Chloro-3,5-dimethylphenol)	88-04-0	428
Chloroxynil	1891-95-8	46
Chlorphenamine Maleate (Chlorpheniramine Maleate)	113-92-8	325
Chlorphenesin	104-29-0	325
Chlorphonium Chloride	115-78-6	46
Chlorphoxim	14816-20-7	117
Chlorpromazine Hydrochloride	69-09-0	325
Chlorpropamide	94-20-2	325
Chlorpropham	101-21-3	47
Chlorpyrifos	2921-88-2	117
Chlorpyrifos (diethyl-D10)	285138-81-0	117, 530
Chlorpyrifos-methyl	5598-13-0	117
Chlorpyrifos-methyl D6	n/a	117, 530
Chlorpyrifos-methyl-oxon	5598-52-7	117
Chlorpyrifos-oxon	5598-15-2	118
Chlorquinaldol	72-80-0	325
Chlorsulfuron	64902-72-3	47
Chlortetracycline Hydrochloride	64-72-2	325
Chlorthal-diacid	2136-79-0	47
Chlorthal-dimethyl	1861-32-1	47
Chlorthiamid	1918-13-4	47

Analyte name	CAS No#	Page No#
Chlorthion	500-28-7	118
Chlorthiophos	60238-56-4	118
Chlorthiophos-I-D10 (diethyl D10)	n/a	118, 531
Chlorthiophos I sulfone	25900-20-3	118
Chlozolate	84332-86-5	5
Chocolate Brown	n/a	405
Cholecalciferol (Vitamin D3)	67-97-0	461
4-Cholesten-3-one	601-57-0	572
Cholesterol	57-88-5	461
Chromafenozone	143807-66-3	118
Chrysanthemic Acid	10453-89-1	118
Chrysene	218-01-9	243
Chrysene D12	1719-03-5	243, 531
Ciclesonide	126544-47-6	325
Ciclosporin (Ciclosporin A; Cyclosporin A)	59865-13-3	325
Cimaterol	54239-37-1	325
Cimaterol D7 (isopropyl D7)	1228182-44-2	326, 531
Cimbuterol	54239-39-3	326
Cinidon (free acid)	175156-71-5	47
Cinidon-ethyl	142891-20-1	47
Cinmethylin	87818-31-3	47
Cinnamaldehyde	104-55-2	572
Cinnamic Acid Benzyl Ester	103-41-3	572
trans-Cinnamyl Alcohol (3-Phenyl-2-propen-1-ol)	4407-36-7	461
Cinosulfuron	94593-91-6	48
Cinoxacin	28657-80-9	326
Ciprofloxacin D8 Hydrochloride	1216659-54-9	326, 531
Ciprofloxacin Hydrochloride	93107-08-5	326
Ciprofloxacin Lactate	97867-33-9	326
Citicoline Sodium	33818-15-4	326
Citral	5392-40-5	572
Citric Acid Monohydrate	5949-29-1	572
(±)-β-Citronellol	106-22-9	572
Clarithromycin	81103-11-9	326
Clenbuterol Hydrochloride	21898-19-1	326

Index

Analyte name	CAS No#	Page No#
Clenbuterol-hydroxymethyl	38339-18-3	327
(±)-Clenbuterol	37148-27-9	327
(±)-Clenbuterol D9 (trimethyl D9)	129138-58-5	327, 531
Clenpenterol Hydrochloride	37158-47-7	327
Clenproperol	38339-11-6	327
Clenproperole D7	1173021-09-4	327, 531
Clethodim	99129-21-2	48
Climbazole	38083-17-9	327
Clindamycin Hydrochloride	21462-39-5	327
Clobetasol Propionate	25122-46-7	327
Clodinafop free acid	114420-56-3	48
Clodinafop-propargyl ester	105512-06-9	48
Cloethocarb	51487-69-5	118
Clofentezine	74115-24-5	118
Clofibrate (Clofibric acid ethyl ester)	637-07-0	328
Clofibric Acid (2-(4-Chlorophenoxy)-2-methylpropionic acid)	882-09-7	328
Clomazone (Command)	81777-89-1	48
Clomeprop	84496-56-0	48
Clomifene Citrate (Clomiphene Citrate)	50-41-9	328
Clophen A 30	55600-34-5	215
Clophen A 40	52306-32-8	215
Clophen A 50	8068-44-8	216
Clophen A 60	11096-99-4	216
Clopidol	2971-90-6	328
Cloprop	101-10-0	177
Clopyralid	1702-17-6	48
Cloquintocet free acid ((5-Chloro-8-quinolinyl)oxyacetic acid)	88349-88-6	177
Cloquintocet-1-methylhexyl ester (Cloquintocet-mexyl)	99607-70-2	177
Cloransulam-methyl	147150-35-4	48
Clorprenaline Hydrochloride	6933-90-0	328
Clorsulon	60200-06-8	328
Closantel	57808-65-8	328
Clostebol	1093-58-9	328

Analyte name	CAS No#	Page No#
Clostebol Acetate	855-19-6	328
Clothianidin	210880-92-5	119
Clothianidin D3 (N'-methyl D3)	1262776-24-8	119, 531
Clothianidin Urea	634192-72-6	119
Clotrimazole	23593-75-1	329
Cloxacillin Sodium Monohydrate	7081-44-9	329
Cloxyfonac	6386-63-6	177
Codlemone	33956-49-9	177
Colistin sulfate	1264-72-8	329
Copper(II) Hydroxide	20427-59-2	5
Copper Oxychloride	1332-40-7	5
Coronene	191-07-1	243
Corticosterone	50-22-6	329
Cortisone	53-06-5	329
Cortisone Acetate	50-04-4	329
Coumachlor	81-82-3	177
Coumafuryl	117-52-2	178
Coumaphos	56-72-4	119
Coumaphos alcohol metabolite	6174-86-3	119
Coumaphos D10 (di(ethyl-D5))	287397-86-8	119, 531
Coumaphos-oxon	321-54-0	119
Coumarin	91-64-5	329
Coumarin 5,6,7,8-D4	185056-83-1	329, 531
Coumatetralyl	5836-29-3	178
Creatine	57-00-1	573
Cresyl Diphenyl Phosphate	26444-49-5	484
Crimidine	535-89-7	178
Crotonaldehyde (Butenal)	123-73-9	268
Crotonaldehyde-2,4-dinitrophenylhydrazone (DNPH)	1527-96-4	269
Crotoxyphos	7700-17-6	119
Crufomate	299-86-5	119
Cuelure	3572-06-3	178
Cumyluron	99485-76-4	48
Curcumin	458-37-7	329
Cyanamide	420-04-2	49
Cyanazine	21725-46-2	49

Index

Analyte name	CAS No#	Page No#
Cyanazine D5 (N-ethyl D5)	1190003-29-2	49, 532
Cyanocobalamin (Vitamin B12)	68-19-9	461
Cyanofenphos	13067-93-1	120
Cyanophos	2636-26-2	120
Cyanuric Acid	108-80-5	573
Cyazofamid	120116-88-3	6
Cyclanilide	113136-77-9	178
Cycloate	1134-23-2	49
Cyclohexane	110-82-7	269
Cyclohexanol-2-(4-tert-butyl-phenoxy)	1942-71-8	178
Cyclohexanone	108-94-1	269
Cyclohexanone-2,4-dinitrophenylhydrazone	1589-62-4	269
Cycloheximide	66-81-9	178
Cyclohexylamine (Cyclohexanamine)	108-91-8	269
N-Cyclohexyldimethylamine	98-94-2	269
2-Cyclohexyl-4,6-dinitrophenol	131-89-5	428
Cyclohexylisocyanate	3173-53-3	573
Cyclopenta(c,d)pyrene	27208-37-3	243
Cyclopentadecanone	502-72-7	269
Cyclopentane	287-92-3	269
Cyclopentanol	96-41-3	269
Cyclopentanone	120-92-3	270
Cyclophosphamide Monohydrate	6055-19-2	330
Cycloprothrin	63935-38-6	120
Cyclosulfamuron	136849-15-5	49
Cycloxaprid	1203791-41-6	120
Cycloxydim	101205-02-1	49
Cycluron	2163-69-1	49
Cyenopyrafen	560121-52-0	120
Cyflufenamid	180409-60-3	6
Cyflumetofen	400882-07-7	120
Cyfluthrin	68359-37-5	120
β-Cyfluthrin	1820573-27-0	120

Analyte name	CAS No#	Page No#
trans-Cyfluthrin D6 (2,2-dimethyl D6)	n/a	120, 532
Cyhalofop	122008-78-0	49
Cyhalofop-butyl	122008-85-9	49
γ-Cyhalothrin	76703-62-3	121
λ-Cyhalothrin	91465-08-6	121
Cyhexatin	13121-70-5	121
Cymiazole	61676-87-7	330
Cymoxanil	57966-95-7	6
Cypermethrin	52315-07-8	121
α-Cypermethrin	67375-30-8	121
β-Cypermethrin	1224510-29-5	121
θ-Cypermethrin	71697-59-1	121
trans-Cypermethrin D6 (dimethyl D6)	82523-65-7	121, 532
ζ-Cypermethrin	1315501-18-8	122
Cyphenothrin	39515-40-7	122
Cyprazine	22936-86-3	50
Cyprazine-desisopropyl	35516-73-5	50
Cyprazine-desisopropyl-2-hydroxy	92510-62-8	50
Cyprazine-2-hydroxy	39095-16-4	50
Cyproconazole	94361-06-5	6
Cyprodinil	121552-61-2	6
Cyprofuram	69581-33-5	6
Cyrosulfamide	221667-31-8	178
Cyromazine	66215-27-8	122
Cyromazine D4	1219804-19-9	122, 532
DL-Cysteine	3374-22-9	573
Cythioate	115-93-5	122
2,4-D (phenyl-13C6)	150907-52-1	50, 532
2,4-D 2-Butoxyethyl Ester ((2,4-Dichlorophenoxy)acetic Acid 2-Butoxyethyl Ester)	1929-73-3	50
2,4-D Butyl Ester ((2,4-Dichlorophenoxy)acetic Acid Butyl Ester)	94-80-4	50

Index

Analyte name	CAS No#	Page No#
2,4-D Dimethylamine Salt	2008-39-1	50
2,4-D Ethyl Ester ((2,4-Dichlorophenoxy)acetic Acid Ethyl Ester)	533-23-3	50
2,4-D 2-Ethylhexyl Ester	1928-43-4	51
2,4-D Isobutyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isobutyl Ester)	1713-15-1	51
2,4-D Isooctyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isooctyl Ester)	25168-26-7	51
2,4-D Isopropyl Ester ((2,4-Dichlorophenoxy)acetic Acid Isopropyl Ester)	94-11-1	51
2,4-D Methyl Ester ((2,4-Dichlorophenoxy)acetic Acid Methyl Ester)	1928-38-7	51
2,4-D Sodium Salt Monohydrate	7084-86-8	51
2,3-D (2,3-Dichlorophenoxyacetic Acid)	2976-74-1	51
2,4-D ((2,4-Dichlorophenoxy) acetic Acid)	94-75-7	51
2,4-D D3 ((2,4-Dichloro-3,5,6-trideuteriophenoxy)acetic Acid)	202480-67-9	51, 532
2,6-D	575-90-6	52
3,4-D ((3,4-Dichlorophenoxy) acetic Acid)	588-22-7	52
Daimuron (Dymron)	42609-52-9	52
Dalapon (2,2-Dichloropropionic Acid)	75-99-0	53
Daminozide	1596-84-5	179
Daminozide D6 (dimethyl D6)	n/a	179, 533
Danofloxacin Mesilate	119478-55-6	330
Dantron (Danthron)	117-10-2	330
Dapsone	80-08-0	330
Dazomet	533-74-4	53
2,4-DB Dimethylamine Salt	2758-42-1	52
2,4-DB-2-ethylhexyl ester	1320-15-6	52
2,4-DB Methyl Ester (4-(2,4-Dichlorophenoxy) butanoic Acid Methyl Ester)	18625-12-2	52
2,4-DB (4-(2,4-Dichlorophenoxy) butanoic Acid)	94-82-6	52
4,4'-DDD (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethane)	72-54-8	122

Analyte name	CAS No#	Page No#
2,4'-DDE (1,1-Dichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl) ethene)	3424-82-6	122
2,4'-DDE D8 (1,1-Dichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl) ethene D8)	1402834-57-4	122, 532
4,4'-DDE (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethene)	72-55-9	122
4,4'-DDE D8 (1,1-Dichloro-2,2-bis(p-chlorophenyl)ethene D8)	93952-19-3	123, 532
4,4'-DDMU (1,1-Bis(p-chlorophenyl)-2-chloroethene)	1022-22-6	123
4,4'-DDNU (1,1-Bis(p-chlorophenyl)ethene)	2642-81-1	123
4,4'-DDT (ring-13C12)	104215-84-1	123, 532
DDT (technical)	8017-34-3	123
2,4'-DDT 13C12 (benzen 13C12)	1396995-26-8	123, 533
2,4'-DDT (1,1,1-Trichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl) ethane)	789-02-6	123
2,4'-DDT D8 (1,1,1-Trichloro-2-(p-chlorophenyl)-2-(o-chlorophenyl) ethane D8 (ring D8))	221899-88-3	123, 533
4,4'-DDT (1,1,1-Trichloro-2,2-bis(p-chlorophenyl)ethane)	50-29-3	123
4,4'-DDT D8 (1,1,1-Trichloro-2,2-bis(p-chlorophenyl)ethane D8 (ring D8))	93952-18-2	124, 533
Deanol (N,N-Dimethylethanolamine)	108-01-0	573
Decabromodiphenyl Ether (PBDE 209)	1163-19-5	573
Decacyclene	191-48-0	243
Decafluorobiphenyl	434-90-2	573
Decafluorotriphenylphosphine (DFTPP)	5074-71-5	495, 574
Decahydronaphthalene	91-17-8	270
cis-Decahydronaphthalene	493-01-6	270
trans-Decahydronaphthalene	493-02-7	270
Decanal-2,4-dinitrophenylhydrazone	1527-95-3	270
n-Decane	124-18-5	270
1,10-Decanedicarboxylic Acid	693-23-2	461
1,10-Decanediol	112-47-0	270
Decanoic Acid (Capric acid)	334-48-5	462

Index

Analyte name	CAS No#	Page No#
Decanoic Acid Ethyl Ester	110-38-3	462
Decanoic Acid Methyl Ester	110-42-9	462
1-Decanol	112-30-1	179
1-Decene	872-05-9	270
Decoquinat	18507-89-6	330
Decoquinat D5 (7-ethoxy D5)	1228182-55-5	330, 533
1-Decylamine	2016-57-1	270
Decylbenzene	104-72-3	271
DEET D7 (methyl D3 phenyl D4)	1219799-37-7	178, 533
Deltamethrin	52918-63-5	124
trans-Deltamethrin D6 (dimethyl D6)	n/a	124, 533
Demeclocycline Hydrochloride	64-73-3	330
Demeton (O+S)	8065-48-3	124
Demeton-O	298-03-3	124
Demeton-S	126-75-0	124
Demeton-S-methyl	919-86-8	124
Demeton-S-methyl D6 (dimethyl D6)	n/a	124, 533
Demeton-S-methyl sulfone	17040-19-6	124
Demeton-S-methyl sulfoxide (Oxydemeton-methyl)	301-12-2	125
Demeton-S-methyl-sulfoxide D6 (di-O-methyl D6)	n/a	125, 534
Denatonium Benzoate (Bitrex (TM))	3734-33-6	330
Desmedipham	13684-56-5	53
Desmetryn	1014-69-3	53
Desoxycortone (21-Hydroxyprogesterone)	64-85-7	331
Dexamethasone	50-02-2	331
Dexamethasone 21-Acetate	1177-87-3	331
Dexpanthenol (D-Panthenol)	81-13-0	462
3,5-Diacetyl-1,4-dihydro-2,6-lutidine	1079-95-4	574
Diafenthuron	80060-09-9	125
Dialifos	10311-84-9	125
Diallate	2303-16-4	53
2,4-Diaminoanisole Dihydrochloride	614-94-8	428

Analyte name	CAS No#	Page No#
4,4'-Diaminobiphenylmethane (Bis-(4-aminophenyl)methane)	101-77-9	405
trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic Acid Monohydrate	125572-95-4	574
4,4'-Diamino-3,3'-dimethyldiphenylmethane	838-88-0	405
4,4'-Diaminodiphenyl Sulfide (4-Aminophenylthioether)	139-65-1	405
1,4-Diamino-2-nitrobenzene	5307-14-2	428
2,6-Diamino-4-nitrotoluene	59229-75-3	428
1,5-Diaminopentane	462-94-2	271
1,2-Diaminopropane-N,N,N',N'-tetraacetic Acid	4408-81-5	574
2,6-Diaminopyridine	141-86-6	428
2,5-Diaminotoluene sulfate (2-Methyl-p-phenylenediamine sulfate)	615-50-9	428
2,3-Diaminotoluene	2687-25-4	428
2,4-Diaminotoluene	95-80-7	428
2,6-Diaminotoluene	823-40-5	429
3,4-Diaminotoluene	496-72-0	429
Diaryl-parl-phenylene-diamines	68953-84-4	574
Diatrizoate Sodium Dhydrate	1185748-84-8	331
Diatrizoic Acid	117-96-4	331
Diatrizoic Acid Meglumine Salt	131-49-7	331
Diaveridine	5355-16-8	331
Diazinon	333-41-5	125
Diazinon D10 (diethyl D10)	100155-47-3	125, 534
Dibenzo[a,c]anthracene	215-58-7	243
Dibenzo(a,e)fluoranthene	5385-75-1	243
Dibenzo(a,e)pyrene	192-65-4	244
Dibenzo(a,h)acridine	226-36-8	244
Dibenzo(a,h)anthracene	53-70-3	244
Dibenzo[a,h]anthracene D14	13250-98-1	244, 534
Dibenzo[a,h]pyrene	189-64-0	244
Dibenzo(a,i)pyrene	189-55-9	244
Dibenzo(a,j)acridine	224-42-0	244
Dibenzo(a,j)anthracene	224-41-9	244
Dibenzo(a,l)pyrene	191-30-0	244

Index

Analyte name	CAS No#	Page No#
Dibenzofuran	132-64-9	245
Dibenzofuran D8	93952-04-6	245, 534
Dibenzothiophene	132-65-0	245
Dibenzothiophene D8	33262-29-2	245, 534
Dibenzylamine	103-49-1	271
Dibromoacetic Acid	631-64-1	574
Dibromoacetonitrile	3252-43-5	386
1,2-Dibromobenzene	583-53-9	429
1,3-Dibromobenzene	108-36-1	429
1,4-Dibromobenzene	106-37-6	429
4,4'-Dibromobenzophenone	3988-03-2	179
Dibromochloroacetic acid-methyl Ester	20428-75-5	574
Dibromochloromethane	124-48-1	386
1,2-Dibromo-3-chloropropane	96-12-8	386
1,2-Dibromo-1,1-dichloroethane	75-81-0	386
Dibromodichloromethane	594-18-3	387
1,2-Dibromoethane	106-93-4	387
1,2-Dibromoethene	540-49-8	387
4',5'-Dibromofluorescein	596-03-2	405
Dibromomethane	74-95-3	387
1,4-Dibromonaphthalene	83-53-4	245
4,4'-Dibromooctafluorobiphenyl	10386-84-2	574
2,4-Dibromophenol	615-58-7	429
2,6-Dibromophenol	608-33-3	429
1,2-Dibromopropane	78-75-1	429
3-(2,2-Dibromovinyl)-2,2-dimethyl-(1-cyclopropane) carboxylic	63597-73-9	125
Dibutyl Phthalate	84-74-2	484
Dibutyl Sulfide	544-40-1	271
Dibutylamine	111-92-2	271
Dibutyldiphenyltin	6452-61-5	574
Di-n-butylether	142-96-1	271
2,5-Di-tert-butylhydroquinone	88-58-4	429
2,6-Di-tert-butyl-4-hydroxymethylphenol	88-26-6	575
2-(3,5-Di-tert-butyl-2-hydroxyphenyl)-2H-benzotriazole	3846-71-7	575

Analyte name	CAS No#	Page No#
2,6-Di-tert-butyl-4-methoxyphenol	489-01-0	430
4,6-Di-tert-butyl-2-methylphenol	616-55-7	430
2,4-Di-tert-butylphenol	96-76-4	430
2,6-Di-tert-butylphenol	128-39-2	430
3,5-Di-tert-butylphenol	1138-52-9	430
Di-n-butylphosphate	107-66-4	484
Dibutyltin Dichloride	683-18-1	575
Dibutyltin Oxide	818-08-6	575
Dicamba	1918-00-9	53
Dicamba 13C6 (ring 13C6)	1173023-06-7	53, 534
Dicamba D3 (methoxy D3)	349553-95-3	53, 534
Dicamba D4 (phenyl D1 methoxy D3)	n/a	53, 534
Dicamba-desmethyl	3401-80-7	54
Dicamba-methyl ester	6597-78-0	125
Dicamba-5-hydroxy (5-Hydroxydicamba)	7600-50-2	54
Dicapthon	2463-84-5	126
Dichlobenil	1194-65-6	54, 126
Dichlofenthion	97-17-6	126
Dichlofluanid	1085-98-9	6
Dichlone	117-80-6	6
Dichlorflurenol-methyl	21634-96-8	179
Dichlormid	37764-25-3	179
Dichloroacetic Acid	79-43-6	575
Dichloroacetic Acid Methyl Ester	116-54-1	575
Dichloroacetonitrile	3018-12-0	387
4-(Dichloroacetyl)-1-oxa-4-azaspiro[4,5]decane	71526-07-3	575
2,3-Dichloroaniline	608-27-5	430
2,4-Dichloroaniline	554-00-7	430
2,5-Dichloroaniline	95-82-9	179
2,6-Dichloroaniline	608-31-1	430
3,4-Dichloroaniline	95-76-1	179
3,5-Dichloroaniline	626-43-7	180
2,3-Dichloroanisole	1984-59-4	430
2,4-Dichloroanisole	553-82-2	431

Index

Analyte name	CAS No#	Page No#
2,6-Dichloroanisole	1984-65-2	431
3,4-Dichloroanisole	36404-30-5	431
3,5-Dichloroanisole	33719-74-3	431
2,6-Dichlorobenzamide	2008-58-4	180
1,2-Dichlorobenzene	95-50-1	387, 431
1,2-Dichlorobenzene D4	2199-69-1	431, 534
1,3-Dichlorobenzene	541-73-1	387, 431
1,4-Dichlorobenzene	106-46-7	387, 431
1,4-Dichlorobenzene D4	3855-82-1	387, 535, 431, 535
3,3'-Dichlorobenzidine Dihydrochloride	612-83-9	405
3,3'-Dichlorobenzidine	91-94-1	406
3,3'-Dichlorobenzidine D6 (ring D6)	93951-91-8	406, 535
2,4-Dichlorobenzoic Acid	50-84-0	575
2,6-Dichlorobenzoic Acid	50-30-6	180
3,5-Dichlorobenzoic Acid	51-36-5	180
2,4'-Dichlorobenzophenone	85-29-0	180
4,4'-Dichlorobenzophenone	90-98-2	180
4,4'-Dichlorobenzophenone D8	1219806-01-5	180, 535
1,4-Dichlorobutane	110-56-5	388
cis-1,4-Dichloro-2-butene	1476-11-5	388
trans-1,4-Dichloro-2-butene	110-57-6	388
3,5-Dichlorocatechol	13673-92-2	432
4,5-Dichlorocatechol	3428-24-8	432
1,1-Dichloroethane	75-34-3	388
1,2-Dichloroethane	107-06-2	388
1,2-Dichloroethane D4	17060-07-0	388, 535
2,2-Dichloroethanol	598-38-9	388
1,1-Dichloroethene	75-35-4	388
cis-1,2-Dichloroethene	156-59-2	388
trans-1,2-Dichloroethene	156-60-5	389
1,1-Dichloro-1-fluoroethane	1717-00-6	389
Dichlorofluoromethane	75-43-4	389
2,4-Dichloro-6-methylphenol	1570-65-6	432
2,6-Dichloro-4-methylphenol	2432-12-4	432
2-(4,6-Dichloro-2-methylphenoxy)propionic Acid	20021-12-9	575

Analyte name	CAS No#	Page No#
1,4-Dichloronaphthalene	1825-31-6	245
1,5-Dichloronaphthalene	1825-30-5	245
1,8-Dichloronaphthalene	2050-74-0	245
2,3-Dichloronitrobenzene	3209-22-1	432
2,4-Dichloronitrobenzene	611-06-3	432
2,5-Dichloronitrobenzene	89-61-2	432
3,4-Dichloronitrobenzene	99-54-7	432
3,5-Dichloronitrobenzene	618-62-2	432
Dichlorophen	97-23-4	6
2,4-Dichlorophenol 13C6	1202864-83-2	180, 535
2,3-Dichlorophenol	576-24-9	433
2,4-Dichlorophenol	120-83-2	180
2,4-Dichlorophenol D3	93951-74-7	181, 535
2,5-Dichlorophenol	583-78-8	433
2,6-Dichlorophenol	87-65-0	433
3,4-Dichlorophenol	95-77-2	433
3,5-Dichlorophenol	591-35-5	433
2,4-Dichlorophenyl Acetic Acid (2,4-DCAA)	19719-28-9	576
2,4-Dichlorophenylacetic Acid Methyl Ester (DCAA Methyl Ester)	55954-23-9	576
2-[(2,6-Dichlorophenyl)amino] benzaldehyde	22121-58-0	331
(1RS)-1-(2,4-Dichlorophenyl)-2-(1H-imidazol-1-yl)ethanol	24155-42-8	331
1-(3,4-Dichlorophenyl)-3-methoxyurea	17356-61-5	54
1-(3,4-Dichlorophenyl)-3-methylurea	3567-62-2	54
1-(3,5-Dichlorophenyl)-3-methylurea	38655-13-9	54
1-(3,4-Dichlorophenyl)urea	2327-02-8	54
1,1-Dichloropropane	78-99-9	389
1,2-Dichloropropane	78-87-5	389
1,3-Dichloropropane	142-28-9	389
2,2-Dichloropropane	594-20-7	389
1,3-Dichloropropan-2-ol	96-23-1	389
2,3-Dichloro-1-propanol	616-23-9	389
1,1-Dichloropropan-2-one	513-88-2	390
1,1-Dichloro-1-propene	563-58-6	390

Index

Analyte name	CAS No#	Page No#
1,3-Dichloropropene	542-75-6	390
2,3-Dichloro-1-propene	78-88-6	390
cis-1,3-Dichloropropene	10061-01-5	390
trans-1,3-Dichloropropene	10061-02-6	390
2,3-Dichloropropionanilide	17839-22-4	576
3,3''-Dichloro-o-terphenyl	2007947-35-3	216
3,3''-Dichloro-p-terphenyl	159392-36-6	216
4,4''-Dichloro-o-terphenyl	21711-56-8	216
4,4''-Dichloro-p-terphenyl	21711-52-4	216
1,2-Dichlorotetrafluoroethane (CFC-114)	76-14-2	390
2,3-Dichlorotoluene	32768-54-0	433
2,4-Dichlorotoluene	95-73-8	433
2,5-Dichlorotoluene	19398-61-9	433
2,6-Dichlorotoluene	118-69-4	433
3,4-Dichlorotoluene	95-75-0	434
α,α-Dichlorotoluene	98-87-3	434
1,2-Dichlorotrifluoroethane	354-23-4	390
(E)-3-(2,2-Dichlorovinyl)-2,2-di(methyl D3)-(1-cyclopropane) carboxylic acid D6	n/a	126, 535
Dichloroxylenol (2,4-Dichloro-3,5-dimethylphenol)	133-53-9	434
Dichlorprop (2,4-DP)	120-36-5	54
Dichlorprop-butoxyethyl ester	53404-31-2	55
Dichlorprop D6 (ring D3, 3,3,3-D3)	n/a	54, 535
Dichlorprop 2-Ethylhexyl Ester	79270-78-3	55
Dichlorprop-methyl Ester	23844-57-7	55
Dichlorprop-P (D-(+)-Dichlorprop)	15165-67-0	55
2,6-Dichlorprop	25140-90-3	55
Dichlorvos	62-73-7	126
Dichlorvos D6 (dimethyl D6)	203645-53-8	126, 535
Diclazuril	101831-37-2	332
Diclazuril-methyl	103337-71-9	332
Diclobutrazol	75736-33-3	7
Diclocymet	139920-32-4	126

Analyte name	CAS No#	Page No#
Diclofenac	15307-86-5	332
Diclofenac D4 (phenyl D4)	153466-65-0	332, 536
Diclofenac Methyl Ester	15307-78-5	332
Diclofenac Sodium	15307-79-6	332
Diclofop (free acid)	40843-25-2	126
Diclofop methyl	51338-27-3	55
Diclofop-methyl-5-hydroxy (2-[4-(2,4-Dichloro-5-hydroxyphenoxy)phenoxy] propanoic Acid Methyl Ester)	124992-48-9	55
Diclomezine	62865-36-5	7
Dicloran (2,6-Dichloro-4-nitroaniline)	99-30-9	7
Diclosulam	145701-21-9	55
Dicloxacillin Sodium Monohydrate	13412-64-1	332
Dicofol	115-32-2	126
Dicofol D8 (ring D8)	n/a	127, 536
2,4'-Dicofol	10606-46-9	127
Dicrotophos	141-66-2	127
Dicyclanil	112636-83-6	127
Dicyclohexylamine (N-Cyclohexylcyclohexanamine)	101-83-7	271
Dicyclopentadiene	77-73-6	181
7,8-Didehydrocholesterol (Provitamin D3)	434-16-2	462
Didodecyl 3,3'-thiodipropionate (3,3'-Thiodipropionic acid didodecyl ester)	123-28-4	462
Dieldrin	60-57-1	127
Dienestrol	84-17-3	332
Dienestrol Diacetate	84-19-5	332
Dienochlor	2227-17-0	127
Diesel Fuel	68334-30-5	271
Bio Diesel Fuel	870530-78-2	272
Diethatyl Ethyl Ester	38727-55-8	55
Diethofencarb	87130-20-9	7
1,1-Diethoxyethane (Acetaldehyde diethylacetal)	105-57-7	272
Diethyl Fumarate	623-91-6	333
Diethyl Phosphorochloridothionate	2524-04-1	181

Index

Analyte name	CAS No#	Page No#
Diethyl Phthalate	84-66-2	484
Diethylamine	109-89-7	272
4-(Diethylamino)azobenzene	2481-94-9	406
2-Diethylaminoethyl Hexanoate	10369-83-2	181
Diethylaminohydroxybenzoyl hexyl benzoate	302776-68-7	576
2-Diethylamino-6-methyl-4-pyrimidinol	42487-72-9	181
2,6-Diethylaniline	579-66-8	181
N,N-Diethylaniline	91-66-7	272
1,2-Diethylbenzene	135-01-3	272
1,3-Diethylbenzene	141-93-5	272
1,4-Diethylbenzene	105-05-5	272
Diethylcarbamazine Citrate	1642-54-2	333
N,N-Diethyldithiocarbamate Sodium Salt Trihydrate	20624-25-3	576
Diethylene Glycol Monobutyl Ether	112-34-5	272
Diethylene glycol monoethyl ether	111-90-0	273
Diethylene Glycol Monohexyl Ether	112-59-4	272
Diethylene Glycol Monomethyl Ether	111-77-3	273
Diethyleneglycol (2,2'-Oxydiethanol)	111-46-6	273
Diethyleneglycoldinitrate (DEGN)	693-21-0	576
Diethylenetriaminepentaacetic Acid (DTPA)	67-43-6	576
Diethylether (Ether)	60-29-7	273
Diethylhexylphthalate (Bis-(2-ethylhexyl) Phthalate)	117-81-7	484
4-N,N-Diethyl-2-methyl-p-phenyldiamine monohydrochloride	2051-79-8	434
N,N-Diethyl-p-phenyldiamine sulfate salt	6283-63-2	576
Diethylstilbestrol	56-53-1	333
Diethyltoluamide (DEET)	134-62-3	181
Difenacoum	56073-07-5	181
Difenoconazole	119446-68-3	7
Difenoconazole D6 (1,1,2,3,3,3-propyl-D6)	n/a	7, 536

Analyte name	CAS No#	Page No#
Difenoconazole-alcohol	117018-19-6	7
Difenoxuron	14214-32-5	56
Difenzoquat Methyl Sulfate (1,2-Dimethyl-3,5-diphenylpyrazolium methyl sulfate)	43222-48-6	56
Difethialone	104653-34-1	181
Difloxacin Hydrochloride	91296-86-5	333
Diflubenzuron	35367-38-5	127
Diflufenican	83164-33-4	56
Diflufenican D3 (3-trifluoromethylphenoxy-2,4,6 D3)	1185009-29-3	56, 536
Diflufenzopyr	109293-97-2	56
Diflufenzopyr sodium salt	109293-98-3	56
Diflumetorim	130339-07-0	7
1,4-Difluorobenzene	540-36-3	434
2,6-Difluorobenzoic Acid	385-00-2	182
2,2'-Difluorobiphenyl (PFB 4)	388-82-9	576
4,4'-Difluorobiphenyl	398-23-2	577
Digitoxin	71-63-6	333
Diglyme (Bis(2-methoxyethyl) ether; Diethylene glycol dimethyl ether)	111-96-6	273
Digoxin	20830-75-5	333
Di-n-heptyltin-dichloride	74340-12-8	577
9,10-Dihydroanthracene	613-31-0	245
Dihydrosafrol	94-58-6	577
Dihydrostreptomycin Sulfate	5490-27-7	333
4,4'-Dihydroxybiphenyl	92-88-6	577
1,4-Dihydroxynaphthalene	571-60-8	246
1,3-Diisopropylbenzene	99-62-7	273
1,4-Diisopropylbenzene	100-18-5	273
Diisopropylether	108-20-3	273
Diisopropylnaphthalene (mixture of isomers)	38640-62-9	182
2,6-Diisopropylnaphthalene	24157-81-1	182
Dikegulac Monohydrate	68539-16-2	182
Dikegulac Na	52508-35-7	182
Diltiazem Hydrochloride	33286-22-5	333
Dimefluthrin	271241-14-6	127

Index

Analyte name	CAS No#	Page No#
Dimefox	115-26-4	127
Dimefuron	34205-21-5	56
Dimepiperate	61432-55-1	56
Dimethachlon	24096-53-5	7
Dimethachlor	50563-36-5	56
Dimethachlor-deschloro	1231710-71-6	57
Dimethachlor ethane sulfonic acid (ESA) sodium salt	1231710-75-0	57
Dimethachlor oxalamic acid (OA)	1086384-49-7	57
Dimethametryn	22936-75-0	57
Dimethenamid	87674-68-8	57
Dimethenamid-P	163515-14-8	57
Dimethipin	55290-64-7	57
Dimethirimol	5221-53-4	8
Dimethoate	60-51-5	128
Dimethoate D6 (O,O dimethyl D6)	1219794-81-6	128, 536
Dimethomorph	110488-70-5	8
Z-Dimethomorph	113210-98-3	8
3,3'-Dimethoxybenzidine	119-90-4	406
1,2-Dimethoxyethane	110-71-4	495
Dimethyl Disulfide	624-92-0	577
Dimethyl Fumarate (Fumaric acid dimethyl ester)	624-49-7	462
Dimethyl sulfoxide	67-68-5	273
Dimethyl Terephthalate (Terephthalic acid-bis-methyl ester)	120-61-6	577
Dimethyl trisulfide	3658-80-8	577
1,3-Dimethyl Urea	96-31-1	577
Dimethylacetamide (Acetic acid dimethylamide)	127-19-5	462
2,4-Dimethylaniline	95-68-1	182
2,6-Dimethylaniline	87-62-7	434
N,N-Dimethylaniline	121-69-7	434
2,5-Dimethylbenzaldehyde-2,4-dinitrophenylhydrazone	152477-96-8	434
1,2-Dimethylbenzene (o-Xylene)	95-47-6	274
1,3-Dimethylbenzene (m-Xylene)	108-38-3	274
1,4-Dimethylbenzene (p-Xylene)	106-42-3	274

Analyte name	CAS No#	Page No#
3,3'-Dimethylbenzidine (o-Tolidine)	119-93-7	406
7,12-Dimethylbenzo[a]anthracene	57-97-6	246
2,2-Dimethylbutane	75-83-2	274
2,3-Dimethylbutane	79-29-8	274
1,3-Dimethylcyclohexane	591-21-9	274
1,4-Dimethylcyclohexane	589-90-2	274
Dimethyldioctylammonium Chloride	5538-94-3	182
N,N-Dimethyldithiocarbamate Nickel Salt	15521-65-0	577
N,N-Dimethyldithiocarbamate Sodium Salt Hydrate	207233-95-2	578
N,N-Dimethylformamide	68-12-2	274
2,6-Dimethyl-4-heptanone	108-83-8	274
Dimethylnaphthalene (technical mixture)	28804-88-8	246
1,2-Dimethylnaphthalene	573-98-8	246
1,3-Dimethylnaphthalene	575-41-7	246
1,4-Dimethylnaphthalene	571-58-4	246
1,5-Dimethylnaphthalene	571-61-9	246
1,6-Dimethylnaphthalene	575-43-9	246
1,7-Dimethylnaphthalene	575-37-1	246
1,8-Dimethylnaphthalene	569-41-5	247
2,3-Dimethylnaphthalene	581-40-8	247
2,6-Dimethylnaphthalene	581-42-0	247
2,7-Dimethylnaphthalene	582-16-1	247
1,2-Dimethyl-3-nitrobenzene	83-41-0	434
1,2-Dimethyl-4-nitrobenzene	99-51-4	435
1,3-Dimethyl-2-nitrobenzene (2-Nitro-m-xylene)	81-20-9	435
2,2-Dimethylpentane	590-35-2	275
2,3-Dimethylpentane	565-59-3	275
2,4-Dimethylpentane	108-08-7	275
3,6-Dimethylphenanthrene	1576-67-6	247
2,3-Dimethylphenol	526-75-0	435
2,4-Dimethylphenol	105-67-9	435
2,5-Dimethylphenol	95-87-4	435
2,6-Dimethylphenol	576-26-1	435
3,4-Dimethylphenol	95-65-8	435

Index

Analyte name	CAS No#	Page No#
3,5-Dimethylphenol	108-68-9	435
N,N-Dimethyl-p-phenylendiamine	99-98-9	435
N-(2,4-Dimethylphenyl) formamide (Form-2',4'-xylidide)	60397-77-5	128
(3,5-Dimethylphenyl) N-Methylcarbamate (XMC)	2655-14-3	128
N-2,4-Dimethylphenyl-N'-methylformamidine (N-Methyl-N'-(2,4-xylyl)formamidine)	33089-74-6	128
O,O-Dimethylphosphorodithioic acid ammonium salt	1066-97-3	578
O,O-Dimethylphosphorothioic acid sodium salt	23754-87-2	578
O,O-Dimethylphosphorothioic acid sodium salt D6	n/a	536, 578
2,2-Dimethylpropane (Neopentane)	463-82-1	275
2,2-Dimethylpropionic Acid	75-98-9	578
N,N-Dimethylpyridin-4-amine (4-Dimethylaminopyridine)	1122-58-3	578
N,N-Dimethylsulfamide	3984-14-3	578
Dimethyltin Dichloride	753-73-1	578
2,4-Dimethyl-3,5,6-trideuteriophenol	93951-75-8	436, 536
Dimethylvinphos (Z type)	67628-93-7	128
Dimeticone (Dimethicone)	9006-65-9	334
Dimetilan	644-64-4	128
Dimetridazole	551-92-8	334
Dimetridazole D3 (N-methyl D3)	64678-69-9	334, 536
Dimetridazole-2-hydroxy D3	1015855-78-3	334, 536
Dimetridazole-2-hydroxy ((1-Methyl-5-nitroimidazol-2-yl) methanol)	936-05-0	334
Diminazene Aceturate	908-54-3	334
Dimoxystrobin	149961-52-4	8
Diniconazole	83657-24-3	8
Dinitolmide (3,5-Dinitro-o-toluamide)	148-01-6	436
Dinitramine	29091-05-2	57
2,4-Dinitroaniline	97-02-9	436
2,6-Dinitroaniline	606-22-4	436
3,5-Dinitroaniline	618-87-1	436

Analyte name	CAS No#	Page No#
2,4-Dinitroanisole	119-27-7	436
1,2-Dinitrobenzene	528-29-0	436
1,3-Dinitrobenzene	99-65-0	436
1,4-Dinitrobenzene	100-25-4	436
2,4-Dinitrobenzoic Acid	610-30-0	437
3,5-Dinitrobenzoic Acid	99-34-3	437
2,4-Dinitrodiphenylamine	961-68-2	437
2,6-Dinitro-4-methylphenol	609-93-8	437
1,3-Dinitronaphthalene	606-37-1	247
1,5-Dinitronaphthalene	605-71-0	247
1,8-Dinitronaphthalene	602-38-0	247
2,6-Dinitrophenol	573-56-8	437
1,3-Dinitropyrene	75321-20-9	247
1,6-Dinitropyrene	42397-64-8	248
1,8-Dinitropyrene	42397-65-9	248
2,3-Dinitrotoluene	602-01-7	437
2,4-Dinitrotoluene	121-14-2	437
2,6-Dinitrotoluene	606-20-2	437
3,4-Dinitrotoluene	610-39-9	437
3,5-Dinitrotoluene	618-85-9	438
Dinobuton	973-21-7	128
Dinocap	39300-45-3	8
Dinoseb	88-85-7	57
Dinoseb acetate	2813-95-8	58
Dinoseb Methyl Ether	6099-79-2	58
Dinotefuran	165252-70-0	129
Dinoterb	1420-07-1	58
Dinoterb Acetate	3204-27-1	58
Di-n-Octyltin Dichloride	3542-36-7	578
Diofenolan	63837-33-2	129
Dioxabenzofos (Salithion)	3811-49-2	129
Dioxacarb	6988-21-2	129
1,4-Dioxane-d8 (Octadeuterodioxane)	17647-74-4	275, 537
1,4-Dioxane	123-91-1	275
Dioxathion	78-34-2	129
Diphenadione (Diphacinone)	82-66-6	182
Diphenamid	957-51-7	58

Index

Analyte name	CAS No#	Page No#
Diphenhydramine Hydrochloride	147-24-0	334
Diphenyl Sulfide	139-66-2	275
Diphenyl Sulfone	127-63-9	182
Diphenylamine	122-39-4	8
Diphenylmercury	587-85-9	579
Diphenylmethan-4,4'-diisocyanate	101-68-8	579
Diphenylmethanol	91-01-0	579
Diphenyltin Dichloride	1135-99-5	579
N,N'-Diphenylurea (Carbanalide)	102-07-8	579
Dipropetryn	4147-51-7	58
Dipropyl Pyridine-2,5-dicarboxylate (MGK 326)	136-45-8	183
Dipropylamine	142-84-7	275
Dipropylene Glycol Monomethyl Ether	34590-94-8	275
Dipropyleneglycol	25265-71-8	276
Di-n-Propyltin Dichloride	867-36-7	579
Diquat Dibromide	85-00-7	58
Diquat dibromide D4	n/a	58, 537
Direct Black 38	1937-37-7	406
Direct Blue 86 (technical product)	1330-38-7	406
Direct Brown 95	16071-86-6	406
Direct Red 28	573-58-0	406
Disperse Black 2	6232-57-1	407
Disperse Blue 1	2475-45-8	407
Disperse Blue 102	69766-79-6	407
Disperse Blue 106	68516-81-4	407
Disperse Blue 124	15141-18-1	407
Disperse Blue 26	3860-63-7	407
Disperse Blue 291	56548-64-2	407
Disperse Blue 3	2475-46-9	407
Disperse Blue 7	3179-90-6	407
Disperse Brown 1	23355-64-8	408
Disperse Orange 1	2581-69-3	408
Disperse Orange 11	82-28-0	408
Disperse Orange 13	6253-10-7	408

Analyte name	CAS No#	Page No#
Disperse Orange 149	85136-74-9	408
Disperse Orange 3	730-40-5	408
Disperse Orange 37	13301-61-6	408
Disperse Orange 61	55281-26-0	408
Disperse Red 1	2872-52-8	408
Disperse Red 11	2872-48-2	409
Disperse Red 13	3180-81-2	409
Disperse Red 15	116-85-8	409
Disperse Red 151	70210-08-1	409
Disperse Red 151 Metabolite 1	58104-55-5	409
Disperse Red 17	3179-89-3	409
Disperse Red 7	4540-00-5	409
Disperse Violet 1	128-95-0	409
Disperse Violet 93	52697-38-8	409
Disperse Violet 93 (Cl derivative)	66557-45-7	410
Disperse Yellow 1	119-15-3	410
Disperse Yellow 23	6250-23-3	410
Disperse Yellow 3	2832-40-8	410
Disperse Yellow 3 D3 (acetyl D3)	947601-96-9	410, 537
Disperse Yellow 39	12236-29-2	410
Disperse Yellow 39 surrogate	56208-37-8	410
Disperse Yellow 49	6858-49-7	410
Disperse Yellow 54	7576-65-0	410
Disperse Yellow 56-methyl	73287-67-9	411
Disperse Yellow 7	6300-37-4	411
Disperse Yellow 9	6373-73-5	411
Disul Sodium	136-78-7	58
Disulfiram	97-77-8	334
Disulfoton	298-04-4	129
Disulfoton D10 (Di-ethyl D10)	n/a	129, 537
Disulfoton-oxon-sulfon	2496-91-5	129
Disulfoton-sulfone	2497-06-5	130
Disulfoton-sulfoxide	2497-07-6	130
Ditalimfos	5131-24-8	8
Dithianon	3347-22-6	8
Dithiopyr	97886-45-8	59
Dithranol	1143-38-0	334

Index

Analyte name	CAS No#	Page No#
Diuron	330-54-1	59
Diuron D6	1007536-67-5	59, 537
N,N'-Divinylethyleneurea	13811-50-2	579
DMSA	4710-17-2	573
DMST	66840-71-9	179
DNOC (2-Methyl-4,6-dinitrophenol)	534-52-1	52
DNOC D5 (ring D2, methyl D3)	n/a	52, 533
n-Docosane	629-97-0	276
Docosanoic Acid (Behenic acid)	112-85-6	462
Docosanoic Acid Ethyl Ester	5908-87-2	463
Docosanoic Acid Methyl Ester	929-77-1	463
1-Docosanol (Behenyl alcohol)	661-19-8	276
1-Docosene	1599-67-3	276
Dodecanal	112-54-9	276
n-Dodecane	112-40-3	276
Dodecanoic Acid (Lauric acid)	143-07-7	463
Dodecanoic Acid Ethyl Ester	106-33-2	463
Dodecanoic Acid Methyl Ester	111-82-0	463
1-Dodecanol	112-53-8	276
2-Dodecanol	10203-28-8	276
1-Dodecene	112-41-4	276
4-Dodecylphenol, mixture of isomers	27193-86-8	438
Dodemorph	1593-77-7	9
trans-Dodemorph	91269-48-6	9
Dodine (Dodecylguanidium Acetate)	2439-10-3	9
Dopamine Hydrochloride	62-31-7	335
Doramectin	117704-25-3	335
n-Dotriacontane	544-85-4	277
Doxepin Hydrochloride	1229-29-4	335
Doxycycline Hyclate	24390-14-5	335
DPMI (Cashmeran)	33704-61-9	573
Drazoxolon	5707-69-7	9
Econazole Nitrate	24169-02-6	335
Edetic Acid (EDTA)	60-00-4	579
Edifenphos	17109-49-8	9

Analyte name	CAS No#	Page No#
EDTA-tetrabutyl ester	14531-12-5	579
11,14-Eicosadienoic Acid Methyl Ester	2463-02-7	463
n-Eicosane	112-95-8	277
Eicosanoic Acid (Arachidic acid)	506-30-9	463
Eicosanoic Acid Methyl Ester	1120-28-1	463
1-Eicosanol (Arachidyl alcohol)	629-96-9	277
1-Eicosene	3452-07-1	277
cis-8-Eicosenoic Acid	76261-96-6	463
Elaidic Acid	112-79-8	464
Elaidic Acid Methyl Ester	1937-62-8	464
Emamectin Benzoate	155569-91-8	130
Embutramide	15687-14-6	335
Empenthrin	54406-48-3	130
Enalapril Maleate	76095-16-4	335
Endosulfan (α + β)	115-29-7	130
Endosulfan-alcohol	2157-19-9	130
Endosulfan-ether	3369-52-6	130
Endosulfan-lacton	3868-61-9	130
Endosulfan-sulfate	1031-07-8	131
α-Endosulfan	959-98-8	131
α-Endosulfan D4	203645-57-2	131, 537
β-Endosulfan	33213-65-9	131
β-Endosulfan D4	203716-99-8	131, 537
Endothal	145-73-3	59
Endrin	72-20-8	131
Endrin-aldehyde	7421-93-4	131
Endrin-ketone	53494-70-5	131
Enilconazole (Imazalil)	35554-44-0	9
Enilconazole-D5 (2-propenyl-D5) (Imazalil-D5)	1398065-91-2	9, 537
Enoxacin	74011-58-8	335
Enoxacin Sesquihydrate	84294-96-2	335
Enoxastrobin	238410-11-2	9
Enrofloxacin	93106-60-6	336
Enrofloxacin D5 Hydrochloride (ethyl d5)	n/a	336, 537
4-Epianhydrotetracycline hydrochloride	4465-65-0	336

Index

Analyte name	CAS No#	Page No#
Epichlorohydrin	106-89-8	9
Epigallocatechin-3-gallate	989-51-5	336
EPN	2104-64-5	130
Epoxiconazole	133855-98-8	10
Eprinomectin	123997-26-2	336
EPTC	759-94-4	59
Equilenin	517-09-9	336
Equilin	474-86-2	336
Erbon	136-25-4	59
Ergocalciferol (Vitamin D2)	50-14-6	464
Ergosterol (Provitamin D2)	57-87-4	464
Erucic acid	112-86-7	464
Erythromycin	114-07-8	336
Erythromycin Ethylsuccinate	1264-62-6	336
Erythromycin Stearate	643-22-1	337
Erythrosine B Disodium Salt (E127)	16423-68-0	411
Esfenvalerate	66230-04-4	132
Esfenvalerate free acid metabolite	55332-38-2	132
Esomeprazole Magnesium Trihydrate	217087-09-7	337
Esprocarb	85785-20-2	59
Estradiol (17β-Estradiol)	50-28-2	337
Estradiol Benzoate (17β-Estradiol 3-benzoate)	50-50-0	337
Estradiol Diacetate (17β-Estradiol 3,17-diacetate)	3434-88-6	337
Estradiol Valerate (17β-Estradiol 17-valerate)	979-32-8	337
17-epi-Estradiol (17α-Estradiol)	57-91-0	337
Estriol	50-27-1	337
Estrone	53-16-7	337
Etaconazole	60207-93-4	10
Ethaboxam	162650-77-3	10
Ethalfuralin	55283-68-6	59
Ethametsulfuron-methyl	97780-06-8	59
Ethanol	64-17-5	277
Ethephon	16672-87-0	183
Ethephon D4 (2-Chloroethyl-1,1,2,2-D4)	1020719-29-2	183, 538

Analyte name	CAS No#	Page No#
Ethephon-hydroxy (2-Hydroxyethanephosphonic Acid)	22987-21-9	183
Ethidimuron	30043-49-3	60
Ethinylestradiol (17α-Ethinylestradiol)	57-63-6	338
Ethiofencarb	29973-13-5	132
Ethiofencarb-phenol-sulfone	53380-27-1	132
Ethiofencarb-phenol-sulfoxide	53380-26-0	132
Ethiofencarb Sulfone	53380-23-7	132
Ethiofencarb-sulfoxide	53380-22-6	132
Ethiolate	2941-55-1	60
Ethion	563-12-2	132
Ethion D20 (tetraethyl D20)	n/a	132, 538
Ethiozin	64529-56-2	60
Ethiprole	181587-01-9	133
Ethirimol	23947-60-6	10
Ethisterone	434-03-7	338
Ethofumesate	26225-79-6	60
Ethofumesate-2-hydroxy	26322-82-7	60
Ethofumesate-2-keto	26244-33-7	60
Ethopabate	59-06-3	338
Ethoprophos (Ethoprop)	13194-48-4	133
Ethoprophos D5 (ethyl D5)	n/a	133, 538
Ethoxyacetic Acid Ethyl Ester	817-95-8	580
2-Ethoxyethanol (Ethyleneglycol monoethyl ether)	110-80-5	277
Ethoxyfen-ethyl	131086-42-5	60
1-Ethoxy-2-propanol	1569-02-4	277
Ethoxyquin	91-53-2	10
Ethoxyquin hydrochloride	3659-01-6	10
Ethoxysulfuron	126801-58-9	60
Ethychlozate	27512-72-7	183
Ethyl acetate	141-78-6	277
Ethyl (1-Naphthyl)acetate (1-Naphthyl acetic acid ethyl ester)	2122-70-5	183
Ethyl Parahydroxybenzoate (4-Hydroxybenzoic acid ethyl ester)	120-47-8	464

Index

Analyte name	CAS No#	Page No#
Ethyl Pentanoate (Valeric acid ethyl ester)	539-82-2	464
N-Ethylaniline	103-69-5	277
Ethylbenzene	100-41-4	278
Ethylbenzene-D10	25837-05-2	278, 538
2-Ethyl-1-butanol	97-95-0	278
2-Ethyl-1-butene	760-21-4	278
2-Ethylbutyric Acid	88-09-5	580
Ethylcyclohexane	1678-91-7	278
4-(1-Ethyl-1,4-dimethylpentyl) phenol	142731-63-3	580
Ethylene Glycol	107-21-1	278
Ethylene Glycol Bis(trichloroacetate)	2514-53-6	60
Ethylene Oxide	75-21-8	580
Ethylene thiourea D4	352431-28-8	183, 538
Ethylenediamine	107-15-3	278
Ethylenediamine-N,N'-diacetic Acid	5657-17-0	580
Ethyleneglycol Monobutyl Ether	111-76-2	278
2-Ethyl-1,3-hexandiol	94-96-2	278
2-Ethylhexanoic Acid	149-57-5	279
2-Ethyl-1-hexanol	104-76-7	279
2-Ethylhexyl Diphenyl Phosphate	1241-94-7	485
2-Ethylhexylnitrate	27247-96-7	580
N-(2-Ethylhexyl)-5-norbornen-2,3-dicarboximide (MGK 264)	113-48-4	183
Ethylmercury Chloride	107-27-7	580
Ethylmethanesulfonate	62-50-0	580
2-Ethyl-6-methylaniline	24549-06-2	183
2-Ethyl-6-methyl-2-chloroacetanilide	32428-71-0	184
4-(2-Ethyl-6-methylphenyl)-5-methyl-3-morpholinone	120375-14-6	580
1-Ethyl-naphthalene	1127-76-0	248
2-Ethyl-naphthalene	939-27-5	248
2-Ethylphenol	90-00-6	438
3-Ethylphenol	620-17-7	438
4-Ethylphenol	123-07-9	438
N-Ethyl-N-phenylurethane	1013-75-8	581

Analyte name	CAS No#	Page No#
2-Ethyltoluene	611-14-3	279
3-Ethyltoluene	620-14-4	279
4-Ethyltoluene	622-96-8	279
Etiocholan-3α-ol-17-one	53-42-9	338
Etobenzanid	79540-50-4	61
Etofenprox	80844-07-1	133
Etofenprox D5 (ethyl D5)	1705649-55-3	133, 538
Etoricoxib	202409-33-4	338
Etoxazole	153233-91-1	133
Etridiazole	2593-15-9	10
Etrimfos	38260-54-7	133
Etrofol (CPMC)	3942-54-9	581
Eucalyptol (Cineole)	470-82-6	581
Eugenol	97-53-0	581
Eulan WA new	55069-01-7	411
Famoxadone	131807-57-3	10
Famphur	52-85-7	133
Farnesol, mixture of isomers	4602-84-0	184
Fast Green FCF	2353-45-9	411
Febantel	58306-30-2	338
Fenamidone	161326-34-7	10
Fenamiosulf	140-56-7	11
Fenamiphos	22224-92-6	133
Fenamiphos-sulfone	31972-44-8	134
Fenamiphos-sulfoxide	31972-43-7	134
Fenarimol	60168-88-9	11
Fenazaflor	14255-88-0	134
Fenazaquin	120928-09-8	134
Fenazox	495-48-7	134
Fenbendazole	43210-67-9	338
Fenbendazole Sulphone (Oxfendazole-sulfone)	54029-20-8	338
Fenbuconazole	114369-43-6	11
Fenbuconazole (phenyl D5)	1398066-06-2	11, 538
Fenbutatin oxide	13356-08-6	134

Index

Analyte name	CAS No#	Page No#
Fenchlorazol-ethyl	103112-35-2	61
Fenchlorphos	299-84-3	134
Fenchlorphos-oxon	3983-45-7	134
Fenclorim	3740-92-9	184
Fenfluramine Hydrochloride	404-82-0	338
Fenfluthrin	75867-00-4	134
cis-Fenfluthrin	74337-23-8	135
Fenfuram	24691-80-3	11
Fenhexamid	126833-17-8	11
Fenitrothion	122-14-5	135
Fenitrothion D6 (O,O-dimethyl D6)	203645-59-4	135, 538
Fenitrothion-oxon	2255-17-6	135
Fenitrothion-oxon D6 (O,O-dimethyl D6)	1185155-54-7	135, 538
Fenobucarb	3766-81-2	135
Fenofibrate	49562-28-9	339
Fenoprofen Calcium Dihydrate	71720-56-4	339
Fenoprop	93-72-1	61
Fenoprop-methyl ester	4841-20-7	61
Fenoterol Hydrobromide	1944-12-3	339
Fenothiocarb	62850-32-2	135
Fenoxanil	115852-48-7	11
Fenoxaprop	95617-09-7	61
Fenoxaprop-ethyl	66441-23-4	61
Fenoxaprop-P (R-enantiomer)	113158-40-0	61
Fenoxaprop-P-ethyl (R-enantiomer)	71283-80-2	61
Fenoxycarb	72490-01-8	135
Fenpiclonil	74738-17-3	11
Fenpropathrin	39515-41-8	136
Fenpropidin	67306-00-7	11
Fenpropidin D10 (piperidine D10)	n/a	12, 539
Fenpropimorph	67564-91-4	12
Fenpyrazamine	473798-59-3	12
(E)-Fenpyroximate	134098-61-6	136
(Z)-Fenpyroximate	149054-53-5	136
Fenridazon	68254-10-4	184
Fenridazon-propyl	78778-15-1	184

Analyte name	CAS No#	Page No#
Fenson	80-38-6	136
Fensulfothion	115-90-2	136
Fensulfothion-oxon	6552-21-2	136
Fensulfothion-oxon-sulfone	6132-17-8	136
Fensulfothion-sulfone	14255-72-2	136
Fenthion	55-38-9	136
Fenthion D6 (dimethoxy D6)	1189662-83-6	137, 539
Fenthion-oxon	6552-12-1	137
Fenthion-oxon-sulfone	14086-35-2	137
Fenthion-oxon-sulfoxide	6552-13-2	137
Fenthion-sulfone	3761-42-0	137
Fenthion-sulfoxide	3761-41-9	137
Fentin-acetate	900-95-8	12
Fentin Chloride (Triphenyltin chloride)	639-58-7	12
Fentin Hydroxide (Triphenyltin hydroxide)	76-87-9	12
Fentin Oxide	1262-21-1	12
Fentrazamide	158237-07-1	61
Fentrazamide Metabolite 1	98377-35-6	62
Fenuron	101-42-8	62
Fenuron D5 (phenyl D5)	1219802-06-8	62, 539
Fenvalerate	51630-58-1	137
Fenvalerate D7 (isopropyl D7)	n/a	137, 539
Fenvalerate free acid metabolite	2012-74-0	138
Ferbam	14484-64-1	12
Ferimzone	89269-64-7	12
Finasteride	98319-26-7	339
Fipronil	120068-37-3	138
Fipronil Amide	205650-69-7	138
Fipronil des F3	154807-27-9	138
Fipronil-desulfinyl	205650-65-3	138
Fipronil Sulfide	120067-83-6	138
Fipronil Sulfone	120068-36-2	138
Flamprop (free acid)	58667-63-3	62
Flamprop-isopropyl	52756-22-6	62

Index

Analyte name	CAS No#	Page No#
Flamprop-M-isopropyl	63782-90-1	62
Flamprop-M-methyl	63729-98-6	62
Flamprop-methyl	52756-25-9	62
Flavomycin	11015-37-5	339
Flazasulfuron	104040-78-0	62
Floxacin	79660-72-3	339
Flocoumafen	90035-08-8	184
Flonicamid	158062-67-0	138
Flonicamid free acid	158063-66-2	138
Flonicamid-carboxylic acid (4-(Trifluoromethyl)nicotinoyl glycine)	207502-65-6	139
Florasulam	145701-23-1	63
Florfenicol	73231-34-2	339
Fluacrypyrim	229977-93-9	139
Fluazifop (free acid)	69335-91-7	139
Fluazifop-butyl	69806-50-4	63
Fluazifop-butyl D9 (n-butyl D9)	n/a	63, 539
Fluazifop-methyl	69335-90-6	63
Fluazifop-P (free acid)	83066-88-0	139
Fluazifop-P-butyl	79241-46-6	63
Fluazinam	79622-59-6	13
Fluazuron	86811-58-7	339
Flubendazole	31430-15-6	339
Flubendazole D3 (methyl D3)	1173021-08-3	340, 539
Flubendiamide	272451-65-7	139
Flubenzimine	37893-02-0	139
Flucarbazone-sodium	181274-17-9	63
Flucetosulfuron	412928-75-7	63
Fluchloralin	33245-39-5	63
Flucloxacillin Sodium	1847-24-1	340
Flucofuron	370-50-3	139
Fluconazole	86386-73-4	340
Fluconazole D4 (bismethylene D4)	1124197-58-5	340, 539
Flucycloxuron	113036-88-7	139
Flucythrinate	70124-77-5	139
Fludioxonil	131341-86-1	13

Analyte name	CAS No#	Page No#
Fluenetil (Lambrol)	4301-50-2	140
Fluensulfone	318290-98-1	184
Flufenacet	142459-58-3	63
Flufenacet-alcohol	54041-17-7	64
Flufenacet ethane sulfonic acid ESA sodium salt	947601-87-8	64
Flufenacet-methylsulfide	n/a	64
Flufenacet-methylsulfone	861391-80-2	64
Flufenacet-oxalat (OA)	201668-31-7	64
Flufenacet-thioglycolate sulfoxide	201668-33-9	64
Flufenoxuron	101463-69-8	140
Flufenpyr-ethyl	188489-07-8	64
Flufenzine	162320-67-4	140
Flufiprole	704886-18-0	140
Flumequine	42835-25-6	340
Flumetasone	2135-17-3	340
Flumethrin	69770-45-2	140
Flumetralin	62924-70-3	184
Flumetsulam	98967-40-9	64
Flumiclorac-pentyl	87546-18-7	64
Flumioxazin	103361-09-7	65
Flumorph	211867-47-9	13
Flunixin	38677-85-9	340
Flunixin D3 (methyl D3)	1015856-60-6	340, 539
Flunixin Meglumine	42461-84-7	340
Fluometuron	2164-17-2	65
Fluometuron-desmethyl (1-Methyl-3-[3-(trifluoromethyl) phenyl]urea)	3032-40-4	65
Fluopicolide	239110-15-7	13
Fluopyram	658066-35-4	13
Fluoranthene	206-44-0	248
Fluoranthene-D10	93951-69-0	248, 539
Fluorene	86-73-7	248
Fluorene-D10	81103-79-9	248, 540
9-Fluorenone	486-25-9	249

Index

Analyte name	CAS No#	Page No#
Fluorescein Sodium	518-47-8	411
Fluoroacetamide	640-19-7	184
Fluoroacetic Acid Sodium Salt	62-74-8	581
Fluorobenzene	462-06-6	438
2-Fluorobiphenyl	321-60-8	581
Fluorodifen	15457-05-3	65
Fluoroglycofen-ethyl	77501-90-7	65
Fluoroimide	41205-21-4	13
1-Fluoronaphthalene	321-38-0	249
2-Fluorophenol	367-12-4	438
4-Fluoro-3-phenoxy benzoic acid	77279-89-1	185
Fluorotrichloromethane (Trichlorofluoromethane)	75-69-4	391
Fluotrimazole	31251-03-3	13
Fluoxastrobin	361377-29-9	13
Fluoxetine Hydrochloride	56296-78-7	341
Flupropanate	756-09-2	65
Flupyrsulfuron-methyl sodium	144740-54-5	65
Fluquinconazole	136426-54-5	14
Flurandrenolide	1524-88-5	341
Flurenol	467-69-6	65
Flurenol Butyl Ester	2314-09-2	65
Flurenol Methyl Ester	1216-44-0	66
Fluridone	59756-60-4	66
Flurochloridone	61213-25-0	66
Fluroxypyr	69377-81-7	66
Fluroxypyr-1-methylheptylester	81406-37-3	66
Flurprimidol	56425-91-3	185
Flurtamone	96525-23-4	66
Flusilazole	85509-19-9	14
Flusulfamide	106917-52-6	14
Fluthiacet-methyl	117337-19-6	66
Fluticasone Propionate	80474-14-2	341
Flutolanil	66332-96-5	14
Flutriafol	76674-21-0	14
Flutrimazole	119006-77-8	14
tau-Fluvalinate	102851-06-9	140
Fluxapyroxad	907204-31-3	14

Analyte name	CAS No#	Page No#
Fluxofenim	88485-37-4	185
Folic Acid	59-30-3	341
Folpet	133-07-3	14
Folpet D4	1327204-12-5	14, 540
Fomesafen	72178-02-0	66
Fonofos	944-22-9	140
Foramsulfuron	173159-57-4	66
Forchlorfenuron	68157-60-8	185
Forchlorfenuron D5 (phenyl D5)	1398065-87-6	185, 540
Formaldehyde-2,4-dinitrophenylhydrazone (DNPH)	1081-15-8	279, 581
Formetanate hydrochloride	23422-53-9	140
Formic Acid (Methanoic acid)	64-18-6	464
Formic Acid Butyl Ester	592-84-7	464
Formic acid-tert-butyl Ester	762-75-4	465
Formic Acid Ethyl Ester	109-94-4	465
Formic Acid Propyl Ester	110-74-7	465
Formothion	2540-82-1	140
4-Formylaminophenazone	1672-58-8	341
Fosamine Ammonium Salt	25954-13-6	67
Fosetyl-aluminium D15	n/a	15, 540
Fosetyl Aluminium Salt	39148-24-8	15
Fosfomycin Calcium	26016-98-8	341
Fosinopril	98048-97-6	341
Fosthiazate	98886-44-3	141
Fructose (D-(-)-Fructose)	57-48-7	465
Fuberidazole	3878-19-1	15
Fuchsin (Basic Violet 14 HCl)	632-99-5	411
Fumaric Acid	110-17-8	465
Furalaxyl	57646-30-7	15
Furalaxyl (free acid)	118597-19-6	15
Furaltadone hydrochloride	3759-92-0	341
Furametpyr	123572-88-3	15
Furan	110-00-9	279
Furan-d4	6142-90-1	279, 540
2,5-Furandicarboxylic acid	3238-40-2	581
Furathiocarb	65907-30-4	141

Index

Analyte name	CAS No#	Page No#
Furazolidone	67-45-8	341
Furfural (Furane-2-carbaldehyde, 2-Furaldehyde)	98-01-1	15
Furfural-2,4-dinitrophenylhydrazone	2074-02-4	279
Furilazole	121776-33-8	185
Furmecyclox	60568-05-0	15
Furosemide	54-31-9	342
Galactose (D-(+)-Galactopyranose)	59-23-4	465
Gallic Acid	149-91-7	465
Gallic Acid Propyl Ester (Propyl Gallate)	121-79-9	465
Gasoline	8032-32-4	280
Natural Gasoline	8006-61-9	280
Gemfibrozil	25812-30-0	342
Genistein	446-72-0	342
Genite	97-16-5	342
Gentamicin Sulfate (Gentamycin sulfate)	1405-41-0	342
(±)-Geosmin	16423-19-1	581
Geraniol	106-24-1	342
Gibberellic Acid	77-06-5	185
Glibenclamide	10238-21-8	342
Glucose Anhydrous (D-(+)-Glucopyranose)	50-99-7	465
Glufosinate-N-acetyl	73634-73-8	67
Glufosinate ammonium	77182-82-2	67
Glufosinate-FMOC	1822429-60-6	67
DL-Glutamic Acid	617-65-2	582
Glutamine (L-Glutamine)	56-85-9	466
Glutaraldehyd-bis(2,4-dinitrophenylhydrazone)	5085-07-4	280
Glutaric Acid (Pentanedioic acid)	110-94-1	466
Glutaric Acid Dimethyl Ester	1119-40-0	466
Glycerol	56-81-5	280
Glyceryl Trinitrate (Nitroglycerin)	55-63-0	438
Glycidol	556-52-5	342
Glycine	56-40-6	466
Glycolic Acid	79-14-1	342

Analyte name	CAS No#	Page No#
Glycopyrronium Bromide	51186-83-5	343
Glyphosate (N-(Phosphonomethyl)glycine)	1071-83-6	67
Glyphosate 1,2-13C2 15N	1185107-63-4	67, 540
Glyphosate-FMOC	1373205-41-4	67
Glyphosate-N-methyl	24569-83-3	67
Glyphosate-N-nitroso mono sodium salt	56516-71-3	67
Glyphosine	2439-99-8	185
Gossypol	303-45-7	343
Griseofulvin	126-07-8	343
Guaifenesin	93-14-1	343
Guanidine Hydrochloride	50-01-1	582
Guanosine-5'-monophosphate	85-32-5	466
Halfenprox	111872-58-3	142
Halofenozide	112226-61-6	142
Halosulfuron-methyl	100784-20-1	68
Halowax 1000	58718-66-4	249
Halowax 1001	58718-67-5	249
Halowax 1013	12616-35-2	249
Halowax 1014	12616-36-3	249
Halowax 1051	70776-03-3	249
Halowax 1099	39450-05-0	249
Haloxon (O,O-Bis(2-chloroethyl) O-(3-chloro-4-methyl-7-coumarinyl) phosphate)	321-55-1	343
Haloxyfop (free acid)	69806-34-4	68
Haloxyfop-ethoxyethyl	87237-48-7	68
Haloxyfop-methyl	69806-40-2	68
Haloxyfop-R-methyl	72619-32-0	68
α-HCH 13C6	222966-66-7	141, 540
α-HCH D6	86194-41-4	141, 540
β-HCH	319-85-7	141
β-HCH D6	86194-42-5	141, 540
δ-HCH	319-86-8	141
δ-HCH D6	n/a	141, 541
ε-HCH	6108-10-7	141

Index

Analyte name	CAS No#	Page No#
γ-HCH D6	60556-82-3	142, 541
α-HCH (alpha-HCH)	319-84-6	142
n-Heneicosane	629-94-7	280
Heneicosanoic Acid	2363-71-5	466
Heneicosanoic acid-ethyl ester	28898-67-1	466
Heneicosanoic Acid Methyl Ester	6064-90-0	466
1-Heneicosanol	15594-90-8	466
n-Hentriacontane	630-04-6	280
Hepachlor-β-dihydro	14168-01-5	142
Heptachlor	76-44-8	142
(±)-cis-Heptachlorepoxide (Heptachlor-exo-epoxide, isomer B)	1024-57-3	142
(±)-trans-Heptachlorepoxide (Heptachlor-endo-epoxide, isomer A)	28044-83-9	142
1,2,4,5,6,9,10-Heptachlorodecane (CP-7)	890302-88-2	582
1,2,5,5,6,9,10-Heptachlorodecane (CP-8)	1278446-65-3	582
2,2'',3',4,4'',5,5'''-Heptachloro-m-terphenyl	n/a	216
n-Heptacosane	593-49-7	280
Heptacosanoic Acid	7138-40-1	467
n-Heptadecane	629-78-7	280
Heptadecanoic Acid	506-12-7	467
Heptadecanoic Acid Methyl Ester	1731-92-6	467
1-Heptadecanol	1454-85-9	281
cis-10-Heptadecenoic Acid	29743-97-3	467
2,2,4,4,6,8,8-Heptamethylnonane	4390-04-9	281
Heptanal	111-71-7	281
Heptanal-2,4-dinitrophenylhydrazone	2074-05-7	281
n-Heptane	142-82-5	281
1-Heptanethiol	1639-09-4	281
Heptanoic Acid (Heptylic acid)	111-14-8	467
Heptanoic Acid Ethyl Ester	106-30-9	467
Heptanoic Acid Methyl Ester	106-73-0	467
1-Heptanol	111-70-6	281
2-Heptanol	543-49-7	281
2-Heptanone	110-43-0	281

Analyte name	CAS No#	Page No#
3-Heptanone	106-35-4	282
4-Heptanone	123-19-3	282
n-Heptatriacontane	7194-84-5	282
trans-2-Heptene	14686-13-6	282
Heptenophos	23560-59-0	143
1-Heptylamine	111-68-2	282
4-Heptylphenol	1987-50-4	582
n-Heptyltin-trichloride	59344-47-7	582
Hesperidin	520-26-3	467
Hexabromobenzene	87-82-1	438
1,2,5,6,9,10-Hexabromocyclododecane	3194-55-6	582
Hexachloroacetone	116-16-5	15
Hexachlorobenzene	118-74-1	16, 439
Hexachlorobenzene 13C6	93952-14-8	16, 541
Hexachlorobutadiene	87-68-3	391
Hexachloro-1,3-butadiene 13C4	93951-70-3	391, 541
1,2,3,4,5,6-Hexachlorocyclohexane (HCH; BHC)	608-73-1	143
Hexachlorocyclopentadiene	77-47-4	582
1,2,4,5,9,10-Hexachlorodecane (CP-6)	890302-87-1	583
1,2,5,6,9,10-Hexachlorodecane (CP-4)	189350-94-5	583
Hexachloroethane	67-72-1	391
1,4,5,6,7,7-Hexachloro-N-(ethylmercury)-5-norbornene-2,3-dicarboximide	2597-93-5	583
Hexachlorophene	70-30-4	16
Hexachloropropene (Perchloropropene)	1888-71-7	391
2,2'',4,4'',5,5'''-Hexachloro-p-terphenyl	n/a	216
3,3'',4,4'',5,5'''-Hexachloro-p-terphenyl	1312994-27-6	216
Hexaconazole	79983-71-4	16
n-Hexacontane	7667-80-3	282
n-Hexacosane	630-01-3	282
Hexacosanoic Acid (Cerotic acid)	506-46-7	467
Hexacosanoic Acid Methyl Ester	5802-82-4	468

Index

Analyte name	CAS No#	Page No#
1-Hexacosanol (Cerotyl alcohol)	506-52-5	282
n-Hexadecane	544-76-3	282
Hexadecanoic Acid Ethyl Ester	628-97-7	468
Hexadecanoic Acid Methyl Ester	112-39-0	468
1-Hexadecanol (Cetyl alcohol)	36653-82-4	283
1-Hexadecene	629-73-2	283
Hexaflumuron	86479-06-3	143
Hexamethylbenzene	87-85-4	283
Hexamethylene-1,6-diisocyanate	822-06-0	583
Hexamethylphosphoramide (Hexametapol)	680-31-9	583
Hexanal (Capronaldehyde)	66-25-1	283
Hexanal-2,4-dinitrophenylhydrazone	1527-97-5	283
n-Hexane	110-54-3	283
1,6-Hexanedicarboxylic Acid (Octanedioic acid, Suberic acid)	505-48-6	468
2,5-Hexanediol	2935-44-6	283
Hexanoic Acid (Caproic acid)	142-62-1	468
Hexanoic Acid Ethyl Ester	123-66-0	468
Hexanoic Acid Methyl Ester	106-70-7	468
1-Hexanol	111-27-3	283
2-Hexanol	626-93-7	283
3-Hexanol	623-37-0	284
2-Hexanone	591-78-6	284
n-Hexatriacontane	630-06-8	284
Hexazinone (Velpar)	51235-04-2	68
Hexazinone D6 (N,N-dimethyl D6)	1219804-22-4	68, 541
1-Hexene	592-41-6	284
cis-3-Hexen-1-ol	928-96-1	583
cis-3-Hexenyl acetate	3681-71-8	583
Hexestrol	84-16-2	343
Hexogen (Hexahydro-1,3,5-trinitro-1,3,5-triazine)	121-82-4	439
1-Hexylamine	111-26-2	284
n-Hexylbenzene	1077-16-3	284
Hexylcinnamal	101-86-0	583
4-Hexylphenol	2446-69-7	439
Hexythiazox	78587-05-0	143

Analyte name	CAS No#	Page No#
HHCB (Galaxolide)	1222-05-5	582
Histamine Dihydrochloride	56-92-8	343
DL-Histidine	4998-57-6	583
Homosalate (3,3,5-Trimethylcyclohexyl salicylate)	118-56-9	584
Hydramethylnon	67485-29-4	143
Hydrochlorothiazide	58-93-5	343
Hydrocortisone (Cortisol)	50-23-7	343
Hydrocortisone Acetate	50-03-3	344
Hydrocortisone 17-Butyrate	13609-67-1	344
Hydroquinone (Benzene-1,4-diol)	123-31-9	439
Hydroquinone Monobenzylether (4-(Benzyloxy)phenol)	103-16-2	344
4-Hydroxybenzoic Acid n-Heptyl Ester	1085-12-7	468
4-Hydroxybenzonitrile	767-00-0	439
6-Hydroxychrysene	37515-51-8	249
7-Hydroxycitronellal (7-Hydroxy-3,7-dimethyloctanal)	107-75-5	584
4-Hydroxycoumarin	1076-38-6	584
7-Hydroxycoumarin	93-35-6	584
2-Hydroxydiphenylmethane	28994-41-4	439
N-(2-Hydroxyethyl)-N-methylper-fluorooctanesulfonamide	24448-09-7	584
1-Hydroxyibuprofen (Ibuprofen-1-hydroxy)	53949-53-4	344
2-Hydroxyibuprofen (Ibuprofen-2-hydroxy)	51146-55-5	344
Hydroxyisohexyl-3-cyclohexene carboxaldehyde	31906-04-4	584
Hydroxymethyl phosphonic acid	2617-47-2	68
2-Hydroxymethylfuran (2-Furfuryl alcohol)	98-00-0	584
4-Hydroxy-4-methyl-2-pentanone	123-42-2	284
2-Hydroxy-2-methylpropiophenone	7473-98-5	584
2-Hydroxy-1,4-naphthoquinone	83-72-7	584
6-Hydroxynicotinic Acid (Nicotinic acid-6-hydroxy)	5006-66-6	468
2-Hydroxy-4-nitroaniline	121-88-0	439
2-Hydroxy-5-nitroaniline	99-57-0	439
1-Hydroxyphenanthrene	2433-56-9	250

Index

Analyte name	CAS No#	Page No#
2-Hydroxyphenanthrene	605-55-0	250
3-Hydroxyphenanthrene	605-87-8	250
4-Hydroxy-phenanthrene	7651-86-7	250
9-Hydroxyphenanthrene	484-17-3	250
Hydroxyprogesterone Acetate (17α-Hydroxyprogesterone 17-acetate)	302-23-8	344
17α-Hydroxyprogesterone	68-96-2	344
1-Hydroxypyrene	5315-79-7	250
2-Hydroxy-5-(trifluoromethyl) pyridine	33252-63-0	585
Hygromycin B	31282-04-9	344
Hymexazol	10004-44-1	16
Hyoscine Hydrobromide Trihydrate (Scopolamine hydrobromide trihydrate)	6533-68-2	344
Ibuprofen	15687-27-1	345
Icaridin (Picaridin)	119515-38-7	185
Ifosfamide	3778-73-2	345
Imazalil sulfate	58594-72-2	16
Imazamethabenz (free acid)	100728-84-5	68
Imazamethabenz-methyl	81405-85-8	69
Imazamox	114311-32-9	69
Imazapic	104098-48-8	69
Imazapyr	81334-34-1	69
Imazaquin	81335-37-7	69
Imazethapyr	81335-77-5	69
Imazosulfuron	122548-33-8	69
Imibenconazole	86598-92-7	16
Imibenconazole-desbenzyl (N-(2,4-Dichlorophenyl)-1H-1,2,4-triazole-1-acetamide)	154221-27-9	16
Imicyafos	140163-89-9	186
Imidacloprid	138261-41-3	143
Imidacloprid D4 (imidazolidin-4,4,5,5 D4)	1015855-75-0	143, 541
Imidacloprid-guanidine Hydrochloride	127202-53-3	143
Imidacloprid Impurity 1	105828-41-9	143
Imidacloprid-olefin	115086-54-9	144

Analyte name	CAS No#	Page No#
Imidacloprid-urea	120868-66-8	144
Imidaclothiz	105843-36-5	144
Imidazole	288-32-4	345
Imidazolidine-2-thione (Ethylene thiourea)	96-45-7	186
2-Imidazolidinone	120-93-4	585
Imidocarb	27885-92-3	345
Imidocarb Dipropionate	55750-06-6	345
Imino bis(methylphosphonic Acid)	17261-34-6	585
Iminoctadine triacetate	57520-17-9	16
Iminoctadine tris(albesilate)	169202-06-6	17
Iminodiacetic Acid	142-73-4	585
Iminostilbene (5H-Dibenzo[b,f] azepine)	256-96-2	345
Imiprothrin	72963-72-5	144
Inabenfide	82211-24-3	69
Indane	496-11-7	284
Indanofan	133220-30-1	69
Indaziflam	950782-86-2	70
Indene	95-13-6	284
Indeno(1,2,3-c,d)fluoranthene	193-43-1	250
Indeno(1,2,3-c,d)pyrene	193-39-5	250
Indeno(1,2,3-c,d)pyrene D12	203578-33-0	250, 541
Indigotine (Indigotindisulfonate Sodium; E132)	860-22-0	411
Indole	120-72-9	251
3-Indolylacetic Acid Ethyl Ester	778-82-5	186
3-Indolylacetic Acid	87-51-4	186
3-Indolylacetoneitrile	771-51-7	585
4-(3-Indolyl)butyric Acid	133-32-4	186
3-(3-Indolyl)propionic Acid	830-96-6	186
Indometacin (Indomethacine)	53-86-1	345
Indoxacarb	144171-61-9	144
Inosine	58-63-9	345
5'-Inosinic Acid Disodium Salt Hydrate	352195-40-5	345
Iodipamide	606-17-7	346

Index

Analyte name	CAS No#	Page No#
Iodixanol	92339-11-2	346
Iodocarb (IPBC)	55406-53-6	17
Iodofenphos	18181-70-9	144
Iodomethane	74-88-4	391
Iodosulfuron-methyl sodium	144550-36-7	70
Iohexol	66108-95-0	346
Iomeprol	78649-41-9	346
Iopamidol	60166-93-0	346
Iopanoic Acid	96-83-3	346
Iopromide	73334-07-3	346
Iothalamate Meglumine (Iotalamic acid meglumine salt)	13087-53-1	346
Ioxynil	1689-83-4	70
Ioxynil-methyl	3336-40-1	70
Ioxynil-octanoate	3861-47-0	70
Ipconazole	125225-28-7	17
Ipfencarbazone	212201-70-2	70
Iprobenfos	26087-47-8	17
Iprodione	36734-19-7	17
Iprodione D7 (isopropyl D7)	n/a	17, 541
Iprodione des-(N-isopropylcarboxamid)	27387-87-7	17
Iprodione isomer 1 (Isoiprodione)	63637-89-8	17
Ipronidazole	14885-29-1	346
Ipronidazole D3 (N-methyl D3)	1015855-83-0	347, 541
Ipronidazole-hydroxy	35175-14-5	347
Ipronidazole-hydroxy D3	1156508-86-9	347, 542
Iprovalicarb	140923-17-7	18
Irgarol 1051	28159-98-0	186
Iron Oxide Pigment Black	1317-61-9	412
Iron Pigment Red (E172)	1309-37-1	412
Iron Pigment Yellow (E172)	51274-00-1	412
Isazofos	42509-80-8	144
D(-)-Isoascorbic Acid	89-65-6	469
Isobenzan (Telodrin)	297-78-9	144
Isobenzofuranone	87-41-2	585
Isobornyl Acetate	125-12-2	585

Analyte name	CAS No#	Page No#
Isobutyl acetate (Acetic acid isobutyl ester)	110-19-0	469
Isobutyl Parahydroxybenzoate (4-Hydroxybenzoic acid isobutyl ester)	4247-02-3	469
Isobutylbenzene	538-93-2	285
2-Isobutyl-3-methoxypyrazine	24683-00-9	585
Isobutyraldehyde	78-84-2	285
Isobutyric Acid	79-31-2	585
Isobutyric Acid Methyl Ester	547-63-7	586
Isocarbamid	30979-48-7	70
Isocarbofos	24353-61-5	144
Isochloridazon	162354-96-3	70
Isodecanol-ethoxylate (technical)	61827-42-7	285
Isododecane	31807-55-3	285
Isodrin	465-73-6	145
Isoeugenol	97-54-1	469
Isofenphos	25311-71-1	145
Isofenphos-des-N-isopropyl	25205-08-7	145
Isofenphos-methyl	99675-03-3	145
Isofenphos-oxon	31120-85-1	145
Isofenphos-des-N-isopropyl-oxon (Isopropyl 2-[Amino(ethoxy) phosphoryl]oxybenzoate)	31120-83-9	145
DL-Isoleucine	443-79-8	586
Isomalathion	3344-12-5	145
Isomethiozin	57052-04-7	70
α-Isomethylionone	127-51-5	586
Isononylphenol-ethoxylate (technical mixture)	37205-87-1	586
Isonoruron	28805-78-9	71
4-iso-Octylphenol-di-ethoxylate	2315-61-9	586
Isophorone	78-59-1	285
Isophthalic Acid	121-91-5	485
Isophthalic acid, bis-n-octyl ester	4654-18-6	485
Isophthalic Acid Dimethyl Ester	1459-93-4	485
Isophthalic Acid Diphenyl Ester	744-45-6	485
Isoprene	78-79-5	285
Isoprobenphos (O,O-Diethyl S-benzyl Phosphorothioate)	13286-32-3	145

Index

Analyte name	CAS No#	Page No#
Isoprocarb	2631-40-5	145
Isopropalin	33820-53-0	71
2-Isopropoxyphenol	4812-20-8	186
Isopropyl Myristate (Myristic acid isopropyl ester)	110-27-0	586
4-Isopropylaniline	99-88-7	186
N-Isopropylaniline	768-52-5	439
Isopropylbenzene (Cumene)	98-82-8	285
Isopropylcyclohexane	696-29-7	285
2-Isopropyl-3-methoxypyrazine	25773-40-4	586
4-Isopropyl-3-methylphenol	3228-02-2	440
2-Isopropyl-6-methyl-4-pyrimidinol	2814-20-2	187
2-Isopropyl-naphthalene	2027-17-0	187
Isopropyl-nitrate	1712-64-7	586
2-Isopropyl-phenol	88-69-7	440
3-Isopropyl-phenol	618-45-1	440
4-Isopropyl-phenol	99-89-8	440
1-(4-Isopropylphenyl)-3-methylurea	34123-57-4	71
1-(4-Isopropylphenyl)urea	56046-17-4	71
2-Isopropylthioxantone (ITX)	5495-84-1	586
2-Isopropyltoluene (o-Cymol)	527-84-4	285
3-Isopropyltoluene (m-Cymol)	535-77-3	286
4-Isopropyltoluene (p-Cymol)	99-87-6	286
Isoprothiolane	50512-35-1	18
Isoproturon	34123-59-6	71
Isoproturon D6	1007461-76-8	71, 542
Isopyrazam	881685-58-1	18
Isoquinoline	119-65-3	251
Isouron	55861-78-4	71
Isovaleraldehyd-2,4-dinitrophenylhydrazone	2256-01-1	286
Isovaleric Acid	503-74-2	587
2-Isovaleryl-1.3-indanedione	83-28-3	587
Isoxaben	82558-50-7	71
Isoxadifen-ethyl	163520-33-0	187
Isoxaflutole	141112-29-0	71
Isoxathion	18854-01-8	146

Analyte name	CAS No#	Page No#
Itraconazole	84625-61-6	347
Ivermectin	70288-86-7	347
Josamycin	16846-24-5	347
Kadethrin	58769-20-3	146
Kanamycin Sulfate	25389-94-0	347
Karanjin	521-88-0	146
Karbutilate	4849-32-5	72
Kasugamycin Hydrochloride	19408-46-9	18
Kelevan	4234-79-1	146
Ketamine Hydrochloride	1867-66-9	347
Ketoconazole	65277-42-1	347
Ketoprofen	22071-15-4	348
(±)-Ketoprofen D3 (methyl D3)	159490-55-8	348, 542
Kinetin (6-Furfurylaminopurine)	525-79-1	587
Kinoprene	42588-37-4	146
Kresoxim-methyl	143390-89-0	18
Lactic acid	50-21-5	348
D-Lactitol Monohydrate	81025-04-9	469
Lactofen	77501-63-4	72
α-D-(+)-Lactose	63-42-3	469
Lamprecide	654-66-0	187
Lansoprazole	103577-45-3	348
Lasalocid A sodium salt	25999-20-6	348
Lead Arsenate	7784-40-9	146
Lenacil	2164-08-1	72
Lenacil (2,2,2',2'-cyclohexyl D4)	n/a	72, 542
Leptophos	21609-90-5	146
Leptophos-oxon	25006-32-0	146
DL-Leucine	328-39-2	587
Leucocrystal Violet	603-48-5	412
Leucocrystal Violet D6	1173023-92-1	412, 542
Leucomalachite Green	129-73-7	412
Leucomalachite Green D6 (phenylmethin D6)	1173021-13-0	412, 542
Leucomycin	1392-21-8	348
Levamisole	14769-73-4	348

Index

Analyte name	CAS No#	Page No#
Levamisole Hydrochloride	16595-80-5	348
Levocarnitine (L-Carnitin)	541-15-1	348
Levomenthol ((-)-Menthol)	2216-51-5	587
Levosulpiride ((S)-(-)-Sulpiride)	23672-07-3	349
Levothyroxine (L-Thyroxine)	51-48-9	349
Lidocaine (Lignocaine)	137-58-6	349
Lidocaine Hydrochloride Monohydrate	6108-05-0	349
Limonene	138-86-3	469
R-(+)-Limonene	5989-27-5	469
S-(-)-Limonene	5989-54-8	469
Linalol (Linalool)	78-70-6	470
Lincomycin Hydrochloride Monohydrate	7179-49-9	349
Lindane (γ-HCH)	58-89-9	146
Linoelaidic Methyl Ester	2566-97-4	470
Linoleic Acid	60-33-3	470
Linoleic Acid Methyl Ester	112-63-0	470
Linolelaidic Acid	506-21-8	470
Linolenic acid	463-40-1	470
Linolenic Acid Methyl Ester	301-00-8	470
Linuron	330-55-2	72
Linuron D6	1219804-76-8	72, 542
Lomefloxacin Hydrochloride	98079-52-8	349
Loperamide Hydrochloride	34552-83-5	349
Lufenuron Anhydrous	103055-07-8	147
DL-Lysine	70-54-2	587
Mabuterol Hydrochloride	54240-36-7	349
Maduramicin ammonium salt	84878-61-5	349
Malachite green oxalate salt	2437-29-8	413
Malachite Green D5 Oxalate Salt	n/a	412, 543
Malachite Green D5 Picrate	1258668-21-1	412, 543
Malaoxon	1634-78-2	147
Malathion	121-75-5	147
Malathion D10 (diethyl D10)	347841-48-9	147, 543
Malathion D6 (dimethyl D6)	1189877-72-2	147, 543
Malathion D7 (dimethyl D6,3-D1)	352438-94-9	147, 543
Malathion dicarboxylic acid	1190-28-9	147

Analyte name	CAS No#	Page No#
Maleic Acid	110-16-7	587
Maleic Acid Dimethyl Ester	624-48-6	587
Maleic Hydrazide	123-33-1	187
Maleic Hydrazide D2	n/a	187, 543
DL-Malic acid	6915-15-7	587
Malonic Acid (Propanedioic acid)	141-82-2	470
Malonic Acid Dimethyl Ester	108-59-8	470
Malonitrile	109-77-3	286
Malonoben (Tyrphostin A9)	10537-47-0	440
Maltol	118-71-8	588
D-(+)-Maltose Monohydrate	6363-53-7	471
Maltotriose Hydrate	207511-08-8	471
Mancozeb	8018-01-7	18
Mandipropamid	374726-62-2	18
Maneb	12427-38-2	18
Mannitol (D-Mannit)	69-65-8	471
D-(+)-Mannose	3458-28-4	471
Mapenterol hydrochloride	54238-51-6	350
Marbofloxacin	115550-35-1	350
MCPA ((4-Chloro-2-methylphenoxy)acetic Acid)	94-74-6	72
MCPA Butoxyethyl Ester	19480-43-4	72
MCPA Butyl Ester	1713-12-8	73
MCPA D3 (phenyl D3)	352431-14-2	73, 542
MCPA D6 (methyl-D3,phenoxy-D3)	n/a	73, 542
MCPA-dimethylammonium	2039-46-5	73
MCPA Ethyl Ester	2698-38-6	73
MCPA-2-ethylhexyl ester	29450-45-1	73
MCPA-methyl ester	2436-73-9	73
MCPA-sodium salt monohydrate	3653-48-3	73
MCPA-thioethyl	25319-90-8	73
MCPB (4-(4-Chloro-2-methylphenoxy)butanoic Acid)	94-81-5	74
MCPB D6 (ring D3, methyl D3) (4-(4-Chloro-2-trideuteriomethyl-3,5,6-trideuteriophenoxy) butanoic Acid)	n/a	74, 543
MCPB Ethyl Ester	10443-70-6	74

Index

Analyte name	CAS No#	Page No#
MCPB-ethylhexylester	94232-74-3	74
MCPB-methyl ester	57153-18-1	74
Mebendazole	31431-39-7	350
Mebendazole D3 (methyl D3)	1173021-87-8	350, 543
Mebendazole-amine ((2-Amino-1H-benzimidazol-5-yl) phenylmethanone)	52329-60-9	350
Mebendazole-5-hydroxy	60254-95-7	350
Mebenil	7055-03-0	18
Mecarbam	2595-54-2	147
Mecarbam D10 (diethoxy D10)	n/a	147, 543
Mecillinam	32887-01-7	350
Meclofenamic Acid	644-62-2	350
Meclofenoxate Hydrochloride	3685-84-5	350
Mecoprop (MCPPE)	93-65-2	74
Mecoprop-2-butoxyethyl ester	23359-62-8	74
Mecoprop-n/iso-butyl ester	n/a	75
Mecoprop D3 (phenyl D3)	352431-15-3	74, 544
Mecoprop D6 (ring D3, methyl D3)	1705649-54-2	74, 544
Mecoprop methyl Ester	23844-56-6	75
Mecoprop-methyl ester D3 (ring D3)	n/a	75, 544
Mecoprop 2-Octyl Ester	28473-03-2	75
Mecoprop-2,4,4-trimethylpentylester	217487-13-3	75
2,6-Mecoprop	35851-12-8	75
D(+)Mecoprop	16484-77-8	75
Mecoprop-2-ethylhexyl ester	71526-69-7	75
Medinoterb Acetate	2487-01-6	75
Medroxyprogesterone	520-85-4	351
Medroxyprogesterone 17-Acetate	71-58-9	351
Mefenacet	73250-68-7	76
Mefenamic Acid	61-68-7	351
Mefenpyr-diethyl	135590-91-9	187
Mefluidide	53780-34-0	76
Megestrol 17-Acetate	595-33-5	351
Melamine	108-78-1	588
Melamine Triamine-15N3	287476-11-3	544, 588

Analyte name	CAS No#	Page No#
Melatonin	73-31-4	351
Melengestrol Acetate	2919-66-6	351
Melezitose Hydrate	10030-67-8	471
α-D(+)-Melibiose Hydrate	585-99-9	471
Meloxicam	71125-38-7	351
Meloxicam D3 (2-methyl D3)	942047-63-4	351, 544
Menadione (Vitamin K3)	58-27-5	471
Menadione Sodium Bisulfite	130-37-0	471
Menbutone	3562-99-0	351
Menthol	1490-04-6	588
Mepanipyrim	110235-47-7	19
Mepanipyrim-2-hydroxypropyl	204571-52-8	19
Mephosfolan	950-10-7	148
Mepiquat Chloride (1,1-Dimethylpiperidinium chloride)	24307-26-4	187
Mepiquat iodide D3 (methyl-d3)	32317-85-4	187, 544
Mepronil	55814-41-0	19
Meptyldinocap	131-72-6	19
Meptyldinocap-phenol	3687-22-7	19
Mepyramine Maleate	59-33-6	352
Merphos (Tributylphosphoro-trithioite)	150-50-5	188
Mesalazine (5-Aminosalicylic acid)	89-57-6	588
Mesosulfuron-methyl	208465-21-8	76
Mesotrione	104206-82-8	76
Mestranol	72-33-3	352
Metacycline Hydrochloride	3963-95-9	352
Metaflumizone	139968-49-3	148
Metalaxyl	57837-19-1	19
Metalaxyl-M	70630-17-0	19
Metaldehyde	108-62-3	188
Metam-sodium dihydrate	6734-80-1	76
(±)-Metamifop	256412-89-2	76
Metamitron	41394-05-2	76
Metamitron-desamino	36993-94-9	76
Metamizole Sodium	68-89-3	352

Index

Analyte name	CAS No#	Page No#
D,L-Metanephrine Hydrochloride	881-95-8	352
Metazachlor	67129-08-2	76
Metazachlor Ethane Sulfonic Acid (2-[(2,6-Dimethylphenyl) (1H-pyrazol-1-yl)methyl]amino]-2-oxoethanesulfonic acid)	172960-62-2	77
Metazachlor oxanilic acid (OA)	1231244-60-2	77
Metconazole	125116-23-6	19
Methabenzthiazuron	18691-97-9	77
Methacrifos	62610-77-9	148
Methacrylaldehyde-2,4-dinitrophenylhydrazone	5077-73-6	286
Methacrylamide	79-39-0	588
Methacrylic Acid Ethyl Ester	97-63-2	286
Methacrylonitrile	126-98-7	286
Methamidophos	10265-92-6	148
Methamidophos (dimethyl D6)	1219799-41-3	148, 544
Methamidophos-N-methyl (O,S-Dimethylmethylphosphoramido thioate)	28167-49-9	148
Methane	74-82-8	286
Methanol	67-56-1	286
Methapyrilene (2-[[2-(Dimethylamino)ethyl]-2-thenylamino]pyridine)	91-80-5	352
Methapyrilene Hydrochloride	135-23-9	352
Methasulfocarb	66952-49-6	19
Methazole	20354-26-1	352
Methfuroxam	28730-17-8	20
Methidathion	950-37-8	148
Methiocarb (Mercaptodimethur)	2032-65-7	148
Methiocarb D3 (N-methyl D3)	1581694-94-1	148, 544
Methiocarb-phenol-sulfone	14763-62-3	149
Methiocarb-phenol-sulfoxide	22454-92-8	149
Methiocarb-sulfone	2179-25-1	149
Methiocarb sulfoxide	2635-10-1	149
dl-Methionine	59-51-8	588
Methomyl	16752-77-5	149

Analyte name	CAS No#	Page No#
Methomyl D3	1398109-07-3	149, 544
Methomyl-oxime	13749-94-5	149
Methomyl-sulfone	55620-24-1	149
Methoprene	40596-69-8	149
trans-S-Methoprene-epoxide	65277-96-5	150
S-Methoprene	65733-16-6	150
Methoprotryne	841-06-5	77
D-Methorphan Hydrobromide	125-69-9	352
Methothrin	34388-29-9	150
Methoxychlor	72-43-5	150
Methoxychlor D14 (bis(4-methoxyphenyl-D7))	1644449-82-0	150, 545
2,4'-Methoxychlor	30667-99-3	150
4,4'-Methoxychlor olefin (1,1-Dichloro-2,2-bis(p-methoxyphenyl)ethene)	2132-70-9	150
2-Methoxy-3,5-dimethylpyrazine	92508-08-2	588
3-Methoxy-2,5-dimethylpyrazine	19846-22-1	588
2-Methoxyethanol (Ethleneglycol monomethyl ether)	109-86-4	287
2-Methoxyethyl Acetate	110-49-6	287
Methoxyfenozone	161050-58-4	150
2-Methoxy-5-methylaniline	120-71-8	413
Methoxyphenone (3,3'-Dimethyl-4-methoxybenzophenone)	41295-28-7	77
4-Methoxyphenylacetone (4'-Methoxyacetophenone)	100-06-1	589
4-Methoxy-1,3-phenylenediamine (2,4-Diaminoanisole)	615-05-4	413
1-Methoxy-2-propanol	107-98-2	287
Methyl Acetate (Acetic acid methyl ester)	79-20-9	471
Methyl Anthranilate (2-Aminobenzoic acid methyl ester)	134-20-3	589
Methyl N-(3,4-Dichlorophenyl) carbamate (Swep)	1918-18-9	77
Methyl Erucate (Erucic acid methyl ester)	1120-34-9	472
Methyl Formate (Formic acid methyl ester)	107-31-3	472

Index

Analyte name	CAS No#	Page No#
Methyl N-(3-Hydroxyphenyl) carbamate	13683-89-1	188
Methyl Methacrylate (Methacrylic acid methyl ester)	80-62-6	287
Methyl Parahydroxybenzoate (4-Hydroxybenzoic acid methyl ester)	99-76-3	472
Methyl Salicylate	119-36-8	589
Methyl Stearate (Octadecanoic acid methyl ester)	112-61-8	472
Methyl Tricosanoate (Tricosanoic acid methyl ester)	2433-97-8	472
N-Methylacetamide	79-16-3	589
3-Methylaniline (m-Toluidine)	108-44-1	440
4-Methylaniline (p-Toluidine)	106-49-0	440
N-Methylaniline	100-61-8	440
1-Methylantracene	610-48-0	251
2-Methylantracene	613-12-7	251
9-Methylantracene	779-02-2	251
Methylarsonic acid mono sodium salt sesquihydrat	2163-80-6	77
3-Methylbenzaldehyde (m-Tolualdehyde)	620-23-5	440
N-Methylbenzenesulfonamide	5183-78-8	589
2-Methylbenzofuran	4265-25-2	251
5-Methylbenzothiophene	14315-14-1	251
5-Methyl-1H-benzotriazole	136-85-6	589
3-(4'-Methyl)benzylidene-bornan-2-one	36861-47-9	589
2-Methyl-1-butanol	137-32-6	287
2-Methyl-2-butanol	75-85-4	287
3-Methyl-1-butanol	123-51-3	287
3-Methyl-2-butanone	563-80-4	287
3-Methylbutylamine	107-85-7	287
Methyl-tert-butylether	1634-04-4	288
3-Methylbutyraldehyde	590-86-3	288
2-Methylbutyric Acid Methyl Ester	868-57-5	288
2-Methylbutyric Acid	116-53-0	589
3-Methylcholanthrene	56-49-5	251
1-Methylchrysene	3351-28-8	251
2-Methylchrysene	3351-32-4	252
4-Methylchrysene	3351-30-2	252

Analyte name	CAS No#	Page No#
5-Methylchrysene	3697-24-3	252
6-Methylchrysene	1705-85-7	252
Methylcyclohexane	108-87-2	288
2-Methylcyclohexanone	583-60-8	288
3-Methylcyclohexanone	591-24-2	288
4-Methylcyclohexanone	589-92-4	288
3-Methylcyclopentadecanone	541-91-3	288
(Methylcyclopentadienyl) manganese Tricarbonyl	12108-13-3	589
Methylcyclopentane	96-37-7	288
N-Methyldiethanolamine	105-59-9	590
Methyldymron	42609-73-4	77
Methylene chloride (Dichloromethane)	75-09-2	391
Methylene Dithiocyanate	6317-18-6	188
4,4'-Methylene-bis(2-chloroaniline)	101-14-4	413
Methyleugenol (4-Allyl-1,2-dimethoxybenzene)	93-15-2	590
1-Methylfluoranthene	25889-60-5	252
2-Methylfluoranthene	33543-31-6	252
3-Methylfuran	930-27-8	289
5-Methyl-2-furfural	620-02-0	590
3-Methylheptane	589-81-1	289
2-Methyl-3-heptanone	13019-20-0	289
5-Methyl-2-hexanone	110-12-3	289
4-Methylimidazol	822-36-6	353
4-Methylimidazolidine-2-thione (Propylene thiourea)	2122-19-2	188
3-Methylindole	83-34-1	252
2-Methylisoborneol	2371-42-8	590
2-Methyl-4-isothiazolin-3-one Hydrochloride	26172-54-3	188
2-Methyl-4-isothiazolin-3-one	2682-20-4	188
Methyl-mercury-chloride	115-09-3	590
Methylmethanesulfonate	66-27-3	590
1-Methylnaphthalene	90-12-0	252
1-Methylnaphthalene D10	38072-94-5	252, 545
2-Methylnaphthalene	91-57-6	253
1-Methyl-3-nitroguanidine	4245-76-5	590

Index

Analyte name	CAS No#	Page No#
2-Methyl-4-nitroimidazole (2-Methyl-5-nitroimidazole)	696-23-1	590
3-Methyl-4-nitroimino-tetrahydro-1,3,5-oxadiazine	153719-38-1	188
2-Methyl-3-nitrophenol	5460-31-1	441
3-Methyl-2-nitrophenol	4920-77-8	441
3-Methyl-4-nitrophenol	2581-34-2	188
4-Methyl-2-nitrophenol	119-33-5	441
5-Methyl-2-nitrophenol	700-38-9	441
Methyl-pentachlorophenyl sulfide	1825-19-0	590
2-Methylpentane	107-83-5	289
3-Methylpentane	96-14-0	289
(2RS)-2-Methylpentanoic Acid (2-Methylvaleric acid)	97-61-0	591
2-Methyl-1-pentanol	105-30-6	289
3-Methyl-3-pentanol	77-74-7	289
4-Methyl-pentan-2-ol (2-Methyl-4-pentanol)	108-11-2	289
4-Methylpentan-2-one (MIBK)	108-10-1	591
4-Methyl-2-pentanone-2,4-dinitrophenylhydrazone	1655-42-1	290
4-Methyl-1-pentene	691-37-2	290
cis-4-Methyl-2-pentene	691-38-3	290
4-Methylpent-3-en-2-one (Mesityloxide)	141-79-7	591
N-Methylperfluorooctanesulfo-namide	31506-32-8	591
1-Methylphenanthrene	832-69-9	253
2-Methylphenanthrene	2531-84-2	253
3-Methylphenanthrene	832-71-3	253
9-Methylphenanthrene	883-20-5	253
2-Methylphenol (o-Cresol)	95-48-7	441
3-Methylphenol (m-Cresol; Metacresol)	108-39-4	189
4-Methylphenol (p-Cresol)	106-44-5	441
2-(4'-Methylphenyl)propionic Acid 3,3,3-D3	n/a	290, 545
2-(4'-Methylphenyl)propionic Acid ((2RS)-2-(4-Methylphenyl) propanoic acid)	938-94-3	290
3-Methylphosphinicopropionic Acid	15090-23-0	150

Analyte name	CAS No#	Page No#
Methylprednisolone Acetate	53-36-1	353
6α-Methylprednisolone	83-43-2	353
2-Methylpropanol (2-Methyl-1-propanol)	78-83-1	290
1-Methylpyrene	2381-21-7	253
2-Methylpyridine (2-Picolin)	109-06-8	591
N-Methylpyrrolidone (1-Methyl-2-pyrrolidone)	872-50-4	353
2-Methylquinoline	91-63-4	253
3-Methyl-quinoxaline-2-carboxylic acid	74003-63-7	189
α-Methylstyrene	98-83-9	290
trans-β-Methylstyrene	873-66-5	290
4-(Methylsulfonyl)-2-nitrobenzoic acid	110964-79-9	591
17α-Methyltestosterone	58-18-4	353
5-Methyl-1,3,4-thiadiazol-2-thiol	29490-19-5	353
2-(Methylthio)benzothiazole	615-22-5	189
Methylthiouracil (6-Methyl-2-thiouracil)	56-04-2	591
Methyltin Trichloride (Trichloromethylstannane)	993-16-8	591
4-Methylvaleric Acid	646-07-1	290
Metiram	9006-42-2	20
Metobromuron	3060-89-7	77
Metobromuron D6 (methyl D3 methoxy D3)	n/a	78, 545
Metoclopramide Hydrochloride	7232-21-5	353
Metofluthrin	240494-70-6	151
Metolachlor	51218-45-2	78
Metolachlor D6 (propyl D6)	1219803-97-0	78, 545
Metolachlor deschloro	126605-22-9	78
Metolachlor Ethane Sulfonic Acid (ESA) Sodium Salt	947601-85-6	78
Metolachlor-2-ethoxy	68544-97-8	78
Metolachlor-2-hydroxy	131068-72-9	78
Metolachlor mercapturate	159956-64-6	78
Metolachlor oxanilic acid (OA)	152019-73-3	78
S-Metolachlor	87392-12-9	79

Index

Analyte name	CAS No#	Page No#
S-Metolachlor CGA 357704	1217465-10-5	79
S-Metolachlor CGA 368208	1173021-76-5	79
S-Metolachlor CGA 37735	97055-05-5	79
S-Metolachlor CGA 50267	82508-03-0	79
S-Metolachlor CGA 50720	152019-74-4	79
S-Metolachlor NOA 413173	n/a	79
Metolcarb	1129-41-5	151
(E)-Metominostrobin	133408-50-1	20
(Z)-Metominostrobin	133408-51-2	20
Metoprolol Tartrate	56392-17-7	353
Metosulam	139528-85-1	79
Metoxuron	19937-59-8	79
Metoxuron-monomethyl	20782-57-4	80
Metrafenone	220899-03-6	20
Metribuzin	21087-64-9	80
Metribuzin-desamino	35045-02-4	80
Metribuzin-desamino-diketo	52236-30-3	80
Metribuzin-diketo	56507-37-0	80
Metrifonate (Trichlorphon)	52-68-6	151
Metronidazole	443-48-1	353
Metronidazole D4 (ethylene D4)	1261392-47-5	354, 545
Metronidazole-hydroxy D2	n/a	354, 545
Metronidazole-hydroxy (1-(2-Hydroxyethyl)-2-hydroxymethyl-5-nitroimidazole)	4812-40-2	354
Metsulfuron-methyl	74223-64-6	80
Metsulfuron-methyl D3 (triazine methoxy D3)	n/a	80, 545
Mevinphos	7786-34-7	151
E-Mevinphos (cis-butenic acid)	26718-65-0	151
Z-Mevinphos	338-45-4	20, 151
Mexacarbate	315-18-4	151
Milbemectin A3	51596-10-2	151
Milbemectin A4	51596-11-3	151
Mineral Oil	8042-47-5	291
Minocycline Hydrochloride Dihydrate	128420-71-3	354

Analyte name	CAS No#	Page No#
Mirex	2385-85-5	152
Mitotane (2,4'-DDD)	53-19-0	152
Molinate	2212-67-1	80
Molinate-sulfoxide	52236-29-0	81
Mometasone Furoate	83919-23-7	354
Momfluorothrin	609346-29-4	152
Monalide	7287-36-7	81
Monascus Red	874807-57-5	472
Monensin Sodium	22373-78-0	354
Monocrotophos	6923-22-4	152
Monolinuron	1746-81-2	81
Monolinuron D6 (methyl D3 methoxy D3)	n/a	81, 545
Monosultap	29547-00-0	152
Monuron	150-68-5	81
Monuron D6 (dimethyl D6)	217488-65-8	81, 546
Monuron Trichloroacetate (Monuron TCA)	140-41-0	81
Morantel Tartrate Monohydrate	26155-31-7	354
Morpholine	110-91-8	591
Moxidectin	113507-06-5	354
Mupirocin	12650-69-0	354
Muscalure	27519-02-4	189
Musk ambrette (Ambrettolide)	83-66-9	592
Musk ketone	81-14-1	592
Musk moskene (1,1,3,3,5-Pentamethyl-4,6-dinitro-indane)	116-66-5	592
Musk NN (Ethylene Brassylate)	105-95-3	592
Musk tibeten (1-tert-Butyl-3,4,5-trimethyl-2,6-dinitrobenzene)	145-39-1	592
Musk xylene	81-15-2	592
Musk Xylene D15	877119-10-3	546, 592
Myclobutanil	88671-89-0	20
Mycrene (β-Myrcene)	123-35-3	592
Myristoleic Acid	544-64-9	472
Myristoleic acid methyl ester	56219-06-8	472, 592
Nabam	142-59-6	20
Nadifloxacin	124858-35-1	355

Index

Analyte name	CAS No#	Page No#
Nadolol	42200-33-9	355
Nafcillin Sodium Monohydrate	7177-50-6	355
Naftalofos	1491-41-4	152
Naled (Dibrom)	300-76-5	152
Nalidixic Acid	389-08-2	356
Nandrolone (19-Nortestosterone)	434-22-0	356
Nandrolone Decanoate (19-Nortestosterone 17-decanoate)	360-70-3	356
Naphthalene	91-20-3	253
Naphthalene D8	1146-65-2	254, 546
Naphthalene-1-aldehyde	66-77-3	254
Naphthenic acids	1338-24-5	291
Naphthenic acids, copper salts (Copper (II) naphthenate)	1338-02-9	20
1-Naphthol	90-15-3	593
2-Naphthol (β-Naphthol)	135-19-3	593
1,4-Naphthoquinone	130-15-4	593
2-Naphthoxyacetic Acid	120-23-0	189
1-Naphthylacetamide (1-Naphthaleneacetamide)	86-86-2	189
1-Naphthylacetic Acid Methyl Ester	2876-78-0	189
1-Naphthylacetic Acid	86-87-3	189
Naproanilide	52570-16-8	81
Napropamide	15299-99-7	82
Naproxen	22204-53-1	356
Naptalam	132-66-1	82
Naringin	10236-47-2	472
Natamycin (Pimaricin)	7681-93-8	356
Navy Blue 018112	118685-33-9	413
Neburon	555-37-3	82
Neomycin Sulfate	1405-10-3	356
Neostigmine Metilsulfate	51-60-5	356
Nequinat	13997-19-8	356
Nereistoxin oxalate	1631-52-3	152
Nereistoxin Oxalate D6 (dimethyl D6)	n/a	152, 546
Nerol	106-25-2	593
Nervonic Acid Methyl Ester	2733-88-2	473

Analyte name	CAS No#	Page No#
New Red	220658-76-4	413
Nicarbazin	330-95-0	356
Niclosamide	50-65-7	189
Nicosulfuron	111991-09-4	82
Nicotinamide	98-92-0	473
Nicotine	54-11-5	357
Nicotinic Acid	59-67-6	473
Niflumic Acid	4394-00-7	357
Nifuroxazide	965-52-6	357
Nifuroxime (5-Nitrofurfural oxime)	555-15-7	357
Nifurpirinol	13411-16-0	357
Nifursol	16915-70-1	357
Nimorazole	6506-37-2	357
Nitarsone	98-72-6	357
(E)-Nitenpyram	150824-47-8	153
Nitralin	4726-14-1	82
Nitrpyrin	1929-82-4	190
Nitrioltriacetic Acid (NTA)	139-13-9	593
5-Nitroacenaphthene	602-87-9	254
2-Nitroaniline	88-74-4	441
3-Nitroaniline	99-09-2	441
4-Nitroaniline	100-01-6	441
2-Nitrobenzaldehyde	552-89-6	593
Nitrobenzene	98-95-3	442
Nitrobenzene D5	4165-60-0	442, 546
6-Nitro-1H-benzimidazole (5-Nitrobenzimidazole)	94-52-0	593
7-Nitrobenzo[a]anthracene	20268-51-3	254
6-Nitrobenzo[a]pyrene	63041-90-7	254
2-Nitrobenzoic Acid	552-16-9	593
4-Nitrobenzoic Acid	62-23-7	593
2-Nitrobiphenyl	86-00-0	594
4-Nitrobiphenyl	92-93-3	594
2-Nitrodiphenylamine	119-75-5	442
4-Nitrodiphenylamine	836-30-6	442
Nitrofen	1836-75-5	82

Index

Analyte name	CAS No#	Page No#
3-Nitrofluoranthene	892-21-7	254
2-Nitrofluorene	607-57-8	254
Nitrofluorfen	42874-01-1	82
Nitrofural (Nitrofurazone)	59-87-0	357
Nitrofurantoin	67-20-9	358
5-Nitroguaiacol Sodium Salt	67233-85-6	594
Nitroguanidine	556-88-7	442
4-Nitroimidazole	3034-38-6	442
Nitromide	121-81-3	442
1-Nitronaphthalene	86-57-7	254
2-Nitronaphthalene	581-89-5	254
Nitropenta (Pentaerythritol tetranitrate)	78-11-5	442
3-Nitroperylene	20589-63-3	255
3-Nitrophenanthrene	17024-19-0	255
9-Nitrophenanthrene	954-46-1	255
4-Nitrophenol Sodium Salt	824-78-2	190
2-Nitrophenol	88-75-5	442
3-Nitrophenol	554-84-7	443
4-Nitrophenol	100-02-7	190, 443
4-Nitrophenol-2,3,5,6-D4	93951-79-2	190, 546, 443, 546
2-Nitropropane	79-46-9	594
1-Nitropyrene	5522-43-0	255
4-Nitropyrene	57835-92-4	255
4-Nitroquinoline N-Oxide	56-57-5	255
(S)-N-Nitrosoanabasine D4 (pyridinyl D4)	1426174-40-4	546, 594
(RS)-N-Nitrosoanabasine	37620-20-5	594
(S)-N-Nitrosoanatabine D4 (pyridinyl D4)	1426174-82-4	547, 594
(RS)-N-Nitrosoanatabine	887407-16-1	595
N-Nitroso-di-n-butylamine	924-16-3	595
N-Nitrosodiethanolamine	1116-54-7	595
N-Nitroso-diethylamine	55-18-5	595
N-Nitrosodiisopropylamine	601-77-4	595
N-Nitroso-dimethylamine	62-75-9	595
N-Nitroso-dimethylamine D6	17829-05-9	547, 595
N-Nitroso-diphenylamine	86-30-6	594

Analyte name	CAS No#	Page No#
N-Nitroso-di-n-propylamine	621-64-7	595
N-Nitroso-di-n-propylamine D14	93951-96-3	547, 595
N-Nitroso-methyl-ethylamine	10595-95-6	594
4-Nitrosomorpholine (N-Nitrosomorpholine)	59-89-2	596
N-Nitrosornicotine ketone	64091-91-4	596
N-Nitrosornicotine-ketone D4 (pyridyl D4)	764661-24-7	547, 596
N-Nitrosornicotine	80508-23-2	596
N-Nitrosopiperidine	100-75-4	596
N-Nitrosopyrrolidine	930-55-2	596
Nitrothal-isopropyl	10552-74-6	21
2-Nitrotoluene	88-72-2	443
2-Nitrotoluene D7	84344-04-7	443, 547
3-Nitrotoluene	99-08-1	443
4-Nitrotoluene	99-99-0	443
4-Nitrotoluene-2-sulfonic Acid	121-03-9	443
Nitrovin Hydrochloride	2315-20-0	358
Nitroxinil (Nitroxynil)	1689-89-0	358
cis-Nonachlor	5103-73-1	153
trans-Nonachlor	39765-80-5	153
1,2,3,4,5,6,7,8,9-Nonachlorodecane (CP-10)	890302-90-6	596
n-Nonacosane	630-03-5	291
n-Nonadecane	629-92-5	291
Nonadecanoic Acid (Nonadecyclic acid)	646-30-0	473
Nonadecanoic Acid Ethyl Ester	18281-04-4	473
Nonadecanoic Acid Methyl Ester	1731-94-8	473
(2E,6Z)-Nona-2,6-dienal	557-48-2	596
Nonanal-2,4-dinitrophenylhydrazone	2348-19-8	291
n-Nonane	111-84-2	291
Nonanoic Acid (Pelargonic acid)	112-05-0	82
Nonanoic Acid Ethyl Ester	123-29-5	473
Nonanoic Acid Methyl Ester	1731-84-6	473
1-Nonanol	143-08-8	291
2-Nonanone	821-55-6	291
n-Nonatriacontane	7194-86-7	291

Index

Analyte name	CAS No#	Page No#
Nonivamide (N-Vanillylnonanamide)	2444-46-4	473
n-Nonylbenzene (Phenylnonane)	1081-77-2	292
Nonylphenol (technical)	25154-52-3	596
4-Nonylphenol Diethoxylate	20427-84-3	597
4-Nonylphenol-ethoxylate (technical)	68412-54-4	597
4-Nonylphenol Monoethoxylate	104-35-8	597
4-Nonylphenol-mono-ethoxylate D2	n/a	547, 597
4-Nonylphenol Triethoxylate	51437-95-7	597
4-n-Nonylphenol	104-40-5	597
4-n-Nonylphenol D8 (ring D4-ethylD4)	n/a	547, 597
4-Nonylphenoxyacetic Acid	3115-49-9	597
Norethisterone Acetate	51-98-9	358
Norfloxacin	70458-96-7	358
Norfloxacin D5 (ethyl D5)	1015856-57-1	358, 547
Norflurane (1,1,1,2-Tetrafluoroethane)	811-97-2	391
Norflurazon	27314-13-2	82
Norflurazon-desmethyl	23576-24-1	83
19-Nortestosterone 17-Propionate	7207-92-3	358
Noruron	18530-56-8	83
Novaluron	116714-46-6	153
Noviflumuron	121451-02-3	153
Novobiocin Sodium	1476-53-5	358
2-NP-AHD	623145-57-3	355
2-NP-AMOZ	183193-59-1	355
2-NP-AMOZ D5	1173097-59-0	355, 546
2-NP-AOZ	19687-73-1	355
2-NP-AOZ D4	1007478-57-0	355, 546
2-NP-SCA	16004-43-6	355
Nuarimol	63284-71-9	21
Nystatin	1400-61-9	358
Octabromodiphenylether	32536-52-0	598
2,3,4,5,6,7,8,9-Octachlorodecane (CP-9)	214327-66-9	598

Analyte name	CAS No#	Page No#
Octachloronaphthalene	2234-13-1	255
Octachlorostyrene	29082-74-4	598
n-Octacosane	630-02-4	292
Octacosanoic Acid (Montanic acid)	506-48-9	474
1-Octacosanol (Montanyl alcohol)	557-61-9	292
n-Octadecane	593-45-3	292
Octadecanedioic Acid	871-70-5	474
Octadecanoic Acid Ethyl Ester	111-61-5	474
trans-11-Octadecanoic Acid Methyl Ester	6198-58-9	474
trans-11-Octadecanoic Acid (Vaccenic acid)	693-72-1	474
1-Octadecanol (Stearyl alcohol)	112-92-5	292
1-Octadecene	112-88-9	292
cis-11-Octadecenoic acid methyl ester	1937-63-9	474
Octafluoronaphthalene	313-72-4	255
Octanal	124-13-0	292
Octanal-2,4-dinitrophenylhydrazone	1726-77-8	292
n-Octane	111-65-9	292
1,8-Octanedicarboxylic Acid Dimethyl Ester	106-79-6	474
Octanoic Acid Ethyl Ester	106-32-1	474
Octanoic Acid Methyl Ester	111-11-5	474
1-Octanol	111-87-5	293
2-Octanol	123-96-6	293
2-Octanone	111-13-7	293
3-Octanone	106-68-3	293
n-Octatriacontane	7194-85-6	293
1-Octene	111-66-0	293
trans-2-Octene	13389-42-9	293
Octhilnone	26530-20-1	21
Octisalate (Salicylic acid 2-ethyl-1-hexyl ester)	118-60-5	359
Octogen	2691-41-0	443
Octopamine Hydrochloride	770-05-8	359
1-Octylamine	111-86-4	293
Octylbenzene	2189-60-8	293

Index

Analyte name	CAS No#	Page No#
4-n-Octylphenol-di-ethoxylate	51437-90-2	598
4-Octylphenol-ethoxylate	9036-19-5	598
4-Octylphenol Polyethoxylate	26636-32-8	598
4-Octylphenol	1806-26-4	598
4-tert-Octylphenol	140-66-9	598
n-Octyltin Trichloride	3091-25-6	598
2-Octynoic Acid Methyl Ester	111-12-6	599
Ofloxacin	82419-36-1	359
Ofloxacin D3 (N-methyl D3)	1173147-91-5	359, 547
Ofloxacin D3 Hydrochloride (N-methyl D3)	1173021-78-7	359, 548
Ofurace	58810-48-3	21
Olaquinox	23696-28-8	359
Oleandomycin phosphate, dihydrate	7060-74-4	359
Oleic acid	112-80-1	599
Oleic Acid Methyl Ester	112-62-9	83
Omeprazole	73590-58-6	359
Omethoate	1113-02-6	83
Omethoate D6 (O-dimethyl D6)	1219804-92-8	83, 548
Orange 2 Sodium Salt	633-96-5	413
Orbencarb	34622-58-7	83
Orbifloxacin	113617-63-3	359
Orciprenaline Sulfate (Metaproterenol hemisulfate salt)	5874-97-5	360
Ormetoprim	6981-18-6	360
Ornidazole	16773-42-5	360
Orotic Acid	65-86-1	475
Orotic acid monohydrate	50887-69-9	475
Orthosulfamuron	213464-77-8	83
Orysastrobin	248593-16-0	21
Oryzalin	19044-88-3	83
Oseltamivir	196618-13-0	360
Oxabetrinil	74782-23-3	190
Oxacillin Sodium Monohydrate	7240-38-2	360
Oxacycloheptadec-8-en-2-one	123-69-3	599
Oxacyclohexadecan-2-one	106-02-5	599
Oxacyclohexadecen-2-one (Globalide)	34902-57-3	599

Analyte name	CAS No#	Page No#
Oxadiargyl	39807-15-3	83
Oxadiazon	19666-30-9	84
Oxadiazon-hydroxy (5-tert-Butyl-3-(2,4-dichloro-5-hydroxyphenyl)-1,3,4-oxadiazol-2-one)	39807-19-7	84
Oxadixyl	77732-09-3	21
Oxalic Acid (Ethanedioic acid)	144-62-7	475
Oxamyl	23135-22-0	153
Oxamyl D3 (2-methyl-D3)	n/a	153, 548
Oxamyl-oxime	30558-43-1	153
Oxandrolone	53-39-4	360
Oxasulfuron	144651-06-9	84
Oxaziclomefone	153197-14-9	84
Oxcarbazepine	28721-07-5	360
Oxfendazole	53716-50-0	360
Oxibendazole	20559-55-1	360
Oxine Copper	10380-28-6	21
Oxolinic Acid	14698-29-4	21
Oxpoconazole fumarate	174212-12-5	21
Oxybenzone (2-Hydroxy-4-methoxybenzophenone; Neo Heliopan BB)	131-57-7	599
Oxycarboxin	5259-88-1	22
Oxyclozanide	2277-92-1	361
Oxyfluorfen	42874-03-3	84
Oxyquinoline (8-Hydroxyquinoline)	148-24-3	22
Oxyquinoline Sulfate (8-Hydroxyquinoline sulfate)	134-31-6	22
Oxytetracycline Dihydrate	6153-64-6	361
Oxytetracycline Hydrochloride	2058-46-0	22
Paclobutrazol	76738-62-0	190
Padimate O	21245-02-3	361
Palmitelaidic acid methyl ester	10030-74-7	475
Palmitic Acid (Hexadecanoic acid)	57-10-3	475
Palmitic Acid n-Butyl Ester	111-06-8	475
Palmitoleic Acid	373-49-9	475
DL-Pantothenic Acid Calcium Salt	6381-63-1	475
D-Pantothenic Acid Sodium Salt	867-81-2	475

Index

Analyte name	CAS No#	Page No#
Para Red	6410-10-2	413
Para Red D6	1014689-16-7	413, 548
Paracetamol	103-90-2	361
Paraformaldehyde	30525-89-4	599
Paraldehyde	123-63-7	361
Paraoxon (Paraoxon-ethyl)	311-45-5	153
Paraoxon-methyl	950-35-6	154
Paraquat Dichloride	1910-42-5	84
Paraquat dichloride D6 (dimethyl)	n/a	84, 548
Paraquat Dichloride D8	347841-45-6	84, 548
Paraquat diiodide D6 (dimethyl)	n/a	84, 548
Parathion-ethyl	56-38-2	154
Parathion-ethyl D10 (diethyl D10)	350820-04-1	154, 548
Parathion-methyl	298-00-0	154
Parathion-methyl D6 (dimethyl D6)	96740-32-8	154, 548
Parbendazole	14255-87-9	361
Paromomycin Sulfate	1263-89-4	294
Parsol 1789 (Avobenzone; 1-[4-(1,1-Dimethylethyl)phenyl]-3-(4-methoxyphenyl)-1,3-propandione)	70356-09-1	600
Parsol 340 (Octocrylene; 2-Ethylhexyl 2-cyano-3,3-diphenylacrylate)	6197-30-4	600
Parsol HS (Ensulizole; 2-Phenyl-1H-benzimidazol-5-sulfonic acid)	27503-81-7	600
Parsol MCX (Octinoxate; 2-Ethylhexyl p-methoxycinnamate)	5466-77-3	600
Parsol SLX	207574-74-1	600
Patent Blue 5 (E131)	20262-76-4	414
PBB 1 (2-Bromobiphenyl)	2052-07-5	309
PBB 2 (3-Bromobiphenyl)	2113-57-7	309
PBB 3 (4-Bromobiphenyl)	92-66-0	309
PBB 4 (2,2'-Dibromobiphenyl)	13029-09-9	309
PBB 9 (2,5-Dibromobiphenyl)	57422-77-2	309
PBB 15 (4,4'-Dibromobiphenyl)	92-86-4	309
PBB 18 (2,2',5-Tribromobiphenyl)	59080-34-1	309
PBB 26 (2,3',5-Tribromobiphenyl)	59080-35-2	309

Analyte name	CAS No#	Page No#
PBB 29 (2,4,5-Tribromobiphenyl)	115245-07-3	309
PBB 30 (2,4,6-Tribromobiphenyl)	59080-33-0	310
PBB 31 (2,4',5-Tribromobiphenyl)	59080-36-3	310
PBB 49 (2,2',4,5'-Tetrabromobiphenyl)	60044-24-8	310
PBB 52 (2,2',5,5'-Tetrabromobiphenyl)	59080-37-4	310
PBB 53 (2,2',5,6'-Tetrabromobiphenyl)	60044-25-9	310
PBB 80 (3,3',5,5'-Tetrabromobiphenyl)	16400-50-3	310
PBB 101 (2,2',4,5,5'-Pentabromobiphenyl)	67888-96-4	310
PBB 103 (2,2',4,5',6-Pentabromobiphenyl)	59080-39-6	310
PBB 153 (2,2',4,4',5,5'-Hexabromobiphenyl)	59080-40-9	310
PBB 155 (2,2',4,4',6,6'-Hexabromobiphenyl)	59261-08-4	311
PBB 209 (Decabromobiphenyl)	13654-09-6	311
PBDE 15 (4,4'-Dibromodiphenyl Ether)	2050-47-7	599
PBDE 47 (2,2',4,4'-Tetrabromodiphenyl Ether)	5436-43-1	599
PCB 1 (2-Chlorobiphenyl)	2051-60-7	217
PCB 2 (3-Chlorobiphenyl)	2051-61-8	217
PCB 3 (4-Chlorobiphenyl)	2051-62-9	217
PCB 4 (2,2'-Dichlorobiphenyl)	13029-08-8	217
PCB 5 (2,3-Dichlorobiphenyl)	16605-91-7	217
PCB 7 (2,4-Dichlorobiphenyl)	33284-50-3	217
PCB 8 (2,4'-Dichlorobiphenyl)	34883-43-7	217
PCB 9 (2,5-Dichlorobiphenyl)	34883-39-1	217
PCB 10 (2,6-Dichlorobiphenyl)	33146-45-1	217
PCB 11 (3,3'-Dichlorobiphenyl)	2050-67-1	218
PCB 13 (3,4'-Dichlorobiphenyl)	2974-90-5	218
PCB 14 (3,5-Dichlorobiphenyl)	34883-41-5	218
PCB 15 (4,4'-Dichlorobiphenyl)	2050-68-2	218
PCB 16 (2,2',3-Trichlorobiphenyl)	38444-78-9	218
PCB 18 (2,2',5-Trichlorobiphenyl)	37680-65-2	218
PCB 19 (2,2',6-Trichlorobiphenyl)	38444-73-4	218
PCB 20 (2,3,3'-Trichlorobiphenyl)	38444-84-7	218

Index

Analyte name	CAS No#	Page No#
PCB 21 (2,3,4-Trichlorobiphenyl)	55702-46-0	218
PCB 22 (2,3,4'-Trichlorobiphenyl)	38444-85-8	219
PCB 23 (2,3,5-Trichlorobiphenyl)	55720-44-0	219
PCB 26 (2,3',5-Trichlorobiphenyl)	38444-81-4	219
PCB 28 (2,4,4'-Trichlorobiphenyl)	7012-37-5	219
PCB 29 (2,4,5-Trichlorobiphenyl)	15862-07-4	219
PCB 30 (2,4,6-Trichlorobiphenyl)	35693-92-6	219
PCB 31 (2,4',5-Trichlorobiphenyl)	16606-02-3	219
PCB 33 (2,3',4'-Trichlorobiphenyl)	38444-86-9	219
PCB 34 (2,3',5'-Trichlorobiphenyl)	37680-68-5	219
PCB 35 (3,3',4'-Trichlorobiphenyl)	37680-69-6	220
PCB 37F (3'-Fluoro-3,4,4'-trichlorobiphenyl)	1191034-39-5	220
PCB 40 (2,2',3,3'-Tetrachlorobiphenyl)	38444-93-8	220
PCB 43 (2,2',3,5-Tetrachlorobiphenyl)	70362-46-8	220
PCB 44 (2,2',3,5'-Tetrachlorobiphenyl)	41464-39-5	220
PCB 46 (2,2',3,6'-Tetrachlorobiphenyl)	41464-47-5	220
PCB 47 (2,2',4,4'-Tetrachlorobiphenyl)	2437-79-8	220
PCB 48 (2,2',4,5-Tetrachlorobiphenyl)	70362-47-9	220
PCB 49 (2,2',4,5'-Tetrachlorobiphenyl)	41464-40-8	220
PCB 50 (2,2',4,6-Tetrachlorobiphenyl)	62796-65-0	221
PCB 51 (2,2',4,6'-Tetrachlorobiphenyl)	68194-04-7	221
PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)	35693-99-3	221
PCB 53 (2,2',5,6'-Tetrachlorobiphenyl)	41464-41-9	221
PCB 54 (2,2',6,6'-Tetrachlorobiphenyl)	15968-05-5	221
PCB 55 (2,3,3',4-Tetrachlorobiphenyl)	74338-24-2	221
PCB 56 (2,3,3',4'-Tetrachlorobiphenyl)	41464-43-1	221
PCB 61 (2,3,4,5-Tetrachlorobiphenyl)	33284-53-6	221
PCB 62 (2,3,4,6-Tetrachlorobiphenyl)	54230-22-7	221

Analyte name	CAS No#	Page No#
PCB 65 (2,3,5,6-Tetrachlorobiphenyl)	33284-54-7	222
PCB 66 (2,3',4,4'-Tetrachlorobiphenyl)	32598-10-0	222
PCB 68 (2,3',4,5'-Tetrachlorobiphenyl)	73575-52-7	222
PCB 69 (2,3',4,6-Tetrachlorobiphenyl)	60233-24-1	222
PCB 70 (2,3',4',5-Tetrachlorobiphenyl)	32598-11-1	222
PCB 73 (2,3',5',6-Tetrachlorobiphenyl)	74338-23-1	222
PCB 74 (2,4,4',5-Tetrachlorobiphenyl)	32690-93-0	222
PCB 75 (2,4,4',6-Tetrachlorobiphenyl)	32598-12-2	222
PCB 77 (3,3',4,4'-Tetrachlorobiphenyl)	32598-13-3	222
PCB 78 (3,3',4,5-Tetrachlorobiphenyl)	70362-49-1	223
PCB 80 (3,3',5,5'-Tetrachlorobiphenyl)	33284-52-5	223
PCB 81 (3,4,4',5-Tetrachlorobiphenyl)	70362-50-4	223
PCB 82 (2,2',3,3',4-Pentachlorobiphenyl)	52663-62-4	223
PCB 85 (2,2',3,4,4'-Pentachlorobiphenyl)	65510-45-4	223
PCB 87 (2,2',3,4,5'-Pentachlorobiphenyl)	38380-02-8	223
PCB 89 (2,2',3,4,6'-Pentachlorobiphenyl)	73575-57-2	223
PCB 90 (2,2',3,4',5-Pentachlorobiphenyl)	68194-07-0	223
PCB 91 (2,2',3,4',6-Pentachlorobiphenyl)	68194-05-8	223
PCB 95 (2,2',3,5',6-Pentachlorobiphenyl)	38379-99-6	224
PCB 97 (2,2',3,4',5'-Pentachlorobiphenyl)	41464-51-1	224
PCB 99 (2,2',4,4',5-Pentachlorobiphenyl)	38380-01-7	224
PCB 100 (2,2',4,4',6-Pentachlorobiphenyl)	39485-83-1	224
PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)	37680-73-2	224
PCB 103 (2,2',4,5',6-Pentachlorobiphenyl)	60145-21-3	224

Index

Analyte name	CAS No#	Page No#
PCB 104 (2,2',4,6,6'-Pentachlorobiphenyl)	56558-16-8	224
PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl)	32598-14-4	224
PCB 106 (2,3,3',4,5-Pentachlorobiphenyl)	70424-69-0	224
PCB 107 (2,3,3',4',5-Pentachlorobiphenyl)	70424-68-9	225
PCB 108 (2,3,3',4,5'-Pentachlorobiphenyl)	70362-41-3	225
PCB 110 (2,3,3',4',6-Pentachlorobiphenyl)	38380-03-9	225
PCB 112 (2,3,3',5,6-Pentachlorobiphenyl)	74472-36-9	225
PCB 113 (2,3,3',5',6-Pentachlorobiphenyl)	68194-10-5	225
PCB 114 (2,3,4,4',5-Pentachlorobiphenyl)	74472-37-0	225
PCB 116 (2,3,4,5,6-Pentachlorobiphenyl)	18259-05-7	225
PCB 118 (2,3',4,4',5-Pentachlorobiphenyl)	31508-00-6	225
PCB 119 (2,3',4,4',6-Pentachlorobiphenyl)	56558-17-9	225
PCB 121 (2,3',4,5',6-Pentachlorobiphenyl)	56558-18-0	226
PCB 123 (2,3',4,4',5'-Pentachlorobiphenyl)	65510-44-3	226
PCB 125 (2,3',4',5',6-Pentachlorobiphenyl)	74472-39-2	226
PCB 126 (3,3',4,4',5-Pentachlorobiphenyl)	57465-28-8	226
PCB 128 (2,2',3,3',4,4'-Hexachlorobiphenyl)	38380-07-3	226
PCB 129 (2,2',3,3',4,5-Hexachlorobiphenyl)	55215-18-4	226
PCB 131 (2,2',3,3',4,6-Hexachlorobiphenyl)	61798-70-7	226
PCB 132 (2,2',3,3',4,6'-Hexachlorobiphenyl)	38380-05-1	226
PCB 136 (2,2',3,3',6,6'-Hexachlorobiphenyl)	38411-22-2	226
PCB 137 (2,2',3,4,4',5-Hexachlorobiphenyl)	35694-06-5	227
PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)	35065-28-2	227
PCB 141 (2,2',3,4,5,5'-Hexachlorobiphenyl)	52712-04-6	227

Analyte name	CAS No#	Page No#
PCB 142 (2,2',3,4,5,6-Hexachlorobiphenyl)	41411-61-4	227
PCB 143 (2,2',3,4,5,6'-Hexachlorobiphenyl)	68194-15-0	227
PCB 145 (2,2',3,4,6,6'-Hexachlorobiphenyl)	74472-40-5	227
PCB 146 (2,2',3,4',5,5'-Hexachlorobiphenyl)	51908-16-8	227
PCB 147 (2,2',3,4',5,6-Hexachlorobiphenyl)	68194-13-8	227
PCB 149 (2,2',3,4',5',6-Hexachlorobiphenyl)	38380-04-0	227
PCB 150 (2,2',3,4',6,6'-Hexachlorobiphenyl)	68194-08-1	228
PCB 151 (2,2',3,5,5',6-Hexachlorobiphenyl)	52663-63-5	228
PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)	35065-27-1	228
PCB 155 (2,2',4,4',6,6'-Hexachlorobiphenyl)	33979-03-2	228
PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl)	38380-08-4	228
PCB 157 (2,3,3',4,4',5'-Hexachlorobiphenyl)	69782-90-7	228
PCB 159 (2,3,3',4,5,5'-Hexachlorobiphenyl)	39635-35-3	228
PCB 160 (2,3,3',4,5,6-Hexachlorobiphenyl)	41411-62-5	228
PCB 161 (2,3,3',4,5',6-Hexachlorobiphenyl)	74472-43-8	228
PCB 163 (2,3,3',4',5,6-Hexachlorobiphenyl)	74472-44-9	229
PCB 164 (2,3,3',4',5',6-Hexachlorobiphenyl)	74472-45-0	229
PCB 166 (2,3,4,4',5,6-Hexachlorobiphenyl)	41411-63-6	229
PCB 167 (2,3',4,4',5,5'-Hexachlorobiphenyl)	52663-72-6	229
PCB 168 (2,3',4,4',5',6-Hexachlorobiphenyl)	59291-65-5	229
PCB 169 (3,3',4,4',5,5'-Hexachlorobiphenyl)	32774-16-6	229
PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl)	35065-30-6	229
PCB 171 (2,2',3,3',4,4',6-Heptachlorobiphenyl)	52663-71-5	229
PCB 172 (2,2',3,3',4,5,5'-Heptachlorobiphenyl)	52663-74-8	229

Index

Analyte name	CAS No#	Page No#
PCB 173 (2,2',3,3',4,5,6-Heptachlorobiphenyl)	68194-16-1	230
PCB 174 (2,2',3,3',4,5,6'-Heptachlorobiphenyl)	38411-25-5	230
PCB 175 (2,2',3,3',4,5',6-Heptachlorobiphenyl)	40186-70-7	230
PCB 177 (2,2',3,3',4,5',6'-Heptachlorobiphenyl)	52663-70-4	230
PCB 178 (2,2',3,3',5,5',6-Heptachlorobiphenyl)	52663-67-9	230
PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)	35065-29-3	230
PCB 181 (2,2',3,4,4',5,6-Heptachlorobiphenyl)	74472-47-2	230
PCB 183 (2,2',3,4,4',5',6-Heptachlorobiphenyl)	52663-69-1	230
PCB 185 (2,2',3,4,5,5',6-Heptachlorobiphenyl)	52712-05-7	230
PCB 187 (2,2'',3,4'',5,5'',6-Heptachlorobiphenyl)	52663-68-0	231
PCB 189 (2,3,3',4,4',5,5'-Heptachlorobiphenyl)	39635-31-9	231
PCB 190 (2,3,3',4,4',5,6-Heptachlorobiphenyl)	41411-64-7	231
PCB 192 (2,3,3',4,5,5',6-Heptachlorobiphenyl)	74472-51-8	231
PCB 193 (2,3,3',4',5,5',6-Heptachlorobiphenyl)	69782-91-8	231
PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)	35694-08-7	231
PCB 195 (2,2',3,3',4,4',5,6-Octachlorobiphenyl)	52663-78-2	231
PCB 196 (2,2',3,3',4,4',5,6'-Octachlorobiphenyl)	42740-50-1	231
PCB 198 (2,2',3,3',4,5,5',6-Octachlorobiphenyl)	68194-17-2	231
PCB 199 (2,2',3,3',4,5,5',6'-Octachlorobiphenyl)	52663-75-9	232
PCB 200 (2,2',3,3',4,5,6,6'-Octachlorobiphenyl)	52663-73-7	232
PCB 202 (2,2',3,3',5,5',6,6'-Octachlorobiphenyl)	2136-99-4	232
PCB 203 (2,2',3,4,4',5,5',6-Octachlorobiphenyl)	52663-76-0	232
PCB 204 (2,2',3,4,4',5,6,6'-Octachlorobiphenyl)	74472-52-9	232

Analyte name	CAS No#	Page No#
PCB 205 (2,3,3',4,4',5,5',6-Octachlorobiphenyl)	74472-53-0	232
PCB 206 (2,2',3,3',4,4',5,5',6-Nonachlorobiphenyl)	40186-72-9	232
PCB 207 (2,2',3,3',4,4',5,6,6'-Nonachlorobiphenyl)	52663-79-3	232
PCB 208 (2,2',3,3',4,5,5',6,6'-Nonachlorobiphenyl)	52663-77-1	232
PCB 209 (Decachlorobiphenyl)	2051-24-3	233
Pebulate	1114-71-2	85
Pefloxacin Mesilate Dihydrate (Pefloxacin methanesulfonate dihydrate)	149676-40-4	361
Penbutolol ((S)-Penbutolol)	38363-40-5	361
Penconazole	66246-88-6	22
Pencycuron	66063-05-6	22
Pendimethalin	40487-42-1	85
Pendimethalin D5 (pent-3-yl (2,2,3,4,4)-D5)	1219803-39-0	85, 549
Penfluron	35367-31-8	154
Penicillin G Benzathine Tetrahydrate	41372-02-5	361
Penoxsulam	219714-96-2	85
Pentabromodiphenylether	32534-81-9	600
2,3,4,5,6-Pentabromoethylbenzene	85-22-3	444
Pentabromophenol	608-71-9	444
Pentachloroaniline	527-20-8	154, 190, 444
Pentachloroanisole	1825-21-4	444
Pentachloroanisole 13C6 (ring 13C)	n/a	444, 549
Pentachlorobenzene	608-93-5	154, 190
1,2,5,6,9-Pentachlorodecane (CP-3)	175801-37-3	600
Pentachloroethane	76-01-7	391
Pentachlorophenol	87-86-5	155
Pentachlorophenol 13C6	85380-74-1	155, 549
Pentachlorophenol Acetate	1441-02-7	155
Pentachlorophenol Sodium Salt Hydrate	123333-54-0	155
3,3',3'',4,4''-Pentachloro-m-terphenyl	1064187-31-0	233

Index

Analyte name	CAS No#	Page No#
3,3',3'',5,5''-Pentachloro-m-terphenyl	n/a	233
Pentachlorothiophenol	133-49-3	444
2,3,4,5,6-Pentachlorotoluene	877-11-2	444
n-Pentacosane	629-99-2	294
Pentacosanoic Acid	506-38-7	476
1-Pentacosanol	26040-98-2	476
n-Pentadecane	629-62-9	294
Pentadecanoic Acid (Pentadecylic acid)	1002-84-2	476
Pentadecanoic Acid Ethyl Ester	41114-00-5	476
Pentadecanoic Acid Methyl Ester	7132-64-1	476
1-Pentadecanol	629-76-5	294
2-Pentadecanone	2345-28-0	294
Pentafluorophenol	771-61-9	444
Pentamethylbenzene	700-12-9	294
2,2,4,6,6-Pentamethylheptane	13475-82-6	294
Pentamidine Diisetonate	140-64-7	362
Pentanal-2,4-dinitrophenylhydrazone	2057-84-3	294
n-Pentane	109-66-0	294
1,5-Pentanedicarboxylic Acid	111-16-0	476
1,5-Pentanediol	111-29-5	295
Pentanochlor	2307-68-8	85
1-Pentanol	71-41-0	295
2-Pentanol	6032-29-7	295
3-Pentanol	584-02-1	295
2-Pentanone	107-87-9	295
3-Pentanone	96-22-0	295
n-Pentatriacontane	630-07-9	295
Penthiopyrad	183675-82-3	22
Pentoxazone	110956-75-7	85
Pentoxifylline	6493-05-6	362
n-Pentylbenzene (n-Amylbenzene)	538-68-1	295
sec-Pentylbenzene (sec-Amylbenzene)	2719-52-0	295
4-Pentylphenol (4-n-Amylphenol)	14938-35-3	444
Perfluidone	37924-13-3	85

Analyte name	CAS No#	Page No#
Perfluorobutanesulfonic Acid	375-73-5	600
Perfluorobutanoic Acid	375-22-4	600
Perfluorododecanoic Acid	307-55-1	601
Perfluoroheptanesulfonic Acid	375-92-8	601
Perfluoroheptanoic Acid	375-85-9	601
Perfluorohexanesulfonic Acid	355-46-4	601
Perfluorohexanoic Acid	307-24-4	601
Perfluorononanoic Acid	375-95-1	601
Perfluorooctane	307-34-6	601
Perfluorooctane sulfonamide PFOSA	754-91-6	601
Perfluorooctane sulfonic acid	1763-23-1	602
Perfluorooctane Sulfonic Acid Potassium Salt (PFOS)	2795-39-3	601
Perfluorooctane-1-sulfonyl Fluoride	307-35-7	602
Perfluorooctanoic Acid	335-67-1	602
Perfluorooctanoic Acid Ammonium Salt (PFOA; POAA)	3825-26-1	602
Perfluoropentanoic acid	2706-90-3	602
Perfluorotributylamine (PFTBA)	311-89-7	602
Perfluorotridecanoic Acid	72629-94-8	602
Perfluoroundecanoic Acid	2058-94-8	602
Pergolide Mesilate	66104-23-2	362
Permethrin	52645-53-1	155
cis-Permethrin	61949-76-6	155
trans-Permethrin	61949-77-7	155
trans-Permethrin D6 (dimethyl D6)	82523-59-9	155, 549
Permethrinic Acid (3-(2,2-Dichlorovinyl)-2,2-dimethylcyclopropanecarboxylic Acid)	55701-05-8	156
Permethrinic Acid Methyl Ester (Methyl 3-(2,2-Dichlorovinyl)-2,2-dimethylcyclopropanecarboxy late)	61898-95-1	156
Perthane	72-56-0	156
Perylene	198-55-0	255
Perylene-d12	1520-96-3	256, 549
Pethoxamid	106700-29-2	85
Phenacetin	62-44-2	362

Index

Analyte name	CAS No#	Page No#
Phenanthrene	85-01-8	256
Phenanthrene D10	1517-22-2	256, 549
Phenazine	92-82-0	256
Phenazone	60-80-0	362
Phenisopham	57375-63-0	85
Phenkapton	2275-14-1	156
Phenmedipham	13684-63-4	86
Phenmedipham D3 (methoxy D3)	1773497-41-8	86, 549
Phenol	108-95-2	445
Phenol D5 (2,3,4,5,6-Pentadeuteriophenol)	4165-62-2	445, 549
Phenol D6	13127-88-3	445, 549
Phenothiazine	92-84-2	362
Phenothrin	26002-80-2	156
D-trans-Phenothrin	26046-85-5	156
Phenoxyacetic Acid	122-59-8	86
2-Phenoxybenzoic Acid	2243-42-7	156
3-Phenoxybenzoic Acid	3739-38-6	156
3-Phenoxybenzyl Alcohol	13826-35-2	156
4-Phenoxybutyric Acid	6303-58-8	602
Phenoxyethanol (Ethyleneglycol monophenyl ether)	122-99-6	296
Phenoxymethylpenicillin Potassium (Penicillin V potassium salt)	132-98-9	362
1-Phenoxy-2-propanol	770-35-4	296
2-Phenoxypropionic Acid	940-31-8	603
Phenthoate	2597-03-7	157
Phenylacetoneitrile (Benzylcyanide)	140-29-4	362
DL-Phenylalanine	150-30-1	603
9-Phenylanthracene	602-55-1	256
Phenylbutazone	50-33-9	362
4-Phenyl-1-cyclohexene	4994-16-5	296
1,3-Phenylenediamine	108-45-2	445
1,4-Phenylenediamine (1,4-Diaminobenzene)	106-50-3	445
2-Phenylethylamine	64-04-0	296
Phenylisocyanate	103-71-9	603

Analyte name	CAS No#	Page No#
Phenylmercury Chloride	100-56-1	22
1-Phenylnaphthalene	605-02-7	256
2-Phenylnaphthalene	612-94-2	256
2-Phenylphenol	90-43-7	22
3-Phenylphenol	580-51-8	445
4-Phenylphenol	92-69-3	445
2-Phenyl-2-propanol	617-94-7	296
N-Phenylthiourea	103-85-5	603
Phloxine B	18472-87-2	414
Phorate	298-02-2	157
Phorate (diethyl-D10)	1219805-45-4	157, 550
Phorate-oxon	2600-69-3	157
Phorate-oxon-sulfone	2588-06-9	157
Phorate-oxon-sulfoxide	2588-05-8	157
Phorate-sulfone	2588-04-7	157
Phorate-sulfoxide	2588-03-6	157
Phosalone	2310-17-0	157
Phosalone D10 (di-ethyl D5)	n/a	158, 550
Phosfolan	947-02-4	158
Phosfolan-methyl	5120-23-0	158
Phosmet	732-11-6	158
Phosmet D6 (dimethoxy D3)	2083623-41-8	158, 550
Phosmet-oxon	3735-33-9	158
Phosphamidon	13171-21-6	158
Phosphonic Acid	13598-36-2	23
N-(Phosphonomethyl) iminodiacetic acid monohydrate	5994-61-6	603
Photoheptachlor	33442-83-0	158
Phoxim	14816-18-3	158
Phthalic Acid Benzyl Butyl Ester	85-68-7	485
Phthalic Acid Benzyl Butyl Ester (3,4,5,6)-D4	93951-88-3	485, 550
Phthalic Acid Bis-dodecyl Ester	2432-90-8	485
Phthalic Acid Bis-hexyl Ester (Mixture of Isomers)	68515-50-4	485
Phthalic Acid Bis-methylglycol Ester	117-82-8	486

Index

Analyte name	CAS No#	Page No#
Phthalic Acid Bis(2-butoxyethyl) Ester	117-83-9	486
Phthalic acid, bis-C6-C10-alkyl ester	68515-51-5	486
Phthalic acid, bis-C6-C8-branched alkyl esters C7-rich	71888-89-6	486
Phthalic Acid Bis(2-ethoxyethyl) Ester	605-54-9	486
Phthalic Acid Bis(2-ethylhexyl) Ester (3,4,5,6)-D4	93951-87-2	486, 550
Phthalic acid, bis-iso-decyl ester	26761-40-0	486
Phthalic acid, bis-iso-heptyl ester	41451-28-9	486
Phthalic Acid, Bis-isononyl Ester	28553-12-0	486
Phthalic Acid Bis(isopentyl) Ester	605-50-5	487
Phthalic Acid Bis(isopentyl) Ester (3,4,5,6)-D4	1346597-80-5	487, 550
Phthalic acid, bis-2-methyloctyl ester	70857-56-6	487
Phthalic Acid Bis(4-methyl-2-pentyl) Ester	84-63-9	487
Phthalic Acid Bis(2-propylheptyl) Ester	53306-54-0	487
Phthalic Acid Bis(undecyl) Ester	3648-20-2	487
Phthalic Acid Butyl 2-Ethylhexyl Ester	85-69-8	487
Phthalic Acid Butyl Octyl Ester	84-78-6	487
Phthalic Acid Diallyl Ester	131-17-9	487
Phthalic Acid Dibenzyl Ester	523-31-9	488
Phthalic Acid Dibutyl Ester (3,4,5,6)-D4	93952-11-5	488, 550
Phthalic Acid Dicyclohexyl Ester	84-61-7	488
Phthalic Acid Dicyclohexyl Ester D4	358731-25-6	488, 550
Phthalic Acid Didecyl Ester	84-77-5	488
Phthalic Acid Diethyl Ester (3,4,5,6)-D4	93952-12-6	488, 550
Phthalic Acid Diheptyl Ester	3648-21-3	488
Phthalic Acid Dihexyl Ester	84-75-3	488
Phthalic Acid Di-n-hexyl Ester (3,4,5,6)-D4	1015854-55-3	488, 551
Phthalic Acid Diisobutyl Ester	84-69-5	489
Phthalic Acid Diisobutyl Ester (3,4,5,6)-D4	358730-88-8	489, 551
Phthalic Acid Di-iso-Octyl Ester	27554-26-3	489

Analyte name	CAS No#	Page No#
Phthalic Acid Diisopropyl Ester	605-45-8	489
Phthalic Acid Dimethyl Ester	131-11-3	489
Phthalic Acid Dimethyl Ester D4	93951-89-4	489, 551
Phthalic Acid Dinonyl Ester	84-76-4	489
Phthalic Acid Di-n-nonyl Ester (3,4,5,6)-D4	1202865-43-7	489, 551
Phthalic Acid Di-n-octyl Ester	117-84-0	489
Phthalic Acid Di-n-octyl Ester (3,4,5,6)-D4	93952-13-7	490, 551
Phthalic Acid Dipentyl Ester	131-18-0	490
Phthalic Acid Di-n-pentyl Ester (3,4,5,6)-D4	358730-89-9	490, 551
Phthalic Acid Diphenyl Ester	84-62-8	490
Phthalic Acid Dipropyl Ester	131-16-8	490
Phthalic Acid Dipropyl Ester (3,4,5,6-D4)	358731-29-0	490, 551
Phthalic Acid 2-Ethylhexyl Hexyl Ester	75673-16-4	490
Phthalic acid, mixed decyl-hexyl-octyl diester	68648-93-1	491
Phthalic Acid Monobutyl Ester	131-70-4	490
Phthalic Acid Mono-2-ethylhexyl Ester	4376-20-9	490
Phthalic acid, mono-2-ethyl-5-oxohexyl ester	40321-98-0	491
Phthalic Acid Monomethyl Ester	4376-18-5	491
Phthalic Acid Octyl Decyl Ester	119-07-3	491
Phthalic Acid n-Pentyl Isopentyl Ester (3,4,5,6)-D4 (mixture of isomers)		491, 551
Phthalic Acid n-Pentyl Isopentyl Ester (mixture of isomers)	84777-06-0	491
Phthalic Acid n-Pentyl Isopentyl Ester	776297-69-9	491
Phthalic Anhydride	85-44-9	491
Phthalide (4,5,6,7-Tetrachlorophthalide)	27355-22-2	23
Phthalimide	85-41-6	414
N4-Phthalylsulfathiazole	85-73-4	363
Phytane	638-36-8	296
Phytomenadione (Vitamin K1)	n/a	476
Picene	213-46-7	256
Picloram	1918-02-1	86

Index

Analyte name	CAS No#	Page No#
Picloram-2-ethylhexyl ester	36374-99-9	86
Picolinafen	137641-05-5	86
Picoxystrobin	117428-22-5	23
Pigment Blue 15	147-14-8	414
Pigment Green 7	1328-53-6	414
Pigment Orange 5	3468-63-1	414
Pigment Red 4	2814-77-9	414
Pigment Red 64	6371-76-2	414
Pigment Violet 3	1325-82-2	414
Pindolol	13523-86-9	363
Pindone	83-26-1	190
α-Pinene	80-56-8	603
β-Pinene	127-91-3	603
Pinoxaden	243973-20-8	86
Piperalin	3478-94-2	23
Piperazine	110-85-0	363
Piperazine Citrate Hydrate	41372-10-5	363
Piperonyl Butoxide	51-03-6	159
Piperophos	24151-93-7	86
Piracetam	7491-74-9	363
Pirimetaphos	31377-69-2	159
Pirimicarb	23103-98-2	159
Pirimicarb D6 (dimethylcarbamate D6)	1015854-66-6	159, 551
Pirimicarb-desamido	40778-16-3	159
Pirimicarb-desamido-desmethyl	78195-30-9	159
Pirimicarb-desmethyl	30614-22-3	159
Pirimicarb-desmethyl-formamido	27218-04-8	159
Pirimiphos-ethyl	23505-41-1	160
Pirimiphos-ethyl D10 (diethoxy D5)	n/a	160, 552
Pirimiphos-ethyl-oxon	36378-61-7	160
Pirimiphos-methyl	29232-93-7	160
Pirimiphos-methyl D6 (dimethoxy D3)	1793055-06-7	160, 552
Pirimiphos-methyl-N-desethyl	67018-59-1	160
Piroxicam	36322-90-4	363
Plifenate	21757-82-4	160

Analyte name	CAS No#	Page No#
Polycarbamate	64440-88-6	23
Poly(hexamethylenediguanide) Hydrochloride	32289-58-0	603
Polymyxin B sulfate	1405-20-5	363
Polyoxin B	19396-06-6	23
Ponazuril	69004-04-2	363
Ponceau 4RC (E124)	2611-82-7	415
Ponceau SX	4548-53-2	415
Potassium Clavulanate (Clavulanic acid potassium salt)	61177-45-5	364
Potassium Cyanate	590-28-3	603
Potassium Sorbate	24634-61-5	476
Prallethrin	23031-36-9	160
Praziquantel	55268-74-1	364
Prednisolone	50-24-8	364
Prednisolone Acetate	52-21-1	364
Prednisone	53-03-2	364
Prednisone Acetate	125-10-0	364
Pregnenolone	145-13-1	364
Pretilachlor	51218-49-6	86
Primidone	125-33-7	364
Primisulfuron-methyl	86209-51-0	87
Pristane	1921-70-6	296
Probenazole	27605-76-1	23
Procaine Hydrochloride	51-05-8	364
Prochloraz	67747-09-5	23
Prochloraz D7 (propyl D7)	n/a	23, 552
Prochloraz desimidazole-amino	139520-94-8	24
Prochloraz desimidazole-formylamino	139542-32-8	24
Prochlorperazine maleate	84-02-6	365
Procyazine	32889-48-8	87
Procymidone	32809-16-8	24
Prodiamine	29091-21-2	87
Profenofos	41198-08-7	160
Profluralin	26399-36-0	87
Profoxydim-lithium	281664-76-4	87
Progesterone	57-83-0	365
Prohexadione-calcium	127277-53-6	191

Index

Analyte name	CAS No#	Page No#
Prohydrojasmon	158474-72-7	191
Proline (L-(-)-Proline)	147-85-3	604
Promazine Hydrochloride	53-60-1	365
Promecarb	2631-37-0	161
Promethazine Hydrochloride	58-33-3	365
Prometon	1610-18-0	87
Prometon D14 (di-iso-propyl D14)	1219803-45-8	87, 552
Prometryn	7287-19-6	87
Prometryn D6 (isopropyl D6)	1705649-52-0	87, 552
Propachlor	1918-16-7	88
Propachlor D5 (phenyl D5)	n/a	88, 552
Propachlor ethane sulfonic acid (ESA) sodium salt	947601-88-9	88
Propachlor-2-hydroxy	42404-06-8	88
Propachlor Oxalamic Acid	70628-36-3	88
Propamocarb free base	24579-73-5	24
Propamocarb free base D7 (O-propyl D7)	1398065-89-8	24, 552
Propamocarb hydrochloride	25606-41-1	24
Propanal (Propionaldehyde)	123-38-6	296
1,3-Propanediol	504-63-2	296
1-Propanethiol	107-03-9	297
2-Propanethiol	75-33-2	297
Propanil	709-98-8	88
Propanoic Acid Ethyl Ester	105-37-3	476
1-Propanol	71-23-8	297
2-Propanol (Isopropyl alcohol)	67-63-0	297
Propaphos	7292-16-2	161
Propaquizafop	111479-05-1	88
Propargite	2312-35-8	161
Propargyl Alcohol	107-19-7	297
Propazine	139-40-2	88
Propazine 13C3 (ring 13C3)	446276-68-2	88, 552
Propazine D6 (isopropyl D6)	1655498-05-7	89, 552
Propazine-2-hydroxy	7374-53-0	89
Propetamphos	31218-83-4	161

Analyte name	CAS No#	Page No#
Propham	122-42-9	89
Propiconazole	60207-90-1	24
Propiconazole D5 (2,2,3,3,3-propyl-D5)	1246818-14-3	24, 553
Propineb	12071-83-9	24
Propionaldehyde-2,4-dinitrophenylhydrazone	725-00-8	297
Propionic Acid (Propanoic acid)	79-09-4	477
Propionic Acid Butyl Ester	590-01-2	477
Propionic Acid Calcium Salt	4075-81-4	477
Propionic Acid Methyl Ester	554-12-1	477
Propionitrile	107-12-0	297
Propionylpromazine hydrochloride	7681-67-6	365
Propiophenone	93-55-0	365
Propisochlor	86763-47-5	89
Propoxur	114-26-1	161
Propoxur D3 (N-methyl D3)	1219798-56-7	161, 553
Propoxycarbazone-sodium	181274-15-7	89
Propranolol Hydrochloride	318-98-9	365
Propyl Acetate (Acetic acid n-propyl ester)	109-60-4	477
Propyl Parahydroxybenzoate (4-Hydroxybenzoic acid propyl ester)	94-13-3	477
1-Propylamine (n-Propylamine)	107-10-8	297
n-Propylbenzene	103-65-1	297
n-Propylcarbamate	627-12-3	604
Propylcyclohexane	1678-92-8	298
1,2-Propylene carbonate	108-32-7	298
Propylene Glycol (Propane-1,2-diol)	57-55-6	298
Propylene Glycol 1-Methyl Ether 2-Acetate	108-65-6	298
Propylene Oxide	75-56-9	298
4-Propylphenol	645-56-7	445
Propylthiouracil (6-Propyl-2-thiouracil)	51-52-5	604
Propylthiourea	927-67-3	604
Propyphenazone	479-92-5	365
Propyzamide	23950-58-5	89

Index

Analyte name	CAS No#	Page No#
Propyzamide D3 (phenyl-2,4,6 D3)	1219805-79-4	89, 553
Proquinazid	189278-12-4	25
Prosulfocarb	52888-80-9	89
Prosulfuron	94125-34-5	89
Prothiocarb hydrochloride	19622-19-6	25
Prothioconazole	178928-70-6	25
Prothioconazole-desthio	120983-64-4	25
Prothiophos (Tokuthion®)	34643-46-4	161
Prothoate	2275-18-5	161
Proximpham	2828-42-4	90
Putrescine Dihydrochloride	333-93-7	365
Pymetrozine	123312-89-0	161
Pyracarbolid	24691-76-7	25
Pyraclofos	89784-60-1	162
Pyraclonil	158353-15-2	90
Pyraclostrobin	175013-18-0	25
Pyraflufen (free acid)	129630-17-7	90
Pyraflufen-ethyl	129630-19-9	90
Pyranocumarin	518-20-7	366
Pyrantel Embonate (Pyrantel Pamoate)	22204-24-6	366
Pyrasulfotole	365400-11-9	90
Pyrasulfotole-desmethyl	936621-05-5	90
Pyrazinamide	98-96-4	366
Pyrazolynate	58011-68-0	90
Pyrazophos	13457-18-6	25
Pyrazosulfuron-ethyl	93697-74-6	90
Pyrazoxyfen	71561-11-0	90
Pyrene	129-00-0	256
Pyrene-d10	1718-52-1	257, 553
Pyrethrin 1	121-21-1	162
Pyrethrins	8003-34-7	162
Pyribencarb	799247-52-2	25
Pyribenzoxim	168088-61-7	91
Pyributicarb	88678-67-5	91
Pyridaben	96489-71-3	162

Analyte name	CAS No#	Page No#
Pyridalyl	179101-81-6	162
Pyridaphenthion	119-12-0	162
Pyridate	55512-33-9	91
Pyridine	110-86-1	298
Pyridine-d5	7291-22-7	298, 553
Pyridoxine Hydrochloride	58-56-0	477
Pyrifenox	88283-41-4	25
E-Pyrifenox	83227-22-9	26
Pyrifluquinazon	337458-27-2	162
Pyriftalid	135186-78-6	91
Pyrimethamine	58-14-0	366
Pyrimethanil	53112-28-0	26
Pyrimidifen	105779-78-0	162
(E)-Pyriminobac-methyl	147411-69-6	91
Pyriminobac-methyl (Z)	147411-70-9	91
Pyrimitate	5221-49-8	162
Pyriproxyfen	95737-68-1	163
Pyrithiobac sodium	123343-16-8	91
Pyroquilon	57369-32-1	26
Pyroxsulam	422556-08-9	91
Pyrrolidine	123-75-1	298
Pyrrolidone (Pyrrolidin-2-one)	616-45-5	366
Quinalphos	13593-03-8	163
Quinchlorac	84087-01-4	91
Quinclorac-methyl ester	84087-33-2	92
Quinine	130-95-0	366
Quinine Hydrochloride Dihydrate	6119-47-7	366
Quinmerac	90717-03-6	92
Quinmerac D4 (methyl(D3)-quinoline-2-D)	n/a	92, 553
Quinocetone (technical product)	81810-66-4	366
Quinoclamine	2797-51-5	92
Quinoline	91-22-5	415
Quinoline-d7 (Quinoline-2,3,4,5,6,7,8-D7)	34071-94-8	415, 553
Quinoline Yellow (E104)	8004-92-0	415
Quinonamid	27541-88-4	92
2-Quinoxalinecarboxylic Acid	879-65-2	366

Index

Analyte name	CAS No#	Page No#
Quinoxifen	124495-18-7	26
Quinthiophos	1776-83-6	163
Quintozene	82-68-8	26
Quizalofop-ethyl	76578-14-8	92
Quizalofop-ethyl D3 (3,3,3-D3)	1398065-84-3	92, 553
Quizalofop free acid	76578-12-6	93
Quizalofop free acid D3 (methyl D3)	n/a	93, 553
Quizalofop-methyl	76578-13-7	93
Quizalofop-P-ethyl	100646-51-3	92
Quizalofop-P-tefuryl	200509-41-7	92
Rabenzazole	40341-04-6	26
Rabeprazole Sodium	117976-90-6	367
Ractopamine Hydrochloride	90274-24-1	367
D-(+)-Raffinose Pentahydrate	17629-30-0	477
Rafoxanide	22662-39-1	367
Reactive Blue 19	2580-78-1	415
Resmethrin	10453-86-8	163
Resorcinol (1,3-Dihydroxybenzene)	108-46-3	367
Retene	483-65-8	257
Retinol Palmitate (Vitamin A palmitate)	79-81-2	477
Retinyl Acetate (Vitamin A acetate)	127-47-9	478
RH 5849	112225-87-3	191
L-(+)-Rhamnose Hydrate	10030-85-0	478
Rhodamine B Chloride	81-88-9	415
Ribavirin	36791-04-5	367
(-)-Riboflavin (Vitamin B2)	83-88-5	478
Riboflavin Sodium Phosphate	130-40-5	478
D-(-)-Ribose	50-69-1	478
Ribostamycin Sulfate Salt	53797-35-6	367
Ricinoleic acid methyl ester	141-24-2	478
Rifampicin	13292-46-1	367
Rifamycin SV Sodium Salt Hexahydrate	n/a	367
Rifaximin	80621-81-4	367
Rimsulfuron	122931-48-0	93

Analyte name	CAS No#	Page No#
Rimsulfuron-desulfon	138724-53-5	93
Rimsulfuron metabolite 2	151331-80-5	93
Robenidine Hydrochloride	25875-50-7	368
Ronidazole	7681-76-7	368
Ronidazole D3	1015855-87-4	368, 554
Rotenone	83-79-4	163
Roxarsone	121-19-7	368
Roxithromycin	80214-83-1	368
Rutoside Trihydrate (Rutin Hydrate)	250249-75-3	368
S 421	127-90-2	191
Saccharin	81-07-2	478
Saccharin Sodium Salt Hydrate	128-44-9	478
Saflufenacil	372137-35-4	93
Salbutamol	18559-94-9	368
Salbutamol-D3	1219798-60-3	368, 554
Salbutamol Sulfate	51022-70-9	368
Salicylic Acid	69-72-7	369
Salicylic Acid Benzyl Ester	118-58-1	369
Salinomycin sodium salt	55721-31-8	369
Salmeterol	89365-50-4	369
Salmeterol Xinafoate	94749-08-3	369
Sarafloxacin Hydrochloride	91296-87-6	369
Sarafloxacin D8 hydrochloride trihydrate	n/a	369, 554
Schradan	152-16-9	163
Scopoletin	92-61-5	604
Sebuthylazine	7286-69-3	93
Sebuthylazine D5 (N-ethyl D5)	1219805-56-7	93, 554
Sebuthylazine-desethyl (6-Chloro-N-sec-butyl-1,3,5-triazine-2,4-diamine)	37019-18-4	94
Sebuthylazine-2-hydroxy	33124-61-7	94
Secbumeton	26259-45-0	94
Secbumeton D5 (N-ethyl D5)	1705649-53-1	94, 554
Secnidazole	3366-95-8	369
Selamectin	220119-17-5	369

Index

Analyte name	CAS No#	Page No#
Semicarbazide Hydrochloride	563-41-7	370
DL-Serine	302-84-1	604
Sethoxydim	74051-80-2	94
Siduron	1982-49-6	94
trans-Siduron	19123-21-8	94
Silafluofen	105024-66-6	163
Silthiofam	175217-20-6	26
Simazine	122-34-9	94
Simazine D10 (diethyl D5)	220621-39-6	94, 554
Simazine D5 (ethyl D5)	220621-41-0	95, 554
Simazine-2-hydroxy	2599-11-3	95
Simazine mercapturate	205764-69-8	95
Simeconazole	149508-90-7	26
Simeton	673-04-1	95
Simetryn	1014-70-6	95
Sirolimus	53123-88-9	370
Sodium Arsenate Tetrahydrate	24992-20-9	604
Sodium Ascorbate (Ascorbic acid sodium salt; Vitamin C sodium salt)	134-03-2	478
Sodium Benzoate (Benzoic acid sodium salt)	532-32-1	479
Sodium Cyclamate (N-Cyclohexylsulfamic acid sodium salt)	139-05-9	479
Sodium Propionate (Propionic acid sodium salt)	137-40-6	479
Solvent Green 3	128-80-3	415
Solvent Green 7	6358-69-6	415
Solvent Red 49	509-34-2	416
Solvent Violet 13	81-48-1	416
Solvent Yellow 33	8003-22-3	416
Solvent Yellow 34 Hydrochloride	2465-27-2	416
Sorbic Acid	110-44-1	479
D-Sorbitol	50-70-4	479
L-(-)-Sorbose	87-79-6	479
Sotalol Hydrochloride	959-24-0	370
Sotalol hydrochloride D6 (isopropyl-1,1,1,3,3,3-D6)	1246820-85-8	370, 554

Analyte name	CAS No#	Page No#
Sparfloxacin	110871-86-8	370
Spectinomycin Dihydrochloride Pentahydrate	22189-32-8	370
Spectinomycin Sulfate Tetrahydrate	64058-48-6	370
Spermidine trihydrochloride	334-50-9	370
Spermine tetrahydrochloride	306-67-2	370
Spinetoram	935545-74-7	163
Spinosad	168316-95-8	163
Spiramycin	8025-81-8	371
Spiramycin adipate base	68880-55-7	371
Spirodiclofen	148477-71-8	164
Spiromesifen	283594-90-1	164
Spiromesifen-alcohol	148476-30-6	164
Spiromesifen-alcohol-4-hydroxymethyl	873423-07-5	164
Spirotetramat	203313-25-1	164
Spirotetramat-enol	203312-38-3	164
Spirotetramat-enol-glucoside	1172614-86-6	164
Spirotetramat-keto-hydroxy	1172134-11-0	164
Spirotetramat-mono-hydroxy	1172134-12-1	164
Spiroxamine	118134-30-8	27
Squalane	111-01-3	298
Squalene	111-02-4	371
Stanozolol	10418-03-8	371
Stearic Acid (Octadecanoic acid)	57-11-4	479
Stearic Acid n-Butyl Ester	123-95-5	604
Stearyl Stearate	2778-96-3	604
Streptomycin Sulfate	3810-74-0	371
Strobane	8001-50-1	165
Strychnine	57-24-9	191
Styrene	100-42-5	299
Styrene D8	19361-62-7	299, 554
Succinic Acid (Butanedioic acid)	110-15-6	479
Succinic Acid Dibutyl Ester	141-03-7	479
Succinic Acid Diethyl Ester	123-25-1	480
Succinic Acid Dimethyl Ester	106-65-0	480
Succinylsulfathiazole	116-43-8	371

Index

Analyte name	CAS No#	Page No#
Sucralose	56038-13-2	480
Sucrose (D-(+)-Saccharose)	57-50-1	480
Sudan 1	842-07-9	416
Sudan 1 D5 (phenyl D5)	752211-63-5	416, 555
Sudan 2	3118-97-6	416
Sudan 3	85-86-9	416
Sudan 4	85-83-6	416
Sudan 4 D6 (naphthyl D6v	1014689-18-9	417, 555
Sudan Black B	4197-25-5	417
Sudan Blue 2	17354-14-2	417
Sudan Orange G	2051-85-6	417
Sudan Red 7B	6368-72-5	417
Sudan Red B	3176-79-2	417
Sudan Red G	1229-55-6	417
Sudan Yellow	60-11-7	417
Sulbactam	68373-14-8	371
Sulcofuron-sodium Monohydrate	1049723-02-5	165
Sulcotrione	99105-77-8	95
Sulfabenzamide	127-71-9	371
Sulfacetamide	144-80-9	371
Sulfacetamide Sodium Monohydrate	6209-17-2	372
Sulfachloropyrazine Sodium	23307-72-4	372
Sulfachloropyrazine sodium monohydrate	73398-14-8	372
Sulfachloropyridazine 13C6 (benzene 13C6)	n/a	372, 555
Sulfachlorpyridazine	80-32-0	372
Sulfaclozine Sodium Monohydrate	1392129-96-2	372
Sulfadiazine	68-35-9	372
Sulfadimethoxine	122-11-2	372
Sulfadimethoxine D6 (2,6-dimethoxy D6)	73068-02-7	372, 555
Sulfadimidine (Sulfamethazine)	57-68-1	373
Sulfadoxine	2447-57-6	373
Sulfafurazole (Sulfisoxazole)	127-69-5	373
Sulfaguanidine Monohydrate	6190-55-2	373

Analyte name	CAS No#	Page No#
Sulfallate	95-06-7	95
Sulfamerazine	127-79-7	373
Sulfamerazine D4	1020719-84-9	373, 555
Sulfameter sodium salt	18179-67-4	373
Sulfamethazine-N4-acetyl	100-90-3	373
Sulfamethizole	144-82-1	373
Sulfamethoxazole	723-46-6	374
Sulfamethoxazole D4 (benzene D4)	1020719-86-1	374, 555
Sulfamethoxypyridazine	80-35-3	374
Sulfametoxydiazine (Sulfameter)	651-06-9	374
Sulfamide	7803-58-9	374
Sulfamonomethoxine	1220-83-3	374
Sulfamonomethoxine sodium	38006-08-5	374
Sulfamoxole	729-99-7	374
Sulfanilamide	63-74-1	374
Sulfanilamide-N-acetyl	121-61-9	375
Sulfanitran	122-16-7	375
Sulfaphenazole	526-08-9	375
Sulfapyridine	144-83-2	375
Sulfaquinoxaline	59-40-5	191
Sulfathiazole	72-14-0	375
Sulfatroxazole	23256-23-7	375
Sulfentrazone	122836-35-5	95
Sulfisomidine Sodium	2462-17-1	375
Sulfisozole sodium	37514-39-9	375
Sulfluramid	4151-50-2	165
Sulfolan	126-33-0	299
Sulfometuron-methyl	74222-97-2	95
Sulfosulfuron	141776-32-1	96
Sulfotep	3689-24-5	165
Sulfoxaflo	946578-00-3	165
Sulfur	7704-34-9	27
Sulprofos	35400-43-2	165
Sunset Yellow (E110)	2783-94-0	417
Syringaldehyde	134-96-3	605
2,4,5-T 2-Butoxyethyl Ester	2545-59-7	96

Analyte name	CAS No#	Page No#
2,4,5-T 1-Butoxy-1-methylethyl Ester	n/a	96
2,4,5-T Butyl Ester	93-79-8	96
2,4,5-T-ester (technical)	n/a	96
2,4,5-T 2-Ethylhexyl Ester	1928-47-8	96
2,4,5-T Isooctyl Ester	25168-15-4	96
2,4,5-T Isopropyl Ester	93-78-7	96
2,4,5-T Methyl Ester ((2,4,5-Trichlorophenoxy)acetic Acid Methyl Ester)	1928-37-6	96
2,4,5-T ((2,4,5-Trichlorophenoxy) acetic Acid)	93-76-5	97
2,4,5-T D4	358731-37-0	97, 555
Tadalafil	171596-29-5	375
Tallowamine polyethoxylated	61791-26-2	605
Tamoxifen Citrate	54965-24-1	376
DL-Tartaric acid	133-37-9	605
Tartrazine (E102)	1934-21-0	418
TCMTB (Busan; 2-(Thiocyanomethylthio) benzothiazole)	21564-17-0	27
Tebuconazole	107534-96-3	27
Tebuconazole D6	n/a	27, 555
Tebufenozide	112410-23-8	165
Tebufenpyrad	119168-77-3	165
Tebupirimfos	96182-53-5	166
Tebutam	35256-85-0	97
Tebuthiuron	34014-18-1	97
Tecloftalam	76280-91-6	27
Tecnazene (2,3,5,6-Tetrachloronitrobenzene)	117-18-0	27
Teflubenzuron	83121-18-0	166
Tefluthrin	79538-32-2	166
Tefuryltrione	473278-76-1	97
Tembotrione	335104-84-2	97
Temephos	3383-96-8	166
Temephos D12 (tetramethyl D12)	1219795-39-7	166, 555
Tenoxicam	59804-37-4	376
O,O-TEPP	107-49-3	165

Analyte name	CAS No#	Page No#
O,S-TEPP	645-78-3	605
Tepraloxydim	149979-41-9	97
Tepraloxydim-5-hydroxy	154749-55-0	97
Terbacil	5902-51-2	98
Terbinafine Hydrochloride	78628-80-5	376
Terbucarb	1918-11-2	98
Terbufos	13071-79-9	166
Terbufos-oxon	56070-14-5	166
Terbufos-oxon-sulfone	56070-15-6	166
Terbufos-oxon-sulfoxide	56165-57-2	166
Terbufos-sulfone	56070-16-7	167
Terbufos-sulfoxide	10548-10-4	167
Terbumeton	33693-04-8	98
Terbumeton-desethyl (N-tert-Butyl-6-methoxy-1,3,5-triazine-2,4-diamine)	30125-64-5	98
Terbutaline Sulfate	23031-32-5	376
Terbuthylazine	5915-41-3	98
Terbuthylazine D5 (ethyl D5)	222986-60-9	98, 556
Terbuthylazine-desethyl (N-tert-Butyl-6-chloro-1,3,5-triazine-2,4-diamine)	30125-63-4	98
Terbuthylazine-desethyl D9 (tert-butyl D9)	1219798-52-3	98, 556
Terbuthylazine-desethyl-2-hydroxy (4-Amino-6-(tert-butylamino)-1,3,5-triazin-2-ol)	66753-06-8	99
Terbuthylazine-2-hydroxy (4-(tert-Butylamino)-6-(ethylamino)-1,3,5-triazin-2-ol)	66753-07-9	99
Terbutryn	886-50-0	99
Terbutryn D5 (ethyl D5)	1219804-47-3	99, 556
Terephthalic Acid	100-21-0	605
Terephthalic Acid Bis(2-ethylhexyl) Ester	6422-86-2	605
m-Terphenyl	92-06-8	299
o-Terphenyl	84-15-1	299
p-Terphenyl	92-94-4	299
p-Terphenyl D14	1718-51-0	299, 556
α-Terpinene	99-86-5	605

Analyte name	CAS No#	Page No#
γ-Terpinene	99-85-4	605
1-Terpinen-4-ol	562-74-3	606
α-Terpineol	98-55-5	606
Terpinolene (δ-Terpinene)	586-62-9	606
Testosterone	58-22-0	376
Testosterone 17-Acetate	1045-69-8	376
Testosterone 17-Propionate	57-85-2	376
1,2,4,5-Tetrabromobenzene	636-28-2	445
Tetrabromobisphenol A-bis(2,3-dibromopropyl Ether)	21850-44-2	606
Tetrabromobisphenol A D10 (dimethyl D3, bisphenol-3,5-D2)	n/a	556, 606
Tetrabromobisphenol A Diallyl Ether	25327-89-3	606
3,3',5,5'-Tetrabromobisphenol A	79-94-7	606
1,1,2,2-Tetrabromoethane	79-27-6	392
Tetrabromomethane	558-13-4	392
Tetrabutyltin	1461-25-2	606
Tetracaine	94-24-6	376
Tetracaine Hydrochloride	136-47-0	376
2,3,4,5-Tetrachloroaniline	634-83-3	446
2,3,5,6-Tetrachloroaniline	3481-20-7	191
2,3,4,5-Tetrachloroanisole	938-86-3	446
2,3,4,6-Tetrachloroanisole	938-22-7	606
2,3,5,6-Tetrachloroanisole	6936-40-9	446
3,3',4,4'-Tetrachloroazobenzene	14047-09-7	418
3,3',4,4'-Tetrachloroazoxybenzene	21232-47-3	418
1,2,3,4-Tetrachlorobenzene	634-66-2	191
1,2,3,5-Tetrachlorobenzene	634-90-2	446
1,2,4,5-Tetrachlorobenzene	95-94-3	446
1,2,4,5-Tetrachlorobenzene-13C6	85380-73-0	446, 556
1,2,9,10-Tetrachlorodecane (CP-2)	205646-11-3	607
2,5,6,9-Tetrachlorodecane (CP-1)	n/a	607
2,3,7,8-Tetrachlorodibenzo-p-dioxine	1746-01-6	607
1,1,1,2-Tetrachloroethane	630-20-6	392
1,1,2,2-Tetrachloroethane	79-34-5	392
Tetrachloroethene	127-18-4	392

Analyte name	CAS No#	Page No#
Tetrachloroguaiacol	2539-17-5	607
Tetrachlorohydroquinone	87-87-6	607
1,2,3,4-Tetrachloronaphthalene	20020-02-4	257
2,3,4,5-Tetrachloronitrobenzene	879-39-0	446
2,3,4,6-Tetrachlorophenol Acetate	5435-60-9	446
2,3,5,6-Tetrachlorophenol Acetate	61925-90-4	446
2,3,4,5-Tetrachlorophenol	4901-51-3	447
2,3,4,6-Tetrachlorophenol	58-90-2	191, 447
2,3,5,6-Tetrachlorophenol	935-95-5	192
2,3,5,6-Tetrachloropyridine	2402-79-1	447
3,3",4,4"-Tetrachloro-o-terphenyl	n/a	233
3,3",4,4"-Tetrachloro-p-terphenyl	1957180-49-2	233
3,3",5,5"-Tetrachloro-o-terphenyl	566197-75-9	233
3,3",5,5"-Tetrachloro-p-terphenyl	500729-81-7	233
α,α,2,6-Tetrachlorotoluene	81-19-6	447
α,α,α,4-Tetrachlorotoluene	5216-25-1	447
2,4,5,6-Tetrachloro-m-xylene	877-09-8	447
Tetrachlorvinphos	22248-79-9	167
Tetraconazole	112281-77-3	27
n-Tetracontane	4181-95-7	299
n-Tetracosane	646-31-1	299
Tetracosanoic Acid (Lignoceric acid)	557-59-5	480
Tetracosanoic Acid Methyl Ester	2442-49-1	480
1-Tetracosanol (Lignoceryl alcohol)	506-51-4	480
Tetracycline Hydrochloride	64-75-5	377
n-Tetradecane	629-59-4	300
Tetradecanoic Acid (Myristic acid)	544-63-8	480
Tetradecanoic Acid Ethyl Ester	124-06-1	480
Tetradecanoic Acid Methyl Ester	124-10-7	481
1-Tetradecanol (Myristyl alcohol)	112-72-1	300
1-Tetradecene	1120-36-1	300
Tetradecylamine	2016-42-4	300
2-Tetradecylcyclobutanone	35493-47-1	300
3,3,4,4-Tetradeuterioadipic Acid Bis(2-ethylhexyl) Ester	n/a	491, 556

Index

Analyte name	CAS No#	Page No#
Tetradifon	116-29-0	167
Tetrahydrofuran (THF)	109-99-9	300
1,2,3,4-Tetrahydronaphthalene	119-64-2	257
cis-1,2,3,6-Tetrahydrophthalimide	1469-48-3	192
Tetramethrin	7696-12-0	167
1,2,3,4-Tetramethylbenzene	488-23-3	300
1,2,3,5-Tetramethylbenzene (Isodurene)	527-53-7	300
1,2,4,5-Tetramethylbenzene (Durene)	95-93-2	300
3,3',5,5'-Tetramethylbenzidine	54827-17-7	418
2,4,7,9-Tetramethyl-5-decyne-4,7-diol	126-86-3	607
Tetramisole Hydrochloride	5086-74-8	377
Tetramisole D5 hydrochloride	1173021-85-6	377, 556
Tetra-n-pentyltin	3765-65-9	607
Tetraphenyltin	595-90-4	607
Tetra-n-propyltin	2176-98-9	607
Tetrasul	2227-13-6	167
n-Tetratetracontane	7098-22-8	301
n-Tetratriacontane	14167-59-0	301
Thenylchlor	96491-05-3	99
Thiabendazole-5-hydroxy (5-Hydroxytiabendazole)	948-71-0	28
Thiacloprid	111988-49-9	167
Thiacloprid-amide	676228-91-4	167
Thiamazole (Methimazole)	60-56-0	377
Thiamethoxam	153719-23-4	167
Thiamethoxam D4 (oxadiazine D4)	1331642-98-8	168, 556
Thiamethoxam urea	902493-06-5	168
Thiamine Hydrochloride	67-03-8	481
Thiamine Nitrate	532-43-4	481
Thiamphenicol	15318-45-3	377
Thiazafluron	25366-23-8	99
Thiazopyr	117718-60-2	99
Thidiazuron	51707-55-2	192
Thiencarbazone-methyl	317815-83-1	99

Analyte name	CAS No#	Page No#
Thifensulfuron-methyl	79277-27-3	99
Thifensulfuron-methyl D3	n/a	100, 557
Thifluzamide	130000-40-7	28
Thiobencarb	28249-77-6	100
Thiocyclam hydrogenoxalate	31895-22-4	168
2,2'-Thiodiacetic Acid	123-93-3	608
Thiodicarb	59669-26-0	168
Thiodiphosphoric Acid Tetrapropyl Ester	3244-90-4	168
Thiofanox	39196-18-4	168
Thiofanox D3 (N-methyl D3)	74829-46-2	168, 557
Thiofanox-sulfone	39184-59-3	168
Thiofanox-sulfoxide	39184-27-5	168
Thiomersal	54-64-8	377
Thiometon	640-15-3	169
Thiometon-sulfone	20301-63-7	169
Thiometon-sulfoxide	2703-37-9	169
Thionazin	297-97-2	169
Thionazin D10 (diethyl D10)	1173019-39-0	169, 557
Thiophanate (Thiophanate-ethyl)	23564-06-9	28
Thiophanate-methyl	23564-05-8	28
Thiophene	110-02-1	301
Thiophenol	108-98-5	447
Thioquinox	93-75-4	28
Thiosultap-sodium	52207-48-4	169
2-Thiouracil	141-90-2	377
Thioxothiamine (Thiothiamine)	299-35-4	481
Thiram	137-26-8	28
Thiram D12	69193-86-8	28, 557
DL-Threonine	80-68-2	608
Thymol	89-83-8	447
Tiabendazole (Thiabendazole)	148-79-8	28
Tiabendazole D6 (Thiabendazole D6)	1262551-89-2	28, 557
Tiadinil	223580-51-6	29
Tiamulin	55297-95-5	377
Tiamulin Hydrogen Fumarate	55297-96-6	377

Index

Analyte name	CAS No#	Page No#
Tiletamine Hydrochloride	14176-50-2	378
Tilmicosin	108050-54-0	378
Timolol Maleate	26921-17-5	378
Tinidazole	19387-91-8	378
Tiocarbazil	36756-79-3	100
Tiopronin	1953-02-2	378
Tobramycin	32986-56-4	378
DL-α-Tocopherol Acetate (Vitamin E acetate)	52225-20-4	481
all-rac-α-Tocopherol (Vitamin E)	10191-41-0	481
Tolbutamide	64-77-7	378
Tolclofos-methyl	57018-04-9	29
Tolfenamic Acid	13710-19-5	378
Tolfenpyrad	129558-76-5	169
Tolnaftate	2398-96-1	378
Tolnifanide	304911-98-6	29
Toltrazuril	69004-03-1	379
Toltrazuril-sulfoxide	69004-15-5	379
m-Tolualdehyd-2,4-dinitrophenylhydrazone	2880-05-9	447
o-Tolualdehyd-2,4-dinitrophenylhydrazone	1773-44-0	448
p-Tolualdehyd-2,4-dinitrophenylhydrazone	2571-00-8	448
Toluene (Methylbenzene)	108-88-3	301
Toluene D8	2037-26-5	301, 557
2,4-Toluenediisocyanate	584-84-9	608
2,6-Toluenediisocyanate	91-08-7	608
p-Toluenesulfonamide	70-55-3	608
Toluidine Red	2425-85-6	418
o-Toluidine	95-53-4	448
Tolyfluanid	731-27-1	29
Tolyltriazole	29385-43-1	608
Topramezone	210631-68-8	100
Toxaphene (Camphechlor)	8001-35-2	169
Toxaphene Parlar 11	165820-10-0	608
Toxaphene Parlar 26	142534-71-2	608
Toxaphene Parlar 32	51775-36-1	608

Analyte name	CAS No#	Page No#
Toxaphene Parlar 40	166021-27-8	609
Toxaphene Parlar 42	187348-02-3	609
Toxaphene Parlar 44	165820-17-7	609
Toxaphene Parlar 50	66860-80-8	609
Toxaphene Parlar 62	154159-06-5	609
Toxaphene Parlar 69	151183-19-6	609
2,4,5-TP butoxyethyl ester	19398-13-1	97
Tralkoxydim	87820-88-0	100
Tralomethrin	66841-25-6	169
Tralopyril	122454-29-9	192
Transfluthrin	118712-89-3	170
D-(+)-Trehalose Dihydrate	6138-23-4	481
Trenbolone	10161-33-8	379
Trenbolone Acetate	10161-34-9	379
Tretinoin	302-79-4	379
Triacetin	102-76-1	481
n-Triacontane	638-68-6	301
Triacontanoic Acid Methyl Ester	629-83-4	481
1-Triacontanol (Melissyl alcohol)	593-50-0	192
Triadimefon	43121-43-3	29
Triadimenol	55219-65-3	29
erythro-Triadimenol (isomer A)	70585-35-2	29
threo-Triadimenol (isomer B)	70585-37-4	29
Triallate	2303-17-5	100
Triamcinolone Acetonide	76-25-5	379
Triamcinolone Acetonide 21-Acetate	3870-07-3	379
Triamcinolone 16,21-Diacetate	67-78-7	379
Triamiphos	1031-47-6	30
Triapenthenol	76608-88-3	192
Triasulfuron	82097-50-5	100
Triazamate	112143-82-5	170
Triaziflam	131475-57-5	100
1H-1,2,4-Triazol	288-88-0	192
Triazophos	24017-47-8	170
Triazoxide	72459-58-6	30
Tribenuron-methyl (technical)	101200-48-0	100
Tribromoacetic Acid	75-96-7	609

Index

Analyte name	CAS No#	Page No#
2,4,6-TribromoanisoIe	607-99-8	448
1,2,4-Tribromobenzene	615-54-3	448
1,3,5-Tribromobenzene	626-39-1	448
2,4,6-Tribromophenol	118-79-6	448
3,5,4'-Tribromosalicylanilide (3,5-Dibromo-N-(4-bromophenyl)-2-hydroxy-benzamide)	87-10-5	609
Tribufos	78-48-8	170
Tri-n-butyl phosphate	126-73-8	492
Tributylamine	102-82-9	301
Tributylethyltin	19411-60-0	609
2,4,6-Tri-tert-butylphenol	732-26-3	448
Tributyltin Chloride (TBTC)	1461-22-9	610
Tributyltin Oxide (TBTO)	56-35-9	610
Tricaine Methanesulfonate	886-86-2	379
Tricaprin	621-71-6	610
Trichlamide (N-(1-Butoxy-2,2,2-trichloroethyl)salicylamide)	70193-21-4	30
Trichloroacetic Acid (TCA)	76-03-9	101
Trichloroacetic Acid Methyl Ester (TCA-methyl ester)	598-99-2	101
Trichloroacetic Acid Sodium Salt	650-51-1	101
Trichloroacetonitrile	545-06-2	392
α,2,4-Trichloroacetophenone	4252-78-2	610
2,3,4-Trichloroaniline	634-67-3	448
2,4,5-Trichloroaniline	636-30-6	449
2,4,6-Trichloroaniline	634-93-5	449
3,4,5-Trichloroaniline	634-91-3	449
2,4,6-TrichloroanisoIe D5	352439-08-8	449, 557
2,3,4-TrichloroanisoIe	54135-80-7	449
2,3,6-TrichloroanisoIe	50375-10-5	449
2,4,5-TrichloroanisoIe	6130-75-2	449
2,4,6-TrichloroanisoIe	87-40-1	449
1,2,3-Trichlorobenzene	87-61-6	449, 450
1,2,4-Trichlorobenzene	120-82-1	450
1,3,5-Trichlorobenzene	108-70-3	450
1,3,5-Trichlorobenzene D3	1198-60-3	450, 557
2,3,6-Trichlorobenzoic Acid	50-31-7	610

Analyte name	CAS No#	Page No#
2,4,6-Trichlorobenzoic Acid	50-43-1	610
1,1,1-Trichloroethane	71-55-6	392
1,1,2-Trichloroethane	79-00-5	392
2,2,2-Trichloroethanol	115-20-8	392
Trichloroethene	79-01-6	393
Trichloronate	327-98-0	170
1,3,5-Trichloro-2-nitrobenzene	18708-70-8	450
2,3,4-Trichloronitrobenzene	17700-09-3	450
2,4,5-Trichloronitrobenzene	89-69-0	450
Trichloronitromethane	76-06-2	170
2,3,4-Trichlorophenol Acetate	61925-89-1	610
2,3,5-Trichlorophenol Acetate	61925-88-0	610
2,3,6-Trichlorophenol Acetate	61925-87-9	610
2,4,5-Trichlorophenol Acetate	5393-75-9	611
2,4,6-Trichlorophenol Acetate	23399-90-8	611
3,4,5-Trichlorophenol Acetate	59190-61-3	611
2,4,5-Trichlorophenol Ethanolamine Salt	56961-45-6	611
2,3,4-Trichlorophenol	15950-66-0	450
2,3,5-Trichlorophenol	933-78-8	450
2,3,6-Trichlorophenol	933-75-5	451
2,4,5-Trichlorophenol	95-95-4	192
2,4,6-Trichlorophenol	88-06-2	451
3,4,5-Trichlorophenol	609-19-8	451
2,4,6-Trichlorophenoxyacetic Acid	575-89-3	611
4-(2,4,5-Trichlorophenoxy) butanoic Acid (2,4,5-TB)	93-80-1	611
1,2,3-Trichloropropane	96-18-4	393
3,5,6-Trichloro-2-pyridinol	6515-38-4	170, 192
3',4,4''-Trichloro-m-terphenyl	n/a	233
2,3,6-Trichlorotoluene	2077-46-5	451
2,4,5-Trichlorotoluene	6639-30-1	451
α,2,4-Trichlorotoluene	94-99-5	451
α,2,6-Trichlorotoluene	2014-83-7	451
α,3,4-Trichlorotoluene	102-47-6	451
α,α,α-Trichlorotoluene	98-07-7	452

Index

Analyte name	CAS No#	Page No#
1,1,1-Trichloro-2,2,2-trifluoroethane	354-58-5	393
1,1,2-Trichloro-1,2,2-trifluoroethane	76-13-1	393
Triclabendazole	68786-66-3	380
Triclabendazole Sulfone	106791-37-1	380
Triclabendazole Sulfoxide	100648-13-3	380
Triclocarban	101-20-2	611
Triclopyr	55335-06-3	101
Triclopyr-butotyl (Triclopyr-2-butoxyethyl ester)	64700-56-7	101
Triclopyr-methyl	60825-26-5	101
Triclopyricarb	902760-40-1	30
Triclosan	3380-34-5	611
Triclosan D3 (2,4-dichlorophenoxy D3)	1020719-98-5	557, 611
Triclosan methyl D3 (methoxy D3)	1020720-00-6	558, 612
Triclosan Methyl ether	4640-01-1	612
n-Tricosane	638-67-5	301
Tricosanoic Acid	2433-96-7	482
Tricosanoic acid ethyl ester	18281-07-7	482
1-Tricosanol	3133-01-5	301
Tri-m-cresyl Phosphate	563-04-2	492
Tri-o-cresyl Phosphate	78-30-8	492
Tri-p-cresyl Phosphate	78-32-0	492
Tricresylphosphate	1330-78-5	492
Tricyclazole	41814-78-2	30
Tricyclohexyltin Chloride	3091-32-5	612
Tridecanal	10486-19-8	482
n-Tridecane	629-50-5	302
Tridecanoic Acid (Tridecylic acid)	638-53-9	482
Tridecanoic Acid Methyl Ester	1731-88-0	482
1-Tridecanol	112-70-9	302
2-Tridecanone	593-08-8	302
Tridemorph	81412-43-3	30
Tridiphane (2-(3,5-Dichlorophenyl)-2-(2,2,2-trichloroethyl)oxirane)	58138-08-2	101
Trietazine	1912-26-1	101

Analyte name	CAS No#	Page No#
Trietazine D5 (ethyl D5)	1397243-73-0	101, 558
Trietazine-desethyl	38902-68-0	102
Trietazine-2-hydroxy	13532-25-7	102
Triethylamine	121-44-8	302
1,3,5-Triethylbenzene	102-25-0	302
Triethyllead Chloride	1067-14-7	612
Triethylphosphate	78-40-0	492
O,O,O-Triethylphosphorothioate	126-68-1	193
Triethyltin Chloride	994-31-0	612
Trifenmorph	1420-06-0	170
Trifloxystrobin	141517-21-7	30
Trifloxysulfuron sodium	199119-58-9	102
Triflumizole	68694-11-1	30
Triflumizole 13C3 (O-propyl 13C3)	n/a	30, 558
Triflumizole-amino	131549-75-2	31
Triflumuron	64628-44-0	170
Trifluoroacetic acid	76-05-1	612
1,1,1-Trifluoroethane	420-46-2	393
2,2,2-Trifluoroethanol	75-89-8	393
4-(Trifluoromethoxy)aniline	461-82-5	193
2-Trifluoromethylaniline	88-17-5	452
3-(Trifluoromethyl)aniline	98-16-8	193
α,α,α-Trifluorotoluene (Trifluoromethylbenzene)	98-08-8	452
Triflupromazine Hydrochloride	1098-60-8	380
Trifluralin	1582-09-8	102
Trifluralin D14 (di-n-propyl D14)	347841-79-6	102, 558
Triflusulfuron-methyl	126535-15-7	102
Triforine	26644-46-2	31
Triforine D8 (piperazine D8)	948595-12-8	31, 558
Trihexyphenidyl Hydrochloride (Benzhexol Hydrochloride)	52-49-3	380
2,3,5-Triiodobenzoic Acid	88-82-4	612
2,3,5-Trimethacarb	2655-15-4	171
3,4,5-Trimethacarb	2686-99-9	171
Trimethoprim	738-70-5	380
Trimethoprim D3 (4-methoxy D3)	1189923-38-3	380, 558

Index

Analyte name	CAS No#	Page No#
1,2,4-Trimethoxybenzene	135-77-3	452
1,3,5-Trimethoxybenzene	621-23-8	452
2,4,5-Trimethylaniline	137-17-7	452
1,2,3-Trimethylbenzene	526-73-8	302
1,2,4-Trimethylbenzene	95-63-6	302
1,3,5-Trimethylbenzene (Mesitylene)	108-67-8	302
O,O,S-Trimethyldithiophosphate	2953-29-9	193
Trimethyllead Chloride	1520-78-1	612
2,3,5-Trimethylnaphthalene	2245-38-7	257
1,1,1-Trimethylolpropane Phosphate	1005-93-2	612
2,2,4-Trimethylpentane (Isooctane)	540-84-1	302
2,2,4-Trimethyl-1,3-pentanediol Diisobutyrate (TXIB)	6846-50-0	492
2,4,4-Trimethyl-1-pentene	107-39-1	303
2,3,5-Trimethylphenol	697-82-5	193
2,3,6-Trimethylphenol	2416-94-6	452
2,4,6-Trimethylphenol	527-60-6	452
3,4,5-Trimethylphenol	527-54-8	193
Trimethylphosphate	512-56-1	613
O,O,O-Trimethylphosphorthioate (O,O,O-Trimethylthiophosphate)	152-18-1	613
O,O,S-Trimethylphosphorthioate (O,O,S-Trimethylthiophosphate)	152-20-5	613
Trimethylsulfonium Iodide	2181-42-2	613
Trimethyltin Chloride	1066-45-1	613
Trinexapac	143294-89-7	193
Trinexapac-ethyl	95266-40-3	193
Trinexapac-methyl	104273-71-4	193
1,3,5-Trinitrobenzene	99-35-4	452
2,4,7-Trinitro-9-fluorenone	129-79-3	257
2,4,6-Trinitrophenol (Picric Acid)	88-89-1	453
2,4,5-Trinitrotoluene	610-25-3	453
2,4,6-Trinitrotoluene (TNT)	118-96-7	453
Trioctyltin Chloride	2587-76-0	613
Tripelennamine Hydrochloride	154-69-8	380
1,3,5-Triphenylbenzene	612-71-5	303
Triphenylene	217-59-4	257

Analyte name	CAS No#	Page No#
Triphenylphosphate	115-86-6	492
Triphenylphosphine	603-35-0	613
Tripropyl Phosphate	513-08-6	492
Tripropylene Glycol Monomethyl Ether (Mixture of Isomers)	25498-49-1	613
Tripropyltin Chloride	2279-76-7	613
Tris-(2-butoxyethyl)phosphate	78-51-3	493
Tris(2-chloro-1-(chloromethyl) ethyl) Phosphate	13674-87-8	493
Tris-(2-chloroethyl)phosphate	115-96-8	493
Tris-(2-chloropropyl)phosphate	13674-84-5	493
Tris(2,3-dibromopropyl)phosphate	126-72-7	493
Tris(2,4-di-tert-butylphenyl) phosphite	31570-04-4	493
Tris-(2-ethylhexyl)phosphate	78-42-2	493
Triticonazole	131983-72-7	31
Tritosulfuron	142469-14-5	102
n-Tritriacontane	630-05-7	303
Troleandomycin (Oleandomycin triacetate)	2751-09-9	380
Trometamol	77-86-1	381
Tryptophan	54-12-6	614
Tulobuterol Hydrochloride	56776-01-3	381
D-(+)-Turannose	547-25-1	482
Tylosin Tartrate	74610-55-2	381
Tyramine Hydrochloride	60-19-5	453
L-Tyrosine disodium salt	69847-45-6	614
DL-Tyrosine	556-03-6	614
Ugilec 141	111483-93-3	234
Undecanal	112-44-7	303
n-Undecane	1120-21-4	303
Undecanoic Acid (Undecylic acid)	112-37-8	482
Undecanoic Acid Ethyl Ester	627-90-7	482
Undecanoic Acid Methyl Ester	1731-86-8	482
2-Undecanone	112-12-9	303
1-Undecene	821-95-4	303
Uniconazole	83657-22-1	102, 194
Uniconazole-P	83657-17-4	102
Urbacide	2445-07-0	31
Urethan (Ethyl carbamate)	51-79-6	614

Analyte name	CAS No#	Page No#
Urotropine (Methenamine)	100-97-0	614
(+)-Valencene (Valencene sesquiterpene)	4630-07-3	614
Valeric Acid Methyl Ester	624-24-8	483
n-Valeric Acid (n-Pentanoic Acid)	109-52-4	483
Validamycin	37248-47-8	31
Valifenalate	283159-90-0	31
DL-Valine	516-06-3	614
Valnemulin hydrochloride	133868-46-9	381
Vamidothion	2275-23-2	171
Vamidothion Sulfone	70898-34-9	171
Vamidothion Sulfoxide	20300-00-9	171
Vancomycin Hydrochloride	1404-93-9	381
Vardenafil Hydrochloride	224789-15-5	381
VAT Red 1	2379-74-0	418
Verbenone	80-57-9	194
Vernolate	1929-77-7	103
Victoria Blue B	2580-56-5	418
Vinclozolin	50471-44-8	31
Vinyl Acetate	108-05-4	303
Vinyl Bromide	593-60-2	614
Vinyl Chloride	75-01-4	393
4-Vinyl-1-cyclohexene	100-40-3	303
Warfarin	81-81-2	194
(±)-Warfarin D5 (phenyl-D5)	75472-93-4	194, 558
Xylazine Hydrochloride	23076-35-9	381
Xylitol	87-99-0	483
D-(+)-Xylose	58-86-6	483
Xylylcarb	2425-10-7	171
Yohimbine Hydrochloride	65-19-0	381
Zeaxanthin	144-68-3	483
Zeranol (racemic)	55331-29-8	381
Zinc neodecanoate	27253-29-8	615
Zinc Phosphate	7779-90-0	614
Zinc Picolinate	17949-65-4	483
Zineb	12122-67-7	31
Ziram	137-30-4	32
Zoxamide	156052-68-5	32



A global leader
in measurement
standards, reference
materials, laboratory
services and
proficiency testing.

2,600

professionals working in

21

countries

We are the UK's designated **National Measurement Institute for chemical and bioanalytical measurement.**

With 2,600 professionals working in 21 countries, our analytical measurement and quality control services are second-to-none.

As a global leader, **we provide the widest range of reference materials available from any single supplier.**

When you make a decision using our resources, you can be sure it's based on precise, robust data. And together, we're creating fairer, safer, more confident societies worldwide.



lgcstandards.com

Science for a safer world



YOUR FIRST CHOICE IN PROFICIENCY TESTING

LGC is an accredited international provider of proficiency testing (PT) schemes. We have over **30 years** of experience in all aspects of providing PT services to laboratories undertaking chemical, clinical, microbiological and physical measurements.

- Food and feed schemes
- Water and environmental schemes
- Clinical and forensic schemes
- Consumer safety schemes
- Beverage schemes
- Industrial schemes
- QC materials

Annually we run over 1,600 PT exercises, and distribute to more than 11,000 laboratories in 160+ countries, produce in excess of 225,000 test materials and provide a turnaround time of 5 weeks from test material despatch to report (typical cycle).



Please see lgcstandards.com/pt for more information



LGC Quality
ISO/IEC 17043
• ISO 9001



Water and environment PT scheme selector

Scheme	Scheme year	Number of distributions per scheme year	Tests	Test material matrix*	Description group*
Air and stack emissions *AIR	Apr to Mar	Six	Chemical and physical	Filters, tubes and impinger solutions.	Gravimetric, organic and elemental analytes for quantitative and qualitative analyses.
Water Chemistry AQUACHECK	Apr to Mar	Twenty	Chemical, ecotoxicology, physical and radiochemical	Clean waters and waste waters, agricultural soils and dry and raw sewage sludge.	Inorganic, organic and elemental analytes for qualitative and quantitative analyses. Determination of radiochemical and ecotoxicological parameters.
Contaminated land CONTEST	Apr to Mar	Five	Chemical and physical	Soil extracts, soil materials, solid waste and standard solutions.	Inorganic, organic and elemental analytes measured in soil, leachates and standard solutions.
Environmental hygiene monitoring HYGIENE	Jan to Dec	Three	Microbiological	Swabs, contact plates and dip slides.	Routine microbiological testing, indicator organisms and pathogens.
Water microbiology QWAS	Jan to Dec	Ten	Microbiological	Waters (e.g. bathing, environmental, mineral, potable, process, recreational, sea, surface, waste) and simulated effluent sludge.	Routine microbiological testing, indicator organisms and pathogens.

*Please see current application form and scheme description for accredited ad non-accredited status of test materials and analytes.

Consumer safety PT scheme selector

Scheme	Scheme year	Number of distributions per scheme year	Tests	Test material matrix*	Description group*
Cosmetics and toiletries COSMETICS	Apr to Mar	Four	Chemical and microbiological	Cream, lipstick, lipgloss, liquids, mouthwash, toothpaste, powder, soap, detergent.	Chemical parameters of relevance to the cosmetics and toiletries testing industries. Microbiological tests including spoilage and indicator organisms.
Nickel Migration NiMS	Jan to Dec	Two	Chemical and physical	Alloy disks, jewellery or other appropriate articles.	Nickel release and surface area.
Pharmaceutical PHARMASSURE	Apr to Mar	Four	Chemical, microbiological and instrument techniques	Pharmaceutical products and standard solutions.	Basic and advanced chemical analysis, microbiological analysis and sterility testing.
Toy safety TOYTEST	Jan to Dec	Four	Chemical, physical, and assessments against toy safety standards	Toys, paper exercises, real materials and standard solutions	Interpretation of toy safety standards, various physical measurements, azo-dyes, metals and phthalates.

*Please see current application form and scheme description for accredited ad non-accredited status of test materials and analytes.

Beverage PT scheme selector


Scheme	Scheme year	Number of distributions per scheme year	Tests	Test material matrix*	Description group*
Brewing analytes BAPS	Jan to Dec	Twelve - Chemistry Six - Microbiology Twelve - Sensory	Chemical, microbiological and sensory	Ales, lagers, and alcohol free/ low alcohol beers.	Chemical testing for quality control and product characterisation.
Alcoholic drinks DAPS	Jan to Dec	Four	Chemical	Distilled spirits, whisky, wort, ciders, wines and fortified wines, liqueurs/ cream liqueurs, and other alcoholic beverages.	Chemical tests including esters relevant for alcoholic beverages.
Malt analytes MAPS	Jan to Dec	Twelve	Chemical and physical	Brewing/distilling malted barley, barley, malt flour and black/crystal malt.	Chemical and physical tests for quality checks and complex analysis, including mycotoxins analysis
Soft drinks and fruit juice QBS	Jan to Dec	Four	Chemical and microbiological	Carbonated and non carbonated soft drinks, dilutable, ready to drink soft drinks, fruit juices and carbonated water.	Chemical testing for quality checks. Comprehensive range of microorganisms of relevance to beverage products, including pathogens, indicator organisms and spoilage organisms.
Sugar SUPS	Jan to Dec	Twelve	Chemical and microbiological	Cane or beet sugar, raw sugar and molasses.	Chemical tests of relevance to the sugar processing, food and beverage industries. Microorganisms of relevance to sugar products, including pathogens and indicator organisms.

*Please see current application form and scheme description for accredited ad non-accredited status of test materials and analytes.

Food and feed PT scheme selector


Scheme	Scheme year	Number of distributions per scheme year	Tests	Test material matrix*	Description group*
Animal feed AFPS	Jan to Dec	Four	Chemical and microbiological	Animal feed: (e.g. broiler, cattle, chicken, pig, sheep), calf replacer, premix, oil and lard.	Nutritional testing and contaminants. Salmonella, other pathogens and indicator organisms
Chocolate QCS	Jan to Dec	Three	Chemical and microbiological	Chocolate and cocoa powder.	Chemical and microbiological parameters relevant to the chocolate and food testing industries including nutritional and elemental analysis.
Dairy chemistry QDCS	Jan to Dec	Four	Chemical	Butter, cheese, cream, milk, milk powder, whey powder, yoghurt and standard solutions.	Nutritional analysis and contaminants, including toxic elements, pesticides and other contaminants.
Food chemistry QFCS	Jan to Dec	Six	Chemical	Bread, cake, cereals, cured meat, flour, fruit/vegetable, hard cheese, nuts, oils, 'ready to eat' products, tea and standard solutions.	Nutritional and quality analysis, toxic elements, pesticides and other contaminants.
Gelatine QGS	Jan to Dec	Two	Chemical and microbiological	Gelatine, gelatine hydrolysate.	Physicochemical testing and microbiological parameters of relevance to gelatine.
Meat and fish QMAS	Jan to Dec	Six	Chemical and microbiological	Meat, fish and shellfish.	Nutritional testing, contaminants and veterinary drug residues, salmonella and other pathogens, indicator organisms.
Food microbiology QMS	Jan to Dec	Twelve	Microbiological	Oatmeal and skimmed milk powder, tea, herb and spice.	Comprehensive range of microorganisms of relevance to food products, including pathogens, indicator organisms, spoilage organisms and probiotics.

*Please see current application form and scheme description for accredited ad non-accredited status of test materials and analytes.



Visit our PT webshop at lgcstandards.com/pt

For further information, or if you require substances or materials not currently listed please contact one of our local sales offices.



1 BRAZIL
T +55 12 3302 5880
E bz@lgcstandards.com

2 BULGARIA
T +359 (0)2 971 4955
E bg@lgcstandards.com

3 CHINA
T +86 400 9216156
E info.china@lgcgroup.com

4 FRANCE
T +33 (0)3 88 04 82 82
E fr@lgcstandards.com

5 GERMANY
T +49 (0)281 9887 0
E de@lgcstandards.com

EXPORT QUERIES
T +49 (0) 281 9887 250
E global.sales@lgcgroup.com

6 HUNGARY
T +40 364 116890
E ro@lgcstandards.com

7 IRELAND
T +44 (0)20 8943 8480
E uksales@lgcstandards.com

8 ITALY
T +39 02 22476412
E it@lgcstandards.com

9 MIDDLE EAST
T +49 (0) 281 9887 0
E fr@lgcstandards.com

10 NETHERLANDS
T +49 (0)281 9887 250
E nl@lgcstandards.com

11 NORDIC COUNTRIES
T +49 (0)281 9887 0
E de@lgcstandards.com

12 POLAND
T +48 22 751 31 40
E pl@lgcstandards.com

13 ROMANIA
T +40 364 116890
E ro@lgcstandards.com

14 RUSSIA
T +7 812 777 04 88
E ru@lgcgroup.com

15 SOUTH AFRICA
T +27 (0)11 466 4321
E sales.za@lgcgroup.com

16 SPAIN
T +34 (0)93 308 4181
E es@lgcstandards.com

17 TURKEY
T +90 530 880 33 90
E tur@lgcstandards.com

18 UK REFERENCE MATERIALS
T +44 (0)20 8943 8480
E uksales@lgcstandards.com

PROFICIENCY TESTING
T +44 (0) 161 762 2500

19 USA + CANADA
T +1 (603) 622 7660
E lgcusa@lgcgroup.com

For full listing of offices and distributors
please see lgcstandards.com





Order online at lgcstandards.com
dr.ehrenstorfer@lgcgroup.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2017. All rights reserved.