



Vials

Content	Page
---------	------

Vials for Agilent 116-119

Wide Opening Crimp Top Vials & Crimp Caps (12 x 32 mm)	116
Crimp Seals With Different Septum Materials	117
Wide Opening Snap It & Snap Caps (12 x 32 mm):	117
Closures Snap Caps for Snap-it Vials	118
Wide Opening Screw Thread Vials, Screw Caps and Septa for Agilent	118
Inserts for Vial of 2 mL Crimp Cap and Screw Cap	119
Kit for Micro Insert Crimp Cap	119
Kit for Micro Insert Screw Cap	119

Vials for Waters Alliance 120-121

Wide Opening Screw Thread Vials (6 mm neck opening, 9 mm Cap) 2 mL (12x32 mm), Caps and Septa	120
Wide Opening Screw Thread Vials, 2 mL (12x32 mm) (10-425) (6 mm neck opening, 10 mm Cap)	121
Wide Opening Screw Thread 4 mL Vial (15x45 mm)	
Screw Thread for Waters 48-Position	
Carrusel WISP Autosampler Vials	121
Shell Vials, Sep Caps, 8x40 mm for Waters-96	
Position Tray	121

Vials for Shimadzu 122

Standard Opening 12x32 mm Screw Thread Vials, Caps and Septa	122
Limited Volume Polypropylene and Glastic Vials with 8-425 Finish	122
Assembled Caps and Septa for 8-425 Convenient standard Opening Vials for Shimadzu	122
Wide Opening Screw Thread Vial, 2 mL, 12x32 for Shimadzu (10-425), and Screw Cap	122

Teknokroma Econo Vials 123-124

Teknokroma Econo Vials 12 x 32 mm	
Crimp Cap Screw Top	123

Teknokroma Cap & Screws Vials 124-125

Teknokroma Econo Head Space Vials	124
Crimpers, Decrimpers, and Decapping Pliers for 8 mm Top Vials	124
Crimpers, Decrimpers, and Decapping Pliers for 11 mm Top Vials	124
Crimpers, Decrimpers, and Decapping Pliers for 20 mm Top Vials	124
Teknokroma 20mm Crimp Top and Screw Sample Vials	125
Teknokroma Screw Thread Vials	125
Teknokroma Butyl Stopper	125
Teknokroma Aluminium Screw Caps	125
Teknokroma Solid Screws (without central hole)	125
Teknokroma Septum	125

Teknokroma Headspace Vials 126

Other Vials 127-137

Target® LoVial™ Wide Opening 12x32 Crimp Top Vials	127
Inserts for Lo Vial™ Crimp-Top & Snap-it Vials	128
Snap-It Vials , Crimp & Snap-It Seals	129
Target Screw Thread Vials, Caps & Septa (10-425)	130
Target Screw Thread Caps, Inserts & Vial Kits	131
Target DP™ Vials & Caps (9-425)	132
Target DP™ Caps & Inserts	133
4 mL 15x45 mm Screw Thread Vials, Caps, Septa & Inserts	134
4 mL Vials Kits, 1 mL Vials & Inserts	135
Sample Storage Vials & Closures	136
Sample Storage Vials Convenience Kits & EPA Vials	137

Vial Comparison Chart 138-139

Autosampler Vials, Closures and Inserts for HPLC and GC for AGILENT

Features:

- Wide opening vial is easier to fill
- Fit most autosamplers that accommodate the 12x32 mm vial format
- Superior quality 33 expansion borosilicate clear glass (Type 1, Class A) or 51A amber glass
- Uniformly flat bottom for security with inserts
- I-D™ vials feature write-on patch with graduations at 0.5, 1.0 and 1.5 mL.
- Target Colorband™ crimp top vials feature color-coded graduated write-on patch
- Vials with an 11 mm crimp finish require aluminium seal closures
- Crimp/snap vials can be used with either 11 mm aluminium seals or Snap-It™ snap caps

Wide Opening Crimp Top Vials & Crimp Caps (12 x 32 mm):



Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-1	5181-3375	Vial 2 ml Clear, Wide Opening Crimp Top	100
CC-C4011-1AP	5181-1210	11 mm Crimp Seal, with PTFE/Red Rubber, standar septum	100



CC-C4011-1W	5182-0543	Vial 2 ml Clear, Wide Opening Crimp Top write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
CC-C4011-2A	5181-1211	11 mm Crimp Seal, with red PTFE/White Silicone/Red PTFE, for aggressive solvents	100



CC-C4011-2		Vial 2 ml Amber, Wide Opening Crimp Top	100
CC-C4011-2W	5181-3376	Vial 2 ml Amber, Wide Opening Crimp Top write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100




CC-C4011-LV1W		MacroVial 250 µl Fused Insert Vial, Clear Glass I-D, Write-on patch.	100
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100



CC-03-FIV		300 µl Glasss insert, fused into a 2 ml Crimp Top Vial, Write-on patch.	500
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100




Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-LV2W		MacroVial 250 ul. Fused Insert Vial, Amber Glass I-D, write-on patch	100
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100



CC-C4011-10		Solid Glass MicroVial 400 µl	12
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100

Crimp Seals With Different Septum Materials



CC-C4011-3A		11 mmCrimp Seal, with Polypropylene, for Corrosive Solvents	100
CC-C4011-4A		11 mmCrimp Seal, with PTFE/White Silicone, for High Purity Septum	100
CC-C4011-5A		11 mm Crimp Seal, with Black Viton Elastomer, for ECD, High Chemical Resistance	100
CC-C4011-1AP		11 mm Crimp Seal, with PTFE/Red Rubber, Standard Septum	100

Wide Opening Snap It & Snap Caps (12 x 32 mm):

Snap-it Vials can be used with Snap Caps or aluminium seals closures



Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-5	5182-0544	Snap-It™ Vial 2 ml Clear, Wide Opening	100
CC-C4011-6		Snap-It™ Vial 2 ml Amber, Wide Opening	100



CC-C4011-5W	5182-0546	Snap-It™ Vial 2 ml Clear, Wide Opening. Write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
CC-C4011-51B	5182-3458	Blue Snap Cap with PTFE/Red Rubber	100



CC-C4011-6W	5182-0545	Snap-It™ Vial 2 ml Amber, Wide Opening. Write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
CC-C4011-54G	5182-3457	Green Snap Cap with red PTFE/white Silicone, High purity septum	100



Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-LV1		MacroVial 250 µl. Fused Insert Snap-It Vial, Clear Glass	100



CC-C4011-LV2		MacroVial 250 µl. Fused Insert Snap-It Vial, Amber Glass	100
--------------	--	--	-----

Closures Snap Caps for Snap-it Vials



CC-C4011-53R		Red Snap Cap with red PTFE/White Silicone/Red PTFE, for aggressive solvents, resists coring.	100
CC-C4011-55B		Blue Snap Cap with blue PTFE/White Silicone Pre-slit , reducing coring, needle damage, allows venting	100
CC-C4011-50		Clear Snap Cap integral molded polypropylene single use economy cap, not resealable	100
CC-C4011-52B		Blue Snap Cap, virgin PTFE, Highly inert, not resealable, single use	100

Wide Opening Screw Thread Vials, Screw Caps and Septa for Agilent



CC-C4000-1	5182-0714	Vial, 2 ml Clear, Wide Opening Screw Top.	100
CC-C4000-51G	5182-0718	Green Screw Cap with Ivory PTFE/Red Rubber.	100



CC-C4000-1W	5182-0715	Vial, 2 ml Clear, Wide Opening Screw Top write-on patch with graduation 0.5, 1.0 and 1.5 ml	100
CC-C4000-55B	5183-2076	Blue, Screw Cap with Blue PTFE/Silicone, Pre-slit	100



CC-C4000-2W	5182-0716	Vial, 2 ml Amber, Wide Opening Screw Top write-on patch with graduation 0.5, 1.0 and 1.5 ml	100
CC-C4000-53G	5182-0724	Green Screw Cap, with Red PTFE/White Silicone/Red PTFE.	100



Cat.No	P/N Manufacturer	Description	Pk
CC-C4000-9	5183-2030	High Recovery Vial, 1,5 ml, Clear Screw Top.	100
CC-C4000-54 R	5182-0722	Red Screw Cap with Red PTFE/White Silicone	100



CC-C4000-LV1	5183-2076	MacroVial , Fused Insert, 350 µl, Screw Top	100
CC-C4000-55B	5183-2076	Blue Screw Cap, with PTFE/White Silicone, Pre-slit	100



CC-C4000-V1	5184-3550	Micro-V™, Vial, Clear, 1.5 ml, Screw Top with 150 µl reservoir	100
CC-C4000-55B	5183-2076	Blue Screw Cap, with Blue PTFE/White Silicone, Pre-slit	100

Inserts for Vial of 2 mL Crimp Cap and Screw Cap



CC-C4010-630	5181-1270	300 µl Polyspring Glass Conical Insert, with graduation and polyspring support	100
CC-C4010-631	5181-3377	250 µl Flat Bottom Glass MicroSert.	500
CC-C4010-630P		300 µl Polyspring Conical Insert Polypropilene	100

Kit for Micro Insert Crimp Cap



TR-400013 + CC-4010-630 + TR-400012 Kit for Micro Insert Crimp Cap

Kit for Micro Insert Screw Cap



TR-400008 + CC-4010-630 + TR-400007B Kit for Micro Insert Screw Cap

Autosampler Vials, Closures and Inserts for HPLC WATERS ALLIANCE

Wide Opening Screw Thread Vials (6 mm neck opening, 9 mm Cap) 2 mL (12x32 mm), Caps and Septa



Cat.No	P/N Manufacturer	Description	Pk
CC-C4000-1	186000273	Vial, 2 ml, Clear, Wide Opening, Screw Top.	100
CC-C4000-54B	WAT058874+ WAT058875	Blue Screw Cap with Red PTFE/White Silicone	100



CC-C4000-1W	186000845	Vial, 2 ml Clear, Wide Opening Screw Top write-on patch with graduation 0.5, 1.0 and 1.5 ml	100
CC-C4000-55B	186000305	Blue Screw Cap, with Blue PTFE/White Silicone, Pre-slit	100



CC-C4000-LV2	186001130	Microvial, Fused Insert Amber 350 µl, Screw Top	100
CC-C4000-LV1	186001128	Macroval , Fused Insert, 350 µl, Screw Top	100
CC-C4000-55B		Blue Screw Cap with Blue PTFE/White Silicone Pre-slit	100
CC-C4000-54B	WAT058874+ WAT058875	Blue Screw Cap with Red PTFE/White Silicone	100



CC-C4000-V1		Micro-V, Vial, Clear, 1.5 ml, Screw Top.	100
CC-C4000-55B	186000305	Blue Screw Cap, with Blue PTFE/White Silicone, Pre-slit	100



CC-C4010-630	WAT094170	300 µl Conical Glass Insert	100
CC-C4010-630P		300 µl Conical Insert, Polypropylene with graduation with polyspring support	100
CC-C4010-627L	WAT072704	350 µl Conical Glass Insert Pulled Point	100
CC-C4011-631	WAT072704	400 µl Glass Flat Bottom Insert	500

CC-C4010-630, 630P, 627L, 631

Wide Opening Screw Thread Vials, 2 mL (12x32 mm) (10-425) (6 mm neck opening, 10 mm Cap)



Cat.No	P/N Manufacturer	Description	Pk
CC-C4010-1	WAT063300	Vial, 2 ml (10-425) 12 x 32 mm, Clear, Wide Opening, Screw thread	100
CC-C4010-55BLK		Black Screw Cap (10-425) with Bue PTFE/Silicone, Pre-slit	100



CC-C4010-1W	186000845	Vial, 2 ml (10-425) 12 x 32 mm Clear, Wide Opening Screw Top write-on patch with graduation 0.5, 1.0 and 1.5 ml	100
CC-C4010-60BLK	WAT058875 + WAT058874	Black Screw Cap (10-425) with Red PTFE/ White Silicone	100




CC-C4010-2W		Vial, 2 ml (10-425) 12 x 32 mm Amber, Wide Opening Screw Top write-on patch with graduation 0.5, 1.0 and 1.5 ml	100
CC-C4010-55BLK		Black Screw Cap (10-425) with Bue PTFE/ White Silicone, Pre-slit	100

Wide Opening Screw Thread 4 mL Vial (15x45 mm) Screw Thread for Waters 48-Position Carrusel WISP Autosampler Vials



CC-C4015-11W	186000840	Vial, 4 ml (15x45 mm) Clear, Screw Top write-on patch with graduation 1.0, 2.0, 3.0 and 4.0 ml	100
CC-C4015-2W	186001135	Vial, 4 ml (15x45 mm) Amber, Screw Top write-on patch	100
CC-C4015-75A	186000841	Black Screw Cap (15x45) with Red PTFE/White Silicone	100
CC-C4015-55BLK	186000842	Black Screw Cap (15x45), Blue PTFE/White Silicone Pre-slit	100

Shell Vials, Sep Caps, 8x40 mm for Waters-96 Position Tray



CC-C4015-96	WAT025054	1 ml (8x40 mm) Clear Glass Shell vial with SepCap	200
CC-C4015-99	WAT025053	1ml (8x40 mm) AmberGlass Shell vial with SepCap	200
CC-C4015-95P		1ml (8x40 mm) Polypropylene Shell vial with SepCap	250

Autosampler Vials, Closures and Inserts for SHIMADZU (8-425)

Standard Opening 12x32 mm Screw Thread Vials, Caps and Septa

- Standard neck vial 8-425 thread finish, 12x32 mm outer profile
- Superior quality 33 expansion borosilicate clear glass (Type 1, Class A) or 51A amber glass
- Uniformity flat bottom for security inserts
- I-D™ vials feature write-on patch with graduations at 0.5, 1.0 and 1.5 mL
- Deactivated vials are Silanized or Kimshield™ pretreated



CC-C4013-12



CC-C4013-2W

CC-C4013-1W

Cat.No	Max. Fill volume	Description	Pk
CC-C4013-1	1.8 mL	2 mL Vial Clear Glass	100
CC-C4013-1W	1.8 mL	2 mL Vial Clear Glass I-D w/patch	100
CC-C4013-2	1.8 mL	2 mL Vial Amber Glass	100
CC-C4013-2W	1.8 mL	2 mL Vial Amber Glass I-D w/patch	100
CC-C4013-12	150 µL	150 µL MicroVial™ Clear Solid Glass	100

Limited Volume Polypropylene and Glastic Vials with 8-425 Finish

- Polypropylene vials are an economical alternative to glass for some applications
- Precision-formed tips are designed to aid sample recovery for limited volume applications
- Glass TPX vial features pre-inserted glass insert



CC-C4013-11,-13



CC-C4013-15GST

Cat.No	Max. Fill volume	Description	Pk
CC-C4013-11	250 µL	250 µL Virgin Polypropylene MicroVial™	100
CC-C4013-13	650 µL	750 µL Virgin Polypropylene MicroVial™	1000
CC-C4013-15GST	100 µL	150 µL Clasic Glass Insert TPX MicroVial™	100

Assembled Caps and Septa for 8-425 Convenient standard Opening Vials for Shimadzu


- Pre-assembled caps and septa are convenient and minimize contamination from handling
- Phenolic caps can be autoclaved
- Kim-Lok cap/septum assemblies feature a molded centering ring to eliminate septum movement during injection
- Storage cap has a solid top and pre-inserted PTFE disk liner



Kim-Lok™ Caps
CC-C4013-61Y

Cat.No	Cap Type	Description	Pk
CC-C4013-63A	Open Top cap w/flange	Black polypropylene cap, Red PTFE/White silicone	100
CC-C4013-64W	Open Top cap w/flange	White polypropylene cap, Blue PTFE/White silicone-Pre-slit	100
CC-C4013-61Y	Kim-Lok cap w/flange	Yellow polypropylene cap, Kim-Lok Red PTFE/White silicone	100

Wide Opening Screw Thread Vial, 2 mL, 12x32 for Shimadzu (10-425), and Screw Cap



CC-C4010-1, CC-C4010-1W



CC-C4015-96

CC-C4015-96A

CC-C4010-LV1

Cat.No	Max. Fill volume	Description	Pk
CC-C4010-1	1.8 mL	2 mL Clear Glass Vial (10-425)	100
CC-C4010-1W	1.8 mL	2 mL Clear Glass Vial w/patch (10-425)	100
CC-C4010-LV1	1.8 mL	350 µL Fused Insert	100
CC-4010-35W		Polypropylene Open Top cap, White cap, Red PTFE/White Silicone for Shimadzu	100
CC-4010-35BLK		Polypropylene Open Top cap, Black cap, Red PTFE/White Silicone for Shimadzu	100
CC-C4015-96A	150 µL	Polyspring Borosilicate Glass Inserts, Conical with pulled-point interior.	100
CC-C4015-96	1 mL	Clear Borosilicate glass vials with SepCaps	200

Teknokroma Vials

- Suitable for almost all automatic injectors in the market
- 2 ml Crimp Top vial and Screw Cap vials
- Crimp Seals and Screw Caps (9-425)
- An economical alternative developed by Teknokroma™

Teknokroma Econo Vials 12 x 32 mm Crimp Cap Screw Top



Cat.No	Description	Pk
TR-400013	Vial, 2 ml Clear, Wide Opening Crimp Top Vial	100
TR-400012	11 mm Crimp Seal, with PTFE/Red Rubber, standard septum	100
TR-400021	11 mm Crimp Seal, with PTFE/Silicone septum	100



TR-400010	Vial 2 ml Clear, Wide Opening Crimp Top write-on patch	100
TR-400012	11 mm Crimp Seal, with PTFE/Red Rubber, standard septum	100
TR-400021	11 mm Crimp Seal, with PTFE/Silicone septum	100



TR-400011	Vial 2 ml Amber, Wide Opening Crimp Top write-on patch	100
TR-400012	11 mm Crimp Seal, with PTFE/Red Rubber, standard septum	100
TR-400021	11 mm Crimp Seal, with PTFE/Silicone septum	100



TR-400014	Vial 2 ml Clear (9-425), Wide Opening Screw Top	100
TR-400007B	Blue screw cap (9-425) with bonded PTFE/Silicone for wide opening vial	100
TR-400020	Blue screw cap (9-425) with bonded Blue PTFE/Silicone for wide opening vial Pre-Slit	100
TR-400025	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone pk	100
TR-400026	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone Pre-Slit	100



TR-400008	Vial 2 ml Clear (9-425), Wide Opening Screw Top write-on patch	100
TR-400007B	Blue screw cap (9-425) with bonded PTFE/Silicone for wide opening vial	100
TR-400020	Blue screw cap (9-425) with bonded Blue PTFE/Silicone for wide opening vial Pre-Slit	100
TR-400025	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone pk	100
TR-400026	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone Pre-Slit	100



TR-400003	Vial 2 ml Amber (9-425), Wide Opening Screw Top write-on patch	100
TR-400007B	Blue screw cap (9-425) with blue PTFE/Silicone for wide opening vial	100
TR-400020	Blue screw cap (9-425) with bonded Blue PTFE/Silicone for wide opening vial Pre-Slit	100
TR-400025	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone pk	100
TR-400026	Blue screw(9-425) Cap with Ring with Blue PTFE/Silicone Pre-Slit	100

Tk Teknokroma Head Space Vials



Teknokroma Econo Head Space Vials


Cat.No	Description	Pk
TR-400016	Rounded Bottom 10 ml 23 x 46 mm	238
TR-400017	Rounded Bottom 20 ml 23 x 75 mm	238
TR-400023	Rounded Bottom 10 ml 23 x 46 mm	100
TR-400024	Rounded Bottom 20 ml 23 x 75 mm	100

Crimpers, Decrimpers, and Decapping Pliers for 8 mm Top Vials



Cat.No	Description	Use	Pk
CC-C4008-100	Manual Crimper	Attaches 8 mm aluminium crimp seals	1
CC-C4008-101	Decapping Pliers	Removes 8 mm aluminium crimp seals	1
CC-C4008-102	Manual Decrimper	Removes 8 mm aluminium seals without vial damage	1

Crimpers, Decrimpers, and Decapping Pliers for 11 mm Top Vials



Cat.No	Description	Use	Pk
CC-C4012-100	Manual Crimper	Attaches 11 mm aluminium crimp seals	1
CC-C4012-101	Decapping Pliers	Removes 11 mm aluminium crimp seals	1
CC-C4012-102	Manual Decrimper	Removes 11 mm aluminium seals without vial damage	1

Crimpers, Decrimpers, and Decapping Pliers for 20 mm Top Vials



Cat.No	Description	Use	Pk
CC-C4020-100	Manual Crimper	Attaches 20 mm aluminium crimp seals	1
CC-C4020-101	Decapping Pliers	Removes 20 mm aluminium crimp seals	1
CC-C4020-102	Manual Decrimper	Removes 20 mm aluminium seals without vial damage	1



Teknokroma 20mm Crimp Top and Screw Sample Vials

Cat.No	Vial Description	Volume/size	Pk
TR-000005	Clear, Crimp Top Vial	8 mL	100
TR-000004	Clear, Crimp Top Vial	12 mL	100
TR-000003	Clear, Crimp Top Vial	25 mL	100
TR-000001	Clear, Crimp Top Vial	60 mL	90
TR-000008	Amber, Crimp Top Vial	8 mL	100
TR-000007	Amber, Crimp Top Vial	12 mL	100
TR-000006	Amber, Crimp Top Vial	25 mL	144



Teknokroma Screw Thread Vials



Cat.No	Vial Description	Volume/size	Pk
TR-000013	Clear, Screw Thread Vial	12 mL	100
TR-000026	Clear, Screw Thread Vial	25 mL	100
TR-000012	Amber, Screw Thread Vial	8 mL	188
TR-000009	Amber, Screw Thread Vial	50 mL	78

Teknokroma Butil Stopper



Cat.No	Vial Description	Volume/size	Pk
TR-000029	Butil Stopper	13 mm	100
TR-000019	Butil Stopper	20 mm	500
TR-000031	Butil Stopper for Liofilized	20 mm	100



Teknokroma Aluminium Screw Caps

Cat.No	Vial Description	Volume/size	Pk
TR-000017	Aluminium Crimp Seal only	20 mm	500
TR-000027	Aluminium Crimp Seal only	13 mm	100




Teknokroma Solid Screws (without central hole)

Cat.No	Vial Description	Volume/size	Pk
TR-000021	Bakelite Screw Cap (without septa)	16 mm	100


Teknokroma Septum

Cat.No	Vial Description	Volume/size	Pk
TR-000018	Teflon silicone	20 mm	100



Cat.No	P/N Manufacturer	Description	Pk
CC-C4020-10	5182-0838	10 ml Vial, Flat Bottom 23 x 46 mm, Beveled Edge	100
CC-C4020-410		10 ml Vial, Flat Bottom 23 x 46 mm, Square Rim	1000
CC-C4020-210	5183-4475	10 ml Vial, Round Bottom 23 x 46 mm, Beveled Edge	100
CC-C4020-20		20 ml Vial, Flat Bottom 23 x 75 mm, Beveled Edge	100
CC-C4020-25		20 ml Vial, Flat Bottom 23 x 75 mm, Square Rim	1000
CC-C4020-2		20 ml Vial, Round Bottom 23 x 75 mm, Beveled Edge	100

20 mm Crimp Seals with Prefitted Septa for Headspace Vials



Cat.No	Seal Type	Septa	Pk
CC-C4020-34A	Regular	Ivory PTFE/Red Rubber, (-40 to +100°C)	100
CC-C4020-34AP	Pressure Release	Ivory PTFE/Red Rubber, (-40 to +100°C)	100
CC-C4020-39A	Regular	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-43A	Magnetic Steel	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-43AP	Pressure Release	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-36A	Regular	Gray PTFE/Black Molded Butyl (-40 to +125°C)	100
CC-C4020-36AP	Pressure Release	Gray PTFE/Black Molded Butyl (-40 to +125°C)	100
CC-C4020-32A	Regular	Tan PTFE/White Silicone (-60 to +200°C)	100
CC-C4020-32AP	Pressure Release	Tan PTFE/White Silicone (-60 to +200°C)	100
CC-C4020-42A	Magnetic Steel	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100
CC-C4020-42AP	Pressure Release	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100
CC-C4020-37AP	Pressure Release	Aluminium Foil/White High Temp. Silicone (-60 to +220°C)	100
CC-C4020-31	Molded Polypropylene Storage Cap	SepCap Integral Molded Polypropylene	250
CC-C4020-42A	Magnetic Steel	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100



Crimpers, Decrimpers, and Decapping Pliers for 20 mm Top Vials

Cat.No	Description	Use	Pk
CC-C4020-100	Manual Crimper	Attaches 20 mm aluminium crimp seals	1
CC-C4020-101	Decapping Pliers	Removes 20 mm aluminium crimp seals	1
CC-C4020-102	Manual Decrimper	Removes 20 mm aluminium seals without vial damage	1

Target® LoVial™ Wide Opening 12x32 Crimp Top Vials Tk



CC-C4011-1, 1W



CC-C4011-2W



CC-C4011-LV1



CC-C4011-1240



CC-C4011-LV1W-LV2W



CC-C4011-10

Target® LoVial™ Crimp Top Vials - Wide Opening 2 mL 12x 32

Cat.No	Vial Description	Maximum Fill Volume	Finish	Pk
CC-C4011-1	2 mL LoVial, 2 mL Clear Glass Vial	1.8 mL	Crimp	100
CC-C4011-1W	2 mL LoVial, 2 mL Clear Glass I-D Vial w/patch	1.8 mL	Crimp	100
CC-C4011-2	2 mL LoVial, 2 mL Amber Glass Vial	1.8 mL	Crimp	100
CC-C4011-2W	2 mL LoVial, 2 mL Amber Glass I-D Vial w/patch	1.8 mL	Crimp	100



CC-C4011-V2-V1

Microsampling Vials for use with 11 mm Aluminium Crimp Seals

Cat.No	Vial Description	Maximum Fill Volume	Finish	Pk
CC-C4011-LV1	MacroVial™ Fused Insert Snap-It Vial, Clear Glass	250 µL	Crimp/Snap	100
CC-C4011-LV1W	MacroVial™ Fused Insert Vial, Clear Glass I-D w/patch	250 µL	Crimp	100
CC-C4011-LV2	MacroVial™ Fused Insert Snap-It Vial, Amber Glass	250 µL	Crimp/Snap	100
CC-C4011-LV2W	MacroVial™ Fused Insert Vial, Amber Glass I-D w/patch	250 µL	Crimp	100
CC-C4011-10	Solid Glass MacroVial™	400 µL	Crimp	12
CC-C4011-V1	MicroV™ Clear Glass Microsampling Vial, 15 µL Reservoir	1.1 mL	Crimp	100
CC-C4011-V2	MicroV™ Amber Glass Microsampling Vial, 15 µL Reservoir	1.1 mL	Crimp	100
CC-C4011-9	1.5 mL High Recovery MicroVial	1.5 mL	Crimp	100
CC-C4011-1.7CB	1.7 mL Conical Bottom Crimp Top LoVial	1.1 mL	Crimp	1000
CC-C4011-1240	2.5 mL Clear Glass Crimp Top Vial 40 mm vial Height	2.3 mL	Crimp	1000
CC-C4011-1PT	900 µL Precision Taper Vial, Clear Glass	700 µL	Crimp	100



CC-C4011-1PT



CC-C4011-14

Deactivated Vials, Crimp Top

Cat.No	Vial Description	Maximum Fill Volume	Finish	Pk
CC-C4011-S1	2 mL Clear Glass Vial, NSC Silanized Crimp	1.8 mL	Crimp	100
CC-C4011-S1W	2 mL Clear Glass I-D Vial w/patch, NSC Silanized Crimp	1.8 mL	Crimp	100
CC-C4011-S2W	2 mL Amber Glass I-D Vial w/patch, NSC Silanized Crimp	1.8 mL	Crimp	100
CC-C4011-S9	1.5 mL High Recovery MicroVial, NSC Silanized Crimp	1.5 mL	Crimp	100
CC-C4011-K1W	2 mL Clear Glass I-D Vial w/patch, Kimshield™ Deactivated	1.8 mL	Crimp	100
CC-C4011-K2W	2 mL Amber Glass I-D Vial w/patch, Kimshield™ Deactivated	1.8 mL	Crimp	100



CC-C4011-13



CC-C4011-16



CC-C4011-11

Limited Volume Plastic Crimp Top Vials

Cat.No	Vial Description	Residual Volume	Finish	Pk
CC-C4011-13	250 µL Virgin Polypropylene MicroVial™, Crimp/Snap	3 µL	Crimp/Snap	100
CC-C4011-16	350 µL Virgin Polypropylene MicroVial™, Crimp/Snap	5 µL	Crimp/Snap	100
CC-C4011-11	600 µL Virgin Polypropylene MicroVial™, Crimp/Snap	6 µL	Crimp/Snap	100
CC-C4011-14	800 µL Virgin Polypropylene Flat Bottom Vial, Crimp/Snap	30 µL	Crimp/Snap	100

Aluminium Crimp Seals with Prefitted Septa for 12x32 mm Vials

Cat.No	Seal Description	Seal Color	Pk
CC-C4011-1AP	Clear PTFE/Red Rubber Standard GC septum	Silver	100
CC-C4011-2A	Red PTFE/White Silicone/Red PTFE for Aggressive Solvents	Silver	100
CC-C4011-3A	Polypropylene for Corrosive Solvents	Silver	100
CC-C4011-4A	Red PTFE/White Silicone, High Purity septum	Silver	100
CC-C4011-5A	Black Viton Elastomer, ECD, High Chemical Resistance	Silver	100



CC-C4011-2A, -3A, -4A, 5A

Tk Inserts for Lo Vial™ Crimp-Top & Snap-it Vials

Target Polyspring™ and Conical Inserts for Wide Opening Target vials



CC-C4010-630 CC-C4010-630P



CC-C4010-629L



CC-C4011-631

- These inserts all have a nominal O.D. of 6mm to fit Target wide opening vials
- Polyspring® inserts feature a polyethylene spring to center the bottom of the insert and to provide a cushion against needle contact.
- Precision point inserts minimize residual sample losses
- Deactivated inserts are Silinized or Kimshield™ pretreated

Cat.No	Insert Description	Pk
CC-C4010-630	Target Polyspring Borosilicate Glass Inserts. Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity.	100
CC-C4010-S630	Silinized Target Polyspring Borosilicate Glass Inserts Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity. Silinized to close off all active sites.	100
CC-C4010-630P	Target Polyspring Polypropylene Inserts. Conical with precision-formed interior and attached Polyspring. 250 µL capacity. Graduated volume marks.	100
CC-C4010-627L	Target Borosilicate Glass Inserts. Conical with pulled-point interior. 300 µL capacity.	
CC-C4010-629L	Target Borosilicate Glass Inserts. Conical with precision-formed mandrel interior. 300 µL capacity.	100
CC-C4010-S629	Silinized Target Borosilicate Glass Inserts. Conical with precision-formed mandrel interior. 300 µL capacity. Silinized to close off all active sites.	100

Target MicroSerts™ - Flat Bottom Inserts

- An economical, limited volume alternative
- Custom fit to Target vials, extremely flat bottom
- Won't push-up, protects syringe needle
- Easy, open-end identification from our exclusive coloured line at top of insert.

Cat.No	Vial Description	Maximum Fill Volume	Pk
CC-C4011-631	Target Borosilicate Glass MicroSerts. Flat Bottom, 400 µL capacity.		500
CC-C4011-S631	Silinized Target Borosilicate Glass MicroSerts. Flat Bottom, 400 µL capacity. Silinized to close off all active sites.		500

11 mm Crimping & Decapping Tools


- Use with any Target or standard-opening 2 mL, 12 x 32 mm crimp top vials



Cat.No	Description	Pk
CC-C4012-100	Manual Crimper, 11 mm	1
CC-C4012-101	Decapping Pliers, 11 mm	1
CC-C4012-102	Manual Decrimper, 11 mm	1

Vial Racks for all 2 mL, 12 x 32 mm Vials

- Polypropylene vial racks are resistant to most solvents
- Racks feature alphanumeric indexing for easier vial identification
- Racks can be stacked for efficient storage



Cat.No 1 pk	Description	Pk
CC-C4011-25	Polystyrene 2 mL, Storage Rack with Clear Lid, 100 holes, (10x10) with alpha-numeric indexing, with cover	1
CC-C4012-25	Polypropylene 2 mL, Vial Rack, 50 holes (5x10), with alpha-numeric indexing	1



CC-C4011-5W



CC-C4011-6W



CC-C4011-LV1



CC-C4011-V5



CC-C4011-4

Target™ Snap-It Glass Vials, Wide Opening 2 mL, 12x32

- Wide-opening 2 mL, 12 x 32 mm vial
- Borosilicate glass
- Meet extremely tight tolerances
- I-D™ vial with graduated write-on patch measures 0.5, 1.0 and 1.5 mL

Cat.No	Target Snap-It Vial Description	Pk
CC-C4011-5	Clear	100
CC-C4011-5W	Clear I-D (graduated write-on patch)	100
CC-C4011-6	Amber I-D	100
CC-C4011-6W	Amber I-D (graduated write-on patch)	100

Microsampling Vials for Use with Snap-It Snap Cap.

- Wide-opening, 12 x 32 mm vials
- Borosilicate glass
- Sample concentrations and injections in one vial, no transferring to limited volume inserts.

Cat.No	Target Snap-It Vial Description	Pk
CC-C4011-4	Snap-It High Recovery Glass Vial, Clear glass, 1.5 mL with 30 µL reservoir	100
CC-C4011-LV1	Macroval™, Snap-It/Crimp, Clear glass, Fused 350 µL, Insert Vial	100
CC-C4011-LV2	Macroval™, Snap-It/Crimp, Amber glass, Fused 350 µL, Insert Vial	100
CC-C4011-V5	Micro-V™, Snap-It/Crimp Top, Clear glass microvial, 1.5 mL capacity, 15 µL reservoir	100
CC-C4011-V6	Micro-V™, Snap-It/Crimp Top, Amber glass microvial, 1.5 mL capacity, 15 µL reservoir	100
CC-C4011-11	Virgin Polypropylene, Microvial, Crimp/Snap-It, 450 µL capacity, 6 µL residual	100
CC-C4011-14	Virgin Polypropylene, Microvial, Crimp/Snap-It, 700 µL capacity, Flat bottom, 30 µl residual	100

Target Crim & Snap-It Seals



Seals 11 mm Crimp Seals with Septa for all 2 mL, 12x32 mm Crimp Top Vials

Cat.No	Seal Colors	Septa	Pk
CC-C4011-1AP	Silver	Standard, Clear PTFE/Red Rubber	100

11 mL Snap-It™ Seals with Septa for all 2 mL, 12x32 mm Snap-It Vials

- Eliminates tedious crimping
- Effective Evaporation resistant closure - apply with simple thumb pressure

Cat.No	Seal Colors	Septa	Pk
CC-C4011-50	Clear	Integral Molded Polyethylene, Single Use Economy	100
CC-C4011-51B	Blue	Natural Rubber Septa, PTFE/Red Rubber	100
CC-C4011-52	Clear	Highly Inert Virgin PTFE	100
CC-C4011-53B	Blue	Red PTFE/White Silicone/Red PTFE for Aggressive Solvents	100
CC-C4011-54G	Green	High Purity Septum, Red PTFE/White Silicone	100
CC-C4011-55B	Blue	Pre-Slit Blue PTFE/White Silicone	100

Tk Target Screw Thread Vials, Caps & Septa (10-425)



CC-C4010-2 CC-C4010-2W

Target Screw Thread Vials (10-425) wide opening 2 mL, 12 x 32 mm

- Wide neck accommodates wide 10-425 Thread cap
- Borosilicate glass
- Uniformly flat bottom for security with inserts
- I-D™ vial with write-on patch measures 0.5, 1.0 and 1.5 mL
- **Vial used for Shimadzu and Varian.**

Cat.No	Target Screw Thread Vial Description	Pk
CC-C4010-1	Clear Glass 2 mL	100
CC-C4010-1W	Clear I-D (graduated write-on patch) Glass 2 mL	100
CC-C4010-2	Amber Glass 2 mL	100
CC-C4010-2W	Amber I-D (graduated write-on patch) Glass 2 mL	100
CC-C4010-SI	Silanized, Clear I-D Glass 2 mL	100
CC-C4010-S1W	Silanized, Clear I-D (graduated write-on patch) Glass 2 mL	100
CC-C4010-S2W	Silanized, Amber I-D (graduated write-on patch) Glass 2 mL	100



CC-C4010-11 CC-C4010-14

Target Microsampling Screw Thread Vials

- Wide-opening, 12 x 32 mm vials. Accommodates wide 10-425 Thread Cap.
- Sample concentrations and injections in one vial - no transferring to limited volume inserts
- Borosilicate glass or virgin polypropylene
- Glastic™ vial has clear plastic TPX body with glass flanged insert molded to interior.



CC-C4010-LV1

Cat.No	Seal Colours	Septa	Pk
CC-C4010-11	Virgin Polypropylene, Microvial 375 µL capacity		100
CC-C4010-14	Virgin Polypropylene, Microvial 550 µL capacity, High Recovery Vial		100
CC-C4010-V1	Micro-V™, Clear glass microvial, 1.5 mL capacity, 150 µL reservoir		100
CC-C4010-LV1	Macroval™, Clear glass, 350 µL fused Insert		100
CC-C4010-V2	MicroV™, 1.5 mL. Tapered Microvial with 150 µL reservoir, Amber glass,		100
CC-C4010-LV2	Macroval™, 350 µL Fused Insert Vial, Amber glass		100



CC-C4010-V1

Target 10 mm Screw with Septa (10-425)



- Pre-assembled caps with septa
- 10-425 thread finish

Cat.No	Cap Colours	Septa	Pk
CC-C4010-30A	Light Blue	Ivory PTFE/Red Rubber	100
CC-C4010-35W*	White	Shimadzu Red PTFE/White Silicone	100
CC-C4010-35BLK*	Black	Shimadzu Red PTFE/White Silicone	100
CC-C4010-40A	Light Blue	Red PTFE/White Silicone/Red PTFE	100
CC-C4010-55A	Light Blue	Pre-Slit Blue PTFE/White Silicone	100
CC-C4010-55BLK	Black	Pre-Slit Blue PTFE/White Silicone	100
CC-C4010-60A	Light blue	Red PTFE/White Silicone	100
CC-C4010-60BLK	Black	Red PTFE/White Silicone	100

* Screw Caps for Shimadzu

Target Screw Thread Caps, Inserts & Vial Kits Tk



CC-C4010-630 CC-C4010-630P CC-C4011-631

Target Polyspring™ and MicroSert™ Inserts

- Eliminates dead volume - with correct autosampler default setting
- Damage-free syringe needle hits
- Allows complete sample extraction
- Shock-absorbing Polyspring or economical MicroSert
- Self-aligning, stands straight
- Easy, open-end identification from our exclusive coloured line at top of MicroSert

Cat.No	Insert Description	Pk
CC-C4010-630	Target Polyspring Borosilicate Glass Inserts, Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity.	100
CC-C4010-S630	Silanized Target Polyspring Borosilicate Glass Inserts, Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity. Silanized to close off all active sites.	100
CC-C4010-630P	Target Polyspring Polypropylene Inserts, Conical with precision-formed interior and attached Polyspring. 250 µL capacity. Graduated volume marks.	100
CC-C4010-631	Target Borosilicate Glass MicroSerts. Flat Bottom 400 µL capacity.	100



Pre-Assembled Kits - Target Screw Thread (10-425) - 12x32 mm in Storage Rack

- Assembled kits include 100 vials with pre-attached caps and septa
- Packed in convenient vial trays with clear covers
- Vial trays are stackable for easy storage
- Vial trays hold vials upright for easy filling and maintaining sample order

Cat.No	Target Screw Thread Vials	Cap Colours	Septa Cat.Nbr.	Septa	Pk
CC-C4010-17	CC-C4010-1, Clear Glass	Light Blue	CC-C4010-60A	PTFE/Silicone	100
CC-C4010-17A	CC-C4010-2, Amber Glass	Light Blue	CC-C4010-60A	PTFE/Silicone	100
CC-C4010-17AW	CC-C4010-2W, Amber Glass I-D™	Light Blue	CC-C4010-60A	PTFE/Silicone	100
CC-C4010-17BLK	CC-C4010-1, Clear Glass	Black	CC-C4010-60BLK	PTFE/Silicone	100
CC-C4010-17ABLK	CC-C4010-2, Amber Glass	Black	CC-C4010-60BLK	PTFE/Silicone	100
CC-C4010-17W	CC-C4010-1W, Clear Glass I-D™	White	CC-C4010-60A	PTFE/Silicone	100
CC-C4010-19	CC-C4010-1, Clear Glass	Light Blue	CC-C4010-40A	PTFE/Silicone/PTFE	100
CC-C4010-19A	CC-C4010-2, Amber Glass	Light Blue	CC-C4010-40A	PTFE/Silicone/PTFE	100
CC-C4010-19AW	CC-C4010-2W, Amber Glass I-D™	Light Blue	CC-C4010-40A	PTFE/Silicone/PTFE	100
CC-C4010-19W	CC-C4010-1W, Clear Glass I-D™	Light Blue	CC-C4010-40A	PTFE/Silicone/PTFE	100
CC-C4010-S19W	CC-C4010-S1W, Clear I-D™ Silinized	Light Blue	CC-C4010-40A	PTFE/Silicone/PTFE	100
CC-C4010-S20W	CC-C4010-S1W, Clear I-D™ Silinized	Light Blue	CC-C4010-55A	PTFE/Silicone-Pre-Slit	100
CC-C4010-21	CC-C4010-1, Clear Glass	Light Blue	CC-C4010-1A/CC-C4010-10	Solid PTFE	100
CC-C4010-21A	CC-C4010-2, Amber Glass	Light Blue	CC-C4010-1A/CC-C4010-10	Solid PTFE	100
CC-C4010-21W	CC-C4010-1W, Clear Glass I-D™	Light Blue	CC-C4010-1A/CC-C4010-10	Solid PTFE	100
CC-C4010-21AW	CC-C4010-2W, Amber Glass	Light Blue	CC-C4010-1A/CC-C4010-10	Solid PTFE	100
CC-C4010-24W	CC-C4010-1W, Clear Glass I-D™	Light Blue	CC-C4010-65A	PTFE/Silicone-Start-Slit	100
CC-C4010-24AW	CC-C4010-2W, Amber Glass I-D™	Light Blue	CC-C4010-65A	PTFE/Silicone-Start-Slit	100
CC-C4010-57	CC-C4010-1, Clear Glass	White	CC-C4010-35W	PTFE/Silicone, Shimadzu®	100
CC-C4010-57A	CC-C4010-2, Amber Glass	White	CC-C4010-35W	PTFE/Silicone, Shimadzu®	100
CC-C4010-57W	CC-C4010-1W, Clear Glass I-D™	White	CC-C4010-35W	PTFE/Silicone, Shimadzu®	100
CC-C4010-57AW	CC-C4010-2W, Amber Glass I-D™	White	CC-C4010-35W	PTFE/Silicone, Shimadzu®	100
CC-C4010-67	CC-C4010-1, Clear Glass	Black	CC-C4010-35BLK	PTFE/Silicone, Shimadzu®	100
CC-C4010-67A	CC-C4010-2, Amber Glass	Black	CC-C4010-35BLK	PTFE/Silicone, Shimadzu®	100
CC-C4010-67W	CC-C4010-1W, Clear Glass I-D™	Black	CC-C4010-35BLK	PTFE/Silicone, Shimadzu®	100
CC-C4010-67AW	CC-C4010-2W, Amber Glass I-D™	Black	CC-C4010-35BLK	PTFE/Silicone, Shimadzu®	100

Tk Target DP™ Vials & Caps (9-425)



CC-C4000-1



CC-C4000-LV1

CC-C4000-2W



CC-C4000-V2



CC-C4000-9



Target DP™ Glass Vials 12x32 mm Screw Thread Glass Vials (9-425)

- Wide Opening 2 mL, 12 x 32 mm vial
- Borosilicate glass
- Fits rotating or robotic arm autosampler trays
- Easy microsample testing with our Microvolume Inserts

Cat.No	Target DP Vial Description	Pk
CC-C4000-1	Clear Glass 2 mL	100
CC-C4000-1W	Clear Glass I-D (graduated write-on patch) 2 mL	100
CC-C4000-2W	Amber Glass I-D (graduated write-on patch) 2 mL	100
CC-C4000-S1	Silanized, Clear Glass 2 mL	100
CC-C4000-S1W	Silanized, Clear Glass I-D (graduated write-on patch) 2 mL	100
CC-C4000-S2W	Silanized, Amber Glass I-D (graduated write-on patch) 2 mL	100

Target DP™ Microsampling Vials

- Wide Opening 12 x 32 mm vials
- Sample concentrations and injections in one vial - no transferring to limited volume inserts
- Borosilicate glass

Cat.No	Target DP Vial Description	Pk
CC-C4000-9	High Recovery, Clear glass, 1.5 mL with 30 µL reservoir	100
CC-C4000-S9	Silanized, High Recovery, Clear glass, 1.5 mL with 30 µL reservoir	100
CC-C4000-9A	High Recovery, Amber I-D glass, 1.5 mL with 30 µL reservoir	100
CC-C4000-LV1	Macroval™, Clear glass, Fused 350 µL Insert vial	100
CC-C4000-LV2	Macroval™, Amber glass, Fused 350 µL Insert vial	100
CC-C4000-V1	Micro-V™, Clear glass microvial, 1.5 mL capacity, 150 µL reservoir	100
CC-C4000-V2	Micro-V™, Amber glass microvial, 1.5 mL capacity, 150 µL reservoir	100

Target DP™ Polypropylene Caps (no septa)

- An economical choice - purchase caps only and replace septa when needed
- Chose from our five colours for separating assays, lots, and standards

Cat.No	Cap Colour	Pk
CC-C4000-98B	Blue (open top)	100

Target DP™ Septa

Cat.No	Septa Description	Pk
CC-C4000-30	Ivory PTFE/Red Rubber	100
CC-C4000-40	Red PTFE/White Silicone/Red PTFE for use with aggressive solvents	100
CC-C4000-55	Pre-Slit Blue PTFE/White Silicone	100
CC-C4000-57	Pre-Slit Ivory PTFE/Red Rubber	100
CC-C4000-60	Red PTFE/White Silicone for high sensibility chromatography	100





Assembled Target DP Caps with Septa (9-425)

- Convenient screw thread cap with profile design of crimp seal
- No crimping or decapping tools needed
- Easy-on, easy-off convenience with just half-a-turn
- Precision positioned septa form evaporation-proof seal
- Full compatibility with all Target DP components

Cat.No	Cap Colours	Septa	Pk
CC-C4000-51G	Green	Ivory PTFE/Red Rubber	100
CC-C4000-53B	Blue	Red PTFE/White Silicone/Red PTFE	100
CC-C4000-54R	Red	Red PTFE/White Silicone	100
CC-C4000-55B	Blue	Pre-Slit Blue PTFE/White Silicone	100

Target Polyspring™ Inserts



CC-C4010-630 CC-C4010-630P

- Eliminates dead volume - with correct autosampler default setting
- Damage-free syringe needle injections
- Permits maximum sample extraction
- Shock-absorbing Polyspring
- Self-aligning, stands straight

Cat.No	Insert Description	Pk
CC-C4010-630	Target Polyspring Borosilicate Glass Inserts, Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity.	100
CC-C4010-S630	Silanized Target Polyspring Borosilicate Glass Inserts, Conical with precision-formed mandrel interior and attached Polyspring. 250 µL capacity. Silanized to close off all active sites.	100
CC-C4010-630P	Target Polyspring Polypropylene Inserts, Conical with precision-formed interior and attached Polyspring. 250 µL capacity. Graduated volume marks.	100



CC-C4011-631

Target MicroSerts™ - Flat Bottom Inserts

- An economical, limited volume alternative
- Custom fit to Target vials - extremely flat bottom
- Easy, open-end identification from our exclusive coloured line at top of insert

Cat.No	Septa Description	Pk
CC-C4011-631	Target Borosilicate Glass MicroSerts. Flat Bottom. 400 µL capacity.	100
CC-C4011-S631	Silanized Target Borosilicate Glass MicroSerts. Flat Bottom. 400 µL capacity.	100

Tk 4 mL 15x45 mm Screw Thread Vials, Caps, Septa & Inserts



CC-C4015-11W CC-C4015-2W

Screw Thread 4 mL Glass Vials 15x45 mm.

- 4 mL, 45 mm vial, 15-425 Thread
- Borosilicate glass

Cat.No	4 mL Screw Thread Vial Description	Pk
CC-C4015-1	15x45mm Clear glass	100
CC-C4015-11W	15x45mm Clear glass I-D™ (graduated write-on patch)	100
CC-C4015-2	15x45mm Amber glass	100
CC-C4015-2W	15x45mm Amber glass I-D™ (graduated write-on patch)	100
CC-C4015-S1	15x45mm Silanized, Clear glass	100
CC-C4015-S2	15x45mm Silanized, Amber glass	100



Screw Caps for 4 mL Thread Vials

- CC-C4015-5A SepCap has integral septa made of 0.25 mm polypropylene - no septa to order
- CC-C4015-30A, CC-C4015-40A, CC-C4015-75A, CC-C4015-45W, are pre-assembled with septa


Cat.No	Cap & Septa Description	Pk
CC-C4015-1A	Black Open Top Screw Cap (no septa)	100
CC-C4015-1W	White Open Top Screw Cap (no septa)	100
CC-C4015-5A	SepCaps with 0.25 mm Polypropylene Septa	100
CC-C4015-30A	Black Screw Cap with Ivory PTFE/Red Rubber	100
CC-C4015-40A	Black Screw Cap with Red PTFE/White Silicone/Red PTFE	100
CC-C4015-75A	Black Screw Cap with Red PTFE/White Silicone, 1.90 mm	100
CC-C4015-45W	White Screw Cap with Red PTFE/White Silicone, 1.15 mm	100
CC-C4015-55BLK	Black Screw Open Top Cap, Polypropylene Septa, Blue PTFE/White Silicone Pre-Slit	100



13 mm Septa for 4 mL Screw Caps

Cat.No	Septa Description	Pk
CC-C4015-10	White Virgin PTFE, 0.25 mm, single use	100
CC-C4015-45	Red PTFE, White Silicone, 1.15 mm for thin gauge needles	100
CC-C4015-30	Ivory PTFE, Red Rubber for general chromatography	100
CC-C4015-40	Red PTFE, White Silicone/Red PTFE for aggressive solvents	100

Polyspring™ Glass Inserts - for 15x45, 4 mL Screw vials



CC-C4015-843

- Eliminates dead volume - with correct autosampler default setting
- Damage-free syringe injections
- Permits maximum sample extractions
- Shock-absorbing Polyspring
- Self-aligning, stands straight

Cat.No	Insert Description	Pk
CC-C4015-638	Polyspring Borosilicate Glass Inserts, Conical with pulled-point interior 700 µL capacity.	100
CC-C4015-641	Polyspring Borosilicate Glass Inserts, Conical with pulled-point mandrel interior. Flange at top of insert. 500 µL capacity.	100
CC-C4015-843	MicroSert™, Glass Flat Bottom Insert 1000 µL capacity.	100
CC-C4015-S641	Silanized Borosilicate Polyspring Glass Inserts, Conical with pulled-point mandrel interior. Flange at top of insert. 500 µL capacity.	100

4 mL Vials Kits, 1 mL Vials & Inserts Tk

Assembled 4 mL Screw Thread Vial Kits



- Assembled kits include 100 vials with pre-attached caps and septa
- Packaged in convenient vial trays with clear covers or in polybags
- Vial trays are stackable for easy storage and feature alphanumeric indexing
- Vial trays hold vials upright for easy filling and maintaining sample order

Cat.No	4 mL Vials	Cap Colour	Septa	Pk
CC-C4015-17	CC-C4015-1, Clear Glass	Black	CC-C4015-75A Red PTFE/White Silicone	100
CC-C4015-17AW	CC-C4015-2W, Amber Glass I-D	Black	CC-C4015-75A Red PTFE/White Silicone	100
CC-C4015-21	CC-C4015-1, Clear Glass	Black	CC-C4015-1A/4015-10 White Virgin PTFE	100
CC-C4015-21AW	CC-C4015-2AW, Amber Glass I-D	Black	CC-C4015-1A/4015-10 White Virgin PTFE	100



Unassembled 4mL Screw Thread Vial Convenience Kits

- Include 100 vials and 100 caps with pre-assembled septa
- All items remain handy and dusty-free
- Caps with pre-inserted septa are not attached to vials - packed in polybags.

Cat.No	4 mL Screw Thread Vials	Cap with Septa	Pk
CC-C4015-88	CC-C4015-1 Clear	CC-C4015-75A Black Cap with Red PTFE/White Silicone	100
CC-C4015-88AW	CC-C4015-2W Amber I-D	CC-C4015-75A Black Cap with Red PTFE/White Silicone	100

1 mL 8x40 mm SepCap Vials for Waters WISP™ autosamplers

- 1 mL, 8 x 40 mm
- For Waters 96 Position Tray
- SepCaps' starburst center design at 0.25 mm eases syringe needle penetration
- Ridge on SepCap underside centers limited volume insert
- SepCap design allows for easy push-in, pull-out from vial

Cat.No	1 mL Vial Description	Pk
CC-C4015-96	1 mL, Clear Borosilicate glass vials with SepCaps	200
*CC-C4015-99	1 mL, Amber Borosilicate glass vials with SepCaps	200
CC-C4015-95P	1 mL, Virgin Polypropylene SepCap vials with SepCaps	250

* Used by Shimadzu AOC-14/1400, LC-10A, SIL6B, SIL-94

Polyspring™ Glass Inserts for 8x40 mm - 1 mL SepCap Vials

- Eliminates dead volume - with correct autosampler default setting
- Damage-free syringe needle injections
- Permits maximum sample extraction
- Shock-absorbing Polyspring
- Self-aligning, stands straight

Cat.No	Insert Description	Pk
*CC-C4015-96A	Polyspring Borosilicate Glass Inserts, Conical with pulled-point interior. 150 µL capacity.	100
CC-C4015-S96A	Silanized Polyspring Borosilicate Glass Inserts, Conical with pulled-point interior. 150 µL capacity. Silanized to close off all active sites.	100
CC-C4015-96PA	Polyspring Polypropylene Inserts, Conical with precision-formed interior and attached Polyspring. 200 µL capacity. Guaranteed volume marks	100

* Used by Shimadzu AOC-14/1400, LC-10A, SIL6B, SIL-94



CC-C4015-96

CC-C4015-99



CC-C4015-96A

CC-C4015-96PA

Tk Sample Storage Vials & Closures



Sample Storage Vials

- Borosilicate glass
- Eliminates leaching of ions
- Provides consistent pH for duration of sample storage life

Cat.No	Glass	Capacity mL	Capacity Drams	Use Cap Size	Vial Size mm	Pk
CC-B7999-1	Clear	2.0	0.5	8-425	12 x 32	100
CC-B7999-2	Clear	4.0	1	13-425	15 x 45	100
CC-B7999-3	Clear	8.0	2	15-425	17 x 60	200
CC-B7999-12	Clear	12.0	3	15-425	19 x 65	200
CC-B7999-4	Clear	16.0	4	18-400	21 x 70	200
CC-B7999-5	Clear	22.0	6	20-400	23 x 85	200
CC-B7999-6	Clear	40.0	10	24-400	28 x 95	100
CC-B7999-1A	Amber	2.0	0.5	8-425	12 x 32	100
CC-B7999-2A	Amber	4.0	1	13-425	15 x 45	100
CC-B7999-3A	Amber	8.0	2	15-425	17 x 60	200
CC-B7999-12A	Amber	12.0	3	15-425	19 x 65	200
CC-B7999-4A	Amber	16.0	4	18-400	21 x 70	200
CC-B7999-6A	Amber	40.0	10	24-400	28 x 95	100



Screw Caps for Sample Storage Vials

- Polypropylene screw cap won't contaminate sample with syringe needle insertion - unlike phenolic caps that can flake
- Open top without septa

Cat.No	Cap Size	Description	Pk
CC-B7807-8	8-425	Black Open top Cap, Polypropylene	100
CC-B7807-13	13-425	Black Open top Cap, Polypropylene	100
CC-B7807-15	15-425	Black Open top Cap, Polypropylene	100
CC-B7807-18	18-400	Black Open top Cap, Polypropylene	100
CC-B7807-20	20-400	Black Open top Cap, Polypropylene	100
CC-B7807-24	24-400	Black Open top Cap, Polypropylene	100
CC-C4013-66	8-425	Black Open top Cap, Phenolic	1000
CC-C4015-66	13-425	Black Open top Cap, Phenolic	1000



PTFE/Silicone Septa for Sample Screw Caps

- Extremely clean, natural-coloured PTFE laminated to clear Silicone

Cat.No	PTFE	Silicone	Fits Cap	Pk
CC-B7995-8	0.25 mm	0.125 mm	8-425	100
CC-B7995-13	0.25 mm	0.125 mm	13-425	100
CC-B7995-15	0.25 mm	2.30 mm	15-425	100
CC-B7995-18	0.25 mm	2.30 mm	18-400	100
CC-B7995-20	0.25 mm	2.30 mm	20-400	100
CC-B7995-24	0.25 mm	2.30 mm	24-400	100
CC-B7995-26	0.125 mm	3.05 mm	24-400	100



PTFE-Lined Sample Storage Caps. Solid Top Cap

- PTFE on film/foam backing offers broad chemical resistance
- Not autoclavable • Solid top Cap

Cat.No	Thread Size	Color	Pk
CC-B7815-8	8-425	White	100
CC-B7815-13	13-425	White	100
CC-B7815-15	15-425	White	100
CC-B7815-18	18-400	White	100
CC-B7815-20	20-400	White	100
CC-B7815-24	24-400	White	100

Sample Storage Vials Convenience Kits & EPA Vials Tk



Non-Assembled Convenience Kits - Sample Vials with PTFE-Lined Storage Caps

- Each carton contains shrink-wrapped vials and PTFE-lined caps in plastic bags

Cat.No	Glass	Capacity mL	Capacity Drams	Use Cap Size	Vial Size mm	Pk
CC-B7800-1	Clear	2	0.5	8-425	12 x 32	100
CC-B-7800-2	Clear	4	1	13-425	15 x 45	100
CC-B7800-3	Clear	8	2	15-425	17 x 60	200
CC-B7800-12	Clear	12	3	15-425	19 x 65	200
CC-B7800-4	Clear	16	4	18-400	21 x 70	200
CC-B7800-20	Clear	20	5	18-400	21 x 70	100
CC-B7800-5	Clear	22	6	20-400	23 x 85	200
CC-B7800-6	Clear	40	10	24-400	28 x 95	100
CC-B7800-1A	Amber	2	0.5	8-425	12 x 32	100
CC-B7800-2A	Amber	4	1	13-425	15 x 45	100
CC-B7800-3A	Amber	8	2	15-425	17 x 60	200
CC-B7800-12A	Amber	12	3	15-425	19 x 65	200
CC-B7800-4A	Amber	16	4	18-400	21 x 70	200
CC-B7800-6A	Amber	40	10	24-400	28 x 95	100



Non-Assembled Convenience Kits - Sample Vials with Caps & Septa

- Each carton contains shrink-wrapped vials and separately packed caps
- Great packaging for storerooms or laboratory cabinets

Cat.No	Glass	Capacity mL	Capacity Drams	Use Cap Size	Vial Size mm	Pk
CC-B7990-1	Clear	2	0.5	8-425	12 x 32	100
CC-B-7990-2	Clear	4	1	13-425	15 x 45	100
CC-B7990-3	Clear	8	2	15-425	17 x 60	200
CC-B7990-12	Clear	12	3	15-425	19 x 65	200
CC-B7990-4	Clear	16	4	18-400	21 x 70	200
CC-B7990-5	Clear	22	6	20-400	23 x 85	200
CC-B7990-6	Clear	40	10	24-400	28 x 95	100
CC-B7990-1A	Amber	2	0.5	8-425	12 x 32	100
CC-B7990-2A	Amber	4	1	13-425	15 x 45	100
CC-B7990-3A	Amber	8	2	15-425	17 x 60	200
CC-B7990-12A	Amber	12	3	15-425	19 x 65	200
CC-B7990-4A	Amber	16	4	18-400	21 x 70	200
CC-B7990-6A	Amber	40	10	24-400	28 x 95	100



EPA Vials with Caps and Septa

- Each carton contains shrink-wrapped vials and separately packed caps
- Great packaging for storerooms or laboratory cabinets

Cat.No	Description	Pk
CC-B7950-B	Clear 40 mL vials, PTFE/Silicone Septa and Caps	100
CC-B7951-B	Amber 40 mL vials, PTFE/Silicone Septa and Caps	100
CC-B7920-B	Clear 20 mL vials, PTFE/Silicone Septa and Caps	100
CC-B7921-B	Amber 20 mL vials, PTFE/Silicone Septa and Caps	100




EPA Vial Components

Cat.No	Description	Pk
CC-B7950-VO	Clear, 40 mL, 28 x 95 mm vials	100
CC-B7951-VO	Amber, 40 mL, 28 x 95 mm vials	100
CC-B7920-VO	Clear, 20 mL, 28 x 57 mm vials	100
CC-B7921-VO	Amber, 20 mL, 28 x 57 mm vials	100
CC-B7950-1A	Caps Only, 24-400	100
CC-B7995-2A	Septa Only, 0.01 PTFE/0.09" Silicone	100

Tk Vial Comparison Chart

Target DP™ Vials - 12 x 32 mm




CC-C4000-1
2 mL

CC-C4000-9
1.5 mL

CC-C4000-LV1
350 µL

CC-C4000-V1
1.5 mL

Target™ Microvolumen Inserts



CC-C4010-630
6 x 31 mm
250 µL


CC-C4010-627L
6 X 31 mm
300 µL

CC-C4010-629L
6 x 31 mm
300 µL

CC-C4011-631
6 X 31 mm
400 µL

CC-C4010-631P
6 X 31 mm
300 µL

Target™ Snap-It-Vials - 12 x 32 mm



CC-C4011-5
2 mL

CC-C44011-24
700 µL

CC-C4011-LV1
350 µL

CC-C4011-V5
1.5 mL


CC-C4011-11
450 µL

CC-C4011-14
700 µL

CC-C4012-15
250 µL

CC-C4011-4
1.5 mL

Target™ Crimp Top Vials - 12 x 32 mm



CC-C4011-1
2 mL

CC-C4011-9
1 mL

CC-C4011-24
700 µL

CC-C4011-10
400 µL


CC-C4011-11
450 µL

CC-C4011-14
700 µL

CC-C4012-15
250 µL

CC-C4011-LV1
350 µL

Target™ Screw Thread Vials - 12 x 32 mm



CC-C4010-1
2 mL

CC-C4010-12
450 µL

CC-C4010-11
375 µL

CC-C4010-14
650 µL


CC-C4010-15
250 µL

CC-C4010-V1
15 mL

CC-C4010-LV1
350 µL


Vial Comparison Chart Tk

Standard Microvolume Inserts




CC-C4012-529L 5 X 31 mm 150 µL
CC-C4012-531L 5 X 31 mm 150 µL
CC-C4012-530 5 X 29 mm 150 µL
CC-C4012-465 5 X 31 mm 200 µL

15x45 4 mL Vial Microvolume Inserts




CC-C4015-638 6 x 38 mm 700 µL
CC-C4015-641 6 X 41 mm 500 µL
CC-C4015-642 6 x 42 mm 500 µL
CC-C4015-643 6 X 42 mm 500 µL

8x40 mm Vials 1 mL Vial Microvolume Inserts



CC-C4015-96A 5 X 36 mm 150 µL
CC-C4015-96PA 6 X 29 mm 200 µL
CC-C4015-536 5 X 36 mm 150 µL
CC-C4015-96P 6 X 29 mm 200 µL

Standard 4mL 15 x 45 mm



CC-C4015-1 4 mL
CC-C4015-48 4 mL
CC-C4015-96 1 mL


Standard Opening Screw Thread Vials 12 x 32 mm



CC-C4013-1 2 mL
CC-C4013-12 150 µL
CC-C4013-11 200 µL

8 mm Crimp Top Vials

Standard Opening Crimp Top Vials 12 x 32 mm



CC-C4012-1 2 mL
CC-C4012-10 150 µL

Headspace Vials - 23 x 75 mm & 23 x 46 mm



CC-C4020-20 20 mL
CC-C4020-10 10 mL
CC-C4020-2 20 mL
CC-C4020-210 10 mL